

# Anesthesia Machines





“Give me to drink mandragora...

That I might sleep out this great gap of time  
my Anthony is away.”

[excerpt from Anthony and Cleopatra]

- Mandrake (*Mandragora officinarum*) is a plant related to the potato family.

# Purpose



- Anesthesia units dispense a mixture of gases and vapors and vary the proportions to control a patient's level of consciousness and/or analgesia during surgical procedures.



# Functions



- Provide oxygen (O<sub>2</sub>) to the patient.
- Blend gas mixtures that can include (besides O<sub>2</sub>) an anesthetic vapor, nitrous oxide (N<sub>2</sub>O), other medical gases, and air.
- Facilitate spontaneous, controlled, or assisted ventilation with these gas mixtures.
- Reduce, if not eliminate, anesthesia-related risks to the patient and clinical staff.

# Anesthesia delivery



- The patient is anesthetized by inspiring a mixture of O<sub>2</sub>, the vapor of a volatile liquid halogenated hydrocarbon anesthetic, and, if necessary, N<sub>2</sub>O and other gases.
- Because normal breathing is routinely depressed by anesthetic agents and by muscle relaxants administered in conjunction with them, respiratory assistance — either with an automatic ventilator or by manual compression of the reservoir bag — is usually necessary to deliver the breathing gas to the patient.

# Principles of operation



- An anesthesia system comprises four basic subsystems:
  - a gas supply and control circuit;
  - a breathing and ventilation circuit;
  - a scavenging system;
  - a set of system function and breathing circuit monitors (e.g., inspired O<sub>2</sub> concentration, breathing circuit integrity).

# Safe practice of anesthesia



- Anesthesia machines incorporate a number of alarms that indicate:
  - levels and variations of several physiologic variables and parameters associated with cardiopulmonary function; and/or
  - gas and agent concentrations in breathed-gas mixtures.

# Safe practice of anesthesia



- Anesthesia machines must monitor:
  - O<sub>2</sub> concentration;
  - airway pressure; and either
  - the volume of expired gas ( $V_{exp}$ ); or
  - the concentration of expired CO<sub>2</sub> (capnography).
- Stand-alone monitors may be used to track other essential variables:
  - electrocardiogram;
  - SpO<sub>2</sub>;
  - blood pressure (invasive / non-invasive);
  - temperature.



# MAJOR COMPONENTS



- Gas Supply
- Pressure Regulators
- Flowmeters
- Vaporizers
- Safety Devices
- Breathing System

# Continuous-flow anesthesia system




Figure . Continuous-flow anesthesia system

Reproduced from Health Care Product Comparison System, ECRI. 2003 – Anesthesia Units

# Breathing circuits used in continuous-flow systems





Figure . Examples of breathing circuits

- Circle systems – advantages:
  - conserve a greater proportion of the anesthetic gases: ↓ cost;
  - conserve body heat and moisture from the patient.
  
- T-piece systems – advantages:
  - lower circuit compliance;
  - easier circuit sterilization;
  - less complex design requiring fewer valves and no CO<sub>2</sub> absorber (although one can be used with it).

Note: T-piece systems are used most often in pediatric anesthesia.


# Gas Supply





- Lack of O<sub>2</sub> delivered to the patient (Hypoxia):
  - can result in brain damage or death.
  
- Administration of O<sub>2</sub> in a concentration of 100%, even for a short duration, may be toxic:
  - resorption atelectasis;
  - particularly acute in neonatal anesthesia; can cause retrolental fibroplasia and bronchopulmonary dysplasia.


# Gas Supply




# Pressure Regulators



Closed




Open





# Flowmeters






N<sub>2</sub>O and O<sub>2</sub> flow controls are interlocked so that the proportion of O<sub>2</sub> to N<sub>2</sub>O can never fall below a minimum value (nominal 0.25) to produce a hypoxic breathing mixture.

# Vaporizers




- Inhaled anesthetic agents, with the exception of N<sub>2</sub>O, exist as liquids at room temperature and sea-level ambient pressure.
- Vaporizers add a controlled amount of anesthetic vapor to the gas mixture.
- Types of vaporizers:
  - variable bypass (conventional);
  - heated blender;
  - Measured flow;
  - draw-over.

# Vaporizers







# Vaporizers



# Safety Devices

## Diameter Index Safety System



# Safety Devices


## Pin Index Safety System

**Gas**  
Oxygen  
Nitrous Oxide  
Air  
Cyclopropane


**Index Pins**  
2-5  
3-5  
3-6  
3-6



Front View




Side View



# Safety Devices

## Oxygen Pressure Failure Devices




# Safety Devices



ORBIS



saving sight worldwide

## Vaporizer Interlock






# Bain System – a common configuration of the T-piece system




Inhalation




Exhalation

# Inhalation-Assisted Manual Ventilation




From Patient

# Exhalation-Assisted Manual Ventilation




To Patient

# Inhalation-Ventilator Controlled Ventilation



From Patient

# Exhalation-Ventilator Controlled Ventilation



To Ventilator

To Scavenger

Fresh Gas

# CO<sub>2</sub> Absorbers



- soda lime
- barium hydroxide lime

# O2 Monitor



- An O2 monitor located on the inspiratory side of the breathing circuit analyzes gas sampled from the Y-piece of the patient's breathing circuit and displays O2 concentration in volume percent.
- O2 monitors sound an alarm if the O2 concentration falls below the preset limit.


# Adjustable Pressure Limiter (APL) Valve



- Pressure imposed on the patient's lungs can cause serious lung damage.
- Either an APL valve or a valve in the ventilator allows excess gas to escape when a preset pressure is exceeded.
- Types:
  - spring-loaded;
  - needle.
- Many APL valves do not have calibrated markings: The anesthetist must adjust them empirically to give a desired peak inspired pressure.
- Circle systems and T-piece systems also include a pressure gauge for monitoring circuit pressure and setting the APL valve.



# Adjustable Pressure Limiter (APL) Valve



# Scavenging System - Rationale



- Captures and exhausts waste gases to minimize the exposure of the operating room staff to occupational risks.
- Exposure to trace levels of anesthetic gases continually present in the operating room can cause adverse health effects in operating room personnel:
  - increased incidence of spontaneous abortion;
  - congenital anomalies in babies.
- Trace gas levels in the air may have a slight anesthetizing effect on the anesthetist and surgeon.

# Scavenging System



- Scavenging systems remove gas by a vacuum, a passive exhaust system, or both.
- Note: Inadequate evacuation of some scavenging systems can cause pressure to build up in the breathing circuit, with the potential for pneumothorax (air in the pleural cavity).

# Scavenging System



## Vacuum (active):

- Vacuum scavengers use the suction from an operating room vacuum wall outlet or a dedicated vacuum system.
- To prevent positive or negative pressure in the vacuum system from affecting the pressure in the patient circuit, manifold-type vacuum scavengers use one or more positive or negative pressure-relief valves in an interface with the anesthesia system.
- Open-type vacuum scavengers have vacuum ports that are open to the atmosphere through some type of reservoir; such units do not require valves for pressure relief.


# Scavenging System



## Passive exhaust system:

- Passive-exhaust scavengers can vent into a hospital ventilation system (if the system is the non-recirculating type) or, preferably, into a dedicated exhaust system.
- The slight pressure of the waste-gas discharge from the anesthesia machine forces gas through large bore tubing and into the disposal system or directly into the atmosphere.

# Scavenging System



# Preventive Maintenance



- Test apparatus and supplies:
  - Lung simulator with adjustable compliance or ventilator tester
  - Pressure gauge or meter with 2 cm H<sub>2</sub>O resolution, from -20 to +120 cm H<sub>2</sub>O
  - Various breathing circuit adapters
  - Leakage current meter or electrical safety analyzer
  - Ground resistance ohmmeter
  - Additional items as required for specific manufacturers' procedures

# Preventive Maintenance



- Qualitative tests:

- Chassis/Housing
- Mount/Fasteners
- Casters/Brakes
- AC Plug
- Line Cord
- Strain Reliefs
- Circuit Breaker/
- Tubes/Hoses
- Cables
- Fittings/Connectors
- Filters
- Controls/Switches
- Fan
- Battery/Charger
- Indicators/Displays
- Alarms/Interlocks
- Labeling
- Accessories
- Bellows



# Preventive Maintenance



- Quantitative tests:
  - Grounding resistance [ $\leq 0.5 \Omega$ ]
  - Leakage current [ $\leq 300 \mu\text{A}$  chassis]
  - Modes and settings [ $\pm 10\%$  accuracy]
  - Monitors and Alarms [ $\pm 10\%$  accuracy]
  - Alarms tested:
    - Airway pressure
    - Tidal volume
    - FIO<sub>2</sub>
    - Others

# Preventive Maintenance



- Others:
  - Gas Supply
  - Pneumatic lines (including air filters)
  - Gas cylinders (and gauges and regulators, if so equipped)
  
  - Patient Circuit  $\Gamma$
  - Breathing circuit (including filters)
  - Humidifiers
  - Pressure-relief mechanism
  - Absorber

# IEC Standards



- International Electrotechnical Commission. Medical electrical equipment — part 1: general requirements for safety [standard]. IEC 60601-1 (1988-12). 1988.
- Medical electrical equipment — part 1: general requirements for safety. Amendment 1 [standard]. IEC 60601-1-am1 (1991-11). 1991.
- Medical electrical equipment — part 1: general requirements for safety. Amendment 2 [standard]. IEC 60601-1-am2 (1995-03). 1995.
- Medical electrical equipment — part 1: general requirements for safety. Section 1. Collateral standard: safety requirements for medical electrical systems. IEC 60601-1-1 (1992-06). 1992.
- Medical electrical equipment — part 1: general requirements for safety. Section 1. Collateral standard: safety requirements for medical electrical systems. Amendment 1. IEC 60601-1-1-am1 (1995-11). 1995.
- Medical electrical equipment — part 1: general requirements for safety. Section 2. Collateral standard: electromagnetic compatibility — requirements and tests. IEC 60601-1-2 (2001-09). 2001.
- Medical electrical equipment — part 2-13: particular requirements for the safety of anesthetic workstations [standard]. IEC 60601-2-13 (1998-05). 1998.

# ISO Standards




- International Organization for Standardization. Anaesthesia and respiratory care alarm signals — part 1: visual alarm signals [standard]. 1st ed. ISO 9703: Part 1:1992. 1992.
- Anaesthesia and respiratory care alarm signals — part 2: auditory alarm signals [standard]. 1st ed. ISO 9703-2:1994. 1994.
- Anaesthetic and respiratory equipment — conical connectors — part 1: cones and sockets [standard]. 2nd ed. ISO 5356:1-1996. 1987 (revised 1996).
- Anaesthetic and respiratory equipment — conical connectors — part 2: screw-threaded weight-bearing connectors [standard]. 1st ed. ISO 5356-2:1987. 1987.
- Anaesthetic and respiratory equipment — heat and moisture exchangers (HMEs) for humidifying respired gases in humans [standard]. ISO 9360:2001. 2001.
- Anaesthetic machines for use with humans [standard]. 2nd ed. ISO 5358:1992. 1992.




ORBIS

saving sight worldwide













ORBIS

saving sight worldwide







ORBIS

saving sight worldwide



