## Ventilog 2

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	Page
<b>For Your Safety and That of Your Patients</b> .....	3
<b>Intended Use</b> .....	4
<b>Design</b> .....	5
<b>Preparation</b> .....	7
Positions for Ventilog 2 .....	7
Connecting to gas supply .....	7
Fitting anaesthetic gas scavenging system .....	8
Connecting Ventilog 2 to Circle System 9 .....	8
<b>Operation</b> .....	10
Set fresh gas flow on anaesthetic machine .....	10
Automatic ventilation .....	10
Manual ventilation .....	13
Spontaneous breathing .....	14
Shut-down .....	14
<b>Fault – Cause – Remedy</b> .....	15
<b>Care</b> .....	16
Stripping down Ventilog 2 .....	16
Disinfecting/cleaning .....	18
Assembling Ventilog 2 .....	19
<b>Testing Readiness for Operation</b> .....	20
Leak test .....	20
Checking manual ventilation .....	21
Checking automatic ventilation .....	22
Checking PEEP .....	23
<b>Maintenance Intervals</b> .....	23
<b>What's What</b> .....	24
<b>Technical Data</b> .....	27
<b>Abbreviations Used/Terminology</b> .....	29
<b>Index</b> .....	30

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## For Your Safety and that of Your Patients

For correct and effective use of the apparatus and to avoid hazards it is essential to read the following recommendations and to act accordingly<sup>1)</sup>:

### Strictly follow the Instructions for Use

Any use of the apparatus requires full understanding and strict observation of these instructions. The apparatus is only to be used for purposes specified here.

### Maintenance

The apparatus must be inspected<sup>2)</sup> and serviced<sup>2)</sup> regularly by trained service personnel at six monthly intervals (and a record kept). Repair<sup>2)</sup> and general overhaul of the apparatus may only be carried out by trained service personnel. We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them. Only authentic Dräger spare parts may be used for maintenance<sup>2)</sup>. Observe chapter "Maintenance Intervals".

### Not for use in areas of explosion hazard

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

1) Insofar as reference is made to laws, regulations or standards, these are based on the legal system of the Federal Republic of Germany.

2) Definitions according to DIN 31 051:  
 Inspection = examination of actual condition  
 Service = measures to maintain specified condition  
 Repair = measures to restore specified condition  
 Maintenance = inspection, service, repair

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**Liability for proper function or damage**

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus is serviced or repaired by personnel not employed or authorized by DrägerService or if the apparatus is used in a manner not conforming to its intended use.

Drägerwerk Aktiengesellschaft cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Drägerwerk Aktiengesellschaft are likewise not modified by the recommendations given above.

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**Intended Use**

Ventilog 2 – an anaesthetic ventilator for use with Dräger inhalation anaesthesia machines with an Circle System 9.

For adults and children – use with bellows E.

For neonates and infants – use with bellows K<sup>1)</sup>.

Because Ventilog 2 is pneumatically controlled and driven, it is particularly suitable for use with flammable gases. However, in situations where electrical monitors and non-conductive hoses are being used with Ventilog 2, flammable gases must not be used.

**In any situation where flammable gas or anaesthetic agents are going to be used, do not use anaesthetic machines with a Ventilog 2 in conjunction with non-conductive hoses, electrical monitors or combustible equipment. Fire risk.**

1) Included in "Infant bellows" modification kit (Option).

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

The switch for I:E ratio also switches the machine on and off.

### In the middle

- Window with scale for setting stroke volume and observing movement of bellows.
- Connection point for control line for ISO switching valve for manual ventilation or spontaneous breathing.

### At the top


- Connector for hose to Circle System 9 and locking system for pressure chamber with detachable patient system.

To open pressure chamber:

- 1 pull knob out – as far as it will go, then
  - 2 turn anti-clockwise:
- the pressure chamber with the patient system is released from its housing.
- 3 Lift handle up and take patient system out – for cleaning and sterilizing or to replace bellows.

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
**At the back**

- 1 Air/O<sub>2</sub> connector for medical gas supply.
- 2 Connector for anaesthetic gas scavenging.
- 3 Connecting cable to monitor switch (**optional**) to switch monitor on automatically.

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**ISO Switching valve**

connects Ventilog 2 to Circle System 9.

This switches over the operating modes which are automatically controlled from Ventilog 2:  
 automatic ventilation  
 or  
 manual ventilation / spontaneous breathing.

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
### Positions for Ventilog 2

#### On Sulla 909 V

- On the base plate – at the left, when viewed from the front. The two bolts on Ventilog 2 fit into the holes on the base plate. Screw Ventilog 2 on firmly from below with tommy screw.

#### On Trajan 808

- On the mounting plate – at the right, when viewed from the front. The two bolts on Ventilog 2 fit into the holes on the mounting plate. Screw Ventilog 2 on firmly from below with tommy screw.

### Ventilog 2 must not be placed above electrically-operated machines. Danger of fire from the oxygen or oxygen-enriched air being discharged.

Electrical machines may be placed on top of Ventilog 2 – using the Dräger latching system.

Assembling Ventilog 2, see page 20.

### Connecting to gas supply

#### on Sulla 909 V with medical air from Sulla 909 V


- Screw medical air hose onto back of Ventilog 2 and push into the coupling on Sulla 909 V – as far as it will go.

When Sulla 909 V is not being supplied with medical air, or if the air supply pressure falls below 2 bar, Ventilog 2 is automatically driven with oxygen from Sulla 909 V.

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
- Preparation
- Connecting to gas supply
- Fitting anaesthetic gas scavenging system
- Connecting Ventillog 2 to Circle System 9

### On Trajan 808 with medical air from wall terminal unit

- Screw medical air hose onto back of Ventillog 2, push probe into medical air wall terminal unit.

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
### Fitting anaesthetic gas scavenging system

- Push scavenging hose onto connector on back of Ventillog 2.
- Fit another scavenging hose to connector below the pressure limiting valve on the Circle System 9.
- Connect the two scavenging hoses with a T-piece.
- Fit scavenging hose with perforated ring onto T-piece.
- Do not connect probe of scavenging hose to the terminal unit on the scavenging system until anaesthetic operation begins.

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
### Connecting Ventillog 2 to Circle System 9

Use ISO switching valve (M 29 900) only, other switching valves may cause faulty operation.

- Screw ISO switching valve onto connector on Circle System 9 – as far as it will go. Secure tightly with knurled nut.
- Screw fresh gas hose onto Circle system 9.

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
- 1 Screw hose onto connectors of  
 Ventilog 2 and ISO switching  
 valve and secure tightly with  
 knurled nuts.
- 2 Fit corrugated hose on breathing  
 bag to connector on ISO switch-  
 ing valve.

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
- 3 Connect control line for ISO  
 switching valve to connector on  
 valve and to front of Ventilog 2 =  
 push probe into socket – until it  
 engages.

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### For neonates and infants

- Use bellows K from "Infant bel-  
 lows" modification kit for stroke  
 volumes of 50 to 150 mL  
 and:
- use both hoses and breathing bag  
 from "ISO Infant Circle System"  
 modification kit.

On ISO switching valve

- Connect 0.5 L breathing bag to  
 the switching valve with hose  
 from the modification kit.

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**Operation**  
Set fresh gas flow on anaesthetic machine  
Automatic ventilation

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
## Operation

### Set fresh gas flow on anaesthetic machine

- 2 to 8 L/min

### Automatic ventilation

1 First set rotary knobs for Frequency PEEP Flow to green dots (initial setting).

Then:

- 2 set stroke volume on the scale to suit the patient using rotary knob. The scale relates to bellows E (adults and children) fresh-gas flow of 4 L/min (with  $f = 10$  /min, I:E = 1:2).
- 3 Switch on Ventilog 2 = set switch to desired I:E ratio. The ISO switching valve connects the bellows on Ventilog 2 to the Circle System 9.

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When using bellows K (for neonates and small children) from the children's bellows modification set, the scale for setting stroke volume that is given on the window is not correct. The scale enclosed with the children's bellows modification set must be used instead.

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For Ventilog 2, the scale given for stroke volume relates to a fresh-gas flow of 4 L/min (with  $f = 10$  /min, I:E = 1:2) flowing continuously in the circle system. At different flows, the resulting stroke volume will deviate from the scale range.

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Therefore:  
correct the stroke volume set for the bellows according to the following tables.

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

Example: with  $f = 10$  /min, I:E = 1:2 and a fresh-gas flow of 2 L/min,  
increase stroke volume by about 65 mL and check tidal volume on monitor.

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
### Automatic pressure limited ventilation

On Circle System 9

- 1 Move lever of pressure limiting valve over – to »MAN«.
- 2 Set pressure limit desired = rotate lever and set value on scale. A pressure reserve of about 70 mbar is available for the ventilator. When setting a low pressure limit – increase the fresh-gas flow.

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
### Rapid venting for instance, for flushing or for releasing pressure in Circle System 9 quickly

- Push lever of pressure limiting valve right down. The pressure in the Circle System 9 is immediately reduced via gas scavenging connector.

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### Adjusting PEEP

- Turn rotary knob for PEEP clockwise; check setting on airway pressure monitor.

The PEEP value stabilizes after about 5 ventilation strokes.

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For PEEP = 0:

- turn rotary knob anti-clockwise – as far as it will go.

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### Hints on adjusting automatic ventilation

#### If bellows are not deflating completely:

increase inspiratory flow  
or  
reduce frequency  
or  
increase inspiratory time with I:E ratio.

#### If bellows do not expand to the lower stop when being filled:

increase fresh gas flow  
or  
increase pressure limit on pressure limiting valve on Circle System 9.

#### If the measured tidal volume is too small in relation to the stroke volume setting:

increase pressure limit on pressure limiting valve on Circle System 9.

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

- 1 turn I:E switch to 0.

On Circle System 9

- 2 Move lever on pressure limiting valve over – to »MAN«.
- 3 Set pressure limit desired = rotate lever and set value on scale.

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
## Spontaneous breathing

- 1 Turn I:E switch to 0:  
Ventilog 2 is switched off.

On Circle System 9

- 2 Move lever of pressure limiting valve over – to »SPONT«. Valve is open – independent of pressure limit set. Pressure cannot build up in Circle System 9.

- Supply an adequate flow of fresh gas to keep breathing bag full.

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## Shut-down

- Turn I:E switch to 0.
- Set fresh gas flow to 0, follow Instructions for Use of anaesthetic machine.
- Disconnect probe of anaesthetic gas scavenging system.

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## Fault – Cause – Remedy

Fault	Cause	Remedy
Bellows do not expand to lower stop	Fresh gas flow too low	Set higher fresh gas flow (at least 2 L/min)
	Leakages in Circle System 9	Remove leakages
Bellows do not deflate	Insp. flow too low	Increase insp. flow using rotary knob "Flow"
	Kinked hose on Circle System 9	Straighten out hose
	Stenosis on tube	Correct stenosis
	Control line for ISO switching valve is kinked or disconnected	Connect control line and straighten out kinks
	Faulty ISO switching valve	Use another ISO switch valve or call DrägerService
No pressure in Circle System 9 during manual ventilation	Lever on pressure limiting valve of Circle System 9 to »SPONT« setting	Lever to »MAN« position and set desired pressure limit
	Leak in Circle System 9	Leak test Circle System 9
	Faulty ISO switching valve	Use another ISO switch valve or call DrägerService
PEEP does not remain constant	Leak in Circle System 9 e.g. leak in conical connector	Leak test Circle System 9

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

## Care

Carry out after each day of anaesthesia.

Care procedures for Circle System 9, anaesthetic machine and accessory machines must be carried out in accordance with relevant Instructions for Use.

### Stripping down Ventillog 2

- Disconnect probe for compressed gas supply.
- 1 Unscrew connecting hose between Ventillog 2 and Circle System 9.
  - 2 Disconnect hose on breathing bag from ISO switching valve.

On the front of Ventillog 2 and on ISO switching valve:


- Detach control line for ISO switching valve, on Ventillog 2: push socket back and at the same time pull out probe, on ISO switching valve: pull the control line off connector.
- Unscrew ISO switching valve from Circle System 9.

Do not dismantle ISO switching valve.

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
**On the back on Ventilog 2**

- Pull anaesthetic gas scavenging hose off the gas scavenging connector.
- 1 Swivel gas scavenging connector into dotted position = unlock, then pull out.

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
**On the front plate of Ventilog 2**

- 2 Pull knob out – as far as it will go and
  - 3 turn anti-clockwise:
- patient system and pressure chamber are released from the housing.
- 4 Lift handle up and take patient system out.

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- Compress bellows completely, and in this compressed condition, pull off from patient system twisting gently to the right and left.

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## Disinfecting/cleaning

For disinfection purposes use preparations from the surface-disinfectant group. For reasons of material compatibility suitable preparations are those based on

- aldehydes
- quaternary ammonium compounds.

Because of the possibility of damage to materials preparations based on the following should not be used

- halogen-releasing compounds,
- strong organic acids,
- oxygen-releasing compounds.

For users in the Federal Republic of Germany we recommend particularly the disinfectants contained in the latest DGHM list (DGHM: German Society for Hygiene and Microbiology).

The DGHM list (mhp-Verlag GmbH, Wiesbaden) also indicates the active substances on which each disinfectant is based. For countries in which the DGHM list is not known, we recommend the above-mentioned substances.

### Wipe disinfecting

Ventilog 2 housing, medical gas connecting hose and control lines of ISO switching valve.

Do not use alcohol-based agents.

- Remove any obvious stains with a disposable cloth.
- Wipe disinfect, with Buraton 10 F (Messrs. Schülke & Mayr, Norderstedt), for instance. Follow manufacturer's instructions.

### Disinfect with moist heat in cleaning and disinfecting machine

Connecting hose, bellows, patient system, ISO switching valve, gas scavenging connector.

- Disinfect with moist heat (93 °C/10 min) in cleaning and disinfecting machine **using cleaning agents only. Do not add disinfectants! Components are subject to corrosion.**
- **Always sterilize the ISO switching valve in hot steam at 134 °C after any moist heat disinfection procedure. Otherwise, it may malfunction due to fluid residues which may remain in the control channels of the switching valve.**

### Sterilizing


Connecting hose, bellows, patient system, ISO switching valve, gas scavenging connector.

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

## Assembling Ventillog 2

- Compress bellows and using a twisting motion push evenly onto the connector on the patient system **as far as it will go**.

Bellows which have not been pushed in completely may move unevenly in the pressure chamber, and may, therefore, not be emptied and filled properly.

- 1 Hang bellows in the pressure chamber and place patient system on the pressure chamber. The two guide bolts on the pressure chamber fit into the appropriate holes on the patient system.
- 2 Lie handle down.
- Push pressure chamber containing patient system into the housing – as far as it will go.


- Move the gas scavenging connector on Ventillog 2 to the dotted position and then swivel into the vertical position = to lock.

- Connect Ventillog 2 to Circle System 9, see page 9.

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## Testing Readiness for Operation

Carry out test after each service and after completely assembling Ventilog 2, anaesthetic machine with Circle System 9 and ISO switching valve.

- 1 Rotary knobs for  
Frequency  
PEEP  
Flow  
I:E  
on green dot = initial setting.
- 2 Stroke volume to about 500 mL.
  - Switch on airway monitor.

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## Leak test

for Ventilog 2, ISO switching valve and Circle System 9.

- Leak test Circle System 9 according to relevant Instructions for Use.
- Fresh gas flow on anaesthetic machine to 0.
- Leave Y-piece open.
- At the beginning of expiration keep Y-piece sealed:
- The bellows should remain hanging freely = the system is sealed.


If not:

- leak test bellows  
or
- replace ISO switching valve  
or
- replace gas scavenging connector on Ventilog 2.
- Leak test hose connecting Ventilog 2 to ISO switching valve.

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## Checking manual ventilation


On Circle System 9

- Pressure limiting valve to »MAN« setting and set scale to 30 mbar.
- Set fresh gas flow on anaesthetic machine to 4 L/min.
- Turn I:E switch to 0.
- Keep Y-piece sealed and compress breathing bag:
- The airway pressure – measured on airway monitor, should increase immediately.

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## Checking automatic ventilation

On Circle System 9

- Pressure limiting valve to »MAN« setting and set scale to 30 mbar.
- Set fresh gas flow on anaesthetic machine to 4 L/min.
- Set rotary knobs:  
Frequency 6/min  
PEEP 0  
Flow to green dot  
I:E 1:1  
 $V_T$  to about 400 mL
- Attach an additional breathing bag, such as a 2.3 L bag to Y-piece as test lung.
- During inspiration the bellows should deflate evenly – up to the upper stop, airway pressure should build up in the Circle System 9.
- The tidal volume  $V_T$  – measured on the airway monitor should match the stroke volume set.
- During expiration, the bellows should fall evenly to the lower stop. The airway pressure in the Circle System 9 should fall to about 1 mbar.

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### If the bellows fill and empty unevenly and jerkily:


check bellows connection on patient system. Push neck of bellows completely and evenly onto the nozzle on the patient system.

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## Checking PEEP

- Turn "PEEP" rotary knob to the right stop.
- The endexpiratory airway pressure should be at least 15 mbar.
- Turn "PEEP" rotary knob back to the left stop.
- Set I:E switch to 0 again.

**The ventilator must not be used if any test requirement has not been met.**

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## Maintenance Intervals

**Clean and disinfect machine and machine parts every time before maintenance is carried out – also before returning for repair.**


Bellows	replace every 2 years
Sealing of pressure chamber	replace every 2 years
Pressure reducer	basic overhaul by trained service personnel every 6 years, we recommend DrägerService
Machine inspection and service <sup>1)</sup>	every six months by trained service personnel
Regular technical safety checks	every six months by trained service personnel

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<sup>1)</sup>  
In accordance with DIN 31 051:  
Inspection = examination of actual condition  
Service = measures to maintain specified condition  
Repair = measures to restore specified condition  
Maintenance = inspection, service and, if applicable, repair

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
- 1 Pressure chamber
- 2 Rotary knob for stroke volume, scale for bellows E next to window
- 3 Rotary knob for ventilation frequency
- 4 Rotary knob for PEEP
- 5 Rotary knob for inspiratory flow
- 6 Switch for I:E ratio and OFF switch (0)
- 7 Rotary knob to unlock pressure chamber
- 8 Patient system
- 9 Connector for hose to Circle System 9
- 10 Coupling for control line to ISO switching valve

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
### Back view

- 11 Connector for scavenging system
- 12 Air/O<sub>2</sub> connector for medical gas supply
- 13 Monitor switch (optional)

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### ISO switching valve

- 14 Connector for hose to Ventilog 2
- 15 Connector for control line on ISO switching valve
- 16 Connector for breathing bag hose
- 17 Connector for Circle System 9

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## Technical Data

### Ambient conditions

During operation	
Temperature	15 to 35 °C
Atmospheric pressure	700 to 1060 hPa
Rel. humidity	20 to 80%
During storage	
Temperature	-10 to 60 °C
Atmospheric pressure	500 to 1060 hPa
Rel. humidity	0 to 80%

<b>Functioning principle</b>	Bellows in Bottle (bellows with primary/ secondary system)
------------------------------	--

<b>Control method</b>	pneumatic, time-cycled, volume-constant
-----------------------	--

### Set values

Ventilation frequency <sup>1)</sup>	6 to 60/min ± 15%
I:E ratio	1:1, 1:2, 1:3 ± 20%
Inspiratory flow	20 to 80 L/min ± 15% at 20 mbar counterpressure
Stroke volume with bellows K <sup>2)</sup> for neonates and infants	50 to 150 mL
with bellows E for adults and children	100 to 1600 mL
Minute volume (according to ISO)	25 L/min at 4 L/min fresh gas flow and I:E=1:2
Fresh gas flow	2 to 8 L/min
PEEP	approx. 1 to >15 mbar

1) The control elements for ventilation frequency have been calibrated with medical air. When oxygen is used as drive gas, there is a reduction in frequency of about 10%.

2) From "Infant bellows" modification kit (Option).

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**Machine data**

Working pressure in primary system	80 mbar (not changeable)
Internal compliance with bellows K <sup>1)</sup> and 1 m connecting hose	0.83 mL/mbar
with bellows E and 1 m connecting hose	1.3 mL/mbar
with bellows E, 1 m connecting hose and Circle System 9 (2 absorbers)	4.7 mL/mbar
Medical gas supply	2 to 6 bar oil-free medical air or oxygen from a medical gas pipeline system or from medical gas cylinders
Connecting thread	M 15 x 1
Medical gas consumption	16 L/min (at I:E = 1:2)
Patient system	as a slide-in unit is detachable and can be autoclaved
Anaesthetic gas scavenging system	via connector at back
Dimensions W x H x D mm	212 x 266 x 300
Weight	12 kg
Noise emission Sound pressure (relative to free field measurement over reflecting surface)	maximum 51 dB (A)

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**Abbreviations Used/Terminology**

Bellows in Bottle	Bellows pressure chamber system
DGHM	German Society for Hygiene and Microbiology
Flow	inspiratory flow applied by Ventillog 2
I:E	ratio inspiration time:expiration time
ISO	International Standards Organization (for medical equipment)
PEEP	positive endexpiratory pressure
V <sub>T</sub>	Tidalvolume

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These Instructions for Use apply only to  
**Ventilog 2**  
with Serial No.:

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If no Serial No. has been filled in by  
Dräger these Instructions for Use are  
provided for general information only and  
are not intended for use with any specific  
machine or device.

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☰ Moislinger Allee 53 – 55  
23542 Lübeck  
☎ (04 51) 8 82 - 0  
☎ 2 68 07 - 0  
FAX (04 51) 8 82-20 80

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Germany  
☰ Moislinger Allee 53 – 55  
D-23542 Lübeck  
☎ (4 51) 8 82 - 0  
☎ 2 68 07 - 0  
FAX (4 51) 8 82-20 80

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