

HMC HIRAYAMA

AUTOCLAVE

HICLAVE

HRG Series

SERVICE MANUAL

株式会社 平山製作所



AUTOCLAVE

HICLAVE

HRG Series

SERVICE MANUAL



IMPORTANT

- Read this manual carefully and follow the instructions to keep the machine in good working condition.

Issued on June 25, 2012, Revised on July 29, 2013

Contents

CONTENTS	1
CHAPTER 1. PROGRAM SPECIFICATION	2
1. SPECIFICATION	2
2. INITIAL SETTING VALUE	3
3. DESCRIPTION OF ERROR MESSAGES	3
4. ERROR MONITORING CHART	4
4.1 FABRIC cycle	4
4.2 LIQUID cycle	5
5. OPERATION	6
5.1 Turning on the power supply	6
5.2 Changing a cycle	7
5.3 Changing a setting	8
5.4 Confirming or setting the auto-start time	9
5.5 Setting the air removal time	10
5.6 Setting the door lock temperature	10
5.7 Setting the floating sensor	11
5.8 Confirmation and correction of the clock	11
5.9 Confirming the floating sensor temperature	13
5.10 Power-saving operation	13
5.11 Error detection	13
5.12 Printer output	14
6. TIMING CHART	15
● Explanation of notes on the timing chart	15
6.1 FABRIC Cycle	16
6.2 SOLID Cycle	17
6.3 LIQ Cycle	18
6.4 ADD DRY Cycle	19
6.5 Vacuum Test Cycle	20
CHAPTER 2. MAINTENANCE AND SERVICE	22
1. REPLACING THE DOOR GASKET	22
2. REPLACING THE SOLID STATE RELAY (SSR)	22
3. REPLACING THE BACKUP BATTERY	23
4. REPLACING THE CONTROL TEMPERATURE SENSOR	23
5. REPLACING THE ROM	24
6. CAUTIONS WHEN REPLACING THE CONTROL BOARD (SWITCHBOARD)	25
CHAPTER 3. REFERENCE DRAWINGS	26
● EXTERNAL VIEW	26
● RIGHT SIDE VIEW	27
● LEFT SIDE VIEW	28
● SWITCHBOARD	29
● WIRING DIAGRAM	30
● PIPING DIAGRAM	31
CHAPTER 4. OPERATION CHECK PROGRAMS	32
1. OUTLINE	32
2. STARTUP METHOD	32
3. ITEMS TO CHECK BEFOREHAND	32
4. OPERATION TEST	33
C1: Checking if LEDs light	33
C2: Checking the switch input	34
C3: Checking the contact input	35
C4: Checking the relay output	35
C5: Checking or adjusting the input of the CN9 control temperature sensor	36
C6: Checking or adjusting the input of the CN10 floating sensor	36
C7: Checking the input of the CN8 pressure sensor	36
C8: Checking the input of the CN4-10 heater temperature sensor	36

Chapter 1. Program Specification

1. Specification

Model	HRG-112	HRG-140	Notes																																
Name of a cycle	FABRIC, SOLID, LIQ(Liquid), ADD DRY																																		
Number of programs	3 programs for each cycle																																		
Sterilization temperature setting range	105~135°C		Least input increment : 1°C																																
Sterilization time setting range	1 min. - 60 min.		Least input increment : 1min.																																
Dry time setting range	1 min. - 99 min.		Least input increment : 1min.																																
Exhaust level setting	4 steps 0: Natural cooling 1:Pulse exhaust (fine exhaust) 2:Pulse exhaust (small amount of exhaust) 3:Exhaust with a valve fully open																																		
Air removal method	•FABRIC, SOLID : Vacuum air removal which is performed twice at -0.085MPa (Vacuum suction starts from -0.075MPa for 2 min.) •LIQ : Vacuum air removal and Gravity displacement air removal. Vacuum air removal is performed at -0.090MPa. (Vacuum suction starts from -0.075MPa for 2 min.), and, when temperature reaches 103°C, the gravity displacement air removal is performed until the preset time elapses.																																		
Air removal time	5 min. - 10 min.		Least input increment : 1min.																																
Door lock temperature	SOLID cycle : 97°C (fixed) LIQ cycle : 60.0 - 95°C (Temperatures for FABRIC and ADD DRY cycles are not specified)		Least input increment : 1°C																																
Over temperature alarm	138°C																																		
Programmable timer	6 day/24 hour programmable timer																																		
S1 dip switch setting	<div style="display: flex; align-items: center;"><div style="margin-right: 10px;">ON</div><div style="text-align: center;"><div>①②③④⑤⑥⑦⑧</div><table border="1" style="border-collapse: collapse; width: 100px; height: 30px;"><tr><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr></table></div></div> <div style="display: flex; align-items: center; margin-top: 5px;"><div style="margin-right: 10px;">OFF</div><table border="1" style="border-collapse: collapse; width: 100px; height: 30px;"><tr><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr><tr><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td><td> </td></tr></table></div>																																		①Model setup ②Model setup (spare) ③Floating sensor ④Option (spare) ⑤Option (spare) ⑥Vacuum test ⑦Safety valve test ⑧Operation check program

2. Initial setting value

Item Cycle	Sterilization temperature	Sterilization time	Exhaust level	Drying time	Air removal time	Lid lock temperature
FABRIC	121°C	20 min.	3 (fixed)	40 min.	-	-
SOLID	121°C	20 min.	3 (fixed)	-	-	97°C (fixed)
LIQUID	121°C	20 min.	0	-	5 min.	60°C
ADD DRY	-	-	-	40 min.	-	-

3. Description of error messages

Error No.	Description of Displayed Error	Action
E 1 (Overheat)	• Heater surface temperature exceeds 270°C.	① Turn off the heater circuit. ② Error buzzer ③ Display of chamber temperature and chamber pressure
E 2 (Disconnection of control temperature sensor wire)	• Disconnection of control temperature sensor wire is detected or temperature is below 0°C	
E 3	(Over-temperature) • The chamber temperature exceeds the upper limit + 3°C of the settable temperature.	
	(Over-temperature during sterilization stage) • The chamber temperature exceeds the preset temperature + 5°C for more than 10 seconds during sterilization stage.	
E 4 (Over-cooling)	• The chamber temperature drops below 102°C for more than 10 seconds during sterilization stage	
E 5	(Over-pressure) • The chamber pressure exceeds the pressure that was converted from the preset temperature by more than 25kPa for more than 15 seconds.	
	(Abnormality in the pressure sensor) • Pressure of more than 0.01MPa or pressure of less than -0.01MPa is indicated when the door is open	
E 6	(Abnormality in the door) • The door lowers during operation.	
E 8 (Abnormality in the water supply and drainage)	• Water is not filled within 10 minutes after startup. • Drainage is not completed within 10 minutes after start of drainage. • Low water is indicated when water level is full. • Full water is indicated when water level is low.	
E 9 (Heater abnormality)	• The chamber temperature does not reach the preset temperature within 300 minutes after startup.	
E P (Vacuum abnormality)	• The degree of vacuum does not reach -0.075MPa within 15 minutes when vacuum air removal is performed. • Air leakage of more than 0.13kPa/min. is detected during a vacuum test.	
E F (Disconnection of floating sensor wire)	• Disconnection of floating sensor wire is detected or temperature is below 0°C. (when a floating sensor is installed)	
E L (Abnormality in the door lock)	• The door lock signal is detected when the door is open. • The door release signal is detected during operation. • The door lock signal is not detected during operation.	

4. Error monitoring chart


4.1 FABRIC cycle

Display	Error Names	Detection	Standby	Auto water supply	Auto start	Vacuum air removal	Heating	Sterilization	Exhaust	Drainage	Dry	To normal pressure	Comp.
E 1	Overheat	CN4-10, CN8-1											
E 2	Disconnection of control temperature sensor wire	CN9											
E 3	Over-temperature	chamber											
	Over-temperature (Sterilization stage)	chamber											
E 4	Over-cool	chamber											
E 5	Over-pressure	chamber											
	Abnormality in the pressure sensor	chamber											
E 6	Abnormality in the door	CN8-3											
E 8	Abnormality in water supply and drainage	CN4-12, CN7-1											
E 9	Abnormality in the heater	chamber											
E P	Vacuum abnormality	chamber											
E L	Abnormality in the door lock	CN8-5,7											
E F	Disconnection of floating sensor wire	CN10											
※ 1	Over-temperature (Set temp. + 3.0°C)	chamber											
※ 1	Over-cool (Set temp. -0.1°C)	chamber											

4.2 LIQUID cycle

Display	Error Names	Detection	Standby	Auto water supply	Auto start	Vacuum air removal	Heating	Gravity displacement air removal	Heating	Sterilization	Exhaust	Drainage	Comp.
E 1	Overheat	CN4-10, CN6-1											
E 2	Disconnection of control temperature sensor wire	CN9											
E 3	Over-temperature	chamber											
	Over-temperature (Sterilization stage)	chamber											
E 4	Over-cool	chamber											
	Over-pressure	chamber											
E 5	Abnormality in pressure sensor	chamber											
E 6	Abnormality in the door	CN6-3											
E 8	Abnormality in water supply and drainage	CN4-12, CN7-1											
E 9	Abnormality in the heater	chamber											
E P	Vacuum abnormality	chamber											
E L	Abnormality in the door lock	CN6-5,7											
E F	Disconnection of floating sensor wire	CN10											
※ 1	Over-temperature (Set temp. + 3.0°C)	chamber											
※ 1	Over-cool (Set temp. -0.1°C)	chamber											

※1 : Over-temperature or overcool mark lights







- If the chamber get overcooled, the sterilization timer stops.
- When the overcooled state is recovered, the sterilization timer resumes operation for the preset time.

※ Others

- Monitoring of "E F" works only when the floating sensor is ON.

5. Operation

5.1 Turning on the power supply

Operation	Display	Action
Power key "ON"	<p>※A model name and program version are displayed for 2 seconds Display of a model name : HRG-112=G112, HRG-140=G140</p> 	
	<p>※If voltage reduction of a backup battery occurs, the following screen is displayed for 2 seconds.</p> 	①A buzzer sounds
	<p>※If the time-keeping of a clock IC stops, the following screen is displayed for 2 seconds.</p> 	①A buzzer sounds ②Initializing the time
(Standby state)		





※When the power supply is turned on by pressing the POWER key, the set value is initialized.

※When the power supply is turned on by pressing the POWER and FUNC keys, the set value and the cycle count (total operation number of the sterilizer) are initialized.

5.2 Changing a cycle

Operation	Display	Action
OFABRIC	<p>"ST-BY FABRIC cycle program 1"</p> <p>PROGM 1 = 2 = 3 =</p> <p>TEMP 121.0 °C TIME 00:20 h:min DRY 40 min PRESS MPa</p> <p>EXHT DRAIN </p> <p> </p>	
OSOLID	<p>"ST-BY SOLID cycle program 1"</p> <p>PROGM 1 = 2 = 3 =</p> <p>TEMP 121.0 °C TIME 00:20 h:min DRY 40 min PRESS MPa</p> <p>EXHT DRAIN </p> <p> </p>	
OLIQ	<p>"ST-BY LIQ cycle program 1"</p> <p>PROGM 1 = 2 = 3 =</p> <p>TEMP 121.0 °C TIME 00:20 h:min DRY min PRESS MPa</p> <p>EXHT DRAIN </p> <p> </p>	
OADD DRY	<p>"ST-BY ADD DRY cycle program 1"</p> <p>PROGM 1 = 2 = 3 =</p> <p>TEMP °C TIME h:min DRY 40 min PRESS MPa</p> <p>EXHT DRAIN </p> <p> </p>	


5.3 Changing a setting

Operation	Display	Action
<p>During standby state</p> 	<p>※The outline characters blink. ※The interval of blinking is 0.5 sec.</p> <p>"Change of the sterilization temperature"</p>  <p>Every time the UP/DOWN key is pressed, the value increases or decreases in increments of 1°C. Holding down the UP/DOWN key increases or decreases the value in increments of 10°C.</p>	
	<p>"Change of the sterilization time and drying time"</p> <p>The set value blinks as described in the change of the sterilization temperature.</p> 	<p>①The changed value is set</p>
<p>A key is not operated for 5 sec.</p>	<p>The display returns to the ST-BY state.</p>	


5.4 Confirming or setting the auto-start time

Operation	Display	Action
<p>Except ADD DRY cycle</p> <p>FUNC.</p>	<p>※The outline characters blink. ※The interval of blinking is 0.5 sec.</p> <p>"The set auto-start time is confirmed."</p>	
<p>SET</p>	<p>"The number of days to a start is set"</p> <p>Every time the UP/DOWN key is pressed the value increases or decreases in increments of 1 day.</p>	
<p>SET</p>	<p>"The hour is set."</p> <p>Every time the UP/DOWN key is pressed the value increases or decreases in increments of 1 hour. Holding down the key increases or decreases the value in increments of 10 hours.</p>	
<p>SET</p>	<p>"The minute is set."</p>	
<p>SET</p>	<p>The display returns to the ST-BY state.</p>	<p>①The auto-start time is set</p>
<p>A key is not operated for 15 seconds when setting the day, hour and minute</p>	<p>The display returns to the ST-BY state.</p>	<p>①The setting is cancelled.</p>
<p>A key is not operated for 5 seconds when confirming the auto-start time</p>	<p>The display returns to the ST-BY state.</p>	


5.5 Setting the air removal time

Operation	Display	Action
<p>LIQ cycle only</p> <p>FUNC.</p> <p>x 2</p> <p>▲ ▼</p>	<p>※The outline characters blink. ※The interval of blinking is 0.5 sec.</p>  <p>Every time the UP/DOWN key is pressed the value increases or decreases in increments of 1 minute.</p>	<p>①The changed value is set.</p>
<p>A key is not operated for 5 sec.</p>	<p>The display returns to the ST-BY state.</p>	



5.6 Setting the door lock temperature

Operation	Display	Action
<p>LIQ cycle only</p> <p>FUNC.</p> <p>x 3</p> <p>▲ ▼</p>	<p>※The outline characters blink. ※The interval of blinking is 0.5 sec.</p>  <p>Every time the UP/DOWN key is pressed the value increases or decreases in increments of 1°C. Holding down the key increases or decreases the value in increments of 10°C.</p>	<p>①The changed value is set.</p>
<p>A key is not operated for 5 sec.</p>	<p>The display returns to the ST-BY state.</p>	










5.7 Setting the floating sensor

Operation	Display	Action
<p>Except ADD DRY cycle</p> <p>FUNC.</p> <p>x 4</p> <p>▲ — ▼</p>	<p>※The outline characters blink. ※The interval of blinking is 0.5 sec.</p>  <p>Every time the UP/DOWN key is pressed, ON and OFF are switched</p>	<p>①The changed value is set</p>
<p>A key is not operated for 5 sec.</p>	<p>The display returns to the ST-BY state.</p>	

5.8 Confirmation and correction of the clock


Operation	Display	Action
<p>FUNC.</p> <p>X5 ※</p>	<p>"Confirming the clock time"</p>  <p>Every time the "UP/DOWN" switch is pressed the value increases or decreases in increments of 1 year. Holding down the switch increases or decreases the value in increments of 10 years.</p>	
<p>SET</p> <p>▲ — ▼</p>	<p>※The correction value display blinks every 0.5 sec.. "Correcting the year"</p>  <p>Every time the "UP/DOWN" switch is pressed the value increases or decreases in increments of 1 year. Holding down the switch increases or decreases the value in increments of 10 years.</p>	

- continued -


	<p>"Correcting the month"</p> 	
	<p>"Correcting the day"</p> 	
	<p>"Correcting the hour"</p> 	
	<p>"Correcting the minute"</p> 	
	<p>Display returns to the ST-BY state.</p>	<p>①The changed value is fixed.</p>
<p>Leaving it untouched for 15 seconds when correcting the time.</p>	<p>Display returns to the ST-BY state.</p>	<p>①The set reservation time is cancelled.</p>
<p>Leaving it untouched for 5 seconds when confirming the clock time.</p>	<p>Display returns to the ST-BY state.</p>	

※ Number of time pressed the "FUNC" switch will be changed according the selected cycle or option installed.


5.9 Confirming the floating sensor temperature

Operation	Display	Action
When process is in progress	<p>※The value is displayed when ON/OFF setting of the floating sensor is ON.</p> 	
	<p>Return to the display of the control sensor temperature.</p>	

5.10 Power-saving operation

Operation	Display	Action
During standby state or at the end of operation	<p>※The outline characters blink. ※The interval of blinking is 0.5 sec.</p> 	①Relay output is turned OFF
Autoclave is not operated for 10 minutes.		
Press any key	The display returns to the ST-BY or COMP.	①Return to the standby or completion stage

5.11 Error detection




Operation	Display	Action
	<p>※The outline characters blink. ※The interval of blinking is 0.5 sec.</p> 	<p>①Relay output is turned OFF ②A buzzer sounds</p>

5.12 Printer output

Operation	Printout
[Example] AGAR cycle	<p>※The printer prints when the printer setting is ON. ※The printer does not print during 104 °C or less of sterilized setting temperature.</p>
Starting the cycle	<div style="display: flex; align-items: flex-start;"> <div style="border: 1px solid black; padding: 10px; margin-right: 20px;"> <pre> HIRAYAMA HICLAVE HRG Version 1.0 Cycle Count = 0001 Date = 23 Mar 2011 Time = 15:30 Program = Fabric 1 Steri Temp = 121 °C Steri Time = 00:20 Dry Time = 40 min Tim C.Tmp F.Tmp Press S 000 24.5 24.5 0.000 H 003 22.3 24.4 -0.085 H 006 32.1 24.4 -0.080 H : 065 122.5 121.1 0.110 S 066 122.0 121.2 0.108 S 067 121.5 121.3 0.106 S : 087 100.5 100.5 0.002 E : 096 96.6 95.5 -0.058 D 099 93.1 95.2 -0.067 D : : Success </pre> </div> <div> <p>←Manufacturer, Model ←Control program version</p> <p>←Cycle count (total operating number of times) ←Date of operation startup ←Time of operation startup</p> <p>←Selected cycle and program ←Set sterilization temperature ←Set sterilization time ←Set dry time</p> </div> </div>
Heating stage	<p>←An elapsed time from a start, control sensor temperature, floating sensor temperature. ("OFF" is printed when a floating sensor is OFF), chamber pressure and the stages (H: Heating, S: Sterilization, E: Exhaust, D: Dry) are printed out every 3 minutes.</p>
Sterilization stage	<p>←An elapsed time from a start, control sensor temperature, floating sensor temperature, chamber pressure and the stages are printed out every 1 minute.</p>
Exhaust stage	<p>←An elapsed time from a start, control sensor temperature floating sensor temperature, chamber pressure and the stages are printed out every 3minutes.</p>
Dry stage	
Completion	<p>←Completion of the cycle is printed out and printing stops.</p>
When changing to the warming stage	<div style="display: flex; align-items: flex-start;"> <div style="border: 1px solid black; padding: 10px; margin-right: 20px;"> <pre> Success </pre> </div> <div> <p>←Completion of the cycle is printed out and printing stops.</p> </div> </div>
When an operation was suspended	<div style="display: flex; align-items: flex-start;"> <div style="border: 1px solid black; padding: 10px; margin-right: 20px;"> <pre> Canceled </pre> </div> <div> <p>←Printing stops.</p> </div> </div>
When an error is detected	<div style="display: flex; align-items: flex-start;"> <div style="border: 1px solid black; padding: 10px; margin-right: 20px;"> <pre> Error # </pre> </div> <div> <p>←Error number is printed and printing stops.</p> </div> </div>

6. Timing chart

●Explanation of notes on the timing chart.

Note Number	Explanation
①	The SSR is turned ON 0.5 sec. after the RY1 heater has been turned ON.
②	The heater is turned OFF 0.5 sec. after the SSR has been turned OFF.
③	The solenoid valve is turned ON when degree of vacuum reaches -0.06MPa.
※1	<p>Completion buzzer : The buzzer rings 3 times</p>  <p>ON OFF</p> <p>0.5 0.5 2 0.5 0.5 1 (sec.)</p>
※2	<p>ON/OFF time of the exhaust valve</p> <p>The exhaust valve is turned on or off as follows according to the exhaust level setting. (Exhaust level 1: the valve is always off. Exhaust level 3: the valve is always on.)</p> <p>•Level 1</p>  <p>•Level 2</p>  <p>ON OFF</p> <p>135 - 130°C 129 - 122°C 121 - 115°C 114 - 108°C 107 - 100°C 0.01 MPa</p> <p>0.1 0.1 0.2 0.2 0.3 0.3 0.5 0.5 1.0 (sec.) 1.0 (sec.)</p> <p>6.5 6.5 6.5 6.5 6.5 (sec.)</p> <p>135 - 130°C 129 - 122°C 121 - 115°C 114 - 108°C 107 - 100°C 0.01 MPa</p> <p>0.2 0.2 0.4 0.4 0.6 0.6 1.0 1.0 2.0 (sec.) 2.0 (sec.)</p> <p>2.0 2.0 2.0 2.0 2.0 (sec.)</p>
※3	When there is no water in the chamber, the solenoid valve is turned on for 20 seconds and turned off for 10 seconds. When there is water in the chamber, the solenoid valve is turned on for 5 seconds and turned off for 10 seconds.
※4	ON/OFF control of the heater is performed at 230°C that is input from the heater temperature sensor. When the floating sensor setting is on, temperature is controlled at the set temperature +1°C.
※5	The ON/OFF control of the heater is performed when temperature input from the heater temperature sensor is 230°C.
※6	The ON/OFF control of the heater is performed when temperature input from the heater temperature sensor is 160°C.
※7	When the vacuum motor valve is open and the degree of vacuum reaches -0.06MPa, the solenoid valve for cooling water supply is turned on.
※8	When pressure is 0.001MPa or more, the solenoid valve for air supply is turned off. When pressure is -0.001MPa or less, the solenoid valve for air supply is turned on.
※9	When pressure is 0.005MPa or more, the solenoid valve for exhaust is turned on. When pressure is 0.002MPa or less, the solenoid valve for exhaust is turned off.
※10	When the door is not locked, or, when water in the chamber is not discharged, the drain motor is turned on.

6.1 FABRIC Cycle

Vacuum air removal → Sterilization → Exhaust/Drainage → Dry

O: open, S: shut

Output	Stage	Standby	Water supply	Auto start	vacuum air removal	Heat	vacuum air removal	Heat	Sterilization	Exhaust	Drain	Dry	Air supply	Comp
SSR (CN3)						① ※5		① ※4	Temp control			① ※6		
Heater (RY1)							②			②			②	
Motor (drainage) (RY2)		※10	Q S								Q S			
Motor (vacuum) (RY3)			S		Q	S	Q	S						
Motor (door lock) (RY5)	Top position	lock												
Air pump (press) SV (air supply) (RY6)	※8													※8
Vacuum pump (CN4-1)						2 s		2 s						2 s
SV (water supply) (CN4-2)			※3											
SV (cleaning) (CN4-3)		30 s x 2												
SV (cooling water) (CN11-1)					③ 30 s		③ 30 s					30 s	※7	
SV (exhaust) (CN11-2)	※9													※9
Status Display LED	ST-BY (LED7)													
	RUN (LED8)													
	COMP (LED9)													
Stage Display LED	Air removal (LED1)					Blink								
	HEAT (LED2)						Blink							
	STERI (LED3)							Blink						
	EXHT (LED4)								Blink					
	DRAIN (LED5)									Blink				
	DRY (LED6)										Blink			
7-seg LED	Left (seg1-4)	Preset temp										Bar (-) display		
	Center (seg5-8)	Preset time										Remaining dry time		
	Right (seg9,10)	Preset drytime		time to a start					Steri. time					
Conditions		Start	Full water	Auto start time	-0.085MPa 2min. from -0.075MPa	-0.01MPa (same)			Sterilization temp.	Stop of timer	0.01MPa	Exhaust	Stop of dry timer	0MPa

ON at -0.05MPa
OFF at -0.07MPa

6.2 SOLID Cycle

Vacuum air removal → Sterilization → Exhaust/Drainage

O : open, S : shut

output	stage	Standby	Water supply	Auto start	vacuum air removal	Heat	vacuum air removal	Heat	Sterilization	Exhaust	Drain	Steam exhaust	Air supply	Comp
SSR (CN3)					① ※5	① ※5	① ※4	① Temp control						
Heater (RY1)						②	②							
Motor (drainage) (RY2)	※10	Q S									Q S	Q S		
Motor (vacuum) (RY3)		S			Q S	Q S	Q S							
Motor (door lock) (RY5)	Top position	lock												
Air pump (press) SV (air supply) (RY6)	※8													※8
Vacuum pump (CN4-1)					2 s	2 s						2 s		
SV (water supply) (CN4-2)			※3											
SV (cleaning) (CN4-3)		30 s x 2												
SV (cooling water) (CN11-1)					30 s	30 s					30 s			
SV (exhaust) (CN11-2)	※9													※9
LED Status Display	ST-BY (LED7)													
	RUN (LED8)													
	COMP (LED9)													
Stage Display LED	Air removal (LED1)				Blink									
	HEAT (LED2)					Blink								
	STERI (LED3)						Blink							
	EXHT (LED4)							Blink						
	DRAIN (LED5)								Blink					
	DRY (LED6)										Blink			
7-seg LED	Left (seg1-4)	Preset temp			Chamber temperature									
	Center (seg5-8)	Preset time			Chamber pressure									
	Right (seg9,10)			time to a start					Steri. time					
Conditions	Start	Full water	Auto start time	-0.085MPa 2min. from -0.075MPa	-0.01MPa (same)			Sterilization temp.	Stop of timer	0.01MPa	Drainage	5min. door lock temp.	0MPa	

6.3 LIQ Cycle

Gravity displacement air removal → Sterilization → Exhaust/Drainage

O: open, S: shut

output		Standby	Water supply	Auto start	vacuum air removal	Heat	Gravity displac.	Heat	Sterilization	Pulse exhaust	Drain	Steam exhaust	Comp
SSR (CN3)						①			Temp				
Heater (RY1)						※5		※4	control				
Motor (drainage) (RY2)		※10	Q S						②				
Motor (vacuum) (RY3)			S		Q	S					Q S		
Motor (door lock) (RY5)		Top position	lock										
Air pump (press) SV (air supply) (RY6)		※8											55°C
Vacuum pump (CN4-1)						2 s							※8
SV (water supply) (CN4-2)				※3									
SV (cleaning) (CN4-3)			30 s x 2										
SV (cooling water) (CN11-1)					③	30 s				※2			
SV (exhaust) (CN11-2)		※9								※2			55°C
Status Display	LED												
	ST-BY (LED7)												
	RUN (LED8)												
Stage Display	LED												
	Air removal (LED1)					Blink							
	HEAT (LED2)						Blink						
	STERI (LED3)							Blink					
	EXHT (LED4)								Blink				
	DRAIN (LED5)									Blink			
7-seg LED	LED												
	Left (seg1-4)	Preset temp											
	Center (seg5-8)	Preset time											
Right (seg9,10)													
Conditions		Start	Full water	Auto start time	-0.09MPa 2min. from -0.075MPa	103°C	Air removal ends	Sterilization temp.	Stop of steri. timer	0.01MPa	Drainage	Door lock temp.	

6.4 ADD DRY Cycle

Dry

O : open, S : shut

stage		Standby	Dry	O : open, S : shut	Air supply	Comp
output						
SSR (CN3)			① ※6			
Heater (RY1)				②		
Motor (drainage) (RY2) ※10			On at -0.05MPa OFF at -0.07MPa		
Motor (vacuum) (RY3)						
Motor (door lock) (RY5) Top position					
Air pump (press) SV (air supply) (RY6) ※8			 ※8	
Vacuum pump (CN4-1)				2 s		
SV (water supply) (CN4-2)						
SV (cleaning) (CN4-3)						
SV (cooling water) (CN11-1)			30 s ※7			
SV (exhaust) (CN11-2) ※9				 ※9
Status Display	LED					
	ST-BY (LED7)					
	RUN (LED8)					
Stage Display	LED					
	Air removal (LED1)					
	HEAT (LED2)					
	STERI (LED3)					
	EXHT (LED4)					
	DRAIN (LED5)					
	DRY (LED6)					
7-seg LED	LED					
	Left (seg1-4)					Temp
	Center (seg5-8)					Pressure

6.5 Vacuum Test Cycle

Vacuum air removal → Pressure holding

O : open, S : shut

stage		Standby	Vacuum air removal	Pressure holding	O : open, S : shut	
output					Air supply	Comp
SSR (CN3)						
Heater (RY1)						
Motor (drainage) (RY2)						
Motor (vacuum) (RY3)			Q	S		
Motor (door lock) (RY5)						
Air pump (press) SV (air supply) (RY6)						
Vacuum pump (CN4-1)				2 s		※8
SV (water supply) (CN4-2)						
SV (cleaning) (CN4-3)						
SV (cooling water) (CN11-1)			③	30 s		
SV (exhaust) (CN11-2)						※9
Status Display LED	ST-BY (LED7)					
	RUN (LED8)					
	COMP (LED9)					
Stage Display LED			All LEDs for stage display are off.			
7-seg LED	Left (seg1-4)		Chamber temperature			
	Center (seg5-8)		Chamber pressure			
	Right (seg9,10)		FO			
Conditions		Start Dip switch S1-6 is ON. The cycle does not begin when temperature is 40°C or more.				
		-0.085MPa 2min. from -0.075MPa				
		Stop of 20 min. timer				
		0MPa				

●Pouring water into the tank

Pressing the UP/DOWN (▲ , ▼) keys simultaneously during standby state pours water into the tank for 3 minutes.

●Inspecting a safety valve

Set the S1-7 of the dip switch to ON position and turn the power supply switch on.
Start the LIQ cycle and close the exhaust valve without performing air removal.

The safety valve does not operate. When the chamber pressure reaches 0.294MPa or more, an over-pressure exhaust is performed and an over-pressure alarm is activated.

●Function number

F0	Vacuum test (Set the S1-6 of the dip switch to ON)
F1	Programmable auto-start
F2	Air removal time
F3	Door lock temperature
F4	Setting of a floating sensor
F5	
F6	Time setting
F7	
F8	
F9	

1. Replacing the door gasket



- ① Open the door.
- ② Pull out the old gasket.
- ③ Using a cloth, wipe off the part where the old gasket was installed.
- ④ Install the new gasket in the chamber

- Insert a little the inner side perimeter of the new gasket into the outer side perimeter of the chamber in the order shown in the figure so that each length becomes almost same.
- Then, gradually and evenly insert the entire gasket while taking care so as not to damage the inner side perimeter of the new gasket.


In addition, touch the gasket with your finger along the inner side perimeter and the outer side perimeter and push it in to eliminate unevenness of the gasket surface.

- If there is unevenness on the gasket surface, it causes air leakage or makes noise when operating the door.

- ⑤ Set the sterilization temperature to the upper limit.
Start an operation and check the gasket for steam leak.




Insert the gasket evenly



2. Replacing the Solid State Relay (SSR)

- ① Pull out the Faston terminals connected to the solid state relay (SSR).
- ② Take out the screws and remove the SSR from the switchboard.
- ③ Wipe off the thermal grease and dust which adhere around the screw holes of the switchboard.
- ④ Clean the mounting surface of the new SSR, then apply thermal grease evenly on it.
- ⑤ Install the SSR on the switchboard and connect the Faston terminals.

- As the IN side of the SSR has polar direction, be sure to connect it as before.




3. Replacing the Backup Battery

● If E_b or E_c is displayed when turning on the main power supply, replace the backup battery as described below.

CAUTION


● Do not throw the used battery into the flames. It can burst.

- ① Slide aside the backup battery mounted on the control board and remove it from the holder.
• Discard the removed backup battery after insulating it with an insulation tape.
- ② Install the backup battery on the holder with the + side facing up
- ③ Referring to the operation manual, correct the clock



4. Replacing the Control Temperature Sensor

- ① Loosen the sensor retainer.
- ② Pull out the temperature sensor from the sensor retainer hole.
- ③ Insert the new sensor through the sensor retainer hole until the sensor tip comes to the same position as the internal surface of the chamber, and firmly tighten the sensor retainer using only fingers. Never use a tool such as an adjustable wrench.




5. Replacing the ROM

! IMPORTANT


- When replacing the ROM, use a specialized tool to avoid damage to the control board or the ROM.
- PCB or ROM can be damaged when it touches the static electricity collected in a person and a cloth. Before replacing the ROM, touch a metal or the like to discharge static electricity.

■ Required tools: ROM puller





- ① Pull the ROM out of the IC socket with a ROM puller as described below.
 - Press the trigger of the ROM puller to open the hooks.
 - Grip the bottom of the ROM with the hooks.
 - Pull the trigger to remove the ROM.
- ② Check if the pins of the new ROM are parallel. If the pins bend outward, correct them on the flat table so that they become parallel.



- ③ Insert the new ROM into the IC socket with care to the direction of ROM.
- ④ Turn on the main power supply switch while pressing the POWER key of the operation panel.
 - The old ROM data is deleted.




IC socket



5. Cautions when replacing the Control Board (Switchboard)

(1) When the control board (switchboard) is replaced, perform the following procedure before turning on the power supply.

- ① After confirming the model name, check if the optional floating sensor is installed, and set the S1 dip switch of the control board as follows.



With a floating sensor	Set S1-3 to ON (to the upper position)
------------------------	--


- The special specification program may be installed. If the version number of the ROM used by the old control board is 3-digit, the program specification needs to be identified because the program is a specially made specification.


② Turn on the main power supply switch while pressing the POWER key of the operation panel.


- The data of the new control board is initialized.



Chapter 3. Reference Drawings


● External view












- | | | | |
|-------|------------------------|------|---|
| CP | Operation panel | H1 | Header (chamber lower part) |
| CPB | Control board | H2 | Header (chamber upper part) |
| RP | Relay board | SSR | SSR (header) |
| ELB | Earth leakage breaker | 7X | Relay (vacuum pump) |
| LF | Noise filter | 8X | Relay (root water supply) |
| VR | Varistor | 9X | Relay (exhaust SV) |
| FUSE1 | Fuse (2A 250V T) | PS | Pressure sensor |
| FUSE2 | Fuse (2A 250V T) | S1 | Control temperature solenoid |
| CPR | Circuit protector | S2 | Floating sensor |
| TR | Transformer | S3 | Header temperature sensor |
| SR | Switching power supply | FSW1 | Flood switch (full water) |
| VP | Vacuum pump | FSW2 | Flood switch (low water) |
| AP | Air pump | EOO | Overheat prevention device for header |
| MV1 | Motor valve (drainage) | LSW1 | Limit switch (door close) |
| MV2 | Motor valve (vacuum) | LSW2 | Limit switch (door lock) |
| M | Motor (door lock) | LSW3 | Limit switch (door unlock) |
| C1~C3 | Condenser for motor | LSW4 | Limit switch (drain motor valve open) |
| SV1 | SV (Soil water supply) | LSW5 | Limit switch (drain motor valve close) |
| SV2 | SV (Exhaust) | LSW6 | Limit switch (vacuum/exhaust motor valve open) |
| SV3 | SV (Soil water supply) | LSW7 | Limit switch (vacuum/exhaust motor valve close) |
| SV4 | SV (air supply) | | |
| SV5 | SV (chamber cleaning) | | |
- SV : Solenoid Valve



Chapter 4. Operation Check Programs

● This chapter describes the way to perform the hardware check of the HRG-112 and HRG-140 autoclaves

1. Outline

C 1	Checking if the LEDs light	All LEDs on the operation panel board are lit and each LED lights up sequentially.
C 2	Checking the switch input	LEDs corresponding to the switches (SW1-16) of the operation panel board or corresponding to the S1 dip switch of the control board are lit.
C 3	Checking the contact input	LEDs corresponding to the input of CN4, CN6 and CN7 connectors light up.
C 4	Checking the relay output	A relay output is activated by pressing an operation switch.
C 5	Checking or adjusting the input of the CN9 control temperature sensor	The input value of the CN9 control temperature sensor is converted into temperature and the temperature is displayed.
C 6	Checking or adjusting the input of the CN10 floating sensor	The input value of the CN10 floating sensor is converted into temperature and the temperature is displayed.
C 7	Checking the input of the CN8 pressure sensor	The input value of the CN8 pressure sensor is converted into pressure and the pressure is displayed.
C 8	Checking the input of the CN4-10 heater temperature sensor	The input value of the CN4-10 heater temperature sensor is converted into temperature and the temperature is displayed.

2. Startup method

- (1) Set the S1-⑧ of the dip switch in the control PCB to ON.
- (2) Turn on the main power supply switch then press the POWER key.

3. Items to check beforehand

● Measure the voltage of the control board.

Measure the voltage between TP1 terminal (Circuit voltage) and TP2 terminal (GND): DC + 5V \pm 0.25V



Operation test

C1: Checking if LEDs light

- "I" is displayed on the 7SEG10 for two seconds. All LEDs on the operation panel board are lit for 4 seconds. Then, each LED is lit in sequence as shown below.

- When all LEDs are lit, the LED 18 and 19 turn to yellow after they light in red for two seconds.
- The 7SEG9 and ⑧ of 7SEG10 do not light up.

Lighting sequence of 7-segment LED





- Pressing the SW14 switch advances the operation test program to C2.

- "2" is displayed on the 7SEG10. A segment of LED which corresponds to the SW1-SW16 switches of the operation panel board or a segment of LED which corresponds to ON/OFF input of the S1 dip switch of the control board is lit.

- The switch input of SW14 is not confirmed with this program, because the SW14 operation switch is used to change the mode.
- When an operation switch is pressed a LED lights up, and when the operation switch is released the LED is turned off.
- When a dip switch is ON, a LED lights up, and when the dip switch is OFF, the LED does not light up.
- After checking the dip switch, restore the ON/OFF setting of the dip switch to a previous state.

Correlation between a switch and a lit segment of LED




- Pressing the SW14 switch advances the operation test program to C3.

C3: Checking the contact input

- "3" is displayed on the 7SEG10. A segment of LED which corresponds to the contact input of the CN4, CN6 or CN7 connectors mounted on the control board is lit.

• When the contact opens, the LED segment lights up, and when the contact closes, it does not light up.

Correlation between the contact input and a lit segment of LED



- Pressing the SW14 switch advances the operation test program to C4.

C4: Checking the relay output

- "4" is displayed on the 7SEG10. When an operation switch is pressed, the relay corresponding to the operation switch operates.

• When an operation switch is pressed a relay operates, and when the operation switch is released the relay stops.

Correlation between a switch and a relay output

SWITCH	OUTPUT CONNECTOR	RELAY	OUTPUT TERMINAL	ACTVATED PARTS
SW5 (Function)	CN3-2	SSR	-	H1,H2 (Heater)
SW6 (Setup)	-	RY1	-	H1,H2 (Heater)
SW7 (▼)	-	RY2	F2	MV1 (Drainage motor valve)
SW8 (▲)	-	RY3	F3	MV2 (Vacuum motor valve)
SW1 (FABRIC)	-	RY5	F5	M (Door lock motor)
SW2 (SOLID)	-	RY6	F6	AP,SV4 (Air pump, Air supply solenoid valve)
SW3 (LIQ)	CN4-2	7X	-	VP (Vacuum pump)
SW4 (ADD DRY)	CN4-4	8X	-	SV1 (Chamber water supply solenoid valve)
SW10 (◀)	CN4-6	9X	-	SV5 (Chamber cleaning solenoid valve)
SW11 (▶)	CN11-2	RL2	-	SV3 (Cooling water supply solenoid valve)
SW15 (Power supply)	CN11-3	RL3	-	SV2 (Exhaust solenoid valve)

- Pressing the SW14 switch advances the operation test program to C5.

C5: Checking or adjusting the input of the CN9 control temperature sensor

- "5" is displayed on the 7SEG10 and the temperature detected by the CN9 control temperature sensor is displayed on the 7SEG1-7SEG4.
 - When temperature is normal, adjust the input of the CN9 control temperature sensor with a VR1 trimmer potentiometer mounted on the control board.
 - When temperature is 121°C, adjust the input of the CN9 control temperature sensor with a VR2 trimmer potentiometer mounted on the control board.
- Pressing the SW14 switch advances the operation test program to C6.

C6: Checking or adjusting the input of the CN10 floating sensor

- "6" is displayed on the 7SEG10 and the temperature detected by the CN10 floating sensor is displayed on the 7SEG1-7SEG4.
 - When temperature is normal, adjust the input of the CN9 control temperature sensor with a VR3 trimmer potentiometer mounted on the control board.
 - When temperature is 121°C, adjust the input of the CN9 control temperature sensor with a VR4 trimmer potentiometer mounted on the control board.
- Pressing the SW14 switch advances the operation test program to C7.

C7: Checking the input of the CN8 pressure sensor

- "7" is displayed on the 7SEG10 and the pressure detected by the CN8 pressure sensor is displayed on the 7SEG5-7SEG8.
- Pressing the SW14 switch advances the operation test program to C8.

C8: Checking the input of the CN4-10 heater temperature sensor

- "8" is displayed on the 7SEG10 and the temperature detected by the CN4-10 heater temperature sensor is displayed on the 7SEG1-7SEG4.
- Pressing the SW14 switch returns the operation test program to C1.



HIRAYAMA MANUFACTURING CORPORATION

Overseas Sales Dept.
2-6-5, Toyono-cho, Kasukabe-shi,
Saitama 344-0014, Japan
Tel. +81-48-735-1242
Fax +81-48-733-2384

HMC HIRAYAMA

株式会社 平山製作所

〒344-0014 埼玉県春日部市豊野町2-6-5
<http://www.hirayama-hmc.co.jp/>

10.06.1TY