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STEAM STERILIZER

IGEA




È un prodotto

0434

Class IIa (93/42/CE)

- ① Istruzioni d'usopag.4
- Ⓜ User's instructionspage 10
- ⓕ Mode d'emploipag.16
- ⓓ Betriebsanweisungenseite 23
- ⓔ Instrucciones para el uso....pag 30


fig. 4


fig. 5

INDICE

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INDEX

Components legend	page 10
Warnings	page 10
General information.....	page 11
Installation	page 12
User's instructions	page 12
Utility functions	page 13
RESET and COUNTER keys (cycle interruption).....	page 13
Cleaning	page 13
Routine maintenance	page 13
First intervention guide	page 14
Guarantee	page 15
WEEE.....	page 15
Sterilizing and drying times.....	page 15
Alarms Table	page 37-38
Safety table.....	page 39
Wiring and hydraulic diagram.....	page 40-42
Technical particulars	page 44
Declaration of conformity.....	page 45

WARNINGS

- The use of this device is limited to qualified staff.
- The sterilizer must be used only with distilled water
- Before using the equipment for the first time, we recommend to check the cleaning of the sterilizing chamber and to effect a loadless sterilizing cycle.
- In case of power failure during the autoclave sterilizing phase when the power comes back the writing ALARM 12 - BLACK OUT appears on the display, the cycle is considered invalid.
Anyway and in any case each cycle is to be considered valid if at the end of the cycle the writing OK - FINISHED CYCLE appears on the display.
- In case of several running cycles the equipment can reach high external temperatures.
- In order not to contaminate the reservoir water, we recommend you not to let any object leaning on the machine and to check that the reservoir cap is closed.
- We recommend to use the autoclave in a clean area.

COMPONENTS LEGEND:

- A) Handwheel
- B) Reservoir lid
- C) Drain tap door
- D) Safety lever
- E) Feet nut
- F) Clutch to unload
- G) Main Switch with fuses
- H) Display
- I) Power green led
- J) Door locked green led
- K) Heating yellow led
- L) Cycle yellow led
- M) START - and ENTER key
- N) CYCLE and RESET key
- O) DRYING and COUNTER key
- P) UTILITY key
- Q) SHIFT key
- R) EMPTY red led
- S) OK green led
- T) FULL yellow led
- U) Manometer
- V) Green Main switch
- W) Serial outlet
- X) -----
- Y) Reset key

- Check that the material can be sterilized with the chosen sterilization cycle.
- In all cases of malfunction empty the chamber connecting the silicone pipe linked to the water drain tap to a container Then open the tap and the water will flow from the chamber into the external container.
- Before screwing out any part of the sterilizer make sure that this one is disconnected from the feeding network.
- This equipment has been designed and manufactured in compliance with the Community regulation 93/42 EEC "Medical devices" and C.B.M. S.r.l. is considered responsible concerning safety, reliability and performances on condition that possible modifications or repairs shall be effected by personnel authorized by C.B.M. , provided that the environment electric equipment is in compliance with IEC regulations and provided that the equipment is used according to these user's instructions.

IMPORTANT: The manufacturer company is relieved of any responsibility for possible damages to

people or things coming from not effecting the earth line connection.

USE SUGGESTION

- We recommend to use first quality disposable bags, in order to avoid little pieces of paper being torn off the envelope during sterilization and obstructing the different parts of the hydraulic system and the autoclave, compromising their working.
- We recommend to clean and dry the instruments to be sterilized before putting them into the autoclave. It is dangerous for the machine to use chloride containing detergents to clean instruments: it may quickly endanger the inner pieces of the autoclave.
- If during the water loading phase, "Alarm 2" appears, it may be a problem of electric conductivity of the water or of earthing of the electric system. If you use bi- or tridistilled water, that is having electric conductivity near zero, the level probe cannot feel its presence. For this reason we suggest to use only demineralized water. On the other hand, if the electric system the autoclave is connected to is provided with an insufficient earthing, the level probes cannot work correctly and they may not feel the presence of water. To avoid malfunctions it is necessary to use an only "Schuko" type socket for the autoclave, to test the earthing in the socket and to avoid the use of reductions and patch cords.
- A frequent cleaning of the exhaust filter placed on the bottom wall of the sterilization chamber is recommended. Before removing the filter itself from its seat, the surrounding area of the chamber is to be cleaned from possible residuum's, in order to avoid the dirt near the filter falling into the exhaust tube while removing the filter itself. Also the exhaust filter shall be manually unscrewed and cleaned with compressed water, if possible.
- Please remember that the alarm display or warning shall be reset through the RESET function, by temporarily pressing the SHIFT and RESET keys. If the machine is switched off when an alarm or a warning are being displayed, that is without resetting them before putting out the machine, on switching on the machine "Alarm 12 Black-Out" is displayed, in order to inform the user that the machine had been put out in "running phase". By resetting this message the machine goes back to its starting conditions, ready for a new cycle.
- The material not intended for immediate use must be packed and sterilised inside the package (envelope and/or rigid container for sterilising by steam).

GENERAL INFORMATION

Main switch

Set on the left-hand side of the equipment and also holds the fuses

During the stage of installing the equipment it is advisable to activate it and leave it switched on permanently. While the equipment is being operated, use only the green main switch

Green Main switch:

Placed near the basement, the switch, when it is pressed, lights up confirming the electric connection and the possible use of the equipment.

writing cycle: (fig.4) Placed on the display, it shows the position where the type of selected cycle is displayed (*time in min./temperature in °C degrees*). Table showing the available cycles at page 15

writing drying: (fig.4) Placed on the display, it shows the position where the selected drying program is displayed (*time in min.*). Table showing available times at page 15

Signalings leds

power green led: The led lights up when the voltage is present on the machine.

Door locked green led: It confirms the right closing of the door.

Heating yellow led: Signalling the working of the resistances.

Cycle yellow led: It shows that the autoclave is effecting a cycle.

RESERVOIR:

The level of the distilled water in the reservoir can be checked through the following luminous indicators:

EMPTY red led:

If it is lighted it shows that the reservoir is empty. For safety reasons when the EMPTY led is lighted the use of the machine is inhibited.

OK green led:

When it is lighted, it shows that the quantity of distilled water contained in the reservoir is enough for at least one sterilization cycle.

FULL yellow led:

It shows the reaching of the maximum water level in the reservoir, after its lighting immediately stop the immission of distilled water in the reservoir.

Functions keys:

START key:

Pressing the **START key** the cycle and the possible selected drying begin.

CYCLE key:

If it is pressed more than one time, it selects the cycle you wish to effect.

DRYING key:

If it is pressed more than one time, it selects the drying time you wish to effect.

UTILITY key:

When it is pressed it enables the access to the utility functions.

SHIFT key:

If it is pressed together with the **RESET** or **COUNTER** keys, it enables their working.

Manometer:

During the sterilization phase a pressure is created in the chamber.

In the Maintenance section controls to be effected are signalled.

Serial outlet:

The writing **printer** identifies the serial outlet SR2 through which it is possible to connect the sterilizer with a PC or with a printer in order to carry out the recording of the sterilization cycles.

Access to the reservoir

In the upper part of the equipment it is possible to enter the reservoir, for its filling or its possible topping up. We recommend to let the reservoir always closed during the use of the autoclave.

We deeply recommend the use of distilled water in order to avoid scalings in the machine that could compromise its functionality.

Access to the chamber

Access to the sterilization chamber is made up of a handwheel and a safety hook. In order to enter inside the chamber keep the safety hook raised, and screw up the handwheel anticlockwise.

Thanks to the safety hook, it is not possible to enter the sterilization chamber when the equipment is working. This in order to prevent the unintentional opening of the autoclave with the consequent risk of not effecting the sterilization or to damage things and/or people. We recommend never to force the lever in order to open the chamber, in this case wait for some seconds and try to open again.

Chamber drain cock

If you open the door, placed on the right part of the sterilizer, you approach the drain cock of the water from the chamber, useful in case of malfunction.

Clutch to unload the water tank

By opening the autoclave door the proper tube can be connected with the clutch for unloading the water tank.

Please remember that the other end of the tube shall be first put into a vessel to collect unloaded water.

INSTALLATION

Place the autoclave on a solid and strong plane. With the help of a level, check the setting on the plane of the sterilizer. If a regulation is necessary release the nuts placed under the autoclave through a 13 mm. key and unscrew or screw the feet until the optimal positioning is reached, then lock the nuts.

Electric system connection

Connect up the power supply cable to the apparatus and turn on the main switch

Connect the equipment to an electric plug equipped with an earthing, with a feeding voltage of 230Volt. and with a minimum capacity of 10A.

On the autoclave safety fuses of the kind 250V - 10A - 5x20 are installed.

USER'S INSTRUCTIONS

Even if it is equipped with multifunction push buttons the autoclave can be used in a very simple way through the following keys:

START key

CYCLE key

DRYING key

ATTENTION: we recommend to clean and dry carefully the material to sterilize.

For a simple and essential use of the autoclave operate in the following way:

- Press the main switch and check the water level.

'
In case the water is missing, remove the reservoir lid and directly pour the distilled water into the opening until the **FULL** yellow led lights up (maximum level) or at least the **OK** green led, then close again the reservoir.

f
Insert the material to sterilize into the chamber and close again, turning the handwheel clockwise. After the lightening of the **door locked** led effect again $\frac{1}{4}$ of a handwheel revolution.

"
Press the **CYCLE** key and select the cycle to effect.

IMPORTANT: Before selecting the cycle check the table at page 15 and verify the compatibility of the material to sterilize with the temperature and the kind of cycle.

...
Press the **DRYING** key in order to select the possible drying (see table at page 15).

†
After checking on the display the exact correspondence of the effected choices press **START** and verify the lighting of the **cycle** led.

The writing CYCLE N°x will appear on the display, and from this moment the autoclave shall begin the sterilization phases in a completely automatic mode.

‡
It is possible to follow all the sterilization phases reading them on the display.

~
At the end of the cycle the writing OK. FINISHED CYCLE will appear on the display with a sound signal.

If the drying has not been selected open the door and take out the sterilized material.

‰

If the drying has been selected the writing DRYING will appear; at this moment half-close the door 2cm., leaving it placed near the closing device and possibly stopped by a turn of screw.

At the end of the drying phase the writing OK FINISHED CYCLE will appear on the display; with a sound signal, then open the door and take out the sterilized material.

UTILITY FUNCTIONS

When the autoclave is in the waiting phase (never during the execution of a sterilizing cycle) it is possible to enter the utility functions through the **UTILITY** key. Pressing this key all the headings listed here below are visualized on the display:

MAINTENANCE:

It indicates maintenance operations to effect and their relating terms.

TEST LEDS:

It effects a functional inspection of the leds.

TEMPERATURE IN THE CHAMBER:

It displays the actual temperature inside the chamber.

LANGUAGE SELECTION:

It makes possible to choose the language of messages visualized on the display among the following languages: ITALIAN, ENGLISH, FRENCH, GERMAN, *SPANISH

(*) if existing in the version.

DATE AND HOUR:

It makes possible to modify the dater and the clock.

Instructions for use:

After selecting the desired heading with the **UTILITY** key operate in the following way:

MAINTENANCE:

If you wish to screen maintenance messages press **ENTER**. The 1st message will appear on the display; in order to screen the following messages press again the **UTILITY key**. If you wish to end the reading before its time press the key **SHIFT** and the **RESET** key at the same time.

TEST LEDS:

In order to check the actual working of the leds press the **ENTER** key. All the working leds will light up for 5 seconds.

TEMPERATURE IN THE CHAMBER:

In order to verify the temperature inside the chamber press

ENTER. Press again the same push button to go out of the function.

LANGUAGE SELECTION:

In case you wish to modify the language press **ENTER**. Choose the desired language by the

UTILITY key.

To confirm your choice press **ENTER**.

DATE AND HOUR:

In order to modify date and hour press the **ENTER** key **twice**, at this moment the **CYCLE** and **DRYING** keys will work respectively from **plus(+)** and **minus (-)**.

Confirm each single modified value through the **ENTER** key.

RESET AND COUNTER KEYS

If they are pressed together with the **SHIFT** key they take on the following functions:

RESET

1. It stops a cycle in progress, whenever any reason makes it necessary.
2. It set alarms to zero resetting the equipment.

COUNTER

1. It displays the total number of cycles executed by the sterilizer, at any moment

CLEANING

In order to do a good cleaning use neutral detergents only and rinse with plenty of water.

Gaskets and all plastic parts can be cleaned with some cotton soaked with alcohol.

Do not use sponges or abrasive materials that can damage the surface.

During the cleaning avoid the use of detergents, disinfectants or other liquids containing chlorides or trichloro-ethylene because they attack the stainless steel; in case of contact an immediate rinsing is recommended with the replacement of the water in the reservoir.

ROUTINE MAINTENANCE

⌘

Every 20 cycles or after two working weeks clean the water drain filter.

ATTENTION: Before unscrewing the filter, be sure you have removed any dirt in the chamber

This filter is placed inside the sterilization chamber, under the bottom grate in the middle of the resistance and can be simply removed screwing it out. It can be cleaned with compressed air or with running water. It is clean when, looking into the light, you can see its wire gauze completely free.

ATTENTION You must effect this operation every time that the manometer indicates a residual pressure in the chamber at the end of the "steam release" phase, when the drying does not produce a good result or it happens that when you open the door some water is coming out of the chamber.

•

Every 50 cycles, or after each working month, change the reservoir water.

In order to execute this operation:

By opening the autoclave door the proper tube can be connected with the clutch for unloading the water tank.

Please remember that the other end of the tube shall be first put into a vessel to collect unloaded water.

When the reservoir is empty clean its bottom and the sensors using some cotton and a pliers long enough.

Ž

Every 50 cycles, or each month, clean the sensor of the water level in the chamber.

Open the door and remove the bottom grate, raising and taking it off. The level sensor is the white/transparent small rod about 3 cm. long coming up vertically from the bottom, in the middle of the resistance.

The small rod and the head nut must be carefully cleaned using a wad of cotton wool soaked with alcohol, making sure not to let any cotton residue.

•

Every 50 cycles, or each month, check that during the sterilization phase the pressure indicated by the manometer is in the band 2~2.3 bars for the temperature 134°C, or in the band 1~1.3 bars for the temperature 121°C. In case this does not take place call the service centre.

•

Every 150 cycles or every 3 months, or everytime you think it is necessary, clean the door gasket. Use a wad of cotton soaked with alcohol to clean the gasket.

Every 300 cycles, or every 6 months check that the stop nut is not wearing out because of some defect.

The stop nut can wear out if the screw mouths are not normal and when you have to force more than necessary closing the door.

You can realize that the nut is wearing out if you find some brass dust in the threads of the closing screw. In order to avoid that the door is kept closed thanks to the safety lever only you must change the stop nut as soon as possible. Each month lubricate the nut and the screw. We recommend to use lubricants for rolling and plain bearings, with a thermic field of -20° ÷ +300°C.

We recommend, in any case, to replace the piece every three years.

Every 300 cycles or every 6 months or everytime you think it is necessary clean the sterilization chamber and the door. In order to clean the chamber and the door use a wad of cotton wool soaked with alcohol, making sure not to let any cotton residue.

“

Every 6 months clean the water load filter, at the same time of the water change in the reservoir, that is when the reservoir is empty. Remove the left panel and screw out with a 19 mm. key the brass nut placed on the pipe fitting under the reservoir. Take out the wire gauze that must be cleaned with compressed air or with running water, then replace it in its seat and retighten the nut.

”

Once a year, at least, it is necessary to verify the integrity of the safety -valve:

The safety-valve is placed inside the carter, right side of the sterilizing chamber.

Unscrew the higher knurled part until you can hear a click, then insert the nose of a screwdriver into the slit (the lower part of the valve), and raise the small piston.

After the execution of this operation retighten the higher part of the valve.

We recommend, in any case, to replace the safety -valve every three years.

FIRST INTERVENTION GUIDE

ATTENTION: Before calling for service assistance check the following guide and the enclosed table of alarms

FAULT: Plug inserted, green main switch on position I but all is inactive (including the green switch pilot light).

CAUSE:

- a) Current is not arriving.
- b) Fuses burnt out

REMEDY:

- a) Check to ensure that power is present
 - b) Check to ensure that the main switch is in posizione I
 - c) Replace the fuses inside the main switch
- Before starting to change the fuses the power cable must be detached from the plug.

If after turning on once more or during the first cycle they burn out again immediately, then technical service assistance must be called in order to rectify the fault causing the protective fuses to jump.

DEFECT: Only the green main switch is lighted, but no pilot lamp is lighted on the control panel.

CAUSE: The overheating cutoff is operating because of an excess temperature in the chamber.

SOLUTION: Take the outlet away and press the button placed in the hole of the left panel with the help of a screwdriver (see Fig. 3 - 5).

Press the button very hard in order to reset the maximum thermostat. If, during the following cycle, the maximum thermostat should interfere again, call the technical assistance in order to remove the trouble causing the intervention.

ATTENTION: In all cases of malfunctioning empty the chamber connecting the silicone pipe, linked to the water drain tap, with a container. Then open the tap and the water in the chamber will flow down in the external container.

ATTENTION: Before unscrewing any part of the sterilizer be sure that this one is disconnected from the feeding network

GUARANTEE:

C.B.M. guarantees its own autoclaves for one (1) year from the recording date of the first cycle, against possible building defects or damages to the components not due according to the irrevocable opinion of C.B.M., to an improper use of the autoclave from the operators.

Warning.

The guarantee does not cover parts subject to wear and repairs or replacements effected by personnel that is not authorized by C.B.M.

During the guarantee period the possible replacement of an autoclave will be decided by C.B.M. only in case it is impossible to effect a repair such to assure the correct running of the machine. Data, descriptions and illustrations are purely indicative.

C.B.M. S.r.l. reserves the right to carry out ameliorative modifications to its own products without any previous notice.

Disposal Of Electric And Electronic Refuse

Waste from Electrical and Electronic Equipment - Directive 2002/96/EC (WEEE)

The symbol (shown above) that appears on this product means that it should not be disposed of together with normal household refuse but should be treated separately.

Electric and electronic equipment may contain substances that are harmful to the environment and to the health of humans and should therefore be disposed of at special waste disposal sites.

There are harsh legal penalties in force against persons found disposing of such equipment in an irresponsible manner.

When the equipment is no longer of any use, it may be taken back to the outlet from which you intend to purchase a replacement.

STERILIZATION TIMES

	Temp.(°C)	Time-
Unpacked instruments-	134°C	5min.
Packed instruments/Glass -	134°C	20min.
Waste in stainless container	134°C	30min.
Porous material/ rubber not packed	121°C	20min.
Porous material/ rubber packed	121°C	30min.
Emergency rapid cycle for unpacked instruments	134°C	3min.

DRYING

Available drying times	0 - 10 - 15 - 20 - 25 - 30min.
-------------------------------	--------------------------------

ATTENTION: In order to have an optimal working of the autoclave follow the instructions below:

1) In case of sterilization of surgical instruments or metallic material do not exceed a quantity of 5 Kg.

2) In case of sterilization of textiles it is a good rule not to obstruct the bottom grate holes completely and not fill the chamber beyond 80% of the available volume.

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ALARMS TABLE

ATTENTION Before contacting your service centre check the following guide

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
wiring diagram

- (A) fuse 10A 5x20
- (B) main switch
- (C) overheating cut off
- (D) resistance 1500 W
- (E) resistance 100 W
- (F) transformer
- (G) unloading solenoid valve
- (H) loading solenoid valve
- (I) minimum level feeler
- (J) chamber level feeler
- (K) maximum level feeler
- (L) thermoresistance PT 100
- (M) door microswitch
- (N) lock microswitch
- (O) cards connection fiat
- (P) power card
- (Q) CPU card
- (R) RS232 fiat
- (S) line filter
- (T) variator
- (U) thermoact.
- (V) fan

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
hydraulikschaltplan

esquema hidraulico

- (A) water drain solenoid valve
- (B) chamber drain tap
- (C) water loading valve
- (D) safety valve
- (E) continuous bleed
- (F) air bised solenoid valve
- (G) manometer
- (I) water reservoir
- (L) sterilization chamber
- (M) chamber level feeler
- (N) minimum water level feeler
- (O) maximum water level feeler
- (P) condenser
- (Q) tank drain tap

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Technical particulars

The machines are manufactured according to the present quality and regulations standards.

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We: C.B.M. S.r.l., Via Castello 10, 26038 Torre de Picenardi (CR) Italy, in the person of its Sole Administrator and Legal Representative Mrs. Trusi Innocenta declare on our own responsibility, that the following electrical equipment:

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meet the provisions of following Directives which apply to them and its modifications:

- Directive CE 93/42 “Medical Devices”
- Directive CE 73/23 “Low Tension”
- Directive CE 89/336 “Electro-Magnetic Compatibility”

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