PRESOCLAVE 30 L 4001135
PRESOCLAVE 75 L 4001136

AUTOCLAVES ELÉCTRICOS PARA ESTERILIZACIÓN A VAPOR
ELECTRIC AUTOCLAVE FOR STERILIZATION

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GENERAL INFORMATION

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be attached to the equipment and it must be available for all users.*
- 4) *If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.*
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY APPARATUS TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*
- 7) *Do not use the apparatus with liquids which can give off vapours capable of making explosive mixtures.*

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The standard equipment consist of the following components:

Descripción / Description	Código / Code	
PRESOCLAVE	4001135	4001136
Cubre-resistencia / Heater protector	5601	5600
Portagomas / Carry Rubbers	34100	
Alargo / Extend		
Manual de instrucciones / Instruction manual	80049	

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ACCESSORIES

(They must be installed in the factory)

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*Safety device that prevents door being opened while the chamber is pressured. Code **4001137** (It must be fitted in the factory).

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* Temperature probe Pt100 for reading the internal temperature of the chamber. A terminal permits the Pt100 to be connected to a printer or temperature measuring device. Code **4001219**.

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*Baskets in perforated plate of AISI 316 stainless steel:

Código cestillo Basket code	Medidas (cm) Dimensions (cm)		Código Autoclave / Autoclave code	Capacidad cestillos/ No of baskets
	Ø	Altura / Depth		
1000495	28	20	4001135	2
1000496	38	28	4001136	2

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Voltage supply 115-230 V 50/60 Hz as indicated on the characteristics plate.

Código Code	Capacidad litros Capacity litres	Medidas útiles (cm) Chamber dimensions (cm)		Medidas exteiores (cm) Overall dimensions (cm)			Rango °C Range °C	Consumo (W) Consump. (W)	Peso (Kg) Weight (Kg)
		Ø	Alto Depth	Alto Height	Ancho Length	Fondo Width			
4001135	30	30	40	90	48	56	115/ 134	2100	55
4001136	75	40	60	105	57	66		3400	85

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EQUIPMENT DESCRIPTION

PRESOCLAVE are autoclaves with semi-automatic control and digital temperature reading.

EQUIPMENT:

** Triple safety system:*

Safety pressure switch:

Prevents the pressure surpassing the set maximum limit.

Safety valve:

Disconnects the heating element in the event of over pressure surpassing the set maximum limit.

Safety thermostat:

Disconnects power to the heating element in the event of over temperature or the lack of water in the reservoir.

Thermal protection:

** Chamber, lid and lock system in AISI 316 stainless steel.*

** Top and outer casing in AISI 304 stainless steel.*

** Automatic air exhaust at start of sterilization cycle.*

** Pressure and vacuum gauge up to 3kg/cm²*

FEATURES:

Sterilization temperature from 115°C up to 134°C.

Sterilization time from 1 minute up to 60 minutes.

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INSTALLATION

Place the autoclave near a suitable electric mains.

Immobilise the autoclave with the stops provided in the front part of the apparatus.

Always use decalcified water.

*Fix a hose in the water/steam exit nipple of the machine, and also fix the other end to a drain or deposit without **OBSTRUCTING THE FLOW**.*

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ATTENTION!!! The **PRESSURE** and **HIGH TEMPERATURE** of the steam exit require the necessary precautions to avoid **SERIOUS BURNS**.

NEVER PRESS THIS HOSE BY HAND.

Do not install the autoclave in places where inflammable liquids are stored or in zones with special protection.

Do not use the autoclave if it is not earthed.

*Before opening the cover **ALWAYS** check that there is **NOT PRESSURE**. The pressure gauge indicator installed in the control panel must read 0kg/cm².*

Be sure that the voltage supply is the same as the one indicated on the characteristics plate.

If you change the plug bear in mind the following:


**Blue cable: Neutral.
Brown cable: Phase.
Yellow/green cable: Earth.**

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STARTING UP:

1º Check that the STEAM valve (23) and the DRAINAGE valve are closed.

2º Fill the tank of the autoclave up to the brim level.

3º Place the material to be sterilized.

4º Close the cover.

5º Turn on the autoclave with the main switch (1).


6º Select the sterilization temperature with the temperature controller switch knob (4) bearing in mind the relation between temperature and pressure of the saturated steam. (See the following table).

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1. Pulsar ▲/▼ para fijar la temperatura deseada.

Temperature set selection:

1. Push ▲/▼ to preset the desired temperature.

<u>Temperatura (°C)</u>	<u>Presión (Kg./cm²)</u>	<u>Temperature (°C)</u>	<u>Pressure (Kg/cm²)</u>
115	0,6	115	0.6
121	1	121	1
126	1,33	126	1.33
130	1,65	130	1.65
134	2	134	2

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7º Select the sterilization time with the TIMER (6) control knob.

At the same time the timer is activated, the indicator lamp (8) will lit up. If the cover is not closed correctly, the indicator lamp (8) will lit off and the (9) will lit up.

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8º *The autoclave has an automatic air exit valve which stays open during heating up to 110-115°C, ejecting the air from the autoclave, and on reaching the previously mentioned temperature, it automatically closes.*

Once the chosen temperature is reached, the sterilisation phase indicator lamp (6) will lit up, starting the discount of the selected time. At the end, the sterilised cycle end indicator lamp (5) will lit up and the apparatus will stop.

All this time, the indicator lamp (8) will work intermitently.

9º *To stop the machine press the main switch (1).*

ATTENTION!!! DO NOT OPEN THE COVER OF THE AUTOCLAVE UNTIL THE PRESSURE GAUGE IS AT «0» AND THE STEAM VALVE (23) IS OPEN

10º *If SOLIDS (gauzes, gloves, cottonwool, instruments, etc...) have been sterilized, steam can be cleared quickly by opening the STEAM valve (23) and emptying the water from the tank with the DRAINAGE valve (28).*

11º *If LIQUIDS (liquid medicines, closed culture medium, bottles, etc...) have been sterilized, the pressure must return to the ambient temperature alone or in extreme case by VERY CAREFULLY opening the STEAM valve (23).*

12º *The RAPID DECOMPRESSION OF LIQUIDS provokes the breakage of flasks, spilling of liquids, etc..., with the resulting obstruction of tubes in the lower part of the autoclave.*

PRACTICAL RECOMMENDATIONS FOR A PERFECT STERILIZATION:

1. *The material to be sterilized has to be perfectly clean, free of any incrustations and impurities, so cleaning with a good detergent and distilled water and rinsing with abundant water is recommended.*
2. *Do not overload the tray, baskets or racks, and leave a space between them of 1 or 2 cms to facilitate the steam circulation and the drying process.*

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3. Sterilization of «UNWRAPPED INSTRUMENTS».

- 3.1 *Place the instruments directly onto water-repellent paper, avoiding any piling up of one upon another.*
- 3.2 *Avoid placing instruments of different metal type in the same tray.*

4. Sterilization of «INSTRUMENTS IN BAGS».

- 4.1 *Put one instrument per bag.*
- 4.2 *Place bags onto water-repellent paper, in vertical position avoiding any piling up of one upon another or the autoclave walls.*

5. Tubes.

- 5.1 *The tubes have to be clean and dried.*
- 5.2 *Place tubes, if it is possible, without obstructing the both ends.*
- 5.3 *Take care that the tubes do not touch the autoclave chamber and avoid folding them because the inner diameter can be obstructed.*

6. Receptacles.

- 6.1 *NEVER put the receptacles hermetically closed.*
- 6.2 *Place the receptacle face down to avoid the water residues.*

7. Liquids.

- 7.1 *Place the liquid to be sterilized in receptacles capable of supporting the sterilization temperature.*
- 7.2 *The receptacles will be filled to about 2/3 of their capacity. They will NEVER be closed hermetically, they will cover with cotton or some type of tap that makes the air easy outlet and it does not make pressure.*
- 7.3 *Please avoid using a receptacle with a narrow neck.*
- 7.4 *At the end of the sterilization cycle, leave the autoclave to cool slowly (until the pressure is 0 Kg/cm²).*

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SAFETY:**Safety Pressure Controller:**

When abnormal overpressure occurs in the sterilization chamber the safety pressure controller acts by stopping the heating and opening the drain electrovalve. The alarm indicator lamp (7) lights up simultaneously.

Safety valve:

If abnormal overpressure occurred the safety valve (24) would open and expel the steam. Said valve is gauged to initiate the opening between 2.2 kg/cm² and 2.5 kg/cm².

Safety Thermostat:

When abnormal increase of the temperature or a considerable decrease in the water level occurs (with the consequent overheating of the heater element), the safety thermostat (26) starts working blocking the heating of the machine, opening the drain electrovalve and the Alarm indicator lamp (7) lights up simultaneously.

To deblock the machine proceed as follows:

- 1º Stop the machine by pressing the switch (1) in OFF position.*
- 2º Release any pressure with valve (23).*
- 3º Open the cover.*
- 4º Fill the tank with water up to the working level.*
- 5º Reset the safety thermostat (26) by pressing its red reset button. The red indicator lamp ALARM (7) will switch off.*


GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damage caused by incorrect use or causes beyond the control of J.P. SELECTA, S.A.

Any manipulation of the apparatus by personnel not authorized by J.P. SELECTA, S.A. automatically cancels the guarantee.

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SPARE PARTS

To guarantee the safety of the equipment, all spare parts must be purchased from J.P.SELECTA, s.a.

Descripción / Description	4001135	4001136
Manómetro / Pressure vacuum gauge		16009
Regulador de temperatura / Temperature controller		16168
Electroválvula de purga / Purge electrovalve		16213
Válvula de vapor / Steam valve		34556
Válvula de seguridad / Safety valve		34200
Resistencia esterilización / Sterilization heater	39158	39001
Rueda / Wheel		40007
Termómetro / Thermometer		42051
Timer 0-60' / Timer 0-60'		37013
Termostato de seguridad / Safety thermostat		43031
Sonda Pt100 / Pt100 probe		43071

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MAINTENANCE

Before removing the cover disconnect the apparatus from the mains.

The manipulation of the internal electronic circuits of the apparatus by unauthorized personnel can cause irreparable damage. Take it to one of the J.P.SELECTA authorized technical services.

CLEANING:

For the cleaning of the different parts use the following products:

Cleaning of stainless steel: Alcohol

Cleaning of plastic: Alcohol with cotton duster.

Cleaning of pipes with lime: Solution of Hydrochloric acid 50% and water 50%.

OTHER RECOMMENDATIONS:

From time to time revise the functioning of the safety valve.

Annual revision of the whole function of the autoclave by qualified personnel is advisable.