

# LIFEPAK® 9P defibrillator/monitor/pacemaker



# Service Manual

Manual No. 805454-01 October 1993

Corporate Headquarters: 11811 Willows Road Northeast

Post Office Box 97006 Redmond, WA 98073–9706 USA
Telephone: 206/867–4000
Toll Free: 800/442–1142 or 800/426–8047
Telex: 990211 D PHYSIO RDMD
Telefax: 206/883–2988

Part No.	
Serial No.	

**About This Manual:** This Service Manual is intended for use by technical service personnel. It describes how to maintain, test, troubleshoot, and repair the LIFEPAK 9P defibrillator/monitor/pacemaker.

A separate publication, the *Operating Instructions*, is intended for use by physicians, clinicians, and emergency care personnel. It provides step-by-step instructions for all operating features of the LIFEPAK 9P defibrillator/monitor/pacemaker as well as operator-level testing and maintenance.

### **Trademarks**

PHYSIO-CONTROL<sup>®</sup>, LIFEPAK<sup>®</sup>, FAST-PATCH<sup>®</sup>, LIFE•PATCH<sup>®</sup>, QUIK-LOOK<sup>®</sup>, DERMA JEL<sup>®</sup>, SAS (SHOCK ADVISORY SYSTEM)™, CODE SUMMARY™, and PARTSLINE™ are trademarks of the Physio-Control Corporation, 11811 Willows Road N.E., PO Box 97006, Redmond, WA 98073–9706.

## Warranty

Refer to the product warranty statement included in the accessory kit shipped with the product. Duplicate copies may be obtained in the USA by calling the Physio-Control PARTSLINE at 1-800-442-1142. Outside the USA, contact your local Physio-Control sales or service office.

# **Table of Contents**

Safety Information	Introduction xii   Terms xii   Warnings xiii   Symbols xiv
General Information	Service Information
Description	Introduction       1-1         Physical Description       1-1         Functional Description       1-7         Theory of Operation       1-9         Power Supply PCB Assembly (803726)       1-9         Power-On       1-9         Regulation       1-9         Overvoltage Protection       1-9         Output Isolation       1-9         Battery       1-9         BATT CHRG Indicator       1-10

Power Conversion PCB Assembly (803724) 1-10
Power Control
Logic Power Monitor
Switching Power Supply1-12
External Control Input
Energy Storage Capacitor Charger
Low Energy Charge Rate
Capacitor Charge Control1-14
Energy Transfer 1-14
Energy Dump
Energy Monitor
Serial Data Interface1-16
Analog Voltage Evaluation
Charge LED Drive
Interconnect Pacing PCB Assembly (805472) 1-17
Interface
Pace Pulse Output Drive
Pace Pulse Wave Shaper1-19
Pace Pulse Amplitude Decoder 1-19
Echo Pulse Width Converter 1-20
Pacing Leads Off Detector
Arming Circuit
Pace Pulse Length Generator 1-22
Isolated Power Supply
Switching Regulator1-22
Power Supply Controller
LED Drive
Main PCB Assembly (805474)
Preamp
System Controller
Display/Recorder/Pacing Interface
Display PCB/CRT Assembly (803706) 1-39
Video Signal Amplifier
Frame Sweep Generator
Raster Scan Oscillator
Raster Voltage Generator
Raster Generator
CRT Rias Generator 1-42

Operation	2	
•		Introduction
		Controls and Indicators
		<b>Option Settings</b>
		Test Menu
		Setup Menu
		Info Menu
		Clock Setting Procedure
Testing	3	
		Introduction
		Performance Inspection Procedure (PIP) 3-1
		PIP-Purpose
		PIP-Scope and Applicability
		PIP-Definitions
		PIP-Requirements 3-2
		PIP-Required Equipment
		PIP-Test Equipment Verification
		PIP-Workstation
		PIP-Personnel
		PIP-General Instructions
		PIP-Physical Inspection
		PIP-AC Power
		PIP-Battery Charge
		PIP-Battery Power
		PIP-Setup Menu Settings
		PIP-Keypad Test
		PIP-Lead Select
		PIP-Recorder Speed/Baseline Centering
		PIP-ECG Display
		PIP-1mV CAL 3-9
		PIP-CRT Gain/Aspect Ratio
		PIP-Heart Rate 3-10
		PIP-Rate Alarms/Alarm-Activated Record
		PIP-Recorder Annotation/Real-Time Clock
		PIP-Test Load
		PIP-Sync/Sync Annotation
		PIP-Recorder/CRT Noise
		PIP-Defib Disarm and CHARGE LEDs Check
		PIP-Paddle Polarity 3-12
		PIP-Charge Time

PIP-Discharge Control
PIP-Energy Dump
PIP-Charge Reduction
PIP-Delivered Energy
PIP-CODE SUMMARY Report
PIP-Pacemaker Power-Up3-16
PIP-Pacing Rate
PIP-Demand Pacing
PIP-Non-Demand Pacing
PIP-Pacing Leads Off
PIP-Leakage Current
PIP-Pacemaker Leads Detect Current
<b>PIP Checklist</b>
Test and Calibration Procedure (TCP) $\ \dots \ 3\text{-}23$
$TCP-Purpose \\ \dots \\ 3-23$
TCP-Scope and Applicability $\hdots 3\text{-}23$
TCP-Definitions
$TCP-Requirements \dots 3-23$
TCP-Equipment Requirements
TCP-Test Equipment Verification 3-24
TCP-Workstation Requirements
TCP-Personnel
TCP-Voltage Checks 3-26
TCP-Low Battery Threshold Checks
Monitoring Check
Defibrillator Charge Check
TCP-Test Menu Access
TCP-Defibrillator Calibration
TCP-Test Load Calibration
TCP-Charge Time
TCP-Output Waveform
TCP-Display Calibration
TCP-Power Supply Calibration
TCP-Recorder Calibration
TCP-Recorder Frequency Response
TCP-CRT Frequency Response
TCP-ECG Output
TCP-Pacemaker Calibration
TCP-Pacemaker Current and Waveform Output 3-43
<b>TCP Checklist</b>
Troubleshooting Aids
Error Codes and SERVICE Indicator

# Service and Maintenance

<b>Introduction</b>
Battery Maintenance
Recharging Batteries
Battery Maintenance Guidelines 4-2
Guidelines for Replacement Batteries $\dots 4-2$
Battery Capacity Check
Recycling Batteries $\dots 4-4$
General Maintenance
Operation from AC or DC Power Source 4-4
Strip Chart Recorder 4-4
Loading Paper
Cleaning
Protecting Printouts
Fuse Replacement
Special Handling Procedures for Static Sensitive
Devices (SSDs)
Look for SSD Symbol
Use Static-Dissipative Mat
Wear Wrist Strap4-6
Transport and Store PCBs Properly 4-6
Keep Work Area Static-Free
Test Work Area Routinely
Disassembly Procedures
Battery Pack Removal
Front Panel Bezel and Keypad Assembly Removal 4-7
Strip Chart Recorder Removal
Case Separation
Power Conversion PCB Removal and Reassembly 4-10
Power Conversion PCB Removal 4-10
Power Conversion PCB Reassembly 4-11
Main PCB Removal
Display PCB/CRT Assembly Removal
Power Supply Removal
Interconnect Pacing PCB Removal 4-12
Inspection Techniques
Exterior Inspection
Interior Inspection
Tools and Materials for Cleaning and Repair 4-14
Cleaning Procedures 4-14
External Cleaning 4-14
Interior Cleaning
PCB Repair Precautions 4-16
Preparation for Storage or Shipping 4-1

······································
ex-1

# List of Figures

Description	1	Figure 1-1	LIFEPAK 9P defibrillator/monitor/pacemaker 1-2
		Figure 1-2	LIFEPAK 9P defibrillator/monitor/pacemaker Functional Block Diagram 1-8
		Figure 1-3	Power Conversion PCB Block Diagram 1-11
		Figure 1-4	Interconnect Pacing PCB Block Diagram 1-18
		Figure 1-5	Main PCB Block Diagram
		Figure 1-6	Preamp Timing
		Figure 1-7	Echo Period
		Figure 1-8	Control Line Timing
		Figure 1-9	Character Cell Locations
		Figure 1-10	Video Intensity 1-35
		Figure 1-11	Printhead Assembly
		Figure 1-12	Display PCB/CRT Assembly Block Diagram 1-40
		Figure 1-13	Display Scanning 1-40
Operation	2	Figure 2-1	Front Panel Controls and Indicators 2-2
operation.	_	Figure 2-2	Rear Panel Controls, Indicators, and Connectors . 2-6
		Figure 2-3	Test Menu
		Figure 2-4	Setup Menu
		Figure 2-5	Info Menu 2-9
		Figure 2-6	Clock Set Display
Testing	3	Figure 3-1	Power-On Display
	•	Figure 3-2	Test Menu
		Figure 3-3	Setup Menu
		Figure 3-4	Display with Flatline ECG Trace 3-9
		Figure 3-5	1mV CAL Pulse
		Figure 3-6	Example of CODE SUMMARY Printout 3-15
		Figure 3-7	Test Connections for Ground Resistance Measurement
		Figure 3-8	Test Setup for Pacemaker Leads Detect Current 3-19
		Figure 3-9	Test Connector J3
		Figure 3-10	Test Menu
		Figure 3-11	Output Waveform
		Figure 3-12	Display PCB Assembly Potentiometer Locations . 3-34
		Figure 3-13	Power Supply Calibration 3-37
		Figure 3-14	Recorder Calibration Adjustments Location 3-38
		Figure 3-15	Test Signal Setup for Function Generator 3-40
		Figure 3-16	Pacemaker Calibration Test Load Setup 3-43
		Figure 3-17	Pacing Pulse Waveform

Service and	4	Figure 4-1	Front Panel Removal 4-8
Maintenance		Figure 4-2	Connector Locations
		Figure 4-3	Connections for Discharging Capacitors on the Power Conversion PCB Assembly 4-11
		Figure 4-4	LIFEPAK 9P defibrillator/monitor/pacemaker Shipping Assembly 4-17
			Simpping Losembly 4 17
Parts Lists/	5	Figure 5-1	LIFEPAK 9P defibrillator/monitor/pacemaker
Schematics			Final Assembly 5-11
		Figure 5-2	Interconnect Diagram
		Figure 5-3	Main PCB Assembly 5-22
		Figure 5-4	Power Supply PCB Assembly 5-31
		Figure 5-5	Power Conversion PCB Assembly 5-40
		Figure 5-6	Interconnect Pacing PCB Assembly 5-50
		Figure 5-7	Keypad Assembly 5-55
		Figure 5-8	Rear Panel Membrane Switch 5-56
		Figure 5-9	Strip Chart Recorder Assembly 5-58
		Figure 5-10	Transfer Relay Assembly 5-60
		Figure 5-11	Paddle Assembly 5-61
		Figure 5-12	Display PCB/CRT Assembly 5-64
		Figure 5-13	Power Supply PCB/AC Receptacle/Ground Wire Harness
		Figure 5-14	Power Conversion PCB/Power Supply PCB/Battery Wire Harness 5-74
		Figure 5-15	Battery Sensor/Power Supply PCB Wire Harness . 5-75
		Figure 5-16	High Voltage Connector Cable Assembly 5-76
		Figure 5-17	Main PCB/Display Cable Assembly 5-77
		Figure 5-18	Main PCB/Recorder Cable Assembly 5-78
		Figure 5-19	Interconnect Pacing PCB/Keypad Cable Assembly
		Figure 5-20	Main PCB/Patient Connector Wire Harness 5-80
		Figure 5-21	AC Receptacle/Ground Wire Harness 5-81
		Figure 5-22	AC Receptacle Wire Harness 5-81
		Figure 5-23	Test Load Contact/Resistor Wire Harness 5-82
		Figure 5-24	Power Conversion PCB/Reed Assy Wire Harness 5-82
		Figure 5-25	Power Conversion PCB/Dump Relay Wire Harness
		Figure 5-26	Power Conversion PCB/Charge Relay Negative Wire Harness
		Figure 5-27	Power Conversion PCB/Charge Relay Positive Wire Harness
		Figure 5-28	Main PCB/Sternum Relay Wire Harness 5-85
		Figure 5-29	Main PCB/Apex Relay Wire Harness 5-85
		Figure 5-30	Storage Capacitor Positive/Relay Wire Harness 5-86
		Figure 5-31	Storage Capacitor Negative/Relay Wire Harness . 5-86
		Figure 5-32	ECG Out Wire Harness

		Figure 5-33 Figure 5-34 Figure 5-35 Figure 5-36 Figure 5-37 Figure 5-38 Figure 5-39 Figure 5-40 Figure 5-41	Main PCB/Speaker Wire Harness5-88Power Supply Bracket/Ground Wire Harness5-89Connector Receptacle/Pacing Cable Assembly5-90Recorder/Bracket Ground Wire Harness5-91Defibrillation Adapter5-93Test Load High Voltage Wire Harness (W1)5-95Test Load High Voltage Wire Harness (W2)5-95Pushbutton Switch Wire Harness5-96High Voltage Connector Cable Assembly5-97
Component	6	Figure 6–1	Buffer
Reference		Figure 6–2	A/D Converter
Diagrams		Figure 6–3	D/A Converter
		Figure 6–4	Dual D/A Converter
		Figure 6–5	Counter
		Figure 6–6	Counter
		Figure 6–7	Font EPROM 6–6
		Figure 6–8	D Flip–Flop 6–6
		Figure 6–9	D Flip–Flop 6–7
		Figure 6–10	D Flip–Flop 6–7
		Figure 6–11	D Flip–Flop
		Figure 6–12	J–K Flip–Flop 6–8
		Figure 6–13	Microprocessor
		Figure 6–14	Multiplexer
		Figure 6–15	Multivibrator 6–10
		Figure 6–16	Multivibrator
		Figure 6–17	Oscillator
		Figure 6–18	Pulse–Width Modulator
		Figure 6–19	Pulse–Width Modulator
		Figure 6–20	Display RAM
		Figure 6–21	Video RAM
		Figure 6–22	Shift Register
		Figure 6–23	Analog Switch

## **List of Tables**

Description	1	Table 1-1	LIFEPAK 9P defibrillator/monitor/pacemaker Specifications 1-3
		Table 1-2	Main PCB Functional Circuits 1-23
		Table 1-3	Preamp Pulse Widths 1-27
		Table 1-4	Display/Recorder/Pacing Microprocessor I/O Lines
Operation	2	Table 2-1	Front Panel Controls and Indicators 2-3
		Table 2-2	Rear Panel Controls, Indicators and Connectors 2-6
Testing	3	Table 3-1	PIP Test Equipment 3-2
		Table 3-2	Delivered Energy Tolerances 3-14
		Table 3-3	Maximum Leakage Current for 120Vac and 60Hz 3-20
		Table 3-4	Maximum Leakage Current for 230Vac and 50Hz 3-20
		Table 3-5	Test and Calibration Test Equipment 3-24
		Table 3-6	Power Conversion PCB Voltage Checks 3-26
		Table 3-7	Low Battery Thresholds
		Table 3-8	Defibrillator Calibration Measurement Ranges 3-31
		Table 3-9	Current and Waveform Levels 3-44
		Table 3-9	Error Codes
Service and	4	Table 4-1	Inspection Techniques 4-13
Maintenance		Table 4-2	Tools and Materials for Cleaning and Repair 4-14
Parts Lists/	5	Table 5-1	LIFEPAK 9P defibrillator/monitor/pacemaker
Schematics			Major Assemblies 5-3
		Table 5-2	Defibrillation Adapter Assemblies 5-4
		Table 5-3	Reference Designator Key 5-4
		Table 5-4	Supplies and Accessories 5-5
Component Reference	6	Table 6-1	Component Reference Diagrams 6-1

# Safety Information

#### Introduction

This safety information includes terms and symbols used in this manual or on the equipment to alert both operating and service personnel of recommended precautions in care, use and handling of this specialized medical equipment.

Refer to NFPA (National Fire Protection Association) 99-1990, *Health Care Facilities*, and NFPA 70-1990, *National Electrical Code*, for specific guidelines on the standards and practices for health-care instruments and environments.

#### **Terms**

Certain terms are used in this manual, or on the equipment. Familiarize yourself with their definitions and significance.

Danger: Immediate hazards which will result in severe personal injury

or death.

Warning: Hazards or unsafe practices which could result in severe

personal injury or death.

Caution: Hazards or unsafe practices which could result in minor

personal injury or product damage.

**Note:** Points of particular interest for more efficient or convenient

instrument operation. Additional information or explanation

concerning the subject under discussion.

## **Warnings**

Following are descriptions of general hazards and unsafe practices that could result in death, severe injury or product damage. Specific warnings and cautions not appearing in this section are found throughout the manual.

#### **Shock Hazard**

When charged and discharged, the LIFEPAK 9P defibrillator/monitor/pacemaker discharges up to 360J of electrical energy through the defibrillator paddles. Unless discharged properly as described in this manual, this electrical energy may cause personal injury or death. Do not attempt to perform this procedure unless you are thoroughly familiar with the operation of the LIFEPAK 9P defibrillator/monitor/pacemaker and all accessories.

# Possible Arcing and Paddle Damage

When discharging the defibrillator into the internal test load, make sure the standard paddles are securely and properly stored in the paddle storage area (STERNUM paddle on the left, APEX paddle on the right). This helps prevent arcing and formation of pits on paddle electrode surfaces. Pitted or damaged paddle electrode surfaces can cause patient skin burns during defibrillation.

### Possible Fire or Explosion

Use care when operating this device close to oxygen sources (such as bag-valve-mask devices or ventilator tubing) and flammable gases and anesthetics.

#### Safety Hazard

Do not mount the product directly above the patient. Place the product in a location where it cannot harm the patient if it should fall from its shelf or other mounting location.

# Electrical Hazard and Possible Equipment Failure

Do not operate equipment using damaged cables and wires. Broken or frayed wires, or loose snap fittings may cause interference or loss of signal and possible equipment failure. Perform frequent electrical and visual inspections on cables and wires. Pay particular attention to the point at which the wires enter the terminals. Repeated flexing at these points eventually causes the wire strands to break.

#### **Shock or Fire Hazard**

Do not immerse any portion of the instrument in water or other fluids. Avoid spilling any fluids on the instrument or accessories.

# Possible Equipment Damage

Do not sterilize this product. Sterilization environments can cause serious damage. Do not autoclave or gas sterilize accessories unless manufacturer instructions clearly approve it.

#### Safety Risk

Do not substitute accessories. Use only recommended accessories listed in Table 5-4 shown on page 5-5. Substitution may cause the instrument to work improperly.

#### **Shock or Fire Hazard**

Make sure that all equipment is interconnected safely in accordance with NFPA (National Fire Protection Association) 70-1990, National Electrical Code.

Within certain governmental jurisdictions, all interconnected accessory equipment must be labeled by an approved testing laboratory. It is important that you verify and observe the required applications in your location. Check leakage current and grounding requirements after interconnecting this instrument with accessory equipment.

## **Symbols**

Any or all of the following symbols may appear in this manual or on the equipment:



Static Sensitive Device (SSD)



Additional information in Section 6, Component References



DANGER — high voltage present



ATTENTION — refer to manual for more information



Defibrillation protected, type CF patient connection



Defibrillation protected, type BF patient connection



Protective ground (earth) terminal



Fusible Link



Equipotentiality connector



Off (Power: disconnection from the mains)

On (Power: connection to the AC mains)



**ECG Output** 

# General Information

# Service Information

Before attempting to clean or repair any assembly in this instrument, the technician should be familiar with the information provided in Section 4, Service and Maintenance.

If the instrument has been dropped, damaged, or abused, a qualified technician should complete the **Performance Inspection Procedure (PIP)** in Section 3, beginning on page 3–1, to confirm whether the instrument is operating within specifications. If calibration or more extensive testing is required, a qualified technician can perform the **Test and Calibration Procedure (TCP)** in Section 3. Component replacement and internal adjustments must be performed only by service personnel qualified by appropriate training and experience.

Use of non-Physio-Control defibrillation electrodes, batteries, accessories, or adapter devices may void Safety Agency Certifications and warranty.

For assistance in servicing the instrument in the US, call Physio-Control at 1-800-442-1142. In other countries contact the local Physio-Control representative.

## Effective Publication Dates

The effective publication date for each page of this manual is listed below.

Title	Page	Date
Trademark and Warranty	ii	October 1993
Table of Contents	iii thru viii	October 1993
List of Figures	ix thru xi	October 1993
List of Tables	xii	October 1993
Safety Information	xiii thru xv	October 1993
Terms	xiii	October 1993
Warnings	xiv	October 1993
Symbols	xv	October 1993
General Information	xvi thru xix	October 1993
Service	xvi	October 1993
Effective Publication Dates	xvii	October 1993
Configuration Information	xviii	October 1993
Content Overview	xix	October 1993
1 Description	1-1 thru 1-42	October 1993
2 Operation	2-1 thru 2-10	October 1993
3 Testing	3-1 thru 3-47	October 1993
4 Service and Maintenance	4-1 thru 4-17	October 1993
5 Part Lists/Schematics	5-1 thru 5-97	October 1993
6 Component Reference	6-1 thru 6-14	October 1993
Index	Index-1 thru Index-7	October 1993

# Configuration Information

This manual is current with the listed revision level of the following part numbers. The assemblies appear in the same order as in Table 5-1.

LIFEPAK 9P defibrillator/monitor/pacemaker	Assembly Name	Part Number	Rev
Power Supply PCB Assembly         803726-01         Orig7           Power Conversion PCB Assembly         803724-03         D6           Interconnect Pacing PCB Assembly         805472-02         B4           Keypad Assembly         805469         B3           Rear Panel Switch Membrane         805464-00         Orig3           Strip Chart Recorder Assembly         804189-00         B2           Transfer Relay Assembly         800240-14         U10           Paddle Assembly         802901-06         F           Display PCB/CRT Assembly         803706         A5           Power Supply PCB/Ac Receptacle/Ground Wire Harness (W1)         803783-30         B5           Power Conversion PCB/Power Supply PCB Water Harness (W2)         803783-23         B5           Battery Sensor/Power Supply PCB Wire Harness (W3)         803783-24         B5           High Voltage Connector Cable Assembly (W4)         803756-13         B5           High Voltage Connector Cable Assembly (W5)         803733-00         Orig1           Main PCB/Monitor PCB Cable Assembly (W6)         803773-0         Orig1           Interconnect Pacing PCB/Keypad Cable Assembly (W7)         805465-00         Orig3           Main PCB/Recorder Cable Assembly (W6)         8037783-00         B5           <	LIFEPAK 9P defibrillator/monitor/pacemaker	805460	E4
Power Conversion PCB Assembly         803724-03         D6           Interconnect Pacing PCB Assembly         805472-02         B4           Keypad Assembly         805469         B3           Rear Panel Switch Membrane         805464-00         Orig3           Strip Chart Recorder Assembly         804189-00         B2           Transfer Relay Assembly         800240-14         U10           Paddle Assembly         802901-06         F           Display PCB/CRT Assembly         803706         A5           Power Supply PCB/Ac Receptacle/Ground Wire Harness (W1)         803783-30         B5           Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2)         803783-23         B5           Battery Sensor/Power Supply PCB Wire Harness (W3)         803783-23         B5           High Voltage Connector Cable Assembly (W4)         803766-13         D4           Main PCB/Monitor PCB Cable Assembly (W5)         803733-00         Orig1           Main PCB/Recorder Cable Assembly (W6)         803773-00         Orig1           Interconnect Pacing PCB/Keypad Cable Assembly (W7)         805465-00         Orig3           Main PCB/Patient Connector Wire Harness (W8)         803783-00         B5           AC Receptacle/Ground Wire Harness (W1)         803783-0         B5	Main PCB Assembly	805474-01	B8
Interconnect Pacing PCB Assembly         805472-02         B4           Keypad Assembly         805469         B3           Rear Panel Switch Membrane         805464-00         Orig3           Strip Chart Recorder Assembly         804189-00         B2           Transfer Relay Assembly         800240-14         U10           Paddle Assembly         802901-06         F           Display PCB/CRT Assembly         803706         A5           Power Supply PCB/AC Receptacle/Ground Wire Harness (W1)         803783-30         B5           Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2)         803783-23         B5           Battery Sensor/Power Supply PCB Wire Harness (W3)         803783-24         B5           Bigh Voltage Connector Cable Assembly (W4)         803756-13         D4           Main PCB/Monitor PCB Cable Assembly (W5)         803733-00         Orig1           Main PCB/Recorder Cable Assembly (W6)         803773-00         Orig1           Interconnect Pacing PCB/Keypad Cable Assembly (W7)         805465-00         Orig3           Main PCB/Patient Connector Wire Harness (W8)         803783-06         B5           AC Receptacle/Ground Wire Harness (W10)         803783-05         B5           Test Load Contact/Resistor Wire Harness (W11)         803783-05         B	Power Supply PCB Assembly	803726-01	Orig7
Keypad Assembly         805469         B3           Rear Panel Switch Membrane         805464-00         Orig3           Strip Chart Recorder Assembly         804189-00         B2           Transfer Relay Assembly         800240-14         U10           Paddle Assembly         802901-06         F           Display PCB/CRT Assembly         803706         A5           Power Supply PCB/AC Receptacle/Ground Wire Harness (W1)         803783-30         B5           Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2)         803783-23         B5           Battery Sensor/Power Supply PCB Wire Harness (W3)         803783-24         B5           High Voltage Connector Cable Assembly (W4)         803756-13         D4           Main PCB/Monitor PCB Cable Assembly (W5)         803733-00         Orig1           Main PCB/Recorder Cable Assembly (W6)         803773-00         Orig1           Interconnect Pacing PCB/Keypad Cable Assembly (W7)         805465-00         Orig3           Main PCB/Patient Connector Wire Harness (W8)         803783-06         B5           AC Receptacle/Ground Wire Harness (W10)         803783-05         B5           Test Load Contact/Resistor Wire Harness (W11)         803783-05         B5           Power Conversion PCB/Charge Relay Negative Harness (W12)         8037	Power Conversion PCB Assembly	803724-03	D6
Rear Panel Switch Membrane         805464-00         Orig3           Strip Chart Recorder Assembly         804189-00         B2           Transfer Relay Assembly         800240-14         U10           Paddle Assembly         802901-06         F           Display PCB/CRT Assembly         803706         A5           Power Supply PCB/AC Receptacle/Ground Wire Harness (W1)         803783-30         B5           Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2)         803783-23         B5           Battery Sensor/Power Supply PCB Wire Harness (W3)         803783-24         B5           High Voltage Connector Cable Assembly (W4)         803756-13         D4           Main PCB/Monitor PCB Cable Assembly (W5)         803733-00         Orig1           Main PCB/Recorder Cable Assembly (W6)         803773-00         Orig1           Interconnect Pacing PCB/Keypad Cable Assembly (W7)         805465-00         Orig3           Main PCB/Patient Connector Wire Harness (W8)         803783-0         B5           AC Receptacle/Ground Wire Harness (W9)         803783-0         B5           AC Receptacle Wire Harness (W10)         803783-0         B5           Test Load Contact/Resistor Wire Harness (W11)         803783-0         B5           Power Conversion PCB/Reed Assembly Wire Harness (W12)	Interconnect Pacing PCB Assembly	805472-02	B4
Strip Chart Recorder Assembly         804189-00         B2           Transfer Relay Assembly         800240-14         U10           Paddle Assembly         8022901-06         F           Display PCB/CRT Assembly         803706         A5           Power Supply PCB/AC Receptacle/Ground Wire Harness (W1)         803783-30         B5           Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2)         803783-23         B5           Battery Sensor/Power Supply PCB Wire Harness (W3)         803783-24         B5           High Voltage Connector Cable Assembly (W4)         803756-13         D4           Main PCB/Monitor PCB Cable Assembly (W5)         803733-00         Orig1           Main PCB/Recorder Cable Assembly (W6)         803773-00         Orig1           Interconnect Pacing PCB/Keypad Cable Assembly (W7)         805465-00         Orig3           Main PCB/Patient Connector Wire Harness (W8)         803783-0         B5           AC Receptacle/Ground Wire Harness (W9)         803783-0         B5           AC Receptacle Wire Harness (W10)         803783-0         B5           Test Load Contact/Resistor Wire Harness (W11)         803783-0         B5           Power Conversion PCB/Reed Assembly Wire Harness (W12)         803783-0         B5           Power Conversion PCB/Charge Relay Wegat	Keypad Assembly	805469	В3
Transfer Relay Assembly         800240-14         U10           Paddle Assembly         802901-06         F           Display PCB/CRT Assembly         803706         A5           Power Supply PCB/AC Receptacle/Ground Wire Harness (W1)         803783-30         B5           Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2)         803783-23         B5           Battery Sensor/Power Supply PCB Wire Harness (W3)         803783-24         B5           High Voltage Connector Cable Assembly (W4)         803756-13         D4           Main PCB/Monitor PCB Cable Assembly (W5)         803733-00         Orig1           Main PCB/Recorder Cable Assembly (W6)         803773-00         Orig1           Interconnect Pacing PCB/Keypad Cable Assembly (W7)         805465-00         Orig3           Main PCB/Patient Connector Wire Harness (W8)         803783-00         B5           AC Receptacle/Ground Wire Harness (W9)         803783-04         B5           AC Receptacle Wire Harness (W10)         803783-05         B5           Test Load Contact/Resistor Wire Harness (W11)         803783-06         B5           Power Conversion PCB/Reed Assembly Wire Harness (W12)         803783-07         B5           Power Conversion PCB/Charge Relay Negative Wire Harness (W13)         803783-08         B5           Pow	Rear Panel Switch Membrane	805464-00	Orig3
Paddle Assembly         802901-06         F           Display PCB/CRT Assembly         803706         A5           Power Supply PCB/AC Receptacle/Ground Wire Harness (W1)         803783-30         B5           Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2)         803783-23         B5           Battery Sensor/Power Supply PCB Wire Harness (W3)         803783-24         B5           High Voltage Connector Cable Assembly (W4)         803756-13         D4           Main PCB/Monitor PCB Cable Assembly (W5)         803733-00         Orig1           Main PCB/Recorder Cable Assembly (W6)         803773-00         Orig1           Interconnect Pacing PCB/Keypad Cable Assembly (W7)         805465-00         Orig3           Main PCB/Patient Connector Wire Harness (W8)         803783-00         B5           AC Receptacle/Ground Wire Harness (W9)         803783-00         B5           AC Receptacle Wire Harness (W10)         803783-05         B5           Test Load Contact/Resistor Wire Harness (W11)         803783-05         B5           Power Conversion PCB/Reed Assembly Wire Harness (W12)         803783-07         B5           Power Conversion PCB/Charge Relay Negative Wire Harness (W14)         803783-08         B5           Power Conversion PCB/Charge Relay Positive Wire Harness (W15)         803783-10         B5 <td>Strip Chart Recorder Assembly</td> <td>804189-00</td> <td>B2</td>	Strip Chart Recorder Assembly	804189-00	B2
Display PCB/CRT Assembly Power Supply PCB/AC Receptacle/Ground Wire Harness (W1) Power Supply PCB/AC Receptacle/Ground Wire Harness (W2) Battery Sensor/Power Supply PCB/Battery Wire Harness (W2) Battery Sensor/Power Supply PCB Wire Harness (W3) Battery Sensor/Power Supply PCB Wire Harness (W5) Battery Sensor/Power Supply PCB Wire Harness (W6) Battery Sensor/Power Supply PCB Wire Harness (W6) Battery Sensor/Power Supply PCB Wire Harness (W6) Battery Sensor/Power Supply PCB/Edb Assembly (W7) Battery Sensor/Power Supply PCB/Edb Assembly (W6) Battery Sensor/Power Supply PCB/Edb Assembly (W6) Battery Sensor/Power Supply PCB/Edb Assembly (W6) Battery Sensor/Power Supply PCB/Edb Assembly Wire Harness (W12) Battery Sensor/Power Supply Bracket/Ground Wire Harness (W16) Battery Sensor/Power Supply Bracket/Ground Wire Harness (W18) Battery Sensor/Power Supply Bracket/Ground Wire Harness (W19) Battery Sensor/Power Supply Bracket/Ground Wire Harness (W22) Battery Sensor/Power Supply Bracket/Ground Wire Harness (W23) Battery Sensor/Power Supply Bracket/Ground Wire Harness (W22) Battery Sensor/Power Supply Bracket/Ground Wire Harness (W23) Battery Sensor/Power Supply Bracket/Ground Wire Harness (W23) Battery	Transfer Relay Assembly	800240-14	<b>U10</b>
Power Supply PCB/AC Receptacle/Ground Wire Harness (W1) 803783-30 B5 Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2) 803783-23 B5 Battery Sensor/Power Supply PCB Wire Harness (W3) 803783-24 B5 High Voltage Connector Cable Assembly (W4) 803756-13 D4 Main PCB/Monitor PCB Cable Assembly (W5) 803733-00 Orig1 Main PCB/Recorder Cable Assembly (W6) 803773-00 Orig1 Interconnect Pacing PCB/Keypad Cable Assembly (W7) 805465-00 Orig3 Main PCB/Patient Connector Wire Harness (W8) 803783-00 B5 AC Receptacle/Ground Wire Harness (W9) 803783-04 B5 AC Receptacle Wire Harness (W10) 803783-05 B5 Test Load Contact/Resistor Wire Harness (W11) 803783-06 B5 Power Conversion PCB/Reed Assembly Wire Harness (W12) 803783-07 B5 Power Conversion PCB/Charge Relay Wire Harness (W13) 803783-08 B5 Power Conversion PCB/Charge Relay Positive Wire Harness (W14) 803783-09 B5 Power Conversion PCB/Charge Relay Positive Wire Harness (W15) 803783-10 B5 Main PCB/Sternum Relay Wire Harness (W16) 803783-11 B5 Main PCB/Apex Relay Wire Harness (W16) 803783-12 B5 Storage Capacitor Positive/Relay Wire Harness (W18) 803783-13 B5 Storage Capacitor Negative/Relay Wire Harness (W19) 803783-15 B5 Main PCB/Speaker Wire Harness (W20) 803783-16 B5 Power Supply Bracket/Ground Wire Harness (W22) 803783-29 B5 Connector Receptacle/Pacing Cable Assembly (W23) 802925-03 C6	Paddle Assembly	802901-06	F
Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2) 803783-23 B5 Battery Sensor/Power Supply PCB Wire Harness (W3) 803783-24 B5 High Voltage Connector Cable Assembly (W4) 803756-13 D4 Main PCB/Monitor PCB Cable Assembly (W5) 803733-00 Orig1 Main PCB/Recorder Cable Assembly (W6) 803773-00 Orig1 Interconnect Pacing PCB/Keypad Cable Assembly (W7) 805465-00 Orig3 Main PCB/Patient Connector Wire Harness (W8) 803783-00 B5 AC Receptacle/Ground Wire Harness (W9) 803783-04 B5 AC Receptacle Wire Harness (W10) 803783-05 B5 Test Load Contact/Resistor Wire Harness (W11) 803783-06 B5 Power Conversion PCB/Reed Assembly Wire Harness (W12) 803783-07 B5 Power Conversion PCB/Dump Relay Wire Harness (W13) 803783-08 B5 Power Conversion PCB/Charge Relay Negative Wire Harness (W14) 803783-10 B5 Main PCB/Sternum Relay Wire Harness (W16) 803783-11 B5 Main PCB/Apex Relay Wire Harness (W17) 803783-12 B5 Storage Capacitor Positive/Relay Wire Harness (W18) 803783-13 B5 Storage Capacitor Positive/Relay Wire Harness (W19) 803783-15 B5 Main PCB/Speaker Wire Harness (W20) 803783-16 B5 Power Supply Bracket/Ground Wire Harness (W22) 803783-29 B5 Connector Receptacle/Pacing Cable Assembly (W23) 802925-03 C6	Display PCB/CRT Assembly	803706	A5
Battery Sensor/Power Supply PCB Wire Harness (W3)  Righ Voltage Connector Cable Assembly (W4)  Main PCB/Monitor PCB Cable Assembly (W5)  Main PCB/Recorder Cable Assembly (W6)  Main PCB/Recorder Cable Assembly (W6)  Main PCB/Recorder Cable Assembly (W6)  Main PCB/Recorder Cable Assembly (W7)  Most465-00  Orig1  Interconnect Pacing PCB/Keypad Cable Assembly (W7)  Main PCB/Patient Connector Wire Harness (W8)  AC Receptacle/Ground Wire Harness (W9)  AC Receptacle Wire Harness (W10)  Test Load Contact/Resistor Wire Harness (W11)  Power Conversion PCB/Reed Assembly Wire Harness (W12)  Power Conversion PCB/Dump Relay Wire Harness (W13)  Power Conversion PCB/Charge Relay Negative Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W15)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Positive/Relay Wire Harness (W19)  B5  Main PCB/Speaker Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Power Supply PCB/AC Receptacle/Ground Wire Harness (W1)	803783-30	B5
High Voltage Connector Cable Assembly (W4)  Main PCB/Monitor PCB Cable Assembly (W5)  Main PCB/Recorder Cable Assembly (W6)  Interconnect Pacing PCB/Keypad Cable Assembly (W7)  Main PCB/Patient Connector Wire Harness (W8)  AC Receptacle/Ground Wire Harness (W9)  AC Receptacle Wire Harness (W10)  Test Load Contact/Resistor Wire Harness (W11)  Power Conversion PCB/Dump Relay Wire Harness (W12)  Power Conversion PCB/Charge Relay Negative Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W15)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Speaker Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  B5  Main PCB/Speaker Wire Harness (W19)  803783-15  B5  Main PCB/Speaker Wire Harness (W20)  M303783-16  B5  Main PCB/Speaker Wire Harness (W21)  803783-15  B5  Connector Receptacle/Pacing Cable Assembly (W23)  803783-29  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Power Conversion PCB/Power SupplyPCB/Battery Wire Harness (W2)	803783-23	B5
Main PCB/Monitor PCB Cable Assembly (W5)803733-00Orig1Main PCB/Recorder Cable Assembly (W6)803773-00Orig1Interconnect Pacing PCB/Keypad Cable Assembly (W7)805465-00Orig3Main PCB/Patient Connector Wire Harness (W8)803783-00B5AC Receptacle/Ground Wire Harness (W9)803783-04B5AC Receptacle Wire Harness (W10)803783-05B5Test Load Contact/Resistor Wire Harness (W11)803783-06B5Power Conversion PCB/Reed Assembly Wire Harness (W12)803783-07B5Power Conversion PCB/Charge Relay Wire Harness (W13)803783-08B5Power Conversion PCB/Charge Relay Negative Wire Harness (W14)803783-09B5Power Conversion PCB/Charge Relay Positive Wire Harness (W15)803783-10B5Main PCB/Sternum Relay Wire Harness (W16)803783-11B5Main PCB/Apex Relay Wire Harness (W17)803783-12B5Storage Capacitor Positive/Relay Wire Harness (W18)803783-13B5Storage Capacitor Negative/Relay Wire Harness (W19)803783-14B5ECG Out Wire Harness (W20)803783-15B5Main PCB/Speaker Wire Harness (W21)803783-16B5Power Supply Bracket/Ground Wire Harness (W22)803783-29B5Connector Receptacle/Pacing Cable Assembly (W23)802925-03C6	Battery Sensor/Power Supply PCB Wire Harness (W3)	803783-24	B5
Main PCB/Recorder Cable Assembly (W6)803773-00Orig1Interconnect Pacing PCB/Keypad Cable Assembly (W7)805465-00Orig3Main PCB/Patient Connector Wire Harness (W8)803783-00B5AC Receptacle/Ground Wire Harness (W9)803783-04B5AC Receptacle Wire Harness (W10)803783-05B5Test Load Contact/Resistor Wire Harness (W11)803783-06B5Power Conversion PCB/Reed Assembly Wire Harness (W12)803783-07B5Power Conversion PCB/Dump Relay Wire Harness (W13)803783-08B5Power Conversion PCB/Charge Relay Negative Wire Harness (W14)803783-09B5Power Conversion PCB/Charge Relay Positive Wire Harness (W15)803783-10B5Main PCB/Sternum Relay Wire Harness (W16)803783-11B5Main PCB/Apex Relay Wire Harness (W17)803783-12B5Storage Capacitor Positive/Relay Wire Harness (W18)803783-13B5Storage Capacitor Negative/Relay Wire Harness (W19)803783-14B5ECG Out Wire Harness (W20)803783-15B5Main PCB/Speaker Wire Harness (W21)803783-16B5Power Supply Bracket/Ground Wire Harness (W22)803783-29B5Connector Receptacle/Pacing Cable Assembly (W23)802925-03C6	High Voltage Connector Cable Assembly (W4)	803756-13	D4
Interconnect Pacing PCB/Keypad Cable Assembly (W7)  Main PCB/Patient Connector Wire Harness (W8)  AC Receptacle/Ground Wire Harness (W9)  AC Receptacle Wire Harness (W10)  Test Load Contact/Resistor Wire Harness (W11)  Power Conversion PCB/Reed Assembly Wire Harness (W12)  Power Conversion PCB/Dump Relay Wire Harness (W13)  Power Conversion PCB/Charge Relay Negative Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W15)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  B5  Main PCB/Speaker Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  B6  Storage Capacitor Receptacle/Pacing Cable Assembly (W23)  803783-29  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Main PCB/Monitor PCB Cable Assembly (W5)	803733-00	Orig1
Main PCB/Patient Connector Wire Harness (W8)  AC Receptacle/Ground Wire Harness (W9)  AC Receptacle Wire Harness (W10)  Rest Load Contact/Resistor Wire Harness (W11)  Power Conversion PCB/Reed Assembly Wire Harness (W12)  Power Conversion PCB/Dump Relay Wire Harness (W13)  Power Conversion PCB/Charge Relay Negative Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W15)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-09  B5  803783-09  B5  803783-10  B5  803783-11  B5  803783-12  B5  803783-13  B5  Storage Capacitor Negative/Relay Wire Harness (W19)  803783-14  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Main PCB/Recorder Cable Assembly (W6)	803773-00	Orig1
AC Receptacle/Ground Wire Harness (W9)  AC Receptacle Wire Harness (W10)  Rest Load Contact/Resistor Wire Harness (W11)  Power Conversion PCB/Reed Assembly Wire Harness (W12)  Power Conversion PCB/Dump Relay Wire Harness (W12)  Power Conversion PCB/Dump Relay Wire Harness (W13)  Power Conversion PCB/Charge Relay Negative Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W15)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-09  B5  803783-09  B5  803783-11  B5  803783-15  B5  803783-15  B5  803783-16  B5  Power Supply Bracket/Ground Wire Harness (W22)  803783-29  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Interconnect Pacing PCB/Keypad Cable Assembly (W7)	805465-00	Orig3
AC Receptacle Wire Harness (W10)  Test Load Contact/Resistor Wire Harness (W11)  Power Conversion PCB/Reed Assembly Wire Harness (W12)  Power Conversion PCB/Dump Relay Wire Harness (W13)  Power Conversion PCB/Charge Relay Negative Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W14)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-05  B5  803783-08  B5  803783-09  B5  803783-10  B5  803783-11  B5  803783-12  B5  803783-13  B5  803783-14  B5  B5  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)	Main PCB/Patient Connector Wire Harness (W8)	803783-00	B5
Test Load Contact/Resistor Wire Harness (W11)  Power Conversion PCB/Reed Assembly Wire Harness (W12)  Power Conversion PCB/Dump Relay Wire Harness (W13)  Power Conversion PCB/Charge Relay Negative Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W15)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-09  B5  B5  B6  B7  B7  B8  B8  B8  B8  B8  B8  B8  B8	AC Receptacle/Ground Wire Harness (W9)	803783-04	B5
Power Conversion PCB/Reed Assembly Wire Harness (W12) 803783-07 B5 Power Conversion PCB/Dump Relay Wire Harness (W13) 803783-08 B5 Power Conversion PCB/Charge Relay Negative Wire Harness (W14) 803783-09 B5 Power Conversion PCB/Charge Relay Positive Wire Harness (W15) 803783-10 B5 Main PCB/Sternum Relay Wire Harness (W16) 803783-11 B5 Main PCB/Apex Relay Wire Harness (W17) 803783-12 B5 Storage Capacitor Positive/Relay Wire Harness (W18) 803783-13 B5 Storage Capacitor Negative/Relay Wire Harness (W19) 803783-14 B5 ECG Out Wire Harness (W20) 803783-15 B5 Main PCB/Speaker Wire Harness (W21) 803783-16 B5 Power Supply Bracket/Ground Wire Harness (W22) 803783-29 B5 Connector Receptacle/Pacing Cable Assembly (W23) 802925-03 C6	AC Receptacle Wire Harness (W10)	803783-05	B5
Power Conversion PCB/Dump Relay Wire Harness (W13)  Power Conversion PCB/Charge Relay Negative Wire Harness (W14)  Power Conversion PCB/Charge Relay Positive Wire Harness (W14)  B03783-09  B5  Power Conversion PCB/Charge Relay Positive Wire Harness (W15)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  B5  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-09  B5  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Test Load Contact/Resistor Wire Harness (W11)	803783-06	B5
Power Conversion PCB/Charge Relay Negative Wire Harness (W14) 803783-09 B5 Power Conversion PCB/Charge Relay Positive Wire Harness (W15) 803783-10 B5 Main PCB/Sternum Relay Wire Harness (W16) 803783-11 B5 Main PCB/Apex Relay Wire Harness (W17) 803783-12 B5 Storage Capacitor Positive/Relay Wire Harness (W18) 803783-13 B5 Storage Capacitor Negative/Relay Wire Harness (W19) 803783-14 B5 ECG Out Wire Harness (W20) 803783-15 B5 Main PCB/Speaker Wire Harness (W21) 803783-16 B5 Power Supply Bracket/Ground Wire Harness (W22) 803783-29 B5 Connector Receptacle/Pacing Cable Assembly (W23) 802925-03 C6	Power Conversion PCB/Reed Assembly Wire Harness (W12)	803783-07	B5
Power Conversion PCB/Charge Relay Positive Wire Harness (W15)  Main PCB/Sternum Relay Wire Harness (W16)  Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-10  B5  803783-12  B5  803783-13  B5  803783-14  B5  B5  C6	Power Conversion PCB/Dump Relay Wire Harness (W13)	803783-08	B5
Main PCB/Sternum Relay Wire Harness (W16)803783-11B5Main PCB/Apex Relay Wire Harness (W17)803783-12B5Storage Capacitor Positive/Relay Wire Harness (W18)803783-13B5Storage Capacitor Negative/Relay Wire Harness (W19)803783-14B5ECG Out Wire Harness (W20)803783-15B5Main PCB/Speaker Wire Harness (W21)803783-16B5Power Supply Bracket/Ground Wire Harness (W22)803783-29B5Connector Receptacle/Pacing Cable Assembly (W23)802925-03C6	Power Conversion PCB/Charge Relay Negative Wire Harness (W14)	803783-09	B5
Main PCB/Apex Relay Wire Harness (W17)  Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-12  B5  803783-13  B5  803783-15  B5  803783-29  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Power Conversion PCB/Charge Relay Positive Wire Harness (W15)	803783-10	B5
Storage Capacitor Positive/Relay Wire Harness (W18)  Storage Capacitor Negative/Relay Wire Harness (W19)  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-13  B5  803783-14  B5  B5  803783-15  B5  803783-29  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Main PCB/Sternum Relay Wire Harness (W16)	803783-11	B5
Storage Capacitor Negative/Relay Wire Harness (W19)  ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-14  B5  803783-15  B5  803783-16  B5  803783-29  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Main PCB/Apex Relay Wire Harness (W17)	803783-12	B5
ECG Out Wire Harness (W20)  Main PCB/Speaker Wire Harness (W21)  Power Supply Bracket/Ground Wire Harness (W22)  Connector Receptacle/Pacing Cable Assembly (W23)  803783-15  B5  803783-16  B5  803783-29  B5  Connector Receptacle/Pacing Cable Assembly (W23)	Storage Capacitor Positive/Relay Wire Harness (W18)	803783-13	B5
Main PCB/Speaker Wire Harness (W21)803783-16B5Power Supply Bracket/Ground Wire Harness (W22)803783-29B5Connector Receptacle/Pacing Cable Assembly (W23)802925-03C6	Storage Capacitor Negative/Relay Wire Harness (W19)	803783-14	B5
Power Supply Bracket/Ground Wire Harness (W22) 803783-29 B5 Connector Receptacle/Pacing Cable Assembly (W23) 802925-03 C6	ECG Out Wire Harness (W20)	803783-15	B5
Connector Receptacle/Pacing Cable Assembly (W23) 802925-03 C6	Main PCB/Speaker Wire Harness (W21)	803783-16	B5
• • • • • • • • • • • • • • • • • • • •	Power Supply Bracket/Ground Wire Harness (W22)	803783-29	B5
Recorder/Bracket Gound Wire Harness (W24) 803783-19 B5	Connector Receptacle/Pacing Cable Assembly (W23)	802925-03	C6
		803783-19	B5

Assembly Name	Part Number	Rev	
Defibrillation Adapter	803747	B2	
Test Load High Voltage Wire Harness	803774-03	<b>A</b> 1	
Test Load High Volage Wire Harness	803774-04	<b>A</b> 1	
Pushbutton Switch Wire Harness	803774-05	<b>A1</b>	
High Voltage Connector Cable Assembly	803756-09	D4	

**Content Overview** This manual contains the following information:

#### Section 1

**Description:** This section details how the instrument works. Input signals, power supplies, PCB functions, and instrument outputs are described. Each PCB circuit description is accompanied by a block diagram to illustrate the major circuits. Circuit names in these block diagrams also appear in the appropriate schematics in Section 5.

#### Section 2

**Operation:** This section familiarizes the user with basic equipment function. It identifies Controls, Indicators, and Connectors, and screen menu options. This section is not intended to instruct the operator in the clinical use of the instrument; such detailed instructions are provided in the separate *Operating* Instructions manual.

#### Section 3

**Testing/Troubleshooting:** This section contains the Performance Inspection Procedure (PIP)—sequential steps to follow when performing an operational closed-case check of the equipment. A PIP checklist is provided which can be duplicated and used during testing. The Test and Calibration Procedure (TCP) in this section describes calibration and more extensive instrument testing; it also includes a checklist. The Troubleshooting Aids includes additional information to support troubleshooting.

#### Section 4

Service and Maintenance: This section provides Disassembly/Assembly Procedures for removing all major subassemblies. The procedures are referenced to numbered parts in the Final Assembly drawing in Section 5. This section includes instructions for inspecting, cleaning, maintaining, and repairing the instrument.

#### Section 5

Parts Lists/Schematics: This section contains a list of Supplies and Accessories, Illustrated Parts Lists, PCB Component Layouts, and Schematic Diagrams for all repairable assemblies.

#### Section 6

Component Reference Diagrams: An aid to troubleshooting, this section consists of Component Reference Diagrams for selected Integrated Circuits (ICs). The ICs in Section 6 are identified in the schematics in Section 5 with a large asterisk \*.

 ·	<del>-</del>	., -	-
•			· .
			•
			*

#### Introduction

This section describes the general features, specifications, functions, and theory of operation of the LIFEPAK 9P defibrillator/monitor/pacemaker. The section is divided into three parts:

- Physical Description describes general features and lists specifications
- Functional Description briefly describes the function of the major assemblies
- Theory of Operation provides circuit descriptions to the major component level.

# Physical Description

The LIFEPAK 9P defibrillator/monitor/pacemaker, shown in Figure 1–1, provides dc defibrillation, synchronous cardioversion, a Cathode Ray Tube (CRT) monitor screen, a strip chart recorder, 3-lead electrocardiogram (ECG) monitoring, heart rate display and alarms, multiple energy selection, external pacing, and CODE SUMMARY critical event record. The ECG signal is detected through patient cable electrodes, QUIK-LOOK defibrillation paddles, or FAST-PATCH disposable defibrillation electrodes.

The instrument operates on either ac power or internal battery power. The sealed lead-acid battery, accessible from the bottom of the instrument, allows temporary operation during ac power failure or disconnection. With the instrument plugged into ac power and the rear panel power switch on, the battery continually recharges, even with the front panel power off.

The front panel pushbuttons provide control for most functions. Charge and discharge pushbuttons are located on the defibrillator paddles. Four rear-panel pushbuttons control the QRS audio volume, a calibration pulse, and the real-time clock.




Figure 1-1 LIFEPAK 9P defibrillator/monitor/pacemaker

The CRT monitor screen provides a no-fade ECG display with a sweep rate of  $25\pm1 \text{mm/second}$ . The monitor screen displays heart rate, QRS indicator, selected lead, ECG gain, heart-rate violations, status messages, defibrillator and pacing information, sync markers, and operation status. Additional menus allow access to power-on default parameters, diagnostic tests, and configuration information.

The defibrillator delivers a controlled monophasic dc defibrillating pulse to the patient with selectable energy levels ranging from 0 to 360 joules. Both the selected energy and the available energy are displayed on the monitor screen. Defibrillation pulses may be delivered through standard paddles, clip-on posterior or pediatric paddles, or optional FAST-PATCH disposable defibrillation/ECG electrodes (with the optional Defibrillation Adapter or Shock Advisory Adapter and cable). When the paddles are stored in the paddle wells, an internal test load is accessible for testing and confirming defibrillation energy delivery.

The recorder can print the ECG trace, time and date, ECG lead and size, heart rate, pacing parameters, and SYNC (if activated). The recorder automatically prints for approximately 16 seconds when a heart rate violation occurs. When the CODE SUMMARY pushbutton is pressed, the recorder prints the stored CODE SUMMARY critical event record data, including event summary and recorded ECG trace data.

Pacing is selectable in either demand or non-demand (asynchronous) mode. Demand mode depends on cardiac activity sensing; non-demand mode paces at a selected rate regardless of cardiac activity. Stimulus markers followed by any cardiac activity are shown on the display along with pacing current and rate.

The ac line power fuses, a ground connector, and an ECG signal output connector are accessible at the rear panel. The LIFEPAK 9P defibrillator/ monitor/pacemaker specifications are listed in Table 1–1.

LIFEPAK 9P defibrillator/monitor/pacemaker Specifications\* Table 1-1

Std, Paddles, I, II, III.  INPUT  Isolated ECG via QUIK-LOOK defibrillator paddles, FAST-PATCH disposable defibrillation/ECG electrodes, or 3-lead patient cable.  PATIENT CABLE LENGTH  COMMON MODE REJECTION  With notch filter engaged, 100dB minimum with respect to chassis ground and 65dB minimum with respect to isolated ground when measured at 60Hz. Common mode range for patient cable input ≥10V peak with respect to isolated ground.  MONITOR DISPLAY SIZE  SWEEP SPEED  FREQUENCY RESPONSE  ECG Leads: Paddles: Paddles:  1.0 to 40Hz (-3dB) 2.2 to 20Hz (-3dB)  ECG SIZE  Adjusts amplitude of ECG trace on monitor, strip chart recorder, and ECG out.  Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display. Heart rate meter is disabled during pacing.
disposable defibrillation/ECG electrodes, or 3-lead patient cable.  PATIENT CABLE LENGTH  COMMON MODE REJECTION  With notch filter engaged, 100dB minimum with respect to chassis ground and 65dB minimum with respect to isolated ground when measured at 60Hz. Common mode range for patient cable input ≥10V peak with respect to isolated ground.  MONITOR DISPLAY SIZE  SWEEP SPEED  FREQUENCY RESPONSE  ECG Leads: Paddles:  1.0 to 40Hz (-3dB) 2.2 to 20Hz (-3dB)  ECG SIZE  Adjusts amplitude of ECG trace on monitor, strip chart recorder, and ECG out.  HEART RATE METER  Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display.
COMMON MODE REJECTION  With notch filter engaged, 100dB minimum with respect to chassis ground and 65dB minimum with respect to isolated ground when measured at 60Hz. Common mode range for patient cable input ≥10V peak with respect to isolated ground.  MONITOR DISPLAY SIZE  SWEEP SPEED  102mm (4in) wide x 76mm (3in) tall, non-fade.  25mm/s  FREQUENCY RESPONSE  ECG Leads: Paddles:  1.0 to 40Hz (-3dB) 2.2 to 20Hz (-3dB)  Adjusts amplitude of ECG trace on monitor, strip chart recorder, and ECG out.  HEART RATE METER  Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display.
ground and 65dB minimum with respect to isolated ground when measured at 60Hz. Common mode range for patient cable input ≥10V peak with respect to isolated ground.  MONITOR DISPLAY SIZE  102mm (4in) wide x 76mm (3in) tall, non-fade.  25mm/s  FREQUENCY RESPONSE  ECG Leads: Paddles:  1.0 to 40Hz (-3dB) 2.2 to 20Hz (-3dB)  ECG SIZE  Adjusts amplitude of ECG trace on monitor, strip chart recorder, and ECG out.  HEART RATE METER  Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display.
SWEEP SPEED  FREQUENCY RESPONSE  ECG Leads: Paddles:  1.0 to 40Hz (-3dB) 2.2 to 20Hz (-3dB)  ECG SIZE  Adjusts amplitude of ECG trace on monitor, strip chart recorder, and ECG out.  HEART RATE METER  Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display.
FREQUENCY RESPONSE  ECG Leads: Paddles: 1.0 to 40Hz (-3dB) 2.2 to 20Hz (-3dB)  ECG SIZE  Adjusts amplitude of ECG trace on monitor, strip chart recorder, and ECG out.  HEART RATE METER  Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display.
ECG Leads: Paddles: 1.0 to 40Hz (-3dB) 2.2 to 20Hz (-3dB)  ECG SIZE Adjusts amplitude of ECG trace on monitor, strip chart recorder, and ECG out.  HEART RATE METER Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display.
Paddles: 2.2 to 20Hz (-3dB)  ECG SIZE Adjusts amplitude of ECG trace on monitor, strip chart recorder, and ECG out.  HEART RATE METER Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display.
ECG out.  HEART RATE METER  Three-digit readout displays rates from 20 to 300bpm. Heart rates outside this range do not yield valid systole tones or heart rate display.
outside this range do not yield valid systole tones or heart rate display.
HEART RATE ALARM  User-selectable alarm limits. Three high/low settings: 150/40, 120/60, 160/90. Other limit options available through the Setup Menu. Heart rate alarms are disabled during pacing.
1 mV Cal Momentary pushbutton on rear panel simulates a 1mV signal pulse to the ECG input.
ECG OUTPUT 1V/mV at x1.0 gain

<sup>&#</sup>x27;All specifications at 20°C unless otherwise stated. Specifications subject to change without notice.

Table 1-1 LIFEPAK 9P defibrillator/monitor/pacemaker Specifications\* (cont.)

#### **DEFIBRILLATOR**

**ENERGY SELECT** 

External Paddles:

1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 20, 30, 50, 100, 200, 300, 360J

Internal Paddles:

1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 20, 30, 50J

**CHARGE CONTROLS** 

Independent momentary pushbutton controls on front panel and APEX

paddle.

PADDLE CORD LENGTH

3m (10ft)

**CHARGE INDICATORS** 

Flashing lamps on paddle and front panel pushbutton along with increasing stored energy display on monitor indicate charge in progress. Upon full charge, energy available is displayed and charge

completed tone sounds.

**CHARGE TIME** 

SYNC INDICATOR

Charge to 360J in less than 10s with a fully-charged battery.

**SYNC** 

Synchronizes defibrillator pulse to patient-generated QRS complex.

Inverted triangle marker on displayed ECG waveform identifies synchronizer trigger point with respect to patient's QRS complex.

DEFIBRILLATOR ENERGY

ACCURACY

 $\pm\,15\%$  at  $50\Omega$ 

 $\pm 30\%$  over the range of 25 to  $100\Omega$ 

I (amps) 80

70


60


50


40

**DEFIBRILLATOR WAVEFORM** 

5ms monophasic pulse (Edmark)







\*All specifications at 20°C unless otherwise stated. Specifications subject to change without notice.

25 Ohm Load 360 Joule

T (msec)

### Table 1-1 LIFEPAK 9P defibrillator/monitor/pacemaker Specifications\* (cont.)

#### THERMAL ARRAY RECORDER

PAPER

Size: 50mm wide x 30m (100ft)

25mm/s Speed:

**DELAY** ECG prints 8 seconds after first appearing on the monitor.

FREQUENCY RESPONSE Non-diagnostic (diagnostic available via Setup Menu).

ECG Leads (non-diagnostic): 1.0 to 40Hz (-3dB). 0.05 to 100Hz (-3dB). ECG Leads (diagnostic):

Paddles: 2.2 to 20Hz (-3dB).

ANNOTATION Time, date, ECG lead, ECG gain, heart rate, defibrillation parameters,

test load discharges, pacemaker parameters.

**CODE SUMMARY Critical Event** Digitally stored record of ECG and device parameters.

Record

#### **NON-INVASIVE PACEMAKER**

**OUTPUT RATE** 40 to 170ppm.

RATE ACCURACY  $\pm 5\%$  over entire range.

**OUTPUT WAVEFORM** Monophasic, truncated exponential current pulse (15 to 25% droop). **OUTPUT CURRENT** 0 to 200mA  $\pm$  5% or 2mA (whichever is greater) at 700 $\Omega$  and  $\pm$  7% or

3mA (whichever is greater) over the range of 100 to  $1500\Omega$ 

#### REFRACTORY PERIOD

Pacing Rates: **Refractory Periods:** 40 to 100ppm 300ms

110 to 120ppm 250ms 130 to 140ppm 220ms 150 to 170ppm 200ms

MODE Demand or non-demand.

#### GENERAL

AC INPUT OPTIONS 120 or 240Vac nominal line voltage without adjustment. 50 or 60Hz

with adjustment in Setup Menu.

**BATTERY TYPE** Sealed lead-acid, 16Vdc nominal, 3Ahr.

**BATTERY CAPACITY** A new, fully-charged battery will provide one of the following prior to

shutdown:

	<u>Typical</u>	<u>Minimum</u>
Number of 360J discharges:	. 75	40
Minutes of monitoring:	90	70
Minutes of pacing:	75	50

\*All specifications at 20°C unless otherwise stated. Specifications subject to change without notice.

### Table 1-1 LIFEPAK 9P defibrillator/monitor/pacemaker Specifications\* (cont.)

#### **GENERAL** (cont.)

LOW BATTERY INDICATOR

Advises operator to connect ac power. Battery capacities before and after warning (LOW BATTERY-CONNECT AC POWER) for a new,

fully charged battery are:

Typical Minimum Number of 360J discharges before warning: 50 20 Number of 360J discharges after warning: 25 5

BATTERY CHARGE INDICATOR Indicator BATT CHRG illuminates when battery is charging.

**BATTERY CHARGE TIME** 

24hrs to full capacity.

SERVICE INDICATOR

Message SERVICE indicates self-diagnostic routines have detected

improper operation requiring service attention.

POWER CONSUMPTION

160W maximum while monitoring with recorder on and defibrillator

charging.

SIZE

Height: 35.2cm (13.9in) Width: 29.7cm (11.7in) Depth: 31.0cm (12.2in) Weight: 13.2kg (29lbs)

POWER CORD LENGTH

3m (10ft)

STANDARD PADDLE

 $82 \text{cm}^2$ 

**ELECTRODE AREA** 

#### **ENVIRONMENTAL**

ATMOSPHERIC PRESSURE

797 to 500mm Hg (-570 to 11,000ft).

RELATIVE HUMIDITY

0 to 95% (non-condensing) at 0 to 34°C (32 to 94°F).

0 to 80% (non-condensing) at 35 to 45°C (95 to 113°F).

TEMPERATURE RANGE

0 to 45°C (32 to 113°F) operating. -30 to 65°C (-22 to 149°F) storage.

\*All specifications at 20°C unless otherwise stated. Specifications subject to change without notice.

## **Functional** Description

The overall function of the LIFEPAK 9P defibrillator/monitor/pacemaker is shown in Figure 1–2. The defibrillator/monitor/pacemaker includes six mechanical assemblies (Keypad, Rear Panel Membrane Switch, Strip Chart Recorder, Transfer Relay, Paddles, and Battery Pack) and five Printed Circuit Board (PCB) assemblies (Power Supply, Power Conversion, Interconnect Pacing, Main, and Display).

The Keypad and Rear Panel Membrane Switch Assemblies receive operator input. The Strip Chart Recorder Assembly provides printed output. The Transfer Relay and Paddle Assemblies help deliver defibrillator pulses. The Battery Pack provides temporary backup power.

The Power Supply PCB functions over the full range of input voltages from 90Vac to 270Vac at 50Hz or 60Hz without input line voltage selection. The Power Supply PCB converts the ac line power to the nominal +18.8Vdc supplied to the Power Conversion PCB.

The Power Conversion PCB performs a variety of power-handling functions: power on/off switching, generation of additional power supply voltages, defibrillation charging, and stored energy transfer. A test load circuit is also present on the Power Conversion PCB.

The Interconnect Pacing PCB provides all the pacing circuitry and the electrical connection between the Main PCB and the Power Conversion PCB. It also contains circuitry to produce the audible tones and to drive the front panel LEDs.

The Main PCB has three major sections of circuitry: the Preamp, the System Controller, and the Display/Recorder/Pacing Interface. The Preamp receives input from both the 3-lead patient cable and QUIK-LOOK paddles inputs. The System Controller directs the overall operation of the other circuits. QRS detection, rate calculation and alarms, and synchronous cardioversion are all implemented in the System Controller software. The current state of the System microprocessor and the ECG data are transmitted from the System microprocessor to the Display/Recorder/Pacing microprocessor, allowing the Display/Recorder/Pacing microprocessor to select and format the appropriate CRT messages and recorder annotations.

The Display PCB/CRT Assembly contains a Monitor PCB that drives the 5-inch diagonal Cathode Ray Tube (CRT). It uses electromagnetic beam deflection and incorporates vertical and horizontal deflection circuitry to produce a raster scan.




Figure 1–2 LIFEPAK 9P defibrillator/monitor/pacemaker Functional Block Diagram

## Theory of Operation

This section contains detailed circuit descriptions of the Power Supply PCB, Power Conversion PCB, Interconnect Pacing PCB, Main PCB, and Display PCB/CRT Assemblies. Refer to schematic diagrams of these electrical circuits in Section 5 while reading the circuit descriptions. Information about selected integrated circuits (indicated by a \* on the schematic) is provided in Section 6.

### **Power Supply PCB Assembly (803726)**

The Power Supply PCB converts input ac voltage to a nominal dc voltage for the Power Conversion PCB and for charging the battery.

Power-On

When ac current flows through J1, it is conditioned by RT1, L1, C2, C3, and CR1. Resistor divider R1, R3, and R4 in series with parallel resistors R11 and R12 turn on Q1. Current builds in T1 primary windings (pins 4 and 6) and is coupled with pins 1 and 2. Capacitor C6 couples the current build-up into Q1, causing saturation.

Increasing current through the T1 primary produces voltage on the base of Q2 and a voltage drop across parallel resistors R11 and R12. Q2 activation overrides Q1 and terminates the primary current. Transformer T2 sustains base drive current to O2.

Regulation

When battery voltage exceeds the voltage threshold set by VR1, Q4 turns on and couples the voltage from T2 pins 4 and 5 to T2 pins 2 and 7. When current subsides to proper levels, Q4 turns off, providing current feedback for regulation.

Overvoltage Protection

If the regulation circuit fails, Silicon Controlled Rectifier (SCR) CR12, together with CR13, R17, and R18, protects the load from damage. Diode CR13 sets a voltage threshold which, when exceeded, triggers SCR CR12 to short the Power Supply PCB circuit output. Removing ac power resets the protection circuit.

Output Isolation

Diode CR11 isolates the Power Supply PCB output from the filter network to protect it from the battery voltage when ac power is removed. This is necessary because the battery is connected in parallel with the Power Supply PCB.

Battery

The 16Vdc, lead-acid battery is wired directly to the Power Supply PCB output. The battery receives a constant charging current whenever the instrument is connected to ac power and the rear panel mains power switch is on. As the battery approaches full charge, the supply operates in constant

voltage mode. Current limiting sets the maximum charge current available to the battery at approximately +4.5A. This limit is implemented when the Power Supply PCB sources current into a dead short, a requirement when charging a severely depleted battery.

#### **BATT CHRG Indicator**

A battery-charging LED signal illuminates the front panel BATT CHRG indicator whenever the ac power is connected to the instrument and the rear panel mains power switch is on. An interlock in the battery harness prevents the indicator from lighting if the battery is not installed.

### Power Conversion PCB Assembly (803724)

The Power Conversion PCB contains circuitry for power on/off switching, additional voltage generation and monitoring, energy storage capacitor charging, and defibrillation energy transfer. Refer to Figure 1–3, page 1–11.

#### **Power Control**

The Power Control circuit regulates battery voltage for the circuits which control instrument operation. The circuits directly affected by the Power Control circuit are: Switching Power Supply, Logic Power Monitor, External Control Input, Low Energy Charge Rate, and Energy Storage Capacitor Charger circuits.

The Front Panel ON momentary pushbutton provides the start-up signal for the Power Control circuit. When the power switch is first pressed (turning power on), V BATT passes to the gate of Q11, causing Q11 to conduct, which in turn causes Q10 to conduct. A portion of V BATT is held on the gate of Q11 by R52, which keeps Q10 conducting. During this time the PWR SW MONITOR line is low and the Power Control circuit provides the SWBATT voltage for the rest of the Power Conversion PCB Assembly circuits.

When the power switch is pressed a second time (turning power off), V BATT is divided by R108 and R109, forward-biasing CR44 and switching the PWR SW MONITOR line high. The System microprocessor reads the low-to-high transition and activates the Logic Power Monitor circuit, disabling the Power Control circuit. The functional integrity of the Power Control circuit may be confirmed by monitoring two pins on the test connector: R60 (connects to the gate of Q10) and R57 (connects to the gate of Q11).




Figure 1-3 Power Conversion PCB Block Diagram

**Logic Power Monitor** 

The Logic Power Monitor checks the logic supply voltage. If the logic voltage is above or below normal limits, it deactivates Q11 and disables the Power Control circuit. Comparator U5 functions as a window detector for high voltages (U5A) and low voltages (U5B). The comparator reference voltages are set by Diode IC U4 at +2.5V. The output of the Power Control circuit, SWBATT, generates the operating voltage for U5, U4, and the Switching Power Supply.

The System microprocessor shuts down the instrument by sending a low signal to the POWER ENABLE line. This switches the open collector output of U5B to ground which turns off Q11, disabling the Power Control circuit. Electrical isolation is provided by CR24 and R39.

A time delay function in the Logic Power Monitor circuit prevents the Power Control circuit from shutting down while operating voltages are still rising during initial power-on. Capacitor C26 accomplishes this for the high-limit detector (U5A) and C18 for the low-limit detector (U5B).

The low- and high-limit detector inputs to the test connector are used to check circuit performance by serving as test points to measure voltage. Voltage measurements check the accuracy of voltage dividers R41, R42, R43, and R56 (relative to the +5V supply).

#### Switching Power Supply

The Switching Power Supply uses a pulse-width regulator (U3) and a flyback-coupled inductor (T2) to produce the +5V to +15V supply voltages used by the other circuits.

The pulse-width regulator U3 contains an oscillator, reference voltages, an error amplifier, and drive circuitry for Q12. Oscillation frequency is set at 90Hz by R30 and C12. The +15V output feedback at U3 pin 2 determines the outgoing pulse-width. The pulse output at U3 pin 6 provides a constant frequency pulse by driving an N-channel FET (Q12) and dumping the V BATT voltage into the primary of T2. Two diodes, CR28 and CR29, protect Q12 from high voltages caused by leakage inductance from the primary.

#### External Control Input

The Energy Storage Capacitor Charger circuit is externally controlled by the N-channel MOSFET, Q9. Two resistors (R58 and R59) provide forward-biasing for Q9 whenever power is applied to the instrument. The System microprocessor on the Main PCB initiates a charge cycle for the energy storage capacitor by setting the CHARGE INHIBIT line at low voltage, turning Q9 off. The inverting input of U2A rises to the required +5V minimum.

# Energy Storage Capacitor Charger

The Energy Storage Capacitor Charger circuit uses a flyback transformer to charge the energy storage capacitor through the Transfer Relay. The core of the circuit consists of the following components: T1, U2A and U2B, Q1, Q2, Q5-Q7, and C17.

Amplitude comparators U2A and U2B determine the drive pulse for the primary of T1. The input to the energy storage capacitor at U2A pin 2 is controlled by the External Control, Low Energy Charge Rate, and Capacitor Charge Control circuits. If any of these incoming signals fall below the +2.5V reference at U2A pin 3, the output switches high, inhibiting the drive to the primary of T1.

Once a charge cycle initiates at U2A pin 2, the active low output simultaneously performs two separate functions: transformer driver Q5 conducts, and Q7 shuts off. Timing for the transformer drive pulse-width is provided by Q7, U2B, C17, and R23. While Q7 conducts, the voltage across C17 is zero. When Q7 shuts off, the voltage at SWBATT starts charging C17. When the voltage at C17 exceeds +2.5V, the output of U2B switches low, forcing U2A output high and preventing current flow to the primary of T1.

The power circuit to the primary of T1 consists of two stages: Q5 and Q6 function as push-pull drivers, and Q1 and Q2 provide improved current handling for the drive pulse.

The secondary winding of T1 supplies the energy stored on the energy storage capacitor for defibrillation. The output is therefore subject to potentials of up to +4.5kV under normal operating conditions. Capacitor C1 limits the single pulse voltage rise to about +1.5kV, providing circuit-response time for control of the output voltage.

Low Energy Charge Rate

This circuit is centered around op amp U17A and is active only during low battery voltage levels, or when an energy level below 30J is selected. When a low battery voltage is detected, the Low Energy Charge Rate circuit slows the charge rate for the energy storage capacitor, thus reducing the amount of power required for charging. Similarly, the charge rate is slowed when a low charge energy setting is selected, thereby making sure the energy storage capacitor will not charge too quickly for the Analog Voltage Evaluation circuit to detect the energy level.

A signal from the System microprocessor at the energy request register, U20, determines the output of U17B. When an energy level of less than 30J is selected, low signals are applied to the inputs of U15D from U20. The U15D output then switches low, forcing CR45 to reverse-bias. This pulls the noninverting input of U17A to +5V, changing the output of U17A.

The signal from U17A controls the charge timing circuit (U2B, C17, and R23) and, in turn, the output of U2A. During normal charging, the effect of this signal is blocked by CR46. When U17A pin 1 is driven high, additional current is supplied to C17. This reduces the on-time of U2A which reduces current to the primary circuit and slows the rate of charge of the energy storage capacitor. This reduces the power consumption of the Energy Storage Capacitor Charger to less than half of normal. At 360J, the charge time is increased to 30 seconds.

The instrument can operate with low battery voltage if the power required for charging the energy storage capacitor is reduced. Op amp U17A detects low battery levels by monitoring the SWBATT voltage at the inverting input and comparing it with the +2.5V at the noninverting input. When the threshold limit set on the noninverting input of U17A exceeds the SWBATT voltage on the inverting input, it indicates a failing battery. The output of U17A then functions as described in the previous paragraph.

#### Capacitor Charge Control

The Capacitor Charge Control circuit controls the Energy Storage Capacitor Charger. It achieves and maintains a selected voltage on the energy storage capacitor and reports the capacitor voltage status to the System microprocessor. The main circuit components consist of op amp U17B and comparators U8A, U8B, and U8C.

Op amp U17B sets and maintains the selected voltage to the energy storage capacitor by controlling an inhibit signal to pin 2 of U2A. When a selected voltage is achieved, U17B shuts down U2A which disables the Energy Storage Capacitor Charger.

To set the gain for a particular charge voltage, digital-to-analog (D/A) converter U18 is connected between the inverting input and the output of U17B. The output of the D/A converter is controlled by selection of a 12-bit word by the System microprocessor. Resistors R2 and R4 function as voltage dividers for the energy storage capacitor output voltage. When this divided voltage is applied to the noninverting input of U17B, the gain of U17B rises accordingly.

Comparator U8C provides switching for the Energy Storage Capacitor Charger in response to input from U17B. When the output of U17B rises above +2.5V, U8C switches low. This low inhibit signal switches the output of U2A high, shutting down the Energy Storage Capacitor Charger.

Comparators evaluate the charge stored on the energy storage capacitor with respect to the defibrillation energy selected. U8B detects values greater than the minimum required for in-tolerance energy; U8A detects values less than the maximum allowed for the same energy setting. The outputs from U8A and U8B form the CHARGE VALID signal, indicating valid energy charge levels. This signal informs the System microprocessor about the status of the selected charge energy.

#### **Energy Transfer**

The Energy Transfer circuit transfers the energy stored in the energy storage capacitor out through the paddles. The main circuit components are U6, Q4, C9, and CR4.

When the sternum and apex pushbuttons on the paddles are pressed, the Sternum Transfer and Apex Transfer lines at J17 produce the Transfer Request signal from the output of U15A which is sent to the System microprocessor. If all energy transfer parameters are correct, the System microprocessor returns a Transfer Enable signal through U11A and U14C.

The Transfer Enable signal enables and maintains a constant current to the Transfer Relay drive. The SWBATT voltage is stored in C9, enabling energy transfer to occur with minimum battery voltage dissipation. Transistor Q4 and SCR CR4 control this voltage and direct it to the Transfer Relay. To enable the relay drive, the signal from U14C turns on CR4. The signal from U11A forces U6A low, causing Q4 to conduct. With both CR4 and Q4 conducting, current starts to flow through R62. To maintain a constant current to the Transfer Relay when current through R62 increases, the noninverting input of U6A goes higher than the reference voltage at the inverting input and the output switches high. This turns off Q4, causing the current to decrease through R62 and switching U6A output low.

The high noise levels involved in the transfer of stored defibrillation energy can jeopardize the function of many types of logic signals, particularly if any are in transition from one state to another. To overcome this effect, the System and Display/Recorder/Pacing microprocessors and other high speed logic components are rendered inactive during transfer. To accomplish this, U6B discharges C35 at the input of comparator U8D and holds U8D low during the time that the Transfer Relay coil is energized. The output of the U8D remains low for an additional period to allow the Transfer Relay to reopen and any stray current to dissipate. This signal, Sleep Inhibit, is used as a logic activity enable.

Sleep Inhibit also clears the shift registers, U19 and U20. This sets all inputs to U17B at 0V, which inhibits the drive to the Energy Storage Capacitor Charger. Inadvertent deactivation of the charge inhibit will then cause no discharge of stored energy from the energy storage capacitor.

**Energy Dump** 

The Energy Dump circuit provides a safe current path for unused energy stored in the energy storage capacitor. If the paddles are charged but not discharged within one minute, the Energy Dump circuit is activated.

The internal dissipation of stored energy is accomplished with the dump relay and a power resistor located within the Transfer Relay Assembly. The dump relay is held in the normally closed condition by a permanent magnet. When a charge is initiated, the contact opens to prevent the stored energy from draining.

The System microprocessor controls Q3, preventing the paddles from being charged for longer than one minute. Prior to any charge activity, the System microprocessor holds the Dump Inhibit line low, keeping Q3 turned off and the dump relay closed. During charging, the Dump Inhibit line is switched high, turning on Q3 and opening the dump relay. If the one-minute charge time-out expires or a different energy setting is chosen, the System microprocessor turns Q3 off again to close the dump relay and dissipate the stored energy.

**Energy Monitor** 

The Energy Monitor circuit confirms the level of energy that is discharged into the internal test load. When an energy transfer is performed as a test, this circuit senses the energy delivered to the test load and informs the System microprocessor. The System microprocessor then displays the information (e.g., 200J DELIVERED).

The value of the defibrillation pulse is acquired by monitoring the delivered energy with a current transformer, T3. A proportional current is generated in the secondary with each transfer pulse. A rectifier bridge, CR8 through CR10, delivers the secondary current to C10. The voltage is indexed in the System microprocessor to a specific delivered energy.

Capacitor C10 must start each transfer cycle with 0V for accurate evaluation of the delivery. FET Q8 is in parallel with C10 and resets C10 before each charge energy pulse. The gate of Q8 is driven by the inverted Charge Inhibit signal. The FET Q8 shorts C10 when the charge is enabled which makes sure that each time the instrument is ready to deliver energy, the Energy Monitor circuit determines the correct energy level. Op amp U1A functions as a noninverting buffer, providing high-impedance isolation from the rest of the signal processing circuitry.

Serial Data Interface

The Serial Data Interface consists of four lines; Converter Data In and Converter Data Out allow communication while Converter Function Select and I/O Clock control data flow. The Analog Voltage Evaluation circuit and the D/A converter U18, used with the Capacitor Charge Control circuit, share this common set of serial lines for communication with the System microprocessor.

A low on the Converter Function Select line at connector P14 pin 7A enables the Analog Voltage Evaluation circuit for data transfer. The data flow rate is set by I/O Clock at P14 pin 6A. A high on the Converter Function Select line at P14 enables 8-bit shift registers U19 and U20, allowing I/O Clock access to the register clock inputs. Information on the Converter Data Out line is stored in the low-order location, and all other bits are advanced one position during the positive-going transitions of I/O Clock. At the same time, U12A through U12C route the output data on the Converter Data In line. A logic high at U11C pin 8 following a byte transfer indicates a communication failure.

Analog Voltage Evaluation

The Analog Voltage Evaluation circuit monitors operating voltages on the Power Conversion PCB. An A/D converter, U9, converts these voltages to digital code for evaluation for correct tolerances by the System microprocessor.

Charge LED Drive

The Charge LED Drive signal from the System microprocessor controls Q14 when the energy storage capacitor is charging. Transistor Q14 drives the front panel CHARGE LED and the paddle CHARGE LED. Resistor R128 on the drain of Q14 limits the current to the Paddle LED. The gate of Q13 is connected to the Charge Valid signal. When the voltage for the selected energy level has been reached, Q13 conducts. Transistor Q14 then clamps the charge LEDs on continuously.

Interconnect Pacing **PCB Assembly** (805472)

The Interconnect Pacing PCB provides the pacing circuitry as well as the electrical connection between the Main PCB and the Power Conversion PCB. Refer to Figure 1–4, page 1–18.

The pacer circuit is an isolated closed-loop programmed current source capable of supplying up to 200mA into a 1.4k $\Omega$  load, or 120mA into a 2.4k $\Omega$ load. When triggered, it automatically applies a 20±1ms pace pulse to the patient load. The pacer circuit returns an echo signal which encodes the current actually applied to the patient so microprocessor U1 on the Main PCB can determine the output current level.

Interface

Control for the pacer circuit operation and monitoring comes from the Display/Recorder/Pacing Interface circuitry on the Main PCB. These controls are: PACER POWER, which turns the pacer on and off, and PACE COMMAND, a single low-going pulse which encodes the pace current amplitude. When the microprocessor (U1, Main PCB) applies a pulse, a pulse is applied to the patient. Signal IP FEEDBACK to the microprocessor remains low when PACER POWER is off and goes high after PACER POWER goes on. LEADS OFF will be low if the impedance between the pacer leads exceeds about  $6k\Omega$ 




Figure 1–4 Interconnect Pacing PCB Block Diagram

Pace Pulse Output Drive

The Pace Pulse Output Drive circuit consists of a high-voltage switch (Q2, Q6, CR5, and R22-24), a current-limiter (Q1, U1, and R4), and a current sink (U5A, O3, R49 current sense resistor).

On the leading edge of PACE TRIGGER the Arming Circuit one-shot biases Q6 on for 25ms. Transistor Q6 biases Q2 on through the voltage drop developed across R23 and R24 (in series to avoid exceeding their maximum voltage ratings). R22 normally holds Q2 off. Diode CR5 protects the gate of Q2 from excessive drive.

The amount of current allowed to flow from the 320V supply into the patient is limited by an independent regulator to less than 227mA. Reference U1 and transistor Q1 prevent the voltage drop developed across R4 from exceeding 2.5V. Q1 remains saturated until the voltage developed by the current flow through R4 exceeds 2.5V. Should more current attempt to flow, Q1 comes out of saturation limiting the total amount of available current by decreasing the gate-source voltage on Q2. This circuit helps protect the patient from excessive current.

Current allowed through the patient is determined by the current sink formed by U5A and Q3. A voltage waveform representing the pace current waveform is applied to U5A noninverting input. U5A controls Q3 so a duplicate of the input voltage waveform also develops across R49, the current sensing resistor. The current required to develop the duplicate voltage can only come from flowing through the patient.

Pace Pulse Wave Shaper

A voltage level proportional to the desired pace pulse amplitude is developed at U2 pin 1. Transistor O7, which is normally on, shunts this voltage to ground through R8. On the trailing edge of PACE TRIGGER a 20ms one-shot is triggered. This one-shot abruptly turns Q7 off. The resulting voltage on the drain of Q7 is rise-time-limited by C3 and R12. The amplitude is determined by the voltage level received at U2 pin 1. This voltage pulse is coupled through high pass filter C2-R7 (t = 90ms) to the current sink circuit described above.

When the 20ms one-shot times out, Q7 turns on and the voltage on the drain immediately drops to ground. This transition is also coupled through C2. Diode CR1 clamps the voltage presented to the current sink to ground, preventing the current-sink loop from opening.

Pace Pulse Amplitude Decoder

Transistor Q4 is an analog switch. When closed, this switch causes the circuit to function as a unity-gain, two-stage, inverting amplifier. The input is the grounded end of R29 and the output is U2 pin 1. The inverting input to U2B is held to 0V through R29. This forces the overall circuit to bias itself through the closed loop so the output of U2A is also 0V. The noninverting

input to U2A is biased to the +5V reference. This forces the overall circuit to bias itself through the closed loop so the inverting input to U2A is also biased at +5V.

With switch Q4 closed, the closed-loop forces U2 pin 5 and U2 pin 2 to be biased at 0V and +5V, respectively. It does not matter whether switch Q8 is open or closed since current drawn through R11 and Q8 can be sourced by U2B through Q4.

Two control waveforms, one at CR4 and the other at CR6, are derived from the falling and rising edges of PACE TRIGGER. When either of these control waveforms are low, switch Q4 is open. When the waveform at CR6 is high, Q8 closes. With Q4 open, the feedback loop around U2A forces the inverting input of U2A to 5V, placing a constant 5V across R11. The only source of current for R11 is through C13. The voltage present on U2 pin 1 must, therefore, be a linear ramp. U2A will maintain this ramp until it either runs out of power supply voltage or until Q8 opens.

Transistor Q8 opens on the trailing edge of PACE TRIGGER. The ramp, therefore, starts from 0V, continues for the duration of PACE TRIGGER, and ends up at a voltage that depends on the low-going width of PACE TRIGGER. With switch Q8 open, the overall circuit functions as a track-and-hold amplifier operating in the hold mode. Since there are no longer any current paths into U2 pin 2, the voltage on U2 pin 1 holds constant.

This output is monitored by comparator U4D. If the ramp exceeds a pre-determined value (approximately 10.8V), this comparator will trip, opening the high voltage switch Q2 and disabling pacemaker output to the patient load. This over-amplitude detector guards against microprocessor upsets which could apply excessive current to the patient. Also, a pace pulse application can be aborted after it is initiated.

Echo Pulse Width Converter

The voltage developed across R49 is a direct indication of current flowing through the patient. This voltage is converted to a pulse width and is returned to the microprocessor as IP FEEDBACK indicating the current actually flowing into the patient.

The Echo Pulse Width Converter circuitry converts an analog voltage to a pulse width. When analog switch Q5 is closed, amplifiers U2C and U2D create a unity gain, inverting, two-stage amplifier. The comparator circuit connected to U2C requires a positive input signal to function properly. Since U2C and U2D are unity gain and inverting, the input must therefore be between 0V and minus the positive power supply rail. Amplifier U5B adds a -13V offset to the voltage developed at R49, multiplied by the gain of U5B. The voltage at the output of U2C will therefore be +13.0V when 0V appears

at R49 and switch Q5 is closed. With switch Q5 closed, current through R21 goes into U2D and the voltage at U2 pin 8 tracks the voltage on R49. On the trailing edge of the pace pulse, Q5 opens and the current through R21 must flow through C5 instead of U2D. This ramps the voltage at U2 pin 1 downward at a constant rate to ground. When it equals ground, comparator U3A trips Q5 closed again. The time required to reach ground is determined by the voltage value just prior to Q5 opening. This pulse width communicates to the microprocessor as IP FEEDBACK. The Echo Pulse Width Converter operates similar to the Preamp Echo circuit (page 1–29).

Pacing Leads Off Detector

Pacing leads status is determined by driving an ac current through the pacing leads. This current develops a voltage across the impedance presented by the patient. The voltage is amplified and if it exceeds a predetermined amount, a comparator is tripped and drives the LEADS OFF optocoupler alerting the microprocessor.

U11D forms a 1kHz oscillator that swings between the positive and negative supply rails supplying the ac carrier voltage. Frequency of operation depends on the supply voltages, R75, and C20. Capacitor C25 couples the ac carrier voltage into resistors R80 and R87 which have impedances high enough to effectively convert the oscillator voltage to current. Capacitors C21 and C22 couple the voltage developed across the patient into diode clamps CR11-CR14. Resistors R81-83 limit the current to the diode clamps.

U11C amplifies the voltage developed across the leads. Eye-close switch Q15 shunts the output of U11C during the 25ms of the PACER ARMED signal (when PACE COMMAND switches the full voltage through Q2). Transistor Q15 conducts as long as the PACER ARMED signal is present. A filter network turns Q15 on quickly, then gradually turns it off which re-enables the signal. The band pass filter is centered around the oscillator frequency and consists of U11B, R86, C23, and C24.

The output of comparator U11A drives optocoupler U12. It compares the peak-detected amplitude of the leads-off carrier to a preset voltage. When the leads impedance exceeds a predetermined amount, U11A trips and alerts the microprocessor.

**Arming Circuit** 

The length of time the pacer is armed (25ms) is determined by one-shot monostable comparator U4C. The one-shot is triggered on the falling edge of PACE TRIGGER. Resistors R31 and R32 bias the inverting input of U4C to +5V, and R33 and R34 bias the noninverting input to 0V. Comparator U4C pin 14 remains at the negative supply rail since its inverting input is 5V higher than its noninverting input. The negative transition of PACE TRIGGER is coupled through C8 to U4C and U4C pin 14 becomes high-impedance. R28 then pulls the output to the positive supply rail. The

time constant formulated by R33, R34, R28, and C9 requires 25ms to discharge to +5V.

Pace Pulse Length Generator The width of the pacer pulse (20ms) is determined by one-shot monostable comparator U3D. The one-shot is triggered on the rising edge of PACE TRIGGER. Resistors R51 and R52 bias the noninverting input of U3D to +1.15V and R37 biases the inverting input to 0V. Comparator U3D pin 13 remains high impedance since its noninverting input is 1.15V higher than its inverting input. Resistor R35 pulls U3D pin 13 to the positive supply rail. The rising transition of PACE TRIGGER is coupled through C10 and U3D, and U3D pin 13 drops to the negative supply rail. Since the negative supply rail is unregulated and the width of the output determines the width of the pace pulse, a stable voltage-supplied reference diode, CR21 (10V), is needed to drive the time constant. Resistor R35 pulls U3D pin 13 to the positive supply rail since ground is blocked by CR22. The time constant formulated by R51, R52, and C11 discharges from -10V (CR21 reference voltage) to 1.15V (derived by R51 and R52). It requires 20ms for the waveform to charge through ground where the circuit resets.

Isolated Power Supply

The Isolated Power Supply is a switching flyback supply regulated on the isolated +15V. It supplies +320V,  $\pm15$ V, and +5V through linear regulator U8. Electrostatic discharge protection is provided by the resistors and capacitors on all the I/O lines that run off the Interconnect PCB.

**Switching Regulator** 

Components U6, Q9, and CR8 form a simple flyback Switching Regulator operating at 17.5kHz. Components R60 and C16 determine the frequency of operation. Transformer T1 is a multi-winding inductor that transfers energy stored in its primary winding to the secondary windings while Q9 is open. Primary winding current is determined by monitoring the voltage developed across R56. The duty cycle of U6 varies to control the energy applied to T1 primary windings. The secondary voltages are proportional to the energy stored.

Power Supply Controller

The Isolated +15V supply is compared to the +5V reference by comparator U4B. When isolated +15V exceeds +15.25V, U4 causes Q11 to conduct which causes the transistor side of U7 to conduct. Component U6 shuts down until the isolated +15V supply has dropped back to 15V. The entire circuit can be shut down by the discrete signal PACER ON. When Q10 conducts, U6 is inhibited. When the microprocessor turns Q10 off, the power supply voltages stabilize within approximately 500ms.

**LED Drive** 

The LED Drive circuitry consists of transistors Q204, Q205, Q213, Q214, and resistors R214, R215, R211, and R225. Static protection is provided by CR201 through CR213.


## Main PCB Assembly (805474)

The Main PCB Assembly, illustrated in Figure 1–5, consists of three major circuits: the Preamp, the System Controller, and the Display/Recorder/ Pacing Interface. The Preamp circuitry converts the ECG from the patient cable to sampled ECG. The converted ECG signal is optically connected to the System microprocessor which performs the initial ECG signal processing. The System microprocessor is central to the System Controller circuitry. The System Controller provides the interface between the operator inputs and the instrument functions. The Display/Recorder/Pacing Interface controls the CRT and recorder displays and the pacing functions according to the information supplied by the System microprocessor.

The three major circuits control the functions listed in Table 1–2. Each circuit is described in the following paragraphs.

Table 1–2 Main PCB Functional Circuits

Preamp	System Controller	Display/R	ecorder/Pacing I	nterface
Isolated Power Supply Front-end Protection Lead Select Optical Interface Instruction Decoding Filtering DC Restore Echo	System Controller Front and Rear Panel Controls Interface Charge Control I/O Protection Pacing Current Adjustment	Display:  Display/ Recorder Interface System Clock Address Decode Display Gate Array Video RAM Font EPROM	Program EPROM Program RAM Recorder Power Motor Speed Control Recorder Gate Array Paper Out Sense Printhead Resistance Evaluation	Pacemaker:  Leads Off Detection  Pacing Current Delivery  ECG and Tone Output



Preamp

**Isolated Power Supply.** The Isolated Power Supply provides the isolated +15Vdc, +12Vdc, -15Vdc, +5Vdc, and -12Vdc supplies required by the Preamp circuitry. The drive for transformer T1 is generated by pulse-width modulator U5. Capacitor C12 and R28 set the frequency between 12kHz and 24kHz. The secondary of T1 is applied to full-wave bridge CR20. The positive half of the ac input waveforms are stored at C87, forming the isolated +15V supply. The negative ac waveforms are stored at C86, forming the isolated -15V supply. Regulators U38 and U39 isolate more sensitive circuitry from the heavier loads. Since the +5Vdc loads are synchronous and small, U37 doubles as a regulated supply and a reference source.

Front-End Protection. To protect against voltage surges, VSP1 through VSP3 limit the input current to 2mA. Isolation for the paddles input is provided by R131 and R132 which guard against excessive charge input to the Preamp. Diodes CR18, CR19, and CR28 clamp the differential paddle voltage when K1 is switched to the leads mode.

**Lead Select.** To select between paddles and the different patient leads, latch U33 is loaded with the decoded data from counter U32. Flip-flop U33 controls relay K1, U43A, and analog switches U42 and U41. To select between patient leads or paddles, U33 causes Q8 or Q9 to conduct, switching the state of K1 to the desired position. Patient leads are selected when Q8 conducts; paddles are selected when Q9 conducts. If patient leads are selected, analog switch U42 selects the lead (I,II,III, or STD) requested by the output of U33.

Optical Interface. Command pulses from the System microprocessor pass through the optical interface to perform three different functions: synchronize the ECG conversion, select 1 of 5 patient cable leads, and induce a +1mV calibration signal. The System microprocessor uses an internal timer to transmit a pulse of variable duration (approximately 80ms to 1.2ms) and a fixed-frequency period (540Hz). The variable-duration pulse from the System microprocessor enables U4 to start counting, creating a 16kHz pulse train at U4 pin 4.

The System microprocessor controls the count of the pulse train by interrupting the output with a CLR signal as shown at TP6 in Figure 1–6. The 16kHz pulse drives the infrared-emitting diode in U7 providing pulse isolation. Comparator U34D performs waveform conditioning for the pulse train as it passes to the Instruction Decoding circuit, U32 and U33.




Figure 1-6 Preamp Timing

Instruction Decoding. The pulse train commands from the System microprocessor are coded by the number of pulses as listed in Table 1–3. The number of pulses represents a binary number which corresponds to a particular command function. Each set of pulses performs two tasks: application of a specific preamp setup, and initiation of ECG signal conversion to a digital equivalent. The rising edge of the first pulse is detected by U31A which responds by sending a pulse to U34A. This initiates the ECG conversion. Monostable multivibrator U31B detects the falling edge of the first pulse and clears U31A so that ECG conversion is not initiated with every pulse on the pulse train.

Flip-flop U36A and counter U32 form a five-stage counter that responds to the pulses from the pulse train (see TP17, Figure 1–6). The first rising edge of the pulse train causes the five-stage counter to load from the preset inputs (pins 6, 5, 3, and 4). Subsequent edges count up from the preset inputs. Because of the timing of RC network R101 and C74, U31B times out at about 1.5ms and sends a positive edge to U33, latching the decoded command signals. These command signals select between leads and paddles.

Table 1–3 Preamp Pulse Widths

Pulses at TP17	Pulse Duration At TP6 (ms)	Result on Preamp Setup
1	64	Not used
2	128	Not used
3	192	Not used
4	256	Lead III
5	320	Lead II
6	384	Lead I
7	448	Std
8	512	Lead III + Cal
9	576	Lead II + Cal
10	640	Lead I + Cal
11	704	Std + Cal
12	768	Lead III + Leads Coil Drive and Switch Smooth
13	832	Lead II + Leads Coil Drive and Switch Smooth
14	896	Lead I + Leads Coil Drive and Switch Smooth
15	960	Std + Leads Coil Drive and Switch Smooth
16	1024	Lead III + Cal + Leads Coil Drive and Switch Smooth
17	1088	Lead II + Cal + Leads Coil Drive and Switch Smooth
18	1152	Lead I + Cal + Leads Coil Drive and Switch Smooth
19	1216	Std + Cal + Leads Coil Drive and Switch Smooth
20	1280	Lead III + Paddles Coil Drive and Switch Smooth
21	1344	Not used
22	1408	Not used

Instrumentation Amplifier. The typical ECG signal at U41 and U42 is no more than 16mV p-p and requires amplification for further signal processing. U44A and U44B form a differential input amplifier with the gain set by RN2 pins 1 and 2, R144, and RN2 pins 3 and 4. This differential stage is followed by a switch network (U48) used to disconnect the other amplifier stages during pacing pulse delivery. EYECLOSE, received from U1 (so pacer-induced spikes are not received as R-waves), opens then closes switch U48 (lasting 40 ±0.1ms) accomplishing pacer rejection. C117 holds the last voltage value so when EYECLOSE terminates and closes U48, the ECG signal may again be sent through U55A and U55B (unity gain buffers). Therefore, the display baseline remains unchanged.

Op amp U44D amplifies the differential signal and rejects a common-mode signal from U44A and U44B. This common-mode signal (mainly 60Hz) is generated from the outputs of U44A and U44B through voltage dividers

R148 and R149. Next, the common-mode signal is inverted and amplified by U35D. This signal then passes through analog switch U41 and relay K1, which routes it to the unused patient lead (LL) cancelling common-mode signals at the patient source.

The +1mV calibration signal is created by an instruction from the system microprocessor which sets U33 pin 7 high. When U43A closes, a -5V offset signal is applied to U44A pin 2 and U44B pin 6. This adds about a +1mV differential offset to the input without changing the gain of the input stages.

Filtering. Filtering provides a uniform signal bandwidth during normal operation. The preamp filtering is performed by two stages centered around op amps U44C, U45D, U45A, and U45B. In the first stage, U44C and U45D function as a high-pass filter for Preamp, DC Restore, and switch smoothing, blocking frequencies below 0.024Hz. In the second stage, low-pass filter U45A removes frequencies above 200Hz. U45B functions as a 540Hz notch filter.

**DC Restore.** The DC Restore circuit compensates for any unwanted dc offset originating from switch smoothing, defibrillation, or motion artifact. This preserves the dc offset level required by the Echo circuit. DC Restore controls any transient signals to prevent the high potential gain of the Instrumentation Amplifier from falsely amplifying the dc level input for conversion. Under normal operation, the A/D conversion dc level at TP19 is from +1V to +11V.

Switch smoothing is the process of limiting the response of the preamp to do offset voltages when switching from one lead (I, II, III, or STD) to another. This function is initiated when the appropriate number of pulses (see Table 1–3, pulses 12 thru 20) turns on Q8 or Q9 and switches U42 to the selected lead. When Q8 switches on, C90 discharges from +12V, causing Q11 to conduct. Q11 changes the high-pass filter by switching R161 in parallel with R165. Because of R161, this filter limits the signal passband to above 40Hz.

Transient signals associated with defibrillation are sensed by op amp U45C. U45C functions as a window comparator and responds with a positive output when the voltage at its input exceeds  $\pm 2$ V. Resistors R172 and R147 apply a negative bias to the noninverting input of U45C; R171 and R143 apply a positive bias to the inverting input. C104 and R151 function as a low-pass filter to prevent very short spikes (such an an internal pacing spike) from inducing DC Restore. The output of U45C drives transistor Q10, routing the signal through high-pass filter U44C and U45D. Resistor R160, with U44C and U45D, limits the signal passband to 5Hz.

When **PADDLES** is selected, further passband limiting is required to eliminate motion artifact introduced by the paddles when used in the QUIK-LOOK mode. A command pulse selects paddles mode, (see Table 1–3, pulse 20), turning Q9 on. When Q9 conducts, flip-flop U36B is loaded so that the next clock pulse at pin 6 switches the Q output high. Analog switch U43B then closes, bringing R155 in parallel with R165 in the high-pass filter, and limiting the ECG signal passband to above 2Hz.

Echo. The ECG signal is converted to a pulse-delay period on command from the System microprocessor. The time from when the convert command from the System microprocessor is received until the "echo" pulse is sent back to the System microprocessor is proportional to a sample of the ECG signal (see Figure 1–7). When the echo pulse is sent back to the System microprocessor (across an opto-isolator), the System microprocessor converts the pulse-delay period to a number.




Figure 1-7 Echo Period

ECG conversion begins when the pulse train at the output of counter U4 passes to U32 pin 2. The rising edge of the first pulse clocks U31A, which sends a trigger pulse every 1.8ms to U34A. The ECG signal is offset to -6.3V at the output of U35A, providing a signal range of -1V to -11V for conversion.

The ECG signal from U35A passes to a track and ramp circuit consisting of U35B, U35C, Q7, and C76. This circuit controls the duration of the echo period. The track phase begins when Q7 conducts, which causes the ECG signal at U35B to charge sampling capacitor C76. The ramp phase starts when comparator U34A receives the trigger pulse (every 1.8ms), driving the output to -12V. This voltage passes to the gate of Q7 turning Q7 off. The +5V line and R122 provide a constant current source to the inverting input of U35C. When Q7 turns off, the constant current source causes a linear discharge from C76. This discharge produces a ramp output from U35C (see Figure 1–7).

While the track and ramp circuit is tracking, comparator U34A receives the trigger pulse from U31A. This drives the output of U34A to -12V, turns off Q6, and allows the +15V supply to charge C72 through voltage dividers R100 and R99. The output of comparator U34C controls the on-time of Q6, which discharges C72 through opto-isolator U6 and creates the echo signal.

The duration of the ramp output of the track and ramp circuit determines how long the output of U34C is low and sets the duration of the echo period. When the ramp output of U35C reaches 0V, the output of comparator U34B switches to -12V, presenting this voltage at the inverting input of U34C. This causes the output of U34C to switch to high impedance which turns on Q6 and allows the echo signal to be coupled through opto-isolator U6.

System Controller

The System Controller circuit receives input signals from the front panel keypad and sends commands to the CRT, recorder, pacing, and defibrillation circuits. The central component of this circuit is System microprocessor U30, together with external RAM and ROM devices. The System microprocessor converts the ECG echo pulse widths to sample values for the CRT and recorder. These ECG samples are then sent to the Display/Recorder/Pacing microprocessor for filtering and display on the CRT.

U30 communicates with the Display/Recorder/Pacing microprocessor (U1) through two asynchronous serial transmit and receive lines at pins 12 and 13. In addition to the three 8-bit I/O ports contained in U30, 16 extra output lines are created by routing D0 through D7 in parallel out to octal latches U19 and U20.

The operating software for the System Controller is contained in U27, a 64k x 8 EPROM. Battery-backed RAM U28A provides data storage, including the calibration constants. DIP (Dual In-line Package) socket U28B contains the lithium energy cell and real-time clock information. Address decoding is performed by 3-to-8 decoder U25 using only the outputs at pins 12, 13, and 14.

Multi-vibrator U26 and flip-flop U16 form the watchdog circuit for the System microprocessor and the Display/Recorder/Pacing microprocessor. The watchdog circuit monitors the 1.8ms preamp trigger pulse and resets both microprocessors if the pulse width is out of tolerance. An RC network consisting of R46 and C31 is used by U26 to detect a pulse width that is too short; R47 and C33 detect a pulse width that is too long.

Controller U46 monitors the +5V supply and causes a system reset as needed. This protects the RAM memory and prevents spurious microprocessor activity. If the controller detects that the +5V supply is out of tolerance, the output at NOT CEO switches high and causes Q5 to conduct, pulling the voltage at R45 low. This low logic level (through U23C, U18A, and U17A) goes to the RESET lines for the System microprocessor and the Display/Recorder/Pacing microprocessor.

Front and Rear Panel Controls Interface. Signals from the front and rear panels are distributed to the System microprocessor which determines what function(s) to activate. Input from the pushbuttons on the front and rear panels are arranged in a matrix of five columns (CAS 1 through 5) and three rows (RAS 1 through 3). The System microprocessor functions as a multiplexer, providing a high logic level for each row and sensing any closed column switch.

Charge Control. The System microprocessor initiates and monitors charge requests to the D/A converter U18 on the Power Conversion PCB. The I/O Clock, Converter Function Select, Converter Data In, and Converter Data Out lines accomplish serial communication with the D/A converter. Pressing CHARGE on the keypad transmits the coded selected energy value from the System microprocessor to the Power Conversion PCB over the Converter Data Out line. Status information on the decoding process on the Power Conversion PCB is sent back to the System microprocessor over Converter Data In. The I/O Clock provides synchronous timing for serial communication.

Following the charge initiation, the System microprocessor performs a set of self-tests which are completed before the selected energy is available. If a self-test fails, the CRT displays the **SERVICE** indicator.

The energy stored in the energy storage capacitor is monitored through a 10-bit A/D converter (U9) located on the Power Conversion PCB. The A/D converter shares the four-line, serial-communication interface with the D/A converter. If the System microprocessor indicates that the energy stored is out of tolerance, an ENERGY FAULT message is displayed.

I/O Protection. Electrostatic discharge protection is provided by the resistors and capacitors on all the I/O lines that connect to the Main PCB. Additional protection is provided by 6.2V zener diodes, CR7 through CR9, on the Video 1, Video 2, and Raster lines directed to the Monitor PCB.

Pacing Current Adjustment. A rotary pulse generator (RPG) adjustment knob on the front panel generates TTL level pulses that the Interconnect PCB sends to flip-flop U54 in digital pulses. System microprocessor U30 can access the knob by reading or writing to a particular address. A read command will put the flip-flop data on the data bus. If U54 pin 4 is high, the knob has been rotated since the last read. If U54 pin 3 is high, it has been rotated clockwise (increased current). A write command will clear both flip-flops and the write data is ignored.

Display/Recorder/Pacing Interface

Display/Recorder/Pacing Controller. The central component of this circuit is the Display/Recorder/Pacing microprocessor U1. It receives the digitized ECG signal data from the System microprocessor and controls the input of this data to the Display PCB Assembly CRT and the recorder. Because the CRT screen resolution is not as high as the recorder resolution, the Display/Recorder/Pacing microprocessor compresses ECG data for the CRT display. The recorder receives all incoming ECG data samples. The Display/Recorder/Pacing microprocessor also formats the ECG data for the CRT so waveforms are not displayed with gaps or overlapping segments.

ECG sample data is transmitted from the System microprocessor to the Display/Recorder/Pacing microprocessor through the Tx serial communication line approximately every 1.8ms. Other data/command bytes from the System microprocessor are also transmitted to the Display/Recorder/Pacing microprocessor when they are required for a particular function. This data is transmitted in packets containing a header byte, two ECG sample bytes, and, if required, a data/command byte.

The first four I/O lines listed in Table 1–4 communicate with the Display Gate Array, U29, driving the input for the Display PCB. The other five I/O lines control the operation of the recorder.

Table 1-4 Display/Recorder/Pacing Microprocessor I/O Lines

I/O Line	Signal	Direction
P(50)	H.RETR	In
P(52)	GO	Out
P(53)	INCMEM	Out
P(54)	TR/CHAR	Out
P(21)	STROBE	Out
P(60)	RECORDER POWER	Out
P(61)	PRINTHEAD CURRENT SENSE	In
P(66)	PRINTHEAD TYPE	In
P(67)	PRINTHEAD LOAD	Out

**System Clock.** The crystal-controlled oscillator U50 supplies the 12MHz system clock to both the System and Display/Recorder/Pacing microprocessors (U30 and U1). Oscillator U10 provides a 16MHz square wave to the Display Gate Array and a 2MHz signal to A/D converter U9 on the Power Conversion PCB.

Address Decode. Programmable logic array U9 controls the reading from and the writing to all external memory and the Recorder Gate Array, U3. U9 functions as an address decoder and is programmed to enable the appropriate devices when the Display/Recorder/Pacing microprocessor addresses different areas in memory.

Buffers U11, U12, and U15 isolate the appropriate RAM/ROM memory devices from interaction with the Display/Recorder/Pacing microprocessor. U11 buffers the data bus. Buffers U12 and U15 control interaction on the address bus between the Display/Recorder/Pacing microprocessor, the Display Gate Array, and the RAM/ROM memory devices.

Display Gate Array. Display Gate Array U29 is a custom IC containing many blocks of standard OR and AND gate circuits (counters, multivibrators, flip-flops, etc.). A metalized layer added on top of the gate circuits connects only those circuits needed for a specific function: Monitor PCB sync timing, ECG memory addressing, ECG display trace generation, video intensity, and additional timing signal generation. Two different internal counters control the timing of the Video RAM (U21) and Font EPROM (U22) addressing during CRT scan periods, and control the Video RAM addressing during CRT retrace periods.

The Display Gate Array uses an 18-bit counter to control the timing of the vertical and horizontal retrace on the CRT. The lower 8 bits correspond to 256 rows, and the upper 10 bits correspond to 768 columns (scan lines). When a count of 256 occurs, U29 initiates a vertical retrace. When a count of 196,608 (256 x 768) occurs, the gate array initiates the horizontal retrace.

Four separate memory loops for two separate trace generators in the Display Gate Array format the ECG signal in memory. Formatting occurs during horizontal retrace. An internal 11-bit counter addresses each memory loop. Each trace generator uses one memory loop for the high point of the scan and another memory loop for the low point of the scan. The trace appearing on the CRT screen is created by turning the video on at the screen location that corresponds to the low point, and turning the video off at the screen location that corresponds to the high point. This saves memory space in Video RAM U21 since the only data loaded in memory is the high point and the low point addresses for each scan.

Serial communication between the Display Gate Array and the Display/ Recorder microprocessor is accomplished by the following control lines: H.RETR, GO, INCMEM, and TR/CHAR. These lines enable the Display Gate Array and the Display/Recorder/Pacing microprocessor to coordinate the storage and retrieval of ECG and character data in the Video RAM during horizontal retrace. Figure 1–8 shows the control line timing sequence.




Figure 1-8 Control Line Timing

Data from the Display/Recorder/Pacing microprocessor loads into the Video RAM memory during horizontal retrace. Six ECG samples update during horizontal retrace to smooth scrolling of the trace displayed on the CRT. Coded character address data used by the Font EPROM also loads into Video RAM at this time.

ECG trace and messages display during raster scan. The Video RAM sends ECG data to the Display Gate Array through EDAT(0-7); character address information passes to the Font EPROM through DD(0-7). Character pattern data stored in the Font EPROM is addressed and routed to the Display Gate Array as CDAT(0-7). The Display Gate Array sends the character data out to the proper address on the Display PCB (see Figure 1–9).




Figure 1-9 Character Cell Locations

The Display Gate Array also functions as a video intensity circuit, producing four levels of brightness for the trace and characters. Two video outputs from the Display Gate Array, Video 1 and Video 2, generate the appropriate brightness level for the digitized display data by transmitting the data in a two digit code to the Video Signal Amplifier circuit on the Monitor PCB. This code determines if the addressed pixel will be at full intensity, normal video, half-tone, or off (see Figure 1-10). Half-tone is used only for the endpoints of the ECG waveform segments; full brightness is used only for the sync marks. The rest of the ECG waveform and all the characters are transmitted with normal brightness.




Figure 1-10 Video Intensity

Character dot information merges with trace information to create the video input signal. Messages have priority over waveforms when the two types of data are mapped to the same point on the CRT.

Video RAM. Trace and character data are stored in a 32k x 8 static RAM, U21. New ECG data is stored during the horizontal retrace and consists of twelve bytes. Horizontal retrace provides an interrupt to the Display/Recorder/Pacing microprocessor which signals permission for a memory access. New data storage and memory address counter incrementing is controlled by the Display/Recorder/Pacing microprocessor.

An additional 512 bytes are used for storage of the character address codes (similar to ASCII codes). The screen can be configured into a matrix of 16 x 32 character cells, with each cell accepting character information from the Font EPROM. The message content is updated as needed by the Display/Recorder/Pacing microprocessor during horizontal retrace.

Font EPROM. A 32k x 8 EPROM, U22, contains character font patterns, with memory for 256 characters. Each character is 16 pixels high and 24 pixels wide to fit in the screen character cells (see Figure 1-9). The displayable character set includes single-sized numerals, punctuation, and

uppercase letters in normal, inverse, and boldface video. Special characters and double-sized numerals are in normal video display.

**Program EPROM.** A 64k x 8 EPROM, U14, contains the operating software for the Display/Recorder/Pacing microprocessor. Included in the software is a list of codes representing all of the words and phrases used in the annotated recorder text.

**Program RAM.** The Program RAM U13 is used by the Display/Recorder/Pacing microprocessor as a scratch-pad while performing program functions. Additionally, all CODE SUMMARY data and the Display/Recorder/Pacing microprocessor bulletin board is stored in the Program RAM.

**Recorder Power.** The Recorder Power circuit provides operating voltage for the Printhead and the Motor Speed Control circuits by responding to a high logic signal from the Display/Recorder/Pacing microprocessor at the base of Q2. This logic signal causes Q2 to conduct which turns on Q4, thereby supplying power for the recorder motor and printhead.

**Motor Speed Control.** Regulation of the recorder speed is necessary because of the variable drag on the recorder motor caused by the thermal printhead and the lack of a spindle on the paper roll. The paper speed of the recorder is regulated by a variable-frequency, pulse-width-modulation circuit that compares the back EMF (electromotive force) on the motor to a reference voltage. The circuit varies the frequency of the motor speed drive pulse enough to maintain the 25mm/s paper speed.

The back EMF produced in the motor when the drive pulse switches off letermines the speed of the recorder motor. The voltage present at voltage dividers R18, R17, and C4 is directly proportional to the speed of the motor when Q3 is off. A frequency-modulated pulse from transistor Q3 drives the recorder motor. The pulse has a fixed on-time and a variable off-time. The off-time is set by comparator U2B and capacitor C4 and determines how much power Q3 conducts to the motor.

The fixed on-time is determined by U2C. When the off-time period is over, C2 has charged to a voltage level approximately 75% of VBATT. At the next on-time cycle, the output of U2B switches to ground, the voltage at U2C pin 10 goes to about 1/3 of VBATT, and C2 starts discharging through R7. U2C changes states and ends the pulse on-time when the voltage on C2 drops below the voltage on pin 10.

The voltage on U2B pin 7 is compared to a calibrated reference voltage at pin 6 to determine the off-time. During the on-time, C4 charges to a voltage proportional to the battery voltage. During the off-time, C4 discharges toward a voltage that is proportional to the back EMF of the motor. Under

normal operating conditions, the voltage on C4 discharges over several time constants before the voltage on U2B pin 7 is equal to the reference voltage on pin 6. When the voltages are equal, the output of U2B switches low and Q3 conducts, supplying voltage to the motor. As the load on the motor increases, the back EMF voltage decreases and the pulse-width off-time gets shorter.

Printhead Control. Control of data to the recorder printhead is accomplished by the Display/Recorder/Pacing microprocessor in conjunction with Recorder Gate Array U3, Program RAM U13, and Program EPROM U14. The Display/Recorder/Pacing microprocessor receives data from the System microprocessor and formats it for transmission to the recorder printhead.

The Recorder Gate Array receives parallel data from the Display/Recorder/Pacing microprocessor and transmits this data in serial form to the recorder printhead. The printhead prints an array of 384 bits from the top to the bottom of the recorder paper. The Recorder Gate Array translates the two ECG signal data bytes from hexadecimal code to the 320 bits required by the printhead. Bytes for generating hexadecimal characters are mapped to the remaining 64 bits (32 for the top, 32 for the bottom) since text for the recorder always appears at fixed locations at the top or bottom of the recorder paper. As in the Display Gate Array, two internal trace generators receive high and low address inputs for the ECG trace and the SYNC marks.

The thermal printhead is an array of 384 switches and dots. An internal 384-bit shift register (see Figure 1–11) receives data input to the printhead. The Recorder Gate Array sends clock pulse CK to the printhead which enters data on the DI line into the shift register and converts the data from serial to parallel. The CK pulse is active only when data is sent to the printhead, thereby preventing the printhead shift register from running continuously. The source of the CK pulses is the E clock used by the System microprocessor and the Display/Recorder/Pacing microprocessor.

The Display/Recorder/Pacing microprocessor moves the data stored in the printhead shift register into the printhead latch via the LD line. Data in the shift register is loaded into a 384-bit latch with a high LD signal from the Display/Recorder/Pacing microprocessor. The output of the latch is gated by 384 AND gates; the drive for the dot heating elements is provided by 384 NPN transistors. After the data is loaded into the latch, the Display/Recorder/Pacing microprocessor Strobe line switches one input of the AND gates high. This causes the gates switched high (data =1) on the other input to switch their outputs high. The selected transistors conduct, dropping a voltage across their respective heating elements and activating the selected printhead dots.




Figure 1–11 Printhead Assembly

**Paper Out Sense.** System microprocessor U30 is informed of the presence or absence of paper by a light emitting diode and photo transistor sensor in the recorder. U47 and R176 provide biasing so the recorder paper reflective properties will be detected.

**Printhead Resistance Evaluation.** The Printhead Resistance Evaluation circuit is active only when the instrument is in motor speed calibration mode. When the motor speed calibration test runs, the Display/Recorder/Pacing microprocessor starts energizing the printhead dots.

The measured printhead resistance used together with printhead temperature and voltage during printing determines the optimum duty cycle of the printhead. The printhead duty cycle is the time period that the printhead dots are energized; this is regulated by the width of the printhead strobe pulse, SB (see Figure 1–11). If the strobe pulse is too short, the print will be light; if it is too long, the excessive heat on the dots will greatly shorten the life of the printhead. The System microprocessor monitors the printhead resistance, battery voltage, and printhead temperature, and calculates the optimum printhead strobe width.

Monitoring the current through the printhead while increasing the number of dots turned on determines the printhead resistance. As the number of dots turned on increases, the combined resistance of the printhead decreases. Comparator U2A monitors the load on the VHD line at pin 5, through current sense resistor R27 and voltage divider R13 and R16. The VBATT line is monitored at pin 4 through voltage dividers R10 and R15. As the load on the printhead increases, the voltage across R16 decreases. When the voltage

across R16 drops lower than the voltage across R15, U2A switches low. This directs the Display/Recorder/Pacing microprocessor to stop sending strobe signals to the Recorder Gate Array at the STROBE input. The Display/Recorder/Pacing microprocessor then determines the printhead resistance from the number of dots turned on and passes this information to the System microprocessor.

Leads Off Detection. The display/recorder/pacing microprocessor U1 monitors a leads off signal from optical coupler U12 on the Interconnect Pacing PCB. A low on LEADS OFF indicates patient impedance exceeds  $6k\Omega$  and U1 sends a signal to the display to alert the operator.

Pacing Current Delivery. The display/recorder/pacing microprocessor U1 controls the delivery of pacing current. PACER POWER activates the pacer power supply. A high on IP FEEDBACK from the Interconnect PCB confirms current delivery through the pacing cable. A low variable length duty cycle PACE COMMAND pulse (from U1 pin 15) determines the output current level.

ECG Out (J29). The ECG signal provided by this circuit is produced for display on an external monitor and is separate from the ECG signal seen on the display or recorder. Data from U1 is converted to an analog signal by DAC U51. The signal is centered around 3.75V and amplified by the noninverting two-pole low-pass filter U52A. U52A also amplifies the signal and provides buffering to drive the high-pass filter C123 and R192. C123 and R192 provide AC coupling and output impedance capable of driving an external device with input impedance of ≥5kΩ Front panel ECG SIZE adjusts ECG out gain. When ECG size is X1.0, ECG out gain is 1V/mV.

Speaker Drive (J35). Tone frequency data is sent from the system microprocessor U30 to DAC U51. The converted output to op-amp U52B, Q14, and Q15 drives a 4W,  $8\Omega$  speaker. An internal oscillator in U30 provides four tone frequencies. Adjust QRS tone volume from the rear panel switch. Tone volume from U1 is sent to DAC U51.

## Display PCB/CRT **Assembly (803706)**

The Display PCB/CRT Assembly (Figure 1–12) is composed of a CRT and the Monitor PCB Assembly. The Monitor PCB assembly processes video information and synchronization signals from the Main PCB and displays the information on the CRT. The incoming signals are Video 1, Video 2, and synchronization signals Frame Sync and Raster Sync.




Figure 1-12 Display PCB/CRT Assembly Block Diagram

Two separate circuits provide horizontal and vertical CRT beam deflection. The CRT raster consists of 768 scan lines in one complete frame with 256 pixels per line. The CRT screen is scanned by a vertical raster, bottom to top and left to right (see Figure 1–13). There is no interlacing during the raster scan.




Figure 1-13 Display Scanning

Video Signal Amplifier

This circuit receives video information at the Video 1 and Video 2 inputs, controlling the voltage to the CRT cathode at the video output (J902 pin 2). The input to Video 1 and Video 2 is supplied from the Display Gate Array on the Main PCB at TTL voltage levels. The outputs of U401A, U401B, and U401C provide biasing for Q401, causing Q401 to conduct and supply a nominal +40V to the video output. Four levels of video intensity are available by switching the outputs of U401A, U401B, and U401C. Brightness is adjustable through Video Gain, R402.

Frame Sweep Generator

A deflection hybrid IC, U101, drives the frame yoke to sweep the CRT during horizontal retrace. The Display Gate Array on the Main PCB supplies a nominal 60Hz synchronizing signal, Frame Sync, for this circuit.

To sweep the CRT accurately, U101 requires biasing by a number of external components. The RC combination, R105, R106, and C102, sets the frequency of the hybrids internal oscillator. Horizontal Hold, R106, controls horizontal rolling. Horizontal Width, R111, controls the width adjustment and Horizontal Linearity, R113, calibrates the CRT beam linearity. R128 senses deflection current and provides feedback to U101 pin 10 to help control both width and linearity. R117 and C112 form a damping circuit to eliminate yoke ringing (oscillation).

Raster Scan Oscillator

The Raster Scan Oscillator sets the sweep frequency for the Raster Generator and the switching frequency for the Raster Voltage Generator power supply. The circuit is designed around voltage regulator U201.

The voltage regulator oscillator frequency is set to about 52kHz by the Raster signal. This signal is sent from the Display Gate Array on the Main PCB to Q201, causing Q201 to conduct and pass the signal to U201 pin 5. Resistor R217 and capacitor C203 function as an RC network, slowing the current at power-on. Resistor R209 decouples the +15V supply voltage so that display size will be minimally affected by changes in the +15V supply.

Input overvoltage protection is provided by Q201, Q202, and D201. If the input voltage rises to 18V or greater, D201 starts to conduct. This signal forward-biases Q203 which routes it to the gate of Q202. When Q202 conducts, it inhibits the internal oscillator in U201 by applying +5V to the input at pin 4.

Vertical Size, R224, controls the width of the pulse from the output stage of U201 which determines the amount of energy delivered to transformer T301.

Raster Voltage Generator

The Raster Voltage Generator supplies the input voltage for the primary of T301 as well as overall CRT bias voltages. The input signal for the circuit comes from the output of voltage regulator U201 pin 10 in the Raster Scan Oscillator circuit.

#### LIFEPAK 9P defibrillator/monitor/pacemaker

The Raster Voltage Generator circuit is designed around transformer T501 and transistors Q501 and Q502. The oscillating output of U201 drives Q501, providing additional current for Q502. The output of Q502 sets the voltage for the primary of T301. The +15V input at T501 pin 4 produces  $\pm15V$  (30V total) output at pin 5. A tap at T501 pin 6 supplies +24V. The secondary winding on T1 adds +14V to boost the +24V to +40V output at T501 pin 8. This voltage provides a ground reference for biasing G1, G2, and G4.

**Raster Generator** 

Vertical deflection for the raster yoke is provided by the Raster Generator circuit. Because of high current and voltage requirements, the vertical deflection circuit consists of discrete components, forming a resonant circuit. The Raster Generator, like the Raster Voltage Generator, is driven from the Raster Scan Oscillator circuit. The Raster Generator circuit includes transformer T301, which is shared with the CRT Bias Generator and provides the drive signal for the 88mH yoke coil. This creates the vertical scan drive for the CRT.

The Display Gate Array on the Main PCB supplies the synchronizing signal, Raster, which is passed to voltage regulator U201. The output of U201 pin 9 causes Q204 to conduct, providing additional current for Q205. The output of Q205 controls the charge and discharge of C208 through C212. The precise values of these capacitors are essential for maintaining the proper resonance while charging and discharging. Transformer T301 provides high voltage for charging the capacitors.

**CRT Bias Generator** 

The CRT Bias Generator supplies half-wave rectified voltage for the CRT anode and grid from the output of transformer T301. It provides up to -150V for bias of G1, up to +700V for bias of G2 and G4, and between +8.5kV and +11.6kV for anode bias. Master Brightness, R302, adjusts the brightness which largely determines the general background brightness of the CRT. Focus, R306, adjusts the focus.

# Introduction

This section provides information about the basic operation of the instrument. It is not intended to instruct an operator in the clinical use of the instrument. For clinical use, refer to the separate *Operating Instructions* manual. This section includes information about the following:

- Front Panel Controls and Indicators, including Paddles Area
- Rear Panel Controls, Indicators, and Connectors
- Option Settings, including self-tests and power-on default settings
- Clock Setting Procedure.

# Controls and Indicators

LIFEPAK 9P defibrillator/monitor/pacemaker controls and indicators are identified in Figures 2–1 and 2–2 and described in the corresponding Tables 2–1 and 2–2.

Note: Pushbuttons are indicated in sans-serif CAPITAL LETTERS and displayed messages are indicated in BOLD CAPITAL LETTERS.




Figure 2-1 Front Panel Controls and Indicators

	Table 2-1	Tiont Fallet Controls and indicators
Item	Control or Indicator	Function
1	ON	Toggles power on and off. Power-on illuminates the LED.
2	▼ LEAD SELECT▲	Selects PADDLES input or patient cable leads (STD, LEAD I, LEAD II, or LEAD III). Switches selected lead display to reverse video and advances to the next lead selection. Defaults to LEAD II at power—on unless PADDLES is specified in the Setup Menu.
3	▼ ENERGY▲ SELECT	Displays the energy selection range incrementally. Select <b>LOW</b> on the high range menu to display the low range selections from 0 to 9 joules. Select <b>HIGH</b> on the low range menu to display the high range selection from 0 to 360 joules. Selecting 0J inhibits energy charging. Pressing during charging or prior to energy transfer causes stored energy to discharge internally.
4	▼ ECG ▲ SIZE	Changes the displayed ECG size to reverse video and increases or decreases the gain by one setting. ECG size is adjustable from 0.2cm/mV to 4.0cm/mV. Defaults to gain of X1.0 at power—on. Affects gain of ECG OUTPUT.
5	CHARGE	Initiates defibrillation charge cycle. Existing defibrillator status message is replaced with XXX JOULES CHARGING message, where XXX displays the stored energy as it ramps up to the selected value. LED flashes during charge. At charge completion XXX JOULES AVAILABLE message displays while LED remains on steadily.
6	HR ALARM	Displays the preset alarm limits in reverse video and lights the LED. Repeated pressing incrementally steps through all three pairs of default limit selections as defined in the Setup Menu (factory preset defaults are 150/40, 120/60, and 160/90).
7	CODE SUMMARY	Starts the recorder which prints a summary of critical event information (defibrillation/cardioversion event, heart rate alarm violations, and selected ECG segments). Pressing pushbutton again stops recorder.
8	SYNC	Activates sync mode message display and LED light. LED flashes with each QRS detection. Pressing SYNC again returns to defibrillate (asynchronous) mode and turns off the LED.
9	RECORD	Toggles the recorder on and off.
10	CRT Display	No-fade CRT display. ECG trace moves from right to left.
11	QRS Indicator	Heart symbol flashes with each QRS detection.
12	Heart Rate	Displays digital QRS rate (20 to 300bpm).
13	Heart Rate Alarm Limits	Display indicators for high/low heart rate alarm settings. Limits can be adjusted through the Setup Menu. Limits flash and alarm sounds if violated.
		$\cdot$

Table 2-1 Front Panel Controls and Indicators (cont.)

Item	Control or Indicator	Function	
15	Lead Select Indicator	Display indicates lead selected.	
16	Calibrated ECG Gain	Display indicates actual ECG gain selected.	
17	Status Message Display	Provides operating status information:	
	SYNC MODE	Indicates the synchronous mode is activated.	
	FOR SYNC: USE LEADS	Flashes if the SYNC pushbutton is pressed and the PADDLES lead is selected with the standard paddles connected (the synchronous mode may not be used with the PADDLES lead selected and standard paddles connected).	
	NO PAPER	Flashes if the recorder is activated but no paper is installed. The recorder will not operate without paper.	
	LOW BATTERY: CONNECT AC POWER	Appears if the defibrillator/monitor is powered by battery and the battery voltage is low. As soon as this message appears, connect the defibrillator/monitor to ac power for continuous operation and to recharge the battery.	
	SERVICE	Indicates the internal self-diagnostic software has detected a fault status condition. Refer to page 3-46 for more information.	
18	Energy Selection Display	Selection from high energy range shown. Press ENERGY SELECT to display the high energy range menu. Select <b>LOW</b> to display the low energy range menu. Energy selection range stays on screen for 10 seconds following the last energy level selected or until pressing another pushbutton.	
19	JOULES SELECTED Indicator	On-screen indication of selected energy level in joules. Number increases up to selected energy level as charging occurs. Defaults to 200J at power-on.	
20	BATT CHRG Indicator	Indicates battery is charging.	
21	Recorder	Records ECG and annotation in delay ECG mode only. Prints date, time, sync mode (if active), ECG lead, ECG gain, and heart rate. Prints selected energy and defibrillation mode during defibrillation. Records automatically for 16 seconds whenever a heart rate alarm violation occurs.	
22	PACER	Pressing pushbutton starts pacing function and lights pacemaker LED. Defaults to <b>DEMAND</b> pacing. Stop pacing by doing one of the following: press PACER again, reduce current to 0mA, press ENERGY SELECT, charge the defibrillator, or press analyze or charge on the LIFEPAK 9 Shock Advisory Adapter. Pacemaker LED flashes with each pacemaker pulse.	

Table 2-1 Front Panel Controls and Indicators (cont.)

Item	Control or Indicator	Function
23	▼ RATE ▲	Selects pacing rate from 40 to 170 pulses per minute (ppm) in 10ppm increments. Pressing the up symbol increases the rate; pressing the down symbol decreases the rate. If a pacing lead disconnects, the ppm setting remains and the pacing current automatically returns to 0mA.
24	CURRENT	Turning rotary control increases or decreases pacemaker current. Current selection is 0mA or 7mA to 200mA in 1mA increments. Turning clockwise increases the current; turning counterclockwise decreases the current. If a pacing electrode disconnects or the pacing is interrupted, the current automatically returns to 0mA.
PADD	LES AREA	
25	Patient Cable Connector	Connection for 6-pin patient cable.
26	STERNUM Paddle	Defibrillation electrode with one discharge pushbutton, usually placed to left of sternum (patient's right). Also serves as negative ECG electrode during QUIK-LOOK paddle monitoring.
27	APEX Paddle	Defibrillation electrode with QUIK-CHARGE control (CHARGE) and one discharge pushbutton, usually placed near cardiac apex.  Also serves as positive ECG electrode during QUIK-LOOK paddle monitoring.
28	CHARGE (QUIK– CHARGE control)	Momentary pushbutton to charge defibrillator from APEX paddle.
29	Charge Indicator	LED flashes during charge and glows steadily when energy has reached selected level.
30	Discharge Pushbuttons	Pushbuttons to discharge energy. <u>Both</u> pushbuttons must be pressed <u>simultaneously</u> to deliver energy to the paddles. ENERGY WILL NOT BE DELIVERED UNLESS THE INSTRUMENT IS FULLY CHARGED TO SELECTED LEVEL.
31	Pacing Cable Connector	Receptacle for connecting pacing cable. Use Physio-Control pacing cable part number 802905.
32	Message Display	Indicates the following: DEMAND PACING, NON-DEMAND PACING, PACING STOPPED, PACER FAULT, PACING LEADS OFF, USE ECG LEADS, CHARGE REMOVED, or ENERGY FAULT.
33	Pacing Pulses Display	Indicates pacing pulses per minute (ppm). Defaults to 60ppm each time the pacemaker is turned on.
34	Pacing Current Display	Indicates the pacing current setting in milliamps (mA). Defaults to 0mA each time the pacemaker is turned on.
35	Test Load (not shown)	$50\Omega$ defibrillator test load. Metal contacts (not shown, in paddles storage area) for receiving defibrillation discharge energy from paddles.




Figure 2-2 Rear Panel Controls, Indicators, and Connectors

Table 2–2 Rear Panel Controls, Indicators, and Connectors

Item		Function
1	TIME/DATE	Pushbuttons for setting the real-time clock.
	MODE	Pressing once selects the clock setting mode. Pressing repeatedly selects minutes, hour, year, month, day fields. Pressing again exits the clock setting mode.
	SET	Pressing advances the value displayed in the selected clock set field.
2	▼ QRS VOL ▲	Pressing the up or down scroll pushbutton increases or decreases the QRS volume. The last setting is retained when power is off.
3	1mV CAL	Pressing displays a 1mV calibration signal on the CRT and recorder.
4	Ground	Equipotential ground tie point.
5	ECG OUTPUT	ECG out connector; 1V/mV output with ECG gain of 1 (press ECG SIZE to set gain).
6	AC Power Connector	AC power connector.
7	Fuse Cover	Access panel for changing fuses.
8	Mains Power	Rocker switch turns ac line power on (1) or off (0). On some models a clip is installed which locks this switch in the on position.
9	NON-DEMAND PACE	Selects non-demand (asynchronous) pacing. Pressing the pushbutton a second time returns the instrument to <b>DEMAND</b> pacing.

# **Option Settings**

The LIFEPAK 9P defibrillator/monitor/pacemaker includes software which allows the user to initiate self-tests, change settings and power-on defaults such as language and heart rate alarm limits, and examine system data. Access for all of these options is gained through the top-level Test Menu which then allows access to the Setup Menu or Info Menu. All menus and options are described below.

### **Test Menu**

The Test Menu is the top-level menu which allows access to the self-tests and the two other menus, the Setup Menu and the Info Menu. To access the Test Menu, do the following:

- Make sure the defibrillator/monitor/pacemaker power is off.
   Simultaneously press and hold LEAD SELECT and HR ALARM, then press ON.
- Continue holding LEAD SELECT and HR ALARM until the message SELECTED TEST: and a list of tests are displayed as shown in Figure 2–3. This information is the Test Menu.




Figure 2-3 Test Menu

The Test Menu allows access to the five self-tests (KEYPAD, DEFIB, DISPLAY, RECORDER and PACER), the Info Menu (INFO), and the Setup Menu (SETUP). The self-tests are executed as part of the Performance Inspection Procedure and the Test and Calibration Procedure and are described in detail in Section 3. The Info Menu and the Setup Menu are described in the following paragraphs.

To exit the Test Menu and return to normal operation, press  $\P$  ECG SIZE to scroll down to EXIT, then press HR ALARM.

## **Setup Menu**

The Setup Menu allows you to examine and change various operating settings such as the language, notch filters, heart rate alarm limits, and the alarm volume. To access the Setup Menu, do the following:

- 1. First, access the Test Menu by turning power off. Then simultaneously press and hold LEAD SELECT and HR ALARM while pressing ON.
- 2. After the Test Menu appears as shown in Figure 2–3, press ▼ ECG SIZE to scroll down to SETUP. Then press HR ALARM to display the Setup Menu as shown in Figure 2–4.




Figure 2-4 Setup Menu

The different option settings are as follows:

The different option settings are as follows.

LANGUAGE: ENG

ENGLISH, FRENCH, GERMAN, or SPANISH

RECORDER:

**MONitor or DIAGnostic** 

NOTCH:

50Hz or 60Hz

**LEADS NOTCH:** 

YES or NO

LEAD:

PADDLES or LEAD II

HR LIMITS:

3 limits can be configured from the following

ranges: 80bpm to 240bpm and 20bpm to 100bpm.

**ALARM VOL:** 

1 to 8 (A short tone sounds with each level change

to indicate the selected volume.)

- 3. Press HR ALARM to scroll through the selections. Press ▼ ECG SIZE ▲ to scroll through the different option settings.
- 4. After you choose the desired option setting, press HR ALARM again to store the setting and move to the next selection. Whatever options are selected become the power-on default settings.
- 5. To exit the Setup Menu and return to the Test Menu, press LEAD SELECT. To exit the Test Menu and return to normal operation, press ▼ ECG SIZE to scroll down to EXIT, then press HR ALARM.

### info Menu

The Info Menu allows you to examine the manufacturing date and EPROM revision numbers, and to examine and clear the number of defibrillator discharges recorded and fault status codes recorded. To access the Info Menu, do the following:

- First, access the Test Menu. Turn power off, then simultaneously press and hold LEAD SELECT and HR ALARM while pressing ON.
- After the Test Menu appears as shown in Figure 2–3, press ▼ ECG SIZE to scroll down to INFO. Then press HR ALARM to display the Info Menu as shown in Figure 2-5.




Figure 2-5 Info Menu

**BUILT:** Manufacturing date (user cannot change).

Revision number for system microprocessor **SYSTEM REV:** 

EPROM; changes automatically with any EPROM

update (user cannot change).

**DISPLAY REV:** Revision number for the display/recorder

> microprocessor EPROM; changes automatically with any EPROM update (user cannot change).

Records the number of defibrillator discharges; user **DEFIB USES:** 

may RETAIN or CLEAR the number.

**FAULT STATUS:** Retains any coded instrument errors using the error

codes listed in Section 3 on page 3-46; user may

**RETAIN** or **CLEAR** codes.

- 3. Press HR ALARM to scroll through the selections. For DEFIB USES: or **FAULT STATUS:**, press ▼ ECG SIZE ▲ to toggle between CLEAR and RETAIN.
- To exit the Info Menu and return to the Test Menu, press LEAD SELECT. To exit the Test Menu and return to normal operation, press ▼ ECG SIZE to scroll down to EXIT, then press HR ALARM.

# Clock Setting Procedure

To set the clock which is displayed on the CRT and recorder, do the following:

- 1. Press TIME/DATE MODE on the rear panel to initiate the clock setting mode. The monitor displays the current clock setting in the lower left corner of the CRT as shown in Figure 2–6. The single minutes field is highlighted.
- Press TIME/DATE SET on the rear panel to change the single minute setting. Each time the pushbutton is pressed, the value of the field increases by one increment.
- 3. Press MODE again to advance to the next field. Press SET to select the desired value. Repeat this to set the entire time and date.
- 4. To exit the clock setting mode, press MODE again after selecting the day field.




Figure 2-6 Clock Set Display

### Introduction

This section describes how to test, calibrate, and troubleshoot the LIFEPAK 9P defibrillator/monitor/pacemaker. Topics cover:

- The Performance Inspection Procedure (PIP)
- The Test and Calibration Procedure (TCP)
- Troubleshooting Aids.

Note: Prefixes are used to distinguish testing steps.

PIP— identifies Performance Inspection Procedure steps.

TCP— identifies Test and Calibration Procedure steps.

### Performance Inspection Procedure (PIP)

This Performance Inspection Procedure (PIP) checks whether a LIFEPAK 9P defibrillator/monitor/pacemaker is performing within specifications. The instrument case does not need to be opened to perform the PIP. The PIP Checklist (page 3–21) may be photocopied and used to record the PIP results. If the instrument fails any part of the PIP, refer to the Test and Calibration Procedure (page 3–23) and the Troubleshooting Aids (page 3–46) for help in locating any failed subassemblies.

#### **PIP-Purpose**

This PIP checks whether the LIFEPAK 9P defibrillator/monitor/pacemaker is performing within specifications. Perform the PIP regularly as a periodic maintenance check and after any repair or calibration.

PIP-Scope and Applicability

This PIP applies to all configurations of the LIFEPAK 9P defibrillator/monitor/pacemaker. This PIP does *not* apply to the LIFEPAK 9, LIFEPAK 9A, or LIFEPAK 9B defibrillator/monitors; refer to their respective service manuals for testing information.

#### **PIP-Definitions**

Acronyms used in this procedure are:

BPM Beats Per Minute

DMM Digital Multimeter

ECG Electrocardiogram

PIP Performance Inspection Procedure

QRS or QRS Complex Refers to portions of the ECG waveform

**NSR** Normal Sinus Rhythm

p-ppeak-peakppmPulse Per Minute

**Note:** Pushbuttons are indicated in CAPITAL LETTERS and displayed messages are indicated in **BOLD CAPITAL LETTERS**.

#### **PIP-Requirements**

The requirements for test equipment, workstation, and personnel are listed below along with important safety information.

PIP-Required Equipment

Performing this PIP requires the equipment listed in Table 3–1. Although specific test instruments are recommended, other test equipment with equivalent specifications may be used.

Table 3-1 PIP Test Equipment

Equipment	Specifications	Manufacturer		
Defibrillator Energy Meter	Power Range: 0-400J Load Resistance: 50Ω±1% Accuracy: ±3% Defib Waveform Amplitude Attenuation: 1000:1 (minimum) with oscilloscope output	Bio-Tek QED-6 or equivalent*		
ECG simulator	Accuracy: 1% of selected rate	Bio-Tek QED-6 or equivalent*		
Leakage Tester/ Safety Analyzer	110Vac Line Voltage: Current Range: 0-2000µA Current Accuracy: ±1% Resistance Accuracy: ±1%	Dynatech Nevada 232, or equivalent*		
	230Vac Line Voltage: Current Range/Accuracy: 200μA/5% of reading 2000μA/1% of reading Resistance Accuracy: ±1%	Dynatech Nevada 432 HD, or equivalent*		

<sup>\*</sup> Equivalent equipment is required to meet or exceed the specifications listed in the specifications column.

Table 3-1 PIP Test Equipment (cont.)

Equipment	Specif	ications	Manufacturer		
Pacing Tester/Simulator	Peak current: Output frequency:	7-250mA ±2% 40-170ppm ±3%	Bio-Tek QED-6 or equivalent*		
Patient Cable	3-lead, AHA or 3-lead, IEC	,	Physio-Control PN 9-10418-02 Physio-Control PN 800947-01		
Pacing Cable	2-conductor		Physio-Control PN 802905		
Digital Multimeter (DMM)	4 1/2 digits DC Voltage Accura	cy: 0.1% of reading	Fluke Model 8012A or equivalent*		
Recorder paper	50mm (use only Ph	ysio-Control paper)	Physio-Control PN 804700		
Metric ruler or recorder paper	Resolution: 1mm divisions		Recorder paper: Physio-Control PN 804700		
Timer	1 second resolution		Aristo or equivalent		
Components for Pacemaker Leads Detect check (construction defined in Figure 3–8, page 3–19)	R1 Resistor, 1.0kΩ, 1/8W, 1% R2 Resistor, 10kΩ, 1/8W, 1% R3 Resistor, 100Ω, 1/8W, 1% C1 Capacitor, 15nF, 50V		R2 Resistor, 10kΩ, 1/8W, 1% R3 Resistor, 100Ω, 1/8W, 1%		Physio-Control PN 200054-193 Physio-Control PN 200054-288 Physio-Control PN 200054-097 Physio-Control PN 201628-007

<sup>\*</sup> Equivalent equipment is required to meet or exceed the specifications listed in the specifications column.

Caution	Possible Equipment Damage. Be sure to use only recorder paper from
	Physio-Control (PN 804700). Other paper may damage the printhead and may not print legibly.

# PIP-Test Equipment Verification

All test equipment used in performing the PIP must have a current calibration label affixed to its exterior. The calibration label must be issued by a certified calibration facility.

#### PIP-Workstation

The ac line power for this workstation must be connected to a grounded power source.

#### PIP-Personnel

Personnel performing this PIP must be thoroughly familiar with the operation of the LIFEPAK 9P defibrillator/monitor/pacemaker. In addition, personnel performing this PIP must meet at least one of the following levels of education or experience:

- Associate of applied science degree with a major emphasis in biomedical electronics
- Certificate of technical training in electronics with a major emphasis in biomedical electronics
- Equivalent biomedical electronics experience.

### Warning

SHOCK HAZARD. When charged and discharged during this PIP, the LIFEPAK 9P defibrillator/monitor/pacemaker discharges up to 360J of electrical energy through the defibrillator paddles. Unless discharged safely as described in this PIP, this electrical energy may cause injury or death. Do not attempt to perform this procedure unless you are thoroughly familiar with the operation of the LIFEPAK 9P defibrillator/monitor/pacemaker.

### **Warning**

POSSIBLE ARCING AND PADDLE DAMAGE. When discharging the defibrillator into the internal test load, make sure the standard paddles are securely and properly stored in the paddle storage area (STERNUM paddle on the left, APEX paddle on the right). This helps prevent arcing and formation of pits on paddle electrode surfaces. Pitted or damaged paddle electrode surfaces can cause patient skin burns during defibrillation.

#### **PIP-General Instructions**

The following sections are the actual PIP steps to be performed. These steps are designed to be executed in sequential order. For instance, if a heart rate input of 30bpm is specified for one test it also applies to the next test unless you are instructed to change it.

#### **PIP-Physical Inspection**

Before beginning the PIP-Physical Inspection, be sure to disconnect the ac power cord from the ac power source.

- Inspect all exterior surfaces for signs of abuse or other damage, including fluid spills, cracks and dents.
- Inspect paddles for pitted electrode plates, presence of dried defibrillation gel, and integrity of cables and switches. Remove gel or other foreign matter from paddles and paddle wells.
- Inspect recorder for ease of door operation and positive latching. Remove dust or foreign matter inside recorder.
- Inspect all external connections for bent, broken, or corroded pins, and debris.
- Inspect the power cord and all cables and accessories for damage, cracks, breaks, bent or broken connector pins, or other signs of wear.
- Open the battery door on the bottom of the instrument and make sure batteries are installed. Close the door and make sure it latches properly.
- Press the ON pushbutton and confirm that the green LED illuminates.
   Press ON again and make sure the green LED is off and power is off (the display is blank). Then press all other pushbuttons and check them for uniform mechanical resistance.

#### **PIP-AC Power**

- 1. Connect the instrument to ac power.
- If the dash number of the instrument is listed below, confirm that the
  power module clip is installed at the bottom of the rear panel power
  switch, locking the switch to the power on (I) position. Then go on to
  step 3.

#### **Dash Numbers with Power Module Clip Installed**

805460-01 805460-03 805460-10

If the dash number of the instrument is *not* listed above, press the bottom of the rear panel power switch to place it in the on (l) position.

- 3. Press ON to apply power to the instrument and confirm that the ON LED illuminates.
- 4. Confirm that © 1990, PCC is displayed for at least 5 seconds.
- 5. Confirm that within 10 seconds after the instrument is turned on, the monitor displays the power-on information as shown in Figure 3–1.




Figure 3-1 Power-On Display

#### **PIP-Battery Charge**

With instrument power on and the ac power source connected, confirm that **BATT CHRG** is displayed in the lower left corner of the display.

#### **PIP-Battery Power**

- With instrument power on, disconnect the power cord from the ac power source. Confirm that the display remains on but BATT CHRG is not displayed.
- 2. Reconnect the ac power cord and press ON to turn off instrument power.

#### PIP-Setup Menu Settings

This procedure describes how to access the Setup Menu and select the proper settings required for the rest of the PIP. Before changing any settings, write down the current settings so you can restore these settings at the completion of the PIP.

- To access the Setup Menu, make sure instrument power is off. Simultaneously press and hold LEAD SELECT and HR ALARM, then press ON.
- Continue holding LEAD SELECT and HR ALARM until the message SELECTED TEST: and a list of tests are displayed as shown in Figure 3–2. This information is the Test Menu.




Figure 3-2 Test Menu

3. Press ▼ ECG SIZE ▲ to scroll through the Test Menu. Scroll to SETUP and then press HR ALARM to display the Setup Menu as shown in Figure 3–3. Before changing any settings, write down the current user settings so you can restore them after completing the PIP.




Figure 3-3 Setup Menu

▼ ECG SIZE ▲ to change an option selection. Select:

LANGUAGE:

Match local requirements

RECORDER:

DIAG

NOTCH:

Match local line requirements

**LEADS NOTCH:** 

NO

LEAD:

**LEAD II** 

HR LIMITS:

150 120 040 060 160 090

**ALARM VOL:** 

Ω

5. To exit and store the changes, press HR ALARM to scroll back to LANGUAGE:, then press LEAD SELECT to return to the Test Menu. Leave the Test Menu displayed for the next test, PIP-Keypad Test.

#### **PIP-Keypad Test**

- The Keypad Test is initiated from the Test Menu (as shown in Figure 3-2) which should be displayed from the end of the previous step. Press ▼ ECG SIZE ▲ to scroll to KEYPAD, then press HR ALARM to start the Keypad Test.
- 2. Press the front panel pushbuttons and then the rear panel pushbuttons in the following order:

#### Front Panel:

PACER

LEAD SELECT

**▼** RATE

RATE A

▼ ECG SIZE

ECG SIZE

**▼ ENERGY SELECT** 

**ENERGY SELECT** 

HR ALARM

CHARGE

**CODE SUMMARY** 

RECORD

SYNC

#### Rear Panel:

**NON-DEMAND PACE** 

▼ QRS VOL

QRS VOL A

MODE

1mV CAL

SET

Press the APEX paddle CHARGE pushbutton.

- 4. Press the APEX paddle discharge pushbutton; confirm no beep sounds.
- 5. Press STERNUM paddle discharge pushbutton; confirm no beep sounds.
- 6. Simultaneously press both paddle discharge pushbuttons; confirm that a beep sounds.
- 7. Confirm that the monitor displays **KEYPAD**: **OK**, then press **LEAD** SELECT to return to the Test Menu. If the pushbuttons are not pressed in the order listed above, or if there is a pushbutton failure, the monitor displays **KEYPAD**: **FAULT**. In this case, press **LEAD** SELECT to return to the Test Menu. Then press **▼** ECG SIZE **△** to scroll to **KEYPAD**, press HR ALARM to restart the Keypad Test, and repeat steps 2 through 7.
- 8. To exit the Keypad Test after completion, press ▼ ECG SIZE ▲ to scroll to EXIT and press HR ALARM to resume normal operation.

#### **PIP-Lead Select**

- 1. Use the patient cable to connect the ECG simulator to the instrument. Set the ECG simulator for a 60 bpm-NSR signal.
- 2. Press LEAD SELECT repeatedly and confirm that the signal displays in all lead (I, II, III) configurations.

#### PIP-Recorder Speed/Baseline Centering

- 1. Confirm there is paper in the recorder.
- 2. Set the ECG simulator for a 60bpm-NSR signal in **LEAD II**.
- 3. Press RECORD, let the recorder run for approximately 10 seconds, then select **STD**. Let the recorder run for an additional 10 seconds, then press RECORD to stop the recorder.
- 4. Measure the distance between the maximum peaks of the recorded ECG trace. Confirm that the distance between consecutive peaks is 25 ±1mm.
- 5. Measure the baseline centering of the recorded **STD** trace. Confirm that the baseline is centered within ±1mm.

**PIP-ECG Display** 

Confirm that the flatline display on the CRT appears level and centered. Confirm that all messages are clearly visible as shown in Figure 3–4.




Figure 3-4 Display with Flatline ECG Trace

PIP-1mV CAL

In order to perform this test the **LEADS NOTCH** must be set to **NO** and the recorder must be set to **DIAG** as described in the PIP-Setup Menu Settings on page 3–6.

- 1. Confirm that the **STD** lead remains selected from the previous PIP test (press LEAD SELECT if necessary to select **STD**).
- 2. Press RECORD, then press the 1mV CAL pushbutton on the rear panel once. Confirm that a 10.0 ±0.5mm p-p calibration pulse is printed out as illustrated in Figure 3–5 (check measurement from the center of the printed trace).




Figure 3-5 1mV CAL Pulse

# PIP-CRT Gain/Aspect Ratio

- 1. Select **LEAD II**.
- 2. Set the ECG simulator for a 60bpm-NSR signal.
- 3. Adjust ▼ ECG SIZE ▲ for a 34 ±2mm p-p ECG display signal.
- Press RECORD.
- 5. Confirm that the signal amplitude printed on the ECG paper is  $20 \pm 1$ mm p-p.

#### **PIP-Heart Rate**

- 1. With **LEAD II** selected, set ECG Size to X1.0.
- 2. Set the ECG simulator for the input rates listed below. Confirm that the systole indicator (the heart-shaped symbol) flashes on the display with each QRS complex and that the displayed HR corresponds:

Input Rates (bpm)	HR
30	30 ±2bpm
60	60 ±2bpm
240	240 ±5bpm

3. Return the ECG simulator selection to 60bpm-NSR.

#### PIP-Rate Alarms/ Alarm-Activated Record

- 1. Press HR ALARM. Confirm that the HR ALARM LED is on and 150/40 is displayed in reverse video.
- While 150/40 remains displayed in reverse video, press HR ALARM twice and confirm that 120/60 and 160/90 display consecutively.
- 3. Press HR ALARM to display **150/40** in reverse video. Confirm that after about 2 seconds **150/40** displays in normal video.
- 4. Set the ECG simulator rate to greater than 150bpm. Confirm that the alarm sounds, the **150** flashes, and the recorder runs for approximately 16 seconds.
- 5. Confirm that the heart rate printed on the printout is 151bpm or higher.
- 6. Press HR ALARM and confirm that the alarm stops.
- Decrease the ECG simulator rate to 30bpm. Press HR ALARM. Confirm that the alarm sounds, the 40 flashes, and the recorder runs for approximately 16 seconds.
- 8. Press HR ALARM and confirm that the alarm stops.

#### PIP-Recorder Annotation/ Real-Time Clock

- 1. Set the ECG simulator to 60bpm.
- Set the instrument to local time. (Press TIME/DATE MODE on the rear panel to display the clock setting; press TIME/DATE SET to increment to the desired value in the displayed field, then press TIME/DATE MODE to advance to the next field.)
- 3. Press RECORD and confirm the recorder begins printing with the date and time annotated.
- 4. Confirm that the time is accurate to within  $\pm 1$  minute.
- 5. Press RECORD to stop the recorder.

#### **PIP-Test Load**

- 1. Make sure the paddles are stored in the paddle storage area of the instrument.
- Press ▼ ENERGY SELECT ▲ as necessary to set the energy level to 200J.
- Press CHARGE.
- 4. After the defibrillator reaches full charge, simultaneously press both paddle discharge pushbuttons to transfer the energy into the test load. Confirm that TEST 200 JOULES DELIVERED displays on the monitor and annotates on the recorder printout.

# PIP-Sync/Sync Annotation

- 1. Press LEAD SELECT to set the instrument to **LEAD II**. Set the ECG simulator for a 60bpm-NSR signal.
- 2. Press SYNC. Confirm that **SYNC MODE** displays on the CRT, the sync marker (the inverted triangle) is clearly visible on each displayed R wave, and the SYNC pushbutton LED flashes.
- 3. Place the paddles in the paddle storage area of the defibrillator. Press ▼ ENERGY SELECT ▲ if necessary to set the energy level to 200J.
- 4. Press CHARGE. After the defibrillator reaches full charge, simultaneously press and hold both paddle discharge pushbuttons until the energy is discharged into the test load.
- 5. Confirm that the recorder prints the message **SYNC TEST 200 JOULES DELIVERED** and the sync markers are clearly visible on each R wave on the printed ECG trace.
- 6. Confirm that the **SYNC MODE** display and SYNC pushbutton LED are off, and no sync markers are visible.

# PIP-Recorder/CRT Noise

- 1. Short together the RA and LL leads of the patient cable.
- 2. Set the ECG gain to X4.0. Select **LEAD II**.
- 3. Press RECORD to start the recorder.
- 4. Confirm that the CRT and recorder baseline noise is less than 3mm p-p.
- 5. Select **PADDLES** electrodes. Short the paddle plates together, wait at least 10 seconds, and then confirm that the CRT and the recorder baseline noise is less than 3mm p-p.

**Note:** The ECG signal that is displayed on the CRT is delayed approximately 10 seconds before being printed by the recorder.

6. Return the paddles to the paddle storage area. Press RECORD to turn off the recorder.

# PIP-Defib Disarm and CHARGE LEDs Check

- Press ▼ ENERGY SELECT ▲ and set the energy level to 0J.
- 2. Press the CHARGE pushbutton.
- 3. Confirm that the CHARGE LED does not light and that the CRT continues to display **0 JOULES SELECTED**.
- 4. Press ▼ ENERGY SELECT ▲ to set the selected energy to 200J.
- 5. Press CHARGE and confirm that both the front panel and the APEX paddle CHARGE LEDs flash during the charge cycle.
- 6. Confirm that the LEDs remain on and a tone sounds when full charge is reached.
- Simultaneously press both paddle discharge pushbuttons to transfer the charge. Confirm that the CHARGE LEDs are off and the recorder prints the message TEST 200 JOULES DELIVERED.

#### **PIP-Paddle Polarity**

- 1. Set the ECG simulator for a 60bpm-NSR signal.
- 2. Press LEAD SELECT to select **PADDLES** lead and remove the paddles from the paddles storage area.
- 3. Hold the paddles on the ECG simulator output connections (APEX on LL, STERNUM on RA) and confirm that a positive-going QRS complex is visible on the CRT.
- 4. Return the paddles to the paddle storage area.

#### **PIP-Charge Time**

- 1. Press ▼ ENERGY SELECT ▲ to set the energy level to 360J. Make sure the paddles are securely stored in the paddle storage area.
- Simultaneously press CHARGE and start the timer. Stop the timer when
  the defibrillator reaches full charge (as indicated by the tone and the
  displayed message 360 JOULES AVAILABLE). Confirm that the charge
  time is less than 10 seconds.
- Press the paddle discharge pushbuttons simultaneously to transfer the energy into the test load. Confirm that TEST USE 200 JOULES displays on the CRT and three tones sound.

#### **PIP-Discharge Control**

1. Press ▼ ENERGY SELECT ▲ to set the energy level to 360J.

### Warning

**SHOCK HAZARD.** Help avoid injury by holding the paddles and coil cords away from your body and away from any other person or object as described in the following steps.

- Remove the paddles from the paddle storage area. Extend the coil cords away from the paddle wells at least three feet and away from your body or any other object. Press CHARGE.
- 3. After the defibrillator reaches full charge (indicated by the tone and the displayed message 360 JOULES AVAILABLE) keep the coil cords extended and away from your body or any other object while shaking the paddles and coil cords with a moderate motion. Confirm that no discharge occurs.
- 4. Return the paddles to the paddle storage area.
- 5. Simultaneously press both discharge pushbuttons and confirm that the energy is discharged into the test load.

#### PIP-Energy Dump

- 1. Place the paddles on the energy meter electrodes.
- 2. Press ▲ ENERGY SELECT ▼ and set the energy level to 360J.
- 3. Press CHARGE.
- 4. After the instrument charges to 360J, change the energy level to 300J.
- 5. Confirm **CHARGE REMOVED** is displayed and tones sound. Return the paddles to the paddle storage area.

#### **PIP-Charge Reduction**

- Leave the paddles stored in the paddle storage area. Press ▲ ENERGY SELECT ▼ and set the energy level to 360J.
- 2. Press CHARGE.
- 3. When full charge is reached, start the timer and immediately change the energy level to 1J and press CHARGE (to perform this test properly, you should scroll to 1J and press CHARGE within a few seconds).
- 4. Confirm that full charge of 1J is reached within 20 seconds or less after you start the timer. Simultaneously press both discharge pushbuttons and confirm that the energy is discharged into the test load.

#### **PIP-Delivered Energy**

Connect the paddles to the energy meter electrodes. Refer to Table 3–2 and confirm that the delivered energy is within tolerances for each energy level.


Table 3–2 Delivered Energy Tolerances


Energy Setting	Tolerance	Energy Setting	Tolerance
<b>1J</b> :	±0.2J	10Ј	±1.0J
<b>2</b> J	±0.2J	20Ј	±2.0J
3Ј	±0.3J	30J	±2.1J
<b>4</b> J	±0.4J	50J	±3.5J
5J	±0.5J	100J	±7.0J
6J	±0.6J	200J	±14.0J
<b>7</b> J	±0.7J	300J	±21.0J
8Ј	±0.8J	360J	±25.2J
9Ј	±0.9J		

# PIP-CODE SUMMARY Report

- 1. Press ON twice to switch instrument power off then on.
- 2. Place the paddles on the energy meter electrodes. Press LEAD SELECT to set the instrument to **LEAD II**.
- 3. Connect the ECG simulator to the patient cable and set for a 60 bpm-NSR signal.
- 4. Press SYNC.
- 5. Press ▲ ENERGY SELECT ▼ and set the energy level to 200J.
- Press CHARGE. When full charge is reached (as indicated by the tone and displayed message 200J AVAILABLE), simultaneously press both paddle discharge pushbuttons to transfer the energy into the energy meter.

- 7. Press RECORD twice to switch the recorder on, then off quickly.
- 8. Press CODE SUMMARY.
- 9. Confirm that the format and appearance of the recorded data is similar to the data shown in Figure 3–6.








Figure 3-6 Example of CODE SUMMARY Printout

#### PIP-Pacemaker Power-Up 1.

- Connect the patient cable and pacing cable to the instrument and to the pacing simulator with pacing simulator power off.
- Press PACER and confirm the PACER LED illuminates and the CRT displays the following: DEMAND PACING, LEAD II, 60 PPM, 00MA.

#### **PIP-Pacing Rate**

- 1. Make sure the pacing cable is connected to the instrument and the pacing simulator with pacing simulator power off.
- Turn the rotary PACER knob to a current level between 50mA and 55mA. Confirm the PACER LED flashes with each new pacer spike on the display.
- 3. Press ▲ RATE ▼ as needed to increase the pacing rate to 150ppm. Push RECORD to print a recording. Examine the printout and confirm that the printout indicates the pacer rate, current level, and pace markers (the arrows at the bottom of the printout) at each pacing spike. Confirm that the distance between pacing spikes is 10±1mm.
- 4. Adjust the rate to 170ppm. Confirm that each rate level is selectable by sequentially stepping down the rate in 10ppm increments to 40ppm.
- 5. Adjust the pacing rate to 60ppm. Print and examine a printout and confirm that the distance between pacing spikes is  $25 \pm 1$ mm.

#### PIP-Demand Pacing

- 1. Make sure the pacing simulator power is off. Turn the CURRENT knob to 130mA and confirm that the instrument remains in the demand mode.
- Set the pacing simulator rate to 40 and turn the pacing simulator power
  on. Make sure the pacemaker pacing rate is set to 60ppm. Confirm that
  the PACER LED flashes off and on in conjunction with the R-waves and
  R-wave markers which are displayed.
- 3. Set the pacing simulator rate to 80. Confirm that normal QRS complexes and R-wave markers are displayed. Confirm that the PACER LED does not repeatedly flash off and on but remains on continuously.

- 1. Make sure the pacing simulator power is on with the pacing simulator rate set to 80.
- 2. Press NON-DEMAND PACE on the rear panel of the instrument and confirm that R-waves are displayed.
- 3. Press RECORD to print a recording. Confirm that the distance between pacer spikes is 25 ±1mm.

#### **PIP-Pacing Leads Off**

1. Make sure the pacing simulator power is on with the pacing simulator rate set to 80.

### Warning


**SHOCK HAZARD.** Hazardous electrical voltages are present at the pacing cable leads when pacing is active. Do not touch the uninsulated portion of the pacing cable lead.

- Turn the CURRENT knob to 7mA. Disconnect one of the pacing cable leads from the pacing simulator and confirm the monitor/defibrillator/ pacemaker emits a continuous beeping sound and displays the message PACING LEADS OFF.
- 3. Press PACER twice to turn the pacemaking function off and on again.
- 4. With the pacing cable still disconnected, turn the CURRENT knob fully clockwise. Confirm that the monitor/defibrillator/pacemaker emits 3 one-second tones and displays the message PACING LEADS OFF.

#### **PIP-Leakage Current**

Check leakage current in accordance with the following standards: AAMI/ANSI (Association for the Advancement of Medical Instrumentation/American National Standards Institute) ES1-1985, AAMI/ANSI DF2-1989, CSA (Canadian Standards Association) C22.2, IEC (International Electrotechnical Commission) 601-1, IEC 601-2-4, and UL (Underwriters Laboratories) 544.

Because of the wide variety of Safety Analyzers that may be used for these tests, specific instructions are not provided in this PIP except for the PIP-Pacemaker Leads Detect Current described below. Refer to your Safety Analyzer operator's manual for instructions about configuring and performing these tests. The maximum allowable leakage current is summarized in Table 3–3 for 120Vac, 60Hz, and in Table 3–4 for 230Vac, 50Hz. When measuring ground resistance for the leakage current tests, connect the ohmmeter between the rear panel ground stud and the ground pin on the ac power connector as shown in Figure 3–7.



**DEFIBRILLATOR/MONITOR REAR PANEL** 

Figure 3-7 Test Connections for Ground Resistance Measurement

PIP-Pacemaker Leads Detect Current To measure the leakage current for the Pacemaker Leads Detect, construct the Pacemaker Load as shown in Figure 3–8. Then perform the following steps:

- 1. Make sure the instrument power is off.
- 2. Connect the Pacemaker Load to the instrument and the DMM as shown in Figure 3–8.
- Select the DMM 200mV ac range.

4. Turn on the instrument power and set the pacemaker current output to 0mA. Confirm that the DMM ac voltage reading is ≤19.6mV (≤19.6mV as measured through the 1kΩ load as shown in Figure 3-8 is ≤20µA current).




Figure 3–8 Test Setup for Pacemaker Leads Detect Current

### Warning

**ELECTRICAL HAZARD.** Make sure you are familiar with your test equipment and how to execute these tests. Failure to properly perform these tests could result in not detecting excessive leakage current.

**Note:** After completing the PIP, be sure to restore the Setup Menu option settings required by the user during normal operation (these settings were noted when you performed PIP-Setup Menu Settings on page 3-6).

Table 3-3 Maximum Leakage Current for 120Vac and 60Hz

Parameter	Closed Safety Ground SOURCE (μΑ)	Open Safety Ground SOURCE (μΑ)	Closed Safety Ground SINK (μΑ)	
Chassis to Ground	≤69	≤69	N/A	
ECG Shorted All Leads	≤9.0	≤9.0	≤19	
Interlead RA-LA	≤9.0	≤9.0	≤19	
Interlead RA-RL	≤9.0	≤9.0	≤19	
Interlead LA-RL	≤9.0	≤9.0	≤19	
Paddles, Shorted	≤69	≤69	≤69	
Paddles, Paired	≤69	<b>≤</b> 69	N/A	
Pacing, All Leads Shorted	≤9.0	≤9.0	≤19	
Pacing, Paired	≤9.0	≤9.0	≤19	

NOTE: 1. All tests are conducted in both hot/neutral normal and reverse polarities.

- 2. Test procedure limits based on AAMI/ANSI Standards with meter inaccuracies factored in.
- 3. N/A = Not Applicable.
- 4. Ground resistance:  $<0.050\Omega$  as measured between the rear panel ground stud and the ground pin on the ac power connector (refer to Figure 3–7).
- 5. Pacemaker Leads Detect Current  $\leq 20\mu A$  ( $\leq 19.6$ mV through 1k $\Omega$  load as per Figure 3–8).

Table 3-4 Maximum Leakage Current for 230Vac and 50Hz

Parameter	Closed Safety Ground SOURCE (μΑ)	Open Safety Ground SOURCE (μΑ)	Closed Safety Ground SINK (μΑ)
Chassis to Ground	≤73	10≤ i* ≤123	N/A
ECG Shorted All Leads	≤6.4	≤24	≤24
interlead RA-LA	≤8.5	≤18	≤18
Interlead RA-RL	≤8.5	≤18	≤18
Interlead LA-RL	≤8.5	≤18	≤18
Paddles, Shorted	≤24	≤24	≤384
Paddles, Paired	≤18	≤18	N/A
Pacing, All Leads Shorted	≤6.4	≤24	≤20
Pacing, Paired	≤8.5	≤18	≤20

NOTE: 1. All tests are conducted in both hot/neutral normal and reverse polarities.

- 2. Test procedure limits based on IEC Standards with meter inaccuracies factored in.
- 3. N/A = Not Applicable.
- 4. Ground resistance:  $<0.100\Omega$  as measured between the rear panel ground stud and the ground pin on the ac power connector (refer to Figure 3–7).
- 5. Pacemaker Leads Detect Current  $\leq 20\mu A$  ( $\leq 19.6$ mV through 1k $\Omega$  load as per Figure 3–8).
- \* i = measured current.

	***************************************	****************	000000000000000000000000000000000000000	***************************************
PIP Checklist	Customer:			
	Dept/Location:_			
	Instrument Type	»:		Model:
	Serial Number:			I act Inon
	Track Name (ID:			Last Insp:
	Tech. Name/ID:			
Performance Inspection	Pass	Fail	Comments	
<b>Physical Inspection</b>				
AC Power				
Battery Charge				
<b>Battery Power</b>				
Setup Menu Settings	User s	LAÑ	GUAGE:	ation (circle or fill blank):  ENG FRE GER SPA  ENG POR ITA SWE
		NOT: LEAI LEAI	OS NOTCH:	MON or DIAG 50HZ or 60HZ NO or YES LEAD II or PADDLES 150 120 160 040 060 090
		ALA	RM VOLUME:	
Keypad Test				
Lead Select				
Recorder Speed/Baseline	_	_		
Centering				
ECG Display	. 🗖			
1mV CAL		Ξ.		
CRT Gain/Aspect Ratio		<u> </u>		
Heart Rate		□ -		
Rate Alarms/Alarm-Activ Record	vated $\Box$			
Recorder Annotation/	_			
Real-Time Clock		Ξ.		
Test Load				
Sync/Sync Annotation Recorder/CRT Noise				
Defib Disarm and CHAR	GE _			
LEDs Check		Π.		
Paddle Polarity				
Charge Time				
Discharge Control				
Energy Dump		Π.		
Charge Reduction		Π.		
Delivered Energy				
CODE SUMMARY Repo				
Pacemaker Power-Up				
Pacing Rate Demand Pacing				
Demano Pacing		$\Box$		

<b>Performance Inspection (Cont.)</b>	Pass	Fail	Comments
Non-Demand Pacing			
Pacing Leads Off			
Leakage Current (see below)			
<b>Pacemaker Leads Detect Current</b>			
Restore user Setup Menu settings			

#### Maximum Leakage Current for 120Vac and 60Hz

Parameter	Closed Safety Ground SOURCE (μA)		Open Safety Ground SOURCE (μΑ)		Closed Safety Ground SINK (μΑ)	
	Limit	Actual	Limit	Actual	Limit	Actual
Chassis to Ground	≤69		≤69		N.	/A
ECG Shorted, All Leads	≤9.0	·	≤9.0		≤19	
Interlead RA-LA	≤9.0		≤9.0		≤19	
Interlead RA-RL	≤9.0		≤9.0		≤19	
Interlead LA-RL	≤9.0		≤9.0		≤19	
Paddles, Shorted	≤69		≤69		≤69	
Paddles, Paired	≤69		≤69		N	/A
Pacing, All Leads Shorted	≤9.0		≤9.0		≤19	
Pacing, Paired	≤9.0		≤9.0		≤19	

NOTE: 1. All tests are conducted in both hot/neutral normal and reverse polarities.

- 2. Test procedure limits based on AAMI/ANSI Standards with meter inaccuracies factored in.
- 3. N/A = Not Applicable.
- 4. Ground resistance:  $<0.050\Omega$  as measured between the rear panel ground stud and the ground pin on the ac power connector (refer to Figure 3–7).
- 5. Pacemaker Leads Detect Current  $\leq 20\mu A$  ( $\leq 19.6$ mV through 1k $\Omega$  load as per Figure 3–8).

#### Maximum Leakage Current for 230Vac and 50Hz

Parameter	Closed Safety Ground SOURCE (μA)		Open Safety Ground SOURCE (μΑ)		Closed Safety Groun SINK (μΑ)	
	Limit	Actual	Limit	Actual	Limit	Actual
Chassis to Ground	≤73		10≤i*≤123		N	/A
ECG Shorted, All Leads	≤6.4		≤24		≤24	
Interlead RA-LA	≤8.5		≤18		≤18	
interlead RA-RL	≤8.5		≤18		≤18	4 E ***
Interlead LA-RL	≤8.5		≤18		≤18	
Paddles, Shorted	≤24		≤24		≤384	
Paddles, Paired	≤18		≤18		N	/A
Pacing, All Leads Shorted	≤6.4		≤24		≤20	
Pacing, Paired	≤8.5		≤18	- ***	≤20	

NOTE: 1. All tests are conducted in both hot/neutral normal and reverse polarities.

- 2. Test procedure limits based on IEC Standards with meter inaccuracies factored in.
- 3. N/A = Not Applicable.
- 4. Ground resistance:  $<0.100\Omega$  as measured between the rear panel ground stud and the ground pin on the ac power connector (refer to Figure 3–7).
- 5. Pacemaker Leads Detect Current  $\leq 20\mu A$  ( $\leq 19.6 \text{mV}$  through  $1 \text{k}\Omega$  load as per Figure 3–8).
- \* i = measured current.

# Test and Calibration Procedure (TCP)

The Test and Calibration Procedure (TCP) contains in-depth testing and calibration instructions for the LIFEPAK 9P defibrillator/monitor/pacemaker. Unlike the PIP (Performance Inspection Procedure), only those tests considered necessary need to be performed, and these may be conducted in any order. It is necessary to separate the case halves for some of the tests.

#### **TCP-Purpose**

These tests may be performed if the defibrillator/monitor/pacemaker fails a step in the PIP or if there is a suspected failure in an assembly. Perform the PIP after completing any of the tests in this section.

**Note:** If you perform any of the tests in this TCP, perform the PIP (page 3–1) afterward to confirm that the defibrillator/monitor/pacemaker is performing within specifications.

# TCP-Scope and Applicability

This TCP applies to the LIFEPAK 9P defibrillator/monitor/pacemaker. This TCP does *not* apply to the LIFEPAK 9, LIFEPAK 9A, or LIFEPAK 9B defibrillator/monitors; refer to their respective service manuals for testing information.

#### **TCP-Definitions**

DMM	Digital Multimeter
ECG	Electrocardiogram
ESD	Electrical Static Discharge
NSR	Normal Sinus Rhythm
PIP	Performance Inspection Procedure
р-р	peak-to-peak
ppm	Pulse Per Minute
QRS or QRS Complex	Refers to portions of the ECG waveform
SSD	Static-Sensitive Device
TCP	Test and Calibration Procedure

#### **TCP-Requirements**

The following are minimum requirements:

# TCP-Equipment Requirements

Test equipment required for maintenance and calibration of the LIFEPAK 9P defibrillator/monitor/pacemaker is listed in Table 3–5. Although specific test equipment is recommended, other test equipment with equivalent specifications may be used.

Table 3–5 Test and Calibration Test Equipment

Equipment Specifications		Manufacturer		
Storage Oscilloscope	Bandwidth: dc to 2kHz Vertical Accuracy: ±5% (5mV-5V/division) Horizontal Time Base Accuracy: ±1%	Tektronix Model 2232 or equivalent*		
Digital Multimeter (DMM)	4-1/2 digits, 2A Range DC Voltage Accuracy: 0.1% of reading Low Resistance: $0.01\Omega$ resolution on $200\Omega$ scale	Fluke Model 8012A or equivalent*		
Function Generator	Output Waveform: Sinusoidal Amplitude/Freq. Accuracy: ±0.5dB/±5% Amplitude/Freq. Range: 10V/11MHz	Krohn-Hite Model 5400A or equivalent*		
ECG Simulator	Sinewave Output: 1.0 ±0.5mVp-p Output Frequency: 10.0 ±0.5Hz	Bio-Tek QED-6 or equivalent*		
Defibrillator Energy Meter	Power Range: 0-400J Load Resistance: 50Ω±1% Accuracy: ±3% Defib Waveform Amplitude Attenuation: 1000:1 (min.) with oscilloscope output	Bio-Tek QED-6 or equivalent*		
Pacing Tester/Simulator	Peak current: 7-250mA ±2% Output frequency: 40-170ppm ±3%	Bio-Tek QED-6 or equivalent*		
Pacing Cable	2-conductor	Physio-Control PN 802905		
DC Power Supply	Output Voltage: 0 to 20Vdc Output Current: 0 to 3A	Topward Model TPS 2000 or equivalent*		
Timer	Resolution: 1 second	Aristo or equivalent		
Components for Recorder Frequency Response test (construction defined in Figure 3–15, page 3–40)	R1, R2 Resistor, 100kΩ, 1/8W, 1% R3, R4 Resistor, 100Ω, 1/8W, 1%	Physio-Control PN 200054-384 Physio-Control PN 200054-097		
Components for Pacemaker Test Load Setup (construction defined in Figure 3–16, page 3–43)	<ul> <li>R1 Resistor, metal-film, 1.0kΩ, 25W, 1%</li> <li>R2 Resistor, metal-film, 10kΩ, 1/8W, 1%</li> <li>R3 Resistor, metal-film, 100Ω, 1/8W, 1%</li> <li>C1 Capacitor, plastic-dielectric, 15nF, 50V</li> </ul>	Ohmite PN 825F1K0 Physio-Control PN 200054-288 Physio-Control PN 200054-097 Physio-Control PN 201628-007		
Phone plug	0.25 inch plug tip			

<sup>\*</sup> Equivalent equipment is required to meet or exceed the specifications listed in the specifications column.

TCP-Test Equipment Verification	All test equipment used in performing the TCP must have a current calibration label affixed to its exterior. The calibration label must be issued by a certified calibration facility.		
TCP-Workstation Requirements	Workstations where the TCP is performed must provide an Electrical Static Discharge (ESD) protected environment. The ac line power for this workstation must be connected to a grounded power source.		

#### **TCP-Personnel**

Personnel performing this TCP must be thoroughly familiar with the operation of the LIFEPAK 9P defibrillator/monitor/pacemaker. In addition, personnel performing this TCP must meet at least one of the following levels of education or experience:

- Associate of applied science degree with a major emphasis in biomedical electronics
- Certificate of technical training in electronics with a major emphasis in biomedical electronics
- Equivalent biomedical electronics experience.

### Warning

SHOCK HAZARD. When charged and discharged during this TCP, the LIFEPAK 9P defibrillator/monitor/pacemaker discharges up to 360J of electrical energy through the defibrillator paddles. Unless discharged safely as described in this TCP, this electrical energy may cause injury or death. Do not attempt to perform this procedure unless you are thoroughly familiar with the operation of the LIFEPAK 9P defibrillator/monitor/pacemaker.

### Warning

POSSIBLE ARCING AND PADDLE DAMAGE. When discharging the defibrillator into the internal test load, make sure the standard paddles are securely and properly stored in the paddle storage area (STERNUM paddle on the left, APEX paddle on the right). This helps prevent arcing and formation of pits on paddle electrode surfaces. Pitted or damaged paddle electrode surfaces can cause patient skin burns during defibrillation.

### Warning

**SHOCK HAZARD.** Terminals and wires carrying high voltages up to +11.6kV are exposed with the instrument cases removed. Do not touch internal circuitry except as indicated in these procedures.

#### Caution

**Possible Component Damage.** Static Sensitive Devices (SSDs) require special handling. Follow the special handling procedures described on page 4–6 to prevent static discharge and damage. Refer to Table 5-1 to identify PCB assemblies containing SSDs.

#### TCP-Voltage Checks

- 1. With the instrument power off, remove the front panel to access the test connector J3 on the Power Conversion PCB (leave the front panel connector cable connected). For disassembly instructions, refer to page 4–7, Front Panel Removal.
- 2. Connect the ground lead of the DMM to pin 7 of J3. Refer to Figure 3–9.
- 3. Turn the instrument power on and check that the voltages are within the limits listed in Table 3–6.




Figure 3-9 Test Connector J3

Table 3-6 Power Conversion PCB Voltage Checks

Test Point	Signal Name	Nominal Voltage	Limits
Pin 4	+15V	+15.0V	±0.3V
Pin 11	+5V	+5.00V	±0.25V
Pin 9	+5V High Limit Detect	+2.25V	±0.15V
Pin 10	+5V Low Limit Detect	+2.625V	±0.125V
Pin 18	SW BATT	+18.0V	+2.0V to -2.4V
Pin 25	V BATT	+18.4V	+1.6V to -2.8V

# TCP-Low Battery Threshold Checks

- 1. With the instrument power off, remove the front panel to access the test connector J3 on the Power Conversion PCB (leave the front panel connector cable connected). For disassembly instructions, refer to page 4–7, Front Panel Removal.
- 2. Remove the battery door on the bottom of the instrument and disconnect the battery. Connect the DC power supply to the battery connector inputs. Connect ground to pin 1 or 2 (black wire) and power to pin 3 or 4 (red wire).

#### **Monitoring Check**

- 1. Use the DMM to monitor the voltage at the J3 pin 25 as shown in Figure 3–9 (connect the DMM ground to J3 pin 7).
- 2. Turn the instrument power on and initially set the power supply output to  $+18.00 \pm 0.05$ Vdc.
- 3. Slowly decrease (by 0.1V steps) the output voltage until the monitor displays the message LOW BATTERY: CONNECT AC POWER (decrease the power supply voltage quickly enough to cause the message to appear within 45 seconds or less). Confirm that the power supply voltage level is within the Low Battery Limits listed for the Monitoring Check in Table 3–7.

Table 3-7 Low Battery Thresholds

Test Point	Function	Low Battery Limits	Shutdown Limits
Pin 25	Monitoring Check	15.50V ±0.23V	14.50V ±0.21V
Pin 25	Defib Charge Check	15.00V ±0.22V	13.50V ±0.20V
Note: Once illuminated, the LOW BATTERY: CONNECT AC			CONNECT AC

POWER message will not reset until the battery voltage rises above 16.5V.

4. Return the power supply voltage to  $+18.00\pm0.05$ Vdc. Slowly reduce the power supply voltage until the CRT display blanks out. Confirm that the power supply voltage is within the Shutdown Limits shown in Table 3–7 for the Monitoring Check.

# Defibrillator Charge Check

- 1. Use the DMM to monitor the voltage at the J3 pin 25 as shown in Figure 3–9 (connect the DMM ground to J3 pin 7).
- 2. Turn the instrument power on and initially set the power supply output to  $+18.00\pm0.05$ Vdc.
- Charge the defibrillator to 360J.
- 4. Gradually reduce the power supply voltage toward the upper Low Battery Limit detection threshold of +15.22Vdc.
- 5. With the defibrillator paddles securely stored in the paddle wells, discharge the defibrillator into the internal test load of the defibrillator.
- 6. Repeat steps 3 through 5 above, lowering the power supply voltage toward the lower Low Battery Limit detection threshold of +14.78Vdc until the LOW BATTERY: CONNECT AC POWER message is displayed. Confirm that the triggering voltage is within the Low Battery Limits shown in Table 3–7 for the Defib Charge Check.
- 7. Return the power supply voltage to  $+18.00 \pm 0.05$ Vdc.
- 8. Charge the defibrillator to 360J.
- 9. Slowly reduce the power supply voltage toward the upper Shutdown Limit detection threshold of +13.70Vdc.
- 10. Discharge the defibrillator into the internal test load of the defibrillator.
- 11. Repeat steps 8 through 10 above, lowering the power supply voltage toward the lower Shutdown Limit of +13.30Vdc until the CRT display blanks out. Confirm that the power supply voltage is within the Shutdown Limits shown in Table 3–7 for the Defib Charge Check.

#### TCP-Test Menu Access

Four of the five tests on the Test Menu (**DEFIB**, **DISPLAY**, **RECORDER** and **PACER**) are used in Test and Calibration Procedures. To access the Test Menu, do the following:

- 1. Turn the instrument power off.
- 2. Press ON while pressing and holding down LEAD SELECT and HR ALARM until the Test Menu is displayed as shown in Figure 3–10.




Figure 3-10 Test Menu

- Press ▼ ECG SIZE ▲ to scroll through the selections.
- Press HR ALARM to start a test.
- Press LEAD SELECT to interrupt a test in progress and return to the Test Menu.
- To exit the test menu, press ▼ ECG SIZE ▲ to scroll to EXIT, then press HR ALARM.

# TCP-Defibrillator Calibration

1. With the instrument power off, place the paddles on the energy meter electrodes.

# Caution

**Possible Energy Meter Damage.** Start each test at the lowest energy level to avoid overheating the energy meter. Test each energy level in successive order.

- Access the Test Menu (turn power off, then press ON while pressing and holding down LEAD SELECT and HR ALARM until the Test Menu is displayed).
- 3. Press ▼ ECG SIZE ▲ to scroll to the **DEFIB** test.
- 4. Press HR ALARM to start the **DEFIB** test.
- 5. Press ▼ ENERGY SELECT ▲ to scroll to LOW. Select 1J.
- 6. Press CHARGE.
- After the instrument reaches full charge as indicated by the tone, simultaneously press both paddle discharge pushbuttons to transfer the energy into the energy meter.
- 8. Confirm that the CRT displays the message **JOULES DELIVERED: XXX**, where XXX is the energy selected.
- 9. Read the value on the energy meter.
- 10. Press ▼ ECG SIZE ▲ as necessary to change the JOULES

  DELIVERED: XXX value to match the value on the energy meter.
- 11. Press HR ALARM to enter the adjusted value.
- 12. Press CHARGE. After charging is complete, simultaneously press both paddle discharge pushbuttons to transfer the energy into the energy meter. Confirm that the energy meter value equals the selected value as indicated in Table 3–8.

Selected Value (J)	Acceptable Range		Selected	Acceptable Range	
	Value (J)	% Tolerance	Value (J)	Value (J)	% Tolerance
1	0.8 to 1.2	±20%	20	19.0 to 21.0	±5%
2	1.8 to 2.2	±10%	30	28.5 to 31.5	±5%
3	2.7 to 3.3	±10%	50	47.5 to 52.5	±5%
4	3.8 to 4.2	±5%	100	95.0 to 105.0	±5%
5	4.8 to 5.2	±5%	200	190.0 to 210.0	±5%
6	5.7 to 6.3	±5%	300	285.0 to 315.0	±5%
7	6.7 to 7.3	±5%	360	342.0 to 378.0	±5%
8	7.6 to 8.4	±5%			
9	8.6 to 9.4	±5%			
10	9.5 to 10.5	±5%			

- 13. Press HR ALARM to reinitiate the **DEFIB** test.
- 14. Repeat steps 5 through 13 for each energy level.
- 15. Return the paddles to the paddle storage area on the defibrillator.

# TCP-Test Load Calibration

**Note:** Complete TCP-Defibrillator Calibration before performing this calibration of the test load.

- 1. Make sure the paddles are securely stored in the paddle storage area of the defibrillator.
- 2. Access the Test Menu (turn power off, then press ON while pressing and holding down LEAD SELECT and HR ALARM until the Test Menu is displayed).
- 3. Press ▼ ECG SIZE ▲ to scroll to the DEFIB test.
- 4. Press HR ALARM to start the **DEFIB** test.
- 5. Press ▼ ENERGY SELECT ▲ to set the selected energy to 200J.
- 6. Press CHARGE.
- 7. After the defibrillator reaches full charge, simultaneously press the paddle discharge pushbuttons to transfer the energy into the test load.
- 8. When the CRT displays **JOULES DELIVERED: 200**, press HR ALARM.

- 9. Press LEAD SELECT to return to the Test Menu.
- Scroll to EXIT and press HR ALARM to resume normal operation. Select STD lead.
- 11. Press CHARGE. After the defibrillator reaches full charge, simultaneously press the paddle discharge pushbuttons to transfer the energy into the test load.
- 12. Confirm that the CRT displays **TEST 200 JOULES DELIVERED**.
- 13. Confirm that the recorder prints the time and date followed by **TEST 200 JOULES DELIVERED**.

#### **TCP-Charge Time**

- 1. Make sure the paddles are securely stored in the paddle storage area of the defibrillator.
- 2. Press ▼ ENERGY SELECT ▲ to set the selected energy to 360J.
- 3. Simultaneously start the timer and press CHARGE.
- 4. Stop the timer when the defibrillator reaches full charge. Confirm that the time required to reach full charge is less than 10 seconds.
- 5. Simultaneously press both paddle discharge pushbuttons to transfer the energy into the test load.

#### TCP-Output Waveform


- 1. Place the paddles on the energy meter electrodes. Place the APEX paddle on the right electrode of the energy meter so that the pulse polarity is positive.
- 2. Connect the storage oscilloscope to the energy meter signal output.
- 3. Set up the oscilloscope for the "single shot" acquisition mode for the pulse waveform capture. Refer to the Operating Instructions for the oscilloscope for assistance. For waveform parameter measurement, place both the left and right measurement markers at both the 50% and 10% levels.
- Press ▼ ENERGY SELECT ▲ to select 360J.
- Press CHARGE.
- After the defibrillator reaches full charge, simultaneously press the paddle discharge pushbuttons to transfer the energy into the energy meter.

7. Inspect the waveform displayed on the oscilloscope and confirm that the waveform meets the specifications listed below and illustrated in Figure 3–11.

Peak current = 46A to 64A peak Time above 50% = 2.10ms to 4.20ms Time above 10% = 3.10ms to 9.20ms

#### **PEAK CURRENT - TIME ABOVE 50%**

#### PEAK CURRENT - TIME ABOVE 10%






Figure 3-11 Output Waveform

# TCP-Display Calibration

The following paragraphs describe the calibration procedure for the Display Assembly PCB. Figure 3–12 shows the location of all Display Assembly PCB potentiometers.

1. Disconnect the ac power cord. Remove front and rear cases to allow access to the potentiometers. For disassembly instructions, refer to page 4–10, Case Separation.

# Warning

**SHOCK HAZARD.** The Display Assembly produces high voltages up to +11.6kV. Do not contact the CRT anode connector or the CRT bias generator circuitry when the instrument power is on.




Figure 3–12 Display PCB Assembly Potentiometer Locations

### Warning

**FLYING GLASS HAZARD.** The CRT is made of glass and contains a vacuum which may implode if broken. Handle with care and wear safety glasses when removing the CRT from the Display Assembly PCB.

2. Reconnect the ac input cord. Reconnect J20 on the front panel.

Note: Reroute J1 ribbon cable underneath C32 on the Power Conversion PCB to prevent strain on the capacitor during the following tests.

# Warning

**SHOCK HAZARD.** High voltages up to 1000V may be present on the leads of R302 and R306. Use an insulated screwdriver for adjustment.

- 3. Press ON. Adjust R302, Master Brightness (CUTOFF), until the vertical raster lines are visible behind the characters.
- 4. Turn R106, Horizontal Hold (FRAME HOLD), approximately 1/8 turn past the point where the horizontal rolling stops.
- Access the Test Menu (turn power off, then press ON while pressing and holding down LEAD SELECT and HR ALARM until the Test Menu is displayed).
- 6. Press ▼ ECG SIZE ▲ to scroll to the **DISPLAY** test, then press HR ALARM to start the test.
- 7. Press ▼ ECG SIZE ▲ to select the PATTERN test.
- 8. Adjust: R224, Vertical Size (SCAN SIZE); L201, Linearity (SCAN LIN); and R111, Horizontal Height (FRAME HEIGHT), to obtain a display pattern that fills most of the screen.
- 9. Loosen the screw holding the yoke clamp.
- 10. Rotate the yoke until the display pattern is level.
- 11. Tighten the screw on the yoke clamp and confirm that the display pattern is still level.
- 12. Adjust the yoke rings to center the display pattern on the CRT.
- 13. Adjust R224 for a display pattern grid height of 73 ±2mm. Adjust L201 to make each box of the display pattern approximately the same height.

**Note:** The adjustments of R224 and L201 are interactive. Repeat the preceding step until both measurements are correct.

- 14. Turn R402, Video Gain (VID GAIN), fully counterclockwise.
- Turn the instrument off and allow the CRT to cool down for approximately 10 minutes. Then, access the Test Menu and select the DISPLAY: PATTERN test.

### Warning

**SHOCK HAZARD.** Potentially dangerous voltages are present on the pins of R302 and R306. Use an insulated screwdriver for adjustment.

- 16. Turn R302 counterclockwise just until the raster disappears. Turn R302 counterclockwise approximately 1/16 turn more.
- 17. Turn R402 clockwise until the normal display intensity is visible.

**Note:** Raster will not reach normal intensity if R302 is turned too far counterclockwise.

18. Adjust R111 for an overall pattern width of 106 ±2mm. Adjust R113, Horizontal Linearity (FRAME LIN) so that each box is approximately the same width.

**Note:** The adjustments of R111 and R113 are interactive. Repeat the above step until both measurements are correct.

- 19. Adjust R402 to obtain a dim display. Adjust R306, Focus, for maximum sharpness of the display characters.
- 20. Adjust the yoke rings to center the test pattern on the CRT within 2mm. Confirm that the difference in height from the right side to the left side of the test pattern is less than or equal to 1mm. Check that the difference in width from the top to the bottom is less than or equal to 2mm.
- 21. Turn R402 clockwise until the normal display intensity is visible.

### TCP-Power Supply Calibration

 To gain access to R23 on the Power Supply PCB, first remove front and rear cases (for disassembly instructions, refer to page 4–10, Case Separation). Then remove the four screws holding the Power Supply PCB and pull it away from the chassis to access R23 as shown in Figure 3–13.




Figure 3–13 Power Supply Calibration

2. Reconnect the battery to the instrument.

**Note:** The Power Supply compensates for temperature fluctuations. The instrument must be powered—on for at least 30 minutes prior to calibration so the temperature around thermistor RT2 stabilizes.

- 3. Connect the DMM between the positive and negative battery terminals.
- 4. Connect the ac power cord and press ON.
- 5. Locate R23. Remove the potting compound from R23 as required to allow adjustment. Adjust R23 for  $\pm 18.40 \pm 0.18$  V.

### TCP-Recorder Calibration

1. Remove the front panel. The two potentiometers (R176 and R8 on the Main PCB) required for Recorder Calibration are accessible with the front panel removed as shown in Figure 3–14.




Figure 3–14 Recorder Calibration Adjustments Location

- Access the Test Menu (turn power off, then press ON while pressing and holding down LEAD SELECT and HR ALARM until the Test Menu is displayed).
- Press ▼ ECG SIZE ▲ to scroll to the RECORDER test, then press HR ALARM. Make sure the RECORDER: PATTERN test is selected.

**Note:** During this calibration procedure the printhead dot intensity is decreased to prevent damage to the printhead and drive roller. Consequently, print legibility may range from light to invisible.

- 4. Turn R176, Paper Out Sense, on the Main PCB fully clockwise.
- Make sure paper is installed in the recorder and press RECORD.
   Confirm that the recorder does not print, or if it does print that the NO PAPER message is displayed.
- 6. Slowly turn R176 counterclockwise while repeatedly pressing RECORD. Adjust R176 just to the point where the recorder is on and the NO PAPER message is *not* displayed. This procedure may need to be repeated to determine this threshold.
- 7. Remove the paper from the recorder and close the recorder doors.
- Counting the number of rotations required, turn R176 counterclockwise while repeatedly pressing RECORD. Adjust R176 just to the point where the recorder turns on. Note the number of rotations required to make the adjustment.
- 9. Adjust R176 to the midpoint of the number of rotations between steps 6 and 8.
- 10. Confirm that the recorder only operates when paper is loaded.
- 11. Press ▼ ECG SIZE ▲ to select the RECORDER: MOTOR SPEED test.
- 12. Press RECORD. Inspect the printout and confirm that the period of the printed sawtooth waveform is  $25 \pm 1$ mm. Adjust R8, Motor Speed Adjust, if necessary to obtain the period of  $25 \pm 1$ mm. Press RECORD to turn off the recorder.
- 13. Reinstall the front panel.

### TCP-Recorder Frequency Response

- Access the Test Menu (turn power off, then press ON while pressing and holding down LEAD SELECT and HR ALARM until the Test Menu is displayed).
- Press ▼ ECG SIZE ▲ to scroll to SETUP, then press HR ALARM to display the Setup Menu. Make sure RECORDER is set to DIAG and LEADS NOTCH is set to NO (if making any changes, press HR ALARM at the end to enter the option changes).
- 3. Press LEAD SELECT to return to the Test Menu, press ▼ ECG SIZE to scroll to EXIT, then press HR ALARM to resume normal operation.
- 4. Press LEAD SELECT to set the instrument to LEAD II. Press ▼ ECG SIZE ▲ as needed to select X1.0 gain.
- 5. Connect the function generator to the patient cable input using the test signal setup shown in Figure 3–15.




Figure 3–15 Test Signal Setup for Function Generator

6. Set the function generator for a 1Vp-p, 10Hz sinewave. Press RECORD to turn on the recorder.

**Note:** Because the ECG trace on the recorder printout lags behind the signal input by approximately 10 seconds, the recorder must run for approximately 10 seconds before a change to the input signal appears on the printout.

- 7. Set the function generator to 0.05Hz. After 15 seconds, examine the printout and confirm that the recorded trace is between 7mm p-p and 10.5mm p-p.
- 8. Set the function generator to 0.5Hz. After 15 seconds, examine the printout and confirm that the recorded trace is between 9mm p-p and 10.5mm p-p.
- 9. Set the function generator to 60Hz. After 15 seconds, examine the printout and confirm that the recorded trace is between 7mm p-p and 10.5mm p-p.
- 10. Set the function generator to 100Hz. After 15 seconds, examine the printout and confirm that the recorded trace is between 7mm p-p and 10.5mm p-p.
- 11. Set the function generator to 110Hz. After 15 seconds, examine the printout and confirm that the recorded trace is less than 7mm p-p.
- 12. Access the Test Menu and then the Setup Menu. Set the **LEADS NOTCH** to **YES**. Press HR ALARM to enter the option change.
- 13. Press LEAD SELECT to return to the Test Menu, press ▼ ECG SIZE to scroll to EXIT, then press HR ALARM to resume normal operation.
- 14. Press LEAD SELECT to select **LEAD II**. Set the function generator output to a 1Vp-p sinewave and match the input line frequency (50Hz or 60Hz).
- 15. Press RECORD. After 15 seconds, examine the printout and confirm that the recorded trace is less than 1mm p-p.
- 16. Access the Test Menu and then the Setup Menu. Set RECORDER to MON and set LEADS NOTCH to NO. Press HR ALARM to enter the option changes, then exit the Setup Menu and Test Menu.
- 17. Set the function generator for a 1Vp-p, 1Hz sinewave.
- 18. Press LEAD SELECT to select **LEAD II**. Press ▼ ECG SIZE ▲ to set the gain to X1.0.
- 19. Press RECORD. Wait 10 seconds, then slowly sweep the function generator frequency from 1Hz to 40Hz. Examine the printout and confirm the recorded trace is 7mm p-p minimum.

### TCP-CRT Frequency Response

- 1. Press  $\nabla$  ECG SIZE  $\triangle$  to increase the ECG gain to X1.4.
- 2. Connect the function generator to the patient cable through the 1000:1 divider (as illustrated in Figure 3–15 on page 3–40).
- 3. Set the function generator for a 1Vp-p, 10Hz sinewave.
- 4. Decrease the function generator signal frequency to 1Hz. Confirm that the trace displayed on the CRT is  $19 \pm 2mm p-p$ .
- 5. Increase the function generator signal frequency to 40Hz. Confirm that the trace displayed on the CRT is 19 ±2mm p-p.
- 6. Disconnect the function generator.

#### **TCP-ECG Output**

- Press LEAD SELECT to select LEAD II.
- 2. Connect the function generator to the patient cable through the 1000:1 divider (as illustrated in Figure 3–15 on page 3–40).
- 3. Set the function generator for a 1Vp-p, 10Hz sinewave at the output of the 1000:1 divider.
- 4. Connect a 0.25-inch phone plug to the ECG output on the rear panel of the instrument. Connect the oscilloscope to the phone plug.
- 5. Press  $\nabla$  ECG SIZE  $\triangle$  to set the gain to X1.0.
- 6. Confirm that the oscilloscope displays a  $1.00 \pm 0.15$  Vp-p sinewave.

### TCP-Pacemaker Calibration

- 1. Connect the pacemaker calibration test load as shown in Figure 3–16 using a X10 oscilloscope probe.
- Access the Test Menu (turn power off, then press ON while pressing and holding down LEAD SELECT and HR ALARM until the Test Menu is displayed).
- 3. Press ▼ ECG SIZE ▲ to scroll to PACER, then press HR ALARM to display the Pacer Menu. Press ▼ ECG SIZE ▲ to select 7 on the Pacer Menu.
- 4. Rotate the CURRENT knob clockwise to start pacer pulses. Confirm that the PACER LED flashes with pacer pulses.
- 5. Adjust the CURRENT knob to set the amplitude of the leading edge of the pacer pulse displayed on the oscilloscope to 7V (7mA output from the pacemaker).

- 6. Press HR ALARM to enter this value.
- 7. Press ▼ ECG SIZE ▲ to select the next value on the Pacer Menu.

  Repeat steps 4 through 6 above for each value listed on the Pacer Menu.

  The pacemaker peak current output values correspond to the following voltages displayed on the oscilloscope:

7mA =	7V	$(\pm 1.75)$
50mA =	50V	(±2.0, 4% tolerance)
100mA =	100V	(±4.0, 4% tolerance)
150mA =	150V	(±6.0, 4% tolerance)
200mA =	200V	$(\pm 8.0, 4\% \text{ tolerance})$

8. Press LEAD SELECT to exit the Pacer Menu and return to the Test Menu. Press ▼ ECG SIZE to scroll to EXIT, then press HR ALARM to resume normal operation. Press ON to turn instrument power off.




Figure 3–16 Pacemaker Calibration Test Load Setup

### TCP-Pacemaker Current 1. and Waveform Output

- Press ON and confirm the instrument powers on. Connect the pacemaker calibration test load as shown in Figure 3–16 using a X10 oscilloscope probe.
- 2. Press PACER to turn pacemaker power on, then press ▼ RATE ▲ to select 170ppm pacemaker output.
- Rotate the CURRENT knob to select the current output levels listed in Table 3–9 and confirm the voltages and pulse widths are within the listed tolerances. Refer to Figure 3–17 for pulse waveform characteristics.

Table 3-9 Current and Waveform Levels

Current	7mA	50mA	150mA	200mA
Vmax (nominal)	7V	50V	150V	200V
Vmax (range)	4-10V	46.5-53.5V	139.5-160.5V	186-214V
Vmin (range)		Within 15 to 2	25% of measured Vi	nax
Tpw (range)	19–21ms	19-21ms	19-21ms	19–21ms




Figure 3-17 Pacing Pulse Waveform

TCP Che	ecklist	Customer:	***************************************	***************************************		Date:	***************************************
	<del> </del>	•				>	,
		Instrument	Type:			Model:	
		Serial Num	ber:			<u> </u>	
		Frequency:				Last Insp:	
		Tech. Name	e/ID:				
Test and C	alibration	Pas	s Fail	Comm	ents	·	
Voltage Che	ecks	. 0					
	y Threshold Cl		_				
M01 Defi	nitoring Check ib Charge Checl	□ k □					
Test Menu	_	<b>к</b> ⊔	<u> </u>				
	r Calibration						
Selected	Acceptabl	-	MEASURE	and the second second			MEASURED
Value (J)	13037 (5)	% Tolerance	VALUE (			% Tolerance	VALUE (J)
1	0.8 to 1.2	±20%		20	19.0 to 21.0	±5%	
2	1.8 to 2.2	±10%		30	28.5 to 31.5	±5%	<del></del>
3	2.7 to 3.3	±10%	·	50	47.5 to 52.5	±5%	
4	3.8 to 4.2	±5%		100	95.0 to 105.0	±5%	
5	4.8 to 5.2	±5%		200	190.0 to 210.0	±5%	
6	5.7 to 6.3	±5%		300	285.0 to 315.0	±5%	<del></del>
7	6.7 to 7.3	±5%		360	342.0 to 378.0	±5%	
8	7.6 to 8.4	±5%	<del></del>	.			
9	8.6 to 9.4	±5%			•		
10	9.5 to 10.5	±5%					
Test Load (	Calibration						
<b>Charge Tin</b>	ne						
<b>Output Wa</b>	veform						
Display Cal	libration						
	ply Calibration	1 🗆	ı 📮				
Recorder C	Calibration		t o				
Recorder F	requency Resp	onse 🗆					
-	ency Response	<b>:</b> □	ı				
<b>ECG Outp</b>					· · · · · · · · · · · · · · · · · · ·		· · · · · · · · · · · · · · · · · · ·
	Calibration					<del></del>	
Pacemaker Waveform	Current and Output						
Curre	<b>ent (Vmax nominal</b>	l) 7mA (	7V)	50mA (50V)	150mA (150V)	200mA (2	00V)
	Vmax						
	Vmin	<u> </u>					<u> </u>

Tpw

## Troubleshooting Aids

The information in this section may help a service technician isolate possible software or circuit problems. It consists of the self-diagnostic error codes. If an instrument problem cannot be isolated and repaired, contact a Physio-Control customer service specialist in the US at 1-800-442-1142; outside the US, contact the local Physio-Control representative.

### Error Codes and SERVICE Indicator

The LIFEPAK 9P defibrillator/monitor/pacemaker contains self-diagnostic software. When a software or circuit problem is detected during operation, the **SERVICE** message is displayed on the monitor. For many problems, an error code is generated and displayed in the Info Menu.

To access the Info Menu, turn power off, then press ON while pressing and holding down LEAD SELECT and HR ALARM until the Test Menu is displayed. Press ▼ ECG SIZE to scroll down to INFO, then press HR ALARM to display the Info Menu. The error code number is listed beside FAULT STATUS. Refer to Table 3–10 for the software or circuit associated with each error code.

To remove or "clear" an error code, press HR ALARM to scroll down to RETAIN (or CLEAR), then press ECG SIZE ▲ to CLEAR the error. If more than one error code has been generated, the next error code will be displayed in sequence. After all error codes are cleared (and the FAULT STATUS displays 00), press LEAD SELECT to exit the Info Menu and return to the Test Menu. Press ▼ ECG SIZE to scroll down to EXIT, then press HR ALARM to resume normal operation.

After the error code is cleared, the **SERVICE** message is removed from the display. If, however, a software or circuit problem still exists or reoccurs, the **SERVICE** message is again displayed. Contact a Physio-Control representative for assistance.

Table 3–10 Error Codes

Error Code	Software or Circuit
1	Watchdog hardware (fast test)
2	Watchdog hardware (slow test)
3	Unknown transfer code
4	Not used
5	Not used
6	Not used
. 7	System external RAM test (U28)
8	System microprocessor RAM (U30)
9	System controller ROM (U27)
10	Critical ROM (U27)

Table 3-10 Error Codes (cont.)

Error Code	Software or Circuit
11	Display interface RAM (U13)
12	Display/recorder microprocessor RAM (U1)
13	Display interface ROM (U14)
14	D/A converter CRC
15	Transfer switching (Q4,CR4)
16	Transfer enable stuck on
17	Charge reference out of tolerance
18	Preamp baseline out of tolerance
19	Clock
20	+5V supply out of tolerance
21	+15V supply out of tolerance
22	Defib cal factors out of tolerance
23	A/D converter (U9) self-test
24	Recorder printhead resistance out of tolerance
25	Pace rate is excessive
26	Pace pulse is too wide
27	Pace current amplitude is too high
28	Pace current amplitude is too low
29–36	Problem with pacing error detection
37–39	Pacing parameter out of range
40	Pacemaker requires calibration, or unsuccessful calibration
41	Display ROM

2 30 0 0 2 30 0 0 3 30 0 0

		•		
				-
	•			

# **Service and Maintenance**

#### Introduction

This section provides the following information:

- Battery Maintenance
- General Maintenance
- Special Handling Procedures for Static Sensitive Devices
- Disassembly Procedures
- Inspection Techniques
- Tools and Materials for Cleaning and Repair
- Cleaning Procedures
- Printed Circuit Board Repair Precautions
- Preparation for Storage or Shipping.

#### Battery Maintenance

The Battery Pack in the LIFEPAK 9P defibrillator/monitor/pacemaker provides short-term interim power while transporting patients or during other brief interruptions of ac power. The defibrillator/monitor/pacemaker automatically switches to battery power if the power cord is disconnected or the ac power source fails.

#### Warning

POSSIBLE DEFIBRILLATOR SHUTDOWN. When operating on battery power, the large current draw required for defibrillator charging may cause the defibrillator to reach shutdown voltage levels with no low battery warning. If the defibrillator shuts down without warning, or if the LOW BATTERY: CONNECT AC POWER message appears on the monitor display, immediately connect the power cord to an ac power source.

#### **Recharging Batteries**

The defibrillator/monitor/pacemaker continually charges the installed Battery Pack whenever the instrument is connected to an ac power source, even with the front panel power switched off (although the rear panel ac mains power

switch must be on). The **BATT CHRG** message appears on the lower left corner of the monitor display during battery charging. A fully depleted Battery Pack is recharged to 90% capacity in three hours, and to full capacity in 24 hours. After full charging, the instrument continues to supply a trickle-charge which maintains the charge on the Battery Pack.

### **Battery Maintenance Guidelines**

The Battery Pack contains sealed lead-acid batteries which must be properly maintained to maximize battery life and performance. To properly maintain Battery Packs, follow these guidelines:

- Whenever ac power is available, keep the defibrillator/monitor/ pacemaker connected to ac power (with the rear panel mains power switch on). This charges the Battery Pack as indicated by the BATT CHRG message on the monitor display.
- When the LOW BATTERY: CONNECT AC POWER message appears on the monitor display, immediately connect the ac power cord to an ac power source (with the rear panel mains power switch on) to continue operation and begin recharging the Battery Pack. Frequent use of the Battery Pack when it is at minimum reserve capacity will reduce battery life. If the low battery message occurs frequently, this may indicate that the Battery Pack needs to be replaced.
- Perform the Battery Capacity Check (p. 4–2) on the installed Battery Pack a minimum of every six months or whenever you suspect the Battery Pack may not be functioning properly.
- Replace the Battery Pack installed in the defibrillator/monitor/ pacemaker a minimum of every two years.

#### Guidelines for Replacement Batteries

End of battery life is inevitable. As batteries age, their charge capacities diminish. The Battery Pack installed in the defibrillator/monitor/pacemaker should be replaced a minimum of every two years as a preventive maintenance practice. To maximize the life and performance of replacement Battery Packs, follow these guidelines:

- Perform the Battery Capacity Check on any replacement Battery Pack at a minimum of every six months. (If it passes the Battery Capacity Check, be sure to recharge the Battery Pack for 24 hours before returning it to storage.)
- Recharge any new replacement Battery Pack as soon as you receive it before storing it or using it.
- Store replacement Battery Packs at 20° to 25°C (68° to 77°F).

#### **Battery Capacity** Check

To check the capacity of a Battery Pack, perform the following steps (this procedure requires the use of a Defibrillator Energy Meter with a power range of 0-400J):

- Install the Battery Pack in a LIFEPAK 9P defibrillator/monitor/ pacemaker and charge the Battery Pack for at least 24 hours (if the Battery Pack is already fully charged, go on to the next step).
- 2. Disconnect the defibrillator/monitor/pacemaker from the ac power source.
- Make sure the instrument power is on and it is operating from battery power (the BATT CHRG message is not displayed). Place the standard defibrillator paddles on the energy meter.

#### Warning

SHOCK HAZARD. When charged and discharged as described in this procedure, the LIFEPAK 9P defibrillator/monitor/pacemaker discharges approximately 360J of electrical energy through the defibrillator paddles. Unless discharged safely as described in this procedure, this electrical energy may cause injury or death. Do not attempt to perform this procedure unless you are thoroughly familiar with the operation of the LIFEPAK 9P defibrillator/monitor/pacemaker.

- Charge the defibrillator to 360J and discharge into the energy meter. Repeat this 4 more times for a total of five (5) 360J discharges within 2 minutes.
- If the instrument shuts down due to power loss or if the LOW BATTERY: CONNECT AC POWER message appears, the Battery Pack has failed this procedure. Discard the Battery Pack and discontinue the procedure.

If the instrument continues to operate after step 4, leave the defibrillator/ monitor/pacemaker operating on battery power for 45 minutes.

At the end of the 45 minutes, repeat step 4.

If the instrument shuts down due to power loss at any time, the Battery Pack has failed this procedure. Discard the Battery Pack.

If the instrument continues to operate on battery power after repeating step 4, the Battery Pack has passed the procedure and has useful life remaining. (It is acceptable for the LOW BATTERY: CONNECT AC POWER message to appear and remain on at this point.) Recharge the Battery Pack for 24 hours before placing into service or storage.

**Note:** Physio-Control Corporation recommends replacing the Battery Pack installed in the LIFEPAK 9P defibrillator/monitor/pacemaker after two years of service as a part of routine maintenance.

#### **Recycling Batteries**

In the United States, recycle Battery Packs locally according to national, state, and local regulations when Battery Packs are no longer useful. If local recycling is not possible, contact Physio-Control customer service specialists at 1–800–442–1142 for information on returning Battery Packs.

Outside the United States, recycle Battery Packs according to local regulations if possible. Otherwise, contact the local Physio-Control representative for information on returning Battery Packs.

#### General **Maintenance**

The following paragraphs contain routine maintenance procedures that can be performed either by an experienced operator or a technician.

#### **Operation From AC or DC Power Source**

To operate from ac power, plug the ac power cord of the instrument into a three-pin, grounded, ac power source. Do not use a free-ground adapter. The defibrillator/monitor/pacemaker operates on ac power whenever it is plugged in as described and the ON LED is on (the rear panel ac mains power switch must be on). Whenever ac power is connected, the Battery Pack is being recharged as indicated by the **BATT CHRG** message in the lower left corner of the monitor display.

To operate the instrument from the internal dc Battery Pack, disconnect the ac power cord from the ac power source and make sure the front panel power switch is pressed on.

#### Strip Chart Recorder

Maintenance on the Strip Chart Recorder consists of loading paper and cleaning. Steps for these procedures follow, parenthetical numbers refer to parts shown in Figure 5-9, page 5-58.

The message NO PAPER will appear when the recorder paper runs out.

#### Caution

Prevent Printhead Damage. Use only Physio-Control recorder paper (PN804700) and do not operate the recorder without paper. Do not use waxed paper as this can damage the printhead.

#### **Loading Paper**

- Pull the latch (14) up and remove the empty paper spool.
- 2. Insert the new paper roll with the grid facing toward the front. Pull up a 2-inch length of paper.
- 3. Push the door (7) forward and the latch down until it latches.

#### Cleaning

Clean the printhead and paper-out sensor periodically to remove paper dust and debris. At a minimum, clean after every 100 rolls of use. Perform these steps for cleaning:

- Pull up the latch (14) and remove the paper roll.
- Gently wipe the surface of the printhead and paper sensor with a cotton swab soaked in isopropyl alcohol. When cleaning the printhead, touch the printhead only with the cotton swab, not with your fingers.

#### **Protecting Printouts**

Protect thermal paper printouts by following these guidelines:

- Do not apply tape or other adhesives over annotations or tracings on the front side of the paper. Apply adhesives only to the back of the paper.
- Store printouts only in paper folders. Do not store printouts in plastic folders.
- Avoid extended exposure to sunlight, temperatures above 27° C (81° F), or humidity above 70 percent.

#### **Fuse Replacement**

The fuses are located in the rear panel of the instrument.

- Place a thin, flat-blade screwdriver in the recessed area under the retaining tab of the fuseholder. Twist the screwdriver to release the retaining tab.
- 2. Pull the fuseholder out. Replace fuses with the same type, voltage rating, and current rating. Refer to the parts list (pg. 5-8).

#### Special Handling Procedures For Static Sensitive Devices (SSDs)

Many electronic semiconductor devices such as MOS ICs, FETs, optical isolators, or film resistors can be damaged by the discharge of static electricity. It is very common for <u>anyone</u> to build up static charges, especially when wearing synthetic clothes, and to transfer this charge to any object touched. Such a discharge can damage or destroy Static Sensitive Devices (SSDs). In most cases, the discharge is not even perceptible to the person at fault.

In order to help prevent damage to SSDs due to static discharge, observe the following precautions when performing any open-case test, maintenance, or repair procedures:

#### **Look for SSD Symbol**

Note that SSDs such as PCBs are indicated in this service manual with the following warning:

### **Use Static-Dissipative Mat**

Always perform repair or maintenance on a static-dissipative mat connected to earth ground.

#### **Wear Wrist Strap**

Always wear a conductive wrist strap connected to the mat and to ground except when working on energized equipment or when discharging high voltage circuits. The strap must be snug enough to make good contact against bare skin.

#### Warning

**SHOCK HAZARD.** Remove the wrist strap when working on energized equipment or when discharging high voltage circuits.

### Transport and Store PCBs Properly

Transport and store PCBs in anti-static racks or inside conductive bags. Label the package containing the PCBs as static-sensitive.

#### Keep Work Area Static-Free

Keep static-generating products such as styrofoam cups or trays away from the work area. Connect all electrical equipment such as soldering irons and test equipment to ground through a three-prong plug.

# Test Work Area Routinely

Test all the anti-static parts of the work area (mat, straps, cables) routinely. Keep a log of the test results.

### Disassembly Procedures

The following procedures provide a logical sequence for disassembling the major components of LIFEPAK 9P defibrillator/monitor/pacemaker. Separate or disassemble only to the extent required. Parenthetical numbers refer to item numbers in the final assembly parts list (beginning on page 5–6) and Figure 5–1 (pg. 5–11). Reassembly procedures are the reverse unless otherwise noted.

#### **Warning**

**SHOCK HAZARD.** Disconnect the AC power cord and remove the Battery Pack before disassembly.

#### **Caution**

Avoid Component Damage. This instrument contains static sensitive devices (SSDs). Use the following special handling procedures for SSDs. SSDs are indicated in Section 5 with this symbol.

#### **Battery Pack Removal**

1. Open and remove the battery door (14) at the bottom of the instrument.


Note: When installing a replacement Battery Pack, be sure to charge the Battery Pack for 24 hours. Refer to Battery Maintenance on page 4-1.

- 2. Remove the Battery Pack (A13).
- 3. Free the Battery Pack by pressing the locking tab on the battery harness and separating the two connectors.

#### Front Panel Bezel and Keypad Assembly Removal

- 1. Press in on the right edge of the front panel bezel (2) as illustrated in Figure 4–1 (arrow A).
- 2. While pressing the edge, push out in the direction of arrow B until the right edge of the bezel releases.
- 3. Slide the bezel to the right slightly (arrow C) until the left tab clears the front case (6).

- 4. Place fingers under the released edge and pull the bezel away from the case (arrow D).
- 5. Press the quick-release levers and disconnect J20 from the Keypad Assembly (A6).
- 6. Remove the two screws (44) that hold the Keypad Assembly and CRT implosion shield (51) to the bezel. Slide out the shield and the Keypad Assembly.
- 7. To separate the Keypad Assembly from the bezel, loosen the allen screw (47) on the shaft of the pacing current switch (34). Next, use a nut driver to remove the nut and washer from the shaft.






Figure 4-1 Front Panel Removal

### **Strip Chart Recorder Removal**

- 1. Pull the upper edge of the top door to open the recorder (A8).
- 2. Remove the paper roll.
- 3. Push down the paper carrier and remove the two screws (44) at the rear of the paper housing.
- 4. Grasp the recorder at the bottom then pull down and away from the instrument.




Figure 4–2 Connector Locations

#### **Case Separation**

- 1. Remove the Battery Pack, front panel bezel, and strip chart recorder.
- 2. Turn the paddles connector counterclockwise, then pull to disconnect. Remove the paddles.
- 3. Remove three screws (44) from the bottom of the instrument, four screws (44) from the rear panel (33), and two screws (44) from the handle.
- 4. Pull the front case (6) forward and remove one of the screws on the pacing connector bracket (3).
- 5. Pull the chassis (8) away from the rear case (7).

# Power Conversion PCB Removal and Reassembly

#### Caution

**Possible Equipment Damage.** Unless capacitors C7, C8, C24 and C25 on the Power Conversion PCB (A3) are discharged before reassembling the Power Conversion PCB, it is possible for surge current from the capacitors to cause damage to U6 on the Interconnect Pacing PCB during reassembly. To help avoid possible component damage, be sure to follow the procedure described for Power Conversion PCB Reassembly.

#### Power Conversion PCB Removal

Complete Case Separation. Refer to Figure 4–2, page 4–9:

- 1. Disconnect J17 and J18, located on the Power Conversion PCB (A3).
- 2. Push the Interconnect PCB (A5) locking tabs up to clear the notches in the Power Conversion PCB.
- Pull the Power Conversion PCB away from the Interconnect PCB to disconnect J14 and lift the Power Conversion PCB out of the board guide.
- 4. To separate the rear panel from the chassis (8), remove the four screws (44) from the rear panel (33) and disconnect J6, J29, J1, and J19.
- 5. Disconnect J32 (wire labeled 2, goes to wave shaping inductor) and J9 (wire labeled 1, goes to Transfer Relay). Cut the tie wrap on the wire bundle to completely separate the PCB.
- 6. Disconnect J10 through J13, then J15 and J16.

Power Conversion PCB Reassembly

To help prevent surge current from capacitors C7, C8, C24 and C25 on the Power Conversion PCB (A3) from damaging U6 on the Interconnect Pacing PCB when reassembling A3, perform the following steps:

- Make sure the AC power cord is disconnected and the Battery Pack is removed.
- Using test clips, connect a 47 $\Omega$ , 1/2W resistor (PN 200471-040) between J13 Pin 5 and the CR2 anode (battery ground) on the Power Conversion PCB as shown in Figure 4-3. This will discharge capacitors C7, C8, C24 and C25 which are connected in parallel. Leave the  $47\Omega$  resistor connected for at least 15 seconds, then disconnect the  $47\Omega$  resistor.
- Connect a digital voltmeter across the terminals of capacitor C7 or C8 on the Power Conversion PCB. Confirm that the measured voltage is less than 0.5V. If the voltage is greater than 0.5V, repeat steps 2 and 3.
- Reassemble the Power Conversion PCB by performing in reverse order the steps for Power Conversion PCB Removal.




Figure 4-3 Connections for Discharging Capacitors on the Power Conversion PCB Assembly

#### **Main PCB Removal**

Complete Case Separation. Refer to Figure 4–2, page 4–9:

- 1. Press down the quick-release levers and disconnect J4 and J5.
- 2. Disconnect J6, J29, and J35 on the Main PCB (A1).
- 3. Push the Interconnect PCB (A5) locking tabs up off of the notches in the Main PCB. Disconnect J22.
- 4. Pull the Main PCB away from the Interconnect PCB to disconnect P7 and lift the Main PCB out of the board guide.
- 5. Disconnect J8 and J9.

#### Display PCB/CRT Assembly Removal

Complete Case Separation. Refer to Figure 4–2, page 4–9:

- 1. Press down the quick–release levers and disconnect J901 on the Display PCB/CRT Assembly (A11).
- 2. Remove four screws (44) to remove the Display PCB/CRT Assembly.

### Power Supply Removal

Complete Case Separation. Refer to Figure 4–2, page 4–9:

- 1. Disconnect J1 on the Power Supply PCB (A2).
- 2. Remove four screws (44) from rear panel (33) and remove rear panel.
- 3. Disconnect J2 on the Power Supply.
- 4. Remove four screws (44) from Power Supply, detach from chassis (8).

### Interconnect Pacing PCB Removal

Complete Case Separation. Refer to Figure 4–2, page 4–9:

1. Remove the Power Conversion PCB, the Main PCB, the Display PCB/CRT Assembly, and the Power Supply PCB.

#### Caution

Possible Equipment Damage. Unless capacitors on the Power Conversion PCB (A3) are discharged before reassembling A3, it is possible for surge current from the capacitors to cause damage to U6 on the Interconnect Pacing PCB (A5) during reassembly. To help avoid possible component damage to A5, be sure to follow the procedure described for Power Conversion PCB Reassembly when reassembling A3.

2. Remove five screws (44) and disconnect J36 to separate the Interconnect Pacing PCB.

#### Inspection **Techniques**

When servicing or repairing the instrument, routinely begin with a visual inspection of the hardware and components for signs of damage. Signs of damage require further inspection of the surrounding area for peripheral breakage or damage.

#### **Exterior Inspection**

Visually inspect the entire instrument for wear, corrosion, deterioration, and damage resulting from extreme temperatures or dropping. Lift and hold the instrument upside down while listening for loose hardware.

#### **Interior Inspection**

Table 4-1, Inspection Techniques, lists the major hardware components of the instrument, possible problems associated with each item, and recommended corrective actions.

**Table 4–1 Inspection Techniques** 

Hardware	Problem	Corrective Action
Chassis, covers, and brackets	Warped, bent, damaged surfaces, or missing hardware	Replace
Components (mechanical and	Loose mountings	Repair
electrical)	broken or damaged leads	Replace
	Deterioration or leakage	Replace
Connector pins	Slightly bent	Straighten
	Badly bent, loose or corroded	Replace
Nameplate, labels, and decals	Not legible	Replace
PCB surfaces	Charred, cracked, or brittle	Replace PCB
	•	PCB surface can be expected due to ratures of some components.
Screws and nuts	Loose or cross-threaded	Tighten or replace
Terminals and connections	Installed incorrectly	Install correctly
	Missing or worn	Replace
	Failed solder connections	Resolder (if not damaged) or replace
Wire insulation and tubing	Deteriorated, worn, pinched, or damaged	Replace

#### Tools and Materials for Cleaning and Repair

Table 4–2 lists recommended tools, materials, and chemicals required for cleaning and repair activities. Although specific items are recommended, tools and materials with specifications equivalent to those listed may be used.

Table 4–2 Tools and Materials for Cleaning and Repair

Product	Description
Static-protected work area	Grounded conductive surface and wrist strap
Electronic Supply Dealers	
X-32B and Xersin <sup>®</sup> solder	Low flux, no cleaning necessary
Multicore	
Westbury, NY 11590	
5-Minute <sup>®</sup> Epoxy	Adhesive
Devcon	•
Danvers, MA 01923	
Tak Pak <sup>®</sup>	Adhesive and 710 accelerator
Loctite	and the second s
Newington, CT 06111	<u>.</u>
Crocus cloth	
Isopropyl alcohol	
Acid brush	
Cotton swabs	
Vacuum cleaner	
Soft-bristle brush	Nonmetallic
Cloth	Clean and lint-free
Compressed air	Clean and dry (60psi, max.)

Cleaning Procedures	Clean the instrument after maintenance, disassembly, or any repair procedure.
External Cleaning	Clean the instrument case, cables, and screen with mild soap and water. Use only a damp sponge or towel to clean.
Warning	SHOCK OR FIRE HAZARD. Do not immerse any portion of the instrument in water. Fluid spills and splashes may damage the instrument's electrical components.

Caution

Possible Instrument Damage. Do not autoclave the instrument.

#### **Interior Cleaning**

Follow disassembly procedures (beginning on page 4–7) to access the interior.

#### Warning

PERSONAL SAFETY HAZARD. Ventilate work area when using solvents. Observe manufacturer warnings regarding personnel safety and emergency first aid. Keep first aid equipment available when using chemicals.

#### **Warning**

CHEMICAL FIRE HAZARD. Observe shop safety and fire precautions. Store solvents and solvent-soaked rags in approved containers. Refer to manufacturers instructions on containers for recommended firefighting procedures. Keep firefighting equipment available.

#### Caution

Possible Component Damage. This instrument contains static sensitive devices (SSDs). Use special handling procedures on page 4-6.

#### Caution

Possible Component Damage. Do not use solvents to clean transformers or inductors.

- Brush surfaces and parts with a nonmetallic soft-bristle brush. Remove loosened dirt and dust using dry low pressure compressed air (60psi) or a vacuum cleaner.
- Wipe metal surfaces with a soft, nonabrasive cloth dampened with isopropyl alcohol.

#### Caution

Prevent Nameplate Damage. Do not use abrasive cleaners or solvents to wipe nameplates and labels—the nameplate information may disappear.

Wipe surfaces of nameplates and labels with a clean, dry cloth.

#### Caution

**Possible Component Damage.** Do not use solvents to clean plastic parts.

- 4. Clean surfaces of plastic parts with a mild soap and water solution. Dry with a clean cloth.
- Clean soldering surfaces with a nonmetallic, soft-bristle brush dipped in isopropyl alcohol then wipe with absorbent cloth. Air dry 10 minutes or use low pressure compressed air (60psi) before soldering.

#### PCB Repair Precautions

The LIFEPAK 9P defibrillator/monitor/pacemaker uses multilayer PCBs which are not easily repaired. They consist of alternate layers of conductive patterns and insulating material bonded together and interconnected with plated-through holes. The internal layers and connection points are quickly damaged by too much heat, rendering the entire PCB useless. Therefore, it is recommended that all multilayer PCBs be returned to Physio-Control for repair or replacement.

#### Caution

**Possible PCB damage.** Improper handing can easily damage the PCBs beyond repair. Plated-through holes connecting the circuitry on two sides of the PCB can be damaged by too much heat.

**Note:** Before attempting any PCB repair, contact Physio-Control. Customer PCB repair may endanger any applicable PCB exchange agreement or other warranty.

#### **Preparation for** Storage or Shipping

Save the original shipping box and packing for the LIFEPAK 9P defibrillator/monitor/pacemaker. If the instrument must be shipped to the service center or factory, the special packaging is required to prevent shipping damage. Refer to Figure 4-4 when repacking.




Figure 4-4 LIFEPAK 9P defibrillator/monitor/pacemaker Shipping Assembly

		 · · · · · · · · · · · · · · · · · · ·	
			٠.,
	,		
`			
•			
			•

#### Introduction

This section includes lists of replaceable parts, component layouts or exploded views, and schematic diagrams for each major assembly of the LIFEPAK 9P defibrillator/monitor/pacemaker and the optional Defibrillation Adapter. This information is provided to assist in troubleshooting and repair. Table 5–1 lists the figure and page number of the major assemblies for the defibrillator/monitor. Table 5–2 lists the figure and page number of the major assemblies for the defibrillator/monitor/pacemaker. Table 5–3 lists standard reference designators for components. Supplies and accessories are listed in Table 5–4.

#### **Parts List**

Use parts lists to identify *replaceable* parts. The format for each list is the same:

Ref: This column contains the reference designators or the assigned item numbers of replaceable parts. Reference designators are abbreviated (see Table 5–3) and are listed alphanumerically by this abbreviation. Assigned item numbers are listed numerically, starting with 2. The number in the heading of this column indicates the figure to which the parts list corresponds.

**Part Number:** Physio-Control Corporation part number is listed in this column. Part ordering information is provided on page 5-2.

**Description:** Descriptive information for a part that is used more than once refers to the initial listing with a "Same as \_\_\_\_" citation. This citation may also list a component value when necessary. Static Sensitive Devices (SSDs) are identified in this column by this symbol:

For information about proper handling of SSDs during repair or replacement, refer to page 4–6.

Use Code: Different hardware configurations and their corresponding letter codes are listed at the beginning of each parts list. Subsequent use of any letter code indicates which part is used in a particular configuration. A blank in this column indicates that the part is used in all configurations.

**Qty:** The quantity column specifies the *total* quantity used for each part number listed. The abbreviation REF indicates a major assembly or subassembly.

### **Component Layouts or Exploded Views**

Use component layouts or exploded views to locate parts. The reference designator or assigned item number on the drawing corresponds to the REF column entry on the parts list.

#### **Schematic Diagrams**

Schematic diagrams depict circuit function and signal flow. Signal names, pin numbers, and polarities are labeled on the schematics as well as reference designators and values from the parts list.

- In general, signal flow is from left to right
- Circuit names on schematics correspond to descriptions in Section 1
- ICs marked with a \* are detailed in Section 6
- To/From labels correspond to connections shown on the interconnect diagram.

This manual documents the revision levels of the assemblies which were in production at the time of publication. Since the time of publication, it is possible that modifications to the assemblies have occurred which are not documented in this manual. For assistance in obtaining updated information in the US, call PARTSLINE at 1-800-442-1142. Outside the US, contact your local Physio-Control representative.

### How to Order Parts

When ordering from Physio-Control Corporation in the US, call the toll free PARTSLINE number (1–800–442–1142). Give the instrument part number and serial number located on the bottom panel of your instrument. Then specify all assembly numbers, part numbers, reference designations, and descriptions. Different parts may be substituted by Physio-Control to reflect modifications and improvements of instrument circuitry. To order parts outside the US, contact your local Physio-Control representative.

Table 5-1 LIFEPAK 9P defibrillator/monitor/pacemaker Major Assemblies

Ref	able 5–1 LII Part Number	FEPAK 9P defibrillator/monitor/pacemaker Major Nomenclature	Fig. No.	Page No.
REF	805460	LIFEPAK 9P defibrillator/monitor/pacemaker Final Assembly	5-1	5-11
KEI	003400	Interconnect Diagram	5-2	5–11 5–12
A1	805474	Main PCB Assembly	5–2 5–3	5–12
A2	803726	Power Supply PCB Assembly	ŀ	ł
A2 A3	803724		5–4 5–5	5-31
		Power Conversion PCB Assembly	5–5	5–40
A4		Not used		
A5	805472	Interconnect Pacing PCB Assembly	5–6	5–50
A6	805469	Keypad Assembly	5–7	5–55
A7	805464	Rear Panel Membrane Switch	5–8	5–56
A8	804189	Strip Chart Recorder Assembly 🛕	5–9	5–58
A9	800240	Transfer Relay Assembly	5–10	560
A10	802901	Paddle Assembly	5–11	5–61
A11	803706	Display PCB/CRT Assembly	5–12	5–64
A12		Chassis/Harness Assembly	N/A	N/A
A13	803704	Battery Pack (shown on Final Assembly)	N/A	N/A
W1	803783–30	Power Supply PCB/AC Receptacle/Ground Wire Harness	5-13	5–73
$ _{W_2}$	803783–23	Power Conversion PCB/Power Supply PCB/Battery Wire Harness	5-14	5-74
W3	803783-24	Battery Sensor/Power Supply PCB Wire Harness	5–15	5-75
W4	803756–13	High Voltage Connector Cable Assembly	5–16	5–76
W5	803733-00	Main PCB/Display Cable Assembly	5–17	5–77
W6	803773-00	Main PCB/Recorder Cable Assembly	5–18	578
W7	805465-00	Interconnect Pacing PCB/Keypad Cable Assembly	5–19	5–79
W8	803783-00	Main PCB/Patient Connector Wire Harness	5–20	5–80
W9	803783-04	AC Receptacle/Ground Wire Harness	5-21	5-81
W10	80378305	AC Receptacle Wire Harness	5–22	5-81
W11	803783-06	Test Load Contact/Resistor Wire Harness	5–23	582
W12	803783–07	Power Conversion PCB/Reed Assembly Wire Harness	5–24	5–82
W13	803783-08	Power Conversion PCB/Dump Relay Wire Harness	5–25	5-83
W14	803783-09	Power Conversion PCB/Charge Relay Negative Wire Harness	5–26	5–84
W15	803783–10	Power Conversion PCB/Charge Relay Positive Wire Harness	5–27	5–84
W16	80378311	Main PCB/Sternum Relay Wire Harness	5–28	5–85
W17	803783–12	Main PCB/Apex Relay Wire Harness	5-29	5–85
W18	803783–13	Storage Capacitor Positive/Relay Wire Harness	5–30	5–86
W19	803783–14	Storage Capacitor Negative/Relay Wire Harness	5–31	5-86
W20	803783–15	ECG Out Wire Harness	5–32	5–87
W21	803783–16	Main PCB/Speaker Wire Harness	5–33	5–88
W22	803783–29	Power Supply Bracket/Ground Wire Harness	5–34	589
W23	802925–03	Connector Receptacle/Pacing Cable Assembly	5–35	5–90
W24	803783–19	Recorder/Bracket Ground Wire Harness	5–36	5-91

Table 5-2 Defibrillation Adapter Assemblies

Ref	Part Number	Nomenclature	Fig. No.	Page No.
REF	803747	Defibrillation Adapter	537	5–93
W1	803774–03	Test Load High Voltage Wire Harness	5–38	5–95
W2	803774-04	Test Load High Voltage Wire Harness	5–39	5–95
W3	803774-05	Pushbutton Switch Wire Harness	5-40	5–96
W4	803756-09	High Voltage Connector Cable Assembly	5–41	5–97

**Table 5–3 Reference Designator Key** 

Designator	Description
A	Main Assemblies
C	Capacitor
CR or D	Diode
F	Fuse
Н	Heat Sink
J	Jack Connector
JP	Jumper Wire
K	Relay
L	Inductor
P	Plug Connector
Q	Transistors
R	Resistor
RN	Resistor Network
RV	Varistor
RT	Thermistor
S	Switch
T	Transformer
TP	Test Point
U	IC (Integrated Circuit)
VSP	Voltage Surge Protector
W	Wire Harness
Y	Crystal
XU	IC Socket

# Supplies and Accessories

Supplies and accessories for the LIFEPAK 9P defibrillator/monitor/ pacemaker are listed in Table 5-4 with the Physio-Control part number. To order an item in the US, call the Physio-Control PARTSLINE at 1-800-442-1142. Outside the US, contact your local Physio-Control representative.

Use of non-Physio-Control defibrillation electrodes, batteries, accessories, or adapter devices may void Safety Agency Certifications and warranty.

Table 5-4 Supplies and Accessories

Description	Quantity	Part Number	
Adapter, Defibrillation		803747	
Adapter, Shock Advisory		803732	
Cable, Defibrillation (for use with FAST-PATCH disposable defibrillation/ECG electrodes)		804089	
Cable, Patient	3-lead, AHA 3-lead, IEC	9–10418 800947	
Cable, Pacing		802905	
Cord, Power		803619	
DERMA JEL Electrode Jelly	4-ounce tube Case (12 tubes)	9–10236–00 9–10236–012	
Electrode, FAST-PATCH disposable defibrillation (use with Defibrillation Adapter for LIFEPAK 9 defibrillator/monitor, P/N 803747)	1 set (2 electrodes) 10 sets (20 electrodes) 50 sets (100 electrodes)	804545-001 804545-010 804545-050	
Electrode, LIFE-PATCH, pre-gelled (disposable)	Box of 30 (10 pkgs/box) Case of 300 (10 boxes/case)	800139-030 800139-300	
Electrode, QUIK-PACE disposable external pacing	One set (two pacing electrodes) Five sets (10 pacing electrodes) Ten sets (50 pacing electrodes)	803377-101 803377-501 803377-251	
<b>Operating Instructions</b>		805455	
Paper, Recorder, 50 mm x 30m (100ft)	Box of 3 rolls Case (150 rolls)	804700–003 804700–150	
Paddle, Pediatric, external, 2 required		800418	
Paddle, Posterior, external, adult		802461	
Service Manual		805454	
Shipping Container for Battery Recycling		805984	

5-1	Part			
Ref	Number	Description	Use Code	Qty
		LIFEPAK 9P defibrillator/monitor/pacemaker Final Assembly:		
	805460-01	English, CSA	A	REF
	805460-02	English, UK, 220V	В	REF
	805460-03	French, CSA, 110V	C	REF
	805460-04	French, 220V	D	REF
	805460-05	German, 220V	E	REF
	805460-06	Spanish, 120V	F	REF
	805460-07	Spanish, 220V	G	REF
	805460-08	Italian, 220V	н	REF
	805460-09	Swedish, 220V	I	REF
	805460-10	English, Intl., CSA, 120V	J	REF
A1	805474-01	Main PCB Assembly		1
		Note: If ordering a Main PCB Assembly, order U14 and U27 separately. Refer to Main PCB Assembly parts list.		
A2	803726-01	Power Supply PCB Assembly		1
<b>A3</b>	803724-03	Power Conversion PCB Assembly		1
<b>A</b> 5	805472-02	Interconnect Pacing PCB Assembly		1
<b>A6</b>	805469-10	Keypad Assembly, front, LP9P, English	A,J	1
	805469-11	Keypad Assembly, front, LP9P, French, CSA	С	1
	805469-12	Keypad Assembly, front, LP9P, English, Intl.	В	1
	805469-13	Keypad Assembly, front, LP9P, French	D	1
	805469-14	Keypad Assembly, front, LP9P, German	E	1
	805469-15	Keypad Assembly, front, LP9P, Spanish	F,G	1
	805469-16	Keypad Assembly, front, LP9P, Italian	н	1
	805469-17	Keypad Assembly, front, LP9P, Swedish	1	1
A7	805464-00	Rear Panel Switch Membrane		1
<b>A8</b>	804189-00	Strip Chart Recorder Assembly		1
<b>A9</b>	800240-14	Transfer Relay Assembly		1
A10	802901-06	Paddle Assembly		1
A11	803706-00	Display PCB/CRT Assembly	A,B,C,J	1
	803706-01	Display PCB/CRT Assembly	D-I	1
A13	803704-02	Battery Pack		1

5-1 Ref	Part Number	Description	Use Code	Qty
W1	803783-30	Power Supply/AC Receptacle/Ground Wire Harness		1
W2	803783-23	Power Conversion PCB/Power Supply/Battery Wire Harness		1
W3	803783-24	Battery Sensor/Power Supply PCB Wire Harness		1
W4	803756-13	High Voltage Connector Cable Assembly		1
W5	803733-00	Main PCB/Display PCB Cable Assembly		1
W6	803773-00	Main PCB/Strip Chart Recorder Cable Assembly		1
W7	805465-00	Interconnect Pacing PCB/Keypad Cable Assembly		1
W8	803783-00	Main PCB/Patient Connector Wire Harness		1
W9	803783-04	AC Receptacle/Ground Wire Harness	·	1
W10	803783-05	AC Receptacle/Ground Wire Harness		1
W11	803783-06	Test Load Contact/Resistor Wire Harness		2
W12	803783-07	Power Conversion PCB/Reed Assembly Wire Harness		1
W13	803783-08	Power Conversion PCB/Dump Relay Wire Harness		1
W14	803783-09	Power Conversion PCB/Charge Relay Negative Wire Harness		1
W15	803783-10	Power Conversion PCB/Charge Relay Positive Wire Harness		1
W16	803783-11	Main PCB/Sternum Relay Wire Harness		1
W17	803783-12	Main PCB/Apex Relay Wire Harness		1
W18	803783-13	Storage Capacitor Positive/Relay Wire Harness		1
W19	803783-14	Storage Capacitor Negative/Relay Wire Harness		1
W20	803783-15	ECG Out Wire Harness		1
W21	803783-16	Main PCB/Speaker Wire Hamess		1
W22	803783-29	Power Supply Bracket/Ground Wire Harness		1
W23	802925-03	Connector Receptacle/Pacing Cable Assembly		1
W24	803783-19	Recorder/Bracket Ground Wire Harness		1
2	805470-00	Bezel, front control panel		1
3	805441-00	Bracket Pacing Connector		1
4	803729-02	Bracket, Strip Chart Recorder		1
5	803705-00	Capacitor, storage		1
6	805467-02	Case, front, with inserts		1
7	803740-01	Case, rear, ECG Out		1
8	803722-06	Chassis		1
9	803737-00	Clip, mounting, test-load resistor		2
10	803718-00	Clip, power entry module	A,C,J	1
11	201393-011	Connector, plug, equipotential		1


5-1 Ref	Part Number	Description	Use Code	Qty
12	202110-000	Connector, receptacle, quick-connect, 6A/250V		1
13	803736-00	Contact, test load		2
14	803711-01	Door, battery		1
15	200363-025	Fastener, plastic rivet, 0.125 D		1
16	802885-00	Foot, mounting		4
17	200256-151	Fuse, slow blow, 4A/250V	A,C,F,J	2
	200619-018	Fuse, slow blow, 2A/250V	B,D,E,G,H,I	2
18	201605-011	Grommet, continuous, serrated (not shown)		1
19	802793-05	Inductor, waveshaping		1
20	802917-04	Knob, pacing current		1
21	805530-01	Label, pacing connector, English	A,J	1
	805530-02	Label, pacing connector, German	E	1
	805530-03	Label, pacing connector, Spanish	F,G	1
	805530-04	Label, pacing connector, Italian	н	1
	805530-05	Label, pacing connector, English, Intl.	В	. 1
	805530-06	Label, pacing connector, French	C,D	1
	805530-07	Label, pacing connector, Swedish	I	1
22	805466-07	Label, front panel, LP9P, English	A,B,J	1
	805466-08	Label, front panel, LP9P, French	C,D	1
	805466-09	Label, front panel, LP9P, German	Е	1
	805466-10	Label, front panel, LP9P, Spanish	F,G	1
	805466-11	Label, front panel, LP9P, Italian	н	1
	805466-12	Label, front panel, LP9P, Swedish	I	1
23	805461-02	Label, instruction, LP9P, French	С	1
	805461-04	Label, instruction, LP9P, French	D	1
	805461-05	Label, instruction, LP9P, German	Е	1
	805461-06	Label, instruction, LP9P, English	В	1
	805461-07	Label, instruction, LP9P, Spanish	F,G	1
	805461-08	Label, instruction, LP9P, Italian	н	1
	805461-09	Label, instruction, LP9P, Swedish	I	1
	805461-11	Label, instruction, LP9P, English, CSA	A,J	1
24	802940-12	Label, paddle, APEX, Italian	Н	1
	802940-11	Label, paddle, APEX, Intl.	A-G,I,J	1
25	802940-13	Label, paddle, STERNUM, Italian (not shown, used on A10)	Н	1
	802940-09	Label, paddle, STERNUM, Spanish (not shown, used on A10)	F,G	1

5-1 Ref	Part Number	Description	Use Code	Qty
	802940-10	Label, paddle, STERNUM, international (not shown, used on A10)	A-E,I,J	1
26	805463-03	Label, rear panel, LP9P, French	C.	1
	805463-04	Label, rear panel, LP9P, German, 220V	E	1
	805463-05	Label, rear panel, LP9P, Spanish	F,G	1
	805463-06	Label, rear panel, LP9P, Italian	Н	1
	805463-07	Label, rear panel, LP9P, Swedish	1	1
	805463-01	Label, rear panel, LP9P, UK, 220V	В	1
	805463-00	Label, rear panel, LP9P, English, CSA	A,J	1
	805463-02	Label, rear panel, LP9P, French, 220V	D	1
27	803727-87	Label, serial number, LP9P, French, CSA	С	1
	803727-84	Label, serial number, LP9P, English, CSA	Α	1
	803727-83	Label, serial number, LP9P, English, UK	В	1
	803727-97	Label, serial number, LP9P, English, CSA	l	1
	803727-86	Label, serial number, LP9P, French, 220V	D	1
	803727-88	Label, serial number, LP9P, German	E	1
	803727-89	Label, serial number, LP9P, Spanish	F	1
	803727-90	Label, serial number, LP9P, Spanish, Intl.	G	1
	803727-93	Label, serial number, LP9P, Italian	н	1
	803727-94	Label, serial number, LP9P, Swedish	I	1
28	801517-04	Label, UL (not shown, next to serial number label)	A	1
29	803456-15	Label, specification, battery (not shown, next to battery pack)		1
30	800943-09	Label, high voltage symbol, International (used on A9 and A11)		4
31	805660-00	Label, warning, hazard, English, Intl.	В	1
	805660-01	Label, warning, hazard, English, CSA	A,J	1
	805660-02	Label, warning, hazard, French	D	1
	805660-03	Label, warning, hazard, German	E	1
	805660-04	Label, warning, hazard, Spanish	F,G	1
	805660-05	Label, warning, hazard, Italian	н	1
	805660-06	Label, warning, hazard, Swedish	I	1
	805660-07	Label, warning, hazard, French, CSA	С	1
32	201508-002	Nut, kep, locking, #6-32		6
33	803713-03	Panel, rear, ECG Out		1
34	805603-00	Potentiometer, rotary		1
35	802935-05	Pushbutton, paddle, charge, domestic	A-D,J	1
	802935-06	Pushbutton, paddle, charge, German	E	1

5-1 Ref	Part Number	Description	Use Code	Qty
nei	802935-07	Pushbutton, paddle, charge, Spanish	F,G	1
	802935-07	Pushbutton, paddle, charge, Italian	H	1
	802935-18	Pushbutton, paddle, charge, Swedish	I	1
36	802935-19	Pushbutton, paddle, discharge, domestic	A,C,J	2
30	802935-08	Pushbutton, paddle, discharge, International	B,D-I	2
37	800516-02	Resistor, defibrillator (not shown, behind inductor)		1
38	802608-00	Resistor, $50\Omega$ , $50W$ , $5\%$	1	1
39	200536-001	Retainer, cable tie, nylon, 0.10 W x 4 L (not shown)		15
40	200536-011	Retainer, cable tie, nylon, 0.35 W x 20.9 L		2
41	200530-011	Retainer Clamp, cable, D-type, plastic		
42	200030-002	Retaining Ring, steel, 1.127 OD		1
43	201874-270	Screw, pan head, #4-40 x 0.312 L		4
44	201874-277	Screw, pan head, #8-32 x 0.375 L	A,C,J	39
77	201874-277	Screw, pan head, #8-32 x 0.375 L	B,D-I	38
45	201874-277	Screw, pan head, #8-32 x 0.500 L	-,	4
46	201874-276	Screw, pan head, #6-32 x 0.375 L		2
47	200543-005	Screw, set, #4-40 x 0.188 L		1 1
48	805557-01	Shield, Display PCB		1
49	803768-00	Shield, recorder bracket		1
50	803759-00	Shield, high voltage, Power Conversion PCB		1
51	803770-02	Shield, implosion, CRT		1
52	803717-01	Shield, preamp, circuit-side		1
53	803716-01	Shield, preamp, component-side		1
54	202250-003	Sleeve, ferrite snap (not shown, near connector		3
		receptacle on rear panel)		
55	803730-00	Spacer, battery pack (not shown, under battery door)		1
56	201625-171	Spring, compression, 0.300 OD, 0.562 L		2
57	803735-01	Spring, paddle retention		2
58	200192-294	Standoff, hex, 0.312 W x 1.875 L		2
59	201501-012	Tape, adhesive foam, 0.125 T x 1 W (not shown, under W12)		A/R
60	200283-005	<b>Tubing</b> , heat shrink, red, 0.250 ID (not shown, used on connector receptacle)		A/R
61	200804-018	Washer, flat, 0.375 OD, 0.152 ID, 0.045 T		2
62	200804-101	Washer, flat, #8		8

Parts Lists and Schematics 5





v

Fa-Sa         Part Number         Description         Use October Of Description	7			Capacitor, 0.1µF/100V, 20%	200049-039	C73
Description   Description	ယ			Capacitor, 0.022µF/50V, 5%	201628-008	C72
Description   Description   Code			***	(Same as C12)		C61-71
Description		·		(Same as C8)		C60
Description   Description   Code	<del> </del>			Capacitor, electrolytic, 100µF/25V, +50%	200205-028	C59
Description   Description   Code				(Same as C12)		C58
Part Number   Description   Code		, f		(Same as C8)		C57
Part Number   Description   Code				(Same as C12)		C53-56
Description   Description   Code		s.		(Same as C8)	434d	C52
Ref   Part Number   Description   Code				(Same as C12)		C43-51
Ref   Part Number   Description   Code			-	(Same as C8)		C38-42
53 Ref         Part Number         Description         Use Code           805474-011         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470µF/10V, 10%         1           201628-006         Capacitor, D01µF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         9           20102-032         Capacitor, electrolytic, 100µF/50V, 10%         1           2010701-006         Capacitor, electrolytic, 33µF/16V, 100%         1           2010701-006         Capacitor, electrolytic, 33µF/16V, 100%         1           2         201628-000         Capacitor, electrolytic, 33µF/16V, 100%         2           2         201628-000         Capacitor, 0.47µF/50V, 5%         2           4         200901-046         Capacitor, 0.01µF/50V, 5%         2           5         200901-046         Capacitor, 100µF/100V, 10%         2           6         200901-046         Capacitor, 100µF/100V, 10%         3           11,22         200803-051         Capacitor, 100µF/50V, 5%         3           12         200901-025         Capacitor, 10µF/50V, 5%         2           13         200901-026         Capacitor, 10µF/50V, 5%         3           14         201628-012	XUENT <b>PA</b>		·	(Same as C12)		C36, 37
53 Ref         Part Number         Description         Use Code           805474-θ1         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, D.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         9           2002102-032         Capacitor, electrolytic, 100μF/50V, 100%         1           200701-006         Capacitor, electrolytic, 100μF/50V, 100%         1           21         200701-042         Capacitor, 0.47μF/50V, 10%         1           22         201628-000         Capacitor, 0.00μF/50V, 5%         1           3-15         Capacitor, 0.00μF/50V, 5%         2           6         200901-046         Capacitor, 0.00μF/50V, 10%         2           8.8         200901-046         Capacitor, 0.00μF/50V, 10%         2           8.8         (Same as C12)         (Same as C2)         3           9         (Same as C8)         3           11, 22         Capacitor, 1.0μF/50V, 5%         2           20         Capacitor, 1.0μF/50V, 5%         4           3.1         Capacitor, 0.00μF/50V, 5%         5           3.2         Capacitor, 0.00μF/50V, 5%         2<				(Same as C8)		C35
53 Ref         Part Number         Description         Use           805473-01         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, Electrolytic, 470μF/10V, 10%         1           200901-026         Capacitor, 0.01μF/50V, 5%         9           202102-032         Capacitor, 220pF/50V, 5%         1           202102-032         Capacitor, 0.00μF/35V, 100%         1           202102-032         Capacitor, 0.00μF/35V, 100%         1           21         200701-006         Capacitor, 0.01μF/50V, 10%         1           22         200701-046         Capacitor, 0.01μF/50V, 5%         1           23         200901-046         Capacitor, 0.01μF/50V, 10%         2           3-15         (Same as CB)         (Same as CB)         2           6         200901-046         Capacitor, 0.01μF/50V, 10%         2           8         (Same as CB)         (Same as CB)         2           9         (Same as CB)         2           11, 22         Capacitor, 0.1μF/50V, 5%         2           12         Capacitor, 0.1μF/50V, 5%         2           13         200901-025         Capacitor, 0.000				(Same as C12)		C34
53 Ref         Part Number         Description         Use Code           Ref         805474-01         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           200501-026         Capacitor, 0.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         1           200701-096         Capacitor, electrolytic, 100μF/35V, 100%         1           200701-096         Capacitor, electrolytic, 33μF/16V, 100%         1           211         200893-042         Capacitor, electrolytic, 33μF/16V, 100%         1           22         201628-000         Capacitor, 0.01μF/50V, 5%         1           3-15         (Same as C8)         (Same as C2)         2           (Same as C2)         (Same as C2)         (Same as C2)         2           (Same as C12)         (Same as C12)         2           11, 22         Capacitor, 10pF/50V, 5%         2           11, 22         Capacitor, 10pF/50V, 5%         2           12, 22         Capacitor, 10pF/50V, 5%         2           13, 29         Capacitor, 10pF/50V, 5%         2           14, 201628-018         Capacitor, 10pF/50V, 5%         2           12, 201628-018         Capacitor				(Same as C31)		C33
53 Ref         Part Number         Description         Use Code           Ref         Part Number         Main PCB Assembly         Code           805474-01         Main PCB Assembly         F           200205-011         Capacitor, D.01µF/50V, 5%         1           201628-006         Capacitor, D.01µF/50V, 5%         1           200901-026         Capacitor, D.01µF/50V, 5%         1           202102-032         Capacitor, D.01µF/50V, 100%         1           202102-042         Capacitor, electrolytic, 100µF/35V, 100%         1           21         200701-06         Capacitor, 0.47µF/50V, 10%         1           2         201628-000         Capacitor, 0.47µF/50V, 10%         2           3-15         Capacitor, 0.01µF/50V, 5%         2           4         200901-046         Capacitor, 0.01µF/50V, 10%         2           5-1         Capacitor, 0.001µF/50V, 10%         2           6         Capacitor, 0.00pF/100V, 10%         3           9         Capacitor, 0.00pF/100V, 10%         2           11, 22         Capacitor, 0.0pF/100V, 10%         2           12         Capacitor, 1.80pF/50V, 5%         2           13         200501-025         Capacitor, 1.0pF/50V, 5%         2				(Same as C16)		C32
53 Ref         Part Number         Description         Use           Ref         805474-01         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470µF/10V, 10%         1           201628-006         Capacitor, 0.01µF/50V, 5%         1           202001-026         Capacitor, 220pF/50V, 5%         9           202102-032         Capacitor, electrolytic, 100µF/35V, 100%         1           1-11         200893-042         Capacitor, electrolytic, 33µF/16V, 100%         1           2         201628-000         Capacitor, 0.01µF/50V, 5%         2           3-15         201628-000         Capacitor, 0.01µF/50V, 10%         2           5-15         200901-046         Capacitor, 0.01µF/50V, 5%         2           6         200901-046         Capacitor, 0.01µF/50V, 10%         2           7         (Same as C2)         (Same as C2)           8         (Same as C3)         (Same as C3)           9         Capacitor, 0.00µF/100V, 10%         9           11, 22         Capacitor, 0.00µF/50V, 5%         9           12         Capacitor, 0.00µF/50V, 5%         1           13         201628-012         Capacitor, 0.00µF/50V, 5%         1           201628-018 <th></th> <th></th> <th></th> <th>(Same as C24)</th> <th>•</th> <th>C31</th>				(Same as C24)	•	C31
53 Ref         Part Number         Description         Use Code           Ref         200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         (Same as C2)         9           200901-026         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 33μF/16V, 100%         1           2         200701-006         Capacitor, 0.47μF/50V, 10%         1           2         201628-000         Capacitor, 0.01μF/50V, 5%         1           3         3-15         201628-000         Capacitor, 0.01μF/50V, 10%         2           6         200901-046         Capacitor, 0.01μF/50V, 10%         3           8         (Same as C2)         (Same as C8)           9         (Same as C8)         9           11, 22         200901-025         Capacitor, 180pF/50V, 5%           12         Capacitor, 190pF/50V, 5%           13         201628-012         Capacitor, 0.1μF/50V, 5%           14         201628-012         Capacitor, 0.01μF/50V, 5%           15         Capacitor, 0.02 per controlled per c	4		-	Capacitor, 1.0µF/50V, 5%	201628-018	C30
53 Ref         Part Number         Description         Use Code           Ref         805474-Φ1         Main PCB Assembly         Part Number           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         Part Number           200505-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           200901-026         Capacitor, 220pF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         1           202102-032         Capacitor, electrolytic, 100μF/35V, 100%         1           -11         200893-042         Capacitor, electrolytic, 33μF/16V, 100%         1           2         201628-000         Capacitor, 0.01μF/50V, 5%         2           3-15         (Same as C8)         (Same as C2)           8         (Same as C2)         (Same as C2)           8         (Same as C12)         9           9         (Same as C8)         9           11, 22         Capacitor, 180pF/100V, 10%         2           14         20628-012         Capacitor, 0.1μF/50V, 5%           20628-012         Capacitor, 180pF/50V, 5%				(Same as C8)		C28, 29
53 Ref         Part Number         Description         Use Code           Ref         805474-Φ1         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         Part Number           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           200901-026         Capacitor, 220pF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         1           202102-032         Capacitor, electrolytic, 100μF/35V, 100%         1           -11         200893-042         Capacitor, electrolytic, 33μF/16V, 100%         1           2         201628-000         Capacitor, 0.47μF/50V, 10%         2           3-15         (Same as C8)         (Same as C2)           8         (Same as C2)         (Same as C2)           8         (Same as C12)         (Same as C8)           9         200893-051         Capacitor, 100pF/100V, 10%         2           11, 22         200901-025         Capacitor, 0.1μF/50V, 5%         2           4         201628-012         Capacitor, 0.1μF/50V, 5%         2				(Same as C12)		C27
53 Ref         Part Number         Description         Use Code           805474-Φ1         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, D.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         9           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           -11         200893-042         Capacitor, electrolytic, 31μF/16V, 100%         1           2         201628-000         Capacitor, 0.01μF/50V, 5%         1           3-15         201628-000         Capacitor, 0.01μF/50V, 10%         2           6         200901-046         Capacitor, 0.01μF/50V, 10%         2           7         Capacitor, 0.01μF/50V, 10%         3           9         Capacitor, 0.00μF/50V, 10%         2           11,22         Capacitor, 0.00μF/50V, 10%         9           12         Capacitor, 100pF/100V, 10%         9           13         200901-025         Capacitor, 100pF/50V, 5%         9	0	r		Capacitor, 0.1µF/50V, 5%	201628-012	C24
53 Ref         Part Number         Description         Use Code           Ref         805474-01         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, Electrolytic, 470μF/10V, 10%         1           200901-026         Capacitor, 220pF/50V, 5%         9           200901-032         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           200701-006         Capacitor, electrolytic, 33μF/16V, 100%         1           2         201628-000         Capacitor, 0.01μF/50V, 5%         1           3-15         (Same as C8)         (Same as C2)           8         (Same as C2)         (Same as C2)           9         (Same as C12)         (Same as C12)           10         Capacitor, 100pF/100V, 10%         2	_			Capacitor, 180pF/50V, 5%	200901-025	C23
Fact         Part Number         Description         Use           805474-01         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, O.01μF/50V, 5%         1           200901-026         Capacitor, O.01μF/50V, 5%         1           202102-032         Capacitor, electrolytic, 100μF/35V, 100%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           211         200893-042         Capacitor, 0.47μF/50V, 10%         1           2         201628-000         Capacitor, 0.01μF/50V, 10%         2           2         201628-000         Capacitor, 0.01μF/50V, 10%         3           3-15         Copolitor, Capacitor, 0.01μF/50V, 10%         3           6         200901-046         Capacitor, 0.01μF/50V, 10%         3           7         (Same as C8)         3           9         (Same as C72)         (Same as C72)           0         (Same as C72)         2           0         Capacitor, 100pF/100V, 10%         2				(Same as C8)		C21, 22
Fact         Part Number         Description         Use           805474-01         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         1           202102-032         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           -11         200893-042         Capacitor, 0.47μF/50V, 10%         1           2         201628-000         Capacitor, 0.001μF/50V, 5%         1           3-15         (Same as C8)         (Same as C8)         3           9         (Same as C12)         (Same as C12)         9	2			Capacitor, 100pF/100V, 10%	200893-051	C20
53 Ref         Part Number         Description         Use           805474-01         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, D.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         9           200901-026         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           2         201628-000         Capacitor, 0.47μF/50V, 10%         1           2         201628-000         Capacitor, 0.001μF/50V, 10%         2           3-15         (Same as C8)         (Same as C8)         3           7         (Same as C8)         9				(Same as C12)		C19
53 Ref         Part Number         Description         Use           805474-01         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         1           200901-026         Capacitor, 0.01μF/50V, 5%         9           202102-032         Capacitor, 220pF/50V, 5%         9           202102-032         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           2-11         200893-042         Capacitor, 0.47μF/50V, 10%         1           2         201628-000         Capacitor, 0.47μF/50V, 10%         2           2         201628-000         Capacitor, 0.01μF/50V, 5%         3           3-15         (Same as C8)         3           7         Capacitor, 0.01μF/50V, 10%         3           6         200901-046         Capacitor, 0.01μF/50V, 10%         3           6         200901-046         Capacitor, 0.01μF/50V, 10%         3				(Same as C8)		C18
S3 Ref         Part Number         Description         Use           805474-01         Main PCB Assembly         P           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         9           202102-032         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           -11         200893-042         Capacitor, electrolytic, 33μF/16V, 100%         1           2         201628-000         Capacitor, 0.01μF/50V, 5%         1           3-15         Capacitor, 0.01μF/50V, 10%         2           6         200901-046         Capacitor, 0.01μF/50V, 10%         9				(Same as C2)		C17
53 Ref         Part Number         Description         Use           805474-01         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         1           202102-032         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           -11         200893-042         Capacitor, 0.47μF/50V, 10%         1           2         201628-000         Capacitor, 0.001μF/50V, 5%         2           3-15         (Same as C8)         3	9			Capacitor, 0.01µF/50V, 10%	200901-046	C16
53 Ref         Part Number         Description         Use           805474-01         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         9           202102-032         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           -11         200893-042         Capacitor, 0.47μF/50V, 10%         1           2         201628-000         Capacitor, 0.001μF/50V, 5%         2				(Same as C8)		C13-15
Feef         Part Number         Description         Use           805474-Φ1         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, Electrolytic, 470μF/10V, 10%         1           200901-026         (Same as C2)         9           202102-032         Capacitor, 220pF/50V, 5%         9           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           200701-006         Capacitor, electrolytic, 33μF/16V, 100%         1           200893-042         Capacitor, 0.47μF/50V, 10%         2	36			Capacitor, 0.001µF/50V, 5%	201628-000	C12
53 Ref         Part Number         Description         Use           805474-01         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         9           200901-026         (Same as C2)         9           202102-032         Capacitor, 220pF/50V, 5%         1           200701-006         Capacitor, electrolytic, 100μF/35V, 100%         1           200701-006         Capacitor, electrolytic, 33μF/16V, 100%         1	24		-	Capacitor, 0.47µF/50V, 10%	200893-042	Ω= 11
53 Ref         Part Number         Description         Use           805474-Φ1         Main PCB Assembly         Code           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         9           200901-026         (Same as C2)         9           202102-032         Capacitor, 220pF/50V, 5%         1           202102-032         Capacitor, electrolytic, 100μF/35V, 100%         1	ш			Capacitor, electrolytic, 33µF/16V, 100%	200701-006	CJ
53 Ref         Part Number         Description         Use           805474-Φ1         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         1           200901-026         Capacitor, 220pF/50V, 5%         9	<b></b> -		esti-si-à	Capacitor, electrolytic, 100µF/35V, 100%	202102-032	Š
53 Ref         Part Number         Description         Use           805474-Φ1         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         9	puret.		<del>oq.,zas</del>	Capacitor, 220pF/50V, 5%	200901-026	C
53 Ref         Part Number         Description         Use           805474-01         Main PCB Assembly         F           200205-011         Capacitor, Electrolytic, 470μF/10V, 10%         1           201628-006         Capacitor, 0.01μF/50V, 5%         9				(Same as C2)		2
53 RefPart NumberDescriptionUse Code805474-01Main PCB AssemblyF200205-011Capacitor, Electrolytic, 470μF/10V, 10%1	9			Capacitor, 0.01µF/50V, 5%	201628-006	3
Ref Part Number Description Code  805474-01 Main PCB Assembly	<b>—</b>			Capacitor, Electrolytic, 470µF/10V, 10%	200205-011	CI
Ref Part Number Description Code  805474-01 Main PCB Assembly				•		
Part Number Description Code	REF	-		Main PCB Assembly	805474-01	Al
	Š	Code		Description	Part Number	75-5 Fef 3

ਜੂ 5 ਦ੍ਰ	Part Number	Description	Use Code	Ş Q
CR3	2000-380002	Diode, Zener, 1N4620, 80mA/3.3V (Same as CR2)		<b>—</b>
CR7-9	200605-010	Diode, Zener, 1N4738A, 550mA/8.2V	:	6
CR10-17		(Same as CR2)		
CR18, 19	200086-019	Diode, Zener, 1N4104, 24.8mA/10V		2
CR31		(Same as CR7)		
CR21-28		(Same as CR2)		, top 1, to 1
CR31-33		(Same as CR7)		
CR34		(Same as CR2)		
CR35	201749-012	Diode, transzorb, SA176		<b></b>
<b>J</b> 4	201191-221	Connector, locking, 20 contact	• • •	
55	201191–222	Connector, locking, 26 contact		
<b>J</b> 6	202128-029	Connector, 6 contact		<b>)—</b>
8	201504-003	Connector, spade 0.187		<b></b>
J9	201504-004	Connector, spade, 0.25	- LIMBURA	<b></b>
J22	202128-026	Connector, locking, 4 contact		<b>,_</b>
J29	202128-024	Connector, locking, 2 contact	\.	<u>,                                    </u>
J35	202128-048	Connector, right angle, 2 contact		-
K1	202045000	Relay, 1500V, 60W		-
jamak Jamak	201409-036	Inductor, 100µH, 5%		-
L2,3	200808-000	Inductor, ferrite bead		2
P7	201678–022	Connector, 96 Contact		jene už
Q1	200585-001	Transistor, VFET, VN0106N3, TO-92		7
Q2	200942-001	Transistor, NPN, PN2222, TO-92		-
Q3	202155-000	Transistor, FET, IRFD9010		<b>)—</b> \
20	200994-001	Transistor, VFET, IRF9531, TO-220		j4
07	200885-000	Chame as Q1)		
08,9		(Same as Q1)		
Q10	200582-004	Transistor, FET, VCR7N, TO-72		<b>—</b>
011	202179-001	Transistor, JFET, J175, TO-92, p-channel		,
Q12,13		(Same as Q1)	1.	- H
Q14	200791-001	Transistor, NPN, 2N3904, TO-92		
Q15	200883000	Transistor, PNP, 2N3905, TO-92	:	<b>j—1</b>

20054-257 Resistor, 4.754Ω, 1/8W, 1%, RN55C 20054-384 Resistor, 20/kQ, 1/8W, 1%, RN55C 20054-413 Resistor, 20/kQ, 1/8W, 1%, RN55C 20054-263 Resistor, 20/kQ, 1/8W, 1%, RN55C 20054-288 Resistor, 1/kQ, 1/8W, 1%, RN55C 20054-222 Resistor, 1/kQ, 1/8W, 1%, RN55C 200054-229 Resistor, 1/kQ, 1/8W, 1%, RN55C 200054-279 Resistor, 20/kQ, 1/8W, 1%, RN55C 200054-173 Resistor, 20/kQ, 1/8W, 1%, RN55C 200054-173 Resistor, 20/kQ, 1/8W, 1%, RN55C 200054-193 Resistor, 20/kQ, 1/8W, 1%, RN55C CSame as R10 200054-194 Resistor, 20/kQ, 1/8W, 1%, RN55C CSame as R18) CSame as R18 CSame as R19 CSAMe as	Fef	Part Number	Description	Use Code	Qty
20054–384 Resistor, 100kQ, 1/8W, 1%, RN55C 20054–413 Resistor, 200kQ, 1/8W, 1%, RN55C 20054–263 Resistor, 200kQ, 1/8W, 1%, RN55C 20054–263 Resistor, 100kQ, 1/8W, 1%, RN55C 20054–288 Resistor, 1,0kQ, 1/8W, 1%, RN55C 20054–279 Resistor, 1,0kQ, 1/8W, 1%, RN55C 20054–279 Resistor, 200Q, 1/8W, 1%, RN55C 20054–279 Resistor, 200Q, 1/8W, 1%, RN55C 20054–279 Resistor, 200Q, 1/8W, 1%, RN55C 20054–117 Resistor, 200Q, 1/8W, 1%, RN55C 20054–126 Resistor, 200Q, 1/8W, 1%, RN55C 20054–127 Resistor, 200Q, 1/8W, 1%, RN55C 20054–193 Resistor, 200Q, 1/8W, 1%, RN55C 20054–193 Resistor, 200Q, 1/8W, 1%, RN55C Same as R10 20054–451 Resistor, 39, 2kQ, 1/8W, 1%, RN55C Same as R10 20054–193 Resistor, 200Q, 1/8W, 1%, RN55C Same as R10 20054–193 Resistor, 200Q, 1/8W, 1%, RN55C Same as R10 20054–193 Resistor, 200Q, 1/8W, 1%, RN55C Same as R10 20054–193 Resistor, 200Q, 1/8W, 1%, RN55C Same as R10 Same as R1) Same as R2) Resistor, 100Q, 1/8W, 1%, RN55C Same as R1) Same as R2) Resistor, 201Q, 1/8W, 1%, RN55C Same as R1) Same as R2) Resistor, 201Q, 1/8W, 1%, RN55C Same as R1) Same as R2) Resistor, 201Q, 1/8W, 1%, RN55C Same as R24) Resistor, 201Q, 1/8W, 1%, RN55C	R1,2	200054–257	1/8W, 1%, RN55C		14
6 (20054-413 (Resistor, 200kQ, 1/8W, 1%, RN55C) (Same as R1) (Same as R1) (20054-363 (Resistor, 60.4kQ, 1/8W, 1%, RN55C) (200527-009 (Potentiometer, 10kQ, 1/2W, 10%, multitum) (20054-480 (Resistor, 10kQ, 1/8W, 1%, RN55C) (20054-288 (Resistor, 10kQ, 1/8W, 1%, RN55C) (20054-222 (Resistor, 10kQ, 1/8W, 1%, RN55C) (20054-259 (Resistor, 20kQ, 1/8W, 1%, RN55C) (Same as R10) (20054-279 (Resistor, 20kQ, 1/8W, 1%, RN55C) (Same as R10) (Same a	R3	200054-384	Resistor, 100kΩ, 1/8W, 1%, RN55C		9,
6 (Same as R1)  20054-263 Resistor, 60.442, 1/8W, 1%, RNSSC 20057-209 Potentiometer, 1042, 1/2W, 10%, mulitum 20054-280 Resistor, 60.442, 1/8W, 1%, RNSSC 20054-222 Resistor, 1042, 1/8W, 1%, RNSSC 200054-229 Resistor, 1042, 1/8W, 1%, RNSSC 200054-279 Resistor, 2942, 1/8W, 1%, RNSSC 200054-279 Resistor, 2042, 1/8W, 1%, RNSSC 200054-279 Resistor, 2042, 1/8W, 1%, RNSSC 200054-216 Resistor, 2042, 1/8W, 1%, RNSSC 200054-126 Resistor, 2042, 1/8W, 1%, RNSSC 200054-126 Resistor, 2042, 1/8W, 1%, RNSSC 200054-126 Resistor, 2042, 1/8W, 1%, RNSSC 200054-193 Resistor, 30.24, 1/8W, 1%, RNSSC 200054-193 Resistor, 102, 1/8W, 1%, RNSSC 200054-194 Resistor, 102, 1/8W, 1%, RNSSC 200054-195 Resistor, 102, 1/8W, 1%, RNSSC 200054-196 Resistor, 102, 1/8W, 1%, RNSSC 200054-197 Resistor, 102, 1/8W, 1%, RNSSC 200054-198 Resistor, 104, 1/8W, 1%, RNSSC 200054-199 Resistor, 104, 1/8W, 1%, RNSSC 200054-297 Resistor, 104, 1/8W, 1%, RNSSC 200054-297 Resistor, 104, 1/8W, 1%, RNSSC 200054-296 Resistor, 104, 1/8W, 1%, RNSSC 200054-296 Resistor, 104, 1/8W, 1%, RNSSC 200054-296 Resistor, 20143, 1/8W, 1%, RNSSC	R4	200054-413	Resistor, 200kΩ, 1/8W, 1%, RN55C		<del></del>
200054-363 Resistor, 60.4kΩ, 1/8W, 1%, RN55C 200637-209 Potentionneter, 10kΩ, 1/2W, 10%, multitum 200054-428 Resistor, 10kΩ, 1/8W, 1%, RN55C 200054-288 Resistor, 10kΩ, 1/8W, 1%, RN55C 200054-259 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-279 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-279 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-317 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-317 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-126 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-173 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-193 Resistor, 10kΩ, 1/8W, 1%, RN55C 200054-193 Resistor, 10kΩ, 1/8W, 1%, RN55C (Same as R10) 200054-194 Resistor, 10kΩ, 1/8W, 1%, RN55C (Same as R1) Same as R1)	R5-6		(Same as R1)		
200527-009  Potentiometer, 10kΩ, 112W, 10%, multitum  200054-480  Resistor, 1MΩ, 18W, 1%, RN55C  200054-222  Resistor, 2kΩ, 1/8W, 1%, RN55C  200054-229  Resistor, 2kΩ, 1/8W, 1%, RN55C  200054-239  Resistor, 20kΩ, 1/8W, 1%, RN55C  200054-279  Resistor, 20kΩ, 1/8W, 1%, RN55C  200054-317  Resistor, 20kΩ, 1/8W, 1%, RN55C  200054-317  Resistor, 20kΩ, 1/8W, 1%, RN55C  200054-126  Resistor, 20kΩ, 1/8W, 1%, RN55C  200054-127  Resistor, 20kΩ, 1/8W, 1%, RN55C  200054-128  Resistor, 20kΩ, 1/8W, 1%, RN55C  200054-193  Resistor, 20kΩ, 1/8W, 1%, RN55C  (Same as R10)  200054-149  Resistor, 20kΩ, 1/8W, 1%, RN55C  (Same as R18)  (Same as R19)  (Same as R1)	R7	200054-363	Resistor, 60.4kQ, 1/8W, 1%, RN55C		w
200054-480     Resistor, 1MΩ, 1/8W, 1%, RN55C 200054-222     Resistor, 2kΩ, 1/8W, 1%, RN55C 200054-259     Resistor, 2kΩ, 1/8W, 1%, RN55C 200054-259     Resistor, 29kΩ, 1/8W, 1%, RN55C 200054-279     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-279     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-279     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-176     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-173     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-193     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-194     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-197     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-286     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-297     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-297     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-297     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-296     Resistor, 20kΩ, 1/8W, 1%, RN55C 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054-296 200054	R8	200527-009	Potentiometer, 10kΩ 1/2W, 10%, multiturn		_
200054–222 Resistor, 2kΩ, 1/8W, 1%, RN55C 200054–288 Resistor, 10kΩ, 1/8W, 1%, RN55C 200054–259 Resistor, 4.99kΩ, 1/8W, 1%, RN55C 200054–279 Resistor, 8.06kΩ, 1/8W, 1%, RN55C 200054–279 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054–279 Resistor, 10Ω, 2W, 5% 200054–126 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054–126 Resistor, 10Ω, 2W, 5% 200054–193 Resistor, 10Ω, 2W, 5% 200054–193 Resistor, 10Ω, 1/8W, 1%, RN55C (Same as R10) 200054–345 Resistor, 0.0Ω, 1/8W, 1%, RN55C (Same as R18) (Same as R19) (Same as R1)	<b>R9</b>	200054-480	Resistor; 1MΩ, 1/8W, 1%, RN55C	٠	w
200054–288 Resistor, 10kΩ 1/8W, 1%, RN55C 200054–259 Resistor, 4.99kΩ, 1/8W, 1%, RN55C 200054–417 Resistor, 221kΩ 1/8W, 1%, RN55C 200054–417 Resistor, 221kΩ 1/8W, 1%, RN55C 200054–317 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054–317 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054–126 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054–193 Resistor, 20kΩ, 1/8W, 1%, RN55C 200054–193 Resistor, 1kΩ 1/8W, 1%, RN55C (Same as R10) 200054–345 Resistor, 10kΩ, 1/8W, 1%, RN55C (Same as R18) (Same as R1) (Same as R18) 200054–164 Resistor, 499Ω, 1/8W, 1%, RN55C (Same as R18) 200054–1097 Resistor, 10kΩ, 1/8W, 1%, RN55C (Same as R18) 200054–266 Resistor, 10kΩ, 1/8W, 1%, RN55C (Same as R1) (Same as R2) (Same as R2) (Same as R3) (Same as R3) (Same as R1) (Same as R4) (Same as R2)	RIO	200054-222	Resistor, 2kΩ 1/8W, 1%, RN55C		4
200054-259   Resistor, 4.99kΩ, 1/8v, 1%, RN55C     200054-417   Resistor, 8.06kΩ, 1/8W, 1%, RN55C     200054-417   Resistor, 8.06kΩ, 1/8W, 1%, RN55C     200054-317   Resistor, 20kΩ, 1/8W, 1%, RN55C     200054-317   Resistor, 20kΩ, 1/8W, 1%, RN55C     200054-317   Resistor, 20kΩ, 1/8W, 1%, RN55C     200054-126   Resistor, 10Ω, 2W, 5%     200054-073   Resistor, 20kΩ, 1/8W, 1%, RN55C     (Same as R10)     200054-193   Resistor, 1kΩ, 1/8W, 1%, RN55C     (Same as R18)     200054-345   Resistor, 0.2\(\hat{L}\) 1W, 1%, RN55C     (Same as R18)     200054-164   Resistor, 39.2\(\hat{L}\) 1W, 1%, RN55C     (Same as R18)     200054-097   Resistor, 100\(\hat{L}\) 1/8W, 1%, RN55C     (Same as R1)     (Same as R1)     (Same as R1)     (Same as R2)     (Same as R1)     (Same as R3)     (Same as R4)     (Same as R2)     (Same as R3)     (Same as R3)     (Same as R4)     (Same as R2)     (Same as R3)     (Same as R3)     (Same as R4)     (Same as R4)     (Same as R2)     (Same as R3)     (Same as R4)     (Same as R2)     (Same as R3)     (Same as R2)     (Same as R2)     (Same as R2)     (Same as R3)     (Same as R2)     (Same as R2)     (Same as R3)     (Same as R3)     (Same as R4)     (Same as R2)     (Same as R2)     (Same as R3)     (Same as R3)     (Same as R3)     (Same as R4)     (Same as R4)     (Same as R3)     (Same as R3)     (Same as R4)     (Same	RI	200054-288	Resistor, 10kQ, 1/8W, 1%, RN55C		16
(Same as R10)  200054-417  Resistor, 221KQ, 1/8W, 1%, RN55C  200054-279  Resistor, 8.06kQ, 1/8W, 1%, RN55C  (Same as R7)  200054-126  Resistor, 200Q, 1/8W, 1%, RN55C  200054-126  Resistor, 200Q, 1/8W, 1%, RN55C  200054-073  Resistor, 10Q, 2W, 5%  200054-193  Resistor, 10Q, 2W, 5%, RN55C  (Same as R10)  200054-345  Resistor, 10Q, 1/8W, 1%, RN55C  (Same as R18)  (Same as R1)  (Same as R2)  (Same as R1)  (Same as R1)  (Same as R1)  (Same as R2)  (Same as R1)  (Same as R1)  (Same as R2)  (Same as R1)  (Same as R2)  (Same as R1)  (Same as R1)  (Same as R2)  (Same as R1)  (Same as R2)  (Same as R1)  (Same as R1)  (Same as R2)  (Same as R1)  (Same as R2)  (Same as R2)  (Same as R2)  (Same as R1)  (Same as R1)  (Same as R2)  (Same as R1)  (Same as R2)	R12	200054-259	Resistor, 4.99kQ, 1/8W, 1%, RN55C	.,	4
200054-417 Resistor, 221KQ, 1/8W, 1%, RN55C 200054-279 Resistor, 8.06kQ, 1/8W, 1%, RN55C (Same as R1) 200054-317 Resistor, 20kQ, 1/8W, 1%, RN55C 200054-126 Resistor, 10Q, 2W, 5% 200054-173 Resistor, 56.2Q, 1/8W, 1%, RN55C (Same as R10) 200054-193 Resistor, 1kQ, 1/8W, 1%, RN55C (Same as R18) (Same as R1) (Same as R2) (Same as R2) (Same as R1) (Same as R2) (Same as R2) (Same as R1) (Same as R2)	R13		(Same as R10)		
200054-279 Resistor, 8.06kΩ, 1/8W, 1%, RN55C (Same as R11) (Same as R11) 200054-317 Resistor, 200Ω, 1/8W, 1%, RN55C 200270-048 Resistor, 200Ω, 1/8W, 1%, RN55C 200054-126 Resistor, 200Ω, 1/8W, 1%, RN55C (Same as R10) 200054-193 Resistor, 10Ω, 1/8W, 1%, RN55C (Same as R18) (Same as R19)	<b>714</b>	200054 417	Resistor, 221kΩ 1/8W, 1%, RN55C		<b>—</b>
(Same as R11) (Same as R7)  200054–317 (Resistor, 20k2, 1/8W, 1%, RN55C 200270–048 (Resistor, 20k2, 1/8W, 1%, RN55C 200054–126 (Resistor, 2002, 1/8W, 1%, RN55C 200054–193 (Resistor, 12A, 1/8W, 1%, RN55C (Same as R10) (Same as R18) (Same as R1) (Same as R18) (Same as R18) (Same as R18) (Same as R18) (Same as R1) (Same as R18) (Same as R18) (Same as R19) (Same as R1) (Same as R24) (Same as R24) (Same as R14)	R15	200054-279	Resistor, 8.06kQ, 1/8W, 1%, RN55C		<del></del>
(Same as R7)  200054-317  Resistor, 20kQ, 1/8W, 1%, RN55C  200070-048  Resistor, 10Q, 2W, 5%  200054-126  Resistor, 200Q, 1/8W, 1%, RN55C  200054-073  Resistor, 5(2Q, 1/8W, 1%, RN55C  (Same as R10)  200054-193  Resistor, 1kQ, 1/8W, 1%, RN55C  (Same as R18)  (Same as R18)  200054-345  Resistor, 20Q, 1/8W, 1%, RN55C  (Same as R18)  200054-097  Resistor, 499Q, 1/8W, 1%, RN55C  (Same as R12)  (Same as R1)	R16		(Same as R11)		
200054-317 Resistor, 20kQ, 1/8W, 1%, RN55C 200270-048 Resistor, 10kQ, 2W, 5% 200054-126 Resistor, 50.2Q, 1/8W, 1%, RN55C 200054-073 Resistor, 56.2Q, 1/8W, 1%, RN55C (Same as R10) 200054-193 Resistor, 1kQ, 1/8W, 1%, RN55C (Same as R18) 200054-345 Resistor, 9.2kQ, 1/8W, 1%, RN55C (Same as R18) (Same as R19) (Same as R19) (Same as R19) (Same as R19) (Same as R1)	R17		(Same as R7)		
200270-048 Resistor, 10Ω, 2W, 5% 200054-126 Resistor, 20ΩΩ, 1/8W, 1%, RN55C 200054-073 Resistor, 56.2Ω, 1/8W, 1%, RN55C (Same as R10) 200054-193 Resistor, 1kΩ, 1/8W, 1%, RN55C (Same as R18) (Same as R18) 200216-002 Resistor, 0.2Ω, 1W, 1% 200054-345 Resistor, 20Ω, 1/8W, 1%, RN55C (Same as R18) (Same as R19) 200471-031 Resistor, 20Ω, 1/2W, 5% 200054-097 Resistor, 10ΩΩ, 1/8W, 1%, RN55C (Same as R1)	R18	200054-317	Resistor, 20kΩ 1/8W, 1%, RN55C		17
200054–126 Resistor, 200Ω 1/8W, 1%, RN55C 200054–193 Resistor, 56.2Ω, 1/8W, 1%, RN55C (Same as R10) 200054–193 Resistor, 1kΩ, 1/8W, 1%, RN55C (Same as R18) (Same as R8) 200216–002 Resistor, 0.2Ω, 1W, 1% 200054–345 Resistor, 39.2kΩ, 1/8W, 1%, RN55C (Same as R1) (Same as R1) (Same as R3) (Same as R3) (Same as R18) 200471–031 Resistor, 20Ω, 1/2W, 5% 200054–297 Resistor, 10ΩΩ, 1/8W, 1%, RN55C (Same as R1) (Same as R1) (Same as R4) (Same as R24) (Same as R24) Resistor, 2.21kΩ, 1/8W, 1%, RN55C (Same as R11) (Same as R11) Resistor, 4.99kΩ, 1/8W, 1%, RN55C	R19	200270-048	Resistor, 10Ω 2W, 5%		-
200054-073   Resistor, 56.2Ω, 1/8W, 1%, RN55C     (Same as R10)   Resistor, 1kΩ, 1/8W, 1%, RN55C     (Same as R18)   (Same as R18)     (Same as R9)     200216-002   Resistor, 0.2Ω, 1W, 1%     200054-345   Resistor, 39.2kΩ, 1/8W, 1%, RN55C     (Same as R1)   (Same as R1)     200054-164   Resistor, 499Ω, 1/8W, 1%, RN55C     (Same as R18)   (Same as R18)     200471-031   Resistor, 20Ω, 1/2W, 5%     200054-097   Resistor, 100Ω, 1/8W, 1%, RN55C     (Same as R1)   (Same as R1)     (Same as R1)   (Same as R24)     (Same as R24)   (Same as R11)     (Same as R11)   (Same as R11)     (Same as R11)   (Same as R11)     (Resistor, 2.21kΩ, 1/8W, 1%, RN55C     (Same as R11)   (Resistor, 4.99kΩ, 1/8W, 1%, RN55C     (Same as R11)   (Resistor, 4.99kΩ, 1/8W, 1%, RN55C     (Resistor, 4.99kΩ, 1/8W, 1/4W, 1/	R21	200054-126	Resistor, 200Q, 1/8W, 1%, RN55C		2
(Same as R10)  200054–193  Resistor, 1kQ, 1/8W, 1%, RN55C  (Same as R18)  (Same as R9)  200216–002  Resistor, 0.2Q, 1W, 1%  200054–345  Resistor, 39.2kQ, 1/8W, 1%, RN55C  (Same as R18)  (Same as R18)  200054–164  Resistor, 499Q, 1/8W, 1%, RN55C  (Same as R18)  (Same as R18)  (Same as R18)  (Same as R11)  (Same as R1)  (Same as R1)  (Same as R24)  (Same as R11)	R22	200054-073	Resistor, 56.2Q 1/8W, 1%, RN55C		<b>,</b>
20054–193 Resistor, 1kΩ, 1/8W, 1%, RN55C  (Same as R18)  200216–002 Resistor, 0.2Ω, 1W, 1%  200054–345 Resistor, 39.2kΩ, 1/8W, 1%, RN55C  (Same as R18)  200054–164 Resistor, 499Ω, 1/8W, 1%, RN55C  (Same as R13)  (Same as R13)  (Same as R12)  (Same as R12)  (Same as R11)  (Same as R14)  (Same as R24)  200054–226 Resistor, 2.21kΩ, 1/8W, 1%, RN55C  (Same as R11)	R23		(Same as R10)		
(Same as R18) (Same as R9)  200216-002  Resistor, 0.2\(\Omega}, 1\W, 1\%  200054-345  Resistor, 39.2\(\Omega}, 1/8\W, 1\%, RN55C  (Same as R18)  (Same as R1)  (Same as R18)  200054-164  Resistor, 499\(\Omega}, 1/8\W, 1\%, RN55C  (Same as R18)  (Same as R18)  (Same as R12)  (Same as R12)  (Same as R1)  (Same as R1)  (Same as R24)  (Same as R24)  (Same as R11)	R24	200054-193	Resistor, 1kΩ, 1/8W, 1%, RN55C		10
(Same as R9)  200216-002  Resistor, 0.2Ω, 1W, 1%  Resistor, 39.2kΩ, 1/8W, 1%, RN55C  (Same as R18)  (Same as R18)  200054-164  Resistor, 499Ω, 1/8W, 1%, RN55C  (Same as R18)  200471-031  Resistor, 20Ω, 1/2W, 5%  Resistor, 100Ω, 1/8W, 1%, RN55C  (Same as R12)  (Same as R1)  (Same as R4)  (Same as R24)  (Same as R21)  (Same as R24)  (Same as R11)  Resistor, 2.21kΩ, 1/8W, 1%, RN55C  (Same as R11)  Resistor, 499kΩ, 1/8W, 1%, RN55C	R25		(Same as R18)		
200216-002 Resistor, 0.2Ω, 1W, 1% 200054-345 Resistor, 39.2kΩ, 1/8W, 1%, RN55C  (Same as R18) 200054-164 Resistor, 499Ω, 1/8W, 1%, RN55C (Same as R3) (Same as R18) 200471-031 Resistor, 20Ω, 1/2W, 5% 200054-097 Resistor, 100Ω, 1/8W, 1%, RN55C (Same as R12) (Same as R12) (Same as R1) (Same as R4) (Same as R24) (Same as R24) (Same as R11) (Resistor, 2.21kΩ, 1/8W, 1%, RN55C (Same as R11) (Resistor, 499kΩ, 1/8W, 1%, RN55C)	R26		(Same as R9)		
200054–345   Resistor, 39.2kΩ, 1/8W, 1%, RN55C     30   (Same as R18)     200054–164   Resistor, 499Ω, 1/8W, 1%, RN55C     200054–164   Resistor, 20Ω, 1/8W, 1%, RN55C     200054–097   Resistor, 20Ω, 1/2W, 5%     200054–097   (Same as R12)     (Same as R1)     (Same as R1)     (Same as R4)     (Same as R24)     (Same as R24)     (Same as R24)     (Same as R1)     (Same as R24)     (Same as R1)     (Same as R1)     (Same as R24)     (Same as R24)     (Same as R1)     (Same as R24)     (Same as R1)     (Same as R24)     (Same as R24)     (Same as R24)     (Same as R1)     (Same as R24)	R27	200216002	Resistor, 0.2Ω, 1W, 1%		<b></b>
, 30 (Same as R18) (Same as R1)  200054–164 Resistor, 499Ω, 1/8W, 1%, RN55C (Same as R3) (Same as R18)  200471–031 Resistor, 20Ω, 1/2W, 5%  200054–097 Resistor, 100Ω, 1/8W, 1%, RN55C (Same as R12) (Same as R1) (Same as R4) (Same as R4) (Same as R24) (Same as R11) (Resistor, 2.21kΩ, 1/8W, 1%, RN55C (Same as R11)	R28	200054-345	Resistor, 39.2kQ, 1/8W, 1%, RN55C		<b></b>
(Same as R1)  200054–164  Resistor, 499\(\Omega\), 1/8W, 1\%, RN55C  (Same as R3)  (Same as R18)  200471–031  Resistor, 20\(\Omega\), 1/2W, 5\%  Resistor, 100\(\Omega\), 1/8W, 1\%, RN55C  (Same as R12)  (Same as R1)  (Same as R4)  (Same as R24)  200054–226  Resistor, 221k\(\Omega\), 1/8W, 1\%, RN55C  (Same as R11)  Resistor, 499k\(\Omega\), 1/8W, 1\%, RN55C	R29, 30		(Same as R18)		
200054–164  Resistor, 499Ω, 1/8W, 1%, RN55C  (Same as R18)  200471–031  Resistor, 20Ω, 1/2W, 5%  200054–097  (Same as R12)  (Same as R1)  (Same as R4)  (Same as R24)  200054–226  Resistor, 2.21kΩ, 1/8W, 1%, RN55C  (Same as R11)  (Same as R11)  (Same as R11)  (Same as R11)	R31, 32		(Same as R1)		
(Same as R3) (Same as R18) 200471–031  Resistor, 20Q, 1/2W, 5% 200054–097  Resistor, 100Q, 1/8W, 1%, RN55C  (Same as R12) (Same as R1) (Same as R4) (Same as R24) (Same as R24) (Same as R11) (Same as R11) (Same as R11) (Resistor, 2.21kQ, 1/8W, 1%, RN55C (Same as R11) (Resistor, 499kQ, 1/8W, 1%, RN55C	R33	200054-164	Resistor, 499Q 1/8W, 1%, RN55C		2
(Same as R18) 200471–031 Resistor, 20Ω, 1/2W, 5% 200054–097 Resistor, 100Ω, 1/8W, 1%, RN55C (Same as R12) (Same as R1) (Same as R4) (Same as R24) (Same as R24) (Same as R21) (Same as R11) (Same as R11) (Same as R11) (Same as R11) (Resistor, 499kΩ, 1/8W, 1%, RN55C	R34		(Same as R3)		
200471–031  Resistor, 20Ω, 1/2W, 5% 200054–097  Resistor, 100Ω, 1/8W, 1%, RN55C  (Same as R12) (Same as R1) (Same as R4) (Same as R24)  200054–226  Resistor, 2.21kΩ, 1/8W, 1%, RN55C  (Same as R11)  Resistor, 499kΩ, 1/8, RN55C	R35		(Same as R18)		
200054-097  Resistor, 100Ω, 1/8W, 1%, RN55C  (Same as R12) (Same as R1) (Same as R4) (Same as R24) (Same as R24) (Same as R24) (Same as R11)	R36	200471-031	Resistor, 20Q, 1/2W, 5%		<u>ь</u>
(Same as R12) (Same as R1) (Same as R1) (Same as R4) (Same as R24) (Same as R24) (Same as R211) (Same as R11) (Same as R11) (Same as R11)	R37	200054-097	Resistor, 100Q, 1/8W, 1%, RN55C	· · · · · · ·	2
(Same as R1) (Same as R1) (Same as R4) (Same as R24) (Same as R24) (Same as R211) (Same as R11) (Same as R11) (Same as R11) (Same as R11)	R38		(Same as R12)		
(Same as R1) (Same as R4) (Same as R24) (Same as R24) (Same as R11)	R39		(Same as R1)		
(Same as R4) (Same as R24) 200054–226	R40		(Same as R1)		
(Same as R24)  200054–226	R41		(Same as R4)		
200054–226 Resistor, 2.21kΩ, 1/8W, 1%, RN55C (Same as R11) 200054–451 Resistor, 499kΩ, 1/8W, 1%, RN55C	R42		(Same as R24)		
(Same as R11) <b>Resistor</b> , 499kQ, 1/8W, 1%, RN55C	R43	200054-226	Resistor, 2.21kQ, 1/8W, 1%, RN55C		31
200054-451 Resistor, 499kQ, 1/8W, 1%, RN55C	R44		(Same as R11)		
	R45	200054 451	Resistor, 499kΩ, 1/8W, 1%, RN55C		3

70 55 Per 33	Part Number	Description	Use	Q V
R46	200054-318	Resistor, 20.5kΩ, 1/8W, 1%, RN55C		<u> </u>
R47	200054-335	Resistor, 30.9kΩ 1/8W, 1%, RN55C		_
R48	200054–321	Resistor, 22.1kQ, 1/8W, 1%, RN55C		5
R49		(Same as R4)		
R50		(Same as R48)		
R51, 52		(Same as R43)		
R53		(Same as R48)		
R54	200054-344	Resistor, 38.3kQ, 1/8W, 1%		_
第55	200054-105	Resistor, 121Q 1/8W, 1%, RN55C		2
R56-58		(Same as R43)	>Ne-Holy	,
R59		(Same as R48)		
R60		(Same as R43)		
R61		(Same as R48)		
R62-66		(Same as R43)		
R67		(Same as R4)		
R68		(Same as R1)		
R69		(Same as R41)		
R70		(Same as R3)		
R71	200054-169	Resistor, 562Q, 1/8W, 1%, RN55C		<del></del> -
R72-76		(Same as R43)		
R77		(Same as R55)		
R78-86		(Same as R43)	/	
R87, 88	200054–139	Resistor, 274Q, 1/8W, 1%, RN55C	:	2
R89, 90		(Same as R24)	٠.	
Ry1		(Same as K1)		
R92, 93	21/17 00	(Same as R43)		
4	4600 471	(Same as KII)		
	<del>personal</del>	(Same as R42)	25 0.2	
R97		(Same as R4)		
R98	201155-003	Resistor, 1000MΩ 1W, 10%		<del></del>
R99	200054-114	Resistor, 150Q, 1/8W, 1%, RN55C		2
R100	A VIDE A	(Same as R18)		
R101	200054-401	Resistor, 150kΩ, 1/8W, 1%, RN55C		<u>ш</u>
R102	250 (cm - 20 cm)	(Same as R18)		
R103, 104		(Same as R11)		
R105		(Same as R18)		
R106	200054–165	Resistor, 511Q, 1/8W, 1%, RN55C		2
R107, 108		(Same as R4)		
KION		(Same as R24)		

5–3 Ref	Part Number	Description		Use Code	Q
R161	200054-267	Resistor, 6.04kΩ, 1/8W, 1%, RN55C			1
R162	200054-373		-		<u></u>
R163	200054-402	Resistor, 154kQ, 1/8W, 1%, RN55C			_
R164	200054-427	Resistor, 280kΩ, 1/8W, 1%, RN55C	ı		2
R165	-	(Same as R145)			
R166		(Same as R164)			••••••
R167		(Same as R11)			
R168	200054-387	Resistor, 107kΩ, 1/8W, 1%, RN55C			<b>–</b>
R169		(Same as R12)			
R170		(Same as R21)	ı		_
R171, 172	200054-525	Resistor, 3.01MΩ1/8W, 1%, RN55C			2
R174, 175		(Same as R18)			<del></del>
R176	200527012	Potentiometer, 50kQ, 1/2W			
R177		(Same as R43)			
R179		(Same as R1)			
R181,182		(Same as R1)			
R184		(Same as R11)			
R185	200054-130	Resistor, 221Q, 1/8W, 1%, RN55C			
R186		(Same as R11)			
R187	200054-259	Resistor, 4.99kQ, 1/8W, 1%, RN55C			<u> </u>
R188		(Same as R33)			
R189		(Same as R11)			
R190		(Same as R18)			
R191	200054-246	Resistor, 3.57kQ, 1/8W, 1%, RN55C			<b>;</b>
R192		(Same as R10)			
R193		(Same as R18)	,		
R194	200054-250	Resistor, 4.02kQ, 1/8W, 1%, RN55C			<b>}4</b>
R195		(Same as R99)			•••
R196		(Same as R43)			•••
R203		(Same as R1)	.*		
R204,205	200054-231	Resistor, 2.94kQ, 1/8W, 1%, RN55C			4
R206-208		(Same as R11)		-	
R209-212		(Same as R24)			- *
R213		(Same as R4)			
R214	200054-001	Resistor, 100, 1/8W, RN55C			-
R215,216		(Same as R204)			
RZI	201425-018	Resistor Network, 4 x 20kQ, 1W, 1%			<u> </u>
RN2	201425-004	Resistor Network, 4 x 20kΩ 1W, 1%			<b>)</b>


5–3 Ref	Part Number	Description Use Code	Oty
	803089-02	Transformer, isolation, 10kV, 290T-CT	
TP1-23	801959-00	Terminal, Test Point	23
C	250015-003	IC, digital, HD63C03Y, microprocessor SMT	2
U2	200669-002		2
U3	804280-01	TMC	
40	201703-000		<del>}</del>
U5	201126-000	IC, linear, 3525AN, pulse-width modulator	) — · ·
U6	201597-002		
U7	201597-000	IC, linear, CNY65, opto-isolator	2
U9	805576-01	IC, digital, programmed gate array logic	<del></del>
U10	202087000		· -
U11, 12	201511-001		. u
	250013-007	IC, digital, 628128 RAM, 128k x 8	
U15			
U16	201588-000	IC, digital, 74HC74, D-type flip-flop	}t
U17	201608-001	IC, digital, 74HCT02, NOR gate	)t
U18	201469-001	IC, digital, 74HCT00, NAND gate	2
U19,20	201638-001	-flop	2
U21	250014-006		j
U22	803710-200	IC, digital, 27C256, EPROM 🛦	<del></del>
U23	201705-000	IC, digital, 74HC14, inverter	2
U24		(Same as U18)	
U25	805575-01	IC, digital, programmed gate array logic	<b></b>
U26	201774-000	IC, digital, 74HC4538N, multivibrator	<b>—</b>
U27	803710-908	IC, EPROM, programmed system 🛕	<b>—</b>
U28A, U28B	202227-000	IC, digital, DS1242Y, smartwatch/RAM, 256k	<b>⊢</b> ∸
U29	803720-01	IC, digital, SCX6218TUA/V4 display gate array SMT	ļt
U30		(Same as U1)	
U31	201611-000	IC, digital, 74HC221, multivibrator	<b></b>
U32	201616-002	IC, digital, 74HC161, counter	<del>}</del>
U33	200165-000	IC, digital, 74C174, D-type flip-flop	<del>}</del>
U34		(Same as U2)	
U35	202211-000	IC, linear, 33174, low power quad op amp	<b></b>
U36	201641-000	IC, digital, 74HC76, J-K Flip-Flop	<b>—</b>
U37	200359-003	IIC, linear, 5MC1404U5, voltage reference	<u></u>

5-3 Ref	Part Number	Description	Use Code	Qty
8£П	200673-007	IC, linear, UA7812ACLP, voltage regulator, +12V, TO-92		1
U39	200706-053	IC, linear, 79L12, voltage regulator, -12V, TO-92		_
U41,42	201743-001	IC, linear, DG509, multiplexer		2
U43	201342-070	IC, linear, DG302, analog switch 🛕		
U44,45	201214-002	IC, linear, LF444, quad op amp		2
U46	202156-000	IC, digital, DS1210, nonvolatile controller		_
U47	200854-006	IC, linear, LM3342, adjustable current source, TO-92		,_
U48	250012-000	IC, digital, DG411DY, quad analog switch SMT		1
U49		(Same as U7)		
U50	202020-013	IC, oscillator, crystal, 12MHz, 5V		_
U51	201654-000	IC, digital, MP7528, dual D/A		_
U52	201444-000	IC, linear, LM358, dual op amp		1
U53	201577-000	IC, digital, 74HC08, 2-input AND		1
U54	202256-000	IC, digital, 74HC173, quad flip-flop		
U55	201486-000	IC, linear, LF442, op amp		1
U56	200933-000	IC, linear, LM336, voltage regulator, -2.5V		-
U <b>57</b>	201596-002	IC, linear, LM385-1, voltage regulator		,
VSP1-3	200324-011	Voltage Surge Protector, 470V, 1pF		υ
VSP4	200324-016	Voltage Surge Protector, 8-10.3kV, 2pF		) mai
VSP5,6		(Same as VSP1)		
XUI	202054-000	Socket, DIP, 64 contact		2
XU9	201525-000	Socket, 28 contact		<b>)—</b>
XU14, 22	200675-026	Socket, DIP, 28 contact		4
XU25	200675-024	Socket, 20 contact		<del></del>
XU27, 28		(Same as XU14)		
XU30		(Same as XU1)		
ы	200536-002	Cable Tie, 8 L x 0.13 W, self-locking (not shown, used on U28)		
ယ	202245-001	RivNut, aluminum, 3/32 D x 0.187 L (not shown, used on P7)		2
4	200624-004	Tubing, teflon, #22		A/R
Un	202412-000	Wire, solid, 22 AWG (not shown, used on L2, L3)		A/R









LIFEPAK 9P defibrillator/monitor/pacemaker

Parts Lists and Schematics

5-22

TO PATIENT CONNECTOR


LIFEPAK 9P defibrillator/monitor/pacemaker

N . <u>1</u> /

....




All Rights Reserved



Main PCB Assembly
Part Number 805474-01
(Sheet 4 of 6)




All Rights Reserved




LIFEPAK 9P defibrillator/monitor/pacemaker

Parts Lists and Schematics

Ŋ





.

Ref	Part Number	Description	Use	2
A2	803726-01	Power Supply PCB Assembly		REF
C1	BOS240005-221	Capacitor, 220pF/250V, 20%		بد
ន	BOS220073-473	Capacitor, 0.047µF/250V, 20%		_ (
CG	BOS200113-227	Capacitor, 200µF/400V		<u> </u>
2	BOS230001-473	Capacitor, 0.047µF/630V		<u></u> ,
S	BOS220019-681	Capacitor, 680pF/1.5kV, 10%		<u> </u>
<u>ς</u>	BOS210031-336	Capacitor, 33µF/50V		<u> </u>
C	BOS240011-221	Capacitor, 220pF/50V, 10%		<u> </u>
Q	BOS220010-104	Capacitor, 0.1µF/63V, 5%		۰ در
C10-13	BOS210046-108	Capacitor, 1000µF/35V		4
C14	BOS220071-105	Capacitor, 1µF/63V, 10%		<b>→</b> .
C15		(Same as C8)		•
C18		(Same as C8)		
C19	BOS220041-332	Capacitor, 3300pF/100V, 5%		
C20	BOS220011-224	Capacitor, 0.22µF/63V, 5%		<u></u>
C21,22		(Same as C1)		
C23	BOS220075-334	Capacitor, 0.33μF/250V, 20%		, <u>.</u>
CRI	BOS100006-600	Diode, bridge, 600 V, 4A		<u>-</u>
CR2	BOS100014-600	Diode, BYV26C, 600V, 1A		- L
CR3	BOS100015-600	Diode, 1N4937, 600V, 1A		<u> </u>
CR4	BOS110005-015	Diode, Zener, 1N4744A, 15V/17mA, 1W, 5%		<u>, , , , , , , , , , , , , , , , , , , </u>
CR5-8	BOS100065-075	Diode, 1N4448, 0.2W		7
CR10	BOS100038-200	Diode, BYV72-200, 200V, 30A		<u> </u>
CR11	BOS100049-045	Diode, 45V, 16A		
CR12	BOS150002-100	Diode, MCR72-3, 100V, 8A		<u></u> ;-
CR13	BOS110005-020	Diode, Zener, 1N4747A, 20V, 1W, 5%		
CR14-16		(Same as CR5)		
<u> </u>	BOS450089-06	Connector, locking, 6 contact	<u> </u>	1
2	BOS450089-12	Connector, locking, 12 contact		<b>→</b> 1
<b>P1-</b> 3	BOS430064-22	Wire, bus, 22 AWG		ω
	BOS800279	Inductor, balanced		
1.2	BOS800280	Inductor, common mode		<del> </del>
Q1	BOS130092	Transistor, HEXFET, IRF840, N-channel	··········	-1
Q2	BOS130060	Transistor, NPN, ZTX649, TO-92		1

E 4	Part Number	Description	Use Code	Qty
Q3	BOS130080	Transistor, FET, 2N7000, N-channel		1
2	BOS120024	Transistor, PNP, MPS4250, TO-92		
2	BOS300040-105	Resistor, 1 MΩ, 1/2W, 5%		<b>)—</b>
R2	BOS300020-182	Resistor, 1.8kΩ, 1/4W, 5%		-
R3	BOS300040-471	Resistor, 470Q, 1/2W, 5%		<b>—</b>
<b>R</b> 4	BOS300020-393	Resistor, 39kΩ, 1/4W, 5%		,
R5	BOS300020-750	Resistor, 75Ω, 1/4W, 5%		2
Ro	BOS300020-4.7	Resistor, 4.7Q, 1/4W, 5%		
R7	BOS300020-102	Resistor, 1kΩ, 1/4W, 5%		w
R8	BOS320010-7680	Resistor, 768Ω, 1/4W, 1%		<b></b>
R9, 10	BOS300020-130	Resistor, 13Ω, 1/4W, 5%		2
R11, 12	BOS33002015	Resistor, 0.15Ω, 1W, 5%		2
R13	BOS330031-331	Resistor, 330Ω 1/4W, 5%		_
<b>7</b> 4	BOS330031-472	Resistor, 4.7kΩ, 5W, 5%		<b>ا</b> سط
R15	****	(Same as R5)		
R16	BOS330031-151	Resistor, 150Ω, 5W, 5%		<b></b>
R17	BOS300020-330	Resistor, 33Ω, 1/4W, 5%		<b>,</b>
R18	BOS300020-100	Resistor, 10Q, 1/4W, 5%		<b>—</b>
R19	BOS300020-103	Resistor, 10kΩ, 1/4W, 5%		<b></b>
R20	BOS300040-151	Resistor, 150Q 1/2W, 5%	_	Н
R21	BOS320010-6341	Resistor, 6.34kQ, 1/4W, 1%		<b>)1</b>
R22	BOS320010-6490	Resistor, 649Q 1/4W, 1%		<b>)—-</b> 4
R23	BOS380002-201	Potentiometer, 20002		<u> </u>
R24		(Same as R7)		
R25	BOS300020-682	Resistor, 6.8kΩ, 1/4W, 5%		نسر
R26	BOS300040-751	Resistor, 750Q, 1/2W, 5%		2
R27, 28	BOS300020-332	Resistor, 3.3kΩ, 1/4W, 5%		2
R29		(Same as R26)		
R30		(Same as R7)	•••	
R31	BOS300020-560	Resistor, 56Ω, 1/4W, 5%		
R32	BOS300020-752	Resistor, 7.5kQ, 1/4W, 5%		<u> </u>
R33	BOS300020-622	Resistor, 6.2kΩ, 1/4W, 5%		
FI	BOS390013	Thermistor, 100, 4A, 15%		<del></del>
RT2	BOS390002	Thermistor, 2200, 5%		<u></u>
RVI	BOS390020	Varistor, 75J, 387Vmin		مستدن برور خصا
-	BOS 800278	Transformer III 544/VDE		<u> </u>
3 :	BO8800026	Transformer miles		
14	DOOOOOO	AT GHENRUR HAVE, PULLOV		-

5-4 Ref	Part Number	Description	Use Code	9 <del>0</del>
VRI	BOS140017	IC, linear, TL431, adjustable shunt regulator		<del>  </del>
2	BOS640007	Adhesive, RTV		A/R
ເມ	BOS500319	Cover (not shown)		<u>-</u>
4.	BOS500318	Heat Sink, diode		<u> </u>
U	BOS500133	Heat Sink, transistor		<u> </u>
6	BOS520058	Insulator (not shown, used on Q1, CR10-12)		4
7	BOS530008	Label, model number (not shown, used on C3)		<u></u>
90	BOS530072	Label, serial number (not shown, used on item 5)		_
9	BOS600006-0375	Screw, pan head, 8-32 x 0.375 L (not shown, used on item 3)		4
10	BOS610011	Spacer, glass, 0.220 x 0.062		2
	BOS500062	Spring Clip (used on Q1, CR10 to 12)		4
12	BOS840021	Tape, 0.0625 x 0.500 L (not shown, used on L2)		A/R
13	BOS600067	Washer, flat, #8 (not shown, used on item 3)		4
14	BOS600083	Washer, lock, #8 (not shown, used on item 3)		4




Figure 5–4 Power Supply PCB Assembly
Part Number 803726-01
(Sheet 1 of 2)

<del>ပ်</del> မ




Figure 5-4 Α2

Power Supply PCB Assembly
Part Number 803726–01
(Sheet 2 of 2)

Parts Lists and Schematics 5

5–5 Ref	Part Number	Description	Use Code	Ωįγ
A3	803724-03	Power Conversion PCB Assembly		REF
CI	202210-000	Capacitor, 4700pF, 9kV		-
S	201628-018	Capacitor, 1.0µF/50V, 5%		ယ
C3,4	200893-046	Capacitor, 1.0µF/50V, 10%		ယ
C	201628-077	Capacitor, 0.68µF/50V, 10%		2
S	200264-010	Capacitor, 0.01µF/50V, 20%		19
C7,8	202102-030	Capacitor, electrolytic, 3300µF/25V, 20%		4
G	201554-068	Capacitor, electrolytic, 22000µF/25V, 100%		<b></b>
C10	202068-001	Capacitor, 20µF/50V, 5%		
C11,12	201628-000	Capacitor, 0.001µF/50V, 5%		ယ
C13		(Same as C2)		
C14	202068-027	Capacitor, 20µF/50V, 20%		2
C15	200901-031	Capacitor, 560pF/50V, 5%		
C16		(Same as C6)		
C17		(Same as C11)	· .	
C18	200205-035	Capacitor, electrolytic, 10µF/35V, +50%		2
C19	201628-072	Capacitor, 0.1µF/50V, 10%	· ·	υş
C20	200901-030	Capacitor, 470pF/50V, 5%		1
C21		(Same as C6)		
C22	200901-058	Capacitor, 0.1µF/50V, 5%		6
C24,25		(Same as C7)		
C26		(Same as C18)		
C27,28		(Same as C6)		
C29		(Same as C22)		
C30		(Same as C14)		
C31	202102-004	Capacitor, electrolytic, 3300µF/6.3V, 20%		<del>ب</del>
C32	202102-021	Capacitor, electrolytic, 2200µF/16V, 20%		2
C33		(Same as C19)		
C34		(Same as C5)		
C35	/	(Same as C3)		
C36,37		(Same as C6)		
C38	201628-006	Capacitor, 0.01µF/50V, 5%	·	4
C39		(Same as C6)		
C40		(Same as C22)		
2		(Same as C38)		
C42		(Same as C6)		
C43		(Same as C22)	-	
C44,45		(Same as C6)		

	Part Number	(Same as C22)	Code
-50		(Same as C6)	
C51 20	201628-060	Capacitor, 0.001µF/50V, 10%	
C52,53		(Same as C6)	
C54		(Same as C38)	
C55		(Same as C22)	
C56		(Same as C38)	
C57,58		(Same as C6)	
C59		(Same as C2)	
C60		(Same as C19)	
C61		(Same as C32)	
CR2 20	200635-053	Diode, 1.5KE82A, 113V	
CR3 20	200971-000	Diode, 1N914B, 8ns, PIV 75	
CR4 20	200959-003	Diode, 2N2324, 100V	
CR5		(Same as CR3)	
CR6 20	200083-002	Diode, 1N5819, 1A/40V	
CR7 20	201528-005	Diode, rectifier, 8A/100V	
CR8		(Same as CR3)	
CR11 20	201527-000	Wire, jumper, 0 Ω, 1/4W, 5%	
CR12 20	201128-000	Diode, 3XF1, 300ns, PIV 100	
CR1315		(Same as CR3)	
CR17-19		(Same as CR3)	
CR20-23 20	200605-010	Diode, Zener, 1N4738A, 55mA/8.2V	
CR24		(Same as CR3)	
CR25 20	200086019	Diode, Zener, 1N4104, 24.8μΑ/10V	
CR26 20	201391-005	Diode, SB320, 3A, PIV 80	
CR27 20	200083-004	Diode, 1N5821, 3A/30V	
CR28 20	200103-015	Diode, Zener, MPTE-22, 40A/22V	
CR29 20	200658-002	Diode, MR812, 200ns, PIV 200	
CR31		(Same as CR3)	
CR32 20	200086-011	Diode, Zener, 1N4625, 55mA/5.1V	
CR33-36		(Same as CR3)	
CR37-41 20	200991-000	Diode, 1N5711, 15mA	
CR42,43		(Same as CR20)	
CR44,45		(Same as CR37)	
CR46		(Same as CR3)	
CR47-53		(Same as CR20)	
	200605018	Diode, Zener, 1N4746A, 250mA/18V	


R12	<b>%</b> 11	R10	R9	R8	R7	₽6	R5	R4	R3	R1,2	(15	2	014	Q13	Q12	Q11	Q10	Q6-9	Q5	Q2	Q3	Q1,2	714	118	<b>§17</b>	J15,16	113	J12	J10,11	Ü	Pany .	CR59	CR58	CR57	CDZX	Ref
200054-250	200054-388	200054-317		200054-091	200054-193	200676-030	200054-288	200054-341	200054-370	801884-00						202076000		200585-001	201738001	200994-001	200942-000	201696-000	201678-015	202128024	201191219		202151-204	202151-200	201504-003	201582-002	202293-000					Part Number
Resistor, 4.02kΩ, 1/8W, 1%, RN55C	Resistor, 110kΩ 1/8W, 1%, RN55C	Resistor, 20kΩ 1/8W, 1%, RN55C	(Same as R7)	Resistor, 86.60, 1/8W, 1%, RN55C	Resistor, 1kQ, 1/8W, 1%, RN55C	Resistor, 120Q 3W, 5%	Resistor, 10kΩ, 1/8W, 1%, RN55C	Resistor, 35.7kΩ, 1/8W, 1%, RN55C	Resistor, 71.5kQ, 1/8W, 1%, RN55C	Resistor, 64MΩ 1W, 1%	(Same as Q11)	(Same as 011)	(Same as Q5)	(Same as Q6)	(Same as Q1)	Transistor, FET, IRFD120, N-channel	(Same as Q4)	Transistor, FET, VN0106N3, TO-92	Transistor, FET, VP0106N3, TO-92	Transistor, FET, IRF9531, TO-220	Transistor, NPN, 2N2222A, TO-18	Transistor, FET, IRF540, TO-220AB	Connector, DIP, 64 contact	Connector, 2 contact	Connector, 14 contact	(Same as J10)	Connector, locking, 6 contact	Connector, locking, 2 contact	Terminal, quick connect, 0.187 W/0.020 T	Connector, receptacle, 25 contact	Fuse, fast blow, 10A/125V, pigtail	(Same as CR37)	(Same as CR20)	(Same as CR3)	(Kome of CD27)	Description
										- '							¥.													٠						
					· · · · · · · · · · · · · · · · · · ·										•														•			٠				Code
2	1	6		<u> </u>	24		31		ω	2						2		<u>ر</u>	2	2	1	ယ	<b>—</b>	<b>–</b>	· <u> </u>		<u> </u>	سر	4	<u> </u>	 —			•		Š

Description 1%, RN55C 1, RN55C RN55C RN55C , RN55C %, RN55C %, RN55C %, RN55C %, RN55C %, RN55C %, RN55C
--

Ref	Part Number 200054-413 200054-265 200054-334 200054-508	(Same as R48) (Same as R36) (Same as R20) (Same as R12) (Same as R20) (Same as R10)
······································	200054–334	Resistor, 30.1kΩ, 1/8W, 1%, RN55C (Same as R15) (Same as R7) (Same as R10)
	200054-508	Resistor, 2MΩ, 1/8W, 1%, RN55C (Same as R20) (Same as R10) (Same as R7) (Same as R5)
- WAS	200054-412	Resistor, 196kQ, 1/8W, 1%, RN55C (Same as R7) (Same as R26) (Same as R7)
R98,99 R100,101	200054-226	(Same as R5) (Same as R7)
R102 R103 R104		(Same as R7) (Same as R98) (Same as R5)
R105		(Same as R5) (Same as R7) (Same as R5)
R107		(Same as R7) (Same as R5)
NIOO-ILI		(Same as R7)


5-5 Ref R114 R115 R116 R118 R119 R120	Part Number 200054-097	Description (Same as R98) (Same as R7) (Same as R98) Resistor, 100Ω 1/8W, 1%, RN55C (Same as R5) (Same as R7) (Same as R7)		Use Code
R123 R124		(Same as R44) (Same as R5)		
R126,127 R128 R129	Manufacture and a second se	(Same as R5) (Same as R7) (Same as R5)		
R129 R130 R131		(Same as R5) (Same as R98) (Same as R73)		
R133 R134	200054-280	Resistor, 8.25kQ, 1/8W, 1%, RN55C (Same as CR11)		
T2 7	804295-01 804294-01 801853-02	Transformer, flyback Transformer, power supply Transformer, test load		
U1 U2	201444-000 200345-000	IC, linear, LM358N, op amp IC, linear, LM393, op amp		
U4 U5,6	200933-000	IC, linear, 2843N, pulse-width modulator IC, linear, LM336BZ, voltage regulator, 2.5V, T0-92 (Same as U2)	92	92
U7 U8	202086000 200669002	IC, linear, AD584JH, voltage regulator, TO-99 IC, linear, LM339, comparator		
U10	202078-003	onve		
U11,12	201469-000	IC, digital, 74HC00, NAND gate	·	
U13	201705-000	1		
U14 U15	202142-000 201608-000	IC, digital, 74HC27, NOR gate IC, digital, 74HC02, NOR gate		
U16	201945-000	HC86, exclusive O		
U17 U18	200361-001	(Same as U1) IC, digital, 7541KN, D/A converter ♠		
U19,20	202085000			


5-5 Ref	Part Number	Description	Use Code	Q
U21		(Same as U4)		:
12	202221-007	Adhesive, hot melt (not shown, used on C7-9, C24-25)		Α/R
Ç.)	202091-035	Heat Sink		2
4	201152002	Insulator, component mount, TO-5 (not shown, used on CR4)	a, 1°	j-d
Si	201797-003	Insulator, silipad, transistor, TO-220		ယ
6	800943-09	Label, symbol, international, high voltage (not shown, used on T1)		_
7	201508-00	Nut, kep, lock, 4-40 x 0.312 L	·	ω
<b>90</b>	201199-100	Nut, 4-40, metal (used on Q1, Q2)		<u>-</u>
9	202245-001	Rivnut, dome head, 3/32 D x 0.187 L (not shown, used on P14)		2
10	200476-761	Screw, pan head, 4-40 x 0.312 L (used on Q12)		<b>-</b>
jasek jerek	200476-762	Screw, pan head, 4-40 x 0.375 L (used on Q1, Q2)		<u> </u>
12	201874-005	Screw, pan head, internal lock washer, 4-40 x 0.312 L		2
3		Sense Leads PCB Assembly (Not repairable)		_
, m.t.	201756-012	Tubing, PVC, clear #10, 105°C (used on Item 17)		Α⁄R
5	202413-000	Washer, lock (used on Q1, Q2)		<u>-</u>
16	200431-849	Washer, nylon (used on Q1, Q2, Q12)		ω
17	803783-31	Wire Harness, diode/jumper		<del>,</del> .



5-4-1

Parts Lists and Schematics 5





Ġ

.

			(Same as CR2)		CR4
14		H	Diode, switching, 1N914B, 8ns, PIV 75, SMT	250003-000	CR2
<b>—</b>	1000110		Diode, switching, IN914B, 8ns, PIV 75	200971-000	CR1
					***************************************
			(Same as C206)		C210
4			Capacitor, 1000pF/50V, 5%	201628-000	C206-208
			(Same as C7)		C41
			(Same as C3)	egge gyfarid	C38
			(Same as C18)		C37
			(Same as C33)		C35,36
			(Same as C7)		C34
w			Capacitor, 1µF/50V, 10%	200893-046	C33
			(Same as C3)		C31
1			Capacitor, electrolytic, 1200µF/25V	202279-183	C30
			(Same as C3)		C28,29
2			Capacitor, electrolytic, 330µ/25V	202102-024	C26,27
<b>)—</b> 4			Capacitor, 0.22µF/50V, 10%, SMT	250017-136	C25
			(Same as C7)		C23,24
2			Capacitor, 200pF/500V, 1%	200274-115	C21,22
			(Same as C7)	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	C20
ယ			Capacitor, 0.01µF/100V, 5%, SMT	250017-012	C18,19
			(Same as C7)		C17
<u> —</u>			Capacitor, 2.7nF/50V, 5%, SMT	250018-052	C16
			(Same as C7)		CLS
<b>,</b>			Capacitor, electrolytic, 150µF/400V	202159-051	C14
			(Same as C9)		C13
			Capacitor, 0.033µF/50V, 1%	200744-240	CI
			(Same as C1)		CIO
2			Capacitor, 0.1µF/50V, 1%	200744-205	C
9			Capacitor, 1000pF/100V, 5%, SMT	250018-036	C7,8
			(Same as C4)		Cø
1			Capacitor, 1µF/50V, 5%	201628-018	C5
2			Capacitor, 150pF/100V, 5%, SMT	250018-026	C4
υ			Capacitor, 0.1µF/50V, 5%, SMT	250017-035	C
<b>-</b>			Capacitor, 0.47µF/50V, 5%	201628-016	S
2			Capacitor, 100pF/100V, 5%, SMT	250018-024	<u>C</u>
<b>7</b>			Interconnect Pacing PCB Assembly	805472-02	A5
			de de la companya de		
Ş Q	Code		Description	Part Number	무선 5-6

5-6 Ref	Part Number	Description	Use	\$
CR5	250022-022	Diode, Zener, MBZ5248, 18V, 5%, 225mW, SMT		
CR6		(Same as CR2)		,
CR8	200177-051	Diode, Zener, P6KE75, tranzorb		<b></b>
CR9-14		(Same as CR2)		
CR16,17	202136-007	Diode, rectifier, UF4004, 400V, 1A		2
CR18	202136-006	Diode, rectifier, 1000V, 1A		
CR19,20		(Same as CR2)		
CR21	200086-019	Diode, Zener, 1N4104, 10V, 8mA, 0.2W, 5%		<b>-</b>
CR22,23	Makada kanan	(Same as CR2)		
CR24	202275-020	Diode, Zener, 1N4698, 11V, 0.25W, 5%		)
CR25		(Same as CR2)		
CR27	200605-009	Diode, Zener, 1N4737A, 7.5V, 1W, 5%		<b></b> -
CR28	200284-005	Diode, silicon, 1N4005, 1A, PIV600	***************************************	<b>,</b>
CR201,202	200605-010	Diode, Zener, 1N4738A, 8.2V, 550mA, 1W, 5%		<b>∞</b>
CR203	200605-018	Diode, Zener, 1N4746A, 18V, 250mA, 1W, 5%		ر. د
CR204		(Same as CR201)		
CR205		(Same as CR203)		
CR206-209		(Same as CR201)		- TPS://www.inc
CR210,211		(Same as CR203)		
CR212		(Same as CR201)		
CR213		(Same as CR203)		Way
<b>J</b> 7	201678-035	Connector, male, 96 contact		<del></del>
44	201678-028	Connector, 64 contact		-1
J19	202258-014	Connector, 26 contact	.,	<del> </del>
J36	202128-198	Connector, right-angle, 2 contact		<del> </del>
	200808001	Inductor, bead (installed on Q3; refer to Figure 5-6)		<b></b>
Q1	250021-001	Transistor, NPN, PN2222A, 40V		<del> </del>
82	202203-001	Transistor, p-channel, power FET, MPT2P50, 500V		<b>James</b>
Q3	201569-000	Transistor, NPN, power FET, IRF830		<u> </u>
Q4,5	200282-004	Transistor, FET, 2N4393		2
Q6	202164-000	Transistor, n-channel, power, MTD2N50, 500V		<del></del> -
Q7,8	200585-001	Transistor, VFET, VN0106N3, 60V		10
Q	201696-000	Transistor, n-channel, power FET, IRF540		<u>, , , , , , , , , , , , , , , , , , , </u>
Q10		(Same as Q7)	······································	
QII	201738-003	Transistor, p-channel, power FET, VP0808L		<b></b>
		The state of the s		

5-6 Ref Q12,13	Part Number	(Same as Q7)	
Q12,13 Q15		(Same as Q7)	
Q204,205		(Same as Q7)	
Q213,214		(Same as Q7)	
<b>F1</b>	201251-157	Varistor, 390V, 4500A	en e
R2,3	200578072	Resistor, 100Q, 1W, 5%	
R4	250001-004	Resistor, 11Ω, 1/8W, 1%, SMT	
RS	250001-164	Resistor, 511Ω, 1/8W, 1%, SMT	
R6	250001-287	Resistor, 10kΩ, 1/8W, 1%, SMT	
₹7	250001-410	Resistor, 191kQ, 1/8W, 1%, SMT	
R8	250001–248	Resistor, 3.92kΩ, 1/8W, 1%, SMT	
R9	250001-096	Resistor, 100Ω, 1/8W, 1%, SMT	
R10	250001-383	Resistor, 100kΩ, 1/8W, 1%, SMT	
R11	250001–258	Resistor, 4.99kΩ, 1/8W, 1%, SMT	in the second of
R12,13	,	(Same as R9)	
R14		(Same as R10)	
R15	250001-420	Resistor, 243kΩ, 1/8W, 1%, SMT	
R16	250001-316	Resistor, 20kΩ, 1/8W, 1%, SMT	- 1 - 2
R17		(Same as R6)	
R18	250001-412	Resistor, 200kQ, 1/8W, 1%, SMT	
R19		(Same as R16)	
R20		(Same as R18)	
R21	250001-390	Resistor, 118kQ, 1/8W, 1%, SMT	
R22	250001-305	Resistor, 15.4kQ, 1/8W, 1%, SMT	
R23,24	200471-124	Resistor, 150kΩ, 1/2W, 5%	
R25	•	(Same as R15)	
R27	250001-192	Resistor, 1kΩ, 1/8W, 1%, SMT	
R28		(Same as R6)	
R29	250001-371	Resistor, 75kΩ, 1/8W, 1%, SMT	
R30		(Same as R6)	
R31	250001-387	Resistor, 110kΩ 1/8W, 1%, SMT	
R32	250001-416	Resistor, 221kQ, 1/8W, 1%, SMT	
R33,34	250001-426	Resistor, 280kΩ, 1/8W, 1%, SMT	
R35	250001-225	Resistor, 2.21kΩ, 1/8W, 1%, SMT	
R37		(Same as R29)	
R38		(Same as R18)	·.
R39		(Same as R16)	

Red         250001-346         Resistor, 1216A, 188V, 1%, SMT         Code         Lity           R41         250001-346         Resistor, 1216A, 188V, 1%, SMT         1           R43         250001-298         (Same as R1)         2           R44         250001-299         Resistor, 134A, 18V, 1%, SMT         1           R46         250001-422         Resistor, 204A, 18W, 1%, SMT         1           R46         250001-423         Resistor, 204A, 18W, 1%, SMT         1           R49         250001-423         Resistor, 204A, 18W, 1%, SMT         1           R45         250001-423         Resistor, 204A, 18W, 1%, SMT         1           R45         250001-423         Resistor, 204A, 18W, 1%, SMT         1           R44         250001-404         Resistor, 204A, 18W, 1%, SMT         1           R43         250001-418         Resistor, 274A, 18W, 1%, SMT         1           R45         250001-418         Resistor, 274A, 18W, 1%, SMT         1           R46         250001-412         Resistor, 274A, 18W, 1%, SMT         1           R47         260470-016         Resistor, 294A, 18W, 1%, SMT         1           R46         250001-312         Resistor, 294A, 18W, 1%, SMT         1           R47	5 G	Dark Nicombag		Use	?
250001–391 Resistor, 121kQ, 1/8W, 1%, SMT  (Same as R11)  250001–298 Resistor, 13kQ, 1/8W, 1%, SMT  (Same as R10)  (Same as R10)  (Same as R10)  (Same as R1)  250001–422 Resistor, 20kQ, 1/8W, 1%, SMT  250001–491 Resistor, 20kQ, 1/8W, 1%, SMT  250001–491 Resistor, 0.301Q, 1W, 1%  (Same as R23)  250001–312 Resistor, 18.2kQ, 1/8W, 1%, SMT  200470–016 Resistor, 18.2kQ, 1/8W, 1%, SMT  200470–016 Resistor, 18.2kQ, 1/8W, 1%, SMT  200001–312 Resistor, 18.2kQ, 1/8W, 1%, SMT  (Same as R6)  250001–400 Resistor, 18.2kQ, 1/8W, 1%, SMT  (Same as R10)  (Same as R10)  (Same as R6)  200054–151 Resistor, 1000MQ, 1/4, 1/9, SMT  (Same as R10)  (Sam	R40	250001-346		1	
(Same as R11)  250001–298 (Same as R1)  250001–029 (Same as R6)  250001–425 (Same as R6)  250001–425 (Same as R7) (Same as R1) (Same as R2) (Same as R3)  200216–004 (Resistor, 2740, 1/8W, 1%, SMT 200470–016 (Same as R6)  250001–312 (Same as R7) (Same as R7)  250001–312 (Same as R7) (Same as R7)  250001–400 (Same as R7) (Same as R7)  250001–400 (Same as R7) (Same as R1) (Same as R	R41	250001–391	Resistor, 121kQ, 1/8W, 1%, SMT		2
250001–298 Resistor, 136.Q. 1/8W, 1%, SMT Same as R 10) Same as R 10) Same as R 10	R43		(Same as R11)		
(Same as R10) (Same as R6) (Same as R7) (Same as R1) (Same as R2) (Same as R2) (Same as R3) (Same as R4) (Same as R4) (Same as R7) (Same as R1)	R44	250001-298	Resistor, 13kQ, 1/8W, 1%, SMT		<u> </u>
Came as R6    Resistor, 20Q, 1/8W, 1%, SMT	R45		(Same as R10)		
250001–029 Resistor, 2002, 1/8W, 1%, SMT 250001–425 Resistor, 274kQ, 1/8W, 1%, SMT 250001–471 Resistor, 274kQ, 1/8W, 1%, SMT 250001–471 Resistor, 274Q, 1/8W, 1%, SMT 250001–138 Resistor, 274Q, 1/8W, 1%, SMT 250001–138 Resistor, 274Q, 1/8W, 1%, SMT 250001–312 Resistor, 18.2kQ, 1/8W, 1%, SMT 250001–400 Resistor, 18.2kQ, 1/8W, 1%, SMT (Same as RO) 200054–151 Resistor, 365Q, 1/8W, 1% (Same as RO) (Same as RO) (Same as RO) 200470–075 Resistor, 137kQ, 1/4W, 5% (Same as RO) 200470–075 Resistor, 1000MQ, 1W, 10% 250001–469 Resistor, 1000MQ, 1W, 10% 250001–479 Resistor, 1000MQ, 1/8W, 1%, SMT 82 250001–479 Resistor, 199kQ, 1/8W, 1%, SMT Resistor, 199kQ, 1/8W, 1%, SMT (Same as R16) (Same as R75) (Resistor, 1000MQ, 1/8W, 1%, SMT Resistor, 221MQ, 1/8W, 1%, SMT Resistor, 221MQ, 1/8W, 1%, SMT Resistor, 1000MQ, 1/8W, 1%, SMT Resistor, 221MQ, 1/8W, 1%, SMT Resistor, 1000MQ, 1/8W, 1%, SMT	R46		(Same as R6)		
250001–425 Resistor, 274kΩ, 1/8W, 1%, SMT 250001–471 Resistor, 825kΩ, 1/8W, 1%, SMT (Same as R1) (Same as R2) 200216-004 Resistor, 0.301Ω, 1W, 1% (Same as R23) 250001–138 Resistor, 1.2kΩ, 1/8W, 1%, SMT 250001–312 Resistor, 1.8 2kΩ, 1/8W, 1%, SMT (Same as R27) 250001–400 Resistor, 1.8 2kΩ, 1/8W, 1%, SMT (Same as R6) (Same as R6) (Same as R6) (Same as R10) (Same as R10) (Same as R6) (Same as R75)	R49	250001-029	Resistor, 200, 1/8W, 1%, SMT		_
250001-471   Resistor, 825kQ, 1/8W, 1%, SMT	R51	250001-425	Resistor, 274kQ, 1/8W, 1%, SMT		<del></del> -
(Same as R1) (Same as R23) (Same as R23) (Same as R23) (Same as R6) (Same as R7) (Same as R8) (Same as R6) (Same as R7) (Same as R6) (Same as R6) (Same as R6) (Same as R6) (Same as R7) (Same as R6) (Same as R7) (Same as R7) (Same as R7) (Same as R75)	R52	250001-471	Resistor, 825kQ, 1/8W, 1%, SMT		<b>,</b>
,55  200216-004  Resistor, 0.301Ω, 1W, 1%  (Same as RC)  250001-138  Resistor, 274Ω, 1/8W, 1%, SMT  200470-016  Resistor, 4.7Ω, 1/4W, 5%  250001-312  Resistor, 1.50kΩ, 1/8W, 1%, SMT  (Same as RC)  200054-151  Resistor, 355Ω, 1/8W, 1%, SMT  (Same as R RO)  200470-075  Resistor, 1.37kΩ, 1/4W, 5%  (Same as R RO)  201155-003  Resistor, 787kΩ, 1/8W, 1%, SMT  (Same as RO)  200054-451  Resistor, 787kΩ, 1/8W, 1%, SMT  (Same as RO)	R53		(Same as R1)		
200216-004 Resistor, 0.301Ω, 1W, 1%  (Same as R6)  250001-138 Resistor, 274Ω, 1/8W, 1%, SMT  200470-016 Resistor, 18, RQ, 1/8W, 1%, SMT  250001-312 Resistor, 150kΩ, 1/8W, 1%, SMT  (Same as R7)  250001-400 Resistor, 150kΩ, 1/8W, 1%, SMT  (Same as R6)  200054-151 Resistor, 355Ω, 1/8W, 1%  (Same as R10)  (Same as R10)  (Same as R10)  (Same as R6)  (Same as R62)  (Same as R62)  (Same as R63)  200155-003 Resistor, 137kΩ, 1/4W, 5%  (Same as R10)  (Same as R75)  (Same as R75)  (Same as R75)  (Same as R75)  Resistor, 1MΩ, 1/8W, 1%, RN55C  (Same as R75)  (Same as R75)  Resistor, 26.7kΩ, 1/8W, 1%, SMT	R54,55		(Same as R23)		سنة سنيت
(Same as R6)  250001–138  Resistor, 274Q, 1/8W, 1%, SMT  200470–016  Resistor, 18.2kQ, 1/8W, 1%, SMT  250001–312  Resistor, 15.0kQ, 1/8W, 1%, SMT  (Same as R27)  250001–400  Resistor, 555Q, 1/8W, 1%, SMT  (Same as R10)  (Same as R10)  (Same as R10)  (Same as R6)  (Same as R6)  (Same as R6)  (Same as R6)  200470–075  Resistor, 1.37kQ, 1/4W, 5%  (Same as R6)  201155–003  Resistor, 1000MQ, 1W, 10%  20001–469  Resistor, 787kQ, 1/8W, 1%, SMT  (Same as R10)  (Same as R10)  (Same as R10)  (Same as R75)  Resistor, 1MQ, 1/8W, 1%, SMT  (Same as R15)  (Same as R75)  Resistor, 1MQ, 1/8W, 1%, RN55C  (Same as R75)  Resistor, 267kQ, 1/8W, 1%, SMT	R56	200216-004	Resistor, 0.301Ω, 1W, 1%		<del></del>
250001–138 Resistor, 274Ω, 1/8W, 1%, SMT 200470–016 Resistor, 4.7Ω, 1/4W, 5% 250001–312 Resistor, 18.2kΩ, 1/8W, 1%, SMT (Same as R27) 250001–400 Resistor, 150kΩ, 1/8W, 1%, SMT (Same as R40) (Same as R10) (Same as R10) (Same as R10) (Same as R4) (Same as R2) (Same as R2) (Same as R6) 250001–469 Resistor, 1000MΩ, 1W, 10% (Same as R10) (Same as R75)	R57		(Same as R6)		
200470-016 Resistor, 4.7Ω, 1/4W, 5% 250001-312 Resistor, 18.2kΩ, 1/8W, 1%, SMT  250001-400 Resistor, 150kΩ, 1/8W, 1%, SMT  250001-400 Resistor, 365Ω, 1/8W, 1%, SMT  (Same as R6) (Same as R10) (Same as R10) (Same as R10) (Same as R6) (Same as R6)  200470-075 Resistor, 1.37kΩ, 1/4W, 5% (Same as R6)  201155-003 Resistor, 1000MΩ, 1W, 10%  250001-469 Resistor, 787kΩ, 1/8W, 1%, SMT  (Same as R10) (Same as R6)  200054-451 Resistor, 2.21MΩ, 1/8W, 1%, SMT  Resistor, 2.21MΩ, 1/8W, 1%, SMT (Same as R16) (Same as R75)	R58	250001-138	Resistor, 274Q, 1/8W, 1%, SMT	-	<b></b>
250001–312   Resistor, 18.2kΩ, 1/8W, 1%, SMT	R59	200470-016	Resistor, 4.7Ω, 1/4W, 5%		<b>)—</b> 4
(Same as R27)  250001-400  Resistor, 150kΩ, 1/8W, 1%, SMT  (Same as R6)  200054-151  Resistor, 365Q, 1/8W, 1%  (Same as R10)  (Same as R10)  (Same as R6)  (Same as R6)  (Same as R62)  (Same as R62)  (Same as R62)  (Same as R6)  201155-003  Resistor, 1000MΩ, 1W, 10%  250001-469  Resistor, 787kQ, 1/8W, 1%, SMT  (Same as R6)  200054-451  Resistor, 499kQ, 1/8W, 1%, RN55C  Resistor, 1MΩ, 1/8W, 1%, SMT  (Same as R16)  (Same as R75)  (Same as R75)  Resistor, 1MΩ, 1/8W, 1%, RN55C  (Same as R75)  Resistor, 1MΩ, 1/8W, 1%, RN55C  (Same as R75)  Resistor, 1MΩ, 1/8W, 1%, RN55C	R60	250001-312	Resistor, 18.2kΩ, 1/8W, 1%, SMT		<u> </u>
250001–400 Resistor, 150kΩ, 1/8W, 1%, SMT (Same as R6) (Same as R10) (Same as R10) (Same as R10) (Same as R10) (Same as R6) (Same as R62) (Same as R62) (Same as R62) (Same as R63) (Same as R62) (Same as R63) (Same as R63) (Same as R63) (Same as R64) (Same as R65) (Same as R10) (Same as R16) (Same as R175) (Same as R75)	R61		(Same as R27)		
(Same as R6)  Resistor, 365Q 1/8W, 1%  (Same as R10)  (Same as R10)  (Same as R6)  201155-003  Resistor, 1000MQ, 1W, 10%  250001-469  Resistor, 787kQ, 1/8W, 1%, 5MT  (Same as R10)  (Same as R10)  (Same as R10)  (Same as R16)  (Same as R16)  (Same as R16)  (Same as R16)  (Same as R75)  Resistor, 1MQ, 1/8W, 1%, SMT  (Same as R75)  Resistor, 1MQ, 1/8W, 1%, RN55C  Resistor, 22/MQ, 1/8W, 1%, RN55C  Resistor, 26/7kQ, 1/8W, 1%, SMT	R62	250001-400	Resistor, 150kQ, 1/8W, 1%, SMT		2
200054-151   Resistor, 365Ω, 1/8W, 1%   (Same as R10)	R63		(Same as R6)		
(Same as R10) (Same as R6) (Same as R6) (Same as R6) (Same as R62) (Same as R6) 201155-003 (Same as R10) (Same as R16) (Same as R16) (Same as R16) (Same as R16) (Same as R75)	R64	200054-151	Resistor, 365Q, 1/8W, 1%		<del></del>
(Same as R6) (Same as R10) (Same as R10) (Same as R62) (Same as R62) (Same as R62) (Same as R6) 201155-003 (Same as R6) 250001-469 (Same as R10) (Same as R16) (Same as R16) (Same as R75)	R66,67		(Same as R10)		*************************
(Same as R10)  Resistor, 1.37kQ, 1/4W, 5%  (Same as R62)  (Same as R62)  (Same as R62)  201155-003  Resistor, 1000MQ, 1W, 10%  250001-469  Resistor, 787kQ, 1/8W, 1%, SMT  (Same as R10)  (Same as R6)  200054-451  Resistor, 499kQ, 1/8W, 1%, RN55C  Resistor, 2.21MQ, 1/8W, 1%, SMT  (Same as R16)  (Same as R75)  Resistor, 1MQ, 1/8W, 1%, RN55C  (Same as R75)  Resistor, 1MQ, 1/8W, 1%, RN55C  (Same as R75)  Resistor, 26.7kQ, 1/8W, 1%, SMT	R68		(Same as R6)		
71 200470-075 Resistor, 1.37kΩ 1/4W, 5%  (Same as R62) (Same as R62) 201155-003 Resistor, 1000MΩ 1W, 10% 250001-469 Resistor, 787kΩ 1/8W, 1%, SMT  (Same as R10) (Same as R6) 200054-451 Resistor, 499kΩ 1/8W, 1%, RN55C 250001-479 Resistor, 1MΩ 1/8W, 1%, SMT (Same as R16) (Same as R75) (Same as R75) (Same as R75) Resistor, 2.7kΩ 1/8W, 1%, RN55C (Same as R75) Resistor, 2.7kΩ 1/8W, 1%, SMT	R69		(Same as R10)		
(Same as R62) (Same as R6) (Same as R6) 201155-003 Resistor, 1000MΩ, 1 W, 10% 250001-469 Resistor, 787kΩ, 1/8W, 1%, SMT (Same as R10) (Same as R6) 200054-451 Resistor, 499kΩ, 1/8W, 1%, RN55C Resistor, 1MΩ, 1/8W, 1%, SMT (Same as R16) (Same as R75) (Same as R75) Resistor, 1MΩ, 1/8W, 1%, RN55C (Same as R75) Resistor, 26.7kΩ, 1/8W, 1%, SMT	R70,71	200470-075	Resistor, 1.37kΩ, 1/4W, 5%		2
(Same as R6) 201155-003 Resistor, 1000MΩ, 1W, 10% 250001-469 Resistor, 787kΩ, 1/8W, 1%, SMT 250001-469 (Same as R10) (Same as R6) 200054-451 Resistor, 499kΩ, 1/8W, 1%, RN55C 250001-479 Resistor, 1MΩ, 1/8W, 1%, SMT (Same as R16) (Same as R16) (Same as R75) (Same as R75) (Same as R75) Resistor, 1MΩ, 1/8W, 1%, RN55C (Same as R75) Resistor, 2.7kΩ, 1/8W, 1%, SMT	R72		(Same as R62)		
201155-003 Resistor, 1000MΩ, 1W, 10% 250001-469 Resistor, 787kΩ, 1/8W, 1%, SMT (Same as R10) (Same as R6) 200054-451 Resistor, 499kΩ, 1/8W, 1%, RN55C 250001-479 Resistor, 1MΩ, 1/8W, 1%, SMT (Same as R16) (Same as R75) (Resistor, 2.7kΩ, 1/8W, 1%, SMT	R73	Į.	(Same as R6)		
250001-469 Resistor, 787kΩ, 1/8W, 1%, SMT  (Same as R10) (Same as R6)  200054-451 Resistor, 499kΩ, 1/8W, 1%, RN55C  250001-479 Resistor, 2.21MΩ, 1/8W, 1%, SMT  (Same as R16) (Same as R16) (Same as R75)	R74	201155-003	Resistor, 1000MΩ, 1W, 10%		<del></del>
-78 (Same as R10) (Same as R6) 200054-451 Resistor, 499kΩ, 1/8W, 1%, RN55C ,82 250001-479 Resistor, 2.21MΩ, 1/8W, 1%, SMT ,84 250001-512 Resistor, 2.21MΩ, 1/8W, 1%, SMT (Same as R16) (Same as R75)	R75	250001-469	Resistor, 787kQ, 1/8W, 1%, SMT		ယ
(Same as R6) 200054-451 Resistor, 499kΩ, 1/8W, 1%, RN55C 250001-479 Resistor, 1MΩ, 1/8W, 1%, SMT 250001-512 (Same as R16) (Same as R16) (Same as R75)	R76-78		(Same as R10)	- "	
200054-451 Resistor, 499kQ, 1/8W, 1%, RN55C  ,82 250001–479 Resistor, 1MQ, 1/8W, 1%, SMT  ,84 250001–512 Resistor, 2.21MQ, 1/8W, 1%, SMT  (Same as R16) (Same as R75)  200054-480 Resistor, 1MQ, 1/8W, 1%, RN55C  (Same as R75)  250001–328 Resistor, 26.7kQ, 1/8W, 1%, SMT	R79		(Same as R6)		
,82       250001–479       Resistor, 1MΩ, 1/8W, 1%, SMT         ,84       250001–512       Resistor, 2.21MΩ, 1/8W, 1%, SMT         (Same as R16)       (Same as R75)         200054–480       Resistor, 1MΩ, 1/8W, 1%, RN55C         (Same as R75)       (Same as R75)         250001–328       Resistor, 26.7kΩ, 1/8W, 1%, SMT	R80	200054-451	Resistor, 499kΩ, 1/8W, 1%, RN55C		<b></b>
,84       250001–512       Resistor, 2.21ΜΩ, 1/8W, 1%, SMT         (Same as R16)       (Same as R75)         200054–480       Resistor, 1ΜΩ, 1/8W, 1%, RN55C         (Same as R75)       (Same as R75)         250001–328       Resistor, 26.7kΩ, 1/8W, 1%, SMT	R81,82	250001-479	Resistor, 1MQ, 1/8W, 1%, SMT		ယ
200054-480 250001-328	R83,84	250001-512	Resistor, 2.21MQ, 1/8W, 1%, SMT		2
200054 <u>4</u> 80 250001 <u>3</u> 28	R85		(Same as R16)		
20005 <del>4-4</del> 80 250001–328	R86		(Same as R75)		
250001–328	R87	200054_480	Resistor, IMQ 1/8W, 1%, RN55C		<del></del>
250001–328	R88		(Same as R75)		
	R89	250001–328	Resistor, 26.7kQ, 1/8W, 1%, SMT		Н.

Fef	Part Number	Description
R90 R91	250001-342	Resistor, 37.4kΩ,1/8W, 1%, SMT (Same as R18)
R92		(Same as R16)
R93	250001-354	Resistor, 49.9kΩ, 1/8W, 1%, SMT
R94		(Same as R16)
R95		(Same as R10)
R96		(Same as R27)
R97		(Same as R29)
R98		(Same as R6)
R99		(Same as R8)
R100,101	-	(Same as R10)
R102	250001-245	Resistor, 3.57kQ, 1/8W, 1%, SMT
R104		(Same as R81)
R105		(Same as R41)
R106-108		(Same as R27)
R109		(Same as R9)
R110		(Same as R7)
R209	200054-226	Resistor, 2.21kQ, 1/8W, 1%, RN55C
R211	200470-071	Resistor, 910Q, 1/4W, 5%
R214,215		(Same as R211)
R220	200850-149	Resistor, 348Ω 1/2W, 1%
R221,222		(Same as R209)
R223	200850-175	Resistor, 649Ω, 1/2W, 1%
R224	***	(Same as R209)
R225		(Same as R211)
T <sub>1</sub>	805534-02	Transformer, pacer
UI	202160-000	IC, linear, TL431CLP, voltage reference
U2	250020-002	IC, linear, 33174, low-power quad op-amp, SMT
U3,4	250023-002	IC, linear, LM339M, comparator, SMT
US	201486-000	IC, linear, LF442, J-FET, op-amp
U6	202050-000	IC, linear, 3843, pulse-width modulator
U7	202116-001	IC, optocoupler-photo-transistor
U8	200359-003	IC, linear, 5MC1404U5, voltage reference
U9,10		(Same as U7)
UII		(Same as U2)
U12		(Same as U7)

5–6 Ref	Part Number	Description	Use Code	Qty
VSP1	200324-016	Voltage Surge Protector, B2-H100-XP, 10kV, 2pF		
2	202091-053	Heatsink, (used on Q9)		<u> </u>
w	202091-033	Heatsink, (used on Q2,3)		2
4	202196-013	Insulator, transistor, white		2
¢η.	201199-100	Nut, hex 440 x 0.188 W X 0.063		
<b>5</b> \	201773-000	Rivet, 11/32 L (not shown, used on J7,14)	· · · · · ·	4
7	200476-762	Screw, pan head, 4-40 x 0.375 L		_
<b>9</b> 0	200283-017	Tubing, heat-shrink, black, 0.125 ID		Α/R
\$	200283-003	Tubing, heat-shrink, red, 0.125 ID		Α/R
Đ	202413-000	Washer, lock #4 x 0.195 OD/0.25 T		
) — i	200431-849	Washer, nylon, shoulder		2





Figure 5-6

∑ טו

Interconnect Pacing PCB Assembly
Part Number 805472-02

(Sheet 1 of 5)





(37)


A5

CIRCUIT SIDE

Figure 5–6 Interconnect Pacing PCB Assembly
Part Number 805472–02
(Sheet 2 of 5)








Parts Lists and Schematics 5


LIFEPAK 9P defibrillator/monitor/pacemaker

5-52


Parts Lists and Schematics








5-7 Ref	Part Number	Description C	Use Code	Qty
A6	805469-10	Keypad Assembly, Front, LP9P, English	·	R H H
	805469-11	Keypad Assembly, Front I DOD Fnolish Intl		
	805469-13	Keypad Assembly, Front, LP9P, French		
	805469-14	Keypad Assembly, Front, LP9P, German		
	805469-15	Keypad Assembly, Front, LP9P, Spanish	···	
	805469-16	Keypad Assembly, Front, LP9P, Italian		
	805469-17	Keypad Assembly, Front, LP9P, Swedish		




80

Figure 5–7 Keypad Assembly
Part Number 805469

* **	
Α7	5–8 Ref
805464-00	Part Number
Rear Panel Membrane Switch	Description
	Use Code
REF	Qty





TO MAIN PCB

Α7

Figure 5–8 Rear Panel Membrane Switch
Part Number 805464–00
(Sheet 1 of 1)

	32	31	30	29	28	27	26	25	24	23	22	<b>1</b> →	20	5	18	17	16	15	4	3	12	)—1 )—1	10	9	∞	7	Ø.	Ŋ	4	w	2	A8	5-9 Ref
	GSI624-11893	GSI615-102331	GSI620-10244	GSI624-10262	GSI625-10263	GSI624-10246	GSIR41-50082	GSIR42-50083	GSIR41-22001	GSIR5E-519	GSIR5E-513	GSI615-200051	GSIR75-00007	GSIR52-41003	GSI420-06175	GSI610-102341	GSI624-10258	GSI624-10257	GSI624-10255	GSI624-10225	GSI624-10224	GSI625-10267	GSI625-10266	GSI625-10268	GSI624-10260	GSI624-10254	GS1624-10270	GSI624-10253	GSI624-10278	GSI624-11759	GSI624-10256	804189-00	Part Number
*Specialized alignment equipment required. Replacement of entire Strip Chart Recorder is recommended.	Tubing (not shown, used on item 17)	Spring Assembly	Spring	Shaft, linkage	Shaft Roller	Shaft	Screw, pan head	Screw, pan head, 4-40 x 0.25	Screw, flat head, 90 M2, OX, 4 x 6mm	Ring, rtng TRU 5133-9	Ring, rtng TRU 5100–18	Printroller Assembly	Printhead, KST 488MPDI (pre-assembled at GSI)*	O-Ring, IB Moore #2-104	Nameplate, label	Motor Assembly/Gear	Link, printhead	Link, latch plate	Latch, chassis	Label, paper feed	Keyshim Gear	Gear, motor 31T	Gear, motor 15T	Gear, idler 108T	Ejector Assembly, paper	Door, chassis	Cover, gear	Bumper, printhead	Bracket, pivot head	Bracket Mod, printhead (pre-assembled at GSI)*	Base, chassis	Strip Chart Recorder Assembly, GSI	Description
															.√				7			,									. '		Use Code
	A/R	نصر	2	2	<u> </u>	2	: 	2	ω	4			<u> </u>	 —		-	<b></b>	<del></del>	<b></b>	H	<del></del>	ш	_	<del></del>	<b></b>	_		2	Н	_	_	REF	Qly






Figure 5–9 Strip Chart Recorder Assembly
Part Number 804189–00
(Sheet 1 of 2)



**№** 

Figure 5-9 Strip Chart Recorder Assembly
Part Number 804189-00
(Sheet 2 of 2)

REF		Transfer Relay Assembly	800240-14	A9
Qty	Use Code	Description	Part Number	5–10 Ref






Figure 5–10 Transfer Relay Assembly
Part Number 800240–14
(Sheet 1 of 1)

A9

5-11 Ref	Part Number	Use Description Code		Oty
A10	802901-06	Paddle Assembly		REF
2	202221-000	Adhesive, hot melt, thermoplastic (not shown, used on item 8)		A/R
w	802930-06	Base, apex paddle	200 P	_
	802930-05	Base, sternum paddle	****	<u></u>
4.	803728-01	Cord Assembly, terminated, high voltage		<u>ш</u>
<b>U</b> I	802934-04	Cover, apex paddle		1
	802934-05	Cover, sternum paddle	U-04	<del></del>
9	802902-02	Electrode Plate Assembly, apex paddle		_
	802902-03	Electrode Plate Assembly, sternum paddle		<b>,</b>
7	802944-01	Gasket, paddle		2
<b>9</b> 0	802931-00	Handle, apex paddle (label not shown, see Figure 5-1)		
	802931-01	Handle, sternum paddle (label not shown, see Figure 5-1)	<u> </u>	<b>-</b>
\$	802932-00	Key, retainer, coil cord		2
10	800819-09	Label, serial number (not shown)		<u> </u>
11	200491-048	LED, panel mount (on apex paddle only)		<u> </u>
12	201508-004	Nut, lock, kep, 8-32 x 0.312/0.100 T		0
13	201103-537	Screw, pan head, black oxide, 4-40 x 0.312 L	**************************************	14
14	80423100	Shield, high voltage		2
15	804697-03	Switch, SPST (pushbutton not shown, see Figure 5-1)		w
16	200283-005	Tubing, heat shrink, red, 0.250 ID (on apex paddle only)		A∕R
	-			




Figure 5–11 Paddle Assembly
Part Number 802901–06
(Sheet 1 of 2)

Figure 5–11 Paddle Assembly
Part Number 802901–06
(Sheet 2 of 2)

A10

Fet 2	Part Number	Description	Use Code	Qty
A11	803706-00	Display PCB/CRT Assembly (DOT 7-038-0101)	Α	REF
	803706-01	Display PCB/CRT Assembly (DOT 7-038-0102)	В	RE T
N	DOT6-002-1678	Monitor PCB Assembly		<b></b>
بن	DOT73027887-913	CRT/Yoke Assembly		1
4	DOT-15005398-003	Cable Tie, 5-1/16in (not shown, used on item 3)	В	<del>-</del>
UI	DOT-048-0110	Cable Tie, nylon (not shown, used in five places)	· · · · ·	5
ō,	DOT15005406-000	Clamp, hose, nylon		نسر
7	DOT6-004-1150	CRT Connector Assembly	,	1
∞	DOT6-004-0452	Deflection Coil Assembly		<u> </u>
\$	DOT5-014-0393	Label, HV caution (not shown, used on item 3)		
10	DOT5-014-0430	Label, model, serial number (not shown, used in item 3)	···	<b>—</b>
jurad. Jurad	DOT32028970-007	Label, X-ray warning (not shown, used in item 3)	В	<b></b>
12	80411700	Mount, CRT, front		<b>—</b>
ن ن	804118-00	Mount, CRT, rear		H
4	DOT3-011-0520	Screw, pan head, 8-32 x 0.357 L		4
S	DOT-29002450-008	Shield, ferrite (not shown, used on item 3)	В	_
<b>1</b>	DOT2-049-0256	Tape, 3/4in adhesive (not shown, used in item 3)		A/R
17	DOT-2-049-0252	Tape, Copper adhesive 1.5in (not shown, used on item 3)	В	A/R
<b>⇔</b>	DOT-37003362-010	Tubing, heat-shrink, 1in black (not shown, used on item 3)	В	A/R
19	DOT3-013-0190	Washer, flat, #8		4
	7073 013 0163	Washer, lock, #8		4




Figure 5–12 Display PCB/CRT Assembly
Part Number 803706
(Sheet 1 of 4)

ΑH


5-12			Use	
2	GICIAGIISC	FOOTING AND ADDRESS OF THE PARTY OF THE PART	0000	
25 64 52 7 7	DOT6-002-1678	Monitor PCB Assembly (part of Display PCB/CRT Assembly)		大艺艺
C101	DOT1-012-2455	Capacitor, 0.1µF, 50V, 20%		13
C102	DOT1-012-2341	Capacitor, 1000μF, 25V		
C103	DOT1-012-2359	Capacitor, 0.15µF, 100V, 5%		<b>—</b>
C104	DOT1-012-2454	Capacitor, 0.01µF, 50V, 20%		2
C105	DOT1-012-2358	Capacitor, 100µF, 50V		2
C106, 107	DOT61220027-000	Capacitor, 0.1µF, 50V		2
C108	DOT1-012-2511	Capacitor, 5600pF, 50V, 5%		_
C109	DOT1-012-0300	Capacitor, 100pF, 500V, 5%		1
C110		(Same as C101)		
C111	DOT61025449-019	Capacitor, electrolytic, 22µF, 16V, 20%		_
C112	DOT1-012-2444	Capacitor, 1000μF, 16V, 20%		
C201	DOT1-012-2403	Capacitor, 1000pF, 100V		1
C202	DOT1-012-2680	Capacitor, 1000pF, 100V, 5%		
C203	DOT1-012-2211	Capacitor, 10µF, 25V		. 2
C204	DOT1-012-2347	Capacitor, 4/ULF, 35 V		-
C200-207	TOTAL 1990013 003	Consistent ONATHE THY SOL		
C209	DOT61220013-001	Capacitor, 0.01µF, 1kV, 5%		-
C210	DOT61220016-010	Capacitor, 2.2μF, 160V, 5%		<u> </u>
C211	DOT61220016-009	Capacitor, IµF, 160V, 5%		1
C212	DOT1-012-2419	Capacitor, 220pF, 500V, 5%		,
C213	-	(Same as C101)		
C214		(Same as C105)		
C301		(Same as C104)		
C302	DOT1-012-2214	Capacitor, 0.01µF, 1000V, 20%		υ J
C303	DOT1-012-2228	Capacitor, 3.3µF, 250V		2
C304, 305		(Same as C302)		
C306	DOT1-012-2375	Capacitor, 0.1µF, 500V, 20%		۲
C307		(Same as C303)		
C401		(Same as C203)		
C402		(Same as C101)		
C406, 407		(Same as C101)		
C408	DOT1-012-2250	Capacitor, 47µF, 100V		<u> </u>
C501, 502		(Same as C101)		
C503	DOT1-012-2627	Capacitor, 100µF, 35V, 20%		2
C504		(Same as C101)		
C505		(Same as C503)		
C506		(Same as C101)		

5–12 Ref	Part Number	Description	Use Code	Qty
D101	DOT1-021-0545	Diode, 1N4002, 100V, 1A		2
D201	DOT1-021-0517	Diode, Zener, Z0517, 18V		<u></u>
D202	DOT1-021-0497	Diode, 1N5282		2
D301	DOT58026488-002	Diode, RGP01-018, 1800V, 100mA		<u>-</u>
D302	DOT1-021-0507	Diode, 1N4937, 600V, 1A		
D303	DOT58005767-019	Diode, Zener, 1N997, 33V, 5%		,t
D401-404	DOT1-021-0550	Diode, 1N916, 100V, 75mA, 2%	·	4
D501		(Same as D101)		
D502, 503	DOT1-021-0543	Diode, UF4004, 400V		2
D506		(Same as D202)	• • •	
H201	DOT1-015-5081	Heat Sink		2
H401	DOT1-015-5134	Heat Sink		<b>-</b>
H501		(Same as H201)	W07-11-1	
1101	DOT30005029-002	Connector, 2 contact		
J201	DOT30005029-003	Connector, 3 contact		_
J901	DOT30257732-020	Connector, 20 contact		نسر
J902	DOT1-034-0416	Connector, 9 contact		<del></del>
L201	DOT1-016-0463	Coil, horizontal linearity	•	<b></b>
L401	DOT1-016-0397	Coil, fixed, 8.2µH	<b>F</b> ithirms	<u>-</u>
L501	DOT29002450-005	Ferrit Bead		<b>j</b> —4
Q201	DOT1-015-1144	Transistor, 2N3904, NPN		2
Q202	DOT1-015-1145	Transistor, 2N3906, PNP		ယ
Q203		(Same as Q201)		
Q204		(Same as Q202)		
Q205	DOT58025451-025	Transistor, FET, BUZ60, 400V, TO 220		<b>—</b>
0501	C071-C10-110G	Transistor, 2N3033A, NFN		<b>,</b>
Cool		(Same as Q202)		
Q502	DOT58025451-026	Transistor, FET, MTP10N15, 10A, 150V		<b>بسر</b>
R102	DOT1-011-2597	Resistor, 330Q, 1/4W, 5%		2
R104	DOT1-011-2551	Resistor, 10Q, 1/4W, 5%		4
R105	DOT1-011-2555	Resistor, 100kΩ, 1/4W, 5%		ω
R106	DOT1-011-5758	Potentiometer, 100kQ 1/8W, 20%	,	2
K107	DOT1-011-2553	Resistor, 1kΩ, 1/4W, 5%	-	10

5–12 Ref	Part Number	Description	Use Code	Q
R108	DOT1-011-2613	Resistor, 4.7kΩ, 1/4W, 5%		2
R109	DOT1-011-2862	Resistor, 300kΩ, 1/4W, 5%		
R110	DOT1-011-2608	Resistor, 390kΩ, 1/4W, 5%		···
R111	DOT1-011-5759	Potentiometer, 250kQ, 1/8W, 20%		_
R112	DOT1-011-2556	Resistor, 1MQ, 1/4W, 5%		2
R113		(Same as R106)		
R114	DOT1-011-2599	Resistor, 33kQ, 1/4W, 5%	18	2
R115	DOT1-011-2614	Resistor, 47kΩ, 1/4W, 5%		1
R116	DOT1-011-2582	Resistor, 220kΩ, 1/4W, 5%		-
R117	DOT1-011-2702	Resistor, 3.3Ω, 1/4W, 5%		<u>.                                    </u>
R118-120	DOT1-011-2620	Resistor, 5.6kΩ, 1/4W, 5%		4
R121	DOT1-011-2704	Resistor, 1Q, 1/2W, 5%		, ,
R122	DOT1-011-2258	Resistor, 330Q, 1/2W, 5%		· -
R201	DOT1-011-2579	Resistor, 220Q 1/4W, 5%		<u> </u>
R202		(Same as R102)		
R203	DOT1-011-2552	Resistor, 100Ω, 1/4W, 5%		ω
R204		(Same as R107)		
R205		(Same as R108)		
R206		(Same as R105)		14
R207		(Same as R107)		
R208	DOT1-011-2580	Resistor, 2.2kQ, 1/4W, 5%		
R209		(Same as R104)		*****
R210		(Same as R203)		
R211	DOT1-011-2619	Resistor, 560Q, 1/4W, 5%		<b>—</b>
R212	180113	(Same as R107)		
R213		(Same as R118)		
R214	DOT1-011-2554	Resistor, 10kΩ 1/4W, 5%		6
R215		(Same as R107)		
R216		(Same as R214)		
R217		(Same as R105)		
R218		(Same as R112)	- 	
R219		(Same as R114)		
R220	DOT69008520-025	Resistor, 24.9k\(\Omega\), 1/4W, 1\%		) mad
R221	DOT1-038-0212	Jumper Wire, 22 AWG		A/R
R223		(Same as R107)		
R224	DOT1-011-5822	Potentiometer, 1kQ, 1/8W, 20%		<b>—</b>
R225		(Same as R107)		•.
R226		(Same as R104)		3
R227		(Same as R107)		
R228	DOT1-011-2878	Resistor, 0.22Ω, 1W, 5%		1


<u>,</u>	W301, 302 DO	W205	W203 DO	U501 DO	U401 DO	U <b>201</b> DO	U101 DO	<b>T501</b> DO		R504	R503	R501, 502 DO	R415 DO	<b>R414</b> DO	R412, 413 DO	<b>R411</b> DO	<b>R410</b> DO	R409	R407 DO	R406 DO	R404	R403 DO	R402 DO	<b>R401</b> DO	R313 DO	R312	<b>R311</b> DO	R310	R309	R308	R307 DC	R306	R305 DC	R303, 304 DC	R302 DO	R301	5–12 Ref
DOT3-012-0140	DOT1-038-0212		DOT1-045-5000	DOT1-025-0143	DOT64200005-053	DOT1-025-0206	DOT64025455-002	DOT1-017-5573	DOT1-017-5574			DOT1-011-2635	DOT1-011-2294	DOT1-011-2596	DOT1-011-2893	DOT1-011-2835	DOT69008520-024		DOT69008520-027	DOT1-011-2567		DOT1-011-2628	DOT1-011-5755	DOT1-011-2713	DOT1-011-2745		DOT1-011-2615				DOT1-011-2639		DOT1-011-2353	DOT1-011-2570	DOT1-011-5783		Part Number
Nut, hex, 4-40	Solid Wire, 22 AWG	(Same as W203)	Jumper Wire, 22 AWG	IC, 78MO5CT, voltage regulator, +5V, TO-2	IC, 74128 quad 2-Input NOR	IC, TL494CN, voltage regulator	IC, TDA1170N, deflection circuit	Transformer, DC–DC, HFS–5	Transformer, flyback, HFS, 32kHz	(Same as R107)	(Same as R203)	Resistor, 82Q, 1/4W, 5%	Resistor, 10kΩ, 1/4W, 5%	Resistor, 33Q, 1/4W, 5%	Resistor, 300Q, 2W, 5%	Resistor, 118Q, 1/4W, 1%	Resistor, 169Q, 1/4W, 1%	(Same as R406)	Resistor, 340Q, 1/4W, 1%	Resistor, 150Q, 1/4W, 5%	(Same as R104)	Resistor, 680Q 1/4W, 5%	Potentiometer, 500Q 1/8W, 20%	Resistor, 56Ω 1/4W, 5%	Resistor, 4.7MQ, 1/4W, 5%	(Same as R214)	Resistor, 470kΩ 1/4W, 5%	(Same as R107)	(Same as R214)	(Same as R214)	Resistor, 820kΩ, 1/4W, 5%	(Same as R302)	Resistor, 3MΩ 1/2W, 5%	Resistor, 150kΩ, 1/4W, 5%	Potentiometer, 7.5MQ, 1W, 20%	(Same as R214)	Description
								1			•		٠																								Use Code
2	Α/R		2	<b> </b>	-	-	<u> </u>	<u> </u>	-			2	_		2	-	H		<u> </u>	2		<b>—</b>	<b></b> 4	<b>;—</b> ↓	<b></b>						_		<b>,</b>	2	2		Qty

5-12 Ref	Part Number	Description	Use Code	Qty
3	DOT3-011-1044	Screw, 4-40 x 5/16		2
4	DOT3-015-0125	Spring, DAG, 1/4 OD		-
U)	DOT3-013-0140	Washer, lock #4		2




AH


Figure 5–12 Monitor PCB
(part of Display PCB/CRT Assembly)
Part Number 803706
(Sheet 2 of 4)




All Rights Reserved



	5–12 Ref
DOT6-004-1150	Part Number
DOT6-004-1150 CRT Connector Assembly (part of Display PCB/CRT Assembly)	Description
	Use Code
REF	Qty





ΑH

Figure 5-12 CRT Connector Assembly (part of Display PCB/CRT Assembly, Part Number 803706) (Sheet 4 of 4)

5–13 Ref	Part Number	Description	Use Code	Qty
IW	803783_30	Power Supply PCB/AC Receptacle/Ground Wire Harness		REF
Z	200267-026	Terminal, receptacle, insulated, 22–18 AWG		2
77	202151-056	Connector, housing, 18 AWG, 6 contact		ш
22	202151-117	Cover, housing, strain relief, 6 contact		<del></del>
<b>د</b> يد	202151-500	Plug, keying		<u> </u>
4	200276-209	Terminal, lug, ring tongue, #10, 22-16 AWG		<del></del>
<b>U</b> I	200283-018	Tubing, heat shrink, black, 0.187 ID		Α/R
Ο.	201756-016	Tubing, PVC, clear, 0.166 ID		Α/R
7	201535-084	Wire, semirigid, PVC, 300V, 18 AWG, blue		Α/R
<b>∞</b>	201535-080	Wire, semirigid, PVC, 300V, 18 AWG, brown		A/R
9	201535-078	Wire, semingid, PVC, 300V, 18 AWG, black		Α⁄R





Figure 5–13 Power Supply PCB/AC Receptacle/ Ground Wire Harness


5-14 Ref	Part Number	Description	Use Code	Qty
W2	803783-23	Power Conversion PCB/Power Supply PCB/Battery Wire Harness		REF
J28	202151–335	Connector, cable, 18 AWG, 6 contact		_
P2-A	202151-425	Connector, housing, feed-thru, 6 contact		<u>ш</u>
P13	202151-056	Connector, housing, 18 AWG, 6 contact		-
2	202151-117	Cover, housing, strain relief, 6 contact		2
w	202151-436	Cover, housing, feed-thru, 6 contact		<u> </u>
4	201535-078	Wire, semirigid, PVC, 300V, 18 AWG, black		A/R
υ <sub>η</sub>	201535-079	Wire, semirigid, PVC, 300V, 18 AWG, red		A/R
0	201535-087	Wire, semirigid, PVC, 300V, 18 AWG, white		A/R



₹ N

Figure 5–14 Power Conversion PCB/Power Supply PCB/Battery Wire Harness


5-15 Ref	Part Number	Description	Use Code	Qty
W3	803783-24	Battery Sensor/Power Supply PCB Wire Harness		REF
P2_B	202151-030	Connector, housing, 22 AWG, 6 contact		<u>, , , , , , , , , , , , , , , , , , , </u>
P19	202125-024	Receptacle Assembly, single row, 2 contact	,,,,,	<u>-</u>
2	202151-117	Cover, housing, strain relief, 6 contact		<u> </u>
Ç.	202151-500	Plug, keying		نسو
4	201535-019	Wire, semirigid, PVC, 300V, 22 AWG, red		A/R




₩З

Figure 5–15 Battery Sensor/Power Supply PCB Wire Harness

REF		High Voltage Connector Cable Assembly	803756-13	W4
Qty	Use Code	Description	Part Number	5–16 Ref




REF		Main PCB/Display Cable Assembly	803733-00	W5
Qty	Use Code	Description	Part Number	5-17 Ref



Ş


Figure 5–17 Main PCB/Display Cable Assembly

REF		Main PCB/Recorder Cable Assembly	803773-00	W6
Qty	Use Code	Description	Part Number	5–18 Ref



W6

Figure 5–18 Main PCB/Recorder Cable Assembly







Figure 5-19 Interconnect Pacing PCB/Keypad Cable Assembly

7.0

5–20 Ref	Part Number	Description	Use Code	Qty
W8	803783-00	Main PCB/Patient Connector Wire Harness		REF
<u></u>	801050-00	Connector: circular chassis	, t. s.	<b>-</b>
P22	202125-026	Receptacle, 4 contact		juani.
2	200734-001	Cable, shield, stranded, 24 AWG, 300V, CSA		A/R
ij.	200283-014	Tubing, heat sbrink, black, 0.046 ID		A∕R
	200283-017	Tubing, heat shrink, black, 0.125 ID		A/R
υı	200283-018	Tubing, heat shrink, black, 0.187 ID		A/R
6.	201756-020	Tubing, PVC, clear, 0.263 ID		A/R



₩8

Figure 5-20 Main PCB/Patient Connector Wire Harness

5–21 Ref	Part Number	Description C	Use Code	Qty
W9	803783-04	AC Receptacle/Ground Wire Harness		REF
2	200276-212	Terminal, lug, 22–16 AWG		<del></del>
ພ	200283-018	Tubing, heat-shrink, black, 0.187 ID		A/R
4	201535-090	Wire, semirigid, PVC, 300V, 18 AWG, green/yellow		A/R





Figure 5-21 AC Receptacle/Ground Wire Harness

5–22 Ref	Part Number	Description	Use Code	Qty
W10	803783-05	AC Receptacle Wire Harness		REF
2	200267-026	Terminal, receptacle, insulated, 22–18 AWG		2
w	201535-084	Wire, semirigid, PVC, 300V, 18 AWG, blue		A/R



**%**10

Figure 5-22 AC Receptacle Wire Harness

ي زر	5-23 Ref	Part Number	Description	Use Code	Qy
	WII	803783-06	Test Load Contact/Resistor Wire Harness		REF
	2	202190-000	Cable, high voltage, UL/CSA, 10kV, 0.0312 OD		A/R
	ဃ	200276-211	Terminal, lug, 22-16 AWG		2




Figure 5–23 Test Load Contact/Resistor
Wire Harness

5–24 Ref	Part Number	Description	Use Code	Qty
W12	803783-07	Power Conversion PCB/Reed Assembly Wire Harness		REF
P18	202125-024	Receptacle, 2 contact, 26-22 AWG		<del></del>
2	803060-00	Reed Assembly, magnetic		1





Figure 5–24 Power Conversion PCB/Reed
Assembly Wire Harness


5–25 Ref	Part Number	Description	Use Code	Qty
W13	803783-08	Power Conversion PCB/Dump Relay Wire Harness		REF
P12	202151-026	Connector, housing, 22 AWG, 2 contact		
Relay 1,4	200267-020	Terminal, receptacle, 22–18 AWG		2
22	202151-113	Cover, housing, strain relief, 2 contact		_
w	200536-001	Retainer, cable tie, nylon, 0.1 W x 4 L		Н.
4	201535-021	Wire, semirigid, PVC, 300V, 22 AWG, orange		A/R
5	201535-027	Wire, semirigid, PVC, 300V, 22 AWG, white		A/R



**₩1**3

Figure 5-25 Power Conversion PCB/Dump
Relay Wire Harness

<b>&gt;</b>	-	Trine, schmight, r vc, 500 v, 22 Avv O, gray	020-020	N
		201525 ONE With comitteed DVC 200V 22 AVVC com	201525 026	<b>.</b>
	•	(Same as P11)		Relay 3, 8
w		200267-020   Terminal, receptacle, 22-18 AWG	200267-020	PII
REF		803783-09 Power Conversion PCB/Charge Relay Negative Wire Harness	803783-09	W14
QIV	Use Code	Description	Part	5-26 Ref



**W14** 

Figure 5-26 Power Conversion PCB/ Charge Relay Negative Wire Harness

5–27 Ref	Part Number	Description	Use Code	Qty
\$1W	803783-10	803783-10 Power Conversion PCB/Charge Relay Positive Wire Harness		REF
P10	200267-020	200267-020 Terminal, receptacle, 22-18 AWG		2
Relay 5	1	(Same as P10)		
2	202190-000	202190-000   Cable, high voltage, UL/CSA, 10kV, 0.0312 OD		Α/R




Figure 5-27 Power Conversion PCB/Charge Relay Positive Wire Harness

5–28 Ref	Part Number	Description	Use Code	Qty
W16	803783—11	Main PCB/Sternum Relay Wire Harness		REF
P8 Relay 11	200267-020	Terminal, receptacle, 22–18 AWG (Same as P8)		2
Kelay II	202190-000	(Same as P8) Cable, high voltage, UL/CSA, 10kV, 0.0312 OD		Α/R





Figure 5–28 Main PCB/Sternum Relay Wire Harness

5–29 Ref	Part Number	Description	Use Code	Qty
W17	803783-12	Main PCB/Apex Relay Wire Harness		REF
<b>P</b> 9	200267-026	Terminal, receptacle, 22–18 AWG		<b></b>
Relay 12	Relay 12   200267-020	Terminal, receptacle, 22–18 AWG		A 7₽
4	202190-000	Capie, figh voltage, UL/CSA, TUKV, U.USTZ UD		AVX.



W17

Figure 5-29 Main PCB/Apex Relay Wire Harness

5–30 Ref	Part Number	Description	Use Code	Qty
W18	803783-13	Storage Capacitor Positive/Relay Wire Harness		REF
Relay 10	200267-020	Terminal, receptacle, 22-18 AWG, 0.187		
CAP+	200267-026	Terminal, receptacle, 22-18 AWG, 0.250		_
2	202190-000	Cable, high voltage, UL/CSA, 10kV, 0.0312 OD		Α/R





Figure 5–30 Storage Capacitor Positive/ Relay Wire Harness


Part Number	Description	Code	Q.V
803783-14	Storage Capacitor Negative/Relay Wire Harness		REF
200267-020	Terminal, receptacle, 22–18 AWG, 0.187		<del> </del>
200267-026	Terminal, receptacle, 22–18 AWG, 0.250		Ъ
202190-000	Cable, high voltage, UL/CSA, 10kV, 0.0312 OD		A/R
	mber	Storage Capacitor Neg Terminal, receptacle, 22 Terminal, receptacle, 22 Cable, high voltage, UL	Imber       Description       Use Code         Storage Capacitor Negative/Relay Wire Harness       Code         Terminal, receptacle, 22–18 AWG, 0.187       Terminal, receptacle, 22–18 AWG, 0.250         Cable, high voltage, UL/CSA, 10kV, 0.0312 OD       Code



₩19

Figure 5–31 Storage Capacitor Negative/ Relay Wire Harness


5-32 Ref	Part Number	Description	Use Code	Qty
W20	803783-15	ECG Out Wire Harness		REF
P29	202125-024	Receptacle, 2 contact		-
2	200317-002	Connector, phone jack		
دي	200283018	Tubing, heat shrink, black, 0.187 ID		A/R
4	201535-021	Wire, 300V, 22 AWG, orange	,	Α/R
Uι	201535-027	Wire, 300V, 22 AWG, white		AR



W20

Figure 5–32 ECG Out Wire Harness


5-33 Ref	Part Number	Description	Use Code	Qy
W21	803783-16	Main PCB/Speaker Wire Harness		REF
P35	202125-024	Receptacle, 2 contact		<del></del>
IX	201148021	Speaker, miniature audio	10 <del>1</del> 111111111111	ped.
0	200283-017	Tubing, heat shrink, black, 0.125 ID		A/R
7	201756-014	Tubing, PVC, clear, 0.133 ID		Α/R
<b>∞</b> ;	201535-021	Wire, 300V, AWG 22, orange		A/R
9	201535-027	Wire, 300V, AWG 22, white		Α/R



W21

Figure 5-33 Main PCB/Speaker Wire Harness


5–34 Ref	Part Number	Description Co.	Use Code	Qty
W22	803783-29	Power Supply Bracket/Ground Wire Harness		REF
2	200276-209	Terminal, lug, ring, #10		_
e,	200276-212	Terminal, lug, ring, #22	·	
4	200283-018	Tubing, black, 0.187 ID	·	A/R
5	201535-090	Wire, semirigid, PVC, 300V, AWG 18, green/yellow		A/R



**₩22** 

Figure 5-34 Power Supply Bracket/ Ground Wire Harness

REF		Connector Receptacle, Pacing Cable Assembly	802925-03	W23
Qty	Use Code	Description	Part Number	5-35 Ref







Figure 5-35 Connector Receptacle, Pacing Cable Assembly

5-36 Ref	Part Number	Use Description Code		νįν
W24	803783-19	Recorder/Bracket Ground Wire Harness	<b>₽</b>	REF
ы	200276–209	Terminal, lug, ring, #10	2	
ယ	201756010	Tubing, PVC, clear, size 12	>	A/R
4	201535-023	Wire, semirigid, PVC, 300V, AWG 22, green	>	A/R




W24


Figure 5–36 Recorder/Bracket Ground
Wire Harness

7	S3 W1 W2 W2 W3 W3 W4		5-37 Ref
803752-37 803752-38 803752-40 803752-42 803731-04 803731-05 803731-07 803731-09 803731-10 803731-11 803731-12	803747-11 803774-01 803774-04 803774-05 803776-09 803738-11 803750-00 803750-00 803752-31 803752-33	803747-00 803747-01 803747-02 803747-03 803747-04 803747-05 803747-06 803747-07 803747-09 803747-09	Part Number
Label, Defibrillation Adapter, Spanish Label, Defibrillation Adapter, Portuguese Label, Defibrillation Adapter, Italian Label, Defibrillation Adapter, Italian Label, Serial Number, English, Domestic (not shown, used on item 2) Label, Serial Number, Spanish, International Label, Serial Number, Spanish, Domestic Label, Serial Number, French, CSA Label, Serial Number, English, CSA Label, Serial Number, English, UK Label, Serial Number, German, International	sh Harness (pg. 5–95) Harness (pg. 5–95) ness (pg. 5–96) e Assembly (pg. 5–97) ede English French	Defibrillation Adapter, English, Domestic Defibrillation Adapter, French, International Defibrillation Adapter, German, International Defibrillation Adapter, Spanish, International Defibrillation Adapter, Spanish, Domestic Defibrillation Adapter, French, Domestic Defibrillation Adapter, French, CSA Defibrillation Adapter, English, CSA Defibrillation Adapter, English, UK Defibrillation Adapter, Portuguese Defibrillation Adapter, Italian	Description
С Н С В В В В В В В В В В В В В В В В В	L A,H,I B,F,G	д н с н в с в в	Use Code
	<b>X</b> EF	REF REF REF REF REF REF REF REF	Qty

5-37 Ref	Part Number	Description	Use Code	Qty
	803731–26	Label, Serial Number, Portuguese	J	
i de la constitución de la const	803731-29	Label, Serial Number, Italian	*	
	803731-32	Label, Serial Number, Swedish	Г	
90	805532-01	Fastener Plate, switch, test load		
9	201508-003	Nut, lock, kep, 6-32	,	
10	803780-02	Pin, connector, electrode		
11	803725-02	Plate, test load		
12	802935-08	Pushbutton, Paddle, Discharge, Domestic	А,Е-Н	
	802935-09	Pushbutton, Paddle, Discharge, International	B-D,I-L	
13	200536-001	Retainer, cable tie, nylon, 0.10 W x 4 L (not shown, used on W3)		
14	200530-002	Retainer, clamp, cable		
5	202039-000	Retainer, wire, adhesive mount (not shown, used on W3)		
16	202253-506	Screw, pan head, 2-56 x 0.375 L		
17	201874-001	Screw, pan head, 6-32 x 0.250 L		
18	201874-236	Screw, pan head, 6-32 x 0.750 L		
19	201625–173	Spring, compression, 0.240 OD, 0.562 L		
20	200283-003	Tubing, heatshrink, red 0.125 ID (not shown, used on S3)		A/R
21	201575-004	Washer, lock, #6 x 0.288 OD, 0.018 T		

Figure 5–37 Defibrillation Adapter
Part Number 803747–00
(Sheet 1 of 2)





LIFEPAK 9P defibrillator/monitor/pacemaker




Figure 5-38 Test Load High Voltage Wire Harness

5–39 Ref	Part Number	Description	Use Code	Qty
W2	803774-04	Test Load High Voltage Wire Harness	·	REF
2	202190-000	Cable, high voltage, UL/CSA, 10kV, 0.0312 OD		A/R
3	200276-211	Terminal, lug, 22-16 AWG		2





Figure 5–39 Test Load High Voltage Wire Harness

1 2	Pushbutton Switch Wire Harness  Receptacle Assembly, polarized, latching, 6 contact  Switch, pushbutton, SPST	<b>803774-05</b> 202125-028 804697-04	W3 P5 S1,S2
Use Code Qty	Description	Part Number	5–40 Ref



₩3

Figure 5–40 Pushbutton Switch
Wire Harness

W4	5–41 Ref
803756-09	Part Number
High Voltage Connector Cable Assembly	Description
	Use Code
REF	Qīy





Figure 5-41 High Voltage Connector Cable Assembly

. 

Түре		IC NUMBER	LOCATION	REF. DES.
Multiplexer	Δ	DG509	Main (A1)	U41, 42







Figure 6–14 Multiplexer

Түре			LOCATION	REF. DES.
Multivibrator		74HC221	Main (A1)	U31



	INPUTS			PUTS
CLEAR	A	8	Q	Q
X X H	X H X L + L	X X H H	7	ᅻᅻᅼᆂᆂᆂ

H = HIGH LEVEL

L = LOW LEVEL

+ = TRANSITION FROM LOW TO HIGH


+ = TRANSITION FROM HIGH TO LOW


T = ONE HIGH LEVEL PULSE
T = ONE LOW LEVEL PULSE

X = IRRELEVANT


FUNCTION TABLE

Figure 6–15 Multivibrator





Түре	IC Number	LOCATION	REF. DES.
Oscillator	NCM031C	Main (A1)	<b>U</b> 10



PIN CONFIGURATION

**OUTPUT FREQUENCY (in MHz)** 

BASE	OUTPUT	DIVIDING	RATIOS		1
FREQUENCY		1/2	1/4	1/8	1/16
16.000	16.000	8.000	4.000	2.000	1.000

FREQUENCY TABLE

Figure 6-17 Oscillator

Туре	IC Number	Location	REF. DES.
Pulse-Width Modulator	2843	Power Conversion (A3)	U3
	3843	Interconnect (A5)	U6




Figure 6-18 Pulse-Width Modulator

Түре	IC NUMBER	Location	REF. DES.
Pulse–Width Modulator	3525A	Main (A1)	U5







Figure 6-19 Pulse-Width Modulator


FUNCTIONAL DIAGRAM







Түре	IC NUMBER	Location	REF. DES.
Shift Register 🛕	74HC164	Power Conversion (A3)	U19, 20




PIN CONFIGURATION


RCK	SCK	SCLR	Ğ	Function
· <b>X</b>	X	×	н	QA thru QH = TRI-STATE
×	X	L	L	Shift Register cleared O'H=0
×	1	н	L	Shift Register clocked  Q <sub>N</sub> = Q <sub>n-1</sub> , Q <sub>0</sub> = SER
1	x	н	L	Contents of Shift Register transferred to output latches


FUNCTION TABLE

Figure 6-22 Shift Register

Түре	IC Number	LOCATION	REF. DES.
Switch, Analog 🛕	DG302	Main (A1)	U43







FUNCTIONAL DIAGRAM

LOGIC	SWITCH
0	OFF
1	ON

FUNCTION TABLE

Figure 6-23 Analog Switch

# Index

## Numbers

1mV CAL pushbutton, 2-6

### A

AC Mains Power switch (rear panel), 2–6
Accessories and Supplies, 5–5
Adapters and Accessories, 5–5
Address Decode circuit, 1–33
Analog Voltage Evaluation circuit, 1–17
APEX Paddle description, 2–5
Arming Circuit, 1–21
Assembly (Theory of Operation)
Display PCB/CRT, 1–39
Interconnect Pacing PCB, 1–17
Main PCB, 1–23
Power Conversion PCB, 1–10
Power Supply PCB, 1–9

## B

BATT CHRG check (PIP), 3-5
BATT CHRG Indicator, 2-4
BATT CHRG Indicator circuit, 1-10
Battery
Battery Pack Removal, 4-7
Low Battery Threshold Checks, 3-27
Maintenance, 4-1
Battery Capacity Check, 4-3
Guidelines, 4-2
Replacement Batteries, 4-2
Recharging, 4-1
Recycling, 4-4
Theory of Operation, 1-9
Battery Power check (PIP), 3-5

### C

Cable Accessories, 5-5 Calibration. See Test and Calibration Procedure (TCP) Capacitor Charge Control circuit, 1-14 Capacity, Battery Check, 4-3 Charge Control circuit, 1-31 Charge LED Drive circuit, 1–17 Charge LED Indicator, 2-5 CHARGE pushbutton, 2-3, 2-5 Charging Batteries, 4–1 Checklist PIP, 3-21 TCP. 3-45 Circuit descriptions. See Theory of Operation Cleaning Procedures, 4–14 External, 4-14 Interior, 4-15 Tools and Materials for, 4-14 Clock Setting Procedure, 2-10 **CODE SUMMARY** PIP-CODE SUMMARY Report, 3-14 Pushbutton description, 2-3 Configuration Information, xviii Controls and Indicators, 2-1 Front Panel, 2-3 Rear Panel, 2-6 CRT Bias Generator circuit, 1–42 **CRT Connector Assy (Part of Display PCB/CRT** Assy), 5-72 CRT Disassembly, 4-12 CRT display, 2-3 CRT Frequency Response test (TCP), 3-42 **CURRENT pushbutton, 2–5** 

	ENERGY SELECT BUS buston 0.0
DO Dantago alimante de co	ENERGY SELECT pushbutton, 2–3
DC Restore circuit, 1–28	Energy Selection Display, 2–4
Defibrillation Adapter	Energy Storage Capacitor Charger circuit, 1–12
Final Assembly, 5–94	Energy Transfer circuit, 1–14
List of Assemblies, 5–4	EPROM
Parts List, 5–92	Font, 1–35
Schematic, 5–93	Program, 1–36
Defibrillator Calibration, 3–30	Error Codes, 3–46
Table of measurement ranges (TCP), 3–31	External Control Input circuit, 1–12
Description	• ,
Functional, 1–7	<b>T</b>
Physical, 1–1	${f F}$
Disassembly Procedures, 4–7	
Battery Pack Removal, 4-7	FAST-PATCH electrodes, ordering, 5-5
Case Separation, 4–10	Fault Status Codes (Error Codes), 3–46
Display PCB/CRT Removal, 4–12	Filtering circuit, 1–28
Front Panel Bezel Removal, 4-7	Final Assembly
Interconnect Pacing PCB Removal, 4–12	Defibrillation Adapter, 5–94
Keypad Asssembly Removal, 4–7	Defibrillator/Monitor/Pacemaker
Main PCB Removal, 4–12	Drawing, 5–11
Power Conversion PCB Reassembly, 4–11	Parts List, 5–6
Power Conversion PCB Removal, 4–10	Font EPROM circuit, 1–35
Power Supply Removal, 4–12	Frame Sweep Generator circuit, 1-41
Strip Chart Recorder Removal, 4–8	Front and Rear Panel Controls Interface circuit,
	1–31
Discharge pushbuttons, 2–5	Front Panel, Controls and Indicators description,
Display Calibration, 3–33	2–3
Display Gate Array circuit, 1–33	Front-End Protection circuit, 1-25
Display PCB/CRT Assembly	Functional Block Diagram, 1–8
Disassembly, 4–12	Functional Description, 1–7
Monitor PCB (part of Display PCB/CRT), Parts List,	Functional Test and Calibration. See Test and
5-65	Calibration Procedure (TCP)
Parts List, 5–63	Fuse Replacement, 4–5
Theory of Operation, 1–39	Fuse Cover (reer penel) 0.6
Display/Recorder/Pacing Controller circuit, 1–32	Fuse Cover (rear panel), 2-6
_	$\mathbf{G}$
${f E}$	
	Ground connection (rear panel), 2-6
ECG Electrodes, Accessories, 5–5	
ECG Gain indication, 2–4	
ECG Out circuit, 1–39	$\mathbf{H}$
ECG Output	
ECG Output (TCP test), 3-42	Heart Rate Alarm Limits, 2–3
ECG OUTPUT rear panel connection, 2–6	Heart Rate display, 2-3
ECG Output specifications, 1–3	HR ALARM pushbutton, 2–3
ECG SIZE pushbutton, 2–3	• <b>,</b>
Echo circuit, 1–29	
	I
Echo Pulse Width Converter circuit, 1–20	

Effective Publication Dates, xvii

Energy Dump circuit, 1–15

I/O Protection circuit, 1–31

Info Menu, 2–9

Inspection Techniques, 4–13	Address Decode, 1–33		
Exterior Inspection, 4–13	Display Gate Array, 1-33		
Interior Inspection, 4–13	Display/Recorder/Pacing Controller, 1–32		
Instruction Decoding circuit, 1–26	ECG Out, 1-39		
Instrumentation Amplifier circuit, 1–27	Font EPROM, 1–35		
Interconnect Diagram, 5–12	Leads Off Detection, 1–39		
	Motor Speed Control, 1–36		
Interconnect Pacing PCB Assembly	Pacing Current Delivery, 1–39		
Disassembly, 4–12	Paper Out Sense, 1–38		
Parts List, 5–44	Printhead Control, 1–37		
PCB Layout, 5–50	Printhead Resistance Evaluation, 1–38		
Schematic, 5–52	Program EPROM, 1–36		
Theory of Operation, 1–17	Program RAM, 1–36 Recorder Power, 1–36		
Isolated Power Supply circuit, 1–22, 1–25	Speaker Drive, 1–39		
	System Clock, 1–32		
	Video RAM, 1–35		
J	Parts List, 5–13		
•	PCB Layout, 5–22		
JOULES SELECTED Indicator, 2–4	Preamp, 1–25		
,	• *		
	DC Restore, 1–28 Echo, 1–29		
$\mathbf{K}$	Filtering, 1–28		
	Front-End Protection, 1–25		
Keypad Assembly	Instruction Decoding, 1–26		
••	Instrumentation Amplifier, 1–27		
Parts List, 5–55	Isolated Power Supply, 1–25		
PCB Layout, 5–55	Lead Select, 1–25		
Schematic, 5–55	Optical Interface, 1–25		
	Schematic, 5–23		
<u>ब</u>	System Controller, 1–30		
${f L}$	Charge Control, 1–31		
	Front and Rear Panel Controls Interface, 1–31		
Lead Select circuit, 1–25	I/O Protection, 1–31		
Lead Select Indicator, 2–4	Pacing Current Adjustment, 1-31		
LEAD SELECT pushbutton, 2–3	Theory of Operation, 1–23		
Leads Off Detection circuit, 1–39	Mains Power switch (rear panel), 2-6		
Leakage Current check (PIP), 3–17	Maintenance		
Pacemaker Leads Detect Current, 3–18	Battery, 4–1		
LED Drive circuit, 1–22	General		
List of Defibrillation Adapter Assemblies, 5–4	Fuse Replacement, 4–5		
List of Major Assemblies, 5–3	Strip Chart Recorder, 4–4		
· · · · · · · · · · · · · · · · · · ·	Cleaning, 4–5		
Logic Power Monitor circuit, 1–11	Loading Paper, 4–5		
Low Battery Threshold Checks (TCP), 3–27	Protecting Printouts, 4–5		
Low Energy Charge Rate circuit, 1–13	Major Assemblies, 5–3		
	· · · · · · · · · · · · · · · · · · ·		
	Message Display, 2–4		
$\mathbf N$	MODE pushbutton. See Setting the Clock		
	Monitor PCB (Part of Display PCB/CRT Assy)		
Main PCB	Parts List, 5–65		
Block Diagram, 1–24	PCB Layout, 5-70		
Disassembly, 4–12	Schematic, 5–71		
Display/Recorder/Pacing Interface, 1-32	Motor Speed Control circuit, 1–36		
· · ·	•		

### N Charge Time, 3–13 Checklist, 3-21 CODE SUMMARY Report, 3-14 NON-DEMAND PACE pushbutton, 2-6 CRT Gain/Aspect Ratio, 3-10 Defib Disarm and CHARGE LEDs Check, 3–12 Definitions, 3-2 Delivered Energy, 3-14 ON pushbutton, 2-3 Demand Pacing, 3–16 Optical Interface circuit, 1-25 Discharge Control, 3-13 ECG Display, 3-9 Option Settings, 2–7 Energy Dump, 3-13 Info Menu, 2–9 Setup Menu, 2-8 General Instructions, 3–4 Heart Rate, 3-10 Test Menu, 2–7 Keypad Test, 3-7 Ordering Replacement Parts, 5-2 Lead Select, 3–8 Output Isolation circuit, 1-9 Leakage Current, 3–17 Overvoltage Protection circuit, 1-9 Pacemaker Leads Detect Current, 3-18 Non-Demand Pacing, 3-17 Pacemaker Power-Up, 3-16 P Pacing Leads Off, 3-17 Pacing Rate, 3-16 Pace Pulse Amplitude Decoder circuit, 1–19 Paddle Polarity, 3-12 Pace Pulse Length Generator circuit, 1-22 Physical Inspection, 3-4 Pace Pulse Output Drive circuit, 1–19 Purpose, 3–1 Pace Pulse Wave Shaper circuit, 1–19 Rate Alarms/Alarm-Activated Record, 3-10 Pacemaker Calibration (TCP), 3-42 Recorder Annotation/Real-Time Clock, 3-11 PACER pushbutton, 2-4 Recorder Speed/Baseline Centering, 3–8 Pacing Cable Connector, 2–5 Recorder/CRT Noise, 3-12 Pacing Current Adjustment circuit, 1–31 Required Test Equipment, 3–2 Pacing Current Delivery circuit, 1-39 Requirements, 3-2 Pacing Current Display, 2-5 Scope and Applicability, 3-2 Pacing Leads Off Detector circuit, 1-21 Setup Menu Settings, 3–6 Pacing Message Display, 2-5 Sync/Sync Annotation, 3-11 Test Equipment Required, 3-2 Pacing Pulse Display, 2-5 Test Load, 3-11 Packing and Storage or Shipping, 4–17 Test Setup, 3-6 Paddle Assembly, 5-61 Training Requirements for Personnel, 3-3 Parts List, 5-61 Workstation Requirements, 3-3 Schematic, 5–61 Power Control circuit, 1–10 Paper Out Sense circuit, 1–38 Power Conversion PCB Assembly, 1-10 Parts Lists, Explanation of, 5–1 Disassembly, 4–10 Parts, How to order, 5-2 Parts List, 5-33 PARTSLINE telephone number, 5-2 PCB Layout, 5-40 Patient Cable Connector, 2–5 Reassembly, 4-11 PCB Repair, Precautions, 4–16 Schematic, 5-41 Performance Inspection Procedure (PIP). See PIP Theory of Operation, 1–10 Physical Description, 1–1 Voltage Checks (TCP test), 3–26 PIP (Performance Inspection Procedure), 3–1 Power Supply Controller circuit, 1-22 1mV CAL, 3-9 Power Supply PCB Assembly, 1-9

AC Power, 3–5

Battery Power, 3–5

Charge Reduction, 3-14

Battery Charge (BATT CHRG) check, 3-5

Calibration (TCP), 3–37

Disassembly, 4–12

PCB Layout, 5–31

Parts List, 5-28

Schematic, 5-32 Defibrillation Adapter, 5–93 Theory of Operation, 1–9 Display PCB/CRT (Monitor PCB), 5-71 Explanation of, 5-2 Power-On circuit, 1-9 Preamp, Pulse Widths, 1-27 Interconnect Pacing PCB, 5-52 Keypad Assembly, 5-55 Preamp circuit, 1-25 Main PCB, 5-23 Printer. See Strip Chart Recorder Paddle Assembly, 5-61 Printhead Control circuit, 1-37 Power Conversion PCB, 5-41 Printhead Resistance Evaluation circuit, 1-38 Power Supply PCB, 5–32 Program EPROM circuit, 1–36 Rear Panel Membrane Switch, 5-56 Program RAM circuit, 1–36 Sense Leads PCB (mounted on Power Conversion Publication Dates, Effective, xvii PCB), 5-41 Strip Chart Recorder, 5-59 Transfer Relay Assembly, 5-60 Sense Leads PCB (on Power Conversion PCB) PCB Layout, 5-40 QRS Indicator, 2-3 Schematic, 5-41 QRS VOL pushbutton, 2-6 Serial Data Interface circuit, 1-16 Service and Maintenance information, 4-1 R **SERVICE Indicator and Error Codes** Error Codes, 3-46 Raster Generator circuit, 1-42 How to clear, 3-46 Raster Scan Oscillator circuit, 1-41 Status Message Display, 2-4 Raster Voltage Generator circuit, 1-41 Service information overview, xvi RATE pushbutton, 2-5 SET pushbutton. See Setting the Clock **Rear Panel Membrane Switch** Setting the Clock, 2-10 PCB Layout, 5-56 Setup Menu Schematic, 5–56 Access, 2-8 Recharging Batteries, 4-1 Settings for PIP, 3–6 RECORD pushbutton, 2-3 Speaker Drive/Tone Control circuit, 1-39 Specifications, 1-3 Recorder See also Strip Chart Recorder **SSD (Static Sensitive Devices)** Brief description, 2-4 Proper Storage, 4-6 Recorder Calibration, 3-38 Special Handling of, 4-6 Static-Dissipative Mat, 4-6 Recorder Frequency Response test (TCP), 3-40 Static-Free Work Area, 4-6 Recorder Power circuit, 1–36 Symbol, 4–6 Recycling batteries, 4-4 Work Area, 4-6 Reference Designator Key, 5-4 Wrist Strap, 4-6 Regulation circuit, 1-9 Static Sensitive Devices. See SSD Removal and Replacement of Assemblies. See Status Message Display, 2-4 **Disassembly Procedures** STERNUM Paddle description, 2-5 Repair and Cleaning, Tools and Materials for, 4-14 Storage or Shipping, 4–17 Replacement Parts. How to order, 5-2 Strip Chart Recorder Assembly Drawing, 5-58 Cleaning, 4-5

### Safety

General Information, xiii Symbols, xv Terms, xiii **Schematic Diagram** 

Disassembly Procedure, 4-8

Loading Paper, 4-5

Protecting Printouts, 4-5

Recorder Calibration (TCP), 3-38

Maintenance, 4-4

Parts List, 5-57

Recorder Frequency Response test (TCP), 3–40	PIP (Performance Inspection Procedure), 3–1		
Schematic, 5–59	TCP (Test and Calibration Procedure), 3-23		
Supplies and Accessories, 5–5	Theory of Operation, 1–9		
Switching Power Supply circuit, 1–12	Display PCB/CRT Assembly, 1-39		
Switching Regulator circuit, 1–22	CRT Bias Generator, 1-42		
	Frame Sweep Generator, 1-41		
Symbols  Find a patient of the	Raster Generator, 1–42		
Explanation of, xv	Raster Scan Oscillator, 1-41		
Used in manual or on equipment, xv	Raster Voltage Generator, 1-41		
SYNC pushbutton, 2–3	Video Signal Amplifier, 1–41		
System Clock circuit, 1–32	Interconnect Pacing PCB, 1–17		
System Controller circuit, 1–30	Arming Circuit, 1-21		
<b>-,</b>	Block Diagram, 1–18		
	Echo Pulse Width Converter, 1–20		
T	Interface, 1–17		
<b>1</b>	Isolated Power Supply, 1–22		
TCP. See Test and Calibration Procedure	LED Drive, 1–22		
	Pace Pulse Amplitude Decoder, 1–19		
Terms, Definition of, xiii	Pace Pulse Length Generator, 1–22		
Test and Calibration Procedure (TCP), 3–23	Pace Pulse Output Drive, 1–19		
Charge Time, 3-32	Pace Pulse Wave Shaper, 1–19		
Checklist, 3–45	Pacing Leads Off Detector, 1–21		
CRT Frequency Response, 3–42	Power Supply Controller, 1–22		
Defibrillator Calibration, 3–30	Switching Regulator, 1–22		
Defibrillator Charge Time, 3–32	Main PCB, 1–23		
Definitions, 3–23	Display/Recorder/Pacing Interface, 1–32		
Display Calibration, 3–33	Address Decode, 1–33		
ECG Output, 3–42	Display Gate Array, 1–33 Display/Recorder/Pacing Controller, 1–32		
Low Battery Threshold Checks, 3–27	ECG Out, 1–39		
	Font EPROM, 1–35		
Output Waveform, 3–32	Leads Off Detection, 1–39		
Pacemaker Calibration, 3–42	Motor Speed Control, 1–36		
Pacemaker Current and Waveform Output, 3–43	Pacing Current Delivery, 1–39		
Personnel training requirements, 3–25	Paper Out Sense, 1–38		
Power Supply Calibration, 3–37	Printhead Control, 1–37		
Purpose, 3–23	Printhead Resistance Evaluation, 1–38		
Recorder Calibration, 3–38	Program EPROM, 1–36		
Recorder Frequency Response, 3–40	Program RAM, 1–36		
Requirements, 3–23	Recorder Power, 1–36		
Scope, 3–23	Speaker Drive, 1–39		
Test Equipment Required, 3-24	System Clock, 1–32		
Test Load Calibration, 3-31	Video RAM, 1–35		
Test Menu Access, 3–29	Preamp, 1–25		
Voltage Checks, 3–26	DC Restore, 1–28		
Workstation Requirements, 3–24	Echo, 1–29		
	Filtering, 1–28		
Test Equipment Required	Front–End Protection, 1–25		
PIP, 3–2	Instruction Decoding, 1–26		
Test and Calibration, 3–24	Instrumentation Amplifier, 1–27		
Test Load, 2-5	Isolated Power Supply, 1–25		
Test Load Calibration (TCP), 3–31	Lead Select, 1–25		
Test Menu	Optical Interface, 1–25 System Controller, 1–30		
Access, 2–7	Charge Control, 1–30		
TCP-Test Menu Access, 3-29	Front and Rear Panel Controls Interface, 1–31		
Testing	I/O Protection, 1–31		
icaung	BO F follotton, 150 t		

Pacing Current Adjustment, 1–31 Power Conversion PCB, 1–10 Analog Voltage Evaluation, 1-17 Capacitor Charge Control, 1-14 Charge LED Drive, 1-17 Energy Dump, 1–15 Energy Monitor, 1–16 Energy Storage Capacitor Charger, 1–12 Energy Transfer, 1-14 External Control Input, 1–12 Logic Power Monitor, 1–11 Low Energy Charge Rate, 1-13 Power Control, 1-10 Serial Data Interface, 1-16 Switching Power Supply, 1–12 Power Supply PCB, 1-9 BATT CHRG Indicator, 1-10 Battery, 1-9 Output Isolation, 1-9 Overvoltage Protection, 1-9

Power-On, 1-9

Regulation, 1–9
TIME/DAY pushbutton. See Setting the Clock
Tools and Materials for Cleaning, 4–14
Tools and Materials for Repair, 4–14
Transfer Relay Assembly, 5–60
Schematic, 5–60
Troubleshooting Aids, 3–46

### ${f V}$

Video RAM circuit, 1–35 Video Signal Amplifier circuit, 1–41

### W

Warnings, xiv General, xiv

		· · · · · · · · · · · · · · · · · · ·	
	•		