

A-dec 500 Adjustments



A-dec 500 Leveling

■ Leveling should be checked and corrected before tension adjustments are made

- Rough leveling is done during assembly
- Final leveling should be done after all equipment is installed on the chair and the systems are in their normal use positions

■ Level from the ground up:

- Link arm and Support arm
- Front mount
- Control head
- Tray


■ A leveling video is available at a-dec.biz



Front Mount Leveling

■ Front Mount Final or Re-Leveling

- A bubble level is fixed to the front mount to check for level
- You will need two $\frac{3}{4}$ " open end wrenches, a $\frac{7}{16}$ " wrench and a $\frac{5}{16}$ " hex wrench
 - Loosen cam bolts, leveling bolt flanged nuts and stabilization screws
 - Rotate cam for left – right level
 - Move leveling studs in or out equally for front to back leveling
- Secure all fasteners to ensure leveling is locked in position



Front Mount Rotation Stops

- **Front Mount Rotation Stops**
- **Stops are positioned in holes to limit arm rotation**
- **If support center is installed, place a stop in forward position on support center side**
- **If converting left to right, the stops can be repositioned**



Front Mount Rotation Tension

■ Front Mount Tension Adjust

- Use 15/16" wrench to adjust rotational tension
- Nut on under side is self locking and captured



Delivery Arm Rotation Tension

■ Delivery System Rotational Tension Adjustment

- Should be adjusted at installation
- Use 3/32" hex wrench to set tension
- Note that this is one adjustment that is not factory set. The clamp must be relaxed to allow the arm to fit over the post.




Delivery Flex Arm Adjustment

■ Delivery System Flex Arm Counter Balance Adjustment

- All instruments and accessories should be installed prior to making this adjustment
- Turn Master Toggle OFF to release air brake
- Remove the 3 covers from rigid arm end of flex arm using a $5/64$ " hex wrench
- Lower control head and slide flex arm cover off
- Use a $7/16$ " open end wrench to adjust spring tension

■ **Note: Apply Lubriplate to the cover grooves to make installation easier**



Delivery System Leveling

■ Delivery System

- Front to Back Leveling Adjustment
 - Remove joint covers with a 5/64" hex wrench
 - Use 3/16" hex wrench to level control head
 - Turn wrench clockwise to raise front of control head
- Position the delivery system centered over the chair when leveling



Control Head Leveling

- **Control Head Leveling Left to Right**
 - Use 5/32" hex wrench
 - Loosen and tighten the two leveling screws on opposite sides of the round support
- **Position the delivery system centered over the chair when leveling**



Control Head Rotation Tension

■ Control Head Rotation Tension Adjustment

- Use 5/32" hex wrench
- Turning wrench clockwise increases brake tension



Tray Holder Leveling

■ Tray Holder Leveling Adjustments

- Use a 5/32" hex wrench
- Left to right
- Front to back



Tray Holder Rotation Tension

■ Tray Holder Adjustments

- Swivel tension in two places
 - Remove cover for access to rotation adjustment for upper arm
- Use 9/16" socket



Support Link Leveling

■ Leveling of Rear Support Link

- Loosen all four $\frac{3}{4}$ " cap screws
- Loosen leveling bar $\frac{3}{4}$ " jamb nut
- Use a $\frac{3}{4}$ " end wrench to turn leveling bar bolt
- Tighten jam nut while holding leveling bolt
- Tighten the four cap screws to secure Support Link

■ You will need to remove the support arm limit switch to gain access.



Lower Support Arm Tension

■ Rotation tension

- Use 3/16" hex wrench to eliminate any free play in the bearing and to increase or decrease rotational tension

■ This adjustment should be done when lower arm is mounted to link arm and prior to mounting modules to the arm.



Assistants Arm Tilt Tension

- **Assistants Arm**
- **Tilt up tension adjustment**
 - Note: this tilt up resistance is factory set and should not need further adjustment
 - Use 3/16" hex wrench to change tilt up resistance



Assistants Arm Rotation Tension

■ Assistants Arm Rotation

- Remove cap from top of joint
- Use 5/32" hex wrench to adjust tension

■ Note position of data cable and the power cord



Assistants Arm Leveling

- **Assistants Arm Leveling**
- **This adjustment is at the top joint**
 - Remove joint cover
 - Use 5/32" hex wrench to loosen the two screws
 - Level top of Assistant's package
 - Tighten the screws to lock the new position
 - Move arm back and forth to observe consistent level



Monitor Mount Tension

- **Monitor Mount Adjustments**
- **Pull tension and push resistance can be adjusted separately**
 - Remove cover from top of mount
 - Use a 5/32" hex wrench
 - Left screw adjusts push resistance
 - Right screw adjusts pull tension



Dental Light Flex Arm Balance

■ Dental Light Adjustments

- To adjust **flex arm** tension
 - Remove cover by removing the Phillips screw and the two 5/64" hex button head screws
 - Using a 1/2" open end wrench adjust tension on the spring so the light head stays in position



Dental Light Adjustments

■ Dental Light Adjustments

- Use 5/32" hex wrench to adjust the tension for:
 - Horizontal
 - Diagonal



Dental Light Adjustments

■ Dental Light Adjustments

- For **vertical** adjustment:
 - Use 5/64" hex wrench to loosen set screw
 - Use flat tip screw driver to adjust tension screw
 - Note: Access to screw is under right side yoke plug
- To **focus** use flat tip screw driver to obtain uniform pattern
 - Focus range is 18 to 31 inches
 - Factory set at 27 inches




Swivel Brake Adjustment

■ Swivel Brake Adjustment, p. 33

- Release brake tension
- Insert 7/64" hex wrench in through-hole and into adjustment screw
- Turn wrench to right to increase brake friction

■ Hinged handle available as of Jan. '06

- P/N TBA



Headrest Adjustment

■ Headrest lever adjustment

- Amount of play in lever actuation can be adjusted
- Use 1/16" hex wrench to control amount of play in lever actuation

■ Headrest clutch assembly is not serviceable.



Glide bar Adjustment

■ Glide bar adjustment

- Use 1/8" hex wrench
- Turn clockwise to increase friction on glide bar



Air Pressure Adjustment

■ Adjust Air Regulator to 80 PSI


- Rotate the pre regulator knob to change system air pressure
 - Clockwise to increase air pressure
 - Counter clockwise to decrease air pressure, note: you will need to release some air from the system to get an accurate reading



Water Pressure

■ Water regulation

- Water bottle has built in regulator set to maintain 40 PSI. Not adjustable.
- City water regulator has built in pressure reducer to maintain 35 to 40 PSI



Cuspidor Water Adjustment

- **Cuspidor Water Intensity Adjustments**
- **Correct intensity adjustment will prevent splash and splatter**
 - Cup fill intensity is the lower of the two top flow controls.
 - Note: the uppermost is reserved for future use
 - Bowl rinse intensity adjust is on bottom
 - Use 1/8" hex wrench to adjust desired water flow



Cuspidor Water Flow Timing

- **Cuspidor Water Cup Level and Bowl Timing Adjustments**
- **Cup fill level is factory set to 2.5 seconds**
- **Bowl rinse timing is factory set to 15 seconds**
- **To change cup level or bowl timing:**
 - Press the program button (1 beep)
 - Press cup or bowl button for as long as flow is desired then release (3 beeps confirms setting)
- **To stop a timed flow, press the button once**



Control Block Adjustments

■ Control block adjustments for handpiece operation

- Begin by closing all control valves: drive air, air coolant and water coolant
- Water coolant adjustment
 - Lift handpiece from holder
 - Select water coolant and press Foot Control
 - Open water coolant flow control until one drop of water per two seconds is observed



Control Block Adjustments


- **Control block adjustments for handpiece operation**
 - Air coolant adjustment
 - Open air coolant flow control until a fine mist is observed



Drive Air Readout


- **Digital readout of Drive Air Pressure on the control PCB**

- Reference only and not to be used to set pressure.
- It will read higher than actual depending on HP, tubing length/material.
- If needed, it can be used to verify HP or control block operation.




Drive Air Pressure Gauge

- **Drive air adjustment using a handpiece pressure gauge**
 - Open control head cover
 - Lift handpiece from holder
 - Press on Foot Control
 - Rotate Drive Air flow control on top of block associated with the selected handpiece position
 - Observe true Drive Air pressure on gauge



Syringe Flow Adjustment

- **Syringe air and water flow adjustment**
- **Delivery System syringe flow control adjustment inside control head**
 - Lift up control head cover for access
- **Assistant's syringe flow control adjustment in lift arm**
 - Remove lift arm cover for access



A-dec 500 Adjustments



A-dec 500 Maintenance



Hydraulic Fluid

■ Chair hydraulic fluid

- Position chair base and back full up
- Remove both the lift arm cover and the safety plate cover
- View fluid level from the rear
- Top up if need, cycle chair, inspect system for leaks

■ Use only A-dec Hydraulic fluid P/N 61.0197.00



Air and Water Filters

■ Air regulator

- Check/replace air filter

■ Water regulator

- Check/replace water filter

■ Notes:

- Turn off supply valves before attempting to remove filter caps.
- Install new filter with stepped edge down against manifold.



Handpiece Oil Collector


■ Oil collector:

- Cap hinges down for access
- Cap pulls off on models after July 2005

■ Gauze pad should be change weekly

■ Note:

- Muffler can be replace if missing or dirty with P/N 77.0503.00



Vacuum Solids Collector

- **Solids Collector pulls off for access to collection screen**
 - Collection basket should be changed or cleaned daily
- **Hot water should be drawn through the vacuum instruments daily**
- **Screens are available for vacuum instruments to prevent large objects from entering the tubing**



Waterline Maintenance

■ ICX Waterline maintenance

- One tablet in empty water bottle
- Fill with treatment water
- Wait one minute
- Connect to unit
- Repeat for each bottle refill

■ Type of Water

- Tap
- Distilled if tap water contains high mineral content

■ Shock treatment

- If test show high bacteria level
- Sterilex Ultra startup
- ICX maintenance



A-dec 500 Adjustments

Questions?

