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




GB To maintain thorough hygiene in your equipment, EURODENT recommends you to replace the parts in drawing hereunder at the timings shown aside.

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The unit can withstand the following weather conditions when packed for shipping and storage:

- Ambient temperature from -40°C to +70°C
- Relative humidity from 10% to 90%; it does not stand condensation.

The unit does not generate electromagnetic- or other forms of interference into other equipments nearby, nor is it affected by electromagnetic interference.

Burrs, as well as suction tips not provided by EURODENT, must be made of bio-compatible materials (for example, compliant with ISO 10993 standards).

Comply with all applicable local laws when disposing of the unit and its parts. Suction circuit components in particular must be considered as medical waste products and therefore they should be treated accordingly.

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OPERATOR'S HANDBOOK

MODES

- 1** - Water / air for spray presetting ● Water selector/Physiological solution for instrument
 - 2** - Water / air for automatic c.b. presetting ● Chair program 1
 - 3** - Water / air for c.b. presetting ● Chair program 2
 - 4** - Power adjustment ● Chair program 3
 - 5** - Cold light adjustment ● Chair program 4
 - 6** - Operating mode ● Chair last position memory
 - 7** - Chair up
 - 8** - Chair down
 - 9** - Forward tilting
 - 10** - Backward tilting
 - 11** - Headrest up (for ISORAMA chair only)
 - 12** - Headrest down (for ISORAMA chair only)
 - 13** - Headrest forward (for ISORAMA chair only)
 - 14** - Headrest backward (for ISORAMA chair only)
 - 15** - Service push-button / nurse bell
 - 16** - Operating lamp preset
 - 17** - Tumbler filler
 - 18** - Bowl rinsing
 - 19** - Alarm mode / timer
 - 20** - Memory ENTER (data storage confirmation)
 - 21** - Displayed data decrease ● Air / water selection
 - 22** - Displayed data increase ● Air / water selection
 - 23** - Display
-
- 24** - Bowl rinsing
 - 25** - Operating lamp preset
 - 26** - Tumbler filler
 - 27** - Chair up
 - 28** - Chair down
 - 29** - Forward tilting
 - 30** - Backward tilting
 - 31** - Program 1
 - 32** - Chair last position memory
-
- 33** - Tumbler filler
 - 34** - Bowl rinsing

CONTENTS

	Page
1 FOREWORD	34
2 TECHNICAL FEATURES	35
3 SWITCHING ON	35
4 POSITIONING THE MODULES ARM	36
5 CLOCK ADJUSTMENT	36
6 PATIENT POSITIONING	37
Driven movements	37
Programmable and non-programmable automatic movements	37
7 ALARM MODE/TIMER	38
8 TUMBLER FILLER	38
9 BOWL RINSING	39
10 LAMP OPERATION	39
11 SERVICE PUSH-BUTTON	39
12 FOOT-CONTROL SWITCHES	40
13 PROGRAMMABLE MODES FOR THE INSTRUMENTS	41
14 MODULES	43
How to remove modules	43
Handpieces pressure adjustments	43
Syringe module and fibre optic syringe module	44
Turbine module and variable speed turbine module	44
Micromotor module	45
Turbine module and variable speed turbine module with fibre optic	45
Fibre optic micromotor module	45
High torque micromotor module with fibre optic	46
Micromotor module for implantology	46
Powder cleaner module	46
Scaler module	47
Scaler module with fibre optic	47
Telecamera module	48
Acta light cure module	48
Mini Led light cure module	49
Electrosurgery module	50
15 STERILIZATION METHODS AND WORKING TIMES	50
16 MESSAGES EXPLANATION	51
17 PARAMETERS PRESETTING	52
18 MAINTENANCE	53
Module-holder arm balancing	53
Water drain from air filter	54
Cleaning of suction filter	54
Filling of medicament tank	55
Filling of physiological solution	55
Disinfectant solution	56
Replacement of sterile solution bottle	56
Emptying the waste collector	57
Surface cleaning	57
19 METASYS AMALGAM SEPARATOR	58
Removing and cleaning the filter	59
Removing the amalgam separator vessel	59
Ecological system of recovery	60

Dear Doctor,

the EURODENT GROUP is pleased with your choice of ISOPLUS and proud of the favour you grant us.

The ISOPLUS unit has been designed according to three main principles:

ASEPSIS CONTROL, OPERATIONAL FLEXIBILITY AND COMFORT.

Surface smoothness, combined with roundish contours and lack of slits and corners, allows an easy external cleaning and prevents from dusts and sprays deposit. Controls are made so that contact with fingers is avoided whenever possible.

All handpieces are fitted with anti-retraction valving, to prevent risks of cross-contamination.

Lightness of movement and breadth of operational displacement of the modules arm, easiness of access to the patient for both Doctor and Assistant and, last but not least, the flexible microprocessor control of the main activities of the apparatus, all make of ISOPLUS a ductile system ready to match and support the Operators' work attitudes.

The soft profile of the unit and handpieces asset, the anatomical shaping of the chair, the reassuring design of control decks, everything cooperates to frame the comfortable sensation induced into the patient, who will more confidently relax himself so making easier the Operators' activity.

We are certain that your choice will prove widely fulfilling both because of the equipment functionality and reliability and thanks to the benefits in terms of prestige as to your Patients and Colleagues.

This booklet will help you to obtain the best from your ISOPLUS. While reminding you that we are at your disposal for any information that you might require, we remain

Sincerely yours,

2**TECHNICAL FEATURES**

Trade-Mark:	EURODENT
Models:	ISOPLUS (Type 22)
	ISOPLUS LITE
Voltage:	220 V (~ alternating)
Frequency:	50 Hz
Power:	1350 VA (550 VA only unit)
Classification:	Class I with attached parts of type B 
Classification headings (EN60 601-1):	
5.3	Common equipments (equipments with no waterproof casing);
5.5	Equipments which cannot be used in presence of an anaesthetic mixture inflammable by air or by oxygen or by nitrous oxide;
5.6	Equipments for continuous working with intermittent load

Electronic bistoury or electrobistoury Attached part of type BF 

0051 The equipment complies with the fundamental requirements of EC Directive 93/42.

Please refer to the Spare Parts catalogue, and apply to EURODENT for further information about what is not included in it.

EURODENT accepts responsibility for the equipment safety and reliability if and only if:

- installation, assembly, extension, adjustment and repair operations have been carried out by authorized people;
- the general and electric installation of the room complies with the equipment pertinent requirements;
- the equipment has been installed in compliance with the "INSTALLATION INSTRUCTIONS" of the equipment;
- the equipment is used in compliance with the "OPERATOR'S HANDBOOK".


EURODENT hereby declares that the equipment is a dental complex intended as an operative tool for therapeutic use on people by personnel highly skilled in dental practices.

Protection earthing

Caution, see enclosed documentation.


3**SWITCHING ON**

The unit sets for operation whenever the knob placed on the front side of the chassis is turned clockwise from position "0" to position "I".

The clock on the operator's keyboard starts functioning.

4

POSITIONING THE MODULES ARM


The air clutch gets released by wrapping the handle with one's palm and fingers. As long as the hand grips the handle, the height of the modules array can be adjusted. When the hand releases the handle, the modules arm stays at the relevant height.

5

CLOCK ADJUSTMENT

The hour is set by the Manufacturer at testing.
For any setting, keeping pushed both keys:

ENTER and **CLOCK**

press key

until the display shows the proper time.

Hours appear in the two left hand cells an

PATIENT POSITIONING

OPERATOR KEYBOARD

DRIVEN MOVEMENTS: keep the blue pertinent keys pressed.

(The grey keys operate with ISORAMA chair only).

Seat up and down movements and backrest frontward and backtilt movements may be controlled also by chair footcontrol:

CHAIR UP

CHAIR DOWN

BACKREST FORWARD

BACKREST BACK

HEADREST UP

HEADREST DOWN

HEADREST FORWARD

HEADREST BACKWARD

The operation of the unit foot-control prevents the driven movements from starting.

PROGRAMMABLE AND NON-PROGRAMMABLE AUTOMATIC MOVEMENTS: press once the relevant violet key:

POSITION MEMORY 1

POSITION MEMORY 2

POSITION MEMORY 3

POSITION MEMORY 4

LAST POSITION MEMORY (NON-PROGRAMMABLE)

Memories **1** and **2** are available also by the chair footcontrol:

Programming procedure: set the chair by footcontrol or pushbuttons at the desired position, then keep pressing the memory pushbutton where the new position should be stored; the confirmation of the occurred recording is got with a blinking of the led placed on ENTER key.

The position previously in memory for such pushbutton is cancelled and replaced by the new one.

Lifting up any instrument but the syringe and the light cure prevents the automatic movements from starting.

A safety device stops the chair downwards movement, when it meets an obstacle.

The ISORAMA chair is also equipped with a similar safety device in order to stop the backrest backwards tilting in such case.

During chair downwards movement (if the chair down key and the safety device on the base are pressed at the same time) a function is stored, which makes the chair move downwards more slowly every time it passes through the point at which that function was stored. To reset this function, move the chair to the lower position and keep the "chair down" key pressed for about 2 s.

ASSISTANT KEYBOARD

DRIVEN MOVEMENTS: keep the blue pertinent keys pressed.

CHAIR UP

CHAIR DOWN

BACKREST FORWARD

BACKREST BACK

The operation of the unit foot-control prevents the driven movements from starting.

PROGRAMMABLE AUTOMATIC MOVEMENTS: press once the relevant violet key:

POSITION MEMORY 1

LAST POSITION MEMORY

7

ALARM MODE/TIMER


Push key:

CLOCK

that shall start blinking.

By pressing key:

obtain on the display the time desired.

The two left cells show the minutes and the two right cells the seconds. To start counting, push rightward the foot-control lever and release it; during countdown any other control of the unit can function normally.

At the end of the count down a sound warns that time has elapsed.

If, after having set the alarm, you do not desire to start it and want to go back to normal operation mode, just push key:

ENTER

8

TUMBLER FILLER

It can be operated either by pressing and releasing the green key:

TUMBLER

or by turning the ring on tip-holder arm or, if the handpieces lie in resting array, by pushing downward the small pedal in the foot-control. (see chap. 12 "Foot-control switches").

TIME ADJUSTMENT:

Push both keys **TUMBLER** and:

ENTER

the display shall show a count-up in tenths of second; when reaching the quantity of water desired in the tumbler, release the keys.

From now until next adjustment, and whenever the key **TUMBLER** is activated, the tumbler filler shall supply water for the programmed time.

9

BOWL RINSING

The bowl flush starts automatically following the tumbler filler stop, and stops automatically when the programmed time has elapsed. To change the programmed time, see the section:

“PARAMETERS PRESETTING”.

To start the bowl flush independently from the tumbler filling, press the green key:

otherwise, turn the ring located on the tip-holder arm.

10

LAMP OPERATION

The lamp can be switched on by photocell or by switch. The unit should be properly programmed just as the lamp installation. The green key:

switches the lamp on. The lamp can also be switched on and off from the Foot-control lever, by pushing it leftward when the instruments are in resting array.


11

SERVICE PUSH-BUTTON


The green key:


closes any low voltage external circuit (“max 24 V - 1 A”) for the assistant call, door opening, and so on.


By shifting rightward the lever any handpiece except the syringe, the electrosurge, and the light cures can be started. In case of adjustable power handpieces, i.e. micromotors, the supplied power stands proportional to the lever stroke. If no instrument has been lifted, this switch activates the ALARM COUNTDOWN (see chap. 7), provided the pertinent key has been pressed.

By shifting the lever leftwards, the "reverse" mode of the micromotor gets preset, if the micromotor is lifted up. Another shift leftwards restores the rightwise mode. The unit stores automatically this switching. If the electrosurge is up, touching leftwards the lever presets one of the functioning modes (cut, clotted cut, clotting). The selected mode gets automatically stored. If the bistoury or any other instrument, but the micromotor is lifted, the switch has no effect. When instruments are not picked up, touching the lever leftwards turns ON or OFF the ISOSUN light.


By pushing the small pedal downwards, with an instrument lifted, except the syringe and the handpieces without spray, the spray can be enabled or disabled. The unit stores automatically this switching. If all the instruments are lying in resting array, such switch operates the tumbler filler. If a handpiece without spray is up, the switch has no effect. This control starts the electrosurge working.


By pushing the small pedal rightwards, the c.b. of the lifted handpiece starts operating. If a handpiece without c.b. or the syringe or all the handpieces lie in resting array, the switch has no effect.

PROGRAMMABLE MODES FOR THE INSTRUMENTS

when all handpieces are laying, allows to choose the spray source (if the unit is fitted with the tank as an optional item for physiological solution): yellow light (1) on means that the physiological solution is preset; no light on means that pipeline water is preset.


Once a handpiece has been lifted, the key allows WATER and AIR selection of the spray of the lifted handpiece. The YELLOW light lights up when WATER is preset, the GREEN light lights up when AIR is preset. If both lights are off, then the handpiece will work in dry conditions; if you desire to program another condition, first press the key so that the green light starts blinking: the blinking means that the microprocessor is waiting for the new instructions.

the yellow light will turn on; by pressing it again the green light will turn on; a further touch shall switch on both yellow and green. Now the handpiece can work in the chosen conditions, but if it is placed back in the resting array the microprocessor restores the previous conditions in memory. If on the contrary the new conditions must be preserved, the key:

must be pressed: one blink shall confirm that the new instruction has been stored and the previous one cancelled.

first AIR, then WATER, and eventually both shall disappear. Therefore: yellow light on means that only water is preset; green light on means that only spray air is preset; yellow and green lights both on at the same time means that complete spray is preset; no light on means that dry operation is preset.

permits selection of WATER and AIR in the automatic purge at Foot-control release.

This function consists of a water or air or spray jet (emitted for an adjustable time) after the handpiece stops, to the purpose of cleaning both the operating zone and the final parts of the hose. If you do not want to use such device, just program the key with yellow and green lights off.

NOTE: the variables of this key (water and air) depend on the ones of the key - 1 - spray (see table A at page 42). When releasing the Foot-control lever, the preset automatic purge persists during a certain time; such time can be adjusted according to the procedure explained in the section "PARAMETERS PRESETTING".


permits selection of WATER and AIR in the Foot-control c.b.. It is therefore possible to preset a normal c.b. of sole air, or a sole water c.b., or eventually a spray c.b. without bur rotation. Introduction and change of modes is obtained in the same way as seen for the SPRAY.

permits to adjust the maximum power provided to the instrument lifted (such max. power will be obtained when the Foot-control lever is at end stroke). The power can be controlled and displayed in N% from 01% to 100%;

to adjust its max., press the key once so that the red light starts blinking; then push:

until the desired N% gets displayed, and then start working.

If you want to store in memory the new max. power, just press once the key:

that shall confirm storage by blinking once; the instrument replaced, when picking it up again the new max. value shall appear on the two rightwards cells of the display.

The violet key:

FIBRE OPTICS ADJUSTMENT

permits to adjust the light intensity of the fibre optics handpiece that has been picked up, from a minimum "0" to a maximum of "9".

When releasing the Foot-control, the light persists during a certain time; such time can be adjusted according to the procedure explained in the section "PARAMETERS PRESETTING" at page 52, and is the same time for each of the fibre optics handpieces.

The light intensity must be adjusted according to the same procedure explained for the key POWER; the intensity lever shall appear on the first cell of the display.


The violet key:

MODE

permits to choose the rotation mode (clockwise/reverse) of the micromotor.

Once the key has been pressed and the pilot light is blinking, by pushing:

the reverse mode gets preset

and the mark "r" (reverse) appears on the third cell of the display;

to the clockwise mode, return

the key:

must be pressed

the key:

MODE

must also be used to select the operation mode of the electronic surgery (see at page 50).

TABLE "A"

SPRAY KEY	AUTO C.B. KEY choice restricted to:
Yellow light = water	Yellow light = water or Ø = nothing
Green light = air	Green light = air or Ø = nothing
Lights off Ø = nothing	Ø = nothing
Yellow and green lights = both water and air	Yellow light = water or Green light = air or Ø = nothing or Yellow and green lights = both water and air

The unit can hold 8 modules plus the Syringe one, at left and/or at right side of the modules holder. These eight modules can be chosen among a wide range of types, and may include two or more modules of the same type (i.e. two or more turbine modules).

TO REMOVE ANY MODULE FOLLOW THE PROCEDURE BELOW:


- 1 Switch off the unit and press the syringe levers up to the complete depressurization.
- 2 Press the release pin of the module (pin located under the module - see following paragraph "HANDPIECES PRESSURE ADJUSTMENT"), holding the module or the group of modules following the one whose pin is pressed.
- 3 Once all air has flown out, shift the module or group of modules aside.
- 4 Press the disconnection pin of the module preceding the one to be taken out, and once all air is out slide the latter aside.
- 5 Place back the module or group of modules previously taken off, or a new module in between, by sliding in the supporting brackets. Taking care of proper matching between modules edges and tightness protection, switch the unit on again.
- 6 If now the handpiece of the module placed in the position previously occupied by another is picked up, a beep is heard and the display shows "**AL04**". In fact the microprocessor, whose memory matches that particular position with the data of the previous module, is asking for the data relevant to the new one. To store them in, first push the key:

to delete AL04 and push again to store the data shown on keyboard and display.

Such data can be varied according to the procedures above explained for the various functions.

If the position of other modules has been changed, this sequence of operations for reprogramming must be repeated for all of them. All modules are fitted with flow rate adjustments for spray/c.b. water, spray/c.b. air and drive/cooling air, provided they include such operative functions. The adjustment devices are located underneath the module, beside the disconnection valve pin (see "HANDPIECES PRESSURE ADJUSTMENT").


HANDPIECES PRESSURE ADJUSTMENT

- 1 Drive air (Turbines, Airmotor, Powder cleaner) or cooling air (Micromotor, Light cure).
- 2 C.B. or spray air.
- 3 Spray water/C.B. water.
- 4 Whip balance adjustment.
- 5 Disconnection valve pin.
- 6 Syringe air adjustment.

SYRINGE MODULE AND SYRINGE MODULE WITH FIBRE OPTIC


MODULES

This module has no programmable functions and it is the only one that cannot be operated through the foot-control.

1. Control push-button for WATER
2. Control push-button for AIR
3. Selector for hot air and water
4. Fibre optics


When pressing simultaneously the push buttons **1** and **2**, spray jets out.

- The outer envelope can be slipped off from the syringe body for sterilization (compare table at page 50).


TURBINE MODULE AND VARIABLE SPEED TURBINE MODULE

MODULES

When picking the handpiece up, spray, automatic c.b. and c.b. light up their pilot leds.

When a variable speed module is selected, the power adjustment lights up as well: on the two right side cells of the display the memorized power value appears (between 01% and 100%).


(only with variable speed turbine module)

- As for the sterilization of the instrument, compare table at page 50.

MICROMOTOR MODULE


MODULES

When picking up the handpiece, spray, automatic c.b., c.b., power adjustment and rotation mode light up their pilot leds.

On the two right side cells of the display the memorized power value appears as well (between 01% and 100%).


If reversed rotation has been preset, the symbol "r" appears on the second cell.

- As for the sterilization of the instrument, compare table at page 50.


TURBINE MODULE WITH FIBRE OPTIC AND VAR. SPEED TURBINE MODULE WITH FIBRE OPTIC

MODULES

When picking up the handpiece, the function keys of the turbine modules without fibre optic light up their pilot leds, as well as the key for the fibre optic: on the first cell of the display a number from 0 to 9 appears.


When the variable speed turbine module has been preset, also the power adjustment key switches on, while on the two right side cells of the display a number appears: it corresponds to the power intensity (in percent between 01% and 100%) memorized.

- As for the sterilization of the instrument, compare table at page 50.

FIBRE OPTIC MICROMOTOR MODULE

MODULES


When picking up the handpiece, the function keys of the micromotor module light up their pilot leds, as well as the key for the fibre optic: on the first cell of the display a number from 0 to 9 appears, indicating the memorized intensity level of the fibre optic.

On the two right side cells of the display a number in percent (between 01% and 100%) appears, corresponding to the memorized power intensity.

If reversed rotation has been preset, the symbol "r" appears on the second cell.


If the max. programmable power is set between 11% and 100%, the module operates just as the MICROMOTOR MODULE WITH FIBRE OPTIC.

If the max. programmable power is set between 01% and 10%, the high torque/slow speed operation mode gets automatically selected.

To obtain best results when operating at slow speed, it is advisable to make use of reducing handpieces.

MICROMOTOR MODULE FOR IMPLANTOLOGY

It operates as the Micromotor, if the maximum power is set between 11% and 100%. If the maximum power is set between 01% and 10% the slow-speed operation with high torque is automatically selected.

Spray can be obtained in two ways:


- I^a function (position **I**): the external spray source supplies only liquid. Make use of the potentiometer under the module to adjust the spray flow rate.
II^a function (position **O**): spray comes from the unit circuit.

- 1 Spray setting switch
- 2 External spray adjustment (position **I**)
- 3 Internal spray adjustment (position **O**).

To obtain best results when operating at slow speed, it is advisable to make use of reducing handpieces.

POWDER CLEANER MODULE


When picking the handpiece up, spray, automatic c.b. and c.b. light their pilot leds up.

- For the handpiece hygiene, use a cold sterilizing liquid.


SCALER MODULE

MODULES

When picking the handpiece up, the keys of spray and power adjustment light up their leds. Notice that spray can only give water, since air is not used in the scaler.

On the two right side cells of the display the memorized power appears (between 01% and 100%).

- The handpiece can be sterilized according to the indications given in the table at page 50.

Ultrasonic oscillations may prevent cardiac pacemakers functioning properly. Therefore, EMS recommend that patients with a cardiac pacemaker should not be treated with scaler.


The scaler instruments oscillate at high frequency during the treatment. In rare cases, these oscillations may cause the instrument to fracture during the treatment. To prevent the patient swallowing or inhaling the fractured instrument, he must be instructed to breathe through the nose during treatment.

Never use the scaler instruments dry. If used dry, the tips will heat immediately. This may cause thermal injury to the tooth. Ensure that adequate liquid coolant is always available.

SCALER MODULE WITH FIBRE OPTIC

MODULES

When picking the handpiece up, the keys of spray, fibre optic and power adjustment light up their leds. Notice that spray can only give water, since air is not used in the scaler.

The display will show:

a number from 0 to 9 in the first left side cell, indicating the stored intensity level of the fibre optic,
and a number in percent (between 01% and 100%) in the two right side cells, corresponding to the stored power intensity.

- The handpiece can be sterilized according to the indications given in the table at page 50.


Ultrasonic oscillations may prevent cardiac pacemakers functioning properly. Therefore, EMS recommend that patients with a cardiac pacemaker should not be treated with scaler.

The scaler instruments oscillate at high frequency during the treatment. In rare cases, these oscillations may cause the instrument to fracture during the treatment. To prevent the patient swallowing or inhaling the fractured instrument, he must be instructed to breathe through the nose during treatment.

Never use the scaler instruments dry. If used dry, the tips will heat immediately. This may cause thermal injury to the tooth. Ensure that adequate liquid coolant is always available.

Lifting the instrument the light on the handpiece is switched on and the camera begins to telecast live. Pressing once the foot-control lever rightwards, the still is got at full screen. For coming back to the live image mode, the lever of the foot-control should be pressed again. For showing up to 4 images in the same page, the lever should be pressed for about 3 seconds. The images are stored and visualised together on the monitor. For coming back to one image mode keep the foot-control pressed for more than 3 seconds.


When the camera is put back, nothing appears on the screen; otherwise if it was stored one or more images, this/these will remain on the screen.

The handpiece cannot be autoclaved and cannot be immersed in a liquid. It is recommended to use the appropriate hygienic protections.

The position of the module cannot be changed.

This module has no programmable functions and is independent from the foot control.

To start operation, press the button located on the instrument itself. The operation time for each start is fixed at 30 s; an acoustic signal gets activated every 10 s.

Press the button again to stop operation.

- The fibre optics tip is autoclavable according to the indications given at page 50.


Do not point the instrument light towards the eyes.

CHECKING THE ACTA LIGHT CURE


Switch the instrument on and point it at the light sensor (No 1): the green light (No 2) must lit on.

If the red light lits on instead, clean the filter and the fibre optic rod.

Replacement of the curing light bulb

- 1) Loosen the fastening screw.
- 2) Remove the front casing, by squeezing it slightly.


- 3) By a slight traction, remove the light bulb from its slot.
- 4) Remove the bulb from the contacts p.c.b.
- 5) To reassemble, follow the above steps from 4) to 2).

Do not touch with fingers the hollow side of the bulb reflector!

MINI LED LIGHT CURE MODULE

MODULES

This module has no programmable functions and is independent from the foot control.

For starting the instrument, press the push-button 5 placed on the side of the optical terminal. The green led placed underneath will light up, enabling the instrument working.

The polymerization cycle begins pressing the key according to the selected mode. The cycle can be interrupted at any time pressing the key again. Cycle start and end are indicated by an acoustic signal.

There are three working modes selectable with the navigation key 1.

Fast Curing

If it is selected, the red led, placed near the symbol 4, will lit.

The lamp will work for 10 s with a power of 1100 W/cm² (110%) with the tip of ø7.5 mm.

Pulse Curing

If it is selected, the orange led, placed near the symbol 3, will lit. The lamp will work at maximum power with 10 emissions alternating with a rest period of 250 ms.

Step Curing

If it is selected, the yellow led, placed near the symbol 2, will lit. The lamp will reach full power progressively in 10 s, then it will remain on at maximum pressure for other 10 s.

The optical tip is autoclavable according to the indications given at page 50.

Do not point the instrument light towards the eyes.

Do not use this instrument if the patient and/or the operator has a heart stimulator or pacemaker.

Do not use on person who suffers, or has suffered in the past, from photo-biological reactions (including urticaria solaris or protoporphyrine erythropoietique), or any person receiving medical treatment and using photosensitive medication (including methoxsalenes or dymethylchlorotetracycline).

Any person who has previously suffered from a retina or crystalline condition or who has undergone eye surgery, in particular cataract treatment, should visit an eye specialist before undergoing treatment with Mini Led. Even with the patient's agreement, it is highly recommended to use protective glasses at all times.

ELECTRO-SURGERY MODULE

MODULES

For detailed instructions on the bistoury module please consult the user's handbook (code 653270020), supplied with the instrument itself.

15

STERILIZATION METHODS AND WORKING TIMES

N°	INSTRUMENT	STERILIZATION		AUTOCLAVE PRESSURE (bar)	NOTES
		Temperature °C	Time (minutes)		
1	Syringe	120 max	20	1	
2	Acta light cure (fibre optics tip only)	134	6	2	The autoclave should be clean and without oil or other chemical substances. Disinfect with hygienizing liquid, alcohol-based but without acetone
3	Mini Led light cure (optical tip only)	134	4	2	
4	Scaler	134	5 max	2,2	Do not use the chemical autoclave
5	Dry scaler	134	5 max	2,2	Do not use the chemical autoclave
6	Micromotor (Bien Air) (removable cap only) or	135 121	3 20	2,2 1	
7	Turbines (Bien Air)	135	3	2,2	
8	Electrobistoury	128/138	7 (20 max)	IN DRY STOVE Temperature °C 160	Time (minutes) 60
9	Lamp handle	Cold sterilization			

Do not leave the instruments in the sterilizer during the night!
The parameters indicated refer to sterilization real time. The preheating and cooling times are excluded.
For further information, see the instructions contained in the instrument cases.

N°	INSTRUMENT	USING TIMES		NOTES
		ON (minutes)	OFF (minutes)	
1	Micromotor (Bien Air)	2	5	
2	Dry handpiece for scaler (EMS)	5	8	
3	6F-Syringe	3	3	
4	Light cure	50 s	3	
5	Telecamera	1	5	
6	Chair	1	14	

The following messages appear on the display in case of operation:

AL01	Clock failure	Inform technical service
AL02	Faulty memory	Inform technical service
AL03	Faulty data transmission between dental unit and chair	Inform technical service
AL04	Module position not recognised	Reprogram relevant module/s
AL05	Faulty data transmission between base and modules holder	Inform technical service
AL06	Module code not recognised	Inform technical service
AL07	Two or more modules calling	Replace instrument; if alarm code continues, please inform technical service
AL11	Disinfectant to be refilled	Fill up
AL13	Separator filled up	Inform technical service
AL15	Disinfectant tank faulty pressure	Inform technical service
AL21	Counteracting movements on the same motor axis	Inform technical service
AL22	Overtime motors functioning (>30 s)	Only manual controls can be used; inform technical service
AL23	General problem on the chair	Inform technical service
AL24	Programming not accepted	Only manual controls can be used; inform technical service

When the message appears, it is accompanied by three beeping sounds.
By pressing once the key:

the alarm (except for **AL07**) goes off.

If the cause of the alarm message is not settled after pressing the ENTER key, the message itself will not appear any more until the next switching on of the unit; if the cause of the alarm message stops and then starts again, the message will again appear.

☞ **ATTENTION:** Before calling the Technical Service, try to switch off the unit and to switch it on again.

17

PARAMETERS PRESETTING

Some programmations are made once in a certain period or once for ever, and do not need continuous variations afterwards. Such programmable parameters are not accessible from the operator keyboard when the unit is in working mode, and must be reached through a separate procedure.

There are 4 parameters:

- persistence of cold light after foot-switch release,
- persistence of automatic c.b. after foot-switch release,
- bowl flush time,
- wet-line suction with amalgam separator.

To access the parameters setting mode, keep the:

key pressed and switch on the unit, then release the key.

The display will show **P008**, where "P" stands for parameters mode, "0" is the first changeable parameter, and the two numbers at right translate the value of such parameter. To change this value, (between 00 and 99), the keys:

must be pressed. To shift to parameter **1**, press the ALARM/CLOCK key; after **P3**, **P0** is shown again.
To get out of the parameters mode and start working, press the key:

THE FOUR PARAMETERS ARE THE FOLLOWING:

P008

End-operation cold light time (in seconds)

P101

End-operation automatic c.b. time (in seconds)


P220

Bowl flush time (in seconds)

P301


Wet-line suction with amalgam separator

MODULE-HOLDER ARM BALANCING


When modules are added or removed, as well as side tray, it is necessary to balance the arm tension and to adjust its horizontality.


To access the mechanism, move the sheath towards the module-holder, so that the carters get uncovered.

Push the sides of the cover so as to release its cogtongues, and pull upwards.


By rotating clockwise the spring Allen head, the arm shall hold a heavier weight; by rotating the same counterclockwise, the arm shall hold a lighter weight.

To settle the module-holder horizontality (longitudinal direction), loosen the two nuts of the hexagon bar. Then, by rotating it clockwise, the keyboard profile shall tilt upwards; by rotating it counterclockwise, the keyboard profile shall tilt downwards. Do not forget to tighten the two reinforcing nuts.

WATER DRAIN FROM AIR FILTER

MAINTENANCE


The unit is fitted with a filter that prevents water vapour condensation from getting to the air circuits.

Such filter is located underneath the cuspidor frame and should be checked weekly; if water is in it, it must be drained by pushing the pin located at the filter cup bottom.

CLEANING OF SUCTION FILTER

MAINTENANCE


The filter is located into the unit front wall.

To get access to the filter, push the wall against the frame (on right top side): the wall is kept by hinges at its bottom, so it shall remain open in tilted position, which makes it very easy to take out the filter. The filter can be handled without touching its dirty parts, by holding the piece as in figure.


The filter must be cleaned daily and replaced at least once a week.

FILLING OF MEDICAMENT TANK

MAINTENANCE


IMPORTANT: preset


SELECT

the yellow pilot light over the key on before starting this operation.
Tilt the unit wall as shown for cleaning of the suction filter.
Turn the lever rightwards as indicated: the side shutter shall swing outwards.
Accompany its movement until it rests in tilted position.

Pull the glass tank out, so as to make it easy to open its plug.
Remove the plug and pour 12 ml. of the appropriate solution by means of the proper dispenser.
Close the plug and turn the tank towards inside.

On the keyboard, press the violet key:

SELECT

The self-filling of the tank will operate automatically and will take about one minute.
To close back the side shutter, push it against the cuspidor side until the click is heard.
A beep shall indicate that the tank is full. Close then the tilting wall.

It is recommended to clean the bottle daily and sterilize it in autoclave. Sterilization must be carried out at max. 135°C, and then the bottle must be left 20 minutes to cool slowly down.

FILLING OF PHYSIOLOGICAL SOLUTION


MAINTENANCE

Proceed as above explained as to opening the side panel and pulling the glass tank off the chassis.
Remove the plug and pour the physiological solution in.

To close back the shutter, push it against the cuspidor until the click is heard.

Close back the cuspidor wall.

Once the tank has been emptied and before filling up with new physiological solution, it is recommended to switch the water supply to mains (yellow pilot light over the key

SELECT

off), and then have sprays and syringe water operated for about 2'.
Such operation can avoid the layering of jellies favourable to bacteria development.

It is recommended to clean the bottle daily and sterilize it in autoclave. Sterilization must be carried out at max. 135°C, and then the bottle must be left 20 minutes to cool slowly down.

IMPORTANT:


ONLY THE USE OF OUR FORMULATION (1 DOSE EVERY 100 LTS. DILUTION) PERMITS TO OBTAIN THE PROPER EFFICACY OF THE AUTOMATIC DOSAGE AND PREVENTS MALFUNCTIONING DUE TO USING IMPROPER PRODUCTS.

Efficacy time: 2 years.

Storage: dry and cool place, repaired from light.


REPLACEMENT OF STERILE SOLUTION BOTTLE

Renew bottle and sterile hose.
Introduce the hose in the outside spray handpiece.


EMPTYING THE WASTE COLLECTOR

MAINTENANCE


To empty the waste collector, use a disposable cup of the type supplied with the unit.


Place a large bowl under the collector and twist it out as indicated from arrow **a**.

Cover the collector with the disposable cup, with the bottom side up.

Pressing the collector and the cup well together to prevent dripping, turn upside down and wait until all the waste has sunk from the collector into the disposable cup.

Close tightly the disposable cup with its cover, and twist the collector back into its seat as indicated from arrow **b** (see fig. 1).


SURFACE CLEANING

MAINTENANCE

Sanitizing liquid for plastics is recommended for surface cleaning.

Handpieces tubings must be cleaned with water and neutre soap, and then must be well dried and lightly spread with talcum powder.

The chair and stool upholsteries must be cleaned with water and neutre soap.

CONTROL PANEL

1. GREEN LIGHT: Power pilot lamp

Normal condition:

When the equipment is connected, the green pilot lamp **1** lights continuously.

Fault:

After the equipment connection, if the green pilot light **1** does not light, check fuse and plate and, if necessary, replace them.

2. ORANGE LIGHT: Centrifuge control

Normal condition:

As a rule, this pilot lamp should not light. In case it lights, turn the dental unit on and off several times (max. 5 times), in order to check whether any fault occurred. If light does not extinguish, a fault is present.

Fault:

A system controls power input to centrifuge's motor. When input is too high, the light **2** starts to operate, to indicate the presence of a faulty condition. At the same time, the suction valve does not open and none operation can be carried out.

- That condition can take place even when filter is removed from its box, because some large particles could have clogged and blocked the centrifuge.
- The same condition can take place when a motor's fault is present.
- In order to eliminate any fault, replace the centrifuge set.


3. YELLOW LIGHT: Control of amalgam level

Normal condition:

When the equipment is switched on, an electronic system controls several activities, among which the amalgam level. When a 100% level is obtained, the yellow pilot lamp lights continuously and the warning signal can not be cut out any more. At the same time, the suction solenoid valve does not open and none operation can be carried out. Thus, it will be necessary to replace the amalgam container, in order to allow the equipment to operate properly.

REMOVING AND CLEANING THE FILTER

METASYS

In order to take out the box filter, follow the easy and short indications hereunder described.

Firstly, remove the threaded cap placed on the centrifuge. Then, take out filter by pulling it outwards.

Now, proceed to clean it, by keeping it under lukewarm water so as to remove all waste. Then, carry out the same operations backwards, by inserting the filter and taking care of placing it properly on its guides. Finally, screw the cap carefully.

REMOVING THE AMALGAM SEPARATOR VESSEL

METASYS


Fig. 1

When the illuminated yellow lamp signals that the complete filling up of the amalgam vessel has been reached, it is necessary to remove the container from its housing and replace the tank, in order to restore the proper working of the equipment.

So as to carry out this last operation, lift the releasing handle placed on the front side of the separator (picture 1).


Fig. 2

After releasing the separator from the plate, take the body out of its sliding guides, simply by pulling it outwards (picture 2/A).

After removing the whole container, you can proceed to dismount the tank, by releasing the fixing hinges placed on the sides of the vessel (picture 2/B). Then, replace tank by a new, supplied one.


Fig. 3

Before forwarding the amalgam tank, carry out the following, simple and precautionary operations.

Cut off one of the corners of the supplied envelope-bag.

Empty the whole content into the container still open (picture 3).


Fig. 4

Close the container carefully!

All tongues have to snap shut loud (picture 4).


Fig. 5

Then, check visually that all tongues are clamped properly.

If necessary, press firmly on the parts indicated by the arrows (picture 5).


Fig. 6

Now, remove the sheet, on which all the addresses of the factories dealing with waste treatments are printed, from the inside of the container packing. The mentioned sheet will be useful at the moment of the sending.

Then, after checking that the container is well closed, place it into the foamed shell (picture 6).

Finally, in order to protect the container, lay the other half-shell on it.

Please, do not forget to enclose your own address before forwarding the packing.


Fig. 7

After packing the container properly by means of the suitable polystyrene shells, put it into the carton box.


Fig. 8

Now, proceed to close the forwarding case in the following way:

- Place the middle tongue in its housing (picture 7)
- Fold the right and left safety tongues as described below:


Fig. 9

- a) Fold the pointed ends upwards and inwards (picture 8).


Fig. 10

- b) By providing a firm pressure, compress the pointed ends (picture 9).
- c) Then, insert tongue into the proper slot, by folding the pointed ends (these last ones have to be turned downwards - picture 10).

ATTENTION!

If you unfold the tongues, the case will be damaged!

Finally, send the case to the factory in charge of waste treatment.

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