C1⁺

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
Operating Instructions

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Deutsch / English / Français / Español

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Deutsch

Operating Instructions

English

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Français

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
Español

C1⁺

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Deutsch

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
General information

Dear customer,

We thank you for purchasing your C1⁺ dental treatment center from the Sirona company.

The technical documentation supplied is also part of the product. You should always keep this documentation within reach.

These **Operating Instructions** describe your treatment center with maximum possible equipment.

For daily care of your treatment center, please observe the instructions provided in "**Care and Maintenance by the Practice Team**".

To safeguard your warranty claims, please complete the attached "**Installation Report / Warranty Passport**" together with the service engineer immediately after the installation of your treatment center.

Separate Operating Instructions with corresponding instructions for care are enclosed with the dental instruments, SIROLUX FANTASTIC, SIROCAM 3, tray and X-ray image viewer.

Prior to start-up, you must read all Operating Instructions to familiarize yourself with the treatment center.

To prevent any personal injury or material damage, pay special attention to any notes printed in bold type or marked in one of the ways indicated below:

Additional information, explanation or supplement

 **NOTE**

 **CAUTION**

The information provided under this keyword directly concerns the functioning of the product and/or warns against possible operating errors.

 **WARNING**

This warning symbol warns against possible hazards for the product or the user.

In case you get stuck despite having thoroughly studied the Operating Instructions, please contact your dental depot.

Your C1⁺ Team

Table of Contents

1	Warning and safety information	7
2	Technical description	11
3	Operating and Functional Elements.....	12
4	Putting the unit into operation	14
5	Foot switch, displays	16
6	Starting the treatment center/instruments.....	17
	6.1 Depositing treatment instruments	18
7	Key functions on the dentist element	19
	7.1 Neutral dialog.....	19
	7.2 Chair program dialog	20
8	Key functions in the assistant dialog.....	22
9	Patient positions, Chair program dialog	26
	9.1 Safety circuit	26
	9.2 Display of safety switches.....	27
	9.3 Chair program dialog	28
10	Patient positions, manual chair setting dialog.....	30
11	Programming automatic positioning.....	32
12	Instrument dialogs	36
	12.1 Instrument dialog: Electric motor	36
	12.2 Instrument dialog: Turbine	39
	12.3 Instrument dialog: SIROSONIC L	41
	12.4 Instrument dialog: SIROTOM	42
	12.5 Safety checks	46
	12.6 Switching over from instrument dialog to Chair program dialog or to manual chair setting dialog ...	46
13	Water unit, assistant element.....	47
	13.1 Swiveling the cuspidor	47
	13.2 Adjusting the water quantity.....	47
	13.3 Amalgam rotor	47
	13.4 Disinfection system.....	48
	13.5 During treatment.....	48
	13.6 Suction handpieces on the assistant element	49
	13.7 SPRAYVIT	50
	13.8 SPRAYVIT setting dialog.....	51

English

14	Service dialog	52
14.1	Function of the keys, left row:	52
14.2	Function of the keys, right row:	53
14.3	Setting the time and date	53
15	Basic setting dialog 1	55
15.1	Function of the keys, left row:	55
15.2	Function of keys, right row:	57
16	Basic setting dialog 2	58
16.1	Function of keys, left row:	58
16.2	Function of the keys, right row:	59
17	SIVISION 3 / SIROCAM 3 / SIROCAM C	61
17.1	SIROCAM 3 Intraoral Camera (additional equipment).....	61
17.2	SIROCAM C Intraoral Camera (additional equipment)	63
17.3	Camera dialog: General notes and SIVISION 3 operating modes.....	65
17.4	Camera dialog: Video mode	65
17.5	Camera dialog: PC mode.....	68
17.6	SIVISION 3 – Second monitor function without SIROCAM 3 / SIROCAM C (PC dialog)	71
18	Simplified dialogs	72
18.1	Simplified Chair program dialog.....	73
18.2	Simplified assistant dialog.....	73
18.3	Simplified Neutral dialog	74
18.4	Simplified Service dialog.....	74
18.5	Simplified instrument dialogs	75
19	Surgery with saline solution (NaCl).....	76
19.1	Preparations on the unit.....	76
19.2	Presetting.....	76
19.3	Operation	77
20	POLYLUX fiber optic light handpiece	78
20.1	The handpiece	78
20.2	Switching on with POLYLUX in the assistant element.....	79
20.3	Switching on with POLYLUX light in the dentist element.....	80
20.4	Programming the exposure time	81
20.5	Application	81
21	Additional equipment and accessories	83
21.1	Media block on the patient chair	83
21.2	Hygienic headrest protection	84
21.3	Seat cushion C.....	84
21.4	Hydrocolloid coolant supply	85
21.5	Tray	85
21.6	X-ray image viewer on the lamp support tube (additional equipment)	86
21.7	X-ray image viewer on the tray (additional equipment).....	86

22

Maintenance.....	87
22.1 Care and cleaning by the practice team	87
22.2 Inspection and maintenance	87
22.3 Safety tests	87
22.4 Maintenance Manual	89

English

1

Warning and safety information

Intended use

This dental treatment center is intended for therapy, diagnosis and dental treatment. in humans.

This unit is **not** intended for **operation in areas subject to explosion hazards**.

General safety information

The ‘**on-site installation**’ must be performed according to our requirements.

As manufacturers of dental medical equipment and in the interest of the operational safety of your system, we stress the importance of having **maintenance and repair** of your treatment center performed only by ourselves or by agencies expressly authorized by us. Furthermore components influencing the safety of the device should always be replaced with **original spare parts** upon failure.

When having such work done, we suggest that you request a certificate stating the type and extent of work performed, if necessary, including statements concerning any modifications of the rated parameters or of the operating ranges, as well as the date, name of organization and signature.

Changes to this unit which could impair its safety for the system owner, patients or third parties are forbidden by legal regulations!

For reasons of product safety, only original Sirona accessories approved for this product, or accessories from third parties approved by Sirona, may be used. The user is responsible for dangers resulting from the use of non-approved accessories.

If any devices not approved by Sirona are connected, they must comply with the applicable standards:

EN 60950 for information technology equipment, and EN 60601-1 for medical electrical equipment.

Combination with other units

Any person who assembles or modifies a medical electrical system complying with the standard EN 60601-1-1 (Safety requirements for electromedical systems) by combining it with other equipment (e.g., by connecting it with a PC) is responsible for ensuring that the requirements of this regulation are met to their full extent for the safety of the patients, operators and environment.

In case of doubt, contact the manufacturer of the system components.

Information on avoiding, recognizing and eliminating unintended electromagnetic effects

Interference of dental medical devices caused by radio telephones

Unit with built-in video camera and monitor

Quality of water / air supply

Water quality compliance

HF surgery

Identification of warning and safety information

Symbols used

Ventilation slots

C1⁺ is a class A device (classified according to CISPR 11, EN 60601-1-2: 1993 based on IEC 60601-1-2). This device may be operated in a residential area provided it is used under the responsibility of a medical specialist.

To guarantee the operational safety of electromedical devices, the **use of mobile radio telephones must be prohibited** in the practice or hospital area.

Electrical disturbances in the environment of the C1⁺ (e.g. lightning) may result in reduced image quality. In such cases, it is advisable to repeat the exposure before saving the image.

Image artifacts may possibly occur in the event of radio-frequency interference caused by the transmitters of e.g. radio services or radio amateurs.

Air and water supply must meet the requirements specified in the Installation Instructions. Use only clean water.

To ensure compliance with the medical and national legal requirements for water from treatment centers, Sirona recommends equipping the treatment center with a disinfection system. As owner of the treatment center, you are responsible for the water quality and may have to take alternative measures to ensure its compliance if you operate the treatment center without a disinfection system.

Please contact your specialized dealer or your relevant dental association for the respective national requirements and measures.

This dental treatment center is available both with and without a high-frequency surgical device (optional).

Only in the Federal Republic of Germany The user is obligated to keep a "medical product log" if a HF surgical device is installed! (see 22.3 "Safety tests")

To avoid personal injury and material damage, you must also observe the warning and safety information provided in the present Operating Instructions. They are specially identified with  **NOTE**,  **CAUTION** or **WARNING**.

Observe the accompanying papers (on the chair nameplate)

Under no circumstances may the ventilation slots on the unit be covered, since otherwise the air circulation will be obstructed.

Do not spray disinfectants or other similar products into the ventilation slots.

Chair and dentist element

Please make sure that no foreign parts (e.g. tray) project into the movement range of the chair and dentist element. This is necessary to avoid collisions during automatic movement sequences!

Dentist and assistant element

In order to avoid damage to the dentist and assistant elements and to the chair upholstery, make sure that these elements do not protrude into the area through which the chair moves. The ventilation slots must be free!

Maximum load

The maximum load on the dentist and assistant element may be 1 to 3kg. If the load is more than 1kg, the spring force must be readjusted by the service engineer (see Installation Instructions).

Vacuum system

The suction of aluminum and other metal oxides from blasting devices via the automatic separator integrated in the treatment center and the amalgam separator is prohibited!

This would cause extreme wear and clogging of the suction and water paths.

A separate vacuum system must be used in connection with metal oxide blasting devices.

Treatment centers equipped with central wet aspiration are generally suitable for suction removal of the above material. However, please be sure to observe the instructions provided by the manufacturer of your vacuum system.

No restrictions apply when using salt blasting devices in connection with Sirona treatment centers. However, in such cases, make sure that the system is subsequently rinsed with an adequate amount of water.

Patient chair

Please observe the maximum load capacity of the chair of **135kg according to EN ISO 6875 (tested with four times safety)**

Maintenance of the treatment center

Despite the outstanding quality of your treatment center and regular maintenance and care by the practice personnel, it is required in the interest of operational safety that preventive maintenance is performed at predetermined intervals.

In order to guarantee the operational safety and functional reliability of your treatment center and to avoid damage due to natural wear, as the system owner it is necessary that you have your dental depot dispatch an authorized service engineer to check your unit according to the appropriate maintenance certificate at regular intervals (at least once a year).

Please contact your dental depot to obtain a maintenance offer.

Dismantling and reassembly

For dismantling and reassembly of the device, proceed according to the Installation Instructions for new installation in order to guarantee the operability and stability of the device.

Disposal

Public legal provisions can contain special regulations regarding the disposal of this device. To prevent environmental damage and personal injuries, we therefore ask you to contact us prior to discontinuing the use of this product, with the aim of disposal.


2 Technical description

Model designation	C1+
Power supply connection	230VAX 50Hz, 115 V~ 50/60 Hz 100 V~ 50/60 Hz
Nominal current	3A at 230 V 9.5A at 115V 9.5A at 100V
Main unit fuse	for 230 V~: T 6,3 A H , 250 V~, Order No. 10 77 452 for 100/115 V~: T 10 A H , 250 V~, Order No. 10 77,460
Operating conditions	Ambient temperature: 10°C – 40°C (50°F – 104°F) Relative humidity: 30% – 75% Air pressure: 700hPa – 1060hPa
Transport and storage conditions	Temperature: -40°C – +70°C (-40°F – 158°F) Relative humidity: 10% – 95% Air pressure: 500hPa – 1060hPa
Protection class	Class I equipment
Degree of protection against electrical shock:	Type B application unit  apart from SIROCAM 3 / SIROCAM C and SIROTOM: Application parts of the type BF 
Degree of protection against ingress of water	Ordinary device (not protected) The foot switch is protected against dripping water IPX 1.
Mode of operation	Continuous operation with intermittent loading corresponding to the dental mode of working. Permanently connected unit.
Year of manufacture	 20XX (on the chair nameplate)
Tests / approvals	This dental treatment center complies with the requirements of EN 60601-1: 1990 + A1: 1993 + A2: 1995 and the EN 60601-1-2: 1993.

DVGW: This unit complies with the technical rules and requirements on safety and hygiene for connection to the drinking water supply.

3 Operating and Functional Elements

- 1 Assistant element
- 2 POLYLUX fiber optic light in the assistant element
- 3 SPRAYVIT in the assistant element
- 4 Suction handpiece
- 5 Saliva ejector
- 6 Display on the assistant element
- 7 Swiveling cuspidor
- 8 Tumbler filler
- 9 Holder for X-ray sensor
- 10 Water unit with amalgam rotor and disinfection system
- 11 X-ray box
- 12 Patient chair
- 13 Receptacle 1: SPRAYVIT in the dentist element

- 14 Receptacle 2, 3, and 4: Electric motor / turbine burr drives
- 15 Receptacle 5: SIROSONIC L scaler or 4th burr drive (electric motor / turbine)
- 16 Holder 6: SIROTOM electrosurgical unit or SIROSONIC L scaler or POLYLUX Fiber Optic Light Handpiece or SIROCAM 3 / SIROCAM C intraoral camera
- 17 Additional instrument holder (only for SIROCAM 3 / SIROCAM C)
- 18 Display on the dentist element
- 19 Dentist element (with motor drive)
- 20 Unit main switch
- 21 Headrest adjustable by motor drive
- 22 Chair foot control
- 23 Unit foot switch

4 Putting the unit into operation

English


Initial start-up

The disinfection system adds a disinfectant to the water (1:100) to prevent the formation of microorganisms in the water system.

Prior to initial start-up of your treatment center, **sanitation** must be performed.

If, on the basis of an agreement with you, sanitation was not performed by the service engineer following installation of your treatment center, please perform sanitation as described in the separate instructions "**Care and Cleaning by the Practice Team**".

Sanitation takes approx. 24 hours.

Main switch

Switch **ON** the main switch (I).

Following power-up, the treatment center automatically performs a self-test.

The displays on the dentist and assistant elements light up; the treatment center is ready for operation.


⚠ CAUTION

For safety reasons, always switch the treatment center OFF (O) after consulting hours.

This cuts off air, water and line voltage.

You can communicate with your treatment center via **dialogs**.

The functions are selected by activating the gray keys **next to** the symbols.

Self-test

The displays on the dentist and assistant elements reach their full brightness step by step.

Simultaneously, the treatment center automatically performs a self-test.

On completion of the self-test, the **neutral dialog** appears.


The chair must be in its neutral position. See Chair program dialog, key **0**.

Ask the patient to take a seat on the chair.

CAUTION

Maximum load capacity of the chair 135 kg according to EN ISO 6875 (tested with six times safety)

The patient's arms and legs must rest on the upholstery of the chair!

Adjusting the headrest

Move headrest out ↗.

Move headrest in ↘.

Hygienic headrest protection see page page 83.

Seat cushion C see page page 83.


Preselecting the user

If this treatment center is used by different persons:

Preselect user **A**, **B** or **C**.

Setting dialog for the dentit's body height page 32.

5 Foot switch, displays

Foot switch

All functions can be selected with the **foot switch** or **by hand** by actuating the gray keys on the displays on the dentist and assistant element.

When using the foot switch for selection, place your entire foot on the surface of the foot switch pedal.

Foot switch functions are only indicated on the display of the dentist element by the cursor points ● .

Preselecting functions

Actuate the foot switch (factory setting)

at **A***: cursor ● moves **upwards** ↑ ,

at **B***: cursor ● moves **downwards** ↓ .

Releasing functions

Swivel


the pedal surface to the left C ← or to the right → D.

The display of the corresponding symbol appears inverted. Symbols of activated functions appear with a **bold frame**.

Displays

⚠ CAUTION

*When selecting by hand,
slightly press the gray keys at the right or left side. Do not
press the symbols next to them!*


NOTE

* On request, foot switch functions A and B can be exchanged by the service engineer.

6

Starting the treatment center/instruments

Starting the chair movement


The Neutral dialog appears.

Brief actuation of the foot switch at **A*** causes the treatment center to travel into the **examination position** and the unit lamp to be switched on.
(freely programmable, see page page 32).

At the same time, the dentist display switches over to the **Chair program dialog** and the **assistant display** switches over to the **assistant dialog**.

The chair starts to move slowly, and then increases its speed. While the chair travels, the second symbol from the top right on the displays changes to a **flashing warning triangle**.


MOVEMENT STOP: see page page 28.

CAUTION

All movements are stopped at once:


- after actuation of the flashing key with the triangle or
- after actuation of any other program key **except the activated key**.

Starting a treatment instrument

A removed instrument is activated by pressing the foot switch at **A***.

6.1 Depositing treatment instruments

A ball stopper for an unoccupied instrument holder is enclosed with the dentist element.

Plug the ball stopper (1) into the unoccupied instrument holder to prevent a treatment instrument from being deposited inadvertently in this holder.

To block further unoccupied instrument holders, more ball stoppers can be ordered under Order No.: 58 99 575.


English

7

Key functions on the dentist element

7.1

Neutral dialog

Function of keys, left row

Move headrest out ".

Switch to the Chair program dialog
(see page page 19).

If this treatment center is used by different persons:
Preselect the corresponding user **A**, **B** or **C**.

Setting dialog for the dentit's body height page 32.

Key for freely selectable function

can be switched **ON / OFF** as required.

Function of keys, right row

Move headrest in ↘.

Switch to the manual chair setting dialog.

X-ray image viewer


can be switched **ON / OFF** as required.

Unit lamp

can be switched **ON / OFF** as required.

7.2 Chair program dialog

English

Function of keys, left row

Program S* with special function.


If this key is pressed a second time, the treatment center returns to the last treatment position (memory function).


Example:

You have programmed this key with the mouth rinsing position. On actuation of this key, the treatment center moves to this position.

The **S**pecial function: If this key is pressed a second time, the treatment center returns to the last treatment position.

Switch to the assistant dialog.

Program 1*, maxilla and mandible

Maxillary molars

Right: palatal and occlusal.

Left: palatal, occlusal and buccal (vestibular).

Maxillary incisors: palatal.

Program 3*, maxilla and mandible

Maxillary molars right:

Buccal (vestibular),

Mandibular molars

Right: lingual.

Left: lingual and buccal (vestibular).

Program 5*, freely assignable,
e.g. collapsing or cleaning position.


Program P* with special function
as with program **S**, e.g. for patient pause position.

i NOTE

* Chair programs: If the keys are actuated longer, the factory settings appear (see page page 27).

Freely selectable function,
e.g. call key (connected by your service engineer).

Function of keys, right row

Program 0*

Neutral position

Switch to the manual chair setting dialog.

Left symbol (triangle in circle) **MOVEMENT STOP** (see page page 28).

Program 2*, mandible.
Mandibular molars right:
Occlusal and buccal (vestibular).
Mandibular molars left: occlusal.
Incisors: lingual

Program 4*, maxilla and mandible.
Incisors: labial.

Program 6*, freely assignable
e.g. collapsing or cleaning position.

i NOTE

* Chair programs: If the keys are actuated longer, the factory settings appear (see page page 27).


Standing or sitting treatment.

Unit lamp ON / OFF.

If this key is actuated longer, the brightness setting dialog appears.

8 Key functions in the assistant dialog


On the assistant element

The **assistant** dialog automatically appears on the assistant display when the treatment center has been switched on and has traveled into the examination position.



In this dialog the assistant can select the assistant functions and, if desired by the dentist, select and set further functions.

For functions of keys, see next page.

On the dentist element

If you are working without an assistant, the **assistant** dialog can also be selected in the Chair program dialog on the display of the dentist element.

**Function of keys, left row:****Automatic tumbler filling.**

As long as the symbol appears with a bold frame, the water is running.

If the key for tumbler filling is actuated longer, the setting dialog for the filling quantity appears.

In this program the filling quantity can be changed in 4 steps with the – / + keys.

The tumbler filling function can also be activated by pressing the **S** key (see page page 34).

Hydrocolloid-operation (see page page 84)

When the **Hydro key (bold frame)** is activated, tumbler filling then switches to continuous operation.

On returning to the assistant dialog (automatically after 15sec or manually with the arrow key), the tumbler filling key is used for activating the hydrocolloid mode.


This mode can be switched off by pressing this key again or with the Hydro key in the settings dialog. After it is switched off, tumbler filling can be switched on again with the programmed tumbler filling time.


**Switch to the camera dialog**

The camera dialog remains disabled as long as another instrument is activated.

**Switch to the PC dialog,**


if no camera is integrated in the dentist element and a PC connection exists (see page page 70).

**Preselection of user
A, B or C****X-ray image viewer ON / OFF.**

**Timer Start / Stop**

Actuate the **Start** key. The time set runs down. At the same time, **Stop** appears. The progress indicator fills up.

An acoustic signal indicates when the time has elapsed. The timing can be stopped by pressing the **Stop** key.

**Switch to the sanitation program**
(see Maintenance Instructions)**Freely selectable function**


e.g. call key (connected by your service engineer)

The function can be preselected as key or switch in the basic setting dialog 1 (see page page 54)
or

SIVISION control

(only as key)


The desired key function can be preselected in the basic setting dialog 1 (see page page 54).

**Function of keys, right row:**

English

**Cuspidor flushing ON / OFF.**

After the key is actuated, the amalgam rotor starts running (symbol flashes). The flushing water starts to flow after a short delay (symbol appears with a bold frame).

**Programming the flushing time:**

Press the **cuspidor flushing** key until the setting dialog appears.

In this dialog the flushing time can be changed in 8 steps ranging from 20 to 300 seconds (shown in the center of the display) with the **-** / **+** keys.

The cuspidor flushing function can also be activated by pressing the **S** key (see page page 34).

**Switch to the**

Chair program dialog(see page page 19).

**Program S with special function.**

On actuation of this key, the treatment center travels into the mouth rinsing position.
If this key is pressed a second time, the treatment center returns to the last treatment position.

Programming the S key (see page page 34).

**Switch**

to the Service dialog


**Preselection of timer setting**

Three individually adjustable timer settings can be selected
(factory setting: 10, 30, 60 sec).


The settings are selected by repeatedly touching the preselection key lightly and briefly.

If this key is actuated longer, the setting dialog appears.
Here the times can be changed separately for seconds and minutes by actuating the - / + keys.

Setting range continuously 1sec – 60min

**Composite key C (for unit lamp)**

This key can be used to set the unit lamp to minimum brightness prior to making a composite filling in order to prevent premature hardening of the filling material.
It can be reset by actuating the same key or the unit lamp ON / OFF key.


**Unit lamp ON / OFF**

If this key is actuated longer, the brightness setting dialog appears. Here you can adjust the brightness of the activated lamp in 9 steps (as long as the C composite key is not actuated).

9 Patient positions, Chair program dialog

9.1 Safety circuit

⚠ CAUTION

The patient's arms and legs must be resting on the chair upholstery during the program run!

ℹ NOTE

A built-in safety circuit stops the chair movement if

- the foot support Δ,
- the backrest Δ or
- the tilting part of the headrest Δ collides with an obstruction, or
- if a swiveling cupid is swiveled in during a chair movement.

The chair then automatically moves upwards a short distance until the path has been cleared (except for the swiveling cupid).

English





9.2 Display of safety switches

When the safety circuit is activated, various symbols under the odontogram will flash.

Cuspidor

These symbols flash alternately:

- if the cuspidor is swiveled in during a chair movement,
- if the cuspidor cannot be swiveled back when a chair movement is started.

Tilting part of the headrest

These symbols flash alternately if the tilting part of the headrest collides with an obstruction during a chair movement.


Backrest

These symbols flash alternately if the backrest collides with an obstruction during a chair movement.

Foot support

These symbols flash alternately if the foot support collides with an obstruction during a chair movement.

9.3 Chair program dialog

The Chair program dialog
is indicated on the display.

The **programs 1–4**, which are assigned to the odontogram, have default settings based on ergonomic aspects. They enable the dentist to look directly at the corresponding area of the jaw from an optimal sitting or standing position.

For this purpose, the dentist **A**, **B** or **C** has to enter his/her height once (for programming see page 32).

All programs can be assigned individually by any user of the treatment center (see page 32).

The factory settings can be restored by pressing and holding down the respective program key. The odontogram then displays the

symbol →←

Select and activate the desired **program**.

NOTE

For program assignment, see “**Chair program dialog** on page 19.”

The chair and the dentist element move into the corresponding treatment position. Then the unit lamp switches on automatically (if programmed).


When the mouth rinsing program **S** or the neutral position program **0** is activated, the lamp switches off automatically (if programmed).


In the manual chair setting dialog

the patient position can be optimized at any time.
(Detailed description see page 29).

Both during manual and automatic chair movement, the second symbol from the top at the right of the displays becomes a **flashing warning triangle**.

MOVEMENT STOP: see next page.

MOVEMENT STOP!


On activation of an automatic chair program, the cursor
• automatically jumps to the second key from the top.

As a **WARNING** during automatic movement, the second key from the top (triangle in a circle) in the right row of keys always flashes together with the selected program key.


!**CAUTION**

All movements are stopped at once:
– after actuation of the flashing key with the triangle or
– after actuation of any other program key except the previously selected program key.

10 Patient positions, manual chair setting dialog

English

Optimizing the patient position

The patient position can be directly optimized for two chair functions with the foot control located on the base of the chair.

This does not require invoking a manual chair setting dialog.

CAUTION

If, in basic setting dialog 2 (see page page 57), the suction hoses and/or the POLYLUX light (in the assistant element) are preselected for being toggled ON / OFF with the foot control, they must be deposited in their receptacles!


Factory default setting

Lever up ↑:

The chair moves up.

Lever down ↓:

The chair moves down.

Lever to the left ←:
Backrest adjustment backward.

Lever to the right →:
Backrest adjustment forward.


Optional setting

You can change the left/right function in **basic setting dialog 2** (see page page 59) to

Lever to the left ←:
Headrest adjustment backward.

Lever to the right →:
Headrest adjustment forward.


11 Programming automatic positioning

Setting the dentist's height


All factory-set programs can be assigned individually by the users **A**, **B** or **C** (overwriting the factory setting).

- In the Chair program dialog, select the Neutral dialog **0** or the assistant dialog.


- In the Neutral dialog **0**


or in the

- Assistant dialog
Preselect the corresponding user **A**, **B** or **C**.
User **B** or **C** is indicated in the center at the top of the display.
- Actuate the key until the **height setting dialog** for the user appears.

- Actuate the – or + key to set your height.
Setting range **150 to 200cm**.
- After a short while, the display switches back to the previous dialog.
You can also switch over with the arrow key.


Programming

Chair and dentist element
(overwriting the factory setting).

Select and activate one of the following programs:

- Programs **1–4** assigned to the odontogram
- Programs **5–6**, freely assignable
- Neutral program **0**
- Mouth rinsing program **S**
(with special **memory function**).
- Freely assignable program **P**
(with special **memory function**).
- Examination position in the Neutral dialog. Actuate foot switch (at **A** or, if it has been exchanged, at **B**,

see page page 15). 

Run the activated program until the end position is reached*!

The selected program symbol appears with a bold frame now (example: program 1).


- Switch the unit lamp **ON** or **OFF**. This state is then also programmed.
- Actuate the foot switch.

The cursor • then jumps to the symbol for the **manual chair setting dialog**.

- Activate the manual chair setting dialog.

The program appears.

- Actuate the different setting keys to move the chair into the desired treatment position (see page page 29).
- Then bring the dentist element into the desired position by hand.

Observe the safety information "Chair and dentist element" on page 9!


- Place the cursor • on the programming symbol (if programming with the foot switch).
- Actuate the programming key until three successive acoustic signals sound.

Programming is completed now.

English

i NOTE

* If the program has not run up to the end position prior to resetting it (i.e. if it is aborted), the new treatment position cannot be programmed! This is indicated by the dimmed programming symbol.


When programming is completed, the **Chair program dialog** automatically appears on the display again. Preselect by activating **sitting** or **standing** dentist (see page 20).

Brief actuation of the desired program causes the treatment center to travel into the treatment position you selected. The unit lamp **then** switches **ON** (if programmed).

If the lamp is programmed **OFF**, e.g. in program **0, S, P** etc., the lamp switches off **before the treatment center moves** into the treatment position.

If the program is activated for a longer period of time, the

appears in the center of the odontogram and the treatment center returns to its factory-set position.

If the factory-set examination position is reprogrammed in the **Neutral dialog**, (foot switch at the front, A, or, if exchanged, at the rear, B) the factory setting is overwritten.

Programming key S




The tumbler filling and cuspidor flushing functions can be programmed at the same time.

- Actuate the **S** key in the Chair program dialog. If required, change the positions of the chair and the dentist element manually as described in the chair setting dialog (see page 29).
- In the Chair program dialog, select the assistant dialog.
- Activate automatic tumbler filling or / and cuspidor flushing.
- Call up the manual chair setting dialog and save it with the programming key.

CAUTION


IMPORTANT: *Tumbler filling and cuspidor flushing must still be activated while the program is saved.*

On selection of this tumbler filling position the also programmed functions of **tumbler filling / cuspidor flushing** switch on automatically.

12 Instrument dialogs

12.1 Instrument dialog: Electric motor


The **electric motor** instrument dialog appears when a motor has been removed from its receptacle.

If the motor is put back into the receptacle, the previously selected dialog appears again.

If several instruments are removed at the same time, only the one removed first will be ready for use.

Motor speed range:
200 to 40000rpm.

English

Function of the keys, left row:

Coolant ON / OFF (example: spray)

i NOTE


The key arrangement

Coolant ⇔ **Reversal of rotation** can be exchanged in the basic setting dialog 2 (see page page 58).

Switch

in Chair program dialog (see page page 19).

Reversal of rotation

Clockwise or counterclockwise rotation.

i NOTE

The key arrangement


Reversal of rotation ⇔ **Chip blower** can be exchanged in the basic setting dialog 2 (see page page 59).

i NOTE

The key arrangement

Reversal of rotation ⇔ **Coolant**

can be exchanged for each motor and for each user **A**, **B** or **C** in the basic setting dialog 2 (see page page 58).


Motor speed reduction to 200 rpm
In steps of 1000 from 40000 – 10000rpm
In steps of 400 from 10000 – 2000rpm
In steps of 200 from 2000 – 200rpm
starting slowly, then increasing rapidly.

Coolant, air or spray,
or **spray or NaCl**.

NaCl (with additional equipment saline solution) can be set alternatively **instead of air**. See basic setting dialog 1 (see page page 56).

Light on the dental burr **ON / OFF**


If this key is actuated longer, the **brightness** setting dialog appears.

Freely selectable function, ON / OFF.
e.g. call key (connected by your service engineer).

The function can be preselected as a key or switch in the basic setting dialog 1 (see page page 56).

Function of the keys, right row:

Chip blower (foot switch to the right)

An air jet comes out of the nozzle on the dental burr.
The light on the dental burr switches on. It remains on until approx. 10sec after the chip blower was switched off.

i NOTE

The key arrangement


Reversal of rotation ⇔ **Chip blower** can be exchanged in the basic setting dialog 2 (see page page 59).


Switch to the manual chair setting dialog (see page page 29).

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
Motor speed increase to 40000rpm.
In steps of 200 from 200 – 200rpm
In steps of 400 from 2000 – 10000rpm
In steps of 1000 from 10000 – 200rpm
starting slowly, then increasing rapidly.

Programming key (save)


Changeover switch for operation as direct starter or **speed controller foot switch**.

Unit lamp OFF / ON

If this key is pressed longer, the **The brightness** setting dialog appears.
(see page page 24).

Motor speed range:
200 to 40000rpm.

The motor is started by actuating the foot switch.


Operation as direct starter (key appears with a normal frame):

When actuating the foot switch, you will be working with the fixed value programmed for the motor.

Operation as speed controller (key appears with a bold frame):

An **expanding triangle** appears. It indicates the regulation of the speed during motor operation.
The foot switch operates as a speed controller within a range of 0 to the programmed speed.

Programming of functions

For each motor, different functions can be programmed for any user (**A**, **B** or **C**; see Neutral dialog **0** or assistant dialog, page 31 and 32) in the three function levels **E1**, **E2** and **E3**:

⚠ CAUTION

Prior to any change of user A, B or C, all instruments must be placed back in their receptacles!

- Remove motor from its receptacle.
The **electric motor** instrument dialog appears.
- Select the desired function level **E1**, **E2** or **E3**.
- Set the desired maximum motor speed.
- Select the desired functions by actuating the corresponding keys.

Example shown for user C:
(symbols with bold frame)

Coolant

Direction of rotation,
clockwise

Function level E1

max. motor speed
24000 rpm

Light on the motor

controller foot switch

- Actuate the programming key until three successive acoustic signals sound. Then put the motor back in its receptacle.

Programming for this motor in function level **E1** is completed now.

After the motor has been removed from the receptacle, the functions programmed in this function level will always be available for treatment.


However, they can be changed manually during a treatment.

In this case, the **E** symbol will appear with a normal frame.

When the functions are reprogrammed, the previous functions will be deleted automatically.

On selection of a key and subsequent actuation of the foot switch, the cursor jumps to the top again.

12.2 Instrument dialog: Turbine

The **turbine** instrument dialog appears when a turbine has been taken out of its receptacle.

When the turbine is put back into the receptacle, the previously selected dialog appears again.

If several turbines/instruments are removed at the same time, only the one removed first will be ready for use.

The turbine is started by actuating the foot switch.

The functions of the keys are identical with those in the **electric motor** instrument dialog.

However, the following functions are not available: Pre-selection of speed, reversal of rotation and changeover to speed controller foot switch.


Programming of functions

For each turbine, different functions can be programmed for any user (**A**, **B** or **C**, see Neutral dialog **0** or assistant dialog, page 31 and 32) in the three function levels **E1**, **E2** and **E3**:

- Select the desired function level **E1**, **E2** or **E3** for the turbine in use.
- Select the desired functions.

Example:

(symbols with bold frame)

Coolant**Function level E1****Light on the turbine**

- Actuate the programming key until three successive acoustic signals sound. Then put the turbine back in its receptacle.

Programming for this turbine in function level E1 is thus completed.


After the turbine has been removed from the receptacle, the functions programmed in this function level will always be available for treatment.

However, they can be changed manually during a treatment.

In this case, the **E** symbol will appear with a normal frame.

When the functions are reprogrammed, the previous functions will be deleted automatically.

12.3 Instrument dialog: SIROSONIC L

The **SIROSONIC L** instrument dialog appears when the ultrasound handpiece has been removed from its receptacle.

When the handpiece is put back into the receptacle, the previously selected dialog appears again.

The handpiece is started by actuating the foot switch.

The cooling water flow is always switched on automatically.

The functions of the displayed keys are identical with those in the other dialogs.

The intensity can be preselected from 1% to 100% with the - and + keys.

The **Endo** key is used to preselect endodontia. Intensity limited to max. 4%.

Programming of functions

For the ultrasound handpiece, different functions can be programmed for any user (**A**, **B** or **C**) in the three function levels **E1**, **E2** and **E3**:

- Select function level **E1**, **E2** or **E3** for the ultrasound handpiece in use.
- Select intensity in %.
- Select the light **ON / OFF** function.

Example: (symbols with bold frame)

Function level E1
Intensity 50%
Light on the hand-


- Actuate the programming key until three successive acoustic signals sound. Then put the handpiece back in its receptacle.

Programming is completed now.

After the handpiece has been removed from the receptacle, the functions programmed in this function level will always be available for treatment.

When the functions are reprogrammed, the previous functions will be deleted automatically.

12.4 Instrument dialog: SIROTOM

The **SIROTOM** instrument dialog appears when the electrosurgical handpiece has been removed from its receptacle.

The functions of the displayed keys are identical with those in the other instrument dialogs.

With the **lower** keys – and + the intensity can be preselected from 1 to 100.

With the **upper** keys – and + the type of current **COAG** can be preselected from **0** to **4**.
(Always press the gray keys!).

Always select COAG 0 for electrotomy!

For **coagulation** select a modulated type of current from **COAG 1** to **4** for the depth of crust based on experience.

i NOTE

The volume of the operating sound can be adjusted by your service engineer.

Programming of functions

For the electrosurgical handpiece, different functions can be programmed for any user (**A**, **B** or **C**) in the three function levels **E1**, **E2** and **E3**:


- Function level **E1**, **E2** or **E3** for the handpiece in use.
- Type of current:
Cutting current **0**
Coagulation current **1–4**
- Intensity.

Example:

Cutting current **0**
Function level **E1**
Intensity 80%

Actuate the programming key until three successive acoustic signals sound. Then put the handpiece back in its receptacle.

Programming is completed now.

Handling

The electrosurgical handpiece is used for electrotomy (cutting), coagulation and desiccation in biterminal technique.

⚠ CAUTION

The patient must always hold the neutral electrode in his hand during this procedure!

⚠ CAUTION

The connection socket is located on the bottom side of the dentist element.

The handpiece is activated by actuating the foot switch.

An intermittent acoustic signal then sounds and a

symbol for non-ionizing radiation is displayed.

⚠ CAUTION

For interference suppression reasons, the handpiece may be used only for a few seconds.


The HF switches off after the handpiece is put back in its receptacle. The programmed values of the previously selected function level appear again the next time the handpiece is removed from its receptacle.

⚠ CAUTION

Incorrect operation and the nonobservance of precautionary measures can cause serious accidents when working with the electrosurgical unit.

⚠ WARNING

There is a risk for patients with cardiac pacemakers that the pacemaker can be disturbed during treatment with the electrosurgical unit. We recommend not using the electrosurgical unit for these patients.


Further WARNINGS for working with the surgical unit

- The patient should not come into contact with metallic parts which are grounded or have substantial capacitance to ground (e.g. operating table, holders).
- The power output should be set to the lowest possible value for the relevant purpose.
- If the surgical unit appears to deliver little power or does not work properly in its normal setting, this may be caused by a poor contact in the supply cable.

- Combustible substances, which are used for instance as cleaning or disinfecting agents, should have evaporated before applying surgery. Cotton wool can be ignited. Endogenous gases can be ignited.
- The HF can impair the function of other electronic devices.
- Check the electrode cable regularly for possible damage to the insulation.
- Prevent contact with metallic parts in or on the body.

SIROTOM technical description


Power characteristic

for cutting and coagulation

Power measured between handpiece and protective ground wire.


English

Power curve

for cutting and coagulation

Power measured between handpiece and protective ground wire.

Output voltage

as a function of the degree of coagulation, the peak voltage increases.

Technical data for SIROTOM

Maximum output voltage peak – peak between hand-piece and protective ground wire

When cutting: 1400V

At max. coagulation: 1600V

Modulation frequency: 1667Hz

Frequency of operating and alarm tone: 1200Hz

The neutral electrode is connected to the protective ground wire via a capacitor.

Connection for the neutral electrode.
Neutral electrode at high frequency relative to ground.

Type **BF** SIROTOM application unit

DAB 25% ED 10s SD 40s
1 MHz / 50 W
300 - 600 Ohm

Intermittent operation 25%

Duty time 10s

Cycle time 40s

HF frequency 1MHz

Output power max. 50W across 300 ohms w/ Koag 0
350 ohms w/ Koag 1
400 ohms w/ Koag 2
500 ohms w/ Koag 3
600 ohms w/ Koag 4


Symbol for non-ionizing radiation

12.5 Safety checks


In Germany, medical products are subject to the Ordinance on the Installation, Operation and Use of Medical Devices (MPBetreibV). Safety checks must be performed and a "Medical Product Log" must be kept. Please refer to Section 22.3 of "Safety Checks" for details.

12.6 Switching over from instrument dialog to Chair program dialog or to manual chair setting dialog


If you have removed an instrument and then switch over from the instrument dialog to the Chair program dialog:

Instead of the chair symbol (right), the instrument symbol for resetting to the instrument dialog appears.


On switching to the manual chair setting dialog:
Instead of the odontogram (left), the instrument symbol for resetting to the instrument dialog appears.

13 Water unit, assistant element

13.1 Swiveling the cupid

The cupid on the water unit can be swiveled manually by approx. 110mm in the direction of the patient chair.


Collision protection

All chair movements are locked with the cupid swiveled in. If a programmed movement is started, then the cupid first moves back into its starting position before the chair moves into the programmed position. When a key for the individual adjustment of the chair is actuated, the cupid also swivels back first and, if the corresponding key is still pressed then, the chair will move. If the cupid is swung in during the chair movement, then the program aborts.

13.2 Adjusting the water quantity

The control for the water amount for the circular flushing is adjusted by the service engineer during installation.

13.3 Amalgam rotor

The indication **Amalgam** appears if the amalgam rotor must be exchanged (see maintenance instructions).


A **continuous sound** indicates that the rotor is full and no longer functioning properly.

An **interval tone** indicates an error. The water supply is interrupted.

Switch the unit **OFF** and then **ON** again. If the error still appears even after the system has been switched **OFF** and **ON** several times, please call your service engineer.

NOTE: The rotor inside the amalgam separator must be replaced at least once a year, even if the **Amalgam** display does not light up; otherwise, unpleasant odors can develop.

13.4 Disinfection system

The disinfection system adds a disinfectant to the water (1:100) to prevent the formation of microorganisms in the water system.


Prior to initial operation

of your treatment center, **sanitation** must be performed.

If, on the basis of an agreement with you, sanitation was not performed by the service engineer following installation of your treatment center, please perform sanitation as described in the separate instructions "**Care and Cleaning by the Practice Team**".

Sanitation takes approx. 24 hours.

13.5 During treatment

If during treatment **Desinf.** flashes on the display, top up with 1 liter of **DENTOSEPT P** before beginning treatment of the next patient:

Remove the cover behind the cuspidor and fill in 1 liter of **DENTOSEPT P**.

The **filler trough** should always be kept **clean**.

Make sure that **no tap water** enters the fill-in trough!

DENTOSEPT P

1 carton = 6x1 liter, Order No. 33 18 156

Microbiological water test

(after a working pause of more than **1 week**)

Sanitation of the unit

Perform sanitation **every 4 weeks**. Only in this way can you ensure that biofilm is effectively controlled in the water paths.


(see instructions "**Care and Cleaning by the Practice Team**").

CAUTION

In Japan the disinfectant **OXYDOL** available there must be used instead of **DENTOSEPT P**.

OXYDOL must be mixed 1:1 with distilled water before it is filled into the disinfection system.

13.6 Suction handpieces on the assistant element

The vacuum system is switched on after a suction hose is removed from the receptacle.

CAUTION

Observe the warning and safety information on page 9 (vacuum system).

The suction handpiece can be angled by turning it.

With removed suction hose, the vacuum system can be switched **OFF-ON** by actuating the foot control of the chair in any direction (if preselected: page 58).

The suction power can be adjusted in the **Basic Setting Dialog 1** (page 56).

This does not apply for water units with wet aspiration.

CAUTION

With wet aspiration, never pour water into the cuspidor with the unit switched off.

CAUTION

Should the vacuum system switch off unintentionally, immediately remove the suction cannula from the patient's mouth.


If the unit has no second saliva ejector hose for the surgical cannula, the thick suction hose must be used for surgical aspiration.

To insert a surgical cannula, please attach the adapter supplied.


If the unit is equipped with a second saliva ejector hose, the surgical cannula can be inserted directly into the saliva ejector handpiece.

CAUTION

Make sure to hold the suction cannula in such a way that the cannula opening cannot attach itself to the oral mucous membranes by accident.

13.7 SPRAYVIT

SPRAYVIT is a multifunctional syringe.

When the pushbuttons on the handpiece are actuated, the light switches on (if preselected).

The light on the removed handpiece can also be switched on without actuating the **SPRAYVIT** pushbuttons by pressing the foot switch in the A (or B) direction.


The light switches off when the handpiece is placed in the receptacle.

A built-in heater heats up the water and air (if preselected).

WARNING

*To avoid excessively high water temperatures causing possible scalds, clean the nozzle regularly. See **SPRAYVIT Operating Instructions**.*

When operating the syringe ensure that the air and water outlets at the nozzle are not blocked in any way.

If you activate the key in the assistant dialog

and then
in the Service dialog the key


the basic setting dialog 1 appears. The **SPRAYVIT** light key is located here on the top left, provided that a **SPRAYVIT** is connected.

The **SPRAYVIT** light can be selected or deselected by briefly actuating this key.

If the light is preselected, the key appears with a bold frame.

The preselection is performed separately for the dentist and assistant **SPRAYVIT** on the corresponding panel.


13.8 SPRAYVIT setting dialog

Each **SPRAYVIT** is preset separately on the corresponding panel for assistant or dentist.

It is not possible to set the assistant **SPRAYVIT** at the dentist panel and vice versa.

After the **SPRAYVIT** with **light** key is actuated for a longer time in the basic settings dialog 1, the **SPRAYVIT** settings dialog appears.

The light intensity of the **SPRAYVIT** can be adjusted with the – and + keys in the top level.

The temperature setting for the **SPRAYVIT** water heater can be adjusted with the – and + keys in the third level,


and for the air heater with the – and + keys in the fourth level.

You can switch the heaters for water and air **ON/OFF** with the **thermometer** key at the bottom left.

The key appears with a bold frame if heating is pre-selected.

When the presettings have been completed, the display changes back to the previous dialog after a short time. You can also switch back with the **arrow** key.

14 Service dialog


The Service dialog appears after actuation of the key in the assistant dialog.

English


14.1 Function of the keys, left row:


- **Contrast setting**
Liquid crystal display
- **Back to previous dialog**
- **Hardware and software status**
of the treatment center
(only for your service engineer)
- **Switch** to the
basic setting dialog 1.
- **Time display ON / OFF**
in the footer of the displays.
If this key is pressed for some time, the setting dialog for the time and date appears, see page page 52.

14.2 Function of the keys, right row:

-  **Contrast setting**
Liquid crystal display
-  **Service information**
(only for your service engineer)
-  **Switch** to the
basic setting dialog **2**.
-  **Configuration** of the treatment center
(only for your service engineer).

14.3 Setting the time and date


The data set are displayed in the top line.

h	-	hours
min	-	minutes
day	-	Day
mon	-	Month
year	-	Year

When the setting dialog has disappeared or has been reset, the time appears in the footer of the display, if the clock symbol was preselected in the service dialog. The date is not displayed. It is required for internal treatment center functions, e.g. for error detection.

15 Basic setting dialog 1


The basic setting dialog 1 appears after actuation of the key in the assistant dialog and subsequent actuation of the key in the service dialog.

15.1 Function of the keys, left row:


SPRAYVIT light ON/ OFF

This symbol appears only if a **SPRAYVIT** is connected. If the key is actuated longer, the **SPRAYVIT setting dialog** (see page 50) appears.

Back to previous dialog


Call key: Function as switch / key. Depends on the activated function.

Connected and set by your service engineer during installation of the unit.

Settings may only be changed by the service engineer!

Only possible when preselected with the "call key". With "SIVISION" preselected, this key disappears.

Preselecting call key or SIVISION 1

With "call key" preselected, the switch/key preselect symbol is displayed and this preselection is possible. With "SIVISION 1" preselected, the switch/key preselect symbol disappears.

The key function is then automatically set.


With C1+ always leave the preselection on the "call key" symbol.

Amalgam display test

Actuate key:


**Amalgamflashes on the display
and a beeper sounds!**

**Restoring the SIROCAM 3 factory color settings**

When the symbol is marked with a bold frame, the factory settings are active. When the symbol frame is normal, the settings are active once the white balance has been made.

***i* NOTE**


The key is displayed only if the camera is removed and active.

**White balance for SIROCAM 3**

A white balance is possible only if the camera is housed in the dentist element or in the additional instrument holder.

***i* NOTE**

The key is displayed only if the camera is removed and active.

The camera handpiece contains data making it possible to restore the factory-set color definitions regardless of the dental treatment center where the camera is used.

This makes it possible to use the camera at different dental treatment centers.

Perform a white balance for adaptations, e.g. for individually wanted color settings.

For this purpose the sliding switch (2) must be set in the middle position (survey exposure).

The camera is then held at a distance of 2 – 5 cm in the live mode with the lens over the color sample sheet A (Order No. 59 12 956) delivered with the dental treatment center or over another color sample sheet.

The complete screen may in this case display **only** the aimed at color surface.

Different color patterns will produce different results with respect to the color settings.

For example, a color pattern with blue parts will reduce the amount of blue color in the video image.

Now activate the "white balance" key.

Once the white balance is finished (may take several seconds), an acoustic signal is emitted and the symbol frame appears bold.


If the new color setting does not correspond with your expectations, then you can restore the factory setting as described or repeat the balance.

15.2 Function of keys, right row:

Adjustment of the unit lamp

Setting of the lamp reflector voltage adaptation value, e.g. 3 (may be set only by the service engineer).

Call key:

Function as switch / button.

Depends on the activated function. Connected and set by your service engineer during installation of the unit.

The setting may only be changed by your service engineer!

Setting the suction power:

The setting can be made separately for the suction hose and saliva ejector hose. If both suction handpieces are used at the same time, the suction power of the suction hose always has priority.



Setting: Remove the suction or saliva ejector hose.

A bold frame around the symbol indicates that **maximum suction power** is set.

If the suction power needs to be **reduced**, actuate the key until the known setting dialog appears.

In this program, the suction power can be adjusted in 3 steps by actuating the **-/+** keys (symbol with normal frame).

When the hose is put back in its receptacle, the suction power setting is programmed.

Surgery (NaCl) preselection keys


(with additional equipment NaCl cooling)

e.g. with instrument in receptacle 3

0 = without saline solution


2 - 6 are receptacles on the dentist element from left to right (see page page 76). All instruments must be deposited in their receptacles!

16 Basic setting dialog 2


The basic setting dialog 2 appears after actuation of the key in the assistant dialog and subsequent actuation of the key in the service dialog.


16.1 Function of keys, left row:

Volume setting "softer" for all signal tones on the dentist panel or on the assistant panel (adjustable separately), e.g. timer signal, programming signal etc.

Back to previous dialog

Cursor • ON / OFF (dentist element only)


CAUTION

When the cursor field is blank, the display cannot be controlled via the foot switch.

Key / cursor “clicking tone” ON / OFF


Acoustic confirmation of key and cursor button actuation.

Key pressure sensitivity

Weaker key pressure

English


Foot control functions

Preselections are made in the following order:

1 Suction and POLYLUX (in the assistant element) can be switched OFF / ON by actuating the 4-way foot control in any direction. The first handpiece removed has priority.

With the handpieces in their receptacles, the chair functions can be selected.

- 2 Only POLYLUX (in the assistant element) and chair functions possible.
- 3 Only chair functions possible.
- 4 Only suction and chair functions possible.


ON / OFF switch for different cursor navigation tones. A different tone pitch is assigned to each level.

i NOTE


The cursor "clicking tone" must be switched ON.

16.2 Function of the keys, right row:

Volume setting "louder"

for all signal tones on the dentist panel or on the assistant panel (adjustable separately), e.g. timer signal, programming signal etc.

Preselection of key position

Reversal of rotation ⇔ Coolant ON / OFF



on the display of the **electric motor instrument dialog**. The symbol displayed bold then appears at the top left in the **electric motor** instrument dialog (page 35).

i NOTE

The symbol in the basic setting dialog 2 does not appear until the motor is removed from its receptacle.

Preselection is thus possible for each individual electric motor.


This applies when the motor is removed for the corresponding preselected user **A**, **B** or **C** and is the same in all three program levels **E1**, **E2** and **E3**.

**Preselection of key position****Reversal of rotation ⇔ Chip blower**on the display of the **electric motor instrument dialog**.
The symbol displayed bold then appears at the top left in the **electric motor** instrument dialog (see page page 35).**NOTE***The symbol in the basic setting dialog 2 does not appear until the motor is removed from its receptacle.*
Preselection is thus possible for each individual electric motor.*This applies when the motor is removed for the corresponding preselected user **A**, **B** or **C** and is the same in all three program levels **E1**, **E2** and **E3**.***Foot control assignment**

for right/left switch position

Preselection of backrest adjustment

or headrest adjustment (see page page 30)

Key
28**Key pressure sensitivity**


Stronger key pressure

**Switching on the PC dialog in the assistant dialog**

only available if no camera is integrated in the dentist element (see page page 70).

17 SIVISION 3 / SIROCAM 3 / SIROCAM C

17.1 SIROCAM 3 Intraoral Camera (additional equipment)

The SIROCAM 3 intraoral camera in combination with SIVISION 3 represents an outstandingly suitable system for effective communication with the patient.

The camera images are not used for diagnosis!

The camera is a sensitive optical instrument and must therefore always be handled with care.

The camera illuminates the parts of the teeth to be exposed with 7 white-light LEDs (1).

These are arranged around the entry prism (2).

The LEDs and optical system are protected by a scratch-proof antiglare glass cover.

CAUTION

* = Intermittent operation of the camera.

The camera heats up during operation.

To avoid high temperatures of the camera this must be set down after an operating time t1 for a period t2.

Please refer to the operating instructions of the SIROCAM 3 for the times for t1 and t2.

Even when set down the camera heats up, because one of the statuses of black image, white image or quad image is always active.

The slide switch (3) is used for presetting the depth of field.


The camera key (4) is intended for controlling different SIVISION 3 functions.

NOTE

Operation of the camera and further information are described in the Operating Instructions for SIROCAM 3.

NOTE

The camera handpiece contains data making it possible to restore the factory-set color definitions regardless of the dental treatment center where the camera is used. This makes it possible to use the camera at different dental treatment centers. To adapt the individually wanted color settings, a white balance can be selected with the camera active (see basic setting dialog 1, chapter 15)

**Hygienic protective covers for SIROCAM 3**

Hygienic protective covers are available for the camera. They can be ordered under **REF 59 38 654**. Packaging unit: 500 pcs.

How to use the hygienic protective covers:

1. Hold the hygienic protective cover with its paper side face-down and insert the camera with the optical system downward into the cover.
2. Pull off the upper protective foil with the blue stripe.
3. Then pull off the backing paper.
4. Finally, stretch the hygienic protective cover over the camera from front to rear.

⚠ CAUTION


In order to achieve optimum image quality, there must not be any folds in the transparent area of the cover above the lens window.

⚠ CAUTION

Use only the hygienic protective covers offered by Sirona.

Prior to each application on a new patient, the camera must be covered with a new hygienic protective cover.

17.2 SIROCAM C Intraoral Camera (additional equipment)

The SIROCAM C intraoral camera in combination with SIVISION 3 represents an outstandingly suitable system for effective communication with the patient.

The camera images are not used for diagnosis!

The SIROCAM C intraoral camera is a sensitive optical instrument and must therefore always be handled with care.

The camera illuminates the parts of the teeth to be exposed with 6 white-light LEDs (1).

These are arranged around the single-lens optical system (2).

The LEDs and optical system are protected by a scratch-proof antiglare glass cover.

CAUTION

* = Intermittent operation of the camera:

t1 = 1 min ON, **t2** = 3min OFF.

The camera heats up during operation.


To avoid high temperatures of the camera this must be set down after an operating time **t1** = 1 min for a period **t2** = 3 min.

Even when set down the camera heats up, because one of the statuses of black image, white image or quad image is always active.

When the camera is in LIVE IMAGE mode, it is possible to toggle between close range and far range by pressing button (3).

Button released – close range **5 – 15mm**

Button pressed – far range **from 15mm** to Smile line.

**Hygienic protective covers for SIROCAM C**

Hygienic protective covers are available for the camera. They can be ordered under **REF 59 14,705**. Packaging unit: 500 pcs.

How to use the hygienic protective covers:

5. Hold the hygienic protective cover with its paper side face-down and insert the camera with the optical system downward into the cover.
6. Pull off the upper protective foil with the blue stripe.
7. Then pull off the backing paper.
8. Finally, stretch the hygienic protective cover over the camera from front to rear.

⚠ CAUTION

In order to achieve optimum image quality, there must not be any folds in the transparent area of the cover above the lens window.

⚠ CAUTION

Use only the hygienic protective covers offered by Sirona.

Prior to each application on a new patient, the camera must be covered with a new hygienic protective cover.

17.3 Camera dialog: General notes and SIVISION 3 operating modes

The camera dialog of the dentist panel opens after the SIROCAM 3 / SIROCAM C3 is removed from instrument holder 6 or from the additional instrument holder.

The camera can be operated by the unit foot switch if no other instrument is active (removed and operable by foot switch) when the camera is removed.

i NOTE

The camera dialog can also be called up:

- with the left dialog change key
of the assistant dialog,*

regardless of whether the camera is removed.

This is possible on the dentist panel and the assistant panel (see page page 22).

The camera dialog does not close when the camera is deposited. Even with the camera deposited, various camera and PC functions can be used.

It is closed only after the foot switch is pressed forward or the dialog is terminated in any way. No image changes occur in this case.


The camera dialog is dependent on the **operating mode**.

- In **Video mode**, the video functions can be used independently of a PC. Up to four images can be saved.
- In the **PC mode**, the video function depends on a connected PC and the software installed on it (e.g. SIDEXIS/VIDEXIS).

17.4 Camera dialog: Video mode

Removing the camera

- If **STILL IMAGE** or **QUAD IMAGE** was displayed with the camera set down in its holder, this display is retained.
- If **BLACK IMAGE** or **WHITE IMAGE** was displayed with the camera set down in its holder, the preceding display status corresponding to the BLACK or WHITE IMAGE appears.
LIVE, STILL or QUAD IMAGE.

Switching from STILL IMAGE to LIVE IMAGE


Actuate foot switch in direction **A*** – cursor upwards.

i NOTE

** The foot switch functions A and B can be reversed on request by the service engineer.*

Depositing the camera


- If a **LIVE IMAGE** was displayed, **BLACK IMAGE** appears.
- **STILL, QUAD, BLACK** and **WHITE IMAGE** remain displayed.

Function of the keys, left row:

English

IMAGE SWITCHING


- In **SINGLE IMAGE** mode, switching to the next image storage location.
If a **STILL IMAGE** is stored there, it will be displayed.
If no **STILL IMAGE** is stored there, a **LIVE IMAGE** will appear with the camera removed.
- In **QUAD IMAGE** mode, switching to storage number **I, II, III** or **IIII** for image selection.
The currently selected image number lights up in the upper right corner of the image.

- Switch to the
Chair program dialog (see page 19).

When the camera is deposited, the camera dialog is also closed.


Clearing storage locations 1–4

- When the camera is deposited, a black image appears.
- When the camera is removed, a **LIVE IMAGE** appears for image **I**.


Switching to PC mode

The image from the PC, e.g. SIDEXIS/VIDEXIS, will then appear on the SIVISION monitor.
The camera switches to the **LIVE IMAGE**.
This mode remains active even if the treatment center is switched OFF and back ON.


Switching to WHITE IMAGE

Freely selectable function, ON / OFF.
e.g. call key (connected by your service engineer).

The function can be preselected as a key or switch in the basic setting dialog 1 (see page page 56).

Function of the keys, right row:

Switching between QUAD IMAGE and SINGLE IMAGE.

When switching to **SINGLE IMAGE**, the image previously selected in the **QUAD IMAGE** is displayed.

Switch to the
manual chair setting dialog

When the camera is deposited, the camera dialog is also closed.

Switch over to BLACK IMAGE


Unit lamp

can be switched **ON / OFF** as required.

i NOTE

In video mode, the monitor is automatically switched off when the electrosurgical handpiece is activated and then switched on again when the handpiece is put back in its holder.

17.5 Camera dialog: PC mode

Removing the camera


The **LIVE IMAGE** window opens.

Switching from STILL IMAGE to LIVE IMAGE

Actuate foot switch in direction **A*** – cursor upwards.

i NOTE

* The foot switch functions **A** and **B** can be reversed on request by the service engineer.

Depositing the camera

The **LIVE IMAGE** window closes.

Displayed **STILL IMAGES** or stored images remain open.

Function keys F1–F12 are available for controlling PC functions.

They can be used e.g. to control SIDEXIS/VIDEXIS functions such as:


- Tile,
- Next Image,
- Full Frame,
- Zoom In,
- Zoom Out,
- Rotate Image.

Other PC applications can also be controlled with the function keys.

i NOTE

The factory settings and the configuration procedure are described in the "**SIVISION 1/2/3 – Operator's Manual for PC Software**".

Function of the keys, left row:

Save image

Switch to the Chair program dialog (see page 19).


When the camera is deposited, the camera dialog is also closed

Function key


Function key

Switching to VIDEO mode

The **SIVISION** monitor then displays the **LIVE IMAGE** of the camera.

This mode remains active even if the treatment center is switched OFF and back ON.

Switching to WHITE IMAGE


Freely selectable function, ON / OFF.

e.g. call key (connected by your service engineer).

The function can be preselected as a key or switch in the basic setting dialog 1 (see page 56).

This key can also be used as a PC function key.

Function of the keys, right row:

Tiled view

All images opened in SIDEXIS/VIDEXIS are displayed side by side.


Switch to the manual chair setting dialog.

When the camera is deposited, the camera dialog is also closed.

Function key

-  **Function key**
-  **Function key**
-  **Switch over to BLACK IMAGE**
-  **Unit lamp**
can be switched **ON / OFF** as required.

i NOTE


In PC mode and with white or black image switched on, the monitor is automatically switched off when the electrosurgical handpiece is activated and then switched on again when the handpiece is put back in its holder. Interferences can occur in the PC image because of the active electrosurgical function.

17.6 SIVISION 3 – Second monitor function without SIROCAM 3 / SIROCAM C (PC dialog)

The PC dialog can be called up with the left switch key of the assistant dialog.

 (Activation via basic setting dialog 2, see page page 57)

The monitor integrated in the treatment center serves as a primary or secondary monitor for the externally connected PC.

Function keys F1–F12 are available for controlling PC functions.

They can be used e.g. to control SIDEXIS/VIDEXIS functions such as:


- Tile,
- Next Image,
- Full Frame,
- Zoom In,
- Zoom Out,
- Rotate Image.

Other PC applications can also be controlled with the function keys.

NOTE

The factory settings and the configuration procedure are described in the “**SIVISION 1/2/3 – Operator’s Manual-for PC Software**”.

18 Simplified dialogs


On the displays of dentist and assistant element there is a possibility for the

- Chair program dialog
- assistant dialog
- Instrument dialog: Electric motor

to switch over to a **reduced display mode**.

i NOTE

The display mode must be set separately for all three dialogs both at the dentist and at the assistant element.

When the **electric motor** instrument dialog is set to reduced display mode, all instrument dialog displays are automatically reduced as well.

With the reduced display mode set in the assistant dialog, the displays in the neutral dialog and in the service dialog are also automatically reduced.

In order to work with reduced dialog displays, call up the corresponding dialog in the **extended display mode**.

In each case, actuate the second key from the top right until the reduced display mode appears (after approx. 10sec).


The presettings for the user **A**, **B** or **C** and the instrument levels **E1**, **E2** or **E3** selected in the extended display mode are used.

Dialogs that are reprogrammed in the reduced display mode and then saved with the corresponding programming key will also be used by the extended display mode when this is reset.

The reduced display mode remains activated when the treatment center is switched **OFF** and **ON** again.


English

18.1 Simplified Chair program dialog

The simplified display mode uses the chair settings programmed for user **A**, **B** or **C** selected in the extended display mode.

The factory settings can be restored also in this display mode.

However, programs 1–4 are not assigned to an odontogram.


If programs are reset in this dialog and saved in the manual chair setting dialog with the programming key, these programs will also be used by the extended display mode when this is activated.

The keys

- **standing / sitting user** and
- **Program P** (memory function)

are not displayed.

18.2 Simplified assistant dialog

The simplified display mode does not have the following keys:

- **Preselection of user A, B or C**,
- **San.** for sanitation and
- composite key **C** for unit lamp.


The user's height can be set only after switching to the extended display mode.

If sanitation of the treatment center is required, the assistant dialog must also be switched over to the extended display mode.

Switching over to the reduced display mode in the assistant dialog also activates the reduced display mode in the

Neutral dialog **0** and in the Service dialog,
see next page.

18.3 Simplified Neutral dialog


Like in the assistant dialog, the simplified display mode also does not have the key for

- **Preselection of user A, B or C.**

The user's height can be selected only after switching to the extended display mode.

English

18.4 Simplified Service dialog


The reduced display mode for the service dialog can be selected in the reduced assistant dialog.

The following keys are missing:

- **Service information**
(only for your service engineer) and
- **Switch to the basic setting dialog.**

To change basic settings, the extended display mode must be selected in the assistant dialog.

18.5 Simplified instrument dialogs

(Example: **Electric motor**)

The user preselected in the extended display mode, i.e. **A, B or C**, is automatically used when switching over to the reduced display mode.

The same applies to the programmed instrument function level selected for this user, i.e. **E1, E2 or E3**.

The following preselection and toggle keys are missing:

- Coolant, spray or air
- Light on the dental burr ON / OFF
- Instrument function level
E1, E2, E3
- Direct starter or speed controller foot switch

Switching over to the reduced display mode in the **electric motor** instrument dialog automatically reduces the display of all other instrument dialogs as well.

The other two instrument dialogs do not have the preselection key for

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nor the following keys, if available in the extended display mode:


- Light on handpiece ON / OFF
- Coolant, spray or air

Using the programming key, the remaining key functions can also be changed in the reduced display mode.

The reprogrammed functions will then also be used in the extended display mode when this is reset.

19 Surgery with saline solution (NaCl)

19.1 Preparations on the unit

- Hook NaCl bottle (1) in place.
- Attach the peristaltic pump set (2).
- Push short end of hose (3) with cannula as shown through stopper and into NaCl bottle.
- The regulator in the hose clip (4) must be in the top position (completely open).
- Run the long hose (5) along the corresponding motor hose up to the right-angle handpiece and fasten with the clips (6).
- Fit coupling (7) onto hose and connect it with the thin silicone hose (8).
- Connect spray clip (9) with the thin hose and attach to angle piece.

The **peristaltic pump set (2)** is a disposable article and can be purchased as consumable material under the article number **F 58707** directly from the manufacturer in packs of 10 pcs.

Ordering address:

Satelec


Industriestrasse 9, D-40822 Mettmann

English

CAUTION

The flow of the pump must be at least 70 ml/min for reasons of safe cooling and to prevent a pressure rise in the hose (risk of bursting). This is not always the case when non-Sirona ultrasound tips are used.

19.2 Presetting


If corresponding NaCl pumps are fitted to the dentist element, the possibility of selecting saline solution NaCl for surgery can be preset for a maximum of 2 burr instrument receptacles (including SIROSONIC L).

For this purpose actuate the key in the assistant dialog.

The Service dialog is displayed.

In the Service dialog, actuate the key.

The basic setting dialog **1** appears.

One or two receptacles for surgery can be preset with the two keys at the bottom right.

All instruments must be deposited in their receptacles!

If only one receptacle should be equipped with NaCl cooling, then "0" must be set on one key and the number of the corresponding receptacle (e.g. 3) on the other key.


With a second NaCl pump, the number of the corresponding receptacle must be set on the lower key.

Receptacles **2–6** and **0** (without NaCl) can be selected.

The basic setting dialog **1** and the Service dialog disappear automatically after some time.

However, you can also manually exit either dialog by pressing the left arrow key.

19.3 Operation


When an **electric motor** with surgery (NaCl) presetting is removed, the **electric motor** instrument dialog appears.

With the key at the top left, the coolant can be switched **ON** (bold frame) or **OFF**.

Spray cooling or NaCl cooling can be preselected with the fifth key from above on the left.

The symbol of the preselected medium then appears on the coolant **ON / OFF** key.

It is not possible to preselect air.


If the coolant **ON / OFF** key with NaCl symbol is actuated longer, then the setting dialog for the intensity adjustment of the saline solution pump appears.

The intensity can be adjusted in 15 steps.

If two receptacles are used, this adjustment applies for both receptacles.

20 POLYLUX fiber optic light handpiece

The **POLYLUX light** can be installed in the assistant or dentist element.


Insert the light guide into the bushing on the bottom of the assistant or dentist element until it engages.

CAUTION

Always grasp the light guide by its anti-kink sleeve (1), and never by its tubing!

*Do not squeeze, twist or kink the light guide!
Observe a minimum bending radius of 80mm.*

20.1 The handpiece

Screw the handle sleeve with fiber optic rod onto the handpiece clockwise.


The handpiece can be turned.

Before making the composite filling, reduce the brightness of the lamp for the working area in order to avoid premature hardening of the filling material.

On the unit model this is done by actuating composite key **C**.

The intensity of the lamp cannot be changed while composite key **C** is activated.

20.2 Switching on with POLYLUX in the assistant element


The assistant display must show the assistant dialog.
Remove the handpiece from its receptacle. (The **POLYLUX light** is switched on via the timer function; automatic time limit: max. 60 seconds).

Set the time of the timer as described below and switch on the **POLYLUX light** with the **Start** timer key.

When switching over to the assistant dialog on the dentist panel, the stored time appears on the timer key as soon as the assistant **POLYLUX** handpiece is removed, and the dentist also has access to the **POLYLUX light** of the assistant.

In addition, the **POLYLUX light** can also be switched ON and OFF with the chair foot control (for configuration, see chapter 16.1).

While the **POLYLUX- light** is on, a short acoustic signal sounds every 20 seconds.

Preselection of timer setting

Three individually adjustable timer settings can be selected
(factory setting: 10, 30, 60 sec).

The settings are selected by repeatedly touching the preselection key lightly and briefly.

If this key is actuated longer, the setting dialog appears.

In this program the times can be changed by pressing the - / + keys.

Setting range continuous 1–60sec.

In the setting dialog, the switch the **soft start function ON / OFF**.

i NOTE

With soft start, the output power is started at 50% and increased to 100% in 10% steps within 8 seconds.


Soft start switched ON

Soft start switched OFF

20.3 Switching on with POLYLUX light in the dentist element

Remove the handpiece from its receptacle.

The **POLYLUX light** instrument dialog appears on the display.

Switch on the **POLYLUX light** with the foot switch.

With the two

- and + an exposure time of up to max. **60s** can be pre-selected.

This is indicated in seconds in the center of the display.

From the assistant panel the assistant also has access to the **POLYLUX light** of the dentist via the instrument dialog.

While the **POLYLUX- light** is on, a short acoustic signal sounds every 20 seconds.

The **POLYLUX light** can be switched **OFF / ON** with the handpiece removed at any time also by actuating the foot control on the patient chair in any direction (if pre-selected, see page 58).

After the time has elapsed, the light goes out automatically and the cooling fan runs on for some time.


Switching soft start ON / OFF

Individually programmable for each user.

i NOTE

With soft start, the output power is started at 50% and increased to 100% in 10% steps within 8 seconds.

20.4 Programming the exposure time


For the **POLYLUX light** in the dentist element, different exposure times can be programmed for any user **A**, **B** or **C** in the three function levels **E1**, **E2** and **E3**:

- Select function level **E1**, **E2** or **E3** for the handpiece removed.
- Set the exposure time.
- Actuate the programming key until three successive acoustic signals sound. Then put the handpiece back in its receptacle.

Programming is completed now.

20.5 Application


The POLYLUX is suitable for the polymerization of composites that harden under the influence of light.

The hardening time recommended by the manufacturer of the material should be observed.

An initial value can be determined by curing a test filling.

The action of light in the oral cavity must be limited to the area requiring clinical treatment.

The light exit surface of the handpiece must be brought as close as possible to the material to be polymerized.

CAUTION

The light exit surface must not contact the tooth filling at the beginning of the hardening process.

A built-in thermal switch protects the reflector lamp against overheating.


This switch may respond following several successive long-term exposure periods.

The lamp can be switched on again as soon as it has cooled down to its operating temperature.

When the filling has hardened, readjust the brightness of the unit lamp with the composite key in the assistant program or with the **ON / OFF** key of the lamp.

WARNING

*Do not direct the light into the eyes of the patient.
Do not look directly into the light.*

We recommend using a suitable glare shield. It can be simply slipped onto the light guide.


Example: glare shield, Order No.: 57311916

Ordering address:

Heraeus Med GmbH
Heraeusstr. 12 – 14
D-63450 Hanau, Germany

21 Additional equipment and accessories

21.1 Media block on the patient chair

Additional equipment

Only approved additional devices (e.g. supplied by EMS) may be connected to a media block on the base of the lower part of the chair (the manufacturer's instructions must be complied with).

The media block is equipped with quick couplings for air and water and an outlet for non-heating appliances with a separate fuse (6.3A).

For treatment centers with a disinfection unit, the integration of a media block complies with the DVGW approval (EN 1717/ DIN 1988). In this case, the additional devices connected are also compliant with the requirements of EN1717 / DIN1988.

- All connected devices must be supplied with disinfected water if a disinfection unit is installed.
0.1% -0.2% hydrogen peroxide are required for regular operation and 1.4% for sanitization.

CAUTION


Before connecting any additional devices, first check whether they are suitable for hydrogen peroxide. If these devices are not suitable for the higher concentration used for sanitization, Sirona recommends separating them from the treatment center during sanitization and then sanitizing them separately according to the manufacturer's instructions.

Sirona shall not be held liable for any damage to connected devices!

Pressure and flow rate specifications

	Pressure	Flow rate
Water	$2,2 \pm 0,2$ bar	max. 300 ml/min
Air	$4,4 \pm 0,5$ bar	max. 70 NI /min

21.2 Hygienic headrest protection


Accessories

Two different hygienic covers are available for the headrest of the patient chair:

- Washable textile headrest cover with Velcro closure and neck pad
Order No.: 59 05 588 (5 pcs.)
- Disposable headrest cover with a self-adhesive fastening
Order No.: 58 72 820 (5 pcs)

Slip headrest cover over headrest and fasten it as illustrated.

21.3 Seat cushion C

Additional equipment

The back part is fastened with a Velcro closure and can be removed.

Covering material: 100% cotton.

The covers can be removed after opening the Velcro strap on the back.

A Teflon coating is applied to the covers.

Instructions for care:


The Teflon coating is dirt-repellent and slightly water-repellent; remove any dirt with lukewarm soap suds (fine washing powder).

CAUTION

Removing stains with stain remover is not suitable, since this damages the Teflon coating!

This comfortable seat is simply placed on top of the upholstery of the patient chair. Thus it is possible to position the head even of small patients on the motor-adjustable headrest.

21.4 Hydrocolloid coolant supply


Additional equipment

If you want to work with a hydrocolloid coolant supply, you can obtain an adapter, Order No: 33 15 814, free of charge with the enclosed voucher.

Please state your address and the serial number of the chair (refer to "Installation Report / Warranty Passport").

The adapter can be connected to the cuspidor as follows:


Pull out the tumbler outlet (1), turning it slightly at the same time, and plug the adapter (2) into the self-closing valve until it snaps into place. Press the locking clip (3) before pulling out the connecting piece.

For water cooling, please activate the **Hydro** key in the setting dialog for the tumbler filling. **Continuous flushing** is then preset (page 22).

Hydrocolloid operation is switched on with the tumbler filling key in the assistant dialog.

21.5 Tray


Additional equipment

i NOTE

Maximum load: 1 kg


21.6 X-ray image viewer on the lamp support tube (additional equipment)

This X-ray image viewer is switched ON or OFF using the shown key on the control panel of the dentist or assistant element.

To view intraoral dental images the anti-glare film (1) must first be attached with two clips. Then fasten the intraoral dental images with the third clip over the cut-out.

21.7 X-ray image viewer on the tray (additional equipment)

This X-ray image viewer is switched ON or OFF using the shown key on the control panel of the dentist or assistant element.

To view intraoral dental images the anti-glare film (1) must first be attached with two clips. Then fasten the intraoral dental images with the third clip over the cut-out.

22 Maintenance

Despite the outstanding quality of your treatment center, it is necessary to perform inspections and preventive maintenance as well as safety tests at predetermined intervals in order to ensure its operational safety.

These events are documented in the Maintenance Manual.

Care and cleaning is done by the practice team.

For details please refer to the following sections.

22.1 Care and cleaning by the practice team

The practice team is responsible for regular cleaning and care of the treatment center, with the aim of minimizing the risk of contamination for patients and users.

The work to be performed is specified in the document **“Care and Cleaning by the Practice Team”**.

22.2 Inspection and maintenance

In order to guarantee the operational safety and reliability of your treatment center and to avoid damage due to natural wear, **yearly** inspection and maintenance must be performed on your treatment center. This is done by an authorized service engineer of your dental depot.

The steps to be performed as well as the parts which must be replaced are specified in the document **“Maintenance Certificate”**.

An overview of the inspection and maintenance work performed is additionally entered by the service engineer in the **“Installation Report / Warranty Passport”**. This document is part of the **“Maintenance Manual”**.

22.3 Safety tests

Medical products are designed in such a way that the first occurrence of a fault does not create a hazard to the safety of patients, users or other persons. Hence it is important to detect such faults before a second fault occurs, which might then lead to safety hazards.

For that reason it is essential to perform safety tests **every 2 years** which aim particularly at detecting electrical faults (e.g. isolation defects). This is done by an authorized service engineer from your dental depot, most practically together with the work to be performed according to section 22.2.

The check includes a visual inspection as well as measurements of the protective ground wire connections and the equivalent leakage currents.

The inspections and measurements to be performed are specified in the “**Maintenance Manual**”. The measured values must be documented there by the service engineer.

Safety tests must also be performed and documented during **initial start-up**, after extensions/upgrades (**conversion**) of your treatment center and after **repair work**.

 **NOTE**

The treatment center must not be operated if it has failed to pass the safety tests!

Safety checks for devices with HF surgery equipment:

In Germany, medical products are subject to the Ordinance on the Installation, Operation and Use of Medical Devices (MPBetreibV) from 06/29/1998.

According to Sec. 6, safety checks for devices with HF surgery equipment must be performed.

According to Sec. 7, a “Medical Product Log” must be kept in which the measurement values and the performance of the tests is documented.

These checks for devices with HF surgery equipment are identical with the checks as described above.

The “Maintenance Manual” is therefore at the same time the “Medical Product Log”.

The operator is obligated to maintain this medical product log.

In order to comply with the provisions of the Ordinance on the Installation, Operation and Use of Medical Devices (MPBetreibV), the following documentation must be maintained for treatment centers with HF surgery equipment in Germany:

- Realization of safety tests
- Personnel trained in the use of the HF surgery device (Sec. 5 of MPBetreibV)
- Repair work on the HF module
- Effects of malfunctions and repeated, similar operator errors
- Reporting of incidents to authorities and manufacturers

The Medical Product Log must be safekept for a period of at least 5 years after putting the system out of service.

Upon request, the Medical Product Log must be made available to the competent authority for inspection purposes at any time.

NOTE

As a user outside of Germany, please observe the legal requirements of your country.

22.4 Maintenance Manual

Keep this document near your treatment center.

Any inspection and maintenance work as well as all safety tests are documented by the service engineer in the Maintenance Manual.

On units with HF surgery equipment, the effects of malfunctions and messages to the appropriate authorities must be documented by the user.

The Maintenance Manual simultaneously acts as the Medical Product Log.

We recommend the user to always keep the documentation in the chapter "**Reporting of incidents to authorities / manufacturers**" up to date, regardless of any legal requirements.

We reserve the right to make any alterations which may be required due to technical improvements.

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

Different light times can be programmed in three levels of function **E1**, **E2** and **E3** for each practitioner **A**, **B** or **C** for the lamp **POLYLUX** on the practitioner element:

- Pre-select the function level **E1**, **E2** or **E3** for the handpiece being used
- Set the light time
- Press the programming button until a triple tone sounds. Then, reposition the handpiece.

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20.5 Application


The POLYLUX allows the polymerization of composites hardening by light.

The curing time is recommended by the manufacturer of the material.

A first empirical value may be determined by applying it to a test filling.

The effect of the light must be limited in the oral cavity to the area to be treated.

The light exit surface of the handpiece must be as close as possible to the material to be polymerized.

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To protect the lamp reflector, a built-in thermostat limits overheating.

It can trigger after several successive long exposures.

Upon return to operating temperature, the lamp can be put back into service.

After curing, increase the intensity of the light of the device by pressing the Composite button in the assistant program or the **MARCHE / ARRET** button of the lamp.

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© Sirona Dental Systems GmbH 2001
Sirona Dental Systems GmbH
Sirona Dental Systems LLC