

4008 HDF

Technical Manual



Fresenius Medical Care

4008 HDF

Technical Manual

4008 HDF Technical Manual

The Technical Manual contains all information necessary for performing maintenance and repair work.

The 4008 HDF option reflects the latest state of technology and complies with the requirements of EN 60601-1

Assembly, extension, adjustment, modification or repair may only be carried out by the manufacturer or persons authorized by him.

Any inquiries should be addressed to:

Head office:

Fresenius Medical Care
Deutschland GmbH
Borkenberg 14
D-61440 Oberursel/Ts., Germany
Tel.: 06171-60-0
Telex: 410805 fres d
Fax: 06171-251-58

Manufacturer:

Fresenius Medical Care
Deutschland GmbH
Plant Schweinfurt
Hafenstraße 9
D-97424 Schweinfurt, Germany
Tel.: 09721-678-0
Fax: 09721-678-200

Local Service:



How to use the Technical Manual


Search and find	<i>What?</i> Tables Contents	<i>Where?</i> Pages 0-5 and at the beginning of each chapter
Purpose	This manual is intended for: <ul style="list-style-type: none">– first studies (to acquire basic knowledge)– reference purposes (for start-up, maintenance and repair)	
Organization	The manual is divided into 6 chapters: <ul style="list-style-type: none">0 General Notes1 Description of machine functions2 Technical safety checks3 Calibration instructions4 Circuit descriptions and circuit diagrams5 Spare parts	
Numbering system	Page number 1-3 is to be interpreted as: Chapter 1, Page 3	
Qualification	This manual is intended for service technicians <ul style="list-style-type: none">– who are familiar with the current Operating Instructions (Operating Instructions relating to this Technical Manual are available under part no. 674 406 1)– who have the necessary background experience in mechanics, electrical and medical engineering– who have been authorized by the manufacturer to perform maintenance and repair work– who have access to the necessary auxiliary and measuring equipment	
Restrictions	The study of this manual does not represent an alternative to the training courses offered by the manufacturer.	
Manual changes	Manual changes will be released as new editions, supplement sheets or product information. Note: Modifications relating to circuit diagrams and component layouts (SP/BP) do not necessarily involve a change of the footer (edition). Refer to the index field of the respective circuit diagram or component layout for the respective state of these diagrams. The identification on the P.C.B. permits the user/technician to verify if the circuit diagram / component layout matches the P.C.B. actually installed. In general, this manual is subject to modification.	

Representation

New circuit symbols are used in the circuit diagrams. Potential data given in the circuit diagrams and setting instructions refer to the respective earth.
For example: \perp 24 means ground for 24 V voltage.

Component marking in circuit diagrams

Example:



This refers to a resistor with a position number 75 with a resistance of 1.5 Ohm.

The decimal point is replaced by a unit symbol (to reduce the possibility of errors).

Resistors:

R1: 0.1 Ω

1R5: 1.5 Ω

1K5: 1.5 k Ω

Capacitors:

μ 1: 0.1 μ F

1 μ 5: 1.5 μ F

1000 μ : 1000 μ F

Note:

When repairing or exchanging replacement parts make sure to take the applicable ESD precautions (e.g. EN 100 015-1).

During repair/troubleshooting in the hydraulic unit, protect the components from dialysate.

Technical data:

The technical data for the 4008 HDF option is to be found in Chapter 1 of the Operating Instructions.

Table of contents

Section		Page
1	Description of machine functions	1-
1.1	Description	1-3
1.2	Description of extended T1 test	1-8
1.3	Error messages.....	1-10
2	Technical safety checks	2-
3	Calibration instructions	3-
3.0	General information on the calibration instructions	3-3
3.1	Calibrating the UF2 pump	3-5
3.2	VDE inspections.....	3-7
3.3	Calibrating the 4008 HDF scale	3-8
3.4	Calibrating the substitute sensor	3-13
3.5	Calibrating the HDF blood pump	3-15
3.6	Repair instructions	3-17
4	Circuit descriptions and circuit diagrams	4-
4.1	LP 625 display board	4-3
4.2	LP 754 control board HDF	4-9
4.3	LP 760 HDF motor control	4-19
4.4	LP 761 HDF drive keyboard	4-27
4.5	LP 762 scale comm	4-33
4.6	LP 763 SSE serial interface extender	4-39
5	Spare parts	5-
5.0	How to use the spare parts catalog	5-3
5.1	P.C.B.s	5-4
5.2	Scales	5-6
5.3	Substitute lift	5-8
5.4	UF2 pump / blood pump (HDF) / valve V 126	5-10

Table of contents

1 Description of machine functions

Section		Page
1.1	Description	1-3
1.1.1	Components.....	1-3
1.1.2	Description	1-3
1.1.3	Component tests	1-3
	<i>Fig.: Block diagram</i>	1-4
	<i>Fig.: Flow diagram</i>	1-5
	<i>Fig.: Flow diagram</i>	1-6
1.2	Description of extended T1 test / error messages in T1 test	1-8
1.2.1	Test UF-Function	1-8
1.3	Error messages during treatment and the HDF test	1-10
1.4	Error messages, substitute pump	1-13

1.1 Description

1.1.1 Components

The components of the 4008 HDF option are permanently connected to the hemodialysis machine. The 4008 HDF option comprises the following components:

- Substitute lift
- Scales (weighing range < 20 kg) with data interface and taring facility
- Substitute pump with bidirectional data interface
- UF2 pump

1.1.2 Description

The scales determine the actual weight of the substitute reserve and signals this to the hemodialysis machine. This calculates the substitute rate depending on the set treatment time.

A substitute sensor mounted on the substitute lift recognises whether:

- Line set is fitted without substitute or
- Line set is fitted with substitute.

The substitute pump delivers the substitute solution to the venous bubble catcher. The filtrate from the dialysis fluid circuit is removed by the UF2 pump.

1.1.3 Component tests

The scales are checked by the operator by means of a plausibility test before treatment. The substitute sensor is tested while the 4008 HDF is filling.

The 4008 HDF test tests:

- The substitute pump
- The UF2 pump (electrically and hydraulically)

(If the 4008 HDF test is not carried out then the connections are to be tested but once).

The UF2 pump is tested cyclically during treatment.

The electric control system of the UF2 pump is checked continuously.

Fig.: Block diagram




Fig.: Flow diagram





Fig.: Flow diagram



Legend

2	Temperature sensor	77	Heat exchanger
3	Temperature sensor	78	Relief valve
5	Float switch	84	Disinfection valve
6	Level sensor	85	Disinfectant connector
7	Conductivity measuring cell	86	Recirculation valve
8	Blood leak detector	87	Discharge valve
9	Pressure transducer	88	Multifunction block
10	Reed contact for concentrate	88a	Degassing chamber
12	Reed contact for bicarbonate	88b	Secondary air separator
20	UF2 pump	88c	Primary air separator
21	Flow pump	89	Degassing orifice
22	UF pump	90a	Rinse chamber concentrate
23	Concentrate pump	90b	Rinse chamber bicarbonate
24	Dialyzer valve 1	91	Rinse valve
24b	Dialyzer valve 2	92	Vent valve
25	Bicarbonate pump	94	Concentrate suction tube
26	Bypass valve	95	Bicarbonate suction tube
29	Degassing pump	97	Ventilation pump
30	Drain valve	99	Rinse valve
31	Balancing chamber valve 1	100	Rinse valve
32	Balancing chamber valve 2	102	Concentrate valve in central delivery system
33	Balancing chamber valve 3	104	Bicarbonate valve in central delivery system
34	Balancing chamber valve 4	109	Temperature sensor
35	Balancing chamber valve 5	111	Hydrophobic filter
36	Balancing chamber valve 6	112	Vent valve
37	Balancing chamber valve 7	114	Dialysate filter
38	Balancing chamber valve 8	115	Disinfection valve sensor
41	Water inlet valve	116	Sampling valve
43	Fill valve	117	Check valve (concentrate)
54	Heater rod	118	Check valve (bicarbonate)
61	Pressure reducing valve	119	Filter (concentrate)
63	Filter/water inlet	120	Filter (bicarbonate)
65	Loading pressure valve	121	Concentrate connector in central delivery system
66	Heater block	122	Bicarbonate connector in central delivery system
66a	Water inflow chamber	123	Pressure switch for V 102
66b	Heater rod chamber	124	Pressure switch for V 104
66c	Float chamber	125	Adapter plate
68	Balancing chamber	126	4008 HDF vent valve
71	Filter/concentrate	148	Filter/rinse valve
72	Filter/bicarbonate	149	Filter/rinse valve
73	Filter/external dialysate	150	Filter
74	Filter/UF	151	Orifice
75	External flow indicator	152	4008 HDF hydrophobic filter
76	Filter/fill valve		

Hydraulic measurement points

- A Reduced water inlet pressure
- B Balancing chamber loading pressure
- C Flow pump pressure
- D Degassing pump pressure

1.2 Description of extended T1 test / error messages in the T1 test

Possible error messages displayed in the T1 test are listed in the Technical Manual for the hemodialysis machine. If the HDF option is installed and activated the following error codes can also appear:

1.2.1 Test UF-Function

Test description:


CPU 1 starts up the UF2 pump at a defined rate.

CPU 2 controls the hydraulic and electrical function of the UF2 pump.

CPU 2 blocks the control line of the UF2 pump and checks for standstill.

UF2 counter check.

Figure:



Error description:

<i>Error message</i>	<i>Description</i>
F11 UF-Function	The interval between the strokes of the UF2 pump was less than 220 ms. Correct delivery of the required volume is not safeguarded because the stroke return travel is too short. <ul style="list-style-type: none">– A pump rate which was too high was signalled by the CPU 1.
F12 UF-Function	The pulse time for the UF2 pump is less than 180 ms. Correct delivery of the required volume is not safeguarded because the discharge time is too short. <ul style="list-style-type: none">– Monoflop on LP 634 defect (IC 42/R65/C45).
F 13 UF-Function	The pulse for the UF2 pump is longer than 500 ms. A maximum rate of 5000 ml/h is not possible. <ul style="list-style-type: none">– Monoflop on LP 634 defect (IC 42/R65/C45).
F14 UF-Function	UF2 pump not active during the test (4 s). <ul style="list-style-type: none">– Acknowledgement (UF_P2, X637/B26) → X632/C7 no LOW pulse– Control line (UF_P2, X634L/A–C24) → X637/B26 no LOW pulse
F15 UF-Function	The UF2 pump cannot be stopped by the CPU2. <ul style="list-style-type: none">– Control line (UF_P_EN, X632/C28) → X634R/A22 no 5 V.– Reset input to IC 42/pin 13 on LP 634 defect.
F16 UF-Function	The UF pump acknowledgement from CPU1 is faulty. <ul style="list-style-type: none">– Acknowledgement (UF_P2, X637/B26) → X633/C23 no LOW pulse
F17 UF-Function	The change in pressure after a UF2 pump stroke less than 20 mmHg. <ul style="list-style-type: none">– UF2 pump mechanically defect.– Control line (UF_P2_CRTL, X632/B24) → X634R/C11 no HIGH pulse.
F20 UF-Function	Difference in pressure between UF pump stroke and UF2 pump stroke more than 20 %. <ul style="list-style-type: none">– Stroke rate for UF pump or UF2 pump not set correctly.

1.3 Error messages during treatment and the HDF test

Message on the display	Cause	Possible error elimination
F321 HDF-failure	The UF pump or UF2 pump test will not be started as long as the level sensor (6) detects air. The level sensor (6) will detect air one minute after the 4008 HDF test has been selected	Check level sensor (6) (defect).
F322 HDF-failure	The UF pump or UF2 pump has carried out more than 50 strokes and still no air has been detected.	Check V126 (does not open, venous bubble catcher collapses), check V43 (whether leaking or open). Check UF pump or UF2 pump, respectively (mech. defect).
F323 HDF-failure	The substitute pump runs at a rate of 10 ml/min. to raise the fluid level by a defined amount in the secondary air separator (the fluid will just be detected). Air will still be detected after 2 minutes.	Check substitute tubing (clamped off, not connected to the venous bubble catcher).
F324 HDF-failure	Since fluid has been recognized in the secondary air separator the substitute pump starts at a rate of 10 ml/min. until a change in weight follows on the scales (i.e. a jump in grammes). There was no difference in weight established after one minute.	Check substitute tubing (clamped off).
F325 HDF-failure	If the expected change in weight (i.e. a jump in grammes) is higher or lower than -1 g, the respective UF pump test will be repeated. The UF pump test was repeated more than 5 times.	Check bag (fluctuations in air draft) Check scales (for drifting).

Message on the display	Cause	Possible error elimination
F326 UF-failure	The fluid level in the secondary air separator is raised again by the substitute pump after 100 strokes of the UF pump or UF2 pump, respectively, until fluid is detected. Plus run-on time at a rate of 10 ml/min. for the gramme jump. Air will still be detected after a specific period of time (depending on the delivery rate of the substitute pump).	Check substitute hose (clamped off, not connected with the venous bubble catcher).
UF1 volume-Error or UF2 volume-Error	One UF pump did not pass the test. The filling volume for the secondary air separator is not within the given tolerance of 100 ml; ± 5 ml. Should the test give a reading of over 105 ml the cause could also be air taken in the flow from a badly vented dialysator.	Check UF pump or UF2 pump, respectively (not calibrated, mechanical defect).
F327 UF-failure	Interval between two strokes of the UF pump less than 220 ms.	Check CPU 1 (defect).
F328 UF-failure	Pulse time of a stroke of UF pump less than 180 ml.	Check LP 634 (regulating monoflop defect).
F329 UF-failure	Pulse time of a stroke of UF pump more than 500 ms.	Check LP 634 (regulating monoflop defect).
F330 UF-failure	Starting time for the UF pump more than 10 seconds.	Check LP 634 (regulating end stage defect).
F331 UF-failure	Difference between desired or actual delivery rate of the UF pump greater than $\pm 10\%$.	Check CPU 1 / CPU 2 (communication problems).
F332 UF-failure	UF pump stops longer than the maximum period time.	Check LP 634 (regulating end stage defect). Check UF pump (interruption, control). Check CPU 1 / CPU 2 (communication problems)
F333 UF-failure	Despite switched off ultrafiltration the change in delivery rate of the UF pump is greater than 10 ml.	Check CPU 1 / CPU 2 (communication problems)
F334 UF-failure	Interval between two strokes of the UF2 pump less than 220ms.	Check CPU 1 (defect).

Message on the display	Cause	Possible error elimination
F335 UF-failure	Pulse time of a stroke of the UF2 pump less than 180 ms.	Check LP 634 (regulating monoflop defect).
F336 UF-failure	Pulse time of a stroke of the UF2 pump more than 500 ms.	Check LP 634 (regulating monoflop defect).
F337 UF-failure	Reaction time of the UF2 pump more than 10 seconds.	Check LP 634 (regulating end stage defect).
F338 UF-failure	Difference in desired or actual delivery rate of UF2 pump greater than $\pm 10\%$.	Check CPU 1 / CPU 2 (communication problems).
F339 UF-failure	UF2 pump stops longer than the maximum period time.	Check LP 634 (regulating end stage defect). Check UF pump (interruption, control). Check CPU 1 / CPU 2 (communication problems).
F340 UF-failure	Despite switched off 4008HDF change in delivery rate of the UF2 pump greater than 10 ml.	Check CPU 1 / CPU 2 (communication problems).
F341 UF-failure	Failure of UF pump.	Check UF pump (spring, screen)
F342 UF-failure	Failure of UF2 pump.	Check UF2 pump (spring, screen)
F343 UF-failure	Volume difference between UF pump and UF2 pump.	Check UF pump and UF2 pump (delivery volume)
F344 HDF-failure	Balance failure recognized by CPU2 greater than ± 500 ml.	Check monitor / scales (communication problems).
F345 HDF-failure	Bolus exceeded by more than + 20 ml, detected by CPU2.	Check monitor / scales (communication problems).
F346 HDF-failure	CPU2 failed to perform the cyclic UF pump test within 5 minutes.	Check V24B and V43 feed-back lines. Check sense of rotation of aspiration pump. Check hydraulics unit for leaks.
F348 HDF-failure	Weight change on the scales during the weighing test greater than 2g.	Check scales Check substitute bag
F349 HDF-failure	The delivery rate correction factor could still not be determined after 20 minutes treatment time.	Check substitute bag (weight fluctuations) Check HDF control system (greater fluctuations)

1.4 Error messages, substitute pump

Code	Description of failure
E.01	Line diameter outside the permitted range
E.02	Non-defined hexaswitch position
E.03	Venous pressure transducer not balanced
E.04	Failure, running time monitoring system, SN operation
E.05	SN stroke volume outside the permitted range
E.06	The SN pressure thresholds outside the value range of the AD converter
E.07	Not defined
E.08	Failure in the AD conversion
E.09	Not defined
E.10	Not defined
E.11	Not defined
E.12	Failure, speed monitoring (hall sensor)
E.13	Failure, monitoring system, current sensing resistors
E.14	Failure, monitoring system, current sensing resistors
E.15	Failure, speed monitoring system

2 Technical safety checks

General Notes

This chapter includes all necessary technical safety checks (TSC).

These inspections must be carried out every 12 months.

The technical safety inspections stipulated for the hemodialysis machine must be carried out in addition to these technical safety checks.

The technical safety checks are to be recorded in the equipment log.

A technical safety checks report is to be found at the page 2-3.

Technical safety checks report

Manufacturer: **Date:**
Machine: **Technician:**
Operating hours:

Extent of checks and time limits for technical safety checks.
TIME LIMIT: every 12 months

The following checks and inspections have to be carried out on the machine at least every 12 months by persons who are capable of carrying out such technical safety checks efficiently as a result of their training, knowledge and experience gained in practice and are not subject to instructions as to their checking activity.
(See also BMA Bulletin dated 02.04.1987)

TSC	WA	No.	Description	Desired value/function	OK
1 UF2 pump					
TSC		1.1	Delivery volume (1 ml/stroke)	60 strokes = 60 ml ± 1 %	<input type="checkbox"/>
2 VDE inspections (values according to EN 60601-1)					
TSC		2.1	Protective ground resistance	max. 0.3 Ω	<input type="checkbox"/>
TSC		2.2	Summarized leakage current	Must fulfill both conditions: 1. No more than 1.5 times the "summarized value first measured" ("Summarized value first measured": Refer to the machine card enclosed with the machine.) 2. No more than 1 mA	<input type="checkbox"/>

The technical safety checks are to be recorded in the equipment log and the results of the checks documented.

If the machines are not safe in function and/or operation they are to be repaired or the operator is to be informed as to the dangers involved when using the machine in its given state.

The correct completion of the listed work and the correctness of statements made are hereby certified.

Signature, technician:

Signature, customer:

....., date

....., date

.....

.....

Table of contents

3 Calibration instructions

Section		Page
3.0	General information on the calibration instructions	3-3
3.1	Calibrating the UF2 pump	3-5
3.2	VDE inspections	3-7
3.3	Calibrating the 4008 HDF scales	3-8
3.4	Calibrating the substitute sensor	3-13
3.5	Calibrating the HDF blood pump	3-15
3.5.1	Pressure transducer	3-15
3.5.2	Checking SN switchover points	3-16
3.6	Repair instructions	3-17

3.0 General information on the calibration instructions

Measuring instruments:

The same measuring instruments as used for the hemodialysis machines are to be employed. Required in addition is a calibration weight PT6 – 5000 g.

Fig.: P.C.B. overview 4008 HDF






Fig.: DIP switch 4008 HDF



SW 1, DIP switch 8:
ON: 4008 HDF test selectable (required/not required)
OFF: 4008 HDF test carried out automatically

3.1 Calibrating the UF2 pump




Switch off hemodialysis machine.
Set service switch to ON (top).
Switch on hemodialysis machine.


Press the **Confirm** key.




Select **ADJ. UF-PUMP VOLUME** by pressing the ▲ and ▼ keys.
Press the **Confirm** key.



Select **UF-PUMP 2** by pressing the ▲ and ▼ keys.
Press the **Confirm** key.








Enter the number of strokes by pressing the ▲ and ▼ keys.



Press the **Override** key.



Text display for approx 3 sec.



Text display

Draw off UF 2 pump line, close T-piece, hang line in a measuring cylinder.

Press the key **UF I/O**.

The **UF** LED emits light.

The UF2 pump delivers.

The remaining UF Pulses are shown in the display.

The number of preselected pulses is shown in the UF Goal display.

Text display after completion.

Check delivered quantity.
If necessary adjust UF 2 pump and repeat process.


Delivery rate: 1 ml/stroke
(60 strokes = 60 ml)
Tolerance: ± 1 %.

3.2 VDE inspections

(Values in accordance with EN 60601-1)

.1 Protective ground resistance:

– max. 0.3 Ω .



.2 Summarized leakage current

Must fulfill both conditions:

1. No more than 1.5 times the “summarized value first measured”
 (“Summarized value first measured”: Refer to the machine card supplied with the machine.)
2. No more than 1 mA

The measurements must be taken in the dialysis mode of operation in the “ON phase” of the heating control system.

The scales must be travelled out to such an extent that neither of the two end switches are actuated (middle position).

The two leakage currents are to be measured each time with different mains connection poling. The two leakage currents are added together and give the summarized leakage current.


3.3 Calibrating the 4008 HDF scales



Note:

The prompts shown in the text display have to be run through correctly by the technician since no *direct* return signal is given by the scales.

Switch off hemodialysis machine.
Set service switch to ON (top).
Switch on hemodialysis machine.



Press the **Confirm** key.



Select **CALIB. HDF-SCALE** by pressing the ▲ and ▼ keys.

Press the **Confirm** key.




Brief text display.


Display on scales is not illuminated.



Text display.




- ① Open housing of scales at the back.
- ② Set calibration switch for the scales in ON position.




Press the **Override** key.

HDF-DISPL shows C?




Text display.

The display on the scales is active.
A "C" appears at the top right.



Press the **Override** key.


HDF-SCALE unloaded ?



Text display.
No weight is to be on the scales.


Press the **Override** key.

HDF-DISPL shows 0 ?



Text display.

Scales display.





Press the **Override** key.



Text display.

Scales display: 5000 g

Hang calibration weight PT6 – 5000 g on scales.





Press the **Override** key.






Text display.

Scales display: 5000 g.








calib. Switch OFF ?



HDF-SCALE unloaded ?



SCALE CALIBRAT. done



CALIB. HDF-SCALE



Press the **Override** key.

Text display.

-
- ① Set calibration switch for the scales to OFF position.
 - ② Close housing.

Press the **Override** key.


Text display.

Remove calibration weight.

Press the **Override** key.

Brief text display.

Text display.



Switch off hemodialysis machine.

Set service switch to OFF (bottom).

Switch on hemodialysis machine.


Hang calibration weight PT6 – 5000 g on scales.

Scales display: 5000 g


Remove calibration weight.

Calibration of 4008 HDF scales completed.


3.4 Calibrating the substitute sensor




Place filled substitute line in substitute sensor.




Set potentiometer P1 on the LP 754 so that the LED D26 just illuminates.



Place empty substitute line in substitute sensor.






Adjust potentiometer P1 on LP 754 until LED D26 just turns off.
Count the number of revolutions the potentiometer is turned. Keep this number in your mind.



Turn back the potentiometer P1 on the LP 754 by half this number.
The LED D26 will illuminate again.

3.5 Calibrating the HDF blood pump

3.5.1 Pressure transducer



Switch the hemodialysis machine off.

Set HEX switch on the LP 754 to position "F".

Switch the hemodialysis machine on.

Display **000**


Open SN pressure transducer to atmosphere.

Press the **Start/Stop** key.


Display **250**

Apply a pressure of 250 mmHg to the SN pressure transducer.

Press the **Start/Stop** key.



Display ---



Switch off hemodialysis machine.

Set HEX switch on the LP 754 to position **1**.

Balancing out has been completed.

3.5.2 Checking SN switchover points

Bottom switchover point: fixed at 75 mmHg (± 6 mmHg)

Top switchover points (± 6 mmHg):

Stroke volumes (ml)	10	15	20	25	30	35	40	45	50
Top switchover point	110	130	150	172	195	219	244	270	299

3.6 Repair instructions



Important:

If the substitute lift has to be dismounted for repair the seal is to be checked and if necessary renewed before fitting the lift again.

Drip water protection!

Table of contents

4 Circuit descriptions and circuit diagrams

Section		Page
	<i>Fig.: Connection diagram 4008 HDF</i>	4-3
	<i>Fig.: Wiring diagram 4008 HDF</i>	4-5
	<i>Fig.: Block circuit diagram, mains voltage 4008 HDF</i>	4-7
	<i>Fig.: Block circuit diagram, low voltage 4008 HDF</i>	4-7
4.1	LP 625 display board	4-9
4.1.1	Description	4-9
4.1.2	Circuit and component layout diagram	4-11
4.2	LP 754 control board HDF	4-15
4.2.1	Description	4-15
4.2.2	Circuit and component layout diagram	4-17
4.3	LP 760 HDF motor control	4-23
4.3.1	Description	4-23
4.3.2	Circuit and component layout diagram	4-25
4.4	LP 761 HDF drive keyboard	4-29
4.4.1	Description	4-29
4.4.2	Circuit and component layout diagram	4-31
4.5	LP 762 scales comm	4-35
4.5.1	Description	4-35
4.5.2	Circuit and component layout diagram	4-37
4.6	LP 763 serial interface extender	4-41
4.6.1	Description	4-41
4.6.2	Circuit and component layout diagram	4-43

Fig.: Connection diagram 4008 HDF




Fig.: Wiring diagram 4008 HDF




Fig.: Block circuit diagram, mains voltage 4008 HDF





Fig.: Block circuit diagram, low voltage 4008 HDF



4.1 LP 625 display board

4.1.1 Description

- **General**

The displays and the keyboard are on this P.C.B..

Plug connections to LP 625:

- X189, connection to LP 760

- **Display**

All necessary information is given to the user through 6 multiplex 7-segment displays and 2 LEDs which are mounted on the socket.


- **Keyboard**

The keyboard comprises 3 single non-illuminated keys and 1 key with integrated LED (HDF). The keys are soldered directly to the P.C.B..

The user receives a acknowledgement through a spring contact when the keys are pressed.

4.1.2 Circuit and component layout diagram
LP 625 display board


LP 625
Circuit diagram



NO.	DATE	DES.	NAME	BLATT
1	24.11.97	ST		1 (11)
2	24.11.97	ST		
3	24.11.97	ST		
4	24.11.97	ST		
5	24.11.97	ST		
6	24.11.97	ST		
7	24.11.97	ST		
8	24.11.97	ST		
9	24.11.97	ST		
10	24.11.97	ST		
11	24.11.97	ST		
12	24.11.97	ST		
13	24.11.97	ST		
14	24.11.97	ST		
15	24.11.97	ST		
16	24.11.97	ST		
17	24.11.97	ST		
18	24.11.97	ST		
19	24.11.97	ST		
20	24.11.97	ST		
21	24.11.97	ST		
22	24.11.97	ST		
23	24.11.97	ST		
24	24.11.97	ST		
25	24.11.97	ST		
26	24.11.97	ST		
27	24.11.97	ST		
28	24.11.97	ST		
29	24.11.97	ST		
30	24.11.97	ST		
31	24.11.97	ST		
32	24.11.97	ST		
33	24.11.97	ST		
34	24.11.97	ST		
35	24.11.97	ST		
36	24.11.97	ST		
37	24.11.97	ST		
38	24.11.97	ST		
39	24.11.97	ST		
40	24.11.97	ST		
41	24.11.97	ST		
42	24.11.97	ST		
43	24.11.97	ST		
44	24.11.97	ST		
45	24.11.97	ST		
46	24.11.97	ST		
47	24.11.97	ST		
48	24.11.97	ST		
49	24.11.97	ST		
50	24.11.97	ST		

LP 625

Component layout diagram



4.2 LP 754 control board HDF

4.2.1 Description

- **General**

The control, section and power pack are on the this P.C.B..

Plug connections to LP 754:

- X186, connection to position sensor
- X188, connection to stepper motor
- X189, connection to LP 748
- X190, connection to cover switch
- X192, connection to pressure transducer
- X348, connection to dialysis monitor

Hexswitches on the LP 754:

- 0 No function
- 1 Single needle blood pump / 4008 HDF / ON LINE HDF
- 2 – A No function
- B Blood pump stop alarm (15, 30 secs.)
- C No function
- D Call-up operating time meter (display x 100 = number of hours)
- E Test operation (only manufacturer)
- F Pressure transducer calibration

- **Voltage generation**

The +24 V and +12 V voltage supply is made available to the blood pump by the monitor. The +5 V supply voltage is generated on the module by the switch controller IC 20 from the +24 V voltage to minimize loss of power.

- **Stepper motor control**

The stepper motor is run in microstep operation to reduce noise. The resolution is 60 microsteps per step. The RISC processor transmits a 8-bit dataword alternately on the pins 3 and 5 of the DA converter from IC 7. Two sine form voltages are available on the output of the converter, phase offset by 90°. They are conducted to the stepper motor controller IC 2 together with the current direction signals. Together with the two SM drivers IC 1/IC 22 and the current sensor resistors R58/R59 as feedback these form a closed control circuit which determine two sine-shaped currents, phase offset by 90°, in the two coils.

- **Microprocessor**

The cycle frequency of the microprocessor is determined by the quartz Q2 between the processor connections 39 and 40. Quartz oscillation is made possible by the capacitors C5 and C6.

The keyboard signals and the signals BSST, BPST are read-in through the port P4.

The IC 9 serves as intermediate memory for the addresses AO/A7.

The ALE signal on pin 50 of the microprocessor is the control line of the data address latch.

The signal from the revolution and direction of turn sensor (position sensor) in the pump bed is read-in through port T1.

The operation data of the pump is saved in the NOVRAM IC 21 through ports 5.1 to 5.3 when the dialysis machine is switched off or there is a voltage failure. Undervoltage detection (power down) is through the comparator IC 23.

With the WR line the data takeover for the display is controlled in the external data latch IC 14.

- **PLL intermediate circuit**

PLL component IC 4 together with the counter IC 19 causes the frequency coming from the processor to control the stepper motor to be multiplied otherwise the processor would be too slow to generate this frequency.

- **RISC processor**

The RISC processor IC 5 is supplied with a cycle signal from the processor. With each cycle an 8-bit data word is read from a look up table by the processor alternately for each phase of the stepper motor, respectively. Included in this data word is both the current direction as well as the current desired value.

A watchdog for the CPU is also incorporated in the RISC processor. When a cycle is missing on pin 8/IC 5 the RISC processor triggers a reset by the CPU through pin 7.

- **Display control**

The data word for controlling the display is stored in IC 14. The multiplex operation of the individual digits is made possible through the decoder IC 18.

- **Speed and flow**

The speed is transmitted by the processor through port P 1.1 and the flow through port P 1.2 to the dialysis monitor.

- **Pressure measurement**

The differential measurement amplifier for the pressure transducer is formed by the IC 6 (1/2/3) and IC 6 (5/6/7). The measuring signal is read-in by the processor through the AD converter input AN 7 and calibrated for zero and steepness by the software. Afterwards the measuring signal is modulated to pulse width through port P1.3 and transmitted through a subsequent DA converter IC 11 (5/6/7) to the monitor.

**4.2.2 Circuit and component layout diagram
LP 754 control board HDF**

LP 754
Circuit diagram 1/2

LP 754
Circuit diagram 2/2

LP 754

Component layout diagram

4.3 LP 760 HDF motor control

4.3.1 Description

- **General**

The LP 760 serves to control the lift drive.

Lift operation presupposes that the scales are swung away. This is recognized by a limit switch. The scales can be moved upwards or downwards by pressing the ▲ or ▼ keys. The drive is switched off automatically, likewise by limit switches, when reaching the uppermost or bottom-most position. When the scales are being travelled upwards or downwards they can be stopped at any time with the **Stop** key.

The motor for the lift drive is safeguarded against overloading by means of an integrated thermoswitch. If this protection facility responds the lift will come to a standstill. The lift can be used again after the motor has cooled down.

- **Circuit structure**

Supply voltage

The control circuit is supplied with +12 V. This voltage is generated from the +24 V voltage of the hemodialysis machine by the linear controller T1 in conjunction with the capacitors C4, C5 and C6.

Motor control

Two types of circuit are implemented to control the motor:

- a) continuous operation (R21 with components, R22 without components)
- b) inching operation for downwards movement (R21 without components, R22 with components).

To Point a):

One D flip-flop is available for each direction of movement (IC 3). The flip-flops are set by a high level on X2/2 (key ▲) or X2/4 (key ▼). The flip-flop for the respective other direction of movement is reset at the same time by pressing the appropriate key. The two flip-flops are reset through high level on X2/3 (**Stop** key). This also takes place should the head of the scales be swung over the machine (high level on X1/6).

To Point b):

There is only one flip-flop available for the upwards movement (IC 3/pin 1/2). The flip-flop for the downward movement is set through IC 5. The control logic for the downward operation is not influenced by this (see Point a). However, the drive is only active as long as the key ▼ is pressed.

The IC 5 works as an 8 to 1 multiplexer. The address (0 to 7) of the input to be interconnected is provided through the pins 11 to 13. Only the address 3 leads to a low level on output IC 5/14 and, as a result, to the motor starting. R 20 serves as a pull-up resistor. A defined output level is obtained with the IC 5/14 in high ohmic state, caused by pressing the **Stop** key (high level on IC 5/15).

- **Other circuit sections**

The motor feed lines are circuited through the relay RL 1 and RL 2. The position of RL 1 determines the travel direction of the lift. The motor is switched on and off through RL 2. Both relays are controlled by the transistors T 2 (RL 1) and T 3 (RL 2). T 2 is controlled directly from the non-inverted output IC 3/1. The output signals of both flip-flops, IC 3/2 and IC 3/12 (inverted) are used to control T 3 and are linked with each other through IC 1/1 and IC 1/2. In the case of low level on one of the two flip-flop outputs IC 3/2 or IC 3/13 (i.e. the respective flip-flop has been set) the IC 1/pin 3 links up high level and T 3 switches RL 2.

- **Motor current recognition**


The motor current is monitored while the lift is operating. The internal thermal-protective switch responds in the case of the motor becoming overheated. Both flip-flops (IC 3) are reset by motor current recognition so that the motor does not start up again automatically after the drive has cooled down. The voltage drop is evaluated through the high load resistor R 18 for motor current recognition. The generated alternate voltage caused by the current flow through R 18 is connected to the capacitor C 14 and C 13 through the diodes D 4 and D 5. The resulting direct voltage through R 13 is approx. double so high as the peak value of a half wave of the alternate voltage on R 18 less the on state voltages of D 4 and D 5. The direct voltage controls the optocoupler IC 4 through R 14. R 13 serves to completely discharge C 13 and C 14 when the motor is switched off.

IC 4 ensures protective isolation of the control circuit from the mains voltage conducting components of the current recognition system.

The output of the optocoupler is circuited with the pull-up resistor R 19 against 12 V. High level on the collector of IC 4 resets both flip-flops through IC 1/pin 12 and 13 and IC 1/8. This would lead to continuous blocking of the control circuit unless further measures were taken. This is why the output of IC 4 is connected up parallel to a further transistor T 4. In switched off state the capacitor C 12 is charged through high level on IC 1/4 through D 3 and R 10. This makes the transistor T 4 conductive. After setting one of the flip-flops (IC 3) the IC 1/4 goes to low. C 12 discharges itself through R 11, R 12 and the B-E path from T 4. T 4 blocks after approx. 500 ms. The output transistor must be linked by IC 4 before this time has elapsed otherwise both flip-flops will be reset again.


4.3.2 Circuit and component layout diagram LP 760 HDF motor control

LP 760 Circuit diagram



LP 760

Component layout diagram



4.4 LP 761 HDF drive keyboard


4.4.1 Description

- **General**

The keys for controlling the substitute lift are found on this P.C.B..


4.4.2 Circuit and component layout diagram
LP 761 HDF drive keyboard

LP 761
Circuit diagram



LP 761

Component layout diagram



20
19
18
17
16
15
14
13
12
11
10
9
8
7
6
5
4
3
2
1
0

				DATUM	NAME	HDF-ANTRIEB-TASTATUR HDF-MOTOR-KEYBOARD	BLATT
				BE 04.08.1994	OTT		
				GE 05.08.1994	ASCH		
						LP761	673 766
a	X1 geändert						FORMAT A4
AUS	AENDERUNG	TAG	NAME				

4.5 LP 762 scales comm

4.5.1 Description

The LP 762 serves as interface between the dialysis machine and the scales.

The serial transmission from the dialysis machine arrives on the plug X1 pin 4 (RXD) and is converted at the interface IC 1 (max. 232) to 5 V level. IC 2 makes a serial parallel conversion and passes the signal on to the scales. The scales are switched on or off, respectively, by the transistor T1.

The serial carry forward (weight) comes from the scales to the plug X3 pin 6 (TXD) and is converted by IC 1 from 5 V to RS 232 level. The signal arrives at the dialysis machine through plug X1 pin 5 (TXD).

SP1 (ML78L05) converts the 12 V into 5 V which serves as the supply voltage for LP 762.

**4.5.2 Circuit and component layout diagram
LP 762 scales comm**

LP 762
Circuit diagram

LP 762
Component layout diagram

4.6 LP 763 SSE serial interface extender

4.6.1 Description

- **General**

The LP 763 is necessary for serial communication of the monitor with the scales and the HDF pump module.

- **Circuit structure**

The circuit of the LP 763 comprises two circuit units independent of each other.

Serial interfaces

Heart of the LP 763 P.C.B. is an 8-fold UART (universal asynchronous receiver/transmitter), i.e. an interface IC (IC 1) with eight independent serial channels. Also to be found on the P.C.B. are two TTL/RS232 level converters (IC 2 and IC 3) and an interrupt logic IC (IC 4), as well as diverse resistors, capacitors and an oscillation quartz Q1 for generating the working cycle for IC 1.

COMMCO connection

The slot 63Z used by the LP 763 P.C.B. is also provided for the installation of P.C.B. LP 752. The LP 752 takes over the connection of the 4008 monitor signals to the COMMCO P.C.B. LP 729. In order to be able to also continue the linking to COMMCO/FINESSE in 4008 HDF machines the circuit section of the LP 752 P.C.B. is taken over in the LP 763 P.C.B..

- **Function**

After switching on the machine the UART is initialised with the transmission parameters.

Transmission parameter channel A:

1200 Baud, 8 databits, 1 stopbit, no paritibit

Transmission parameter channel B:

9600 Baud, 8 databits, 1 stopbit, no paritibit

Channel A is responsible for communication with the scales, channel B for communication with the HDF pump module.

After going through the watchdog test the scales are switched on through the transmission line from channel A and the weight requirements of the 4008 monitor are communicated to the scales. The weight value of the scales is read by the 4008 monitor on the reception line.

Commands and data from the 4008 monitor are sent through the transmission line from channel B to the HDF pump modules. Data which comes from the HDF pump module is received by the 4008 monitor through the receiving line.

The remaining six communication channels are not used in 4008 operation.

If UART data is ready for collecting or if the transmission buffer is ready to accept new digits during a transmission the component causes an interrupt at CPU 1. CPU 1 processes this interrupt demand (in so far as nothing more important is to be done) and reads a digit from the reception buffer or writes a digit on the transmission buffer.

**4.6.2 Circuit and component layout diagram
LP 763 SSE serial interface extender**

LP 763
Circuit diagram 1/2

LP 763
Circuit diagram 2/2

LP 763
Component layout diagram

Table of contents

5 Spare parts

Section		Page
5.0	How to use the spare parts catalog	5-3
5.1	P.C.B.s	5-4
5.2	Scales	5-6
5.3	Substitute lift	5-8
5.4	UF2 pump / blood pump (HDF) / valve V126	5-10

5.0 How to use the spare parts catalog

Objective:

Spare parts can be defined and ordered.

Organization:

This spare parts catalogue contains 4 assemblies.

5.1 P.C.B.

5.2 Scales

5.3 Substitute Lift

5.4 UF2 pump / blood pump (HDF) / valve V126

Each assembly has its assembly number. 5.X.0 refers to the basic version.

Major changes (modifications) within an assembly can be identified by an increasing decimal (e.g. 5.X.1, 5.X.2, ...)

In the spare part lists it is referred to:

Built-in

From Equipment Code

To Equipment Code

Update service:

Updates to this spare parts catalog will be released as:

- Replacement pages
- Supplementary pages
- Technical information


Generally we reserve the right to make changes.

Requirements for correct spare part orders:

The equipment code must be entered in the machine record.

All modifications must be entered in the machine record. Update the equipment code, if necessary, to ensure that the correct spare parts will be ordered.

5.1 P.C.B.s





Built-in
From/until Equipment Code

000 – 000

Item	Part number	Description
1	673 790 1	P.C.B. LP 762 scales comm
2	674 182 1	P.C.B. LP 763 SEE serial interface extender
3	673 766 1	P.C.B. LP 761 HDF drive keyboard
4	673 872 1	P.C.B. LP 625 display board
5	673 720 1	P.C.B. LP 760 HDF motor control
6	673 041 1	P.C.B. LP 754 control board HDF

5.2 Scales




Built-in
From/until Equipment Code

000 – 000

Item	Part number	Description
1	673 983 1	Sliding bush
2	552 244 1	Screw
3	646 352 1	Screw
4	674 409 1	Microswitch
5	552 312 1	Screw
6	553 330 1	Screw
7	673 635 1	System unit, scales
8	643 161 1	Screw
9	673 988 1	Hook
10	673 987 1	Latch
11	642 861 1	Screw
12	553 621 1	Nut
13	673 986 1	Adapter
14	563 151 1	Retaining washer
15	674 106 1	Line holder
16	552 275 1	Screw
17	673 971 1	Fastening plate
18	642 693 1	Screw
19	673 985 1	Bearing
20	672 530 1	Cover
21	552 276 1	Screw
22	674 551 1	Caution label: Stability limit
23	673 977 1	Display plate
24	673 922 1	Back shell
25	673 921 1	Front shell
26	553 620 1	Nut
27	559 690 1	Washer
28	642 648 1	Spacer sleeve
29	640 143 1	Spacer
30	673 923 1	Front foil display
	674 068 1	Scales, complete


5.3 Substitute lift






Item	Part number	Description
1	674 399 1	Connection cable
2	644 205 1	Screw
3	673 990 1	Locking ring
4	674 106 1	Line holder
5	653 099 1	Sensor
6	552 244 1	Screw
7	674 334 1	Pipe
8	673 975 1	Clamping plate
9	552 275 1	Screw
10	553 733 1	Toothed washer
11	674 395 1	Drive
12	673 973 1	Dog
14	673 755 1	Screw
16	553 621 1	Nut
17	537 203 1	Fastening clip
18	646 274 1	Tolerance sleeve
19	674 078 1	Stripper
20	673 974 1	Collar bush
21	673 920 1	Housing, lift
22	673 639 1	Front foil, lift
23	559 690 1	Washer
24	674 400 1	Connection cable
25	674 358 1	Connection cable
26	673 918 1	Base plate
27	673 972 1	Funnel
28	674 255 1	P.C.B. cover
29	646 606 1	Spacer bolt (M 3 x 8.5)
30	650 439 1	Spacer bolt (M 3 x 12.5)
31	559 690 1	Washer
32	574 617 1	Spacer sleeve
33	644 040 1	Spacer bolt (M 3 x 13)
34	673 982 1	Bush
35	674 487 1	Sealing tape
36	675 052 1	Split pin
37	553 694 1	Washer
38	674 333 1	Locating pin
39	552 312 1	Screw M5x16
40	675 222 1	Fork
41	675 051 1	Bolt with head
	673 915 1	Substitute lift, complete

5.4 UF2 pump / blood pump (HDF) / valve V126

Hydraulic section 4008 / 4008 E



Hydraulic section 4008 B



Built-in
From/until Equipment Code

000 – 000

Item	Part number	Description
-------------	--------------------	--------------------

UF2 pump

The UF2 pump is identical with the ultrafiltration pump
(see appropriate Spare Parts Catalogue 4008)

	670 706 1	Ultrafiltration pump, complete
--	-----------	--------------------------------

Blood pump (HDF)

The blood pump (HDF) is similar to the blood pump (single needle)
(see appropriate Spare Parts Catalogue 4008)

40	671 726 1	Front foil
	673 924 1	Blood pump module HDF, complete

Valve V126

1	645 881 1	Filter
2	545 361 1	Silicon tube (0.12 m)
3	674 349 1	Solenoid valve MV 126 / solenoid valve MV43
4	674 347 1	Valve block
5	552 275 1	Screw
6	579 072 1	O-ring
	674 941 1	Valve block, complete

