
# Serviceanleitung Servicing Instructions

## CardioServ

Version 3

**Achtung:**

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Version 3

**Caution:**

During repairs/service interventions,  
observe the protective measures against  
damage due to ESD.

- \* MARQUETTE HELLIGE considers itself responsible for the effects on safety, reliability, and performance of the equipment, only if:
  - assembly operations, extensions, readjustments, modifications, or repairs are carried out by MARQUETTE HELLIGE or by persons authorized by MARQUETTE HELLIGE,
  - the electrical installation of the relevant room complies with the applicable national and local requirements, and
  - the instrument is used in accordance with the instructions for use.
- \* This manual contains service information; operating instructions are provided in the user manual of the instrument.
- \* This manual is in conformity with the instrument at printing date.
- \* All rights are reserved for instruments, circuits, techniques, and names appearing in the manual.

**Contents**

1.	Functional Description . . . . .	page 4 to 10
2.	Signal names . . . . .	page 11 to 18
3.	Block Diagram . . . . .	page 19
4.	Assembly Instructions . . . . .	page 20
5.	Troubleshooting Hints . . . . .	page 21 to 23
6.	Adjustment Instructions . . . . .	page 24 to 30
7.	PCB Application in different CardioServ Versions . . . . .	page 31
8.	Spareparts List . . . . .	page 33 bis 36
9.	Technical Specifications . . . . .	page 37 bis 42
10.	Pcb. Digital 388 028 09 . . . . .	Appendix 1
11.	Pcb. Digital 388 031 06 . . . . .	Appendix 2
12.	Pcb. Digital 388 031 05 . . . . .	Appendix 3
13.	Pcb. Analog 388 028 10 . . . . .	Appendix 4
14.	Pcb. Analog 388 031 10 . . . . .	Appendix 5
15.	Pcb. Power Supply 930 117 18 . . . . .	Appendix 6
16.	Pcb. Pacemaker partly assembled . . . . .	Appendix 7
17.	Pcb. Pacemaker 388 028 74 . . . . .	Appendix 8
18.	Pcb. Potential Separation 388 030 88 . . . . .	Appendix 9
19.	Wiring Diagram 101 117 01 / 02 / 09 . . . . .	Appendix 10
20.	Wiring Diagram 101 117 03 / 04 / 05 . . . . .	Appendix 11
21.	Car-Holder for CardioServ 202 307 01 . . . . .	Appendix 12
22.	Master Record Index . . . . .	Appendix 13

## 1. Functional Description

The Therapy System CardioServ comprises a defibrillator, a monitor ( LC display 240 \* 320 pixels ), a printer and a power supply unit. These components are arranged in a casing. It is a compact unit, which can also be implemented using a 220-V power supply system, even without a battery. The unit has a 140-mA standard charge for the battery. Additional versions available have an external pacemaker and an SpO<sub>2</sub> module, separately or combined.

The Defi is switched on by turning the rotary switch to position ECG or energy levels 2 ... 360 J.

Following a self-test, which tests

- RAM data store
- ROM program memory
- Watchdog / Reset
- EEPROM permanent memory
- Display all Pixels on / off
- Time base of internal quartz crystal
- Energy charge from last shock
- Defective charge or Defib keys
- Position of energy selector switch
- Connected electrodes
- Pacemaker
- SpO<sub>2</sub>

the Defi is ready for operation. If an error arises, depending on the severity of the error, the Defi is either on standby or not. Serious errors are watchdog, ROM or RAM errors. In the case of minor errors, an error message appears, but the instrument can be rendered operative by pressing a key.

The Defi shock is delivered via exchangeable paddles. This permits the use of different electrodes, i.e., in addition to conventional electrodes, large-surface adhesive electrodes or internal electrodes, for example. Adhesive electrodes are used for the pacemaker. The CardioServ also has Start and Trigger buttons on the instrument.

The charge operation is actuated with the Charge key on the apex electrode (for internal and adhesive electrodes on the instrument). During the charge operation, the accumulation of energy is registered and when the selected energy level has been reached the energy to be released in the 50-ohm resistor appears on the display. An acoustic and visual display on the monitor indicates the configured energy level has been reached. When the battery is fully charged or with line-power operation, charging at 360 J should take about 8 s.

The buttons on both electrodes are pressed to trigger the defibrillation pulse (safety two-hand operation). This means that the button on the apex electrode (or the corresponding one on the instrument) has two functions:

- start the charge operation
- trigger defibrillation (together with the Trigger button)

The type of electrodes connected are recognized and, in the case of internal electrodes, the energy is automatically limited to 50 joules.

Following defibrillation, the energy released is displayed on the screen for 6 seconds. Documentation output is by the printer.

The defibrillator has two operating modes:

- \* synchronized defibrillation (cardioversion)
- \* unsynchronized defibrillation (emergency therapy)

Cardioversion (operating mode Sync) can proceed either via the paddles or a patient cable. The operating mode Sync appears on the display and Sync marks are superimposed on the display. Following cardioversion, for safety reasons, the CardioServ switches automatically over to the unsynchronized operating mode, i.e., before each cardioversion the instrument must be switched over to Sync.

ECG signal pick-up is possible via the paddles or ECG cable. ECG is output on the LCD display or by the built-in 50-mm printer.

A safety device ensures the discharge of the storage capacitor in the following situations:

- when defibrillation is not triggered 30 seconds after storage operation is complete,
- 0.2 seconds after triggering defibrillation,
- immediately after switching back the energy selector switch during or after storage,
- when the battery voltage drops below a certain level,
- when the selected energy level still has not been reached after 32 seconds.

When switching to a higher energy level, the energy deficit is supplied by recharging.

The pacemaker is a transthoracic pacemaker operating in the Demand and Fix mode (over 180 P/min overpace).

- The frequency is between 0...300 P/min
- The current is between 0...200 mA
- Adjustment is in 5-unit intervals.

The pacemaker is operated using 6 keys on the control panel. One key switches the pacemaker on, one key changes the operating mode, the other keys are used to adjust the frequency and current. The defibrillator cannot be disabled as long as the pacemaker is switched on. The adhesive electrodes of the Defi are used to release a pulse. This ensures a rapid switch from defibrillator to pacemaker and back again.

SpO<sub>2</sub> is based on one of the modules manufactured by Nellcor (MP204). This module is connected to the PCB Pacemaker via the PCB Potential Separation and controlled by the latter. Measurement of the oxygen saturation level and pulse frequency is made by this PCB.

Pulse frequency from 0...250 P/min

Saturation from 0...100%

We display the plethysmogram, saturation level and, when selected, the pulse frequency instead of the HR.

#### **Control options:**

A battery charge check is performed under load application during the charge operation. The message "charge battery" on the display means the battery needs recharging. At least five discharges at maximum energy are still possible.

Performance check of defibrillation is performed by discharging into the integrated 50-ohm resistor.

Monitoring battery charge via the connector for power cord by the yellow LED.

The monitoring function of the CardioServ records and monitors the ECG and oxygen saturation level. The cardiac-action potentials are fed to the ECG preamplifier via a 5-wire patient cable. Selection can be made from 12 leads. The ECG curve is displayed on the LC screen with an illuminated background. The heart rate obtained from the ECG signal is continuously displayed on the screen.

Communication between the user and the instrument is via the keyboard and the LC screen. These also serve to configure the default settings and functions (see User's Manual for further details).

**The PCBs have the following functions:**

**PCB Power Supply (930 117 18)**

This PCB supplies the 16-V instrument power supply. 10 A can be supplied for short-term capacitor charging. It is a primary clock-rated wide-range power supply system. It operates from 64 .. 260 V, at 45 ... 120 Hz. It is protected by 2 fuses in the primary circuit. One on the PCB and one in the feed line.

**PCB Analog (388 028 10)**

This PCB has several functional units.

The first unit accepts the ECG signals (patient input) received via the patient cable from the patient. In accordance with its task, it is connected to the patient by a cable. It is thus constructed as an application component isolated from other circuit components and the circuit reference (floating). The tasks of this circuit component are:

- amplification of ECG signals
- lead selection (lead selector)
- identification and reporting of electrode defects
- identification and reporting of patient lead connection
- digitization and serialization of ECG data
- potential separation of the application component from the remaining circuit components (floating)
- floating component power supply
- ECG trigger
- pacing-pulse detection

There are two of these circuit components on this PCB. Once for the ECG via paddles, once for the ECG via patient cable.

A further circuit component controls high-voltage generation for the Defi shock.

There is a high-voltage capacitor on this PCB. The functions are as follows :

- safety control during charging and shock release (1st path : hardware)
- controlling storage operation and energy release, independent of the CPU
- controlling the high-voltage generator
- elimination of interference in the operating elements in the Defi electrodes
- conversion of the voltage and current signals into a frequency to measure energy
- activation and control of the high-voltage transformer
- generation of the emergency cutoff signal in the event of inadmissibly high voltage
- measurement of the high voltage on the primary winding of the high-voltage transformer

- rectification of the high voltage
- protection of the high-voltage diode
- discharging the storage capacitor as a safety precaution
- measurement of the high voltage on the storage capacitor
- energy delivery via the high-voltage relay

A further unit organizes the entire instrument's power supply. This is provided by the built-in battery, the power supply unit or, optionally, from an external battery (e.g., ambulance power supply system). It has the following tasks:

- switching instrument on/off
- automatic cutoff when battery voltage is too low
- monitoring the "charge battery" threshold
- power source selection: battery or power supply system
- charging the built-in battery

### **PCB Digital (388 028 09)**

The PCB Digital is equipped with the microcontroller, which controls the entire instrument. It comprises a Motorola 68332 microcontroller, RAM, EPROM, EEPROM, watchdog and I/O ports. It has the following functions:

- safety control during charging and shock release (2nd path : software)
- displaying the ECG and text on the LC screen (delete bar procedure)
- processing the curves and texts for the printer controller for output onto paper
- 50 / 60-Hz filter (software)
- alarm management

Circuit component display :

- graphics controller to control the display, which ensures screen refresh
- display RAM
- generation of voltages for the LCD ( LCD power supply )
- generation of high voltage for background illumination
- Pal for chip select generation for display
- contrast adjustment

Circuit component real-time clock :

- real-time clock with lithium battery buffer and integrated quartz crystal
- Pal for chip select generation for clock

Circuit component 1-channel recorder :

- single-chip microcontroller to output data at thermal-array printhead
- monitoring printhead temperature
- motor control for paper transport
- switches off motor in the event of overload or underload

- generation of power supply voltages for the recorder
- processing of all signals for the printer microcontroller

Further circuit components on the PCB are :

- generator of 5 V from 12 V via an integrated switching controller, short-circuit-proof and protected against overvoltage by Transzorb diodes
- sound generator, frequency controlled by SW
- RS-232 interface, protected against overvoltage by Transzorb diodes. Interface integrated into microcontroller, external conversion of level from 5 V to +- 12 V.
- energy selector switch, binary coded, delivers its information to PCB Analog and to the microcontroller
- keyboard interface to connect the membrane keypad
- control of LEDs on the keypad
- watchdog
- reset generation, when +5-V supply drops below 4,65 V or exceeds 5,5 V
- RAM buffering by battery when instrument is switched off
- options interface for pacemaker and SpO2.

#### **PCB Pacemaker (38802874)**

The CONTROLLER block acts as a microcontroller system. STIMCON is equipped with the parallel-to-serial converter, the analog protective circuits and some glue logic. Control of the optoelectronic coupler, protective circuit against DC stimulation with integrated powerdown of the floating part, DC/AC conversion and the floating part are hidden behind the FLOATING part. The function of each connection is shown in a table.

The chips required for the microcomputer system are listed on the CONTROLLER sheet. Excluded is the address decoder which, for clarity's sake, is included in the STIMCON plan.

The main component is the microcontroller 8032 (Z501). Contrary to the 8052, it requires an external program memory store, here integrated into the circuit in the form of a 32-kbyte EPROM (Z502). An external, static RAM (Z8) is envisaged, also with a capacity of 32 kbytes. This provides sufficient memory space to save graphics displays, such as the plethysmogram, generated by the SpO2 module.

Integrated into the system is a watchdog, MAX695 (Z1). This must be reset at least every 1.6 s since, otherwise, a processor reset of 50 ms duration is released by pins 15 and 16 (RESET, RESET\_). A reset is also released when the supply voltage at pin 9 (PFI) drops below 4.75 V, since under this level the integrated circuits no longer operate within their specific dynamic ranges, thus failing to ensure reliable function.

Every pacemaker circuit input is equipped with a  $10k\Omega$  resistor after +5V supply. This is a protective measure if the feed-in circuits can adopt high-impedance states. This would lead to uncertain levels on the input cables, which could result in the succeeding components showing unwanted behavior.

## 2. Signal names

Interface Signals PCB DIGITAL <-> PCB ANALOG Connector AA

Note: I/O from processor 68332 viewpoint

Signal Name	I/O	Meaning	Level	Polarity
ENSELD-ENSELA	I/O	Selection of energy level or device On/Off or position ECG DCBA HHHH -> locked mechanically, reaction to error HHHL -> Device off HHLH -> ECG HHLL -> 2J HLHH -> 5J HLHL -> 7J HLLH -> 10J HLLL-> 20J LHHH -> 30J LHHL -> 50J LHLH -> 70J LHLL -> 100J LLHH -> 150J LLHL-> 200J LLLH -> 300J LLLL -> 360J	<del>XXX</del>	----
TACHARGE	----	Charging signal from keypad	+5V	high-active
TADISCHARGE	----	Discharge signal from keypad	+5V	high-active
TAENABLE	O	Release signal for shock release via keypad	+5V	high-active
CHARGE_	I	Charging OR signal from keypad and paddle	+5V	low-active
DISCHARGE_	I	Discharge OR signal from keypad and paddle	+5V	low-active
PAPRINT	I	Triggers printout	+5V	high-active
ENERGY-SELECTED	O	L=enables HV generator start H=switches HV generator and HV off Reset position or lead interrupt = high	+5V	low-active
Release	I	Message D relay on	+5V	high-active
SYNC_	O	Triggers Defi shock in SYNC mode off. Is low in normal mode. Prevents shock release when high	+5V	low-active
RESDEF_	O	Reset Defi	+5V	low-active
AIMESS	O	Enables pick-up of C relay, toggles between current and voltage measurement L=C relay off, voltage reading H=C relay on, current reading	+5V	high-active
RELE	O	Switches E relay on, thereby enabling pick-up of C relay L=E relay not picked up H=E relay picked up	+5V	high-active
LOADC	O	Processor via LOADC HV generator L=Generator off H=Generator on	5V	high-active

Signal Name	I/O	Meaning	Level	Polarity
GEST	I	With a L->H change shows when the generator was started and with H->L change when the selected energy level is reached. (Detection by oscillator reverse voltage)	+5V	high-active
ELEC1 - ELEC3	I	Defi Electrode Identification 321 HHH -> External Electrode HLL -> Adhesive Electrode LHL -> Internal Electrode Other	+5V	----
U_F	I	Current or Voltage Measurement Converted into Frequency	+5V	----
PWON_	----	Stop Signal for Electronic Relay (Device-On/Off)		low-active
BATLOW_	I	Battery Undervoltage	+5V	low-active
BATLAD_	I	Battery Charge Display	+5V	low-active
12VBAT	----	12-V Standby Power Supply for Static RAMs	+12V	----
LINE	----	Line-Power Operation Display	+12V	high-active
MISOEXT	I	Master-In Slave-Out (Buffered)	+5V	----
MOSIEXT	O	Master-Out Slave-In (Buffered)	+5V	----
SCKEXT	O	QSPI Serial Clock (Buffered)	+5V	----
PCS1EXT_	O	Peripheral Chip Select1 (Buffered)	+5V	low-active
PCS2EXT_	O	Peripheral Chip Select2 (Buffered)	+5V	low-active
GND	----	Ground Reference Level 5-V Supply	0V	----
+5V	----	+%V Power Supply	+5V	----

## Interface Signals PCB DIGITAL &lt;-&gt; PCB ANALOG Connector AB

Signal Name	I/O	Meaning	Level	Polarity
GND_ANA	----	Ground	0V	----
+12V_ANA	----	+12V Power Supply	+12V	----
+5V_ANA	----	+5V Power Supply	+5V	----

## Interface Signals PCB DIGITAL &lt;-&gt; Printhead/Motor Interface Connector AC

Note: I/O from processor 80C51 or motor control viewpoint

Signal Name	I/O	Meaning	Level	Polarity
KLA,	O	Printhead Latch Signal	+5V	high-active
KCLK	O	Printhead Clock Signal	+5V	L->H
KDAT	O	Printhead Data Signal	+5V	----
KSTRB	O	Printhead Strobe Signal	+5V	high-active
TERM1	I	Printhead Thermistor Connector1	+5V	----
TERM2	I	Printhead Thermistor Connector2	+5V	----
GNDJ	----	Ground Thermal-Array Printhead	0V	----
+5V	----	++5V Power Supply Thermal-Array Printhead	+5V	----
+22V	----	+22V Power Supply Thermal-Array Printhead	+22V	----
VDD1	O	Release Signal for Thermal-Array Printhead Overvoltage	+22V	high-active
MOTOR+	O	Motor Connector +	var	----
MOTOR-	O	Motor Connector -	var	----
SENSE	I	Motor Current Sensor Signal	var	----

Interface Signals PCB DIGITAL <-> LCD Connector AD (HOSIDEN HLM6323)
--

Note: I/O from Graphics Controller viewpoint

Signal Name	I/O	Meaning	Level	Polarity
V0	O	Contrast Voltage	var.	----
VEE	----	Negative Power Supply for LCD	-23V	----
D3-D0	O	Serial Pixel Data	+5V	----
M	O	LCD Drive Signal (AC Signal)	+5V	----
VSS	----	Ground	+0V	----
VDD	----	+5V Power Supply	+5V	----
CP2	----	Display Data Shift Clock	+5V	H->L
CP1	----	Display Data Latch Signal	+5V	H->L
S	O	Frame Signal	+5V	high-active
DISP_OFF_	O	Display ON/OFF H=Display ON L=Display OFF	+5V	low-active
FG	----	Frame Ground	0V	----

Interface Signals PCB DIGITAL <-> LCD Connector AE (HOSIDEN HLM8619)
--

Note: I/O from Inverter viewpoint

Signal Name	I/O	Meaning	Level	Polarity
FL HOT	O	High Voltage for CCFT	---	----
FL GND	O	High-Voltage Ground	---	----

Interface Signals PCB DIGITAL<->LCD Connector AF (OPTREX DMF50174NF- FW)

Note: I/O from Graphics Controller viewpoint

Signal Name	I/O	Meaning	Level	Polarity
FLM	O	Frame Signal	+5V	high-active
LP	----	Data Latch Signal	+5V	H->L
CP	O	Clock Signal for Shifting Serial Ddata	+5V	H->L
M	O	Alternate Signal for LCD Drive	+5V	----
VADJ	----	Contrast Signal	var	----
VCC_	----	+5V Power Supply	+5V	----
VSS	----	Ground	0V	----
VEE	----	-23V Power Supply	-23V	----
D3-D0	O	Display Data	+5V	----
DISP_OFF__	O	Display ON/OFF H=Display ON L=Display OFF	+5V	low-active

Interface Signals PCB DIGITAL<->LCD Connector AG (OPTREX DMF50174NF- FW)

Note: I/O from Inverter viewpoint

Signal Name	I/O	Meaning	Level	Polarity
FL HOT	O	High Voltage for CCFT	----	----
FL GND	O	High-Voltage Ground	----	----

**Interface Signals PCB DIGITAL <-> Keypad Interface Connector AH**

Note: I/O from processor 68332 viewpoint

Signal Name	I/O	Meaning	Level	Polarity
KBROW0_ KBROW1_ KBROW2_ KBROW3_	O	Keypad Matrix Columns	+5V	low-active
KBCOL0_ KBCOL1_ KBCOL2_ KBCOL3_	I	Keypad Matrix Lines	+5V	low-active
TAENABLE	O	Release Signal for Key Charge/Defib + Defib	+5V	high-active
TACHARGE	I	Charge/Defib Key	+5V	high-active
TADISCHARGE	I	Defib Key	+5V	high-active
ON_OFF_	O	On/Off Key	+5V	----
GND	----	Ground	0V	----
LINE	----	Line LED	+12V	high-active
BAT	O	Battery Undervoltage LED	+5V	high-active

Interface Signals PCB DIGITAL <-> RS-232 Interface Connector AI
---

Note: I/O from processor 68332 viewpoint

Signal Name	I/O	Meaning	Level	Polarity
TxD_E	O	Transmit Signal (RS232))	+/-12V	----
RxD_E_	I	Receive Signal (RS232)	+/-12V	----
CONNECT1	----			
CONNECT2	I	Identifies RS232 Lead is Plugged in	+5V	high-active
GND_	----		0V	----
12VEX	----	+12V Power Supply (Optional)	+12V	----

Interface Signals PCB DIGITAL<->PCB Pacemaker AJ
--

Note: I/O from processor 68332 viewpoint

Signal Name	I/O	Meaning	Level	Polarity
PACEEN	O	Stimulator Release	+5V	high-active
SPO2_	O	Chip Select for SPO2 Module	+5V	low-active
RWE_	O	Read/Write Control Signal for Options Board L=Write H=Read	+5V	----
ERRSPO2	I	Error Message SPO2	+5V	high-active
ERRSTIM	I	Error Message Stimulator	+5V	high-active
EKTRIG_	O	ECG Trigger	+5V	low-active
AE0-AE11	O	Buffered Addresses 68332	+5V	----
DE_0-DE_7	I/O	Buffered Data 68332	+5V	----
GND	----	Ground	0V	----
+5V	----	+5V Power Supply	+5V	----
+12V	----	+12V Power Supply	+12V	----


## Interface PCB Pacemaker&lt;-&gt;PCB Potential Separation

Signal Name	I/O	Meaning	Level	Polarity
TX	O	Transmit	+5V	high-active
RX	O	Receive	+5V	high-active
EKG_TRIG	1	ECG Trigger Signal	+5V	high-active
RESOXY_	0	Reset	+5V	low-active
GNDF	----	Ground	0V	----
+5VF	----	+5V Supply	+5V	----
12V	----	12V Supply	12V	----

## Interface Signals PCB Potential Separation &lt;-&gt;Nellcor MP204

Signal Name	I/O	Meaning	Level	Polarity
TX	O	Transmit	+5V	high-active
RX	O	Receive	+5V	high-active
EK_SYN	1	ECG Trigger Signal	+5V	high-active
RESET_	0	Reset	+5V	low-active
GNDF	----	Ground	0V	----
+8.5V	----	+8.5V Supply	+8.5V	----
+5VF	----	+5V Supply	+5V	----
-8.5V	----	-8.5V Supply	-8.5V	----

### 3. Block Diagram CardioServ

## 4. Assembly Instructions

### **General Information:**

Before opening the CardioServ switch the instrument off, disconnect power cord plug, remove battery and implement ESD protection measures. Return replacement PCB in ESD packaging only

### **Opening the instrument:**

Undo the 9 screws on the bottom of the device, place upper case shell to the left next to the instrument.

### **Removing the PCB Digital:**

The PCB Digital is located in the upper shell. Disconnect the 40-pin flat strip cable and the 2-pin 12-V power supply cord. The recorder, plugged into the lower case shell, can be removed and the flat strip cable connected to the thermal-array printhead and motor disconnected.

Lay the upper shell "face down", remove the screening plate fixed in place by a screw and 5 plastic catches. Remove the insulating card. Undo the 8 PCB Digital fastening screws (leave the 4 display screws in place) and remove the rotating knob of the energy selector switch. **Carefully** lift up the PCB and pull out the membrane keypad connector on the left. The PCB Digital can now be removed along with the display screen.

### **Removing the display screen:**

- Disconnect plug-in connector to display illumination.
- **Carefully** pull the zero power connector safety catch of the display-membrane keypad connector upwards and disconnect membrane keypad lead.
- Remove the 4 display fastening screws

### **Removing the PCB Analog:**

Disconnect 40-pin flat strip cable, connection cable to the battery, power supply unit and PCB Digital . Detach connector to the ECG socket and disconnect connections to Defi socket. Remove ECG amplifier cover with 3 screws and 2 plastic catches. Unscrew the 10 screws on the PCB Analog (4 capacitors, 2 high-voltage oscillators) and 3 hexagonal columns.

## 5. Troubleshooting Hints

### Device cannot be switched on during line-power operation:

- 220-V power supply via power supply unit OK ? Measure at connector CA.
- Connector 220-V plug to power supply unit OK ?
- Fuse in feed line and on power supply unit OK ?
- Lead line-power supply unit -- PCB Analog OK ?
- Lead PCB Analog -- PCB Digital (2-pin power supply 16V) OK ?
- Lead PCB Analog -- PCB Digital (40-pin control) plugged in correctly?
- 16 V on PCB Analog at connector ? no --> power supply unit defective or short-circuit ( power supply unit tackeder )
- Fuse SI500 (16V) on PCB Digital OK ?
- Fuse SI502 (5V) on PCB Digital OK ?
- "Device On Signal" PWON\_ present and 16V available on PCB Digital : yes --> Digital defective
- PCB Analog or PCB Digital defective ?

### Device cannot be switched on during battery operation:

- Battery charged?
- Lead battery contacts -- PCB Analog OK ?
- Fuse for battery on PCB Analog OK ?
- Connector battery contact and connector power supply PCB Digital confused
- Fuse SI500 (16V) , SI502 (5V) on PCB Digital OK ?
- PCB Analog defective ?

### Device does not print out:

- Paper inserted correctly? (thermosensitive side facing up)
- Flat strip cable to printer OK ?
- Mechanical failure in paper transport, paper feed ?
- Without paper platin runs approx. 1 sec
- Flap not completely closed
- Print key on membrane keypad does not respond when pressed (perform key test)
- Does motor run ? or printhead voltage there ? or movement on the data lines to printhead ? no --> PCB Digital defective ?

**No display on LCD**

- Contrast poorly adjusted ? (change contrast)
- Contrast keys defective ? (perform key test)
- Connector for illumination or data OK ?
- Illumination defective ?
- Fuse SI500 (16V) , SI501(background illumination) , SI502 (5V) on PCB Digital OK ?
- Display defective ?
- PCB Digital defective ?

**No ECG on display via ECG cable or paddles**

- Connector to ECG socket OK ?
- Connector to paddles OK ?
- Lead PCB Analog -- PCB Digital (40-pin control ) OK ?
- ECG via paddles or ECG cable OK : --> PCB Analog defective
- PCB Analog or PCB Digital defective ?

**Battery does not charge:**

- Yellow LED must glow when battery is being charged
- Fuse SI501 on PCB Analog defective ?
- Check connector to battery
- Voltage on power supply unit under 16 V

**Self-test:**

F1 + F5 simultaneously to go to configuration and with Print print out the current error messages

In configuration with F3 + F4 key test and carry out printer adjustment (Print key)

In configuration with F3 + F5 restore default (factory setting)

Display test chessboard pattern at power up

Error in time base : The clock must be set once or PCB Digital defective

Error in energy storage : Adjust energy storage / measurement or PCB Analog defective

Error RAM, ROM, Display-RAM, permanent memory, watchdog : PCB Digital defective.

Error Charge or Discharge actuated (key sticking).

**Troubleshooting hints for PCB Pacemaker**


No adjustments can be made on the PCB Pacemaker.

When inserting this PCB make sure it is plugged in correctly; especially in the case of the red pacemaker cable from connector C it is important to observe the correct polarity on connectors L and M. Mixing these up will result in the pacemaker pulse not being transmitted to the paddles.

**Troubleshooting hints for PCB SPO2**

When the PCB SPO2 is inoperative or operates defectively check the 5V floating voltage. The voltage must lie between 4.8-4.9V.

## 6. Adjustment Instructions

Lpl Digital 388 028 09 R  
PCB Digital 388 028 09 R

## Motor Adjustment

What to adjust or to check?	What to measure with?	How to activate testing mode?	Where to turn?	How much and exact?	What else to note?
Adjust motor speed under small load (recorder flap open) on stationary motor pinion.	Flash at motor pinion with stroboscope (LED).	1) Call up configuration (Press F1 + F5 simultaneously) 2) Call up key test (Press F3+F4 keys simultaneously) 3) Press Print key	Speed selector R521	Stroboscope frequency = 181.07 Hz  When a square wave generator activated by LED is used, a low key ratio should be observed.	Motor should be at room temperature. The testing mode terminates automatically after a few seconds. This can be suppressed by pressing any key.
Adjust motor speed under normal load (recorder flap closed) and small load (recorder flap open) so that there is no difference in RPM.	"	"	Compensation adjuster R520		If the RPM changes after compensation adjustment, this should be corrected once more only with R521.

**No Stroboscope Available:**

What to adjust or to check?	What to measure with?	How to activate testing mode?	Where to turn?	How much and exact?	What else to note?
Adjust motor speed under small load (recorder flap open wide enough so that paper is just transported)	Remeasure grid printout. min. 200 mm	1) Call up configuration (Press F1 + F5 simultaneously) 2) Call up key test (Press F3+F4 keys simultaneously) 3) Press Print key	Speed selector R521	< 1%	Motor should be at room temperature. The testing mode terminates automatically after a few seconds. This can be suppressed by pressing any key.
Adjust motor speed under normal load (recorder flap closed)	"	"	Compensation adjuster R520	"	The adjust-procedure is complete when the RPM is 25 mm/s under small load and normal load.

**Heating Time Adjustment**


What to adjust or to check?	What to measure with?	How to activate testing mode?	Where to turn?	How much and exact?	What else to note?
Heating time	Connect oscilloscope to MP2/1 and MP1/1 (GND)	Press Print key in main menu	Heating time adjuster R522	500us +/- 1% referred to 22°C printhead temperature	Adjustment ensues at 25 mm/s. Important: printhead warms up during printing

**5-V Overvoltage Cutoff Threshold**

What to adjust or to check?	What to measure with?	Where to measure?	Where to turn?	How much and exact?	What else to note?
5-V power supply monitoring	Voltmeter. Input impedance >= 1MΩ	Between MP3/1 and GND	R523	2,720V <1%	

**Control of the Power Supply Voltages on PCB Digital**

What to adjust or to check?	What to measure with?	Where to measure?	How much and exact?	What else to note?
5-V power supply	Voltmeter.	Between +5V and GND	5,0V +/-2%	
+22-V power supply	Voltmeter.	Between +22V and GND (C504)	+22V +/-5%	Recorder stationary
-23-V power supply	Voltmeter.	between -23V and GND	-23V +/-5%	

Lpl Analog    388 028 10 R  
PCB Analog    388 028 10 R

## Energy Adjustment

What to adjust or to check?	What to measure with?	Where to measure?	Where to turn?	How much and exact?	What else to note?
Energy released in 50-ohm resistor	Defi tester		R509	.100J +-2J	Turn R508 left to minimum. Selector switch to 100J Defibrillate in Defi tester
Adjust charging time	Timer		R512	< 8s	Starting from the left adjust R512 Selector switch to 360J Defibrillate in Defi tester
Adjust display of energy to be released in 50-ohm resistor	Defi tester	Check screen display	R508	Energy displayed in Defi tester +-2J	Selector switch to 100J Defibrillate in Defi tester
Adjust display of energy released	Defi tester	Check screen display (50 ohms)	R507	Energy displayed in Defi tester +-2J	Selector switch to 100J Defibrillate in Defi tester
Check all energy levels	Defi tester		R509, R508, R507	Selected energy level +-5% or 2J Higher value permissible	Adjust all energy levels on energy selector switch. If required, correct previous settings

**EC Amplifier Adjustment**

What to adjust or to check?	What to measure with?	Where to measure?	Where to turn?	How much and exact?	What else to note?
Offset EC amplifier cable	Optical scale	Display	R502	Baseline offset < 1 mm	Toggle amplification between 0.5, 1 and 2
Amplification EC amplifier cable	Optical scale	Printer	R501	10 mm +-0.5 mm	ECG pickup 1 mV
Offset EC amplifier paddles	Optical scale	Display	R510	Baseline offset < 1 mm	Toggle amplification between 0.5, 1 and 2
Amplification EC amplifier paddles	Optical scale	Printer	R511	10 mm +- 0.5 mm	Connect ECG pickup with paddles Amplitude 1 mV

## 7. PCB Application in Different CardioServ Versions

	CardioServ 101 117 01  Version 2.0/2.1	CardioServ 101 117 05  Version 3.x With Pacemaker and SPO2	CardioServ 101 117 03  Version 3.x With Pacemaker without SPO2	CardioServ 101 117 04  Version 3x Without Pacemaker with SPO2	CardioServ 101 117 01  Version 3x Without Pacemaker without SPO2
PCB Digital 388 028 09 Exch. 389 003 96	X				
PCB Digital 388 031 06 Exch. 389 004 06					X
PCB Digital 388 031 05 Exch. 389 004 08		X	X	X	
PCB Analog 388 028 10 Exch. 389 003 97	X				
PCB Analog 388 031 10 Exch. 389 004 07		X	X	X	X
PCB Power Supply 930 117 18	X	X	X	X	X
PCB Pacemaker partially equipped 388 031 14				X	
PCB Pacemaker 388 028 74		X	X		
PCB Potential Separation SPO2 388 030 88		X		X	
PCB SPO2 930 117 26 AT 930 117 36		X		X	

### Pcb. usage on different CardioServ

	CardioServ 101 117 01  Version 2.0/2.1	CardioServ 101 117 05  Version 3.0/3.01 With Pacemaker and SPO2	CardioServ 101 117 03  Version 3.0/3.01 With Pacemaker without SPO2	CardioServ 101 117 04  Version 3.0/3.01 Without Pacemaker with SPO2	CardioServ 101 117 01  Version 3.0/3.01 Without Pacemaker without SPO2	CardioServ 101 117 09 USA Version  Version 3.0/3.01 Without Pacemaker without SPO2
Pcb. Digital 388 028 09 Exch389 003 96	X					
Pcb. Digital 388 031 06 Exch389 004 06					X	X
Pcb. Digital 388 031 05 Exch389 004 08		X	X	X		
Pcb. Analog 388 028 10 Exch389 003 97	X					
Pcb. Analog 388 031 10 Exch389 004 07		X	X	X	X	X
Pcb. Power Supply 930 117 18	X	X	X	X	X	X
Pcb. Pacemaker partly assembled 388 031 14				X		
Pcb. Pacemaker 388 028 74		X	X			
Pcb. Potential Separation Spo2 388 030 88		X		X		
Pcb. SPO2 930 117 26 Exch930 117 36		X		X		

### Pcb. usage on different CardioServ Version

	CardioServ Esaote 101 117 02  Version 2.0 / 2.1	CardioServ Esaote 101 117 08  Version 3.x With Pacemaker and SPO <sub>2</sub>	CardioServ Esaote 101 117 06  Version 3.x With Pacemaker without SPO <sub>2</sub>	CardioServ Esaote 101 117 07  Version 3.x Without Pacemaker with SPO <sub>2</sub>	CardioServ Esaote 101 117 02  Version 3.x Without Pacemaker without SPO <sub>2</sub>
Pcb. Digital 388 028 09 Exch.389 003 96	X				
Pcb. Digital 388 031 06 Exch.389 004 06					X
Pcb. Digital 388 031 05 Exch.389 004 08		X	X	X	
Pcb. Analog 388 028 10 Exch.389 003 97	X				
Pcb. Analog 388 031 10 Exch.389 004 07		X	X	X	X
Pcb. Power Supply 930 117 18	X	X	X	X	X
Pcb. Pacemaker partly assembled 388 031 14				X	
Pcb. Pacemaker 388 028 74		X	X		
Pcb.Potential Separation Spo <sub>2</sub> 388 030 88		X		X	
Pcb. SPO <sub>2</sub> 930 117 26 Exch.930 117 36		X		X	

## 8. Spareparts List

### Partnumber      Description

101 117 ..	CardioServ SCP 910
101 117 80	Exchange CardfioServ SCP 910
227 432 05	Service manual CardioServ SCP 910
227 432 01	Operation Manual Version 2 ( german )
227 446 01	User Manual Version 3 ( german )
227 446 04	User Manual Version 3 ( USA )

### Defi paddles

A	217 304 03	Defibrillation Electrode, external, with recorder start, for CardioServ
A	217 308 01	Electrode Lead f.Defi-El.“internal“ without contact insert, for SCP 913
	926 074 97	Handle

### Recorder

	384 015 81	Exchange Recorder preman. f. CardioServ
B	303 439 76	Thermal printhead
	303 434 78	Transportation roller
	918 083 49	Motor 18 V-
	480 159 46	Gear 10 mm diameter
	480 159 47	Gear 24 mm diameter
	480 159 48	Gear 40 mm diameter
	388 010 56	Pcb. Motor
A	432 522 30	Flap, Recorder
A	430 517 84	Sticker in Flap

### Printed circuit boards for Version 2.0 / 2.1

B	388 028 09	Pcb. Digital      (CardioServ SCP 910)
B	388 028 10	Pcb. Analog      (CardioServ SCP 910) Info: w/o pcb HF-protection (433 671 65) (Pls. order additional two spacer!)
	927 229 71	Spacer 6 mm
	927 229 72	Spacer 17,7 mm
	903 448 29	Capacitor 32 uF +-5% 5,2 KV-
B	390 001 07	Keypad without Pacemaker
	389 003 96	Exch. pcb. Digital CardioServ SCP 910
	389 003 97	Exch. pcb. Analog CardioServ SCP 910

**Housing**

C	303 442 85	Upper shell
C	303 443 61	Upper shell ( Esaote ) Version
C	303 444 46	Upper shell USA Version
	430 517 93	Sticker "operation information", english
	430 517 94	Sticker „operation information“, french
	430 518 09	Sticker „operation information“, italia
	430 518 10	Sticker „operation information“ spanish
	430 518 11	Sticker „operation information“ russian
	432 522 08	Lower Shell
	430 518 05	Sticker „Attention, high voltage...“ d
	430 518 06	Sticker „Attention, high voltage...“ e
	430 518 07	Sticker „Attention, high voltage...“ f
	430 518 12	Sticker „Attention high voltage...“ it
	430 518 13	Sticker „Attention high voltage...“ sp
	430 518 14	Sticker „Attention high voltage...“ ru
	504 657 54	Contact —> Station
	432 522 27	latch
	432 522 35	Filler Plug
	924 017 22	Foot
	482 035 19	Knob for energy selector
	923 096 69	Mounting Base
	829 074 11	SCOTCH MOUNT, WHITE

**In-/Output connector**

C	303 442 99	ECG Inputconnector completely wired for Cardio Serv_
C	303 441 25	Connector for defi paddles (Employment only by Version 2.0)
	915 417 95	Mains plug, 3 wire
	929 166 33	Ferrit Core for Round Cable 7,8-8,5mm

**Power supply**

C	930 117 18	SWITCHING POWER SUPPLIES
B	303 440 30	Battery Pack for SCP 910/913/915/922 12 V / 1,2 Ah
	908 122 28	Fuse block
A	912 084 53	FUSE FAST 5.0A UL

**LCD- Display**

C    930 117 17    LCD Graphic Display

**Accessories**

- 223 287 01    Patient Cable,3-w.2.2 m,not caut. proof clips,IEC (4th gen.a.VICOM a CardioServ)
- 223 288 02    Patient Cable,3-wire, 2.2 m, caut. proof clips, USA (4.G.a. VICOM .a. CardioServ)
- 226 130 02    Recording Chart,T,55mm, roll 40 m for CardioServ
- 931 098 68    Accessory Bag for CardioServ
- 931 098 88    Instrument Bag for CardioServ, nylon fabric, orange, with belt

**Software**

- 443 295 02    EPROM Cardioserv V 2.1
- 443 295 05    EPROM Cardioserv V 3.02

**Printed Circuit boards for version 3.x**

- B    388 031 06    Pcb.Digital  
(usable without pacemaker and SPO2 )  
389 004 06    Exchange pcb. digital
- B    388 031 05    Pcb. Digital with RAM Extended  
(usable with Pacemaker and SPO2 )\_  
389 004 08    Exchange pcb Digital with RAM Extended
- B    388 031 10    Pcb. Analog  
389 004 07    Exchange pcb. Analog

**Pacemaker**

- 388 028 74 P.C. Pacemaker  
388 031 14 pcb. partly assembled

930 117 26	Pcb. SPO2
930 117 36	Exchange pcb.SPO2
388 030 88	Pcb. Potential Separation for SPO2
303 443 58	Adapter Cable for SpO2 Sensors (NELLCOR) to CardioServ + EAGLE10xx
303 443 57	SPO2 Inputconnector kompletely wired for Cardio Serv

**Keypad**

390 001 07	Keypad without Pacemaker
390 001 08	Keypad with Pacemaker
390 001 31	Keypad with Pacemaker Esaote
390 001 28	Keypad witout Pacemaker Esaote

**Auto Holder for CadioServ**

504 658 18	Case
504 658 21	Insulating Device
416 118 28	CONTACT SPRING
504 658 23	Recording Upper
427 338 70	BOLT
504 658 22	Recording Lower
504 658 20	Handle
402 236 14	SOCKET
914 326 43	Pushbutton switch
427 338 69	ROD
388 011 38	Pcb. 12V- option
908 092 00	FUSE HOLDER
908 093 00	CAP
912 049 00	Fuse F 10 A
912 084 43	Fuse T 1,6 A
912 084 44	Fuse T 2,5 A (TR 5)
504 658 19	upper Foot
504 658 56	Lower Foot / red

## 9. Technical Specifications

### Operating modes

- \* non-synchronized (defibrillation at any time)
- \* synchronized (cardioversion)
- \* HR monitoring

### Energy selection

adjustable, energy to be delivered into 50 ohms  
displayed digitally

- \* selectable energy levels, energy to be delivered into a 50-ohm resistance (max. energy of 50 joules for internal defibrillation):

2  
5  
7  
10  
20  
30  
50  
70  
100  
150  
200  
300  
360 joules

- \* possible deviation from selected energy below values specified by IEC

### Energy storage

by means of capacitor, capacitor is charged from battery, from a 12-volt power source (emergency vehicle) or from the power line (95 to 240 V; 49 to 65 Hz); when capacitor charging is terminated buzzer sounds:

- \* capacitor charging time for energy setting of 360 J:
  - from power line: typically 8 s
  - from fully charged battery: typically 8 s
  - from partially discharged battery: typically 10 s (15 s max.), measured at least 5 minutes after 15 shocks of 360 joules each (for 200 joules typically 4 s)

### Defibrillation pulse

capacitor discharge via induction coil (damped serial resonant circuit), pulse shape resembling a sinusoidal halfwave with decay period:

- \* pulse duration for an external resistance of 50 ohms approx. 4 ms, measured from the beginning of the pulse to the intersection of the zero line and the inflection point of the trailing pulse edge
- \* in synchronized mode the defibrillation pulse is released approx. 40 ms after the R-wave trigger

### Discharge circuit

serial oscillating circuit in series with external resistance (patient):

- \* capacitance 32  $\mu$ F
- \* inductance 26 mH
- \* equivalent resistance 6.5 ohms

### Pulse output

isolated, no galvanic connection with enclosure, open-circuit and short-circuit-proof

- \* isolation test voltage 8 kV DC, type CF according to IEC requirements

## Safety discharge

capacitor discharge via internal load resistance:

- \* when the defibrillation shock is not triggered within 30 s after charging
- \* when the defibrillation shock is triggered, but the discharge circuit is interrupted, after approx. 0.2 s
- \* immediately when reducing the selected energy during or after charging
- \* when the selected energy is not reached, after 32 s
- \* in the event of technical malfunctions

## Test features

- \* pilot lamp for battery charging
- \* defibrillator test by discharging the stored energy into the internal 50-ohm load resistance; 3-digit display of the delivered energy
- \* warning on LCD when discharge circuit is interrupted (e.g., defibrillation paddle not applied)
- \* automatic defibrillator test on power up with display of error message, if applicable

## Synchronization

with ECG signal of either polarity:

- \* min. ECG amplitude for reliable triggering approx. 0.5 mV with QRS duration of 80 ms

## ECG signal input via paddles

ECG signal acquisition via defibrillator paddles, ECG trace displayed on LCD, automatic switching to ECG electrodes when patient cable is connected; cardioversion both via ECG patient cable and via paddles; differential input; isolated, class CF according to IEC, with overvoltage protection

- \* input voltage range  $\pm 4$  mV
- \* input impedance  $>1.5$  Mohm
- \* max. polarization voltage  $\pm 1$  V
- \* frequency response 2.2 Hz to 20 Hz (-3 dB)
- \* common-mode rejection  $>80$  dB
- \* patient leakage current: in normal condition  $<10$   $\mu$ A, in single-fault condition  $<50$   $\mu$ A
- \* detection of pacing pulses
  - pulse duration  $d_p >$  approx. 0.1 ms  $<2.0$  ms
  - pacemaker independent of polarity
  - pulse amplitude  $a_0 \pm 20$  to  $\pm 700$  mV
  - reverse-current pulse  $a_0 \pm 1$  mV
  - time constant  $t_0 = 25$  to 100 ms

## ECG signal input via patient cable

via ECG electrodes, automatic switching to ECG electrodes when patient cable is connected; cardioversion both via ECG electrodes and via paddles; differential input, symmetrically referred to N, isolated, class CF according to IEC; 7 standard leads selectable via lead selector; input with overvoltage protection (defibrillation-proof):

- \* input voltage range  $\pm 4.4$  mV for recorder,  $\pm 4$  mV for display
- \* input impedance  $>2.5$  Mohm for 10 Hz
- \* common-mode dynamic range  $\pm 3$  V

- \* differential DC voltage compatibility  $\pm 1$  V
- \* common-mode rejection (CMRR) R,L referred to N 65 dB, N referred to chassis >110 dB
- \* trigger sensitivity 250  $\mu$ V to 300  $\mu$ V
- \* bandwidth .5 to 100 Hz
- \* patient leakage current: in normal condition <10  $\mu$ A, in single-fault condition <50  $\mu$ A
- \* ground leakage current: in normal condition .5 mA, in single-fault condition 1 mA
- \* voltage resistance referred to circuit reference 4 kV
- \* detection of pacing pulses
  - pulse duration  $d_p$  > approx. 0.1 ms <2.0 ms
  - pacemaker independent of polarity
  - pulse amplitude  $a_0$   $\pm 20$  to  $\pm 700$  mV
  - reverse-current pulse  $a_0$   $\pm 1$  mV
  - time constant  $t_0$  = 25 to 100 ms

## Signal display

backlit LCD, 2-channel erase bar mode, calibration pulse in left-hand corner (for ECG)

alphanumeric presentation of alarm messages, sensitivity, leads, systole blinker, alarm limits, heart rate, energy and softkey labels

ECG freeze with simultaneous display of the current ECG and (for modules with SpO<sub>2</sub> function) plethysmogram at a smaller scale

ECG trace 1.5 times larger than on recording: with a sensitivity of 1 cm/mV a 1-mV signal has an amplitude of 1.5 cm on the display

- \* erase bar sweep speed 25 mm/s
- \* trace length in real-time mode 4.6 s
- \* display dimensions: 115 mm wide, 86 mm high

- \* resolution 320 x 240 pixels (pitch of .36 x .36 mm)
- \* displayed image can be rotated 180 °

## Signal transmission

signal input —> amplification —> signal sampling —> AD conversion —> digital processing —> LCD or recorder

- \* selectable sensitivity: .5-1-2 cm/mV (with max. sensitivity of 2 cm/mV a 1-mV input signal is 2 cm in amplitude on the recorder and about 3 cm on the LCD), amplitude limited to approx.  $\pm 2$  cm on the recorder and approx.  $\pm 3$  cm on the LCD
- \* signal sampling rate 1000 Hz

## Systole check

- \* heart symbol flashing on the LCD
- \* QRS beep (can be disabled)

## Heart-rate measurement

derivation of trigger pulses from the ECG of either polarity, adaptive trigger threshold, calculation of the average rate, storage of the result, 3-digit display on LCD, alarm limits to the left of the reading:

- \* measuring range 15 to 300 BPM
- \* digit height of HR reading 7.5 mm
- \* digit height of alarm limits 2.5 mm
- \* min. amplitude for reliable triggering >.25 mV for ECG signal with a QRS duration of 80 ms

## Alarm system

electronic release of alarm

- \* when the HR violates one of the set limits for at least 4 s: alarm tone sounds (can be disabled), message "Alarm, high HR" or "Alarm, low HR", recorder starts (if configured)
- \* when asystole is present for more than 4 s: alarm tone sounds, "Alarm, Asystole" is shown on LCD, recorder starts (if configured)
- \* when at least one of the selected electrodes drops off: audible signal (if patient cable is plugged in), message "Alarm, Electrode" on display
- \* alarm limits adjustment range 15 to 300 BPM (not overlapping)
- \* digit height of alarm limits 2.5 mm
- \* keys to cancel alarm and to silence alarm tone
- \* softkeys to adjust alarm limits

## Recording

delayed recording of the ECG stored in the signal memory (strip length 16 s, incl. 4-s history) in the event of an alarm plus alphanumeric annotations on the paper margin:

- heart rate
- lead
- filter
- date
- time
- paper speed
- cause of recording (defib, alarm, manual)
- selected energy [J into 50 ohms]
- delivered energy [J into 50 ohms]
- sync mark
- text (name of user/hospital/practice)

After the ECG recording, a patient ID sheet is printed indicating name, date of birth, user, comments, date, time, energy and alarm limits.

Direct writing with rectangular coordinates using thermorecording technology (printhead with electronically controlled thermal elements records on thermosensitive paper), baseline fixed at the center of the space available for recording of the ECG trace, grid imprint, roll paper, paper transport by electronically controlled DC motor, limited duration of transport

- \* number of recording channels 2
- \* paper width 55 mm
- \* roll diameter 60 mm max. (roll with 40 m of HELLIGE CONTRAST®)
- \* printhead resolution vertical 6 dots/mm, horizontal 24 dots/mm
- \* paper speed 25 mm/s ±5%
- \* paper transport after both manual and automatic start 16 s (incl. history of 4 s)

To prevent damage to the printhead, use only the original HELLIGE CONTRAST® chart paper.

## Memory

- \* storage of 10 or 40 ECG strips initiated by defibrillation or alarm with a length of 16 seconds and a history of 4 seconds each, incl. a full report
- \* storage of the 40 or 80 most recent actions (e.g., power on, power off, alarms, defibrillation energy) incl. date and time

**SpO<sub>2</sub>**

- \* saturation: 0 to 100 %, in 1-% increments
- \* frequency: 0 to 250 BPM, in 1-BPM increments
- \* alarm limit: off, 15 to 100 %
- \* display of plethysmogram
- \* C-LOCK ECG synchronization
- \* integration time: 4, 8 and 12 seconds
- \* measuring accuracy:      70 to 100%  $\pm 2$  digits  
                                50 to 69%  $\pm 3$  digits  
                                pulse display 1.2% or  
                                 $\pm 1$  BPM

**Pacemaker**

- \* operating modes: demand, fixed-rate, overpace
- \* pacing rate: 30 to 150 BPM (overpace up to 300 BPM)
- \* pacing current:      0 to 200 mA (for 500 Ohms), voltage up to 120 V
- \* pulse width (demand, fixed-rate): 40 ms  
(overpace):            20 ms
- \* pulse shape: monophase square-wave pulse

**Power**

- from the power line
- \* 95 V to 240 V, 49 Hz to 65 Hz
- \* from a 12-Volt power source of the emergency vehicle
- \* from exchangeable, rechargeable NiCd batteries
- \* rated voltage 12 V
- \* rated capacity 1.4 Ah
- \* battery is charged while inserted in the device
- \* charging time for depleted battery approx. 16 hours
- \* operating time with a fully charged battery approx. 35 defibrillation shocks of 360 joules each (into 50 ohms) or 2 hours of monitor operation (1.2 hours with pacemaker and SpO<sub>2</sub> measuring system)

**Operational readiness**

4 s after power up (incl. automatic selftest)

**Operating position**

any

## Environment

### Operation

under the following conditions regarded as normal:

- \* temperature between 0 and +40 °C
- \* rel. humidity between 30 and 95%, no condensation
- \* atmospheric pressure between 700 and 1060 hPa

### Storage and transport

- \* temperature between -20 and +60 °C
- \* rel. humidity between 10 and 95%, no condensation
- \* atmospheric pressure between 500 and 1060 hPa

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


## Installation Instructions for the CardioServ Wall Mount System 202 307 02

- \* The mounting system has 4 bore holes both at the top and bottom edges, allowing it to be mounted to the ambulance in a vertical or horizontal position (take care that the mounting system is flush against the wall of ambulance; in the vertical position the wider part with the cable aperture must face down!) (see Fig. 1).


Figure 1

- \* Replace the 4 rubber feet of **CardioServ** with the metal feet enclosed. When screwing down the metal feet, please note that they have a recess on one side for the screw head.
- \* Use the appropriate type of screws for attachment of the metal feet: If the **CardioServ** unit is equipped with metal bushings, please use the fillister head screws DIN 912 M 4 x 12; if there is no metal bushing but only a hole, please use Plastite screws 4 x 16.

- \* Please be sure to use the red feet (angular groove) at the bottom, i.e., on the side with the contacts. (See Fig. 2).



Figure 2

- \* If required, you may attach the enclosed rubber feet on top of the metal feet.
- \* Mount the **CardioServ** by first inserting the feet at the bottom, then by introducing the feet at the top.
- \* To remove the device press on the latches (top center).

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

A PART OF THE  
MARQUETTE  
HELLIGE  
UNITY  
NETWORK®

The Marquette Hellige Unity Network® is our response to the critical challenge of managing patient information. The Unity icon represents more than a network. It symbolizes our commitment to integrating information and devices throughout a healthcare system.

European Headquarters  
Marquette Hellige GmbH  
Postfach 728  
79007 Freiburg • Germany  
Tel. +49 761 45 43 - 0  
Fax +49 761 45 43 - 233

World Headquarters  
Marquette Medical Systems, Inc.  
8200 West Tower Avenue  
Milwaukee, WI 53223 • USA  
Tel. +1 414 355 5000  
Fax +1 414 355 3790

Asia Pacific  
Marquette Medical Systems  
26/F, Catic Plaza  
8 Causeway Road, Hong Kong  
Tel. +852 2804 2320  
Fax +852 2804 1776

**marquette**  
**HELLIGE**  
MEDICAL SYSTEMS

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L. D. LAMBERT ©