

CardioSmart/CardioSmart ST

Servicing Instructions

Version 1.2 and 1.3
227 435 11 SA(e) Revision D



marquette

A GE Medical Systems Company

Caution:

During repairs/service interventions, observe the protective measures against damage due to ESD.

- * Marquette Hellige GmbH considers itself responsible for the effects on safety, reliability, and performance of the equipment, only if:
 - assembly operations, extensions, readjustments, modifications, or repairs are carried out by Marquette Hellige GmbH or by persons authorized by Marquette Hellige GmbH,
 - the electrical installation of the relevant room complies with the applicable national and local requirements, and
 - the instrument is used in accordance with the instructions for use.
- * This manual contains service information; operating instructions are provided in the user manual of the instrument.
- * This manual is in conformity with the instrument at printing date.
- * All rights are reserved for instruments, circuits, techniques, and names appearing in the manual.

Contents

| | |
|--|-----------|
| 1. Device | 6 |
| 1.1 Block circuit diagram of entire instrument | 9 |
| 1.2 Mechanical Components | 10 |
| 2. Functional description | 11 |
| 2.1 PCB Power Supply_CS | 11 |
| 2.1.1 Switch on/off circuit | 11 |
| 2.1.2 Charging circuit for NC Battery | 11 |
| 2.1.3 Device behavoir depending on the state of battery charge | 12 |
| 2.2 PCB Control CS_S | 13 |
| 2.2.1 Generation of the logic power supply +5V | 13 |
| 2.2.2 ECG recording and front-end processing | 13 |
| 2.2.3 Controller core | 14 |
| 2.2.4 Real-time clock | 14 |
| 2.2.5 Memory | 14 |
| 2.2.6 Printhead control | 15 |
| 2.2.7 Motor controller | 15 |
| 2.2.8 RS-232 interface | 16 |
| 2.2.9 Keypad interface | 16 |
| 2.3 PCB Graphics Display_CS | 17 |
| 2.3.1 Power supply | 17 |
| 2.3.2 Background illumination | 17 |
| 2.3.3 Contrast control | 17 |
| 2.4 Internal interfaces | 18 |
| 2.4.1 Mechanical interfaces | 18 |
| 2.4.2 Electronic interfaces | 19 |
| 2.4.2.1 Interface Character Display | 19 |
| 2.4.2.2 Interface Graphics Display | 21 |
| 2.4.2.3 Interface PCB Power Supply_CS | 21 |
| 2.4.2.4 Option interfaces | 23 |
| 2.5 Interfaces to peripherals | 24 |
| 2.5.1 Electronic interfaces | 25 |
| 2.5.1.1 RS-232 interface | 26 |
| 2.5.1.2 Patient input | 27 |
| 2.6 Limitations | 27 |
| 3. System test functions | 30 |
| 3.1 System test functions with graphics display | 30 |
| 3.1.1. General information | 30 |
| 3.1.2 Test start | 30 |
| 3.1.3 Display test (1) | 31 |
| 3.1.4 Key test (2) | 32 |
| 3.1.5 Motor test (3) | 32 |
| 3.1.6 Test results (4) | 33 |
| 3.1.7 Self-test | 34 |
| 3.1.8 Recording test (5) | 34 |
| 3.1.9 Interface test (6) | 35 |
| 3.1.10 Time and date (7) | 36 |
| 3.1.11 Electrode test (8) | 37 |
| 3.1.12 Test time constants (9) | 38 |
| 3.1.13 Setting the Device model (V) | 38 |

| | |
|--|-----------|
| 3.1.14 Switching over the program for Interpretation (C) | 40 |
| 3.2 Device test functions with character display | 42 |
| 3.2.1 General information | 42 |
| 3.2.2 Test start | 42 |
| 3.2.3 Display test (1) | 42 |
| 3.2.4 Key test (2) | 43 |
| 3.2.5 Motor test (3) | 43 |
| 3.2.6 Test results (4) | 44 |
| 3.2.7 Self-test | 45 |
| 3.2.8 Recording test (5) | 45 |
| 3.2.9 Interface test (6) | 46 |
| 3.2.10 Time and date (7) | 47 |
| 3.2.11 Electrode test (8) | 48 |
| 3.2.12 Test time constants (9) | 49 |
| 3.2.13 Setting the Device model (V) | 49 |
| 3.2.14 Switching over the program for Interpretation (C) | 51 |
| 4. Repair instructions | 52 |
| 4.1. Safety instructions | 52 |
| 4.2 Replacing components | 52 |
| 5. Troubleshooting tips | 56 |
| 6. Adjustment instructions | 59 |
| 7. Servicing and maintenance | 61 |
| 7.1 Technical inspection | 61 |
| 7.1.1 Visual check | 61 |
| 7.1.2 Test functions | 62 |
| 7.1.2.1 Recommended testing instruments and accessories | 62 |
| 7.1.2.2 Test preparations | 62 |
| 7.1.2.3 Operating and display unit performance tests | 62 |
| 7.1.2.4 Test for recording speeds 25 and 50 mm/s | 63 |
| 7.1.2.5 Device test result check | 63 |
| 7.1.2.6 RS232 interface test | 63 |
| 7.1.2.7 Analysis of the ECG signals and HR value | 63 |
| 7.1.2.8 Pacemaker identification test | 64 |
| 7.1.2.9 Identification of disconnected electrodes | 65 |
| 7.1.2.10 Checking the charge status of the NC battery | 65 |
| 7.1.2.11 Pump leakage | 65 |
| 7.1.3 Safety Analysis Tests | 66 |
| 7.1.3.1 General information | 66 |
| 7.1.3.2 Protective Earth Resistance Test | 66 |
| 7.1.3.3 Measurement of Leakage Current | 66 |
| 7.1.3.3.1 Enclosure Leakage Current Test | 67 |
| 7.1.3.3.2 Patient Leakage Current Test | 68 |
| 7.2 Maintenance, cleaning, disinfection | 69 |
| 8. Parts Lists | 70 |
| 9. Specifications | 73 |
| 9.1 Specifications CardioSmart | 73 |
| 9.2 Specifications CardioSmart ST | 79 |

| | |
|-----------------------------------|-----------|
| 10. Device documents | 85 |
| Annexe | 86 |

Revision History

| | | |
|-----------|---------------|-----------------|
| V 1.0 | May 1995 | Initial Release |
| V 1.1 | March 1996 | Update 1 |
| V 1.2 | February 1997 | Update 2 |
| V 1.2/1.3 | August 1997 | Update 3 |

1. Device description

This Service Manual for device versions V1.2 and V1.3 describes both CardioSmart and CardioSmart ST. Unless indicated specifically, the description applies to both CardioSmart and CardioSmart ST.

The device version V1.3 hardware components are identical to those of device version V1.2. The new combinations of these components in V1.3 result in new models, so-called regional models. Differences between device versions V1.2 and V1.3 are indicated clearly in the relevant sections.

The CardioSmart is a portable electrocardiograph with an integrated printing unit. It is used to acquire, record and process ECG signals. It is designed for line-power and battery operation. Operation without a battery is not possible. Power supply unit and battery are integrated into the instrument.

CardioSmart ST has an additional ergometry mode option with its own ergometry keypad.

The patient cable for the acquisition of ECG signals is connected by means of a 15-pin connector as used in earlier ECG recorders. This means that the patient cables used in the past can still be used.

Stored patient data and ECGs can be transmitted to a PC via an RS-232 interface. This requires a CardioProm version with additional patient data memory.

CardioSmart and CardioSmart ST are both based on the same hardware platform V1.2 or V1.3 and represent different types of the basic unit 101 116 xx.

With device version V1.2, the following models are available:

| | | | |
|------------|---------------------------|-------|----------------------|
| 101 116 01 | CardioSmart international | 230V~ | |
| 101 116 02 | CardioSmart US Version | 115V~ | |
| 101 116 03 | CardioSmart Russian | 230V~ | |
| 101 116 04 | CardioSmart international | 230V~ | with integrated pump |
| 101 116 05 | CardioSmart Russian | 230V~ | with integrated pump |

| | | | |
|------------|------------------------------|-------|----------------------|
| 101 116 10 | CardioSmart ST international | 230V~ | |
| 101 116 11 | CardioSmart ST international | 230V~ | with integrated pump |
| 101 116 12 | CardioSmart ST US Version | 115V~ | |
| 101 116 13 | CardioSmart ST US Version | 115V~ | with integrated pump |

With device version V1.3, the following models are additionally available:

| | | | |
|------------|----------------------|-------|----------------------|
| 101 116 06 | CardioSmart EUROPE 2 | 230V~ | |
| 101 116 07 | CardioSmart EUROPE 2 | 230V~ | with integrated pump |

| | | | |
|------------|-------------------------|------------|----------------------|
| 101 116 14 | CardioSmart ST Russian | 230V~ | |
| 101 116 15 | CardioSmart ST Russian | 230V~ | with integrated pump |
| 101 116 16 | CardioSmart ST EUROPE 2 | 230V~ | |
| 101 116 17 | CardioSmart ST EUROPE 2 | 230V~ | with integrated pump |
| 101 116 18 | CardioSmart ST ASIA | 115V/230V~ | |
| 101 116 19 | CardioSmart ST ASIA | 115V/230V~ | with integrated pump |

The CardioSmart entire instrument comprises the following components:

- **CardioSmart Basic Instrument**, without display and CardioProm
- **CardioProm Storage Module**, there are several modules with different features, with or without ECG storage facility
- **CardioVision Display**, there are 2 different displays available:
 - CardioVision Text, LCD module 2-line
 - CardioVision Graphics, graphics LCD

The CardioSmart ST entire instrument comprises the following components:

- **CardioSmart ST Basic Instrument**, without display and CardioProm
 - **CardioProm ST Storage Module**, there are several modules with different features, with or without ECG storage facility
 - **CardioVision Graphics**, graphics LCD
- CardioSmart ST is not enabled for CardioVision Text

This concept permits the user to configure an instrument himself to meet his own application specifications. It is also possible to upgrade the basic version.

The hardware comprises the following functional blocks:


- PCB Control CS_S (basic instrument)
- Power supply unit with battery (basic instrument)
- Keypad (basic instrument)
- Printing unit (basic instrument)
- PCB Energy Storage_CS (storage module)
- Display (character or graphics display)

The following functional blocks are constructed directly as PCBs:

- PCB Control_CS
- PCB Energy Storage_CS
- PCB Power Supply_CS with battery charging circuit
- PCB Graphics Display_CS

The applications, functions and operation of the CardioSmart are described in the User's Manual.

1.1 Block circuit diagram of entire instrument



1.2 Mechanical Components

The main mechanical components comprise the **upper and lower shell** of the **CardioSmart Basic Instrument**. The lower shell serves as a basic unit to receive the following assemblies:

- Mains input module with instrument On/Off switch
- PCB Power Supply_CS
- Battery
- Thermal printing unit
- Paper container
- PCB Control_CS

The upper shell accommodates the keypad, which is connected to the PCB Control_CS via a 22-pin connector. The display (character or graphics display) is also integrated into the upper shell.

The **CardioProm** (program module) also has its own housing, comprising an upper and lower assembly section. It contains the PCB Energy Storage_CS, which is connected to the PCB Control_CS via a multipoint connector when inserted into the CardioProm slot.

The two available **displays** are also located in a housing, comprising an upper and lower shell. The graphics display contains the graphics LC display itself, as well as the associated control electronics, it can be swivelled within the range from 0° to 45° with 5 lock-in positions. The character display only contains the character LCD and cannot be swivelled. The two displays can be connected to the PCB Control_CS by insertion into the depression in the upper shell of the basic instrument provided with a female multipoint connector.

The 15-pin input connector to **connect the patient cable** and the 9-pin connector for the **RS-232 interface** are located directly on the PCB Control_CS.

Models 101 116 04/05 have an added integrated suction pump, comprising the pump module 303 442 91 and the PCB Pump 388 029 18.

2. Functional description

The **block circuit diagram** of the entire instrument in **Section 1.1** and the **functional blocks of the P-plans** describe the individual functional blocks.

2.1 PCB Power Supply_CS

The PCB Power Supply provides the following functions:

- Secondary clock-rated power supply unit to generate the fundamental device voltage UVERS
- Device switch on/off circuit
- Charging circuit for NC battery
- Battery almost empty identification
- Battery empty instrument switch-off

2.1.1 Switch on/off circuit

The device on/off switch (momentary toggle switch) supplies the starting pulse. The actual stop signal comes from the PCB Control_CS via the DOFF cable.

In the line-power operating mode CardioSmart can be switched off by:

- Device on/off switch
- Processor (DOFF)

In the battery operating mode CardioSmart can be switched off by:

- Device on/off switch
- Processor (DOFF)
- Battery empty identification

In the line-power operating mode recording can be switched off by:

- Battery empty identification

2.1.2 Charging circuit for NC Battery

A charging IC is used to charge the battery. During charging it monitors the battery voltage function and switches over from rapid charging to trickle charging in conformity with the charging current.

The charging circuit is activated by applying the power supply voltage or battery. If a recording is made whilst the device is in the operating mode the charging circuit is deactivated.

2.1.3 Device behaviour depending on the state of battery charge

As a general rule the following applies:

No applied line-power voltage:

Operation of the device without a battery is not possible

Applied line-power voltage, battery is charged:

Device completely operational

Applied line-power voltage, battery **depleted**:

Device is only partially operational (everything except recording)

Battery charging is triggered by:

Application of line-power voltage

Switching off a current recording, manually with the Stop key
or automatically when the cutoff voltage has been attained

Device behaviour:

- State 1 Battery depleted, power supply system is connected up:
Battery charged in accelerated mode; when full, switchover to
conservation charging, no time limit
- State 2 Battery full, power supply system is connected up:
Battery charged in accelerated mode; when full, switchover to
conservation charging, no time limit
- State 3 Battery depleted, no line-power supply, device is switched on:
Device inoperative
- State 4 Battery full, no line-power supply, device is switched on:
Device fully operative, when battery almost empty message appears
LEDBAT signal activated, when battery empty message appears the
device becomes inoperative
- State 5 Battery depleted, power system supply is connected, charging begins,
device is switched on:
All functions except recording, as long as the battery cannot supply the
printhead with sufficient power.
- State 6 Battery full, power supply system is connected, device switched on:
Device fully operative. During recording charging is interrupted, at stop
recharging begins until the charging IC switches to conservation
charging.

2.2 PCB Control CS_S

2.2.1 Generation of the logic power supply +5V

The logic power supply +5V is generated directly on the PCB Control_CS.
The stabilized voltage UVERS, which is generated by the PCB Power Supply_CS, serves as the input voltage.

The +5V power supply is generated by the clock-rated **voltage controller**.
The "**adjustable version**" of the controller is used in order to be able to use it later for applications requiring other voltages, the output voltage is determined by appropriate dimensioning of the voltage divider at the feedback input of the voltage controller.

2.2.2 ECG recording and front-end processing

The main element of this operational unit is a **set of chips comprising 4 ASICs** located on the PCB.

The **patient lead** is connected with a **15-pin connector**. Overvoltage arresters at the electrode inputs provide high-voltage protection. A hybrid chip acts as a microfuse for each electrode, it serves to identify a disconnected electrode.

The patient input is classified as cardiac floating and is defibrillation-proof.

Each electrode input is preamplified after the hybrid chip with a low-noise operational amplifier with a small input voltage, before it is fed into the two ASICs with the sigma-delta resistors.

The filter time constants can be adjusted by the controller core by a command. 8 time constants are available for selection (system test).

N common-mode compensation ensures suppression of interference, at the same time serving to improve the in-phase suppression of the input electrodes.

To protect patients the ECG recorder and conditioning processor are assembled as floating components. Digital signals from and to the controller core are transmitted via optoelectronic couplers, the floating power supply +/-5V is generated by a flyback converter from the logic power supply +5V.

2.2.3 Controller core

The actual core comprises the **Motorola Controller 68332**, containing the following integrated components:

- CPU32, computer core, internal 32-bit register, external 16-bit processor
- TPU, an independently operative clock generator processor
- QSM with SCI to implement a simple RS-232 interface and a serial QSPI interface with up to 16 channels
- SIM with chipselect generation, system monitoring, clock synthesizer

The **watchdog/reset generation** is implemented with an integrated system monitoring chip. It has the following functions:

- Power-up reset for the 68332 when the instrument is switched on
- Voltage monitoring with reset generation
- Watchdog function
- Switchover to battery for buffered RAMs
- Access protection for buffered RAMs

2.2.4 Real-time clock

This provides the time of day and date. When the device is operative it is powered by the logic power supply, when the device is switched off a 3-V lithium cell is switched in to keep the clock running.

The control signals for the clock (chipselect read/write signal) are generated in the **PAL RTC_Control**, which also contains the release sequence for the character display.

2.2.5 Memory

Static RAMs are used as **RAM memory** with an internal organization 128kbit x 8. They are equipped with 4 chips, providing an organisation of 256kbit x 16. It is possible to add 4 further chips to give a RAM memory area of 1MByte.

The access time is 70 ns, i.e., access **without wait states**.

2.2.6 Printhead control

The printhead controller takes on the complete control of a 216-mm thermal printhead with a line width of 200 mm. It also contains the printhead monitor, which reduces the heating time with increasing printhead temperature, terminating the printhead power supply when the printhead temperature reaches 55°C. The printhead power supply is reactivated only when the printhead temperature drops below 50°C.

The core is a **gate array**, which takes on the heating time management for each dot. The speed-dependent heating parameters are transmitted serially.

For the **printhead monitoring** the temperature of the thermal array printhead is measured by a thermistor located on the printhead. A constant current source effects a temperature-dependent voltage drop at a comparator input, the switching threshold being adjusted by hysteresis at the other comparator input.

If a printhead temperature of 55°C is exceeded, the message PRINTHEADTEMP is activated and the current supply to the printhead prevented by the STROBE signal. Only after the temperature drops below 50°C is the PRINTHEADTEMP message disabled and the current supply via STROBE reenabled.

A second comparator generates a VTEMP signal when the temperature reaches 45°C, it is disabled again when the temperature drops below 40°C. It represents a prewarning signal, so that the quantity of data at the printhead is reduced to prevent it being switched off.

In addition, a continuous reduction in the current activation time occurs with increasing temperature, resulting in a regular type face throughout the entire temperature range.

The **printhead is operated directly at the battery**, i.e., the printhead voltage is dependent on the state of battery charge. Compensation is achieved by reducing or extending the printhead current activation time accordingly, ensuring that the type face is maintained within the fluctuation range of the battery voltage.

2.2.7 Motor controller

This activates the CardioSmart's DC motor. This motor does not have a speedo, control is achieved according to the principle of current compensation.

The main components are the speed selector, motor controller, and motor current monitoring.

The resistor for the **control** according to the principle of current compensation is attached directly to the motor casing.

The voltage of +/- 12 V required to send control signals is generated from the supply voltage as follows: An in-phase regulator generates + 12 V, due to the small current requirements the voltage -12 V is generated by an inverting voltage pump. The supply voltage +/-12 V is also used for the temperature monitoring system of the printhead controller.

It is only enabled during printing!

When the motor current is raised > 140 mA (e.g., as a result of paper jam, transmission problems) the **motor current monitoring system** terminates paper transport and the message "Please insert paper and press return key" appears on the display.

2.2.8 RS-232 interface

The CardioSmart has an RS-232 interface which, except for the RS-232 driver chip, is directly integrated into the controller.

It has the following attributes:

- Software handshake with XON/XOFF
- No HW handshake
- Transmission speed 1200 19200 Baud

To achieve **protection against ESD up to 10kV**, the interface is equipped with additional suppressor diodes.

2.2.9 Keypad interface

The interface to the keypad is via the **22-pin connector AS**. The keypad itself is a membrane keypad with touch-sensitive check-back.

The keypad comprises a **matrix with 7 columns and 7 rows**.

The CardioSmart ST keypad has additional keys for ergometer/treadmill control. Depending on the keypad used (CardioSmart or CardioSmart ST) the PC receives a corresponding signal, identifying the device used.

Identification of the key pressed is as follows:

The controller activates a column, activation is via low-level, simultaneously it reads out the register for the lines and by a low-level signal identifies whether one or more keys had been pressed. Then the next column is activated. The run procedures are repeated automatically until all columns have been activated.

In addition, the LEDs for start and stop, for the power supply system display and the battery status control are located directly on the keypad.

2.3 PCB Graphics Display_CS

This PCB is only used in connection with the component parts of the CardioSmart graphics display unit. It comprises the following functional units:

- Voltage generation
- CPU interface, display control
- Background illumination
- Contrast control

2.3.1 Power supply

A **negative voltage of -23 V** is required for the power supply of the display. It is generated by an inverting switching controller from the 5-V power supply. The output voltage -23 V is enabled by the additional external wiring of the switching controller. The voltage -23 V also serves as a voltage source for the generation of the **negative contrast voltage**.

The **functional unit LCD BLOCK** ensures that the logic power supply +5 V is active at the graphics display before the negative power supply and the contrast voltage. When the instrument is switched off the negative power supply and the contrast voltage are switched off first.

2.3.2 Background illumination

The background illumination of the graphics display is effected by a **CFL tube**, which requires an alternating voltage of approx. 800 Vss.

The alternating voltage is generated by a voltage inverter, driven by 5 V. Since the CardioSmart is a portable, battery-operated device, the background illumination can be **enabled and disabled** (configuration).

2.3.3 Contrast control

The contrast control has an additional **temperature compensation control**. However, when there are rapid fluctuations in temperature it takes a certain time until the compensation is complete, since the display itself has a greater thermal time constant.

2.4 Internal interfaces

2.4.1 Mechanical interfaces

The 3 following mechanical interfaces are important within the CardioSmart:

- Interface to CardioProm
- Interface to display CardioVision Text
- Interface to display CardioVision Graphics

Interface to CardioProm

The CardioProm has its own housing, which is conducting to prevent damage by ESD. To shunt any possible static charge it has 3 rails on the underside which, on the PCB Control_CS in the CardioSmart basic device, when inserted make contact with the 3 rails on GND, thereby shunting the charge.

This mechanical insertion of the rails ensures that a fail-safe contact is made with the PCB Control_CS.

Furthermore, mechanical coding ensures that the non-compatible CardioProms of the EK56/EK512 series cannot be inserted into the slot. Similarly, the insertion of CardioSmart CardioProms into the EK56/512 slot is not possible.

Interface to display CardioVision Text

The display CardioVision Text is placed into the depression provided in the upper shell of the CardioSmart basic device. The mechanical rail on the outside of the display ensures the electrical contact is automatic and fail-safe. The display CardioVision Text can be replaced, e.g., by a CardioVision Graphics using the appropriate tools without opening the device.

Interface to display CardioVision Graphics

The display CardioVision Graphics is also placed in the above-mentioned depression, the mechanical rail on the outside also ensures a fail-safe contact. In contrast to the CardioVision Text, the CardioVision Graphics can only be replaced by opening the upper shell of the basic device.

2.4.2 Electronic interfaces

This section describes the pinning, function and significance of the signals of the internal interfaces of the functional components.

2.4.2.1 Interface Character Display

The interface to the character display provides all the signals to connect the display CardioVision Text. It comprises the upper 8 bits of the data bus, control signals and the signal for the contrast control.

Connection denotation: **AG/**

Female multipoint connector **2 x 7-pin, upright** (180 °), reverse terminal protection achieved mechanically.

The function of the individual pins is given in the following table. The definition as an Input/Output is seen with reference to PCB Control_CS.

| Pin Number | Input/Output | Denotation | Function |
|------------|------------------------|-------------------------|-------------------------------|
| 1 | Output | GND | System ground |
| 2 | Output | + 5 V | Voltage supply |
| 3 | Output | CHAR_CONTR 0 5V | Contrast character display |
| 4 | Output | AD_12 | Controller address bus bit 12 |
| 5 | Output | RW_ | Read/write signal |
| 6 | Output | ENCHARD | Enable character display |
| 7 | Bi-Direct, tristate | DB_8 | Controller data bus bit 8 |
| 8 | Bi-Direct, tristate | DB_9 | Controller data bus bit 9 |
| 9 | Bi-Direct, tristate | DB_10 | Controller data bus bit 10 |

| | | | |
|----|------------------------|-------|----------------------------|
| 10 | Bi-Direct, tristate | DB_11 | Controller data bus bit 11 |
| 11 | Bi-Direct, tristate | DB_12 | Controller data bus bit 12 |
| 12 | Bi-Direct, tristate | DB_13 | Controller data bus bit 13 |
| 13 | Bi-Direct, tristate | DB_14 | Controller data bus bit 14 |
| 14 | Bi-Direct, tristate | DB_15 | Controller data bus bit 15 |

2.4.2.2 Interface Graphics Display

The interface to the graphics display provides all the signals to connect the display CardioVision Graphics. It comprises the upper 8 bits of the data bus, address bus, control signals and the signal for the contrast control.

Connector denotation: **AH/**

2.4.2.3 Interface PCB Power Supply_CS

This has the voltages from the power supply unit to generate the + 5 V and the other voltages, the control of the LEDs for the power lamp and battery monitoring, the control signals for battery charging and switching on the printhead and motor power supply as well as the signal to switch off the device on the controller's part.

Connector denotation: **AE/**

Male multipoint connector **2 x 10-pin, upright** (180 °), reverse terminal protection and coding with coding pin 12.

The function of the individual pins is given in the following table. The definition as an Input/Output is seen with reference to PCB Control_CS.

| Pin Number | Input/Output | Denotation | Function |
|----------------|--------------|------------|---|
| 1 | Input | LEDLIN | Control of LED for power supply system "- " = mains not connected "AC" = mains connected |
| 2 | Output | REG_ON | Printhead power supply, motor enabled by controller: "0" = Enable battery charging "1" = Enable printhead power supply |
| 3 | Output | GOFF_ | Switch on device with controller "0" = Device off "1" = Device on |
| 4, 5 | Input | UVERS | Input voltage to generate +5 V |
| 6 | Output | + 5 V | Logic power supply |
| 7,8,9, 10 | Input | COMMON_IN | Input voltage to generate motor- printhead voltage and +/- 12V |
| 11 | Input | REG_OFF | Printhead power supply, motor enabled by power-supply unit: "0" = Enable battery charging "1" = Enable printhead power supply |
| 12 | | CODE | Coding pin |
| 13, 14, 15, 16 | Input | GND | System ground |
| 17 | Input | VRAM | Buffered voltage to generate the RAM standby power supply |
| 18 | Output | COMMON_S | Voltage switched in from COMMON_IN |
| 19 | | NC | NC |
| 20 | Input | LEDBAT | Charging status display for BATTERY: "0" = battery needs charging "1" = battery charged, normal operat. |

2.4.2.4 Option interfaces

Option Power Check

This constitutes a service connector to enable checking of all the voltages active within the system. This **is not an extension port**, the voltages must not be used to supply power to other components!

Connector denotation: **AF/**

Male multipoint connector **1 x 6-pin, upright (180 °), for servicing purposes only!**

The function of the individual pins is given in the following table. The definition as an Input/Output is seen with reference to PCB Control_CS.

| Pin Number | Input/Output | Denotation | Function |
|------------|--------------|------------|--|
| 1 | Output | + 12 V | Testing point power supply + 12 V |
| 2 | Output | GND | System ground |
| 3 | Output | - 12 V | Testing point power supply - 12 V |
| 4 | Output | COMMON_S | Testing point power supply motor voltage |
| 5 | Output | GND | System ground |
| 6 | Output | + 5 V | Testing point logic power supply + 5 V |

2.5 Interfaces to peripherals

The CardioSmart has only 3 interfaces for peripherals:

- Mains input
- Patient Input
- RS-232 interface

The **mains input** interface on the device is a 3-pin standard cold appliance socket connection, which is integrated into the mains input module. Connection to the mains is effected via a 3-pin power cord with a non-fused earth conductor.

Although the mains input is not designed as a "wide range input", adjustment to the two mains voltage ranges 110V~ ... 120V~ and 220V~ ... 240V~ is easily achieved without having to open the device or a wiring change of the bridges:

- Check voltage setting in the inspection window of the interrupter of the mains input module
- If the voltage value is not the same as the mains voltage:
- **Disconnect mains plug!**
- Open interrupter with appropriate tool
- Take out fuse holder with fuses
- Insert fuses rated according to the rating plate specifications
- Reinsert holder so that the required voltage is visible
- Close interrupter
- The required mains voltage should now appear in the inspection window

This readjustment is to be carried out by the manufacturer or servicing agent only!

No further description of the mechanical and electronic parameters of the mains input is given in Section 2.5.

The mechanics and the interface configuration as used in the EK53/EK56 and in the EK512 are implemented in the **patient input**. This means that all the patient cables used up to date can still be used, this also being true for the suction system with its integrated pump.

A 9-pin sub-D connector with a standard configuration of the signals TXD, RXD and GND is implemented in the construction of the **RS-232 interface**. In addition, it is attached to the connected plug housing by means of threaded inserts. The EK53/56 and EK 512 RS-232 power cables can thus no longer be used.

The cable 223 362 03 is used as an RS-232 cable to the PC.

2.5.1 Electronic interfaces

2.5.1.1 RS-232 interface

This comprises the transmitter signal TXD, the incoming signal RXD and the reference ground GND. The pin configuration is the same as that used for a standard PC 9-pin sub-D socket. The interface does not operate with a hardware handshake.

Furthermore, a REMOTE START input is provided at pin 8.


Transmission is in accordance with the standard V-24 protocol, input and output voltages are defined as follows:

- maximum input voltage at RXD: $\pm 15\text{V}$
- minimum output voltage at TXD: $\pm 5\text{V}$
- additional input protection against ESD effects: $\pm 10\text{kV}$

In addition, there is a REMOTE START input at pin 8 with the following specifications:

- external short-circuiting contact via circuit reference referred to ground
- max. source resistor switch: $R_i < 300\Omega$
- make time $> 100\text{ ms}$
- continuous load limits $\pm 10\text{ V}$
- ESD protection to max. $\pm 10\text{ kV}$

Connection principle to PC



The connector is denoted: **AC/**

Male multipoint connector **sub-D socket housing, 9-pin, 90° offset**

The function of the individual pins is given in the following table. The definition as an Input/Output is seen with reference to PCB Control_CS.

| Pin Number | Input/Output | Denotation | Function |
|------------|--------------|--------------|----------------------|
| 1 | | NC | NC |
| 2 | Input | RXDE | Receive RS-232 data |
| 3 | Output | TXDE | Transmit RS-232 data |
| 4 | | NC | NC |
| 5 | Output | GND | System ground |
| 6,7,9 | | NC | NC |
| 8 | Input | REMOTE START | Remote Start Input |

2.5.1.2 Patient input

Below is a list of the most important features of the patient input:

- Classification: cardiac floating
- Defibrillation-proof
- Simultaneous data acquisition from all leads
- Individual electrode monitoring
- Voltage feed for floating power supply via floating transformer
- Control and data signals from and to floating part via optoelectronic coupler

The connector is denoted: **AA/**

Male multipoint connector **sub-D socket housing, 15-pin**

The function of the individual pins is given in the following table. The definition as an Input/Output is seen with reference to PCB Control_CS.

| Pin Number | Input/Output | Denotation | Function |
|------------|--------------|------------|---|
| 1 | Input | C2 | Input electrode C2 |
| 2 | Input | C3 | Input electrode C3 |
| 3 | Input | C4 | Input electrode C4 (Nap for Nehb cables) |
| 4 | Input | C5 | Input electrode C5 |
| 5 | Input | C6 | Input electrode C6 |
| 6 | Output | S | Shielding |
| 7 | Output | PLGND | Identification of patient cable used |
| 8 | Input | PL | Identification of patient cable used |
| 9 | Input | R | Input electrode R |
| 10 | Input | L | Input electrode L |
| 11 | Input | F | Input electrode F |
| 12 | Input | C1 | Input electrode C1 |
| 13 | Input | Nst | Input electrode Nst (for Nehb) |
| 14 | Output | N | Output electrode N (common mode compensation) |
| 15 | Input | Nax | Input electrode Nax (for Nehb) |

2.6 Limitations

The following operating modes are not implemented in the CardioSmart:

- No ergometry
- No SAO2
- No spirometry
- No late potentials, no RR variability
- No phono, no US Doppler

A scope output is not available.

No analog inputs.

The existing RS-232 interface does not provide an ECG trigger output.

Recording without or with an empty NC battery is not possible.

Lead accumulators cannot be used.

Primary cells cannot be used.

3. System test functions

With country default setting German all messages displayed during the self-test are in German, in all other country settings they appear in English.

3.1 System test functions with graphics display

3.1.1. General information

The functions for the system test are mostly menu-guided.

For complete execution some tests require special auxiliary resources. These include interface testers or PC, connection cables, signal generators, etc. The tests that are necessary are described in the various test descriptions.

3.1.2 Test start

After simultaneously pressing the keys

Shift + Auto (CardioSmart) or **Shift + Mode** (CardioSmart ST)

the initial display menu appears to enable selection of a specific test. (see below)
Keys used for menu selection appear in parenthesis.

System test

Test functions:

- (1) Display test
 - (2) Key test
 - (3) Motor test
 - (4) Test results
 - (5) Recording test
 - (6) Interface test
 - (7) Time and date
 - (8) Electrode test
 - (9) Time constant test
 - (V) Device model (from SW V4.21)
 - (C) 12SL / HEART (from V4.21)
- Terminate with SH_MAN or SH_MODE key

3.1.3 Display test (1)

After pressing the "1" key the following menu appears:

Display test

- (1) Test pattern
- (2) Character set
- (3) Delete display
- (4) Inverse
- (5) Display illumination
- Contrast control
 - Shift_Cursor_Down
 - Shift_Cursor_Up
- Terminate with any key

Pressing any other key than one of those listed in the menu above leads back to the initial display menu.

(1) Test pattern

Pressing this key generates a chessboard pattern. Each subsequent pressing generates the inverse display of its predecessor.

(2) Character set

The character set used for the display is output.

(3) Delete display

All display pixels are inactive.

(4) Inverse

The contents of the existing display appear in reverse video.

(5) Display illumination

The illumination is switched on or off by pressing this key.

Contrast control

The contrast setting can always be adjusted as described above (cursor keys).

3.1.4 Key test (2)

By pressing the "2" key the following menu appears:

Key T e s t

- (P) Bleeper test
- (K) Cover for paper compartment
- (E) End

When a key is pressed this key or its function is displayed. This can either occur by a triple character display, e.g., "AAA" or as a text, e.g., "CURS_UP". Moreover, pressing the "P" key tests the bleeper (audible sound) and pressing the "K" key tests a micro key to identify whether the cover for the paper compartment is open or closed. The "E" key terminates this test and simultaneously undergoes a self-test.

3.1.5 Motor test (3)

When the "3" key is pressed the following menu appears:

Motor/Speed Test

- Speed select key
- Start/Stop
- (A) Measurement marking pulse ON
- (S) Measurement marking pulse OFF

The speed select key is used to set the required speed and the motor set into motion with Start/Stop.

Whilst running the motor speed can be controlled by using the "A" key to activate the setting of the measurement marking pulse and to deactivate it with the "S" key. A making pulse is set once every second. The running speed can then be calculated or its accuracy checked from the distance between marking pulses. Any other key pressed not included in this menu leads to the initial display menu.

3.1.6 Test results (4)

Pressing the "4" key triggers the output of the test results.

The output of the test protocol supplies data on the software integrated in the device (Part No., Version No., production date of the firmware), printhead voltage and the test results on the memory test performed during power-up.

The areas listed under RAM have the following meaning:

CardioSmart

| | |
|-----------|---|
| 0 | - 080000h |
| 512 kByte | RAM area on the PCB Control CS_S |
| 100000 | - 140000h |
| 256 kByte | RAM area on the CardioProm (depending on CardioProm type) |
| 140000 | - 1C0000h |
| 512 kByte | RAM area buffered on the CardioProm (dep. on CardioProm type) |

CardioSmart ST

| | |
|-----------|---|
| 0 | - 080000h |
| 512 kByte | RAM area on the PCB Control CS_S |
| 100000 | - 180000h |
| 256 kByte | RAM area on the CardioProm (depending on CardioProm type) |
| 180000 | - 200000h |
| 512 kByte | RAM area buffered on the CardioProm (dep. on CardioProm type) |

Values may deviate from those given above, depending on the device or software version used.

A specification in the form 100000 - 100000 permits the conclusion that this RAM area does not exist in the device. However, if it does exist in conformity with the software version implemented, this means that this RAM cell at address 100000H is defective. The memory could also be defective when the initial and final addresses of the memory area differ and do not correspond to the above-mentioned limits. This then means, for example, that the specification 100000 - 123000 on the defective cell at address 123000H.

Following output the initial display menu reappears automatically.

3.1.7 Self-test

The test results elucidated in Section 3.1.6 are ascertained during the self-test, which is always performed on power-up. Should errors be detected, a message appears on the display after the test, indicating the possible errors. The following error codes are used to identify the error.

Error codes

The following error codes appear on the display together with the message "Self-test failed":

ERROR_CODES: 0 - reserved
 1 VEKT - error in vector table
 2 ERAM - RAM error on PCB Control
 3 LCDDR - LCDRAM error (Graphics Display)
 4 - reserved
 5 - reserved
 6 RAMB - RAM error on the CardioProm
 7 GRAM - error in buffered RAM on the CardioProm
 8 ROMB - ROM error (checksum) on the CardioProm

 9 - 10 not used

3.1.8 Recording test (5)

The recording test is activated by pressing the "5" key and interrupted with the "Start/Stop" key.

Immediately after pressing the "5" key the printer begins to record the first channels at 25mm/s and 20mm/mV.

In addition, the message "Recording test" appears on the display. In this test only the first two channels are passed on to the printer. No other channels can be configured.

3.1.9 Interface test (6)

Several possibilities are available to test the serial interface. Moreover, on the one hand, the signal transmission and receiving of the interface can be tested by an internal feedback from TXD and RXD and, on the other hand, the signal transmission and receiving with a remote station.

For the test with a remote station a standard PC cable is required for a serial interface (zero modem cable). In addition, the following transmission protocol should be adjusted at the remote station:

1 start bit, 8 data bits, parity even
1 stop bit
Baud rate 9600

Pressing the "6" key calls up the following menu:

Interface test

- (1) Transmitting and receiving
- (2) Transmit test string to device
- (3) Transmitting and receiving with device

(1) Transmitting and receiving

This test enables complete testing of the RS232 signal path for transmitting and receiving, including RS232 driver and connector.

This test requires an RS232 connector with an internal bridge from pin 2 (RXDE) to pin (TXDE). Depending on the result of the test, the following message appears:

Test result : OK
or
Test result : Error

This procedure requires a software version as from V4.1. Earlier software versions only permit internal testing of the RS232 interface without including the RS232 driver.

(2) Transmit test string to device

In this test a remote station, e.g., PC must be connected up and have a terminal program which can receive signals and be configured for the above-mentioned protocol. If the remote station is on receive, then every time the "2" key is pressed the test string "Interface test" is transmitted to the remote station. Simultaneously, the following message appears on the display:

Transmitting : interface test

(3) Transmitting and receiving with device

Transmitting and receiving can be tested with a remote device by pressing the "3" key. The device should be connected as described in (2).

After pressing the "3" key a test string is sent the remote station. The following message appears on the display:

Transmitting: interface test

Receiving:

A subsequent input at the remote station is sent back to the device and displayed there. Moreover, it should be noted that an input must take place within 10 s and the input terminated with "return". If the key is pressed a second time, then the string "CardioSmart Test" is transmitted. Altogether there are three test strings available, which are sent in the following sequence with repeated key activation:

"Schnittstellenest", "CardioSmart-Test", "Sieht doch gut aus, oder?"
("Interface Test", "CardioSmart Test", "Looks good, doesn't it?")

The last string corresponds to the expected test result.

This test is terminated with any other key not used for a test.

3.1.10 Time and date (7)

This section deals with the quick checking and setting of the time and date.

Pressing the "7" key calls up the following display:

```
Time    /    Date
10:12:03    10.11.1994
(S) Set
End with any key
```

From this screen the clock can be adjusted by pressing the "S" key, the following mask appearing:

```
Time      Date
[  ]      [  ]
```

The first mask is for the time and the second for the date.

Example: [1012] [10119] for 10:12 o'clock and 10.11.94

Pressing the return key terminates the entry. The adjusted data are subsequently accepted and seconds set to zero.

The test is terminated with any other key not used for a test.

3.1.11 Electrode test (8)

The electrode test is started by pressing the "8" key.

The electrode test is terminated by pressing any key.

To evaluate the electrode states the following status WORD is available.

Status = 0 -> Electrode OK
Status = 1 -> Electrode disconnected

| Bit15 | | | | | | | | | | Bit0 | | | | | | | | | | | | | | | | | | | | | | |
|-------|---|--|---|--|-----|--|-----|--|-----|------|-----|--|-----|--|----|--|----|--|----|--|----|--|----|--|----|--|---|--|---|--|---|--|
| | x | | x | | DAT | | ELN | | NAP | | NST | | NAX | | C6 | | C5 | | C4 | | C3 | | C2 | | C1 | | F | | L | | R | |

And precisely this status WORD extended to LONG is displayed on the screen during the electrode test to be performed. It is important to note, however, that this status WORD is not displayed as a binary number but as a hexadecimal number.

However, to make the status of a particular electrode quickly recognizable, defective electrodes are represented by an ID code.

EX_R, for example, stands for the defective electrode R and C1 for electrode C1.

This permits the identification of every individual electrode. Naturally, since several electrodes can become disconnected simultaneously, it was necessary to set a certain priority in the display of disconnected electrodes.

The following display priority was programmed:

EX_R, EL_N, EX_L, EX_F, C1, C2, C3, C4, C5, C6, NAX, NST, NAP

For example, if electrodes R and N have become disconnected, then first R is reported with EX_R as having become disconnected. Only after R is OK is EL_N displayed.

3.1.12 Test time constants (9)

The test is started by pressing the "9" key and interrupted with the "Start/Stop" key.

Immediately after the "9" key has been actuated the printer begins to record the first two channels at 25 mm/s and 20 mm/mV.

Moreover, the message "Testing time constants" appears on the display.

In this test only the first two channels send signals to the printer. No other channels can be configured. The default time constant is output on the printout as a control.

Setting the time constants

As a general rule, the time constants are adjusted to 2.04 s.

This setting can be changed with the CORR_CONFIG key.

The following settings are available:

Time constants:

- (1) 0.260s (2) 0.130s (3) 0.033s
- (4) 0.004s (5) 4.080s (6) 2.040s
- (7) 1.020s (8) 0.510s

The figures given in parentheses indicate the key used to select the corresponding time constant.

Verification of the currently active time constant is only possible by printing out the self-test results (see Section 3.1.6) or by testing the time constants.

3.1.13 Setting the Device models (V)

The setting of the Device model is activated by pressing the "V" key, the following display image appears:

D e v i c e m o d e l

- (1) INT
- (2) USA
- (3) EU_2
- (4) ASIA
- (5) RUSS

Terminate with any key.

The currently selected device model appears in reverse video.

Warning!

Changing the device model affects the selection of the output formats and languages available, can switch the program for Interpretation from HEART to 12 SL, or vice versa. Also affected are the default configuration settings, e.g., override function yes/no, line filter 50/60 Hz.

The device model configured during manufacture should only be changed when deemed really necessary.

The following table shows the most important combinations:

| Device model | Interpretation | Output formats | Languages | Default settings | |
|--------------|----------------|----------------|-----------|---------------------------|-------------|
| | | | | Override function enabled | Line filter |
| 01 INT | HEART | International | g,e,f,i,s | no | 50 Hz |
| 02 USA | 12SL | USA | e,f,s | yes | 60 Hz |
| 03 EU_2 | 12SL | Int. +H1 +H2 | e,g,f,sw | yes | 50 Hz |
| 04 ASIA | HEART | International | g,e,f | no | 50 Hz |
| 05 RUSS | HEART | International | g,e,f | no | 50 Hz |

Output formats: - International: 12_FS, 12_F1, 12_F2, 6_F1, 6_F2, 3_F1

- USA: STND1, STND3, 1x10, 2x5, 2x5_ex, 3RHY1

The required device model is selected by pressing the appropriate key, "1" to "5".
Quit selection menu by pressing any key.

Quit the self-test with "Shift" + "Man" or "Shift" + "Mode".

CardioSmart configures the appropriate items and initiates a cold start automatically.
Thus, when the device is rebooted automatically all the new settings are adopted.

Selecting a particular device model leads to the configuration of the program for Interpretation as indicated in the table above, even when the interpretation program was configured differently beforehand as described in section 3.1.14.

Note on saved ECGs when configuring the device model:

Saved ECGs are not lost.

If when a new device model is selected, the program for Interpretation is switched over, however, when printing out the saved ECGs the display in the status line (12SL, or without in the case of HEART) is related to the currently enabled program for

Interpretation, although the results are based on the program previously configured. In this case, the saved ECGs should be printed out or transferred to a PC before configuration of the new model.

3.1.14 Switching over the program for Interpretation(C)

Switching over the program for Interpretation is activated by pressing the “C” key, the following display image appears:

I n t e r p r e t a t i o n

- (1) HEART
- (2) 12SL

Terminate with any key.

The currently enabled interpretation program appears in reverse video.

Warning!

Switching over to another program for Interpretation affects the measurement results and the interpretation!

The program for Interpretation configured during manufacture should only be changed when deemed really necessary!

The required program for Interpretation is selected by pressing the appropriate key, “1” or “2”.

Quit the selection menu by pressing any key.

Quit the self-test with “Shift” + “Man” or “Shift” + “Mode”.

CardioSmart configures the appropriate items and initiates a cold start automatically.

Thus, when the device is rebooted automatically all the new settings are adopted.

Note on saved ECGs when switching over the program for Interpretation:

Saved ECGs are not lost.

After switching over to a different program for Interpretation when printing the saved ECGs the display in the status line (12SL, or without in the case of HEART) is related to the newly selected interpretation program, although the results are based on the program previously configured. In this case the saved ECGs should be printed out or transferred to a PC before switching over to the new program for Interpretation.

3.2 Device test functions with character display (only CardioSmart)

3.2.1 General information

In contrast to the graphics display, at the current time the device test functions all require menu-guided execution. This is because the rendition of the character display of max. 80 characters is very limited.

To be performed in their entirety some tests require special auxiliary aids, e.g., interface tester or PC, connection cables, signal generators, etc. Just what is required for which test are given in the various test descriptions.

3.2.2 Test start

After simultaneous pressing of the keys

Shift + Auto (CardioSmart) or **Shift + Mode** (CardioSmart ST)

the following initial display image appears:

S Y S T E M T E S T
END TEST: with Shift MAN (CardioSmart) or Shift + Mode
(CardioSmart ST)

Key used to select menu options are given in parentheses.

3.2.3 Display test (1)

After pressing the "1" key the following menu appears:

Display test
(3) Delete display (4) Inverse

Pressing any other key than those given above leads back to the initial display image.

(3) Delete display

All the display pixels are disabled.

(4) Inverse

Inverse here does not mean inverse display of the rendered text, but inverse to delete display, i.e., all the display pixels are enabled.

Contrast control

Furthermore, the contrast control can also be adjusted under this option "Display test". Moreover, there are 7 levels available which can be run though with the CURSOR_UP and CURSOR_DOWN keys.

3.2.4 Key test (2)

Pressing the "2" key calls up the following screen:

K e y T e s t
(E) End

When a key is pressed this or its meaning is displayed. This can either occur by a triple character display, e.g., "AAA" or as a text, e.g., "CURS_UP". Moreover, pressing the "P" key tests the bleeper (audible sound) and pressing the "K" key tests a micro key to identify whether the cover for the paper compartment is open or closed. The "E" key terminates this test and simultaneously undergoes a self-test. The initial display image appears.

3.2.5 Motor test (3)

When the "3" key is pressed the following menu appears:

Motor/Speed Test

The speed select key is used to set the required speed and the motor set into motion with Start/Stop.

Whilst running the motor speed can be controlled by using the "A" key to activate the setting of the measurement marking pulse and to deactivate it with the "S" key. A marking pulse is set once every second. The running speed can then be calculated or its accuracy checked from the distance between marking pulses.

Any other key pressed not included in this menu leads to the initial display menu.

3.2.6 Test results (4)

Pressing the "4" key triggers the output of the test results.

The output of the test protocol supplies data on the software integrated in the device (Part No., Version No., production date of the firmware), printhead voltage and the test results on the memory test performed during power-up.

The areas listed under RAM have the following meaning:

CardioSmart

| | |
|-----------|---|
| 0 | - 080000h |
| 512 kByte | RAM area on the PCB Control CS_S |
| 100000 | - 140000h |
| 256 kByte | RAM area on the CardioProm (depending on CardioProm type) |
| 140000 | - 1C0000h |
| 512 kByte | RAM area buffered on the CardioProm (dep. on CardioProm type) |

CardioSmart ST

| | |
|-----------|---|
| 0 | - 080000h |
| 512 kByte | RAM area on the PCB Control CS_S |
| 100000 | - 180000h |
| 256 kByte | RAM area on the CardioProm (depending on CardioProm type) |
| 180000 | - 200000h |
| 512 kByte | RAM area buffered on the CardioProm (dep. on CardioProm type) |

Values may deviate from those given above depending on the device or software version used.

A specification in the form 100000 - 100000 permits the conclusion that this RAM area does not exist in the device. However, if it does exist in conformity with the software version implemented, this means that this RAM cell at address 100000H is defective. The memory could also be defective when the initial and final addresses of the memory area differ and do not correspond to the above-mentioned limits. This then means, for example, that the specification 100000 - 123000 on the defective cell at address 123000H.

The output can be interrupted with the Start/Stop key. Following output the initial

display menu reappears automatically.

3.2.7 Self-test

The test results elucidated in Section 3.2.6 are ascertained during the self-test, which is always performed on power-up. Should errors be detected, a message appears on the display after the test indicating the possible errors. The following error codes are used to identify the error.

Error codes

The following error codes appear on the display together with the message "Self-test failed":

ERROR_CODES: 0 - reserved
 1 VEKT - error in vector table
 2 ERAM - RAM error on PCB Control
 3 LCDR - LCDRAM error (Graphics Display)
 4 - reserved
 5 - reserved
 6 RAMB - RAM error on the CardioProm
 7 GRAM - error in buffered RAM on the CardioProm
 8 ROMB - ROM error (checksum) on the CardioProm

 9 - 10 not used

3.2.8 Recording test (5)

The recording test is activated by pressing the "5" key and interrupted with the "Start/Stop" key.

Immediately after pressing the "5" key the printer begins to record the first channels at 25mm/s and 20mm/mV.

In addition, the message "Recording test" appears on the display. In this test only the first two channels are passed on to the printer. No other channels can be configured.

3.2.9 Interface test (6)

Several possibilities are available to test the serial interface. Moreover, on the one hand, the signal transmission and receiving of the interface can be tested by an internal feedback from TXD and RXD and, on the other hand, the signal transmission and receiving with a remote station.

For the test with a remote station a standard PC cable is required for a serial interface (zero modem cable). In addition, the following transmission protocol should be adjusted at the remote station:

- 1 start bit, 8 data bits, parity even
- 1 stop bit
- Baud rate 9600

Pressing the "6" key calls up the following menu:

Interface test Keys (1), (2), (3)
'1' key tests transmission and receiving

(1) Transmitting and receiving

This test enables complete testing of the RS232 signal path for transmitting and receiving, including RS232 driver and connector.

This test requires an RS232 connector with an internal bridge from pin 2 (RXDE) to pin (TXDE). Depending on the result of the test, the following message appears:

Test result : OK
or
Test result : Error

This procedure requires a software version as from V4.1. Earlier software versions only permit internal testing of the RS232 interface without including the RS232 driver.

(2) Transmit test string to device

In this test a remote station, e.g., PC must be connected up and have a terminal program which can receive signals and be configured for the above-mentioned protocol. If the remote station is on receive, then every time the "2" key is actuated the test string "Interface test" is transmitted to the remote station. Simultaneously, the following message appears on the display:

Transmitting : interface test

(3) Transmitting and receiving with device

Transmitting and receiving can be tested with a remote device by pressing the "3" key. The device should be connected as described in (2).

After pressing the "3" key a test string is sent the remote station. The following message appears on the display:

Interface test Keys (1), (2), (3)

A subsequent input at the remote station is sent back to the device and displayed there. A received text would then be displayed on the second line of the screen.

Moreover, it should be noted that an input must take place within 10 s and the input terminated with "return". If the key is actuated a second time then the string "CardioSmart Test" is transmitted. Altogether there are three test strings available, which are sent in the following sequence with repeated key actuation:

"Interface Test", "CardioSmart Test", "Looks good, doesn't it?"

The last string corresponds to the expected test result.

This test is terminated with any other key not used to start a test.

3.2.10 Time and date (7)

This section deals with the quick checking and setting of the time and date. Pressing the "7" key calls up the following display:

| Time | / | Date |
|----------|---|------------|
| 10:12:03 | | 10.11.1994 |

From this screen the clock can be adjusted by pressing the "S" key the following mask appearing:

| Time | Date |
|-------|-------|
| [] | [] |

The first mask is for the time and the second for the date.

Example: [1012] [10119] for 10:12 o'clock and 10.11.94

Pressing the return key terminates the entry. The adjusted data are subsequently accepted and seconds set to zero.

The test is terminated with any other key not used for a test.

3.2.11 Electrode test (8)

The electrode test is started by pressing the "8" key.

The electrode test is terminated by pressing any key.

To evaluate the electrode states the following status WORD is available.

Status = 0 -> Electrode OK
Status = 1 -> Electrode disconnected

| Bit15 | | | | | | | | | | | Bit0 | | | | | | | | | | | | | | | | | | | | | |
|-------|---|--|---|--|-----|--|-----|--|-----|--|------|--|-----|--|----|--|----|--|----|--|----|--|----|--|----|--|---|--|---|--|---|--|
| | x | | x | | DAT | | ELN | | NAP | | NST | | NAX | | C6 | | C5 | | C4 | | C3 | | C2 | | C1 | | F | | L | | R | |

And precisely this status WORD extended to LONG is displayed on the screen during the electrode test to be performed. It is important to note, however, that this status WORD is not displayed as a binary number but as a hexadecimal number.

However, to make the status of a particular electrode quickly recognizable, defective electrodes are represented by an ID code.

EX_R, for example, stands for the defective electrode R and C1 for electrode C1.

This permits the identification of every individual electrode. Naturally, since several electrodes can become disconnected simultaneously, it was necessary to set a certain priority in the display of disconnected electrodes.

The following display priority was programmed:

EX_R, EL_N, EX_L, EX_F, C1, C2, C3, C4, C5, C6, NAX, NST, NAP

For example, if electrodes R and N have become disconnected, then first R is reported with EX_R as having become disconnected. Only after R is OK is EL_N displayed.

3.2.12 Test time constants (9)

The test is started by pressing the "9" key and interrupted with the "Start/Stop" key.

Immediately after the "9" key has been actuated the printer begins to record the first two channels at 25 mm/s and 20 mm/mV.

Moreover, the message "Testing time constants" appears on the display.

In this test only the first two channels send signals to the printer. No other channels can be configured. The default time constant is output on the printout as a control.

Setting the time constants

As a general rule, the time constants are adjusted to 2.04 s.
This setting can be changed with the CORR_CONFIG key.
The following settings are available:

Time constants:

(1) 0.260s (2) 0.130s (3) 0.033s
(4) 0.004s (5) 4.080s (6) 2.040s
(7) 1.020s (8) 0.510s

The figures given in parentheses indicate the key used to select the corresponding time constant.

Verification of the currently active time constant is only possible by printing out the self-test results (see Section 3.2.6) or by testing the time constants (see Section 5.2.8).

3.2.13 Setting the Device model (V)

This setting is only available from Software Release V4.21.

Setting the Device model is activated by pressing the “V” key, the following display image appears:

D e v i c e M o d e l --> (1) INT

(1) INT (2) USA (3) EU_2 (4) ASIA (5) RUSS

Terminate with any key.

The currently enabled device model appears in line 1.

Warning!

Changing the device model affects the selection of the output formats and languages available, can switch the program for Interpretation from HEART to 12 SL, or vice versa. Also affected are the default configuration settings, e.g., override function yes/no, line filter 50/60 Hz.

The device model configured during manufacture should only be changed when deemed really necessary.

The following table shows the most important combinations:

| Device model | Interpretation | Output formats | Languages | Default settings | |
|--------------|----------------|----------------|-----------|---------------------------|-------------|
| | | | | Override function enabled | Line filter |
| 01 INT | HEART | International | g,e,f,i,s | no | 50 Hz |
| 02 USA | 12SL | USA | e,f,s | yes | 60 Hz |
| 03 EU_2 | 12SL | Int. +H1 +H2 | e,g,f,sw | yes | 50 Hz |
| 04 ASIA | HEART | International | g,e,f | no | 50 Hz |
| 05 RUSS | HEART | International | g,e,f | no | 50 Hz |

Output formats: - International: 12_FS, 12_F1, 12_F2, 6_1F, 6_F2, 3_F1

- USA: STND1, STND3, 1x10, 2x5, 2x5_ex, 3RHY1

The required device model is selected by pressing the appropriate key, "1" to "5".
Quit selection menu by pressing any key.

Quit the self-test with "Shift" + "Man" or "Shift" + "Mode".

CardioSmart configures the appropriate items and initiates a cold start automatically.
Thus, when the device is rebooted automatically all the new settings are adopted.

Selecting a particular device model leads to the configuration of the program for Interpretation as indicated in the table above, even when the interpretation program was configured differently beforehand as described in section 3.2.14.

Note on saved ECGs when configuring the device model:

Saved ECGs are not lost.

If when a new device model is selected, the program for Interpretation is switched over, however, when printing out the saved ECGs the display in the status line (12SL, or without in the case of HEART) is related to the currently enabled program for Interpretation, although the results are based on the program previously configured. In this case, the saved ECGs should be printed out or transferred to a PC before configuration of the new model.

3.2.14 Switching over the program for Interpretation (C)

Switching over the program for interpretation is activated by pressing the “C” key, the following display image appears:

I n t e r p r e t a t i o n --> (1) HEART

Terminate with any key.

The currently enabled device model appears in line 1.

Warning!

Switching over to another program for Interpretation affects the measurement results and the interpretation!

The program for Interpretation configured during manufacture should only be changed when deemed really necessary!

The required program for Interpretation is selected by pressing the appropriate key, “1” or “2”:

“1” for HEART, or “2” for 12SL.

Quit the selection menu by pressing any key.

Quit the self-test with “Shift” + “Man” or “Shift” + “Mode”.

CardioSmart configures the appropriate items and initiates a cold start automatically. Thus, when the device is rebooted automatically all the new settings are adopted.

Note on saved ECGs when switching over the program for Interpretation:

Saved ECGs are not lost.

After switching over to a different program for Interpretation when printing the saved ECGs the display in the status line (12SL, or without in the case of HEART) is related to the newly selected interpretation program, although the results are based on the program previously configured. In this case the saved ECGs should be printed out or transferred to a PC before switching over to the new program for Interpretation.

4. Repair instructions

4.1. Safety instructions

Repairs may only be conducted by approved specialist personnel.
Before opening the CardioSmart switch off device and disconnect mains plug!
Never connect mains plug when the device is open!

Before replacing the primary fuses in the power input module or when resetting the voltage selector switch to a different voltage range the device should also be switched off and the mains plug disconnected!

When replacing electronic components implement **ESD protection**.
Return replaced PCBs in **ESD packaging** only.

Defective batteries, empty batteries should be **disposed of in accordance with the applicable legal stipulations** or returned to the factory.
Batteries returned to the factory should be labeled **'For disposal'**.

4.2 Replacing components

For all the following points the safety instructions in Section 4.1. are to be observed!

The CardioProm (program module) can be unlocked with a coin and then pulled out.
The character display CardioVision Text can be disengaged without opening the device by applying a screwdriver three times left, right and in the middle between the display housing and the upper housing shell of the CardioSmart.

Warning: The device must be switched off!

Opening the device

To open the CardioSmart undo the 5 fastening screws on the underside of the device, open the paper feed flap, carefully raise the upper shell of the housing, disconnect keypad. The display remains in the upper shell of the casing.

During reassembly ensure that no cables are pinched. Exert gentle pressure in the vicinity of the display to ensure that the display is in contact with the PCB and that it remains snapped into position in the upper shell.

Replacing the mains input module

After opening the device undo the screw located on the underside of the housing used to secure the holding angle bracket for the mains input module and the device power-on switch. Disconnect the primary connector DC and switch on/off connector DA from the PCB Power Supply.

Replace mains input module (snap-in connection).

Warning: Only use fuses indicated on the rating plate! (see Chapter 8, Parts List)

For safety reasons insert the fuse holder such that the correct mains voltage range is adjusted!

Visual check at the inspection window of the mains input module. (see Section 2.4)

Replacing the battery

Disconnect and lift out the battery compartment on the underside of the device by raising the middle fastening catch. Remove the battery and disconnect the battery plug DD from the PCB Power Supply.

Warning: Only original Marquette Hellige batteries may be used! (see Chapter 8, Parts List)

Defective batteries, empty batteries should be **disposed of in accordance with the applicable legal stipulations** or returned to the factory.

Batteries returned to the factory should be labeled **“For disposal”**.

Before putting in the new battery insert plug DD. Push in battery compartment, secure by pressing on the fastening catch.

Replacing the PCB Power Supply

After opening the device first of all disconnect the primary connector DC and the switch on/off connector DA from the PCB Power Supply. Open the battery compartment and disconnect the battery connector DD (see "Replacing the battery").

Disconnect connector AE from the PCB Control. Remove PCB Control (see "Replacing the PCB Control").

Undo the 2 fastening screws next to the transformer and fuse Si500.

PCB Control has 2 soldered fuses:

SI 500: Battery fuse

SI 501: Secondary circuit fuse

No adjustments on the PCB Control are necessary.

Replacing the PCB Control

With the transition to device version V1.2 the PCB Control has a new part number:

PCB Control CS_S equipped: 388 032 13, Spare Part No: 389 004 16

Before replacing the PCB, if still possible, printout the settings configured by the user.

After opening the device first of all disconnect the plug-in connector AE to the power supply unit, then disconnect the remaining connectors to the recording unit. Undo the 6 PCBs fastening screws.

Remove shielding plate under the patient input section (1 screw, 3 x snap-in catches) and assemble it on the new PCB.

Check whether the lithium battery for the real-time clock (BA 500) is correctly soldered. Insert new PCB and plug in connectors.

Since the motor controller is located on the PCB Control, the motor adjustment for speed and compensation must be carried out. Motor adjustment is described in Chapter 6.

In addition, disengage the retaining clips by pressing them together, then remove from the upper shell and plug into the connector AG (character display) or AH (graphics display). However, simpler is a character LCD module connected to connector AG via a flat strip cable.

Connect keypad in the upper shell using connector AS.

Setting the time (Section 3.1.10 / 3.2.10).

If known, adopt the user-programmed configuration; otherwise, default setting.

If known, reset the user-programmed time constants; otherwise by outputting the self-test results (Section 3.1.6 / 3.2.6).

Check default setting is 2.04 s, and adjust to 2.04 s if necessary (Section 3.1.12 / 3.2.12)

Adjust contrast with the appropriate display.

Replacing the motor

After replacing the motor the motor adjustment of speed and compensation must be carried out as described in Chapter 6.

Replacing the graphics display

The CardioSmart must be opened to be able to replace the graphics display. Disengage the graphics display from the upper shell of the housing by pressing together the retaining clips and remove it. Reassemble the upper shell. Then carefully insert the graphics display into the basic device until the retaining clips snap securely into place.

Adjust the contrast with the new display.

Replacing the character display

The character display CardioVision Text can be disengaged without opening the device by applying a screwdriver three times to the left, right and in the middle between the display housing and the upper shell of the CardioSmart housing. Carefully insert the new display (character or graphics) into the basic device until the retaining clips snap securely into place. Adjust the contrast with the new display.

Replacing pump module

Open device, take care in removing the upper shell to avoid damage to vacuum tubing. Detach the vacuum tubing from the upper shell. Disconnect power plug SA from the pump module. The module can be replaced after unscrewing the two fastening screws and detaching the evacuated air tubes. After replacement perform a performance check as described in section 7.1.2.11.

Take the negative pressure reading at PTT after approx. 10 s. This value must lie within the range 225...265 mmHg (300...350 mbar).

5. Troubleshooting tips

Device cannot be switched on even though power plug is plugged in

- green power lamp LED is off and device cannot be switched on:
 - power cable defective or not plugged in correctly?
 - primary fuses in the mains input module defective?
 - connector DC on PCB Power Supply plugged in correctly?
 - mains input module defective?
 - keypad via connector AS on PCB Control plugged in correctly?
 - are all above-mentioned points OK: ==> PCB Power Supply defective (transformer or rectifier)
- green power lamp LED is on, but device still cannot be switched on:
 - device on/off switch or its lead defective?
 - connector DA on PCB Power Supply plugged in correctly?
 - keep device on switch depressed and measure voltage UVERS on connector AE pins 4 + 5 on PCB Control: UVERS greater than 25V-?
 - if not: - fuse SI501 on PCB Power Supply defective?
 - if SI501 OK ==> PCB Power Supply defective.
 - if so: ==> PCB Power Supply defective.

Device cannot be switched on when being battery-operated only

- Is the battery empty?
- Plug in power plug, green power lamp LED should glow and the device can be activated, if not, refer to "Device cannot be switched on even though power plug is plugged in".
- By connecting the power plug rapid charging is activated, LED1 on PCB Power Supply should glow, via the PTC resistor R500 a power supply of approx 1 V should be available on PCB Power Supply.
(measured with a multimeter):
 - if so: battery charging circuit OK.
 - if not: battery disconnected or defective?
 - if not: ==> error on PCB Power Supply
- After charging for 10 minutes disconnect power plug, can the device be switched on and trigger a recording?
 - if not: - fuse SI500 on PCB Power Supply defective? (Device does not switch on)
 - if SI500 OK ==> battery defective or capacity too small.

No display on the screen

- does the yellow Stop LED on the keypad glow after the device is switched on?
- if not: - refer to "Device cannot be switched on"
 (if the green Start LED glows CardioProm can also be defective)
- contrast badly adjusted (adjust contrast)
- graphics display: can the background illumination be activated by pressing a key?
- does a beep sound approx. 10s after switching on the CardioSmart? (indicates that the self-test is over)
- if so: - display correctly engaged?
- if so: ====> display defective?
 ====> or driver on the PCB Control defective
- if not: - error in basic device:
 ====> PCB Control defective?
 ====> or CardioProm defective

Error in self-test identified

When an error is detected during the self-test, in addition to the message "Self-test failed" the error code number and a short description also appear on the display. The meaning of the error codes is described in Section 3.1.7 / 3.2.7. The error codes refer to the PCB Control, graphics display and the CardioProm:

- error code 1, 2: PCB Control
- error code 3: PCB Graphics Display
- error code 6,7,8: CardioProm

CardioSmart fails to give printout, no paper transport

Perform the following test in the manual operating mode:

- Paper available? Paper correctly inserted? Paper transport defective?
- Paper feed flap correctly engaged on both sides?
- after pressing the Start key the green Start LED must glow
 - if not: ====> CardioProm defective?
 ====> or PCB Control defective
- Power supply +/-12V at test connector AF/pin 1 and pin 3 of PCB Control available?
 - if not: - Signal REG_ON on high-level? (Measure at connector AE/pin 2)
 - if not: ====> CardioProm defective?
 ====> or PCB Control defective
 - if not: - Signal REG_OFF on high-level (Measure at connector AE/pin 11)
 - if not: ====> Battery depleted or defective?
 ====> or PCB Power Supply defective
- Motor feed line defective? Connector AK plugged in correctly?
- Motor blocked? (check roller, transmission)

- Measuring voltage COMMON_S at test connector AF/pin 4 of the PCB Control available in the range 16V ... 23V?
 - if so: ==> Motor defective?
 - ==> or PCB Control defective (Motor controller)
 - if not: - Battery voltage too low, recharge

Paper transport functions, no printout

- Paper feed flap correctly engaged on both sides?
- Printhead impact pressure switch OK? In addition, in the system test Sections 3.1.4/3.2.4 use the "K" key to check whether opening/closing of the switch is identified.
- if not: - Is connector AJ on PCB Control inserted correctly?
- Good feed line contact on printhead impact pressure switch?
- Printhead impact pressure switch defective?
- Printhead power supply on? Measurement at plus terminal of C500 on PCB Control, voltage in range 16V ... 23V
- if not: - Signal REG_ON on high-level? (Measure at connector AE/pin 2)
- if not: ==> CardioProm defective?
 - ==> or PCB Control defective
- Signal REG_OFF on high-level? (Measure at connector AE/pin 11)
- if not: ==> Battery depleted or defective?
 - ==> or PCB Power Supply defective
- if both signals are on high:
 - ==> PCB Control defective (common/VDD1 switch)
- Printhead power supply cable (connector AO) defective, or not plugged in correctly?
- Printhead data feed line (connector AL) defective, or not plugged in correctly?
- After an extensive printjob, printhead still overheated? Allow to cool.
- none of the above-mentioned errors apply:
 - ==> PCB Control defective?
 - ==> or thermal array printhead defective

CardioSmart only prints on the upper or lower section of the printout

- Paper feed flap is only engaged on one side.

CardioSmart prints, but only baselines are printed out

- Electrodes applied correctly?
- Electrode cables plugged into the patient trunk cable terminal box correctly?
- Patient trunk cable defective (e.g., N core defective)
- Contact problems at the patient input of the CardioSmart?
- none of the above-mentioned errors apply: ==> PCB Control defective

6. Adjustment instructions

CardioSmart only requires motor adjustment. Adjustment is necessary when a new motor (replacement of paper feed flap) or a new PCB Control is implemented.

Motor Adjustment

| What to adjust or to check? | What to measure with? | How to adjust test mode? | Where to turn? | How much and exact? | What else to note? |
|---|--|---|-------------------------|--|---|
| Motor speed under small load, slow down driving roller by applying a finger (recorder flap open), adjust on stationary motor pinion. | Flash at motor pinion with a stroboscope (LED) | 1. Shift + Auto 2. "3" key 3. "A" key 4. With mm/s key adjust to 25mm/s 5. START key | Speed adjuster R505 | Stroboscope frequency = 283.6 Hz When an LED driven by a square-wave generator is used, a small key ratio should be observed. | Motor should be at room temperature be switched on for only 5 s for each adjustment. |
| Adjust motor speed on the slowed down driving roller, so that there is no difference in the RPM. Slow down the driving roller with a finger sufficiently, so that the LED500 just remains illuminated. (Recorder flap open) | see above | see above | Speed compensation R506 | see above | If the RPM has changed after compensation adjustment, it should be corrected again with R505. |

When a stroboscope is not available:

| What to adjust or to check? | What to measure with? | How to adjust test mode? | Where to turn? | How much and exact? | What else to note? |
|--|--|--|-------------------------|---------------------|--|
| Adjust motor speed with the recorder flap closed. | Check the measurements of the grid with a ruler. min. 200mm | 1. Shift + Auto 2. "3" key 3. "A" key 4. With mm/s key adjust to 25mm/s 5. START key | Speed adjuster R505 | < 0,5% | Motor should be at room temperature be switched on for only 5 s for each adjustment. |
| With an increased load adjust motor speed (additional speed reduction, finger on the driving roller until LED 500 just goes out) | see above | see above | Speed compensation R506 | see above | The adjustment procedure is complete when the RPM is 25 mm/s both under normal and increased load. |

7. Servicing and maintenance

7.1 Technical inspection

A technical inspection is to be performed once a year. The following items, including the accessories used, are to be performed:

- Check device and accessories for mechanical defects which impair their function.
- Perform a function check as detailed in Section 3 “System test functions”.
- Check labels and inscriptions on the device relating to safety are clearly legible.
- Measure the resistance of the non-fused earth conductor as per measuring circuit from VDE 0751:1990
- Measure the device leakage current as per measuring circuit from VDE 0751:1990
- Measure the patient leakage current as per measuring circuit from VDE 0751:1990

Warning!

The following checks may only be performed by persons whose training, knowledge and practical experience enable them to carry out such checks reliably and correctly.

Notes:

The operational and functional reliability of the device is checked using the following checklists.

They serve the experienced technician when checking the device.

A knowledge of device operation as detailed in the “User’s Manual” is assumed.

The checklist items are based on the testing instruments given below.

The tests should be carried out using the customer’s accessories, so that defective accessories are also detected automatically.

If other testing instruments are used besides those mentioned, the items on the check list and tolerance specifications may need to be modified.

7.1.1 Visual check

Device and accessories are to be checked to ensure that

- fuse cartridges comply with vendor’s specifications;
- labels and inscriptions on the device relating to safety are clearly legible;
- the mechanical state of the device permits its further use;
- there is no fouling which could cause any reduction in safety.

7.1.2 Test functions

7.1.2.1 Recommended testing instruments and accessories

| | |
|-------------|--|
| 1x | Multi-parameter simulator Lionheart |
| 1x | RS232 interface connector with internal connection between pins 2 and 3 (TXD and RXD) |
| 1x | Patient trunk cable and customer electrode leads, or 1X patient trunk cable, 10-pin 233 402 04 with 1 set of electrode leads (10 leads) with 4 mm plugs 38401129 |
| 1x mmHg) | Pneumatic Transducer Tester PTT, for example: X - Caliber (Display in |
| 1x | 1 x Y-Adapter 303 444 89 |

7.1.2.2 Test preparations

The following descriptions apply to device version V1.2 and software version V4.2.

In general, the device test functions implemented in CardioSmart are used for the tests. These are described in Section “3. Device test functions”.

Connect CardioSmart up to the mains, the green LED for standby must glow.
Switch the device on, the self-test runs automatically, no error message should appear.
When the self-test has finished the device is in the automatic mode, the yellow LED for still disabled operating mode must glow.

Modifications in the user-programmed configuration may need to be made in order to carry out the test. Should such a change need to be made to enable testing, make a printout of all the modified configuration lists and mark the changes made.

Important: After completing the test the original user-programmed configuration is to be retrieved and activated.

7.1.2.3 Operating and display unit performance tests

- Carry out the “Display test (1)” as detailed in Sections 3.1.3 and 3.2.3, respectively.
- Carry out the “Keypad test (2)” as detailed in Sections 3.1.4 and 3.2.4, respectively.

7.1.2.4 Test for recording speeds 25 and 50 mm/s

- Carry out the “Motor test (3)” as detailed in Sections 3.1.5 and 3.2.5, respectively.
- Enable measurement marks by pressing the “A” key.
- Should there be feed speed deviations of > 2%, adjust motor as detailed in Section “6 Adjustment instructions”.

7.1.2.5 Device test result check

- Output “Test results 4)” as detailed in Sections 3.1.6 and 3.2.6, respectively.
 - Main parameters:
 - all memory stores free from error?
 - ASIC test O.K?
 - Sample rate 1000?
 - Selected time constants: 2.04 s (or 4.08 s)?
 - Printhead voltage > 18V, battery charge O.K?
 - Printout clearly legible and without any lapses or

interference?

7.1.2.6 RS232 interface test

- Carry out the “Interface test (6)” test item “(1) Transmitting and receiving”.

Note: For devices with CardioProms V3.0x the RS232 driver is not yet tested in this test; therefore, in this case test item “(3) Transmitting and receiving with device” should be carried out.
- Remote start input test:
 - Adapt push button (NO contact) between RS232 interface pin 8 and pin 5
 - Select manual mode
 - Press push button once: initiates recording in manual mode
 - Press again: terminates recording in manual mode

7.1.2.7 Analysis of the ECG signals and HR value

Carry out the following settings on the ECG simulator:

- Amplitude 1 mV
- Heart rate (RATE) 60 bpm

Connect the electrode leads as indicated below:

| | | | |
|----|--------------|--------|----|
| R | red | -----> | RA |
| L | yellow | -----> | LA |
| F | green | -----> | LL |
| N | black | -----> | RL |
| C1 | white/red | -----> | V1 |
| : | : | : | : |
| C6 | white/violet | -----> | V6 |

Switch in **manual operating mode** and start recording by pressing the Start key.

By pressing the lead scrolling key check whether all leads **are being recorded**.
The ECG traces must be **noise-free**.

Record one page in the “manual” operating mode. The following annotations must be present:

- Heart rate (top right)
- Lower status line - date and time
 - recording speed
 - sensitivity
 - active filter, e.g., 50/60 Hz, 35/20 Hz, ADS
 - frequency range of the recording (not at SW level V3.0x)

The **heart rate** of 60 bpm +/- 2 bpm appears on the display and is printed out on the recording.

Activate the square-wave function on the ECG transmitter at 1 mV.

Using the lead scrolling key select lead II. The square-wave pulse trace must correspond in **amplitude** with the displayed 1 mV reference pulse (applicable to named transmitter only).

Switch back to ECG signal on the ECG transmitter.

Start the recording in manual mode.

Increase the **heart rate to 200 bpm** on the ECG transmitter. The acoustic **warning signal** must sound for about 1s. Reduce the heart rate back to 60 bpm, the warning signal no longer sounds.

7.1.2.8 Pacemaker identification test

Make the following settings on the multifunction simulator:

- pace setting
- pace amplitude 6 mV
- pace duration 0.2 ms

Adjust manual mode on CardioSmart device to be tested and select lead groups I, II and III.

Start the recording. The pace pulses must be visible as needles on the recording output.

7.1.2.9 Identification of disconnected electrodes

Reset the simulator to ECG signal at a heart rate of 60 bpm. Remove one electrode after the other from the ECG transmitter.

Activate the automatic mode in the CardioSmart device to be tested without activating it by pressing the Start key.

Check to ensure that each disconnected electrode is displayed correctly and that an acoustic alarm signal sounds.

7.1.2.10 Checking the charge status of the NC battery

The NC battery can, among other things, be checked as follows:

discharging the battery, followed by charging up completely (duration 4 h), followed by discharging in standby mode without recording.

If the operating time for this procedure is under 2.5 h, the battery should be replaced.

7.1.2.11 Pump leakage

Start the pump by< blowing against the suction pipe still open and then closing the pipe with your finger as shown in the test setup (within 3 s). Take the negative pressure reading at PTT after approx. 10 s. This value must lie within the range 225...265 mmHg (300...350 mbar). Then switch off the negative pressure pump (ECG recorder or power supply unit). The drop in negative pressure must not exceed max. 1 mmHg/s (1.33 mbar/s).

If this value is not observed, the pump in the CardioSmart must be replaced .

7.1.3. Safety Analysis Test

7.1.3.1 General Information

The suggested Safety Analysis Test refer to the international Standard IEC 601-1. The tests are generally performed with Safety Testers, on most of them, the measuring circuits according IEC 601 are already implemented. The tests which have to be performed are described generally, for the handling of your Safety Tester follow the user manual.

The tests may be performed under normal ambient conditions of temperature, humidity and pressure and with line voltage.

The leakage currents correspond to 110 % of rated voltage for the tested unit. Most Safety Testers take this into account, otherwise the measured values have to be calculated.

Recommendet test equipment

- Safety Tester for measurements according IEC 601.
- Testing connector according the following description.

Testing connector for measuring patient leakage current.

For testing the ECG input,a patient cable, with all leads connected together,is used.


7.1.3.2 Protective earth resistance test

The protective earth resistance test has to be performed including its power cord. This test determins whether the device has a power ground fault.

- The protective earth resistance from power connector to any protective earth connected exposed conductive part is measured.
- Specs. of test circuit: AC current source 50 Hz/60 Hz of at least 10 A up to 25 A with limited output voltage of 6 V.
- If resistance is greater than 100 mOhm

7.1.3.3 Measuring of leakage current

To proceed the suggested measurements, the unit under test has to be separated from any interconnection to a system. If the unit is part of a system, extended tests according IEC 601-1-1 have to be performed. The following diagram shows the needed Measuring Circuit [M] for leakage current. The reading in mV corresponds to μA (leakage current). The Safety Testers generally work with this Measuring Circuit [M] and the displayed values are already converted to leakage current.



7.1.3.3.1 Enclosure Leakage Current Test

This test is performed to measure leakage current from chassis to ground during normal conditions (N.C.) and single fault conditions (S.F.C.).

In any case, the leakage current is measured from any exposed conductive parts to ground, the unit under test has to be switched on and off.

Connect the unit under test to your Safety Tester.

- During normal conditions (N.C.), referring to the electrical diagram, measurements have to be done under following conditions:

- * Polarity switch Norm and RVS
- * GND switch GND closed
- * S1 closed and open


- During single fault conditions (S.F.C.), referring to the electrical diagram, the measurements have to be done under following conditions:

- * Polarity switch NORM and RVS
- * GND switch GND open
- * S1 closed

Test has failed if the measured values are greater than:

| | |
|-------------|---------------------------------|
| N.C. | S.F.C |
| 100 μ A | 500 μ A |
| | 300 μ A (U.L. requirements) |

Electrical diagram for Enclosure Leakage Current Test



7.1.3.3.2 Patient Leakage Current Test

This test performs a leakage current test under single fault conditions (S.F.C.) dependent of domestic power outlet with 115 or 230 V AC as source into the floating inputs.

In any case, the leakage current is measured from Input Jack, of unit under test, to ground.


For testing the ECG input, a patient cable, with all leads connected together, is used.

Connect the unit under test to your Safety Tester.

- Referring to the electrical diagram, measurements have to be done under following conditions
 - * Polarity switch NORM and RVS
 - * GND switch GND closed
 - * S1 closed

Test has failed if the measured values are greater than 50 μ A

Electrical Diagram for Patient Leakage Current Test



For protection of the test person the following values of resistor R may be used:

| | |
|--------|-------------------------|
| Typ BF | 22 kOhm (120 to 130 V) |
| | 47 kOhm (220 to 240 V) |
| Typ CF | 100 kOhm (220 to 240 V) |

7.2 Maintenance, cleaning, disinfection

CardioSmart maintenance, cleaning, disinfection is performed in accordance with the CardioSmart or CardioSmart ST User's Manual, Chapter 11, "Cleaning, disinfection and maintenance", as applicable.

Applications requiring extensive recordings may result in deposits on the thermal array printhead which normally do not, however, have any adverse effect on the printing quality. This can be removed with a soft, fluff-free cloth soaked in an alcohol-based cleaning agent (e.g., surgical spirit).

8. Parts Lists

| Part Number | Description CardioSmart | | |
|-------------|----------------------------|-------------|-------------------------|
| 101 116 01 | CardioSmart | (int.) | Basic unit without Pump |
| 101 116 02 | CardioSmart | (USA) | Basic unit without Pump |
| 101 116 03 | CardioSmart | (russ.) | Basic unit without Pump |
| 101 116 04 | CardioSmart | (int), | Basic unit with Pump |
| 101 116 05 | CardioSmart | (russ), | Basic unit with Pump |
| 101 116 06 | CardioSmart | (Europe2) | Basic unit without Pump |
| 101 116 07 | CardioSmart | (Europe2) | Basic unit with Pump |

CardioSmart ST

| | | | |
|------------|----------------|-------------|--|
| 101 116 10 | CardioSmart ST | (int.) | Basic unit Stress Test without Pump |
| 101 116 11 | CardioSmart ST | (int.) | Basic unit Stress Test with Pump |
| 101 116 12 | CardioSmart ST | (USA.) | Basic unit Stress Test without Pump |
| 101 116 13 | CardioSmart ST | (USA.) | Basic unit Stress Test with Pump |
| 101 116 14 | CardioSmart ST | (russ.) | Basic unit Stress Test without Pump |
| 101 116 15 | CardioSmart ST | (russ) | Basic unit Stress Test with Pump |
| 101 116 16 | CardioSmart ST | (Europe 2.) | Basic unit Stress Test without Pump |
| 101 116 17 | CardioSmart ST | (Europe 2.) | Basic unit Stress Test with Pump |
| 101 116 18 | CardioSmart ST | (ASIA) | Basic unit Stress Test without Pump |
| 101 116 19 | CardioSmart ST | (ASIA) | Basic unit Stress Test with Pump |

User Manual CardioSmart:

| | |
|------------|-----------------------------------|
| 227 435 15 | User Manual CardioSmart (d) V 1.1 |
| 227 435 16 | User Manual CardioSmart (e) V 1.1 |
| 227 435 17 | User Manual CardioSmart (f) V 1.1 |
| 227 435 41 | User Manual CardioSmart (d) V 1.3 |
| 227 435 42 | User Manual CardioSmart (e) V 1.3 |
| 227 435 43 | User Manual CardioSmart (f) V 1.3 |

| | |
|------------|--|
| 227 435 45 | User Manual CardioSmart (i) V 1.3 |
| 227 435 46 | User Manual CardioSmart (sp) V 1.3 |
| 227 43518 | User Manual CardioSmart (USA) V 1.3 |
| 227 435 52 | User Manual CardioSmart (e) for UK V 1.3 |
| 227 435 53 | User Manual CardioSmart (f) for BeneluxV 1.3 |
| 227 435 58 | User Manual CardioSmart (sw) V 1.3 |

User Manual CardioSmart ST:

| | |
|------------|--|
| 227 465 01 | User Manual Cardio Smart ST (d) V. 1.2 |
| 227 465 02 | User Manual Cardio Smart ST (e) V. 1.2 |
| 227 465 03 | User Manual Cardio Smart ST (f) V. 1.2 |
| 227 465 11 | User Manual Cardio Smart ST (d) V. 1.3 |
| 227 465 12 | User Manual Cardio Smart ST (e) V. 1.3 |
| 227 465 13 | User Manual Cardio Smart ST (f) V. 1.3 |
| 227 465 15 | User Manual Cardio Smart ST (i) V. 1.3 |
| 227 465 16 | User Manual Cardio Smart ST (sp) V. 1.3 |
| 227 465 24 | User Manual Cardio Smart ST (US) V. 1.3 |
| 227 465 17 | User Manual Cardio Smart ST (russ) V. 1.3 |
| 227 465 32 | User Manual Cardio Smart ST (e) for UK V. 1.3 |
| 227 465 33 | User Manual Cardio Smart ST (f) for Benelux V. 1.3 |
| 227 465 38 | User Manual Cardio Smart ST (sw) V. 1.3 |
| 227 465 42 | User Manual Cardio Smart ST (asia) V. 1.3 |
| 227 435 11 | Service Manual .d/e, 101116.. Cardiosmart |

Housing / Driving unit

| | |
|------------|-------------------------------------|
| 432 521 90 | Unter -case shell |
| 924 017 14 | Rubber foot |
| 432 521 98 | Battery cover |
| 303 445 33 | Upper case with keypad (russ) |
| 303 443 18 | Upper case with keypad (Int) |
| 303 444 50 | Upper case with keypad (USA) |
| 303 446 14 | Upper case with keypad (Europe 2) |
| 303 445 34 | Upper case with keypad ST (Int) |
| 303 446 16 | Upper case with keypad ST (Russ) |
| 303 44617 | Upper case with keypad ST (Europe2) |
| 303 44615 | Upper case with keypad ST (US) |
| 303 446 13 | Upper case with keypa ST (ASIA) |
| 303 443 19 | Handle |
| 303 443 20 | Drive Flap |
| 303 438 29 | Code reader |

| | |
|------------|---------------------------------|
| 303 442 79 | Motor |
| 914 326 69 | Microswitch |
| 432 523 54 | Stop for Microswitch |
| 801 777 62 | Screw for Stop |
| 432 522 05 | Plate for Driver Flap |
| 504 657 60 | Tinplate |
| 929 166 45 | Luer- connection kit (female) |
| 929 166 54 | Luer-cap |

Electronic components / pcb's

| | |
|------------|--|
| 929 166 40 | Line filter with fuse holder |
| 914 325 94 | Rocker switch |
| 303 442 71 | Comb |
| 303 442 70 | Battery |
| 388 029 10 | pcb.Printhead connection |
| 388 028 14 | pcb. Power Supply |
| 912 084 37 | FuseT 4 A TR5 |
| 912 084 35 | FuseT 1,25 A |
| 912 084 04 | FuseT 0,25 A |
| 912 084 14 | FuseT 0,5 A |
| 388 028 12 | pcb. Control (usable up to V 1.1) |
| 389 004 05 | Exchange pcb. Control (usable up to V 1.1) |
| 303 446 18 | pcb. Control ST (compatible for all CardioSmarts) |
| 389 004 16 | Exchange pcb. Control ST (compatible for all CardioSmarts) |
| 929 165 73 | Battery 3 V 0,8 Ah |
| 303 442 91 | Pump 12 V (Standard) |

Display-Modul Graphics

| | |
|------------|--|
| 235 066 01 | CardioVision Graphics -Module for CardioSmart |
| 235 066 02 | CardioVision Graphics -Module for CardioSmart(USA) |
| 388 028 87 | pcb Mother board |
| 388 028 11 | pcb digital Display |
| 930 117 17 | LCD-Graphics Display |
| 303 443 1 | Upper-Case Part |
| 303 442 80 | Joint |
| 432 522 01 | Display Bottom Part |
| 423 203 04 | Insulation |

Display-Modul Text

| | | |
|------------|--|---------|
| 235 06 701 | CardioVision Text Display-Modul for CardioSmart | |
| 235 067 02 | CardioVision Text Display-Modul for CardioSmart | (USA) |
| 388 028 88 | pcb. Mother board | |
| 303 442 87 | LCD Modul | |
| 432 522 03 | Display Top Part | |
| 432 522 04 | Display Bottom Part | |
| 931 098 80 | Instrument Bag for CardioSmart, nylon material, with pockets for accessories | |

Connectioncable (for CardioSmart)

| | |
|------------|---|
| 223 330 04 | Connection cable M(L)700 - C'Smart(ST) 5m lang |
| 223 362 03 | Connection cable 5 m long for CardioSmart to PC |
| 223 366 04 | Connection cable EC1200-CardioSmart(ST) 5m long |
| 223 378 01 | Connection cable for ECG- Transmission CardioSmart to Modem ELSA,9 pin.,3m long |
| 223 378 02 | Connection cable for ECG- Transmission CardioSmart to Modem Multitech 25 pin.,3m long |

Connectioncable (Only for CardioSmart-ST)

| | |
|------------|---|
| 223 362 04 | Connection cable EC(B)560/561-C'Smart ST 5m long |
| 223 380 02 | Connection cable ERG900(L) - C'Smart ST 5m long |
| 223 368 02 | Connection cable TM400AC - C'Smart ST 5m long |
| 223 390 03 | Connection cable T2000-CardioSmart ST 5m long, (RS232 to RS 232) |
| 384 017 67 | Upgrade Interface Kit T 2000 (RS 422 to RS 232) —> CardioSmart for 230 V |
| 384 017 68 | Upgrade Interface Kit T 2000 (RS 422 to RS 232) —> CardioSmart für 115 V |

9. Technical Specifications

For operation of the CardioSmart it is necessary to plug in one CardioProm module and one of the two available displays. The following technical specifications apply to all CardioSmart versions.

Recording

direct recording of waveforms and alphanumeric characters with rectangular coordinates by means of thermal-array printhead printing on thermosensitive paper.

- * recording channels: 3 or 6, 12 in automatic mode; overlapping
- * space between baselines3 channels: 62 mm (arrhythmia)
6 channels: 31 mm (manual)
12 channels: 16 mm (auto.)
- * writing width 200 mm max.
- * annotation of device settings, date, time of day and patient name, if entered, in margin of recording strip
- * with appropriate software module documentation of analysis results and report in each mode on separate pages
- * resolution of the recording:
vertically 8 dots/mm
horizontally 25 µm at 25 mm/s

Chart paper

Marquette Hellige CONTRAST® Z-fold pad, 150 pages per pad, equivalent to a chart length of 45 m

paper width 210 mm or 8.5 in. (US format)
sheet length 300 mm or 11 in. (US format)

To prevent debris from collecting on the printhead, use the Marquette Hellige CONTRAST chart paper only!

Chart transport

- * paper speed
5 – 25 – 50 mm/s (key-selectable)
error limits at 25 and 50 mm/s typically $\pm 1\%$
at 5 mm/s $\pm 10\%$ max.
- * when supply of chart paper is depleted, recorder switches off and sounds an alarm; the last pages of the pad bear a red stripe in the margin

Membrane keypad

membrane keys with tactile feedback

- * function keys for all important routine operations
- * alphanumeric keypad for entry of texts

Display

alphanumeric LCD, 2 x 40 characters, contrast adjustable, or

graphics display for 24 x 40 characters, contrast adjustable

resolution of 320 x 240 pixels, backlit

Indicators (LEDs)

for mains supply, battery charge level and start/stop function

Lead selection

manual selection of individual lead groups or automatic sequencing of lead groups

- * lead programs (c):

EINTHOVEN, GOLDBERGER, WILSON, NEHB, CABRERA
lead sequence, custom lead programs in manual
mode

Automatic functions

assisting and facilitating device operation:

- * automatic amplifier blocking
- * automatic lead sequencing, paper feed, calibration (c)
- * report formatting (c)
- * automatic baseline adjustment
- * anti-drift system compensating for polarization voltage fluctuations (c)

Detection of pacing pulses

- * pulse duration between 0.1 and 2.5 ms
- * marks irrespective of pulse polarity
- * pulse amplitude between ± 5 mV and ± 700 mV

Heart-rate display

derivation of the heart rate from all ECG signals

- * display range 30 to 300 BPM
- * display update with every heart beat, max. every 2 s

Signal inputs

isolated patient signal input, type CF according to IEC, high-voltage protection for all lead connections and N-electrode connection, interference compensation via neutral electrode (N), monitoring for detection of open leads

- * electrode connections for R, L, F, N, C1 to C6, Nax, Nst and Nap (=C4)
- * input impedance for differential signals applied between any two electrode connections $>10\text{ M}\Omega$ for 10 Hz
- * input impedance for common-mode signals referred to N $>50\text{ M}\Omega$ for up to 60 Hz
- * dynamic range for differential signals between any two electrode connections for AC voltage ± 10 mV, for superimposed DC voltage (polarization voltage) ± 600 mV
- * dynamic range for common-mode signals referred to N ± 1 V, referred to chassis 263 V AC voltage (rms)
- * quiescent input current via any electrode connection for 1-k Ω termination referred to N <50 nA
- * patient leakage current (rms values) according to IEC class CF: under normal conditions $<10\text{ }\mu\text{A}$, in single-fault condition (e.g., patient in contact with line voltage) $<20\text{ }\mu\text{A}$
- * non-destructive range for lead-electrode connections and the N-connection referred to N ± 50 V, referred to chassis ± 1500 V
- * pulse voltage resistance of all lead-electrode connections and of the N-connection referred to chassis (e.g., defibrillation) 5000 V
- * monitoring of each electrode for disconnection: R, L, F, N, C1, C2, C3, C4, C5, C6, Nap, Nax, Nst audible alarm signal upon program start

Data interface

one serial RS 232 interface for exchange of data with adapted peripherals and software handshake

RS 232 interface (standard V.24 interface):

- * input voltage range ± 15 V max.
- * output voltage range ± 5 V min.
- * ESD interface protection for up to ± 10 kV

Remote start

remote start connection for paper feed (irrespective of selected operating mode), external make contact referred to chassis via circuit reference:

- * source impedance $R_i < 300 \Omega$
- * dwell of contact > 100 ms
- * non-destructive load ± 10 V
- * ESD interface protection up to ± 10 kV max.

Pin assignment of data port

- 2 RXDE
- 3 TXDE
- 5 circuit reference
- 8 remote start

Signal transmission

Patient signal input to recording

after lead formation and digitization simultaneous transmission of all electrode signals to the digital processing system; muscle filter, AC filter, pacing pulse identification, automatic or manual sensitivity adjustment, automatic baseline adjustment and drift compensation by means of the anti-drift system (A.D.S.) can be enabled or disabled simultaneously for all channels; digital output of processed signals via thermal-array printhead

- * low cut-off frequency (-3-dB limits) 0.08 Hz, equivalent to a time constant of 2.04 s
- * high cut-off frequency (-3-dB limits) in Auto and manual modes: 150 Hz (IEC/AHA) in Arrhy mode: 100 Hz (IEC)
- * signal sampling rate 1000/s
- * resolution, referred to the input, 5 μ V
- * output rate to recorder 2000/s
- * for all leads adjustment of sensitivity in four steps: 40-20-10-5 mm/mV
- * with muscle filter (low-pass characteristic) switched on, 3-dB drop of the amplitude frequency response at approx. 35 Hz or 20 Hz
- * with AC filter switched on, identification and compensation of periodic 50 or 60-Hz frequency components (depending on recorder model: attenuation > 40 dB)
- * non-linear distortion below values specified in IEC and AHA recommendations
- * coincidence error limits between any two channels ± 0.5 mm
- * identification of pacing pulses in C2 or other C-leads and marking in all channels for signals referred to patient input: duration $+ 0.1$ ms, amplitude > 5 mV
- * noise in the signal-transmission path below values specified in IEC and AHA recommendations: $\leq 2.5 \mu$ V rms
- * common-mode rejection for 50 or 60-Hz signals (depending on CardioSmart model) with AC filter switched on > 140 dB

ECG calibration

automatic recording of a defined voltage step, valid for all channels

- * calibration voltage of 1 mV, referred to ECG signal input

pulse width depending on paper speed

25 mm/s 5 mm

50 mm/s 10 mm

5 mm/s 1 mm

Automatic sensitivity adjustment of ECG signals

automatic adaption of the signal gain in dependence of the incoming signal; highest amplitude of a lead group or of all leads determines gain

- * automatic adjustment range 5 to 40 mm/mV
- * amplitude range (6 channels) 18 to 31 mm

Baseline

automatic adjustment of the baseline to the optimal recording range, in dependence of the signal amplitude

Anti-Drift System (A.D.S.)

automatic compensation of baseline fluctuations caused by polarization voltage fluctuations at the lead electrodes (recording delay 4.2 s)

ECG storage

- * in automatic mode storage of up to 45 ECGs
- * stored ECGs can be deleted (individually or all in one pass), printed, transferred or the patient data can be modified
- * when memory is full user is informed of the possible actions

Blocking

rapid charge reversal of the coupling capacitors in the preamplifiers after electrode application

Electrode monitoring

audible and visual indication on the LCD of disconnected electrodes or line break; each single electrode is monitored

Text input

patient and user data as well as comments can be entered via the panel keyboard and are annotated on the recording strip

Copy function

in the automatic mode, after ECG recording, copies of the ECG can be printed from memory and/or transferred to a connected PC (c)

Test

automatic performance test upon power up, including verification of the signal path starting at the signal input

stored test data for demonstration of the device functions

Power

from the power line or from a built-in rechargeable battery, automatic switchover; automatic battery charging during line-power operation; no recording when battery depleted or missing

Line-power operation

- * instrument design in protection class I according to IEC 601-1, internally powered equipment

| | | |
|---|------------------------------|--|
| rated voltage range | 220 to 240 V | Operational readiness |
| * operating voltage range | 198 to 264 V; 49 to 65 Hz | after successful selftest within 10 s of power-up |
| * rated current fuse 2 x T 0.25 A, 5 x 20 | 0.16 A | Operating position |
| rated voltage range | 110 to 120 V | horizontal |
| * operating voltage range | 99 to 132 V; 49 to 65 Hz | Environment |
| * rated current fuse 2 x D 500 mA, 3 A-G | 0.32 A | Operation |
| * power consumption typically: during battery charging | 15 W | * ambient temperature between +10 and +40 °C |
| device turned on | 13 W | * rel. humidity between 25 and 95% |
| * max. power consumption | 28 W | * atmospheric pressure between 700 and 1060 hPa |
| Battery-power operation | | Storage and transport |
| * battery type nickel cadmium | | * ambient temperature between -30 and +60 °C (also with battery) |
| * rated battery voltage | 18 V | * rel. humidity between 25 and 95% |
| * rated battery capacity | 1.4 Ah | * atmospheric pressure between 500 and 1060 hPa |
| * fully charged battery sufficient for up to 50 automatic ECG recordings, if device is switched on only for recording | | Case dimensions |
| * battery charging time approx. 4 hours | | * width 370 mm |
| * (min. charging time for recording of one automatic ECG 10 min) | | * height 95 mm |
| | | * depth 320 mm (incl. handle) |
| * battery life approx. 2 to 3 years, replacement by service technician only | | Weight |
| * battery for built-in clock: lithium battery, life of approx. 5 years, replacement by service technician | | CardioSmart with battery and text display approx. 5.3 kg CardioSmart with battery and graphics display approx. 5.6 kg |

CardioSmart ST

For operation of the CardioSmart ST it is necessary to plug in one CardioProm module and one of the two available displays. The following technical specifications apply to all CardioSmart ST versions.

Recording

direct recording of waveforms and alphanumeric characters with rectangular coordinates by means of thermal-array printhead printing on thermosensitive paper.

- * recording channels: 6, 12 in automatic mode; overlapping
- * space between baselines 6 channels: 31 mm (manual)
12 channels: 16 mm (auto.)
- * writing width 200 mm max.
- * annotation of device settings, date, time of day and patient name, if entered, in margins of recording strip
- * with appropriate software module documentation of analysis results and report in each mode on separate pages
- * resolution of the recording:
vertically 8 dots/mm
horizontally 25 µm at 25 mm/s

Chart paper

HELLIGE CONTRAST® Z-fold pad, 150 pages per pad, equivalent to a chart length of 45 m

paper width 210 mm or 8.5 in. (US format)
sheet length 300 mm or 11 in. (US format)

To prevent debris from collecting on the printhead, use the HELLIGE CONTRAST chart paper only!

Chart transport

- * paper speed
5 – 25 – 50 mm/s (key-selectable)
error limits at 25 and 50 mm/s typically $\pm 1\%$
at 5 mm/s $\pm 10\%$ max.
- * when supply of chart paper is depleted, recorder switches off and sounds an alarm; the last pages of the pad bear a red stripe in the margin

Membrane keypad

membrane keys with tactile feedback

- * function keys for all important routine operations
- * alphanumeric keypad for entry of texts

Display

graphics display for 24 x 40 characters, contrast adjustable

resolution of 320 x 240 pixels, backlit

Indicators (LEDs)

for mains supply, battery condition and start/stop function

Lead selection

manual selection of individual lead groups or automatic sequencing of lead groups

- * lead programs (c):

EINTHOVEN, GOLDBERGER, WILSON, NEHB, CABRERA
lead sequence and custom lead programs in manual and stress test mode

Automatic functions

assisting and facilitating device operation:

- * automatic amplifier blocking
- * automatic lead sequencing, paper feed, calibration (c)
- * report formatting (c)
- * automatic baseline adjustment
- * anti-drift system compensating for polarization voltage fluctuations (c)

Detection of pacing pulses

- * pulse duration between 0.1 and 2.5 ms
- * marks irrespective of pulse polarity
- * pulse amplitude between ± 5 mV and ± 700 mV

Heart-rate display

derivation of the heart rate from all ECG signals

- * display range 30 to 300 BPM
- * display update with every heart beat, max. every 2 s

Signal inputs

isolated patient signal input, type CF according to IEC, high-voltage protection for all lead connections and N-electrode connection, interference compensation via neutral electrode (N), monitoring for detection of open leads

- * electrode connections for R, L, F, N, C1 to C6, Nax, Nst and Nap (=C4)

- * input impedance for differential signals applied between any two electrode connections $>10\text{ M}\Omega$ for 10 Hz
- * input impedance for common-mode signals referred to N $>50\text{ M}\Omega$ for up to 60 Hz
- * dynamic range for differential signals between any two electrode connections for AC voltage ± 10 mV, for superimposed DC voltage (polarization voltage) ± 600 mV
- * dynamic range for common-mode signals referred to N ± 1 V, referred to chassis 263 V AC voltage (rms)
- * quiescent input current via any electrode connection for 1-k Ω termination referred to N <50 nA
- * patient leakage current (rms values) according to IEC class CF: under normal conditions $<10\text{ }\mu\text{A}$, in single-fault condition (e.g., patient in contact with line voltage) $<20\text{ }\mu\text{A}$
- * non-destructive range for lead-electrode connections and the N-connection referred to N ± 50 V, referred to chassis ± 1500 V
- * pulse voltage resistance of all lead-electrode connections and of the N-connection referred to chassis (e.g., defibrillation) 5000 V
- * monitoring of each electrode for disconnection: R, L, F, N, C1, C2, C3, C4, C5, C6, Nap, Nax, Nst audible alarm signal upon program start

Data interface

one serial RS 232 interface for exchange of data with adapted peripherals and software handshake

RS 232 interface (standard V.24 interface):

- * input voltage range ± 15 V max.

- * output voltage range ± 5 V min.
- * ESD interface protection for up to ± 10 kV

Remote start (hardware)

remote start connection for paper feed (depending on selected operating mode), external make contact referred to chassis via circuit reference:

- * source impedance $R_i < 300 \Omega$
- * dwell of contact > 100 ms
- * non-destructive load ± 10 V
- * ESD interface protection up to ± 10 kV max.

Pin assignment of data port

- 2 RXDE
- 3 TXDE
- 5 circuit reference
- 8 remote start

Signal transmission

Patient signal input to recording

after lead formation and digitization simultaneous transmission of all electrode signals to the digital processing system; muscle filter, AC filter, pacing pulse identification, automatic or manual sensitivity adjustment, automatic baseline adjustment and drift compensation by means of the anti-drift system (A.D.S.) can be enabled or disabled simultaneously for all channels; digital output of processed signals via thermal-array printhead

- * low cut-off frequency (-3-dB limits) 0.08 Hz, equivalent to a time constant of 2.04 s

- * high cut-off frequency (-3-dB limits)
in Auto and manual modes: 150 Hz (IEC/AHA)
in stress test mode: 75 Hz (IEC)
- * signal sampling rate 1000/s
- * resolution, referred to the input, 5 μ V
- * output rate to recorder 2000/s
- * for all leads adjustment of sensitivity in four steps: 40-20-10-5 mm/mV
- * with muscle filter (low-pass characteristic) switched on, 3-dB drop of the amplitude frequency response at approx. 35 Hz or 20 Hz
- * with AC filter switched on, identification and compensation of periodic 50 or 60-Hz frequency components (depending on recorder model: attenuation > 40 dB)
- * non-linear distortion below values specified in IEC and AHA recommendations
- * coincidence error limits between any two channels ± 0.5 mm
- * identification of pacing pulses in C2 or other C-leads and marking in all channels for signals referred to patient input: duration + 0.1 ms, amplitude > 5 mV
- * noise in the signal-transmission path below values specified in IEC and AHA recommendations: $\leq 2.5 \mu$ V rms
- * common-mode rejection for 50 or 60-Hz signals (depending on CardioSmart ST model) with AC filter switched on > 140 dB

ECG calibration

automatic recording of a defined voltage step, valid for all channels

- * calibration voltage of 1 mV, referred to ECG signal input

pulse width depending on paper speed

25 mm/s 5 mm

50 mm/s 10 mm

5 mm/s 1 mm

Automatic sensitivity adjustment of ECG signals

automatic adaption of the signal gain in dependence of the incoming signal; highest amplitude of a lead group or of all leads determines gain

- * automatic adjustment range 5 to 40 mm/mV
- * amplitude range (6 channels) 18 to 31 mm

Baseline

automatic adjustment of the baseline to the optimal recording range, in dependence of the signal amplitude

Anti-Drift System (A.D.S.)

automatic compensation of baseline fluctuations caused by polarization voltage fluctuations at the lead electrodes (delay in recording: 4.2 s)

ECG storage

- * in automatic mode storage of up to 45 ECGs
- * stored ECGs can be deleted (individually or all in one pass), printed, transferred or the patient data can be modified

- * when memory is full user is informed of the possible actions

Blocking

rapid charge reversal of the coupling capacitors in the preamplifiers after electrode application

Electrode monitoring

audible and visual indication on the LCD of disconnected electrodes or line break; each single electrode is monitored

Text input

patient and user data as well as comments can be entered via the panel keyboard and are annotated on the recording strip

Copy function

in the automatic mode, after ECG recording, copies of the ECG can be printed from memory and/or transferred to a connected PC (c)

Test

automatic performance test upon power up, including verification of the signal path starting at the signal input

stored test data for demonstration of the device functions

Power

from the power line or from a built-in rechargeable battery, automatic switchover; automatic battery charging during line-power operation; no recording when battery depleted or missing

Line-power operation

- * instrument design in protection class I according to IEC 601-1
- * **rated voltage range** **220 to 240 V**
- * operating voltage range 198 to 264 V;
49 to 65 Hz
- * rated current 0.16 A
fuse 2x slow-blow 0.25 A, 5 x 20
- * **rated voltage range** **110 to 120 V**
- * operating voltage range 99 to 132 V;
49 to 65 Hz
- * rated current 0.32 A
fuse 2 x D 500 mA, 3AG
- * power consumption typically:
during battery charging 15 W
device turned on 13 W
- * max. power consumption 28 W

Battery-power operation

- * rated battery voltage 18 V
- * rated battery capacity 1.4 Ah
- * fully charged battery sufficient for up to 50 automatic ECG recordings, if device is switched on only for recording
- * battery charging time approx. 4 hours
- * (min. charging time for recording of one automatic ECG 10 min)
- * battery life approx. 2 to 3 years, replacement by service technician only

- * battery for built-in clock: lithium battery, life of approx. 5 years, replacement by service technician

Operational readiness

after successful selftest within 10 s of power-up

Operating position

horizontal

Environment

Operation

- * ambient temperature between +10 and +40 °C
- * rel. humidity between 25 and 95%
- * atmospheric pressure between 700 and 1060 hPa

Storage and transport

- * ambient temperature between -30 and +60 °C (also with battery)
- * rel. humidity between 25 and 95%
- * atmospheric pressure between 500 and 1060 hPa

Case dimensions

- * width 370 mm
- * height 95 mm
- * depth 320 mm (incl. handle)

Weight

CardioSmart ST with battery and graphics display
approx. 5.6 kg

10. Device Documents

The following documents are enclosed:

- Entire instrument mechanical specifications:

| | | |
|--|----------|--|
| - Assembly drawing: | | 101 116 01...19 |
| - Entire instrument wiring specifications: | | 101 116 01...07/10...19 S Sheet 1 |
| - Master Record Index: | MRI V1.0 | 101 116 01...03 S Sh. 2 |
| | MRI V1.1 | 101 116 01...05 S Sh. 3 |
| | MRI V1.2 | 101 116 01...05/10...13 S Sheet 4 |
| | MRI V1.3 | 101 116 01...07/10...19 S Sheet 5,6 |
| - PCB Control CS: | | 388 028 12 P Sh. 1/8 and Sheet 8/8 388 028 12 R Sh. 1/2 and Sheet 2/2 |
| - PCB Control CS_S: | | 388 032 13 P Sh. 1/8 and Sheet 8/8 388 032 13 R Sh. 1/2 and Sheet 2/2 |
| - PCB Power Supply CS: | | 388 028 14 P 388 028 14 R |

Note on Master Record Index (MRI)

| MRI | modified/new functions | PCB Control used |
|-----------|---|------------------------------|
| MRI V1.0 | CardioSmart w/o remote start input Device types 101 116 01...03 Device status V1.0 | 388 028 12, PCB Control CS |
| MRI V.1.1 | CardioSmart with remote start input and optional suction pump Device types 101 116 01...05 Device status V1.1 | 388 032 12, PCB Control CS_S |
| MRI V.1.1 | CardioSmart with remote start input and optional suction pump Device types 101 116 01...05 Device status V1.1 | 388 032 13, PCB Control CS_S |
| MRI V.1.2 | CardioSmart ST with remote start input and optional suction pump Device types 101 116 10...13 Device status V1.2 | 388 032 13, PCB Control CS_S |
| MRI V.1.3 | CardioSmart with remote start input and optional suction pump, additional Regional models Device types 101 116 10...13 Device status V1.2 | 388 032 13, PCB Control CS_S |

and

| | | |
|-----------|--|------------------------------|
| MRI V.1.3 | CardioSmart ST with remote start input and optional suction pump, additional Regional models Device types 101 116 10...13 Device status V1.2 | 388 032 13, PCB Control CS_S |
|-----------|--|------------------------------|

For reference numbers, compatibilities and correspondence with the respective versions of the device, please refer to the appropriate Master Record Index.

11. Appendix

Schutzvermerk nach DIN 34 beachten

Sicherungen:
fuses:

für Variante 01/03_07/10/11/14_19 220V
for version 01/03_07/10/11/14_19 220V
912 084 04 10.25A
für Variante 02/12/13/18/19 15V
for version 02/12/13/18/19 15V
912 084 14 0500mA
für Variante 18/19 mit Hinweisblatt 440 083 01
verpackt in Plastikbeutel 931 098 31

416 118 22
504 657 19
M3 x 10 DIN7985
914 325 94
432 522 13
3 x 10 Plastile 45 Senk
2x SicherRing 6x0,7 DIN471

In Variante 01/03_07/10/11/14_19
for version 01/03_07/10/11/14_19
432 522 29
3 x 14 Plastile5 Senk

2x 925 096 31
432 522 09
432 523 54
3 x 25 Plastile 45
Nase oben, Anschlag
des Mikroschalterhebels

In Variante 01/03_07/10/11/14_19
for version 01/03_07/10/11/14_19
432 522 12
2x 2,2 x 6 Plastile Unsenk

432 522 34
3 x 10 Plastile 45 Senk
504 657 74
2x M3 x 6 DIN7985

388 029 10
915 414 86 6-pol. auf Lpt.
6-pin on pcb
915 414 82 6-pol. auf Kamm
6-pin on printhead

915 416 21 6-pol. auf Kamm
6-pin on printhead

416 118 18
2x 416 118 19

504 657 25
3x 923 081 00
432 522 10
915 415 34 3-pol.
915 414 82

432 521 99
303 442 72

432 522 11
gekittet
glued

405 174 01
M2,5 x 8 DIN965
432 521 90
4x M4 x 40 DIN7985 mit Schaft
1x M4 x 20 DIN7985

303 438 29
303 442 71
6x 415 154 39
3 links montiert
3 rechts montiert

504 657 26
2x M3 x 8 DIN965
504 657 23
9x 3 x 12 Plastile 45

927 230 20
In Variante 02/12/13/18/19
for version 02/12/13/18/19
Siehe Beschriftungsblatt
see labels catalog

430 518 20
In Variante 02/12/13
for version 02/12/13

430 501 00
In Variante 02/12/13
for version 02/12/13

430 515 47
über Falscheneingang geklebt
nur in Variante 02
for version 02

430 515 49
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 59
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 67
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 75
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 83
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 91
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 99
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

929 166 40
923 096 69
923 081 00
914 326 07
2,2 x 12 Plastile 45

915 417 89
M3 x 6 DIN7985
B3 DIN127 (unter Schraubenkopf)
(under screw head)
A3,2 DIN6798
(zwischen Kabelschuh und Druckbalken)
(between cable lug and press bar)

915 415 34 3-pol.
915 414 82
388 028 14
2x M3 x 6 DIN7985

In Variante 02/12/13
for version 02/12/13
927 229 01
In Variante 01/03_07/10/11/14_19
for version 01/03_07/10/11/14_19
504 657 24

432 522 06
M 2,1
A4,3 DIN6798
M4 x 10 DIN7985
B4 DIN127
916 118 00

an beiden Enden montiert
clipped to both ends

In Variante 04/05/07/11/13/15/17/19
for version 04/05/07/11/13/15/17/19

303 442 91
923 096 84
388 031 25
383 273 27
923 096 65
927 228 48

388 032 13
388-028-12
7x M3 x 6 DIN7985

432 521 91
2x ZyLSH ISO 2338-B-4x24
805 124 16

2x 432 521 92
2x 421 032 54

927 229 55
In Batteriefach geklebt
glued into battery compartment

Beschriftung siehe BeschBlatt
nach Variante
mit Etikettendrucker bedruckt

829 074 08 2x 70mm
In Batteriefach geklebt
glued into battery compartment
ca. 190mm Abstand, mittig

433 671 66
404 029 00
in Lpt. Abschirmung gekittet
riveted to pcb shielding

927 229 76 3x
924 017 14 4x

430 515 47
über Falscheneingang geklebt
nur in Variante 02
for version 02

430 515 49
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 59
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 67
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 75
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 83
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 91
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 99
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 107
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 115
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 123
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 131
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 139
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 147
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

929 166 40
923 096 69
923 081 00
914 326 07
2,2 x 12 Plastile 45

915 417 89
M3 x 6 DIN7985
B3 DIN127 (unter Schraubenkopf)
(under screw head)
A3,2 DIN6798
(zwischen Kabelschuh und Druckbalken)
(between cable lug and press bar)

915 415 34 3-pol.
915 414 82
388 028 14
2x M3 x 6 DIN7985

In Variante 02/12/13
for version 02/12/13
927 229 01
In Variante 01/03_07/10/11/14_19
for version 01/03_07/10/11/14_19
504 657 24

432 522 06
M 2,1
A4,3 DIN6798
M4 x 10 DIN7985
B4 DIN127
916 118 00

an beiden Enden montiert
clipped to both ends

In Variante 04/05/07/11/13/15/17/19
for version 04/05/07/11/13/15/17/19

303 442 91
923 096 84
388 031 25
383 273 27
923 096 65
927 228 48

388 032 13
388-028-12
7x M3 x 6 DIN7985

432 521 91
2x ZyLSH ISO 2338-B-4x24
805 124 16

2x 432 521 92
2x 421 032 54

927 229 55
In Batteriefach geklebt
glued into battery compartment

Beschriftung siehe BeschBlatt
nach Variante
mit Etikettendrucker bedruckt

829 074 08 2x 70mm
In Batteriefach geklebt
glued into battery compartment
ca. 190mm Abstand, mittig

433 671 66
404 029 00
in Lpt. Abschirmung gekittet
riveted to pcb shielding

927 229 76 3x
924 017 14 4x

430 515 47
über Falscheneingang geklebt
nur in Variante 02
for version 02

430 515 49
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 59
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 67
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 75
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 83
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 91
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 99
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 107
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 115
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 123
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 131
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 139
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 147
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

929 166 40
923 096 69
923 081 00
914 326 07
2,2 x 12 Plastile 45

915 417 89
M3 x 6 DIN7985
B3 DIN127 (unter Schraubenkopf)
(under screw head)
A3,2 DIN6798
(zwischen Kabelschuh und Druckbalken)
(between cable lug and press bar)

915 415 34 3-pol.
915 414 82
388 028 14
2x M3 x 6 DIN7985

In Variante 02/12/13
for version 02/12/13
927 229 01
In Variante 01/03_07/10/11/14_19
for version 01/03_07/10/11/14_19
504 657 24

432 522 06
M 2,1
A4,3 DIN6798
M4 x 10 DIN7985
B4 DIN127
916 118 00

an beiden Enden montiert
clipped to both ends

In Variante 04/05/07/11/13/15/17/19
for version 04/05/07/11/13/15/17/19

303 442 91
923 096 84
388 031 25
383 273 27
923 096 65
927 228 48

388 032 13
388-028-12
7x M3 x 6 DIN7985

432 521 91
2x ZyLSH ISO 2338-B-4x24
805 124 16

2x 432 521 92
2x 421 032 54

927 229 55
In Batteriefach geklebt
glued into battery compartment

Beschriftung siehe BeschBlatt
nach Variante
mit Etikettendrucker bedruckt

829 074 08 2x 70mm
In Batteriefach geklebt
glued into battery compartment
ca. 190mm Abstand, mittig

433 671 66
404 029 00
in Lpt. Abschirmung gekittet
riveted to pcb shielding

927 229 76 3x
924 017 14 4x

430 515 47
über Falscheneingang geklebt
nur in Variante 02
for version 02

430 515 49
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 59
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 67
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 75
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 83
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 91
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 99
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 107
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 115
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 123
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 131
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 139
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 147
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

929 166 40
923 096 69
923 081 00
914 326 07
2,2 x 12 Plastile 45

915 417 89
M3 x 6 DIN7985
B3 DIN127 (unter Schraubenkopf)
(under screw head)
A3,2 DIN6798
(zwischen Kabelschuh und Druckbalken)
(between cable lug and press bar)

915 415 34 3-pol.
915 414 82
388 028 14
2x M3 x 6 DIN7985

In Variante 02/12/13
for version 02/12/13
927 229 01
In Variante 01/03_07/10/11/14_19
for version 01/03_07/10/11/14_19
504 657 24

432 522 06
M 2,1
A4,3 DIN6798
M4 x 10 DIN7985
B4 DIN127
916 118 00

an beiden Enden montiert
clipped to both ends

In Variante 04/05/07/11/13/15/17/19
for version 04/05/07/11/13/15/17/19

303 442 91
923 096 84
388 031 25
383 273 27
923 096 65
927 228 48

388 032 13
388-028-12
7x M3 x 6 DIN7985

432 521 91
2x ZyLSH ISO 2338-B-4x24
805 124 16

2x 432 521 92
2x 421 032 54

927 229 55
In Batteriefach geklebt
glued into battery compartment

Beschriftung siehe BeschBlatt
nach Variante
mit Etikettendrucker bedruckt

829 074 08 2x 70mm
In Batteriefach geklebt
glued into battery compartment
ca. 190mm Abstand, mittig

433 671 66
404 029 00
in Lpt. Abschirmung gekittet
riveted to pcb shielding

927 229 76 3x
924 017 14 4x

430 515 47
über Falscheneingang geklebt
nur in Variante 02
for version 02

430 515 49
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 59
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 67
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 75
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 83
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 91
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 99
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 107
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 115
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 123
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 131
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 139
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

430 515 147
In Gehäusemitte aufgeklebt
In Variante 02/12/13
for version 02/12/13

929 166 40
923 096 69
923 081 00
914 326 07
2,2 x 12 Plastile 45

915 417 89
M3 x 6 DIN7985
B3 DIN127 (unter Schraubenkopf)
(under screw head)
A3,2 DIN6798
(zwischen Kabelschuh und Druckbalken)
(between cable lug and press bar)


915 415 34 3-pol.
915 414 82
388 028 14
2x M3 x 6 DIN7985

In Variante 02/12/13
for version 02/12/13
927 229 01
In Variante 01/03_07/10/11/14_19
for version 01/03_07/10/11/14_19
504 657 24

432 522 06
M 2,1
A4,3 DIN6798
M4 x 10 DIN7985
B4 DIN127
916 118 00

an beiden Enden montiert
clipped to both ends

In Variante 04/05/07/11/13/15/17/19
for version 04/05/07/11/13/15/17/19









Legende Sicherungen
Fuses and ratings

| | | |
|------|------------|------------|
| 230V | 2 x T0,25A | 912 084 04 |
| 115V | 2 x T0,5A | 912 084 14 |

Option Graphics Display CS
235 066 01 oder/or 235 066 02

(L) Display 930 117 17



| | | | | | | |
|---------------------|-------------------|----------|------------|---|-----------------|------------|
| Werkstoff | f | - | 202295 | neue Varianten 06.07.19 | 27.03.97 | GKaltentb. |
| | e | - | 202294 | neue Varianten 10.13 | 26.06.96 | GKaltentb. |
| | d | - | 202308 | Massenverbindung auf G/2 geändert. | 02.05.96 | GKaltentb. |
| | c | - | 202087 | Hinweis Drahtsatz 383 273 27 nur in 04.05 | 06.12.95 | GKaltentb. |
| | b | - | 202019 | neue Varianten 04.05 | 10.10.95 | GKaltentb. |
| | a | - | 207782 | neue Variante 03 | 11.05.95 | GKaltentb. |
| Änd.-komf. Index | Änd.Nr. | | | Änderungstext | Datum | Name |
| Maßstab | Entwerfen | 25.10.94 | GKaltentb. | Sach-Nr. | 01...07 | |
| Normierung | Geprüft | 25.10.94 | NSchnurr | | 101 116 10...19 | S Bl.1 |
| Palmaff. Oberfläche | Freimachfreieranz | | | Format | Benennung | |
| ACAD-Nr. | 01847S11 | | | A1 | CardioSmart | |

▲

Diese Tabelle definiert die gültigen Konfigurationen des CardioSmart Version 1.0
This table defines the valid configurations of the CardioSmart Version 1.0

in Variante
in Device
Variation

| | | | | | | | |
|--------------|-------------------|-----------------|------------|---|---|------------|--------|
| Steuerung CS | Lpl. Steuerung CS | PCB Control CS | 388 028 12 | I | A | E , F1 , H | 01..03 |
| | Ersatzteil-Nr. | Spare Part No. | - | | | | |
| | Austausch Nr. | Replacement No. | 389 003 99 | | | | |

| | | | | | | |
|-------------------------|------------------------------|------------------------|------------|---|---|--------|
| Stromver- sorgung CS | Lpl. Stromvers- orgung CS | PCB Power Supply CS | 388 028 14 | B | D | 01..03 |
| | Ersatzteil-Nr. | Spare Part No. | - | | | |
| | Austausch Nr. | Replacement No. | - | | | |

| | | | | | | |
|-------------------|------------------------|-------------------------|------------|---|---|---|
| Displays: | | Cardio Vision Graphics | 235 066 01 | | | |
| Grafik Display CS | Lpl. Grafik Display CS | PCB Graphics Display CS | 388 028 11 | D | C | C |
| | Ersatzteil-Nr. | Spare Part No. | - | | | |
| | Austausch Nr. | Replacement No. | - | | | |

| | | |
|-----------|--------------------|------------|
| oder / or | Cardio Vision Text | 235 067 01 |
| | LCD Modul-CS | 303 442 87 |

Master Record Index
CardioSmart Version 1.0

| | | | | | | | | |
|---|--|--|--|--|-------------------------------|-----------------------|--|---------------------------------------|
| A | 1 | | 2 | | 3 | | 4 | |
| | Diese Tabelle definiert die gültigen Konfigurationen des CardioSmart Version 1.1 This table defines the valid configurations of the CardioSmart Version 1.1 | | | | | | | |
| | Komponente Component | | Sach-Nr. Part No. | | Index (Prod.) | Referenz Reference | kompatibel(Service) Compatible(Service) | in Variante in Device Variation |
| | Steuerung CS | | Lpl. Steuerung CS Ersatzteil-Nr. Austausch Nr. | PCB Control CS Spare Part No. Replacement No. | 388 028 12 - 389 004 05 | J A | H , I | 01..05 |
| B | oder / or | | | | | | | |
| | Steuerung CS_S | | Lpl. Steuerung CS_S Ersatzteil-Nr. Austausch Nr. | PCB Control CS_S Spare Part No. Replacement No. | 388 032 13 - 389 004 16 | B A | | 01..05 |
| C | Stromver-sorgung CS | | Lpl. Stromvers-sorgung CS Ersatzteil-Nr. Austausch Nr. | PCB Power Supply CS Spare Part No. Replacement No. | 388 028 14 - - | C D | | 01..05 |
| | Displays: | | Cardio Vision Graphics | | 235 066 01 | | | |
| D | Grafik Display CS | | Lpl. Grafik Display CS Ersatzteil-Nr. Austausch Nr. | PCB Graphics Display CS Spare Part No. Replacement No. | 388 028 11 - - | D C | C | |
| | oder / or | | Cardio Vision Text | | 235 067 01 | | | |
| E | | | LCD Modul-CS | | 303 442 87 | | | |
| | Pumpenmodul: | | | | 303 442 91 | | 04,05 | |
| | | | Lpl. Pumpe | PCB Pumpe | 388 029 18 | B | | 04,05 |
| | | | Ersatzteil-Nr. Austausch Nr. | Spare Part No. Replacement No. | - - | | | |

| | | | | | | | | | | |
|------------|-------------------|--|------------------|---|-----------|--|---------------|--|----------|-----------|
| | | Werkstoff | | | | | | | | |
| | | | e | -- | 202498 | Index 388 028 12 von "I" nach "J" | | | 03.09.96 | GKaltenb. |
| | | | d | -- | 202354 | Lpl Steuerung neue Sach-Nr. 388 032 13 | | | 25.06.96 | GKaltenb. |
| | | | c | -- | 202019 | Varianten 04.05 hinzu. | | | 09.10.95 | GKaltenb. |
| | | | b | -- | 201999 | Index v. Lpl Steuerung auf "I" | | | 09.08.95 | GKaltenb. |
| | | Rohteil Entw. | a | -- | 201901 | Index 388 028 14 von "B" nach "C". | | | 09.08.95 | GKaltenb. |
| | | | N.Schnurr | Änd.-Index | kommt vor | Änd.Nr. | Änderungstext | | | Datum |
| | | HELLIGE GMBH Munzinger Str. 3 D-79111 Freiburg | Maßstab — | | Datum | Name | | Sach-Nr. 101 116 01...05 S Bl.3 | | |
| | | | | Entworfen | 11.05.95 | GKaltenbach | | | | |
| | | | | Geprüft | 11.05.95 | N.Schnurr | | | | |
| | | | | Normgepr. | | | | | | |
| Paßmaß | Abmaße | | | | | | | | | |
| Oberfläche | Freimaßtoleranz | | Format A4 | Benennung Master Record Index CardioSmart Version 1.1 | | | | | | |
| | ACAD-Nr. 01847S3e | | | | | | | | | |

| | | | | | | | |
|---|--|---|---|-------------------------------|-----------------------|--|---------------------------------------|
| | 1 | | 2 | | 3 | | 4 |
| A | <div>Diese Tabelle definiert die gültigen Konfigurationen des CardioSmart Version 1.2 This table defines the valid configurations of the CardioSmart Version 1.2</div> | | | | | | |
| | Komponente Component | | Sach-Nr. Part No. | Index (Prod.) | Referenz Reference | kompatibel(Service) Compatible(Service) | in Variante in Device Variation |
| | Steuerung CS_S | Lpl. Steuerung CS_S Ersatzteil-Nr. Austausch Nr. | PCB Control CS_S Spare Part No. Replacement No. | 388 032 13 - 389 004 16 | C | A B | 01...05,10..13 |
| | | | | | | | |
| B | Stromver- sorgung CS | Lpl. Stromvers- orgung CS Ersatzteil-Nr. Austausch Nr. | PCB Power Supply CS Spare Part No. Replacement No. | 388 028 14 - - | C | D | 01...05,10..13 |
| C | Displays: | | | | | | |
| | | Cardio Vision Graphics | | 235 066 01 | | | 01...05,10..13 |
| | Grafik Display CS | Lpl. Grafik Display CS Ersatzteil-Nr. Austausch Nr. | PCB Graphics Display CS Spare Part No. Replacement No. | 388 028 11 - - | D | C C | 01...05,10..13 |
| | | | | | | | |
| D | oder / or | | | | | | |
| | | Cardio Vision Text | | 235 067 01 | | | nur/only 01...05 |
| | | LCD Modul-CS | | 303 442 87 | | | nur/only 01...05 |
| | | | | | | | |
| | Pumpenmodul: | | | | | | |
| | | | | 303 442 91 | | | 04,05,11,13 |
| | | Lpl. Pumpe | PCB Pumpe | 388 029 18 | B | | 04,05,11,13 |
| | | Ersatzteil-Nr. Austausch Nr. | Spare Part No. Replacement No. | - - | | | |
| E | | | | | | | |

| | | | | | | | | | | |
|------------|--------|--|------------------|--|----------|-----------------------------------|-------------------------|------------------|----------|--------------|
| | | Werkstoff | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | Rohteil Entw. | | | | | | | | |
| | | | a | - | 202713 | 388 032 13 v. Index "B" auf "C" . | | | 04.03.97 | G.Kaltenbach |
| | | N.Schnurr | Änd.-Index | kommt vor | Änd.Nr. | Änderungstext | | | Datum | Name |
| | | | | | | | | | | |
| | | HELLIGE GMBH Munzinger Str. 3 D-79111 Freiburg | Maßstab — | | Datum | Name | Sach-Nr. 101 116 | 01..05 10..13 | S | BL.4 |
| | | | | Entworfen | 26.06.96 | G.Kaltenbach | | | | |
| | | | | Geprüft | 26.06.96 | N.Schnurr | | | | |
| | | | | Normgepr. | | | | | | |
| Paßmaß | Abmaße | | | | | | | | | |
| Oberfläche | | Freimaßtoleranz | Format | Benennung | | | | | | |
| | | ACAD-Nr. 01847S4a | A4 | Master Record Index CardioSmart Version 1.2 | | | | | | |

| | | | | | | | |
|---|--|--|---|-------------------------------|-------------------------|--|---|
| | 1 | 2 | 3 | 4 | | | |
| A | Diese Tabelle definiert die gültigen Konfigurationen des CardioSmart Version 1.3 This table defines the valid configurations of the CardioSmart Version 1.3 | | | | | | |
| | Komponente Component | | Sach-Nr. Part No. | Index (Prod.) | Referenz Reference | kompatibel(Service) Compatible(Service) | in Variante in Device Variation |
| | Steuerung CS_S | Lpl. Steuerung CS_S Ersatzteil-Nr. Austausch Nr. | PCB Control CS_S Spare Part No. Replacement No. | 388 032 13 - 389 004 16 | D | A | B , C 01...07,10..19 |
| B | Stromver- sorgung CS | Lpl. Stromvers- sorgung CS Ersatzteil-Nr. Austausch Nr. | PCB Power Supply CS Spare Part No. Replacement No. | 388 028 14 - - | C | D | 01...07,10..19 |
| C | Graphics Displays: Cardio Vision Graphics oder/or Cardio Vision Graphics | | | | | | 235 066 01 (US) 235 066 02 01,03...07,10,11,14...17 02,12,13,18,19 |
| | Grafik Display CS | Lpl. Grafik Display CS Ersatzteil-Nr. Austausch Nr. | PCB Graphics Display CS Spare Part No. Replacement No. | 388 028 11 - - | D | C | C 01...07,10..19 |
| D | oder/or Text Displays : Cardio Vision Text oder/or Cardio Vision Text | | | | | | 235 067 01 (US) 235 067 02 01,04,06,07 02 |
| | LCD Modul-CS | | | 303 442 87 | | 01,02,04,06,07 | |
| E | Pumpenmodul: | | | | | | 303 442 91 04,05,07,11,13,15,17,19 |
| | Lpl. Pumpe | PCB Pumpe | 388 029 18 | B | 04,05,07,11,13,15,17,19 | | |
| | Ersatzteil-Nr. Austausch Nr. | Spare Part No. Replacement No. | - - | | | | |

| | | | | | | | | |
|------------|------------------------|-------------------|----------------|--------------|--|---------------------------|--------------|----------|
| | Werkstoff | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | Rohteil Entw. | | | | | | | |
| | | | a | - | 202936 | neuer Index Lpl Steuerung | | 17.06.97 |
| | N.Schnurr | | Änd.- Index | kommt vor | Änd.Nr. | Änderungstext | | Datum |
| | | | | | | | | Name |
| | Marquette Hellige GmbH | | Maßstab | | Datum | Name | Sach-Nr. | |
| | Munzinger Str. 3 | | - | | Entworfen | 01.04.97 | G.Kaltenbach | 01..07 |
| | D-79111 Freiburg | | | | Geprüft | 01.04.97 | N.Schnurr | 10..19 |
| | | | | | Normgepr. | | | S Bl.5 |
| Paßmaß | Abmaße | | | | | | | |
| | | | | | | | | |
| Oberfläche | | Freimaßtoleranz | | Format | Benennung | | | |
| | | ACAD-Nr. 01847S5a | | A4 | Master Record Index CardioSmart Version 1.3 | | | |

Diese Tabelle definiert die gültigen Konfigurationen des CardioSmart Version 1.3
This table defines the valid configurations of the CardioSmart Version 1.3

Referenztabelle der Geräte Versionen zu den zugehörigen CardioProm Software-Versionen

| Geräte-Versionen | Software-Versionen | | | | | |
|------------------|--------------------|-------|------|-------|------|-------|
| | V3.0...3.0c | V3.0d | V4.1 | V4.11 | V4.2 | V4.21 |
| V1.0 | X | | | | | |
| V1.1 | | X | X | X | | |
| V1.2 | | X | | X | X | |
| V1.3 | | | | | | X |

X = vorgesehene Verwendungen,vollkomptibel,alle Funktionen entsprechend Gebrauchsanweisung verfügbar.


Generell können alle Software-Versionen in allen Geräte-Versionen verwendet werden,jedoch stehen dann bestimmte Funktionen nicht,oder nur eingeschränkt zur Verfügung.Neue Features der Geräte Versionen werden nur mit den korrespondierenden Software Versionen unterstützt.

Verlauf der Geräte- und Software-Versionen


| Geräte-Version | Software-Version | Neue Geräte Funktionen | Neue Software Funktionen | Datum |
|----------------|------------------|---|---|-------|
| V1.0 | V3.0 | | CS-MI | 12/94 |
| | V3.0a | | CS, CS-E, CS-I | 12/94 |
| | V3.0b | | CS-M, CS-ME | 12/94 |
| | V3.0c | | italienisch spanisch | 02/95 |
| | V3.0d | | Fernstart russisch | 05/95 |
| V1.1 | | Fernstart Pumpe | | 09/95 |
| | V4.1 | | 12SL-Programm Modem Speicherung 45EKGs CSI-Protokoll amerikanische Variante | 08/96 |
| V1.2 | V4.2 | CardioSmart ST mit Ergometrie-Tastatur | Ergometrie | 12/96 |
| | V4.11 | | Override Funktion | 02/97 |
| V1.3 | V4.21 | Gebietsvarianten | schwedisch | 04/97 |
| | | | | |


Schutzvermerk nach DIN 34 beachten

| | | | | | | | | | |
|------------|-------------------|------------------------|----------------|--|---------|---------------|----------|--------------------------|--------|
| | | Werkstoff | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | Rohteil Entw. | | | | | | | |
| | | N.Schnurr | Änd.- Index | kommt vor | Änd.Nr. | Änderungstext | | Datum | Name |
| | | Marquette Hellige GmbH | Maßstab | - | Datum | Name | Sach-Nr. | 01..07 101 116 10..19 | S Bl.6 |
| | | Munzinger Str. 3 | | | | | | | |
| | | D-79111 Freiburg | | | | | | | |
| Paßmaß | Abmaße | | | | | | | | |
| Oberfläche | Freimaßtoleranz | | Format | Benennung | | | | | |
| | ACAD-Nr. 01847S6_ | | | | | | | | |
| | | | | Master Record Index CardioSmart Version 1.3 | | | | | |





| | | | | | | | | | |
|--------------|-------|--------------|--|----------------------------------|--------------|---------------------|--|--|--|
| HELLIGE GMBH | | | | Postfach 728 D-79007 Freiburg | | 388 028 12 P BL 1/8 | | | |
| ÄNDERUNGEN | | | | A03 | F67139 | | | | |
| Antrag vom | Index | Datum/Name | | | Datum/Name | | | | |
| ----- | A | 20.04.94/MSG | | | 14.12.93/NTS | | | | |
| 201357 | B | 24.08.94/NTS | | gez. | 14.12.93/NTS | | | | |
| | | | | gepr. | 14.12.93/NTS | | | | |
| | | | | Entw. | SCHNURR | | | | |


\$USER_CAE_D/cardiosmart/controller/ecg_core




| | | | | | | | | | |
|--------------|-------|--------------|--------|----------------------------------|--|-------------------------------------|--|--|--|
| HELLIGE GMBH | | | | Postfach 728 D-79007 Freiburg | | 388 028 12 P BL 3/8 | | | |
| ÄNDERUNGEN | | | A03 | F67139 | | LPL. STEUERUNG CS PCB CONTROL CS | | | |
| Antrag vom | Index | Datum/Name | | Datum/Name | | | | | |
| ----- | A | 20.04.94/NTS | | | | | | | |
| 201357 | B | 24.08.94/NTS | gez. | 30.03.94/HTW | | | | | |
| 201449 | C | 10.11.94/JOK | gespr. | 14.12.93/HTW | | | | | |
| | | | Entw. | WEBER | | | | | |









HELLIGE GMBH

Postfach 728
D-79007 Freiburg

388 028 12 P BL. 8/8


| ÄNDERUNGEN | | | A03 | F67139 |
|------------|-------|--------------|-------|--------------|
| Antrag vom | Index | Datum/Name | gez. | Datum/Name |
| 201357 | A | 28.04.94/NTS | gez. | 14.12.93/NTS |
| 201512 | B | 24.08.94/NTS | gepr. | 14.12.93/NTS |
| | C | 13.12.94/JOK | | |
| | | | Entw. | SCHNURR |

LPL. STEUERUNG CS
PCB CONTROL CS

MOTOR CONTROL




| | | | | | |
|--------------|-----|----------------------------------|-------|---------------------|--|
| HELLIGE GMBH | | Postfach 728 D-79007 Freiburg | | 388 028 12 R BL 1/2 | |
| ÄNDERUNGEN | | | | A02 F67139 | |
| ANTR. V | IND | DAT. | NAME | | |
| 201761 | C | 13.08.95 | GE2 | 18.05.94/JOK | |
| 201229 | C | 13.08.94/JOK | GE2 | 18.05.94/JOK | |
| 201601 | C | 14.12.94/JOK | GEPR. | 18.05.94/JOK | |
| 201579 | C | 10.02.95 | GEPR. | | |
| 201989 | E | 13.10.95/MNS | ENTR. | SCHNURR | |




* = NICHT BESTUECKT
* = NOT EQUIPPED


| | | | | | |
|--------------|-----|----------------------------------|--|----------------------|--|
| HELLIGE CMBH | | Postfach 728 D-79007 Freiburg | | 388 028 12 R BL2/2 | |
| ÄNDERUNGEN | | | | A02 F67139 | |
| ANTR. V | IND | DAT. / NAME | | A02 F67139 | |
| 201810 | D | 14. 05. 95/JOK | | 1994 18. 05. 94/JOK | |
| 201780 | F | 15. 05. 95/JOK | | 18. 05. 94/JOK | |
| 201999 | F | 13. 10. 95/MNS | | GFPR. 18. 05. 94/JOK | |
| 202080 | G | 06. 05. 94/JOK | | GFPR. 18. 05. 94/JOK | |
| 202498 | H | 04. 09. 96/JOK | | ENTR. SCHNURR | |
| | | | | LPL. STEUERUNG CS | |
| | | | | PCB CONTROL CS | |





| | | | | | | | |
|--------------|-------|------------|--|----------------------------------|--------------|--|--|
| HELLIGE GMBH | | | | Postfach 728 D-79007 Freiburg | | 388 032 13 P BL 1/8 | |
| ÄNDERUNGEN | | | | A03 | F67139 | LPL . STEUERUNG CS_S PCB CONTROL CS_S | |
| Antrag vom | Index | Datum/Name | | | Datum/Name | | |
| | | | | gez. | 24.06.96/NTS | | |
| | | | | gepr. | 24.06.96/NTS | | |
| | | | | Entw. | SCHNURR | | |





[1993, 7, 30, 6]
\$USER_CAE_D/cardiosmart/controller/ecg_core


















| | | | | | |
|---|--|--|---|---------------------|--|
| Marquette Hellige GmbH D-79007 Freiburg | | | | 388 032 13 R BL 1/2 | |
| REVISIONS Revision-No. Index Date/Name 202 713 A 06.02.97/JOK | | | A02 F67139 Date/Name 26.06.96/JOK Approved 26.06.96/JOK Issued EKH/IRB | | |
| | | | | LPL.STEUERUNG CS_S | |
| | | | | PCB CONTROL CS_S | |

| Marquette Hellige GmbH D-79007 Freiburg | | | | 388 032 13 R BL2/2 | | | | | | | | | | | | | | | | | | | | | |
|---|--------------|--------------|-------------|--------------------|-----------|---------|---|--------------|---------|---|--------------|--|--|-----|--------|-----------|--|-------|--------------|----------|--------------|--------|--|--------------|--|
| REVISIONS <table border="1"> <thead> <tr> <th>Revision-No</th> <th>Index</th> <th>Date/Name</th> </tr> </thead> <tbody> <tr> <td>202 498</td> <td>A</td> <td>04.09.96/JOK</td> </tr> <tr> <td>202 713</td> <td>B</td> <td>06.02.97/JOK</td> </tr> </tbody> </table> | | | Revision-No | Index | Date/Name | 202 498 | A | 04.09.96/JOK | 202 713 | B | 06.02.97/JOK | <table border="1"> <thead> <tr> <th>A02</th> <th>F67139</th> </tr> </thead> <tbody> <tr> <td colspan="2">Date/Name</td> </tr> <tr> <td>Drawn</td> <td>26.06.96/JOK</td> </tr> <tr> <td>Approved</td> <td>26.06.96/JOK</td> </tr> <tr> <td colspan="2">Issued</td> </tr> <tr> <td colspan="2">CONTROL CS_S</td> </tr> </tbody> </table> | | A02 | F67139 | Date/Name | | Drawn | 26.06.96/JOK | Approved | 26.06.96/JOK | Issued | | CONTROL CS_S | |
| Revision-No | Index | Date/Name | | | | | | | | | | | | | | | | | | | | | | | |
| 202 498 | A | 04.09.96/JOK | | | | | | | | | | | | | | | | | | | | | | | |
| 202 713 | B | 06.02.97/JOK | | | | | | | | | | | | | | | | | | | | | | | |
| A02 | F67139 | | | | | | | | | | | | | | | | | | | | | | | | |
| Date/Name | | | | | | | | | | | | | | | | | | | | | | | | | |
| Drawn | 26.06.96/JOK | | | | | | | | | | | | | | | | | | | | | | | | |
| Approved | 26.06.96/JOK | | | | | | | | | | | | | | | | | | | | | | | | |
| Issued | | | | | | | | | | | | | | | | | | | | | | | | | |
| CONTROL CS_S | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | LPL.STEUERUNG CS_S | | | | | | | | | | | | | | | | | | | | | |
| | | | | PCB CONTROL CS_S | | | | | | | | | | | | | | | | | | | | | |

Schutzvermerk nach DIN 34 beachten.



433 671 41



Kodierung PIN 12

NICHT INS ENGLISCHE UEBERSETZTE HIN-
WEISE NUR FUER FERTIGUNG WICHTIG!
TERMS NOT TRANSLATED INTO ENGLISH
ONLY RELEVANT FOR PRODUCTION!

BT-AS = 82
BT-BS = 17
LOET = 259

| | | | | | | | | | |
|--------------|-----|--------------|--|----------------------------------|-------|--------------|--------------|---------------|--|
| HELLIGE GMBH | | | | Postfach 728 D-79007 Freiburg | | 388 028 14 R | | | |
| ÄNDERUNGEN | | | | | | M1.5:1 | | A03 F67141 | |
| ANTR. V | IND | DAT./NAME | | | | | DAT./NAME | | |
| 201357 | A | 30.03.94/bus | | | 1994 | | 12.01.94/JOK | | |
| 201449 | B | 17.08.04/JOK | | | GEZ. | | 12.01.94/JOK | | |
| 201449 | C | 07.11.94/JOK | | | GEPR. | | 12.01.94/JOK | | |
| 201901 | D | 17.07.95/JOK | | | GEPR. | | 12.01.94/JOK | | |
| | | | | | ENTW. | | M&F/SCHULER | | |



marquette

A GE Medical Systems Company

World Headquarters
GE Marquette Medical Systems
8200 West Tower Avenue
Milwaukee, WI 53223 • USA
Tel. +1 414 355 5000
800.558.5120 (US only)
Fax +1 414 355 3790

Europe Region
Marquette Hellige GmbH
A GE Medical Systems Company
Postfach 60 02 65
D-79032 Freiburg • Germany
Tel. +49 761 45 43 - 0
Fax +49 761 45 43 - 233

Asia Region
26/F, Catic Plaza
8 Causeway Road
Causeway Bay
Hong Kong
Tel. +852 2804 2320
Fax +852 2804 1776