BSM- 9510J BSM- 9510K MU- 950RJ MU- 950RK RY- 002PA AY- 900PA AY- 910PA AA- 900PA WS- 920PA EK- 900P

Life Scope M BEDSIDE MONITOR

BSM-9510

0634-001352

Contents

Conventions Used in this Manual and Instrument	i
Warnings, Cautions and Notes	i
Explanations of the Symbols in this Manual and Instrument	ii

General	1C.1
Introduction	1.1
General Information on Servicing	
Service Policy, Service Parts and Patient Safety Checks	1.4
Service Policy	1.4
Service Parts	1.4
Patient Safety Checks	1.5
Composition	1.6
Bedside Monitor Main Unit	1.6
Specifications	1.7
Bedside Monitor Main Unit, MU-950RJ/RK	1.7
Sound	1.7
Alarm	1.7
Display	1.7
Module Slots	1.7
Multi Parameter Module, AY-900PA	1.8
ECG	
Respiration (Transthoracic impedance pneumography)	1.9
SpO, and Pulse Wave (Arterial Plethysmographic Waveform)	1.9
Non-invasive Blood Pressure, NIBP	1.10
Multi Amplifiers	1.10
Invasive Blood Pressure, IBP	1.10
Temperature	1.11
Cardiac Output, CO	1.11
Respiration (Thermistor probe pneumography)	1.12
Inspired Oxygen Fractional Concentration, FiO	
Expired Carbon Dioxide Tension, CO ₂	
ECG/BP Output	1.13
Multi Parameter Module, AY-910PA	1.13
SpO, and Pulse Wave (Arterial Plethysmographic Waveform)	1.13
Recorder Module, WS-920PA	1.14
Power Requirement	1.14
Environment for All Units and Modules	1.14
Dimensions and Weight (Approximate)	1.14
Electromagnetic Compatibility	1.15
Safety Standard	1.15
- Panel Description	
Control Panel	1.16
Remote Control	1.17

Main Unit	1.18
Multi Parameter Module	1.21
Smart Module	
Recorder Module	
Connection Diagram	
Block Diagram	

Section 2	Troubleshooting	
	Troubleshooting Table	2.1
	How to Use the Troubleshooting Table	2.1
	Power-Related Problem	2.2
	Display Problems	2.2
	Sound Problem	2.2
	Key Operation Problems	2.3
	Recorder Problem	2.3
	Other Module-Related Problem	2.3

Section 3	Diagnostic Check	3C.1
	Introduction	
	Power On Self Check	
	Calling up the Diagnostic Check and System Setup Screen	
	MU Manual Check	
	Memory Check	
	Flash ROM (program) Check	
	Flash ROM (data) Check	
	SRAM Check	
	DRAM Check	
	Com Check	
	Network I/F Check	
	Serial I/F Check	
	JA I/F Check	
	Display Check	
	Frame Mem Check	
	Graphic Check	
	Waveform Check	
	Backlight Check	
	Key LED Check	
	Key Check	
	Remote Check	
	Alarm Indicator Check	
	Alarm Pole Check	
	Other Check	3.17
	Sound Check	3.18
	Power Check	
	Card I/F Check	
	Timer IC Check	

Section 4	Board/Unit Description	4C.1
	Signal Flow	4.1
	Vital Sign Signals from Patients	4.1
	Display Data	4.1
	Power Control Signal by Front Power Switch or Power Button on	
	Optional Remote Control	4.1
	MAIN Board UR-3485	4.2
	EXT JA Board UR-3489	4.2
	JA Motherboard UR-3486	4.2
	IR DETECT Board UR-3487	4.3
	LCD JUNC Board UR-3504	4.3
	LED Board UR-3393	4.3
	OPERATION Board UR-3506	4.3
	Power Supply Unit SC-036R	4.4

Section 5	Disassembly and Assembly	
	Before You Begin	5.1
	Warnings and Cautions	5.1
	Required Tools	5.1
	Replacing MAIN Board	5.2
	Replacing Fuse	5.5
	Replacing Power Supply Unit	5.7
	Replacing JA Motherboard	5.10
	Replacing IR DETECT Board	5.13
	Replacing LCD Unit	5.16
	Replacing DC-AC Inverter	
	Replacing Backlight Lamps	
	Replacing LCD Filter	5.26
	Replacing OPERATION Board	5.29
	Replacing Lithium Battery	

Section 6	Maintenance	
	Maintenance Check Items	6.1
	External	6.1
	Safety	
	Modules	
	Display	
	Measuring Parameters	
	Recorder	
	Backup	
	Others	6.4

Section 7	Adjustment	7C.1
-----------	------------	------

Section 8	Replaceable Parts List	. 8C.1
	Main Unit Parts	8.1

Section 9	Connector Pin Assignment	
	Input/Output Connector Pin Assignment	9.1
	Alarm Output Socket, ALARM	9.1
	General Serial Socket, SERIAL	9.1
	Network Socket, NETWORK	
	JA Output Socket	
	Connector on MAIN Board	
	Connector on EXT JA Board	
	Memory Card Connector	
	OPERATION Connector	
	LED Connector	
	LCD Connector	
	JA Motherboard Connector	
	Connectors on the Power Supply Unit	
	BDM Connector	
	DEBUG Connector	9.7

Conventions Used in this Manual and Instrument

Warnings, Cautions and Notes

Warnings, cautions and notes are used in this manual to alert or signal the reader to specific information.

WARNING

A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.

CAUTION

A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

NOTE

A note provides specific information, in the form of recommendations, prerequirements, alternative methods or supplemental information.

Explanations of the Symbols in this Manual and Instrument

The following symbols found in this manual/instrument bear the respective descriptions as given.

On panels

Symbol	Description	Symbol	Description
Ċ	Standby lamp	Ŕ	Type BF applied part
\odot	"On" only for a part of equipment	IPX4	Splash-proof equipment
Ò	"Off" only for a part of equipment	IPX7	Watertight equipment
\sim	Alternating current	~~	Year of manufacture
Å	Equipotential terminal	SN	Serial number
Â	Attention, consult operator's manual	⇔	Fuse
×	Alarm off/suspend	Â	High voltage
\ominus	Output terminal	Ð	Protective earth
ł₩ŀ	Defibrillation-proof type CF applied part		The CE mark is a protected
	Type CF applied part	0086	conformity mark of the European Community. The products herewith comply with the requirements of the Medical Device Directive 93/42/EEC.
۱ <u>۴</u> ۲	Defibrillation-proof type BF applied part		

On screen

Symbol	Description	Symbol	Description
	Open pulldown menu	×	Alarm off/suspend
	Close window button	溪2	Alarm off/suspend with remaining time
•	QRS sync mark	ARR	Arrhythmia analysis off

Section 1 General

Introduction	1.1
General Information on Servicing	1.2
Service Policy, Service Parts and Patient Safety Checks	1.4
Service Policy	1.4
Service Parts	1.4
Patient Safety Checks	1.5
Composition	1.6
Bedside Monitor Main Unit	1.6
Specifications	1.7
Bedside Monitor Main Unit, MU-950RJ/RK	1.7
Sound	1.7
Alarm	1.7
Display	1.7
Module Slots	1.7
Multi Parameter Module, AY-900PA	1.8
ECG	1.8
Respiration (Transthoracic impedance pneumography)	1.9
SpO_2 and Pulse Wave (Arterial Plethysmographic Waveform)	1.9
Non-invasive Blood Pressure, NIBP	1.10
Multi Amplifiers	1.10
Invasive Blood Pressure, IBP	1.10
Temperature	1.11
Cardiac Output, CO	1.11
Respiration (Thermistor probe pneumography)	1.12
Inspired Oxygen Fractional Concentration, FiO ₂	1.12
Expired Carbon Dioxide Tension, CO ₂	1.12
ECG/BP Output	1.13
Multi Parameter Module, AY-910PA	1.13
SpO_2 and Pulse Wave (Arterial Plethysmographic Waveform)	1.13
Recorder Module, WS-920PA	1.14
Power Requirement	1.14
Environment for All Units and Modules	1.14
Dimensions and Weight (Approximate)	1.14
Electromagnetic Compatibility	1.15

Safety Standard	
Panel Description	1.16
Control Panel	
Remote Control	
Main Unit	
Multi Parameter Module	
Smart Module	
Recorder Module	
Connection Diagram	
Block Diagram	

Introduction

This service manual provides useful information to qualified personnel to understand, troubleshoot, service, maintain and repair the BSM-9510J/K Bedside Monitor (referred to as "the instrument" in this service manual).

All replaceable parts or units of this instrument are clearly listed with exploded illustration to help you locate the parts quickly.

The "Maintenance" section in this service manual describes the maintenance that should be performed by qualified service personnel. The "Maintenance" section in the operator's manual describes the maintenance that can be performed by the user.

The information in the operator's manual is primarily for the user. However, it is important for service personnel to thoroughly read the operator's manual and service manual before starting to troubleshoot, service, maintain or repair this instrument. This is because service personnel need to understand the operation of the instrument in order to effectively use the information in the service manual.

General Information on Servicing

Note the following information when servicing the instrument.

Safety

CAUTIONS

• There is the possibility that the outside surface of the instrument, such as the operation keys, could be contaminated by contagious germs, so disinfect and clean the instrument before servicing it. When servicing the instrument, wear rubber gloves to protect yourself from infection.

• There is the possibility that when the lithium battery is broken, a solvent inside the lithium battery could flow out or a toxic substance inside it could come out. If the solvent or toxic substance touches your skin or gets into your eye or mouth, immediately wash it with a lot of water and see a physician.

Liquid ingress

The instrument is not drip-proof, so do not install the instrument where water or liquid can get into or fall on the instrument. If liquid accidentally gets into the instrument or the instrument accidentally drops into liquid, disassemble the instrument, clean it with clean water and dry it completely. After reassembling, verify that there is nothing wrong with the patient safety checks and function/performance checks. If there is something wrong with the instrument, contact your Nihon Kohden representative to repair.

Environmental Safeguards

Depending on the local laws in your community, it may be illegal to dispose of the lithium battery and CRT unit in the regular waste collection. Check with your local officials for proper disposal procedures.

Disinfection and cleaning

To disinfect the outside surface of the instrument, wipe it with a non-abrasive cloth moistened with any of the disinfectants listed below. Do not use any other disinfectants or ultraviolet rays to disinfect the instrument.

 Chlorohexidine gluconate solution: 	0.5%
- Benzethonium chloride solution:	0.2%
- Glutaraldehyde solution:	2.0%
- Benzalkonium chloride:	0.2%

- Hydrochloric alkyl diaminoethylglycine: 0.5%

Transport

Use the specified shipment container and packing material to transport the instrument. If necessary, double pack the instrument. Also, put the instrument into the shipment container after packing so that the buffer material does not get into the inside of the instrument.
When transporting the board or unit of the instrument, be sure to use a conductive bag. Never use an aluminum bag when transporting the power board, power unit or board on which a lithium battery is mounted. Also, never use a styrene foam or plastic bag which generates static electricity to wrap the board or unit of the instrument.

Handling the instrument

• Because the outside surface of the instrument is made of resin, the outside surface of the instrument is easily damaged. So when handling the instrument, remove clutter from around the instrument and be careful to not damage the instrument or get it dirty.

• Because most of the boards in the instrument are multilayer boards with surface mounted electrical devices (SMD), when removing and soldering the electrical devices, a special tool is required. To avoid damaging other electrical components, do not remove and solder SMD components yourself.

Measuring and Test Equipment

Maintain the accuracy of the measuring and test equipment by checking and calibrating it according to the check and calibration procedures.

Service Policy, Service Parts and Patient Safety Checks

Service Policy

Our technical service policy for this instrument is to replace the faulty unit, board or part or damaged mechanical part with a new one. Do not perform electrical device or component level repair of the multilayer board or unit. We do not support component level repair outside the factory for the following reasons:

- Most of the boards are multilayer boards with surface mounted electrical devices, so the mounting density of the board is too high.
- A special tool or high degree of repair skill is required to repair the multilayer boards with surface mounted electrical devices.

Disassemble the instrument or replace a board or unit in an environment where the instrument is protected against static electricity.

As background knowledge for repair, pay special attention to the following:

- You can reduce the repair time by considering the problem before starting repair.
- You can clarify the source of most of the troubles using the information from the diagnostic check function of the instrument. Refer to "Diagnostic Check " of this manual.

Service Parts

Refer to "Replaceable Parts List" of this manual for the service parts for technical service that we provide.

NOTE

When ordering parts or accessories from your Nihon Kohden representative, please quote the NK code number and part name which is listed in this service manual, and the name or model of the unit in which the required part is located. This will help us to promptly attend to your needs. Always use parts and accessories recommended or supplied by Nihon Kohden Corporation to assure maximum performance from your instrument.

Patient Safety Checks

Periodic maintenance procedures and diagnostic check procedures are provided in this manual to ensure that the instrument is operating in accordance with its design and production specifications. To verify that the instrument is working in a safe manner with regard to patient safety, patient safety checks should be performed on the instrument before it is first installed, periodically after installation, and after any repair is made on the instrument.

For patient safety checks, perform the following checks as described in the International Electrotechnical Commission's standard, IEC60601-1 (1988):

- Protective earth resistance check
- Earth leakage current check
- Enclosure leakage current check
- Patient leakage current check
- Withstanding voltage check

Composition

Bedside Monitor Main Unit



Specifications

Bedside Monitor Main Unit, MU-950RJ/RK

For other details, refer to the specifications of the respective module and unit.

Sound		
	Sound type:	Alarm, synchronization, click
	Alarm sound:	Volume variable
	Synchronization sound:	Volume variable, pitch variable for SpO_2 or BP
	Click sound:	Volume variable
Alarm		
	Alarm items:	Upper/lower limits alarm, apnea alarm, arrhythmia alarm, module alarm, external instrument alarm, electrode check alarm, faulty instrument alarm, connector disconnection alarm, operating environment alarm
	Alarm types:	Crisis (red, blinking), Warning (yellow, blinking), Advisory (yellow, lighting), Message, System Guidance
	Alarm indication:	Alarm indicator, highlighted numerical display, numerical display color, alarm sound, highlighted message for arrhythmia
	Alarm suspend:	Provided (For 1, 2, 3, 4, 5, 10 min)
	Alarm silence:	Provided (For 1, 2, 3, 4, 5, 10 min)
	Alarm setting:	Individual upper/lower limits setting for each parameter
Display		
	Sweep speed:	Respiration wave (respiration, CO_2): 6 or 25 mm/s
		Others: 25 mm/s
	Display waveforms:	ECG, IBP, ICP, SpO_2 , Respiration wave, CO_2 and other parameters
		depending on the module and unit
	Numerical data display:	Heart rate, VPC rate, arrhythmia message, ST level, IBP (systolic,
		diastolic, mean), NIBP (systolic, diastolic, mean), respiration rate, pulse rate, temperature, SpO ₂ , CO ₂ , O ₂ , and other parameters depending on the module and unit.
	Synchronization mark:	Heart rate sync mark, pulse rate sync mark
	Trendgraph display time:	30 min, 1, 2, 4, 6, 8, 12 or 24 h
	Display size:	10.4 inch
	Viewing area:	211×158 mm
	Resolution:	$640 \times 480 \text{ dots}$
	Max. No. of waveform trace	es: 6 traces
	Sweep time:	about 6 s (When sweep speed is 25 mm/s)
	Waveform display mode:	Moving mode

Module Slots

No. of	slots:
--------	--------

6

Multi Parameter Module, AY-900PA ECG

• Electrode offset potential tolerance:		$\pm 500 \text{ m}$	$\pm 500 \text{ mV}$		
• Input dynamic range:		$\pm 10mV$	$\pm 10 \text{ mV}$		
• Internal noise:		$\leq 20 \ \mu V$	\leq 20 µVp-p, referred to input		
• Input impedance:		$\geq 5 M\Omega$	(at 10 Hz)		
• Common Mode	e Rejection Ratio:	≥95 dB	(with a 51 l	$\Omega/47$ nF imbalance)	
• Input bias curre	ent:	$\leq 100 \text{ nA}$	4		
• Heart rate coun	t				
	Calculation method:	8-beat m beats av	noving aver reraged) (Se	age/Instantaneous beat-to-beat (consecutive 2 lectable)	
	Counting range:	15 to 30	0 beats/min	ı	
• Arrhythmia ana	llysis				
	Analysis method:	Multi-te	emplate mat	ching method	
	No. of channels:	2 chann	els		
	VPC counting range:	0 to 99	DI E V FIR	VF/VT V TACHY VPC RUN COUPI ET FARI V	
	r ming uning undry 515.	V, BIGE	V, BIGEMINY, FREQ VPC, TACHY, BRADY, PROLONGED RR		
• Arrhythmia rec	all				
	Number of recall files:	32			
	Storage time per file:	8 s			
• ST level measu	rement				
	No. of measurement channels:	3-electro	ode: 1 ch		
		6-electro	ode: max. 8	ch	
	Magguring range	10-elect	rode: max.	12 ch	
• Decometer pul	measuring range.			1002 competible	
Defibrillation a			ECG input protected against 400 L discharge		
			De chal		
• ESU interference	ce filter:	Provideo	Provided		
• AC hum filter:		Provide	d		
• Lead					
	3-electrode cable:	I, II, III			
	6-electrode cable:	I, II, III,	aVr, aVl, aV	VF, Va, Vb (Va and Vb: any 2 leads from chest	
	10 alactroda cabla:	leads)	aVr aVr aV	Va V1 to V6	
• FCG frequency	range	1, 11, 111,	а v к, а v L, а v	r, v 1 to v 0	
• ECO nequency	Signal handwidth:	DC to 9	0 Hz		
	Display bandwidth:	Filter	Drift free		
	Ţ	Off	Off	0.3 to 70 Hz	
		On	Off	0.3 to 20 Hz	
		Off	On	1.0 to 70 Hz	
		On	On	1.0 to 20 Hz	

• Waveform display

Display sensitivity:	10 mm/mV \pm 5% (at \times 1 sensitivity)
Sensitivity control:	$\times 1/4$, $\times 1/2$, $\times 1$, $\times 2$, $\times 4$ or AUTO
Drift rejection filter:	Available
Pacing spike display:	Available
Auto positioning:	Available
• Heart rate display update cycle:	Every 3 s or when alarm is generated

• Alarm

Upper limit range:	20 to 300 beats/min in 5 beats/min steps, OFF
Lower limit range:	OFF, 15 to 295 beats/min in 5 beats/min steps
Alarm items:	TACHY, BRADY

Respiration (Transthoracic impedance pneumography)

• Measuring imp	edance range:	$\geq 2 k\Omega$	
• Internal noise:		$\leq 0.1 \ \Omega$	
• Excitor current:		$30 \pm 10 \ \mu Arms$ at 40 kHz	
• Frequency resp	onse:	0.1 to 3 Hz (-3 dB)	
• Respiration cou	inter		
	Counting range:	0 to 150 breaths/min	
	Apnea setting range:	5 to 40 s (set on the main unit)	
• Defibrillation p	proof:	Respiration input protected against 400 J discharge	
• Waveform display			
	Display sensitivity:	$10 \text{ mm}/\Omega \text{ (at} \times 1 \text{ sensitivity)}$	
	Sensitivity control:	$\times 1/4, \times 1/2, \times 1, \times 2, \times 4$ or AUTO	
• Respiration rate display update cycle:		Every 3 s or when alarm is generated	
• Alarm			
	Upper limit range:	2 to 150 breaths/min in 2 breaths/min steps, OFF	
	Lower limit range:	OFF, 0 to 148 breaths/min in 2 breaths/min steps	
	Apnea time:	5 to 40 s in 5 s steps, OFF	

SpO₂ and Pulse Wave (Arterial Plethysmographic Waveform)

With Nihon Kohden probe

• Measuring range

	SpO ₂ :	1 to 100% SpO ₂ in 1% SpO ₂ steps
	Pulse rate:	30 to 300 beats/min
• SpO ₂ accuracy		
	80 to 100% SpO ₂ :	$\pm 2\%$ SpO ₂
	50 to 79% SpO ₂ :	$\pm 3\%$ SpO ₂
• SpO ₂ display		
	Pulse rate display update cycle:	Every 3 s or when alarm is generated
	Sync tone modulation:	Changes in 20 steps at 81 to 100% SpO_2
• Alarm		
	Upper limit range:	OFF, 51 to 100% SpO_2 in 1% SpO_2 steps
	Lower limit range:	OFF, 50 to 99% SpO_2 in 1% SpO_2 steps
Service Manual B	SM-9510	

NOII-IIIVASIVE DIOOU FIESSUIE, NIDI	Non-invasive	Blood	Pressure,	NIBF
-------------------------------------	--------------	-------	-----------	------

• Measuring method	1:	Oscillometric	
• Cuff pressure disp	lay range:	0 to 300 mmHg (0 to	o 40 kPa)
• Accuracy			
C	to 200 mmHg:	±3 mmHg	
2	200 to 300 mmHg:	±4 mmHg	
• Safety			
Ν	Maximum pressurization value cu	iff inflation limiter:	Adult 300 mmHg
(Suff inflation time limiter		Adult 150 s
			Neonate 80 s
• Measurement mod	e:	Manual (Single mea	asurement)
		Continuous (Succes	sive repetition or at 1 min interval for 15 min
		period)	
		(In ICU m	ode at 2, 2.5, 5, 10, 15, 30 min interval) ode at 5, 10, 15, 30 min, 1, 2, 4, 8 h interval)
• NIBP data display	update cycle:	Updated every measured	surement
• Measurement end	sound:	Generated when me	asurement ends (Set on the System Setup screen)
• Alarm (Systolic, D	iastolic, Mean)		
τ	Jpper limit range:	22 to 250 mmHg in	2 mmHg steps, OFF
Ι	Lower limit range:	OFF, 20 to 248 mmI	Hg in 2 mmHg steps
Multi Amplifiers			
	Measuring parameters:	IBP, Temp, CO, Res	p (thermistor), FiO,, CO, (main stream)
Ι	nput impedance:	1 MΩ	
Invasive Blood Pre	ssure, IBP		
• Measuring range:		-50 to 300 mmHg	
• Measuring accurate	су	±1% (≥100 mmHg)), ±1 mmHg (< 100 mmHg)
• Auto zero balancin	ng range:	±200 mmHg	
• Auto zero balancin	ng accuracy:	±1 mmHg	
• Pulse rate countin	g range:	30 to 300 beats/min	n
• Pulse rate countin	g accuracy:	±1 beat/min	
• Noise:		$\leq 0.25 \text{ mmHg}$	
• Temperature zero d	lrift:	±0.1 mmHg/°C	
• Frequency response	e:	DC to 10/20 Hz, set software)	on the main unit (Digital filter processing by
• Blood pressure dis	play range:	0 to 300 mmHg	
• Display update cy	cle:	Every 3 s or when a	larm is generated
• BP tone:		Provided, systolic v 5 mmHg	alue 20 to 120 mmHg, changes in 20 steps every

• Alarm (Syst	olic, Diastolic, Mean)	
	Upper limit range:	2 to 300 mmHg in 2 mmHg steps, OFF
	Lower limit range:	OFF, 0 to 298 mmHg in 2 mmHg steps
Temperature		
• Measuring r	ange:	0 to 45°C (32 to 113°F)
• Measuring a	accuracy	
	Sensor:	±0.1°C
	Module:	±0.1°C (25 to 45°C, 77 to 113°F)
		±0.2°C (Other range)
• Noise:		$\leq 0.014^{\circ}$ C at 37°C (99°F)
• Temperature	e drift:	±0.005°C/°C
• Temperature	e display	
	Display range:	0 to 45°C (32 to 113°F)
	Display update cycle:	Every 3 s
• Alarm		
	Upper limit range:	0.1 to 45°C (33 to 113°F) in 0.1°C (0.5°F) steps, OFF
	Lower limit range:	OFF, 0 to 44.9°C (32 to 112°F) in 0.1°C (0.5°F) steps
Cardiac Outpu	ut, CO	
• Measuring 1	method: Thermodilution method	
 Measuring r 	ange	

- Injectate temperature (Ti):0 to $27^{\circ}C$ (32 to $81^{\circ}F$)Blood temperature (Tb):15 to $45^{\circ}C$ (59 to $113^{\circ}F$)Thermodilution curve (Δ Tb):0 to $2.5^{\circ}C$ (32 to $36.5^{\circ}F$)Cardiac output (CO):0.5 to 20 L/min
- Measuring accuracy (catheter sensor accuracy is not included)

	Ti:	±0.1°C
	Tb 25 to 45°C:	± 0.1 °C
	Tb 15 to 25°C:	±0.2°C
	CO:	$\pm 5\%$
• Noise		
	Ti:	≤ 0.025°C
	Tb:	$\leq 0.016^{\circ}$ C at 37°C (99°F)
	$\Delta Tb:$	$\leq 0.005^{\circ}$ C at 37°C (99°F)
• Temperatu	re drift	
	Ti:	±0.005°C/°C
	Tb:	$\pm 0.005^{\circ}C/^{\circ}C$
• Frequency	response (ΔTb):	0 to 3 Hz (Digital filter processing)
• Catheter s	ize:	5F, 7F or 7.5F
• Injectate v	volume range:	3, 5, 10 cc

• Cardiac output display

CO value display update cycle: Every measurement Thermodilution curve display time: 45 s

- 1. GENERAL
 - Alarm limits

	CO:	None
	Blood temperature:	15 to 45°C in 0.1°C steps, OFF
Respiration (Th	ermistor probe pneumogra	phy)
• Respiration ra	te counting range:	0 to 150 resp/min
		Apnea, 5 to 40 s
• Accuracy:		± 1 resp/min
• Temperature n	neasuring range:	10 to 40°C (50 to 104°F)
• Maximum det	ection resistance range:	1 kΩ
• Recorder sens	itivity:	$\Delta 100 \ \Omega/400 \text{ digits } \pm 10\%$
		(400 digits is equivalent to 1 cm on paper at \times 1 recorder sensitivity)
• Noise:		$\leq 1.0 \Omega$ (referred to input)
• Frequency res	ponse:	0.1 to 3 Hz (Digital filter processing)
• Waveform disp	play	
	Display sensitivity:	$10 \text{ mm}/100 \Omega (at \times 1 \text{ sensitivity})$
	Sensitivity control:	$\times 1/4, \times 1/2, \times 1, \times 2, \times 4$ or AUTO
• Respiration ra	te display update cycle:	Every 3 s or when alarm is generated
• Alarm		
	Upper limit range:	2 to 150 breaths/min in 2 breaths/min steps, OFF
	Lower limit range:	OFF, 0 to 148 breaths/min in 2 breaths/min steps
	Apnea time:	5 to 40 s in 5 s steps, OFF
Inspired Oxyge	n Fractional Concentration,	FiO ₂
• Measuring ran	ge:	10 to 100% O_2 in 1% steps
• Amplifier accu	iracy:	±1% full scale
• Accuracy incl	uding sensor	
	21% O_2 calibration:	±3% full scale
	100% O_2 calibration:	±2% full scale
• Noise:		$\leq 0.12\% \text{ O}_2$
• Temperature d	rift:	±0.12% O ₂ /°C
• FiO ₂ display u	pdate cycle:	Every 3 s or when alarm is generated
• Alarm		
	Lower limit range:	OFF, 18 to 100% in 1% steps, OFF
Expired Carbon	Dioxide Tension, CO ₂	
• Measurement	method:	Main stream method
• Measuring ran	ge:	0 to 76 mmHg
• Warm-up time	:	Main stream method: none
• Response time	::	200 ms (typ.) for step from 10 to 90%

• Detectable respiration rate:

- Main stream: 3 to 60 breaths/min
- Measuring accuracy (When 1 atmospheric pressure, air inspiration, non condensation)

	0 to 40 mmHg:	±4 mmHg
	41 to 76 mmHg:	±10% reading
	O ₂ gas effects:	Approx10% reading (When 100% oxygen is inspired)
	CO_2 gas, N_2O anesthetic gas effects:	Effects when 1 mmHg CO_2 gas is inspired, approx. 10% reading
		Accuracy in using N_2O anesthetic gas is not guaranteed
• CO_2 value display	y update cycle:	Every 3 s or when alarm is generated
• Alarm		

EtCO ₂ upper limit range:	2 to 99 mmHg in 1 mmHg steps, OFF
EtCO ₂ lower limit range:	OFF, 1 to 98 mmHg in 1 mmHg steps, OFF
Apnea time:	5 to 40 s in 5 s steps, OFF

ECG/BP Output

• Output impedance

ECG:	100 Ω
BP:	100 Ω

• Output waveform

ECG:	$\geq \pm 4.096 \text{ V} (\text{at } 1 \text{ mV/V})$
BP:	-0.64 to +4.48 V (at 100 mmHg/V)
HT:	Open collector output (Maximum sink current 3 mA, pulse width 100
	ms)

• Frequency response

•

	ECG:	DC to \geq 35 Hz
	BP:	DC to ≥15 Hz
Sensitivity	y accuracy	
	ECG:	±5%
	BP:	$\pm 1\%$
Delay		
	ECG:	≤ 35 ms
	BP:	\leq 35 ms
	HT:	≤ 100 ms

Multi Parameter Module, AY-910PA

Same as AY-900PA Multi Parameter Module, except for SpO₂.

SpO₂ and Pulse Wave (Arterial Plethysmographic Waveform)

With Nellcor probe

• Measuring range

SpO ₂ : 1	to 100% SpO_2 in 1% SpO_2 steps
Pulse rate: 2	20 to 250 beats/min

• SpO₂ accuracy

70 to 100% SpO₂:

 $\pm 2\%$ SpO₂

- 1. GENERAL
 - SpO_2 display

ery 3 s or when alarm is generated
anges in 20 steps at 81 to 100%
F, 51 to 100% in 1% steps
F, 50 to 99% in 1% steps

Recorder Module, WS-920PA

• Recorder

Recording method:	Thermal array recording
Number of channels:	2 channels (maximum)
Maximum recording width:	\geq 45 mm
Paper speed:	5, 25 or 50 mm/s, $\pm 2\%$
Recording paper:	FQW50-3-100

Power Requirement

Line voltage:	MU-950RJ: AC 100 to 127 V ±10%
	MU-950RK: AC 220 to 240 V $\pm 10\%$
Line frequency:	50 or 60 Hz
Power consumption:	120 VA maximum, including the color display unit

Environment for All Units and Modules

• Operating environment

Temperature:		10 to 40°C (50 to 104°F) excluding recording paper
	Humidity:	30 to 90% RH (0 to 40°C, 32 to 104°F non-condensing)
	Atmospheric pressure:	70 to 106 kPa
• Storage environment		
	Temperature:	-20 to 65°C (-4 to 149°F) excluding recording paper
		Recording paper: -15 to 55°C (5 to 131°F)
	Humidity:	15 to 90% RH (non-condensing)
	Atmospheric pressure:	70 to 106 kPa

Dimensions and Weight (Approximate)

MU-950RJ/RK Main Unit:	345 W \times 300 H \times 205 D, 6.5 kg
AY-900PA/910PA Multi Parameter Module:	
	113 W \times 117 H \times 160 D, 1.5 kg
AA-900PA Smart Module:	37 W \times 117 H \times 160 D, 0.4 kg
RY-002PA Remote Control:	$50~\mathrm{W}\times22~\mathrm{H}\times162~\mathrm{D},$ 0.09kg
WS-920PA Recorder Module:	113 W \times 117 H \times 160 D, 0.9 kg
EK-900P Blank Module:	$37~\mathrm{W} \times 117~\mathrm{H} \times 160~\mathrm{D}, 0.17~\mathrm{kg}$

Electromagnetic Compatibility

IEC 60601-1-2 (1993)

Safety Standard

Safety Standard: IEC 60601-1 (1988)		IEC 60601-1 (1988)
		IEC 60601-1 (1991) Amendment 1
		IEC 60601-1 (1995) Amendment 2
		IEC 60601-1-1 (1992)
		IEC 60601-1-1 (1995) Amendment 1
		IEC 60601-2-27 (1994) Particular requirements for the safety of
		electrocardiographic monitoring
		IEC 60601-2-30 (1995) Particular requirements for the safety of
		automatic cycling in direct blood pressure monitoring equipment
		IEC 60601-2-34 (1994) Particular requirements for the safety of
		direct blood pressure monitoring equipment
	According to the type of protection against electrical shock:	
		CLASSIEQUIPMENT
	According to the degree of protection against electrical shock:	
	AY-900PA/910PA:	Defibrillator-proof type CF applied part: ECG, Resp (impedance),
		SpO ₂ , NIBP, IBP, Temp, Resp (thermistor), FiO ₂
		CF applied part: CO
		BF applied part: CO_2 (main stream)
	AA-900PA:	Defibrillator-proof type CF applied part: IBP, Temp, Resp
		(thermistor), FiO ₂
		CF applied part: CO
		BF applied part: CO_2 (main stream)
	According to the degree of protection against harmful ingress of water:	
		IPX0 (Ordinary EQUIPMENT)
	According to the degree of safety of application in the presence of a FLAMMABLE ANAESTHETIC	
	MIXTURE WITH AIR, OR OXYGEN OR NITROUS OXIDE:	
		EQUIPMENT not suitable for use in the presence of FLAMMABLE
		ANAESTHETIC MIXTURE WITH AIR, OR WITH OXYGEN OR
		NITROUS OXIDE
	According to the mode of operation:	

CONTINUOUS OPERATION

Panel Description

Control Panel



Remote Control



Attach a string or strap to hang the remote control.



Main Unit

Front Panel



Left Side Panel



Insert EK-900P blank modules in empty slots to prevent dust from accumulating on the connectors.

Service Manual BSM-9510

Rear Panel



Connects to the optional alarm pole.

Cover the NETWORK socket when not in use.

Multi Parameter Module

AY-900PA is for Nihon Kohden SpO_2 probes. AY-910PA is for Nellcor SpO_2 probes.

CAUTION

When pressing the keys on the module, use your finger. Do not press the keys with a sharp object. Otherwise keys may be broken.



No.	Name	Description
1	ECG/RESP key	Press to open the ECG window.
2	SpO ₂ key	Press to open the SpO_2 window.
3	NIBP key	Press to open the NIBP window and change the NIBP measurement mode.
4	START/STOP key	Measures the NIBP in the selected mode. Pressing again during measurement stops measurement.
5	MULTI 1, MULTI 2 keys	Press to open the windows of the parameters connected to the MULTI 1 or MULTI 2 socket.
6	ECG RESP socket	Connects to the ECG connection cord.
7	SpO ₂ socket	Connects to the SpO_2 connection cord.
8	NIBP socket	Connects to the air hose.
9	MULTI 1, MULTI 2 sockets	Connects to the connection cord of the parameter to be measured (IBP, temperature, CO, CO_2 , FiO_2 or respiration by thermistor method). The type of parameter is automatically recognized. When using the BP output signal, use the MULTI 1 socket.
10	ECG/BP OUT socket	Outputs 100 mmHg/V IBP waveform of the pressure connected to the MULTI 1 socket, 1 mV/V ECG waveform of the first trace and heart rate trigger. Refer to the below warning.

Using the Output Signal from the ECG/BP OUT Socket

WARNING

When using the output signal from the module as the synchronization signal for other equipment such as IABP (intra-aortic balloon pump) or defibrillator:

- Set the timing of the other equipment by checking the waveform on the monitoring screen.
- Check the condition of the bedside monitor at all times. The output signal may become unstable.
- The hum filter is always set to on for the ECG output. ECG output differs from the ECG on the monitor screen regardless of the filter setting.
- Check that the delay time of the output signal (heart rate trigger 100 ms maximum) is within the range of the connected equipment. Refer to Section 11 of the operator's manual.
- Do not use the heart rate trigger as the synchronization signal for defibrillator.

NOTE

The output signal from the ECG/BP OUT socket may become unstable in the following conditions.

- Electrode is dry or detached.
- Electrode lead is damaged or disconnected from the electrode.
- Electrode lead is pulled.
- AC interference or EMG noise superimposed.
- Air bubbles or blood clog in the circuit for monitoring IBP.
- Any cord or cable disconnected or damaged.

Using MULTI Sockets for CO and CO, Monitoring

WARNING

- When performing defibrillation during CO monitoring, never touch the CO connection cord. Otherwise the discharged energy may cause serious electrical burn, shock or other injury.
- When performing defibrillation during CO₂ monitoring with the TG-900P CO₂ sensor kit, remove the sensor from the patient. When the sensor cannot be removed, do not touch the sensor cable because the discharged energy may cause serious electrical burn, shock or other injury.

NOTE

- CO monitoring using the MULTI socket does not comply with the Defibrillator proof type CF.
- CO₂ monitoring using the MULTI socket does not comply with the Defibrillator proof type BF.

Smart Module

CAUTION

When pressing the keys on the module, use your finger. Do not press the keys with a sharp object. Otherwise keys may be broken.



Using MULTI Sockets for CO and CO, Monitoring

WARNING

- When performing defibrillation during CO monitoring, never touch the CO connection cord. Otherwise the discharged energy may cause serious electrical burn, shock or other injury.
- When performing defibrillation during CO₂ monitoring with the TG-900P CO₂ sensor kit, remove the sensor from the patient. When the sensor cannot be removed, do not touch the sensor cable because the discharged energy may cause serious electrical burn, shock or other injury.

NOTE

- CO monitoring using the MULTI socket does not comply with the Defibrillator proof type CF.
- CO₂ monitoring using the MULTI socket does not comply with the Defibrillator proof type BF.

Recorder Module



Connection Diagram



Block Diagram


Section 2 Troubleshooting

Troubleshooting Table	2.1
How to Use the Troubleshooting Table	2.1
Power-Related Problem	2.2
Display Problems	2.2
Sound Problem	2.2
Key Operation Problems	2.3
Recorder Problem	2.3
Other Module-Related Problem	2.3

Troubleshooting Table

Use the troubleshooting table to locate, identify and solve a problem in the instrument. The problems are divided into general problem areas. Each category has its own troubleshooting table for fast and easy troubleshooting.

- Power-Related Problems
- Display Problems
- Sound Problems
- Key Operation Problems
- Recorder Problems
- Other Module-Related Problems

Refer to Section 21 "Maintenance" in the operator's manual.

How to Use the Troubleshooting Table

- 1. Determine which troubleshooting table to use.
- 2. In the "Problem" column, find the trouble item that matches the problem.
- 3. Do the action recommended in the "Action" column.
- 4. If the problem is not solved, do the action for the next possible cause or criteria.
- 5. If none of the actions solve the problem, contact your Nihon Kohden representative.

2. TROUBLESHOOTING

Power-Related Problem

Problem	Possible	e Cause/Criteria	Action
The power of the	The AC power switch on the r	ear panel is set to off.	Set the AC power switch to on.
main unit does not	When the standby lamp does	No AC power input.	Check the AC power input.
	not light even if the AC power switch is set to on.	Faulty power cord.	Replace the power cord.
		One or both of the AC inlet fuses are	Determine and correct the cause of the
		blown.	blown fuse, then replace the fuse.
		Faulty power supply unit.	Replace the power supply unit.
	It happens after the main unit	Poor internal connection.	Check the continuity of cables and boards.
	is assembled.		

Display Problems

Problem	Possible Cause	Action	
No display.	NIBP START/STOP key on the	Faulty LCD unit.	Replace the LCD unit.
	control panel is operational with AY-	Faulty DC-AC inverter.	Replace the DC-AC inverter.
	900PA/910PA.	Faulty MAIN board.	Replace the MAIN board.
	The power lamp does not light.	Faulty MAIN board.	Replace the MAIN board.
The screen is distorted or	The same part* of the LCD has an abnormal on any screen.	Faulty LCD unit.	Replace the LCD unit.
partially abnormal.	The LCD unit has an abnormal on specified Display Check screens.	Faulty MAIN board	Replace the MAIN board.
Diagnostic Check and System	The error message is shown at the Power ON Check Result on the		Replace the board considered faulty
Setup screen is displayed.	screen.		with the error message.
System error is displayed.	Faulty MAIN board.		Replace the MAIN board.

* For TFT LCD screen, it is considered normal if some pixels have abnormal color or do not light.

Sound Problem

Problem	Possible Cause/Criteria	Action
No sound.	Poor contact between the speaker and power supply	Check the continuity of the connection cable and socket at
	unit.	the power supply unit.
	Faulty power supply unit.	Replace the power supply unit.
	Faulty MAIN board.	Replace the MAIN board.

Key Operation Problems

Problem	Possi	Action	
No key operation.	Any key is not operational.	Poor internal connection.	Check the continuity of the cables and boards in between.
		Faulty OPERATION board.	Replace the OPERATION board.
		Faulty MAIN board.	Replace the MAIN board.
	Specified key is not	The key on the IR DETECT board.	Replace the IR DETECT board.
	operational.	The key on the MAIN board.	Replace the MAIN board.
No remote control operation.	Any key is not operational.	Remote control channel is set to a specified cahnnel.	Confirm which channel is set on the Setup Menu window.
		The two batteries are weak.	Replace the two batteries.
		Faulty IR DETECT board.	Replace the IR DETECT board.
		Faulty remote control.	Replace the remote control.
	Specified key is not operational.	Faulty key switch on the remote control.	Replace the remote control.

Recorder Problem

Problem	Possible Cause/Criteria		Action
No operation of the	There is an error at	Faulty JA motherboard.	Replace the JA motherboard.
recorder module. JA I/F Check.		Faulty MAIN board.	Replace the MAIN board.
	There is no error at	Deen incertion of the needed module	Completely insert the recorder moule into the
	JA I/F Check.		slot.
Faulty		Faulty recorder module.	Replace the recorder module.

Other Module-Related Problem

Problem	I	Action	
No waveform display.	Any waveform is not displayed on the screen.	Faulty JA motherboard.	Replace the JA motherboard.
	Waveform of the specified parameter is not	Breaking of the input cable or faulty connector, sensor, transducer.	Replace the corresponding cable.
	displayed.	Faulty module.	Replace the module.
		Faulty slot.	Replace the JA motherboard.

Section 3 Diagnostic Check

Introduction	3.1
Power On Self Check	
Calling up the Diagnostic Check and System Setup Screen	
MU Manual Check	
Memory Check	
Flash ROM (program) Check	
Flash ROM (data) Check	
SRAM Check	
DRAM Check	
Com Check	
Network I/F Check	
Serial I/F Check	
JA I/F Check	
Display Check	
Frame Mem Check	
Graphic Check	
Waveform Check	3.14
Backlight Check	
Key LED Check	
Key Check	
Remote Check	
Alarm Indicator Check	
Alarm Pole Check	
Other Check	
Sound Check	
Power Check	
Card I/F Check	
Timer IC Check	

Introduction

The instrument has an automatic power on self check as well as a complete set of diagnostic checks that you can perform at any time.

All errors detected during the power on self check, diagnostic checks, and any time during operation are stored in the error history table.

The diagnostic checks, error history, system setup, and initialization are accessed from the Diagnostic Check and System Setup screen.

In this section, functions which are displayed on the screen are indicated by brackets, for example, the [Monitor Mode] function on the Diagnostic Check and System Setup screen.

Power On Self Check

This self check is performed every time the power switch on the front panel is turned on. "Check Program Running" message appears during the power on self check. If no error is detected, the normal operating mode begins and the patient monitoring display appears. If an error is detected, the screen displays one of the following items according to the check items.

Check Item	How to Check	Action at Error Occurrence
Flash ROM Check	The sum of the stored data from the beginning	The check is interrupted. An error message is
for system program	address to the last second address is compared with	displayed.
	the prestored check sum at the last address.	
DRAM Check	The test patterns written to the DRAM are compared	The check is interrupted. An error is displayed.
	with the test patterns which were read out from the	
	DRAM.	
Flash ROM Check	The data which is read out from the flash ROM is	The flash ROM is initialized and its message is
for Setting Conditions	compared with the fixed data.	displayed. Next check item follows.
JA RAM Check	The test patterns written on the JA RAM are	An error message is displayed. Next check item
	compared with the test patterns which were read out	follows. After the power on self check, the
	from the JA RAM.	Diagnostic Check and System Setup screen appears.
SRAM Backup	The data read out from the SRAM is compared with	The SRAM is initialized and its message is
Check	the previous data which was written on the SRAM.	displayed. Next check item follows.
Timer IC Memory	The data read out from the real time clock IC built-	The real time clock IC built-in memory is initialized
Backup Check	in memory is checked.	and its message is displayed. Next check item
		follows.

Calling up the Diagnostic Check and System Setup Screen

 With the power off, press the power switch on the front panel while pressing the SILENCE ALARMS key on the control panel. Continue pressing the SILENCE ALARMS key until the Diagnostic Check and System Setup screen appears.

Power DN Ch	*** Di Ieck Result	agnostic Chec OK	k and System	Setup ***	
MU	Modu	lle	EXT JA	Module	(EXT JA)
MU-950R					
(MAIN) Ver 00-53					
(800T) Ver 00-22					
MU Manual Check	JA Manual Check	System Setup		System Initialize	Monitor Mode

Power ON Check Result: When no error is found during the power on self check, OK is displayed. If an error is found during the power on self check, ERROR and the check item at which the error is found are displayed. MU: Shows model number of the main unit, system program version, and boot program version Module: Shows model numbers and software versions of the modules inserted into the slots of the main unit. EXT JA: Shows description of the optional unit connected to the EXT JA socket. Module (EXT JA): Shows model number and software version of any optional unit connected to the EXT JA socket of the main unit. While this screen is displayed, the connection check results of Module, EXT JA, Module (EXT JA) are updated in real time.

2. To exit the Diagnostic Check and System Setup screen and return to the patient monitoring mode, select the [Monitor Mode] using the left or right arrow key on the control panel and pressing the center key on the control panel.

To perform the diagnostic check, system setup, or initialization or to view the error history, select one of the functions at the bottom of the screen using the left or right arrow key on the control panel and pressing the center key on the control panel.



MU Manual Check

Select [MU Manual Check] using the left or right arrow key on the control panel and pressing the center key on the control panel. The MU Manual Check screen appears.

On the screen, using the up or down arrow key on the control panel and pressing the center key on the control panel allows you to select the various function check screens. To return to the Diagnostic Check and System Setup screen, select [Return].



Memory Check

On the MU Manual Check screen, select [Memory Check]. The Memory Check screen, which checks the memory on the MAIN board, appears. On the Memory Check screen, the following checks are available. To return to the MU Manual Check screen, select [Return].

Memory Check	1
Flash ROM (program) Check	1
Flash ROM (data) Check	
SRAM Check	
DRAM Check	
Return	
	-

Flash ROM (program) Check

This screen allows you to check if the program stored in the flash memory has an error.

To start the check, select [Start]. "Checking..." message is displayed during the check. After the check, OK or Error appears.

To return to the Memory Check screen, select [Return].

*** Flash ROM (program) Check ***	
Checking	
Start	
Return	

Flash ROM (data) Check

This screen allows you to check if the writable area of the flash memory has an error by writing test data and reading it after deleting the original data. To start the check, select [Start]. A "Erasing..." message is displayed while the data is being deleted. After the data is deleted, OK appears. A "Checking..." message is displayed during the check. After the check, OK or Error appears. To return to the Memory Check screen, select [Return].

*** Flash ROM (data) Check ***	
Checking	
Start	
Return	

SRAM Check

This screen allows you to check if the SRAM has an error by writing test data and reading it after deleting the original data.

To start the check, select [Start]. A "Checking..." message is displayed during the check. After the check, OK or Error appears. If an error occurs, the Error is displayed with the address.

To return to the Memory Check screen, select [Return].

*** SRAM Check ***	
Checking	
Start	
Return	

DRAM Check

This screen allows you to check if the DRAM has an error by writing test data and reading it after deleting the original data.

To start the check, select [Start]. A "Checking..." message is displayed during the check. After the check, OK or Error appears. If an error occurs, the Error is displayed with the address.

To return to the Memory Check screen, select [Return].

*** DRAM Check ***	
Checking	
Start	
Return	

Com Check

On the MU Manual Check screen, select [Com Check]. The Com Check screen, which checks each communication interface on the MAIN board, appears. On the Com Check screen, the following checks are available. To return to the MU Manual Check screen, select [Return].

*** Com Check ***	
Network I/F Check	
Serial I/F Check	
JA I/F Check	
Return	

Network I/F Check

This screen allows you to check the network interface, interface loopback, and network connection.

	*** Networ	k I/F Check ***	
	MAC Address	00A02Axxxxx	
	Reply		
	Traffic		
	Return		

Reply Check

This screen displays the MAC addresses of which the other bedside monitors (up to 48 monitors) in the network transmit the acknowledge signals to this bedside monitor.

		*** Reply Check **	*		
00A0	2Axxxxxx	00A02Axxxxxx	00A02Axxxxxx		
00A0	2Axxxxxx	00A02Axxxxxx	00A02Axxxxxx		
00A0	2Axxxxxx	00A02Axxxxxx	00A02Axxxxxx		
	•	•	•		
	•	•	•		
	·	·	·		
	Return				

Traffic Check

This screen displays the traffic status in the network when test data is uploaded to the network.

*** Traffic Check ***
Return

MAC Address Write

Since this function is hidden, press the SILENCE ALARMS key on the control panel to display the [Adr WR] in addition to the [Reply] and [Traffic] functions on the Network I/F Check screen. Select the [Adr WR] to display the MAC Address Write screen which allows you to change the MAC address of the bedside monitor. Use the up, down, left, or right arrow key on the control panel to move the cursor to the following items in the table. Press the center key on the control panel to execute the item.

Network I/F Check screen when the SILENCE ALARMS key is pressed

***	Netwo	ork l	/F (Cheo	:k **	*			
MAC Address 00A02Axxxxx									
							_		
4	dr W	R							
F	Reply								
	Traffic								
	Returr	۱							
			ļ						
*** N	ЛАС А	٩dd	ress	Wr	ite *	**			
	00A0	09A	xxx	xx					
0 1 2	3	4	5	6	7	8 9	,		
A B C	; D	Е	F	÷	→	Save Can	sel .		

- 0 to F: Hexadecimal numbers for input of the MAC address
- \leftarrow , \rightarrow : Moves the MAC address input position left or right.
- Save: Stores the entered MAC address and returns the screen to the Network I/F Check screen.
- Cancel: Leaves the previous MAC address and returns the screen to the Network I/ F Check screen.

Serial I/F Check

This screen allows you to perform the loopback test of the serial interface RS-232C under the following conditions by connecting the loopback tester* to the RS-232C port.

Parity:NoneData length:8 bitsStop bit:1 bitBaud rate:9600 bps* Make the tester locally.Refer to figure xxx.

*** Serial I/F Check ***
ок
Start
Return

To start the check, select [Start]. A "Checking..." message is displayed during the check. After the check, OK or Error appears.

To return to the Com Check screen, select [Return].

JA I/F Check

This screen allows you to check the communication interfaces for any modules and external units.

The following checks are available.

*** JA I/F Check ***					
	RAM Check 1				
	Ram Check 2				
	Internal Loopback Check				
	External Loopback Check				
	CTL Check				
	Return				

RAM Check 1

This screen allows you to check the RAM, which stores the communication data between the modules and MAIN board, by writing the test data into the RAM and reading the stored data from the RAM. The addresses 0000H to FFFFH (Bank 1) and 0000H to 07FFH (Bank 2) are used for this check.

To start the check, select [Start]. "Checking..." message is displayed during the check. After the check, OK or Error appears. If an error occurs, the Error is displayed with the bank and address.

To return to the JA I/F Check screen, select [Return].

*** RAM Check 1 ***	
Checking	
Start	
Return	

RAM Check 2

This screen allows you to check the RAM, which stores the communication data between the modules and MAIN board, by writing the test data into the RAM and reading the stored data from the RAM. The addresses 0000H to 03FFH (Bank A) and 0400H to 07FFH (Bank B) are used for this check.

To start the check, select [Start]. A "Checking..." message is displayed during the check. After the check, OK or Error appears. If an error occurs, the Error is displayed with the bank and address.

To return to the JA I/F Check screen, select [Return].

*** RAM Check 2 ***					
	Checking				
	Start				
	Return				

Internal Loopback Check

NOTE

Remove all the modules from the main unit before starting the internal loopback check.

This screen allows you to check the loopback from the JA motherboard by transmitting 512 byte packet (all the data is either zero or one) to the JA motherboard. When all the data is set to one, the source address changes to FEH. To start the check, select [Start]. A "Checking..." message is displayed during the check. When the data returns to the MAIN board within 5 ms after the data is transmitted to the JA motherboard, OK appears. Otherwise, Error appears. To return to the JA I/F Check screen, select [Return].

*** Internal Loopback Check ***	
ОК	
Start	
Return	

External Loopback Check

NOTE

Connect the Loopback tester to the JA socket of the main unit before starting the external loopback check.

This screen allows you to check the loopback from the EXT JA board by transmitting a 512 byte packet (all the data is either zero or one) to the EXT JA board. When all the data is set to one, the source address changes to FEH. To start the check, select [Start]. A "Checking..." message is displayed during the check. When the data returns to the MAIN board within 5 ms after the data is transmitted to the EXT JA board, OK appears. Otherwise, Error appears. To return to the JA I/F Check screen, select [Return].

*** External Loopback Check ***
ок
Start
Return

CTL Check

NOTE

Connect the loopback tester to the JA socket of the main unit before starting the CTL check.

This screen allows you to check the control signal by transmitting it to the EXT JA board every 16 ms.

To return to the JA I/F Check screen, select [Return].

*** CTL Check ***	
Return	

Display Check

On the MU Manual Check screen, select [Display Check]. The Display Check screen, which checks the LCD functions, appears. On the Display Check screen, the following checks are available. To return to the MU Manual Check screen, select [Return].

*** Display Check ***	
Frame Mem Check	
Graphic Check	
Waveform Check	
Backlight Check	
Return	

Frame Mem Check

This screen allows you to check the graphic frame memory (character generation area included) and waveform frame memory (waveform data buffer included) by writing the test data into the memories and reading the stored data from the memories.

To start the check, select [Start]. A "Checking..." message is displayed during the check. After the check, OK or Error appears. If an error occurs, the Error is displayed with the address.

To return to the Display Check screen, select [Return].

*** Frame Mem Check ***	
ОК	
Start	
Return	

Graphic Check

This screen allows you to check the graphic screen display by displaying the following screens.

To select them, use the up, down, left, or right arrow key on the control panel. White whole screen \rightarrow Red whole screen \rightarrow Green whole screen \rightarrow Blue whole screen \rightarrow "H" pattern screen \rightarrow White whole screen

To return to the Display Check screen, press the center key on the control panel.

Waveform Check

This screen allows you to check the waveform screen display by displaying the following screen.

There are 12 rectangular waveforms which show the waveform traces 1 to 6 on the W1 screen and the waveform traces 1 to 6 on the W2 screen.

To return to the Display Check screen, select [Return].



- A: 10 pixels
- B: 100 pixels
- C: 20 pixels

Waveform: white Smoothing bit S1: red Smoothing bit S2: blue Smoothing bit S3: green



Backlight Check

This screen allows you to check the backlight brightness control by displaying the white whole screen.

Use the up or down arrow key on the control panel to increase or decrease the brightness.

To return to the Display Check screen, press the center key on the control panel.

Key LED Check

On the MU Manual Check screen, select [Key LED Check]. The Key LED Check screen, which checks the keys and LED functions, appears. On the Key LED Check screen, the following checks are available. To return to the MU Manual Check screen, select [Return].



Key Check

This screen allows you to check the front power switch and the keys on the control panel.

All the keys are displayed on the screen. Press one of the keys and confirm that the corresponding key display is reverse shaded in the color on the screen. If not, the key may have a failure.

If the center key has no failure, the screen returns to the Key LED Check screen after the corresponding key display is reverse shaded when the center key on the control panel is pressed.

NOTE

Pressing the front power switch does not turn off the main unit in this check.



Remote Check

This screen allows you to check the keys on the remote control. All the keys are displayed on the screen. Press one of the keys and confirm that the corresponding key display is reverse shaded in the color on the screen. If not, the key may have a failure.

If the selection knob has no failure, the screen returns to the Key LED Check screen after the corresponding key display is reverse shaded when the selection knob is pressed at the center of its position.

Pressing the center key on the control panel is another way to return to the Key LED Check screen.

NOTE

Pressing the power button does not turn off the main unit in this check.



Alarm Indicator Check

This screen allows you to check the LEDs of the alarm indicator. To select the LED color (red or yellow), use the up, down, left, or right arrow key on the control panel. The selected color LEDs blink at the alarm indicator. To return to the Key LED Check screen, press the center key on the control panel so that [Return] is executed

*** Alarm Indicator Check ***
YELLOW
Return

Alarm Pole Check

This screen allows you to check the LEDs of the optional alarm pole.

Connect the alarm pole to the main unit and use the up, down, left, or right arrow key on the control panel to select the color (red, yellow, or green). The selected color LEDs blink at the alarm pole.

To return to the Key LED Check screen, press the center key on the control panel so that [Return] is executed.

	*** Alarm Pole Check ***
	YELLOW
[Return

Other Check

On the MU Manual Check screen, select [Other Check]. The Other Check screen appears. On the Other Check screen, the following checks are available. To return to the MU Manual Check screen, select [Return].

*** Other Check ***	
Sound Check	
Power Check	
Card I/F Check	
Timer IC Check	
Return	

Sound Check

This screen allows you to check the FM sound source and audio amplifier. The following checks are available.



[Sound 1]:	Generates each alarm sound continuously.			
[Sound 2] to [Sound 5]:	Generate a musical scale of one octave. Sound 2 is			
	lowest pitch. The difference between Sound 2 and Sound			
	3, between Sound 3 and Sound 4, and between Sound 4			
	and Sound 5 is one octave.			
[Return]:	Returns to the Other Check screen.			

Power Check

This screen allows you to check the status signals from the power supply unit and 24 V for an external unit in future. The following checks are available. This screen checks each status (AC fail status, AC line frequency automatic detection, +24 V over current fail status) in real time.

To return to the Other Check screen, press the center key on the control panel so that [Return] is executed.

*** Power Check ***		
ACFAIL Status	ОК	
AC Frequency	50 Hz	
+24 V Over Current	OK	
Return		

Card I/F Check

This screen allows you to check the memory card interface by writing the test data into the memory card and reading the stored data from the memory card.

To start the check, insert the optional memory card into the slot of the main unit and select [Start]. A "Checking..." message is displayed during the check. After the check, OK or Error appears. If an error occurs, the Error is displayed with the address.

To return to the Other Check screen, select [Return].

*** Card I/F Check ***
Insert QM-010P (1 MB SRAM)
ОК
Start
Return

Timer IC Check

This screen allows you to check the memory of the real time clock IC by writing the test data into the memory and reading the stored data from the memory. The current time is displayed on the screen. Before the check, the previous data stored into the memory is backed up. The backed-up data is written into the memory again when the screen returns to the Other Check screen after the check. To start the check, select [Start]. A "Checking…" message is displayed during the check. After the check, OK or Error appears. If an error occurs, the Error is displayed with the address.

To return to the Other Check screen, select [Return].

*** Timer IC Check ***						
	Timer	1999/06/14	17:24:46			
	Memory	ОК				
Γ	Start					
	Retur	n				

Section 4 Board/Unit Description

Signal Flow	4.1
Vital Sign Signals from Patients	4.1
Display Data	4.1
Power Control Signal by Front Power Switch or Power Button on	
Optional Remote Control	4.1
MAIN Board UR-3485	4.2
EXT JA Board UR-3489	4.2
JA Motherboard UR-3486	4.2
IR DETECT Board UR-3487	4.3
LCD JUNC Board UR-3504	4.3
LED Board UR-3393	4.3
OPERATION Board UR-3506	4.3
Power Supply Unit SC-036R	4.4

Signal Flow

Vital Sign Signals from
PatientsThe vital sign signals such as ECG/respiration wave, plethysmographic pulse
wave, blood pressure waveform, and NIBP go to the modules through the electrode
leads, sensor probe, transducer, and cuff and tube. These signals are digitized in
the modules. The digitized data is transferred to the MAIN board through the JA
motherboard.Display DataOn the MAIN board, the CPU writes the binary code waveform data into the WAVE
RAM area of the CRT controller (gate array IC). The CRT controller converts the
waveform data to graphic data and writes it into the WAVE RAM area. The CPU
writes alphanumeric data and review data such as trendgraph into the GRAPH
RAM area of the CRT controller. The CRT controller mixes the data in the WAVE
RAM and GRAPH RAM area to create the display data for screen. The CRT
controller converts the display data to RGB signals which go to the LCD unit.

Power Control Signal by Front Power Switch or Power Button on Optional Remote Control When the AC power switch on the rear panel is set to on, the BSM-9510 is on standby. The MAIN board outputs the power control signal to the power supply unit when the power switch on the front panel or power button on the optional remote control is pressed during standby. The power control signal allows the BSM-9510 to turn on or be on standby.

MAIN Board UR-3485

The MAIN board has one 68EN360 which is used as the CPU of the BSM-9510. The CPU clock is 25 MHz.

The following components are mounted on this board.

• Flash memory for system program	h: Has the ROM and RAM area on which the
	system program runs.
• Flash E2PROM:	Included in the flash memory. Stores the
	conditions set on the System Setup screen.
• Real time clock (RTC):	Used for date and time.
• Key operation controller:	Reads the status of the front power switch and
	each key on the control panel
• Ethernet driver:	Used for local area network LAN.
• RS232C driver:	Used for communication with an external unit
	in future or a personal computer.
• CRT controller:	Converts the screen data to RGB signals and
	outputs them to the LCD unit.

EXT JA Board UR-3489

This daughter board has a socket for connection of external units in the future and is mounted on the MAIN board.

JA Motherboard UR-3486

This board has 6 sockets for modules. This board is connected to the MAIN board. The MAIN board communicates with the modules through this board. External units in the future will also communicate with the MAIN board through the EXT JA board and this board.

IR DETECT Board UR-3487

This board has the infrared sensor for remote control and outputs the remote control signals to the MAIN board.

LCD JUNC Board UR-3504

This daughter board is mounted on the LCD unit and converts the signals from the MAIN board to the LCD unit.

LED Board UR-3393

This board has LEDs for alarm indicator and is mounted inside the top of the display.

OPERATION Board UR-3506

This board has switches for key operation and is mounted into the control panel.

Power Supply Unit SC-036R

The power supply unit has the following outputs.

- $+5V_{D}$ constant voltage output: Used for all the digital circuits.
- $+5V_s$ constant voltage output: Used for the front power switch activation in standby mode.
- +12 V_U constant voltage output: Used for the LCD backlight and alarm indicator LEDs
- $+V_{IA}$ constant voltage output: Used for the JA motherboard and modules.
- $+24V_{c}$ constant voltage output: Will be used for future external units.

The power switch on the front panel turns on or off the outputs except $+5V_s$ when the $+5V_s$ output is turned on. If the $+5V_s$ output is turned off, the power switch has no function. The AC power switch on the rear panel turns on the $+5V_s$ output only or turns it off together with the other outputs.

When an external unit (future) is connected to the socket, the $+24V_{c}$ output is turned on.

This power supply unit allows BSM-9510 to resume monitoring after power recovery if the BSM-9510 has momentary power loss during monitoring.

Section 5 Disassembly and Assembly

Before You Begin	5.1
Warnings and Cautions	5.1
Required Tools	5.1
Replacing MAIN Board	5.2
Replacing Fuse	5.5
Replacing Power Supply Unit	5.7
Replacing JA Motherboard	5.10
Replacing IR DETECT Board	5.13
Replacing LCD Unit	5.16
Replacing DC-AC Inverter	5.19
Replacing Backlight Lamps	
Replacing LCD Filter	5.26
Replacing OPERATION Board	5.29
Replacing Lithium Battery	5.32

The procedures in this section tell how to remove, replace and install major components in the instrument.

Before You Begin

Removing, replacing and installing major components should be done by qualified service personnel.

Warnings and Cautions

WARNINGS

- To avoid the possibility of injury to yourself or damage to the instrument, do not install or remove any component or change switch settings while the power is on and wait 10 minutes after the power is off before installing or removing any component from the instrument.
- To avoid accidental discharge of static electricity which could damage the instrument components, use a wrist ground strap when installing or removing any component of the instrument.

CAUTIONS

- Before connecting or disconnecting any cables, turn off the instrument and unplug the AC power cord from the instrument.
- Fuses cut off the power when an abnormality occurs in the instrument. Eliminate the malfunction before replacing the fuse. Use the correct fuse* only. The fuse rating is shown on the printed circuit board of the power supply unit.
- Removal and replacement of any component in the instrument should be done by qualified service personnel.
- Use only parts recommended by Nihon Kohden to assure maximum performance from your instrument.

* Refer to Section 8 "Replaceable Parts List".

Required Tools

- Anti-static bench mat
- Wrist ground strap
- M3 Phillips screwdriver (insulated type)
- M4 Phillips screwdriver (insulated type)
- Small flat-blade screwdriver (insulated type)

Replacing MAIN Board



1. Loosen the 2 screws (M4 binding head screws) at the bottom of the main unit.

2. Remove the 2 screws (M4 binding head screws) at the top of the rear cover and the screw (M3 pan head screw with spring washer) at the rear of the control panel.



3. Slightly pull out the front cover assy until the 2 loosened screws are in the clearance holes at the bottom of the front cover assy. Move the front cover assy downward until the screws come out of the clearance holes. Put the front cover assy on a smooth, soft surface to avoid scratching the screen.



- 4. Loosen the screw (M3 pan head screw with spring and flat washers) which holds the LCD JUNC board with ribbon cable. Remove it from the front cover assy.
- 5. Disconnect the cable between the MAIN board and IR DETECT board at the IR DETECT board to separate the front cover assy from the main unit.
- Lay the main unit on its back. Remove the 2 screws (M3 binding head screws) which hold the ferrite core on its holder. Disconnect all 6 cables at the MAIN board. The 6 cables are as follows:
 - Cable between JA Motherboard and MAIN boards
 - Cable between LCD JUNC and MAIN boards
 - Cable between IR DETECT and MAIN boards
 - Cable between LED and MAIN boards
 - 2 cables between power supply unit and MAIN board

Remove the 12 screws (M3 pan head screws with spring washer) which hold the shield case 2 and MAIN board. Remove them from the main unit.



7. Replace the MAIN board with a new one and assemble the main unit by reversing the above procedure.

CAUTION

Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing Fuse

- 1. Remove the MAIN board according to the "Replacing MAIN board" section.
- 2. Remove the 12 screws (M3 pan head screws with spring washer) which hold shield case 1. Remove the shield case from the main unit.



There is a square hole as shown below. You can see the two fuses through the square hole.



3. Replace the blown fuse(s) with a new one and assemble the main unit by reversing the above procedure.

CAUTION

Be careful not to pinch or stretch the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of the front cover assy and the alarm indicator of the main unit.


Replacing Power Supply Unit

- 1. Remove the MAIN board according to the "Replacing MAIN Board" section.
- 2. Remove the 12 screws (M3 pan head screws with spring washer) which hold shield case 1. Remove it from the main unit.



3. Remove the 2 screws (M3 pan head screws with spring washer) which hold the power supply unit to the main unit. Remove it from the main unit.



4. Remove the 2 screws (M3 pan head screws with spring and flat washers) which hold the speaker. Remove it from the power supply unit. Disconnect the speaker cable from the power supply unit.



5. Remove the 2 screws (M3 pan head screws with spring washer) which hold the fan with holder. Remove it from the power supply unit. Disconnect the fan cable from the power supply unit.



6. Replace the power supply unit with a new one and assemble the main unit by reversing the above procedure.

CAUTION

Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing JA Motherboard

- 1. Remove the MAIN board according to the "Replacing MAIN Board" section.
- 2. Remove the 12 screws (M3 pan head screws with spring washer) which hold the shield case 1. Remove it from the main unit.



3. Remove the 2 screws (M3 binding head screws) which hold the JA socket holder. Remove it from the main unit.



4. Remove the 2 screws (M3 binding head screws) which hold the module case and JA motherboard to the main unit. Remove them from the main unit.



 Remove the 5 screws (M3 pan head screws with spring washer) which holds the module case shield cover. Remove it from the module case. Remove the 3 screws (M3 pan head screws with spring and flat washers) which hold the JA motherboard. Remove it from the module case.



6. Replace the JA motherboard with a new one and assemble the main unit by reversing the above procedure.

CAUTION

Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing IR DETECT Board



1. Loosen the 2 screws (M4 binding head screws) at the bottom of the main unit.

2. Remove the 2 screws (M4 binding head screws) at the top of the rear cover and the screw (M3 pan head screw with spring washer) at the rear of the control panel.



3. Slightly pull out the front cover assy until the 2 loosened screws are in the clearance holes at the bottom of the front cover assy. Move the front cover assy downward until the screws come out of the clearance holes. Put the front cover assy on a smooth, soft surface to avoid scratching the screen.



- 4. Loosen the screw (M3 pan head screw with spring and flat washers) which holds the LCD JUNC board with ribbon cable. Remove it from the front cover assy.
- 5. Disconnect the cable between the MAIN board and IR DETECT board at the IR DETECT board to separate the front cover assy from the main unit.
- Disconnect the cable between the IR DETECT board and DC-AC inverter at the IR DETECT board. Remove the 2 screws (M3 pan head screws with spring washer) which hold the IR DETECT board. Remove it from the front cover assy.



7. Replace the IR DETECT board with a new one and assemble the main unit by reversing the above procedure.

CAUTION Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing LCD Unit



1. Loosen the 2 screws (M4 binding head screws) at the bottom of the main unit.

2. Remove the 2 screws (M4 binding head screws) at the top of the rear cover and the screw (M3 pan head screw with spring washer) at the rear of the control panel.



3. Slightly pull out the front cover assy until the 2 loosened screws are in the clearance holes at the bottom of the front cover assy. Move the front cover assy downward until the screws come out of the clearance holes. Put the front cover assy on a smooth, soft surface to avoid scratching the screen.



- 4. Loosen the screw (M3 pan head screw with spring and flat washers) which holds the LCD JUNC board with ribbon cable. Remove it from the front cover assy.
- 5. Disconnect the cable between the MAIN board and IR DETECT board at the IR DETECT board to separate the front cover assy from the main unit.
- 6. Disconnect the cable between the IR DETECT board and DC-AC inverter at the DC-AC inverter. Remove the 4 screws (M3 pan head screw with spring and flat washers) and screw (M3 pan head screw with spring washer) which hold the LCD unit and DC-AC inverter. Remove them from the front cover assy.



7. Disconnect the 2 cables from the DC-AC inverter to separate the LCD unit and DC-AC inverter.



8. Replace the LCD unit with a new one and assemble the main unit by reversing the above procedure.

CAUTION Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing DC-AC Inverter

- Screws
- 1. Loosen the 2 screws (M4 binding head screws) at the bottom of the main unit.

2. Remove the 2 screws (M4 binding head screws) at the top of the rear cover and the screw (M3 pan head screw with spring washer) at the rear of the control panel.



3. Slightly pull out the front cover assy until the 2 loosened screws are in the clearance holes at the bottom of the front cover assy. Move the front cover assy downward until the screws come out of the clearance holes. Put the front cover assy on a smooth, soft surface to avoid scratching the screen.



- 4. Loosen the screw (M3 pan head screw with spring and flat washers) which holds the LCD JUNC board with ribbon cable. Remove it from the front cover assy.
- 5. Disconnect the cable between the MAIN board and IR DETECT board at the IR DETECT board to separate the front cover assy from the main unit.
- 6. Disconnect the cable between the IR DETECT board and DC-AC inverter at the DC-AC inverter. Remove the 4 screws (M3 pan head screw with spring and flat washers) and screw (M3 pan head screw with spring washer) which hold the LCD unit and DC-AC inverter. Remove them from the front cover assy.



7. Disconnect the 2 cables from the DC-AC inverter to separate the LCD unit and DC-AC inverter.



8. Replace the DC-AC inverter with a new one and assemble the main unit by reversing the above procedure.

CAUTION

Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing Backlight Lamps

NOTE

Always replace both backlight lamps together even if one works well.

1. Loosen the 2 screws (M4 binding head screws) at the bottom of the main unit.



2. Remove the 2 screws (M4 binding head screws) at the top of the rear cover and the screw (M3 pan head screw with spring washer) at the rear of the control panel.



3. Slightly pull out the front cover assy until the 2 loosened screws are in the clearance holes at the bottom of the front cover assy. Move the front cover assy downward until the screws come out of the clearance holes. Put the front cover assy on a smooth, soft surface to avoid scratching the screen.



- 4. Loosen the screw (M3 pan head screw with spring and flat washers) which holds the LCD JUNC board with ribbon cable. Remove it from the front cover assy.
- 5. Disconnect the cable between the MAIN board and IR DETECT board at the IR DETECT board to separate the front cover assy from the main unit.
- 6. Disconnect the cable between the IR DETECT board and DC-AC inverter at the DC-AC inverter. Remove the 4 screws (M3 pan head screw with spring and flat washers) and screw (M3 pan head screw with spring washer) which hold the LCD unit and DC-AC inverter. Remove them from the front cover assy.



7. Disconnect the 2 cables from the DC-AC inverter to separate the LCD unit and DC-AC inverter.



8. Remove the 2 backlight lamps from the LCD unit by moving the black hook with a small flat-blade screwdriver in the direction of the arrow.



5. DISASSEMBLY AND ASSEMBLY

9. Replace the 2 backlight lamps with the new backlight lamp set as shown below. The backlight lamp set consists of one backlight lamp marked with
→ at its side and one backlight lamp without a
→ mark at its side. The lamps are fragile and must be kept clean. Gently insert the backlight lamp marked with
→ into the slot marked with
△ until the black hook clicks. Gently insert the backlight lamp without the
→ mark into the slot without the
→ mark until the black hook clicks. Assemble the main unit by reversing the above procedure.



CAUTION Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing LCD Filter



1. Loosen the 2 screws (M4 binding head screws) at the bottom of the main unit.

2. Remove the 2 screws (M4 binding head screws) at the top of the rear cover and the screw (M3 pan head screw with spring washer) at the rear of the control panel.



3. Slightly pull out the front cover assy until the 2 loosened screws are in the clearance holes at the bottom of the front cover assy. Move the front cover assy downward until the screws come out of the clearance holes. Put the front cover assy on a smooth, soft surface to avoid scratching the screen.



- 4. Loosen the screw (M3 pan head screw with spring and flat washers) which holds the LCD JUNC board with ribbon cable. Remove it from the front cover assy.
- 5. Disconnect the cable between the MAIN board and IR DETECT board at the IR DETECT board to separate the front cover assy from the main unit.
- 6. Disconnect the cable between the IR DETECT board and DC-AC inverter at the DC-AC inverter. Remove the 4 screws (M3 pan head screw with spring and flat washers) and screw (M3 pan head screw with spring washer) which hold the LCD unit and DC-AC inverter. Remove them from the front cover assy.



- 7. Remove the LCD filter together with the packing. Remove the packing from the LCD unit.
- 8. Replace the LCD filter with a new one and assemble the main unit by reversing the above procedure.

NOTE

When attaching the LCD filter with packing into the original position, be careful not to leave an edge of the LCD filter with packing on a rib of the front cover.



CAUTION Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing OPERATION Board

- Screws
- 1. Loosen the 2 screws (M4 binding head screws) at the bottom of the main unit.

2. Remove the 2 screws (M4 binding head screws) at the top of the rear cover and the screw (M3 pan head screw with spring washer) at the rear of the control panel.



3. Slightly pull out the front cover assy until the 2 loosened screws are in the clearance holes at the bottom of the front cover assy. Move the front cover assy downward until the screws come out of the clearance holes. Put the front cover assy on a smooth, soft surface to avoid scratching the screen.



- 4. Loosen the screw (M3 pan head screw with spring and flat washers) which holds the LCD JUNC board with ribbon cable. Remove it from the front cover assy.
- 5. Disconnect the cable between the MAIN board and IR DETECT board at the IR DETECT board to separate the front cover assy from the main unit.
- 6. Remove the screw (M3 binding head screw) and spacer bolt which hold the control panel assy 2. Remove it from the front cover assy.



7. Remove the 3 screws (M3 pan head screws with spring washer) which hold the OPERATION board. Remove it from the front cover assy.



8. Replace the OPERATION board with a new one and assemble the main unit by reversing the above procedure.

CAUTION

Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Replacing Lithium Battery

1. Remove the MAIN board according to the "Replacing MAIN Board" section.



2. Remove the lithium battery on the MAIN board with a soldering iron.

1

ŀ..

- 3. Solder the new lithium battery onto the original position of the MAIN board.
- 4. Assemble the main unit according to the reverse order of the "Replacing MAIN Board" section.

CAUTION

Be careful not to pinch or strain the cables.

NOTE

When attaching the front cover assy to the main unit, minimize the gap between the top of front cover assy and alarm indicator of the main unit.



Section 6 Maintenance

Maintenance Check Items	6.1
External	6.1
Safety	6.2
Modules	6.2
Display	6.3
Measuring Parameters	6.3
Recorder	6.3
Backup	6.3
Others	6.4

Maintenance Check Items

Maintenance check sheets are provided at the end of this section. Make a copy of this check sheet before using it. Check all the items on this check sheet when doing any service or maintenance. The check sheet contains the check items grouped as follows:

- External
- Safety
- Modules
- Display
- Measuring Parameters
- Recorder
- Backup
- Others

The rest of this section describes how to check each of the above items.

External

Item	Check Procedure	Action at Error Occurrence	
Dirt	Check that the outside of the instrument is not dirty.	If the outside of the instrument is dirty, clean it with a cloth moistened with neutral soap or alcohol.	
Loose screws	Check that there are no loose screws.	If any screw is loose, tighten it.	
Damaged or bent parts	Check that there are no physically damaged or bent parts. This includes switch, key top cover, cable, cord and pins on the connector or socket.	If any part is damaged or bent, replace it.	
Abnormal vibration or resonant sound	Check that the instrument has no abnormal vibration and resonant sound.	Remove the cause of abnormal vibration or resonant sound.	
Warning and caution labels	Check that the labels are clearly readable.	If a label cannot be read, replace it with a new one.	
Safety standard labels such as CSA mark and CE mark	Check that the labels are clearly readable.	Replace damaged or worn labels.	

6. MAINTENANCE

Safety

ltem	Check Procedure	Action at Error Occurrence	
Power cord and ground lead	Check that the power cord has no physical damage and no AC cutoff, no heat, no sound and no smell while the power cord is bent at each part.	Replace the power cord if it has any damage.	
	Check that the ground lead has no physical damage and no cutoff of the continuity while the ground lead is bent at each part.	Replace the ground lead if it has any damage.	
Grounding	Check that the instrument is grounded to the dedicated grounding terminal in the facility.	Use the provided ground lead to ground the instrument to a dedicated grounding terminal if the instrument is not grounded.	
Fuse	Check that the fuse is not blown.	Remove the cause and replace the fuse if it is blown.	
Fuse type	Check that the fuse is the correct fuse type (T 4 A/250 V).	Replace the fuse with the correct type if the fuse is not correct.	
AC power	Check that the input AC power is within the correct range.	Use only AC power within the correct range.	
AC operation	Check that the AC power switch, standby lamp, front power switch and power lamp work normally.	If one of them has a malfunction, remove its cause.	
Protective earth impedance (Refer to IEC601-1 18.(f).)	Check that the impedance between the protective earth contact and any accessible metal part does not exceed 0.1Ω .	Remove the cause if the impedance exceeds 0.1 Ω .	
Earth leakage current (Refer to IEC601-1 19.)	Check that the earth leakage current does not exceed 0.5 mArms under normal condition and 1.0 mArms under each single fault condition.	Remove the cause if the earth leakage current exceeds one of the maximum values.	
Enclosure leakage current (Refer to IEC601-1 19.)	Check that the enclosure leakage current does not exceed 0.1 mArms under normal condition and 0.5 mArms under each single fault condition.	Remove the cause if the enclosure leakage current exceeds one of the maximum values.	
Patient leakage current (Refer to IEC601-1 19.)	Check that the patient leakage current to type CF or defibrillation-proof type CF applied part does not exceed 0.01 mArms under normal condition and 0.05 mArms under each single fault condition.	Remove the cause if the patient leakage current exceeds one of the maximum values.	
Withstand voltage (Refer to IEC601-1 20.)	Check that the instrument has the following withstand voltage. - A-a1: 1500 VAC for one minute - A-f: 1500 VAC for one minute - B-a: 4000 VAC for one minute* - B-d: 1500 VAC for one minute* * Perform the check with the module inserted.	Remove the cause if the instrument does not have one of the withstand voltages.	

Modules

ltem	Check Procedure	Action at Error Occurrence
Input parts on modules	Check that no input part on the modules has poor contact by connecting the corresponding input cable and slightly pulling or pushing the connector.	Replace the input cable or socket of the module if the corresponding parameter disappears during the check.

Display

Item	Check Procedure	Action at Error Occurrence
LCD screen	Check that the LCD screen is neither dirty nor damaged.	Clean the LCD or replace it with a new one.
Brightness	Check that the brightness is correctly adjusted by controlling the brightness on the Display window.	Remove the cause if the brightness cannot be adjusted.
Diagnostic check	Check that the Display check has no error.	Remove the cause if there is an error.

Measuring Parameters

ltem	Check Procedure	Action at Error Occurrence		
Waveform and numeric data	Check that the waveform and numeric data are correctly displayed by using the AX-800P, AX-110G or calibration gas.	Remove the cause if the waveform and numeric da are not displayed correctly.		
Alarm function	Check that the alarm function works correctly.	Remove the cause if the alarm function has a problem.		
BP zeroing and comparison with mercury manometer	Check that the zeroing and comparison with mercury manometer has no error.	Remove the cause if there is an error.		
NIBP air leak	Check that the cuff and hose has no air leak.	Replace the cuff or hose if there is any air leak.		

Recorder

ltem	Check Procedure	Action at Error Occurrence	
Paper feeding	Check that the recorder feeds the paper correctly.	Remove the cause if there is any problem.	
Printing	Check that the recorder prints the waveforms and alphanumeric data correctly.	Remove the cause if there is any problem.	
Abnormal sound	Check that there is no abnormal sound in the recorder.	Remove the cause if there is any abnormal sound.	

Backup

ltem	Check Procedure	Action at Error Occurrence
Backup of system setup setting	Check that the setting data is correctly stored after the power off for 2 hours or more.	Contact Nihon Kohden.
Backup of alarm setting	Check that the setting data is correctly stored after turning on the instrument for 30 minutes or more and turning off for 2 hours or more.	Contact Nihon Kohden
Backup of real time clock	Check that the date and time is updated after the power off for 10 minutes or more.	Contact Nihon Kohden.

6. MAINTENANCE

Others

ltem	Check Procedure	Action at Error Occurrence
External device	Check that the external device is correctly connected to the instrument according to the operator's manual.	Follow the operator's manual.
Installation location	Check that the instrument is correctly installed in a suitable location according to the operator's manual.	Follow the operator's manual.

Maintenance Check Sheet

Fill out and save this check sheet each time you do maintenance or service. Circle each item.

	Date:		
Customer:			
Customer Address:			
Service Personnel:	Service Company:		
Instrument Name: Bedsid	e Monitor Instrument Model: BSM-9510J	K	
Instrument Serial Number	: Hardware Revision:		
Software Version:			
External	Outside of instrument is clean.	Yes	No
	All screws are tight.	Yes	No
	Instrument is not physically damaged and has no bent parts.	Yes	No
	The warning and caution labels are clearly readable	Yes	No
Safety	Power cord and ground lead are not damaged.	Yes	No
	Instrument is grounded to the dedicated grounding terminal.	Yes	No
	Fuse is not blown.	Yes	No
	Undamaged and correct fuse is used.	Yes	No
	Correct AC power is used.	Yes	No
	Protective earth resistance is within prescribed range.	Yes	No
	Earth leakage current is within prescribed range.	Yes	No
	Enclosure leakage current is within prescribed range.	Yes	No
	Patient leakage current is within prescribed range.	Yes	No
	Instrument can withstand prescribed withstand voltage.	Yes	No
Modules	Input parts on the modules have no poor contact.	Yes	No
Display	Brightness is correctly adjusted.	Yes	No
	LCD display is not dirty or damaged.	Yes	No
Measuring Parameters	Waveforms and numeric data are correctly displayed.	Yes	No
	Alarm functions correctly.	Yes	No
	BP zeroing is correct.	Yes	No
	NIBP cuff and hose has no air leak.	Yes	No
Recorder	Specified recording paper is used.	Yes	No
	Recorder does not make abnormal sounds or record		
	at an incorrect speed, and paper does not skew.	Yes	No
	There are no missing dots on the recording.	Yes	No
	The date and time are correctly recorded.	Yes	No
Backup	System settings, alarm settings, and date and time are correct	ly	
	backed up.	Yes	No
Others	External device is correctly connected.	Yes	No
	Instrument is correctly installed in a suitable location.	Yes	No

Cut here

Section 7 Adjustment

Only the following adjustment is necessary. This adjustment can be done locally. A frequency counter is necessary for this adjustment.

VC561 on the MAIN board

Adjust the variable capacitor VC561 so that the frequency counter reads 2.000000 \pm 0.000022 Hz between the test points TP561 (SQW) and TP851 (DGND).



Section 8 Replaceable Parts List
When ordering parts or accessories from your Nihon Kohden Corporation distributor, please quote the NK code number and part name which are listed in this service manual, and the name or model of the unit in which the required part is located. This will help us to promptly attend to your needs. Always use Nihon Kohden parts and accessories to assure maximum performance from your instrument.

Main Unit Parts

<u>Index</u>	NK Code No.	<u>Qty</u>	Description	
1	6141-000665A	1	Front cover assy	フロントカバASSY
2	6124-025928	1	NIHON KOHDEN plate	シャメイプレート56(721)
3	6124-028533	1	Model number label	BSM-9510 カタシキラベル
4	6124-027409B	1	Power off attention panel	デンゲンチュウイパネルA
5	6114-073556B	1	Key top cover	キートップ
6	6143-009278	1	Control panel assy 2	MU-950RA ソウサブASSY2
7	6114-084936	1	IR sensor holder	IRDAコテイカナグ
8	6114-085677	4	Alarm indicator packing	アラームインジケータパッキン
9	6113-026387D	1	Packing	フィルタパッキン
10	6114-081208A	1	LCD filter	LCDフィルタ
11	127997	1	Spacer bolt	カンカクボルトUN18-2102-0005 (L9)
12	132777	1	Flat cable clip	ミニフラットケーブルクリップFCA-10
13	6112-009906B	1	Shield case 2	シールドケース 2
14	6114-099617	3	Earth spring	ウルトラソフトフィンガ
15	6114-097869A	1	Ferrite core holder	コアホルダ
16	6114-081235B	1	Connector holder	コネクタホルダ
17	1114-176166	2	Spacer nut	カンカクナットL=13.2(C1)
18	6114-094132	1	D-sub socket spring C	D-SUB コネクタバネC
19	6114-096282	2	Heat radiation sheet (38×38)	ホウネツシート(38×38)
20	6113-031274A	1	Heat sink	ホウネツバン
21	6112-009898E	1	Shield case 1	シールドケース 1
22	6114-080869	1	Earth spring	クリップフィンガ
23	6114-080708A	1	Indicator packing	インジケータパッキン
24	6114-084963	1	Alarm indicator reflection cheet	アラームインジケータハンシャシート
25	6112-008649B	1	Alarm indicator	アラームインジケータ
26	6113-026378B	1	Socket cover	コネクタカバ
27	6113-026369B	1	JA socket holder	JA コネクタホルダ
28	6111-005075B	1	Rear cover	リアカバジョウチャクシジズ1
29	6114-099626	2	Earth spring	ウルトラソフトフィンガ(4 レン)
30	6114-096371A	1	Front cover fixing metal fittings	フロントケツゴウカナグ
31	1114-167087B	4	Rubber foot	ゴムアシ
32	2229-000798A	2	Attention mark label	チュウイラベル
33	6124-027392A	1	Memory card attention panel	メモリカードパネルA
34	6114-100919	2	Collar	カラ
35	369085	2	Grommet	グロメットゴムブッシュ φ9
36	6114-100901	1	Speaker holder	スピーカホルダ

8. REPLACEABLE PARTS LIST

37	107243	4	O-ring	Oリング B-2401 (P-3) ダイニチ
38	6114-100759	2	Shaft	シャフト
39	6114-098779A	1	Fan holder	ファンホルダ
40	6112-011609A	1	Module case shield cover	JAシールドケース
41	6114-096273	1	Heat radiation sheet (10×10)	ホウネツシート(10×10)
42	6113-031292	1	Insulation sheet	マザーゼツエンシート
43	6111-005066A	1	Module case	モジュールケースジョウチャクシジズ
44	6114-073476	1	FET metal fittings	FETトリツケカナグ
А	UR-3485	1	MAIN Board	メインボード
В	UR-3489	1	EXT JA Board	EXT JA ボード
С	SC-036R	1	Power Supply Unit	電源ユニット
	371215	2	Fuse T4 A/250 V	ヒューズ 218 004
D	UR-3486	1	JA Motherboard	JAマザーボード
Е	525121	1	LCD Unit	LCDユニットLQ10D41
			(Two backlight lamps are include	d)
	540078	1	Backlight Lamp Set (2 pcs/set)	交換用ランプユニットLQ0B192
F	524826	1	DC-AC Inverter	DC-ACインバータ CXA-0217
G	UR-3487	1	IR DETECT Board	IR DETECT ボード
Н	UR-3506	1	OPERATION Board	OPERATIONボード
Ι	UR-3393	1	LED Board	LEDボード
J	073662	1	Speaker	スピーカ マツシタ EAS-65P34S 160 Ω
Κ	558273	1	Fan	FAN F412R-12MA
L	UR-3504	1	LCD JUNC Board	LCD JUNCボード



Section 9 Connector Pin Assignment

Input/Output Connector Pin Assignment	9.1
Alarm Output Socket, ALARM	9.1
General Serial Socket, SERIAL	9.1
Network Socket, NETWORK	9.2
JA Output Socket	9.2
Connector on MAIN Board	9.2
Connector on EXT JA Board	9.3
Memory Card Connector	9.3
OPERATION Connector	9.4
LED Connector	9.4
LCD Connector	9.5
JA Motherboard Connector	9.5
Connectors on the Power Supply Unit	9.6
BDM Connector	9.6
DEBUG Connector	9.7

Input/Output Connector Pin Assignment

CAUTION

- When connecting an external instrument to the output jacks, ensure that the external instrument complies with the IEC60601-1 safety standard for medical equipment or IEC60601-2-27 particular requirements for the safety of electrocardiographic monitoring equipment. When the instrument does not comply with IEC60601-1 or IEC60601-2-27, use locally available medical use isolation transformer unit between the instrument and the AC outlet.
- The monitor should only be connected to an external instrument which complies with the IEC60601-1-2 or CISPR 11 Second Edition 1990-09, Group 1 and Class B standard.

Alarm Output Socket, ALARM



<connector></connector>

Round 6-pin connector	Model	HR12-10R-5SDL
<pin assignment=""></pin>		

Pin No.	Signal Name	Pin No.	Signal Name
1	YLRED	4	+12V
2	YLYELLOW	5	+12V
3	YLGREEN		

♦ General Serial Socket, SERIAL



<connector></connector>	
D sub 9-pin connector	Model JES-9S-2A3B
<connection cable=""></connection>	
9P-9P serial cable	Code no. YS-061P2
9P-25P serial cable	Code no. YS-061P3
<pin assignment=""></pin>	

Pin No.	Signal Name	Pin No.	Signal Name
1	CD	6	DSR
2	RXD	7	RTS
3	TXD	8	CTS
4	DTR	9	Not used
5	GND		

Network Socket, NETWORK



<Connector>

8-pin modular connector Model TM11R-3C-88

<Connection Cable>

4 twisted pair cable with shield

<Pin Assignment>

Pin No.	Signal Name	Pin No.	Signal Name
1	TX+	5	Not used
2	TX-	6	RX-
3	RX+	7	Not used
4	Not used	8	Not used

♦ JA Output Socket



<connectors></connectors>		
Washer male 15-pin right angle	Model	JAY-15S-1A2G13-N
D sub 15-pin connector		
<connection cables=""></connection>		
JA connection cable	Code no	o. YJ-900P (0.4 m)
JA connection cable	Code no	o. YJ-901P (2.5 m)
JA connection cable	Code no	o. YJ-902P (5.0 m)
<pin assignment=""></pin>		

Pin No.	Signal Name	Pin No.	Signal Name
1	+24VOUT	9	Not used
2	E2	10	CNNTIN
3	RXIN+	11	Not used
4	RXIN-	12	Not used
5	TXOUT+	13	Not used
6	TXOUT-	14	Not used
7	CTLOUT+	15	Not used
8	CTLOUT-		

Connector on MAIN Board



<Connector>

Model XH3A2041-2A

<pin assignment=""></pin>					
Pin No.	Signal Name	Pin No.	Signal Name		
a10	+24VOUT	b10	+24VOUT		
a9	JGND	b9	JGND		
a8	RXIN+	b8	RXIN-		
a7	TXOUT+	b7	TXOUT-		
a6	CTLOUT+	b6	CTLOUT-		
a5	Not used	b5	CNNTIN		
a4	Not used	b4	Not used		
a3	Not used	b3	Not used		
a2	Not used	b2	Not used		
a1	Not used	b1	Not used		

♦ Connector on EXT JA Board

<Connector>

Model XH3B2041-A



Pin No.	Signal Name	Pin No.	Signal Name
a10	+24VOUT	b10	+24VOUT
a9	JGND	b9	JGND
a8	RXIN+	b8	RXIN-
a7	TXOUT+	b7	TXOUT-
a6	CTLOUT+	b6	CTLOUT-
a5	Not used	b5	CNNTIN
a4	Not used	b4	Not used
a3	Not used	b3	Not used
a2	Not used	b2	Not used
a1	Not used	b1	Not used

Memory Card Connector

<Connector>

Model IC7-68PD-1.27SF-EJR-C

<Pin Assignment>



Pin No.	Signal Name	Pin No.	Signal Name	Pin No.	Signal Name	Pin No.	Signal Name
1	DGND	18	+5VD (VPP1)	35	DGND	52	+5VD (VPP2)
2	MEMD03	19	MEMA16	36	XCD1	53	MEMA22
3	MEMD04	20	MEMA15 <gdispout></gdispout>	37	MEMD11	54	MEMA23
4	MEMD05	21	MEMA12 <gb1></gb1>	38	MEMD12	55	MEMA24
5	MEMD06	22	MEMA07 <gg0></gg0>	39	MEMD13	56	MEMA25
6	MEMD07	23	MEMA06 <gr3></gr3>	40	MEMD14	57	NC (RFU)
7	XMEMCE1	24	MEMA05 <gr2></gr2>	41	MEMD15	58	XMEMREST
8	MEMA10 <gg3></gg3>	25	MEMA04 <gr1></gr1>	42	XMEMCE2	59	XWAIT
9	XMEM0E	26	MEMA03 <gr0></gr0>	43	NC (RFSH)	60	NC (RFU)
10	MEMA11 <gb0></gb0>	27	MEMA02 <xgvsync></xgvsync>	44	NC (RFU)	61	REG
11	MEMA09 <gg2></gg2>	28	MEMA01 <xghsync></xghsync>	45	NC (RFU)	62	XBVD2
12	MEMA08 <gg1></gg1>	29	MEMA00 <gdclk></gdclk>	46	MEMA17	63	XBVD1
13	MEMA13 <gb2></gb2>	30	MEMD00	47	MEMA18	64	MEMD08
14	MEMA14 <gb3></gb3>	31	MEMD01	48	MEMA19	65	MEMD09
15	XMEMWE (XMEMPGM)	32	MEMD02	49	MEMA20	66	MEMD10
16	MEMRDY (XMEMBSY)	33	XWP	50	MEMA21	67	XCD2
17	+5VD	34	DGND	51	+5VD	68	DGND

♦ OPERATION Connector



<Connector>

<Pin Assignment>

Model S30B-PHDSS

Pin No.	Signal Name	Pin No.	Signal Name
1	XRESET	16	+12VU
2	IRRXD	17	+12VU
3	IRBSY	18	+12VU
4	IRTXD	19	UGND
5	IRDTR	20	UGND
6	REMOTE	21	UGND
7	Not used	22	RETURN3
8	Not used	23	RETURN4
9	+5VD	24	RETURN5
10	+5VD	25	RETURN6
11	+5VS	26	SCAN0
12	Not used	27	SCAN1
13	DGND	28	SCAN2
14	DGND	29	SCAN3
15	+5VU	30	BLCNTL

♦ LED Connector

<Connector>

<Pin Assignment>

Pin No.	Signal Name	Pin No.	Signal Name
1	+12VU	8	LED06
2	+12VU	9	LED07
3	LED01	10	LED08
4	LED02	11	Not used
5	LED03	12	Not used
6	LED04	13	Not used
7	LED05	14	Not used

Model DF11-14DP-2DS



♦ LCD Connector

<Connector>

1

Model DF13-40PD-1.25V (59)

<Pin Assignment>

40

Pin No.	Signal Name	Pin No.	Signal Name		
1	DGND	21	DGND		
2	DGND 22 B0				
3	DCLK	23	B0		
4	XHSYNC	24	B0		
5	XVSYNC	25	B1		
6	DGND	26	B2		
7	R0	27	B3		
8	R0	28	DGND		
9	R0	29	DISPOUT		
10	R1	30	+5VD		
11	R2	31	+5VD		
12	R3	32	+5VD		
13	DGND	33	DGND		
14	G0	34	DGND		
15	G0	35	Not used		
16	G0	36	Not used		
17	G1	37	Not used		
18	G2	38	Not used		
19	G3	39 Not used			
20	DGND	40	Not used		

♦ JA Motherboard Connector



<Connector>

<Pin Assignment>

Model S22B-PHDSS

Pin No.	Signal Name	Pin No.	Signal Name
1	JRXIN	12	JCTLEN
2	MTXIN	13	JPCLKEN
3	MCTLIN	14	Not used
4	JNUMIN	15	VJA
5	JCNNTIN	16	VJA
6	MCLKIN	17	VJA
7	XJRESET	18	JGND
8	MRXOUT	19	JGND
9	JCTLOUT	20	JGND
10	MNUMOUT	21	+5VD
11	JTXEN	22	DGND

♦ Connectors on the Power Supply Unit

Power supply line connector





<Pin Assignment> Pin No. Signal Name Pin No. Signal Name 1 +5VD 7 Not used 2 +5VD 8 Not used 3 DGND 9 +24VC 4 DGND CGND 10 5 +12VU 11 +VJA6 UGND 12 JGND

Power supply status connector



<Connector>

<Pin Assignment>

Model S16B-PHDSS

Model 1-171826-2

Pin No.	Signal Name	Pin No.	Signal Name
1	Not used	9	XC24CNT
2	Not used	10	XC240C
3	Not used	11	XC24RES
4	ACFSP	12	SGND
5	XACFAIL	13	+5VS
6	Not used	14	SGND
7	XPWSTT	15	SNDSIG
8	XPWDWN	16	SNDGND

♦ BDM Connector

<Connector>

Model FFC-10BSM1



<Pin Assignment>

Pin No.	Signal Name	Pin No.	Signal Name
1	DS	2	BERR
3	DGND	4	BKPT/DSCLK
5	DGND	6	FREEZE
7	RESET	8	IFETCH/DSI
9	+5VD	10	IPIPE/DS0

♦ DEBUG Connector

Connector 1

<Connector>

<Pin Assignment>

Model TFM-140-12-S-D-LC

Pin No.	Signal Name	Pin No.	Signal Name	Pin No.	Signal Name	Pin No.	Signal Name	Pin No.	Signal Name
1	A0	17	A15	33	D5	49	XCS2	65	SIZE0
2	A1	18	DGND	34	D6	50	XCS3	66	SIZE1
3	A2	19	A16	35	D7	51	XCS4	67	XDSACK0
4	A3	20	A17	36	D8	52	DGND	68	XDSACK1
5	DGND	21	A18	37	D9	53	XCS5	69	XAS
6	A4	22	A19	38	D10	54	XCS6	70	XDS
7	A5	23	A20	39	D11	55	Not used	71	XIPIPE0
8	A6	24	A21	40	DGND	56	FREEZE	72	XIFETCH
9	A7	25	A22	41	DGND	57	Not used	73	XBERR
10	A8	26	A23	42	D12	58	Not used	74	Not used
11	A9	27	D0	43	D13	59	Not used	75	XRESETH
12	A10	28	D1	44	D14	60	Not used	76	DGND
13	A11	29	DGND	45	D15	61	Not used	77	XIFETCH
14	A12	30	D2	46	XBERR	62	CKL01	78	XBKPT
15	A13	31	D3	47	XCS0	63	DGND	79	XIPIPE0
16	A14	32	D4	48	XCS1	64	RXW	80	FREEZE
	Cor	nnector 1		80	79		Connector 2		1
>	>•••••••••••••••••••••••••••••••••••••								
79 80 2									

Connector 2

<Connector>

Model TFM-140-12-S-D-LC

	<pin assignment=""></pin>								
Pin No.	Signal Name	Pin No.	Signal Name	Pin No.	Signal Name	Pin No.	Signal Name	Pin No.	Signal Name
1	D16	17	D31	33	FC0	49	NC (ex11)	65	+5VD (VCC)
2	D17	18	DGND	34	FC1	50	NC (ex12)	66	+5VD (VCC)
3	D18	19	A24	35	FC2	51	NC (ex13)	67	+5VD (VCC)
4	D19	20	A25	36	FC3	52	DGND	68	+5VD (VCC)
5	DGND	21	A26	37	NC (ex1)	53	NC (ex14)	69	Not used
6	D20	22	A27	38	NC (ex2)	54	NC (ex15)	70	Not used
7	D21	23	A28	39	NC (ex3)	55	NC (ex16)	71	Not used
8	D22	24	A29	40	DGND	56	Not used	72	Not used
9	D23	25	A30	41	DGND	57	Not used	73	Not used
10	D24	26	A31	42	NC (ex4)	58	Not used	74	Not used
11	D25	27	XHALT	43	NC (ex5)	59	+5VD (VCC)	75	Not used
12	D26	28	XIPIPE1	44	NC (ex6)	60	+5VD (VCC)	76	DGND
13	D27	29	DGND	45	NC (ex7)	61	+5VD (VCC)	77	Not used
14	D28	30	XRESETH	46	NC (ex8)	62	+5VD (VCC)	78	Not used
15	D29	31	XAVEC	47	NC (ex9)	63	DGND	79	Not used
16	D30	32	XCS7	48	NC (ex10)	64	+5VD (VCC)	80	Not used