


HP Component Monitoring System Anesthesia/Standard/Neonatal



Notice

Hewlett-Packard makes no warranty of any kind with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. Hewlett-Packard shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.



Declaration

The M1165A/66A Component Monitoring System, and the plug-in modules used in this system that carry the CE mark, comply with the requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

©Copyright Hewlett-Packard GmbH. 1995.




INOP MESSAGES
or PATIENT NAME

YELLOW ALARM
or DATE + TIME

OWN BED LABEL

RED ALARM



ARRHYTHMIA
MESSAGES

To return to the main screen:

Press **Main Screen**.

To stop patient monitoring, but keep all settings and patient data information:

Press **Alarms/Volume** → **Monitor Standby**.

To resume patient monitoring:

Press any control panel key, except **Suspend**.

CMS/ACMS Main Screen

This User's Guide contains information on all the most commonly performed tasks for the Release F HP Component Monitoring System (CMS), HP Anesthesia Component Monitoring System (ACMS) and HP Neonatal Component Monitoring System (NCMS). Please refer to the Release F User's Reference Manual, Volume I and II for full details.

Keys on the control panel are all printed on a white background, for example

Procedures.

Keys on the screen are all printed on a shaded background, for example

Print Calc.

Cards inside the first row of tabs contain general information about the monitor.

Cards inside the second row of tabs contain information about specific parameters.

Where relevant, troubleshooting information for these parameters is also given.

Basic Troubleshooting

Problem	Possible Cause	Corrective Action
Wave is not on the screen	<p>Transducer is not plugged in. Wave is not switched on</p> <p>The monitor is not showing the display required.</p>	<p>Plug in transducer.</p> <p>Switch on wave:</p> <ol style="list-style-type: none"> 1. Press Module Setup → Parameter On/Off 2. Press Select Parameter (to highlight the parameter). 3. Press Parameter On/Off (to switch wave on) <p>Change to the required screen display:</p> <ol style="list-style-type: none"> 1. Press Change Screen. 2. Press Screen A, Screen B or Screen C.

Basic Troubleshooting (cont'd)

Problem	Possible Cause	Corrective Action
Wave is not on the screen	The wave has not been assigned to a channel.	<p>Assign wave to a channel:</p> <ol style="list-style-type: none"> 1. Press Monitor Setup → Display 1 Setup. 2. Press Select Channel (to highlight required channel). 3. Press Assign Wave (to highlight wave). <p>(Change speed and # of waves on screen or wave layout as required.)</p>
Numeric is not on the screen.	Parameter/Numeric is not switched on.	<p>Switch on Parameter/Numeric:</p> <ol style="list-style-type: none"> 1. Press Module Setup → Parameter On/Off. 2. Press Select Parameter (to highlight the required parameter). 3. Press Monitor Setup → Numeric 1 On/Off.

To turn ECG, Respiration and NBP on:

1. Press **Module Setup** → **Parameter On/Off**
2. Press **Select Parameter** (to highlight the required parameter).
3. Press **Parameter On/Off** (to turn the required parameter on).

To turn all other parameters on:

Plug the transducer into the parameter module.
(Parameter is turned on automatically.)

or when the transducer is unplugged

1. Press **Module Setup** → **Parameter On/Off**
2. Press **On/Off Setup** (to turn the required parameter on)

To enter a parameter setup window:

Press key on module which has parameter name (e.g. **ECG**)

or

1. Press **Module Setup**.
2. Press softkey with parameter name (e.g. **ECG**).

To change a setting in a Task Window:

Move the highlighting by pressing:

the softkey for the action you wish to perform
(e.g. **Change Lead**)

or

the highlighted arrow keys:



Parameter Setup Reminders



To switch Split Screen Trending on:

1. Press **Change Screen**
2. Press **Split Screen**
(surrounded by a highlighted frame, if selected)

To switch Split Screen Trending off:

1. Press **Change Screen**
2. Press **Split Screen**
(surrounded by a highlighted frame, if selected)

To Switch Split Screen Trending On or Off

Not available on the ACMS

To display the oxyCRG:

1. Press **Change Screen**
2. Press **oxyCRG Screen**
(surrounded by a highlighted frame, if selected)

To change oxygen channel (on the display and recording) and recorder speed:

1. Press **Monitor Setup** (Press **More Choices** if required)
2. Press **oxyCRG Setup**
 - Press **Change Channel** to change the oxygen channel for recording and displaying SpO₂ or tcpO₂.
 - Press **Change Speed** to change the recorder speed

To make an oxyCRG recording:

1. Press **Realtime Record**
2. Press **Preset Recording** (if available)
3. Press **oxyCRG PLUG-IN** to produce the recording

To configure oxyCRG to be the alarm recording:

1. Press **Monitor Setup**
2. Press **Recording Setup**. (Press **More Choices** if required)
3. Press **Delayed & Alarm**
4. Press **Change AlRecTyp** to highlight "oxyCRG"

To Adjust oxyCRG Measurements



Alarms/Volume - to turn *all* alarms off/on and *individual* alarms on/off, adjust alarm limits, show the alarm messages and go to standby mode.

Other Patients - to view data and alarms between beds.

Monitor Setup - to change settings such as screen displays and patient type, setup OxyCRG¹ and recording and make general configuration changes.

Procedures - to perform Cardiac Output and Wedge Pressure measurements, make ST segment monitoring adjustments, review ST wave segments, admit and discharge patients, or end a particular patient case and transfer patient data.

Trends/Calcs - to view data in graphs and tables, do physiological calculations in graphs, print reports and mark events to view in graphs.

Module Setup - to set up and change settings for all modules that are plugged into the CMS.

¹ Only available in Standard CMS.

Middle Row of Keys on the Control Panel

Silence/Reset - to silence an alarm and, if cause has been eliminated, reset it.

Suspend - to switch all alarms on/off.

Main Screen - to return to the standard monitoring screen.

Change Screen - to freeze waves on the screen. To change between screen layouts; combinations of OxyCRG¹ or the Split Screen Trending, (the selected screens are surrounded by a highlighted rectangular frame).

Realtime Record - to start a realtime preset recording or a realtime selected recording.

Delay Record - to start a delayed recording.

Airway Gases/Ventilation - to measure airway gases or show ventilator data.

Arrow Keys -     - to move highlighting and make selections in Task Windows - only active when lit.

1 Only available in Standard CMS.

Bottom Row of Keys on the Control Panel





A DELAYED recording records waves that occurred 10 to 15 seconds before you requested the recording.

A REALTIME recording records waves that occur at the time you request the recording.

To make a DELAYED recording (normally ECG):

(On Plug-In recorder ONLY, press **RUN/CONT**)

1. Press **Delayed Record**

To make a REALTIME (preset) recording:

1. Press **Realtime Record**
2. Press **Preset Recording** (if available).
3. Press **Mode A/B/C PLUG-IN**¹

4. Press **Stop Recording** to stop
(or **STOP** on recorder, if key available).

To make a REALTIME recording you select:

1. Press **Realtime Record**
2. Press **Select RecWaves** (if available)
3. Select the waves
4. Press **Start Recording**
5. Press **Stop Recording** to stop
(or **STOP** on recorder, if key available).

¹ The name of the actual recorder appears on the key. This can be PLUG-IN, Central, 4CH BS-1, or 4CH BS-2.

To Make Recordings



To turn alarm recordings on/off:

1. Press **Monitor Setup**
2. Press **Recording Setup** (Press **More Choices** if required)
3. Press **AlarmRec On/Off**
4. Configure alarm recordings for each parameter.

The availability of the following softkeys depends on the configuration of the monitor.

- Press **Select Parametr**
- Press **Change AlarmRec** to select sensitivity of alarm
- Press **All Rec On** to enable red or red and yellow conditions
- Press **All Rec Off** for no alarm recordings.

To change the recorder:

1. Press **Monitor Setup**
2. Press **Recording Setup**. Press **More Choices** if required)
3. Press **Procedur Recording**
4. Press **Change Recorder** to select the appropriate recorder

To Make Recordings (cont.)

To make a single Vital Signs recording:

1. Press **Realtime Record**.
2. Select **Preset Recordng**.
3. Select **RecVital PLUG-IN**.

To start/stop a sequence of Vital Signs recordings:

1. Press **Monitor Setup**.
2. Select **Recording Setup**, then **VitSigns Recording**.
3. To start:
 - a. Turn VitSigns recordings on with **Vitals Rec**.
 - b. Press **Start Rec On** to select the trigger source (NBP or Timer). If Timer, select **Repeat Time** to choose the frequency.
4. To stop:
 - a. Press **Vitals Rec** to select the Off position.

To Make Recordings (cont.)



Delayed

Realtime/Preset

Realtime/Your Selection

<p><i>To Start</i></p>	<p>▶ Delayed Record</p>	<p>▶ Realtime Record</p> <p>▶ Preset Recording (if available)</p> <p>▶ Mode X PLUG-IN¹</p>	<p>▶ Realtime Record</p> <p>▶ Select RecWaves (if available)</p> <p>▶ Select the wave(s).</p> <p>▶ Start Recording</p>
<p><i>To continue indefinitely</i></p>	<p>▶ RUN/CONT on recorder</p>	<p>▶ RUN/CONT on recorder if preset to timed</p>	<p>Realtime/Your Selection is preset to continued</p>
<p><i>To extend</i></p>	<p>▶ Delayed Record</p>	<p>▶ Mode X PLUG-IN¹</p>	<p>Realtime/Your Selection is preset to continued</p>

1 The name of the actual recorder appears on the key. This can be PLUG- IN, Central, 4CH BS-1, or 4CH BS-2.

Recordings Summary

Delayed

Realtime/Preset

Realtime/Your Selection

<p><i>To Stop</i></p>	<p>Stops automatically or</p> <ul style="list-style-type: none"> ▶ STOP on recorder <i>or</i> ▶ Record ▶ Preset Recordng (if available) ▶ Stop Recordng while highlighted. (Highlighted for local delayed recordings). 	<ul style="list-style-type: none"> ▶ STOP on recorder <i>or</i> ▶ Record ▶ Preset Recordng (if available) ▶ Stop Recordng while highlighted 	<ul style="list-style-type: none"> ▶ STOP on recorder <i>or</i> ▶ Record ▶ Preset Recordng (if available) ▶ Stop Recordng
<p><i>To Change Speed while Recording is running</i></p>	<ul style="list-style-type: none"> ▶ Realtime Record ▶ Preset Recordng (if available) ▶ Change RecSpeed 	<ul style="list-style-type: none"> ▶ Realtime Record ▶ Preset Recordng (if available) ▶ Change RecSpeed 	<ul style="list-style-type: none"> ▶ Realtime Record ▶ Preset Recordng (if available) ▶ Change RecSpeed

Recordings Summary (cont.)





To silence and reset alarms that are sounding,


Press **Silence/Reset**.

To suspend all alarms (turn alarms off):

Press **Suspend**

or

Press **Alarms/Volume** → **Suspend Alarms**

Either "Alarms Suspended" or "Alarms Suspended 1, 2 or 3Min" message, depending on how the monitor is configured; Large  appears on control panel.

To set/change alarm limits:

1. Press **Alarms/Volume** → **Alarm Limits**
2. Press **Select Parameter** (to highlight the required parameter)
3. Press **Low Limits↑** or **High Limits↑** (to adjust limits).

To turn individual parameter alarms on/off:

via alarm limits Task Window

1. Press **Alarms/Volume** → **Alarm Limits**.

Then either:

2. Press **Select Parameter** (to highlight required parameter).
 3. Press **On/Off Alarms** (to turn alarms on or off).
- or**
4. Press **Confirm** to switch all alarms off *then proceed as in 2. and 3. above to switch specific alarms on.*

via parameter Task Window

1. Press key on module.
2. Press **Adjust Alarms** → **On/Off Alarms** (to turn alarms on or off).

(If Off, an  appears next to the parameter numeric).

Alarms/Volume



To turn monitor alarms on:

Press **Suspend**.

or via Alarm Task Window

Press **Alarms/Volume** → **SwitchOn Alarms**

NOTE: Any individual alarms previously turned off will still be off (X next to numeric).

To view Current Alarms/INOPs:

Press **Alarms/Volume** → **Alarm Messages**.

To stop patient monitoring

(but keep all settings and patient data):

Press **Alarms/Volume** → **Monitor Standby**.

To resume patient monitoring:

Press any key on the control panel, except **Suspend**.

ACMS only

To Accept Ohmeda Ventilator Alarms:

1. Press **Airway Gases/Ventilation**
2. Press **VENT** to bring up the Ventilator Setup Task Window.
3. If the message **Ventilator alarms ignored** is displayed, Press the **Ventilatr Alarms** softkey.

To Ignore Ohmeda Ventilator Alarms:

1. Press **Airway Gases/Ventilation**
2. Press **VENT** to bring up the Ventilator Setup Task Window.
3. If the message **Ventilator alarms accepted** is displayed, Press the **Ventilatr Alarms** softkey.

To Accept or Ignore Ohmeda Ventilator Alarms





To view waveforms and vital signs on another patient:

Press **Other Patients** and move highlighting to line with bed or group label you want.

Press the softkey for that bed or group.

For groups, select bed, using **More Beds** if necessary.

To view other waves:

Press **Show NextWave**

To view another bed:

Press **Show Next Bed** or group label.

To send alarms to/receive alarms from another bed:

1. Press **Other Patients**

2. Press **Controls**

3. Press **Show Roster** to see which group the bed you're sending to/receiving from is in.

If needed, Press **Select Group** and highlight the desired group. Press **Send Alarms** and/or

Receive Alarms and highlight YES.

To Use Overview



To assign Arrhythmia:

1. Press either **Trends/Calcs** or **Alarms/Volume**
2. Press **arrhyth functns**

"Arrhyth assigned" on screen shows that you have been assigned arrhythmia. If arrhythmia has not been assigned to your bed and there are no channels available, you must UNASSIGN one of the other beds using the arrhythmia channel so as to assign arrhythmia to this bed. You cannot UNASSIGN the other bed from your monitor, you must do so at the Central Station or arrhythmia monitor. Assign your bed to the Central Station. (The QRS detection mode switches to "Automatic". The softkey **Auto/Man & Pacing** does not function).

To set up Arrhythmia:

1. Press **Paced Yes/No** to indicate the patient's pacer status.
2. The heartrate shown in any window is calculated by the arrhythmia computer, not the heartrate counter of the monitor.

To display the delayed Arrhythmia wave:

1. Press **Module Setup**
2. Press **Delayed ArrhWave**
(This key will only work if arrhythmia is assigned)

To turn off the autmatic pop up of the Arrhythmia Validation task window:

1. Press **Trends/Calcs**
2. Press **Rcv Arrhyth -> No**

Press **Rcv Arrhyth -> Yes** to turn pop up back on.





To admit or edit patient information:

1. Press **Procedures**
2. Press **Admit Patient** to admit new patient or edit existing patient data
3. If necessary, press **Select Line** to select the field
4. Enter or change patient information using the keypad
5. Press **Confirm**

To print patient information:

Press **Print Admit**

To Admit/Edit Patient Information

To end a case :

1. If arrhythmia is assigned through an HP Arrhythmia Monitor, press **Trends/Calcs** → **Arrhyth Functns** and unassign arrhythmia.
2. Press **Procedures**
3. Press **End Case**
4. Press **Confirm** to:
 - print a scheduled report (if previously selected)
 - discharge the patient
 - erase the database
 - reset all settings to user defaults
 - enter Monitor Standby

or press **Main Screen** to cancel.



or press **Confirm** again to cancel an ongoing report and continue with End Case.

5. If the monitor is connected to an HP Component Central Monitor, discharge the patient at the Central.

To discharge a Patient:

1. If arrhythmia is assigned through an HP Arrhythmia System, press **Trends/Calcs** → **Arrhyth Functns** and unassign arrhythmia.
2. Press **Monitoring Procedures**
3. Press **Discharg Patient**
4. Press **Confirm** to delete previous patient's data, or **Main Screen** to cancel.
5. If Monitor is connected to an HP Component Central Monitor, discharge the patient at the Central.





To perform calculations:

1. Press **Trends/Calcs**
2. Press **Hemo Calc**, **Ventil Calc** or **Oxygen Calc** to display the Calculations Task Window. (Displayed values may have been automatically or manually entered).
3. Optionally, press **Resample Vitals** and check that the values are representative of the patient state at the current time. If values are not representative press **Resample Vitals** again.

Press **Confirm** to store values without calculation.

4. Press **Perform Calc** to store any resampled values and perform calculation with displayed values.

To change or enter an input value:

- a. Press **▲** and **▼** to highlight value
- b. Enter new value using handheld keypad
- c. Press **Confirm**

To print the displayed calculations:

Press **Print Calc**

To view collection times of non-calculated values:

Press **On/Off DataTime**

To specify a calculation time:

1. Press **Change Time**
2. Press **◀** and **▶** to move highlighting to hours/minutes
3. Change hours/minutes if needed
4. Press **Confirm**
5. Press **Perform Calc.**

To view ranges of calculated values:

Press **On/Off Ranges**

To view previous calculations:

Press **Review Calc.**

To Perform Calculations

To review calculations:

1. Press **Trends/Calcs**
2. Press **Hemo Review**, **Ventil Review** or **Oxygen Review** to display Calculation Review Task Window.

To view older or newer data:

Press **◀** or **▶**

To modify input values used to perform calculations:

- a. highlight the time of the calculation to be modified using the **◀** and **▶** keys
- b. Press **Original Calc**
- c. Enter modifications, as needed, using the handheld keypad
- d. Press **Confirm**
- e. To recalculate, press **Perform Calc**

To show normal ranges for calculated parameters:

Press **On/Off Ranges**

To print the displayed calculations:

Press **Print Calc**

To print scheduled reports:

1. Press **Trends/Calcs**
2. Press **Print SchedRep**

To cancel a report:

Press **Cancel Report**

To mark an event:

1. Press **Trends/Calcs**
2. Press **Mark Event**
3. Press an event key (labelled A-D).
The message "Mark recorded" appears.

To review an event:

1. Press **Trends/Calcs**
2. Press **Graph Trends**
3. move the time bar along the Events line by pressing **◀** and **▶** until it rests on the Events marker you want to examine
4. Press **Zoom in Time**

To Mark and Review Events



To make drug calculations:

1. Press **Trends/Calcs**
2. Press **Drug Calc**
3. Press **Select Drug** or **▲** and **▼** to highlight desired drug
4. Press **Perform Calc** to access the Drug Calculation Task Window

To modify an input field:

1. Press **Select Item**
2. Enter the value using the keypad
3. Press **Confirm** after each entry.

To assign a field to be the unknown entity in the calculation:

1. Press **Choose Unknown**
2. Press **▶** and **◀** to select the unknown field

To calculate rate, dose, amount or volume when all inputs are complete:

Press **Perform Calc**

To display the Titration Table:

1. Press **Display Table**
2. Use the **Scroll** keys to view nondisplayed parts of the table.

To print the Titration Table:

Press **Print Table**





CAUTION:



- **Be sure that the patient in the monitor matches the one in the module before performing a transfer to the monitor.** To prevent accidental overwriting of patient data or name, the module will ask you to choose a name/MRN if the monitor name and module name don't match or if there is no name/MRN.
- **Be sure to discharge the previous patient at central before performing a transfer,** or the name moved into CMS from the module will be overwritten by the name from the central monitor.

Notes:

- After insertion, the Data Transfer Module requires approximately 8 seconds to synchronize with the monitor.
- If you are transferring data to the Component Transport System, erase any existing data in the transport monitor before performing the transfer.

To transfer quickly:

1. Plug the Data Transfer module into the rack.

2. Press  or  to select transfer direction.

3. Press **Confirm**.



Remove the module when the transfer is complete.

Quick Transfer



To Transfer to the Module:

1. Plug the Data Transfer module into the rack.

2. 
Press 

3. Select the type of transfer.

- **All Data:**

Press **Transfer All Data**

Any existing data in the module will be erased.

The monitor data will be transferred, followed by monitor updates once a minute.

- **Updates:**

Press **Collect New Data** to *send only monitor updates* once a minute starting at the current time.

Press **Clear Yes/No**

- *Yes* erases any data existing in the module before the transfer begins.

- *No* appends new data to data currently existing in the module.

If the patient name and ID do not match or if the monitor name is blank, the Combine in Module Task Window will be displayed.

Select the patient name to be used in the module.

4. Press  to start the transfer.




Blinking

- Rapid blinking of module light indicates a transfer in progress.
- Slow blinking indicates that the module is collecting new data.

To Stop a Transfer to the Module:

Unplug the Data Transfer module from the rack.


- or -

1. 
Press 
2. Press 

When the module is unplugged, it retains data for at least one hour.

Transfer to Module

To Transfer to the Monitor:

1. Plug the Data Transfer module into the rack.
2.  Press
3. Select whether or not to discharge the current patient data in the monitor before starting the transfer.
 - a. **If Names Match**, press **Confirm** to start the transfer without discharging the patient in the monitor.
 - b. **If Names Don't Match** (or either is blank) **and Discharge = No**, the Combine in Monitor Task Window will be displayed.
 - i. Select the patient name to be used in the monitor.
 - ii. Press **Confirm** to start the transfer and incorporate the module data into the monitor's existing patient record.
 - c. **If Discharge = Yes**, the Discharge Task Window will be displayed.
 - i. Press **Confirm** to discharge the monitor patient first and then start the transfer of data from the module.

To Stop a Transfer to the Monitor:

- Unplug the Data Transfer module from the rack.
- After a transfer is completed, the data in the module is erased. If the transfer is interrupted, the module data will remain intact.

Transfer to Monitor





To adjust ECG measurements:

Press **ECG** on module

To select the required lead:

Press **Change Lead**

To adjust the size of the QRS complex:

Press **Adjust Size**

To select the filter setting:

Press **Filter Mon/Diag**

FILTER -use only if excessive electrical interference

(can impede ECG analysis)

MONITOR - use for normal monitoring

DIAG - use when diagnostic quality required

(may increase false alarms).

To view multi-lead ECG:

Press **Multi-Ld ECG/ST** → **View ECG**

To record multi-lead ECG:

Press **Multi-Ld ECG/ST** → **Record ECG**

To select the QRS detection method:

Press **Auto/Man & Pacing**

PACED or **NONPACED** / **MANUAL** or **AUTO**

Use **AUTO** except in rare cases where the monitor miscounts - if change to **MANUAL**, adjust detection level so only QRS complexes cross dotted threshold line.

To turn off channels 2 and 3:

1. Press **Setup Next ECG**
2. Press **On/Off ECG-CH2** to highlight **off**
(Channel 3 is turned off automatically).

To Adjust ECG Measurements



To Adjust Respiration Measurements:

Press **RESP** on module

To adjust the size of the RESP wave:

Press **Adjust Size**


To select RESP detection method:

Press **Auto/Manual**

AUTO - monitor counts respiration and adjusts detection level automatically

MANUAL - you adjust the detection level by moving the dotted threshold line just below peak of curve) using the **▲** and **▼** keys.

To Adjust Respiration Measurements



The ST measurement for each beat complex is the vertical difference between the points where the wave intersects two measurement points:

- The **isoelectric point** provides the baseline for the measurement
- The **ST point** provides the other measurement point

You can adjust both measurement points. Your monitor is configured so that you set the ST Point in reference to:

- *either* the R-wave peak
- *or* the J Point, which is the transition between the QRS and the ST segment. The ST point is referenced from the j point at either +60 or +80 msec.

To Adjust the ST Measurement Points



To setup for ST monitoring:

1. Ensure that the ECG channel(s) to be used for ST measurement are turned on. (Channel does not have to be displayed).
2. If necessary, turn ST on in the Parameters On/Off Task Window.
3. Press **ECG** on the module (or **Procedures**), then **ST Analysis** to display the ST Analysis Task Window.
4. If necessary, turn individual ST channels on using **On/Off ST**. Then press **ST Analysis** to return to ST Analysis Task Window.

To adjust the measurement points: *(from the ST Analysis Task Window)*

1. Press **Adjust Meas Pts**
2. Press **On/Off Freeze** to freeze the waveform
3. adjust the isoelectric point and J/ST points.
4. Press ST Analysis to return to ST Analysis Task Window.

To change the reference wave:

(from the ST Analysis Task Window)

1. Move the cursor to the point of interest in the trend. (This displays the trended beat in the cursor recall spot)
2. Press **Temp Reference** to display the beat in the reference spot.

To make the temporary reference beat permanent,

Press **Change Reference** → **Confirm**

To redisplay the original saved reference beat,

Press **Saved Reference**

3. Press **ST Analysis** to return to the ST Analysis Task Window.

To adjust ST alarms: (from the ST Analysis Task Window)

1. Press **Adjust Alarms**
2. Press **Select Channel**

To adjust the alarm limits:

Press **Low Limit** or **High Limit** softkeys

To turn alarms on or off:

Press **On/Off Alarms** .

Press **ST Analysis** to return to ST Analysis Task Window.

To Adjust ST Segment Monitoring (cont.)



Problem	Possible Cause	Corrective Action
<p>No display of ECG</p> <p>LEADS OFF message. ? instead of HR/RESP numerics. INOP tone.</p> <p>ECG waveform is present but the RESP channel shows a flat line.</p>	<p>System not plugged in or switched on. ECG or ECG/RESP module is not plugged in. Brightness/contrast not adjusted. The video cable is not plugged in.</p> <p>ECG cable is not connected. One or more electrodes are detached. No gel on electrodes. Strong solvent used for cleaning skin. Extremely large pace pulses.</p> <p>One of the RESP electrodes is loose. No gel on the electrodes. Patient's breathing is too shallow to be detected in auto mode.</p>	<p>Connect and switch on system. Check ECG or ECG/RESP Module is plugged in correctly. Check that brightness and contrast controls are correctly adjusted. Check the video cable connection.</p> <p>Check ECG cable is connected. Check that electrodes are not detached. Check that electrodes have enough gel. Change to a suitable skin cleaning solvent. Ensure that channels 2 and 3 are switched off if a 3-lead set is used. Reduce size or select lead with lower pacer amplitude.</p> <p>Check that RESP electrodes are not loose. Check that electrodes have enough gel. Reposition electrodes to try and optimize signal. Change RESP detection mode from auto to manual if necessary.</p>

Problem	Possible Cause	Corrective Action
<p>NO PULSE SOURCE message (only if PULSE is the active alarm parameter).</p> <p>Lots of noise on the ECG wave signal.</p>	<p>The pulse source (PRESS or SpO₂/PLETH Module) is unplugged, or the transducer is not connected, or equipment malfunction.</p> <p>Bad electrode contact, high frequency disturbance.</p>	<p>Check that module is plugged in. Check that transducer is connected.</p> <p>Remove noise source.</p> <p>Increase the filter on the waveform:</p> <p>DIAG - No filter; all interference left in</p> <p>MONITOR - Small filter; usually enough to clear ECG</p> <p>FILTER - Stronger filter; especially good in conditions of excessive muscle movement.</p>

ECG Troubleshooting (cont.)





To adjust pressure measurements:

Press **PRESS** on module

To change the size of the waveform:

Press **Change Scale**

(to automatically optimize the scale to suit the amplitude of the incoming wave, select "optimum" → **confirm**).

To change the pressure label:

Press **Change Label**

IMPORTANT: For information on Pressure labeling, refer to the HP Component Monitoring System User's Reference Manual Volume 2.

To adjust the alarm limit:

1. Press **Adjust Alarms**
2. change alarm parameter (sys, dia, mean), if needed.
3. Press **[Low Limit]** or **[High Limit]**
(to adjust alarm limits).

To zero transducer:

from the module

1. open transducer to air (closing it to the patient)
2. Press **ZERO** on module for one second - you should see the message "Zero in Progress"

from task window

1. Press **Zero Xducer**
2. open transducer to air (closing it to the patient)
3. Press **Zero Xducer** again -

Wait for screen message saying that procedure is complete, then close transducer to air and open it to the patient.

To Adjust Pressure Measurements



Remember: In order to do a Wedge measurement, the pressure label **must** be PAP.

To measure Pulmonary Artery Wedge Pressure:

Press **Procedures** → **Wedge**

To measure wedge pressure:

1. start the Wedge procedure
Wave freezes 12 seconds after wedge detected.
2. Press **Stop Trace** to freeze trace at any time.
3. If satisfied with the trace, deflate the balloon

To adjust where the monitor is taking the reading from:

Press **Edit Wedge** and move the cursor up and down using the **↓ Cursor ↑** keys.

To record the curve:

Press **Record Curve**

To store in patient data:

Press **Confirm**



Note: ssCO_2 cannot be used as a reference waveform.


To change the pressure arterial source:

1. Press **Module Setup** → **CPP**
2. Press **Change Source** (to select arterial source).

Remember: In order to measure CPP, both ICP and Mean Arterial Pressure must be monitored.

To Adjust Cerebral Perfusion Pressure Measurements

Problem	Possible Cause	Corrective Action
<p>No pressure wave on screen</p> <p><i>Noise or fling in the pressure waveform</i></p>  <p><i>Resonating Arterial Waveform</i></p>  <p><i>or Abnormally High Systolic Pressure</i></p>	<p>Selected pressure scale is incorrect.</p> <p>Movement of P.A. catheter tip. Motion of pressure tubing.</p> <p>Small air bubbles in, or very close to the transducer.</p> <p>Tubing too long.</p>	<p>Change the size of the waveform by re-optimizing or changing pressure scale.</p> <p>Reposition catheter (in accordance with your Hospital Procedures Policy). Make sure patient or ventilator tubing is not touching pressure tubing.</p> <p>Carefully flush the transducer and tubing when setting up system. Debubble carefully after each flush.</p> <p>Shorten tubing.</p>

Problem	Possible Cause	Corrective Action
<p><i>Dampened Arterial Waveform</i></p>  <p><i>Abnormally Low Pressure</i></p> <p><i>No Pressure</i></p>	<p>Thrombus formation or blood left in catheter following blood sampling.</p> <p>Large air bubble in tubing.</p> <p>Catheter tip against the vessel wall. Catheter kinking or arteriospasm.</p> <p>Tubing too long or too compliant.</p> <p>Improper sequence of stopcock operation.</p> <p>Defective transducer and/or amplifier.</p> <p>Transducer level higher than heart.</p> <p>Loose connections.</p> <p>If monitor turned on and there is no pressure. The transducer or connector wire could be broken.</p>	<p>Use syringe to withdraw air or particles in catheter, then flush the line with fresh solution.</p> <p>Use syringe to withdraw air in tubing, then flush line with fresh solution.</p> <p>Reposition the catheter to relieve the spasm.</p> <p>Shorten the tubing (recommended length 3-4 ft. or replace with large diameter stiff tubing.</p> <p>Flush line, re- zero and recalibrate.</p> <p>Replace transducer and have the monitor checked.</p> <p>Check patient and transducer positions.</p> <p>Check and tighten connections.</p> <p>Replace transducer and cable. Check stopcocks.</p>

Pressure Troubleshooting (cont.)



To adjust NBP measurements:

Press **(NBP)** on module

To change mode:

Press **Auto/Manual**

To change repetition time (AUTO mode only):

Press **Change Rep Time**

To do a single manual measurement or initiate the AUTO cycle:

Press **Start NBP**

or

(START) on module

To do a STAT measurement (as many times as possible in a 5-minute period):

Press **Stat NBP**

or

(STAT) on module

To STOP a measurement:

Press **Stop NBP**

or

(STOP) on module

Note: Use appropriate cuff size for NBP measurements.

To Adjust NBP Measurements



The following factors can make the NBP measurement unreliable or impossible to detect

Patient movement: Movement, shivering and convulsions may all interfere with the measurement of the arterial pressure pulses. Measurement time will also be increased.

Heart-lung machine: Measurements will not be possible if the patient is connected to a heart-lung machine.

Severe shock: Reduced blood flow to the peripheries will cause reduced pulsation to the arteries in conditions of severe shock or hypothermia. This will make measurements unreliable.

Incorrect cuff size: Using the incorrect cuff size will give erroneous readings.

Cardiac arrhythmias: An irregular heart beat will make the measurement unreliable or impossible. The measurement time will be increased.

Pressure changes: Rapid changes in blood pressure, while the arterial pressure pulses are being analyzed, may make the measurement unreliable or impossible.

Heart rate extremes: Measurements can not be made at a heart rate of less than 40 bpm or greater than 300 bpm.

Obese patients: A thick layer of fat surrounding a limb tends to dampen oscillations coming from the artery, and stops them from reaching the cuff. Accuracy may be lower than normal.





To adjust SpO₂/Pleth measurements:

Press **SpO₂** on module

To switch Tone Modulation on/off:

Press **On/Off Tone Mod**

To adjust the QRS volume from the SpO₂ window:

Press **Volume Control** → **QRS Tone ↑ ↓**

To select Pleth or one of the pressures as the pulse source:

1. Press **Module Setup** → **HR/PULSE**
2. Press **Change Source** (to highlight pulse source)

Setting the high SpO₂ alarm limit to 100% is the equivalent to switching off the high alarm. High oxygen levels may predispose a premature infant to retrolental fibroplasia. Therefore the upper alarm limit for oxygen saturation must be carefully selected in accordance with accepted clinical practices.

To Adjust SpO₂/Pleth Measurements



The following factors can cause errors in the measurement of the SpO₂

Reduced arterial blood flow can be caused by:

- shock
- hypothermia
- use of vasoconstricting drugs
- anemia.

A "SpO₂ NON PULSATILE" message appears on the display, if the pulsations are undetectable.

High levels of ambient light

An "SpO₂ LIGHT INTERF" message will appear on the display. The measurement quality can be improved by covering the transducer with a suitable non see-through material.

Absorption of light by other chemicals in the blood

for example:

- carboxyhemoglobin
- methemoglobin
- methylene blue
- indocyanine green
- indigocarmine.

Problem	Possible Cause	Corrective Action
Light interference	Surgical lamps, bilirubin lights, overhead fluorescent lights, infra-red heater/lamp, direct sunlight.	The sensor site must be covered with opaque, non-reflecting material.
Erratic readings, or screen messages: SpO₂ non pulsatile or Weak Signals	<p>Using inappropriate sensor for patient or clinical condition.</p> <p>Sensor applied to inappropriate site.</p> <p>Sensor applied incorrectly.</p> <p>Poor sensor adhesion.</p> <p>Using too many adhesive dots.</p> <p>Tissue between LED and detector too thick.</p>	<p>Select correct sensor based on patient weight, application site, level of activity.</p> <p>Apply sensor to recommended site.</p> <p>Apply sensor to recommended site, check LED/detector are lined up, check sensor still sticks.</p> <p>Use new sensor if necessary. Do not use additional tape.</p> <p>Do not use more than 3 dots for each LED and detector.</p> <p>Do not use dots with sensors for which they were not intended.</p> <p>Apply sensor to thinner tissue site.</p>

SpO₂ Troubleshooting (cont.)

Problem	Possible Cause	Corrective Action
Different values on different extremities	Cardiac shunts	Check with blood sample. Accuracy of value questionable, depending on shunt size.
No SpO₂ numeric, or Weak Signal message	<p>Reduced arterial flow due to shock, hypothermia, vasoconstricting drugs.</p> <p>Nail polish</p> <p>Hematoma</p>	<p>Change transducer site. In neonates try moving sensor to thinner part of extremity, for instance, towards toes.</p> <p>Remove nail polish.</p> <p>Change transducer site</p>
SpO₂ waveform looks flattened and monitor is questioning the reading	<p>Poor application of probe</p> <p>Inappropriate waveform size</p>	<p>Check probe application</p> <p>Adjust waveform size (Perfusion Mode only): Press Pleth on front of module → Optimize Size</p>





Attaching the sensor:

After calibrating the sensor, attach it to either:

mainstream - the patient's breathing circuit

sidestream - the sidestream module

The following waveform represents a normal capnogram:



To adjust CO₂ measurements:

Press **CO₂** on module

To change size of waveform,

Press **Change Scale**

To set anesthetic gas ON (if it contains N₂O)
or OFF (if it doesn't),

Press **N₂O Corr On/Off**

To adjust alarm limits,

Press **ETCO₂ Alarms**, **AWRR Alarms** or
IMCO₂ Alarms, and use **[Low Limit↑]** or
[High Limit↓] softkeys.

- ETCO₂: set low and high limits
- IMCO₂: set high limit
- AWRR: set apnea delay time, low and high limit.

To Adjust CO₂ Measurements



For accurate measurements, allow sensor to warm up until "CO₂ Warmup" message disappears (up to 20 minutes if it has not been plugged into a module).

When to perform an accuracy check . . .

- Once a day
- If module is moved with sensor plugged in to another monitor or used with another patient.

When to calibrate . . .

- If a new (or different) sensor is attached to the module
- When the *Accuracy Check* values displayed do not coincide with the calstick value
- Once every 1-2 weeks.

To perform an accuracy check

(this should be done once a day)

Press **CAL** on module, and check calstick value on the screen.

If calstick and screen values are the same ...

1. Place sensor on cell labelled 0.0 mmHg.
Reading on display should be within +/- 0.3mmHg within 1 minute.
Place sensor on other cell.
Reading on display should be +/- 1mmHg within 1 minute.
If the values match, you don't need to calibrate; Press

Main Screen.

2. If values don't match, you need to calibrate. Press **Start Calibr.**
3. Go to **Perform the calibration.**

If calstick and screen values are different ...

1. Press **Change Cal Value** to set the correct value.
2. Press **Confirm**
3. go to **Perform the calibration.**

To calibrate CO₂:

1. Place sensor on calstick cell.
Press **Start Calibr.**
2. Repeat process with other cell.
(For both cells, values on screen should match values on calstick).
3. Press **Main Screen**

A more detailed version of this procedure can be found in the HP Component Monitoring System User's Reference Manual Volume 2.

To Calibrate CO₂

To turn sidestream pump on (if pump turns off automatically):

Press **CO₂** on module

Press **Restart Sidestrm** key:

ONCE if key is already highlighted,

TWICE if key is not highlighted.

To Turn Sidestream Pump On

Problem	Possible Cause	Corrective Action
CO₂ wave line disappears	Module not plugged in, transducer not connected, Equipment malfunction, sidestream pump off.	Check that module and accessories are plugged in. Sidestream - check pump switched on.
CO₂ wave set to zero, other numerics displayed as -?-	Not in monitoring mode, calibration taking place or failed, barometric pressure out of range.	Change filter (may be blocked), check sample tube (may be kinked). Check whether monitor in calibration mode. Check that transducer is placed on correct cal cell. Call service engineer to check barometric pressure.
Wave Display Clipped	Incorrect scale selected - complete wave display cannot be seen	Change wave scale size, by pressing (CO₂) on module then Change Scale .

CO₂ and Sidestream Mode Troubleshooting



FiO₂ Stabilization Times:

No more than 1 min. is needed if sensor is left connected to module for more than 15 min.(power off up to 4 hr.)

15 min. is needed if disconnected from module for less than 16 hours.

When to calibrate FiO₂:

The FiO₂ transducer should be calibrated:

every 8 hours

when sensor is moved to another module

when sensor is moved to another monitor.


To perform a 21% calibration:

1. remove sensor from T-Piece and expose to room air.

2. Press **CAL** on module for at least one second.

(For more information on FiO₂ stabilization, the 21% or 100% calibration procedure, please refer to the User's Reference Manual Volume 2.)





Problem	Possible Cause	Corrective Action
Readings constantly drift	Transducer improperly applied. Air leaks. Old electrolyte not fully removed during remembraning procedure.	Reapply fixation ring. Reattach transducer. Place transducer in cal chamber and turn on gas. Do not calibrate. If transducer does not give stable reading after 5 minutes, remembrane, taking care to clean out old electrolyte. Check groove. Calibrate.
Readings not stable 20 minutes after application.	Patient status unstable Inadequate vasodilation	None Increase transducer temperature if not contra-indicated by patient's condition. (Reduce the transducer site time if appropriate).
tcpO₂ values too low or too high	Capillary samples taken (Astrups, heelsticks) Blood taken while patient agitated (crying, feeding, suctioning, positioning).	Take arterial samples. Capillary samples do not correlate well with paO ₂ . Take another sample when patient is quiet.

tcpO₂/tcpCO₂ Troubleshooting (cont.)

Problem	Possible Cause	Corrective Action
tcpO₂ values too low	Patient has severe respiratory distress	No corrective action possible. tcpO ₂ levels will underestimate paO ₂ .
tcpO₂ values too high	Incorrect handling of blood sample Transducer improperly applied - air leaks.	Take another sample. Ensure that atmospheric oxygen does not enter sample. Reapply fixation ring. Reattach transducer.
tcpCO₂ values too high	Correction factor not ON. Hypoperfusion Shock Leakage in calibration tubing, or non-original tubing used	<i>Either</i> ask biomed or service to configure correction ON, <i>or</i> take the rise of pCO ₂ due to heating and metabolic CO ₂ production into account. None possible None possible Replace tubing. Use correct tubing (M2205A).
tcpCO₂ values too low	Transducer improperly applied - air leaks.	Reapply fixation ring. Reattach transducer.

When tcpO₂/CO₂ Values Do Not Match Blood Gas Values





To adjust Cardiac Output measurements:

Press **C.O.** on module

To adjust the computation constant:

Press **Adjust Constant** or **↓ Cursor ↑** keys
(see catheter insert for comp. constant appropriate to catheter type, injectate temperature and volume)

Press **Confirm**

To measure Cardiac Output:

1. Press **Measure C.O.**
2. Press **Start C.O.** or **START** on module
3. when "Inject Now" appears, inject.
(For accuracy start within 15 seconds - complete in 4 seconds).
4. Repeat Step 2 for each injection. Wait at least one minute (or longer, depending upon patient's clinical condition) between injections.
5. Press **Record Curve** (if you want to record the curve)
6. Press **Edit C.O.** (when enough readings have been taken).

To delete any unwanted values:

1. Press **Select Curve** and highlight curve to delete
2. Press **Delete Yes/No**

To store average in patient data:

Press **Confirm**

To perform hemodynamic calculations:

Press **Hemo Calc**

To Make Cardiac Output Measurements



The following factors can cause errors in the measurement of the Cardiac Output

Physiological reasons:

Patient movement during the procedure
Anxious patient
Variations in cardiac rate and rhythm
Cardiac abnormalities, for example,
incompetant valves
Shock.

Catheter related errors:

Balloon inflated during measurement
Catheter not positioned properly
Damaged catheter.

Injection errors:

Use of the wrong catheter injection port
Poor timing of injection
Incorrect volume of injectate
Inaccurate injectate temperature.

Instrument errors:

Incorrect computation constant
Instrument failure.

Screen Message	Possible Cause	Corrective Action
1. Noisy Baseline	Mechanical Ventilation Variations in venous return (coughing, shivering, abnormal respiration)	Synchronise injections to end-expiration Maximize signal amplitude by: 1. Increasing injectate volume to 10cc (if not contra- indicated) 2. Using iced injectate (change comp constant) 3. Increasing number of measurements to be averaged (if not contra-indicated)
2. Thermal Baseline Drift	Administration of I.V. fluids; Rewarming after hypothermia, fever, shock	Try to stabilize patient temperature; Maximize signal amplitude (as above); Increase number of measurements to be averaged.
3. Unsteady Baseline	Combination of messages 1 and 2	see above
4. Multiple Peaks	Injection not smooth Cardiac shunt	Repeat smoothly and rapidly. Check with blood samples. Accuracy of values questionable, depending on shunt size.
5. Abnormal Decay Time	Low Cardiac Output	Maximize signal output

Cardiac Output Troubleshooting (cont.)

Screen Message	Possible Cause	Corrective Action
6. Very Short Curve	Injection too fast	Inject again as soon as possible after the prompt.
7. Very Long Curve	Injection too slow	Inject again as soon as possible after the prompt.
8. Irregular Curve	Combination of messages 4-7	See above.
9. Delayed injection	Low stroke volume, Low flow rate Delayed injection Kinked catheter	Inject as soon as possible after prompt. Low flow rates need 15-20 seconds to pass by thermistor. Repeat Check for kinks below dressing. Check chest x-rays for interior kinks.
10. Injectate temperature too high	Incorrect positioning of injectate probe Hypothermic patient <34°C, 93.2°F Room temperature injectate Room temperature > 25°C (77.8°F)	Make sure that probe is positioned in a solution identical to the one injected. Use iced injectate with hypothermic patients. (Remember to change comp. constant.) Cool injectate. (Remember to change computation constant.)
11. Excessive Curve Height	Injectate too cold	Heat up injectate and start measurement again
12. Curve Below Baseline	Rewarming after hypothermia, fever, shock	Try to stabilize patient temperature and start measurement again.

To adjust temperature measurements:

Press **Temp** on module

To select temperature label:

Press **Change Label**


To measure temperature differences when monitoring more than one temperature:

1. Press **Module Setup**
2. Press **Diff**
3. Press **Select First**
(to select the first temperature)
4. Press **Select Second**
(to select the second temperature)

IMPORTANT: For information on the Temperature labeling, refer to the HP Component Monitoring System User's Reference Manual Volume 2.

To Adjust Temperature Measurements





Problem	Possible Cause	Corrective Action
SvO ₂ CONFIGURATION message appears.	Computation mode was set to SaO ₂ when the optical module was connected to an Abbott monitor.	Enter SvO ₂ Task Window and return to SvO ₂ mode, by pressing Confirm .
SvO ₂ CONNCT OPTMOD message appears.	The optical module was disconnected during storage of data.	Reconnect optical module to SvO ₂ module until storage of data is complete.
SvO ₂ CAL REQUIRED message appears.	No valid calibration data, usually due to connection of a new optical module.	If the catheter is still in the tray, perform a pre-insertion calibration. If the catheter is already inserted into the patient, perform an in-vivo calibration.
SvO ₂ UNABL TO MEAS message appears (No SvO ₂ value calculated).	Signal received is outside the normal range.	Perform an In-vivo calibration. If problems persists, contact Abbott.

SvO₂ Troubleshooting



Problem	Possible Cause	Corrective Action
<p>SvO₂ LOW LIGHT message appears. (Unusually low intensity of light received from the catheter).</p>	<p>If the catheter is still in the tray:- Catheter is not properly connected.</p> <p>The catheter or optical module is faulty.</p> <p>Catheter has been removed from optical reference, but not yet inserted into patient.</p>	<p>Check connections and ensure that the catheter tip is fully inserted into the optical reference.</p> <p>Disconnect the catheter and position optical module out of direct light:-</p> <ul style="list-style-type: none"> ● If INOP message disappears, the catheter is faulty. Manipulate or use a new catheter. ● If INOP message persists, the optical module is faulty. Check optical module by substituting with a replacement module. <p>Insert catheter into patient.</p>
<p>SvO₂ LIGHT INTENS message appears (Change in light intensity).</p>	<p>Catheter tip against blood vessel wall.</p> <p>Inadequate blood flow past the catheter tip.</p> <p>Damaged optical fibers.</p>	<p>Ensure distal lumen is patent. Reposition the catheter.</p> <p>Straighten any kinks in the catheter. Check connections at the SvO₂ and optical modules.</p>





Only for use with the ACMS

To set up the Anesthetic Gas Module:

1. Press the power-on switch
2. Fit the bacterial filter (1) to the gas sample tube (2).
3. Connect the gas sample tube to the Anesthetic Gas Module:
 - a. Move the cover (4) to reveal the inlet (3).
 - b. Insert the tube into the slot
 - c. In a clockwise direction, rotate the tube through 90°.
4. Wait until the Anesthetic Gas Module has finished its self test and has warmed up.

5. Connect the other end of the gas sample tube to the airway adapter (5).
6. Connect the airway adapter to the patient.

When the module is ready for use, the waveforms and numerics are shown on the Main Screen.

Setting up the Anesthetic Gas Module





Only for use with the ACMS

To set up the Anesthetic Gas Module:

1. Press the power-on switch
2. Fit the bacterial filter (1) to the gas sample tube (2).
3. Connect the gas sample tube to the Anesthetic Gas Module:
 - a. Move the cover (4) to reveal the inlet (3).
 - b. Insert the tube into the slot
 - c. In a clockwise direction, rotate the tube through 90°.
4. Wait until the Anesthetic Gas Module has finished its self test and has warmed up.

5. Connect the other end of the gas sample tube to the airway adapter (5).
6. Connect the airway adapter to the patient.

When the module is ready for use, the waveforms and numerics are shown on the Main Screen.

Setting up the Anesthetic Gas Module

To manually select the anesthetic agent:

1. Press the **Airway Gases** set-up key on the module to bring up the Overview Task Window, or the **Airway Gases/Ventilation** hardkey.
2. Press **Airway Agent** softkey
3. Press **Select Agent** softkey
4. Use **←** or **→** to select agent.

The measured agents are as follows:

Anesthetic Agents	Display Abbreviation
Halothane	HAL
Isoflurane	ISO
Enflurane	ENF
Sevoflurane	SEV
Desflurane	DES

5. If no agent is to be measured select "None".
6. Press **Main Screen** when finished.

To Manually Select the Anesthetic Gas

To adjust the alarms:

1. Press the **Airway Gases** set-up key on the module to bring up the Overview Task Window, or the **Airway Gases/Ventilation** hardkey.
2. Press the softkey corresponding to the gas to adjust (for example, **Airway CO2**)
3. Select the **et** or **in** value as required, using the arrow keys.
4. Press **CO2 Alarms**
5. Use **↓ Low Limit ↑** or **↓ High Limit ↑** softkeys to adjust alarm limits.
6. Press **Main Screen** when finished.

To Adjust Anesthetic Gas Module Alarms



To Select Automatic Agent Identification:

This feature automatically detects a set of anesthetic gases.
(This feature is optional.)

1. Press the **Airway Gases** set-up key on the module to bring up the Overview Task Window, or the **Airway Gases/Ventilation** hardkey.
2. Press **Airway Agent**
3. Press **Select Agent**

4. Select **Auto** for automatic agent identification

NOTE: **Auto** is only available for selection if you have the agent identification option

5. If no agent is to be measured, select **None**
6. Press **Main Screen** when finished.

Calibration


Zero Calibration:

Done automatically when the module is switched on at the following intervals: 8, 15, 30 and 90 minutes, then every 8 hours. The calibration takes 15 seconds. During zero calibration, the patient's AWRR is not measured.

Span Calibration:

The Anesthetic Gas Module should be checked at least once per year to determine if span calibration is necessary.







HEWLETT®
PACKARD

M1046-9041F

Copyright © 1995
Hewlett-Packard Company
Printed in U.S.A. 12/95

