

User's Guide

# ECG Interpretation Program

ECAPS 12C

For Bedside Monitors

If you have any comments or suggestions on this manual, please contact us at: [www.nihonkohden.com](http://www.nihonkohden.com)

---

**Copyright Notice**

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved. No part of this document may be reproduced, stored, or transmitted in any form or by any means (electronic, mechanical, photocopied, recorded, or otherwise) without the prior written permission of Nihon Kohden.

# Contents

<b>Section 1</b>	<b>Introduction .....</b>	<b>1.1</b>	<b>1</b>
<b>Section 2</b>	<b>Precautions (Discrepancies between Physician's and ECAPS 12C's Findings) .....</b>	<b>2.1</b>	<b>2</b>
<b>Section 3</b>	<b>Outline of ECG Automatic Recording .....</b>	<b>3.1</b>	<b>3</b>
	Entering Patient Data .....	3.3	<b>4</b>
	Recording ECG Data .....	3.4	<b>5</b>
	Improving Waveform Quality .....	3.5	<b>6</b>
	Measuring ECG Waveform .....	3.6	<b>6-1</b>
	General Measurement Values to Be Printed .....	3.10	<b>6-2</b>
	Classification and Printing ECG Data .....	3.11	<b>6-3</b>
<b>Section 4</b>	<b>Recorded Analysis Results .....</b>	<b>4.1</b>	<b>6-4</b>
<b>Section 5</b>	<b>How to Read Analysis Results .....</b>	<b>5.1</b>	<b>6-5</b>
	ECG Findings .....	5.3	<b>6-6</b>
	Priority of ECG Findings .....	5.3	<b>A</b>
	Criteria .....	5.4	
	Overall Judgement .....	5.4	
	Analysis Result When ECG Cannot Be Analyzed .....	5.5	
<b>Section 6</b>	<b>Criteria of Findings .....</b>	<b>6.0.1</b>	
	Indication of Analysis Criteria .....	6.0.2	
<b>Section 6-1</b>	<b>Arrhythmias .....</b>	<b>6.1.1</b>	
	Rhythm analysis .....	6.1.2	
	1002 Marked rhythm irregularity, possible non-conducted PAC, SA block, AV block, or sinus pause .....	6.1.13	
	1100 Sinus rhythm .....	6.1.3	
	1102 Sinus arrhythmia .....	6.1.12	
	1108 Marked sinus arrhythmia .....	6.1.13	
	1120 Sinus tachycardia .....	6.1.3	
	1130 Sinus bradycardia .....	6.1.3	
	1200 Atrial rhythm .....	6.1.4	
	1210 Atrial fibrillation .....	6.1.7	
	12101 Atrial fibrillation with rapid ventricular response .....	6.1.7	
	12102 Atrial fibrillation with slow ventricular response .....	6.1.7	
	12103 Atrial fibrillation with aberrant conduction, or ventricular premature complexes .....	6.1.18	

12108	Atrial fibrillation with rapid ventricular response with aberrant conduction, or ventricular premature complexes.....	6.1.18
12109	Atrial fibrillation with slow ventricular response with aberrant conduction, or ventricular premature complexes.....	6.1.18
1220	Rapid atrial rhythm .....	6.1.4
1250	Atrial flutter .....	6.1.8
12503	Atrial flutter with aberrant conduction, or ventricular premature complexes ...	6.1.19
12505	Cannot rule out atrial flutter .....	6.1.8
1300	Junctional rhythm.....	6.1.5
1320	Rapid junctional rhythm .....	6.1.5
1400	Undetermined rhythm (Possible supraventricular rhythm).....	6.1.6
1420	Undetermined rhythm (Possible supraventricular tachycardia).....	6.1.6
1430	Undetermined rhythm (Possible supraventricular bradycardia) .....	6.1.6
1470	with occasional supraventricular premature complexes .....	6.1.14
1474	with frequent supraventricular premature complexes .....	6.1.14
1475	with frequent supraventricular premature complexes in a pattern of bigeminy .....	6.1.14
1570	with occasional ventricular premature complexes .....	6.1.15
15708	with occasional ventricular premature complexes (Unreliable analysis due to noise).....	6.1.15
1574	with frequent ventricular premature complexes .....	6.1.15
15748	with frequent ventricular premature complexes (Unreliable analysis due to noise).....	6.1.15
1575	with frequent ventricular premature complexes in a pattern of bigeminy.....	6.1.15
15758	with frequent ventricular premature complexes in a pattern of bigeminy (Unreliable analysis due to noise).....	6.1.15
1577	with couplet ventricular premature complexes.....	6.1.15
15778	with couplet ventricular premature complexes (Unreliable analysis due to noise).....	6.1.15
16006	Electronic atrial pacemaker .....	6.1.9
16007	Electronic ventricular pacemaker.....	6.1.9
16008	Electronic atrial pacemaker (Unreliable analysis due to noise) .....	6.1.9
16009	Electronic ventricular pacemaker (Unreliable analysis due to noise) .....	6.1.9
1901	Undetermined regular rhythm.....	6.1.10
1902	Undetermined rhythm .....	6.1.10
1921	Undetermined regular rhythm (tachycardia) .....	6.1.10
1922	Undetermined rhythm (tachycardia) .....	6.1.10
1931	Undetermined regular rhythm (bradycardia).....	6.1.10
1932	Undetermined rhythm (bradycardia).....	6.1.10
1938	Extreme bradycardia .....	6.1.11
1970	with occasional ectopic premature complexes .....	6.1.17
19708	with occasional ectopic premature complexes (Unreliable analysis due to noise).....	6.1.17
1974	with frequent ectopic premature complexes .....	6.1.17
19748	with frequent ectopic premature complexes (Unreliable analysis due to noise).....	6.1.17
1975	with frequent ectopic premature complexes in a pattern of bigeminy.....	6.1.17

<b>Section 6-2</b>	<b>Conductive Defect .....</b>	<b>6.2.1</b>	<b>1</b>
	2210 Short PR interval.....	6.2.2	
	2216 Type-A Wolff-Parkinson-White syndrome.....	6.2.3	<b>2</b>
	2217 Type-B Wolff-Parkinson-White syndrome.....	6.2.3	
	2218 Atypical Wolff-Parkinson-White syndrome .....	6.2.3	<b>3</b>
	2219 Intermittent Wolff-Parkinson-White syndrome.....	6.2.3	
	2231 First degree AV block .....	6.2.5	<b>4</b>
	2232 2nd degree AV block, Mobitz type I.....	6.2.5	
	2233 2nd degree AV block, Mobitz type II.....	6.2.5	<b>5</b>
	2234 Possible 3rd degree AV block.....	6.2.5	
	2320 Nonspecific intraventricular conduction delay .....	6.2.10	<b>6</b>
	2330 Nonspecific intraventricular conduction block .....	6.2.10	<b>6-1</b>
	2420 RSR (QR) in lead V1/V2, consistent with right ventricular conduction delay ...	6.2.6	
	2440 Incomplete right bundle branch block.....	6.2.6	<b>6-2</b>
	2450 Right bundle branch block.....	6.2.6	
	24501 Right bundle branch block, plus possible RVH.....	6.2.6	<b>6-3</b>
	2540 Incomplete left bundle branch block.....	6.2.8	
	2550 Left bundle branch block .....	6.2.8	<b>6-4</b>
	2630 Left anterior fascicular block.....	6.2.9	
	2730 Left posterior fascicular block.....	6.2.9	<b>6-5</b>
<b>Section 6-3</b>	<b>Myocardial Infarction.....</b>	<b>6.3.1</b>	<b>6-6</b>
	Analysis Criteria.....	6.3.2	
	Anterior Myocardial Infarction .....	6.3.5	<b>A</b>
	Septal Myocardial Infarction .....	6.3.8	
	Lateral Myocardial Infarction.....	6.3.10	
	Inferior Myocardial Infarction.....	6.3.12	
	Children .....	6.3.14	
<b>Section 6-4</b>	<b>ST-T Abnormality .....</b>	<b>6.4.1</b>	
	ST Depression.....	6.4.2	
	Injury.....	6.4.4	
	Subendocardial Ischemia .....	6.4.8	
	Early Repolarization .....	6.4.11	
	Pericarditis .....	6.4.12	
	T Wave Abnormality.....	6.4.13	
	Nonspecific ST Elevation.....	6.4.15	
<b>Section 6-5</b>	<b>Ventricular Hypertrophy.....</b>	<b>6.5.1</b>	
	Point Score System .....	6.5.2	
	Analysis Criteria for RVH.....	6.5.2	
	5120 Possible right ventricular hypertrophy.....	6.5.4	
	5130 Right ventricular hypertrophy.....	6.5.4	
	5134 Right ventricular hypertrophy, probably repolarization abnormality .....	6.5.5	
	Analysis Criteria for LVH.....	6.5.6	
	5211 Minimal voltage criteria for LVH, may be normal variant .....	6.5.7	
	5220 Possible left ventricular hypertrophy.....	6.5.7	
	5222 Moderate voltage criteria for LVH, may be normal variant.....	6.5.7	
	5233 Voltage criteria for LVH .....	6.5.7	

5234 Left ventricular hypertrophy with repolarization abnormality ..... 6.5.7

**Section 6-6 Atrial Enlargement, Abnormal Axis Deviation and Others..... 6.6.1**

6120 Possible right atrial enlargement ..... 6.6.2  
6130 Right atrial enlargement ..... 6.6.2  
6220 Possible left atrial enlargement..... 6.6.3  
6230 Left atrial enlargement..... 6.6.3  
7100 Abnormal right axis deviation ..... 6.6.4  
7102 Moderate right axis deviation..... 6.6.4  
7200 Abnormal left axis deviation..... 6.6.4  
7202 Moderate left axis deviation ..... 6.6.4  
7300 Indeterminate axis ..... 6.6.4  
7400 S1-S2-S3 pattern, consistent with pulmonary disease, RVH, or normal  
variant..... 6.6.6  
7500 Abnormal QRS-T angle ..... 6.6.7  
8003 Consistent with pulmonary disease ..... 6.6.8  
8100 Low QRS voltage..... 6.6.9  
8101 Low QRS voltage in limb leads..... 6.6.9  
8102 Low QRS voltage in chest leads..... 6.6.9  
8200 Dextrocardia ..... 6.6.10  
8304 Long QTc interval ..... 6.6.11  
8305 Short QTc interval ..... 6.6.11  
0101 Possible arm leads reversed, check lead requested ..... 6.6.12  
0102 ARTIFACT PRESENT..... 6.6.13  
0103 CANNOT ANALYZE ECG..... 6.6.13  
0104 ELECTRODE(S) FAILURE...Repeat ECG is required..... 6.6.13  
0201 ...Analysis based on intrinsic rhythm ..... 6.6.13

**Appendix Modified Minnesota Code ..... A.1**

General ..... A.2  
Code List..... A.3  
Priority of Code Printing..... A.7  
Detailed Criteria..... A.8

# *Section 1 Introduction*

## 1. INTRODUCTION

ECAPS 12C is the ECG analysis program for the Nihon Kohden's instruments, such as electrocardiographs. A computer analysis program is merely a collection of ECG evaluation criteria created by physicians. It is not possible for a computer program to correctly judge every unique ECG, so sometimes it makes wrong interpretations where a physician could very easily read and interpret the waveforms. The final decision can only be made by the qualified physicians. Use this system only as a diagnostic aid, based on proper understanding of its features and limitations.

This manual describes the criteria of the analysis results of the ECAPS 12C output data.

For the detailed operation procedure of the system, refer to the operator's manual of the instrument.

### **NOTE**

The contents of this manual are subject to change without prior notice for improvement of analysis precision.



**ADVISORY**

The ECAPS 12C analysis program is applicable to ages 3 and older. Ages below 3 years are treated according to the criteria for the age of 3 years.

Final determination of overall interpretation judgement, diagnosis and treatment must be made by a qualified physician.

**Patient Age and the Age Used for Analysis**

- When classifying as adult or child:

<b>Age</b>	<b>Age which is used for analysis</b>
Child (3-15)	12
Adult (>=16)	35

- When classifying in age range:

<b>Age</b>	<b>Age which is used for analysis</b>
<=5	3
6-9	7
10-14	12
15-34	25
>=34	35

- When patient's actual age is entered:

The patient's actual age is used for analysis.

## ***Section 2 Precautions (Discrepancies between Physician's and ECAPS 12C's Findings)***

## 2. PRECAUTIONS (DISCREPANCIES BETWEEN PHYSICIAN'S AND ECAPS 12C'S FINDINGS)

The causes for discrepancy between computer analysis and physician's findings and the countermeasures to be taken for these causes are given below.

<b>Causes</b>	<b>Countermeasures</b>
Difference in judgement criteria between physician and computer program, or lack of applicable findings by computer program.	Refer to Section 6 "Criteria of Findings" which list all the judgement criteria used by ECAPS 12C.
The ECG waveform is on the borderline of the judgement criteria.	Compare the measured data with the data in Section 6 "Criteria of Findings". This may be the limit of computer analysis. Patient data other than ECG should be taken into consideration.
Artifact (EMG, AC interference, baseline wandering, etc.) not recognized, leading to wrong interpretation.	Try to record ECG with as little artifact as possible.
Arrhythmia, etc. which are intrinsically difficult for the computer to analyze.	This is a limit for computer analysis. The advice of a physician should be obtained.

## 2. PRECAUTIONS (DISCREPANCIES BETWEEN PHYSICIAN'S AND ECAPS 12C'S FINDINGS)

Although various means of eliminating errors and discrepancies between computer analysis and physician's findings are employed in the ECAPS 12C, not all ECG diagnostic cases are incorporated, therefore, consider the following points when using the computer's analysis results.


- (1) ECAPS 12C is not programmed to compensate for the influence of medicine. Check the dosing record when reading the ECG. However, the influence of digitalis is noted on the recording paper according to the presence of atrial fibrillation.
- (2) ECAPS 12C is not programmed to compensate for fluid balance such as abnormal electrolytes. When interpreting ECG, read the QTc interval, T waveform and U waveform, along with other relevant test results.
- (3) ECAPS 12C does not further classify premature complex into trigeminy, short-run, etc.
- (4) ECAPS 12C is not programmed to identify escaped beats and pararrhythmia. These may be interpreted as "undetermined rhythms".
- (5) ECAPS 12C is not programmed to identify LGL syndrome. Judge this syndrome from "short PR interval".
- (6) ECAPS 12C does not identify wandering pacemakers.

# *Section 3 Outline of ECG Automatic Recording*

Entering Patient Data .....	3.3
Recording ECG Data.....	3.4
Improving Waveform Quality.....	3.5
Measuring ECG Waveform .....	3.6
General Measurement Values to Be Printed .....	3.9
Classification and Printing ECG Data.....	3.10

### 3. OUTLINE OF ECG AUTOMATIC RECORDING

#### ECG Measurement Flowchart



## Entering Patient Data

Enter the following data before acquiring ECG data. With some instruments, some of these items cannot be entered. Refer to the operator's manual of the instrument for details of entering the patient data.

- Patient identification number (ID)
- Name
- Gender
- Birth Date
- Age
- Height
- Weight
- Blood pressure
- Medication
- Date\*
- Hour\*

\* Date and hour are automatically set.

### NOTE

- Among these items, only age and gender affect the analysis. Other items have no effect on the analysis.
- If no age is input, the factory default setting of 35 is used. If no gender is specified, the factory default setting of male is used.
- For accurate analysis results, input gender and age.

## Recording ECG Data

The ECGs of all 12 standard leads are acquired simultaneously for 10 seconds at an accuracy of 1.25  $\mu\text{V}/\text{bit}$  and 500 samples/s.

Refer to the operator's manual of the instrument for details of recording ECG.




## Improving Waveform Quality

During ECG data acquisition, the quality of the ECG waveforms is improved with digital filters and the adverse influence of baseline wandering due to electrode potential drifting and AC interference is minimized.


- AC interference

A digital filter is used in the system to reduce the AC frequency components. A digital filter eliminates adverse influence on the ECG waveform more than a conventional analog filter.



50/60<sub>Hz</sub>

Conventional analog filter




50/60<sub>Hz</sub>

Digital filter

Baseline wandering

Shortening the time constant distorts the ST segments which reduces diagnostic accuracy. A digital filter is used to remove the components which cause baseline wandering.




- High frequency noise

High cut filters of different frequencies are incorporated in the instruments. The high cut filter cuts the high frequency components of ECG but reduces the effect of EMG. The high cut filter attenuates the QRS amplitude although the influence on judgement of LVH (left ventricular hypertrophy) is reduced as much as possible. However, the ECAPS 12C program always analyzes the ECG waveform acquired by 150 Hz filter. Therefore, some differences may occur between the recorded waveform and the analysis result.



## Measuring ECG Waveform

ECG waveforms are measured as shown below.



Waveforms are measured as shown on the next page. The reference point in measuring waveforms is the starting point of the QRS wave.

(a) Waveform measurement parameters




Total QRS amplitude = R max. amp. + S max. amp.


Net QRS amplitude = R max. amp. - S max. amp.

### 3. OUTLINE OF ECG AUTOMATIC RECORDING

(b) QRS area



(c) Upward oriented T



T is upward oriented:  $C > 0.05 \text{ mV} + (1.5 \times \text{STJ})$

(d) Modified T amplitude (T amp. (mod))

To simplify the treatment of biphasic T wave, and to explain the T amplitude when ST and T differ from the QRS starting point, T amplitude is modified as below.

When T' is present:

T amp. (mod) = (T amp. or T' amp., whichever is smaller) - (STE or T end, whichever is larger)


When T' is not present:

T amp. (mod) = T amp. - (STE or T end, whichever is larger)

**General Measurement Values to Be Printed**

Printout	Meaning
HR (Vent. rate) -- bpm	Heart rate
PR int. -- ms	PR interval
QRS dur. -- ms	QRS duration
QT int. -- ms	QT interval
QTc int. -- ms	QTc interval (See p. 3-7(a))
P axis -- °	P axis deviation
QRS axis -- °	QRS axis deviation
T axis -- °	T axis deviation
RV5 amp -- mV	R amplitude in V5
SV1 amp -- mV	S amplitude in V1

When a measurement value cannot be measured, “\*\*\*” is displayed.



## Classification and Printing ECG Data

The ECGs are classified into the findings. For the details of the classification method, refer to Section 6 “Criteria of Findings”.

The findings are printed out as “Analysis Result”. Refer to the next section.

# *Section 4 Recorded Analysis Results*

## 4. RECORDED ANALYSIS RESULTS

The analysis results of the system are printed out in the format selected by the operator. The outline of the formats and the operation are described in the operator's manual of the instrument. The explanation of the common features of the printout with typical examples are given below. These are factory default settings.

Not all of the following items are printed with some instruments. For details, refer to the operator's manual of the instrument.

### (1) Patient's data and recording conditions

- ① Data and Time
- ② Name
- ③ ID No.
- ④ Gender
- ⑤ Birth date and Age (in some formats, birth date is omitted): years
- ⑥ Vent. rate (Heart Rate)
- ⑦ Paper speed
- ⑧ Lead name

#### **NOTE**

①, ⑥, ⑦ and ⑧ are automatically entered. Items ② to ⑤ are printed as entered by the operator.

### (2) ECG waveform printout (In some formats, these are omitted.)

- ⑨ Dominant ECG waveform (averaged)
- ⑩ Rhythm lead
- ⑪ Calibration wave


### (3) Analysis result

- ⑫ Overall judgement
- ⑬ ECG findings name [criteria]
- ⑭ General measurement values
- ⑮ Physician's signature (Reviewed)

#### **NOTE**

- The findings [criteria] are printed as a supplementary which is one of the judgement criteria. For details, refer to Section 5 "How to Read Analysis Results".
- In some findings, [criteria] are not printed.






# *Section 5 How to Read Analysis Results*

ECG Findings .....	5.3
Priority of ECG Findings.....	5.3
Criteria.....	5.4
Overall Judgement.....	5.4
Analysis Result When ECG Cannot Be Analyzed .....	5.5

## 5. HOW TO READ ANALYSIS RESULTS

The analysis result printout, as shown below, contains the following data.

1. ECG findings [Criteria]
2. Overall judgement
3. General measurement values



## ECG Findings

The ECAPS 12C classifies the ECGs into about 200 findings by comparing the features of the ECGs with the analysis criteria specified for each finding. The analysis criteria are all taken from the judgement criteria used by the physicians and arranged for computer processing. For details, refer to Section 6 “Criteria of Findings”.

For the analysis process, refer to Section 3 “Outline of ECG Automatic Recording”.

When the obtained ECG exactly conforms to the analysis criteria, the finding name is classified as “determined”. However, there are findings which cannot be determined definitely and findings which are on the borderline between conformance and nonconformance, liable to be influenced by slight noises. Those findings which are close to “determined” but not definite are classified as “possible”, and those findings which may or might not conform but cannot be completely ruled out are classified as “Cannot rule out”.

Some arrhythmia ECGs are difficult to classify into finding names; these are classified as “undetermined rhythm”.

### NOTE

When the ECG is classified as “undetermined rhythm” (1901, 1902, 1921, 1922, 1931, 1932), there is no further rhythm analysis.

## Priority of ECG Findings

When the obtained ECG conforms to two or more findings criteria, only the most important finding name is printed. For example, if “ischemia” is found to apply, less important findings such as “nonspecific ST & T wave abnormality” or “ST elevation” are not printed.

Further details of analysis are given as NOTE under the respective findings in Section 6 “Criteria of Findings”.

## Criteria

The analysis criteria can be printed out on the recording paper after the findings, as shown on p. 5.2. They should serve as an aid in observing the ECG.

For example, when the instrument judges the ECG as abnormal junctional ST depression because there is junctional ST depression of more than 0.1 mV in V5 and V6, the criteria is printed as “4023 Abnormal junctional ST depression [0.1 + mV junctional ST depression (V5, V6)].

For reading convenience, the analysis findings and criteria are printed in simplified form.

For further details of the analysis criteria, refer to Section 6 “Criteria of Findings”.

## Overall Judgement

The findings are classified into one of five overall judgements. The list of all findings for each overall judgement is shown in Section 6 “Criteria of Findings”.

	<b>Overall Judgement</b>	<b>Code</b>
1	abnormal ECG	9150
2	abnormal rhythm ECG	9140
3	borderline ECG	9130
4	normal ECG	9110
5	atypical ECG	9120

Where two or more findings are output, only the highest priority finding is printed.

Although only the highest priority finding is selected for judgement, the findings requiring immediate treatment, such as “Myocardial infarction, possible acute” are printed as “abnormal ECG” instead of “borderline ECG”. This is to raise an alarm in case of emergency although there is actually little possibility of danger. When an unusual ECG pattern is recognized as neither “normal ECG” nor “abnormal ECG”, it is judged as “atypical ECG”. For example, “Low QRS voltage”, “undetermined axis”, and “undetermined rhythm” are judged as “atypical ECG”.

## Analysis Result When ECG Cannot Be Analyzed

When a measurement value cannot be measured, “\*\*\*” is displayed. When no measurement values can be measured, the ECG cannot be analyzed and a “0103 CANNOT ANALYZE ECG” message is displayed. For recording example, refer to “General Measurement Values to Be Printed” in Section 3.

# *Section 6 Criteria of Findings*

Section 6-1	Arrhythmias .....	6.1.1
Section 6-2	Conductive Defect .....	6.2.1
Section 6-3	Myocardial Infarction .....	6.3.1
Section 6-4	ST-T Abnormality .....	6.4.1
Section 6-5	Ventricular Hypertrophy .....	6.5.1
Section 6-6	Atrial Enlargement, Abnormal Axis Deviation and Others .....	6.6.1

## Indication of Analysis Criteria

### 1. How to read analysis criteria

The analysis criteria are marked with (1), (2) and ①, ②.

(1) ...indicates AND

① ...indicates OR

[ Example 1 ]

(1) condition A

(2) condition B

(3) ① condition C

② condition D

Conditions A **and** B are satisfied, **and** either condition C or condition D are satisfied.

[ Example 2 ]

(1) ① condition A

(2) ② • condition B

• condition C

} Either condition A is satisfied, **or** both conditions B and C are satisfied.

### 2. Age reference tables

The following tables show the values used for certain analysis parameters which vary according to age. The applicable items are marked with “\*” in the criteria column.

#### NOTE

- When no age is input in the patient information, the factory default setting of age 35 is used.
- When the age is 2 or younger, the computer's results may not be accurate.

[ PR interval ] (%)

Age	PR ratio
<=5	76
<=9	83
<=11	85
<=14	88
<=18	91
>18	100

(Example)

PR interval = 170 ms

Age = 6 years old

PR criteria = PR int × PR ratio/100

= 170 × 83/100

= 141 ms



[ QRS duration ] (%)

Age	QRS duration
<=5	82
<=9	85
<=11	85
<=14	92
<=18	95
>18	100

(Example)

QRS duration: 0.12 s

Age = 10 years old

$$\begin{aligned} \text{QRS criteria} &= \text{QRS dur.} \times \text{QRS dur (\%)/100} \\ &= 0.12 \times 85/100 \\ &= 0.102 \text{ s} \end{aligned}$$


[ S duration: Lateral lead (I, aVL, V4, V5, V6) ] (s)

Age	S duration	S duration
	0.040 (s)	0.060 (s)
<=5	0.031	0.047
<=9	0.034	0.051
<=11	0.037	0.056
<=14	0.039	0.059
<=18	0.039	0.059
>18	0.039	0.059

 shows the parameters used in criteria for each finding.

[ R duration (V1, V2) ] (s)

Age	R duration	R duration
	0.020 (s)	0.030 (s)
<=5	0.016	0.024
<=9	0.017	0.026
<=11	0.019	0.028
<=14	0.020	0.030
<=18	0.020	0.030
>18	0.020	0.030

 shows the parameters used in criteria for each finding.

## 6. CRITERIA OF FINDINGS

[ R duration (I, aVL, V4, V5, V6) ] (s)

Age	R duration	R duration	R duration
	0.060 (s)	0.100 (s)	0.250 (s)
<=5	0.049	0.082	0.207
<=9	0.051	0.086	0.218
<=11	0.054	0.090	0.228
<=14	0.057	0.095	0.240
<=18	0.059	0.099	0.250
>18	0.059	0.099	0.250

■ shows the parameters used in criteria for each finding.

[ Electrical axis (left axis) ] (degree)

Age	LAXD1		LAXD2	
	Male	Female	Male	Female
<=5	15°	19°	5°	9°
<=9	9°	9°	-1°	-1°
<=11	4°	23°	-6°	13°
<=14	5°	20°	-5°	10°
<=18	-14°	13°	-24°	3°
>18	-20°	-20°	-30°	-30°

[ Electrical axis (right axis) ] (degree)

Age	RAXD1		RAXD2	
	Male	Female	Male	Female
<=5	97°	101°	107°	111°
<=9	97°	97°	107°	107°
<=11	92°	100°	102°	110°
<=14	97°	97°	107°	107°
<=18	99°	100°	109°	110°
>18	90°	90°	100°	100°

[ R amplitude (V1) ] (mV)

<b>Age</b>	<b>Male</b>	<b>Female</b>
<=5	2.20	1.75
<=9	1.65	1.60
<=11	1.30	1.40
<=14	1.60	1.20
<=18	1.40	1.15
>18	1.00	1.00

[ Average R amplitude (V1) ] (mV)

<b>Age</b>	<b>Male</b>	<b>Female</b>
<=5	0.90	0.85
<=9	0.75	0.65
<=11	0.60	0.55
<=14	0.70	0.45
<=18	0.60	0.40
>18	0.40	0.30

[R amplitude (V6) ] (mV)

<b>Age</b>	<b>Male</b>	<b>Female</b>
<=5	1.45	1.45
<=9	1.70	1.65
<=11	1.65	1.50
<=14	1.60	1.35
<=18	1.50	1.20
>18	1.20	1.00

## 6. CRITERIA OF FINDINGS

[ S depth (V1) ] (mV)

<b>Age</b>	<b>Male</b>	<b>Female</b>
<=5	1.20	1.25
<=9	1.20	1.25
<=11	1.35	1.30
<=14	1.50	1.15
<=18	1.65	1.15
>18	1.20	0.95

[ S depth (V6) ] (mV)

<b>Age</b>	<b>Male</b>	<b>Female</b>
<=5	0.54	0.54
<=9	0.65	0.60
<=11	0.61	0.61
<=14	0.47	0.35
<=18	0.50	0.36
>18	0.40	0.30

[R/S (V1) ratio]

<b>Age</b>	<b>Male</b>	<b>Female</b>
<=5	2.7	2.0
<=9	1.8	1.8
<=11	2.3	1.3
<=14	2.2	1.6
<=18	1.3	1.7
>18	1.0	1.0

[ VAT (V1) ] (s)

<b>Age</b>	<b>Male</b>	<b>Female</b>
<=5	0.037	0.039
<=9	0.033	0.029
<=11	0.045	0.028
<=14	0.041	0.032
<=18	0.039	0.034
>18	0.050	0.050

[ Heart Rate ] (BPM) - For arrhythmia analysis -

<b>Age</b>	<b>Bradycardia base</b>	<b>Tachycardia base</b>
<=5	65	140
<=8	60	135
<=12	55	130
<=16	50	120
>16	50	100

BPM : Beats per minute

# Section 6-1 Arrhythmias

Rhythm analysis .....	6.1.2
1002 Marked rhythm irregularity, possible non-conducted PAC, SA block, AV block, or sinus pause .....	6.1.13
1100 Sinus rhythm .....	6.1.3
1102 Sinus arrhythmia .....	6.1.12
1108 Marked sinus arrhythmia .....	6.1.13
1120 Sinus tachycardia .....	6.1.3
1130 Sinus bradycardia.....	6.1.3
1200 Atrial rhythm .....	6.1.4
1210 Atrial fibrillation.....	6.1.7
12101 Atrial fibrillation with rapid ventricular response .....	6.1.7
12102 Atrial fibrillation with slow ventricular response.....	6.1.7
12103 Atrial fibrillation with aberrant conduction, or ventricular premature complexes.....	6.1.18
12108 Atrial fibrillation with rapid ventricular response with aberrant conduction, or ventricular premature complexes .....	6.1.18
12109 Atrial fibrillation with slow ventricular response with aberrant conduction, or ventricular premature complexes .....	6.1.18
1220 Rapid atrial rhythm .....	6.1.4
1250 Atrial flutter .....	6.1.8
12503 Atrial flutter with aberrant conduction, or ventricular premature complexes.....	6.1.19
12505 Cannot rule out atrial flutter.....	6.1.8
1300 Junctional rhythm .....	6.1.5
1320 Rapid junctional rhythm.....	6.1.5
1400 Undetermined rhythm (Possible supraventricular rhythm).....	6.1.6
1420 Undetermined rhythm (Possible supraventricular tachycardia) .....	6.1.6
1430 Undetermined rhythm (Possible supraventricular bradycardia).....	6.1.6
1470 with occasional supraventricular premature complexes .....	6.1.14
1474 with frequent supraventricular premature complexes .....	6.1.14
1475 with frequent supraventricular premature complexes in a pattern of bigeminy .....	6.1.14
1570 with occasional ventricular premature complexes .....	6.1.15
15708 with occasional ventricular premature complexes (Unreliable analysis due to noise).....	6.1.15
1574 with frequent ventricular premature complexes.....	6.1.15
15748 with frequent ventricular premature complexes (Unreliable analysis due to noise).....	6.1.15
1575 with frequent ventricular premature complexes in a pattern of bigeminy .....	6.1.15
15758 with frequent ventricular premature complexes in a pattern of bigeminy (Unreliable analysis due to noise).....	6.1.15
1577 with couplet ventricular premature complexes .....	6.1.15
15778 with couplet ventricular premature complexes (Unreliable analysis due to noise) .....	6.1.15
16006 Electronic atrial pacemaker.....	6.1.9
16007 Electronic ventricular pacemaker .....	6.1.9
16008 Electronic atrial pacemaker (Unreliable analysis due to noise).....	6.1.9
16009 Electronic ventricular pacemaker (Unreliable analysis due to noise) .....	6.1.9
1901 Undetermined regular rhythm.....	6.1.10
1902 Undetermined rhythm.....	6.1.10
User's Guide ECAPS 12C .....	6.1.1

## 6. CRITERIA OF FINDINGS

1921	Undetermined regular rhythm (tachycardia).....	6.1.10
1922	Undetermined rhythm (tachycardia).....	6.1.10
1931	Undetermined regular rhythm (bradycardia).....	6.1.10
1932	Undetermined rhythm (bradycardia).....	6.1.10
1938	Extreme bradycardia.....	6.1.11
1970	with occasional ectopic premature complexes.....	6.1.17
19708	with occasional ectopic premature complexes (Unreliable analysis due to noise).....	6.1.17
1974	with frequent ectopic premature complexes.....	6.1.17
19748	with frequent ectopic premature complexes (Unreliable analysis due to noise).....	6.1.17
1975	with frequent ectopic premature complexes in a pattern of bigeminy.....	6.1.17

### **Rhythm analysis**

Arrhythmia analysis is divided into three major categories; “basic rhythm analysis”, “basic rhythm fluctuation analysis”, and “premature complex analysis”. In the basic rhythm analysis, the presence of P wave and P wave axis are important factors. The connection between P waves and QRS waves is an important factor in classifying the waveform into basic rhythm, basic rhythm fluctuation and premature complex.

## 1. Basic rhythm analysis

## (1) With P wave

- Sinus

Code	Findings [Criteria]	Judgement
1100	Sinus rhythm	normal ECG
1120	Sinus tachycardia	abnormal rhythm ECG
1130	Sinus bradycardia	abnormal ECG

**Analysis criteria**

Findings	Criteria
Sinus rhythm	(1) Electronic atrial pacemaker is not used. (2) P waveform; constant, PR interval; regular (3) $-30^{\circ} \leq P \text{ axis} < 120^{\circ}$ (4) $50 \leq \text{heart rate} < 100^*$
Sinus tachycardia	(1) Electronic atrial pacemaker is not used. (2) P waveform; constant, PR interval; regular (3) $-30^{\circ} \leq P \text{ axis} < 120^{\circ}$ (4) $100 \leq \text{heart rate}^*$
Sinus bradycardia	(1) Electronic atrial pacemaker is not used. (2) P waveform; constant, PR interval; regular (3) $-30^{\circ} \leq P \text{ axis} < 120^{\circ}$ (4) $50 > \text{heart rate}^*$

**NOTE**

The values marked with “\*” vary with age. For details, refer to p.6.0.7.

	50/minute	100/minute
Bradycardia	Sinus rhythm	Tachycardia



## 6. CRITERIA OF FINDINGS

- Atrial

Code	Findings [Criteria]	Judgement
1200	Atrial rhythm	abnormal rhythm ECG
1220	Rapid atrial rhythm	abnormal rhythm ECG

### Analysis criteria

Findings	Criteria
Atrial rhythm	(1) Electronic atrial pacemaker is not used. (2) P waveform; constant, PR interval; regular (3) ① ① $120^\circ \leq \text{P axis} \leq 240^\circ$ ② $-30^\circ > \text{P axis} \geq -60^\circ$ ② (1) PR interval $> 0.14\text{s}$ (2) ① $120^\circ \leq \text{P axis} \leq 270^\circ$ ② $-30^\circ > \text{P axis} \geq -90^\circ$ (4) Heart rate $\leq 70$
Rapid atrial rhythm	(1) Electronic atrial pacemaker is not used. (2) P waveform; constant, PR interval; regular (3) ① ① $120^\circ \leq \text{P axis} \leq 240^\circ$ ② $-30^\circ > \text{P axis} \geq -60^\circ$ ② (1) PR interval $> 0.14\text{ s}$ (2) ① $120^\circ \leq \text{P axis} \leq 270^\circ$ ② $-30^\circ > \text{P axis} \geq -90^\circ$ (4) $70 < \text{heart rate}$

70/minute

Atrial rhythm	Rapid
---------------	-------

- AV junction

Code	Findings [Criteria]	Judgement
1300	Junctional rhythm	abnormal rhythm ECG
1320	Rapid junctional rhythm	abnormal rhythm ECG

#### Analysis criteria

Findings	Criteria
Junctional rhythm	(1) Electronic atrial pacemaker is not used. (2) P waveform; constant, PR interval; regular PR $\leq 0.14$ s (3) ① $-60^\circ > \text{P axis} > -90^\circ$ ② $240^\circ < \text{P axis} \leq 270^\circ$ (4) Heart rate $\leq 70$ (age $> 3$ years) Heart rate $\leq 80$ (age $\leq 3$ years)
Rapid junctional rhythm	(1) Electronic atrial pacemaker is not used. (2) P waveform; constant, PR interval; regular PR $\leq 0.14$ s (3) ① $-60^\circ > \text{P axis} > -90^\circ$ ② $240^\circ < \text{P axis} \leq 270^\circ$ (4) $70 < \text{heart rate}$ (age $> 3$ years) $80 < \text{heart rate}$ (age $\leq 3$ years)

70/minute

Junctional rhythm	Rapid
-------------------	-------

## 6. CRITERIA OF FINDINGS

- (2) Without P wave
- Without atrial fibrillation or without atrial flutter

Code	Findings [Criteria]	Judgement
1400	Undetermined rhythm (Possible supraventricular rhythm)	abnormal rhythm ECG
1420	Undetermined rhythm (Possible supraventricular tachycardia)	abnormal rhythm ECG
1430	Undetermined rhythm (Possible supraventricular bradycardia)	abnormal rhythm ECG

### Analysis criteria

When the P wave is not joined to the QRS, the following criteria are used for classification.

Findings	Criteria
Undetermined rhythm (Possible supraventricular rhythm)	(1) Electronic atrial pacemaker is not used. (2) Atrial fibrillation and atrial flutter are not present. (3) QRS duration < 120 ms* (4) 50 ≤ heart rate < 100*
Undetermined rhythm (Possible supraventricular tachycardia)	(1) Electronic atrial pacemaker is not used. (2) Atrial fibrillation and atrial flutter are not present. (3) QRS duration < 120 ms* (4) heart rate ≥ 100*
Undetermined rhythm (Possible supraventricular bradycardia)	(1) Electronic atrial pacemaker is not used. (2) Atrial fibrillation and atrial flutter are not present. (3) QRS duration < 120 ms* (4) heart rate < 50*

### NOTE

- The values marked with “\*” vary with age. For details, refer to p.6.0.3 and 6.0.7.
- When these findings are recognized, other less important findings are not printed on the recording paper.

50/minute	100/minute
Bradycardia	Supraventricular rhythm
	Tachycardia

- With atrial fibrillation or with atrial flutter

Code	Findings [Criteria]	Judgement
1210	Atrial fibrillation	abnormal rhythm ECG
12101	Atrial fibrillation with rapid ventricular response	abnormal rhythm ECG
12102	Atrial fibrillation with slow ventricular response	abnormal rhythm ECG

#### Analysis criteria

Findings	Criteria
Atrial fibrillation	(1) ① (1) P wave is not present (2) RR interval deviation $> 0.125 \times$ mean RR interval (3) Random RR interval ② Fibrillation wave (f wave) is detected. (2) $50 \leq$ heart rate $< 100^*$
Atrial fibrillation with rapid ventricular response	(1) ① (1) P wave is not present (2) RR interval deviation $> 0.125 \times$ mean RR interval (3) Random RR interval ② Fibrillation wave (f wave) is detected (2) $100 \leq$ heart rate $^*$
Atrial fibrillation with slow ventricular response	(1) ① (1) P wave is not present (2) RR interval deviation $> 0.125 \times$ mean RR interval (3) Random RR interval ② Fibrillation wave (f wave) is detected (2) $50 >$ heart rate $^*$

#### NOTE

The values marked with “\*” vary with age. For details, refer to p.6.0.7.

	50/minute	100/minute
Bradycardia	Atrial fibrillation	Tachycardia

## 6. CRITERIA OF FINDINGS

Code	Findings [Criteria]	Judgement
1250	Atrial flutter	abnormal rhythm ECG
12505	Cannot rule out atrial flutter	abnormal rhythm ECG

### Analysis criteria

Findings	Criteria
Atrial flutter	(1) Flutter wave (including fibrillation or flutter wave) is present.
Atrial flutter with aberrant conduction, or ventricular premature complexes	(1) Atrial flutter is present. (2) Ectopic QRS is not pacemaker waveform. (3) Ectopic QRS with duration > 0.12 s*

### NOTE

The values marked with “\*” vary with age. For details refer to p. 6.0.3.

## (3) With pacemaker

Code	Findings [Criteria]	Judgement
16006	Electronic atrial pacemaker	atypical ECG
16007	Electronic ventricular pacemaker	atypical ECG
16008	Electronic atrial pacemaker (Unreliable analysis due to noise)	atypical ECG
16009	Electronic ventricular pacemaker (Unreliable analysis due to noise)	atypical ECG

## Analysis criteria

Findings	Criteria
Electronic atrial pacemaker rhythm	(1) QRS is dominant waveform. (2) Pacemaker pulse recognized within 0.08 s before and after the beginning of the P wave (3) Three or more heartbeats satisfying above two conditions are present.
Electronic ventricular pacemaker rhythm	Pacemaker pulse and dominant QRS are joined.

## NOTE

- When there is any detached electrode or noise during analysis, “16008” or “16009” appear in the findings column even if the above conditions are satisfied.
- If spontaneous systoles occur during ECG observation, the ECG system analyzes the spontaneous systoles. In this case, “0201 Analysis based on intrinsic rhythm” is additionally printed out in the findings column.

## 6. CRITERIA OF FINDINGS

### 2. Basic rhythm fluctuation analysis

Code	Findings [Criteria]	Judgement
1901	Undetermined regular rhythm	atypical ECG
1902	Undetermined rhythm	atypical ECG
1921	Undetermined regular rhythm (tachycardia)	atypical ECG
1922	Undetermined rhythm (tachycardia)	atypical ECG
1931	Undetermined regular rhythm (bradycardia)	atypical ECG
1932	Undetermined rhythm (bradycardia)	atypical ECG

#### Analysis criteria

Findings	Criteria
Undetermined regular rhythm	(1) Normal P wave is not present. (2) QRS duration $\geq 0.12$ s* (3) Atrial fibrillation is not present. (4) Atrial flutter is not present. (5) Electronic pacemaker is not used. (6) All QRSs are dominant type. (7) Maximum RR interval - Minimum RR interval $< 1/8$ mean RR interval
Undetermined rhythm	(1) Normal P wave is not present. (2) QRS duration $\geq 0.12$ s* (3) Atrial fibrillation is not present. (4) Atrial flutter is not present. (5) Electronic pacemaker is not used. (6) Regular rhythm is not present.
Undetermined regular rhythm (tachycardia)	(1) Cannot determine the rhythm (PR interval; regular) (2) $100^* \leq$ Heart rate
Undetermined rhythm (tachycardia)	(1) Cannot determine the rhythm (2) $100^* \leq$ Heart rate
Undetermined regular rhythm (bradycardia)	(1) Cannot determine the rhythm (PR interval; regular) (2) $50^* >$ Heart rate
Undetermined rhythm (bradycardia)	(1) Cannot determine the rhythm (2) $50^* >$ Heart rate

#### NOTE

The values marked with "\*" vary with age. For details refer to p. 6.0.3 and 6.0.7.

<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
1938	Extreme bradycardia	abnormal rhythm ECG

**Analysis criteria**

<b>Findings</b>	<b>Criteria</b>
Extreme bradycardia	(1) $40 >$ Heart rate (2) 2° AV block (Mobitz type II) is not present.



## 6. CRITERIA OF FINDINGS

- With abnormal rhythm

<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
1002	Marked rhythm irregularity, possible non-conducted PAC, SA block, AV block, or sinus pause	abnormal rhythm ECG

### Analysis criteria

<b>Findings</b>	<b>Criteria</b>
Marked rhythm irregularity, possible non-conducted PAC, SA block, AV block or sinus pause	(1) 2° AV block is not present. (2) Heart rate < 100 (3) Random RR interval

Code	Findings [Criteria]	Judgement
1102	Sinus arrhythmia	normal ECG
1108	Marked sinus arrhythmia	normal ECG

### Analysis criteria

When sinus rhythm, sinus tachycardia and sinus bradycardia are judged and satisfy the below conditions, then the classification is as below.

Findings	Criteria
Sinus arrhythmia	(1) Premature complexes are not present. (2) Marked rhythm irregularity is not present. (3) RR interval deviation $> 0.2 \times$ mean RR interval
Marked sinus arrhythmia	(1) Premature complexes are not present. (2) Marked rhythm irregularity is not present. (3) RR interval deviation $> 0.4 \times$ mean RR interval

### NOTE

“Premature complexes: absent” and “Marked rhythm irregularity: absent” are analyzed according to their analysis criteria on p. 6.1.13 to 6.1.16.

## 6. CRITERIA OF FINDINGS

### 3. Premature complex analysis

#### (1) Supraventricular

Code	Findings [Criteria]	Judgement
1470	with occasional supraventricular premature complexes	abnormal rhythm ECG
1474	with frequent supraventricular premature complexes	abnormal rhythm ECG
1475	with frequent supraventricular premature complexes in a pattern of bigeminy	abnormal rhythm ECG

#### Analysis criteria

Findings	Criteria
with occasional supraventricular premature complexes	(1) Intermittent WPW is not present. (2) 2° AV block is not present. (3) Marked rhythm irregularity is not present. (4) ① RR interval < mean RR interval × 3/4 ② • P wave is not sinus-induced • RR interval < mean RR interval – mean RR interval × 1/10 (maximum 100 ms) • The second heartbeat on the recording paper. Or, when the heartbeat is the third or later on the recording paper; the RR interval of the previous heartbeat ≥ RR interval + 10 ms. Or, when the heartbeat is the third or later on the recording paper; the previous heartbeat is premature complexes.
with frequent supraventricular premature complexes	Three or more supraventricular premature complexes given above (code 1470) and/or ectopic premature complexes are present.
with frequent supraventricular premature complexes in a pattern of bigeminy	Supraventricular premature complexes given above (code 1470) and dominant waveform appear alternately.

## (2) Ventricular

<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
1570	with occasional ventricular premature complexes	abnormal rhythm ECG
1574	with frequent ventricular premature complexes	abnormal rhythm ECG
1575	with frequent ventricular premature complexes in a pattern of bigeminy	abnormal rhythm ECG
1577	with couplet ventricular premature complexes	abnormal rhythm ECG
15708	with occasional ventricular premature complexes (Unreliable analysis due to noise)	abnormal rhythm ECG
15748	with frequent ventricular premature complexes (Unreliable analysis due to noise)	abnormal rhythm ECG
15758	with frequent ventricular premature complexes in a pattern of bigeminy (Unreliable analysis due to noise)	abnormal rhythm ECG
15778	with couplet ventricular premature complexes (Unreliable analysis due to noise)	abnormal rhythm ECG

## 6. CRITERIA OF FINDINGS

### Analysis criteria

Findings	Criteria
with occasional ventricular premature complexes	(1) Intermittent WPW is not present. (2) Ectopic QRS duration > 0.12 s* (3) ① For the first heartbeat on the recording paper: <ul style="list-style-type: none"> <li>• P wave does not precede QRS-complexes</li> <li>• next RR interval – 40 ms &gt; mean RR interval</li> </ul> ② For the first heartbeat on the recording paper: <ul style="list-style-type: none"> <li>• P wave precedes QRS-complex</li> <li>• next RR interval – 100 ms &gt; mean RR interval</li> </ul> ③ When the heartbeat is the second or later: <ul style="list-style-type: none"> <li>• P wave does not precede QRS-complexes</li> <li>• RR interval + 40 ms &lt; mean RR interval</li> </ul> ④ When the heartbeat is the second or later: <ul style="list-style-type: none"> <li>• P wave precedes QRS-complexes</li> <li>• RR interval &lt; mean RR interval – mean RR interval × 1/10 (maximum 100 ms)</li> <li>• The second heartbeat on the recording paper. Or, when the heartbeat is the third or later on the recording paper; the RR interval of the previous heartbeat ≥ RR interval + 10 ms. Or, when the heartbeat is the third or later on the recording paper; the previous heartbeat is premature complexes.</li> </ul>
with frequent ventricular premature complexes	Three or more ventricular premature complexes given above (code 1570) are present.
with frequent ventricular premature complexes in a pattern of bigeminy	The ventricular premature complexes given above (code 1570) and dominant waveform appear alternately.
with couplet ventricular premature complexes	More than two ventricular premature complexes given above (code 1570) appear consecutively.
with occasional ventricular premature complexes (Unreliable analysis due to noise)	Analysis criteria for code 1570 is satisfied and there is electrode detachment or noise during analysis.
with frequent ventricular premature complexes (Unreliable analysis due to noise)	Analysis criteria for code 1574 is satisfied and there is electrode detachment or noise during analysis.
with frequent ventricular premature complexes in a pattern of bigeminy (Unreliable analysis due to noise)	Analysis criteria for code 1575 is satisfied and there is electrode detachment or noise during analysis.
with couplet ventricular premature complexes (Unreliable analysis due to noise)	Analysis criteria for code 1577 is satisfied and there is electrode detachment or noise during analysis.

#### NOTE

- When any of above conditions is satisfied and there is any electrode detached or noise exists during analysis, “15708”, “15748”, “15758”, or “15778” appears in the finding column.
- The value marked with “\*” varies with age. For details, refer to p. 6.0.3.

## (3) Ectopic

Code	Findings [Criteria]	Judgement
1970	with occasional ectopic premature complexes	abnormal rhythm ECG
1974	with frequent ectopic premature complexes	abnormal rhythm ECG
1975	with frequent ectopic premature complexes in a pattern of bigeminy	abnormal rhythm ECG
19708	with occasional ectopic premature complexes (Unreliable analysis due to noise)	abnormal rhythm ECG
19748	with frequent ectopic premature complexes (Unreliable analysis due to noise)	abnormal rhythm ECG

## Analysis criteria

Findings	Criteria
with occasional ectopic premature complexes	(1) Satisfies all criteria for ventricular premature complexes other than ectopic QRS (2) $0.06\text{ s} < \text{ectopic QRS duration} \leq 0.12\text{ s}$
with frequent ectopic premature complexes	Three or more ectopic premature complexes given above are present.
with frequent ectopic premature complexes in a pattern of bigeminy	Above occasional ectopic premature complexes and dominant waveform appear alternately.

## NOTE

- When either the conditions 1970 or 1974 is satisfied and there is any electrode detached or noise exists during analysis, “19708” or “19748” appears in the finding column.
- Ectopic premature complexes do not appear in the finding column when it is judged together with the occasional supraventricular premature complexes. This is because the instrument judges it to be occasional supraventricular premature complexes with aberrant conduction.

## 6. CRITERIA OF FINDINGS

### (4) With atrial fibrillation or with atrial flutter

Code	Findings [Criteria]	Judgement
12103	Atrial fibrillation with aberrant conduction, or ventricular premature complexes	abnormal rhythm ECG
12108	Atrial fibrillation with rapid ventricular response with aberrant conduction, or ventricular premature complexes	abnormal rhythm ECG
12109	Atrial fibrillation with slow ventricular response with aberrant conduction, or ventricular premature complexes	abnormal rhythm ECG

### Analysis criteria

Findings	Criteria
Atrial fibrillation with aberrant conduction, or ventricular premature complexes	(1) Atrial fibrillation is present. (2) Ectopic QRS is not pacemaker waveform (3) Ectopic QRS with interval $\geq 0.12$ s*
Atrial fibrillation with rapid ventricular response with aberrant conduction, or ventricular premature complexes	(1) Atrial fibrillation with rapid ventricular response is present. (2) Ectopic QRS is not pacemaker waveform (3) Ectopic QRS with interval $\geq 0.12$ s*
Atrial fibrillation with slow ventricular response with aberrant conduction, or ventricular premature complexes	(1) Atrial fibrillation with slow ventricular response is present. (2) Ectopic QRS is not pacemaker waveform (3) Ectopic QRS with interval $\geq 0.12$ s*

### NOTE

- For “Atrial fibrillation”, “Atrial fibrillation with rapid ventricular response” and “Atrial fibrillation with slow ventricular response”, the same analysis criteria given on p. 6.1.7 are used.
- The values marked with “\*” vary with age. For details, refer to p. 6.0.3.

Code	Findings [Criteria]	Judgement
12503	Atrial flutter with aberrant conduction, or ventricular premature complexes	abnormal rhythm ECG

#### Analysis criteria

Findings	Criteria
Atrial flutter with aberrant conduction, or ventricular premature complexes	(1) Atrial flutter is present. (2) Ectopic QRS is not pacemaker waveform (3) Ectopic QRS with interval $\geq 0.12$ s*

#### NOTE

- For “Atrial flutter”, the same analysis criteria given on p. 6.1.8 is used.
- The values marked with “\*” varies with age. For details, refer to p. 6.0.3.



## *Section 6-2 Conductive Defect*

2210	Short PR interval .....	6.2.2
2216	Type-A Wolff-Parkinson-White syndrome .....	6.2.3
2217	Type-B Wolff-Parkinson-White syndrome .....	6.2.3
2218	Atypical Wolff-Parkinson-White syndrome .....	6.2.3
2219	Intermittent Wolff-Parkinson-White syndrome .....	6.2.3
2231	First degree AV block .....	6.2.5
2232	2nd degree AV block, Mobitz type I .....	6.2.5
2233	2nd degree AV block, Mobitz type II .....	6.2.5
2234	Possible 3rd degree AV block.....	6.2.5
2320	Nonspecific intraventricular conduction delay .....	6.2.10
2330	Nonspecific intraventricular conduction block.....	6.2.10
2420	RSR (QR) in lead V1/V2, consistent with right ventricular conduction delay.....	6.2.6
2440	Incomplete right bundle branch block.....	6.2.6
2450	Right bundle branch block .....	6.2.6
24501	Right bundle branch block, plus possible RVH.....	6.2.6
2540	Incomplete left bundle branch block.....	6.2.8
2550	Left bundle branch block .....	6.2.8
2630	Left anterior fascicular block.....	6.2.9
2730	Left posterior fascicular block.....	6.2.9

## 6. CRITERIA OF FINDINGS

### 1. A-V conductive defect


Code	Findings [Criteria]	Judgement
2210	Short PR interval	atypical ECG

### Analysis criteria

Findings	Criteria
Short PR interval	(1) Pacemaker is not used. (2) P waveform and PR interval are both constant. (3) $-60^{\circ} \leq \text{P axis} \leq 240^{\circ}$ (4) PR interval $< 0.12 \text{ s}^*$

### NOTE

The value marked with "\*" varies with age. Refer to p. 6.0.2.



- Wolff-Parkinson-White syndrome

Code	Findings [Criteria]	Judgement
2216	Type-A Wolff-Parkinson-White syndrome	abnormal ECG
2217	Type-B Wolff-Parkinson-White syndrome	abnormal ECG
2218	Atypical Wolff-Parkinson-White syndrome	abnormal ECG
2219	Intermittent Wolff-Parkinson-White syndrome	abnormal ECG

### Analysis criteria

Findings	Criteria
Type-A WPW syndrome	(1) ① • PR interval $\leq 0.12$ s* • Delta waves are recognized in at least two leads. ② • PR interval $\leq 0.14$ s* • Delta waves are recognized in at least three leads. ③ Delta waves are recognized in at least five leads. ④ • PR interval $\leq 0.12$ s* • Q wave is not present and VAT $> 0.08$ s in at least two leads (2) Maximum R amplitude $>$ Maximum S amplitude in V1
Type-B WPW syndrome	(1) ① • QRS area ratio $> 0.4$ in at least two leads among I, V5 and V6 • R duration $> 0.03$ s in V2 • PR interval $\leq 0.14$ s* ② • PR interval $\leq 0.12$ s* • Delta waves are recognized in at least two leads. ③ • PR interval $\leq 0.14$ s* • Delta waves are recognized in at least three leads. ④ Delta waves are recognized in at least five leads. ⑤ • PR interval $\leq 0.12$ s* • Q wave is not present and VAT $> 0.08$ s in at least two leads (2) Maximum R amplitude $\leq$ Maximum S amplitude in V1
Atypical WPW syndrome	① • PR interval $\leq 0.12$ s* • Delta waves are recognized in at least two leads. ② • PR interval $\leq 0.14$ s* • Delta waves are recognized in at least three leads. ③ Delta waves are recognized in at least five leads. ④ • PR interval $\leq 0.12$ s* • Q wave is not present and VAT $> 0.08$ s in at least two leads
Intermittent WPW syndrome	(1) Four or more ectopic QRS without pacemaker pulses are present. (2) Heart rate $< 120$ (3) RR interval of the ectopic beat $+ 0.16$ s $>$ RR interval of dominant beat (4) Delta waves are recognized in at least two leads (5) PR interval of the ectopic beat $< 0.14$ s* (6) PR interval of the ectopic beat $<$ Mean PR interval of the dominant QRS- 0.02 s (7) PJ interval of the ectopic beat $>$ PJ interval of the dominant QRS-0.02 s

## 6. CRITERIA OF FINDINGS

### NOTE

- (1) If any of the following is satisfied, WPW check is not done.
  - ① P wave is not the same type as the dominant beat's.
  - ② PR interval  $> 0.17$  s\*
  - ③ QRS duration  $< 0.10$  s\*
  - ④ QRS duration  $> 0.20$  s\*
  - ⑤ Heart rate  $> 120$ \*
- (2) If WPW is determined by the above analysis criteria, other waveform analysis will be omitted.
- (3) The values marked with "\*" vary with age. Refer to p. 6.0.2, 6.0.3 and 6.0.7.
- (4) PJ interval is the time between the starting point of P wave to the end point of QRS wave (STJ point).

- AV block

Code	Findings [Criteria]	Judgement
2231	First degree AV block	abnormal ECG
2232	2nd degree AV block, Mobitz type I	abnormal ECG
2233	2nd degree AV block, Mobitz type II	abnormal ECG
2234	Possible 3rd degree AV block	abnormal ECG

#### Analysis criteria

Findings	Criteria
First degree AV block	(1) No pacemaker is used. (2) P waveform and PR interval are both constant. (3) $-60^{\circ} \leq \text{P axis} \leq 240^{\circ}$ (4) PR interval $\geq 0.21 \text{ s}^*$
2nd degree AV block, Mobitz type I	(1) 2nd degree AV block (Mobitz type II) is not present. (2) Both the preceding and current heartbeats are in dominant waveform. (3) QRS dropout (which is characteristic of Mobitz type I) exists which is calculated from RR interval.
2nd degree AV block, Mobitz type II	QRS dropout (which is characteristic of Mobitz type II) exists which is calculated from RR interval.
Possible 3rd degree AV block	(1) P wave is not present. (2) Heart rate $< 50$ (3) Differences between RR intervals are less than 2% of the mean RR interval.

#### NOTE

The values marked with “\*” vary with age. Refer to p. 6.0.2.

## 6. CRITERIA OF FINDINGS

### 2. Intra-ventricular conductive defect

- Right bundle branch block

Code	Findings [Criteria]	Judgement
2420	RSR (QR) in lead V1/V2 consistent with right ventricular conduction delay	borderline ECG
2440	Incomplete right bundle branch block	borderline ECG
2450	Right bundle branch block	abnormal ECG
24501	Right bundle branch block, plus possible RVH	abnormal ECG

### Analysis criteria

Findings	Criteria
RSR (QR) in lead V1/V2 consistent with right ventricular conduction delay	In V1 or V2 ① • R amplitude > 0.1 mV • R duration > 0.02 s* • S wave is not present. ② • R' amplitude > 0.1 mV • R' duration > 0.02 s*
Incomplete right bundle branch block	(1) 0.09 s < QRS duration < 0.12 s* (2) In two leads among I, aVL, V4, V5 and V6, S duration >= 0.04 s* (3) Right ventricular conduction delay is present.
Right bundle branch block	(1) QRS duration >= 0.12 s* (2) QRS area > 0 in V1 (3) S duration > 0.04 s* in 2 or more leads among I, aVL, V4, V5, V6 (4) R duration < 0.10 s* in 4 or more leads among I, aVL, V4, V5, V6 (5) Does not end with S or S' wave in V1 Or (1) QRS duration >= 0.105 s* (2) S duration > 0.06 s* in 3 or more leads among I, aVL, V4, V5, V6 (3) R or R' duration > 0.06 s* in V1 (4) QRS area > 0 in V1

Findings	Criteria
Right bundle branch block, plus possible RVH	(1) RBBB is present. (2) ① • Age $\geq$ 1 year old • R or R' amplitude $>$ 1.5 mV in V1 ② R or R' amplitude $>$ 2.0 mV in V1 (3) $110^\circ <$ QRS axis $\leq 270^\circ$ ( $>$ 14 years old) $120^\circ <$ QRS axis $\leq 270^\circ$ ( $\leq$ 14 years old)

**NOTE**

- The values marked with “\*” vary with age. Refer to p. 6.0.3 and 6.0.4.
- When RBBB is recognized, no right axis deviation is judged.
- With right bundle branch conduction defects, the terminating vector is directed towards the anterior and right and it is extended. In the ECAPS 12C criteria, not only QRS duration but also the presence of R at V1, and of wide S in at least two lateral leads are required.
- Previously, QRS duration over 0.12 second was a criterion for bundle branch blocks.

However, since QRS duration over 0.105 second with wide S in the lateral leads and wide R in V1 appearing is also judged as RBBB, this case has been included in the ECAPS 12C program.

## 6. CRITERIA OF FINDINGS

- Left bundle branch block

Code	Findings [Criteria]	Judgement
2540	Incomplete left bundle branch block	abnormal ECG
2550	Left bundle branch block	abnormal ECG

### Analysis criteria

Findings	Criteria
Incomplete left bundle branch block	<ol style="list-style-type: none"> <li>(1) QRS duration <math>&gt; 0.105 \text{ s}^*</math></li> <li>(2) Net QRS amplitude <math>&lt; 0</math> in V1, V2</li> <li>(3) Q/S duration <math>\geq 0.080 \text{ s}^*</math> in V1, V2</li> <li>(4) Q wave is not present in at least 2 leads among I, V5, V6</li> <li>(5) R duration <math>\geq 0.060 \text{ s}^*</math> in at least 2 leads among I, aVL, V5, V6</li> </ol>
Left bundle branch block	<ol style="list-style-type: none"> <li>(1) Incomplete LBBB is present.</li> <li>(2) R + R' duration <math>\geq 0.1 \text{ s}^*</math> in any leads among I, aVL, V6</li> <li>(3) QRS area ratio <math>&gt; 0.25</math> in I or V6</li> <li>(4) <ol style="list-style-type: none"> <li>① QRS duration <math>\geq 0.16 \text{ s}^*</math></li> <li>② <ul style="list-style-type: none"> <li>• QRS duration <math>\geq 0.14 \text{ s}^*</math></li> <li>• total of each R + R' duration <math>\geq 0.25 \text{ s}^*</math> in I, aVL, V6</li> </ul> </li> <li>③ QRS duration <math>\geq 0.12 \text{ s}^*</math> <ul style="list-style-type: none"> <li>• total of each R + R' duration <math>\geq 0.25 \text{ s}^*</math> in I, aVL, V6</li> <li>• QRS area ratio <math>&gt; 0.4</math> in at least 2 leads among I, aVL, V6</li> </ul> </li> </ol> </li> </ol>

### NOTE

- The values marked with “\*” vary with age. Refer to p. 6.0.3 and 6.0.4.
- When “incomplete LBBB” is judged, “Moderate level left axis deviation” is not printed out.
- When “LBBB” is judged, “Moderate level left axis deviation”, “Left anterior fascicular block”, and “Left posterior fascicular block” are not printed out.
- QRS area in the analysis criteria for LBBB means the area from the start to the end of QRS. Refer to p. 3.8 (b). This area increases through R type, expansion, notch, etc. Therefore, the boundary values are determined on the basis of typical LBBB cases. This value is used instead of R wave to classify between the true LBBB and R pattern where R wave is expanded by nonspecific slur at the end of the wave.
- There is no specific definition for incomplete LBBB in ECG analysis. In the ECAPS 12C program, incomplete LBBB is defined very narrowly, and whenever other findings are applicable, such as left posterior fascicular block, then one of those findings is taken.



- Fascicular block

Code	Findings [Criteria]	Judgement
2630	Left anterior fascicular block	abnormal ECG
2730	Left posterior fascicular block	abnormal ECG

#### Analysis criteria

Findings	Criteria
Left anterior fascicular block	(1) $-90^{\circ} < \text{QRS axis} \leq -45^{\circ}$ (2) R amplitude > Q amplitude in I and aVL (3) Q wave is present in I (4) S or S' amplitude > maximum R amplitude in II
Left posterior fascicular block	(1) Age $\geq 1$ year old (2) S pattern is not present. (3) Right atrial enlargement is not present. (4) Lung disease is not recognized. (5) $110^{\circ} \leq \text{QRS axis} \leq 270^{\circ}$ (> 14 years old) $120^{\circ} \leq \text{QRS axis} \leq 270^{\circ}$ ( $\leq 14$ years old) (6) R amplitude > Q amplitude (III and aVF) (7) Q wave is present (III and aVF).

#### NOTE

- When “Left anterior fascicular block” is judged, “Moderate left axis deviation” and “Left axis deviation” are not printed out.
- When “Left posterior fascicular block” is judged, “Moderate right axis deviation” and “Right axis deviation” are not printed out.

## 6. CRITERIA OF FINDINGS

- Nonspecific

<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
2320	Nonspecific intraventricular conduction delay	borderline ECG
2330	Nonspecific intraventricular conduction block	abnormal ECG

### Analysis criteria

<b>Findings</b>	<b>Criteria</b>
Nonspecific intraventricular conduction delay	(1) Block-related findings are not present. (2) QRS duration > 0.11 s*
Nonspecific intraventricular conduction block	(1) Criteria for RBBB and LBBB are not satisfied. (2) QRS duration > 0.13 s*

### NOTE

- The values marked with “\*” vary with age. Refer to p. 6.0.3.
- Intraventricular conduction delay is judged when the criteria on p. 6.2.5 to 6.2.8 are not satisfied, and QRS duration is not large enough to be judged as blocks.

## *Section 6-3 Myocardial Infarction*

Analysis Criteria .....	6.3.2
Anterior Myocardial Infarction.....	6.3.5
Septal Myocardial Infarction .....	6.3.8
Lateral Myocardial Infarction .....	6.3.10
Inferior Myocardial Infarction .....	6.3.12
Children.....	6.3.14

### Analysis Criteria

Myocardial infarction arises from the stricture and obstruction of coronary arteries and death of heart muscles. It is mainly caused by coronary arteriosclerosis. It is said that myocardial infarction shows a characteristic ECG waveform.

#### 1. Features and analysis method of ECG of myocardial infarction

The Q wave duration is commonly taken as the main factor for judging the presence of myocardial infarction and likewise, ECAPS 12C also checks the Q wave duration. In addition to Q wave duration, repolarization abnormality that is sometimes reported to accompany “acute” or “recent” cases of myocardial infarction is also considered to be useful in detecting myocardial infarction. For example, elevated ST junctions and negative T waves are very clear signs of myocardial infarction, even when significant Q waves are absent.

When abnormal repolarization is taken into consideration, both the sensitivity and discrimination level of analysis are improved for acute and recent cases. For old cases, observing the QRS amplitude and duration improves the analysis precision.

With this program, these factors are integrated into “equivalent Q duration” for analysis processing.

The use of this new factor does not differ greatly from conventional analysis. It measures and processes age, gender, Q duration, Q amplitude QRS duration and QRS amplitude as a whole to improve the precision level of analysis. This “equivalent Q duration” is used in the analysis criteria unless stated otherwise.

The ECAPS 12C analysis program classifies ECGs for myocardial infarction diagnostic purposes as follows:

Item	Division	Criteria	Description
Age classification	Acute (?)	ST elevation	Varies with the extent of ST elevation
	Recent (?)		
	Old (?)	No ST elevation	
	Age undetermined	Age cannot be determined.	
Level classification	Determined	Equivalent Q duration: 40 ms ~	Classified by eq. Q duration (Eq. Q duration is one corrected for Q amp., R amp. and QRS duration.)
	Possible	35 ~ 39 ms	
	Cannot rule out	30 ~ 34 ms	

#### (1) Age classification

It is said that after the occurrence of myocardial infarction, the T wave increases height and ST elevates within several hours, and then, with some delay, abnormal Q waves appear. It is also said that the coronary T wave starts to appear while the ST elevation is improving (2 days to 1 week) and remains for a long period of time. Although the time division, “acute (?)”, “recent (?)”, “age undetermined” and “old (?)”, is classified according to the ST elevation level, this is a classification of ECG diagnosis, not a clinical diagnosis.

In atypical instances, ST elevation may be observed for more than a few months.

#### (2) Level classification

Abnormal Q wave is the most important factor in identifying myocardial infarction. The divisions “determined”, “possible” and “cannot rule out” given here according to the durations of abnormal Q waves are the terms expressing the level of definiteness of the ECG features or the level of computer analysis, and they are not directly related to the level of seriousness in clinical diagnosis.

## 6. CRITERIA OF FINDINGS

### 2. Diagnosis of infarct portion by abnormal Q wave


The locations of infarction are classified according to the leads showing abnormal Q waves:

	I	II	III	aVR	aVL	aVF	V1	V2	V3	V4	V5	V6
Anterior								V	V	V	V	
Septal							$\Delta^1$	V				
Lateral	V				V						V	V
Inferior		V				V						
Posterior inferior		V				V	$\Delta^2$	$\Delta^2$				
Anterolateral	V				V			V	V	V	V	
Anteroseptal							$\Delta^1$	V	V	V	V	

1: Q wave is present

2: By R duration

The program analyzes the infarct portions in the following sequence:



### 3. Other features

- (1) Abnormal Q waves may appear in cases of “LBBB”, “WPW syndrome”, “pulmonary embolism”, “LVH”, and “RVH”.
- (2) Sometimes the characteristic ECG of infarction does not appear at all depending on the time after the nosogenesis, the size and portion of the infarction. Clinical reviews including chemical blood tests (GOT, GPT, LDH, CPK) are recommended.
- (3) An intermixed EMG or AC noise is sometimes mistaken for small R waves. This program is designed to use the typical heartbeat waves that contain Q waves with priority among all typical heartbeat waves, but record favorable waveforms by eliminating artifact as much as possible.

## Anterior Myocardial Infarction

When the following criteria are satisfied, anterior myocardial infarction is not analyzed.

- LBBB
- QRS duration > 140 ms, negative net QRS amplitude at V1

### Analysis criteria

Age classification	Criteria	Other
Acute (?)	STM > 0.2 mV in V3 and V4 and STE > 0.2 mV in V3 and V4 and Modified T amplitude $\geq 0^*$ in V3 and V4	
Recent (?)	STM > 0.05 mV in V3 or V4 and STE > 0.05 mV in same leads as above and Modified T amplitude < $0^*$ in V3 or V4	
Old (?)	STM < 0.03 mV in V3 and V4 and Modified T amplitude $\geq 0^*$ in V3 and V4	Acute (?) Recent (?) are not satisfied.
Age undetermined	Acute (?), recent (?) and Old (?) are not satisfied.	

\* Refer to p. 3.8 (d) for measurement method.

## 6. CRITERIA OF FINDINGS

Level classification	Criteria	Other
Cannot rule out	Q duration $\geq$ 30 ms in V2 and V3, or V3 and V4, or V4 and V5 or R amplitude $<$ 0.2 mV in V4	
Possible	Q duration $\geq$ 30 ms in V2 or V4 and Q duration $\geq$ 35 ms in V3 or Q duration $\geq$ 30 ms in V3 or V5 and Q duration $\geq$ 35 ms in V4	LVH is not present.
(Determined)	Q duration $\geq$ 30 ms in V2 or V4 and Q duration $\geq$ 40 ms in V3 or Q duration $\geq$ 30 ms in V3 or V5 and Q duration $\geq$ 40 ms in V4 or “Cannot rule out” is satisfied	LVH is not present. Chest lead low voltage is not present. Nonspecific intraventricular conduction block is not present.  Acute (?) or Recent (?) is satisfied.

### Code and findings

Code	Findings [Criteria]	Judgement
3113	Cannot rule out anterior myocardial infarction, probably old	abnormal ECG
3114	Cannot rule out anterior myocardial infarction, age undetermined	abnormal ECG
3121	Possible anterior myocardial infarction, possibly acute	abnormal ECG
3122	Possible anterior myocardial infarction, probably recent	abnormal ECG
3123	Possible anterior myocardial infarction, probably old	abnormal ECG
3124	Possible anterior myocardial infarction, age undetermined	abnormal ECG
3131	Anterior myocardial infarction, possibly acute	abnormal ECG
3132	Anterior myocardial infarction, probably recent	abnormal ECG
3133	Anterior myocardial infarction, probably old	abnormal ECG
3134	Anterior myocardial infarction, age undetermined	abnormal ECG
3213	Cannot rule out anteroseptal myocardial infarction, probably old	abnormal ECG
3214	Cannot rule out anteroseptal myocardial infarction, age undetermined	abnormal ECG
3221	Possible anteroseptal myocardial infarction, possibly acute	abnormal ECG
3222	Possible anteroseptal myocardial infarction, probably recent	abnormal ECG
3223	Possible anteroseptal myocardial infarction, probably old	abnormal ECG
3224	Possible anteroseptal myocardial infarction, age undetermined	abnormal ECG
3231	Anteroseptal myocardial infarction, possibly acute	abnormal ECG
3232	Anteroseptal myocardial infarction, probably recent	abnormal ECG
3233	Anteroseptal myocardial infarction, probably old	abnormal ECG
3234	Anteroseptal myocardial infarction, age undetermined	abnormal ECG
3313	Cannot rule out anterolateral myocardial infarction, probably old	abnormal ECG
3314	Cannot rule out anterolateral myocardial infarction, age undetermined	abnormal ECG
3321	Possible anterolateral myocardial infarction, possibly acute	abnormal ECG
3322	Possible anterolateral myocardial infarction, probably recent	abnormal ECG
3323	Possible anterolateral myocardial infarction, probably old	abnormal ECG



<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
3324	Possible anterolateral myocardial infarction, age undetermined	abnormal ECG
3331	Anterolateral myocardial infarction, possibly acute	abnormal ECG
3332	Anterolateral myocardial infarction, probably recent	abnormal ECG
3333	Anterolateral myocardial infarction, probably old	abnormal ECG
3334	Anterolateral myocardial infarction, age undetermined	abnormal ECG

## 6. CRITERIA OF FINDINGS

### Septal Myocardial Infarction

When the following criteria are satisfied, septal myocardial infarction is not analyzed.

- LBBB
- QRS duration > 140 ms and net QRS amplitude is negative in V1.
- Cannot rule out anterior infarction and Q wave is not present in V1.

#### Analysis criteria

Age classification	Criteria	Other
Acute (?)	STM and STE > 0.2 mV in V2 and Modified T amplitude $\geq 0^*$ in V2	
Recent (?)	STM and STE > 0.05 mV in V2 and Modified T amplitude < $0^*$ in V2	
Old (?)	STM < 0.05 mV in V2 and Modified T amplitude $\geq 0^*$ in V2	Acute (?) and Recent (?) is satisfied.
Age undetermined	Acute (?), Recent (?) and Old (?) are not satisfied.	

Level classification	Criteria	Other
Cannot rule out	Q duration $\geq 30$ ms in V2 or Q duration > 20 ms in V2	RBBB is present.
Possible	Q duration $\geq 35$ ms in V2	LVH is not present.
(Determined)	Q duration $\geq 40$ ms in V2	LVH is not present.

\* Refer to p. 3.8 (d) for measurement method.

**Code and findings**

<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
3411	Cannot rule out septal myocardial infarction, possibly acute	abnormal ECG
3412	Cannot rule out septal myocardial infarction, probably recent	abnormal ECG
3413	Cannot rule out septal myocardial infarction, probably old	abnormal ECG
3414	Cannot rule out septal myocardial infarction, age undetermined	abnormal ECG
3421	Possible septal myocardial infarction, possibly acute	abnormal ECG
3422	Possible septal myocardial infarction, probably recent	abnormal ECG
3423	Possible septal myocardial infarction, probably old	abnormal ECG
3424	Possible septal myocardial infarction, age undetermined	abnormal ECG
3431	Septal myocardial infarction, possibly acute	abnormal ECG
3432	Septal myocardial infarction, probably recent	abnormal ECG
3433	Septal myocardial infarction, probably old	abnormal ECG
3434	Septal myocardial infarction, age undetermined	abnormal ECG

## 6. CRITERIA OF FINDINGS

### Lateral Myocardial Infarction

#### Analysis criteria

Age classification	Criteria	Other
Acute (?)	STM and STE > 0.2 mV in V5 and V6 and STM and STE > 0.1 mV in I and aVL and Modified T amplitude $\geq 0^*$ in I, aVL, V5 and V6	
Recent (?)	STM and STE > 0.05 mV in one lead among I, aVL, V5 and V6 and Modified T amplitude < $0^*$ in one lead among I, aVL, V5 and V6	
Old (?)	STM < 0.03 mV in I, aVL, V5 and V6 and Modified T amplitude $\geq 0^*$ in I, aVL, V5 and V6	Acute (?) and Recent (?) are not satisfied.
Age undetermined	Acute (?), Recent (?) and Old (?) are not satisfied.	

Level classification	Criteria	Other
Cannot rule out	Q duration $\geq 30$ ms in 2 leads among I, aVL, V5 and V6	
Possible	Q duration $\geq 35$ ms in one lead among I, V5 and V6	Cannot rule out is satisfied.
(Determined)	Q duration $\geq 40$ ms in one lead among I, V5 and V6 or Cannot rule out	Cannot rule out is satisfied. Acute (?) and Recent (?) are satisfied.

\* Refer to p. 3.8 (d) for measurement method.

**Code and findings**

<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
3513	Cannot rule out lateral myocardial infarction, probably old	abnormal ECG
3514	Cannot rule out lateral myocardial infarction, age undetermined	abnormal ECG
3521	Possible lateral myocardial infarction, possibly acute	abnormal ECG
3522	Possible lateral myocardial infarction, probably recent	abnormal ECG
3523	Possible lateral myocardial infarction, probably old	abnormal ECG
3524	Possible lateral myocardial infarction, age undetermined	abnormal ECG
3531	Lateral myocardial infarction, possibly acute	abnormal ECG
3532	Lateral myocardial infarction, probably recent	abnormal ECG
3533	Lateral myocardial infarction, probably old	abnormal ECG
3534	Lateral myocardial infarction, age undetermined	abnormal ECG

6. CRITERIA OF FINDINGS

**Inferior Myocardial Infarction**

**Analysis criteria**

<b>Age classification</b>	<b>Criteria</b>	<b>Other</b>
Acute (?)	STM and STE > 0.1 mV in II and aVF and Modified T amplitude $\geq 0^*$ in II and aVF	
Recent (?)	STM and STE > 0.05 mV in II or aVF and Modified T amplitude < $0^*$ in II or aVF	
Old (?)	STM < 0.03 mV in II and aVF and Modified T amplitude $\geq 0^*$ in II and aVF	Acute (?) and Recent (?) are not satisfied.
Age undetermined	Acute (?), Recent (?) and Old (?) are not satisfied.	

<b>Level classification</b>	<b>Criteria</b>	<b>Other</b>
Cannot rule out	Q duration $\geq 30$ ms in II or aVF and Q amplitude (I) < Q amplitude (II) or Q duration $\geq 30$ ms in II or aVF and Q amplitude (I) < Q amplitude (aVF)	
Possible	Q duration $\geq 35$ ms in II or aVF	Cannot rule out is satisfied.
(Determined)	Cannot rule out is satisfied and Q duration $\geq 40$ ms in II or aVF or Cannot rule out is satisfied and Recent (?) is satisfied or Cannot rule out is satisfied and Acute (?) is satisfied.	

<b>Findings</b>	<b>Criteria</b>
Posterior Extension	Any of the inferior myocardial infarction is satisfied. and complete RBBB is not present. and Q amplitude = 0 mV (in V1 and V2) and R duration $\geq 40$ ms (in V1 and V2) or R duration $\geq 35$ ms and Net QRS amplitude > 0 mV (in V1 or V2) or R duration $\geq 30$ ms and Net QRS > 0 (in V1 and V2)

\* Refer to p. 3.8 (c) and (d) for measurement method.

## Code and findings

Code	Findings [Criteria]	Judgement
3613	Cannot rule out inferior myocardial infarction, probably old	abnormal ECG
36132	Cannot rule out inferior myocardial infarction with posterior extension, probably old	abnormal ECG
3614	Cannot rule out inferior myocardial infarction, age undetermined	abnormal ECG
36142	Cannot rule out inferior myocardial infarction with posterior extension, age undetermined	abnormal ECG
3621	Possible inferior myocardial infarction, possibly acute	abnormal ECG
36212	Possible inferior myocardial infarction with posterior extension, possibly acute	abnormal ECG
3622	Possible inferior myocardial infarction, probably recent	abnormal ECG
36222	Possible inferior myocardial infarction with posterior extension, probably recent	abnormal ECG
3623	Possible inferior myocardial infarction, probably old	abnormal ECG
36232	Possible inferior myocardial infarction with posterior extension, probably old	abnormal ECG
3624	Possible inferior myocardial infarction, age undetermined	abnormal ECG
36242	Possible inferior myocardial infarction with posterior extension, age undetermined	abnormal ECG
3631	Inferior myocardial infarction, possibly acute	abnormal ECG
36312	Inferior myocardial infarction with posterior extension, possibly acute	abnormal ECG
3632	Inferior myocardial infarction, probably recent	abnormal ECG
36322	Inferior myocardial infarction with posterior extension, probably recent	abnormal ECG
3633	Inferior myocardial infarction, probably old	abnormal ECG
36332	Inferior myocardial infarction with posterior extension, probably old	abnormal ECG
3634	Inferior myocardial infarction, age undetermined	abnormal ECG
36342	Inferior myocardial infarction with posterior extension, age undetermined	abnormal ECG

## 6. CRITERIA OF FINDINGS

### Children

With children under 18 years old, the following analysis are executed. The criteria are the same as the myocardial infarction. When the last number of the code (right) is 1, see the criteria for the code with the same first 4 numbers. When the last number is 3, see the criteria for the code with the same first 4 numbers but with the last (fifth) number 2.

For example, for the criteria of 31211 see the code 3121. For 36213, see the code 36212.

#### Code and findings

Code	Findings [Criteria]	Judgement
31211	Abnormal Q wave in lead V3/V4, cannot rule out cardiomyopathy	abnormal ECG
31221	Abnormal Q wave in lead V3/V4, cannot rule out cardiomyopathy	abnormal ECG
31311	Abnormal Q wave in lead V3/V4, cannot rule out cardiomyopathy	abnormal ECG
31321	Abnormal Q wave in lead V3/V4, cannot rule out cardiomyopathy	abnormal ECG
32211	Abnormal Q wave in lead V2 + V3/V4, cannot rule out cardiomyopathy	abnormal ECG
32221	Abnormal Q wave in lead V2 + V3/V4, cannot rule out cardiomyopathy	abnormal ECG
32311	Abnormal Q wave in lead V2 + V3/V4, cannot rule out cardiomyopathy	abnormal ECG
32321	Abnormal Q wave in lead V2 + V3/V4, cannot rule out cardiomyopathy	abnormal ECG
33211	Abnormal Q wave in lead I/aVL/V3-V6, cannot rule out cardiomyopathy	abnormal ECG
33221	Abnormal Q wave in lead I/aVL/V3-V6, cannot rule out cardiomyopathy	abnormal ECG
33311	Abnormal Q wave in lead I/aVL/V3-V6, cannot rule out cardiomyopathy	abnormal ECG
33321	Abnormal Q wave in lead I/aVL/V3-V6, cannot rule out cardiomyopathy	abnormal ECG
34111	Abnormal Q wave in lead V2, cannot rule out cardiomyopathy	abnormal ECG
34121	Abnormal Q wave in lead V2, cannot rule out cardiomyopathy	abnormal ECG
34211	Abnormal Q wave in lead V2, cannot rule out cardiomyopathy	abnormal ECG
34221	Abnormal Q wave in lead V2, cannot rule out cardiomyopathy	abnormal ECG
34311	Abnormal Q wave in lead V2, cannot rule out cardiomyopathy	abnormal ECG
34321	Abnormal Q wave in lead V2, cannot rule out cardiomyopathy	abnormal ECG
35211	Abnormal Q wave in lead I/aVL/V5/V6, cannot rule out cardiomyopathy	abnormal ECG
35221	Abnormal Q wave in lead I/aVL/V5/V6, cannot rule out cardiomyopathy	abnormal ECG
35311	Abnormal Q wave in lead I/aVL/V5/V6, cannot rule out cardiomyopathy	abnormal ECG
35321	Abnormal Q wave in lead I/aVL/V5/V6, cannot rule out cardiomyopathy	abnormal ECG
36211	Abnormal Q wave in lead II/aVF, cannot rule out cardiomyopathy	abnormal ECG
36221	Abnormal Q wave in lead II/aVF, cannot rule out cardiomyopathy	abnormal ECG
36311	Abnormal Q wave in lead II/aVF, cannot rule out cardiomyopathy	abnormal ECG
36321	Abnormal Q wave in lead II/aVF, cannot rule out cardiomyopathy	abnormal ECG
36213	Abnormal Q wave in lead II/aVF, cannot rule out cardiomyopathy	abnormal ECG
36223	Abnormal Q wave in lead II/aVF, cannot rule out cardiomyopathy	abnormal ECG
36313	Abnormal Q wave in lead II/aVF, cannot rule out cardiomyopathy	abnormal ECG
36323	Abnormal Q wave in lead II/aVF, cannot rule out cardiomyopathy	abnormal ECG



# *Section 6-4 ST-T Abnormality*

ST Depression.....	6.4.2
Injury .....	6.4.4
Subendocardial Ischemia .....	6.4.8
Early Repolarization .....	6.4.11
Pericarditis.....	6.4.12
T Wave Abnormality .....	6.4.13
Nonspecific ST Elevation.....	6.4.15

**ST Depression**

**Analysis criteria**

- (1) • Nonspecific intraventricular conduction block
  - RBBB\*
  - LBBB
  - ST elevation\*
  - Possible acute pericarditis
  - RVH (with repolarization abnormality)\*
  - LVH (with repolarization abnormality)

When any of these findings is recognized, the process of “(2) ST depression classification” is omitted. However, when the findings marked with “\*” are recognized, ST depression classification is performed, but for RBBB and RVH (repolarization abnormality), leads V1 and V2 are not analyzed.

(2) ST depression classification

Processing the data of all the leads except for a VR and III, the system adopts the finding uniformly judged by the data of two or more leads.

When the age is less than 16, V1 and V2 leads are not considered for “ST depression, possible digitalis effect” and “Minimal ST depression”.

Classification	STJ	STM	STE
Junctional ST depression, probably normal	< -0.1 mV		>= 0 mV
Abnormal junctional ST depression	< -0.1 mV		< 0 mV and >= STJ/2
ST depression, possible digitalis effect		< STJ and < -0.05 mV	
			< STJ and < -0.05 mV
Minimal ST depression	< -0.025 mV	< -0.025 mV	< -0.025 mV
Moderate ST depression		< -0.05 mV	< 0 mV
			< STJ and < STM and < -0.05 mV
*Marked ST depression, possible subendocardial injury	< -0.1 mV	< -0.1 mV	< -0.1 mV
Marked ST depression consistent with subendocardial injury	< -0.2 mV	< -0.2 mV	< -0.2 mV

**NOTE**

- When the finding marked with “\*”, i.e. “marked ST depression, possible subendocardial injury” is judged, and at the same time “atrial fibrillation” is found, “possible digitalis effect” is added to the findings.
- When any of the ST depression findings is recognized, “nonspecific ST elevation” is not analyzed.

**Code and findings**

<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
40106	ST depression, possible digitalis effect	abnormal ECG
4011	Minimal ST depression	borderline ECG
40116	Minimal ST depression, probably digitalis effect	borderline ECG
4012	Moderate ST depression	abnormal ECG
40126	Moderate ST depression, probably digitalis effect	abnormal ECG
4016	Marked ST depression, possible subendocardial injury	abnormal ECG
40166	Marked ST depression, possible subendocardial injury or digitalis effect	abnormal ECG
4017	Marked ST depression, consistent with subendocardial injury	abnormal ECG
4021	Junctional ST depression, probably normal	borderline ECG
4023	Abnormal junctional ST depression	borderline ECG

6. CRITERIA OF FINDINGS

**Injury**

**Analysis criteria**

- (1) • Left bundle branch block
- Right bundle branch block
- Nonspecific ventricular conductive defect
- Possible acute percarditis

When any of these findings is recognized, subendocardial injury is not analyzed.

(2) Injury Classification

Portion	Finding	Judgement Criteria	Lead
Anterior	Possible injury	All of the following conditions are satisfied: (1) $STJ \geq J$ threshold value. Refer to the table 1. (2) $STE - STJ \leq JE$ threshold value. Refer to the table 2. (3) Inclination from STJ to STM is $-1.0 \mu V/ms$ to $3.5 \mu V/ms$ . (4) The difference between the inclination from STJ to STM and the inclination from STE to Tpeak is $-0.8 \mu V/ms$ to $2.5 \mu V/ms$ . (5) T wave is not upward oriented. (6) The findings of possible left ventricular hypertrophy or left ventricular hypertrophy are not present. (7) The findings of anterior myocardial infarction are not present.	2 leads consecutive in V2, V3, V4, and V5
		All of the following conditions are satisfied: (1) $STJ \geq 450 \mu V$ and inclination from STJ to STM is positive. (2) $R \text{ amplitude} \times 1.8 < S \text{ amplitude}$ , or $S \text{ amplitude} = 0$	
	Injury	All of the following conditions are satisfied: (1) The criteria of possible anterior injury is met. (2) $STJ > \text{Total QRS amplitude}/4$	2 leads or more among V2, V3, V4, V5
Anteroseptal	Possible injury	All of the following conditions are satisfied: (1) The criteria of possible anterior injury is met. (2) The criteria of possible septal injury is met. (3) The criteria of possible anterolateral injury is not met.	-
	Injury	All of the following conditions are satisfied: (1) The criteria of possible anteroseptal injury is met. (2) The criteria of septal injury is met.	
Anterolateral	Possible injury	All of the following conditions are satisfied: (1) The criteria of possible anterior injury is met. (2) The criteria of possible lateral injury is met.	-
	Injury	All of the following conditions are satisfied: (1) The criteria of possible anteroseptal injury is met. (2) The criteria of anterior injury is met.	
		All of the following conditions are satisfied: (1) The criteria of possible anteroseptal injury is met. (2) The criteria of lateral injury is met.	

Portion	Finding	Judgement Criteria	Lead
Septal	Possible injury	All of the following conditions are satisfied: (1) $STJ \geq J$ threshold value. Refer to the table 1. (2) $STE - STJ \leq JE$ threshold value. Refer to the table 2. (3) Inclination from STJ to STM is $-1.0 \mu V/ms$ to $3.5 \mu V/ms$ . (4) The difference between the inclination from STJ to STM and the inclination from STE to Tpeak is $-0.8 \mu V/ms$ to $2.5 \mu V/ms$ . (5) T wave is not upward oriented. (6) The findings of possible left ventricular hypertrophy or left ventricular hypertrophy are not present. (7) The findings of septal myocardial infarction are not present.	V1 and V2
		All of the following conditions are satisfied: (1) $STJ \geq 450 \mu V$ and inclination from STJ to STM is positive. (2) $R \text{ amplitude} \times 1.8 < S \text{ amplitude}$ , or $S \text{ amplitude} = 0$	
	Injury	All of the following conditions are satisfied: (1) The criteria of possible septal injury is met. (2) $STJ > \text{Total QRS amplitude}/4$	
Lateral	Possible injury	All of the following conditions are satisfied: (1) $STJ \geq J$ threshold value. Refer to the table 1. (2) $STE - STJ \leq JE$ threshold value. Refer to the table 2. (3) Inclination from STJ to STM is $-1.0 \mu V/ms$ to $3.5 \mu V/ms$ . (4) The difference between the inclination from STJ to STM and the inclination from STE to Tpeak is $-0.8 \mu V/ms$ to $2.5 \mu V/ms$ . (5) T wave is not upward oriented. (6) The findings of lateral myocardial infarction are not present.	V5 and V6
		All of the following conditions are satisfied: (1) $STJ \geq 450 \mu V$ and inclination from STJ to STM is positive. (2) $R \text{ amplitude} \times 1.8 < S \text{ amplitude}$ , or $S \text{ amplitude} = 0$	
		All of the following conditions are satisfied: (1) $STJ > STLMT/2$ (2) T wave is not upward oriented. (3) The findings of lateral myocardial infarction are not present.	
	All of the following conditions are satisfied: (1) $STJ > \text{Total QRS amplitude}/6$ (2) The findings of lateral myocardial infarction are not present.		
	Injury	All of the following conditions are satisfied: (1) The criteria of possible lateral injury is met. (2) $STJ > \text{Total QRS amplitude}/4$	2 leads or more among I, aVL, V5, V6
Inferior	Possible injury	All of the following conditions are satisfied: (1) $STJ > STLMT/2$ (2) T wave is not upward oriented. (3) The findings of inferior myocardial infarction are not present.	II and aVF
		All of the following conditions are satisfied: (1) $STJ > \text{Total QRS amplitude}/6$ (2) The findings of inferior myocardial infarction are not present.	
	Injury	All of the following conditions are satisfied: (1) The criteria of possible inferior injury is met. (2) $STJ > \text{Total QRS amplitude}/4$	

## 6. CRITERIA OF FINDINGS

The threshold values of STJ and the threshold of potential difference between STJ and STM are as follows.

**Table 1: J threshold values**

Age	Lead	Threshold value ( $\mu\text{V}$ )
18 years or less	No Judgement	
19 to 29 years	V1	195
	V2	208
	V3, V4, V5 and V6	221
30 to 39 years	V1	180
	V2	192
	V3, V4, V5 and V6	204
40 years or more	V1	150
	V2	160
	V3, V4, V5 and V6	170

**Table 2: JE threshold values**

Lead	Threshold value ( $\mu\text{V}$ )
V1	200
V2	300
V3, V4, V5 and V6	320

STLMT is used for comparison of STJ. The values are below  
limb leads: STLMT = 0.3 mV

### NOTE

- (1) When injury or possible injury is found, the following findings are not printed out:
  - ST elevation, probably early repolarization
  - Early repolarization
  - ST elevation consistent with epicardial injury, pericarditis, or early repolarization
  - Nonspecific ST elevation
- (2) Relationship between injury portion and leads  
The portion of the injury is judged by the relationship between the standard 12 leads and heart portions as shown in the table below.
- (3) For the method of measuring upward oriented T, refer to p. 3.8 (c) of the ECAPS12C user's guide.
- (4) For the method of measuring total QRS amplitude and net QRS amplitude, refer to p. 3.7 of the ECAPS12C user's guide.

Portion	I	II	III	aVR	aVL	aVF	V1	V2	V3	V4	V5	V6
Anterior								*	*	*	*	
Anteroseptal							*	*	*	*	*	
Anterolateral	*				*			*	*	*	*	*
Septal							*	*				
Lateral	*				*						*	*
Inferior		*				*						

6

## Code and findings

Code	Findings [Criteria]	Judgement
4136	Possible anterior injury or acute infarct	abnormal ECG
4137	Anterior injury or acute infarct	abnormal ECG
4236	Possible anteroseptal injury or acute infarct	abnormal ECG
4237	Anteroseptal injury or acute infarct	abnormal ECG
4336	Possible anterolateral injury or acute infarct	abnormal ECG
4337	Anterolateral injury or acute infarct	abnormal ECG
4436	Possible septal injury or acute infarct	abnormal ECG
4437	Septal injury or acute infarct	abnormal ECG
4536	Possible lateral injury or acute infarct	abnormal ECG
4537	Lateral injury or acute infarct	abnormal ECG
4636	Possible inferior injury or acute infarct	abnormal ECG
4637	Inferior injury or acute infarct	abnormal ECG

6-4

## Subendocardial Ischemia

### Analysis criteria

(1) The new version has some changes in the criteria to withhold judgement of T wave abnormality (possible subendocardial ischemia).

The criteria changes are underlined>.

When any of the following findings is present, T wave abnormality (possible subendocardial ischemia) is not judged.

- Nonspecific intraventricular conduction block
- LBBB
- Myocardia infarction in relevant portion
- Subendocardial injury in relevant portion
- Possible subendocardial injury in relevant portion
- Possible acute pericarditis
- Marked ST depression consistent with subendocardial injury
- RVH (repolarization abnormality)
- LVH (repolarization abnormality)

(2) T wave abnormality and ST-T wave abnormality classification

Portion	Finding	Criteria	Lead	Judgement criteria
Anterior	T wave abnormality (Possible subendocardial ischemia)	modified T amplitude < -0.1 mV	V3 and V4	Cannot rule out anterior myocardial infarction is not satisfied.
			two or more leads among V2, V3 and V4	Cannot rule out anterior myocardial infarction is not satisfied. Complete RBBB is not satisfied. Age $\geq$ 16 years old
	T wave abnormality (subendocardial ischemia)	modified T amplitude < -0.5 mV	V3 or V4	Possible anterior subendocardial ischemia is satisfied
			and of V2, V3 and V4	Possible anterior subendocardial ischemia is satisfied Complete RBBB is not satisfied. Age $\geq$ 16 years old
Lateral	T wave abnormality (Possible subendocardial ischemia)	modified T amplitude < -0.1 mV	two or more leads among I, V4, V5 and V6	Cannot rule out lateral myocardial infarction is not satisfied.
		maximum R amplitude > 0.5 mV	aVL	Cannot rule out lateral myocardial infarction is not satisfied.
		modified T amplitude < -0.1 mV	at least two in I, aVL, V4, V5 and V6	



Portion	Finding	Criteria	Lead	Judgement criteria
Lateral	T wave abnormality (subendocardial ischemia)	modified T amplitude < -0.5 mV	any of I, V5 and V6	Possible lateral subendocardial ischemia is satisfied
		maximum R amplitude > 0.5 mV	aVL	Possible lateral subendocardial ischemia is satisfied
Inferior	T wave abnormality (Possible subendocardial ischemia)	modified T amplitude < 0 mV	II and aVF	Cannot rule out inferior myocardial infarction is not satisfied.
		modified T amplitude < -0.1 mV	II	
		modified T amplitude < 0 mV	II and aVF	Cannot rule out inferior myocardial infarction is not satisfied.
		modified T amplitude < -0.1 mV	aVF	
	Net QRS amplitude > 0 mV	aVF		
	T wave abnormality (subendocardial ischemia)	modified T amplitude < -0.5 mV	II or aVF	Possible inferior subendocardial ischemia is satisfied
Net QRS amplitude > 0 mV		aVF	Possible inferior subendocardial ischemia is satisfied	
	modified T amplitude < -0.5 mV	aVF		

**NOTE**

- When atrial fibrillation is judged together with subendocardial ischemia, “possible digitalis effect” is added to the findings.
- When both anterior ischemia and lateral ischemia are judged, anterolateral ischemia is printed out.

## 6. CRITERIA OF FINDINGS

### Code and findings

<b>Code</b>	<b>Findings [Criteria]</b>	<b>Judgement</b>
4164	T wave abnormality, possible anterior ischemia	abnormal ECG
41646	T wave abnormality, possible anterior ischemia or digitalis effect	abnormal ECG
4165	T wave abnormality, consistent with anterior ischemia	abnormal ECG
4364	T wave abnormality, possible anterolateral ischemia	abnormal ECG
43646	T wave abnormality, possible anterolateral ischemia or digitalis effect	abnormal ECG
4365	T wave abnormality, consistent with anterolateral ischemia	abnormal ECG
4564	T wave abnormality, possible lateral ischemia	abnormal ECG
45646	T wave abnormality, possible lateral ischemia or digitalis effect	abnormal ECG
4565	T wave abnormality, consistent with lateral ischemia	abnormal ECG
4664	T wave abnormality, possible inferior ischemia	abnormal ECG
46646	T wave abnormality, possible inferior ischemia or digitalis effect	abnormal ECG
4665	T wave abnormality, consistent with inferior ischemia	abnormal ECG

## Early Repolarization

### Analysis criteria

- (1) • QTc interval > 450 ms
- Nonspecific intraventricular conduction block
  - RBBB
  - LBBB
  - Myocardial infarction
  - LVH
- When any of these findings is recognized, the analysis processes of (2) and (3) are not executed.
- (2) When following two conditions are satisfied, the total number of the leads and the total amplitude of STJ are used to classify the ECG data as shown in (3).
- Chest lead: STJ and STM amplitude > 0.075 mV (age ≥ 20 years old)  
STJ and STM amplitude > 0.150 mV (age < 20 years old)
  - Limb lead: STJ and STM amplitude > 0.049 mV (age ≥ 20 years old)  
STJ and STM amplitude > 0.099 mV (age < 20 years old)

### (3) Classification

Findings	No. of leads	Total STJ amplitude	T waveform*
ST elevation consistent with subepicardial, pericarditis, or early repolarization	3 or more	≥ 0.45 mV (Age ≥ 20 years old) ≥ 0.90 mV (Age < 20 years old)	not upward
ST elevation, probably early repolarization	3 or more	≥ 0.45 mV (Age ≥ 20 years old) ≥ 0.90 mV (Age < 20 years old)	upward in more than half the number of leads satisfied (2)
Early repolarization	6 or more	≥ 0.45 mV	

\* For the method of measuring upward T waveform, refer to p. 3.8 (c).

### Code and findings

Code	Findings [Criteria]	Judgement
40302	ST elevation, probably early repolarization	borderline ECG
40303	Early repolarization	normal ECG
40371	ST elevation, consistent with subepicardial injury, pericarditis, or early repolarization	abnormal ECG

**Pericarditis****Analysis criteria**

- (1) • Nonspecific intraventricular conduction block
- RBBB
  - LBBB
  - Myocardial infarction
  - LVH
- When any of these findings is recognized, the analysis processes of (2), (3) and (4) are not executed.
- (2) The number of leads that satisfy the following conditions is counted.
- STJ and STM amplitude > 0.075 mV (in I, II, aVF)
  - STJ and STM amplitude > 0.09 mV (in V2 ~ V6)
- (3) The number of leads that satisfy the following conditions is counted.
- STJ and STM amplitude > 0.09 mV (in I, II, aVF)
  - STJ and STM amplitude > 0.11 mV (in V2 ~ V6)

## (4) Classification

Findings	STJ amplitude × 4	STJ and STM amplitude	No. of leads
Possible acute pericarditis (40304)	> T amplitude > 0 mV 4 or more leads among I, II, V4, V5, V6	> -0.1 mV All leads other than aVR	5 or more leads in (2)
Possible acute pericarditis (40305)	> T amplitude > 0 mV 4 or more leads among I, II, V4 V5, V6	> -0.1 mV All leads other than aVR	5 or more leads in (3)

**NOTE**

When “possible acute pericarditis” is found, the following findings are not printed out.

- ST elevation, probably early repolarization
- Early repolarization
- ST elevation, consistent with subendocardial injury, pericarditis, or early repolarization
- Nonspecific ST elevation

**Code and findings**

Code	Findings [Criteria]	Judgement
40304	Possible acute pericarditis	abnormal ECG
40305	Possible acute pericarditis	abnormal ECG

## T Wave Abnormality

### Analysis criteria

- (1) • Nonspecific intraventricular conduction block
- RBBB
  - LBBB
  - Myocardial infarction
  - Possible acute pericarditis
  - ST elevation
  - Subendocardial injury
  - Subendocardial ischemia
  - RVH (repolarization abnormality)
  - LVH (repolarization abnormality)

When any of these findings is recognized, the analysis processes of (2) and (3) are not executed.

### (2) Classification 1

Findings	Modified T amplitude*	Maximum R amplitude	ST abnormality	Tall T wave
Nonspecific T wave abnormality	< T min 2 or more leads among I, II, aVL, aVF, V3 to V6	> 0.5 mV In same leads as left		
Nonspecific ST&T abnormality	< T min 2 or more leads among I, II, aVL, aVF, V3 to V6	> 0.5 mV In same leads as left	Present	Absent

Here, T min is defined as follows:

- 1: Net QRS amplitude > 0 mV:  $T \text{ min} = 0.025 \text{ mV} + \text{Net QRS amplitude}/20$
- 2: Net QRS amplitude < 0 mV:  $T \text{ min} = 0.025 \text{ mV}$

\* For the method of measuring modified T amplitude, refer to p. 3.8 (d).

## 6. CRITERIA OF FINDINGS

### (3) Classification 2

When  $QTc \leq 0.45$  s, the following analysis is made.

Findings	T amplitude
Tall T wave, possible hyperkalemia	(1) $> 1.0 \text{ mV} > R \text{ amplitude}/2$ 3 or more leads among I, II, V1 to V6
	(2) $> 1.5 \text{ mV} > R \text{ amplitude}/2$ any lead among I, II, V1 to V6

### NOTE

When “nonspecific T wave abnormality” or “nonspecific ST&T wave abnormality” is recognized with “atrial fibrillation”, “probably digitalis effect” is added to the findings.

### Code and findings

Code	Findings [Criteria]	Judgement
4048	Nonspecific ST&T wave abnormality	normal ECG
40486	Nonspecific ST&T wave abnormality, probably digitalis effect	borderline ECG
4050	Tall T waves, possible hyperkalemia	abnormal ECG
4068	Nonspecific T wave abnormality	borderline ECG
40686	Nonspecific T wave abnormality, probably digitalis effect	borderline ECG

## Nonspecific ST Elevation

### Analysis criteria

- (1) • Nonspecific intraventricular conduction block
- RBBB
  - LBBB
  - Myocardial infarction
  - LVH (repolarization abnormality)
- When any of these findings is recognized, the analysis processes for (2), is not executed.
- (2) STJ  $\geq$  0.049 mV (Age  $\geq$  20 years old):
- STJ  $\geq$  0.099 mV (Age < 20 years old):
- 2 or more leads in I, II, III, aVF, V3, V4, V5 and V6
- STM  $\geq$  0.049 mV (Age  $\geq$  20 years old):
- STM  $\geq$  0.099 mV (Age < 20 years old):
- In same leads as above.
- STE  $\geq$  0.049 mV (Age  $\geq$  20 years old):
- STE  $\geq$  0.099 mV (Age < 20 years old):
- In same leads as above.

T wave is not upward\*: In same leads as above.

\* For measurement method, refer to p. 3.8 (c).

### Code and findings

Code	Findings [Criteria]	Judgement
4038	Nonspecific ST elevation	normal ECG

# *Section 6-5 Ventricular Hypertrophy*

Point Score System .....	6.5.2
Analysis Criteria for RVH .....	6.5.2
5120 Possible right ventricular hypertrophy .....	6.5.4
5130 Right ventricular hypertrophy .....	6.5.4
5134 Right ventricular hypertrophy, probably repolarization abnormality .....	6.5.5
Analysis Criteria for LVH.....	6.5.6
5211 Minimal voltage criteria for LVH, may be normal variant .....	6.5.7
5220 Possible left ventricular hypertrophy.....	6.5.7
5222 Moderate voltage criteria for LVH, may be normal variant .....	6.5.7
5233 Voltage criteria for LVH.....	6.5.7
5234 Left ventricular hypertrophy with repolarization abnormality .....	6.5.7



## 6. CRITERIA OF FINDINGS

### Point Score System

For judging cardiac hypertrophy, the major identifying features such as amplitude, QRS duration and repolarization abnormality are given points that vary with age and gender, and their total is used as the criteria.

### Analysis Criteria for RVH

#### Analysis criteria

- ① Complete RBBB
- ② Complete LBBB
- ③ QRS duration > 140 ms, net QRS duration < 0 mV (V1)

When any of these findings is recognized, the analysis of RVH is not executed.

- Judgement criteria for ages 17 years old and below

No.	Judgement	Points
1	In V1 : Q amplitude > 0.04 mV R amplitude > 0.7 mV S amplitude < 0.5 mV	3
2	In V1 : R amplitude > 1.0 mV*	3
3	In V6 : S amplitude > 0.4 mV*	3
4	Right ventricular conduction delay is present and S amplitude (V1) < 0.5 mV, and R' amplitude > 1.0 mV (V1) (Age >= 1 year old) or R' amplitude > 1.5 mV (V1)	2
5	T amplitude (V1) >= 0, Age <= 9 years old	2
6	(RV1 + SV5) / (SV1 + RV5) > 1.0 (Age >= 3 years old)	3
7	Right axis deviation is present	2
8	Moderate right axis deviation is present	1
9	Moderate right axis deviation is present or Right axis deviation is present, and Incomplete right bundle branch block is present.	1

#### NOTE

When the below condition is satisfied, right ventricular hypertrophy is not analyzed:

maximum S depth > 2 × maximum R height (V1)

The values marked with "\*" vary with age. Refer to p. 6.0.5 and 6.0.6.

- Judgement criteria for 18 years and over

No.	Judgement	Point
1	R or R' amplitude > 0.5 mV in V1	1
2	Net QRS amplitude > 0 mV in V1	1
3	Net QRS amplitude > 0.5 mV in V1	1
4	Net QRS amplitude < 0 mV in V5 or V6 and S amplitude > 0.5 mV in V5 or V6	1
5	QRS axis $\geq 90^\circ$ (no unidentified axis)	1
6	QRS axis $\geq 100^\circ$ (no unidentified axis)	1
7	QRS axis $\geq 110^\circ$ (no unidentified axis)	1
8	Possible Right atrial enlargement is present	1
9	S pattern is present	1
10	Age > 30 years old	1
11	Unidentified axis is present	1

However, when the following conditions are satisfied, right ventricular hypertrophy is not analyzed.

- ① Q, S, S' amplitude < 0.25 mV in I
- ② QRS axis <  $60^\circ$
- ③ Maximum S amplitude > 1 mV (in V1)

## 6. CRITERIA OF FINDINGS

Code	Findings [Criteria]	Judgement
5120	Possible right ventricular hypertrophy [Some/all of; prominent R in V1, late transition, RAD, RAE, SSS]	borderline ECG
5130	Right ventricular hypertrophy [Some/all of; prominent R in V1, late transition, RAD, RAE, SSS]	abnormal ECG

### Analysis criteria

Findings	Criteria
Possible RVH	Points $\geq$ 4
RVH	Points $\geq$ 6

### NOTE

When possible RVH is found, the following findings are not analyzed.

- Nonspecific intraventricular conduction delay
- RSR (QR) in lead V1/V2 consistent with right ventricular conduction delay
- Left posterior fascicular block
- Moderate right axis deviation
- Abnormal right axis deviation
- S1-S2-S3 pattern, consistent with pulmonary disease, RVH, or normal variant
- Low QRS voltage
- Low QRS voltage in limb leads
- Low QRS voltage in chest leads

Code	Findings [Criteria]	Judgement
5134	Right ventricular hypertrophy with repolarization abnormality [Some/all of; prominent R in V1, late transition, RAD, RAE, SSS, ST depression in lead V1, V2, V3]	abnormal ECG

#### Analysis criteria

Findings	Criteria
Right ventricular hypertrophy with repolarization abnormality	(1) RVH or possible RVH : present
	(2) ① STJ > STM > STE : in V1, V2, V3
	② • STM < -0.1 mV : in V1, V2, V3
	• STE < -0.1 mV
③ T amplitude < -0.1 mV : in V1, V2, V3	
(3) QRS duration < 120 ms*	

#### NOTE

The value marked with “\*” varies with age. Refer to p. 6.0.3.

## 6. CRITERIA OF FINDINGS

### Analysis Criteria for LVH

- (1) ① LBBB is present  
 ② QRS duration > 0.14 s\* and  
 Net QRS amplitude < 0 mV in V1
- ] When either of these findings is recognized, LVH is not analyzed.

\*The value marked with \* varies with age. Refer to p. 6.0.3.

- (2) Point score calculation and classification

The following four items are analyzed and the points are calculated.

Obtained point scores and the respective finding analysis criteria are used to identify the findings given on p. 6.5.7.

Item	Lead	Judgement	Points
1. R or R' amplitude	aVL	> 1.1 mV	2 points, 1 point added at every + 0.1 mV
2. S or S' amplitude	V1	> Threshold value***	2 points, 1 point added at every + 0.2 mV*
3. R or R' amplitude	V5	> Threshold value***	2 points, 1 point added at every + 0.2 mV*
4. R or R' amplitude + S or S' amplitude	V5 or V6 V1	> Threshold value***	2 points, 1 point added at every + 0.3 mV**

\* The values marked with “\*” are used for 17 years old and over.

1 point is added at every another 0.3 mV when 16 years old and under.

\*\* The values marked with “\*\*” is used for 17 years old and over.

1 point is added at every another 0.45 mV when 16 years old and under.

\*\*\* Threshold values for age groups are given below.

YEARS	R (V5)		S (V1)		R (V5 / V6) +S (V1)	
	MALE	FEMALE	MALE	FEMALE	MALE	FEMALE
3 to 5	3.50 mV	3.50 mV	2.50 mV	2.50 mV	6.00 mV	6.00 mV
6 to 8	3.50 mV	3.50 mV	2.75 mV	2.75 mV	6.00 mV	6.00 mV
9 to 12	3.75 mV	3.50 mV	3.00 mV	3.00 mV	6.00 mV	5.30 mV
13 to 16	3.75 mV	3.00 mV	3.75 mV	3.00 mV	6.00 mV	4.80 mV
17 to 19	3.75 mV	3.00 mV	3.75 mV	3.00 mV	6.00 mV	4.50 mV
20 to 23	3.75 mV	3.00 mV	3.75 mV	3.00 mV	5.55 mV	4.50 mV
24 to 29	3.50 mV	3.00 mV	3.50 mV	3.00 mV	5.00 mV	4.50 mV
>= 30	3.00 mV	3.00 mV	3.00 mV	3.00 mV	4.50 mV	4.50 mV

Code	Findings [Criteria]	Judgement
5211	Minimal voltage criteria for LVH, may be normal variant	borderline ECG
5220	Possible left ventricular hypertrophy	abnormal ECG
5222	Moderate voltage criteria for LVH, may be normal variant	borderline ECG
5233	Voltage criteria for LVH	abnormal ECG
5234	Left ventricular hypertrophy with repolarization abnormality	abnormal ECG

#### Analysis criteria

Findings	Criteria
Minimal voltage criteria for LVH	Point $\geq$ 2
Moderate voltage criteria for LVH	Point $\geq$ 3
Voltage criteria for LVH	Point $\geq$ 5
Possible left ventricular hypertrophy	(1) Point $\geq$ 2 (2) ① Maximum change point of inclination – QRS start point > 68 ms (V5)* ② Left atrial enlargement or Possible right atrial enlargement.
Left ventricular hypertrophy with repolarization abnormality	(1) Point $\geq$ 2 (2) Atrial fibrillation is not present (3) ① T amplitude (V1) > T amplitude (V6) + 0.2 mV ② • STE height < STJ height in any lead among I, aVL, V4, V5, V6 • STE height < -0.05 mV in same lead as above • R amplitude > 1.1 mV in same lead as above

#### NOTE

- When left ventricular hypertrophy with repolarization abnormality is found, “Nonspecific intraventricular conduction delay”, “Incomplete left bundle branch block” and “ST-T Abnormality” are not analyzed.
- The value marked with “\*” varies with age. Refer to p. 6.0.3.

# *Section 6-6 Atrial Enlargement, Abnormal Axis Deviation and Others*

6120	Possible right atrial enlargement .....	6.6.2
6130	Right atrial enlargement .....	6.6.2
6220	Possible left atrial enlargement .....	6.6.3
6230	Left atrial enlargement .....	6.6.3
7100	Abnormal right axis deviation .....	6.6.4
7102	Moderate right axis deviation .....	6.6.4
7200	Abnormal left axis deviation .....	6.6.4
7202	Moderate left axis deviation.....	6.6.4
7300	Indeterminate axis.....	6.6.4
7400	S1-S2-S3 pattern, consistent with pulmonary disease, RVH, or normal variant .....	6.6.6
7500	Abnormal QRS-T angle .....	6.6.7
8003	Consistent with pulmonary disease.....	6.6.8
8100	Low QRS voltage .....	6.6.9
8101	Low QRS voltage in limb leads .....	6.6.9
8102	Low QRS voltage in chest leads .....	6.6.9
8200	Dextrocardia .....	6.6.10
8304	Long QTc interval.....	6.6.11
8305	Short QTc interval .....	6.6.11
0101	Possible arm leads reversed, check lead requested .....	6.6.12
0102	ARTIFACT PRESENT .....	6.6.13
0103	CANNOT ANALYZE ECG.....	6.6.13
0104	ELECTRODE(S) FAILURE...Repeat ECG is required.....	6.6.13
0201	...Analysis based on intrinsic rhythm .....	6.6.13

## Atrial Enlargement

Code	Findings [Criteria]	Judgement
6120	Possible right atrial enlargement	borderline ECG
6130	Right atrial enlargement	abnormal ECG

### Analysis criteria

Findings	Criteria
Possible right atrial enlargement	(1) Heart rate < 120/minute (2) P amplitude > 0.25 mV in any lead among II, III, aVF, V1, V2
Right atrial enlargement	(1) Heart rate < 120/minute (2) P amplitude > 0.3 mV in any lead among II, III, aVF, V1, V2

### NOTE

- The reason for not analyzing when heart rate is over 120/minute is that at such a high heart rate, the increased P wave height may not definitely be caused by atrial enlargement.
- As shown in (2) in the table, when the P wave amplitude is 0.25 to 0.3 mV, the finding is marked “possible”.



Code	Findings [Criteria]	Judgement
6220	Possible left atrial enlargement	borderline ECG
6230	Left atrial enlargement	abnormal ECG

**Analysis criteria**

Findings	Criteria
Possible left atrial enlargement	(1) Negative P amplitude $< -0.1$ mV in V1 (2) Negative P area $\geq 4.0$ mV $\times$ ms in same lead as above
Left atrial enlargement	(1) Negative P amplitude $< -0.15$ mV in V1 (2) Negative P area $\geq 6.0$ mV $\times$ ms in same lead as above

**NOTE**

As shown in (1) in the table, when the negative P amplitude is between  $-0.15$  and  $-0.1$  mV, the finding is marked "possible".

## Abnormal Axis Deviation


Code	Findings [Criteria]	Judgement
7100	Abnormal right axis deviation	borderline ECG
7102	Moderate right axis deviation	normal ECG
7200	Abnormal left axis deviation	borderline ECG
7202	Moderate left axis deviation	normal ECG
7300	Indeterminate axis	atypical ECG

### Analysis criteria

Findings	Criteria
Moderate right axis deviation	$90^\circ < \text{QRS axis} \leq 100^\circ$ (RAXD1)*
Abnormal right axis deviation	$100^\circ < \text{QRS axis} \leq 270^\circ$ (RAXD2)*
Moderate left axis deviation	$-30^\circ \leq \text{QRS axis} < -20^\circ$ (LAXD1)*
Abnormal left axis deviation	$-90^\circ \leq \text{QRS axis} < -30^\circ$ (LAXD2)*
Indeterminate axis	Net QRS amplitude < 33% of total QRS amplitude in I, II, III

### NOTE

- The values marked with "\*" vary with age. Refer to p. 6.0.4.



## 6-6. ATRIAL ENLARGEMENT, ABNORMAL AXIS DEVIATION AND OTHERS

- When the measured values are at the boundary, the expression “moderate” is added to the findings.
- Since measuring the axis is irrelevant when the net QRS amplitude in I, II and III is smaller than 1/3 of the total QRS amplitude, the expression “indeterminate” is used.  
(For net QRS amplitude and total QRS amplitude, refer to p. 3.7.)

## 6. CRITERIA OF FINDINGS

Code	Findings [Criteria]	Judgement
7400	S1-S2-S3 pattern, consistent with pulmonary disease, RVH, or normal variant	borderline ECG

### Analysis criteria

When QRS duration  $\geq 120$  ms\*, S pattern is not analyzed.

The value marked with \* varies with age. Refer to p. 6.0.3.

Findings	Criteria
S1-S2-S3 pattern	(1) ① R < S (amplitude) in I, II, III ② • S > 0.3 mV in I • S > 0.4 mV in II • S > 0.7 mV in III (2) R' wave is not present in I, II, III (3) S > 0.2 mV in I, II, III (4) Age $\geq 16$ years old

Code	Findings [Criteria]	Judgement
7500	Abnormal QRS-T angle	borderline ECG

**Analysis criteria**

- Nonspecific intraventricular conduction block
- RBBB
- LBBB
- Myocardial infarction
- ST elevation
- Subendocardial injury
- Acute pericarditis
- RVH (with repolarization abnormality)
- LVH (with repolarization abnormality)
- Possible marked ST depression consistent with subendocardial injury
- Age < 1 year old

When any of these findings are recognized, the abnormal QRS-T angle is not analyzed.

Findings	Criteria	
	QRS axis – T axis	T axis
Abnormal QRS-T angle 1	> 60°	< 0°
2	< -60°	> 90°

## Others

Code	Findings [Criteria]	Judgement
8003	Consistent with pulmonary disease	abnormal ECG

### Analysis criteria

Findings	Criteria
Consistent with pulmonary disease	(1) QRS duration < 120 ms* (2) Total points >= point 4

### NOTE

Pulmonary disease is judged from the points representing the features of pulmonary diseases. The point scores are calculated as shown in the table below. This logic is not enough to identify pulmonary diseases, but if 4 or more of these features in the table are present in the ECGs, the probability of pulmonary diseases can be said to be fairly high.

Features	Points
1. Right atrial enlargement or possible right atrial enlargement	1
2. $-90^\circ \leq \text{QRS axis} < \text{LAXD2}^*$	1
3. $\text{RAXD1}^* < \text{QRS axis} \leq 270^\circ$	1
4. Indeterminate axis	1
5. S1-S2-S3 pattern	1
6. Low QRS voltage in limb leads	1
7. Low QRS voltage in chest leads	1
8. (1) Net QRS amplitude < 0 mV in V5 (2) R (and R') amplitude < 0.5 mV in V6	3

- The value marked with “\*” varies with age. Refer to p. 6.0.4.

Code	Findings [Criteria]	Judgement
8100	Low QRS voltage	abnormal ECG
8101	Low QRS voltage in limb leads	atypical ECG
8102	Low QRS voltage in chest leads	atypical ECG

#### Analysis criteria

Findings	Criteria
Low QRS voltage	(1) QRS duration < 0.12 s* (2) Total QRS amplitude < 0.5 mV in all limb leads (3) Total QRS amplitude < 1.0 mV in all chest leads
Low QRS voltage in limb leads	(1) QRS duration < 0.12 s* (2) Total QRS amplitude < 0.5 mV in all limb leads
Low QRS voltage in chest leads	(1) QRS duration < 0.12 s* (2) Total QRS amplitude < 1.0 mV in all chest leads

#### NOTE

The values marked with "\*" vary with age. Refer to p. 6.0.3.

## 6. CRITERIA OF FINDINGS

Code	Findings [Criteria]	Judgement
8200	Dextrocardia	atypical ECG

### Analysis criteria

Findings	Criteria
Dextrocardia	(1) $90^\circ < \text{QRS axis} \leq 270^\circ$ (2) $90^\circ < \text{P axis} \leq 270^\circ$ (3) PR interval $\geq 110 \text{ ms}^*$ (4) ① Q wave is present in I ② • Q wave is not present in I • R amplitude $< 0.15 \text{ mV}$ in I (5) R amplitude $< 0.5 \text{ mV}$ in V6 (6) Net QRS amplitude $\leq 0 \text{ mV}$ in V6 (7) P amplitude $< 0.02 \text{ mV}$ in V6 (8) Negative P amplitude $< -0.02 \text{ mV}$ in V6

### NOTE

- Refer to “0101 Possible arm leads reversed, check lead requested”.
- The values marked with “\*” varies with age. Refer to p. 6.0.2.



Code	Findings [Criteria]	Judgement
8304	Long QTc interval	abnormal ECG
8305	Short QTc interval	abnormal ECG

**Analysis criteria**

- Nonspecific intraventricular conduction block
- RBBB
- LBBB
- Myocardial infarction
- ST elevation (Subendocardial injury)
- Subendocardial ischemia
- Possible acute pericarditis
- Marked ST depression (Subendocardial injury)
- RVH (with repolarization abnormality)
- LVH (with repolarization abnormality)
- Age < 1 year old

When any of these findings are recognized, long QTc interval is not analyzed.

Findings	Criteria
Long QTc interval	QTc interval > 0.45 s
Short QTc interval	(1) QTc interval < 0.36 s (2) heart rate < 140/minute

## 6. CRITERIA OF FINDINGS

Code	Findings [Criteria]	Judgement
0101	Possible arm leads reversed, check lead requested	---

### Analysis criteria

Findings	Criteria
ARM LEADS REVERSED	(1) $90^\circ < \text{QRS axis} \leq 270^\circ$ (2) $90^\circ < \text{P axis} \leq 270^\circ$ (3) PR interval $\geq 110 \text{ ms}^*$ (4) ① Q wave is present in I ② • Q wave is not present in I • R amplitude $< 0.15 \text{ mV}$ in I (5) ① R amplitude $\geq 0.5 \text{ mV}$ at V6 ② R amplitude $>$ S amplitude in V6 ③ P amplitude $\geq 0.02 \text{ mV}$ in V6 ④ Negative P amplitude $\geq -0.02 \text{ mV}$ in V6

### NOTE

- In properly recorded ECGs, the P and QRS waves of Lead I are not expected to appear in negative simultaneously. If the type of QRS wave is Qr (or rSr'), the most probable cause is either dextrocardia or reversed arm lead electrodes. If the V6 lead shows the typical upward-oriented waveform, the probability of reversed electrode is high, otherwise, the probability for dextrocardia is high.
- Pulmonary disease tends to show right axis deviation in both P wave and QRS wave, and the rS type occurs. The same is true with other diseases which show right axis deviation. A Qr type rarely appears in myocardial infarction, but an inverted P is not expected to appear at the same time.
- Both "Arm leads reversed" and "Dextrocardia" have "inverted P and QRS wave" as criteria, but in the actual decision, the Qr or rSr type of QRS wave is taken as a significant factor.
- The value marked with "\*" varies with age. Refer to p. 6.0.2.

6-6. ATRIAL ENLARGEMENT, ABNORMAL AXIS DEVIATION AND OTHERS

Code	Findings [Criteria]	Judgement
0102	ARTIFACT PRESENT	---

Artifact on waveform.

Artifact may cause incorrect analysis. Remove artifacts and record ECG again.

Code	Findings [Criteria]	Judgement
0103	CANNOT ANALYZE ECG	---

ECG could not be analyzed because

- Artifact on waveform, or
- Heart beats are too small, or
- The ECG has less than 3 heart beats which the instrument can measure.

Code	Findings [Criteria]	Judgement
0104	ELECTRODE(S) DETACHED...Repeat ECG is required.	---

Electrode detachment

Attach electrodes correctly and record ECG again.

Code	Findings [Criteria]	Judgement
0201	...Analysis based on intrinsic rhythm	---

Both artificial pacemaker rhythm and intrinsic heart beats are in ECG recording and the instrument analyzed intrinsic heart beats.

Refer to p. 6.1.9 “(3) With pacemaker”.

# *Appendix Modified Minnesota Code*

General.....	A.2
Code List .....	A.3
Priority of Code Printing .....	A.7
Detailed Criteria.....	A.8

## General

The Minnesota code is a classification of adult ECG waveforms according to a certain criteria for the purpose of disease research. This criteria is adopted by WHO.

In Japan, the Japanese Association for Cerebro-Cardiovascular Disease Control arranged classification criteria of Minnesota code under consideration for the Japanese body size (modified Minnesota code).

The electrocardiograph uses a modified Minnesota code which has been adopted by the Japanese Association for Cerebro-Cardiovascular Disease Control.

The modified Minnesota code after exercise is coded by comparing the code for rest ECG with the code for After Exercise ECG. Therefore, to print out the modified Minnesota code after exercise, be sure to analyze the rest ECG before attempting exercise tests.

Up to 8 codes can be printed at the same time. Modified Minnesota code for rest ECG is 1-n to 9-n, modified Minnesota code after exercise is 11-n to 16-n. For classification criteria, refer to page 3.2 and later.

Some instruments do not print out modified Minnesota codes. For the procedure to print out modified Minnesota code, refer to the electrocardiograph operator's manual.

### NOTE

- The modified Minnesota code classifies ECG waveform according to a criteria different from that of the ECAPS12C. Therefore, the analysis result of ECAPS12C and the classification by the modified Minnesota code may differ.
- The ECAPS12C classifies averaged waveforms according to the modified Minnesota code.

## Code List

- |  |     |        |
|--|-----|--------|
|  | 1-0 | Normal |
|--|-----|--------|
1. Q and QS patterns
 

1-1 Class1	1-1-1	1-1-5
	1-1-2	1-1-6
	1-1-3	1-1-7
	1-1-4	
1-2 Class2	1-2-1	1-2-5
	1-2-2	1-2-6
	1-2-3	1-2-7
	1-2-4	1-2-8
1-3 Class3	1-3-1	1-3-4
	1-3-2	1-3-5
	1-3-3	1-3-6
  
  2. QRS Axis Deviation
 

2-1	Left axis deviation
2-2	Right axis deviation
2-3	Right axis deviation
2-4	Extreme axis deviation
2-5	Indetermination axis
  
  3. High Amplitude R Waves
 

3-1	Left: High amplitude R waves
3-2	Right: High amplitude R waves*
3-3	Left: Moderate high amplitude R waves**

\* S amplitude > R amplitude in either V2, V3, V4, V5 or V6  
R amplitude > 0.5 mV and R amplitude > S amplitude in V1  
(If criteria for 3-2 is met, 7-3 is not coded.)

\*\* For easy understanding, 3-3 is divided into 3-3-1, 3-3-2 and 3-3-3.
  
  4. ST Junction (J) and Segment Depression
 

4-1*	4-3
4-2	4-4

\* For easy understanding, 4-1 is divided into 4-1-1 and 4-1-2 according to degree of phenomenon.
  
  5. T-Wave Items
 

5-1	5-4
5-2	5-5
5-3	
  
  6. A-V Conduction Defect
 

6-1	Complete (third degree) A-V block
6-2*	Partial (second degree) A-V block
6-3	First degree A-V block

- 6-4\* Wolff-Parkinson-White Pattern (WPW)
  - 6-5 Short P-R interval
  - 6-8 Artificial pacemaker
  - \* 6-2 is divided into 6-2-1 (Mobitz type II) and 6-2-3 (Wenckebach's phenomenon)
  - \* 6-4 is divided into 6-4-1 (WPW pattern, persistent) and 6-4-2 (WPW pattern, intermittent)
7. Ventricular Conduction Defect
- 7-1 Complete left bundle branch block
  - 7-2 Complete right bundle branch block
  - 7-3 Incomplete right bundle branch block
  - 7-4 Intraventricular block
  - 7-5 R-R' pattern
  - 7-6 Incomplete left bundle branch block
8. Arrhythmias
- 8-1\* with frequent premature complexes
  - 8-2 Ventricular tachycardia
  - 8-3\* Atrial fibrillation or atrial flutter
  - 8-4 Supraventricular tachycardia
  - 8-5 Ventricular rhythm
  - 8-6 Atrioventricular (A-V) nodal rhythm
  - 8-7 Sinus tachycardia
  - 8-8 Sinus bradycardia
  - 8-9 other arrhythmias
  - \* 8-1 is divided into 8-1-1 (supraventricular) and 8-1-2 (ventricular)
  - \* 8-3 is divided into 8-3-1 (Atrial fibrillation) and 8-3-2 (Atrial flutter).
9. Miscellaneous Items
- 9-1\* Low QRS amplitude
  - 9-2 ST elevation
  - 9-3-1 Tall P waves
  - 9-3-2 Widened P waves
  - 9-4-1 Transition zone
  - 9-4-2 Transition zone
  - 9-5 Tall T waves
  - 9-6 Dextrocardia
  - 9-8 Measurement failure because of technical problems (electrode detached, arm leads reversed)
  - \* 9-1 is divided into 9-1-1 (limb leads and chest leads), 9-1-2 (limb leads), 9-1-3 (chest leads)
10. ST Items after Exercise
- 11-1 No ST Junction (J) and Segment Depression code (4-x) is present in rest ECG and 4-1 appears after exercise.
  - 11-2 No ST Junction (J) and Segment Depression code (4-x) is present in rest ECG and 4-2 appears after exercise.

- 11-3 No ST Junction (J) and Segment Depression code (4-x) is present in rest ECG and 4-3 appears after exercise.
  - 11-4 No ST Junction (J) and Segment Depression code (4-x) is present in rest ECG and 4-4 appears after exercise.
  - 11-5 Any ST Junction (J) and Segment Depression code present in rest ECG changes to a lower ST Junction (J) and Segment Depression code after exercise.
  - 11-6 Any ST Junction (J) and Segment Depression code present in rest ECG changes to a higher ST Junction (J) and Segment Depression code after exercise.
  - 11-7 Any ST Junction (J) and Segment Depression code is present in rest ECG and the same ST Junction (J) and Segment Depression code appears after exercise.
  - 11-8 Any ST Junction (J) and Segment Depression code is present in rest ECG and no ST Junction (J) and Segment Depression code appears after exercise.
  - 11-9 both 11-x and 9-8
11. T Items after Exercise
- 12-1 No T-wave Items code (5-x) is present in rest ECG and 5-1 appears after exercise.
  - 12-2 No T-wave Items code (5-x) is present in rest ECG and 5-2 appears after exercise.
  - 12-3 No T-wave Items code (5-x) is present in rest ECG and 5-3 appears after exercise.
  - 12-4 Any T-wave Items code present in rest ECG changes to a lower T-wave Items code after exercise.
  - 12-5 Any T-wave Items code present in rest ECG changes to a higher T-wave Items code after exercise.
  - 12-6 Any T-wave Items code is present in rest ECG and the same T-wave Items code appears after exercise.
  - 12-7 Any T-wave Items code is present in rest ECG and no T-wave Items code appears after exercise.
  - 12-8 both 12-x and 9-8
12. A-V Conduction after Exercise
- 13-1 No A-V Conduction Defect code (6-x) is present in rest ECG and Complete (third degree) A-V block appears after exercise.
  - 13-2 No A-V Conduction Defect code (6-x) is present in rest ECG and Partial (second degree) A-V block appears after exercise.
  - 13-3 No A-V Conduction Defect code (6-x) is present in rest ECG and First degree A-V block (P-R interval  $\geq 0.22$  sec) appears after exercise.
  - 13-4 No A-V Conduction Defect code (6-x) is present in rest ECG and 6-4-1 (WPW pattern, persistent) appears after exercise.
  - 13-5 6-3 or 6-2 present in rest ECG changes to any other A-V Conduction Defect code after exercise.
  - 13-6 Any A-V Conduction Defect code is present in rest ECG and the same code appears after exercise.



- 13-7 Any A-V Conduction Defect code is present in rest ECG and no A-V Conduction Defect code appears after exercise.
13. Ventricular Conduction after Exercise
- 14-1 No Ventricular Conduction Defect code (7-x) is present in rest ECG and complete left bundle branch blocks (7-1) appears after exercise.
- 14-2 No Ventricular Conduction Defect code (7-x) is present in rest ECG and complete right bundle branch block (7-2) appears after exercise.
- 14-3 No Ventricular Conduction Defect code (7-x) is present in rest ECG and incomplete right bundle branch block (7-3) appears after exercise.
- 14-4 No Ventricular Conduction Defect code (7-x) is present in rest ECG and intraventricular block appears after exercise.
- 14-5 7-1, 7-2, 7-3, 7-4, 7-5 or 7-6 present in rest ECG changes to another Ventricular Conduction Defect code after exercise.
- 14-6 Any Ventricular Conduction Defect code is present in rest ECG and the same code appears after exercise.
- 14-7 7-1, 7-2, 7-3, 7-4, 7-5 or 7-6 is present in rest ECG and no Ventricular Conduction Defect code appears after exercise.
14. Arrhythmia after Exercise
- Sinus arrhythmia, sinus tachycardia (8-7) and sinus bradycardia (8-8) are excluded from judgement, but when sinus bradycardia (8-8) in rest ECG changes to atrioventricular (A-V) nodal rhythm (8-6) after exercise, this is included regardless transient or persistent.
- 15-1 No Arrhythmias code (8-x) in rest ECG and any of Arrhythmias code appears after exercise.
- 15-2 Any Arrhythmias code (8-x) present in rest ECG changes to another Arrhythmias code after exercise.
- 15-3 Any Arrhythmias code (8-x) is present in rest ECG and the same code appears after exercise.
- 15-4 Any Arrhythmias code (8-x) is present in rest ECG and no Arrhythmias code appears after exercise.
15. Miscellaneous Items after Exercise
- 16-1 ST elevation code (9-2) is not present in rest ECG and 9-2 (ST elevation) appears after exercise.
- 16-2 ST elevation code (9-2) is present in rest ECG and it appears after exercise too.
- 16-3 ST elevation code (9-2) is present in rest ECG and it does not appear after exercise.

#### NOTE

The judgement criteria in this section are extracted from the handbook by the Japan Association for Cerebro-Cardiovascular Disease Control. These criteria were not written for computer analysis. Where the computer needs more detailed criteria for analysis, Nihon Kohden has added some criteria.

## Priority of Code Printing

- 1) For each top level group 1 to 5, only the highest priority item is coded.  
Example (1) 1-1-1 suppresses 1-2-4 and 1-3-2 In this case only 1-1-1 is coded and 1-2-4 and 1-3-2 are ignored.  
Example (2) When 1-1-1, 1-2-4, 2-1 and 3-1 are present in analysis, only 1-1-1, 2-1 and 3-1 are coded.
- 2) For 6 to 9, two or more codes are coded together among each 6 to 9. All the present codes in analysis result are coded.  
Example (1) When 8-1 and 8-3 are present, both 8-1 and 8-3 are coded. For 6, 7, 8 and 9, each group is coded separately.
- 3) Q and QS patterns (1-x)  
Do not code if 6-4 or 7-1 is present. Q wave < 0.1 mV is not coded.
- 4) QRS axis deviation (2-x)  
Do not code if 6-4, 7-1, 7-2, 7-4 or 9-1 is present.
- 5) High amplitude R waves (3-x)  
Do not code if 6-4, 7-1, 7-2 or 7-4 is present.
- 6) ST Junction (J) and Segment Depression (4-x)  
Do not code if 6-4, 7-1, 7-2 or 7-4 is present.
- 7) T-wave Items (5-x)  
Do not code if 6-4, 7-1, 7-2 or 7-4 is present.

## Detailed Criteria

### 1. Q and QS patterns

#### 1-1 Class1

- 1-1-1 Q/R amplitude ratio  $\geq 1/3$ , plus Q duration  $\geq 0.03$  sec in any of leads I, II, V2, V3, V4, V5, V6.
- 1-1-2 Q duration  $\geq 0.04$  sec in any of lead I, II, V1, V2, V3, V4, V5, V6.
- 1-1-3 Q duration  $\geq 0.04$  sec, plus R amplitude  $\geq 0.3$  mV in lead aVL.
- 1-1-4 Q duration  $\geq 0.05$  sec in lead III, plus Q amplitude  $\geq 0.1$  mV in lead aVF.
- 1-1-5 Q duration  $\geq 0.05$  sec in lead aVF.
- 1-1-6 QS pattern when initial R-wave is present in adjacent lead to the right on the chest wall, in any of leads V2, V3, V4, V5, V6.
- 1-1-7 QS pattern in all of leads V1-V4, V1-V5 or V1-V6.

#### 1-2 Class2

- 1-2-1 Q/R amplitude ratio  $\geq 1/3$ , plus Q duration  $\geq 0.02$  sec and  $< 0.03$  sec in any of leads I, II, V2, V3, V4, V5, V6.
- 1-2-2 Q duration  $\geq 0.03$  sec and  $< 0.04$  sec in any leads of I, II, V2, V3, V4, V5, V6.
- 1-2-3 QS pattern in lead II.
- 1-2-4 Q duration  $\geq 0.04$  sec and  $< 0.05$  sec in lead III, plus a Q-wave  $\geq 0.1$  mV in lead aVF.
- 1-2-5 Q duration  $\geq 0.04$  sec and  $< 0.05$  sec in lead aVF.
- 1-2-6 Q amplitude  $\geq 0.5$  mV in leads III or aVF.
- 1-2-7 QS pattern in all of leads V1, V2, and V3.
- 1-2-8 Initial R amplitude decreasing to 0.2 mV or less in every beat (and absence of codes 3-2, 7-2, 7-3) between any of leads V2 and V3, V3 and V4, V4 and V5, V5 and V6. All beats in the lead immediately to the right on the chest must have an initial R  $> 0.2$  mV.

#### 1-3 Class3

- 1-3-1 Q/R amplitude ratio  $\geq 1/5$  and  $< 1/3$  plus Q duration  $\geq 0.02$  sec and  $< 0.03$  sec in any of leads I, II, V2, V3, V4, V5, V6.
- 1-3-2 QS pattern in lead V1 and V2. (Do not code in the presence of 3-1, 6-4, 7-1.)
- 1-3-3 Q duration  $\geq 0.03$  sec and  $< 0.04$  sec, plus R amplitude  $\geq 0.3$  mV in lead aVL.
- 1-3-4 Q duration  $\geq 0.03$  sec and  $< 0.04$  sec in lead III, plus a Q-wave  $\geq 0.1$  mV in lead aVF.
- 1-3-5 Q duration  $\geq 0.03$  sec and  $< 0.04$  sec in lead aVF.
- 1-3-6 QS pattern in each of leads III and aVF.

2. QRS Axis Deviation
  - 2-1 Left axis deviation:  
QRS axis from  $-30^{\circ}$  through  $-90^{\circ}$  in leads I, II, III. (The algebraic sum of major positive and major negative QRS waves must be positive in I, zero or negative in II and negative in III.)
  - 2-2 Right axis deviation:  
QRS axis from  $+120^{\circ}$  through  $-150^{\circ}$ . (the algebraic sum of major positive and major negative QRS waves must be negative in I, zero or positive in III, and the sum in I must be one-half or more of the sum in III.)
  - 2-3 Right axis deviation:  
QRS axis from  $+90^{\circ}$  through  $+119^{\circ}$ . (The algebraic sum of major positive and major negative QRS waves must be zero or negative in I and positive in II and III.)
  - 2-4 Extreme axis deviation (usually S1-S2-S3 pattern)  
QRS axis from  $-91^{\circ}$  through  $-149^{\circ}$ . (The algebraic sum of major positive and major negative QRS waves must be negative in each of leads I, II and III.)
  - 2-5 Indeterminate axis:  
(The algebraic sum of major positive and major negative QRS waves is zero in each of leads I, II and III.)
  
3. High amplitude R waves
  - 3-1 Left: High amplitude R waves  
R amplitude  $> 2.6$  mV in either V5 or V6, or R amplitude  $> 2.0$  mV in either I, II, III, aVF, or R amplitude  $> 1.2$  mV in lead aVL.
  - 3-2 Right: High amplitude R waves  
R amplitude  $\geq 0.5$  mV and R amplitude  $\geq$  S amplitude in lead V1, when S amplitude is  $>$  R amplitude somewhere to the left on the chest of V1 (Code as 3-2 if criteria for 3-2 and 7-3 are met. Do not code 1-2-8 in the presence of 3-2)
  - 3-3-1 Left: Moderate high amplitude R waves  
R amplitude  $\geq 1.5$  mV but  $\leq 2.0$  mV in lead I.
  - 3-3-2 Left: Moderate high amplitude R waves  
R amplitude in lead V5 plus S amplitude in lead V1  $\geq 3.5$  mV
  - 3-3-3 Left: Moderate high amplitude R waves  
R amplitude in lead V6 plus S amplitude in lead V1  $\geq 3.5$  mV.
  
4. ST Junction (J) and Segment Depression
  - 4-1-1 STJ depression  $\geq 0.2$  mV and ST segment is horizontal or downward sloping in any of leads I, II, aVL, aVF, V1, V2, V3, V4, V5, V6.
  - 4-1-2 STJ depression  $\geq 0.1$  mV but  $< 0.2$  mV, and ST segment is horizontal or downward sloping in any of leads I, II, aVL, V1, V2, V3, V4, V5, V6.
  - 4-2 STJ depression  $\geq 0.05$  mV and  $< 0.1$  mV and ST segment is horizontal or downward sloping in any of leads I, II, aVL, V2, V3, V4, V5, V6.

- 4-3 No STJ depression as much as 0.05 mV but ST segment downward sloping and segment or T-wave nadir  $\geq 0.05$  mV below P-R baseline, in any of leads I, II, aVL, V2, V3, V4, V5, V6.
  - 4-4 STJ depression  $\geq 0.1$  mV, and ST segment upward sloping or U-shaped, in any of leads I, II, aVL, aVF, V1, V2, V3, V4, V5, V6.
5. T Wave Items
- 5-1 T amplitude negative 0.5 mV or more in any of leads I, II, V2, V3, V4, V5, V6, or in lead aVL when R amplitude is  $\geq 0.5$  mV, or in lead aVF when QRS is mainly upright.
  - 5-2 T amplitude negative or diphasic (negative-positive or positive-negative type) with negative phase at least 0.1 mV but not as deep as 0.5 mV in any of leads I, II, V2, V3, V4, V5, V6, or in lead aVL when R amplitude is  $\geq 0.5$  mV, or in lead aVF when QRS is mainly upright.
  - 5-3 T amplitude zero (flat) or negative or diphasic (negative-positive type) with less than 0.1 mV negative phase, in any leads of I, II, V3, V4, V5, V6, or in lead aVL when R amplitude is  $\geq 0.5$  mV. (Do not code in lead aVF.)
  - 5-4 T amplitude positive and T/R amplitude ratio  $< 1/20$  in any of leads I, II, aVL, V3, V4, V5, V6: R wave amplitude must be  $\geq 1.0$  mV.
  - 5-5 T amplitude positive and T/R amplitude ratio  $< 1/10$  and  $\geq 1/20$  in any of leads I, II, aVL, V3, V4, V5, V6: R wave amplitude must be  $\geq 1.0$  mV.
6. A-V Conduction Defect
- 6-1 Complete (third degree) A-V block (permanent or intermittent) in any lead.
  - 6-2 Partial (second degree) A-V block (permanent or intermittent, 2:1 or 3:1 block, Wenckebach's phenomenon) in any lead.
    - 6-2-1 Mobitz Type II
    - 6-2-3 Wenckebach's phenomenon
  - 6-3 P-R (P-Q) interval  $\geq 0.22$  sec in any of leads I, II, III, aVL, aVF.
  - 6-4-1 Wolff-Parkinson-White Pattern (WPW)
 

P-R interval  $< 0.12$  sec, plus QRS duration  $\geq 0.12$  sec, plus R peak duration  $\geq 0.06$  sec, coexisting in the same beat and present in all beats in any of leads I, II, aVL, V4, V5, V6.
  - 6-4-2 Wolff-Parkinson-White Pattern (WPW), intermittent
 

P-R interval  $< 0.12$  sec, plus QRS duration  $\geq 0.12$  sec, plus R peak duration  $\geq 0.06$  sec, coexisting in the same beat and present in some beats in any of leads I, II, aVL, V4, V5, V6.
  - 6-5 Short P-R interval:
 

P-R (P-Q) interval  $< 0.12$  sec in all beats of any two of leads I, II, III, aVL, aVF. (in the absence of 8-6, 8-7)
  - 6-8 Artificial pacemaker:
 

Artificial pacemaker pulse is present.

7. Ventricular Conduction Defect
  - 7-1 Complete left bundle branch block (Do not code in presence of 6-4)  
QRS duration  $\geq 0.12$  sec in any of leads I, II, III, aVL, aVF, plus R peak duration  $\geq 0.06$  sec and a codable Q-wave is not present in any of leads I, II, aVL, V5, V6.
  - 7-2 Complete right bundle branch block (Do not code in the presence of 6-4)  
QRS duration  $\geq 0.12$  sec in any of leads I, II, III, aVL, aVF, plus: R'  $> R$  or R peak duration  $\geq 0.06$  sec in V1 or V2.
  - 7-3 Incomplete right bundle branch block:  
QRS duration  $< 0.12$  sec in each of leads I, II, III, aVL, aVF, and R'  $> R$  in either of leads V1, V2. (Code as 3-2 if criteria for 3-2 is met.)
  - 7-4 Intraventricular block (Do not code in presence of 6-4, 7-1, 7-2)  
QRS duration  $\geq 0.12$  sec.
  - 7-5 R-R' pattern:  
R-R' pattern which do not meet the criteria of 7-2 and 7-3 in V1 or V2
  - 7-6 Incomplete left bundle branch block:  
QRS duration  $\geq 0.10$  sec but  $< 0.12$  sec and codable Q-wave is not present in leads I, aVL and either V5 or V6.
  
8. Arrhythmias
  - 8-1 with frequent atrial, junctional or ventricular premature complexes (10% or more of recorded complexes)
    - 8-1-1 with frequent supraventricular premature complexes (10% or more of recorded complexes)
    - 8-1-2 with frequent ventricular premature complexes (10% or more of recorded complexes)  
If supraventricular or ventricular is undetermined, code as 8-1.
  - 8-2 Ventricular tachycardia ( $\geq 100$ /min)
  - 8-3-1 Atrial fibrillation
  - 8-3-2 Atrial flutter
  - 8-4 Supraventricular tachycardia ( $\geq 100$ /min)
  - 8-5 Ventricular rhythm ( $\leq 100$ /min)
  - 8-6 atrioventricular (A-V) nodal rhythm ( $\leq 100$ /min)  
Negative P in aVF, and P-R interval  $\leq 0.12$  sec in any of leads I, II, III, aVL, aVF.
  - 8-7 Sinus tachycardia ( $\geq 100$ /min)
  - 8-8 Sinus bradycardia ( $\leq 50$ /min)
  - 8-9 other arrhythmias
  
9. Miscellaneous Items
  - 9-1 Low QRS amplitude:
    - 9-1-1 Low QRS amplitude: QRS peak-to-peak amplitude  $< 0.5$  mV in all beats in each of leads I, II, III, and  $< 1.0$ mV in all beats in each of leads V1, V2, V3, V4, V5, V6.
    - 9-1-2 Low QRS amplitude: QRS peak-to-peak amplitude  $< 0.5$  mV in all beats in each of leads I, II, III.

## APPENDIX MODIFIED MINNESOTA CODE

- 9-1-3 Low QRS amplitude: QRS peak-to-peak amplitude  $< 1.0$  mV in all beats in each of leads V1, V2, V3, V4, V5, V6.
- 9-2 ST segment elevation  $\geq 0.1$  mV in any of leads I, II, III, aVL, aVF, V5, V6, or  $\geq 0.2$  mV in any of leads V1, V2, V3, V4. (Do not code in the presence of 6-4, 7-1, 7-2 or 7-4)
- 9-3-1 P-wave amplitude  $\geq 0.25$  mV in any of leads II, III, aVF.
- 9-3-2 P-wave duration  $\geq 0.10$  sec in any of leads I, II, aVL.
- 9-4-1 QRS Transition zone to the right of V3 on the chest wall. (Do not code in the presence of 6-4, 7-1, 7-2 or 7-4)
- 9-4-2 QRS Transition zone at V4 or to the left of V4 on the chest wall. (Do not code in the presence of 6-4, 7-1, 7-2 or 7-4)
- 9-5 T wave amplitude  $> 1.2$  mV in any of leads I, II, III, aVL, aVF, V1, V2, V3, V4, V5, V6. (Do not code in the presence of 6-4, 7-1, 7-2 or 7-4)
- 9-6 Dextrocardia
- 9-8 Measurement failure because of electrode detachment or reversed arm leads.



Manufacturer

## NIHON KOHDEN CORPORATION

1-31-4 Nishiochiai, Shinjuku-ku

Tokyo 161-8560, Japan

Phone +81 (3) 5996-8036

Fax +81 (3) 5996-8100

### Sales

---

#### USA


##### NIHON KOHDEN AMERICA, INC.

90 Icon Street, Foothill Ranch, CA 92610, USA

Phone +1 (949) 580-1555

Fax +1 (949) 580-1550

#### Europe

 European Representative

##### NIHON KOHDEN EUROPE GmbH

Raiffeisenstrasse 10

D-61191 Rosbach v.d.H., Germany

Phone +49 6003 827-0

Fax +49 6003 827-599

##### NIHON KOHDEN ITALIA S.r.l.

Via Fratelli Bronzetti 28

I-24124 Bergamo, Italy

Phone +39 035 219 543

Fax +39 035 232 546

##### NIHON KOHDEN FRANCE SARL

8, rue Francois Delage, 94230 Cachan, France

Phone +33 1 49 08 05 50

Fax +33 1 49 08 93 32

##### NIHON KOHDEN IBERICA S.L.

C/Ulises 75A

E-28043 Madrid, Spain

Phone +34 917 161 080

Fax +34 913 004 676

##### NIHON KOHDEN UK LTD

Tolworth Tower, Ewell Road, Surbiton

Surrey KT6 7EL, UK

Phone +44 20-8390-8622

Fax +44 20-8390-4675

#### Asia

##### SHANGHAI KOHDEN MEDICAL ELECTRONIC INSTRUMENT CORPORATION

7th Floor, Dawning Centre Tower A No.500

Hongbaoshi Road, Changning District

Shanghai 201103, China

Phone +86 (21) 6270-0909

Fax +86 (21) 6270-9700

##### Beijing Branch

Room 1701, East Ocean Centre

No. 24A JianGuoMenWai Street, Beijing, 100004

Phone 010-6515-5750

Fax 010-6515-5758

##### Guangzhou Branch

Room 2514, Yian Plaza

No. 33 Jian She Liu Ma Road, Guangzhou, 510060

Phone 020-8363-3737

Fax 020-8363-3807

##### NIHON KOHDEN SINGAPORE PTE LTD

1 Maritime Square, #10-34 (Lobby C), Harbour Front Centre

Singapore 099253

Phone +65 6376-2210

Fax +65 6376-2264

##### NIHON KOHDEN KOREA, INC.

Hannam Tower Annex Bldg. Suite 203

730 Hannam-dong, Yongsan-gu, Seoul, Korea 141-210

Phone +82 (2) 3273-2310

Fax +82 (2) 3273-2352

### Production

---

Reagent Production

##### NIHON KOHDEN FIRENZE S.r.l.

Via Torta 72/74

I-50019 Sesto Fiorentino Firenze, Italy

Phone +39 055 3045 1

Fax +39 055 308548

Contact information is accurate as of May 2013. Visit [www.nihonkohden.com](http://www.nihonkohden.com) for the latest information.

The model and serial number of your instrument are identified on the rear or bottom of the unit. Write the model and serial number in the spaces provided below. Whenever you call your representative concerning this instrument, mention these two pieces of information for quick and accurate service.

Model \_\_\_\_\_

Serial number \_\_\_\_\_

Your Representative