
Dräger

ENGLISH TEXT
P.T.O.

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Babytherm® 4200

Von Dräger: Babytherm® 4200

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
OPERATING INSTRUCTIONS

Important Notice

For correct and effective use of the device, and to avoid hazards, we wish to point out the following:

- 1 Any use of the device requires precise knowledge and observation of these operating instructions.
- 2 The device is intended only for the purposes specified in the operating manual or for purposes specified in writing by Drägerwerk AG.
- 3 The device must be serviced by experts at regular intervals. An official record of inspections and maintenance services must be kept. We recommend that inspections and maintenance services be carried out by the Technical Customer Service of your Dräger subsidiary or distributor. Regular inspection is best ensured by entering into an Inspection Service Contract with the Customer Service

of your Dräger subsidiary or distributor.

- 4 As regards medical-engineering equipment with an electrical connection, particular attention is to be paid to VDE Regulation 0750, which states that the equipment may only be repaired by us or by an agency expressly authorized by us to do so.
- 5 Only original Dräger spare parts may be employed in maintenance and repair work.
- 6 For reasons of safety, pressure reducers should be overhauled at least every 6 years.
- 7 This device is not approved for use in explosion-hazard areas.
- 8 The reliable function of the device becomes the responsibility of the owner or operator in all cases where

the device has been improperly maintained or repaired by persons not employed by the Dräger subsidiary or distributor in question or when it has been used in a manner which does not conform to the specified conditions of use.

Drägerwerk AG, its subsidiaries and distributors are not liable for damage which is a result of not observing this Important Notice. Warranty and liability stated or implied in the sales and delivery conditions of Drägerwerk AG, its subsidiaries or distributors are not extended by the above Important Notice.

We also wish to point out that the national recommendations, regulations and laws governing the use of technical equipment should be observed.

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Contents

| | Page | | Page |
|--|------|---|------|
| Important Notice | 2 | 12 Infusion stand | 9 |
| Babytherm 4200, a modern pre-fabricated system | 3 | 13 Restrainer | 9 |
| Intended Use | 3 | 14 Terminal box | 9 |
| Composition and Instructions for Use | 4 | 15 Swivel table | 9 |
| 1 Basis Babytherm 4200 unit | 4 | 16 Small operating lamp | 9 |
| 2 Radiant heater | 5 | 17 Ventilation in the Babytherm | 10 |
| 3 Skin temperature regulator | 5 | 18 Bronchial aspirator | 10 |
| 4 Hinged cover | 7 | 19 Oxygen distributor with cylinder connection | 10 |
| 5 Injector | 7 | 20 Oxygen distributor (connection to a central supply system) | 10 |
| 6 Water nebulizer | 7 | 21 Pillow to facilitate intubation .. | 10 |
| 7 Medicament nebulizer | 8 | 22 Timer | 10 |
| 8 Oxygen inhalation device | 8 | 23 Accessories for gas supply | 10 |
| 9 Hood | 8 | Care | 11 |
| 10 Rail | 8 | Order List | 13 |
| 11 Cabinet | 8 | | |

Titel: Babytherm 4200
used for intensive care

Babytherm® 4200 — an advanced modular system

Babytherm 4200 makes it possible to carry out

- simple thermal treatment of infants
- oxygen or aerosol therapy
- resuscitation
- blood exchange
- intensive care

In addition, Dräger offers an extensive range of equipment for accommodation and treatment of premature and sick babies:

- **PP 800** for initial care resuscitation in the labour room
- **Intensive Care Transport Incubator 5400** for the transport of premature and sick babies; also available in standard version.
- **Incubator 7310** for nursing premature babies.

- **Intensive Care Incubator 7510** for accommodation of premature babies and infants requiring special intensive care.
- And the heated cot **Babytherm 4200**

another important link in the equipment range. By adding various additional equipment the simple heated cot can be extended to a very versatile intensive-care unit.

There are considerable differences in construction, function and, consequently, in the uses of the incubators and the Babytherm 4200, which are compared in the following table:

| Incubator 7310 and 7510 | Babytherm 4200 |
|---|--|
| Closed-circuit air circulating system | No air circulating system |
| Bacterial filter for the fresh air supplied | No filter |
| Precise regulation of the incubator air temperature or skin temperature | Regulation of the mattress temperature; Adjustable power of the radiant heater |
| Temperature alarms | No temperature alarms |
| Restricted access to the baby (hand ports) | Total access (all side panels hinged or removeable) |
| Dimensions of the resting surface: 37 x 67 cm | Dimensions of the resting surface: 46 x 78 cm, i. e. it is suitable for babies up to six months of age |

Intended use

The components mentioned in the proposals composition are described on pages 4 to 10. The Babytherm can be used as:

a) Heated cot

The mattress heater protects newborn babies and infants from loss of heat. The baby is easily accessible for nursing and care when the Plexiglas side panels are swung down. The resting surface can be easily swivelled, so that the baby can be positioned with the head up or down. For incubator babies the unit is used for the transition stage.

Suggestion for composition: Basic unit 1 with mattress heater (Fig. 1).

b) For oxygen and aerosol therapy

When the Plexiglas cover is placed over the unit, it is an almost sealed compartment, into which an oxygen/air mixture can be administered via an injector, or aerosols can be applied by way of a nebulizer. Frequent cleaning which is required when sputum-dissolving agents (e. g. Tacholiquin) are nebulized is facilitated since the cot consists of easily accessible smooth panels.

Suggestion for composition: 1, 4, 5, 6, 8, 9, 10, 23.

c) For intensive care

The manifold accessories make the Babytherm 4200 particularly suitable for intensive care of infants.

Suggestion for composition: Basic unit 1 and all the accessories from 2 to 23, (cf. cover photo).


Fig. 1 Basic unit with mattress heater

40-430

Composition and instructions for use

1 Basic unit Babytherm 4200 with mattress heater (Fig. 1)

The basic unit – equipped with mattress heater – consists of a trolley with four castors (two of which can be locked) and the resting surface, which can be tilted forwards and backwards and arrested in four different inclined positions (max. 20°). The angle can be adjusted with one hand by pulling the lever behind the handle and then tilting the resting surface (Fig. 1a).

The two Plexiglas side panels can be swung down for feeding and nursing the baby. Care must be taken that the side panels are later properly re-adjusted, whereby the locks must noticeably engage. The front panels which are also made of Plexiglas can be removed from their holders and exchanged against each other if required. The padding on the front panels can be detached.

Bores in the four corner posts are intended as ducts for infusion hoses and other hoses. Ventilation hoses can be inserted through the nebulizer or injector supporting flange. The trolley holds two stainless steel instrument trays.

The mattress heater is used for heating the Babytherm mattress.

Check that type of current and voltage of the unit correspond to the ratings of line voltage, prior to taking the heater into operation.

Once the heater had been removed it must be pushed into the compartment provided for in the Babytherm until it engages (Fig. 1b).

The mains cable is then plugged in and the switch set to setting »I«. The white pilot light goes on, indicating that the apparatus is ready for use. The mattress temperature required can be set between 25°C and 39°C on the control knob.

If the mattress is covered with a pillow, the temperature set is reached after a warming-up period of 1 hour. This value is kept constant by a built-in thermostat.

The yellow pilot light of the heater goes on and off periodically. When the maximum adjustable temperature is exceeded the temperature is cut-off and a red excess-temperature alarm lamp as well as an alarm buzzer are switched on by a temperature limiting switch which operates independent of the thermostat. The temperature limiting switch is not self-resetting, e. g. the heater remains cut off even after cooling down.

Owing to the long warming-up period it


Fig. 1a Basic unit dimensions


Fig. 1b Mattres heater

Technical data of the mattress heater

Type of current: alternating current, 50/60 cps Power: 180 W
Mains voltage: as ordered Classification: Type B 

The device is suppressed in accordance with VDE regulations 0875 and carries the interference suppression symbol

is recommended to keep the heating switched on when the unit is kept at the ready.

The air in the Babytherm is warmed only to a negligible extent by the mattress heater, even when it is closed by putting on the hinged cover 4. The baby should, therefore, be covered with a quilt.

The bores, threaded connections and pins necessary for fitting the accessories described in the following

are fitted to the basic unit. Almost all the accessories can be fitted at different points on the basic unit, so that the composition suited to the purpose envisaged can be achieved.

The accessories can be transferred from Babytherm to Babytherm, so that it is not necessary to equip each Babytherm 4200 with a complete set of accessories on wards which have several Babytherms.

2

Radiant heater

The radiation of heat from the radiant heater can protect the baby against loss of heat when it is lying naked and uncovered in the Babytherm, for example during routine treatment and in resuscitation after birth or during blood exchange, X-rays and in many instances of intensive care.

The radiant heater can be used alone or together with the mattress heater.

The two supporting rods of the radiant heater are pushed into two corner posts on the Babytherm. The heater itself can be swivelled through 90° on both sides. This is particularly necessary when X-rays are being taken.

Before putting the heater into operation, it must be ensured that the type of current and voltage available in the Department are suitable for the heater.

The mains cable should be plugged in, the switch set to position I. The power of the heater can be set between 80 and 300 W using a control knob. The heating power set is kept constant by means of a regulator even when the mains voltage fluctuates considerably. The heating power to be chosen depends on the ambient temperature and on the baby's body temperature, which should be taken at frequent intervals.

Experience has shown that the settings 9 and 10 which are marked in colour do only have to be used when the device must be rapidly warmed up or when an undercooled baby has to be brought to normal temperature.

The electric socket which is arranged at the side is envisaged for connection of a skin-temperature regulator, which has to be ordered separately.

The radiant heater can be used both when the top of the Babytherm is open and when it is covered by the hinged cover 4.

a) Babytherm open at the top

Where the treatment of the baby takes only a short time, the Babytherm is left open. The heat radiated from the radiant heater acts directly on the baby. The ambient air is heated only to a very slight extent.

b) Babytherm closed

If the naked baby has to remain for some time in the Babytherm, the apparatus should be closed by the hinged cover 4. The heat radiated from the heater no longer acts directly on the baby, but heats the Plexiglas panels of the hinged


Fig. 2 Radiant heater

cover, which in turn give up their heat to the air in the Babytherm. The temperature must be watched on the thermometer, which is fitted on the hinged cover.

Technical data of the radiant heater

Type of current: alternating current, 50/60 cps

Mains voltage: as requested

Power consumption: 320 W

Max. heating capacity: 300 W

Test symbol: 

The apparatus is suppressed in accordance with VDE regulations 0875 and carries the interference suppression symbol

Classification: 

Type B

Symbol  means:

The electric socket arranged at the side of the device is solely intended for connection of the skin-temperature regulator (item 3)

3

Skin-temperature regulator

The skin-temperature regulator permits regulation of the infant's skin temperature. Given certain conditions it is also possible to regulate the temperature within the Babytherm even with the hinged cover closed. In addition, the regulator can be used for measuring the skin temperature.

The regulator can be supplementary screwed onto any existing radiant heater of the Babytherm 4200 (Fig. 3).

Usage

Plug mains cable of skin-temperature regulator into connecting socket of radiant heater. Switch on both units regulator and heater.

a) Regulation of skin temperature

(Fig. 3)

Switch setting »measure + regulate«.

Note!

The skin-temperature regulator may **not** be used on infants who are under **shock**. They may have fever (rectal more than 37°C) in spite of low skin temperature. Subjected to skin-temperature regulation they are liable to get too hot. In such a case the heating output of the radiant heater has to be manually adjusted.

The rotary knob for the heating output of the radiant heater is set to stage 8 (190 Watt). This makes for sufficient warmth on the one hand, while it prevents overheating of the infant on the other.


Fig. 3 Radiant heater with skin-temperature regulator

The desired skin-temperature value is adjusted on the screw marked »set value« by using the screwdriver supplied; the value is displayed on the upper scale of the measuring instrument. Attach temperature sensor to regulator and tighten by means of the lock nut. The sleeve at the end of the cable contains the temperature measuring element which is fixed in place between the navel and the sternum of the baby, by means of adhesive tape¹⁾. Care must be taken that the measuring point is at all times directly subjected to the radiant heater and that the Babytherm is open at the top (hinged cover removed).

¹⁾ e. g. »Leukoflex« and »Leukoplast waterproof« manufactured by the Beiersdorf Company; »Micropore« and »Blenderm« manufactured by the 3M Company.

The actual skin temperature of the infant, which is indicated on the lower »actual value« scale of the instrument is brought to the set desired value by the regulator (quasi-proportional regulation).

The pilot light of the radiant heater lights up during the heating process. The light flashes in the proportional range of regulation (set value equal to or similar to actual value).

The lamp lights continuously if the actual value is smaller than the set value, whereas the lamp does not light if the actual value is greater than the set value.

Sudden changes in the infant's skin temperature are liable to occur when the infant is fed, falling asleep or waking up, or being attended to.

Although the regulating system reacts immediately to any disturbance by switching the heater on or off, it takes some time for the infant's temperature to return to the desired value. Deviations of a few tenths of a degree between the set value and indicated skin temperature, are liable to occur any time and can be regarded as normal.

The temperature sensor must be firmly attached to the infant's skin. Inadvertent loosening of the sensor will, as a result, not keep the infant's skin temperature to the desired value, but, for example, the temperature of the colder mattress, onto which the sensor might have dropped. In this case, the heat supplied would rise to a value higher than that required by the infant, but it will not exceed the heating level set on the radiant heater.

The settings 9 and 10 marked in colour (full heating power) should only be switched when an infant suffering from hypothermia is to be brought to normal body temperature or when the ambient air is not sufficiently warm.

b) Air temperature regulation

Switch setting »measure + regulate«.

The skin-temperature regulator can also be used for regulating the air temperature. However, the hinged covers of the Babytherm 4200 must be closed for this purpose.

Set rotary knob of radiant heater to position 10.

Push skin-temperature sensor from above through a rubber duct in the hinged cover (near thermometer) roughly 2 cm deep. The desired air temperature is to be adjusted at the »set-value« screw. The radiant heat from the heater no longer acts directly on the infant, but heats the Plexiglas of the hinged cover, which, in turn, transmits its heat to the air inside the Babytherm.

It should be noted that it takes a comparatively long time for the Babytherm to warm up. The ambient temperature may not be more than 8°C max. below the desired Babytherm temperature.

c) Measuring the skin temperature

If the skin-temperature regulator is used for this purpose, the switch on the front panel of the device is set to »measure«. In this case the radiant heater need not be switched on. The device may even be completely removed from the radiant heater. It is then connected to the mains by means of a cable which has to be ordered separately.

Temperatures within the range of 32°C and 38°C can be read from the actual value scale of the measuring instrument.

Technical data of the skin-temperature regulator

Temperature measuring range and desired-value adjusting range: 32...38°C

Type of current: alternating current, 50/60 cps

Mains voltage: 100 to 240 V, according to request

Power: approx. 5 W

Test symbol:  The apparatus is suppressed in accordance with VDE regulations 0875 and bears the interference suppression symbol

Classification:  Type B

Symbol  means:

- 1) For power supply the device is connected to the respective socket of the radiant heater (item 2)
- 2) The socket underneath the display instrument is solely intended for connection of the skin-temperature sensor or the test plug

Dimensions: Height x width x depth = 90 x 180 x 140 mm Weight: 1.4 kg

4

Hinged cover

The Babytherm, which is normally open at the top, can be closed with the hinged cover. This is necessary, for example, when carrying out oxygen or aerosol therapy. Uniform heating of the air in the Babytherm is possible using the radiant heater.

The cover consists of a frame and two Plexiglas panels which can be swung upwards. Even with the cover closed, the expired air does not lead to an inadmissibly high CO₂ concentration inside the Babytherm. The bores and clearances in the unit ensure adequate air exchange. A thermometer makes it possible to measure the air temperature inside the Babytherm.


Fig 4 Hinged cover

5

Injector

One front panel of the Babytherm features a supporting flange, into which the injector for oxygen therapy can be pushed. An oxygen/air mixture, the concentration of which is adjustable to 30, 40 or 50%, can be administered into the unit through the oxygen-fed injector. In this case the Babytherm has to be closed with the hinged cover 4.

The administered oxygen is humidified

by means of a nebulizer, so that, depending on the temperature conditions, a relative humidity of approx. 50–70% is obtained inside the Babytherm. See item 23 for the required connecting parts.

It is recommended to set the following values at the oxygen flowmeter:

| Injector setting (% O ₂) | 30 | 40 | 50 |
|--------------------------------------|----|----|----|
| O ₂ flow (L/min) | 5 | 10 | 15 |


Fig. 5 Injector

6

Water nebulizer

To carry out aerosol therapy, the Babytherm is closed with the hinged cover 4. The water nebulizer is either pushed into the supporting flange at the front panel (Fig. 6), or fixed to the rail 10, whereby the aerosol is routed into the inside of the Babytherm (Fig. 6a) by way of a corrugated hose (cf. 23f). The advantage of fixing the nebulizer at the rail is, that there is no noise inside the Babytherm.

Water or, if required, aqueous medicaments, are nebulized by means

of compressed gas (compressed air or oxygen). When using oxygen as compressed gas, the oxygen concentration of the gas administered to the Babytherm can be adjusted at the nebulizer to 100, 70, 50 or 40%.

Up to 65 g/h water can be nebulized depending on the operating conditions of the nebulizer. Double or triple the nebulizer capacity can be achieved by using the heater rod 2M 13125, which is pushed into the nebulizer in place of the plug (Fig. 6b) thus also heating the gas. The water in the reservoir of the nebulizer is heated by the rod heater.

The desired water temperature is adjustable and kept constant by a thermostat.

This may be a rough guide for the temperature adjustment: The greater the rate of the set flow, and the greater the distance between nebulizer and patient (great distance means a hose between patient and nebulizer) the higher can the temperature be set on the rod heater.

For connecting parts required for oxygen or compressed air supply cf. item 23.


Fig. 6 Water nebulizer

28 796


Fig. 6a Water nebulizer fitted to suspension rail

29 630


Fig. 6b Water nebulizer with rod heater

29 831

7

Medicaments nebulizer

The medicaments nebulizer which consists of the nebulizer element M 22282 and an inhalation mask, can, by using the holder 2M 11170, be placed in front of the infant's face such, that a great part of nebulized medicaments is inhaled.

For connecting parts required for oxygen or compressed air supply cf. item 23.


Fig. 7 Medicaments nebulizer

40 393

8

Oxygen inhalation device

This inhalation device which consists of a hose and a funnel makes it possible for oxygen to be rapidly administered to the baby in case of acute asphyxia. A flow of approximately 4 L/min is set on the oxygen flowmeter.

For connecting parts required for oxygen supply cf. item 23.


Fig. 8 Oxygen inhalation device

40 395

9

Hood

As with the oxygen inhalation device for infants described under 8, humidified oxygen can be rapidly administered to the baby by using the transparent Plexiglas hood. Since the baby rests in an almost sealed compartment, the oxygen content of inspired air can be brought to nearly 100%.

For connecting parts required for oxygen supply cf. item 23. An oxygen flow of approx. 10 L/min has to be set at the flowmeter.

If a particularly high oxygen concentration is desired, the leak between the hood and the infant's neck has to be sealed by using a nappy.

Dimensions of the hood (without handle): length 21 cm, width 20 cm, height 16.5 cm.


Fig. 9 Oxygen inhalation device with hood

40 396

10

Rail

Both ends of the trolley can be fitted with a rail for suspension of various auxiliary equipment e. g. flowmeter or bronchial aspirator.


Fig. 10 Rail

40 392

11

Cabinet 2M 14780, small

One or two small cabinets can be fitted between the rods of the trolley. The cabinet 2M 14780 contains one large and two small drawers, which can be pulled out completely on smoothly running telescope tracks. The drawers can easily be removed for cleaning.

If work is to be carried out sitting at the front of the Babytherm (for example, during blood exchange), only one cabinet should be fitted (Fig. 11), to leave room for the doctor's legs. If two cabinets are provided, it may be advantageous to fit them such that the drawers open from opposite sides of the Babytherm. If two doctors are working on opposite sides of the Babytherm, each of them then has a cabinet available for his instruments.


Fig. 11 Cabinet, small

40 424

Cabinet 2M 12910, large (Fig. 11a)

This cabinet, which consists of a single element, occupies the entire space between the rods of the trolley.

It is divided into a cabinet with door and two large drawers, which can be pulled out completely on smoothly running telescope tracks.


Fig. 11a Cabinet, large

40 423

**12
Infusion Stand**

To carry out blood exchange and infusions, a rod, the height of which is adjustable, is pushed into one of the corner posts of the Babytherm.


Fig. 12 Infusion stand

**13
Restrainer**

During resuscitation, blood exchange and other treatment, the baby's arms and legs must be restrained. For this purpose, the baby is fastened to two rods, which are fitted between the corner posts of the Babytherm by means of cushioned straps with Velcro fastening. Eight straps (Order No. 2M 13898) are supplied with the restrainer. The straps can be washed in the same way as nylon.


Fig. 13 Restrainer

Fitting the restrainer rods

At one end of the rods is a retaining ring, which can be removed with a small screw-driver. This rod end is pushed right into the bore in the left-hand post. The two pins on the other end of the rod are then pushed into the small holes in the right-hand post. The retaining ring is then forced into the groove in the rod again (Fig. 13a).


Fig. 13a Fitting supporting rods

**14
Terminal Box**

A terminal box equipped with six sockets permits simultaneous connection of various pieces of electrical equipment (e. g. heater, infusion pump, monitoring instruments). The box has a pocket, into which any excess cable can be pushed. The terminal box is suspended from the trolley.

It is for reasons of maximum safety for both patient and operator, that the current conducted to the earth by way of the nonfused earth conductor of the cable lead may not exceed the value of 0.5 mA at all the auxiliary sockets (IEC 601/1, chart 4). This conducted earth current is the total of all earth conducted currents from all attached units.


Fig. 14 Terminal box

**15
Swivel table**

When treating the baby (e. g. in blood exchange, etc.), the doctor needs a »working-table«. One, or, if required, two hinged tables, adjustable in height, can be fitted to the Babytherm trolley. The table has a stainless steel top. When not in use, it can be swung under the baby's resting surface.


Fig. 15 Swivel table

**16
Small Operating Lamp**

The small operating lamp gives a delimited, concentrated area of illumination with cool, colour-corrected achromatised light, which is useful for diagnosis and treatment of the baby. The lamp is pushed into one of the corner posts of the Babytherm. It can be swivelled and tilted, so that any point in the Babytherm can be illuminated.


Fig. 16 Small operation lamp

Bulb replacement

There are two screws on the outer rim of the lamp. These have to be unscrewed, using a screwdriver. Now remove complete filter insert, whereupon the bulb is accessible. It can be removed by pressing it in and turning it anticlockwise (Swan-holder). Upon assembly of the lamp care must be taken that the rubber ring behind the chrome ring of the filter insert is seated on the Plexiglas disc. There is only one position in which the filter insert can be fitted into the lamp body.

Technical data of the operating lamp

Mains voltage: 220 V
Power consumption: 50 W
Brightness: approx. 3000lux
Classification: type B

17

Ventilation in the Babytherm 4200

There are two units available for this purpose: either Babylog 1 or Babylog 2 (cf. respective operating manuals).

18

Bronchial aspirator

The aspirator is suspended from the rail 10 which is attached to the Babytherm.

The required vacuum is generated by an adjustable ejector, driven either by compressed air or oxygen. Secretion is extracted from the infant's air passages by using a thin catheter. The rinsing jar is used to suck water or a detergent through the catheters.

For connecting parts required for oxygen or compressed air supply cf. item 23.

19

Oxygen distributor with cylinder connection

This device can be suspended from the rail 10 which is attached to the Babytherm. Three built-in oxygen outlets are ready to take over oxygen supply of various units simultaneously (e. g. ventilator, aspirator, oxygen inhalation device).

Oxygen for the distributor is provided for either by connection to a central supply unit or via an oxygen cylinder, which has to be ordered separately. Pressure gauge and pressure reducer are integrated in the distributor. For oxygen connecting hoses cf. item 23.

20

Oxygen distributor

Essentially, this device corresponds to the oxygen distributor with cylinder connection 19; however, cylinder connection is dispensed with; oxygen is solely supplied from the central system.

Take great care in the handling of oxygen

- Increased fire risk!
- Smoking, naked flames and fire are prohibited!
- Prevent sparks!
- Keep oxygen fittings and seals free from oil or grease!
- Open cylinder valves slowly!
- Valve handwheel may only be operated manually!

The cylinder valve is a precision instrument which can be easily damaged if force is applied. A leaking cylinder valve must be expertly repaired.


Fig. 17 Babylog 2 attached to rail


Fig. 18 Bronchial aspirator


Fig. 19 O₂ distributor with cylinder connection


Fig. 20 O₂ distributor (connection to a central supply system only)

- Electrical apparatus with the exception of »intrinsically safe« units may not be used in the interior of the Babytherm! Combustible liquids such as alcohol or ether may neither be used nor stored in the Babytherm.

21

Pillow to facilitate intubation

A foam rubber pillow is pushed under the baby's head for intubation. A hollow for the head and a neck support facilitate the desirable backwards bend of the head.


Fig. 21 Pillow to facilitate intubation

22

Timer with holder

This timer is designed for temporal monitoring of resuscitation measures that have to be carried out. It is fitted with a minute and centre second hand. Scale diameter: 95 mm; green starting knob and red stop-knob, white knob for resetting.

The spring mechanism has a running time of roughly 20 h.

The timer can be clamped to the supporting rod of the infusion stand (Fig. 22), to the small operating lamp (Fig. 16) or to one of the supporting rods of the radiant heater.


Fig. 22 Timer (clamped to infusion stand)

23

Gas supply

Fig. 23 (Page 12) illustrates the required accessories (hoses, flowmeter, humidifier etc.) for the gas supply of various auxiliary equipment.

Care

Cleaning and disinfection of the Babytherm should be carried out at least once a week, and, of course, upon every patient change. For auxiliary equipment cf. the respective operating manuals.

Cleaning

Empty all the jars containing liquid. The mattress and both head pads have to be removed from the Babytherm for cleaning purposes. If applicable, disassemble water nebulizer, humidifier nebulizer, medicaments nebulizer and bronchial aspirator. Remove all water absorbing parts such as nappies, cellulose and the like. Dirt on the device and the connected parts is removed using a cloth soaked in detergent¹⁾. The parts are then to be thoroughly dried to prevent the growth of bacteria.

¹⁾ Please pay attention to the attached compatibility list for Plexiglas 2M 8016.

Disinfection in the Dräger Aseptor[®]

The device has to be prepared for disinfection as described under Section »Cleaning«. All the parts must be dry since otherwise unpleasant odours are liable to occur following disinfection. Prior to disinfection the Babytherm must have been out of use for approx. 1½ hours, to allow for cooling down, or else, efficient disinfection may be questionable. The Babytherm is wheeled into the Aseptor. Keep cabinet doors and drawers open as well as the Plexiglas panels and the hinged cover to grant free access of the disinfecting gas. The mattress and the two head pads are placed onto the depositing tray of the Aseptor.

Following disinfection the mattress heater should be switched on for approx. 1 hour with the hinged covers open in order to prevent even the

slightest touch of unpleasant odour. Prior to the next operation the unit has to be subjected to a functional check.

Disinfection with liquid disinfectants

The device has to be prepared for disinfection as described under Section »Cleaning«. It is recommended to place small parts e. g. injector, as well as small components of auxiliary equipment, into a disinfecting solution. Following disinfection, the parts have to be thoroughly rinsed with water and then dried.

Wipe Babytherm with a disinfecting solution¹⁾. Prior to re-use the disinfectant must have been completely evaporated.

Check device for proper function.

¹⁾ Cf. Plexiglas compatibility list 2M 8016 which is supplied with the unit.


Fig. 23 Accessories for gas supply

Order List

| Item | Name and description | Order No. | Item | Name and description | Order No. |
|------|--|--|------|---|--|
| 1 | Babytherm 4200 , basic unit with mattress heater 220 V (temperature adjustment up to 39°C and 2nd safety thermostat) Supplement for previous Babytherm 4200 versions: Mattress heater, 220 V with temperature adjustment up to 39°C and 2nd safety thermostat | 2M 18702 2M 18280 | 23 | Accessories for gas supply | |
| 2 | Radiant heater | 2M 13500 | a | O₂ pressure reducer with manual connection O₂ cylinder, 11 litres (R ^{3/4"} , charged, 200 bar) Trolley for 11-litre cylinder Cylinder jacket | D 12340 B 2710 M 20315 M 3691 |
| 3 | Skin-temperature regulator (not calibrated) to be placed onto radiant heater Mains cable only required if skin-temperature regulator is used without radiant heater, e. g. for measurement purposes | 82 00 008 83 01 037 | c | O ₂ connecting hose for cylinder operation: O₂ connecting hose, 3 m O₂ connecting hose, 5 m | M 14565 M 17263 |
| 4 | Hinged cover | 2M 13369 | b | O ₂ connecting hose for operation from central supply system: O₂ connecting hose, 1 m O₂ connecting hose, 3 m O₂ connecting hose, 5 m | M 22478 M 22346 M 22347 |
| 5 | Injector | 2M 14190 | d | O ₂ -compressed-air connecting hose for operation from central supply system: O₂-compressed-air connecting hose, 2 m O₂-compressed-air connecting hose, 3 m O₂-compressed-air connecting hose, 5 m | M 22360 M 22361 M 22362 |
| 6 | Water nebulizer Holder for water nebulizer Rod heater Thermometer Sleeve for thermometer | 2M 12836 2M 15915 2M 13125 2M 13259 2M 15666 | e | O₂-compressed-air connecting hose, 1 m O₂-compressed-air connecting hose, 1.5 m O₂-compressed-air connecting hose, 3 m O₂-compressed-air connecting hose, 5 m | M 17615 M 17616 M 17617 M 17618 |
| 7 | Medicaments nebulizer, comprising Nebulizer element Set of inhalation masks (2 ea) Holder for nebulizer element | M 22282 M 23456 2M 11170 | f | Hose, 1 m (set of 5 ea) | 2M 17478 |
| 8 | Inhalation device, infants | M 11950 | g | O₂ Flowmeter 0...15 L/min | M 21535 |
| 9 | Hood | 2M 14790 | h | O₂ Flowmeter 0...15 L/min with rail clamp | M 24669 |
| 10 | Rail (2 ea can be fitted) | 2M 13390 | i | Humidifier nebulizer | M 19995 |
| 11 | Cabinet, small Cabinet, large | 2M 14780 2M 12910 | k | Adapter | 2M 16843 |
| 12 | Infusion stand | 2M 12940 | l | Socket (set of 4 ea) Dust cover | M 22199 2M 6835 |
| 13 | Restrainer | 2M 13190 | | Spare parts | |
| 14 | Terminal box | 2M 13300 | | Mattress | 2M 12925 |
| 15 | Swivel table | 2M 13380 | | Head pad | 2M 12955 |
| 16 | Small operating lamp | 2M 13400 | | Skin-temperature sensor | 2M 18553 |
| 18 | Aspirator, pediatrics (ejector operation) | M 27174 | | Screwdriver 0.5 x 2.3 mm | 13 22 923 |
| 19 | O₂ distributor with cylinder connection (R ^{3/4"} connecting thread) O₂ cylinder, 3 litres (R ^{3/4"} , charged, 200 bar) | 2M 11150 B 2533 | | Retaining strap | 2M 13898 |
| 20 | O₂ distributor (supplied from central system) | 2M 11200 | | Bulb 220 V for operating lamp | 2M 14413 |
| 21 | Pillow to facilitate intubation | 2M 15298 | | | |
| 22 | Timer with holder | 2M 15333 | | | |

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DRÄGERWERK AG LÜBECK