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These instructions for use apply only to Incubator 8000 IC with Serial No.:

Without entry of the Serial No. by Dräger these instructions for use are provided for general information only and are not destined to be used with a specific device!

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
M 110 03

90 27 494 - GA 6141.222 d/e  
© Drägerwerk AG  
3rd edition - January 1992  
Subject to alteration

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ED 212891/02 / Printed in the FRG / Serien-Nr. AREA - F001 - Printed on chlorine free paper

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## For Your Safety and That of Your Patients

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For correct and effective use of the apparatus and to avoid hazards it is essential to read the following recommendations and to act accordingly<sup>1)</sup>:

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### Strictly follow the instructions for use Intended Use

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For the intensive care of premature babies (including those below 1500 g in weight) and sick newborn babies, with controlled temperature (optionally alternatively skin temperature), humidity and oxygen concentration.

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Skin temperature control must not be used on babies in shock (page 28)

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Correct physiological concentrations of oxygen must be maintained.  
(Hypoxaemia / Hyperoxaemia, see page 32)

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The incubator should only be used by experienced nursing staff and under the supervision of qualified medical staff, who are familiar with the presently known risks and advantages of incubator application!

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## 3) Module faults

Warning for faults in the control modules for

- skin temperature and
- air humidity

produce the following warning signals:

Modul operation status	Actual value display	Alarm tone	"Inop"-LED
Off	"Err" continuous	mute	off
On	"Err" flashing	continuous warning tone non-suppressible	lights continuously

The operator may continue to use the incubator after the faulty module has been switched off.

The functions of the other modules which are operating correctly are retained.

Call DrägerService.

## 4) Total shut-down

Faults in the central microprocessor control system and in the air temperature control module result in the complete shut-down of the incubator.

Warning signals:

- o Red "Inop" LED lights.
- o Continuous warning tone sounds (non-suppressible).

Individual checks should be carried out to determine whether the fault is temporary.

If the incubator shows the "Inop" signal after it is switched on again, it is inoperative.

The incubator must not be used and must be clearly marked faulty.

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


## Before Using for the First Time

- Remove all packaging material.
- Remove transport safety lock from the fan (polystyrene block). (See "Stripping Down", page 37).

- Fit the canopy at a slight angle.

- Fit the front door in such a way that the lower section rests on the hinge, almost vertically, then push this lower section down firmly until the door engages in the hinge.

- Closing front door: Press the latch together with thumb and indexfinger, thereby pressing front door against the canopy. After engaging the front door the red index shall no longer be visible.


## Fitting ventilator hose support

- Open front door.
- With door open raise the bed.
- Insert the cone on the support arm into the holes on the right or left on the far side of the bed and secure from below with knurled screw.

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

**Preparation**

**Fitting mounting plate for ancillary equipment**

- May be attached on either side.
- Maximum load 25 kg.

**Fitting ancillary equipment to plate:**

- Remove both bottom strips.
- Tilt equipment forwards by about 45°.
- Insert front latches into the slots in the mounting plate.
- Lower the equipment at the rear. Insert the rear latches into the slots in the mounting plate. Lock with knurled screws.

**Fitting drawer unit**

Units may be fitted to the left and/or right of the pedestal.

- Lift drawers at the front and remove.
- Place unit onto the trolley and secure according to the instructions. (Assembly instructions in unit).
- Replace drawers.
- Fill the drawers with whatever is required.

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**Appendix Description**

- The actual value display of the relevant control flashes.
- The specific warning LED flashes.
- The general warning LED flashes.
- Intermittent warning tone sounds (can be muted for 10 minutes)

The warning tone for the difference between set and actual values of air temperature of more than 1.5 °C is suppressed automatically for 30 minutes when the incubator is switched on to avoid a continuous warning signal sounding during the warming-up phase.

If a warning has been suppressed, the red alarm LED, which has previously flashed, stays alight but remains constant, and a yellow LED indicator above the suppression button indicates the suppressed mode.

The «Reset» warning button for over-temperature has to be re-set by the operator so that he/she is fully aware of the fault.

All other alarms are re-set automatically when the particular conditions have returned to normal.

If another warning occurs when warnings are suppressed, the audible warning is re-activated immediately. Depending on the type of fault, this later warning can also be suppressed again. The period for automatic re-activation of the first warning tone is then prolonged by the time interval between the two warnings.

**2) Sensor warnings, fan failure warnings**

Faults in the sensors for

- oxygen,
- skin temperature,
- air humidity
- air temperature

as well as fan failure produce the following alarm signals:

- o Three dashes (" - - - ") flash on the relevant actual value display.
- o The actual air temperature value flashes in the event of fan failure.
- o The specific warning LED flashes.
- o The general warning LED flashes.
- o Continuous warning tone (non-suppressible).

When there have been sensor warnings for O<sub>2</sub> oxygen concentration and relative humidity and the appropriate module has been switched off, the incubator can continue to operate on the other functions.

This is also the situation when there is a skin temperature control warning, though in this case the faulty sensor must be removed.

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The function of actuators, acknowledgement signals and displays is checked by switching them on and off.  
This test is also repeated every 10 minutes during operation.

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All operations activated and present in the unit are tested. An "Error" message may be indicated for a faulty module even if it is switched off.

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When the incubator is switched on again, or when a module is switched on, the relevant required value or set value is set to "starting value".

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The set value display then flashes requesting confirmation or modification of this value. The actual value indication is displayed alternately with the prompt "SET".

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**Warnings**

The incubator has a hierarchical system of warning signals. Faults which occur are signalled in order of importance. If a non-essential option fails, the functions which are of vital importance remain in operation.

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**Continuous warning tone:**

Cannot be silenced for faults which entail the greatest potential danger.

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**Intermittent warning tone:**

Can be silenced for up to 10 minutes for set value deviation, water shortage and over-temperature.

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In addition, the relevant warning LEDs flash to enable the operator to locate the fault rapidly. Each individual warning also activates an indicator on the central warning LED so that a visual signal can still be seen even when the flap which shows the brief instructions for use is closed.

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**Warning Signals**

**1) Warnings for set value deviations and water shortage:**


- Deviations outside the set limits of more than
- ± 5 % O<sub>2</sub> by vol. oxygen concentration
- ± 0.5 °C skin temperature
- ± 1.5 °C air temperature,

in addition to overtemperature (over 38° or 40°C) and water shortage result in the following warning signals:

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
**Accessories for rail mounting**

- Fit rail clamp for accessories to side handles and tighten knurled screw.


To attach:  
O<sub>2</sub> flowmeter  
Infusion stand  
Bronchial aspirator and  
Instrument tray

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



**Oxygen**

The oxygen concentration required inside the canopy can be pre-selected and controlled automatically by the oxygen control module.

If the actual value deviates from the set value, warnings are given.

The control system will provide an oxygen concentration of 40% in 1.5 minutes or of 60% in less than 7 minutes.

Rapid action is, therefore, possible in emergencies and similarly the selected oxygen concentration is rapidly re-tained after the front door or the hand ports have been opened.

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**Air Humidity**

The humidity control gives the air humidity required.

Sterile water bottles supply the water for evaporation, thus guaranteeing hygienic air humidification.

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**Hygiene**

All parts of the incubator which come into contact with the gas breathed by the baby can be removed from the basic unit.

The incubator 8000 IC can be thermally disinfected in the Dräger Aseptor 9000.

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**Control**

The incubator is controlled by a microprocessor and carries out a general self-test whenever it is switched on and then automatically every 10 minutes thereafter. This test checks the various elements which are particularly important for safety.

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The operator is automatically requested to re-calibrate the oxygen sensors after switching on the O<sub>2</sub> control system and, thereafter, after every 24 hours of operation.

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Measured value deviations or system failures are indicated visually and audibly.

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The heating system and valves are switched off automatically if operating conditions are outside permitted limits.

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An additional air cooler produces rapid cooling as soon as the actual value of the air temperature is a few tenths of a degree above the set value.

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**Safety Features**

After switching on, the incubator carries out a self-test. This involves a test of all memories in the microprocessor control system and a check to establish whether the various programme segments are running correctly.

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
## Appendix Description

### Design

Heated air is channelled between the double walls of the entire front face (front door), over the roof of the canopy and then drawn down past the rear panel by means of an extraction system.

Hence the baby lies in a still area with very low air flow, and so convective heat loss is prevented.

When the hinged front door and the hand port doors are open, there is still an effective hot air curtain which reduces cooling inside the canopy.

The mattress is made of a soft foam plastic, hygienically encased in film, which provides a high level of comfort and very low heat loss.

### Ergonomics

The working height of the incubator may be adjusted to give the operator optimum working conditions, e.g. prolonged operations can be carried out in a sitting position.

The incubator has rotating castors so that it can be moved around easily though two of the castors may be locked in position when necessary.

Ancillary equipment can be attached to the side handles.

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

## Preparing for Operation

Equipment already cleaned, disinfected and fitted.

### Checking air filter, filling water supply

1 Open the rear panel and lower.

2 Open filter cover plate downwards.

When filter is already fitted:

- Check date of filter.
- If the filter has been in use for more than two months, replace with a new filter.
- Write the fitting date on the filter label.

**When fitting: the arrow on the edge of the filter must point towards the incubator.**

• Close door.

• Remove all empty water bottles (three), from the socket and replace with full bottles.


• Use only infusion bottles which conform to IL 500 - DIN 58363 - glass 1. These may be refilled with 500 mL clear, distilled or demineralised water from the pharmacy. Additives must not be used.

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

Fitting environmental sensor

This measures temperature, relative humidity and O<sub>2</sub> concentration.

- Fit two O<sub>2</sub> sensors with similar expiry dates:
- Unscrew the retaining ring. Fit O<sub>2</sub> sensor. Replace retaining ring.

- 1 Fit the environmental sensor onto the spindle and swing it into place.
- 2 Open the rear panel and lower.
- 3 Slide the sensor connector from below through the slot between incubator housing and rear panel. Fit the probe into the socket and secure.

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Sensor : Only for earthed skin temperature sensors produced by Drägerwerk AG (see Order List page 57)

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	Display	Skin temperature sensor
Calibration error limits:	± 0.2 °C	± 0.1 °C
Calibration error limits of the entire measuring system over the measuring range: ± 0.3 °C		

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**Gewicht** 89 kg

**Weight** : 89 kg

**Additional informations for use of incubators specified for 100 V**

Heating power (air): 303 W

Heating power (evaporation of water): 76 W

Warming-up time: approx. 45 minutes from 20°C to 31°C with or without air humidification


Humidification: If the relative humidity of the ambient atmosphere is very low, the maximum relative humidity of the incubator air may not exceed 70%.

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
- Route hose (or cable) through the flexible hose seals.

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
- Guide hose (or cable) around the support hook on the inside surface.
- Snap the ventilator hoses and cable into the clip on the end of the ventilator hose support arm.

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

**Testing Readiness for Operation**

**Electrical supply**

- The mains voltage must correspond to the mains voltage specified on the rating plate of the incubator.
- Lower the flap to see the brief instructions for use.

**O<sub>2</sub> connection**

**With O<sub>2</sub> control**

- 1 Use a connecting hose with an M 12 x 1 threaded connector and an oxygen probe. The probe should be in the «parking» position in the terminal unit.

**Without O<sub>2</sub> control**

- 2 Use a flowmeter with an oxygen probe. The probe should be in the «parking» position in the terminal unit.

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**Technical Data**

Tilt angle of bed: Continuously adjustable right and left up to a 10° tilt angle.

**Display and adjustment ranges:**

	Adjustment range	Display range
Air temperature:	+28 to +39°C	+10 to +50°C
Skin temperature:	+35 to +37°C	+33 to +38°C
Air humidity:	35 to 85% rel. humidity	0 to 99% rel. h.
Oxygen concentration:	21 to 75% by vol. at 50 Hz 21 to 65% by vol. at 60 Hz	0 to 99% by vol.

Measuring principles of sensors:  
 Air temperature : NTC (2 x)  
 Air humidity : capacitative  
 Oxygen content : galvanic cell (2 x)  
 Skin temperature: NTC

Noise level (inside the canopy) : < 49 dB (A) at 50 Hz operation < 56 dB (A) at 60 Hz operation

Fresh air filter : Particle class P2 in acc. with DIN 3181 factor 2 %

Classification : Class I Equipment Type B  
 The incubator complies with IEC 601, part 1 and 2-19  
 The skin temperature sensor is not defibrillation resistant.

RFI suppression : in accordance with relevant regulation 1046 of the German Post Office

Identification : on the rating plate: Serial No. Part No. Registration No. The equipment type is defined by the order number. (see also Order List page 59)

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## Technical Data

### Ambient conditions

Ambient temperature : +20 °C to +30 °C (operation)  
0 °C to +70 °C (storage)

Ambient humidity : 15 % to 95 % relative humidity

Ambient air pressure : 900 hPa to 1100 hPa

### Performance Data

Operating voltages: 100/120/127/230/240 VAC, to order; 50/60 Hz

Current consumption : see rating plate

Heating (air) : 400 W

Heating (boiler) : 100 W

### Operating Data

Warming-up time : 35 minutes from 20 °C to 31 °C with or without air humidification \*

Temperature drop when opening two hand ports : < 1 °C

Opening front door: < 5 °C in each case at 25 °C room temperature and 36 °C incubator inside temperature

Oxygen increase : < 7 min. from 21 to 60 % O<sub>2</sub> by vol.

Humidification : Evaporation of distilled water by heating

Air velocity : about 8 cm/sec over the bed

Fresh air supply rate : 30 L/min at 50 Hz operation  
37 L/min at 60 Hz operation

CO<sub>2</sub>-Flushing : (accordingly, IEC 601-2-19/105.1): The maximum CO<sub>2</sub>-Concentration present inside the baby compartment is 0,2 %.

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

## Testing Readiness for Operation

### Switching on

- Lower the flap to see the brief instructions for use.

- 1 Switch ON
- 2 Green LED is lit.

Whilst the microprocessor and its peripherals are self-tested, dashes will be shown on the actual value displays of the modules for about 35 seconds.

### Function test


- 3 Press button - all displays except the power failure LED will light up briefly, a warning tone will sound briefly, and all digital displays will indicate "88.8.". All displays will then go blank until the previous settings are shown.

This function test for all lamps, display segments and the audible warning signal can be carried out during use.

**The function test must be carried out at least once a day.**

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
# Testing Readiness for Operation

- 1 Central warning:  
The red warning LED flashes at the same time as individual control warnings are given.
- 2 "inop" warning:  
Red warning LED flashes: fault in the central control system of the incubator, or in a control module which will also display "Err".  
All functions of the incubator then shut down.  
As soon as the faulty module has been switched off, all other functions will again be operational.  
Call DrägerService.

# Suppressing Warnings

- 3 Press button.  
The acoustic warning is suppressed for 10 minutes (does not apply to "inop", mains failure warnings, fan failure and sensor failure).
- 4 The yellow pilot LED will light to show that warnings are suppressed

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**Front panel**

- 1 Humidity control
- 2 Oxygen concentration control
- 3 Skin temperature control (Option)
- 4 Air temperature control
- 5 On/Off switch and self-test

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

Front view

- 1 Environmental sensor housing
- 2 Front door latch
- 3 Hand port door latch
- 4 Bed
- 5 Side handle
- 6 Tilt adjustment for bed
- 7 Castors
- 8 Flexible hose seals
- 9 Front control panel
- 10 Brief instructions for use
- 11 Drawer unit

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
Rear view

- 1 Rear panel latch
- 2 O<sub>2</sub> terminal unit
- 3 O<sub>2</sub> threaded connection
- 4 Foot pedal for height adjustment
- 5 Air filter
- 6 Water supply
- 7 Plug for environmental sensor

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
Power failure warning and battery testing

- Do not connect the incubator to the mains power supply.

- 1 Switch ON.  
An intermittent tone sounds and maintains a constant volume for at least 30 seconds.

- 2 Red LED for mains failure warning is lit.

- 1 Switch OFF.  
Mains failure warning is cancelled.

If the volume of the warning tone drops prematurely, the internal battery must be charged for at least 24 hours:

- Connect incubator to mains supply.

- Switch ON.  
Leave switched on at the mains for 24 hours to recharge.  
The battery is recharged.

- Disconnect from the mains power supply:  
The volume of the intermittent tone must remain constant for at least 30 seconds.


- Reconnection to mains or switch OFF = intermittent tone should cease.

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Adjusting working height

Keep the space above the incubator clear of obstructions.

Clear all objects from the top of the drawer unit.

Hoses and cables must be long enough to allow for adjusting the height of the incubator so that kinking, forced disconnection or squashing cannot occur.

- 1 Press foot pedal. The incubator can be raised.
- 2 Press foot pedal. The incubator can be lowered.

• Adjust incubator to comfortable working height.

Humidity Control Module			
Fault Symbol	Description	Cause	Remedy
"H <sub>2</sub> O"	Red LED and intermittent tone.	Water shortage. Fault in heater	Replace empty water bottles with full bottles Switch off control module. Call DrägerService!
"Sensor"	Red LED and intermittent tone.	Faulty sensor. Sensor plug not properly inserted.	Switch off control module. Call DrägerService! Check connections and correct, if necessary.
Display "Err"	Continuous tone.	Fault in control module.	Switch off control module. All other functions remain operational. Call DrägerService!

Skin Temperature Control Module			
Fault Symbol	Description	Cause	Remedy
"± 0.5°C"	Red LED and intermittent tone.	Temperature deviates set value by more than ± 0.5°C.	Below set value: Check that the sensor is properly fitted to the baby. Above set value: measure the core temperature of the baby rectally. Inform the doctor responsible for the baby immediately.
"Sensor"	Red LED and intermittent tone.	Sensor plug not properly inserted Faulty sensor.	Check connection and correct, if necessary Replace sensor.
Display "Err"	Continuous tone.	Fault in control module.	Switch off control module. All other functions remain operational. Call DrägerService!
"Check" button display <b>NOT</b> 36°C ± 0,1°C	Temperature reading outside tolerance levels.	Temperature sensor not operating correctly.	Call DrägerService! Have equipment re-calibrated!

Oxygen Control Module			
Fault Symbol	Description	Cause	Remedy
"± 5 % by vol."	Red LED and intermittent tone.	Oxygen concentration deviates from set value more than ±5 % by vol.	Below set value: Close openings on canopy. Check oxygen connecting hose to piped medical gas system. Above set value: Switch off other O <sub>2</sub> sources. If O <sub>2</sub> concentration does not drop, disconnect from piped medical gas system and continue operating with manual control (see page 33). Call DrägerService!
"Sensor"	Red LED and intermittent tone.	Sensor not connected. Sensor swung out. Plug not properly connected. Faulty sensor.	Connect Sensor. Swing sensor back into place. Check connection and correct, if necessary. Remove faulty sensor and replace with a new one. Calibrate the new sensor.
"Sensor"	Red LED and intermittent tone. "CAL" flashes.	Calibration not started.	Calibrate sensor (see page 32).
"Cal" Display "CAL"	Yellow LED and intermittent tone. Display flashing.	Requesting re-calibration of O <sub>2</sub> sensors. During continuous operation every 24 hours and during start-up.	Calibrate sensors (see page 32).
Display "Err"	Continuous tone.	Fault in control module.	Switch off control module. All other functions remain operational. Call DrägerService!
"Inop" O <sub>2</sub> display	Red LED and continuous tone. Flashing O <sub>2</sub> actual value	O <sub>2</sub> increase after O <sub>2</sub> control module is switched off.	Disconnect both from piped medical gas system and other O <sub>2</sub> sources. Switch incubator off and then on again. If fault persists, call DrägerService

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## Testing Readiness for Operation

### Visual inspection

- No gaps, cracks or holes in the incubator canopy.

The hinges and pivots of the canopy must be in good working order.


Opening the canopy:

- Press semicircular button.  
Check that hand ports stay properly engaged and securely locked.

If not, call DrägerService.

- The incubator must be cleared and assembled in accordance with these Instructions for Use.  
See also «Care» on page 37.

Check that all the ancillary equipment fitted is ready for operation, in accordance with the relevant instructions for Use.

The unit is ready for operation when all tests have been successfully completed.

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**Precautions**

**Warming-up time**

Allow adequate time for warming-up before placing the baby in the incubator (about 35 minutes).

**Avoiding overheating**

**Incubator temperature control**

- Desired temperature increases occur very rapidly (high heating power).
- Temperature drops slowly due to good thermal insulation.

Additional external heat sources such as sunshine, heat lamps, spotlamps, electric cushions etc. must be avoided.

These increase the air temperature inside the incubator in an uncontrolled manner.

**Incubator air temperature setting**

The baby experiences minimal

- heat loss by radiation as a result of double wall insulation
- heat loss by convection due to low air speed
- heat loss by conduction through the mattress
- heat loss by evaporation because of the high humidity setting in the incubator.

Therefore:

in comparison with other incubators such as the Incubator 6000 or 7000, a lower incubator air temperature can be set.

The core temperature of the baby must be monitored continually - particularly during the first few hours of incubator care.

Mains Switch - Front Panel			
Fault Symbol	Description	Cause	Remedy
Warning	General warning: red LED and intermittent tone.	A module has signalled a fault.	Lower flap to see brief instructions for use and locate fault in appropriate module. Refer to "symbol" column for further details.
"Inop"	Red LED and continuous tone.	Central fault in the electronic control system.	Switch incubator off and then on again. Renewed "Inop" message shows incubator not operational. Call DrägerService! If no "Inop" message recurs, set required values, press "reset" button.
"Inop" Display "Err"	Red LED and continuous tone. "Err" on actual value display of a control module.	Fault in the relevant control module.	Switch off the relevant control module. All other functions remain operational. Call DrägerService!
Power symbol 	Red LED and continuous tone.	Power failure.	Connect to mains power.  Inform hospital technician. Call DrägerService!

Height Adjustment			
Fault Symbol	Description	Cause	Remedy
Arrow symbols  	No adjustment when the switch is operated	Thermal cutout has been activated.  The unit is already in final position.	Wait until the motor has cooled down.  No further height adjustment is possible.

Air Temperature Control Module			
Fault Symbol	Description	Cause	Remedy
"± 1.5°C"	Red LED and intermittent tone.	Temperature deviates from set value by more than ± 1.5°C.	Below set value: Close opening in canopy. Wait for incubator to warm up. Above set value: Avoid additional heat sources such as radiators, sunshine, lamps. Wait until incubator has cooled down
Thermometer symbol 	Red LED and intermittent tone.	Temperature above 38 °C (or above 40 °C *)	Avoid additional heat sources such as radiators, sunshine, lamps etc. Wait until the incubator has cooled down. <b>Press "reset" button when the temperature has dropped below the warning temperature limit.</b>
Fan symbol 	Red LED and continuous tone.	Fan blade not fitted.  Fan blade not rotating or rotating too slowly.	Fit fan blade.  Check that fan blade and drive shaft are firmly coupled.
"Sensor"	Red LED and continuous tone.	Environmental sensor disconnected. Sensor not swung back into place.	Push the environmental sensor firmly into its socket.  Swing back under canopy.
		Defective temperature sensor	Call DrägerService!
		Incubator temperature below 5°C.	Wait for incubator to warm up fully

\*) For set values up to 37°C switch-off occurs at 38°C; for set values of 37.1°C to 39°C switch-off occurs at 40°C.

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**Reducing the internal temperature of the incubator**

The cooling time is dependent on the design and can be accelerated by:

- reducing the outside temperature,
- reducing the air-humidity setting,

The rate of cooling is **not** accelerated by reducing the air temperature setting below the actual value required.

In urgent cases; open front door or hand ports.

Then it is absolutely essential that there is continuous supervision so that the baby cannot fall out of the incubator whilst the door is open.

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**Phototherapy in the incubator**

Absorption of light through the baby's skin will supply heat which may increase the baby's core temperature.

Therefore:

- Reduce setting of incubator air temperature by about 2°C 15 minutes before phototherapy.
- Reduce incubator air-humidity setting.
- Reduce the room temperature to at least 3°C below the incubator air temperature.

(This value applies to the 800/8000 type Phototherapy Unit, see pages 10 and 36).

(Other phototherapy units, particularly those without a built-in fan, may cause even greater heating of the incubator; the room temperature must then be reduced still further).

**The core temperature of the baby must be monitored particularly carefully during phototherapy.**

The supply of fluids to the baby, e.g. by parenteral infusion, must be increased in order to compensate for increased water loss during phototherapy.

The phototherapy lamp and incubator canopy must not be covered with cloths, aluminium foil, or other materials, to boost the photo-therapeutic effect.

**Risk of heat build-up.**

The incubator cannot be adequately cooled by the ambient air; therefore there is a danger of overheating the patient.

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**Monitoring the breathing gas temperature**

The temperature of the breathing gas must be monitored on ventilated or CPAP patients because:

breathing gas hoses may be warmed by the heated air circulating in the incubator.

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**Preventing increased noise**

Unacceptably high noise levels may be caused by:

- using head boxes to deliver pressurised gas.
- wear on the bearings of the air circulation motor.

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**Electrical safety of additional equipment**

Particular care must be taken to ensure that additional equipment connected to the baby is electrically safe.

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**Fire risks from oxygen**

- No naked lights and burning matches.

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- Textiles, oils and plastics can very easily catch fire and burn rapidly in an oxygen-enriched atmosphere.
- All oxygen fittings and seals must be kept free of oil and grease.
- Open valves slowly.
- Do not use the incubator where there are flammable anaesthetic gases or disinfection agents. Risk of explosion.
- Flammable liquids, such as alcohol, ether and acetone, must not be used or stored in the incubator.
- Do not use any electrical equipment inside the incubator.

Except for: equipment expressly designed for use in rooms where there is a risk of explosion.

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**Physiological risks from oxygen**

The air in the incubator should only be enriched with oxygen when prescribed by a doctor.

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It is absolutely essential that such oxygen enrichment is controlled on the basis of the arterially-measured oxygen partial pressure in the blood of the patient.

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Doing this is the only way of avoiding both hyperoxemia (damage to the eyes) and hypoxemia (damage to the brain).

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

Operation

- Connect to power supply.
  - Connect to gas supply, if necessary.  
If standby oxygen is required, put the probe of the oxygen supply line in the «parking» position at the terminal unit.
- Lock the castors in position, if necessary.
  - Lower the flap so that the brief instructions for use can be seen.
  - Open the front door.
- Pull the bed right out, if necessary.
  - Place a suitable nappy/diaper on the bed and slide it back in.
  - Close the front door.

Switching on

- Switch ON.
  - Green LED lights.
- The functioning of the microprocessor and its peripherals is self-tested. Dashes are displayed on the actual value display of the module for about 35 seconds.
- Press button.  
Function test of all lamps, display segments and warning tone.
- All displays except the mains failure LED are briefly lit, a warning tone sounds and all digital displays show «88.8». All displays will then go blank until the measured values and settings are shown.
- Carry out this test at least once a day.
- This test may be carried out during operation. It does not change or cancel any settings.

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- Adjust incubator to comfortable working height.

## Standby-operation

To make use of the device without delay, the following settings are suggested:

Air temperature: 33°C; see page 26

Relative humidity: 60%; see page 34

Oxygen: Oxygen control module switched off; see page 32.

## Placing the baby in the incubator

After a warming-up time of 35 minutes.

- Put the baby in the incubator.
- Slide the bed back in position.
- Close the front door and check that the latches are **properly engaged.**  
After engaging the front door the red index shall no longer be visible.

## Positioning the baby

- Tilt the bed using the handwheels (left or right) so that the baby's head is either raised or lowered as necessary.

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Clean and disinfect unit and parts before any maintenance - also before sending for repair.

**O<sub>2</sub> sensor capsule:** In the event of a malfunction, always replace both capsules. Ensure that expiry dates on both are similar.  
("Fault - Cause - Remedy" page 45)

**Air filter:** Replace after two months. (See "Preparing for Operation", page 11)

**Hand port seals:** The seals must be replaced when the material becomes brittle or tacky, or when a proper seal can no longer be obtained.

Maintenance works to be carried out by trained servicemen only:

Inspection of incubator Every six months

The fan motor shall be lubricated with 10 drops of oil 51524-HLP 32 Every six months

NiCd battery for mains failure warning After 1 year

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
**Sterilisation**

up to a maximum of 120°C:

- water hoses
- fan
- hand port seals
- hose seals

up to a maximum of 134°C:

- Water bottles and bottle holder

**Before re-use**

- Re-assemble the incubator, as for »Stripping down«, page 37.

in addition:

- Mounting the seals on hand port:
  - seal marked with an L on nose is mounted in the left hand port.
  - seal marked with an R on nose is mounted in the right hand port.
  - nose shall rest on hinge.
  - sealing lip (this is the thin part of the seal) should show outwards.

- open rear panel.
- Open filter cover plate downwards
- re-insert air filter, or replace with new one.

- When incubator has been wipe disinfected, switch on the assembled incubator and operate for a few hours unoccupied to eliminate all traces of disinfectant.

- Running the water heater dry results in:
  - disinfection of heater
  - periodic display of »water shortage« and function check of warning at the same time.

- Immediately before baby is placed in the incubator:
  - re-fill with distilled water
  - check for readiness of operation, page 14.

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
In the event of a warning, do not alter the settings, since this would alter the control values; instead suppress warning tone for 10 minutes:

- Press warning suppression button.

It should be possible to correct the fault in this time.

**Setting controls and limit values**

The controls for air temperature, skin temperature and oxygen concentration are given limit values automatically:

- For air temperature control:
  - Upper limit value:  
Value more than 1.5°C above actual setting
  - Lower limit value:  
Value more than 1.5°C below actual setting

- For skin temperature control:
  - Upper limit value:  
Value more than 0.5°C above actual setting
  - Lower limit value:  
Value more than 0.5°C below actual setting

- For oxygen concentration control:
  - Upper limit value:  
Value more than 5% by vol. above actual setting
  - Lower limit value:  
Value more than 5% by vol. below actual setting

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
## Controlling Air Temperature

Range 28°C to 37°C

Check the baby's core temperature at regular intervals.

Correct the air temperature setting, if necessary.

The »air temperature control system« is automatically activated when the incubator is switched on.

If the incubator is set in the »skin temperature control« mode, it can be switched over to the »air temperature control« mode (Option. Only skin-temperature-control where installed).

- 1 Press »Control« button. (Option. Only skin-temperature-control where installed.)
- 2 Green LED lights. (Option. Only skin-temperature-control where installed.)

It is still possible to measure and display the skin temperature. (Option. Only skin-temperature-control where installed.)

In the same way it is possible to switch from air temperature control to skin temperature control.

- 3 The display of the air temperature inside the incubator alternately flashes with the display »SEt« (request to confirm or reset the value).

- 4 Pre-set value flashes:

- A pre-set value of 33°C is shown when switching on the incubator.
  - After switching from skin temperature control to air temperature control the actual air temperature is shown as the set air temperature value.
- Option. Only skin-temperature-control where installed.

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## Cleaning

Use only the cleaning materials and disinfectants specified

Otherwise there is the risk of damaging the plexiglass and macrolon as, for example, when using alcohol.

Wipe all surfaces with a damp cloth soaked in detergent. Lower the inner panels, wipe the inside carefully.

Use a disinfectant for cleaning if there is any risk of infection.

## Wipe disinfection

Wipe over all cleaned surfaces of the canopy, dismantled parts, and the unit itself using a suitable disinfectant. We recommend the following disinfectants. Please note the maximum concentrations specified:

Bacillotox®	required concentration up to 2 %
Buraton 10F®	as specified by the manufacturer
Dismozon® pur	as specified by the manufacturer
Frekanol®	as specified by the manufacturer
Gevisol®	required concentration up to 5 %

Allow time for the disinfection process, then wipe with a clean, damp cloth and dry well.


The sensor housing must not be treated with a spray disinfectant or disinfected in a bath.

Do not use UV-radiation in incubator as this may damage the plexiglass.

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

Additional cleaning, after two months use:

- Swing the sensor housing out from the canopy and remove it.

- Tilt the incubator canopy slightly backwards, though not quite as far as it will go. Grip from the rear. Lift canopy off in the direction of the arrow.

- Strip down, as described on page 37, in addition:
- Remove the trough.

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
- 1 Press the button briefly to confirm the set value required.
- 2 The incubator acknowledges: flashing "SET" ceases; only actual temperature is shown. The set value stops flashing.


**Setting air temperature value**

- Increasing or reducing the value (max. 37°C min. 28°C):
- 1 Press button on left = value decreases; press button on right = value increases; keep button pressed until the desired value is shown.
- 3 The green LED lights: the set value has not yet been reached. The green LED is not lit: set value has been reached, or exceeded. The green LED flashes: set temperature has been reached. The incubator heating responds to the adjustment of the control system.

**Increasing the temperature range from 37 °C to 39 °C**

- Only if prescribed by a doctor. The baby's core temperature must be monitored very carefully.
- 4 Press button.
- 5 Yellow LED is lit. The warning limit for over-temperature is increased from 38 °C to 40 °C.

Set value control as for "Controlling air temperature", see page 27.

A subsequent reduction below 37°C automatically cancels the increase in the range of set values. When the incubator temperature is still above 38°C, the over-temperature alarm is triggered.

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
Recommendation: first adjust set value to 37.1 °C and allow incubator to cool down.

As soon as the incubator temperature drops below 38,0°C, set value below 37°C, when applicable.

6 Warning LED's and reset button, see pages 45; 49 and 56 for explanation of symbols.

### Controlling skin temperature

#### (Option)

Not to be used on babies in shock when skin temperature lies well below the core body temperature.

(Using a skin temperature control could lead to overheating.)

In these circumstances operate incubator with air temperature only.

The core temperature must be monitored carefully.

### Preparation for skin temperature control

- Establish an appropriate temperature inside the incubator using the "air temperature control" mode.

- Plug the skin temperature sensor into the socket in the environmental sensor housing.

When switched to operating mode "air temperature control", the skin temperature is merely indicated but not controlled.

- Fit a double-sided adhesive ring to the skin temperature sensor:

- Remove the backing from the adhesive ring.

- Stick the adhesive ring around the centre of the metal part of the sensor.

- Remove the protective foil from the adhesive ring.

- Attach the sensor to the skin where the temperature is to be measured with the adhesive ring. Keep sensor cable in place with adhesive plaster.

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
- after each change of baby
- at least once a week

Clean and disinfect the incubator thoroughly.

Clean and disinfect the accessories, such as the bronchial aspirator:

See the relevant instructions for Use.

### Stripping Down

- Switch off incubator. Disconnect from mains and disconnect the oxygen supply.

- Remove all accessories.

- Open front door.

- Tilt the canopy backwards as far as it will go.

1 Swing out the sensor housing.

2 Disconnect skin temperature sensor.

- Remove flexible hose seals.

- Remove seals from hand ports.

3 Remove the mattress from the bed.

4 Lift off the bed. Undo the screw on the ventilator hose support. Remove the hose support.

5 Pull the metal brackets supporting the bed up and out.

6 Remove the base.

Caution. Be careful not to burn yourself.

This risk is particularly great when the incubator is closed and even exists long after the unit has been switched off. It takes about one hour for the temperature to fall to 70°C.

7 Remove fan.

- Open the rear panel.

- Remove water bottles (see page 11 or 35)

- Disconnect hose. Unscrew the bottle holder.

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
## Using vacuum mattress (Accessory)

The contours of the mattress can be formed into any shape, and will keep this shape after the mattress has been evacuated.

In this way babies can be supported in extreme positions for special treatment.

The standard mattress may also remain inside the incubator.

- Open the front door and lower it.
- Position the mattress inside and form it into the shape required.
- Position the baby on the mattress and «mould» the mattress around the baby.
- Evacuate the vacuum mattress with suction unit.
- Close valve and remove hose.
- Close front door.

## Using phototherapy unit 800/8000 (Accessory)

- Push the trolley carrying the phototherapy unit under the incubator.
- Adjust the height of the incubator: The underside of the phototherapy unit must be at least 0.4 m away from the patient.
- Operate the phototherapy unit in accordance with the relevant instructions for Use.

If the working height or the height of the bed is changed later, the clearance of 0.4 m must be maintained.

Note the safety advice on page 20-21 about avoiding overheating.

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If the baby is lying on its back:  
attach sensor to the abdomen between the navel and the sternum.

If the baby is lying on its stomach:  
Attach the sensor to the back in the region of the kidneys.

The skin temperature sensor must never be under the baby since it would then measure core temperature.

Result:  
Incubator temperature too low, and rapid temperature fluctuations on the skin cannot be controlled.

Check regularly that the sensor is properly fixed to the baby's skin.

Loose fitting may lead to faulty measurement and result in the baby becoming overheated.

(The incubator air temperature does not, however, exceed 37°C).

Do not use the skin temperature sensor as a rectal temperature sensor!

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



### Controlling skin temperature

Allow the skin temperature sensor at least 5 minutes to reach the baby's temperature:

When displayed value remains constant,

1 Press button = ON

The air temperature control is replaced by the skin temperature control. The air temperature is still displayed.

(Press button again = switch back to air temperature control).

2 Green LED lights.

3 The skin temperature is shown alternately with the display "SET": request to confirm or adjust the set value.

4 Set value flashes initially. It then shows the actual value of the skin temperature control system.

(Max. 37°C, min. 35°C)

Otherwise set value as for "Setting air temperature value" as described on page 27.

5 Press button. Display must show 36°C ± 0.1°C.

6 Warning LED's, see pages 45; 49 and 56 for details.

If the skin temperature sensor is going to be temporarily removed, the air temperature control must be switched on first.

Display for a measured value outside the measuring range: below 33 and above 38°C.

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

### Refilling the water supply

● Open the rear panel and lower it. Wait a few seconds for the water to flow back from the vaporiser unit.

● Pull the empty water bottles out of the socket and fit full water bottles. Check that the water bottles are properly secured.

● Only use infusion bottle IL 500 - DIN 58363 - glass 1. Refill them with a maximum of 500 mL clear, distilled and de-mineralised water from the hospital pharmacy.

**Do not use any additives.**

● Lift the rear panel and close it.

### Using baby scales (Accessory)

1 Adjust the bed to maximum height with both handwheels.

● Open and then lower the front door.

2 Slide the electronic scales under the bed. Guide the scales along the ridges on the baseplate.

● Close the front door.

● Lower the bed using both handwheels until it comes to rest on the scales.

● Weighing: in accordance with the Instruction for Use for the scales.

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
### Controlling humidity

High relative humidity (over 70 %) may lead to condensation inside the incubator canopy.

Condensation will impair the view of the baby.

Condensation should also be avoided for reasons of hygiene.


- 1 Press button = ON. (press again = OFF).
- 2 Green LED lights.
- 3 Relative humidity inside the incubator is shown alternately with «SEI» (request: to confirm or reset set value).
- 4 Pre-set value (60 %) flashes. Possible values: Minimum 35 % relative humidity and maximum 85 % relative humidity.
- 5 Set value as for «air temperature», page 27.
- 5 Warning indicators (LED), see pages 45 and 49 for details.

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
### Control characteristics

The skin temperature is the measurement parameter

The temperature in the incubator is controlled by the difference between set value and actual value, within the limits of a minimum of 28°C and a maximum of 37°C.

- If the set skin temperature is higher than the actual temperature (skin is cooler than desired), the temperature in the incubator will be increased by this very value.

Which means correspondingly:

- If the set skin temperature is lower than the actual temperature (skin is warmer than desired), the temperature in the incubator will be reduced by this very value.

- It is also the duration of temperature difference between set and actual skin temperature which changes the air temperature in the incubator correspondingly.

- The skin temperature of the patient is subject to frequent changes, e.g. as a result of feeding or nursing procedures. Deviations of a few tenths of a degree between actual value and set value will therefore exist until temperature equilibrium is re-established.

Therefore:

The set value for skin temperature should only be changed if the core temperature needs to be adjusted.

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

### Oxygen supply with O<sub>2</sub> concentration control

Remove probe on O<sub>2</sub> connecting hose from parking position and push fully into O<sub>2</sub> terminal unit in the pipeline system.

- 1 Press button = ON. (press again = OFF).
- 2 Green LED lights.
- 3 Yellow LED and "CAL" display flash to request calibration of the sensor.

Swing the environmental sensor housing out of the incubator as far as it will go for not more than one minute to avoid a sensor warning.

- 4 Press "CAL" button,
- 5 Yellow LED remains continuously lit.
- 6 During calibration, "CAL" and "-- --" flash alternately on the display.

A sensor warning is given if calibration has not commenced within one minute of the housing being swung out. (Press "Cal" button)

After about 45 seconds the calibration value of 21 % O<sub>2</sub> by vol. (= oxygen content of ambient air) is shown alternately with "SET" (request to confirm or reset value).

After calibration:

Swing the housing back into place within one minute, otherwise the sensor warning will be given.

The oxygen concentration inside the incubator is measured, and the actual value is shown.

The pre-set value flashes (21 % O<sub>2</sub> by vol.).


The set value can be adjusted between a maximum of 40% and a minimum of 21%.

Set oxygen concentration value as for "setting air temperature" see page 27.

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### Oxygen concentration above 40 % O<sub>2</sub> by vol.

- 1 Press button Values of up to 75 % O<sub>2</sub> by vol. can be set.
- 2 Yellow LED is lit.
- 3 Press button to set required value. The set value is shown.
- 4 If the set value is then reduced to below 40 % O<sub>2</sub> by vol., the extension of the range is automatically cancelled.
- 4 Warning indicators (LED), see page 45 for details.

If the O<sub>2</sub>-Module is switched off:

Remove O<sub>2</sub> probe from terminal unit of pipeline system and replace in parking position.

If any malfunction occurs in the oxygen control module (display "Err") or there is a sensor warning which cannot be cancelled, switch off oxygen control module.

However, if further urgent oxygen treatment is required:

please refer to the next section.

### Oxygen supply without O<sub>2</sub> concentration control

Supply the oxygen from the pipeline system via a flowmeter.

- 1 Connect hose to socket.
- 2 Set the O<sub>2</sub> flow rate on the basis of the recommended values in the table below.

Always use an oxygen meter when supplying oxygen!