



A  
HALMA GROUP  
C O M P A N Y

# ***Pulsair***

## **Non-contact Tonometer**

### **Service Instructions**



ALWAYS READ THE INSTRUCTIONS



## **CONTENTS**

1. Service Equipment Checklist and Technical Data.....1
2. Pulsair Spare parts list .....2
3. Software installation guide .....3-5
4. Hardware set-up guide (Electronic unit).....6-7
5. Dis-assembly and Re-assembly of the Pulsair2000/3000.....8
6. Dis-assembly and Re-assembly of the PulsairEasyEye.....9
7. Calibration andVerification Software Procedure.....10-22
8. Cleaning of the Optics.....23
9. Replacing the Air Filter.....24
10. Replacing the Umbilical Cord Assembly Pulsair2000/3000...25-26
11. Replacing the Umbilical Cord Assembly PulsairEasyEye.....27-28

## **CONTENTS**

12. Replacing the Mainframe and pcb Assemblies.....	29
13. Replacing the Pump and Pump Assembly Pulsair2000/3000..	30
14. Replacing the Pump and Pump Assembly PulsairEasyEye....	31
15. Set-up of Mainframe Mirror Assembly.....	32
16. Upgrade Power Supply pcb Pulsair2000.....	33
17. General Troubleshooting .....	34-38
18. Software Troubleshooting .....	39-44
19. Replacement of Case Halves Pulsair2000/3000.....	45-46
20. Replacement of Case Halves PulsairEasyEye.....	47-48
21. Replacement of thermal printer PulsairEasyEye.....	49

## ***1. SERVICE EQUIPMENT CHECKLIST AND TECHNICAL DATA***

1.	2401-P-5100	TUBE & CONNECTOR ASSEMBLY
2.	2401-P-5119	RS232 COMMUNICATION CABLE
3.	2401-P-5127	PC SERIAL PORT CABLE ASSEMBLY
4.	2401-P-6736	CALIBRATION ELECTRONIC UNIT
5.	2401-P-6787	CALIBRATION JIG Pulsair 2000/3000
6.	2401-P-6005	CALIBRATION JIG Pulsair (EasyEye)
7.	EP79-09146	14/15mm SPANNER
8.	EP79-09154	DAQ CABLE ASSEMBLY
9.	EP79-09162	DAQ CARD 1200
10.	EP59-10901	25 WAY DONGLE
11.	EP84-01130	SOFTWARE INSTALLATION DISC V2.7
12.	RP18-00094	LOCTITE 601 ADHESIVE
13.	RP18-00318	BOSTIC CLEAR ADHESIVE

### **CALIBRATION OF TEST JIG**

This calibration function can only be carried out at KEELER LTD

The test unit must be returned to Keeler for recalibration periodically {once a year} as per the due date on the back of electronic unit.

The items to be returned are as follows

1. CALIBRATION ELECTRONIC UNIT AND CABLE ASSEMBLES
2. TRANSDUCER ASSEMBLY
3. CALIBRATION JIG & DUMMY EYE
4. DAQ CABLE ASSEMBLY AND CARD 1200

PLEASE RETURN TO  
KEELER LIMITED  
CLEWER HILL ROAD WINDSOR  
BERKSHIRE SL4 4AA  
ATTN DAVE GRAHAM

## 2. Pulsair Spare parts list

### **Parts common to all models**

2401-P-6744	AIR FILTER ASSEMBLY
EP39-05869	PUMP ADAPTOR
EP39-06896	BSP\NPT ADAPTOR
EP39-05500	BULB CONTACT
2110-P-7018	WIPES (PACK OF 24)
PP02-82056	PTFE TAPE
EP79-40899	CABLE-TIE
RP18-00094	LOCTITE 601 RETAINER
EP39-04639	NEEDLE VALUE

### **Pulsair model 2000 only**

2401-P-6664	UMBILICAL CORD ASSEMBLY (2000)
2401-P-6955	P.C.B (2000)
1024-P-5062	BULB (2000)
EP18-00203	CASE LEFT (2000)
EP18-00211	CASE RIGHT (2000)
2401-P-6437	MAINFRAME ASSEMBLY
EP79-03334	SPIRAL PIN
EP79-05161	M3 INSERT

### **Pulsair model 3000 only**

2401-P-6963	UMBILICAL CORD ASSEMBLY (3000)
2401-P-6998	P.C.B (3000)
1024-P-5118	BULB (3000)
EP19-01848	CASE LEFT (3000)
EP19-01856	CASE RIGHT (3000)
2401-P-6437	MAINFRAME ASSEMBLY
EP79-03334	SPIRAL PIN
EP79-05161	M3 INSERT

### **Pulsair model EasyEye only**

2414-P-5002	UMBILICAL CORD ASSEMBLY (EasyEye)
2414-P-6002	MAINFRAME ASSEMBLY
EP29-50860	P.C.B (EasyEye)
EP39-50643	DISPLAY P.C.B
EP39-50651	SWITCH P.C.B
EP09-00501	CASE RIGHT (EasyEye)
EP09-00528	CASE LEFT (EasyEye)
1024-P-5118	BULB
EP39-50387	SCREW PLUG




### 3. SOFTWARE INSTALLATION FOR PULSAIR CALIBRATION

#### ITEMS REQUIRED

1. NATIONAL INSTRUMENTS PCMCIA ANALOGUE TO DIGITAL CONVERTER
2. NATIONAL INSTRUMENTS SOFTWARE DISC (6.1V)
3. PULSAIR CALIBRATION SOFTWARE DISC (2.7.0V)

#### LOADING PROCEDURE

It is **IMPORTANT** that the National Instruments Software for the **A.D.C** is installed first to minimise possible problems with the PCM CIA card

1. Place the **NATIONAL INSTRUMENTS CD ROM DISC** into CD-ROM drive bay.
2. **DOUBLE CLICK** My Computer on your Desktop.
3. **DOUBLE CLICK** the CD ROM icon /Drive Letter in My Computer.
4. **DOUBLE CLICK** onto the set-up icon.
5. Follow the on screen instructions. **DOUBLE CLICK** onto the install **NI-DAQ** option **FOR Windows 95 ONLY** select the language setting (this should be highlighted in blue). **DOUBLE CLICK** ok. Click next at each screen instruction until the install screen appears, then click ok.
6. When ok box has been clicked, the computer will automatically re-start.
7. After the computer has restarted, close the computer down completely and then switch off.
8. Connect the ribbon cable EP79-09154 to the ADC PCMCIA card carefully.
9. Place the national instruments PCMCIA **DAQ-1200** card carefully into the PCMCIA slot in the computer, ensuring that the card is in the correct position.(as the instruction written on top of the DAQ –1200 card) *insert card* 
10. Turn the computer back on.
11. Windows 95/98 will automatically detect the new ADC card and install the appropriate software for it.
12. **CLICK ON TO START** in main menu, then **PROGRAMS** and then select **NATIONAL INSTRUMENTS DAQ MENU**. Run the NI-DAQ Configuration Utility programme from the National Instruments DAQ program group, accept new device in socket 1 or 2.

### **3. SOFTWARE INSTALLATION FOR PULSAIR CALIBRATION**

13. Single click Configure, which will find the ADC PCMCIA card. Follow the configuration wizard. **DO NOT** change the **AI Polarity/Range or Mode Settings** to there default settings. Change the **AO Polarity** to **UNIPOLAR**. Leave **Accessory** setting at none and **OPC** to disabled, click ok to continue.
14. Then click X on the top right hand corner of the window for exit from configuration utility. Click save if requested.
15. Remove NI-DAQ CD disc from CD bay
16. Place CD-ROM DISC for Pulsair Calibration (EP84-01130) into CD-ROM drive bay.
17. **Double Click** My Computer on main desktop.
18. **Double Click CD-ROM** icon.(CAL SW-V2-7) or (KEELER)
19. **Double Click Calibration Software V2.7.0** folder(**Calibration S**)
20. **Double Click SET-UP.** Follow all on screen instructions.
21. Allow Windows to install the calibration software. Click OK to Tonometer set up, and then click the *install icon* to install software. When completed click OK. Exit all open files by clicking on the **X** in the right hand corner

### **DONGLE SECURITY DEVICE SOFTWARE INSTALLATION INSTRUCTIONS**

1. **Double Click** My Computer on main desktop
2. **Double Click CD-ROM** icon.(CAL SW-V2-7) or (KEELER)
3. Select and **Double Click** Hasp Command line Installer Latest version 120599 v3.83
4. **Double Click** Install Dongle
5. When Aladdin device Driver Installation Utility for Win32 screen instruction appears **Single Click OK**
6. Exit all open files by clicking on the **X** in the right hand corner.



### ***3. SOFTWARE INSTALLATION FOR PULSAIR CALIBRATION***

#### **SHORT CUT INSTRUCTIONS**

1. To obtain a **Short Cut** to the Calibration Software, first right click the mouse when in the main desktop.
2. Select new then Select short cut (single click)
3. At command line Select browse (single click). At browse, select folder programme files (double click).
4. At Programme files, select Tonometer, then (double click).
5. At tonometer, select eye symbol with Tonometer, then (double click).
6. At create shortcut, the command line will read "C:\ Programme Files.  
\Tonometer\Tonometer. exe" select next then (single click).
7. When select a name for shortcut appears, it will read Tonometer. Single click **Finish**.
8. A short cut symbol of an eye will now appear. When you double click the eye symbol it will take you into main Pulsair calibration software.

#### 4.

### **CALIBRATION JIG ELECTRONIC UNIT**

Take cable assembly 2401-P-5127 (PC SERIAL) and connect into fifteen way A.D.C female socket on the front of calibration electronic unit (as per fig 1)

Take the other end of cable assembly and connect into the nine way male socket on the back of laptop computer.

Plug transducer assembly into nine way female socket on the front of calibration unit, which is engraved "Pressure Transducer" (as per fig 1)

Connect RS232 communication cable (2401-P-5127) into nine way male socket on front of calibration unit, which is engraved "Tonometer" (as per fig 1)

Take the other end of the cable and connect into the back of the Pulsair tonometer mother unit.




FIG 1

#### **4. CALIBRATION JIG ELECTRONIC UNIT**

The last cable assembly to be connected is the national instruments **ADC** ribbon cable (EP79-09154)

This is connected first to the multi way socket on the front of calibration unit engraved “**ADC**” (as per fig 1).

Ensuring that the arrow on blue multi way **ADC** sockets aligns with the arrow on multi way, plug on ribbon cable.



**NOTE! WHEN CONNECTING DAQCARD INTO THE LAPTOP COMPUTER ENSURE LAPTOP IS SWITCHED OFF, AS THIS MAY CAUSE DAMAGE TO THE DAQCARD DEVICE**

Once connected, lock into position with locking clips. The smaller end is connected to the **DAQCard-1200** (EP79-09162) aligning first the white dots together as plug is pushed into socket of **DAQCARD**.

Plug blue air tube assembly (2401-P-5100) into the front of the electronic unit marked AIR INPUT

Connect dongle software security device to back of laptop computer. Locate into printer port LP1.



**NOTE! IF THE DONGLE IS NOT CONNECTED THE SOFTWARE WILL NOT BE ACCESSIBLE**

## 5. ***DIS-ASSEMBLY AND RE-ASSEMBLE OF PULSAIR***

***PULSAIR 2000/3000***

1. Remove the four M4 screws and shakeproof washers holding the base moulding to the Pulsair mother unit (as per fig 2)




FIG 2

2. While in same position remove the two M4 transit screws and shakeproof washers holding the pump assembly in position at location **B** (as per fig 3)




FIG 3

3. Remove the four M5 screws and shakeproof washers holding the top moulding to the mother unit as seen in fig 3 position **A**


*Ensure when removing top moulding that the white microswitch plunger is not mislaid*

4. Re-assemble in reverse order fitting white plunger first

## 6. *DIS-ASSEMBLY AND RE-ASSEMBLE OF PULSAIR*

### *EasyEye*

1. Remove the six M4 screws and shakeproof washers holding the base moulding to the Pulsair top unit removing the wall/desk support **A** first to reveal the two hidden screws under the support (as per fig 2A)



**FIG 2A**

2. While in same position remove the pump transit screw holding the pump assembly in position at location **B** (as per fig 2A)

*Ensure when removing top moulding that the white microswitch plunger is not mislaid*

3. Re-assemble in reverse order fitting white plunger first

## **7. CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000- EasyEye**

1. Remove the 8 screws holding the top and bottom mouldings from the Pulsair 2000/3000 mother unit (as in fig 2). Remove the 6 screws for the Pulsair EasyEye refer to fig 2A
2. Remove the 2 transit screws if fitted from the pump assembly (as in fig 3) or fig 2A for Pulsair EasyEye
3. Disconnect the blue pipe from the umbilical cord assembly, which is connected to the filter, which is located in the inner part of the Pulsair mother unit.
4. Connect the blue pipe from umbilical cord assembly and the connector assembly (2401-P-5100) from the calibration pressure gauge together and then into filter.
5. Place the Pulsair mother unit back into bottom moulding (**do not replace the screws at this stage**) as this is only a method of holding the Pulsair while calibration and repair work is taking place.
6. Place the hand unit assembly onto the calibration jig, first removing puff tube shroud and trim ring from front of the hand unit assembly. Procedure 19 and 20 (fig 34 and 34A).
7. Locate the hand unit assembly from the Pulsair **2000/3000** onto the 3 pillars **A - B** and **C** and slide into the central locking device **D** (as per fig 4) **LOCK INTO POSITION WITH LOCKING KNOB**. When calibrating Pulsair **Easyeye** use the jig shown in fig 4a
8. Connect RS232 communication cable (2401-P-5119) from calibration electronic unit to the back of Pulsair mother unit.




### **ELECTRIC SHOCK HAZARD**

Electrical power is required for the next operations


9. Fit the mains lead to underside of the Pulsair **2000/3000** mother unit (**DO NOT SWITCH ON**). Fit the lead from the power supply for the Pulsair **EasyEye** to the front or back of the mother unit

***IF MAINFRAME ASSEMBLY OR MAIN P.C.B ASSEMBLY HAS BEEN  
REPLACED REFER TO MAINFRAME MIRROR SETTING PROCEDURE  
(SECTION 15)***

## 7. ***CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000-EasyEye***



The X Y Z indecking head is now interchangeable between the two calibration jigs undo the two locking screws and remove the head from the two dowels.



10. Once the Pulsair is in position on calibration jig you must then follow the on screen instructions from the calibration software.

## 6. *CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000-EasyEye*

Select the software icon, which is identified by an eye symbol on main screen menu. Double click and this will take you into the main Pulsair menu (as per fig 5)




FIG 5

You are required to select a Pulsair instrument type; to achieve this you must first click on **select instrument type** box. A tonometer menu will appear (as per fig 6)




FIG 6

You must now select a Pulsair type to be calibrated by clicking onto appropriate Pulsair instrument type. This is verified by a black dot in the circular window



## 6. ***CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000-EasyEye***

***WARNING IF INCORRECT PULSAIR TYPE IS SELECTED THE WRONG CALIBRATION DETAILS WILL BE INSTALLED INTO THE PULSAIR EPROM AFTER CALIBRATION***

Once selected press continue box to continue, this will take you back to main Pulsair menu (as per fig 5)

You then have a choice of six instruction boxes. If calibration is required, select “Calibrate Tonometer” and follow on screen instructions.

When screen instruction has been carried out, press continue to proceed to the next screen menu. This is for setting of pump pressure follow on screen instructions (as per fig 7)



***THE PULSAIR TYPE IS DISPLAYED AT THE TOP OF ALL DATA AND INFORMATION MENUS ENSURE CORRECT TYPE HAS BEEN SELECTED***




FIG 7

Press continue box and this will take you to ***DATA ENTRY***. This information data is required for ***KEELER*** records. Fill in data were required (as per fig 8)

Operator Name:	CHRIS.W
Date:	23 September 1998
Serial Number	23131251
Product Number:	4007
PCB Serial Number:	MB8049560069
Calibration Jig number:	105
Pump Serial Number	1166504
Voltage	240

FIG 8

## **6. CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000-EasyEye**

When completed press the continue box. A switch on screen message will appear (as per fig 9)




FIG 9

**ENSURE THAT THIS INSTRUCTION IS FOLLOWED AS STATED.**

If the instructions are not followed as stated an error message will appear (as per fig 10)

The error message will take 16 seconds to appear. You must then click on the OK button, which will take you back to main Pulsair menu to start procedure again.

**ENSURE WHEN COMMUNICATION IS MADE WITH THE PULSAIR THAT THE PULSAIR IS NOT SWITCHED OFF DURING THE CALIBRATION PROCEDURE.**

**AS THIS WILL CAUSE AN ERROR BETWEEN THE CALIBRATION SOFTWARE AND THE PULSAIR AND WILL RESULT IN YOU HAVING TO EXIT THE PROGRAMME AND RESTART.**

## 7. ***CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000-EasyEye***




FIG 10

Fit dummy eye block to calibration jig in position E (as per fig 4 and 4A)

Follow on screen instructions for setting of **K FACTOR VALUE** (as per fig 11) the position is achieved by using the **X Y Z** adjusters on the calibration jig (as per fig 4 and 4A)

***REFER TO SETTING OF MAINFRAME MIRROR ASSEMBLY PROCEDURE TO ENSURE CORRECT METHOD FOR SETTING IS OBTAINED (SECTION 15)***




FIG 11

Once correct K VALUE has been obtained the **K VALUE** box will change from red to green to indicate that setting is correct.

Click onto continue box to proceed to next on screen instruction as seen in fig11A




FIG 11A

## 7. ***CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000-EasyEye***

Remove dummy eye from calibration jig position **E** and re-fit with pressure transducer block into the same position

Click onto continue box to proceed to solenoid time graph as indicated in fig 11A. This will test valve response time for Pulsair (as per fig 12)




FIG 12

The time period in which the solenoid valve will open is set automatically when the Pulsair type is selected at beginning of calibration procedure.

This time period will change from 3.1ms normal time and 3.4ms turbo time on the **Pulsair 3000** and the **Pulsair EasyEye** then to a higher time period of 4ms normal time and 4.3ms turbo time on **Pulsair 2000 mod state 1to10** and **Pulsair 2000 mod state 11**.

The time period is tested in steps of .3ms until correct time period is achieved and the Pressure difference is not > than the acceptable tolerance set in the software,

Each level is verified with 8 test puffs, (as seen in fig 12).

Once completed click continue box to proceed to next instructions.

## 6. CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000- EasyEye

You must now click onto the **OK** box below **Calibration eeprom Values Written** this will now transfer calibration data into eeprom (as per fig 15).

Once completed the main data information box will appear again (as per fig14)

A **tick** will now be displayed in the **EEPROM Status** box. This will now verify that the eeprom data has been successfully transferred as seen next to the EEPROM Data Written statement




FIG 15




FIG 15A

## 7. ***CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000- EasyEye***

The next graph to appear is the Align Timing Regress graph (as per fig 13)

The Pulsair will fire 32 times automatically to produce two graphs, the plot of external pressure versus the plot of internal pressure (as seen in fig 13)

Once completed the plot of internal versus external ADC VALUES will appear on right hand side of main graph.

All other calibration details will appear at bottom of graph.

***ENSURE CORRECT PULSAIR TYPE IS SELECTED AT THE TOP OF GRAPH*** before the *continue* button is pressed




FIG 13

## 7. ***CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000- EasyEye***

The tonometer's Calibration data details will appear next (as per fig 14)

Ensure all details are correct before clicking onto **WRITE EEPROM VALUES**,

At this stage all calibration details are ready to be transferred to the Pulsair eeprom,

Parameter	Value
Instrument Type	T2000 (Mod 11)
Operator name	CHRIS.W
Calibration Date	23/09/98
Serial Number	23131251
Product Number	4007
PCB Serial Number	MB8049560069
Calibration Jig Number	105
Pump Serial Number	1166504
Operating Voltage	240
Best K Value	0.122
Photodiode Values (A,B,C)	1790, 217, 1770
Normal Pulse Time (mS)	4.0
Turbo Pulse Time (mS)	4.3
Average Regression Slope	14.167
Average Regression Intercept	4.437
Skew	4.0
Skew Range	0.1
Software Version	2.4
Calibration Ae Value	-2.15
Calibration Be Value	0.68

Parameter	Value
Best Contrast	3
KIEE normal	49
KIEE Subflex	38
Normal Pulse Time	66
Turbo Pulse Time	71
aa Value	-152
ba Value	12582
Skew Value	6

EEPROM Status

EEPROM Data Written ☐

Print Tables Save Tables Write eeprom Values Cancel Continue

FIG 14

A verification box will now appear to ensure you that the details will now be transferred (as per fig 15)

**ENSURE WHEN CALIBRATION IS COMPLETED THAT THE PUMP ADAPTOR LOCK NUT IS SEALED WITH LOCTITE 601 ADHESIVE**





## 6. CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000- EasyEye

You must now click onto the **OK** box below **Calibration eeprom Values Written** this will now transfer calibration data into eeprom (as per fig 15).

Once completed the main data information box will appear again (as per fig14)

A **tick** will now be displayed in the **EEPROM Status** box. This will now verify that the eeprom data has been successfully transferred as seen next to the EEPROM Data Written statement




FIG 15




FIG 15A

## 7. CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000- EasyEye

At this stage if the calibration data has not been transferred, the continue button will not be highlighted and will not allow you to proceed,

If you were to click onto Print tables or Save tables, a warning message will appear (as per fig 16) as a reminder that you must write to the eeprom before the completion of the calibration.

When calibration is completed, click continue box to return to main menu.




FIG16

## 7. ***CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000- EasyEye***

To verify calibration you must start at main Pulsair menu as before

Then **Select Instrument Type** where you must then select correct Pulsair to be checked i.e. Pulsair 2000 mode state 0-10...Pulsair 2000 mode state 11...Pulsair 3000... Pulsair *EasyEye*

Then select **CHECK CALIBRATION** by clicking onto the highlighted Check Calibration box.

All on screen instructions except for the setting of pump pressure (**see note**) are the same as the previous calibration instructions, but will change after the Align Timing Regress graph as (per fig 13).

At this stage a new graph will appear, which is the comparison graph (as per fig17). This will verify the calibration

**NOTE! The pump pressure on screen instruction will not appear in this part of the procedure, as the pump pressure must not be adjusted in the verification of the Pulsair calibration.**




FIG 17

The software will compare old eeprom values against the new eeprom values, which will then be calculated to see if they correspond

## 7. CALIBRATION AND VERIFICATION PROCEDURE FOR PULSAIR 2000-3000-EasyEye

The pressure errors of the Pulsair will now be displayed on the graph (as seen in fig17)

The errors at 10 and 20mmHg should not exceed more than **plus or minus** 1 mmHg. If this occurs the Pulsair must be recalibrated.




FIG 17A

If when checking the old eeprom values against the new values (as seen in fig 17A) the Best Contrast is set at a value of 255 and the KIEE normal is set at a value of 255 it would indicate that the calibration data is corrupt.

The Pulsair must now be recalibrated, as the unit will not auto fire when focussed on the eye.

Once completed, click continue box to return to main Pulsair menu

When finished, **DISCONNECT MAINS LEAD** from mother unit. Remove RS232 communication cable from back of Pulsair.

Refit blue pipe from umbilical cord assembly to filter assembly and add tested label (EP39-05372) to top of pump assembly diaphragm block

Remove hand unit from calibration jig

Re-assemble puff tube shroud and trim ring to hand unit and locking screw. Refer to fig34 and fig34A

Re-fit transit screws to pump and 8 screws holding the top and bottom moulding (**Pulsair 2000/3000**)  
 Re-fit the 6 screws holding the top and bottom moulding to the **Pulsair Easyeye** in reverse order.  
 Refer to **procedure 5 and 6**.

## 8.


## ***CLEANING OF THE PULSAIR OPTICS***

1. First remove puff tube shroud and eyepiece shrouds before cleaning procedure can take place.



**Ensure when cleaning the lenses that they are not touched with the fingers as these are coated lenses and will easily smear**

2. Using the cleaning wipes (2110-P-7018) in a circular motion remove any dirt deposits from front of puff Tube assembly and viewing lens assembly (as per fig18)



## 9. ***REPLACEMENT OF AIR FILTER***



### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug  
Should be removed before starting next operation.

1. Remove top and bottom mouldings as described in procedure 5 and 6
2. Disconnect blue tube from both sides of filter assembly (as in fig 19)
3. Remove filter from bracket and replace with new filter assembly (2401-P-6744)  
Reconnect blue tubes
4. Refer to calibration procedure 7 for recalibration instructions.
5. Refer to re-assemble of Pulsair procedure 5 and 6




FIG 19



## **10. UMBILICAL CORD ASSEMBLY REPLACEMENT**

### **PULSAIR 2000/3000**



#### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug Should be removed before starting next operation.

1. Remove the Pulsair top and bottom mouldings as described in procedure **5 and 6**
  2. Disconnect blue tube from filter in mother unit.
  3. Remove either M4 nut from earth stud on microswitch bracket or M4 screw from side of transformer
  4. Remove cable ties holding earth wiring to wire loom, unplug 8 way Molex connector from power supply P.C.B.
  5. Remove M3 screw holding “P” clip and blue pipe from power supply P.C.B
  6. Remove umbilical cord assy from mother unit.
  7. At hand unit end, remove puff tube shroud, trim ring, eyepiece shroud and M3 locking screw at bottom of hand moulding (as per fig 34).
  8. Remove top hand unit moulding
  9. Disconnect Molex connector containing 8 wires from P.C.B.(as per fig 20)
  10. Remove the M3 screw from centre of mainframe holding earth wires (as per fig 20).
  11. Remove the screw holding yellow wire and bulb contact to top of mainframe. Cut the blue tube from hand unit value to mainframe and then remove the umbilical cord assembly.
- NOTE!** The blue tube removed from hand unit value must be replaced before reassemble can take place.
12. Re-assemble new umbilical cord assy in reverse order. Refer to calibration procedure 7. For re-calibration.



**10. UMBILICAL CORD ASSEMBLY REPLACEMENT**  
**PULSAIR 2000/3000**




FIG 20

## 11. **UMBILICAL CORD ASSEMBLY REPLACEMENT** *EasyEye*



### **ELECTRIC SHOCK HAZARD**


Electrical power is not required for this operation and the mains supply plug should be removed before starting next operation.

1. Remove the Pulsair top and bottom mouldings as described in procedure **5 and 6**
2. Disconnect blue tube from filter in the mother unit.
3. Unplug the 8 way Molex connector from power supply P.C.B.
4. Remove the umbilical and cable gland assemble from the mother unit.
5. At the hand unit end, unscrew the puff tube shroud and the end cap from the hand unit mouldings refer to Fig 34A, remove the four locking screws holding the top hand unit moulding to the bottom moulding.
6. Remove the top hand unit moulding and the front panel assembly
7. Unplug the 8-way molex connector from the P.C.B and remove the clip from the coil on the handunit value (as per fig 20A).
8. Remove the M3 screw from centre of mainframe holding it to the bottom hand unit moulding (as per fig 20A).
9. Cut the blue tube from hand unit value to the mainframe, disconnect the blue tube from the reservoir you now remove the umbilical cord assembly.

**NOTE!** The blue tube removed from hand unit value must be replaced before reassemble can take place.

12. Re-assemble new umbilical cord assy in reverse order. Refer to calibration procedure 7. For re-calibration.

# 11. *UMBILICAL CORD ASSEMBLY REPLACEMENT* *EasyEye*



**FIG 20A**

## 12. REPLACEMENT OF MAINFRAME & P.C.B. ASSEMBLIES



### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug Should be removed before starting next operation.

1. Disconnect the umbilical cord assembly from hand unit as described in the previous procedure 10 and 11. (from operation 7 onwards is required)
2. Remove the four screws holding P.C.B. to mainframe and earth wire; disconnect blue tube to transducer. Replace defective part (mainframe assembly or main hand unit P.C.B)
3. Ensure when re-assembling that the serial number label is removed from the P.C.B and re attached to the top of the mainframe. The insulation tape must be fitted to the P.C.B with 7mm cut out to clear solder earth pad (as per fig 21) and blue tube to transducer is reconnected from mainframe assembly with brass restrictior.
4. Fit M2 washer between bulb contact and P.C.B prior to the fitting of the mainframe; ensure the second bulb contact is fitted to the mainframe and the earth lead is attached.
5. Re-assemble umbilical cord assembly to bottom hand unit moulding.
6. Fit the dummy handunit case (EP29-04207) Pulsair 2000/3000(EP39-50272) Pulsair EasyEye prior to calibration so setting of mainframe mirror assembly can be carried out. (See fig 24 procedure 15.) and fig 24A for Pulsair Easy Eye ensure the calibration ring is attached
7. Refer to calibration procedure 7 for recalibration instructions.

(CUT OUT 7MM SQUARE IN RP99-99325 TAPE TO CLEAR SOLDER EARTH PAD)





FIG 21

### 13.

## ***REPLACEMENT OF PUMP PARTS PULSAIR 2000/3000***



### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug should be removed before starting next operation.

1. Disconnect blue tube fitting from elbow connection on side of pump adaptor as in (fig 22A)
2. Remove the three **mains** wires connected to the coil and valve assembly (**BLACK, WHITE, GREEN/YELLOW**) as in (fig 22B).
3. Remove coil and valve assembly from side of pump adaptor as in (fig 22C), using 14 and 15mm spanners.
4. Remove filter silencer from underside of pump adaptor fitted to pump block (as per fig 22D).
5. Remove pump adaptor (EP39-05869) and BSP/NPT adaptors (EP39-06896) from pump (as per fig 22E) using 14mm spanners  
*ENSURE WHEN REFITTING BSP/NPT adaptors a small amount of PTFE tape (PP02-82056) is secured around threads*
6. Replace faulty components and re-assemble in reverse order.
7. Refer to calibration procedure 7 for re-calibration instructions.




FIG 22

## 14.


### ***REPLACEMENT OF PUMP PARTS*** ***PULSAIR EasyEye***



#### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug should be removed before starting next operation.

1. Disconnect blue tube fitting from elbow connection on side of pump adaptor as in (fig 22A/A)
2. Remove the two wires connected to the coil and valve assembly (**BLACK, WHITE**) as in (fig 22A/B)
3. Remove coil and valve assembly from side of pump adaptor as in (fig 22A/C) using a 14 and 15mm spanners.
4. Remove filter silencer from underside of pump adaptor fitted to pump block (as per fig ).
5. Remove pump adaptor (EP39-05869) and BSP/NPT adaptors (EP39-06896) from pump (fig 22A/D) using 14mm spanners
6. *ENSURE WHEN REFITTING BSP/NPT adaptors* a small amount of PTFE tape (PP02-82056) is secured around threads
7. Replace faulty components and re-assemble in reverse order.
8. Refer to calibration procedure 7 for re-calibration instructions.



**FIG 22A**

## 15. **SETTING OF MAINFRAME MIRROR ASSEMBLY AND TARGET FOR PULSAIR EasyEye**

1. Place the hand unit assembly onto calibration jig as shown in (FIG 4 or 4a). **NOTE!** At this stage ensure that the dummy hand unit case is fitted.
2. Refer to calibration procedure 7 and follow on screen instructions for calibration until the setting of 'K' value is required (as per fig 23)
3. Loosen the two screws holding the mirror assembly in position on mainframe (as per fig 24) and adjust until the best reading is obtained on 'K' value. Refer to example shown below.





FIG 23

4. A and C photodiode value should be set equally first before obtaining the best 'K' value. Once achieved, lock mirror assembly into position as per fig 24.
5. Only on the Pulsair EasyEye will you see the bow tie target when you view through the Eye lens, the target must be adjusted to the central position as shown in fig 24A.



## 16. UPGRADE OF POWER SUPPLY (P.C.B.)

### THIS PROCEDURE APPLIES TO PULSAIR 2000 ONLY



1. Verify voltage at diode D6 (anode side of diode) and PL4/1 connectors on power supply P.C.B. as per figure 25. If voltage is 17 volts, the power supply P.C.B. must be modified.



#### ELECTRIC SHOCK HAZARD

Electrical power is not required for this operation and the mains supply plug Should be removed before starting next operation.

2. First remove P.C.B. from the mother unit by disconnecting wire looms PL1-PL4 and PL3 From the P.C.B plug connections.
3. Remove the five-m3 screws holding the power supply P.C.B. in position.
4. The P.C.B. can now be removed.
5. The next part of this procedure should only be carried out on **Pulsair 2000** as the Pulsair 3000 has a higher drive voltage of 23.8 volts on the solenoid valve
6. REPLACE resistor R4 from 1k5 to 1k8 OHMS (**on Pulsair 2000 only**) this will increase voltage at D6 from 17 volts to 19.95 volts which enables the solenoid valve to be more responsive.
7. Re-assemble power supply P.C.B. in reverse order.
8. Refer to calibration procedure 7. For re-calibration.





## **17. GENERAL TROUBLE SHOOTING**

### ***PULSAIR WILL NOT AUTO FIRE***

Envelope of bulb discoloured - *Requires replacement of bulb*


Front puff tube lens contaminated with splash back from tear film - *Requires cleaning off the puff tube lens*

Puff tube lens orientation incorrect or broken (as per fig 26) - *Requires realigning or new mainframe is required*

Calibration data erased in eeprom - *Requires recalibration*

No red light from hand unit - *Check 6volt line to bulb and bulb contacts are not broken*

**FIG 26**



### ***PULSAIR HAND UNIT FAILURE***



#### **ELECTRIC SHOCK HAZARD**

Electrical power is required for the next test operations and the mains supply plug should be inserted into the *Pulsair* mother-unit before starting.

## **17. GENERAL TROUBLE SHOOTING**

### **Hand unit display will not function Pulsair 2000/3000**

Verify the T2amp fuse inside the *Pulsair* mother unit  
*If faulty remove and replace.*

Verify crimps on fuse holder  
*If faulty remove crimps and solder wires to fuse holder.*

Verify the 10volt AC supply from transformer to power supply. The red and orange wires are required for this test (Edge connector PL1 pins 3&4)  
*If there is low or no voltage verified at this point, the transformer must be replaced.*

Verify the 5volts DC supply from power supply P.C.B in mother unit, to hand unit. The black and red wires are required for this test (PL4 edge connector on power supply PCB pins 1&4)  
*If there is low or no voltage verified at this point, the power supply P.C.B must be replaced*

Verify the 5volts DC supply from power supply P.C.B to the inside of the hand unit main P.C.B. The black and red wires are required for this test (Edge connector PL2 on main P.C.B 1&4)  
*If there is low or no voltage verified at this point, the umbilical cord assembly must be replaced*

### **Hand unit will not puff on demo Pulsair 2000/3000**

Verify the 400ma fuse inside mother unit  
*If faulty remove and replace*

Verify crimps on fuse holder  
*If faulty replace crimps and solder wires to fuse holder*

Verify the 20volts AC supply from transformer to power supply P.C.B. The Grey and blue wires are required for this test (Edge connector PL1 on power supply P.C.B pins 1&2)  
*If there is low or no voltage verified at this point, the transformer must be replaced*

## **17. GENERAL TROUBLE SHOOTING**

Verify the 19.9volts DC supply on the Pulsair 2000 and 23.8volts DC on the Pulsair 3000 from the power supply P.C.B in the mother unit to the hand unit solenoid valve. The black and blue wires are required for this test

(Edge connector PL4 on power supply PCB pins 1&2)

*If there is low or no voltage verified at this point, the power supply P.C.B must be replaced*

Verify the 19.9volts DC supply on the Pulsair 2000 and 23.8volts DC on the Pulsair 3000 in side the hand unit on the main P.C.B. The black and blue wires are required for this test (Edge connector PL2 on main P.C.B pins 1&3)

*If there is low or no voltage verified at this point, the umbilical cord assembly must be replaced*

Verify inside Pulsair mother unit the inline fuse to the pump (T250ma)

*If faulty remove and replace*

Verify airflow from pump.

*If no airflow from the pump the pump must be replaced*

### **PULSAIR MOTHER UNIT NO POWER Pulsair 2000/3000**

Verify fuses on under side of Pulsair mother unit

*If faulty remove and replace*

Verify mainswitch on side of Pulsair mother unit

*If faulty remove and replace*

Verify Pulsair mains power cord lead

*If faulty replace*

### **INCONSISTENT AIR FLOW**

Check pump adaptor for air flow

*If faulty replace*

Check pump assembly for airflow

*If faulty replace*

## **17. GENERAL TROUBLE SHOOTING**

### ***BATTERY not charging Pulsair easyeye***

Check 5amp fuse inside Pulsair  
*If faulty replace*

### ***PRINTER will not function Pulsair easyeye***

Check paper is present inside printer  
Check lid is shut  
Check printer is switched on

*If Printer is faulty refer to section 21*

### ***NO POWER Pulsair easyeye***

Check voltage from wall mounted P.S.U  
12volts should be present.

*If faulty replace*

### ***NO POWER Pulsair easyeye MOTHER UNIT***

Verify the voltages from the mother unit power supply P.C.B  
Check the voltages at the molex connector SK1

SK1 pin8 (orange) = common  
SK1 pin7 (blue) = 24volts  
SK1 pin6 (yellow) = 6volts  
SK1 pin5 (red) = 5volts  
SK1 pin4 (black) = common  
SK1 pin3 (white) = -12volts

*If faulty replace the P.S.U*

## ***17. GENERAL TROUBLE SHOOTING***

### ***NO POWER Pulsair easyeye HAND UNIT***

Verify the voltages at the hand unit P.C.B

Check the voltages at the Molex connector SK2

SK2 pin8 (orange) = common

SK2 pin7 (yellow) = 6volts

SK2 pin6 (blue) = 24volts

SK2 pin5 (red) = 5volts

SK2 pin4 (black) = common

SK2 pin3 (white) = -12volts

*If faulty replace the pulsair umbilical cord assembly (Procedure 11)*

## 18. *TROUBLE SHOOTING THE SOFTWARE*

When calibrating or verifying the calibration of the *Pulsair* the software has been designed to show warning messages or graphs to give a indication that there is a problem with the *Pulsair* at this stage of the procedure (as seen in fig 27)



FIG 27

**Communication Link Error** (as shown in fig 27) is an indication that the RS232 communication cable has not been connected between the *Pulsair* and calibration electronic unit, or when switching on the *Pulsair* after being instructed by the software, has not taken place in the time period requested.




FIG 28


The graph above - **Puff Pressure Error** (fig 28) indicates that there is no pressure seen by the pressure transducer on calibration jig. This would suggest that there is a pressure problem from pump (no air pressure!) or that there may be a connection problem with air pipes from air filter! If the correct pump pressure is not set to 21psi this warning graph will then also appear.

## 18. *TROUBLE SHOOTING THE SOFTWARE*

The graph below (Pressure versus Time) shows how the Pulsair calibration software has tried to set a Pulsair 3000 valve time out setting to a time period of 3.1ms (as shown in fig 29) this has not been achieved as illustrated.

Once this has occurred the software is designed to go back above the point of failure to try and reset the time period by reducing the time period's adjustment in step of .06ms. This has also failed (as seen in fig 29).

At this stage the umbilical cord assembly must now be replaced, as the valve is faulty. Refer to procedure 9 and 9A for replacement of umbilical cord assembly.




**FIG 29**

The graph reference (fig 30) **Puff PressureWarning** is an indication that the pump pressure has fallen out of the recommended calibration specification pressure setting of 21psi.

The reason for this could be as follows! Low output pressure from pump assembly, Pump adaptor faulty or requires internal cleaning. Air filter blocked requires replacing.  
Leak in main air system (mother unit or hand unit)!

## 18. ***TROUBLE SHOOTING THE SOFTWARE***

When fault of pressure problem has been identified and then rectified the Pulsair must be recalibrated.



**FIG 30**

The graph reference (fig31 and 31A) **INTERNAL PUFF WARNING** and **Run time error 6 overflow** indicates that there is no internal pressure profile or a low signal from the hand unit's P.C.B pressure transducer (as seen in fig 31 and 31A).

**Run time error 6 overflow will only appear on Pulsair mod state 0-10 units**

A faulty internal transducer, the air pipe not connected to transducer, blockage in air pipe restrictor could be the cause of this problem.



### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug should be removed before starting next operation.

If found to be the transducer on the hand unit Printed circuit board, the printed circuit board must be replaced

Reference procedure 12. Replacement of Mainframe and P.C.B assemblies.

After replacement the Pulsair must be recalibrated - refer to recalibration procedure 7.



18. *TROUBLE SHOOTING THE SOFTWARE*




FIG 31




FIG 31A

## 18. *TROUBLE SHOOTING THE SOFTWARE*

In reference to (fig 32) **Alignment of K value** there is no reading shown in the K value box, only 0.000 in red.

This is a warning that the bulb is not functioning correctly, which may be caused by a loose or broken bulb contact, a blown bulb or a faulty bulb

Verify the voltage to bulb from power supply (check 6-volt power line),

Verify the main frame optical system; check fig 24 or 24a of procedure15. (**Setting of mainframe mirror assembly**).

Ensure transducer has not been left on calibration jig when setting **K value**, this would give a similar fault




FIG 32

If the dongle security device is not connected into the back of laptop computer which is located into printer port LP1, a Error message will appear (as per fig 33)




Fig 33

## 18. TROUBLE SHOOTING THE SOFTWARE

When calibrating a Pulsair you are requested to select an instrument type a Pulsair 3000 Pulsair *EasyEye* or a Pulsair 2000 mod state 11 or a Pulsair 2000 mod state 0 to10 (refer to fig 6).

If the incorrect Pulsair type has been selected when calibrating a Pulsair 2000 mod state 0 to10 the Pulsair calibration software will allow you to calibrate the Pulsair until you are requested to write calibration values to the eeprom.

The software has been designed to detect at this stage that there has been a selection error and will reject your request to write calibration to the eeprom.

The same effect will be seen if trying to calibrate a Pulsair *EasyEye*/ 3000/2000 mod state 11 and selecting a 2000 mod state 0-10 instrument type

A warning will now be displayed as seen as **Comms Error!** Error writing values to eeprom (fig 33A)




FIG 33A

To correct this error you must first click onto the **OK** on the **Comms error** box and then click onto the continue box to return to the start of the calibration procedure 7.

You must then re-select the correct Pulsair to be calibrated and re-start the calibration.

The software will not be able to distinguish between the Pulsair *EasyEye* Pulsair 3000 and the Pulsair 2000 mod state 11 after calibration if the incorrect Pulsair type has been selected

To verify that the correct type after calibration has been selected you must auto fire the Pulsair on to the glass eye.

The Pulsair 3000 display will show a flashing Er and the 2000 mod state 11 will show a zero and two flashing dots.

## 19.

### ***REPLACEMENT OF CASE HALVES PULSAIR 2000/3000***




#### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug should be removed before starting next operation.

1. Unscrew puff tube shroud and remove (as per fig 34)
2. Remove trim ring by sliding over case halves
3. Unscrew eye piece shroud and remove
4. Remove M3 locking screw from bottom of top case half
5. Remove top case half
6. Unscrew M3 centre locking screw holding mainframe and remove.
7. Remove bottom case half.
8. Replace bottom case halves with new adding spirol pin(EP79-03334) and M3 insert (EP79-01561) , add bostic clear adhesive to slot in case half to retain display p.c.b
9. Re-assemble in reverse order

19.

**REPLACEMENT OF CASE HALVES**  
**PULSAIR 2000/3000**



## ***REPLACEMENT OF CASE HALVES*** ***PULSAIR EasyEye***




### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug Should be removed before starting next operation.


1. Unscrew the puff tube shroud from the case halves and remove (as per fig 34A)
2. Unscrew end cap on unbilical cord assembly and remove from case halves  
(As per fig 34A)
3. Push the Allen key supplied through the four screw plugs to access the cap headed screws  
Unscrew the four screws holding the two cases halves together.
4. Remove the top case half
5. Remove the front panel assembly.
6. Unscrew the M3 centre locking screw holding the mainframe and P.C.B.
7. Remove the bottom case half
8. Replace top and bottom case half adding new spiral pin (EP79-03334) to top case half and led covers (EP39-50387) if required
9. Re-assemble in reverse order.

20.

## ***REPLACEMENT OF CASE HALVES*** ***PULSAIR EasyEye***



**FIG 34A**








## ***21. REPLACEMENT OF THERMAL PRINTER (PULSAIR EasyEye)***



### **ELECTRIC SHOCK HAZARD**

Electrical power is not required for this operation and the mains supply plug should be removed before starting next operation.



1. Unscrew the top assembly from the base as described in section 6
2. Unplug the printer switch wiring loom from the powers supply P.C.B
3. Unplug the two ribbon cables from the lid which are connected to the power supply P.C.B
4. Remove the three screws from inside of the lid holding the printer in position
5. Remove the printer p.c.b located on the side of the top moulding.
6. Open the printers lid and then extract the printer.
7. Before fitting the new printer remove the lid attached, this is not required and can be discarded. Attach P.C.B from old printer prior to re-fitting.
8. Relocate the printer into the top assembly, before tightening the screws shut the printer's lid to allow the printer roller and gearing to be aligned correctly.
9. Re assemble in reverse order