



AESCULAP®

MEDITEC



## QL 10 DIRECTIONS FOR ALIGNMENT AND CHECKS

### GENERAL

- Layout of parts see sketch.
- Before adjusting the QL 10 attachment, check if slit lamp is adjusted properly i.e. slit projection and focal plane must be in the rotation center (focal stick) of the slit lamp.

### OPTICAL ADJUSTMENT OF THE QL 10 ATTACHMENT


#### HeNe Aiming beam adjustment

1. The HeNe is put centric through the HeNe beam expander (3) together with the beam doubler (2).
2. Lens of the HeNe beam expander in central setting.
3. The two subbeams are adjusted symmetrically with the HeNe passive reflector (4) towards the center of the expanding lens (5). The test is performed with a piece of transparent graph paper. While being illuminated adequately (from underneath) the lens holder is drawn as a dark rim and permits an exact setting of the two HeNe beams.
4. After putting the transparent graph paper away, both HeNe subbeams have to hit the achromat (6) symmetrically in the center. The testing is performed with an adequate piece of transparent paper, again.
5. By reciprocal displacement of beam doubler (2) and passive reflector (4) it must be achieved that the beam is directed perfectly. The two HeNe beams should leave the QL 10 housing mainly parallel.
6. By using the HeNe YAG passive reflector (11) the aiming beam is directed onto the slit lamp center. The HeNe has to be within the smallest aperture of the slit lamp.
7. By shifting the expanding lens (5) the HeNe aiming beam is to be directed onto the focal stick.

. To be continued.../2

YAG Adjustment

1. Adjust beam with YAG mirror (8) through the filter wheel on the beam divider (10), and check the centric pass on the filter wheel with burn paper.
2. The YAG beam has to hit the expander lens exactly in the center between the two HeNe spots. Check with burn paper.
3. Position of the expander beam has to be controlled in front of the achromat (6) by aid of IR Phosphor card.
4. At inclined pass, correct with YAG mirror (8) and beam divider (10).
5. Check the YAG focus with burn paper at min. energy level. The YAG focus must coincide with HeNe focus and focal plane. Minor YAG variations (to the rear) are permissible.
6. Micro correction of the YAG/HeNe focus concurrence is performed by displacing the HeNe expander lenses (3) to each other.
7. If YAG and HeNe are centricly adjusted no x/y variations must be visible on the burn paper.
8. At the max. energy level an optical breakdown must be perceptible in the air.
9. Control of the energy distribution 1 cm in front of the focus with burn paper.  
  
Energy distribution has to be cyclic and symmetrical.



 Institut für Laser- und Strahlentechnik Obere Kirschgasse 3 8501 Hechingen		Projekt: <b>Strahlenverlauf v. MQL 10          (TOPCON SL-2E)</b>
Zeichnung: <b>10</b>	Blatt: <b>1</b>	MQL 330 4379
Datum:	Änderung:	Folie: <b>3</b>



AESCULAP<sup>®</sup>

MEDITEC




- 4 -

## ELECTRONIC ADJUSTMENT OF MQL 10

ATTENTION when opening the MQL 10:

At the HeNe cavity there may still be high voltage. It must be discharged, but not directly by short circuit, because this may damage the power supply. Discharge it by resistor ( $> 100 \text{ k}\Omega$ ).

### 1. Adjustment of the voltages for burst 1, 2, 3 Circuitry for burst measurements:



Connect the voltmeter to the Iskra-card PIN 26c (+) and 32a (-). (Attention: range 1000 V =).

Arrange the measuring circuitry in the beam path of the laser. (Even a Joulemeter can be used, because at 2 bursts the energy will be double, at 3 bursts three times as much.) Connect The foot switch, put the energy selector (handwheel) in position 10. Switch on the device and activate the "ready" switch.

#### a) Adjustment of voltage at burst 1

Adjust the voltage between PIN 26 c and 32a by potentiometer P4, so that one burst is produced (note this voltage). After that determine the upper limit of voltage at which just two shots per burst are produced (note down this voltage, too). Adjust the voltage midway between both determined limits by the help of P4.

To be continued.../5

If you want to ascertain the number of shots per burst, you have to trigger one laser pulse by pressing the foot switch.

**b) Adjustment of voltage at burst 2**

Press the "burst" button once. (burst indicator LED 2 is now illuminated).

Adjust the voltage between PIN 26c and 32a by potentiometer P3 so that just three shots per burst are emitted (note this voltage). The lower limit was already determined at burst 1.


Set the voltage midway between both limits by the help of potentiometer P3.

**c) Adjustment of voltage at burst 3**

Press the "burst" button once again. (Burst indicator LED 3 is now illuminated).

Adjust maximum voltage (1000 V) between PIN 26c and 32a by potentiometer P2. If 4 shots per burst can be measured now, reduce voltage until only 3 shots are emitted. The lower limit was already determined at burst 2. Adjust voltage midway between these limits by P2.

If potentiometer P2 stands at stop (max. voltage) and you can only measure 3 shots per burst, adjust voltage by potentiometer P2 midway between maximum and voltage determined at burst 2.




## 2. Adjustment of energy check indicator

Connect footswitch, put the laser energy selector in position 10, switch on the instrument and deactivate the "ready" switch.

### a) Adjustment of the preamplifier V101

Connect oscilloscope at JC 101 PIN 6. Push energy check button. Adjust output voltage to 13 V by R 106. (Only necessary at factory.)



### b) Adjustment of the energy check indicator for 10 mJ

Put the laser energy selector in position 10. Push energy check button. Align threshold by R 117 so that LED for 10 mJ is surely illuminated. (Press energy check button three times; LED for 10 mJ has to glow each time).

### c) Adjustment of the energy check indicator for 0.5 mJ

Put the laser energy selector in position 0.5. Push energy check button. Adjust threshold by R 128 so that LED for 0.5 mJ is surely illuminated. (Same check as under b).

### d) Adjustment of the energy check indicator for 1 mJ

Put the laser energy selector in position 1. Push energy check button. Adjust threshold by R 127 so that LED for 1 mJ is surely illuminated. (Same check as under b).

### e) Adjustment of the energy check indicator for 2 mJ

Put the laser energy selector in position 2. Push energy check button. Adjust threshold by R 126 so that LED for 2 mJ is surely illuminated. (Same check as under b).

To be continued.../7

f) Adjustment of the energy check indicator for 4 mJ


Put the laser energy selector in position 4. Push energy check button. Adjust threshold by R 125 so that LED for 4 mJ is surely illuminated. (Same check as under b).

g) Adjustment of the energy check indicator for 7 mJ

Put the laser energy selector in position 7. Push energy check button. Adjust threshold by R 124 so that LED for 7 mJ is surely illuminated. (Same check as under b).


3. Test of the entire machine

- a) After adjustment of the energy check indicator, the energy check has to be tested once more. If you turn the energy selector to position 0.5 and you push the energy check button, only the LED for 0.5 mJ must be illuminated. Otherwise you have to adjust the thresholds again.
- b) Energy selector in position 10. Push energy check button. After approx. 8 sec. all LEDs have to switch off.
- c) Energy selector in position 10. Push energy check button. Activate "ready" switch. All LEDs should switch off immediately.
- d) Device is in Ready position. Energy selector is in position 10. Push energy check button - none of the energy check LEDs must glow.




		Oskar-Brückner-Str. 3 85041 Herrsching	8 MOL 10
GOMTEC OPTISCHE MESSTECHNIK		Blockschaltbild	MOL 2104.577
Name: _____ Matrikel-Nr.: _____ Datum: _____	Blatt: _____ von _____	Gezeichnet: _____ Geprüft: _____	Gezeichnet: _____ Geprüft: _____






0 +5V @ +15V @ -15V

MEDITEC		Ober- & Bergstraße 3		83041 Heroldsberg	
RENHARDT TRITZEL		QM-Karte:		Schallbild	
MQL 210 3127		MQL 10		g	





1)	St. 508/3	-----	QS a 28
2)	St. 508/4	-----	QS c 28
3)	St. 508/5	-----	QS a 30
4)	St. 508/6	-----	QS c 30
5)	St. 508/7	-----	QS c 32
6)	St. 508/8	-----	QS a 32
7)	St. 501/1	-----	QS c 16
8)	St. 501/2	-----	QS a 10
9)	St. 501/4	-----	QS a 8
10)	St. 505/1	-----	QS c 24
11)	St. 505/3	-----	QS c 20
12)	St. 502/1	-----	QS c 18
13)	St. 502/3	-----	QS a 26
14)	St. 502/5	-----	QS a 20
15)	St. 506/1	-----	QS c 2
16)	St. 504/1	-----	QS c 12
17)	St. 504/2	-----	QS c 22
18)	St. 504/4	-----	QS c 26
19)	St. 503/1	-----	QS a 24
20)	St. 503/3	-----	QS c 10
21)	St. 503/5	-----	QS c 14
22)	St. 507/1	-----	NE a 10
23)	St. 507/2	-----	NE b 10
24)	St. 507/3	-----	NE c 10
25)	St. 507/4	-----	NE a 8
26)	St. 507/5	-----	NE b 8
27)	St. 507/6	-----	NE c 8
28)	24V-Trafo oben	-----	ISKRA c 8
29)	24V-Trafo unten	-----	ISKRA c 6
30)	V 501	-----	QS a 22
31)	V 501	-----	ISKRA a 10
32)	R 501	-----	ISKRA c 12
33)	R 501	-----	ISKRA c 10
34)	C 502	-----	QS c 4
35)	YAG Trigger	-----	ISKRA c 26
36)	YAG	-----	ISKRA a 32
37)	YAG	-----	ISKRA c 18
38)	YAG	-----	ISKRA c 20
39)	ISKRA c 2	-----	QS a 14
40)	ISKRA c 4	-----	QS a 2
41)	ISKRA a 2	-----	QS a 4
42)	ISKRA a 4	-----	QS a 12



Alle Dioden 1N4148

Medibtec Obere Bergstraße 3 8501 Heerlosberg		NE-Karte: Schaltbild		Meditec
EIL 310 2426				
Tag	Name	Tag	Name	
1		1		
2		2		
3		3		
4		4		
5		5		
6		6		
7		7		
8		8		
9		9		
10		10		
11		11		
12		12		
13		13		
14		14		
15		15		
16		16		
17		17		
18		18		
19		19		
20		20		
21		21		
22		22		
23		23		
24		24		
25		25		
26		26		
27		27		
28		28		
29		29		
30		30		
31		31		
32		32		
33		33		
34		34		
35		35		
36		36		
37		37		
38		38		
39		39		
40		40		
41		41		
42		42		
43		43		
44		44		
45		45		
46		46		
47		47		
48		48		
49		49		
50		50		







Part No.	QTY	Part No.	QTY
SI 50112	1	SI 50113	1
SI 50114	1	SI 50115	1
SI 50116	1	SI 50117	1
SI 50118	1	SI 50119	1
SI 50120	1	SI 50121	1
SI 50122	1	SI 50123	1
SI 50124	1	SI 50125	1
SI 50126	1	SI 50127	1
SI 50128	1	SI 50129	1
SI 50130	1	SI 50131	1
SI 50132	1	SI 50133	1
SI 50134	1	SI 50135	1
SI 50136	1	SI 50137	1
SI 50138	1	SI 50139	1
SI 50140	1	SI 50141	1
SI 50142	1	SI 50143	1
SI 50144	1	SI 50145	1
SI 50146	1	SI 50147	1
SI 50148	1	SI 50149	1
SI 50150	1		

Part No.	QTY	Part No.	QTY
SI 50112	1	SI 50113	1
SI 50114	1	SI 50115	1
SI 50116	1	SI 50117	1
SI 50118	1	SI 50119	1
SI 50120	1	SI 50121	1
SI 50122	1	SI 50123	1
SI 50124	1	SI 50125	1
SI 50126	1	SI 50127	1
SI 50128	1	SI 50129	1
SI 50130	1	SI 50131	1
SI 50132	1	SI 50133	1
SI 50134	1	SI 50135	1
SI 50136	1	SI 50137	1
SI 50138	1	SI 50139	1
SI 50140	1	SI 50141	1
SI 50142	1	SI 50143	1
SI 50144	1	SI 50145	1
SI 50146	1	SI 50147	1
SI 50148	1	SI 50149	1
SI 50150	1		

Part No.	QTY	Part No.	QTY
SI 50112	1	SI 50113	1
SI 50114	1	SI 50115	1
SI 50116	1	SI 50117	1
SI 50118	1	SI 50119	1
SI 50120	1	SI 50121	1
SI 50122	1	SI 50123	1
SI 50124	1	SI 50125	1
SI 50126	1	SI 50127	1
SI 50128	1	SI 50129	1
SI 50130	1	SI 50131	1
SI 50132	1	SI 50133	1
SI 50134	1	SI 50135	1
SI 50136	1	SI 50137	1
SI 50138	1	SI 50139	1
SI 50140	1	SI 50141	1
SI 50142	1	SI 50143	1
SI 50144	1	SI 50145	1
SI 50146	1	SI 50147	1
SI 50148	1	SI 50149	1
SI 50150	1		

Part No.	QTY	Part No.	QTY
SI 50112	1	SI 50113	1
SI 50114	1	SI 50115	1
SI 50116	1	SI 50117	1
SI 50118	1	SI 50119	1
SI 50120	1	SI 50121	1
SI 50122	1	SI 50123	1
SI 50124	1	SI 50125	1
SI 50126	1	SI 50127	1
SI 50128	1	SI 50129	1
SI 50130	1	SI 50131	1
SI 50132	1	SI 50133	1
SI 50134	1	SI 50135	1
SI 50136	1	SI 50137	1
SI 50138	1	SI 50139	1
SI 50140	1	SI 50141	1
SI 50142	1	SI 50143	1
SI 50144	1	SI 50145	1
SI 50146	1	SI 50147	1
SI 50148	1	SI 50149	1
SI 50150	1		






Gradivo		Meh. obdelava		Toplotna obdelava		Odstal. del. dela		Povrhl. zaključa		NASTY	
Konstr.	MAHTIGALI, Ljubljana	Iskila	1	2							Koda
Projekt.	MAHTIGALI, Ljubljana	Znak	HR	NI							700 110 099
Prejela.	MERCINA, Ljubljana	St. obr.	L-100								Prilozitev
St. izved.	MERCINA, Ljubljana	Datum	19.2.85	12.8.85							700 110 000
Yehning		Podpis									Merilo
											Form
											ANV

# NAPAJALNIK ZA OKULISTIČNI LASER

YAG 02




Iskra  
Center za  
elektrooptiko  
Ljubljana


Trasno brailom celesno in podroben in sprejeto v redovno delo in izvedeno


11







		Freimaßtoleranzen		 Meditec Obere Bergstraße 3 8501 Heroldsberg		
		Tag	Name			QF-Karte Schaltbild
		Bearb.	13.5.86	<i>M...</i>	MQL 410 2972	
		Gepr.	20.7.86	<i>J...</i>		18
		Norm.				
1QL	10	Ausgabe	Änderung	Tca	Name	






MERZU  
 Layout B5 MQL 210 2486  
 Layout L6 MQL 210 2487  
 Layoutnachb. MQL 210 2489


MQL 210 2488/A	
QF-Karte	
Bauelementenliste	
1:1	
MQL 210 2488	

# Spallampentrafo


6V ~ 0V





		Freimadtoleranzen		 Meditec Obere Bergstraße 3 8501 Heroldsberg	Maßstab
		Tag	Name		
		Bearb.	11.7.86	QI-Karte Schaltbild	
		Gepr.	29.7.86		
		Norm.			
		MQL 410 3104			20
1QL 10	Aus- gabe	Änderung	Tag	Name	




HIERZU:  
 Layout Ls. MQL 210 3015  
 Lotstopplack MQL 210 3015







		Freimastlötlötungen		 <b>REINHARDT THYZEL</b> <small>Germany</small>		Obere Bergstraße 3 8501 Heroldsberg	
		Tag	Name			QL-Karte: Schaltbild	
		Bearb.	16.5.86 Meuser				
		Gepr.	27.7.86 i.e.				
		Norm.					
MQL 10		A SFH 202 ohne 76.04.11 AT 334		MQL 410 2983			
Änderung Tag Name							




HIERZU:

Layout Ls. MQL 410 2483 / B


Lotslopplock MQL 410 2485 / B



Sicht auf Bestückungsseite

		Freimaßtoleranzen		 Meditec Obere Bergstraße 3 8501 Heroldsberg	Maßstab /
		Tag	Name		
		Bearb.	5.6.86 <i>Mueller</i>	QT-Karte: Schaltbild	/
		Gepr.	19.7.86 <i>lue</i>		
		Norm.			
9 <i>Ansicht</i> <i>24.3.88</i> <i>Mueller</i>		MQL 410 3031		2.4	





HIERZU :  
 Layout Ls MQL 310 2436 / B  
 Lötlstopplack MQL 310 2438 / B