BD 900®

Spaltlampe

Lampe à fente

Slit lamp

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Mode d'emploi

Instruction manual

Vorwort

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Introduction

We would like to thank you for your decision to purchase this Haag-Streit product.

If the instructions in this manual are carefully followed we are confident that this product will give you reliable and trouble-free usage.

Only properly trained and authorized persons are to operate the **Slit Lamp BD 900®**.

We will be forced to disclaim any warranty and liability if the instrument is altered in any way or if routine maintenance is neglected or not carried out according to factory specifications.

Only original Haag-Streit spare parts are to be used to ensure continued reliability of the instrument.

Purpose of use

The **Slit Lamp BD 900®** is used, at room temperature, in the examination, diagnosis and documentation of the human eye. It is usually used by Ophthalmologists, Optometrists or Opticians in their consulting rooms, clinics, hospitals or teaching facilities.

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Contents

Summary operating instruction	page 7
1 Overview	
1.1 Headrest	9
1.2 Slit lamp BD 900 [®]	9
2 Operating safety	11
3 Description of the equipment	
3.1 Illumination unit	15
3.2 Stereomicroscope	15
3.3 Accessories	16
3.4 Power supply	19
3.5 Adjustable focusing target	21
3.6 Video equipment (option)	21
4 Operating the equipment	
4.1 Illumination unit	23
4.2 Adjusting the eye-pieces and the eyeguards	25
4.3 Possibilities for regulating the brightness	25
4.4 Focussing with the video	25
Appendix A Installation of the equipment	
A.1 Attaching the headrest and the power supply to the instrument table HSM-801 or -901	27
A.2 Attaching the headrest and the power supply to tables and units of other manufacturers	27
A.3 Installation of background illumination	29
A.4 Installation of camera cables	33
A.5 Installation of video camera with C-mount interface	33
A.6 Handling the slit lamp when fitted with video equipment	33
A.7 Installation of video camera with Lipstick interface	35
A.8 Mounting of headrest and power supply	35
Appendix B Equipment maintenance routines	
B.1 Replacement of the halogen bulb	37
B.2 Replacement of the mirror	39
B.3 Cleaning of the gliding plate and the rails	39
B.4 Cleaning of the optical parts	39
B.5 Cleaning of the axle	39
B.6 Plastic dust cover	39
Appendix C Technical specifications	41

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
Summary operating instructions

1. Move the chinrest (6) up or down with the adjustment screw (7) until the eyes of the patient are level with the black mark on the headrest column (4).
2. Tighten the fixing screw of the eye-piece (25). Focus the eye-pieces to suit your refraction by turning the knurled rings and then set to your interpupillary distance.
3. Turn on illumination with the green toggle-switch on the power supply.
4. The brightness can be progressively set and adjusted by means of the potentiometer (18).
5. Rotate the control lever (31) until the light beam is at eye level.
6. Holding the control lever (31), inclined towards yourself, move the instrument base until the slit appears to be approximately in focus on the cornea. This coarse setting is achieved with the naked eye. Fine adjustments are obtained by rotating and tilting the lever while observing the image through the microscope (24).
7. The slit width is adjusted to the left or right by means of the knob (17), which also adjusts the angle between the microscope and the illumination unit.
8. The slit image can be shown in the vertical, horizontal or slanted positions, as required, by turning the illumination unit (lock engagement at 90°; stop engagements at 0° and 180°). The number of degrees can be read-off on the scale (11).
9. The magnification of the microscope in the Slit Lamp BD 900® is altered by changing the lens with the lever (26) or by exchanging the eye-pieces.
10. For general observations over a wide field the grey filter is placed in position by means of lever (9), and the slit is opened wide.
11. The video output facility (14) can be used for documentation.
12. The contrast enhancing filter, yellow (22) can be engaged to obtain the fluo-image display.

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
1 Overview

1.1 Headrest

- 1 Lamp cable
- 2 Headrest
- 3 Forehead band
- 4 Headrest height marker (patient's eye)
- 5 Adjustable focusing target
- 6 Chinrest
- 7 Chinrest height adjustment

1.2 Slit Lamp BD 900®

- 8 Lamp cover
- 9 Lever for filter grey and redfree
- 10 Lever for slit length, blue filter
- 11 Scale for angled position of the slit image
- 12 Interchangeable illumination mirror
- 13 Diffusor
- 14 Video connection (C-mount)
- 15 Protective cover
- 16 Illumination unit / microscope angle scale
- 17 Slit width control
- 18 Brightness control
- 19 Plug for power supply
- 20 Registered trade mark BD 900®
- 21 Slit diaphragm scale and blue filter
- 22 Yellow filter
- 23 Cap for accessories base
- 24 Stereo microscope with eye-pieces
- 25 Fixing screw of the eye-piece
- 26 Lens-changing lever
- 27 Knurled screw for fixation of
- 28 Breath shield
- 29 Joy stick base locking screw
- 30 Rail covers
- 31 Control lever
- 32 Gliding plate

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
2 Safety

Ambient conditions

• Transportation	Temperature	-40 °C to +70 °C
	Air pressure	500 hPa to 1060 hPa
	Relative humidity	10% to 95%
• Storage	Temperature	-10 °C to +55 °C
	Air pressure	700 hPa to 1060 hPa
	Relative humidity	10% to 95%
• Working	Temperature	+10 °C to +35 °C
	Air pressure	800 hPa to 1060 hPa
	Relative humidity	30% to 75%

Installation

- To avoid condensation, allow the instrument to adjust to room temperature for several hours before unpacking.
- Always mount the slit lamp and the headrest on an electric insulated and fire resistant table top.
- The rail covers (30) prevent from a tip over of the slit lamp.
- Check: do the connection pieces of the accessories fit (like screwed connections and bayonet fitting)?

Operation and surrounding

- Only qualified and trained personnel should operate the equipment, the training is at the owner's responsibility.
- Use only Haag-Streit accessories.
- Shut down after every use. In case the dust cover is used: risk of overheating.
- Never use the ocular to look at the sun.
- Never operate the instrument in the same room with combustible gases, volatile solvents (alcohol, benzol) or flammable anesthetic agents.
- Avoid humidity.

Changing the light bulbs

- Switch off the main switch - disconnect the mains connector.
- Allow the bulb to cool down.
- Do not open the slit lamp for any other reason than to change the light bulb.

Electrical

- Only a hospital grade 3-conductor electrical power supply cable must be used. *For USA and Canada: Detachable Power Supply Cord Set, UL Listed, type SJE, SJT or SJO, 3-conductor, not smaller than 18 AWG.*
- Plug, cable and ground lead connection of the socket have to be in perfect condition.

Cleaning

- Clean the slit lamp housing only with a slightly water dampened cloth.
- No liquids or corrosive agents.
- Clean only the exterior surfaces of the optical parts.

Light toxicity

- Long and intensive illumination could damage the retina. The brightness of illumination should be kept at a minimum, to view the structure of the target clearly. To avoid potential eye damage, the examination should not last longer than necessary.

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Instrument Table HSM-901 type M

Observe the following safety precautions

Höhenverstellung nur mit aufgesetztem Gerät betätigen
Do not use height adjustment without instrument in place
Ne pas utiliser le réglage en hauteur quand l'instrument n'est pas en place

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- The retinal dose for potential photo chemical damage depends on ray density and examination time. When the ray density is cut in half, the examination time doubles until the border value is reached.
- Until today, there is no proof established that slit lamp examination can endanger the eye. Nevertheless, we recommend to use the minimum amount of light possible for retinal examination.
- Children, Aphakics and people with eye problems are especially endangered. If the retina is exposed to the same or another light source in 24 hours, there is also a higher risk. This is especially true when a retinal picture is taken with flash.

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Warranty / product liability

- The instrument should be operated in accordance with the chapter 'Safety'. Incorrect operation can damage the instrument. Thus no warranty claims can be accepted.
- Continued use of an instrument that has been damaged by incorrect operation can lead to personal injury. The manufacturer cannot accept liability in this case.
- Repairs and alterations on this equipment should only be carried out by Haag-Streit service technicians or by authorized persons.

Statutory requirements

- The slit lamp BD 900® has been designed to conform with the IEC / EN 60 601-1 and ISO 10939 standards. Manufacturing procedures, testing, commissioning, maintenance and repair are conducted under the observance of Swiss and international regulations.
- When combining different medical and/or nonmedical electrical equipment standard IEC / EN 60 601-1 applies.
- The 'CE' marking confirms compliance of the slit lamp BD 900® with the directives 93/42/EEC and 89/336/EEC and the conformity module A.
- A copy of the declaration of conformity and information concerning EMC (Electromagnetic compatibility) of the present instrument can be requested at any time at Haag-Streit company.
- All statutory accident prevention regulations are to be observed.
- Classification

IEC / EN Standard 60 601-1 Slit Lamp BD 900® equipment safety class I

*Application part Type B
Operation mode: continuous operation*

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Environment

- Electric and electronic equipment must be separated from house waste! This equipment has been introduced into the market after the 13th of August 2005.
- Disposal via your local collecting point or your Haag-Streit distributor.
- Thus it is ensured, that no harmful materials get into the environment and that valuable raw materials can be used again.

a


b

c

d

Pictograms

- Read very carefully (in the instruction manual)!
- Warning for hot surfaces!
- Disconnect the electrical supply plug!
- Safety class II
- Manufacturer
- Serial number
- HS-Part Number
- Type B equipment
- Date of manufacture

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3.1 Illumination unit

Halogen bulb


The halogen bulb is centered ex-works to obtain an optimal illumination.

The light beam from the pre-centered halogen bulb passes through the condenser, the slit diaphragm and the illumination lens, before it is directed to the patient's eye by the oblique mirror.

The following can be mounted in the path of the light beam of the illumination unit

Filters operating by lever (34)

- 36 Open aperture
- 37 Redfree filter
- 38 Grey filter (10%)

36 37 38

Due to standard ISO 10939 the heat absorption filter and UV filter are permanently mounted!

Apertures at window (33) operating by lever (35)

- 39 Aperture producing area of 14 mm Ø
- 40 Wedge-shaped diaphragm for continuous slit length adjustment
- 41 Blue filter


39 40 41

3.2 Stereomicroscope

- 42 Two objective lenses
- 43 Lever for yellow filter
- 44 Cap for accessories base
- 45 Serial number of microscope
- 46 Eye-piece 10x or 25x
- 47 Connection for the video camera (C-mount)
- 48 Lens-changing lever
- 49 Eye-piece fixing screw

Lens	1x	1.6x	1x	1.6x
Eye-piece	10x	10x	25x	25x
Total magnification	10x	16x	25x	40x
Object image section Ø mm	18	11.3	8	5

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(model type designation: LC-SLT)


CAUTION

The power supply is fitted with a voltage selection switch (50) 115 V / 230 V. This must be adjusted to correspond with the local electric power supply before first use.

The main switch (52) will illuminate when the appliance is switched on.

Possible fault causes

1. The main switch is on, but fails to illuminate:
 - primary fuse defective
2. The main switch is on and illuminated, but the slit lamp fails to illuminate:
 - secondary fuse or halogene bulb defective


To change the primary voltage fuse

Open the fuse box (53) with the point of a screwdriver and pull out. Change the defective fuse cartridge. Push in the fuse box until it is heard to engage.

Fuse 2 x T 800 mA

HS-Part Number 1001330

To change the secondary voltage fuse

Turn the slit of the fuse holder (51) to the left with a screwdriver. Change the fuse cartridge.

Fuse F 3.15 A

HS-Part Number 1001331

The degree of brightness can be regulated with the brightness control (18).

The illumination intensity as per Appendix C - Technical specifications.

19


3.5 Adjustable focusing target

If the lamp cable (54) is connected, the fixation lamp burns as soon as the power supply switch is activated, then the focusing target (59) is visible.

The fixation target can be brought into any position by moving the bearer's arm (55). This permits uninterrupted observation through the microscope and simultaneous guidance of the patient's eye to any desired position. This is especially advantageous in examining the fundus.

It's a condition, that the patient's eye is correctly positioned at the headrest's height marker (58).

Ametropia of the patient from +10 D to -15 D (56) can be corrected (57). Steady fixation of the eye is maintained as accommodation and convergence are avoided.

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3.6 Video equipment (optional)

Dimension of object area

Magnification microscope	Sensor video camera 1/3 "	1/2 "
1x	16 x 12	21.3 x 16
1.6x	10 x 7.5	13.3 x 10

The video connection consists of the following components:

- Either **Video camera with C-mount**

Maximum dimensions of camera

depth 110 mm
height 50 mm
width 55 mm


- Or **Lipstick video camera**

needs an adapter from the relevant video camera manufacturer so that it can be mounted to the C-mount interface of the BD 900®.

- **Cable lead guides**

consisting of

- Cable junction box for camera cable
- Cable guide for optical light lead
- Rail cover cable lead guide

Maximale Grösse

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Maximum dimensions

Moreover Haag-Streit recommends:

- **Background illumination with light source**

consisting of

- Cold light source
- Optical light lead
- Background illumination EcoLite EL 02

21

CAUTION

The following components must comply with
EN 60601-1-1: Power supply, monitor, video
recorder and video printer.

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4 Operating the equipment

CAUTION

It is imperative to read the chapter "Operating safety" before operating the equipment and to observe its precautions.

- The height of the table must be selected to enable the patient to obtain a firm rest for the forehead and chin.
- In order to ensure that only that part of the eye to be examined is illuminated and to avoid disturbing glare, the height of the light beam should be appropriately adjusted.

NOTICE

- Only use as much light intensity as required.
- The light intensity should be reduced to a minimum on completion of each examination. This extends the life of the light bulb and protects the patient from severe glare at the beginning of an examination.
- Those parts of the equipment coming into contact with patients should be cleaned with a dry cloth before each examination.

4.1 Illumination unit

Slit width

The slit width is set with either of the control knobs (63). The division on the left hand knob enables any given slit width setting to be subsequently recovered; but provides no information on the actual width of the slit light.


Slit length

The slit length is set with the lever (61).


Slit positioning

The slit image can be positioned by turning the upper part of the illumination unit (60) from the vertical, through a range of any given oblique positions, to the horizontal.

A scale (62) and (64) indicates the angle of the slit.

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4.2 Adjusting the eye-pieces and the eyeguards


Eye-pieces

Before the first examination, the eye-pieces should be adjusted singly to correspond with the refraction divergence of the examiner.

Insert the provided focusing rod (65) in place of the protecting cover (66) with the flat projection surface facing the microscope. Line up the illumination unit and the microscope in the central position (0°). Set the eye-piece with the fixing screw until it no longer continues to turn when adjusting the diopter settings. Adjust each eye-piece separately, by turning the knurled ring with the diopter scale separately, until the projected slit is seen in focus.

Adjust from (+) to (-) at low magnification.

Younger examiners are recommended to make an additional correction on both eye-pieces of -1 to -2 diopters each, to compensate for the accommodation that occurs as a result of convergence when conducting binocular examinations.

Eyeguard

The adjustable eyeguard (71) serves to set up the correct working distance of the examiner to the eye-piece.

For examiners **without** spectacles:

pull out the eyeguard fully to the extended position.

For examiners **wearing** spectacles:

push in the eyeguard fully to the retracted position.

4.3 Possibilities for regulating the brightness

The brightness of the slit light illumination can be varied by means of

- the knob on the instrument base (18) which adjusts the intensity of the slit illumination
- the 10% grey filter on the upper part of the illumination unit of the slit lamp

Regulation of the brightness for the video

The brightness of the background illumination (optional) on the cold light source housing can be varied by means of the knob for a continuous regulation of the brightness.

4.4 Focussing with the video

Focussing onto an object during an examination should always be carried out as a rule by moving the microscope towards the object until the area of interest comes into focus.

If the microscope is moved further towards the object, it may continue to be physically observed due to the users accommodation, but no further compatibility will exist with the video camera.

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Eye-piece with prism housing

- 67 Prism housing
- 68 Eye-piece mounting
- 69 Index (white spot)
- 70 Knurled ring with diopter scale for setting the refraction of the examiner (± 6 dpt)
- 71 Sliding eyeguard
- 72 Eye-piece fixing screw

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Appendix A Installation of the equipment

A.1 Attaching the headrest and the power supply to the instrument table HSM-801 or HSM-901

The operating instructions of the concerning instrument table should be observed!


A.2 Attaching the headrest and the power supply to tables and units of other manufacturers

1. Fix the power supply with 4 screws (73).
2. Remove the knurled screws (70) and the cable plate (68).
3. Fix the ground lead (67) with the screw (69) to the unit.
Tables of other manufacturers: connect the ground lead only if necessary!
4. Connect the four pole electrical plug connector (77) between the headrest and the power supply.
5. Install the headrest in the required position on the adapter piece of the table top by means of both Allan screws (82).
6. The cable plate (74) should be installed in such a manner that the plug connection and cable are affixed.
7. Mount the handles (83) and connect the main lamp lead (81) to the lamp housing.
8. Screw the angled plug (80) up tight to the plug socket (84).
9. Attach the potentiometer cable with the cable junction box (85) on the right side of the rail cover by using the center groove.
10. Connect the power cable (78).

CAUTION

- Only qualified and trained personnel may work on cables and components being connected to the power supply.
- To prevent short circuits, do not crimp or squeeze the electric wires!

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A.3 Installation of background illumination

a) It is possible to take images and video using only the inbuilt **illumination** of the slit lamp.


Streuscheibe
Diffuseur
Diffusor
1004744

b) Diffusor

Results can be improved by use of an optional diffusor, but this excludes simultaneous slit illumination.

Installation

Insert the diffusor from below over the rounded feature into the foreseen reception slots (82) where it will snap in to secure.

Haag-Streit recommends, for a good illumination of the anterior part of the eye, use an optional cold light source and background illumination through the optional holder for fiberoptic light cable.

c) Optical light lead

The variable output cold light source with light lead holder is the optimal solution for video images.

Installation

The cold light source should be positioned according to the instrument table.


Slacken the set screw (86) until the light lead adapter (87) can be removed downwards.

Insert the light lead (90) for the background illumination into the female receiver piece boring of the light lead adapter (87) and secure slightly by means of the set screw (89).

The light lead adapter (87) can now be re-inserted in the foreseen slot under the mirror. The set screw (86) can then again be re-inserted and gently tightened up.

Insert the diffusor from below over the rounded feature into the foreseen reception slots (88) where it will snap in to secure.

The light lead is fixed with the cable guides on the illumination arm (91) and on the rail cover (92), and then passed through the cable cover (93) to the cold light source.

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
EcoLite EL 02

d) Background illumination EcoLite EL 02

The background illumination EcoLite EL 02 allows greater flexibility simply by pivoting it.

Installation


- Remove protection cover (95)
- Install the EcoLite EL 02 (94)
- Insert the angled plug (96) into EcoLite EL 02
- Connect the power supply (97)
- Turn ON / OFF the light by moving up and down the switch (98)
- Pivoting
- Filter can be removed

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
A.4 Installation of camera cables

Secure the cable to the camera (100), lead it through the table (102) and connect it to the camera control unit. On the slit lamp the cable is fixed on the microscope arm (99) and additionally on the rail cover with the cable guide (101).

A.5 Installation of video camera with C-mount interface

Unscrew the cover (103) on the C-mount interface. First screw the adapter (104) on the camera and then affix it with 3 set screws (105).

Image position on the monitor

Note the position of the camera! The camera must be oriented parallel to the microscope arm and must be mounted in the way, that the picture is correctly displayed on the screen.


In the case of misaligned images on the screen, the camera should be rotated until the correct position image is displayed.

Image focus on the monitor

Needs not to be adjusted, the interface has been designed so that optimal image focus is set automatically!

Breath shield


Depending on camera model, a distance screw must be built in (106), thus the breath shield can be fixed correctly.

A.6 Handling the slit lamp when fitted with video equipment

To achieve the same focus through both eye-pieces and monitor, the eye-pieces must be exactly adjusted with the focus test rod.

When coupling the accommodation facility with convergent viewing, the normal refraction often has to be corrected by -1 to -2 dpt on the eye-piece.

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
A.7 Installation of video camera with Lipstick interface

Installation of video camera head

Affix the adapter for Lipstick to C-mount. Only tighten the 3 set screws sufficiently enough to enable the camera adapter still to be turned.

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Image position on the monitor

Turn the camera head until the image on the monitor is correctly margin-adjusted.



Image focus on the monitor

Use the focus test rod. Focus a slit of 8 mm height and 4 mm width onto the focus test rod. Shift the illumination to 45° from the microscope.

Set a 16x magnification. Adjust the video camera depth. The area, in which the structure of the focus rod surface is seen sharply in focus, should lie in the centre of the light slit (107), whereas the structure to the left and right is shown out of focus (108).

The 3 set screws are still slightly loose and the camera height can be adjusted vertically by turning the camera adapter - but not the camera which can be held with one hand so that it cannot turn. When the image on the video monitor screen appears sharply in focus, the lock nut of the camera can be firmly tightened.

The position of the image on the monitor as to margin correctness should be checked again and correctly adjusted. The camera adapter can then be secured to the microscope housing by tightening the 3 set screws.

CAUTION

Only switch on the camera control unit when all components have been correctly installed.
Observe the instruction manual for the video camera when making these connections.

A.8 Mounting of headrest and power supply

The corresponding Haag-Streit operating instructions should be observed!

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Appendix B Equipment maintenance routines

B.1 Replacement of the halogen bulb

**CAUTION**

1. Before changing the lamp:
turn off the switch on the power supply and disconnect the power cord!

**CAUTION**




Increased temperatures can occur on the upper surface of the lamp housing! Allow the light bulb to cool down before removing it or use protective gloves to remove.
Use only original HS halogen bulbs!

2. Hold the lamp bulb cover at the sides and pull upwards.
3. Unplug the contact base (112).
4. Unhook the retaining spring (110) on the left side and swing out sideways.

**CAUTION**

The lamp bulb should not be touched with the fingers!

5. Replace the halogen lamp bulb (111). The groove on the lamp holder (109) should point upwards.
6. Hook up the retaining spring, plug in the contact base, push on the lamp bulb cover until it engages and connect the lamp cable.

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HS-Part number 1001434

B.2 Replacement of the mirror

The mirror need not be removed for cleaning. It should only be replaced when mechanically damaged.

CAUTION


Do not touch the mirror surface with the fingers!

Dismantling

- Depress the mirror in the direction of the arrow. Care should be taken to ensure that it does not jam.

Assembly

- Depress the spring retaining plate downwards
- Slide the mirror half way in with the other hand
- Release the retaining plate
- Continue to slide the mirror into position until it engages

B.3 Cleaning of the gliding plate and the rails

If it becomes difficult to move the slit lamp around on the gliding plate, the plate should be cleaned with a slightly oily rag.


After removing the rail covers (30), the rails should be cleaned with a stiff brush.

B.4 Cleaning of the optical parts

A dust brush is to be employed for dusting the exposed glass surfaces.

Surfaces which have become really dirty should first be dusted, then wiped carefully with a soft dry clean cloth, washed linen, chamois leather, or some such material which will not scratch the surface.

Brush for cleaning the optics
HS-Part number 1001398

39

B.5 Cleaning of the axle

Clean only with dry, lint-free cloths!

B.6 Plastic dust cover

It is recommended to protect the slit lamp with a dust cover in case it is not in use.

Plastic dust cover, small (for slit lamp)
HS-Part number 1001395
Plastic dust cover, large (for several instruments)
HS-Part number 1001434

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Appendix C Technical specifications

Slit illumination

Slit image width	0 - 14 mm continuous
Slit image length	1 - 14 mm continuous
illumination field diameter	14 mm
Slit image radial range	0° - 90°
Radial movement of the slit illumination relative to the microscope axis	Horizontal ± 90°
Filter	Blue, grey (10%) and red. The UV filter and the heat absorption filter are permanently mounted.
Light source	Halogen lamp bulb 12 V / 2.5 A
Illumination intensity	0 V / 0 Lux at position 0 7 V / 75 000 Lux at position 1/4 9 V / 150 000 Lux at position 1/2 11.5 V / 300 000 Lux at position 1

Irradiance depending on light wave length
see **diagram** at right hand side

Irradiance at max. intensity (11.5 V lamp) and max. aperture*)
0.01 mW / cm² (305 nm - 400 nm)
90 mW / cm² (380 nm - 700 nm)
5 mW / cm² (700 nm - 1100 nm)

Spectrally-weighted photochemical aphakic output source radiance L_A at maximum intensity (11.5 V lamp) and maximum aperture: 270 mW / (cm² sr).*)

Spectrally-weighted photochemical phakic output source radiance L_B at maximum intensity (11.5 V lamp) and maximum aperture: 230 mW / (cm² sr).*)


Spectrally-weighted photochemical radiance L_A and L_B give a measure of the potential hazard that exists for a beam of light to cause photochemical damage to the retina. L_A gives the measure either for eyes in which the crystalline lens has been removed (aphakes) and has not been replaced by a UV-blocking lens or for eyes of very young children. L_B gives the measure for eyes in which a crystalline lens is in place. Values of L_A and L_B over 80 mW / (cm² sr) are considered high for beams which completely fill a dilated pupil. *)

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Diagram Irradiance at positions 1, 1/2 and 1/4

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Stereo-microscope

Stereo angle	13°
Magnification changer	1x and 1.6x
Magnification	
with ocular 10x (standard)	10x and 16x
with ocular 25x	25x and 40x
Range of adjusting eye-pieces	+ 6 to -6 diopters
Inter pupillary distance	54 - 94 mm
Interface video camera	C-mount

Instrument base

Operation	Single handed 3-dimensional opera- tion of the guide lever
Adjustment of the instrument base	80 mm (length) 30 mm (height) 100 mm (side)
Electric supply	Class II, type B 230 V~, 0.13 - 0.21 A or 115 V~, 0.25 - 0.42 A 50 - 60 Hz

Miscellaneous

Weight	8.1 kg (without power supply, headrest and options)
--------	---

Classification

Against electric shock	Class I, type B
Against water	IPX0
Sterilisation and disinfection	
• for the slit lamp	not necessary
• for the headrest	use chin paper and clean the forehead with cotton soaked in alcohol
Operation mode	continuous

Certified Quality Management System

Power supply

Input	AC 115V, ±10% / 0.25 - 0.42A AC 230V, ±10% / 0.13 - 0.21A 50 - 60Hz
Output 1	DC 7 - 14.4V / max. 36W
Output 2	DC 5V / 1W

Classification	<ul style="list-style-type: none"> • Class II • IPX0 • continuous operation • equipment not suitable for use in the presence of a flammable anaesthetic mixture
----------------	--

CLASSIFIED BY UNDERWRITERS LA-
BORATORIES INC. WITH RESPECT
TO ELECTRIC SHOCK, FIRE AND
MECHANICAL HAZARDS ONLY IN
ACCORDANCE WITH UL 60 601-1

Haag-Streit AG
Gartenstadtstrasse 10
CH-3098 Koeniz, Switzerland
Phone ++ 41 31 978 01 11
Fax ++ 41 31 978 02 82
eMail info@haag-streit.ch

Haag-Streit AG
Gartenstadtstrasse 10
CH-3098 Koeniz, Switzerland
Phone ++ 41 31 978 01 11
Fax ++ 41 31 978 02 82
eMail info@haag-streit.ch
internet www.haag-streit.com

01.09-.5

1500.7200631.04040