BP 900®

Spaltlampe

Lampe à fente

Slit lamp

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Mode d'emploi

Instruction manual

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Informations supplémentaires

Additional informations

Vorwort

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Introduction

We would like to thank you for your decision to purchase this Haag-Streit product.

If the instructions in this manual are carefully followed we are confident that this product will give you reliable and trouble-free usage.

Purpose of use

The **BP 900®** is an AC-powered slitlamp biomicroscope intended for use in eye examination of the anterior eye segment, from the cornea epithelium to the posterior capsule. It is used to aid in the diagnosis of diseases or trauma which affect the structural properties of the anterior eye segment.

Contents

1 Safety	5
2 Overview	
2.1 Headrest	8
2.2 Slit lamp.....	8
3 Operating	
3.1 Adjust the eyepieces	11
3.2 Prepare the patient	11
3.3 Operating the instrument.....	11
4 Equipment maintenance routines	
4.1 Replacement of the bulb	15
4.2 Check electrical contacts of bulb and lamp housing	15
4.3 Utilisation du moyen de contact.....	15
5 Technical specifications	15
Additional informations	21

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1 Safety

Ambient conditions

• Transportation	Temperature	-40 °C to +70 °C
	Air pressure	500 hPa to 1060 hPa
	Relative humidity	10% to 95%
• Storage	Temperature	-10 °C to +55 °C
	Air pressure	700 hPa to 1060 hPa
	Relative humidity	10% to 95%
• Working	Temperature	+10 °C to +35 °C
	Air pressure	800 hPa to 1060 hPa
	Relative humidity	30% to 75%

Installation

- To avoid condensation, allow the instrument to adjust to room temperature for several hours before unpacking.
- Always mount the slit lamp and the headrest on an electric insulated and fire resistant table top.
- The rail covers (33) prevent from a tip over of the slit lamp.
- Check: do the connection pieces of the accessories fit tightly (like screwed connections and quick fit)?

Operation and surrounding

- Only qualified and trained personnel should operate the equipment, the training is at the owner's responsibility.
- Use only Haag-Streit accessories.
- Shut down after every use. In case the dust cover is used: risk of overheating.
- Never use the ocular to look at the sun.
- Never operate the instrument in the same room with combustible gases, volatile solvents (alcohol, benzol) or flammable anaesthetic agents.
- Avoid humidity.

Changing the light bulbs (not for LI 900®)

- Switch off the main switch - disconnect the mains connector.
- Allow the bulb to cool down.
- Do not open the slit lamp for any other reason than to change the light bulb.

Electrical


- Only a hospital grade 3-conductor electrical power supply cable must be used. *For USA and Canada: Detachable Power Supply Cord Set, UL Listed, type SJE, SJT or SJO, 3-conductor, not smaller than 18 AWG.*
- Plug, cable and ground lead connection of the socket have to be in perfect condition.

Cleaning

- Clean the slit lamp housing only with a slightly water dampened cloth.
- No liquids, no alcohol or corrosive agents.
- Clean only the exterior surfaces of the optical parts.
- Use chin paper and clean the forehead band with a cotton pad soaked in alcohol.

Light toxicity

- Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures. This device should be used with filters that eliminate UV radiation (< 400 nm) and, whenever possible, filters that eliminate short-wavelength blue light (< 420 nm).
- The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

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Accessory on the auxiliary power

- If the 5V auxiliary power (20) is not used, the socket must be covered with the cover cap.
- If an accessory is connected, check plug!

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- While no acute optical radiation hazards have been identified for slit lamps, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. Infants, aphakes and persons with diseased eyes will be at greater risk.

- The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

- See caution notes on page 17

Warranty / product liability

- The instrument should be operated in accordance with the chapter 'Safety'. Incorrect operation can damage the instrument. Thus no warranty claims can be accepted.

- Continued use of an instrument that has been damaged by incorrect operation can lead to personal injury. The manufacturer cannot accept liability in this case.

- Repairs and alterations on this equipment should only be carried out by Haag-Streit service technicians or by authorized persons.

Statutory requirements

- The slit lamp BP 900® has been designed and constructed to conform with the IEC / EN 60 601-1-2 and EN ISO 15004-2 standards. Manufacturing procedures, testing, commissioning, maintenance and repair are conducted under the observance of Swiss and international regulations.

- When combining different medical and/or nonmedical electrical equipment standard IEC / EN 60 601-1 applies.

- The 'CE' marking confirms compliance of the slit lamp BP 900® with the directive 93/42/EEC.

- The slit lamp BP 900® meets the electromagnetic compatibility requirements of IEC / EN 60601-1-2. The appliance is designed to keep the generation and emission of electromagnetic interference to a level that does not affect the normal operation of other items of equipment. The appliance is likewise designed to provide a certain level of resistance to outside sources of electromagnetic interference. • All statutory accident prevention regulations are to be observed.

Classification

Standard IEC / EN 60 601-1 Slit Lamp BP 900® equipment safety class I.

Application part Type B.

Operation mode: continuous operation

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Environment

- Electric and electronic equipment must be separated from house waste! This equipment has been introduced into the market after the 13th of August 2005.

- Disposal via your local collecting point or your Haag-Streit distributor.

- Thus it is ensured, that no harmful materials get into the environment and that valuable raw materials can be used again.

a


b

c

Pictograms

- a Read very carefully (in the instruction manual)!
- b Warning for hot surfaces!
- c Disconnect the electrical supply plug!
- d Manufacturer
- e Serial number
- f HS-Part Number
- g Application part Type B
- h Date of manufacture

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
2 Overview

- 2.1 Headrest**
- 1 Lamp cable
 - 2 Headrest
 - 3 Forehead band
 - 4 Headrest height marker (patient's eye)
 - 5 Adjustable focusing target
 - 6 Chin rest
 - 7 Chin rest height adjustment
- 2.2 Slit Lamp BP 900®**
- 8 Lamp cover tungstene and halogen (LED LI 900® see separate instruction manual)
 - 9 Lever for grey - and redfree filter
 - 10 Scale for angled position of the slit image
 - 11 Illumination mirror
 - 12 Diffusor
 - 13 Magnification changer
 - 14 Mounting screw for the stereo microscope
 - 15 Protective cover
 - 16 Illumination unit / microscope angle scale
 - 17 Illumination arm locking screw
 - 18 Microscope arm locking screw
 - 19 Slit width controls
 - 20 5V auxiliary power (eg. EcoLite)
 - 21 Interface power supply PS-B
 - 22 Registered trade mark BP 900®
 - 23 Slit length scale / diaphragms
 - 24 Slit length, slit rotation, blue filter and fixation star control
 - 25 Cap screw (accessories pin)
 - 26 Quick fit connector accessory
 - 27 Stereo microscope with eyepieces
 - 28 Mounting hole for tonometer AT 900® Mod. BP
 - 29 Mounting screw for breath shield
 - 30 Centering screw
 - 31 5° to 20° inclination angle latch
 - 32 Breath shield
 - 33 Rail covers
 - 34 Axle
 - 35 Screw to block horizontal movement
 - 36 Push-button for slit illumination level control
 - 37 Control lever
 - 38 Gliding plate

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3 Operating

3.1 Adjusting the eyepieces

Before the first examination, the eyepieces should be adjusted singly to correspond with the refraction divergence of the examiner. Insert the provided focusing rod (39) in place of the protective cover (40) with the flat projection surface facing the microscope. Line up the illumination unit and the microscope in the central position (0°).

Adjust each eyepiece separately, by turning the knurled ring with the diopter scale separately, until the projected slit is seen in focus. Adjust from (+) to (-) at low magnification.

Younger examiners are recommended to make an additional correction on both eyepieces of - 1 to - 2 diopters each, to compensate for the accommodation that occurs as a result of convergence when conducting binocular examinations. The requisite pupil distance on the microscope can then be adjusted.

The adjustable **eyecup** (41) serves to set up the correct working distance of the examiner to the eyepiece.

For examiners **without spectacles**: pull out the eyecup fully to the extended position.

For examiners **wearing spectacles**: push in the eyecup fully to the retracted position.


3.2 Prepare the patient

- The height of the table must be selected to enable the patient to obtain a firm rest for the forehead and chin.
- In order to ensure that only that part of the eye to be examined is illuminated and to avoid disturbing glare, the height of the light beam should be appropriately adjusted.
- Those parts of the equipment coming into contact with patients should be cleaned with a dry cloth before each examination.
- The light intensity should be reduced to a minimum on completion of each examination. This extends the life of the light bulb, slit lamp and protects the patient from severe glare at the beginning of an examination (not for LI 900®).

3.3 Operating the instrument

1. Move the chin rest (43) up or down with the adjustment screw (44) until the **eyes** of the patient are **level** with the black mark on the headrest column (42).

2. Focus the eyepieces (45) to suit **your refraction** by turning the knurled rings and then set to your **interpupillary distance**.

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3. Switch on the illumination unit by activating the power supply switch. Control level of slit **illumination** by pushing button (46) on control lever (56)

- 1x click: 1 position brighter
- 2x click: 1 position less bright
- Press and hold down the button: goes to maximum brightness and back after release.

LI 900® see separate instruction manual

4. Rotate the control lever (56) until the **light beam** is at **eye level**.

5. Microscope **magnification** is changed by rotating the magnification changer (52).

6. Holding the control lever (56), inclined towards yourself, move the instrument base until the **slit** appears to be approximately in focus on the cornea. This coarse setting is achieved with the naked eye. Fine adjustments are obtained by rotating and tilting the lever while observing the image through the microscope (50).

7. The **slit width** is adjusted to the left or right by means of the knob (55), which also adjusts the **angle** between the microscope and the illumination unit.

8. The **slit image** is made vertical, horizontal, or given any desired **angle** by means of control knob (49) (notches at 45°, 90° and 135°; stops at 0° and 180°, markers in 10° steps).

9. In order to ensure the possibility of conducting a **binocular fundus examination** in optical cross section without difficulty from side angles between 3 and 10°, the shorter mirror (51) should be installed, the illumination unit is rotated through 90° with the knurled control knob (49) and with the latch (54) in stages of 5°, the illumination unit and the microscope are turned to the center position (0°).


10. Loosening the centering screw (53) will allow the slit image to be moved off center for **scleral illumination**. Tightening the screw will recenter the slit image in the visual field of the microscope.

11. For **general observations** over a wide field the grey filter is placed in position by means of lever (48), and the slit is opened wide.

12. For examination of the **posterior of the eye**, preset lenses or diagnostic contact lenses have to be used.

Filters (48) / apertures (47) operating (49)

- 57 Open aperture
- 58 Grey filter (10%)
- 59 Redfree filter
- 60 Empty space for accommodating an add. filter
- 61 Fixation star (*particularly suitable for examining the fixation in amblyopia*)
- 62 Apertures producing areas of 8, 5, 3, 2, 1 and 0.2 mm Ø
- 63 Indication of wedge-shaped diaphragm for continuous slit length adjustment in mm
- 64 Blue filter

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4 Equipment maintenance routines

The slit lamp BP 900® in combination with the illumination set LI 900® may be used without any maintenance during the entire life cycle.

4.1 Replacement of the bulb (not for LI 900®)

The lamp cover can be removed by depressing both pressure buttons (65) to expose the bulb.

CAUTION



Before changing the bulb, turn off the switch on the power supply and disconnect the power cord!

CAUTION

Increased temperatures can occur on the upper surface of the lamp housing! Allow the light bulb to cool down before removing it or use protective gloves to remove. Use only original HS (66) lamp bulbs!

When installing a new bulb, care should be taken to place the notch in the centering base into the positioning catch. Replace the lamp cover and push it down until the two pressure buttons engage securely.

4.2 Check electrical contacts of bulb and lamp housing

With suspicion of bad electrical contacts (flickering or loss of light):

- Clean contact surface (with sandpaper, screwdriver etc.) of bulb (67), the contact springs in the lamp housing (69) and (70)
- If necessary bend slightly the contact surfaces in the lamp housing (69) and (70) in direction of the bulb
- Bend slightly both of the retaining clips (68) outwards


4.3 Use of contact medium on the electrical contacts of the bulb and lamp housing

At each change of light bulb, lubricate the contact points (A) and (B) with enough contact medium. The contacts of the bulb and the lid of the lamp should not show any scorch marks.

- Fit the lamp cover and remove it to check if there is enough contact medium on the points of contact (C) and (D). Reapply some medium if needed.
- The product does not require labeling. (Perfluorocarbons, Polyetheroil PTFE).

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5 Technical specifications

Slit illumination

Slit image width	0 - 8 mm continuous
Slit image length	1 - 8 mm continuous
Illumination fields	Circle ø: 8 / 5 / 3 / 2 / 1 / 0.2 mm Test mark with fixation star
Slit image radial range	± 90°
Radial movement of the slit illumination relative to the microscope axis	Horizontal ± 90°, vertical 0 - 20°
Filters	Blue, redfree (green), grey (10%). The UV filter and the heat absorption filter are permanently mounted.
Light source	Tungsten / halogen incandescent light bulb 6 V, 4.5 A / 7.5 V, 38 W
Illumination intensity	0 V 0 Lux no power 5 V ~150 000 Lux level 1/2 6 V ~300 000 Lux level 1 7.5 V ~600 000 Lux level 2

LED illumination LI 900®, see separate instruction manual


Slit illumination tungsten

Maximum examination times according to ISO 15004-2 and ISO 10939

CAUTION – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum output will exceed the safety guideline after 3.2 minutes with slit illumination or 2.5 minutes with both slit and peripheral illumination.

Spektrum Wolfram Beleuchtung / Spectre éclairage tungstène / Spectrum tungsten illumination

X= Wellenlänge / Longueur d'onde / Wavelength [nm]
Y= Rel. Intensität / Intensité relative / Rel. intensity


Slit illumination tungsten

Maximum examination times according to ISO 15004-2 and ISO 10939

CAUTION – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum output will exceed the safety guideline after 5.5 minutes with slit illumination or 3.5 minutes with both slit and peripheral illumination.

Spektrum Halogen Beleuchtung / Spectre éclairage halogène / Spectrum halogen illumination

X= Wellenlänge / Longueur d'onde / Wavelength [nm]
Y= Rel. Intensität / Intensité relative / Rel. intensity

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Stereo-microscope

Stereo angle	13°
Magnification changer	10x / 16x / 25x
Ocular magnification	12.5x
Range of adjusting eyepieces	+ 8 to - 8 diopters
Inter pupillary distance	53 - 78 mm

Instrument base

Operation	Single handed 3-dimensional opera- tion of the control lever
Adjustment of the instrument base	80 mm (length) 90 mm (side) 30 mm (height)
Weight	11.7 kg (without power supply, headrest and options)

Classification

Against electric shock	Class I, type B
Against water	IPX0
Sterilisation and disinfection	
• for the slit lamp	not necessary
• for the headrest	use chin paper and clean the forehead band with a cotton pad soaked in alcohol
Operation mode	continuous

Certified Quality Management System

EN ISO 13485
Reg. Nr. 11956PRODUCTS CERTIFIED FOR BOTH THE U.S.
AND CANADIAN MARKETS, TO THE APPLI-
CABLE U.S. AND CANADIAN STANDARDS**Haag-Streit AG**Gartenstadtstrasse 10
3098 Koeniz, Switzerland
Phone ++ 41 31 978 01 11
Fax ++ 41 31 978 02 82
eMail info@haag-streit.com

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Additional informations

A Indications perating the equipment

A.1 Lamp bulb	21
A.2 Regulating the brightness	21
A.3 Microscope and eyepiece	23
A.4 Attaching an accessory in the beam path.....	23
A.5 Adjustable focusing target.....	23

B Mounting of the equipment

B.1 Microscope and illumination unit	25
B.2 Headrest and power supply.....	25
B.3 Digital photo adapter DP 01.....	27
B.4 Video adapter VP 01	29
B.5 Background illumination.....	31

C Equipment maintenance routines

C.1 Replacement of the mirror	33
C.2 Cleaning of the optical parts	33
C.3 Adjustment of friction of the slit width control	33
C.4 Cleaning of the gliding plate and the rails	33
C.5 Cleaning of the axle.....	33
C.6 Plastic dust cover.....	33


D Accessories	34
---------------------	----

A Indications instrument and operating

A.1 Lamp bulb (not for LI 900[®])

The tungsten / halogen lamp bulb is centered ex-works to obtain an optimal illumination. The 'H-S' marking on the centering base-cap indicates a genuine Haag-Streit spare part.

The light beam from the pre-centered bulb passes through the condenser, the slit diaphragm and the illumination lens, before it is directed to the patient's eye by the oblique mirror.

A.2 Regulating the brightness (not for LI 900[®])

The brightness of the slit light illumination can be varied by means of

- the three step push-button on the control lever
- the 10% grey filter on the upper part of the illumination unit of the slit lamp

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
A.3 Microscope and eyepiece

- 71 Front lens f 103 mm
- 72 Turning knob with display of selected magnification
- 73 Magnification changer 3 step
(Galilean system)
- 74 Knurled ring for quick fit
- 75 Binocular tube f 130 for convergent viewing
range of pupillary distance 53 - 78 mm
- 76 Eyepiece 12.5x
Diameter of field 16 mm
- 77 Index (white spot)
- 78 Knurled ring with diopter scale for setting the refraction
of the examiner (± 8 D)
- 79 Sliding eyecup (for spectacle wearers)

Total magnification	10x	16x	25x
Object image section ø in mm	20.0	12.7	8.0


A.4 Attaching an accessory in the beam path

Align the positioning points (80) on the upper side of the parts to be joined.

Close

Push together and tighten firmly the joint in direction of the arrow by turning it.


A.5 Adjustable focusing target

If the lamp cable (81) is connected, the fixation lamp burns as soon as the power supply switch is activated, then the focusing target (86) is visible.

The fixation target can be brought into any position by moving the bearer's arm (82). This permits uninterrupted observation through the microscope and simultaneous guidance of the patient's eye to any desired position. This is especially advantageous in examining the fundus.

It's a condition, that the patient's eye is correctly positioned at the headrest's height marker (85).

Ametropia of the patient from +10 D to -15 D (83) can be corrected (84). Steady fixation of the eye is maintained as accommodation and convergence are avoided.

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
B Installation of the equipment

B.1 Microscope and illumination unit


Place the stereo microscope on the saddle bearing of the microscope arm. Move it forward against the stop piece and secure with the screw (87) appearing from obliquely below (use the Allen key 5 mm).

Secure the breath shield (89) by means of the knurled screw (88) on the inner side of the bearer arm.


Before attaching the illumination unit to the column of the instrument base, the set screw (91) in the hub of the bearer arm is to be sufficiently loosened so that its tip no longer protrudes into the drill hole.


Grasp the illumination unit arm (90) as shown and place in position.

Tighten the set screw (91) in the hub of the bearer arm with the small screwdriver.

Small screwdriver
Haag-Streit part number 1001436

Allen key 5 mm
Haag-Streit part number 1001602

B.2 Power supply

The corresponding Haag-Streit operating instructions should be observed!

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HS-Part Number 1007160

Beam splitter 70/30

Rotating knob in inclined position, as shown on (98)

70% of the light go to the camera and

30% to the observer

Rotating knob in horizontal position

100% of the light go to the observer

Interface

Usually a Haag-Streit adapter ring for digital camera and a camera specific Step down ring are needed. A list of recommended cameras and requested components can be found in the internet under: www.haag-streit.com

Installation

- Remove protective cover of the video adapter.
- Remove binocular tube, attach the digital photo adapter (97) in the beam path and mount again the binocular tube.
- Loosen the fixation screw (92) of the digital photo adapter, the connecting ring (93) can now easily be removed.
- Screw the Haag-Streit adapter ring for digital camera (94) to the connecting ring.
- If necessary: screw the 'Step Down Ring' (95) to the Haag-Streit adapter ring for digital camera.
- Fix camera (96) onto the 'Step Down Ring'.
- Position the complete set on the digital photo adapter, align it horizontal and fix it with fixation screw (92).
- Check connection!

Camera cables

Haag-Streit recommends the use of a cable guide with angle (99)

(see accessories) which can be fixed on the headrest (102).

Thus all camera cables (100) for data transfer and for power can be stored in an elegant way.

• The angle is fixed with a screw (101) on the standardized tripod thread of the camera. For details see separate sheet attached to the cable guide.

Camera settings




- The camera's focus must be set at '**infinity**' (do not work in the autofocus mode!)
- **Zoom** the desired image section, observation on the camera's display.
- Experience shows, that exposure times exceeding **1/60** second do not produce good results (eye movements).
- The integrated **flash** must be **disabled**, otherwise it could shade the object.

The instruction manual of the camera is to be observed!

- Haag-Streit recommends the use of the eyepiece with double crosshair reticule (103).

The images can be viewed on the computer by using any picture viewer program. For patient administration, the photographs can also be imported into an advanced program such as **EyeCap**.

Please consult your Haag-Streit Distributor for further information.

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Beam splitter 70/30

Lever in position up (104)


70% of the light go to the camera and
30% to the observer

Lever in position down

100% of the light go to the observer

Interface

C-mount

Focal length

f = 60 mm

Installation


- Remove protective cover of the video adapter.
- Remove binocular tube, attach the video adapter (111) in the beam path and mount again the binocular tube.
- Loosen the fixation screw (107), the knurled ring (106) can now be twisted around.
- Screw the video camera onto the C-mount interface (105).
- Connect cables (109).
- Rotate the camera until the correct image position is displayed on the monitor.
- Set fixation screw again (107).

Image focus on the monitor

Needs not to be adjusted, the interface has been designed so that optimal image focus is set automatically!

Camera cables

Haag-Streit recommends the use of a cable guide (108) (see accessories) which can be fixed (110) on the headrest. Thus all camera cables for data transfer and for power can be stored in an elegant way.

Dimension of object area (in mm)

Magnification microscope	with sensor video camera	
	1/2 "	1/3 "
10x	17.3 x 13.0	13.0 x 9.7
16x	11.0 x 8.2	8.2 x 6.2
25x	6.9 x 5.2	5.2 x 3.9

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
B.5 Background illumination

a) Slit illumination and diffusor

It is possible to take images and video using the inbuilt illumination of the slit lamp and the diffusor. But this excludes simultaneous slit projection.

Installation

- Insert the diffusor into the foreseen reception slots (112) where it will snap in to secure.

b) EcoLite

For video images good results can be obtained by using the LED background illumination EcoLite.

Installation


- Remove protection cover (115) and protecting cap (118).
- Install the EcoLite (113).
- Insert the angled plug of the delivered connection cable (116) into EcoLite (114) and connect the auxiliary power (117).

c) Optical light lead and pivoting mirror

Haag-Streit recommends as optimal solution for digital photography the use of the cold light source CLS 01 with background illumination pivoting mirror in addition to the slit illumination.

Installation

- The cold light source (123) should be positioned according to the instrument table.
- Remove protection cover (120).
- Install the pivoting mirror (119)
- The light lead is fixed with the cable guides on the rail cover (122), and then passed through the cable cover (121) to the cold light source.

d) LED illumination set LI 900®

Integrated LED slit and background illumination

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HS-Part number 1001434

C Additional maintenance routines

C.1 Replacement of the mirror


Removal of the mirror is easiest when the microscope and illumination unit are well separated and the latter inclined by 10° or so.

C.2 Cleaning of the optical parts

A dust brush is to be employed for dusting the exposed glass surfaces.


Not for LI 900®: From time to time it is necessary to clean the upper surface of the condenser. The condenser surface is now visible and can be cleaned with a brush for cleaning the optics. Surfaces which have become really dirty should first be dusted, then wiped carefully with a soft dry clean cloth, washed linen, chamois leather, or some such material which will not scratch the surface.

Brush for cleaning the optics
HS-Part number 1001398

C.3 Adjustment of friction of the slit width control


The friction of the slit width control (125) can be regulated by means of the small screw in the centre of the right hand control. The movement can be stiffened by turning the screw clockwise or loosened by turning it counter-clockwise. The screw should be at least tight enough to ensure that the slit does not close spontaneously.

C.4 Cleaning of the gliding plate and the rails

If it becomes difficult to move the slit lamp around on the gliding plate (129), the plate should be cleaned with a slightly oily rag.

After removing the rail covers (128), the rails (127) should be cleaned with a stiff brush.

C.5 Cleaning of the axle (126)

Clean only with dry, lint-free cloths!

C.6 Plastic dust cover

It is recommended to protect the slit lamp with a dust cover in case it is not in use.

Plastic dust cover, small (for slit lamp)
HS-Part number 1001395

Plastic dust cover, large (for several instruments)
HS-Part number 1001434

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