LED Illumination LI 900®

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Introduction

We would like to thank you for your decision to purchase this Haag-Streit product.

If the instructions in this manual are carefully followed we are confident that this product will give you reliable and trouble-free usage.

Purpose on use

The LED illumination LI 900® have been designed exclusively for use of the slit lamps BQ 900®, BM 900® and BP 900®.

- Slit and background illumination are continuously variable controlled by a control panel.
- This set may not be used for other purposes. Slit lamps are usually used by Ophthalmologists, Optometrists or Opticians in their consulting rooms, clinics, hospitals or teaching facilities.

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1 Safety

Ambient conditions

• Transportation	Temperature	-40 °C to +70 °C
	Air pressure	500 hPa to 1060 hPa
	Relative humidity	10% to 95%
• Storage	Temperature -	10 °C to +55 °C
	Air pressure	700 hPa to 1060 hPa
	Relative humidity	10% to 95%
• Working	Temperature	+10 °C to +35 °C
	Air pressure	800 hPa to 1060 hPa
	Relative humidity	30% to 75%

Installation

- To avoid condensation, allow the instrument to adjust to room temperature for several hours before unpacking.
- Always mount the slit lamp and the headrest on an electric insulated and fire resistant table top.

Operation and surrounding

- Only qualified and trained personnel should operate the equipment, the training is at the owner's responsibility.
- Use only Haag-Streit accessories.
- Shut down after every use.
- Never operate the instrument in the same room with combustible gases, volatile solvents (alcohol, benzol) or flammable anesthetic agents.
- Avoid humidity.

Electrical

- Only a hospital grade 3-conductor electrical power supply cable must be used. *For USA and Canada: Detachable Power Supply Cord Set, UL Listed, type SJE, SJT or SJO, 3-conductor, not smaller than 18 AWG.*
- Plug, cable and ground lead connection of the socket have to be in perfect condition.

To avoid the risk of an electric shock, these devices should only be connected to a supply network with an earth wire.

5

Cleaning

- Clean the housing with a water dampened cloth only.
- No liquids or corrosive agents.

Warranty / product liability

- This instrument should be handled in accordance with the 'Safety' chapter. Improper handling may damage the instrument. Any resultant claim made under the guarantee will not be accepted.
- Continued use of an instrument that has been damaged by incorrect handling may lead to personal injury. The manufacturer cannot accept liability in this case.
- Repairs and alterations to this equipment should only be carried out by Haag-Streit service technicians or by authorized persons.

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Statutory requirements

- The LI 900® illumination set has been designed to conform with the IEC / EN 60 601-1-2 and ISO 15004-2 standards.
- Manufacturing procedures, testing, commissioning, maintenance and repair are conducted in compliance with Swiss and international regulations.
- The 'CE' mark confirms that the LI 900® illumination set complies with directive 93/42/EEC.
- The lighting set should only be used when operating type B application parts (device offering special protection against electric shock, especially with regard to reliable discharge current and reliability of the ground wire connection).
- When combining different medical and/or non medical electrical equipment standard IEC/EN 60601-1 applies.
- A copy of the declaration of conformity for the present instrument may be requested at any time from Haag-Streit.
- All statutory accident prevention regulations are to be observed.
- Classification
*IEC/EN Standard 60601-1 Equipment safety class I
Operation mode:
continuous operation*

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Environment

- Electric and electronic equipment must be disposed off separately from house waste! This device first brought to market after the 13th of August 2005.
- Disposal via your local collecting point or your Haag-Streit distributor.
- Thus ensures that no harmful materials will be released into the environment and that valuable raw materials can be used again.

EMC

The LI 900® illumination set meets the electromagnetic compatibility requirements of IEC / EN 60601-1-2. The appliance is designed to keep the generation and emission of electromagnetic interferences to a level that does not affect the normal operation of other items of equipment. The appliance is likewise designed to provide a certain level of resistance to outside sources of electromagnetic interference. Also see Appendix A.

Maintenance

The LI 900® illumination set can be operated without maintenance throughout its lifetime.

Pictograms

- 1 Read very carefully (in the instruction manual)
- 2 Refer to instruction manual
- 3 Always disconnect the electrical supply before opening the equipment!
- 4 Functional earth ground
- 5 Non protected against ingress of solid foreign objects
- 6 Direct current
- 7 Alternating current

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
Label

- 9 Serial number
- 10 Manufacturer
- 11 Date of manufacture
- 12 HS part number
- 13 Connection of non-HS illumination controller
- 14 IC 01, connection of Haag-Streit illumination control
- 15 Functional earth connection
- 16 L1, LED indicator
- 17 S1, pressure switch
- 18 Environnement

In case of LED lighting LI 900®, the slit lamps BQ 900® (a), BM 900® (b) and BP 900® (c) are powered by the power supply PS-LED or PS-LED HSM901. See separate instruction manuals.

2 Device description

The PS-LED power supply can only be used in combination with an illumination controller and an LI 01 type LED light.

2.1 Functions

Standby mode

- If the rotary button on the IC 01 illumination control is set to '0', the light source will be switched off (19).
- After running for one minute, light sources will switch off automatically and the controller will idle in energy saving mode. Any slight movement of the cross-slide or adjustment of brightness on the illumination control will cause the light sources to switch on again.

Normal operation

- The rotary button on the illumination control is set to a position between '1' and '10'. The light source is switched on.

Brightness levels, slit illumination:

Level 1 (min.) Approx. 20'000 Lux

Level 10 (max.) Approx. 450'000 Lux

Brightness levels, fill illumination FI 01f:

Level 1 (min.) Approx. 700 Lux

Level 10 (max.) Approx. 15'000 Lux




Brightness levels, fill illumination FI 01p:

Level 1 (min.) Approx. 1'500 Lux

Level 10 (max.) Approx. 30'000 Lux

Reduced operation

In order to be able to ensure a long life time for the light source, fill lighting power will be slightly reduced once the maximum operating temperature has been reached. Full power can be used again after a couple of minutes. This operating state only arises when both light sources have remained switched on simultaneously for a considerable

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


2.2 PS-LED power supply unit

2.2.1 Variant types

- 20 Power supply unit PS-LED, for mounting on non-HS tables and units
- 21 Power supply unit PS-LED HSM 901 Slide-in type for HSM-901 instrument table

20

21

2.2.2 Description

- 22 Mains switch
Switch in position 0 = 'OFF'
Power supply is disconnected from mains
Switch in position 1 = 'ON'
Control lamp shines green
- 23 3-pole mains connection
- 24 Connection for main and fixation lamps
- 25 Type plate (on housing)
- 26 Connection for Haag-Streit illumination controller IC 01 (USB mini B)
- 27 Connection for non-HS illumination controller (power jack 2.1 mm)
- 28 Protective earth connection (M3)
- 29 L1, LED indicator (green - red)
- 30 S1, Push button

2.2.3 LED indicator

Operating state LED indicator
Normal operation: Green

2.3 LED light LI 01x

4.3.1 Variant types



- 31 LED light LI 01-plus with fill illumination
- 32 LED light LI 01 without fill illumination

2.3.2 Description


- 33 24 VDC connector cable
- 34 Filter wheel for blue filter
- 35 Fill illumination connection
- 36 Cover

2.3.3 Blue filter

- 37 The rotating wheel allows the blue filter to be swivelled into position. Marker spots at same level = blue filter engaged.

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

39		40	
LED Verpolung LED tension inversée LED voltage reversal		LED Spaltbeleuchtung LED illumination fente LED Slit illuminationl	
X X X	Grün kurz blinkend Clignote brièvement vert Green, short flashes		

LED Verpolung	LED Spaltbeleuchtung
LED tension inversée	LED illumination fente
LED voltage reversal	LED Slit illuminationl
X X X	Grün Vert Green
X X X	Grün Vert Green
X X X	Grün kurz blinkend Clignote brièvement vert Green, short flashes

LI 01

LED Verpolung	LED Spaltbeleuchtung
LED tension inversée	LED illumination fente
LED voltage reversal	LED Slit illuminationl
X X X	Grün blinkend Clignote vert Green flashes

13

2.4 Illumination control IC 01

2.4.1 Variant types


- 41 Illumination control IC 01
- 42 Illumination control IC 01-1

- 43 Illumination control IC 01T
- 44 Illumination control IC 01T-1


Power cable 2 m
Power cable 5 m


- Power cables may only be used to connect an IC 01 illumination control with a PS-LED power supply unit. Do not connect the illumination control to any other devices than those intended!

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

2.4.2 Description

45 Rotary button for fill illumination
46 Rotary button for slit illumination

Optional power cable connection
47 For use to left of cross-slide
48 For use to right of cross-slide
49 LED indicateur (green - red)

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

- Observez les instructions de la table d'instruments 'HSM-901'!

2.6 Fill illumination

Fill illumination can only be used with an LI 01-plus LED illumination.

- 50 Fill illumination FI 01f
- 51 Fill illumination FI 01p

2.7 Mains connection accessories


- 52 Mains cable CH 2500 mm
- 53 Mains cable USA 3100 mm
- 54 Mains cable USA 760 mm

3 Device installation

3.1 Power supply unit on HSM-901 instrument table

- Ground wire must not be connected!

- 55 Connection for main and fixation lamps

- 56 Mains connection

- 57 Type plate on housing

3.2 Head support on instrument table




- Cable must never be clamped in, as this may lead to the possibility of a short circuit!
- Any work on cables or device components that come into connection with the power network should only be undertaken by qualified personnel.

- Connect multi-pole plug-type connector between head support and instrument table according to the illustration (58).
- Screw head support tightly into correct position on table with the two cylinder screws. Mount handgrips and connect main lamp cable to lamp housing.

- See User guide: 'HSM-901' instrument table!

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3.3 Upper part of light on slit lamp

- If the middle LED shines red during operation, the 24 V plug (60) has been connected in reverse.
- Turn off the device with the mains switch.
- Remove cover from upper part of light by releasing the screw (59).
- Rotate two-pole 24 V plug connector through 180°.

3.4 Power supply and head support on units and non-HS tables, including HSM-801

- Connect all cables. If necessary, relieve strain on mains cable with cable tie.
- Mount power supply firmly with 4 screws.
- Connect multi-pole electric plug connector between head support and instrument table.
- HSM-801: Do not connect ground wire (61)!
- Units / non-HS tables: Only connect ground wire (61) if necessary.
- Screw head support tightly into correct position on table with the two cylinder screws.
- Mount handgrips and connect main lamp cable (62) to the upper part of light.

- When assembling the guide rail and head support, make sure no electrical connection exists to any metal part of the table whatsoever. Otherwise it will be necessary to mount an earth wire on the head support.
- The power supply unit must always be connected with an earth wire on the mains side!
- Cable must never be clamped in, as this may lead to the possibility of a short circuit!
- Any work on cables or device components that come into connection with the power network should only be undertaken by qualified personnel.

- See User guide: 'HSM-801' instrument table'!




Directives for mounting on non-HS tables and HSM-801

If a power supply unit is mounted on a non-HS table the following directive will apply:

- 93/42 EEC

or the appropriate national standards

- IEC/EN 60601-1 Medical electrical equipment - General requirements for basic safety and essential performance
- IEC/EN 60601-1-2 Electromagnetic Compatibility

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3.5 IC 01T illumination controller on non-HS table




- 63 Recess in table top for IC 01T or IC 01T-1 illumination control
- 64 sufficient space must be available underneath the table top.
- 65 The illumination control is fastened to the table top with self-adhesive tape.

3.6 Non-HS illumination controllers


- An illumination controller from a foreign supplier can be plugged into the connection (66).
- Prior calibration to the appropriate potentiometer value is necessary.
- Controllers can be calibrated between 5 and 50 kΩ.

3.6.1 Calibrating non-HS illumination controllers

- Remove protective cap to access connector (59).
- Connect external illumination controller (power jack 2.1 mm) and set to maximum illumination (turn adjustment knob to right until limit stop).
- Switch on PS-LED power supply unit.
- Press button S1 for longer than 1 second. The device will change over to calibration mode. The orange LED (L1) indicates that the illumination controller is being calibrated.
- After approx. 1 to 3 seconds (depending on the type of illumination controller) the calibration process will be complete. This will be indicated by two short blinks of the green LED (L1).

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

Y = Wavelength (nm)

Power supply unit PS-LED / PS-LED HSM901

Input	100 to 240 VAC ($\pm 10\%$) 50 – 60 Hz 60 VA
Output 1	24 VDC / 1 A
Output 2	5 VDC/ 200 mA
Classification	Protection class I • IPX0 • Continuous operation • Not suitable for operation in vicinity of flammable anaesthetic products!

Dimensions:

PS-LED	PS-LED HSM901
Length	164 mm
Width	140 mm
Height	60 mm
Weight	450 g
IC 01	IC 01T
Length	75 mm
Width	35 mm
Height	33 mm
Weight	32 g
Length	316 mm
Width	146 mm
Height	69 mm
Weight	750g

IC 01	IC 01T
Length	75 mm
Width	35 mm
Height	33 mm
Weight	32 g
Length	90 mm
Width	51 mm
Height	33 mm
Weight	32 g

Maximum testing times according to
ISO 15004-2 and ISO 10939

CAUTION – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum output will exceed the safety guideline after 2.9 minutes with slit illumination or 2.5 minutes with both slit and peripheral

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
Appendix A

Guidance and manufacturer's declaration for electromagnetic compatibility (EMC)

General

- Electrical medical appliances and systems are subject to special measures regarding EMC and should be installed in accordance with these accompanying instructions.
- Portable and mobile high frequency communications equipment can negatively influence electrical medical appliances.

BQ 900°

System Components:

The Slit Lamp with LED Illumination System includes the following equipment:

- Slit Lamp BQ 900°, BM 900° or BP 900°
- LED Illumination LI 900°
- Power Supply Unit PS-LED HSM 901


Connecting Cables

The following connecting cables are required:

- Power Supply Unit CH: maximum length = 2500 mm
- Power Supply Unit USA: maximum length = 3100 mm


- Any connecting cables or accessories used to run the appliance other than those indicated can lead to increased emissions or to a reduced interference factor in the Slit Lamp with LED Illumination System.

BM 900°

Main features of performance:

The Slit Lamp with LED Illumination System serves to examine and diagnose the human eye. This is achieved by projecting onto the eye a very narrow beam of light, which can be adjusted to vary its length and width.

BP 900°

Guidelines and Manufacturer's Declaration on Electromagnetic Emissions

The Slit Lamp with LED Illumination System is intended for use in the following type of environment.

The customer or user of the Slit Lamp with LED Illumination System should ensure that it is only operated in such an environment. - See page 36.

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Table 1: Emission

Guidance and manufacturer's declaration – electromagnetic emissions		
The Slit Lamp with LED Illumination System is intended for use in the electromagnetic environment specified below. The customer or the user of the Slit Lamp with LED Illumination System should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Slit Lamp with LED Illumination System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2: Immunity (all devices)

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Table 3: Immunity (not life-supporting devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
Electromagnetic environment - guidance			
Immunity test standard	IEC 60601 test level	Compliance level	Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz - 80 MHz	3 Vrms 150 kHz - 80 MHz	$D = 1.2 \sqrt{P} \text{ 150 kHz - 80 MHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz - 800 MHz	3 V/m 80 MHz - 800 MHz	$D = 1.2 \sqrt{P} \text{ 80 MHz - 800 MHz}$
Radiated RF IEC 61000-4-3	3 V/m 800 MHz - 2,5 GHz	3 V/m 800 MHz - 2,7 GHz	$D = 2.3 \sqrt{P} \text{ 800 MHz - 2.7 GHz}$
<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and D is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^{a)} should be less than the compliance level in each frequency range ^{b)}.</p>			
Interference may occur in the vicinity of equipment marked with the following symbol:			
<p>Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Slit Lamp with LED Illumination System is used exceeds the applicable RF compliance level above, the Slit Lamp with LED Illumination System or the device, which contains it should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device containing the Slit Lamp with LED Illumination System.		
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.		
c	Possible shorter distances outside ISM bands are not considered to have a better applicability of this table.		

Table 4: Recommended separation distances (not life-supporting devices)

Recommended separation distances between portable and mobile RF communications equipment and the Slit Lamp with LED Illumination System			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz - 80 MHz $D = 1.2 \sqrt{P}$	80 MHz - 800 MHz $D = 1.2 \sqrt{P}$	800 MHz - 25 GHz $D = 2.3 \sqrt{P}$
0,01	0.12 m	0.12 m	0.23 m
0,1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (**m**) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (**W**) according to the transmitter manufacturer.


Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.


Note 3: An additional factor of 10/3 is used in calculating the recommended separation distance to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

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






PRODUCTS CERTIFIED FOR BOTH THE U.S.
AND CANADIAN MARKETS, TO THE APPLI-
CABLE U.S. AND CANADIAN STANDARDS

Haag-Streit AG
Gartenstadtstrasse 10
3098 Koeniz, Switzerland
Phone ++ 41 31 978 01 11
Fax ++ 41 31 978 02 82
eMail info@haag-streit.com
Internet www.haag-streit.com

Haag-Streit AG
Gartenstadtstrasse 10
3098 Koeniz, Switzerland
Phone ++ 41 31 978 01 11
Fax ++ 41 31 978 02 82
eMail info@haag-streit.com
Internet www.haag-streit.com

10.09 - 1

1500 7220243.04000