## HEINE OMEGA 500®

CE 93/42 EWG/CEE/EEC

## HEINE

HEINE Optotechnik GmbH & Co. KG  
 Kientalstr. 7 · 82211 Herrsching · Germany  
 Tel. +49(0)8152/38-0 · Fax +49(0)8152/38-202  
 E-mail: info@heine.com · www.heine.com  
 med 2505 2/207

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## Warranty

We guarantee this device (except bulbs) for five years from date of purchase. We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product, the buyer must prove that any such defect was present at the time the product was delivered. This guarantee excludes any defects which occur because of wear and tear, incorrect use, the use of non-original HEINE parts, (especially bulbs, as these are designed with emphasis on the following criteria: Color temperature, life expectancy, safety, optical quality and performance). Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.

For repairs, please contact your supplier, who will return the goods to us or our authorized agent.

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
position of the headband or to choose an alternative working position by pressing your view and screw the adjuster knob in (2). (It may be necessary to correct the lever (5) and swivelling the optics overband into a different position. If the light spot needs horizontal adjustment, simply rotate the headband slightly to the side. Adjust up the PD setting by closing first one then the other eye and observing an object in the middle of the light spot while adjusting the appropriate eyepiece. Repeat until the object is in the middle of the field of view, and a single image is obtained. Take off the instrument and check that the PD is symmetrically-adjusted. If not, repeat the selection procedure as above. Correct adjustment of PD is particularly important when examining through a small pupil. Each user should adjust the setting to match his own PD.)

## Instructions for use

## ENGLISH

Please read these instructions carefully before using the HEINE OMEGA 500® and keep them in a safe place for future reference.

## Safety Information

**Range of application:** The HEINE OMEGA 500® headband ophthalmoscope is designed solely for the examination of the retina.

**General Safety:** Use as little light as possible for the examination and ensure that the light is switched off after each examination.

## Safety in use:

- Use only power sources, which are approved for use in medical applications.
- The HEINE OMEGA 500® is compatible with existing 6 V HEINE ophthalmoscope power sources with the following exceptions:
  - HEINE EN 15®
  - HEINE OMEGA 500® fitted with the HC50 headband rheostat is not compatible with the power sources HEINE E10®, Accubox II, EN20-1, EN30, E8.
- We can only guarantee the performance of the instrument when fitted with original HEINE accessories and spare parts.

## Setting up and Use

## 1 HEINE OMEGA 500®

## 1.1 HEINE OMEGA 500® Ophthalmoscope

## 1.1.1 Control Locks

The aperture and filter controls can be locked into position by means of the lock (8) to avoid unintentional operation.

## 1.1.2 Aperture Control

3 different aperture sizes and a diffuser aperture can be selected by means of the aperture control (7). The choice of aperture depends mainly on the size of the patient's pupil. The diffuser is useful for photophobic patients and the examination of the periphery.

## 1.1.3 Filter Control

The filter control (6) can be used to select an interference red-free filter, a blue or a yellow filter which can be introduced into the illumination beam.

## 1.1.4 Synchronized adjustment of Convergence and Parallax

The HEINE OMEGA 500® includes HEINE's patented synchronized (US Pat. 4,684,227) adjustment of Convergence and Parallax. This feature ensures the best possible stereopsis and image quality for any pupil size by simultaneous adjustment of the observation and illumination beams. For a **dilated pupil**, move the control (11) underneath the instrument to the large circle position. For a **small pupil**, move the control to the small circle position. The control (11) can be positioned anywhere between these two settings to select the optimum parallax and stereopsis (convergence) for any pupil size. The illumination beam can be adjusted vertically by the control (9).

## 1.1.5 Teaching Mirror

The optional teaching mirror can be slid onto the fixture (10) of the optics.

## 1.1.6 Changing the bulb

- First allow the instrument to cool down.
- Disconnect the instrument from the power source.
- Unscrew the control (2) so that the optics unit (3) is free to move.
- Slide the bulb cover (12) upwards. Do not pull on the cord.
- Swivel the optics unit (3) to enable the bulb (13) to be pulled out of the bulb holder.
- Clean the glass dome of the new bulb carefully with a soft cloth.
- Insert the new bulb so that the locating pin on the bulb body fits into the guide in the bulb holder.
- Re-fit the bulb cover (12) into the optics unit (3) and slide the bulb cover back down into position.

## 1.2 HEINE OMEGA 500® Headband

## 1.2.1 Selecting the position of the optics overband

The overband (4) swivels and can be fixed into the following positions by means of the position lock (5) (upwards for the rest position, down for working positions 1 or 2).

To release the overband, press the position lock control (5) and rotate it.

Do not try to force the overband up or down.

With practice a comfortable working position can be selected and, when properly adjusted, the optics unit will always return to the desired working position.

## 1.2.2 Fitting the headband rheostat HC50

When fitting, ensure that the guide pin (15) of the headband rheostat HC 50 (14) locates correctly in the opening (16) of the overband (4). Locate the cord between the power source and the headband rheostat in the cord clip (21) at the back of the headband.

## 1.2.3 Fitting the position control

Before fitting the position control (5) ensure that the positioning plate (19) is swivelled upwards as far as possible. To do this, press the positioning plate (19) and swivel it upwards simultaneously. Ensure that the optics over-band is in the lower position and that the guide pin (20) locates in the opening (16). Now press the position control (5) on to the guide (17) until it clicks into placed.

## 1.2.4 Removing the headband rheostat HC50 and position control.

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