
SL 120

Spaltlampe Slit Lamp

Lampe à fente

Lámpara de rendija

Gebrauchsanweisung
User manual
Mode d'emploi
Instrucciones de uso

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Knowledge of this manual is required for the operation of the instrument. Therefore please make yourself familiar with the contents of this manual and pay special attention to hints concerning the safe operation of the instrument.

The specifications are subject to change; the manual is not covered by an update service.

- © Unless expressly authorized, forwarding and duplication of this document, and the utilization and communication of its contents are not permitted.
Violations will entail an obligation to pay compensation. All rights reserved in the event of granting of patents or registration of a utility model.

	Page
Copyright.....	1
Safety Instructions	4
Warning and instruction labels	4
Standards and regulations.....	4
Instructions for installation and usage	5
Reliable functioning	6
Intended use	7
Description.....	8
General construction.....	8
Construction of subassemblies/Controls	9
Accessories	12
Assembly.....	16
Assembly of the SL 120.....	16
Assembly of accessories	18
Operation.....	23
Preparations.....	23
General notes on instrument operation	24
Illumination equipment	25
Tonometry	26
Gonioscopy.....	26
Maintenance.....	27
Replacing the bulb	27
Cleaning	28
Safety inspections	28
Disposal	29
Specification	30
Slit Lamp SL120	30
Power supply (optional) for DigiCam illuminator	31
Declaration of the manufacturer	32

Important safety notes:	Page
Safety Instructions	4
General notes on instrument operation	24
Maintenance	27

Warning

Correct operation of the SL 120 Slit Lamp is imperative for its safe function. Therefore, please read these operating instructions thoroughly before using the instrument the first time. If you need supplementary information, please contact our service or authorized agents.

Warning and instruction labels

Warning

Danger for the user.

Caution

Danger for the instrument.

Instrument type B according to EN 60601-1.

Unplug the system before servicing.

Note

Information and notes for a better understanding of the operating instructions

Standards and regulations

- The European Medical Device Directive (MDD) classifies this device as a Class I device.
- The device complies with the EC Medical Device Directive 93/42/EEC and its national equivalent in form of the German Medical Product Act (MPA).

Instructions for installation and usage

- ❑ Do not use the instrument in potentially explosive atmospheres. Do not operate it in the presence of flammable anesthetics and volatile solvents such as alcohol, benzine or similar agents.
- ❑ Temperature range for intended use: +10 ... +35 °C.
- ❑ Clean optical components only externally and in accordance with the instructions given in this manual.
- ❑ The instrument may be modified and repaired only by our service or other authorized persons. The manufacturer is not liable for damages caused by unauthorized tampering of the equipment. Such tampering will also forfeit any rights to claim under warranty.

Warning

Never look into the sun through the binocular tube and eyepieces.

Warning

Let the bulb cool down before changing.

- ❑ In the vicinity of the system, do not use any mobile phones and other instruments that do not comply with EMC Class B, because their signals may cause malfunction of the equipment. The effects of radio signals on medical instruments depend on various factors and are therefore unpredictable. To avoid electromagnetic disturbances, the system has to be installed and started up in accordance with the operating instructions exclusively using the components delivered by Carl Zeiss Meditec.

Safety instructions

Reliable functioning

Caution

This instrument is a high-quality technical product. To ensure perfect and reliable function, have it checked at least every 12 months by authorized service technicians.

Electrical safety

- The user may only open the SL 120 Slit Lamp for changing the bulb, since any other operation in the instrument may endanger both operator and patient.
- Switch off the instrument with the mains switch before changing the bulb.
- The instrument's power supply is short-circuit-proof.
- Opening of the power supply is only allowed to the service staff.

Requirements for operation

Please take care that the following requirements are met before using the instrument each time:

- Check all detachable connections (e.g. binocular tube, **1**, Fig. 2) prior to each use.
- The instrument is electrically connected by means of the instrument's power cable.
- The power cable is connected to a socket with perfect earthing contact.
- All cables and plugs are in perfect condition.
- Additional movable multiple sockets or extensions lines must not be connected to the system. Do not interrupt the earthing contact intentionally either inside or outside the instrument, as this may cause the instrument casing to carry dangerous voltages. If any cables are damaged, discontinue use until repairs are carried out.

For every use of the instrument

- Avoid looking directly into the light emitting prism.

After every use

- To switch off the instrument, use the mains switch.
- Always switch off after use.

Mechanical system of slit lamp

- On slit lamp and instrument base various mechanical units can be moved relative to each other. Avoid that your fingers are caught between them.
- The patients should rest their hands on their lap to avoid contact with the microscope. This also assists in obtaining the correct posture on the headrest.

Intended use

The SL 120 Slit Lamp is a universal instrument system for observation, diagnosis, measurement and documentation of the anterior eye segment.

- Use only accessories obtained from Carl Zeiss Meditec. If you intend to use other accessories, absolutely consult the Carl Zeiss Meditec Service or Sales Department prior to use.
- The manufacturer shall not be liable for damages resulting from unauthorized tampering with the device or from the use of unapproved accessories. Besides, this will void all warranties.
- Any use differing from the one described is in the responsibility of the user.
- Read and observe the highlighted safety instructions and notes of this manual with particular attention and care.
- The instrument may only be operated by instructed and trained staff.
- It is the responsibility of the owner of the equipment to have staff trained and instructed in its use.
- Please always keep this manual at hand.

Description

General construction


Fig. 1 General construction of the SL 120 Slit Lamp

The SL 120 Slit Lamp consists of the following main subassemblies:

1 Stereomicroscope

mounted at swivel arm, comprising 5x magnification changer and binocular tube (optionally convergent or parallel tube)

2 Slit projector

with lamp housing

3 Headrest

4 Instrument base

Construction of subassemblies/Controls

Stereomicroscope

1 Hexagon socket screw

for fastening the binocular tube or the accessories for observation or documentation

2 Mounting plane

for accessories

3 Objective, f = 110 mm

4 Galilean magnification changer, 5-steps,

Magnifications of 5/8/12/20 und 32x (in combination with 110 mm objective and 10x eyepiece). The set value is facing the operator. In position 12x, which is recommended for tonometry, the rotary knob is horizontal.

Magnifications of 6/10/16/25 und 40x (in combination with 110 mm objective and 12.5x eyepiece).

5 Binocular tube

Two different types can be used:

- convergent tube or
- parallel tube

6 Plug-in eyepieces

super high eyepoint, with easily moved sliding eyecups, suitable for spectacle wearers


Fig. 2 Stereomicroscope

Slit projector with lamp housing

1 Lamp housing

2 Slit width control,

0 ... 14 mm, continuously adjustable, marks at 0.3; 1; 2 mm

3 Decentering controls for retroillumination.

4 Diffusor

for diffuse illumination, flip up cover upward onto the prism head

5 Scale for horizontally inclined illumination

(only in the equipment with lowerable prism head)

6 Knurled control

for setting 0 ... 20° horizontally inclined illumination
(only in the equipment with lowerable prism head)

7 Filter setting wheel

with three positions: red-free, free aperture, fluorescence exciter filter.

8 Scale for slit rotation

(Rotation is via the thermally isolated protective cover of the lamp housing.)

9 Slit length control

variable in steps of 0.3, 3.5, 8, and 14 mm, or 1 ... 6 mm continuously with scale


Fig. 3 Slit projector

Description

Instrument base with headrest


Fig. 4 Instrument base with headrest

- 1 3D joystick control**
for lateral and height adjustment (the latter by turning the control)
- 2 Fast action brake**
if moved off the 3D joystick control, the control clamps the slit lamp in X and Y
- 3 Fixation lamp mounting**
- 4 Head band**
- 5 Red ring marks**
indicating the middle eye level of the patient for optimum position of the patient's
- 6 Chin rest**
- 7 Hole**
for mounting the focusing rod
- 8 Lock screw**
for difference angel setting
- 9 Heigh control**
of headrest
- 10 Cable for slit illumination**
- 11 Tension control**
adjusting the tension of lateral motion

- 12 **Cable for fixation lamp**
to be connected to the power supply
- 13 **Two threaded holes**
for fastening screws for fitting the base plate to instrument table
- 14 **Red mark**
indicating the mid-position of the height adjusting range
- 15 **Brightness control**
continuously adjustable, no slit light in the left final position, in the right final position corresponds to maximum illumination intensity
- 16 **Hole for mounting the fixation lamp**
right and left
- 17 **Connection for electrical connection of fixation lamp**
right and left

Power supply


Fig. 6 Power supply

The power supply NAG SL 6 V provides the voltages for the slit lamp (illumination and fixation lamp). **It is a short-circuit-proof unit.**

Nominal line voltages: 100 ... 240 V at 50 / 60 Hz.

Power consumption: max. 75 VA.

Front panel:

- 1 **Power ON/OFF rocker switch**
- 2 **Fastening plate** with four threaded holes

Rear panel:

- 3 **Connector for the cable for slit illumination**
- 4 **Connector for fixation lamp cable**
- 5 **Connector for DigiCam Illuminator**
- 6 **Connector for mains cable**


Fig. 5 Power supply rearview

Caution

You may use only a mains cable with straight cable connector.

Description

Accessories

Optional focusing rod


Fig. 7 Focusing rod

The focusing rod can be installed in hole (7, Fig. 4). The rod can be used for determining the individual diopter setting on the eyepieces to compensate for ametropia or instrument myopia. The exact adjustment of the eyepieces is a prerequisite for good examination results and is necessary when using a video or photographic camera on a slit lamp or when using a laser slit lamp.

- To use the rod, install it in the hole at the axis with the black observation plane towards the microscope. Push the rod in all the way.
- Adjust the slit width to approx. 1 mm and set the illumination angle to 0° (i.e. slit projector in middle position).
- Set the magnification to 12x.
- Set the eyepieces to a high + diopter setting. Look through the microscope and slowly turn the diopter ring in the - direction until the image is sharply focused.

Note

The eyepieces can be adjusted even more exactly using the optional eyepiece reticles.

Fixation lamp


Fig. 8 Fixation lamp

Caution

Only use the appropriate special incandescent bulbs/blingking diodes.

AT 020/AT 030 Applanation Tonometer

This instrument for intraocular pressure measurement easily mounts to the microscope body. It is moved into the operating position by swinging its measuring arm down.

The operation of this instrument is described in a separate manual for AT 020 or AT 030. Observe the manuals of these instruments.

Beam splitter

The beam splitter is inserted in the observation path between magnification changer and convergent or parallel tube. It is mounted by means of annular dovetail mounts and hexagon head screw (**1**, Fig. 9). The beam splitter is required for the use of further accessories for documentation and co-observation (TV adapter, camera attachment, co-observation tube).

With push-pull rod (**2**, Fig. 9) a prism slide is moved to send either 100 % of light to the binocular tube (rod pushed in) or divide it into 50 % each to binocular tube and accessory attachments (rod pulled out).


Fig. 9 Beam splitter

TV-Compact Camera 020

As the beam splitter, the Compact Camera 020 is mounted between microscope body and tube.

The TV-Compact Camera delivers a video signal. For this, 50 % of the light energy of the right-hand light path are used. This, however, will not affect viewing through the eyepiece.


Fig. 10 TV-Compact Camera 020

Filter changer with barrier filter (yellow)

The filter serves as barrier filter for the work with fluorescein and with a blue filter in the observation path. You can swing the yellow filter into the light path by means of the lever shown beside. In the drawn position of the filter control lever, the filter is not in the light path.


Fig. 11 Filter changer with barrier filter (yellow)

Special eyepieces

12.5x eyepieces and a 10x micrometer eyepiece with linear and angular scales are available.

Description

10x micrometer eyepiece with linear and angular scales

This eyepiece has a linear scale of 15 mm length graduated to 0.2 mm. Use of this eyepiece requires the 8x position of the magnification changer.

A TABO angle scale graduated to 2° provides determination of the angle of inclination in fitting toric soft contact lenses. The angular measurement requires not a definite position of the magnification changer.


Fig. 12 10x micrometer eyepiece with linear and angular scales

The contact lenses (fitted and prescription lenses) carry an index mark on their periphery outside the "optical zone".

The contact lens will quickly take a preferred position by the natural movement of the eye. The index mark on the lens will no longer agree with the horizontal. This deviation, called inclination, is to be determined.

Subtractive inclination:

The index mark, the so-called DS axis (DS = dynamic stabilization axis) lies in the 1st quadrant of the TABO scheme, i.e. between 0° and 90° .

Additive inclination:

The DS axis lies in the 2nd quadrant (between 90° and 180°).

The special eyepiece has an angular scale graduated to 2° for the measurement of the inclination angle.


A gravity ball (1, Fig. 13) produces the artificial horizon required for the angular measurement.

Measurement process

- Grasping the knurled ring, turn the entire eyepiece in the tube socket until the continuous line of the linear scale coincides with the index mark on the contact lens.
- Read the value from the measuring edge (2, Fig. 13) of the gravity ball. The inclination angle is displayed according to the TABO scheme.

The reproducibility of the measurement is $\pm 0.5^\circ$. When taking the measurement, allow for a short settling time of the ball.

The image scale in the eyepiece plane is insignificant for the measurement of the inclination angle. Select the magnification changer setting of the slit lamp microscope only so as to obtain a sufficiently large object field.

1 Gravity ball
2 Measuring edge

Fig. 13 Reading the inclination

Fender plate

The fender plate is used to prevent collision of slit lamp and headrest, when the slit lamp is mounted to instrument setups. Fix the fender plate below the pivot of stereomicroscope and lamp holder by means of two M4 cheese-head screws.


Fig. 14 Fender plate

DigiCam adapter

The DigiCam Adapter allows standard digital cameras to be connected for the purpose of photographic documentation. There is a separate user's manual for connecting and operating the adapter.


Fig. 15 DigiCam Adapter

DigiCam Illuminator

The DigiCam Illuminator is a separate illumination system for illuminating the area surrounding the projected slit for the purpose of photographic or video documentation.


Fig. 16 DigiCam Illuminator

Other accessories


Breathing shield

The breathing shield may be mounted to the microscope to protect patient and physician mutually from the other's breath.

Paper pads

They provide a more comfortable chin rest for the patient and improve hygienic conditions.

Assembly of the SL 120

- 1 Screw holes in the instrument base
- 2 Fastening screws
- 3 Plastic caps
- 4 Eyepieces
- 5 Binocular tube
- 6 Hexagon socket screw at the stereomicroscope

Fig. 17 Assembly of the SL 120

- After you have unpacked the slit lamp, remove the transport locks and keep them for any future transportation. The transport locks are mainly marked by red color.

Caution

After unpacking remove the cable clip that fastens the cable of the fixation lamp to the headrest.

- Place the slit lamp onto the chosen instrument table and fasten it by means of two M4 screws (2, Fig. 17) to be screwed through the holes (1, Fig. 17). Afterwards cover the screws with the provided plastic caps (3, Fig. 17).
- If you use the slit lamp without base plate, the service staff must drill the required holes with the drilling template. In this case the further installation is also done through the service staff.
- Check if binocular tube (5, Fig. 17) is firmly mounted. Tighten up the hexagon socket screw (6, Fig. 17), if necessary.
- Make sure the eyepieces (4, Fig. 17) have been plugged into the binocular tube (5, Fig. 17) as far as it will go.
- Taking the cable length into consideration the power supply can be placed anywhere, however it is best under the table top.
- Afterwards, connect the cables of the slit lamp and fixation lamp to the power supply and, finally, connect the power cable of the power supply to a standard socket with earthing contact. The cables have different plugs to preclude misconnection.

Assembly of accessories

Mounting a beam splitter



Fig. 18 Mounting a beam splitter

- Loosen hexagon socket screw (**3**, Fig. 18) and remove the binocular tube.
- Attach the beam splitter with its annular dovetail mount to the microscope and fasten socket screw (**3**, Fig. 18) again.
- Reattach the binocular tube to the beam splitter and fasten it by means of hexagon socket screw (**2**, Fig. 18).
- The chosen accessories (e.g. TV adapter, co-observation tube) can then be attached to the ports (**1** and **4**, Fig. 18) of the beam splitter.

Mounting filter changer with barrier filter

Mount these accessories like the beam splitter between microscope and binocular tube by means of hexagon socket screws.

Assembling the DigiCam Illuminator


- | | |
|------------------------------------|---|
| 1 Plastic cover | 7 Centring hole |
| 2 Connector | 8 Rear of the power supply |
| 3 LED | 9 Socket for fixation lamp cable |
| 4 DigiCam Illuminator | 10 Socket for the DigiCam Illuminator |
| 5 Centring bolt | 11 Separate power supply (6 V) |
| 6 Slit projector swivel arm | 12 Connector for the DigiCam Illuminator |

Fig. 19 Assembling the DigiCam Illuminator

- Remove the plastic cover (**1**, Fig. 19) on the swivel arm (**6**, Fig. 19) of the slit projector.
- Insert the DigiCam Illuminator's (**4**, Fig. 19) centring bolt (**5**, Fig. 19) into the centring hole (**7**, Fig. 19), such that the LED (**3**, Fig. 19) faces towards the patient side.
- Insert the DigiCam Illuminator's power cable (**2**, Fig. 19) into the relevant socket (**10**, Fig. 19) on the slit lamp's power supply (**8**, Fig. 19).
- If the slit lamp does not have a separate power supply (e.g. examination units), a separate 6 V power supply (**11**, Fig. 19) is required for the DigiCam Illuminator. Insert the DigiCam Illuminator power cable (**2**, Fig. 19) into the relevant socket (**12**, Fig. 19) on the power supply.

Assembly

Assembling the compact video attachment

- 1** Binocular tube
- 2** Locking screw for the miniature camera
- 3** Locking screw for the binocular tube
- 4** Compact video attachment
- 5** Locking screw for the compact video attachment
- 6** Slit lamp microscope base body

Fig. 20 Assembling the compact video attachment

- Remove the binocular tube by releasing the locking screw (**5**, Fig. 20).
- Attach the compact video attachment (**4**, Fig. 20) to the microscope using the annular dovetail mount and lock it in place using the locking screw (**5**, Fig. 20).
- Now attach the binocular tube (**1**, Fig. 20) to the compact video attachment (**4**, Fig. 20) and lock it in place using the locking screw (**3**, Fig. 20).

 **Note**

Ensure that all screws are properly tightened and that the compact video attachment and binocular tube are not misaligned.

Attaching the miniature camera


Fig. 21 Attaching the miniature camera

- Bring the slit lamp sharply into focus on a test object (business card, focusing rod, etc.).

Note

The eyepieces must be precisely adjusted beforehand.
12x magnification is recommended.

- Rotate the camera head's focus lock ring (**5**, Fig. 21) until it is as far back as it will go.
- Screw the adapter ring (**4**, Fig. 21) for the camera being used onto the thread of the miniature camera (**6**, Fig. 21). Secure the adapter ring carefully using the focus lock ring (**5**, Fig. 21).

- Connect the camera to the connecting cable without turning on the camera's power unit.

Caution

Do not turn the miniature camera's power unit on until the camera cable is firmly connected at both ends!

When disassembling, the procedure should be carried out in the reverse order:

The power unit must be turned off before removing the camera cable.

- Insert the adapter ring into the opening (**2**, Fig. 21) in the compact video attachment (**1**, Fig. 21) and observe the monitor image.

Note

The miniature camera can be focused by moving it along the axis. The image can be adjusted so that it is not reversed by rotating about the axis.

- When the focus and image position are correct, tighten the locking screw for the miniature camera (**3**, Fig. 21) in order to lock the camera.

Caution

Please observe the user's manual for the camera you are using.

Warning

Ensure that the miniature camera set is approved for use as a medical product accessory.

Caution

Before setting up the instrument for use, absolutely read and observe chapter "Safety instructions".

Preparations

Slit lamp

- Set interpupillary distance of binocular tube.
- Adjusting the eyepieces:
If you wear glasses, push the eyecups in and set dioptre setting ring to 0 D. If necessary, consider any instrument myopia.
Other operators should pull the eyecups out and set their refractive powers at the dioptre scale of the eyepieces.
- Setting the diffusor:
Normally, keep the diffusor in the bottom position. Flip it up only, if you need diffuse illumination for external examination of the eye.

Warning

Before moving the slit lamp towards the patient, make sure that the diffusor (**1**, Fig. 22) on the prism head is in one of its vertical positions (**E**) to avoid any eye injuries.


Fig. 22 Possible positions of the diffusor on the prism head

Power supply

- Actuate the power switch (1, Fig. 5) of the power supply to switch it on.
- Afterwards, turn on the brightness control (2, Fig. 23) on the instrument base. No slit light in the left final position, in the right final position corresponds to the maximum illumination intensity.

Patienten positionieren

- Adjust the patient's seat height as required.
- If possible, adjust the table height in such way, that the patient is slightly leaning forward. Thus, a stable rest of the patient's chin and forehead is achieved. For the same reason, one may move the patient's seat slightly away from the instrument or set it higher.
- Adjust the chin rest (6, Fig. 23) vertically until the patient's eyes are approximately level with the red ring marks (7, Fig. 23).
- Turn the joystick (1, Fig. 23) to set the slit lamp to mid-position of the vertical adjustment range (red mark (3, Fig. 23).


Fig. 23 Adjust slit lamp

General notes on instrument operation

- Choose a slit length that only illuminates that section of the eye which is to be examined. Otherwise, sections of high reflectance (e.g. sclera) may be illuminated, thus causing disturbing glare.
- Turn the filter wheel (4, Fig. 23) to choose a filter as required.
- For retroillumination, set the illumination angle as small as possible and slightly decentre the illumination by means of control (5, Fig. 23).
- When using the tonometer, we recommend one looks with the right eye through the left eyepiece while observing the patient with the left eye.
- The use of either convergent or parallel viewing path depends on your choice.

Illumination equipment

So far, no damage to the patient's eyes by optical irradiation has been observed. Nevertheless, we recommend that you reduce the light intensity directly incident on the patient's eye to the minimum level required for the diagnosis.

Principle

Apply as little light as possible and only as much as necessary.

The light dose depends on illuminance and exposure time.

Very long exposure times at the same area of the retina are not typical, as slit lamp examination is a dynamic process.

You may reduce the illuminance on the fundus as follows:

- Reduce the light brightness.
- Do not dilate the patient's pupil unless it is really necessary.
- Only adjust the size of the illuminated area (slit width, slit length) as large as possible.
- For observation, adjust the illumination angles as large as possible.
- Use colour filter or, in special cases, a diffusing screen.

Warning

Take particular care, when examining infants, aphakic patients.

Recommendation

Turn down the lamp before every examination. This will protect the patient against excessive brightness at the beginning of the examination and increase the lamp life.

Tonometry

- The measuring prisms must be disinfected and undamaged.
Therefore, we recommend the use of several measuring prisms.
- The instrument base must be horizontal to avoid movement of the slit lamp with attached measuring prism moves towards the patient. This is of particular importance, when you take tonometric measurements without locking the instrument base.
Therefore, we recommend that you always clamp the instrument base as soon as the measuring prism touches the cornea.
- Observe the manual of the instrument AT 020 or AT 030 Applanation Tonometer.

Gonioscopy


Fig. 24 Adjusting the illumination angle on the prism head


- Contact lenses must be disinfected and undamaged.
- For gonioscopy, the three-mirror contact lens is used.
- To cover different angular ranges, you can vary the illumination angle relative to the horizontal between 0° and 20°. For this, turn the knurled ring on the prism head. A scale permits the inclination angle to be set reproducibly.

Caution

Make sure not to turn on forcibly when you reach the mechanical stop while adjusting the tilting prism head.

Replacing the bulb

- 1** Lamp holder
- 2** Bulb (6 V/20 W halogen)
- 3** Spring-loaded catch
- 4** Knurled screw

Fig. 25 Replacing the bulb

Warning

Before replacing the bulb, switch the instrument off and disconnect the power cord.

Let the bulb cool down or grasp it with protective gloves.

- Loosen knurled screw (**4**, Fig. 25) to remove the cover of the lamp housing.
- Withdraw the lamp holder (**1**, Fig. 25) of the bulb (**2**, Fig. 25).
- Press spring-loaded catch (**3**, Fig. 25) to unlock the bulb (**2**, Fig. 25).
- Then, pull out the defective bulb while tilting it slightly, if necessary. Without touching the glass bulb, push new bulb (**2**, Fig. 25) (6 V, 20 W) into the opening until it snaps in. Observe the alignment of the locating pins.
- Slide the lamp holder (**1**, Fig. 25) onto the contact pins of the bulb (**2**, Fig. 25).
- Let the catch (**3**, Fig. 25) engage again.
- Finally, attach the cover and fix it by means of the knurled screw (**4**, Fig. 25).

Cleaning

Warning

Before cleaning, turn off the instrument and disconnect the power cable!

- ❑ Due to its simple and covered construction, the SL 120 Slit Lamp is almost maintenance-free. It is, however, necessary to clean the optical components (objective, eyepiece) occasionally. To remove dust from accessible surfaces, use a soft brush.
- ❑ Clean very dirty paint surfaces with a cloth moistened (not dripping) with weak detergent.
- ❑ Clean all contaminated parts the patient gets in touch with during the examination (chin rest, forehead rest) using a disinfectant approved for this application. These parts are resistant to agents of the following categories: "Low" (e.g. wash liquids, quaternary ammonium compounds) und "Intermediate" (e.g. alcohols, Javel water, iodine; Categorization acc. to Disinfectants and range of effectiveness in compliance with the Center of Disease Control and Prevention; Atlanta/USA).
- ❑ Paper pads should be changed after every patient.
- ❑ If necessary, clean the front surface of the objective carefully with an ether/spirit mixture (1:1) applied to a cotton swab. For cleaning, switch on the illumination so that you can check the cleaning effect. Move the cotton swab or an equivalent appropriate tool of the optics cleaning kit circularly from the lens center to the edge.
- ❑ When not in use cover the instrument, to protect it from dust using the dust cover provided. Store all accessories in suitable storage boxes.

Safety inspections

To ensure proper function of the instrument, it should be subjected to a safety inspection once a year (visual inspection, measurement of the resistance of the protective conductor, and discharge current measurement). The technical safety inspection is to be performed by an authorized specialist.

Consider the applicable national safety requirements.

Disposal

This instrument contains electronic components. At the end of its life-time, the instrument should be disposed of in accordance with the relevant national regulations.

Disposal of the product within the EU

EU directives and national regulations currently in force at the time of marketing prohibit the disposal of the product specified on the delivery note in domestic waste or by municipal waste disposal companies.

For more information on disposal of the product, please contact your local dealer, the manufacturer or their legal successors. Please also take note of the manufacturer's current online information.

If the product or its components are resold, the seller has the duty to notify the buyer that the product must be disposed of in accordance with currently valid national regulations.

Slit Lamp SL120

Line voltage	100 V -10 % ... 240 V +10 %, 50 / 60 Hz U.S. and Canada only: single-phase connection 120 V AC with plug type NEMA 5-15P
Power consumption	75 VA
Protection class	I
Instrument type	B (acc. to EN 60601-1)
Enclosure protection	IP 20
Environment conditions for intended use	Temperature: +10 °C ...+35 °C Humidity: 30 ... 75 % (no condensation) Air pressure: 800 ... 1060 hPa
Storage conditions (in original packing)	Temperature: -10 °C ...+55 °C Humidity: 10 ... 95 % Air pressure: 700 ... 1060 hPa
Shipment conditions (in original packing)	Temperature: -40 °C ...+70 °C Humidity: 10 ... 95 % Air pressure: 500 ... 1060 hPa
Mode of operation	continuous operation
Light source for projection	6 V / 20 W, halogen lamp
Lamp voltage	continuously adjustable
Fixation lamp	5 V / 1.5 W
Magnification	5x, 8x, 12x, 20x, 32x (with high eyepoints 10x) 6x, 10x, 16x, 25x, 40x (with high eyepoints 12.5x)
Field of view	40 mm ... 6 mm (with high eyepoints 10x) 31 mm ... 5 mm (with high eyepoints 12.5x)
Eyepiece magnification	either 10x or 12.5x high eyepoint, ±8 D compensation of ametropia
Slit length	0 to 14 mm continuously progressive
Slit width	in steps at 0.3 / 3.5 / 8 / 14 mm; 1 ... 6 mm continuously
Slit rotation	±90° continuously
Slit decentration	±4° horizontal, click-stop at 0°
Swivel range of slit prism	180°, scale for angular difference click-stops at -10°, 0°, +10°
Angle of incidence	0° or 0° ... 20° with tild prism head
Filters	blue, green (red-free); heat absorbing filter (fixed); swing-in diffusor in front of the prism
Free working distance	65 mm
Travel of instrument base	30 mm (vertical), 110 mm (X), 90 mm (Y)
Vertical travel of chin rest	58 mm
Mass	9.75 kg (basic unit), 1.25 kg (headrest)

Power supply (optional) for DigiCam illuminator

Power supply voltage	100 V -10 % - 240 V +10 %, 50 / 60 Hz
Rated output voltage	6 V DC
Dimensions (W x H x D)	52 mm x 34 mm x 88 mm
Weight	140 g
Power consumption	400 mA
Electrical protection class	II
Protection type	IP 20

Subject to changes of design and scope of supply due to continuous product improvement.

Declaration of the manufacturer

The device SL 120 Slit Lamp meets the requirements of the EC Medical Device Directive 93/42/EEC and of their national equivalent in form of the German Medical Product Act.

Instrument class according to the German Medical Product Act: I

UMDNS-No.: 12-281

This declaration ceases to be valid if the product was modified without the manufacturer's authorization.

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**CARL ZEISS MEDITEC AG**