
**INSTRUCTION AND SERVICE
MANUAL**

MODELS 8500 & 8500M

HAND HELD PULSE OXIMETERS

English
Français
Deutsch
Italiano
Español


NONIN[®]

**INSTRUCTION AND SERVICE
MANUAL**

MODELS 8500 & 8500M

HAND HELD PULSE OXIMETERS

Declaration of Conformity

Manufacturer : NONIN Medical, Inc.
Address : 2605 Fernbrook Lane North
: Plymouth, Minnesota 55447-4755
USA
Model Numbers : 8500, 8500M and NONIN Pulse
Oximeter Sensors

We herewith declare that the above mentioned pulse oximeters and sensors are classified as Class IIb and comply with EC Directive No. 93/42/EEC concerning medical devices.

This declaration is based on:


**Certification of a Quality System to ISO 9001 : 1994 /
EN 46001 : 1993**

Certificate Number: Q1 95 08 24497 001
Issued by: TUV Product Service GmbH
Date: August 30, 1995

**Certification according to Annex II, Clause 3 of EC
Directive No. 93/42/EEC concerning medical devices.**


Certificate Number: GI 95 08 24497 002
Issued by: TUV Product Service GmbH
Date: August 30, 1995

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Authorized EC Representative:
MPS, Medical Product Service GmbH
Borggasse 20
D-35619 Braunfels, Germany

Telephone: +49-(06442)-962073
FAX: +49-(06442)-962074

ABOUT THE MANUAL

Boxes draw attention to contraindications, warnings, and cautions.

The boxed messages are important. Be sure to read them.

TRADEMARKS

NONIN is a registered trademark of NONIN Medical, Inc.

REFERENCES

References to "NONIN" in this manual shall imply NONIN Medical, Inc.

References to "8500" in this manual shall imply Models 8500 and 8500M.

The information in this manual has been carefully checked and is believed to be accurate. In the interest of continued product development, NONIN reserves the right to make changes and improvements to this manual and the products it describes at any time, without notice or obligation.

COPYRIGHT 1997

NONIN Medical, Inc.
2605 Fernbrook Lane North
Plymouth, MN 55447-4755 USA

• (612) 553-9968 • (800) 356-8874 • FAX (612) 553-7807


TABLE OF CONTENTS

I. CONTRAINDICATIONS, WARNINGS, AND CAUTIONS.....	1
II. UNPACKING YOUR PULSE OXIMETER	3
III. INTRODUCTION	4
A. INDICATIONS FOR USE.....	4
B. GENERAL	4
IV. OPERATING INSTRUCTIONS.....	5
A. BATTERIES.....	5
B. GENERAL	7
1. Connect Sensors.....	7
2. Turn on the Pulse Oximeter.....	7
3. Verify Operation	9
4. Cleaning the Pulse Oximeter.....	9
V. FEATURES.....	10
A. CONTROLS.....	10
1. Power	10
2. Display Brightness.....	10
3. Setup Mode	11
a. Printer Settings.....	11
b. Calendar Settings.....	13
c. Clock Settings	14
B. DISPLAYS	15
1. SpO ₂ Display.....	15
2. ♥ (Pulse Rate) Display	15
3.  (Perfusion) Indicator	15
4. Flashing Display	16
5. Dash in SpO ₂ Display.....	16
C. PRINTER/SERIAL OUTPUT.....	17
D. MEMORY OPTION (8500M ONLY)	18
1. Recording Sessions	18
2. Memory Output Mode.....	19
VI. PULSE OXIMETER SENSORS.....	21
A. GENERAL SENSOR DESCRIPTION	21
B. FINGER CLIP SENSORS.....	24
1. Adult Finger Clip Sensor	24
2. Adult Articulated Finger Clip Sensor	24
3. Pediatric Finger Clip Sensor.....	25

C.	FLEX SENSORS	26
1.	Adult Flex Sensor	26
2.	Infant Flex Sensor	27
3.	Neonatal Flex Sensor	27
D.	EAR CLIP AND REFLECTANCE SENSORS	28
1.	Ear Clip Sensor	28
2.	Reflectance Sensor	29
E.	ADULT AND PEDIATRIC DISPOSABLE SENSORS.....	29
F.	INFANT AND NEONATAL DISPOSABLE SENSORS	31
1.	Infant Flexi-Form® Sensors	31
2.	Neonatal Flexi-Form® Sensor.....	32
G.	CLEANING THE SENSORS.....	34
H.	SENSOR COMPATIBILITY	34
VII.	THEORY OF OPERATION	35
VIII.	SPECIFICATIONS	37
IX.	SERVICE	39
X.	WARRANTY	43
XI.	ACCESSORIES	44
XII.	TROUBLESHOOTING GUIDE.....	45

TABLES


Table I:	Printer, Clock, and Calendar Mode Parameters	11
Table II:	8500P Printer Modes	12
Table III:	Printer/Sensor Interface Assignments	17
Table IV:	8500M Memory Data Format	20
Table V:	Reusable Sensor Reference	22
Table VI:	Disposable Sensor Reference	23
Table VII:	Sensor Accessory Reference	23

FIGURES

Figure I:	Replacing Batteries in the 8500	6
Figure II:	Connecting Sensors to the 8500	7
Figure III:	Front View of the 8500	8
Figure IV:	Rear View of the 8500	8
Figure V:	Flow Chart for Setting Printer Mode	12
Figure VI:	Flow Chart for Setting Calendar	13
Figure VII:	Flow Chart for Setting Clock	14
Figure VIII:	Finger Clip Sensor	24
Figure IX:	Adult Articulated Finger Clip Sensor	25
Figure X:	Flex Sensor	26
Figure XI:	Flex Sensor with Tape Strain Relief	26
Figure XII:	Infant Flex Sensor on Big Toe	27
Figure XIII:	Neonatal Flex Sensor on Foot	28
Figure XIV:	Adult/Pediatric Flexi-Form® Sensor	30
Figure XV:	Infant Flexi-Form® Sensor	32
Figure XVI:	Neonatal Flexi-Form® Sensor	33

GLOSSARY OF SYMBOLS

SYMBOL	DEFINITION OF SYMBOL
	ATTENTION Consult Accompanying Documents
	Type BF Equipment Type B equipment with an F-type applied part (patient isolation from electric shock)
	Canadian Standards Association
UL	Underwriters Laboratories
	CE Marking indicating conformance to EC directive No. 93/42/EEC concerning medical devices

I. CONTRAINDICATIONS, WARNINGS, AND CAUTIONS

CONTRAINDICATIONS

- Do not operate in an explosive atmosphere.
- Do not operate the NONIN 8500 in an MRI environment.
- Do not use the 8500 in a situation where alarms are required. The 8500 and 8500M models have no audible alarms.

WARNING

- The 8500 is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.
- Use only NONIN manufactured sensors. These sensors are manufactured to meet the calibration requirements for NONIN pulse oximeters. Use of other manufacturer's sensors may cause improper pulse oximeter performance.
- Check application site **frequently** to determine circulation, positioning, and skin sensitivity of the patient. Each patient's sensitivity to NONIN sensors may vary depending on their medical status or the condition of their skin.
- Use of NONIN double-backed adhesive strips or the Hydrogel tape strips should be discontinued if the patient exhibits allergic reactions to the adhesive material.
- Do not stretch the adhesive tape while applying the sensors. This may cause inaccurate readings or skin blisters.

 **CAUTION**

- Federal law (USA) restricts this device to sale by or on the order of a physician.
- Read this manual carefully before use of the 8500.
- Carefully read the instructional insert provided with the sensor before use.
- The 8500 is intended for spot checking or continuous monitoring by an attending health care professional. Because the 8500 has no audible alarms, international labeling requirements dictate it be labeled "Not for continuous monitoring."
- The 8500 must be able to measure the pulse properly to obtain accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO₂ measurement.
- Fingernail polish may reduce light transmission and thereby affect SpO₂ accuracy.
- The 8500 may not work on all patients. If you are unable to achieve stable readings, discontinue use.
- The 8500 is sensitive and must be repaired by knowledgeable and specially trained personnel only.
- The 8500 may interpret motion artifact of sufficient amplitude and regularity as good perfusion (green).
- The 8500 is calibrated to determine the percentage of arterial oxygen saturation of functional hemoglobin. Significant levels of dysfunctional hemoglobin such as carboxyhemoglobin or methemoglobin may affect the accuracy of the measurement.
- Cardiogreen and other intravascular dyes, depending on the concentration, may affect the accuracy of the SpO₂ measurement.
- Do not use different types of batteries at the same time. Do not mix fully charged and partially charged cells at the same time. These actions may cause the batteries to leak.
- Do not immerse the 8500 or NONIN sensors in liquid to clean.
- Do not use caustic or abrasive cleaning agents.
- Ear Clip and Reflectance sensors are not recommended for pediatric or neonatal use. The accuracy of the sensors has not been established for pediatric or neonatal use.
- Do not remove any covers other than the battery cover when battery replacement is necessary. There are no user serviceable parts inside other than the replaceable batteries.
- Alkaline batteries may leak or explode if used or disposed of improperly.

II. UNPACKING YOUR PULSE OXIMETER

Contact the carrier immediately if the shipping carton for the 8500 is damaged. Carefully unpack the instrument and its accessories. Confirm that the items listed below are packed with the 8500 Hand Held Pulse Oximeter. The 8500/8500M shipment includes:

- 8500 or 8500M Hand Held Pulse Oximeter
- Instruction and Service Manual for Models 8500 and 8500M
- Six AA batteries
- 8000AA Sensor (Adult Articulated Finger Clip) or 8000K2 Sensor (Adult Finger Clip)
- 8500P brochure

If any item on this list is missing or damaged, do not use the pulse oximeter. Contact your local distributor or, if you do not know your local distributor, contact NONIN's Technical Service department at (800) 356-8874.

III. INTRODUCTION

A. INDICATIONS FOR USE

The 8500 Hand Held Pulse Oximeter is intended to be used for monitoring oxygen saturation and pulse rate for adult, pediatric, and neonatal patients in hospital, ambulatory, home, and EMS environments. The 8500 may be used for spot checking and/or continuous monitoring when attended by a healthcare professional. The variety of individual sensors available must be checked frequently to ensure proper circulation and application.

B. GENERAL

The 8500 Hand Held Pulse Oximeter is small and light weight. The 8500 has visual alarms for tracking patient status. It typically will operate for 100 hours continuously between battery replacements. **The 8500 requires no routine calibration or maintenance.**

The 8500 determines arterial oxyhemoglobin saturation (%SpO₂) by measuring the absorption of red and infrared light passed through tissue. Changes in absorption caused by pulsation of blood in the vascular bed are used to determine arterial saturation and pulse rate.

Oxygen saturation and pulse rate are displayed on light emitting diode (LED) digital displays. On each detected pulse, the perfusion LED flashes. Patient perfusion signals are graded as good, marginal or inadequate and are indicated as such by the LED flashing green, yellow, or red. This simple method gives the user a pulse-by-pulse visual indication of waveform signal quality without requiring the user to perform complex waveform analysis during critical patient care situations.

Sensor disconnect or malfunction is indicated by lack of good perfusion flashes and/or a dash to the left of the SpO₂ value on the LED display. Ultimately, if adequate perfusion pulses are not received, the SpO₂ and Pulse Rate (♥) numerical values will be replaced by dashes. When the batteries are low, the digital displays will blink.

The 8500 Pulse Oximeter may be used with all NONIN pulse oximeter sensors (except fiber optic) enabling proper operation for nearly every patient.

☛ **NOTE:** *NONIN sensors do not contain natural rubber latex. Natural rubber latex has been found to cause allergic reactions in some individuals.*

IV. OPERATING INSTRUCTIONS

A. BATTERIES

 CAUTION	<ul style="list-style-type: none">• Do not remove any covers other than the battery cover when battery replacement is necessary. There are no user serviceable parts inside other than the replaceable batteries.• Alkaline batteries may leak or explode if used or disposed of improperly.• Do not use different types of batteries at the same time. In addition, do not mix fully charged and partially charged cells at the same time. These actions may cause the batteries to leak.
--	--

The 8500 Pulse Oximeter is powered by 6 AA Alkaline cells that will typically provide 100 hours of continuous operation. **The 8500 indicates when the batteries are low by flashing the digital displays once each second.** When the displays begin flashing, the batteries should be replaced as soon as possible. Replace the batteries by removing the battery door on the bottom of the 8500. Be sure to follow the polarity markings on the rear label of the pulse oximeter when installing new batteries. Refer to Figure 1 for an illustration of battery replacement.

Rechargeable Nickel Cadmium batteries may be used in the 8500 if desired. Since NiCad batteries have less than half the capacity of alkaline batteries, the batteries will have to be recharged more often than every 100 hours.

- **NOTE:** *Reducing the LED display brightness can extend battery life up to 100%! Refer to Section V (A. 2).*
- **NOTE:** *Replacing batteries erases the clock settings of the 8500. The memory of the 8500M will also be erased when the batteries are replaced. These settings should be reset after the batteries are replaced. To do this, refer to Section V (A. 3).*
- **NOTE:** *Batteries should be removed if the 8500 is going to be stored for more than 30 days. Batteries may leak if left in the device for a long period of time.*


Figure I: Replacing Batteries in the 8500

B. GENERAL

The 8500 Hand Held Pulse Oximeter is portable and is intended for attended patient monitoring by trained personnel. It displays numerical values for oxygen saturation and pulse rate. Because the 8500 Pulse Oximeter has no patient alarms, the user must frequently observe the SpO₂ and Pulse Rate displays.

1. Connect Sensors

Connect the sensor to its 9-pin mating jack on the top of the 8500 as shown in Figure II. To select a sensor type, refer to Section VI (A). If additional cable length is necessary, connect the Model 8500I Patient Cable between the sensor and the 8500 Pulse Oximeter. Position the appropriate sensor on the patient. Refer to Section VI for sensor positioning information.


Figure II: Connecting Sensors to the 8500

2. Turn on the Pulse Oximeter

Turn on the 8500 Hand Held Pulse Oximeter by depressing the " | " button on the front of the pulse oximeter. Refer to Figure III.

When the 8500 is powered on, the SpO₂ and ♥ displays will cycle through the following sequence before displaying valid data values:

- "888 888"
- time saved in memory in hours and minutes
- software revision number
- " - - "


Figure III: Front View of the 8500



Figure IV: Rear View of the 8500

3. Verify Operation

-  **CAUTION**
- The 8500 must be able to measure the pulse properly to obtain accurate SpO₂ measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO₂ measurement.

Verify that the sensor is properly positioned. Ensure the system is sensing adequate perfusion by observing that the  (perfusion) indicator is blinking green, and the blinking is correlated to the pulse rate for 10 seconds. Should the perfusion light be red or yellow or flashing erratically, reposition the sensor or try a different sensor.

4. Cleaning the Pulse Oximeter

-  **CAUTION**
- Do not immerse the 8500 in liquid to clean.
 - Do not use caustic or abrasive cleaning agents.

The 8500 Pulse Oximeter may be cleaned with a mild detergent and a soft cloth or with an isopropyl alcohol wipe. Allow enough time for the 8500 to dry thoroughly before reusing.

V. FEATURES

A. CONTROLS

All functions of the 8500 are controlled by switches found on the front of the unit. Refer to Figure III for an illustration of these switches.

1. Power

Pressing the ON switch ("|") causes power to be applied to all internal circuitry. Pressing the OFF switch ("Ø") causes power to be removed from the displays and puts the pulse oximetry circuitry into a low power standby mode.

In order to conserve battery life, the 8500 will automatically power off after 10 minutes of inactivity. Inactivity is indicated by dashes on the displays and is caused by:

- no sensor connected to the pulse oximeter
- patient pulse too low
- sensor not attached to a patient

Each time a reading is displayed, the 10 minute timer is restarted.

The "|" switch has additional clock and printer mode setting functions when used in conjunction with the "↑" switch. Refer to Section 3 of this section for Setup Mode instructions.

2. Display Brightness

The arrow switch ("↑") causes the brightness of the digital displays to change. When powered up, the digital display defaults to the maximum brightness. Pressing the "↑" will advance the brightness to the lowest setting, and each subsequent press will advance the brightness through 8 different settings. The function is circular which means it will cycle through the entire brightness range and then start at the beginning again. The lower light brightness may be used to preserve battery life. The higher light intensity may be used to view the displays from a distance.

➤ **NOTE:** *Reducing the LED display brightness can extend battery life up to 100%!*

The "↑" switch has additional clock mode setting functions when used in conjunction with the "|" switch. Refer to Section 3 of this section for Setup Mode instructions.

3. Setup Mode

Setup mode is used to control the internal time-of-day clock and the external real-time printer (purchased separately). The setup mode is initiated by holding the "↑" switch when the unit is turned on by pressing the "↓" switch. In setup mode, the "↓" switch and the "↑" switch are used to make the selections.

Advance to the next sequential mode by pressing the "↓" switch. Each time the "↑" switch is pressed, the number on the ♥ display will increment. It starts with the current value stored in memory for the parameter designated in the SpO₂ display. When the correct value appears in the ♥ display, pressing the "↓" switch will advance the SpO₂ display to the next sequential parameter as listed in Table I. This process is continued until all parameters are set. The settings can be easily checked, since the first value displayed for each parameter represents the current setting. When the setting sequence has been completed, the 8500 exits the setup mode and begins normal operation.

<u>Sequence</u>	<u>Appears in SpO₂ Display</u>	<u>Range of Values</u>	
		<u>From</u>	<u>To</u>
Printer	P r n	00	15
Year	y	00	99
Month	n n	00	12
Day	d	01	31
Hours	h h	00	23
Minutes	n n	00	59

Table I: Printer, Clock, and Calendar Mode Parameters

a. Printer Settings

➤ **NOTE:** *The 8500P prints out real-time data only. Data stored in the 8500M can not be downloaded to the 8500P.*

"P r n" will appear in the SpO₂ display indicating print setup mode. There are 16 options for the printer mode: 00 through 15. Each printer mode is explained in Table II. The modes determine how often data is written out to the printer and the format of data written to the printer. Refer to Figure V for a flow chart of setting the printer mode. When the printer setting sequence has been completed, the 8500 continues to the calendar settings (refer to section b of this section).


Figure V: Flow Chart for Setting Printer Mode

<u>Printer Mode</u>	<u>Seconds per data point</u>	<u>Minimum SpO₂ printed?*</u>	<u>Touch print capability?*</u>
00	10	No	Yes
01	30	No	Yes
02	120	No	Yes
03	10	Yes	Yes
04	30	Yes	Yes
05	120	Yes	Yes
15	---	No	Yes

* For modes where the minimum SpO₂ data is written, there are two lines of data written for each data output. The first line contains the minimum value for SpO₂ since the last printout, and the second line contains the current data values.

** The touch print mode enables the user to print out data at any time. This is activated by pressing the " | " switch.

Table II: 8500P Printer Modes

☛ **NOTE:** Printer modes 06 – 14 are not used at this time. They are reserved for future development.

b. Calendar Settings

After the printer setting has been determined in the setup mode, "y" will appear in the SpO₂ display indicating calendar setup mode for the year. The year may be set to "00" through "99". After selecting the year, the display will show "n n" indicating the setup mode for the month. The month may be set to "00" through "12". After selecting the month, the display will show "d d" indicating the setup mode for the day of the month. The day may be set to "01" through "31". Refer to Figure VI for a flow chart of setting the calendar. When the calendar setting sequence has been completed, the 8500 continues to the clock settings (refer to section c. of this section).

➤ **NOTE:** Setting the month to "00" disables the clock function and helps conserve battery life.


Figure VI: Flow Chart for Setting Calendar

c. Clock Settings

After the calendar settings have been determined in the setup mode, "h" will appear in the SpO₂ display indicating clock setup mode for the hour. The hour may be set to "00" through "23". After selecting the hour, the display will show "n n" indicating the setup mode for the minutes. The minutes may be set to "00" through "59". After selecting the minutes, the display will return to normal operation. Refer to Figure VIII for a flow chart of setting the clock.


Figure VII: Flow Chart for Setting Clock

The printer modes are explained in Table II. The modes determine how often data is written out to the printer and the format of data written to the printer. For modes where the minimum SpO₂ data is written, there are two lines of data written for each data output. The first line contains the minimum value for SpO₂ since the last printout, and the second line contains the current data values. The touch print mode enables the user to print out data at any time. This is activated by pressing the " | " switch.

B. DISPLAYS

1. SpO₂ Display

The left digital display is a 3-digit light emitting diode (LED) digital display that indicates oxygen saturation percentage. Refer to Figure III. Oxygen saturation is calculated by the methods outlined in Section VII.

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
The right digital display is a 3-digit LED digital display that indicates pulse rate in pulses per minute. Refer to Figure III.

3. (Perfusion) Indicator

The perfusion indicator (identified by the waveform symbol ) will flash once for each pulse while measuring oxygen saturation. The perfusion indicator changes color to indicate changes in the pulse waveform signal that may affect the SpO₂ data.

The perfusion indicator may blink one of three colors: green, yellow, or red. These colors are similar to the colors of a stoplight such that:

Red indicates the pulse amplitude is too small. During red perfusion, SpO₂ and pulse rate values are not updated. After ten seconds, the values are replaced with dashes indicating SpO₂ measurement is not possible.

Yellow indicates the pulse waveform amplitude is marginal or the pulse oximeter has detected artifact. Although SpO₂ data is acceptable, corrective measures should be considered to improve sensor placement, change sensor type, or reduce patient movement.

Green indicates the pulse waveform signal is of good quality and SpO₂ data is accurate.

CAUTION

- The 8500 may interpret motion artifact of sufficient amplitude and regularity as good perfusion (green).

4. Flashing Display

The numerical displays will flash once each second if the 8500 determines that a battery low condition exists. Replace all six batteries immediately.

• **NOTE:** *Inaccurate SpO₂ and/or pulse rate measurement may result if the 8500 is operated in a low battery condition.*

5. Dash in SpO₂ Display

If the 8500 determines that a sensor fault exists (sensor disconnect, sensor failure, or sensor dislodgment), a dash (-) appears in the left-most digit of the SpO₂ display. The readings that are being displayed will remain unchanged while the sensor fault exists. If the sensor fault is not corrected, dashes will be displayed 10 seconds after the minus sign appears.

C. PRINTER/SERIAL OUTPUT

Both the 8500 and 8500M Hand Held Pulse Oximeters provide output capability to a custom printer via the 9-pin Sub-D connector. This connector serves as a sensor input connector as well as a printer interconnect device. The 9-pin Sub-D connector pin assignments are listed in Table III.

<u>Pin Number</u>	<u>Assignment</u>
1	Battery Voltage
2	Infrared Anode, Red Cathode
3	Infrared Cathode, Red Anode
4	Serial Data, TTL Levels
5	Detector Anode
6	Logic Level
7	Cable Shield
8	Coaxial Shield
9	Detector Cathode, +5 V

Table III: Printer/Sensor Interface Assignments

The information from the 8500 in the real-time mode is sent in an ASCII serial format at 9600 baud with 9 data bits, 1 start bit, and 1 stop bit. The data is output at a rate of once per second.

• **NOTE:** *The 9th data bit is used for parity in memory playback mode. In real-time mode, it is always set to the mark condition. Therefore the real-time data may be read as 8 data bits, no parity.*

The data printed by the 8500P printer is in the following format:

HH:MM:SS SPO₂=XXX HR=YYY

where "HH" represents the hour the real-time clock is set to, "MM" represents the minutes, "SS" represents the seconds, "XXX" represents the SpO₂ value, and "YYY" represents the heart rate. The SpO₂ and heart rate will be displayed as "---" if there is no data available for the data reading.

D. MEMORY OPTION (8500M ONLY)

The memory option is identified by the "M" in the model number (i.e. 8500M as opposed to the 8500). This model number is located just above the serial number on the back of the unit. The 8500M Pulse Oximeter can collect and store up to eighteen hours of SpO₂ and pulse rate information. An infrared optical link or 8500YC cable transfers the data to the 8586PI Printer Interface unit where it can be printed on an Epson®-compatible graphics printer (parallel version).

The solid-state memory in the 8500M functions much like the "endless loop" tapes that are used in some telephone answering machines. When the memory fills up, the unit begins overwriting the oldest locations with the latest data.

Each time the 8500M is powered up, the current time/date information (if the clock is set properly) is stored in memory to allow quick differentiation of recording sessions. Patient SpO₂ and pulse rate are sampled and stored every four seconds. The stored resolution of the oxygen saturation is in 1% increments in the range of 0 to 100 %. The stored pulse rate ranges from 18 to 300 BPM. The stored values have a resolution of 1 BPM from 18 to 200 and a resolution of 2 BPM from 201 to 300.

During the printing of the data, the last data recorded is the first data printed. For example, the last four minutes of data recorded would be the first four minutes of printout.

1. Recording Sessions

Each time the 8500M is turned on (except while setting the clock) data is automatically collected.

• **NOTE:** *Only recording sessions greater than one minute in length are kept in memory for later printing.*

2. Memory Output Mode

To output the data stored in the memory of the 8500M, start with the unit OFF and then:

- 1) Hold the "↑" switch while pressing the "↓" switch;
- 2) Release the "↑" and "↓" switches when "000 000" is displayed in the SpO₂ and pulse rate LEDs;
- 3) Observe "Prn XX" will be displayed in the SpO₂ and ♥ LEDs;
- 4) Data is automatically transferred from the memory in approximately 8 to 10 minutes.

Data is transferred at a rate of 20 minutes of collected data per second. An 18 hour recording session (the maximum memory saved) is transferred in approximately 1 minute. After all the data is transferred, the 8500M should be shut off prior to collecting new patient data. The patient information is held in memory as long as the batteries are good, so if the memory has to be cleared, remove the battery door for a period of 60 seconds or longer. Outputting the memory does not clear any data from the memory.

☛ **NOTE:** *The memory is cleared each time the batteries are changed.*

Additional information on the operation of the 8500M, and the 8586PI Printer Interface system can be found in the manual included with the 8586PI Printer Interface.

The format of the data transferred is given in Table IV. The size of this file will depend on the amount of data saved in the memory. The most recent data is transferred first. The memory data format is in binary. Bad data is represented by FF (hexadecimal) or 255 (decimal). If the memory "wrapped around" (the recording time exceeded 18 hours) and the final (i.e. the oldest) file of data has been truncated, the final start time will be represented by zeroes and the start times for that file will then not match up.

8500M MEMORY DATA FORMAT

DATA 1	DATA 2	CHECK SUM**	NOTE
\$FE	\$FE	\$FC	Header of first file
month	Day	XX	Current time
year	Minute	XX	
second	Hour	XX	
month	Day	XX	Stop time of first file saved in memory
year	Minute	XX	
second	Hour	XX	
month	Day	XX	Start time of first file saved in memory
year	Minute	XX	
second	Hour	XX	
pulse rate	SpO2	XX	Data in first file
pulse rate	SpO2	XX	
pulse rate	SpO2	XX	
.	.	.	
.	.	.	
pulse rate	SpO2	XX	
\$FE	\$FE	\$FC	Header of second file
month	Day	XX	Start time of first file saved in memory
year	Minute	XX	
second	Hour	XX	
month	Day	XX	Stop time of file <i>N</i> saved in memory
year	Minute	XX	
second	Hour	XX	
month	Day	XX	Start time of file <i>N</i> saved in memory
year	Minute	XX	
second	Hour	XX	
pulse rate	SpO2	XX	Data in file <i>N</i>
pulse rate	SpO2	XX	
pulse rate	SpO2	XX	
.	.	.	
.	.	.	
pulse rate	SpO2	XX	
\$FE	\$FE	\$FC	Header of file <i>N+1</i>
month	Day	XX	Start time of file <i>N</i> saved in memory
year	Minute	XX	
second	Hour	XX	
0	0	0	End of data marker (18 zeroes)
0	0	0	
0	0	0	
0	0	0	
0	0	0	
0	0	0	

NOTE: ** CHECKSUM = DATA1 + DATA2

The data represented by the shaded area is repeated for additional data blocks that are saved in the 8500M's memory. *N* represents the "next file" saved in the memory.

Table IV: 8500M Memory Data Format

VI. PULSE OXIMETER SENSORS

WARNING

- Use only NONIN manufactured sensors. These sensors are manufactured to meet the calibration requirements for NONIN pulse oximeters. Use of other manufacturer's sensors may cause improper pulse oximeter performance.
- Use of NONIN double-backed adhesive strips or the Hydrogel tape strips should be discontinued if the patient exhibits allergic reactions to the adhesive material.
- Check application site **frequently** to determine circulation, positioning, and skin sensitivity of the patient. Each patient's sensitivity to NONIN sensors may vary depending on their medical status or the condition of their skin.

CAUTION

- Before use, carefully read the instructional insert provided with the sensor.

➤ **NOTE:** *NONIN sensors do not contain natural rubber latex. Natural rubber latex has been found to cause allergic reactions in some individuals.*

A. GENERAL SENSOR DESCRIPTION

The 8500 features sensors which are ideal for nearly every application. Refer to Tables V and VI for a complete list of NONIN sensors and specific information on each sensor. Each sensor is designed for a specific site application and specific patient size and weight range. Reusable sensors are designed for multiple patient use, while disposable sensors are designed for single patient use.

Factors that may degrade performance:

- excessive ambient light
- excessive motion
- electrosurgical interference
- arterial catheters, blood pressure cuffs, infusion lines, etc.
- moisture in the sensor
- improperly attached sensor
- incorrect sensor type
- poor perfusion
- venous pulsations
- anemia or low hemoglobin concentrations
- cardiovascular dyes
- sensor not at heart level

NONIN Reusable Sensor Reference

Sensor Description	NONIN Model Number	Approximate Patient Weight Range	Preferred Site	Application
Adult Finger Clip	8000K2	greater than 30 kg (greater than 66 lbs)	Index, middle, or ring finger	<ul style="list-style-type: none"> • spot checking • short-term continuous use • for adults or adolescents
Adult Articulated Finger Clip	8000AA	greater than 30 kg (greater than 66 lbs)	Index, middle, or ring finger	<ul style="list-style-type: none"> • spot checking • short-term continuous use • for adults or adolescents
Pediatric Finger Clip	8000AP	greater than 8 kg (greater than 18 lbs)	Index, middle, or ring finger	<ul style="list-style-type: none"> • spot checking • short-term continuous use • for children or small adults
Ear Clip 1	8000Q	greater than 40 kg (greater than 88 lbs)	Ear lobe	<ul style="list-style-type: none"> • spot checking • short-term continuous use • stress testing
Reflectance 1	8000R	greater than 30 kg (greater than 66 lbs)	Middle of forehead	<ul style="list-style-type: none"> • continuous monitoring • stress testing
Adult Flex	8000J	greater than 20 kg (greater than 44 lbs)	Index, middle, or ring finger	<ul style="list-style-type: none"> • extended duration • motion situations • for adults
Infant Flex	8008J	2 to 20 kg (4 to 44 lbs)	Great toe	<ul style="list-style-type: none"> • extended duration • motion situations • for infants
Neonatal Flex	8001J	less than 2 kg (less than 4 lbs)	Lateral aspect of foot	<ul style="list-style-type: none"> • extended duration • motion situations • for neonates

Table V: Reusable Sensor Reference

¹ **Reflectance and Ear Clip** sensors generally do not perform as well as sensors located on the fingers or toes. They are not recommended for applications where the best possible SpO₂ accuracy is important.

➤ **NOTE:** *The reusable sensors must be secured to the patient using either Velcro®-type straps or medical tapes. NONIN's **Tape Strips** or **Hydrogel Tape Strips** are recommended to stabilize the reusable sensors.*

NONIN Disposable Sensor Reference

Sensor Description	NONIN Model Number	Patient Weight Range	Preferred Application	Application
Adult Flexi-Form®	9000A	greater than 30 kg (greater than 66 lbs)	Index, middle, or ring finger	<ul style="list-style-type: none"> • extended monitoring • self-adhesive • motion situations • for adults
Pediatric Flexi-Form®	9000P	10 to 40 kg (22 to 88 lbs)	Index, middle, or ring finger	<ul style="list-style-type: none"> • extended duration • self-adhesive • motion situations • for children
Infant Flexi-Form®	9000I	2 to 20 kg (4 to 44 lbs)	Great toe	<ul style="list-style-type: none"> • extended duration • self-adhesive • motion situations • for infants
Neonatal Flexi-Form®	9000N	less than 2 kg (less than 4 lbs)	Lateral aspect of foot	<ul style="list-style-type: none"> • extended duration • self-adhesive • motion situations • for neonates

Table VI: Disposable Sensor Reference

NONIN Sensor Accessory Reference

Sensor Accessory Description	NONIN Model Number	For Use With	Packaging Information
Tape Strips	8000T	Flex sensors	Packaged with 8000J
Hydrogel Tape Strips	8000TH	Flex sensors	Packaged with 8001J and 8008J
Sensor Holder	8000H	Reflectance sensor	Packaged with 8000R
3 foot extension cable	8500I	All sensors for 8500	---

Table VII: Sensor Accessory Reference

B. FINGER CLIP SENSORS

- **NOTE:** For the best results, secure the sensor cable independently from the sensor, preferably around the base of the finger. Make sure that the tape securing the cable does not restrict the blood flow.
- **NOTE:** If the patient's finger is too small or the sensor is not adjusted properly, light may bypass the tissue and result in SpO₂ inaccuracies.

1. Adult Finger Clip Sensor

Insert finger into the sensor until the end of the finger reaches the finger stop. If the patient has very long nails preventing the finger tip from reaching the stop, the nail must be trimmed. The thumb is specifically not recommended for use with the Adult Finger Clip Sensor.


Figure VIII: Finger Clip Sensor

2. Adult Articulated Finger Clip Sensor

Insert finger into the sensor until the end of the finger reaches the finger stop. Ensure that long fingernails are not interfering with proper finger position. The thumb is specifically not recommended for use with the Adult Articulated Finger Clip Sensor.


Figure IX: Adult Articulated Finger Clip Sensor

3. Pediatric Finger Clip Sensor

Insert finger into the sensor until the end of the finger reaches the finger stop. Ensure that long fingernails are not interfering with proper finger position. The thumb is specifically not recommended for use with the Pediatric Finger Clip Sensor.

C. FLEX SENSORS

- **NOTE:** For optimum light transmission, attach the sensor on the finger or toe. For best results, secure the cable independently from the sensor. Make sure that the tape securing the cable does not restrict the blood flow.
- **NOTE:** Ensure sensor placement is positioned such that no stray light bypasses the tissue or SpO₂ inaccuracies may result.

1. Adult Flex Sensor

The Adult Flex Sensor is designed for monitoring of adult and pediatric (greater than 20 kg) patients in which patient movement may occur and/or long term monitoring is required. Apply the double stick tape to the side of the sensor which is applied to the patient's finger.

Position the sensor on the top and bottom of the finger or toe as shown in Figure X. Place the light emitter portion on the finger/toe nail side and the detector on the side opposite of the nail. In all sensor placement applications, align the windows (detector and emitter portions of the sensor) over the tissue as illustrated in Figure X. Attach the sensor by wrapping the tape or wrap over the sensor assembly. Wrap the sensor snug, but not so tight as to restrict the blood flow.


Figure X: Flex Sensor


Figure XI: Flex Sensor with Tape Strain Relief

2. Infant Flex Sensor

- **NOTE:** *The preferred sensor site for infants is on the big toe. This placement provides increased light transmission compared to the foot placement. Refer to Table V for approximate size requirements for infant sensors.*

The Infant Flex Sensor is designed for continuous monitoring of infants where finger tip monitoring is impractical. It is designed for application on the large toe of infants. Apply the Hydrogel Tape Strips to the side of the sensor which is applied to the patient's toe.

Position Infant Flex Sensors on the big toe as illustrated in Figure XII. Make sure that the emitter portion of the sensor is aligned with the detector portion.

Attach the sensor by wrapping the tape or sensor wrap over the sensor assembly snugly, but not so tight as to restrict blood flow. For best results, secure cable independently from the sensor, preferably around patient's ankle or lower leg. Make sure that the tape securing the cable does not restrict the blood flow.

Figure XII: Infant Flex Sensor on Big Toe

3. Neonatal Flex Sensor

The Neonatal Flex Sensor is designed for continuous monitoring of neonates where finger tip monitoring is impractical. The Neonatal Flex Sensor is designed for application on the foot of infants. Apply the Hydrogel Tape Strips to the side of the sensor which is applied to the patient's foot.

Position Neonatal Flex Sensors on the foot as far forward near the toes as practical, ensuring sensor light does not pass between the toes as illustrated in Figure XIII. Make sure the emitter portion of the sensor is aligned with the detector portion.

Attach the sensor by wrapping the tape or sensor wrap over the sensor assembly snugly, but not so tight as to restrict blood flow. For best results, secure cable independently from the sensor, preferably around patient's ankle or lower leg. Make sure the tape securing the cable does not restrict blood flow.

Figure XIII: Neonatal Flex Sensor on Foot

D. EAR CLIP AND REFLECTANCE SENSORS

CAUTION

- The Ear Clip and Reflectance sensors are not recommended for pediatric or neonatal use. The accuracy of these sensors has not been established for pediatric or neonatal use.

➤ **NOTE:** *Reflectance and Ear Clip sensors generally do not perform as well as sensors located on the fingers or toes. They are not recommended for applications where the best possible SpO₂ accuracy is important.*

1. Ear Clip Sensor

The Ear Clip Sensor is designed for adults where fingertip monitoring is impractical. Rub the ear lobe vigorously for 5 seconds and then apply the Ear Clip Sensor to the lobe of the ear. Make sure the Ear Clip Sensor is positioned so the LED emitter and detector are completely covered by the earlobe. This ensures no stray light bypasses the earlobe, which can lead to SpO₂ inaccuracies.

2. Reflectance Sensor

The Reflectance Sensor is for use on well-vascularized skin surfaces. For adults, this will usually be the center of the forehead slightly above and between the eyebrows. The Sensor Holder provides the precise pressure this sensor needs against the skin.

Remove the backing from one side of the double-back tape and apply to the flange of the Sensor Holder. Then remove the back from the other side of the tape. Press the Reflectance Sensor into the foam with the sensor window out, and apply to the patient. Use additional tape to secure the sensor secondarily to the patient to avoid pulling or tipping the sensor.

E. ADULT AND PEDIATRIC DISPOSABLE SENSORS

WARNING

- Do not stretch the tape while applying the sensor. This may cause inaccurate readings or skin blisters.


The Flexi-Form[®] Sensors are designed for monitoring adult and pediatric patients as a single patient use sensor and are intended for use where moderate patient movement is expected or cross contamination is a concern.

The preferred application site is the index, middle, or ring finger. However, the thumb or toes may be used where the tissue thickness is between 5 and 21 millimeters. Other sites may not give acceptable results because of inadequate perfusion or inadequate light transmission. The application of these sensors is the same for either the adult or the pediatric patient. The difference is in the size of the sensor. For best results, secure the cable independently from the sensor. Make sure that the tape securing the cable does not restrict the blood flow.


- **NOTE:** *Ensure sensor placement is positioned such that no stray light bypasses the tissue or SpO₂ inaccuracies may result.*

Step 1: Grip the tab on the adhesive backing on the bottom of the sensor and carefully peel it off.

Step 2: Place the patient's finger or toe in the sensor with the nail side up and the tip of the finger or toe centered against the center line marked on the sensor as illustrated. This will ensure vertical alignment between the emitter and detector. Wrap the sensor around the finger. The finger nail should not be covered with tape during this step.

Step 3: Grip the tab on the adhesive backing of the top of the sensor and carefully peel it off.

Step 4: Fold the sensor's top over the top of the finger. Ensure that the detector and emitter are vertically aligned as illustrated by the dotted axis line.


Figure XIV: Adult/Pediatric Flexi-Form[®] Sensor

F. INFANT AND NEONATAL DISPOSABLE SENSORS


WARNING

- Do not stretch the tape while applying the sensor. This may cause inaccurate readings or skin blisters.
- Check application site frequently to determine circulation, positioning, and skin sensitivity of the patient. Each patient's sensitivity to NONIN sensors may vary depending on their medical status or the condition of their skin.


1. Infant Flexi-Form® Sensor

The preferred application site for the Infant Flexi-Form Sensor is the large toe of infants. Other sites may not give acceptable results because of inadequate perfusion or inadequate light transmission.

For best results, secure cable independently from the sensor, preferably around the patient's ankle or lower leg. Make sure the tape securing the cable does not restrict blood flow.

Step 1: Grip the sensor's adhesive backing and peel the adhesive off carefully, from top to bottom, as illustrated.

Step 2: Place the sensor on the large toe with the center line indicated on the sensor centered on the side of the large toe, taking care to have the detector centered on the bottom of the toe and the emitter located directly opposite on the toe nail, as illustrated.



Figure XV: Infant Flexi-Form[®] Sensor


Step 3: Wrap the tape around the toe taking care to ensure that the detector and emitter are vertically aligned as illustrated by the dotted axis line.

2. Neonatal Flexi-Form[®] Sensor


The preferred application site for the Neonatal Flexi-Form Sensor is the outstep of the foot, close to the toes for neonates. Other sites may not give acceptable results because of inadequate perfusion or inadequate light transmission. For best results, secure the cable independently from the sensor. Make sure that the tape securing the cable does not restrict the blood flow.

Step 1: Grip the sensor's adhesive backing and peel the adhesive cover off.

Step 2: Place the sensor on the outstep of the foot just behind the toes (as close to the toes as possible but not over the toes). Align the dashed line mark to the center of the outstep edge of the foot. The light emitter should be on the top and the detector on the bottom.

Step 3: Wrap the tape on top and bottom as illustrated.


Figure XVI: Neonatal Flexi-Form[®] Sensor

Step 4: Ensure that the detector and emitter are vertically aligned as illustrated by the dotted axis line.

G. CLEANING THE SENSORS

CAUTION

- Do not immerse the sensors in liquid to clean.
- Do not use caustic or abrasive cleaning agents.

Carefully clean the NONIN reusable sensors with an isopropyl alcohol wipe and ensure that all tape residue is removed. Allow enough time for the sensor to dry thoroughly before reusing. The sensors may be sterilized using ethylene oxide (EtO) (cold cycle).

H. SENSOR COMPATIBILITY

WARNING

- Use only NONIN manufactured sensors. These sensors are manufactured to meet the calibration requirements for NONIN Pulse Oximeters.

The 8500 is compatible with all NONIN manufactured sensors (except the Fiber Optic sensors). Refer to Tables V, VI, and VII for a list of the sensors, available sensor accessories, and their model numbers.

VII. THEORY OF OPERATION

The pulse oximeter shines light (red and infrared) through the tissue and detects the fluctuating signals caused by arterial blood pressure pulses. Well oxygenated blood is bright red, while poorly oxygenated blood is dark red. The pulse oximeter determines oxygen saturation from this color difference by measuring the ratio of absorbed red and infrared light. Since steady conditions (steady venous blood flow, skin thickness, bone, finger nails, etc.) do not cause fluctuations, they do not affect the saturation reading.

Mathematically:

$$SpO_2 = f \left[\frac{\ln \left(\frac{\text{min}}{\text{max}} \right)_{\text{red}}}{\ln \left(\frac{\text{min}}{\text{max}} \right)_{\text{infrared}}} \right]$$

Note that anything that affects the intensity of the light such as thick or colored skin will affect the maximum and minimum proportionally and thus the ratio max/min does not change. However, if too little light gets through, the pulse oximeter will not work and will signal sensor alarm.

CAUTION

- The NONIN 8500 Pulse Oximeter is calibrated to determine the percentage of arterial oxygen saturation of functional hemoglobin. Significant levels of dysfunctional hemoglobin such as carboxyhemoglobin or methemoglobin may affect the accuracy of the measurement.
- Cardiogreen and other intravascular dyes, depending on the concentration, may affect the accuracy of the SpO₂ measurement.

Pulse oximeters use two different wavelengths of light (colors) and thus have the ability to determine one component of blood. The NONIN 8500 Hand Held Pulse Oximeter is calibrated to closely approximate functional oxygen saturation values. The 8500 oxygen saturation values will closely approximate laboratory instrument fractional saturation values if the dysfunctional hemoglobin saturation levels are negligible. If the dysfunctional hemoglobin is carboxyhemoglobin or methemoglobin, then the difference between the oxygen saturation value displayed by the 8500 and the oxygen saturation values determined by the laboratory instrument will be greater as the dysfunctional hemoglobin levels rise approximately in accordance with the following formulas:

$$SpO_2 = O_2Hb + COHb + MetHb$$

$$SaO_2 = 100 \times O_2Hb / (100 - COHb - MetHb)$$

Where:

- SpO₂ = NONIN 8500 determined and displayed oxygen saturation in percent
- O₂Hb = Fractional oxyhemoglobin saturation in percent
- COHb = Carboxyhemoglobin saturation in percent
- MetHb = Methemoglobin saturation in percent
- SaO₂ = Functional oxygen saturation in percent

Example 1

O₂Hb = 96
 COHb = 0.5
 MetHb = 0.6

SpO₂ = 97
 SaO₂ = 97.07

Example 2

O₂Hb = 88
 COHb = 8
 MetHb = 2

SpO₂ = 98
 SaO₂ = 97.78

The mathematics are fixed in the pulse oximeter hardware and software. Thus no field calibrations are needed or are they possible. There are no adjustable parts within the pulse oximeter that affect the calibration.

The function, *f*, does depend on the color of the LEDs producing the light. These colors are fixed by the manufacturing processes and material. The sensors are checked for correct operation before shipping so no adjustment or calibration is needed or possible.

In summary, because the pulse oximeter does all critical computations in software and there are no critical parts to drift; no re-calibration is needed. Also, since the pulse oximeter reading does not depend on the absolute light intensity, but only on the fluctuations, then extremely thick skin, finger nails, colored skin, etc. do not affect the readings.

VIII. SPECIFICATIONS

- | | |
|---|---|
| 1. Oxygen Saturation Range (SpO₂) | 0 to 100% |
| 2. Pulse Rate Range | 18 to 300 Pulses Per Minute |
| 3. Displays | |
| Patient Indicator | Perfusion LED |
| Digital Displays | 3-digit 7-segment LEDs |
| 4. Measurement Wavelengths | |
| Red | 660 nanometers |
| Infrared | 910 nanometers |
| 5. Accuracy | |
| SpO ₂ | 70 - 95% ± 3 digits for neonates using infant or neonatal sensors |
| (± 1 Standard Deviation) ♦ | |
| | 70 - 100% ± 2 digits for adults using the Finger Clip sensor |
| | 70 - 100% ± 3 digits for adults using Flex or Reflectance sensors |
| | 70 - 100% ± 4 digits using Ear Clip sensor |
| | Below 70% is not specified for all sensors |
| Pulse Rate | ± 3% ± 1 digit |

♦ Standard Deviation is a statistical measure: up to 32% of the readings may fall outside these limits.

6. Temperature	
Operating	-20 to +50 °C
Non-operating	-30 to +50 °C
7. Humidity	
Operating	10 to 90% non-condensing
Non-operating	10 to 95% non-condensing
8. Power Requirements	6 AA alkaline batteries 100 hours typical operation (maximum display brightness) 200 hours typical operation (minimum display brightness)
9. Patient Isolation	> 12 MΩ
10. Leakage Current	Not applicable
11. Dimensions	3" wide x 6" high x 1" deep 8 cm x 15 cm x 2 cm
12. Weight	10 oz. (280 g) with batteries

IX. SERVICE

CAUTION

- 8500 Hand Held Pulse Oximeters are sensitive and must be repaired by knowledgeable and specially trained personnel only. Any sign or evidence of opening the system, field service by non-NONIN personnel, tampering, or any kind of misuse or abuse of the system, shall void the warranty in its entirety.

The solid state circuitry within the 8500 requires no periodic maintenance or calibration other than battery replacement.

NONIN does not recommend field repair of the 8500 Hand Held Pulse Oximeter. The circuit board in the 8500 series is a multi-layer board using traces 0.01" wide. Due to the very small trace size, extreme care must be used when replacing components to prevent permanent non-repairable damage to the circuit board. Most components are surface-mounted and require special hot air jet soldering and desoldering equipment. After any repairs are made, the pulse oximeter must be tested to ensure correct operation.

- **NOTE:** *All repair work on the 8500 should be done by trained NONIN personnel. For NONIN Technical Service contact:*

NONIN Medical, Inc.
2605 Fernbrook Lane North
Plymouth, Minnesota 55447 - 4755 USA

(612) 553-9968
(800) 356-8874 (USA and Canada Only)
FAX: (612) 553-7807


All non-warranty work shall be done according to NONIN standard rates and charges in effect at the time of delivery to NONIN. All repairs include a 24 hour burn-in and a complete retest of the pulse oximeter using factory test fixtures.

NONIN MEDICAL, INC.
PARTS LIST MODEL 8500 SERIES

ASSY ID	PART NUMBER	DESCRIPTION	U.M.	QTY
U1	0563-000	IC, ANALOG VOLTAGE REG (SO)	EA	1
U2,9	0235-103	IC, ANALOG	EA	2
U3	0562-000	IC, ANALOG	EA	1
U4	0679-000	IC, CONTROLLER	EA	1
U5	0427-003	IC, ANALOG	EA	1
U6	0430-000	IC, ANALOG	EA	1
U7	0427-103	IC, ANALOG	EA	1
U8	0429-001	IC, ANALOG	EA	1
Q1,2	0432-000	TRANSISTOR (SOT23)	EA	2
D1,2,3	0431-000	DIODE, SIGNAL (SOT23)	EA	3
R17,21	0423-121	RESISTOR, 1/8W,5%,100K Ω	EA	5
R23,28,31				
R3	0423-025	RESISTOR,1/8W,5%,10 Ω	EA	1
R4,8	0423-097	RESISTOR,1/8W,5%,10K Ω	EA	2
R5,6	0423-089	RESISTOR,1/8W,5%,4.7K Ω	EA	2
R7	0423-077	RESISTOR,1/8W,5%,1.5K Ω	EA	1
R9	0423-080	RESISTOR,1/8W,5%,2.0K Ω	EA	1
R10	0424-189	RESISTOR,1/8W,1%,909 Ω	EA	1
R11	0424-381	RESISTOR,1/8W,1%,90.0K Ω (1206)	EA	1
R12,JP1	0423-018	RESISTOR,1/8W,5%,5.1 Ω	EA	1
R13	0423-029	RESISTOR,1/8W,5%,15 Ω	EA	1
R14	0424-319	RESISTOR,1/8W,1%,20.5K Ω	EA	1
R15,16,25, R26,27	0424-481	RESISTOR,1/8W,1%,1.00M Ω	EA	5
R18,19,20	0423-061	RESISTOR,1/8W,5%,330 Ω	EA	3
R24	0423-464	RESISTOR,1/8W,1%,665K Ω	EA	1
R29	0423-169	RESISTOR,1/8W,5%,10M Ω	EA	1
C1,8,10,13C 18,19,20	0422-033	CAP, MONO CHIP,10%, 0.1 μ F(1206)	EA	7
C2,3,5	1902-106	CAP, TANTALUM CHIP,10%, 16V,10 μ F	EA	3
C4	0588-007	CAP,AL.ELECTROLYTIC,10 V, 470 μ F	EA	1
C7	0422-009	CAP, MONO CHIP,10%, 1000pF (1206)	EA	1
C9	0421-038	CAP, MONO CHIP,10%, 1200pF NPO (1206)	EA	1
C11,12	1902-684	CAP, TANTALUM CHIP,5%, 35V,3.3 μ F	EA	2
C14,15	0421-017	CAP, MONO CHIP,10%, 22pF (1206)	EA	2
C16,17	1902-079	CAP, TANTALUM CHIP,10%, 10V,47 μ F	EA	2

NONIN MEDICAL, INC.
PARTS LIST MODEL 8500 SERIES (CONTINUED)

ASSY ID	PART NUMBER	DESCRIPTION	U.M.	QTY
R1,2,22	0423-073	RESISTOR,1/8W,5%,1KΩ	EA	3
Y1	0568-001	CRYSTAL,4.0MHZ LOW PROFILE	EA	1
DS1,2,3,4, DS5,6	0569-006	LED, 7-SEGMENT,LOW CURRENT	EA	6
DS8	0577-002	LED, BI-COLOR	EA	1
J2	1725-000	CONNECTOR,9-PIN SUB-D (SM)	EA	1
J3	0587-004	CONNECTOR,SINGLE ROW 4-PIN	EA	1
J4,5,6,7,8, J9	0570-401	SOCKET, 10-PIN DIP (SM)	EA	6
J10	1299-002	CONNECTOR,RIGHT ANGLE LOCKING	EA	1
J11,12,13, J14,15,16	0602-201	SOCKET, 10-PIN DIP	EA	6
BC1	0593-001	CONTACT,POSITIVE BATTERY	EA	1
1	0829-000	PC BOARD, SM 8500	EA	1
2	0671-000	SUB-ASSY., NEGATIVE BATTERY CONTACT	EA	1
4	0148-000	HUMISEAL	A.R.	
5	1519-000	SCREW,#4-40X3/8 BUTTON SOCKET CAP	EA	2
6	1520-000	NUT,NYLON INSERT STOP	EA	2
	0498-250	SPACER, NYLON	EA	1
	0592-003	RED LENS 8500 & 8500M	EA	1
	0592-001	CASE, TOP 8500 & 8500M	EA	1
	0592-002	CASE, BOTTOM 8500 & 8500M	EA	1
	0592-005	COLLAR 8500 & 8500M	EA	1
	0670-000	BATTERY DOOR w/CONTACT 8500 & 8500M	EA	1
	1282-000	MEMBRANE SWITCH ASSY.	EA	1

X. WARRANTY

NONIN MEDICAL, INCORPORATED, (NONIN) warrants to the purchaser, for a period of three years from the date of delivery, each system exclusive of sensors, cables, and batteries. NONIN shall repair all systems found to be defective in accordance with this warranty, free of charge, for which NONIN has been notified by the purchaser by serial number that there is a defect, provided said notification occurs within the applicable warranty period. This warranty shall be the sole and exclusive remedy by the purchaser hereunder for any systems or accessories delivered to the purchaser which are found to be defective in any manner whether such remedies be in contract, tort or by law.

This warranty excludes cost of delivery to and from NONIN. All repaired units shall be received by the purchaser at NONIN's place of business. For any system or accessory sent to NONIN for warranty repair which is found to be within specification, the purchaser agrees to pay \$100.00.

These systems are sensitive and must be repaired by knowledgeable and specially trained personnel only. Accordingly, any sign or evidence of opening the system, field service by non-NONIN personnel, tampering, or any kind of misuse or abuse of the system, shall void the warranty in its entirety.

All non-warranty work shall be done according to NONIN standard rates and charges in effect at the time of delivery to NONIN.

DISCLAIMER/EXCLUSIVITY OF WARRANTY

THE EXPRESS WARRANTIES SET FORTH IN THIS MANUAL (SECTION X) ARE EXCLUSIVE AND NO OTHER WARRANTIES OF ANY KIND, WHETHER STATUTORY, WRITTEN, ORAL, OR IMPLIED INCLUDING WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE OR MERCHANTABILITY SHALL APPLY.

XI. ACCESSORIES

The following NONIN accessories function with the 8500 Hand Held Pulse Oximeter:

<u>Model</u>	<u>Description</u>
<i>Reusable Sensors</i>	
8000AA	Adult Articulated Finger Clip Sensor
8000AP	Pediatric Finger Clip Sensor
8000J	Adult Flex Sensor
8000K2	Finger Clip Sensor
8000Q	Ear Clip Sensor
8000R	Reflectance Sensor
8001J	Neonatal Flex Sensor
8008J	Infant Flex Sensor
<i>Disposable Sensors</i>	
9000A	Adult Finger Flexi-Form® Sensor, 10 per box
9000P	Pediatric Finger Flexi-Form® Sensor, 10 per box
9000I	Infant Toe Flexi-Form® Sensor, 10 per box
9000N	Neonatal Foot Flexi-Form® Sensor, 10 per box
9000D	Flexi-Form® Sensor Assortment Pack, 10 per box
<i>Accessories</i>	
8000H	Reflectance Sensor Holder
8000S	Simulator
8000T	Sensor Attachment Tape, clear, 100 per box
8000TH	Hydrogel Tape Strips, 25 per bag
8500CC-B	Carrying Case – Black
8500CC-Y	Carrying Case – Yellow
8500I	Patient Cable (3 foot extension)
8500MB	Mounting Bracket (Wall or Pole Mount System)
8500MC	Memory Transfer Cable (for use with PC)
8500RTC	Real-Time Transfer Cable (for use with PC)
8500RB	Rubber Bumper
8500YC	Real-Time Cable
8586PI	Printer Interface
9440	Finger Phantom Pulse Oximeter Tester

For further information about NONIN parts and accessories contact your local distributor.

If you do not have your local distributor information, call NONIN at (800) 356-8874, and ask for your local distributor's name and phone number.

XII. TROUBLESHOOTING GUIDE

Symptom	Possible Cause	Possible Solution
8500 will not power up	Batteries are completely depleted	Replace all six batteries of the 8500
	Incorrect battery installation	Verify battery sequence. Refer to Figure I: Replacing Batteries in the 8500
	Battery door is missing metal contact	Replace the battery door
Numeric displays are blinking at once per second	Battery voltage is low	Replace all six batteries of the 8500
	Incorrect battery installation	Verify battery sequence. Refer to Figure I: Replacing Batteries in the 8500
Minus sign appears in SpO ₂ display	Sensor fault exists. Sensor may have become dislodged from 8500 or from the patient	Verify sensor is connected to the 8500 and the patient correctly; Try a new sensor if condition persists
Segments of SpO ₂ or ♥ display are missing	Defective LED displays	Displayed values may not be reliable; discontinue use of 8500
Displayed pulse rate does not correlate to pulse rate displayed on ECG monitor	Excessive motion at sensor site may be prohibiting the 8500 from acquiring a consistent pulse signal	Eliminate or reduce cause of motion artifact or reposition sensor to new sensor site where motion is not present
	Patient may have an arrhythmia resulting in some heart beats that do not yield a perfusion signal at sensor sight	Examine the patient: Condition may persist even though both monitors are functioning properly if patient's arrhythmia persists
	Non-NONIN sensor is being used	Replace sensor with a NONIN sensor
	ECG monitor may not be functioning properly	Examine the patient: Replace ECG monitor or refer to operator's manual for ECG monitor

TROUBLESHOOTING GUIDE (continued)

Symptom	Possible Cause	Possible Solution
Erratic ♥ display and/or yellow perfusion LED during concurrent use of electrosurgical equipment (ESU)	ESU may be interfering with oximeter performance	Examine the patient: move oximeter, cables, and sensor as far away from ESU as possible <u>or</u> refer to the ESU operator's manual
Perfusion is blinking yellow with each pulse	Perfusion signal at sensor site is marginal	Examine the patient: reposition sensor <u>or</u> select alternate sensor site
Unable to obtain green perfusion	Low patient pulse strength	Reposition sensor on patient
	Sensor site poorly perfused	
	Sensor not correctly positioned	
	Sensor attached too tightly or tape or other items are restricting perfusion at sensor site	Reapply sensor, select alternate sensor site, or remove restrictive material from sensor site
	Circulation reduced due to excess pressure between the sensor and a hard surface	Allow sensor and finger to rest comfortably on surface
	Excessive ambient light	Reduce ambient light
	Excessive patient motion	Reduce patient motion
	Sensor applied to a polished fingernail	Remove fingernail polish
	Interference from: <ul style="list-style-type: none"> • arterial catheter • blood pressure cuff • electrosurgical procedure • infusion line 	Reduce or eliminate interference

TROUBLESHOOTING GUIDE (continued)

Symptom	Possible Cause	Possible Solution
Perfusion LED is blinking red and SpO ₂ and ♥ displays show dashes	Inadequate perfusion signal at sensor site	Examine the patient: reposition sensor <u>or</u> select alternate sensor site
	Excessive motion at sensor site may be prohibiting 8500 from acquiring a consistent pulse signal	Eliminate or reduce cause of motion artifact or reposition sensor to sensor site where motion is not present
Printer not printing out after changing Prn mode	Printer mode did not get updated internally	Reset pulse oximeter by turning it off and then on

If any of these solutions do not correct the problem with your 8500, please contact NONIN Technical Service at (800) 356-8874.

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