TECHNICAL ASSISTANCE

SUCTION UNITS SERIES VORTECO

The documents of quality registration filled in our Factory at the moment of final checking and testing of the units are kept in the record office of the Company ALSA APPARECCHI MEDICALI S.R.L. and can be requested, if considered as necessary (specify the serial number of the unit).

REQUIRED INSTRUMENTS/EQUIPMENT AND REFERENCE RULES

1. Vacuum tester (tol. $\pm 5\%$) for depression from 0 to 250 mmHg (1 mmHg = 0.0013 atm or bar, 0.1333 kpa, 1.359 cmH₂O)
The instrument can be even for a vacuum higher than 250 mmHg, but in this case it must be able to measure a low vacuum.
2. Free air flow tester (tol. $\pm 5\%$) for a free air flow from 1 to 10/15 l/min (or 600/1000 l/h).
3. Multimeter (tol. $\pm 1\%$).
4. Tester for L.F. leakage currents measurements and protection hearth conducting resistance (the instruments must be able to make the measurements according to the International Safety Rules for L.F. leakage currents and protection hearth conducting).
5. There are a lot of Firms that manufacture this kind of items, for example BENDER, METRON, BIOTECH and so on.
6. The reference Rules are: IEC 601-1 (EN 60601-1) and ISO 10079-1.

GENERAL DESCRIPTION

A medical suction unit is an apparatus with 3 general sections:

- a. Electrical vacuum pump (AS 130 is provided also with an electronic board which allows to fix the activation times “active times/passive times”)
- b. Vacuum system (vacuum tubing, internal connections and vacuum controls)
- c. Protection filters (when used), internal tubing and external tubing, including also connections/bottles/overflow devices. This system, so, must be airtight (any leakage prevent to reach the desired maximum vacuum/flow) and free from obstructions (an obstruction allows to the unit to reach the maximum vacuum, but it can decrease the aspiration flow).

ORDINARY MAINTENANCE

For each model the ordinary maintenance is detailed in the specific instructions to use the unit.

MAINTENANCE IN CASE OF MALFUNCTIONING

There are 3 possibilities:

- 1) The unit doesn't work because of a breakdown of the electric part.
- 2) The unit hasn't the intended of the aspiration capacity (the vacuum data are those indicated in the instructions to use the units).
According to the International Rules, a suction unit is with low vacuum/low flow when it reaches a vacuum of at least 60 kpa and a free air flow of at least 20 l/min measured to the joint of aspiration of the bottle.
- 3) The unit doesn't work because suction group is partially obstructed or completely broken.

-
- 1) Using the multimeter, check:
 - Supply and mains fuses (see the technical data on the label)
 - General switch and working switches (The correct working of the electronic board - only mod AS130)
 - Internal electric connections and suction group pump
 - 2) Using the vacuum tester and/or the flow tester, check that the circuit has not any leak or obstruction (you have to check every single component if necessary excluding it from the circuit). A good method is to get through directly to the suction group to verify its suction capacity, and so, after having reconnected it to the circuit, make the same test on all the points of the circuit.
Remember that with an obstruction (for example an obstructed bacterial filter) it seems the unit works regularly, but really it has a very low suction capacity (or even nothing).
 - 3) If the problem is a presumed obstruction of the suction group, act in the following way:
Unscrew the screws that block the suction part, disassembling all the internal parts, cleaning them and if necessary changing them. Reassembling the group, please note that the screws must be uniformly clasped.
 - 4) If the suction group is completely broken, it must be disconnected from the suction circuit, disassembled, changed and reconnected to the suction circuit.
 - 5) At the end of the technical intervention (connecting to the suction tubing first the vacuum tester, then the flow tester) check the normal working of the unit (intended vacuum and or flow) and so, in order to be sure that in disassembling and reassembling conditions of electrical risk haven't been created, make a test on the L.F. leakage and of the protection hearth conducting using a specific tester (the procedures of execution of this test must be those specified in the relative International Safety Rules ISO 10079-1).

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