
Aesculap Power Systems


-
- GB** Service Manual
Oscillating depth saw handpiece GB125R/Micro pendulum
saw GB129
 - D** Service-Manual
Oszillierendes Tiefensägehandstück GB125R/Mikro-Pendel-
säge GB129

GB129R

GB125

Oscillating depth saw handpiece GB125R/Micro pendulum saw GB129

Copyright information

Copyright© 2006
 Aesculap AG & Co. KG®
 All rights reserved 09/06

This service information is copyrighted. This service information must not be copied or reproduced by other means either in whole or in part without the express permission of Aesculap AG & Co. KG.

This manual is intended for informational purposes only. Ownership of this manual alone does not constitute or imply authorization to service the product.

The technical information, illustrations and dimensions contained in this manual are non-binding. No claims may be made on the basis of the information contained herein. We reserve the right to make improvements without altering this documentation. We reserve the right to make technical changes without prior notice.

Main address for servicing:

Aesculap Technischer Service (ATS)
 Am Aesculap-Platz
 78532 Tuttlingen / Germany
 Phone: +49 7461 95-2700
 Fax: +49 7461 16-2887
 E-mail: ats@aesculap.de

Other service addresses can be obtained through the address indicated above.

Manufacturer's liability

We expressly point out that we can accept responsibility for any effects on the safety, reliability, and performance of our medical products, if and only if:

- any assembly, extensions, readjustments, modifications, or repairs are performed by technically experienced, knowledgeable and trained personnel and
- the medical products are used as set forth in the instructions for use.

It is possible to learn to service the products through an appropriate course of instruction given on the premises of Aesculap about the relevant medical products. To arrange for such a course, contact Aesculap Technical Service (ATS).

To ensure that your Aesculap warranty remains valid, we recommend that only spare parts that have been factory tested by Aesculap be used for repairs. Spare parts, as well as the relevant tools, can be ordered from Aesculap Technical Service (ATS). Any unauthorized opening and/or alterations of the medical product by third parties lead to the exclusion of our liability, as far as a fault is attributable to such unauthorized opening and/or alteration of the product. Aesculap cannot accept responsibility for the use of unsuitable spare parts, tools or devices. After repairs, dropping, severe damage or misuse, the product should be inspected by a qualified person.

This manual refers to the product as it was when the manual was prepared. Technical changes may be made at any time, particularly in software.

General advisory

This service manual contains illustrations and explanations. Each explanation covers the following: modes of operation, spare parts list, assembly, disassembly, functionality test, and the section tools/devices (if applicable).

For further information about a specific product (e.g., troubleshooting list, accessories) consult the relevant user instructions.

Literature	Art. no.
Instructions for use of GB125	TA005327
Instructions for use of GB129	TA005320

It is crucial that extreme cleanliness be maintained when the products described below are repaired. All products must undergo a functionality test after being repaired. If, in the course of a repair, a labeled component is replaced, the labeling must be transferred to the spare part. Either the original label should be used or the spare part bearing the relevant label should be ordered. To this end, be sure to include the applicable information when ordering.

Unless otherwise indicated, all medical products should be inspected annually.

General information about motors and handpieces

It is critically important that ball bearings be inspected annually. Even the slightest defect or soiling can lead to overheating during use, which would render the product unusable. If for the repair of, e.g., a motor a ball bearing has to be replaced, it is best in most cases to replace all other ball bearings as well.

If it is necessary to heat up adhesive joints during disassembly, a hot-air blower is the appropriate tool. The smallest possible nozzle should be used so as to avoid damaging other components. When gluing parts together, make certain that the surfaces to be joined are absolutely clean and free of grease. For cleaning we recommend Loctite Rapid Cleaner 7063 (to be ordered under WS. no. 520001750). Only use adhesives listed in this manual. Follow the glue manufacturer's recommendations.

Lubricate components with the recommended Aesculap product only.

Following disassembly, clean all components thoroughly and inspect them to ensure that they are undamaged. Do not install any component of whose status you are uncertain. Products must always be test-run for several minutes after being repaired.

Contents

1. Safe handling	3
2. Tools, auxiliary materials, supplies	3
2.1 Tools	3
2.2 Auxiliary materials	3
2.3 Supplies	3
3. Expendable parts/replacement parts	4
3.1 GB125R	4
3.2 GB129	4
4. Basic preparations and maintenance	4
4.1 Preparations	4
4.2 Maintenance	4
5. Troubleshooting list	5
5.1 Oscillating depth saw handpiece GB125R	5
5.2 Micro pendulum saw GB129	6
6. Disassembling	8
6.1 Disassembling the saw head (GB125R)	8
6.2 Removing the drive unit from the housing (GB125R)	8
6.3 Disassembling the drive unit (GB125R)	8
6.4 Disassembling the saw head (GB129)	8
6.5 Removing the drive unit from the housing (GB129)	8
6.6 Disassembling the drive unit (GB129)	8
7. Assembling	9
7.1 Preliminary work	9
7.2 Assembling the drive unit (GB125R)	9
7.3 Assembling the housing (GB125R)	9
7.4 Assembling the saw head (GB125R)	9
7.5 Assembling the drive unit (GB129)	9
7.6 Installing the drive unit in the housing (GB129)	9
7.7 Assembling the saw head (GB129)	9
8. Postmaintenance	10
8.1 Surface inspection	10
8.2 Function tests	10
9. Parts list	10
9.1 GB125R	10
9.2 GB129	11

1. Safe handling

- Operate oscillating depth saw handpiece GB125 and micro pendulum saw GB129 with authentic Aesculap accessories only.
- To avoid damage to the products, do not knock depth saw handpiece GB125 and micro pendulum saw GB129 against hard objects.

2. Tools, auxiliary materials, supplies

2.1 Tools

Designation	Art. no.
Wrench WW3 DIN 911	TA004801
Box wrench, angled, WW14	TA010100
Press-in device, saw handpiece	W-81000180
Press-in punch for axial shaft seal TA004527	W-90000005
Press-in punch for bearing ring seal TA004527	W-90000026
Punch	W-90000026 no. 72
Auxiliary tool	W-90000030 no. 66
Wrench for bearing bush GB125211	W-99010115

2.2 Auxiliary materials

Designation	Art. no.
Instructions for use of GB125R	TA005327
Instructions for use of GB129	TA005320

2.3 Supplies

Designation	WS. no.
Barrierta L25DL	537001568
Special grease MI-setral-FKR 2	537001602
Loctite 275	516001814
Loctite 572	560001817
Loctite 640	560001820
Eloxa Cleaner	JG601

Oscillating depth saw handpiece GB125R/Micro pendulum saw GB129

3. Expendable parts/Replacement parts

3.1 GB125R

Designation	Art. no.
Saw blade axle	GB125213
Rocker arm	GB125302
Calotte bearing	GB125306
Eccentric shaft	GB125307
Nut with pressure plate	GB125804
Housing, complete	GB125805
Locking washer	TA003131
Simmerring	TA003257
Locking washer (x2)	TA003904
Radial ball bearing	TA003907
Radial ball bearing (x2)	TA003916
Radial ball bearing	TA004110
Shaft seal, axial	TA004527
Bearing ring seal	TA005307

3.2 GB129

Designation	Art. no.
Saw blade axle	GB125213
Calotte bearing	GB125306
Nut with pressure plate	GB125804
Radial ball bearing (x2)	TA003907
Radial ball bearing	TA003983
Radial ball bearing	TA004110
Bearing ring seal	TA005307

4. Basic preparations and maintenance

4.1 Preliminary work

Function tests

- Switch on the oscillating depth saw handpiece GB125R/micro pendulum saw GB129.
- Make certain that the motor runs smoothly.
- Couple saw blade GC650R.
- Test the power of oscillating depth saw handpiece GB125R/micro pendulum saw GB129 by sawing into a block of wood.

4.2 Maintenance

In order to ensure reliable operation, Aesculap recommends maintenance after every 300 reprocessing cycles or at least once a year.

- Lubricate the gears and all ball bearings.

5. Troubleshooting list

5.1 Oscillating depth saw handpiece GB125R

Malfunction	Finding	Cause	Remedy
Loud running noise	Major slack of saw blade axle GB125213 with the drive shaft fixated	Worn gearing	Replace ball bearing TA004110/ rocker arm GB125302
	Major slack of eccentric shaft GB125307	Worn or defective bearing	Replace bearing of eccentric shaft GB125307
Excessive heat-up	Hot saw blade axle GB125213	Bearing ring seal TA005307 overheated	Replace bearing ring seal TA005307 Sterilize oscillating depth saw handpiece GB125R Allow bearing ring seal TA005307 to run in
	Hot oscillating depth saw handpiece GB125R	Oscillating depth saw handpiece GB125R overloaded	Move oscillating depth saw handpiece GB125R back and forth during sawing Allow oscillating depth saw handpiece GB125R to cut free Do not cant the saw blade
Saw blade fails to move	Major slack of saw blade axle GB125213 with the drive shaft fixated	Worn saw blade coupling	Replace ball bearing TA004110/ rocker arm GB125302
	Driver pivot broken	Drive unit defective	Repair drive unit (ELAN-E/EC-flexible shaft/pneumatic motor GA200)
Deficient cutting performance	Traces of wear on the saw blades	Blunt saw blades	Insert a new tool Check instrument sets
Oscillating depth saw handpiece GB125R cannot be coupled to drive unit	Dirty or damaged drive pivots/ drive hole of the saw	Dirty or damaged drive pivots/ drive hole of the saw	Clean or smoothen the surfaces of the pivots/hole
	Bent, damaged or broken eccentric shaft GB125307	Drive pinion damaged	Repair drive unit (ELAN-E/EC-flexible shaft/pneumatic motor GA200)
Oscillating depth saw handpiece GB125R cannot be fixed on the drive unit	Dirty or damaged drive pivots/ drive hole of the saw	Worn drive	Replace drive unit/housing adapter

Oscillating depth saw handpiece GB125R/Micro pendulum saw GB129

Malfunction	Finding	Cause	Remedy
Tool cannot be coupled	Dirty tool/ adapter chuck	Dirt	Couple a new tool Take off nut GB125804 and clean saw blade axle GB125213
	Defective or broken driver (star)	Deformation	Replace saw blade axle GB125213
	Wrong tool/ coupling shape mismatch	Wrong tool	Check art. no. and compare to the following information: Catalog no. O17599, Burrs&Blades

5.2 Micro pendulum saw GB129

Malfunction	Finding	Cause	Remedy
Loud running noise	Major slack of saw blade axle GB125213 with the drive shaft fixated	Worn gearing	Replace ball bearing TA004110/ rocker arm GB125302
	Major slack of eccentric shaft GB125307	Worn or defective bearing	Replace bearing of eccentric shaft GB125301
Excessive heat-up	Hot saw blade axle GB125213	Bearing ring seal TA005307 overheated	Replace bearing ring seal TA005307 Sterilize micro pendulum saw GB129 Allow bearing ring seal TA005307 to run in
	Hot micro pendulum saw GB129	Micro pendulum saw GB129 overloaded	Move micro pendulum saw GB129 back and forth during sawing Allow micro pendulum saw GB129 to cut free Do not cant the saw blade
Saw blade fails to move	Grooved pin TA005387 broken	Overloaded/worn grooved pin	Replace grooved pin TA005387/ coupling piece TA006659
	Broken drive pin	Drive unit defective	Repair drive unit (ELAN-E/EC-flexible shaft/pneumatic motor GA200)
Deficient cutting performance	Traces of wear on the saw blades	Blunt saw blades	Insert a new tool Check instrument sets
Micro pendulum saw GB129 cannot be coupled to the drive unit	Dirty or damaged drive pivots/ drive hole of the saw	Dirty or damaged drive pivots/ drive hole of the saw	Clean or smoothen the surfaces of pivots/hole
Micro pendulum saw GB129 cannot be fixed on the drive unit	Locking on the drive pivot (INTRA) defective	Worn drive	Replace locking mechanism

Malfunction	Finding	Cause	Remedy
Tool cannot be coupled	Dirty tool/adapter chuck	Dirt	Take off nut GB125804 and clean saw blade axle GB125213 Couple a new tool
	Defective or broken driver (star)	Deformation	Replace saw blade axle GB125213
	Wrong tool/ coupling shape mismatch	Wrong tool	Check art. no. and compare to the following information: Catalog no. 017599, Burrs&Blades

Oscillating depth saw handpiece GB125R/Micro pendulum saw GB129

6. Disassembling

6.1 Disassembling the saw head (GB125R)

- Heat up housing GB125805 if necessary and unscrew complete saw head.
 - Unscrew nut with pressure plate GB125804 from the saw head.
 - Clamp in the saw head between clamping jaws.
 - Unscrew gear cover GB125305 from the saw head.
 - Unscrew flat-headed screw GB125304 on the saw head.
 - Expel saw blade axle GB125213 from the saw head, using an appropriate punch.
 - Heat up bearing bush GB125211 if necessary and unscrew it, using wrench W-99010115.
 - Remove rocker arm GB125302.
 - Remove ball bearing TA004110 from bearing bush GB125211.
 - Remove bearing ring seal TA005307 from bearing bush GB125211.
 - Remove ball bearing TA003907 from the saw head.
- This completes the disassembly of the saw head.

6.2 Removing the drive unit from the housing (GB125R)

- Remove Seeger circlip ring TA003909 from housing GB125805.
 - Push out, from back to front, the complete drive unit from housing GB125805.
- This completes the removing of the drive unit from the housing.

6.3 Disassembling the drive unit (GB125R)

- Remove 2 locking washer TA003904 and locking washer TA003131.
 - Remove ball bearing TA003907, 2 radial ball bearings TA003916 and calotte bearing GB125306 from eccentric shaft GB125307.
 - Remove axial shaft seal TA004527 from housing GB125805.
 - Unscrew clamping screw GB115801 from housing GB125805.
 - Unscrew cylinder head screw TA004018 from the rear of housing GB125805.
 - Remove spring washer TA004692.
 - Slide off sleeve GB131202 from housing GB125805.
 - Remove spring bracket GB125312 from housing GB125805.
- This completes the disassembly of the drive unit.

6.4 Disassembling the saw head (GB129)

- Heat up housing GB125800 if necessary and unscrew complete saw head.
 - Unscrew nut with pressure plate GB125804 from the saw head.
 - Clamp in the saw head between clamping jaws.
 - Unscrew gear cover GB125305 from the saw head.
 - Unscrew flat-headed screw GB125304 from the saw head.
 - Expel saw blade axle GB125213 from the saw head, using punch W-90000026 no. 72.
 - Heat up bearing bush GB125211 if necessary and unscrew it, using wrench W-99010115.
 - Remove rocker arm GB125302.
 - Remove ball bearing TA004110 from bearing bush GB125211.
 - Remove bearing ring seal TA005307 from bearing bush GB125211.
 - Remove ball bearing TA003907 from the saw head.
- This completes the disassembly of the saw head.

6.5 Removing the drive unit from the housing (GB129)

- Remove lock ring TA003976 from housing GB129800.
 - Push out, from back to front, the complete drive unit from housing GB125800.
- This completes the removing of the drive unit from the housing.

6.6 Disassembling the drive unit (GB129)

- Expel pin TA005387 from eccentric shaft GB129301.
 - Remove coupling piece TA006659 and coupling spring TA005385 from eccentric shaft GB129301.
 - Remove locking washers TA003131 and TA003222 from eccentric shaft GB129301.
 - Press off calotte bearing GB125306 and ball bearings TA003907 and TA003983 from eccentric shaft GB129301.
- This completes the disassembly of the drive unit.

7. Assembling

7.1 Preliminary work

- Rethread every thread, using a taper.
- Clean all components, remove all silicone and Loctite residues, and degrease.
- Clean saw head surfaces with Eloxal Cleaner JG601.
- Lubricate open ball bearings TA003907, TA003983, TA003916, TA004110 with special grease MI-setral-FKR 2, WS.-no. 537001602. Externally lubricate calotte bearing GB125306 with Barrierta L25DL, WS.-no. 537001568.

7.2 Assembling the drive unit (GB125R)

- Press calotte bearing GB125306 and 2 radial ball bearings TA003916 onto eccentric shaft GB125307.
- Install locking washers TA003131 and TA003222 on eccentric shaft GB129301.
This completes the assembly of the drive unit.

7.3 Assembling the housing (GB125R)

- Install spring bracket GB125312 in housing GB125805. Take care to position it correctly.
- Install sleeve GB131202 on housing GB125805.
- Insert spring washer TA004692 in sleeve GB131202.
- Screw cylinder head screw TA004018 into the rear part of housing GB125805.
- Screw in clamping screw GB115801 in housing GB125805.
- Press axial shaft seal TA004527 into housing GB125805, using press-in punch W-90000005.
- Slide in, from front to back, the complete drive unit into housing GB125805.
- Insert Seeger circlip ring TA003909 in housing GB125805.
This completes the assembling of the housing.

7.4 Assembling the saw head (GB125R)

- Press bearing ring seal TA005307 into bearing bush GB125211, using press-in punch W-90000026.
- Press ball bearing TA004110 into bearing bush GB125211.
- Press ball bearing TA003907 into the saw head.
- Apply Loctite 640 on bearing bush GB125211 and screw the bush onto the saw head, using wrench W-99010115.
- Insert rocker arm GB125302 in the saw head in such a way that the chamfer at the slotted hole is oriented upwards.
- Position rocker arm GB125302 in such a way that the slotted hole is aligned with the drill hole.

- From the outside, insert saw blade axle GB125213 in the saw head in such a way that it sits flush in rocker arm GB125302.
- Press saw blade axle GB125213 into rocker arm GB125302, using press-in device W-81000180.
- Apply Loctite 572 on flat-headed screw GB125304 and insert the screw in the saw head.
- Apply Loctite 572 on gear cover GB125305 and screw the gear cover to the saw head.
- Screw nut with pressure plate GB125804 into the saw head.
- Apply Loctite 572 on the saw head and screw the saw head into housing GB125805.
This completes the mounting of the saw head.

7.5 Assembling the drive unit (GB129)

- Press calotte bearing GB125306 and ball bearings TA003907 and TA003983 onto eccentric shaft GB129301.
- Install locking washers TA003131 and TA003222 on eccentric shaft GB129301.
- Mount coupling spring TA005385 on eccentric shaft GB129301.
- Mount coupling piece TA006659 on eccentric shaft GB129301 in such a way that the slotted hole is aligned with the drill hole.
- Impel pin TA005387 in eccentric shaft GB129301 and check for smooth movement.
This completes the assembly of the drive unit.

7.6 Installing the drive unit in the housing (GB129)

- Slide the drive unit into housing GB129800 as far as it will go.
- Install lock ring TA003976 in housing GB129800.
This completes the installation of the drive unit in the housing.

7.7 Assembling the saw head (GB129)

- Press bearing ring seal TA005307 into bearing bush GB125211, using press-in punch W-90000026.
- Press ball bearing TA004110 into bearing bush GB125211.
- Press ball bearing TA003907 into the saw head.
- Apply Loctite 640 on bearing bush GB125211 and screw the bush onto the saw head, using wrench W-99010115.
- Insert rocker arm GB125302 in the saw head in such a way that the chamfer at the slotted hole is oriented upwards.
- Position rocker arm GB125302 in such a way that the slotted hole is aligned with the drill hole.
- From the outside, insert saw blade axle GB125213 in the saw head in such a way that it sits flush in rocker arm GB125302.
- Press saw blade axle GB125213 into rocker arm GB125302, using press-in device W-81000180.
- Apply Loctite 640 on flat-headed screw GB125304 and insert the screw in the saw head.
- Screw gear cover GB125305 onto the saw head.

Oscillating depth saw handpiece GB125R/Micro pendulum saw GB129

- Screw nut with pressure plate GB125804 into the saw head.
 - Apply Loctite 275 on the saw head and screw the saw head into housing GB125800.
- This completes the mounting of the saw head.

8. Postmaintenance

8.1 Surface inspection

- Check surfaces for:
 - dents
 - stains
 - sharp edges
 - glue residues

8.2 Function tests

- Switch on the oscillating depth saw handpiece GB125R/micro pendulum saw GB129.
- Couple saw blade GC650R.
- Test the power of oscillating depth saw handpiece GB125R/micro pendulum saw GB129 by sawing into a block of wood.
- Check oscillating depth saw handpiece GB125 R/micro pendulum saw GB129 for heat-up.

After 3 minutes of running, the temperature must not exceed 41 °C.
- Check oscillating depth saw handpiece GB125R/micro pendulum saw GB129 for untypical running noise, which would indicate the presence of defective bearings, gears or connections.

9. Parts list

9.1 GB125R

Art. no.	Designation
GB115305	Lever lower part
GB115801	Clamping screw
GB125110	Lever upper part (blank)
GB125211	Bearing bush
GB125212	Nut M 7 x 0.75
GB125213	Saw blade axle
GB125214	Pressure plate
GB125215	Saucer spring
GB125302	Rocker arm
GB125304	Flat-headed screw DIN 921 M 4.0 x 4.0
GB125305	Gear cover
GB125306	Calotte bearing
GB125307	Eccentric shaft
GB125308	Housing
GB125310	Lever upper part
GB125312	Spring bracket
GB125400	Gearhead
GB125804	Nut with pressure plate
GB125805	Housing, complete
GB131202	Sleeve
TA003131	Locking washer DIN 6799
TA003904	Locking washer (x2)
TA003907	Ball bearing, radial 7 x 19 x 6, stainless steel
TA003909	Seeger circlip ring 22 x 1 DIN 472
TA003916	Radial ball bearing 8 x 22 x 7 stainless ISOFLEX (x2)
TA004018	Cylinder head screw DIN 84 M 2.6 x 6.0
TA004110	Ball bearing, radial 10 x 22 x 6
TA004527	Axial shaft seal
TA004692	Spring washer DIN 7980
TA004797	Cylinder head screw DIN 912 M 4.0 x 10.0
TA005307	Bearing ring seal

9.2 GB129

Art. no.	Designation
GB125211	Bearing bush
GB125211	Bearing
GB125212	Nut M 7 x 0.75
GB125212	Nut M 7 x 0.75
GB125213	Saw blade axle
GB125214	Pressure plate
GB125215	Saucer spring
GB125302	Rocker arm
GB125304	Flat-headed screw DIN 921 M 4.0 x 4.0
GB125305	Gear cover
GB125306	Calotte bearing
GB125804	Nut with pressure plate
GB125804	Nut with pressure plate
GB129200	Housing
GB129301	Eccentric shaft
GB129302	Gearhead
GB129800	Housing, complete
GB601202	Parallel pin 2.0 R 6 x 7.0
TA003131	Locking washer DIN 6799
TA003222	Locking washer DIN 6799
TA003905	Retaining ring DIN 472 19.0 x 1.0
TA003907	Ball bearing, radial 7 x 19 x 6 (x2)
TA003983	Ball bearing, radial 5 x 16 x 5
TA004110	Ball bearing, radial 10 x 22 x 6
TA005307	Bearing ring seal
TA005385	Coupling spring for GB130
TA005387	Grooved pin DIN 1473 1.2 x 5 NR
TA006659	Coupling piece

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