-ebro-

MEASUREMENTS FOR LIFE

ebro® Electronic GmbH & Co. KG


Peringerstraße 10 • D-85055 Ingolstadt

Tel. +49 (0) 841 - 9 54 78-0

Fax +49 (0) 841 - 9 54 78-80

Internet: www.ebro.com

E-Mail: info@ebro.com

-ebro-

MEASUREMENTS FOR LIFE

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
Operating Instructions

Mode d'emploi

Datenlogger
Data Logger
Enregistreur
Serie EBI 10

-ebro-[®]

Datenlogger
Data Logger
Enregistreur

Serie EBI 10

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Table of contents

| | |
|---|-----------|
| Overview..... | 5 |
| The Interface (write/read device)..... | 7 |
| Safety Instructions | 9 |
| Tips for the data logger with flexible probe..... | 11 |
| Unpacking/What is Included in the Delivery | 13 |
| Operating the Data Logger..... | 15 |
| Data Transfer Mode..... | 15 |
| Swapping Out the Spacer Ring..... | 19 |
| Replacing the Battery | 23 |
| What to Do, if... | 29 |
| Cleaning and Maintenance..... | 33 |
| Calibration service | 33 |
| Disposal..... | 33 |
| Technical Data | 35 |
| Approvals | 37 |
| Declaration of Conformity | |

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
Overview

EBI 10 series data loggers are small, handy devices used to measure and record temperature and/or pressure. They are powered by a lithium battery and programmed from a PC.

To program the data logger, you will need a suitable interface (write/read device, which can be supplied as an accessory) along with the appropriate read/write program on the PC (e. g. "Winlog.pro", also available as an accessory). The interface is connected to the PC via a universal serial bus (USB) port.

EBI 10 series data loggers are available in a range of models to suit a variety of uses:

- channels 1, 2, 3 or 4
- temperature or temperature/pressure
- temperature control versions flexible or stiff
Ø 1.5, 3.0 or 6.0 mm
Lengths from 10 to 500 mm
Position radial or axial
Available with sharp or blunt point
- pressure sensor connection
Standard
M10 internal thread for hose, flange or glass adapter
Luer lock

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
These different models are inserted in the interface and used for programming the data loggers and extracting the recorded data.

For more information about how to program the data logger and to extract the data, please refer to the user instructions for your read/write program.

The data logger does not have an on/off switch. It starts operating as soon as it is supplied with power by the battery and continues until the battery runs out of power.

The Interface

To program the data logger you will need a suitable interface (available as an accessory). This interface is connected to the PC via the universal serial bus (USB), with a special cable that can be supplied as an accessory, and either powered from there or from its own power supply unit (which can also be supplied as an accessory).

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Safety instructions

Do not use the data logger in areas in which there is a danger of explosions. Danger of death!

Do not introduce the data logger into an environment that is hotter than 150°C. The lithium battery may explode!

Do not grasp the logger with bare hands immediately after use! Danger of burning!

You must only use the data logger with batteries supplied by Ebro. These are the only batteries that have adequate heat resistance!

Do not subject the device to any microwave radiation. The lithium battery may explode!

Data loggers emit wireless signals in the 2.4 GHz frequency range. These may cause significant interference to other technical devices! For this reason, only use data loggers in areas where these wireless emissions do not represent a risk!

Data loggers are currently only approved for use in Europe. Approvals for other, non-European states have been applied for.

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
Tips for the data logger with flexible probe

! Do not bend these parts of the probe: the first 10 mm that protrude from its housing, or the last 50 mm at its tip. Bending it in these areas could result in damage to the probe. The edge of the housing will place too much mechanical stress on the probe if you bend it too close to the housing. The sensor is located in the last 50 mm before the tip of the probe.

Do not use pincers to bend the probe. To bend the probe, hold its lower part securely with two fingers and then use your other hand to bend it into the required position.

The radius of the bend must be at least 10 mm. A smaller radius will lead to material fatigue and cause the probe to break.

Do not remove data loggers that are supplied in a protective silicone box from that box. Removing and manipulating the flexible probe will cause it unnecessary stress and might damage it.

To fetch data from, or to program, data loggers that are supplied in a protective silicone box, simply pull down the lower part of this box. After you have fetched data from, or programmed the data logger, close the protective silicone box.

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Operating the Data Logger

The data logger does not have an on/off switch. It starts operating as soon as it is supplied with power by the battery and continues until the battery runs out of power. However, you must program it before you can use it. To do this, use a specially designed software program (e.g. "Winlog.pro") on a PC or laptop to which a special read/write device (interface) is connected.

To program an EBI 10 data logger, place it in one of the interface ports so that the plastic side faces upwards. Then use the programming software to program it. For more information, please refer to the user instructions for this software.

Always allow the data logger to cool down slowly from high to lower temperatures. Fast cooling may cause the data logger to leak.

Data Transfer Mode

Wireless mode:

The logger sends its data to the interface remotely at a programmed time interval.

Interface mode:

The logger is inserted in the interface for programming and extracting the recorded data.

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If you want measurement values to be transferred, you ensure that there is an uninterrupted wireless connection between the data logger and the interface (for more details, please refer to the interface's user instructions) and that the evaluation program on the computer can accept this constant stream of data (for more details, please refer to the software's operating instructions).

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



(3)

Swapping Out the Spacer Ring

EBI 10 series data loggers are fitted either with a simple metal spacer ring or with a spacer ring that has a fixing loop (Figure 1). If necessary, you can swap the existing spacer ring for a replacement. The appropriate spacer ring set is available as an accessory.

To exchange the spacer ring:

- Use the opening tool to open the data logger (see Figure 2).
- Remove the elastic, silicone O-ring.
- Remove the metal spacer ring.
- Clean the slot in which the silicone O-ring was located (Figure 3).
- Insert the new metal spacer ring.

! Take a new silicone O-ring, cover it with a thin layer of grease and place it in the slot.

A used silicone O-ring must always be disposed of. Under no circumstances should it be reused! This is because it is deformed when the two parts of the housing are screwed together. If it is used again it will not make a proper seal!

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


Take care not to overstretch the elastic ring and not to damage it!

- Screw the two parts of the logger housing together so they are hand-tight.
- Wipe away any grease that may be squeezed out.

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Replacing the Battery

Once a battery is spent, it must be replaced with another one.

 **Whilst changing the battery they must protect themselves against electrostatic discharges using suitable precautionary measures such as an earthed wristband! Electrostatic discharges may destroy the device!**

To swap out a battery:

- Use the opening tool to open the housing of the data logger (see Figure 1).
- Remove the battery retaining clip (Figure 2). (Loosen the screw, lift the retaining clip and remove it carefully.)
- Pull the battery plug out of the circuit board (Figure 3).
- Remove the empty battery.
- Place a new battery in the battery compartment in exactly the same position as the one you have just removed.

Dispose of the used battery in an environmentally friendly way at a suitable re-cycling point!

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



- Insert the retaining clip carefully in the two slots so that the red battery wire lies below it, so that it cannot be crushed or damaged (Figure 4).
- Plug the battery plug carefully into the socket on the circuit board.
- Tighten the retaining screw again. (Figure 5).

! Take great care throughout this procedure and do not use excessive force at any point. The plug is polarized and can only fit into the socket on the circuit board in one direction.

To the right of the battery compartment, you will see a LED. This LED will flash as soon as the circuit board is receiving power from the battery again (Figure 6).

- Remove the elastic silicone O-ring from the red plastic part of the logger housing.
- Cleaning the slot in which the O-ring was located (Figure 7).
- Take a new silicone O-ring, cover it with a thin layer of grease and place it in the slot.

! The layer of grease should cover the O-ring completely (so that it glistens) but not leave lumps of residue on it.

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 Take care not to overstretch the elastic ring and not to damage it!

- Screw the two parts of the logger housing together so they are hand-tight.
- Wipe away any grease that may be squeezed out.
- Use the programming software to reset the battery counter. For more information, please refer to your programming software's operating instructions.

No data is lost when you change a battery or if the battery fails. However, if a longer power failure occurs, the built-in clock will no longer show the correct time. In this situation, extract the data from the data logger and then reprogram it. This resynchronizes the date and time.

Recording continues immediately after you insert a new battery.

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
**What to Do, if...**

... there is no wireless transmission to the data logger:

- Check that the data logger's wireless mode is active. To do this, place the data logger close to the interface and call the appropriate function in the programming software.
- If necessary, use the programming software to switch the data logger to wireless mode.
- If the data logger does not react to this, check if the battery is flat and needs to be replaced (see below).

... the data logger cannot be contacted from the interface:

- Check that the interface has been connected correctly and that the evaluation software has recognized it.
- Open the data logger (see also the Figures in the "Replacing the Battery" section on page 22 ff.).
- Check to see whether condensation has penetrated the data logger housing.
- Either wipe the condensation from the housing or allow it to air dry by leaving the data logger housing open.

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
- Even if no condensation could be found in the data logger housing: pull the battery plug from the socket and wait 30 seconds.
 - Then reconnect the battery and check that the LED to the right of the battery compartment is flashing. If the LED does not flash, the battery may be flat and needs to be replaced (see “Replacing the Battery” on page 22 ff.)
 - As a test, insert a fresh battery in the socket on the circuit board.
- If the LED still does not flash (despite the fresh battery), the problem lies with the data logger.**
- Before you close the data logger housing again, replace the silicone O-ring (see page 24/26).

! A used silicone O-ring must always be disposed of. Under no circumstances should it be reused! This is because it is deformed when the two parts of the housing are screwed together. If it is used again it will not make a proper seal!

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Cleaning and Maintenance

Use a damp cloth to clean the data logger's external surfaces.

Do not use a solvent, such as acetone, because this damages the plastic.

Calibration service

The device takes extremely accurate measurements. To ensure it continues to work to this high standard, it should be calibrated annually. For this reason, we offer a calibration service.

- To take advantage of this service, please complete the accompanying service card and send it to ebro Electronic (address on the reverse of these instructions), Fax +49 841 95478-80.

Every year we shall then recall your device for calibration.

We shall return the fully calibrated device to you within one week.

Disposal

If the device becomes no longer fit for purpose, it must be disposed of in a suitable, environmentally-friendly manner.

Do not, under any circumstances, simply dispose of it in domestic garbage. You must return it to the manufacturer.

Dispose of the battery at the appropriate recycling point.

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Technical Data

| | |
|---------------------------------------|---|
| Measurands | |
| Pressure: | (kPa), 1 pressure channel (10 mbar = 1 kPa) |
| Temperature: | (°C), 1-4 temperature channels |
| Measuring ranges | |
| Pressure: | 0 to 4000 mbar |
| Temperature: | |
| Series EBI 10-TP | 0 to +150 °C |
| Series EBI 10-T probe length >= 100mm | -85 to +400°C |
| Series EBI 10-T remaining | -85 to +150°C |
| Series EBI 10-T model 101 | -85 to +85°C |
| Accuracy | |
| Pressure: | ±10mbar (50 to 150 mbar) ±10mbar (2050 to 2250 mbar) ±10mbar (3000 to 3250 mbar) ±15mbar (for the remaining) |
| Temperature: | ±0,5 °C (-85 to -40°C) ±0,2 °C (-40 to 0°C) ±0,1 °C (0 to +140°C) ±0,2 °C (+140 to +250°C) ±0,5 °C (+250 to +400°C) |
| Measurement signal resolution | |
| Pressure: | 1 mbar (100 Pa) |
| Temperature: | 0,025 °C |
| Data memory | 100,000 measurement values |
| Sensor | PT1000, piezo-resistive pressure sensor |
| Interface | wireless 2.4 GHz / IEEE 802.15.4 |
| Operating temperature | |
| Record mode: | -85 bis +150 °C |
| Radio mode: | -30 bis +150 °C |
| Storage | -40 to +125 °C |
| Protection type | IP68/NEMA 6P |
| Battery: | 3.6 V lithium battery exchangeable |
| Battery life | up to 2 years depending on use |
| Dimensions (D x H) | 46 X 24 mm |
| Housing material | stainless steel (V4A), PEEK |
| Weight | approximately 45 g |

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Approvals

The conformity certificate confirms that this product fulfills CE guidelines.

The product **EBI 10** fulfills the guidelines in accordance with EN 12830.

Suitability:

S (Storage), T (Transport)

Location:

A, B, C, D

Accuracy classification:

1, for measurement of the air temperature and of the internal product temperature.

Measurement range:

-40 ... +150°C

In accordance with EN 12830, this device is subject to regular (yearly) inspections as per EN 13486

Wireless:

This device has been tested against ETSI EN:

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This device has been tested against CFR 47 Part 15 FCC Rules and Regulations

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This device has been tested against RSS Rules RSS Rules and Regulations

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Operation is subject to the following two conditions:
 (1) This device may not cause harmful interference, and
 (2) This device must accept any interference received including interference that may cause undesired operation.

This device has been designed to operate with an internal antenna.
 Other Antennas are strictly prohibited for use with this device.

Changes or modifications of the device could void the user's authority to operate the equipment.

(Exception: Replacing the battery as described in the user manual.

Users manual or instruction manual may downloaded at www.ebro.com over the internet.

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



Table of Contents

| | |
|--|----|
| Overview..... | 5 |
| Safety Instructions | 7 |
| Unpacking/What is Included in the Delivery | 9 |
| Operating the Interface..... | 11 |
| Connecting the Interface to the PC | 11 |
| Connecting a suitable antenna | 13 |
| If required: connecting a power supply unit..... | 15 |
| Installing the Interface Driver | 15 |
| Meaning of the colored status display lights | 17 |
| Interface operating types (modes)..... | 19 |
| Operation in standard mode..... | 21 |
| Operation in wireless mode..... | 25 |
| What to Do, if | 29 |
| ...the status display lights remain dark? .. | 29 |
| ...the PC responds with “interface not found”? | 29 |
| ...no connection is established to the data logger? | 29 |
| Cleaning and Maintenance | 33 |
| Disposal..... | 33 |
| Technical Data | 35 |
| Approvals | 35 |
| Conformity declaration | 36 |

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Overview

The EBI IF-100/150/200/300 interface is a programming and reading device for EBI 10/15/11 series that also acts as a wireless receiver for EBI 10 series data loggers. These data loggers are programmed, and with EBI 10 data is extracted from them, by wireless in the 2.4 GHz frequency range.

The interface is connected by a cable to a PC, from which it is controlled, using suitable software (e. g. “Winlog.pro”), and usually also supplied with power. If the PC cannot supply enough power for the interface, you can purchase a separate power supply unit as an accessory.

To program the data logger and extract the data recorded on it in standard mode (see Page 21), the device must be placed in the interface ports. To receive data in wireless mode (see Page 25) a powerful antenna must first be connected to the interface.

Structure of the interface:

- 1..... Antenna connection
- 2..... USB connection
- 3..... Power supply unit connection
- 4..... Data logger ports

For more information about how to program the data logger, and how to extract the data, please refer to the user instructions for your read/write program.

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Safety Instructions

Do not use the device in areas in which there is a danger of explosions. Danger of death!

Be sure to install the device in such a way so that the disconnection of power supply unit easily is accessible.

Do not introduce the device into an environment that is hotter than 60°C!

To avoid damaging the device, protect it from direct contact with water. Please note the device's protection class (IP 20).

The device emits wireless signals in the 2.4 GHz frequency range. These may cause significant interference to other technical devices! For this reason, only use the device in areas where these wireless emissions do not represent a risk!

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Unpacking/What is Included in the Delivery

Check the device packaging and contents to ensure they are complete and undamaged. You must also check that the contents of the packaging match your order.

What is Included in the Delivery

The delivery includes these parts:

- EBI IF-100/150/200/300 interface
- USB cable
- rod antenna 5 dbi (Part Name AL 111), only EBI IF-100/200 and rod antenna 2 dbi (Part Name AL 110) with EBI IF-150!
- these operating instructions
- CD-ROM with driver software

The delivery may also include various accessories:

| <u>Description</u> | <u>Name</u> |
|---------------------|------------------|
| Wire antenna | Part Name AL 112 |
| Evaluation software | Winlog.x |
| Power supply unit | Part Name AL 120 |

If you have any reason for complaint, do not hesitate to contact us. Our contact details are printed on the reverse of the envelope.

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Operating the Interface

Bringing the interface into operation involves four steps:

- 1.) Connecting the interface to the PC with a USB cable
- 2.) Connecting a suitable antenna
- 3.) Connecting a power supply unit (only if the PC cannot supply enough power)
- 4.) Installing the interface driver on the PC

Connecting the Interface to the PC

The interface is supplied along with a USB cable. Insert its plug (type A, with rectangular cross-section) into the PC. The square cross-section plug (type B) fits the socket on the interface.

- Insert the ends of the cable plug into the appropriate socket on the PC and on the interface. We recommend that you use a high-power USB socket on the PC, if present.

Take great care throughout this procedure and do not use excessive force at any point. The plugs are polarized and only fit into the socket in one particular alignment.

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
**Connecting a suitable antenna**

The interface is supplied along with two antennae: one rod antenna and one wire antenna.

The rod antenna is designed for normal operations (in both standard mode and wireless mode). The wire antenna is designed to record measurements in processes that have critical wireless transmission attributes (e.g. in steam sterilizers, autoclaves, ovens etc.). Each antenna has a thread that fits into the interface's antenna socket. Alternatively, they can be screwed into the interface socket.

 **Before you change an antenna, pull the power supply cable or the USB cable out of the interface so that the status display lights go out.**

To ensure correct operation, you must only connect antennae supplied by ebro to the interface.

- Either screw the rod antenna, or the wire antenna's plug, into the interface socket.

Take great care throughout this procedure and do not use excessive force at any point.

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





If required, connect a power supply unit

If the PC cannot supply enough power along the USB cable to operate the interface, you must use a separate power supply to run the interface (this is available as an accessory). In this case, a message appears on screen immediately after you plug the interface into the PC's USB port.

You will only need to attach a power supply unit if this message appears.

Only use the power supply unit supplied by ebro to operate the interface. This is the only unit that has the required technical functionality.

- If necessary, insert the mains plug into the mains socket and then insert the low voltage plug into the socket on the interface.

As soon as the interface is supplied with power, all status display indicators light up blue.

Installing the Interface Driver

The interface is supplied together with a CD-ROM that holds the required driver software. This CD-ROM also contains precise instructions about how to install the software.

- The instructions on the CD-ROM will tell you how to install the driver software on your PC.

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Meaning of the colored status display lights

The status display lights on the data logger ports come on as soon as the interface is supplied with electricity. These lights change color, depending on which function the interface is performing. The colors have these meanings:

| | |
|-----------------|--|
| blue | interface ready for use |
| yellow | program/extract data from the data logger in the particular port |
| green | programming/data extraction successful |
| red | programming/data extraction failed |
| flashing violet | wireless mode |

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
Interface operating types (modes)

The interface can be used in two different operating modes:

- in standard mode and
- in wireless mode.

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In wireless mode the interface is used as a wireless receiver. In this situation, the interface continuously receives signals from data loggers that are also running in wireless mode and sending measurement data as soon as they record the measurements. When operating in wireless mode the interface can receive and process radio signals from up to 40 data loggers at the same time.

When the interface is running in wireless mode, the violet status display lights flash on and off.

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



Operation in standard mode

In standard mode the interface is used either to program the data logger or to extract all the data collected in the data logger at once. In this way a different data logger can be programmed, or data can be extracted from it, on each one of the interface's ports.

To program a data logger, or extract data from it, follow these steps:

- If the data logger is enclosed in a protective silicon box, first remove the lower part of this box.
- Place the data logger, with its plastic side facing downwards, in a free port on the Interface.
- On the PC that is running evaluation software, start the programming or extraction process and follow the on-screen instructions.

If a port's status display light is yellow, this means the programming or extraction process is running on that particular data logger.

If the port's status display light is green, this means the programming or extraction process was successful. If this process was not successful, this light goes to red.

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

- Once you have completed the programming or extraction process, replace the lower part of the data logger's protective silicon box.
- Repeat the entire process as required for the next data logger.

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

Operation in wireless mode

In wireless mode, the interface operates as a radio receiver. It continuously receives and processes signals from data loggers that are also running in wireless mode and sending their measurement data as soon as they record the measurements. When operating in wireless mode the interface can receive and process radio signals from up to 40 data loggers at the same time.

Radio waves behave like light, and metal and concrete objects, in particular, prevent the waves from being diffused. To ensure clear reception, the space between the data loggers and the interface antenna should be kept as clear as possible of objects of this type. The best reception is achieved when there is a direct line of sight between the data logger and the interface antenna.

The data logger should also be sited fairly close to the interface to ensure the radio signals can be received clearly and easily. The data logger has a **maximum range** of approximately 30 m from the interface if there are no obstacles in the way.

Antennae receive the clearest radio signals from directions that are vertical to their axis. For this reason we recommend that the rod antenna is sited vertically.

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 To ensure correct operation, and to avoid interference, we recommend that only antennas supplied by ebro are connected to the interface.

The wire antenna supplied with the interface is designed for devices which are to be used in wireless mode.

Many devices have metal cases which obstruct the data logger's radio signals. Despite this, these signals can still be sent to the interface if the wire antenna is fitted into the device.

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What to do, if...**...the status display lights remain dark?**

In this case, it is clear that no power is reaching the interface.

- Check that the PC is switched on, that the USB cable plugs are inserted correctly and whether the USB socket on the PC is a high-power socket (see Page 15).
- If an external power supply unit is being used, check that its plug is correctly connected to the socket and the interface.

...the PC responds with “Interface not found”

- Check that the interface’s status display lights are on. If not, carry out the same checks as you would for if the status display remains dark.
- Check the USB connection between the interface and the PC.
- Check you have installed the correct drivers (see the instructions for installing the drivers on the CD-ROM supplied with the interface).

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...no connection is established to the data logger?

- Check whether the data loggers are running in wireless mode and you have selected data extraction from the data loggers in standard mode. To communicate with each other, the data logger and interface must be running in the same operating mode.

Please refer to the programming software instructions on the PC to find out how you check the data logger's operating mode.

- Check that the interface is using the correct radio channel for its attempts to establish contact (refer to the software's programming log for the data logger from which data is to be extracted). The data logger and interface will not be able to contact each other if they are running on different radio channels.

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Technical Data

| | |
|-----------------------------|--|
| Power supply | USB (high power USB 500 mA) or external power supply unit Part Name AL 120 |
| Wireless interface | 2.4 GHz IEEE 802.15.4 |
| USB port | Type B 500 MA |
| Antenna port | RP-SMA |
| Operating temperature | -30 to +60°C |
| Storage | -40 to +60°C |
| Protection type | IP20 |
| Dimensions Plan (L x W x H) | 150 x 180 x 45 mm |
| Housing material | ABS |
| Weight | approximately 358 g |

Approvals

The EBI IF-100/200 Interface has these approvals:

Wireless:

| | |
|--|-----------|
| This device has been tested against ETSI EN: | |
| This device complies with ETSI EN: | |
| EN 301 489-1 V1.6.1 | : 2005 |
| EN 301 489-3 V1.4.1 | : 2002 |
| EN 301 489-17 V1.2.1 | : 2002 |
| ETSI EN 300 328 V1.7.1 | : 10/2002 |
| ETSI EN 300 330-2 V1.3.1 | : 04/2006 |
| ETSI EN 300 330-1 V1.5.1 | : 04/2006 |

This device has been tested against CFR 47 Part 15 FCC Rules and Regulations
CFR 47 Part 15 Subpart B and C May 2007

FCC ID: VQ5-EBIIF

This device complies with part 15 of the FCC Rules.

This device has been tested against RSS Rules RSS Rules and Regulations
RSS-Gen Issue 2, RSS-102 Issue 2, RSS-210 Issue 7

IC: 7412A-EBIIF

This device complies with RSS Rules.

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TDS PAL

SPECIFICATIONS

| Model | Resolution | Range |
|-------|-------------|---------------|
| TDS-1 | 1 x 10 ppm | 10 to 1990 |
| TDS-2 | 1 x 100 ppm | 100 to 10,000 |
| TDS-3 | 1 x 10 uS | 10 to 1990 |
| TDS-4 | 1 x 100 uS | 100 to 19900 |

Accuracy : $\pm 2\%$ (full scale)

Battery : 4 x 1.5V
button cell
(Alkaline A76 or
equivalent)

Battery life : Approx. 100 hours


Operating
temperature : 0° to 50°C

Temperature
compensation : Automatic

Size (LxWxH) : 147 x 32 x 15mm

Weight : Approx. 50 gm

PRODUCT LAYOUT

OPERATION

- Remove protective cap from bottom (See product layout)
- To switch on the unit, slide the 'ON/OFF' switch located on top of the tester to 'ON'.
- Dip tester into sample solution up to the immersion level as shown.
- Stir gently to remove bubbles and wait for 2 minutes so as to allow full temperature compensation before readout.
- Multiply reading according to the display resolution.
- Always rinse with distilled water, blot dry before and after each test.
- Switch off the tester and replace protective cap before storing away.

NOTES:

In the presence of certain radio transmitters, this product may produce erroneous readings. If this occurs then measurements should be repeated at another location.

CALIBRATIONS

It is recommended to check meter every six months for accuracy.

1. To perform this check, dip the tester in Standard Calibration solution.

Order Code : 2020 for TDS-1
Order Code : 2100 for TDS-2
Order Code : 2010 for TDS-3
Order Code : 2100 for TDS-4


2. Stir gently to remove bubbles and wait for reading to stabilize.
3. With the use of a small screw driver, locate the "Calibration trimmer" at the back of the tester (see product layout) and trim it to read as below:

TDS-1 = 138
TDS-2 = 64
TDS-3 = 141
TDS-4 = 129

4. Rinse tester with distilled water.
5. Calibration is now complete.

PRECAUTIONS IN HANDLING


Do not touch, rub or scratch the sensor. It might lose its sensitivity and accuracy.

Do not drop unit. The tester is fragile not made to withstand impact. It will crack or break under impact, which is beyond repair.


Do not drop unit in water or dip unit beyond the immersion level. The unit is not watertight and is beyond repair if water get in unit.

Do not store unit under high temperature or direct sunlight. This will shorten the life span of the product.

Do not clean unit with thinner or solvents. This will damage the unit. Use only a damp cloth to clean unit if needed.

MAINTENANCE

- To improve performance of tester, clean the electrode periodically every 6 or 12 months by rinsing it in alcohol for few minutes.
- If the displays become faint or disappear, this could be due to run-down batteries. Pull out the battery case from top of unit (see layout) and replace all four batteries accordingly.
- Check batteries voltage should be 4.0 volts and above. Reading accuracy will deteriorate with weak batteries.
- If unit is not used or stored for a long period of time, it is a good practice to remove the batteries in case of battery leaks and cause damage to the contacts.

APPLICATION

- Home water purifier Maintenance
- Pools and Spas maintenance
- Photograph development chemical test
- Waste water, sewage treatment
- Aquariums , Aqua-culture
- Agriculture, gardening, Hydroponic
- Clinics & health care
- Cooling & Boilers water
- Chemical industry
- Food industry
- Wood and pulp industry

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(GB)

5 TAKING A MEASUREMENT

- 1 Remove protective cap from bottom (See product layout.)
- 2 Press the 'ON-OFF / Hold' button once to switch on.
- 3 Shake the tester in the same way you would use a mercury thermometer, before each test.
- 4 Dip tester into the sample solution up to the immersion level. Shake the sensor area in solution to remove bubbles and wait for about a minute or when reading stabilizes.
- 5 Pressing the ON-OFF/Read-Hold button once will freeze the reading. Pressing again will release it for another readout.
- 6 If the sensor is dry, a slow response will result. Dip the sensor area in a cup of water or preferably Buffer 7 calibration solution for 1-2 hours before testing again.
- 7 Always rinse the sensor area with water and blot it dry before and after each test.

AUTO-LOCK MEASUREMENT

- 1 The Auto-Lock feature is user friendly, allowing the tester to take a stable reading without guesswork.
- 2 To enable auto-lock measurement, press both the ON-OFF/Read-Hold button and the CAL button and hold till a blinking small 'A' sign appears on the right side of display.
- 3 From now on, everytime you dip the tester into the solution, press the ONOFF button once, the display digit will blink while sensing for a stable reading. Once the reading is stabilised, the display digit will stop. You can now take the reading.
- 4 Press the ON-OFF button again for another reading.
- 5 To disable the Auto-Lock feature, press both the ON-OFF button and the CAL button till the small 'A' sign disappears. Now you can resume standard measurement as before.
- 6 To Switch off the tester, press the 'ON-OFF/Read-Hold' button for 3 seconds.
- 7 Replace protective cap before storing away.

3 From now on, everytime you dip the tester into the solution, press the ONOFF button once, the display digit will blink while sensing for a stable reading. Once the reading is stabilised, the display digit will stop. You can now take the reading.

- 4 Press the ON-OFF button again for another reading.
- 5 To disable the Auto-Lock feature, press both the ON-OFF button and the CAL button till the small 'A' sign disappears. Now you can resume standard measurement as before.
- 6 To Switch off the tester, press the 'ON-OFF/Read-Hold' button for 3 seconds.
- 7 Replace protective cap before storing away.

6 MAINTENANCE

- When the battery symbol appears on the display, this indicates a low battery and only 2 hours of continuous use remain. Though the unit may continue to function, the accuracy of the unit will be affected beyond the 2 hours.
- To change batteries, turn to the back of the unit and lift the latch on the battery case from the top of unit, (see layout) slide the clip upward to open. Replace all four batteries according to polarity. Push cap back till latch clicks in place.
- If the unit is not used or stored for a long period of time, the sensor will become dry. This will result in a slow response to a stable reading. Soaking the sensor area in a cup of tap water or preferably Buffer 7 solution for 1-2 hours will restore sensitivity to the sensor.
- Keep in mind that all pH sensors age with time and usage. Therefore, recalibration is necessary to maintain accurate reading.
- Note that the pH sensor has a limited life span of about 365 tests. When the unit fails to calibrate or responds very slowly, it means that the unit should be replaced. It is not possible to repair broken sensor, defective or expired unit.

7 SPECIFICATION

| | |
|------------------------|------------------------|
| RANGE: | 0.0 to 14.0 pH |
| RESOLUTION: | 0.1 pH |
| ACCURACY: | ±0.2 pH |
| BATTERY: | 4 x 1.5V A76/LR44 |
| BATTERY LIFE: | approx. 150 hours |
| AUTO SHUT-OFF: | approx. 15 minutes |
| PRODUCT LIFE: | 365 tests |
| CASING MATERIAL: | ABS plastic |
| OPERATING TEMPERATURE: | 0° to 50°C/32 to 122°F |
| SIZE (LxWxH): | 170 x 32 x 15mm |
| WEIGHT: | approx. 70 gm |

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



2 INSTALLING BATTERY CAP

The unit is shipped with the battery cap open. Close the battery cap by pressing the cap on a hard surface until the latch clicks, indicating a secure lock.

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2 REPLACING THE BATTERIES

1. Lift the latch on the battery cap with a small flat object, such as a small screwdriver. Use your thumb to push the cap forwards as shown. Firmly grasp the cap and carefully separate it from the pH Meter body.

2. Ensure that the batteries are orientated correctly (correct orientation is marked inside the battery cap).

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

UNDERSTAND YOUR PRODUCT

Change battery according to polarity
 Lift latch here then push cap up to open


Installing battery cap:

The unit is shipped with the battery cap open. Close the battery cap by pressing cap on table till the latch „clicks“ for a secure lock.


How to open battery cap:

- 1 Use a mini-screwdriver to lift latch till it pops up. DO NOT PULL latch out completely.
- 2 Use the thumb to push cap forward as shown. Turn over to the front and pull cap out completely.


PRECAUTIONS IN HANDLING


Do not touch, rub or scratch the sensor. It is very delicate and might break or lose its sensitivity.


Do not submerge the unit under water. Although the unit is splash-proof and water resistant, it cannot come under high pressure underwater and is beyond repair if water gets into the unit. If it is dropped into water, retrieve it immediately and wipe dry with a cloth.

Do not store the unit under high temperature or direct sunlight. This will shorten the life span of the unit.

Do not store unit without the protective cap. Chemical in the unit will expire faster, thus shortening its life span.

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Applications

- Industries
- Process water
- Washer & Disinfectant rinsing water
- Laboratories
- Brewery process
- Food production
- Pool, Aquarium

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-ebra-

pH Digital Tester

Laboratories
 Industries
 Pool
 Aquarium
 Spa

water resistant - float on water - drop shock - simple to use

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Attributes:

- Rugged
- Replaceable battery
- Waterproof
- Calibrateable

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Technical data:

| | |
|-----------------------|-------------------------------|
| Measurement range | 0.0 to 14.0 pH |
| Resolution | 0.1 pH |
| Accuracy | ± 0.2 pH |
| Battery | 4 x A76/LR44 |
| Operating temperature | 0°C to +50°C 32°F to 122°F |
| Housing material | ABS plastic |
| Weight | approx. 70g |

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pH is the most frequent test parameter in any applications involving water. Each industry has its own regulations for the monitoring and controlling of pH value. Maintaining pH value is important in the quality control of various production processes.

Countless applications include leather production, chemical processing, food production, water production, clinical inspection and wastewater disposal inspection. Only proper control and inspection of pH value can assure higher production yields and quality products.

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ebro Electronic GmbH & Co. KG