

# Rafale EV 30

Mobile Radiographic Unit

## Service Manual

**APELEM**

D M S G R O U P

Parc Scientifique G.Besse, 175 allée Von Neumann ♦ 30035 Nimes cedex 1

+33 (0)4 66 29 09 07 ♦ FAX +33 (0)4 66 29 71 23

e-mail [export@apelem.com](mailto:export@apelem.com)

Site web: [www.apelem.com](http://www.apelem.com)

Original version of this manual is in Italian language therefore, for further information and clarification, please refer to the Italian one.

## **DOCUMENT REVISION STATUS**

Document code: 04.02.007.005

Revision 1 – Valid from 01<sup>st</sup> April 2011

- Doc. 10551 - Valid from 01<sup>st</sup> April 2011
- Doc. 7819 - Valid from 01<sup>st</sup> July 2004

# Index

# Service Manual

---

INDEX

## General index (doc. 10551)

Valid from 01<sup>st</sup> April 2011

### Section 0 - Index

#### Section 1 – Installation

**Page**

1	Introduction	6
2	Unpacking	7
	2.1 Packing with cellophane	7
	2.2 Packing with carton	9
	2.3 Packing with wood crate	10
3	Electrical connections	11
4	Acceptance test	11
	4.1 Checking cables and connectors	11
	4.2 Mechanical checks	11
	4.3 Electrical checks	12

#### Section 2 – Maintenance

1	Routine maintenance	16
	1.1 General recommendations	16
	1.2 Frequent checks and inspections	17
	1.3 General checks and inspections	17
	1.4 Cleaning and disinfection	18
2	Special maintenance	19
	2.1 Dismantling the Unit Casings	19
	2.2 Extraordinary Mechanical Maintenance	21
	2.3 Mechanical Adjustments	22
	2.3.1 Parking brake	22
	2.3.2 Anti-tilting System	22
	2.3.3 Arm-X-ray tube head Balancing System	24

Section 0

2.3.4	Locking the Arm in the Transport Position	26
2.3.5	X-ray tube head-Collimator Assembly Rotation	26
2.3.6	Centring the X-ray tube head – Collimator Assembly	28
2.4	Electrical troubleshooting	29
2.5	Thermomagnetic Circuit-breaker	30
2.6	Replacement of Components	31
2.7	List potentiometers	32
2.8	List leds	33
2.9	List fuses	34
3	Spare parts	34

## Section 3 – X-ray settings

1	Introduction	36
1.1	Dip switch	37
2	Rotating anode stator power supply adjustment	37
3	Set-up procedure	38
3.1	Service menu	39
4	kV adjustment	39
4.1	Adjustment “set-up” voltage	39
4.2	Adjustment of maximum kV setting safety	40
4.3	Checking maximum kV setting safety	40
5	Adjustment of the filament current	41
5.1	Adjustment of filament current	42
6	mAs adjustment	44
7	Adjustment of the battery voltage	45
8	Checking the radiographic parameter	46
8.1	Checking kV and mAs during radiography	46

## Section 4 – Drawings

### Annex – 115Vac Power supply (optional)

# Service Manual

---

INDEX

White Page

# Installation


# Service Manual

---

## INSTALLATION

### 1 INTRODUCTION

The installation procedure consists of different operations to be carried out in the order shown in the *flow chart* below.



The various operations indicated are described in detail in the following paragraphs.

*NOTE: the radiological system is normally pre-installed and configured in the factory according to the specific requirements of the customer. Certain adjustments may, however, be necessary when components are replaced.*

*The whole machine parameter adjustment procedure is described in detail in this manual.*



*After 30 minutes of non-use, the equipment automatically switches into stand-by mode, that means the capacitors battery is unloaded; the words “Stand By” appears on the display of control console and the value of the battery voltage is zero. Pressing any key of control console, the battery reloads to the maximum value and it is possible to make X-Ray exposure. During capacitors battery loading the words “Waiting” will appear on control console display.*

## 2 UNPACKING

The radiological system is packed in a single case containing all the parts of the apparatus. To unpack the unit, proceed as described on the following paragraph, according to the type of packing:

- simple pallet with cellophane – SEE § 2.1
- pallet with carton – SEE § 2.2
- wood packing crate – SEE § 2.3

If necessary, follow the unpacking procedure backwards to pack the unit again.

### 2.1 PACKING WITH CELLOPHANE

1. Cut the tape (4) and remove the slide (5) – SEE **FIGURE 1**
2. remove the protective bag (6)
3. unscrew the four self-threading screws (7) which fix the pedal fixing bracket (8) to the bottom of the crate
4. unscrew the four screws (9) which fix the pedal fixing bracket (8) to the mobile unit
5. mount the tilting pedal (10), located inside the cassette holder (11), by means of the screws and washers previously disassembled from the bracket (8)
6. free the X-ray Tube Head fixing brackets (12) from the tie-rod (17) - SEE **FIGURE 2** - working on nuts (13)
7. unscrew the self-threading screws (14) and remove the X-ray Tube Head fixing brackets (12)
8. put the slide (5) on the floor, as shown in **FIGURE 1**
9. slowly move the mobile unit towards the slide (5), keeping the "Dead-Man-Control" brake lever (18) pressed. Note that the front wheel (15) of the unit must come out of its packing housing (16)
10. let the unit move slowly from the slide (5) till its complete positioning on the floor and release the brake lever (18).

# Service Manual

## INSTALLATION




FIGURE 1 - LATERAL VIEW





FIGURE 2 - VIEW FROM ABOVE

## 2.2 PACKING WITH CARTON

1. Remove the screws (1) - SEE **FIGURE 3** - which fix the packing carton to the crate pallet (2)
2. cut the band (3)
3. remove the carton packing from its housing; the mobile unit is as shown in **FIGURE 1**
4. proceed according to the packing with cellophane procedure (§ 2.1).



**FIGURE 3** - PALLET WITH CARTON

# Service Manual

## INSTALLATION

### 2.3 PACKING WITH WOOD CRATE

1. Dismantle the crate cover (1) – SEE FIGURE 4
2. remove the crate walls: be careful not to ruin the wall (5) that has to be used as slide; the mobile unit is as shown in FIGURE 1

*Note: in the wood packing case, the slide (5) is not mounted as shown in FIGURE 1, because it is substituted by the wall (5), previously mentioned*

3. proceed according to the packing with cellophane procedure (§ 2.1).




FIGURE 4 - WOOD CRATE

## 3 ELECTRICAL CONNECTIONS

The power supply single-phase alternate voltage and maximum absorbed current values are reported both on the unit label and on the technical data in the user manual.



*Make sure that the power supply socket is approved for the values reported on the unit label*

## 4 ACCEPTANCE TEST

The **ACCEPTANCE TEST** includes three main kinds of operations:

- ➡ checking cables and connectors
- ➡ mechanical checks
- ➡ electrical checks

### 4.1 CHECKING CABLES AND CONNECTORS

Visually control each cable and each connector, checking that there are no interruptions and/or crushing.

- A.1** Power supply cable of the unit
- A.2** Radiography control pushbutton cable
- A.3** Connection cable between unit/X-ray Tube Head

### 4.2 MECHANICAL CHECKS

The mechanical checks include the following operations:

- ➡ **MOVEMENTS:** all the movements foreseen must be possible without excessive force or jerks
- ➡ **BRAKES:** all the movement locking brakes must be operating and easily applied

# Service Manual

---

## INSTALLATION

- ➡ **SOUNDNESS:** check that there are no breakages and/or damage which might have jeopardised operation and safety of the apparatus.

Follow the procedure below:

- B.1** Parking brake of the radiological unit: push the unit without activating the parking brake on “dead-man control” (21) and check that the unit does not move
- B.2** Turning the wheels of the radiological unit: release the parking brakes of the unit by pressing the “dead-man” lever (21). Move it a few metres and check that there are no movement imperfections and that there are no noises due to the wheels turning.
- B.3** Correct operation of the anti-tilting assemblies:
  - Tilting: place the unit on an horizontal surface in the transport position, carry out the tilting operation (SEE THE USER’S MANUAL) and check that the unit lowers over the rear but with slight friction due to the anti-tilting system.
  - Flag movement: loosen the arm brake by means of the handle (4), lift the arm/X-ray tube head assembly into the horizontal position and turn it 90° onto one side. Then check that the unit does not tilt onto the side in question. Repeat the procedure on the opposite side.
- B.4** Soundness of the control console: check that there are no scratches or breakages on the surface
- B.5** Soundness and fixing of the casings: check that there are no scratches or breakages on the casings and that these are coupled correctly without any visible openings.
- B.6** Balancing the pantograph arm: move the arm/X-ray Tube Head assembly vertically and check that it remains blocked for release by the operator in all the positions.

## 4.3 ELECTRICAL CHECKS

With regard to the electrical part, correct operation of the following must be checked:

- ➡ **SAFETY DEVICES**
- ➡ **SIGNALS**

*NOTE: Remove the carters to carry out the electrical checks (SEE SECTION 2 - §2.1).*

The procedure is described below:

- C.1** POWER SUPPLY: connect the unit to the mains power supply;
- C.2** INITIAL TEST: check that the initial automatic test is correct of the unit;
- C.3** SOFTWARE: check that the visualization of installed software;
- C.4** LOADING: check that the switching on the unit appears the increment of voltage on display until 320 V. The complete loading of the capacitors battery is ready in about 40/45 sec;
- C.5** DISCHARGING: switch the unit off and check the complete discharging of the capacitors battery by means the switching off the **LD3** led mounted on **S25** board;



*Do not carry out switching the unit On and Off*

- C.6** SET-UP DATA IN MEMORY: enter on SET-UP (SEE SECTION 3 - § 3) and check that the SET-UP data are in according with data selected (SEE SECTION 3 - § 3 – **TABLE 5**);
- C.7** X-RAY TUBE FORMATION: when not in use the unit for a long time, before to use the unit normally, carry out the “formation” of the X-ray tube with a serial of exposures using the following data (carry out for each value of **kV** and **mAs** 10 exposures and keep a pause of 1 minute between an exposure and another one):

40 kV	➤	100 mAs
50 kV	➤	80 mAs
60 kV	➤	63 mAs
70 kV	➤	63 mAs
80 kV	➤	50 mAs
90 kV	➤	50 mAs
100 kV	➤	40 mAs
110 kV	➤	40 mAs
120 kV	➤	40 mAs

# Service Manual

---

## INSTALLATION

During this formation check with oscilloscope the value of **kV** [**S19** board **Tp11** (-); **Tp3** (+) ] and **mA** [**S19** board: **Tp11** (-); **Tp6** (+) ] on inverter. Select the oscilloscope in the following way:

CHANNEL A: 1 V/DIV ( 1V read is equal to 50 mA)

CHANNEL B: 1 V/DIV ( 1V read is equal to 20 KV)

BASE TIME: 25ms/DIV

In case kV wave form presents high voltage discharge or you hear discharge inside X-ray tube interrupt the sequence, wait few minutes before starting again; repeat the initial formation checking that the fuse **F7** is not interrupted (100A) (SEE. ELECTRICAL DRAW);



**Warning: presence of 320 Vdc voltage!**

- C.8** EXPOSURE TIME: check that during exposure time the X-ray passage led until the end of exposure the display shows the real exposure time in seconds.
- C.9** STAND BY: let the equipment switched on in “Stand by” mode, that means without make X-Ray exposure, and verify that after 30 minutes the capacitors battery is unloaded, as explained in § 1. Press any key of control console to restore the capacitors loading and repeat the whole procedure for a total time of at least 4 hour.

# Maintenance

# Service Manual

---

## MAINTENANCE

### 1 ROUTINE MAINTENANCE

#### 1.1 GENERAL RECOMMENDATIONS

The radiological system requires regular checks and maintenance. The following recommendations have the aim of helping the operator to keep the apparatus in good working and safe conditions during service.

The system contains mechanical parts subject to wear according to use: following prolonged use, wear on parts may decrease safety during use. For this reason, it is essential for the checking and maintenance operations indicated below to be carried out consistently to protect the operators and patients against any damage caused by mechanical breakdowns.

Correct adjustment of the electrical and electronic systems has a direct influence on the operation of the system, on the quality of the image and on the electrical safety of the system, as well as on the level of exposure to radiation the operators and patients are subjected to.

The **MAINTENANCE PROGRAMME**, described in the following paragraphs, consists of controls and interventions to be carried out by specialised personnel authorised by the manufacturer. All maintenance operations are charged to the owner of the apparatus.



*Should it be necessary to replace components or parts which may in any way condition the safety of the machine, only use original spare parts.*

## 1.2 FREQUENT CHECKS AND INSPECTIONS

The operating personnel must be suitably trained to be able to carry out the daily and weekly checks indicated in **TABLE 1**.

The other 6-monthly and annual controls described in this chapter and the interventions described in the following chapters are reserved for qualified and authorised personnel of the technical assistance service.

**TABLE 1**

INTERVAL	CHECK	REFERENCE
<b>DAILY CHECKS</b>	Operation of the signals, displays and LEDs	
	Operation of the stationary brake	SERVICE MAN. – SECTION 2 - § 4.2
	Check that the warning and danger labels are intact	USER'S MAN. – SECTION 2 - § 1.7.1
<b>WEEKLY CHECKS</b>	Absence of oil leaks from the X-ray tube head	
	Absence of unusual noises in the X-ray tube head during X-ray emission	
<b>6-MONTHLY CHECKS</b>	General checks and inspections	SERVICE MAN. – SECTION 2 - § 1.3
<b>ANNUAL CHECKS</b>	Efficiency of the unit as described in the extraordinary mechanical maintenance	SERVICE MAN. – SECTION 2 - § 2.2

## 1.3 GENERAL CHECKS AND INSPECTIONS

Every six months, and in any case **according to local legislation in force**, the procedures indicated in the **ACCEPTANCE TEST** must be carried out completely.

Moreover the checks listed below must be added:

- D.1** AC and DC power supply
- D.2** Protection earth
- D.3** Fixing and general state (dust and corrosion) of the boards and of the connectors
- D.4** Centering of the X-ray tube head/collimator assembly (SEE SECTION 2 - § 2.3.6)

# Service Manual

---

## MAINTENANCE

### 1.4 CLEANING AND DISINFECTION

Products with an high content of alcohol, corrosive and/or abrasive detergents or solvents must not be used to clean the surfaces of the apparatus.

To disinfect the system, only use methods in compliance with the laws in force regarding disinfection and protection procedures against explosion.

To carry out the cleaning and disinfection operations, take the following precautions:

- ▣ turn the system off and disconnect the mains power supply cable
- ▣ make sure that no liquid gets into the apparatus so as to avoid any short-circuits or corrosion of the electrical and electromechanical parts.



*The unit has not to be used in presence of anaesthetic and/or inflammable disinfectant and cleaning products.*

*If, producing explosive gaseous, mixture, are used, make sure that gases are dispersed before switching on the unit.*

## 2 SPECIAL MAINTENANCE

### 2.1 DISMANTLING THE UNIT CASINGS

Any maintenance operation requires the external casings of the unit to be dismantled as indicated below (SEE FIGURE 5):

- ➡ dismantle the cassette-holder basket (1):
  - open the basket, releasing it from the relative locking jacks (2)
  - release the respective hinges (3) located under the basket, by means of the pins on them
  - to completely remove the basket, it is necessary to release the earthing cable fixed to the sheet of the basket
- ➡ dismantle the plastic nuts (4) which lock the assembly casing to the keyboard (5).
- ➡ loosen and completely dismantle the release handle (6) of the arm/X-ray tube head rotation



**ATTENTION:** *from this moment, the arm can rotate freely so do not move the unit because a sudden movement of the arm might occur, endangering people in the vicinity*

- ➡ release the two casings (7) and (8) of the arm
- ➡ unscrew the screws (9) and (10) and dismantle the relative shaft casings (11) and (12)
- ➡ dismantle the tilting handle (13), if provided
- ➡ release the assembly casing with the keyboard (5), which was released previously (see instruction n°2). To remove the assembly completely, disconnect the cables connected to the keyboard
- ➡ dismantle the casing with the cable support (14). To remove the casing completely, it is necessary to release the earthing cable fixed to the aluminium cable support
- ➡ dismantle the remaining casings (15) and (16) from the rubber strip section of the base.

At this point the unit is completely accessible for any maintenance operations.

# Service Manual

## MAINTENANCE

To remount the casings, repeat the procedure described above in reverse order, taking care to couple the casings correctly in their relative seats and grooves.




FIGURE 5 – EXPLODED VIEW OF COMPACT

## 2.2 EXTRAORDINARY MECHANICAL MAINTENANCE

Apart from the ordinary maintenance procedures (SEE § 1), it is advisable to carry out the checks indicated below (SEE FIGURE 6) annually, or following heavy use of the apparatus and, if necessary, adjust/replace the component (SEE § 2.3 and § 2.6).

Before starting, dismantle the unit casings as explained in § 2.1.

➡ Pantograph arm: check that there is grease in the moving parts. When there is any dirt or dust, remove the remaining grease and add some new grease as explained below:

- in the rotation points of the tie-rods (1) and (2) of the pantograph system: Loctite 8150 (A) type grease
- in the chain supporting part (3) on the cam (4): common C2 type grease (B)

Also check that there are no worn or broken parts which need replacement

- ➡ X-ray tube head cables: check that the connection cables to the X-ray tube head have no cuts or stripping, particularly in the points corresponding with rotation of the shaft and with vertical arm movement
- ➡ vulcanized wheels: clean the internal surface of the 2 wheels (5) and (6) and the contact surface of the pads (7) and (8) with degreasers to guarantee better locking of the unit
- ➡ general check: check correct fixing of the parts mechanically assembled and screwed up.




FIGURE 6 - SHAFT WITH ARM – BASE WHEELS

# Service Manual

---

## MAINTENANCE


### 2.3 MECHANICAL ADJUSTMENTS

#### 2.3.1 PARKING BRAKE (SEE FIGURE 7)

- ➡ The mobile unit remains braked:
  - the tension of the rope (1) is not sufficient: adjust it by means of the tensiometer (2)
  - the connection rod (3) to the pad (13) is not released or adjusted correctly: unscrew/screw up the rod, position the pad correctly (the space between the pad and wheel must be between 0.5 and 1 mm) and then lock using the relative counter nut (4)
- ➡ the mobile unit does not remain braked:
  - the two thrust dowels (14) of the springs (15) are not adjusted correctly and the pad (13) does not press correctly against the wheel (10): adjust by means of the dowels; if the spring (15) is lacking grease (common C2 type), add the amount needed
- ➡ when the wheels of the unit turn, there are noises or intermittent jerks:
  - carry out the adjustments indicated in the two previous points.

#### 2.3.2 ANTI-TILTING SYSTEM (SEE FIGURE 7)

- ➡ In the flag position 90° to the side, the gas spring is not released and the machine is unstable, attempting to tilt over:
  - the tension of the spring control ropes (11) is insufficient: adjust by means of the tie-rods (9)
  - the cam/bearing system (5) controlled by rotation of the shaft (8) does not move correctly: grease the parts if necessary, or replace the return spring (6) or the bearing (7) activated by the cam anchored to the shaft
- ➡ it is not possible to carry out the Tilting function:
  - the gas spring control rope (11) is broken: replace it
  - the gas spring (12) is faulty: replace it.



**FIGURE 7 - EXPLODED VIEW OF MACHINE**

# Service Manual

---

## MAINTENANCE

### 2.3.3 ARM-X-RAY TUBE HEAD BALANCING SYSTEM (SEE FIGURE 8)

#### ▣ Adjustment of the arm clutch


The balanced suspension system of the mobile unit allows smooth vertical movement of the arm and therefore of the x-ray assembly. This system has been studied and realised to reduce the need for maintenance to a minimum.

However, as time passes and due to continual stresses, the balancing system may lose a minimum amount of its efficiency. It is usually sufficient to adjust the boom pin (2) which acts directly on the clutch plates (1)

#### ▣ Adjustment of the counterweight spring force

If the previous operation is not sufficient to restore the optimal balancing condition, it is necessary to carry out an adjustment of the spring thrust (5) as described below:

- with a setscrew spanner, loosen the boom pin (2), watching out for sudden and undesired movements of the arm and of the x-ray assembly which takes the arm to the highest possible position
- just for device with 115V power supply, disassemble the 115V power supply group (SEE SECTION 4 - EXPLODED DRAWINGS)
- with the help of a second person, lower the arm until the nut (7) in the window (6) is accessible
- using a spanner (mm 12), screw up or loosen the nut (7) through the window (6) present on the shaft with the spring (4). This operations modifies the pre-load applied to the compression spring (5)
- to check correct balancing of the spring, make sure that the arm can stay in the horizontal position without any support
- on completion of balancing, tighten the boom pin (2) to take the clutch back to a compression value which, although keeping the arm stopped in any position in a stable way, allows its movement without difficulty (balancing to be checked with the arm casing mounted)
- on completion of the operations, put the casings back into their original position.



**FIGURE 8 - EXPLODED VIEW OF THE PANTOGRAPH ARM**

# Service Manual

---

## MAINTENANCE

### 2.3.4 LOCKING THE ARM IN THE TRANSPORT POSITION (SEE FIGURE 8)


- ➡ If no lock is noticed or the arm tends to rise during the movement stage when putting the unit into the transport position (SEE USER'S MANUAL) with the arm completely retracted downwards, the arm stops (3) are worn, and must therefore be replaced.

### 2.3.5 X-RAY TUBE HEAD-COLLIMATOR ASSEMBLY ROTATION (SEE FIGURE 9)



*ATTENTION – In cases where dismantling components of the X-ray tube head assembly or the assembly itself are foreseen, make sure the arm is locked so as to prevent any sudden and involuntary movements of the arm which would endanger the technician.*

- ➡ 360° rotation is difficult or not correctly adjusted:
  - inside the hole where the pin of the chrome-plated bracket (1) is housed dirt has accumulated due to wear on the parts in contact with the system of clamps (2) and (3): dismantle, clean carefully and lightly grease. (If this is not sufficient, replace the clamps). Then remount and lock using the correct force – just enough to make rotation of the bracket possible.
- ➡ X-ray tube head inclination ( -15° ÷ +105° ) is difficult or not correctly adjusted, or the X-ray tube head does not remain blocked in the position for operator release:
  - dismantle the lateral goniometer (4) and clean it of any dirt caused by wear on the clutch plate (5), then remount it locking the screws (6) with the correct force to prevent rotation movement blocking.
  - to guarantee correct pressure of the Belleville washers (7) which press on the clutch plate, it is also possible to increase or decrease the number of Belleville washers (7)
  - if cleaning is not sufficient, the clutch plate (5) may be worn, and must therefore be replaced.



**FIGURE 9 - X-RAY TUBE HEAD**

# Service Manual

---

## MAINTENANCE

### 2.3.6 CENTRING THE X-RAY TUBE HEAD – COLLIMATOR ASSEMBLY

#### ▣ Preliminary centring:

- turn the x-ray assembly so that the x-ray output is turned upwards
- position the adapter ring of the collimator and screw it onto the x-ray assembly, without tightening the fixing screws completely so that the ring itself can make small movements
- position the adapter and the centring device inside the ring
- position the x-ray assembly behind a glass sheet sealed with lead so that the graded part of the centring device can be seen
- carry out an exposure and check the position of the luminous dot which appears on the graded circular scale of the centring device
- adjust the adapter ring of the collimator so that the luminous dot is exactly in the centre of the graded circular scale of the centring device
- lock the ring in question, fully tightening the four screws which fix it to the x-ray assembly.

#### ▣ Final centring:

- position a washer at the centre of the collimator and, by means of the collimator light, position the “Fluorescent Screen” so that its centre coincides with the centre of the washer
- carry out an exposure and check that the washer has been impressed on the fluorescent screen.

***Note: For any other centring operations of the collimator, please refer to the TECHNICAL MANUAL of the collimator.***

## 2.4 ELECTRICAL TROUBLESHOOTING

Troubleshooting can be facilitated by the indications given in **TABLE 2**, which illustrates the most common incorrect operating conditions, showing their probable causes. In the presence of the following alarms visualized on the display, the X-ray exposition is disqualified and the console's alarm red led light up. In case it is also present after the execution of the suggested intervention, please contact the service assistance department.

**TABLE 2**

N°	NOTED FAULT	LIKELY CAUSE	RECOMMENDED INTERVENTION
1	The unit doesn't switch on (on key LED switched on)	No power supply on input	Check / Replace the S22-F2 or TR2-F24V fuse
	The unit doesn't switch on (on key LED switched off)	No power supply on input	Check / Replace the S22-F1 or F1 fuse (on input)
2	The apparatus works but does not emit rays – No alarm on the unit display	Radiography pushbutton fault	Check the Radiography pushbutton cable and connector
3	Alarm: "SUPPLY FAULT"	+15V or -15V fault of the S22 board supply	Check power supply (+15VA/-15VA/+24VA voltage on S22 voltage)
4	Alarm: "FILAMENT FAULT"	400 Hz power supply absent	Check / Replace TR2-F130V0 fuse
		X-ray tube filament interrupted	Replace the X-ray tube head
		S15 board fault	Replace the board
5	Alarm: "INVERTER FAULT"	If the indication appears on stand-by: P1 Set-Up on S22 board faulty	Carry out set-up as in Section 3 - § 4.2
		If the indication appears while carrying out X-ray: faulty on X-ray power circuit	Switching the unit off, switching it on and repeat the X-ray
		F2 and F3 fuses blown	Replace F2 and F3 fuses
		F6 and F7 fuses blown on the Power circuit	Replace F6 and F7 fuses
6	Alarm: "kV FAULT"	During a radiograph the effective kV are less than 85% of those set: fault on the power circuit	Switching the unit off, switching it on and repeat the X-ray
7	Alarm: "mA FAULT"	During an exposure the mA value is lower than the allowed limit	Switching the unit off, switching it on and repeat the X-ray
8	Alarm: "THERMIC SAFETY"	X-ray tube head overheated	Do not carry out any X-rays and wait for the X-ray tube head to cool down
		Breakage of the thermal sensor	Replace the thermal sensor located on the X-ray tube head and/or check the relevant wiring on S22 board
9	Alarm: "MAN STOP RX"	During a radiograph with cassette, the control pushbutton for X-ray command has been released early	Repeat the radiography
10	Alarm: "X-RAY LACKING"	Error in the high voltage generation circuit	Switching the unit off, switching it on and repeat exposure
11	Alarm: "MAX TIME"	The radiography time has exceeded the concurred limit	Check that the mA are correct Check the mA reading on board S22-Tp1 or S19-Tp6
12	Alarm: "STARTER FAULT"	During the anode preparation was not put into rotation	Check the fuses on board S21. Check the stator winding of the X-ray tube head

# Service Manual

## MAINTENANCE

N°	NOTED FAULT	LIKELY CAUSE	RECOMMENDED INTERVENTION
		Fault on board S21 Fault on board S20	Replace board S21 Replace board S20
13	Alarm: "SWITCH OFF FOR 1 MIN"	Capacitors bank still loads	With the unit off wait 1 minute before switching it on
14	Alarm: "BATTERY FAULT"	Power Circuit Fault	Contact technical service
15	Alarm: "OVERVOLTAGE BATTERY"	Battery Circuit Fault	Contact technical service
16	Alarm: "X-RAY COMMAND ACTIVE"	The operator has pressed the radiography command before the system had finished the initial control stage	Release the radiography pushbutton and wait until the system is ready
17	Alarm: "X-RAY TUBE TOO HOT"	It is not possible to begin exposure since the remaining thermal units available are too few	Wait for the X-ray tube to cool down.
18	Alarm: "WAIT CONNECTION"	The keyboard does not communicate with the unit	Switching the unit off, switching it on and repeat exposure
19	Alarm: "POTTER FAULT" (only with Potter installed)	Potter Fault	Check the F4 - F5 fuses of the Potter power supply and the K3 relè

## 2.5 THERMOMAGNETIC CIRCUIT-BREAKER (SEE FIGURE 10)

The system is protected on the power supply circuit by a thermomagnetic circuit-breaker located inside the machine. Should this circuit-breaker trip, dismantle the cassette-holder basket (1) as explained in the first point of § 2.1 and reset the circuit-breaker (2) which is accessible at the back of the right-hand side of the machine.




*Take special care not to touch the resistances placed above the thermomagnetic circuit-breaker because they may be very hot.*



*Any trip of the thermomagnetic circuit-breaker immediately after resetting can be an indication of incorrect operation of the system: contact the technical service assistance immediately.*



*In case of device with 115V power supply, another thermomagnetic circuit-breaker is provided in the 115V power supply assembly (SEE ANNEX - § 1).*



**FIGURE 10 - THERMOMAGNETIC CIRCUIT-BREAKER**

## 2.6 REPLACEMENT OF COMPONENTS

To dismantle the unit and for the list of spare parts, please refer to SECTION 4 – EXPLODED DRAWINGS and to § 2.1. If necessary, repeat the checks for the Acceptance Test regarding the component replaced (SEE SECTION 1 - § 4). Furthermore, remember that replacement of a component may require checks and adjustments to restore correct operation of the radiological unit. For this purpose, please see § 2.3 for the mechanical parts, and **TABLE 3** and SECTION 3 of this manual for the electrical parts.

**TABLE 3**

COMPONENT REPLACED	ADJUSTMENT AND CHECKS	REFERENCE
15 BOARD	Verify <b>mA</b> set up in SET UP	SECTION 3 - § 5
S19 BOARD	Check of <b>kV e mA</b>	SECTION 3 - § 5.1 and 8
S20 BOARD	Check Run voltage Check start 800 ms	SECTION 3 - § 2
S22 BOARD	Check max Set kV (P1) Check max exposure time (P2)	SECTION 3 - § 4.3 SECTION 3 - § 4.1 and 4.2
S25 BOARD	Check max voltage of battery	SECTION 3 - § 7
S62 BOARD	Check LCD Display contrast (P1)	

# Service Manual

## MAINTENANCE



*Before replacing any component it is necessary to make the unit safe in order to prevent dismantling of the parts causing sudden and involuntary movements, which endanger personnel in the vicinity. In particular, pay attention to the following:*

- dismantling the X-ray tube head-collimator assembly
  - add the counter-weight on the pivoting wheel of the unit or under the tilting pedal to prevent the unit itself tilting over
  - secure the pantograph arm to the base of the unit (using the ropes) to prevent the arm from making undesired movements caused by the lack of weight which occurs when the X-ray tube head is removed
- anti-tilting assembly system
- rotation of shaft with arm–X-ray tube head
- vertical movement of the pantograph arm
- parking brake system.

## 2.7 LIST POTENTIOMETERS

BOARD	POTENTIOMETER	VALUE	NAME	NOTE
S15	P1	5K	MANUAL SET I FIL (ONLY FOR TEST)	Only factory adjustment
	P2	5K	MAX IFG	Only factory adjustment
	P3	5K	MAX IFP	Only factory adjustment
S17	P1	10K	FEEDBACK KV+	Only factory adjustment
	P2	10K		Only factory adjustment
	P3	10K	FEEDBACK KV-	Only factory adjustment
	P4	10K		Only factory adjustment
S19	P1	2K	MAX CURRENT	Only factory adjustment
	P2	2K	FREQUENCY SET	Only factory adjustment
	P3	2K	MA FLUORO OFFSET ADJ.	N.U.
	R66		FEEDBACK MA ADJ	Only factory adjustment
S20	P1	1K	RUN VOLTAGE	Only factory adjustment
	P2	1M	START 800MS	Only factory adjustment
S22	P1	2K	MAX SET KV	Only factory adjustment
	P2	100K	MAX TIME RAD	Only factory adjustment
S25	P1	5K	MAX V CAPACITORS BANK	Only factory adjustment
S62	P1	20K	LCD CONTRAST	Only factory adjustment

# Service Manual

MAINTENANCE

## 2.8 LIST LEDS

BOARD	LED	COLOR	NAME	NOTE
S15	LD1	YELLOW	+24V	
	LD2	GREEN	+15VA	
	LD3	RED	SAFETY MAX I FIL	
S19	LD1	YELLOW	COM-RX	
	LD2	RED	kV MAX	
	LD3	RED	kV > 110%	
	LD4	RED	kV MIN	
	LD5	RED	I MAX	
	LD6	YELLOW	COM. 1	
	LD7	YELLOW	COM. 2	
	LD8	YELLOW	PRE-RX	
	LD9	GREEN	+15V	
	LD10	GREEN	-15V	
S20	LD1	YELLOW	STARTER ON	
	LD2	GREEN	STARTER OK	
	LD3	YELLOW	FLUORO PREP (N.U.)	
S22	DL1	YELLOW	N.U.	
	DL2	YELLOW	N.U.	
	DL3	YELLOW	N.U.	
	DL4	YELLOW	N.U.	
	DL5	YELLOW	N.U.	
	DL6	YELLOW	N.U.	
	DL7	YELLOW	N.U.	
	DL8	YELLOW	N.U.	
	DL9	YELLOW	STARTER OK	
	DL10	YELLOW	PREP K1A	
	DL11	YELLOW	RAD	
	DL12	YELLOW	kV 85%	
	DL13	YELLOW	POTTER OK	
	DL14	GREEN	THERMIC	
	DL15	YELLOW	kV OK	
	DL16	GREEN	FIL OK	
	DL17	YELLOW	PRE RX	
	DL18	YELLOW	RX	
	DL19	YELLOW	PREP K1A	
	DL20	YELLOW	START POTTER	
	DL21	GREEN	+24VA	
	DL22	GREEN	+24VA	
	DL23	YELLOW	ON/OFF	
	DL25	GREEN	+24VA	
	DL26	YELLOW	CAPACITORS BANK CHANGE	
	DL27	YELLOW	CAPACITORS BANK CHANGE	
	DL28	YELLOW	OUT_WD	
	DL29	YELLOW	PREP STARTER	
	DL30	YELLOW	FP	
	DL31	YELLOW	CORR MA	

Section 2

# Service Manual

## MAINTENANCE

### 2.9 LIST FUSES

BOARD	NAME	TYPE	NOTE
LOADING BATTERY	F2	16A	Input power supply protection
	F3	16A	
INVERTER	F6	63A	Inverter protection
	F7	63A	
TR1 TRANSFORMER	F1	200mA-T	TR1 Transformer protection
TR2 TRANSFORMER	F230V	2A-T	TR2 Transformer protection
	F115V	3.15A-T	
	F135V	1A-T	
	F24V	2A-T	
	F12V	10A-T	
POTTER BUCKY	F4V	2A-T	Potter Bucky protection(Optional)
	F5V	2A-T	
AUTOTRANSFORMER	F10	2A-T	115Vac power supply protection(Optional)
S21	F1	10A-T	S21 Board Input power supply protection
	F2	10A-T	
	F3	500mA-T	
S22	F1	1A-T	Bridge Rectifier protection (24Vac)
	F2	2A-T	
S25	F1	80mA-T	Input power supply protection S25 Board
S81	F1	200mA-T	Group 1 Bridge Rectifier protection (24Vac)
	F2	200mA-T	
	F3	200mA-T	Group 2 Bridge Rectifier protection (24Vac)
	F4	200mA-T	
	F5	200mA-T	Group 3 Bridge Rectifier protection (24Vac)
	F6	200mA-T	
	F7	200mA-T	Group 4 Bridge Rectifier protection(24Vac)
	F8	200mA-T	
	F9	315mA-T	S81 Board Input power supply protection

## 3 SPARE PARTS

In case of replacing of parts that can be negative for the safety of the machine, use only original parts. For the list of the spare parts and relative codes make reference to the SECTION 4 – EXPLODED DRAWINGS.

**NOTE:** When you require spare parts, it is necessary to tell the code of the piece and the serial number of the unit too.

# X-ray settings

# Service Manual

## X-RAY SETTINGS

### 1 INTRODUCTION

The radiological system is fully regulated in the factory: it is therefore only necessary to install it in order to use it.

All the information contained in this section and SECTION 2 of this manual must be used when extraordinary maintenance is necessary.



*In the case of any intervention, it is essential for the adjustment operations to be carried out in the sequence indicated.*



*In Table 4 indicates the potentiometers which do not require any further adjustments since they have already been regulated and sealed in the factory.*

**TABLE 4**

BOARD	POTENTIOMETER	DESCRIPTION
S15	P1	Manual Setting of the filament current
	P2	Maximum current of large focus
	P3	Maximum current of small focus
S17 (X-ray tube head)	P1	Return of the kV+
	P2	
	P3	Return of the kV-
	P4	
S19	P1	Maximum current
	P2	Frequency setting
S22	P1	Maximum kV setting
	P2	Maximum exposure time
S25	P1	Maximum battery voltage

**NOTE : Jumper S19-J2A and J2B ON: safety  $kV_{max}$  set up for rotating anode**

**Jumper S19-J2A and J2B OFF: safety  $kV_{max}$  set up for stationary anode**


## 1.1 DIP SWITCH

Here follows the function of the dip-switch **S22-SW1** pins:

- |   |                             |
|---|-----------------------------|
| 1 – Only factory adjustment                           | 6 – Only factory adjustment |
| 2 – Only factory adjustment                           | 7 – Only factory adjustment |
| 3 – ON – It doesn't display the name when switched on | 8 – Only factory adjustment |
| 4 – Only factory adjustment                           |                             |
| 5 – ON – Small Focus when switched on                 |                             |

## 2 ROTATING ANODE STATOR POWER SUPPLY ADJUSTMENT

- ➡ Connect an **AC** voltmeter (*true RMS* if possible) with full-scale set at **300 V** between the **TB3-2** and **TB3-3** test points of board **S21** (power supply of the main winding);
- ➡ connect the probe of an oscilloscope between **IC2-8 (+)** and **Dz1 Anode (-)** of board **S20** and prepare the instrument to read a time of **800 ms**;
- ➡ take the **RADIOGRAPHY COMMAND PUSHBUTTON**.



*For maximum operator safety, at this point it is advisable to temporarily disable the X-ray command signal by disconnecting the wire coming from S22-X4-3.*

- ➡ only press the **first click** of the **RADIOGRAPHY COMMAND PUSHBUTTON** (*preparation pushbutton*);
- ➡ check that the stator is supplied with a voltage of about **230 Vac** for approximately **0,8 seconds**. If necessary, adjust this time interval using potentiometer **P2** of board **S20**;
- ➡ check that once the 0,8 seconds have passed the stator is supplied with a voltage of about **70 Vac**. If necessary, adjust this value by using potentiometer **P1** of board **S20**.

# Service Manual

## X-RAY SETTINGS


### 3 SET-UP PROCEDURE

To access the modality of SET-UP procedure to follow what underneath reported:


- press the keys of the console indicated to side at the same time



- the wording "SET-UP PROCEDURE ?" will appear on the display; to confirm selection, press the key indicated at the side marked with "YES", whereas to cancel the procedure, press the key marked "NO".



**NOTE:** In the Set-up procedure, the functions for editing the content of the fields shown on the display have been assigned to the keys as indicated below.



**TABLE 5 – MINIMUM, MAXIMUM AND OF DEFAULT VALUES OF SET-UP ADJUSTMENT**

FUNCTION		VALUE			
DIRECTORY	FILE	MINIMUM	MAXIMUM	DEFAULT	
kV ADJUSTMENT	KV40	300	600	380	
	kV120	1000	1500	1190	
SMALL FIL. ADJ	7.5 kW	150 mA (40 kV)	1000	2000	1600
		75 mA (100 kV)	900	2000	1360
LARGE FIL. ADJ	30 kW	425 mA (40 kV)	1000	2000	1750
		300 mA (100 kV)	900	2000	1560
MAs ADJ	kV 050 (50 mAs)	250	1000	530	
BATTERY	V = 320				

## 3.1 SERVICE MENU

When the unit is off, pushing both the keys underneath you enter the SERVICE MENU of the unit which subsequently visualizes the following data:

- name of the device and software version
- chosen anatomical program
- X-ray tube model
- total number of radiographic exposures
- total mAs accumulated.





## 4 kV ADJUSTMENT

### 4.1 ADJUSTMENT “SET-UP” VOLTAGE

The SET KV voltage is supplied by the microcomputer and must be adjusted according to the relationship  $kV/20$  (i.e.:  $1V = 20kV$ ).

- First of all, enter the SET-UP procedure, following the instructions indicated in the previous paragraph;
- On board S19 connect the probes of a Vdc digital multimeter between CP1-10 (+) and TP11 (-);
- in Set-up select "kV ADJ at 40 kV";
- adjust the value until a reading of 2,00 V is obtained on the voltmeter;
- in SET-UP select "kV ADJ at 120 kV";
- adjust the value until a reading of 6,00 V is obtained on the voltmeter;
- exit the SET-UP procedure by using the key;
- then select 100 kV, 80 kV and 60 kV, checking that the voltage read on the voltmeter is  $V_{read} = (Kv_{SET}/20) \pm 25 mV$ .




# Service Manual

## X-RAY SETTINGS

### 4.2 ADJUSTMENT OF MAXIMUM kV SETTING SAFETY



The maximum parameter value of **SET kV** under operating conditions must be **6 V**, corresponding to **120 kV**. The **MAX SET kV** safety device intervention is normally set to **130 kV**.

- ➡ connect the probes of a **Vdc** digital voltmeter between test points **TP10 (+)** and **TP9 (-)** of board **S22**;
- ➡ turn the system on and use potentiometer **P1** until a reading of **6,5 V ± 0,1 V** is obtained.



### 4.3 CHECKING MAXIMUM kV SETTING SAFETY

- ➡ Connect the voltmeter to test points **TP10 (+)** and **TP9 (-)** of board **S22**
- ➡ lower the **MAX SET kV** safety intervention to **100 kV** adjusting potentiometer **P1** until a reading of **5 V ± 0,1 V** is obtained on the voltmeter;
- ➡ checking that by setting the **kV** to a higher value than **100 kV** the "**kV FAULT**" alarm goes off;
- ➡ put the **MAX SET kV** safety intervention back to **130 kV** using potentiometer **P1** again until a reading of **6.5 kV ± 0,1 V** is obtained.



## 5 ADJUSTMENT OF THE FILAMENT CURRENT

The following tables indicate the **mA** anode values associated with the **kV** and the typical **FILAMENT CURRENTS** determined on a sample **X-RAY TUBE HEAD** for **30 kW** and **7.5 kW** power.

**TABLE 6**

		RADIOGRAPHY – 30 kW	
Set kV	kV	mA	I <sub>Filament</sub>
S22: TP10 – TP9 [V]	on display [kV]	S19: TP6 - TP11 [1V=10mA]	S15: CP6-2 [mA]
2	40	425	870
2,5	50	400	817
3	60	375	822
3,5	70	355	803
4	80	337	814
4,5	90	320	806
5	100	300	795
5,5	110	204	769
6	120	141	753
		Pre-switching:	466

**TABLE 7**

		RADIOGRAPHY – 7.5 kW	
Set kV	kV	mA	I <sub>Filament</sub>
S22: TP10 – TP9 [V]	on display [kV]	S19: TP6 - TP11 [1V=10mA]	S15: CP6-2 [mA]
2	40	150	826
2,5	50	150	777
3	60	125	757
3,5	70	107	747
4	80	93	737
4,5	90	83	727
5	100	75	716
5,5	110	68	706
6	120	62	696
		Pre-switching:	471

# Service Manual




## X-RAY SETTINGS

### 5.1 ADJUSTMENT OF FILAMENT CURRENT

**NOTE:** by *Filament Current* the current which runs through the primary winding of the transformer supplying the power to the filament itself is intended (*Coolidge Transformer*).

#### MEASUREMENTS

- To measure the Filament Current, connect an ammeter with *true RMS* reading in series with terminal **CP6-2** of board **S15**. Set the full-scale to **1A** in **AC**.
- to measure the anode **mA** in **RADIOGRAPHY**, connect the probe of an oscilloscope between test points **TP6 (+)** and **TP11 (-)** of board **S19** to display the wave form of the mA (reading correspondence is **1V = 50 mA**) and between test points **TP3 (+)** and **TP11 (-)** of board **S19** to display the wave form of the **kV** (reading correspondence is **1V = 20 kV**).



#### PROCEDURE (example with Small focus LP 7,5 kW)

- activate the **SET-UP** procedure;
- select "**SMALL FIL ADJ**";
- adjust the "**DAC**" value set until a filament current, displayed on the digital voltmeter, equal to the one indicated in **TABLE 7** for **40 kV (826)** is obtained – Typical **DAC 1600** – **FIGURE 11**;
- give the **X-ray** command and check that the **anode mA** are **150**. If necessary, adjust the **SET-UP** until the desired value is obtained;
- proceed with **SET-UP**;
- adjust the set value of the "**DAC**" until a **filament current**, displayed on the digital voltmeter, equal to the one indicated in **TABLE 7** for **100 kV (716)** is obtained – Typical **DAC 1360** – **FIGURE 12**;

- give the **X-ray** command and check that the **anode mA** are **75**;  
if necessary, adjust the **SET-UP** until the desired value is obtained;
- proceed with **SET-UP** by selecting "**LARGE FIL ADJ**";
- Repeat the same operations described for "**SMALL FIL ADJ**"  
with **40 kV** and **100 kV**. For the anode and filament current mA  
values, refer to **TABLE 6**.
- Exit **SET-UP**.




FIGURE 11 – DAC (40 kV) ADJUSTMENT






FIGURE 12 – DAC (100 kV) ADJUSTMENT



# Service Manual

---

## X-RAY SETTINGS


### 6 mAs ADJUSTMENT

#### ⇒ MEASUREMENTS

- Remove the bridge between terminals **CP3-1** and **CP3-3** of board **S19**
- connect a mAsmeter between the above-mentioned terminals.

#### ⇒ PROCEDURE


- activate the **SET-UP** procedure
- select "**mAs ADJ**" and give the **X-ray** command
- check that the **mAs** value on the instrument is really the one set (**50 mAs** – Typical **DAC 0530**); If necessary, adjust the **SET-UP** until the desired value is obtained
- exit **SET-UP**.



## 7 ADJUSTMENT OF THE BATTERY VOLTAGE

### ➡ MEASUREMENTS

- To measure the battery voltage, connect a voltmeter with true RMS reading between terminals **CP2-1 (+)** and **CP2-3 (-)** of board **S25**.



*Warning: presence of 320 Vdc voltage*

### ➡ PROCEDURE

- activate the **SET-UP** procedure;
- select "**V BATTERI ADJ**";
- Adjust the Set-up until the point in which the battery voltage measure value, reaches the battery voltage set value (tip. **320V**);
- exit **SET-UP**.

# Service Manual

## X-RAY SETTINGS

### 8 CHECKING THE RADIOGRAPHIC PARAMETER


After having completed the operations indicated in the previous paragraphs (§ 2, § 4, § 5), it is necessary to check the values of the signals proportional to the **kV** and to the **mA** which come from the **INVERTER**. The tolerance for the **kV** signal is  $\pm 5\%$  whereas for the **mA** it is  $\pm 15\%$ .

*NOTE: If the kV value is not correct, it is necessary to repeat the SET kV adjustment.*

*NOTE: An oscilloscope with two channels is required for the checks below.*

#### 8.1 CHECKING kV AND mAs DURING RADIOGRAPHY

- ➡ Connect the probe of the oscilloscope as indicated below:  
 CHANNEL A: **TP3 (+) – TP11 (-)** of board **S19 (REAL kV)**;  
 TIME BASE: 10 ms
- ➡ If necessary, connect the probe of the oscilloscope as indicated below:  
 CHANNEL B: **TP6 (+) – TP11 (-)** of board **S19 (mA)**;  
 TIME BASE: 10 ms
- ➡ set the minimum value of the **kV (40 kV – 425 mA)**
- ➡ give the X-RAY command and check that the values read on the oscilloscope correspond with those indicated in **TABLE 8**; (**FIGURE 13** shows an example of a 100 kV wave form)
- ➡ repeat the previous operations, setting the **kV** values to **80 kV** and **120 kV** for **7.5 kW** power.




**TABLE 8 – 30Kw**


kV		mAs		mA	Time (s)
Set	Acceptance range S19: TP3 - TP11 [1V = 20 kV]	Set	Acceptance range ( $\pm 10\%$ )	Theoretical	Theoretical
40	38 ÷ 42	16	14.5 ÷ 17.5	425	0.037
80	76 ÷ 84	16	14.5 ÷ 17.5	337	0.047
120	114 ÷ 126	16	14.5 ÷ 17.5	141	0.113

**TABLE 9 – 7.5KW**

kV		mAs		mA	Time (s)
<i>Set</i>	<i>Acceptance range</i> S19: TP3 - TP11 [1V = 20 kV]	<i>Set</i>	<i>Acceptance range</i> (±10%)	<i>Theoretical</i>	<i>Theoretical</i>
40	38 ÷ 42	8	7.2 ÷ 8.8	150	0.053
80	76 ÷ 84	8	7.2 ÷ 8.8	93	0.086
120	114 ÷ 126	8	7.2 ÷ 8.8	62	0.129



**FIGURE 13 – ACQUISITION OF THE kV WAVE FORMS DURING RADIOGRAPHY [100 kV – 16 MAS – 30 kW]**



**FIGURE 14 – ACQUISITION OF THE mA WAVE FORMS DURING RADIOGRAPHY [100 kV – 16 MAS – 30 kW]**

# Service Manual

---

X-RAY SETTINGS

White Page

# Drawings

# Service Manual

---

DRAWINGS

## Index

**Document n°**

### Exploded views


- Mechanical 10551
- Electrical 10551

### Electrical diagrams

- Electronic group 100 10465
- S15/B board – filament power supply 9678
- S17/A board – X-ray tube head interface 7548
- S19/C board – inverter control 9686
- S20/A board - X-ray tube stator starter 9674
- S21/B board – X-ray tube stator rotation 10192
- S22 board – microprocessor 9706
- S25 board – capacitors bank control 9680
- S62/A board – keyboard interface 8202
- S81 board – inverter control 9331

### Topographic diagrams

- S14 board – line filter 5040
- S15/B board – filament power supply 9677
- S17/A board – X-ray tube head interface 7547
- S19/C board – inverter control 9685
- S20/A board - X-ray tube stator starter 9669
- S21/B board – X-ray tube stator rotation 9670
- S22 board – microprocessor 9705
- S25 board – capacitors bank control 7613
- S62/A board – keyboard interface 8203
- S81 board – inverter control 9574



N°	DESCRIPTION	CODE	NOTE
<b>MOBILE UNIT</b>			
1	Handle for tilting polished	grey	01.02.052.006
		green	01.02.052.004
2	Carter with flat band stock	01.02.010.143	
3	Carter without flat band stock	01.02.010.144	
4	Carter with flat band arm	01.02.018.014	
5	Carter without flat band arm	01.02.018.015	
6	Carter prof. inf. articulated Arm	01.02.010.085	
7	Carter back	01.02.010.140	
8	Carter basement	01.02.010.141	
9	Carter front	01.02.010.139	
10	Cassette holder basket	01.02.043.003	
11	Pedal for tilting	01.02.010.052	
12	Carter control panel	01.02.010.142	

Drawing N° 01.02.043.001 - 007681

# Service Manual

---

MECHANICAL EXPLODED VIEW


Section 4

Page 52 of 69

Doc: 10551

# Service Manual

MECHANICAL EXPLODED VIEW



Drawing N° 01.02.008.012 - 003507

N°	DESCRIPTION	CODE	NOTE
<b>X-RAY TUBE HEAD SUPPORT GROUP</b>			
13	X-ray tube head 100KHz	01.02.008.012	
14	S17//A X-ray tube head board	01.17.008.025	Depends on the customer
15	Sup. carter with flat band	01.02.008.002	
16	Inf. carter without flat band	01.02.008.003	
17	Goniometer	light blue	01.02.014.009
		light green	01.02.010.148
18	Brake disk Ø 62.5	01.02.001.108	
19	Small X-ray tube head handle	grey	01.02.008.006
		green	01.02.025.058
20	Distance ring for rotating group	01.02.025.065	
21	Rotating adaptor Ralco mm. 18	49.18.002.001	
22	Lamp for centering 12V 100W	49.05.012.001	
23	Collimator R105 12V	49.18.006.001	

# Service Manual



---

MECHANICAL EXPLODED VIEW

Section 4

Page 54 of 69

Doc: 10551



N°	DESCRIPTION	CODE	NOTE
	<b>CASSETTE HOLDER BASKET</b>	01.02.043.003	
24	Cassette holder pushlock system	49.07.038.001	
	<b>KEYBOARD GROUP</b>	01.24.012.015	
25	Control Panel	Depends on the customer	
26	S62/A Serial keyboard card	01.17.037.002	

# Service Manual


---

MECHANICAL EXPLODED VIEW

Section 4

Page 56 of 69

Doc: 10551



Drawing N° 01.02.043.002 - 007680

N°	DESCRIPTION	CODE	NOTE
	MECHANICS/ELECTRONIC GROUP	01.02.043.002	
27	Double-click pushbutton	01.02.010.229	
28	Holdfast with thread	01.02.001.170	
29	Holdfast with through hole	01.02.001.169	

# Service Manual


---

MECHANICAL EXPLODED VIEW


Section 4

Page 58 of 69

Doc: 10551



Drawing N° 01.02.043.004 - 007671



Drawing N° 01.02.043.006 - 007669

N°	DESCRIPTION	CODE	NOTE
<b>BASEMENT/CASING/COLUMN</b>			
30	Rubber buffer	01.02.001.111	
31	Column move lock pull knob	49.07.018.005	
<b>CASING GROUP</b>			
32	Miniplast foot MP 68-10	49.13.018.001	
33	Pull spring T31620	49.07.021.007	
34	Radial bearing D.8x22x7 608-2Z	49.07.019.005	
35	Steel brakes rope with pawl	49.07.033.001	
36	Needle cage 55x78x3 AXK-5578	49.07.019.006	

# Service Manual


---

MECHANICAL EXPLODED VIEW


Section 4

Page 60 of 69

Doc: 10551



Drawing N° 01.02.010.029 - 000469



Drawing N° 01.02.010.028 - 000241

N°	DESCRIPTION	CODE	NOTE
<b>ARM GROUP</b>			
37	Lock for arm-transport mode	01.02.010.123	
38	Plate brake Ø 55	01.02.001.106	
<b>COLUMN GROUP</b>			
		01.02.010.028	

# Service Manual


---

MECHANICAL EXPLODED VIEW


Section 4

Page 62 of 69

Doc: 10551




Drawing N° 01.02.010.206 - 000640




Drawing N° 01.02.006.083 - 000765

N°	DESCRIPTION	CODE	NOTE
<b>SLIDING BLOCK GROUP</b>			
39	Pin with sliding block	01.02.010.204	
40	Brake pin spring D12400	49.07.021.004	
<b>BASEMENT/WHEEL GROUP</b>			
41	Pivoting wheel 80X30 64-4551	49.07.016.023	
42	Rubber protection	light-blue	01.02.014.007
		light-green	01.02.025.060
<b>COMPLETE ANTI-OVERTURN ASSEMBLY</b>			
43	Spring head-attachment	49.13.019.002	
44	Compression spring D12550	49.07.021.008	
45	Block spring KO-B1-K-3-030	49.13.019.001	
46	1 mt. cable	49.13.019.003	
47	Pivoting wheel 2470 UOK 050 S70	49.07.016.007	



Drawing N° 01.02.043.005 - 007668



Drawing N° 01.02.006.116 - 000742

# Service Manual

---


MECHANICAL EXPLODED VIEW

Section 4


Page 64 of 69

Doc: 10551

Drawing N° 007823



N°	DESCRIPTION	CODE	NOTE
<b>ELECTRONICS GROUP</b>			
BR1	Graetz bridge Semikron SKB 30/16A1	01.02.047.055	
C1/Starter	30 µF 450 Vdc Power capacitor	49.06.018.015	
FAN	Fan 80x80x25 230Vac 50Hz	49.06.023.008	
F1 - Fvent	Fuse 5x20 200mA 250V	49.06.022.025	
F2-F3	Fuse 10x38 16A 500V GL	49.06.022.003	
F6	Fuse Brush 63FE (63A 660V)	49.06.022.001	
F12V	Fuse 5x20 10A 250V	49.06.022.022	
F24V-F230V	Fuse 5x20 2A 250V	46.06.022.013	
F135V	Fuse 5x20 1A 250V	49.06.022.011	
K1B-K2-K0	BC6-30-01 24Vdc 4KW Contactor	49.06.015.017	
K2/A	CA6-11M Auxiliary contactor	49.06.015.021	
L1	2,5 mm <sup>2</sup> Filter inductor	01.02.010.195	
LF1	EPCOS B84112-B-B120 Filter	49.06.030.006	
R1-R2	ATE RB 106 25E Resister	49.06.011.080	
S14	S14 Filter board	01.17.005.005	
S15/B	S15/B Filament board	01.17.001.029	
S19/C	S19/C Inverter board	01.17.008.040	
S20/A	S20/A Starter board	01.17.006.003	
S21/B	S21/B Rotation board	01.17.003.006	
S22	S22 Micro board	01.17.010.004	
S25	S25 Capacitors bank Loading board	01.17.005.009	
SW1	Int. Aut. Multi 9 C60N C16A 2P	49.05.008.002	
TH1	Powerblock TT61 N12 KOF	49.06.021.013	
TR1	TM 15/24 Module transformer	01.05.001.017	
TR2	300VA Tor. transformer	01.05.001.026	
TR3	230-250/330 600VA Tor. transformer	01.05.001.016	



Drawing N° 01.02.047.012 - 007824

# Service Manual


---

ELECTRICAL EXPLODED VIEW

Section 4


Page 66 of 69

Doc: 10551



Drawing N° 01.02.047.015 - 007830

N°	DESCRIPTION	CODE	NOTE
C36-C27	1500 µF 400Vdc Elect. capacitor	49.06.012.079	
C29	(80A) 2 µF 1700 Vdc mk V capacitor	49.06.012.078	
F7	Fuse Brush 100 FE (100A 660V)	49.06.022.010	
L2	5,1 µH 45112 Tor. inductor	01.05.003.006	
Q1-Q2	Tr. IGBT Siemens BSM 300GB 120 DLC	49.06.017.025	
S81	S81 Driver board	01.17.081.001	
TR4	100KHz Driver for. transformer	01.05.001.018	
TA1	ENCO 10701 1/500 Amper transformer	01.05.004.002	



Drawing N° 01.02.025.018 - 003302

# Service Manual



---

ELECTRICAL EXPLODED VIEW


Section 4

Page 68 of 69


Doc: 10551



Drawing N° 01.02.047.016 - 007831



Drawing N° 01.02.047.017 - 007832




Drawing N° 01.02.025.027 - 003259

N°	DESCRIPTION	CODE	NOTE
C1-C20B	V 16000 µF 350Vdc Itecond Elect. capacitor	49.06.012.076	
R4-R4A-R5-R5A	IRE RCC 47 Ω Resister	49.06.011.088	










The callbacks refer the same document

Item No.	Description	Quantity	Unit	Remarks
1	RESISTOR	100	PCB	
2	CAPACITOR	50	PCB	
3	CONNECTOR	20	PCB	
4	WIRE	1000	PCB	
5	WIRE	1000	PCB	
6	WIRE	1000	PCB	
7	WIRE	1000	PCB	
8	WIRE	1000	PCB	
9	WIRE	1000	PCB	
10	WIRE	1000	PCB	
11	WIRE	1000	PCB	
12	WIRE	1000	PCB	
13	WIRE	1000	PCB	
14	WIRE	1000	PCB	
15	WIRE	1000	PCB	
16	WIRE	1000	PCB	
17	WIRE	1000	PCB	
18	WIRE	1000	PCB	
19	WIRE	1000	PCB	
20	WIRE	1000	PCB	
21	WIRE	1000	PCB	
22	WIRE	1000	PCB	
23	WIRE	1000	PCB	
24	WIRE	1000	PCB	
25	WIRE	1000	PCB	
26	WIRE	1000	PCB	
27	WIRE	1000	PCB	
28	WIRE	1000	PCB	
29	WIRE	1000	PCB	
30	WIRE	1000	PCB	
31	WIRE	1000	PCB	
32	WIRE	1000	PCB	
33	WIRE	1000	PCB	
34	WIRE	1000	PCB	
35	WIRE	1000	PCB	
36	WIRE	1000	PCB	
37	WIRE	1000	PCB	
38	WIRE	1000	PCB	
39	WIRE	1000	PCB	
40	WIRE	1000	PCB	
41	WIRE	1000	PCB	
42	WIRE	1000	PCB	
43	WIRE	1000	PCB	
44	WIRE	1000	PCB	
45	WIRE	1000	PCB	
46	WIRE	1000	PCB	
47	WIRE	1000	PCB	
48	WIRE	1000	PCB	
49	WIRE	1000	PCB	
50	WIRE	1000	PCB	
51	WIRE	1000	PCB	
52	WIRE	1000	PCB	
53	WIRE	1000	PCB	
54	WIRE	1000	PCB	
55	WIRE	1000	PCB	
56	WIRE	1000	PCB	
57	WIRE	1000	PCB	
58	WIRE	1000	PCB	
59	WIRE	1000	PCB	
60	WIRE	1000	PCB	
61	WIRE	1000	PCB	
62	WIRE	1000	PCB	
63	WIRE	1000	PCB	
64	WIRE	1000	PCB	
65	WIRE	1000	PCB	
66	WIRE	1000	PCB	
67	WIRE	1000	PCB	
68	WIRE	1000	PCB	
69	WIRE	1000	PCB	
70	WIRE	1000	PCB	
71	WIRE	1000	PCB	
72	WIRE	1000	PCB	
73	WIRE	1000	PCB	
74	WIRE	1000	PCB	
75	WIRE	1000	PCB	
76	WIRE	1000	PCB	
77	WIRE	1000	PCB	
78	WIRE	1000	PCB	
79	WIRE	1000	PCB	
80	WIRE	1000	PCB	
81	WIRE	1000	PCB	
82	WIRE	1000	PCB	
83	WIRE	1000	PCB	
84	WIRE	1000	PCB	
85	WIRE	1000	PCB	
86	WIRE	1000	PCB	
87	WIRE	1000	PCB	
88	WIRE	1000	PCB	
89	WIRE	1000	PCB	
90	WIRE	1000	PCB	
91	WIRE	1000	PCB	
92	WIRE	1000	PCB	
93	WIRE	1000	PCB	
94	WIRE	1000	PCB	
95	WIRE	1000	PCB	
96	WIRE	1000	PCB	
97	WIRE	1000	PCB	
98	WIRE	1000	PCB	
99	WIRE	1000	PCB	
100	WIRE	1000	PCB	



**PREMIERE IL PULSANTE S25-SW1 PER SCARICARE I CONDENSATORI PRIMA DI RIMUOVERE LA PROTEZIONE. PUSHA BUTTON S25-SW1 TO DISCHARGE CAPACITORS BEFORE REMOVING PROTECTION.**



\* ONLY FOR CONTROL TABLE UNIT

X-RAY TUBE ASSEMBLY

The callbacks refer the same document

Modification Description	Coordinates	Date	Signature	Replaced N°	Replaced by N°
"TR2 MF24V" RISE IMPORTANT: RISE-NOTE	pag.1 - 4/B	23/11/10	V. Perico	10189	
	pag.1 - 6/D				

Drawn:	Date:	Checked:	Date:	Approved:	Date:
V. Perico	30/09/09	R. Doldi	23/11/10	R. Tringali	23/11/10

References	Date:
I-15	
I-30	


Code	Draw
01.02.047.018	10465

PG3/3









Drawn V. Perico	Checked R. Solda	Approved /	Weight (kg) /	Scale /	Code Art. 01.17.008.025
DATE 30-01-2004	DATE 30-01-2004	DATE /			

Questo disegno non può essere riprodotto o reso noto a terzi senza autorizzazione scritta del costruttore. Reproduction or Presentation to Third Persons Prohibited Without Manufacturer Written Authorization

Doc. N° 7548	Foglio N° 1/1
Doc. 7548	Sheet 1/1

Needed quantities for drawn groups	CARTOUCHE A3-E
	NEXT CODE /
	PREVIOUS CODE /
	UNI CEI IEC

Eliminato rinvi e ridotto dimensioni	30-01-04	5684
MODIFICATION DESCRIPTION	Coordinates	Date
1	2	3
4	5	6
7	8	9
10	11	12
13	14	15
16	17	18
19	20	21
22	23	24
25	26	27
28	29	30
31	32	33
34	35	36
37	38	39
40	41	42
43	44	45
46	47	48
49	50	51
52	53	54
55	56	57
58	59	60
61	62	63
64	65	66
67	68	69
70	71	72
73	74	75
76	77	78
79	80	81
82	83	84
85	86	87
88	89	90
91	92	93
94	95	96
97	98	99
100	101	102
103	104	105
106	107	108
109	110	111
112	113	114
115	116	117
118	119	120
121	122	123
124	125	126
127	128	129
130	131	132
133	134	135
136	137	138
139	140	141
142	143	144
145	146	147
148	149	150
151	152	153
154	155	156
157	158	159
160	161	162
163	164	165
166	167	168
169	170	171
172	173	174
175	176	177
178	179	180
181	182	183
184	185	186
187	188	189
190	191	192
193	194	195
196	197	198
199	200	201
202	203	204
205	206	207
208	209	210
211	212	213
214	215	216
217	218	219
220	221	222
223	224	225
226	227	228
229	230	231
232	233	234
235	236	237
238	239	240
241	242	243
244	245	246
247	248	249
250	251	252
253	254	255
256	257	258
259	260	261
262	263	264
265	266	267
268	269	270
271	272	273
274	275	276
277	278	279
280	281	282
283	284	285
286	287	288
289	290	291
292	293	294
295	296	297
298	299	300


Eliminato rinvi e ridotto dimensioni	30-01-04	5684
MODIFICATION DESCRIPTION	Coordinates	Date
1	2	3
4	5	6
7	8	9
10	11	12
13	14	15
16	17	18
19	20	21
22	23	24
25	26	27
28	29	30
31	32	33
34	35	36
37	38	39
40	41	42
43	44	45
46	47	48
49	50	51
52	53	54
55	56	57
58	59	60
61	62	63
64	65	66
67	68	69
70	71	72
73	74	75
76	77	78
79	80	81
82	83	84
85	86	87
88	89	90
91	92	93
94	95	96
97	98	99
100	101	102
103	104	105
106	107	108
109	110	111
112	113	114
115	116	117
118	119	120
121	122	123
124	125	126
127	128	129
130	131	132
133	134	135
136	137	138
139	140	141
142	143	144
145	146	147
148	149	150
151	152	153
154	155	156
157	158	159
160	161	162
163	164	165
166	167	168
169	170	171
172	173	174
175	176	177
178	179	180
181	182	183
184	185	186
187	188	189
190	191	192
193	194	195
196	197	198
199	200	201

Eliminato rinvi e ridotto dimensioni	30-01-04	5684
MODIFICATION DESCRIPTION	Coordinates	Date
1	2	3
4	5	6
7	8	9
10	11	12
13	14	15
16	17	18
19	20	21
22	23	24
25	26	27
28	29	30
31	32	33
34	35	36
37	38	39
40	41	42
43	44	45
46	47	48
49	50	51
52	53	54
55	56	57
58	59	60
61	62	63
64	65	66
67	68	69
70	71	72
73	74	75
76	77	78
79	80	81
82	83	84
85	86	87
88	89	90
91	92	93
94	95	96
97	98	99
100	101	102
103	104	105
106	107	108
109	110	111
112	113	114
115	116	117
118	119	120
121	122	123
124	125	126
127	128	129
130	131	132
133	134	135
136	137	138
139	140	141
142	143	144
145	146	147
148	149	150
151	152	153
154	155	156
157	158	159
160	161	162
163	164	165
166	167	168
169	170	171
172	173	174
175	176	177
178	179	180
181	182	183
184	185	186
187	188	189
190	191	192
193	194	195
196	197	198
199	200	201

Eliminato rinvi e ridotto dimensioni	30-01-04	5684
MODIFICATION DESCRIPTION	Coordinates	Date
1	2	3
4	5	6
7	8	9
10	11	12
13	14	15
16	17	18
19	20	21
22	23	24
25	26	27
28	29	30
31	32	33
34	35	36
37	38	39
40	41	42
43	44	45
46	47	48
49	50	51
52	53	54
55	56	57
58	59	60
61	62	63
64	65	66
67	68	69
70	71	72
73	74	75
76	77	78
79	80	81
82	83	84
85	86	87
88	89	90
91	92	93
94	95	96
97	98	99
100	101	102
103	104	105
106	107	108
109	110	111
112	113	114
115	116	117
118	119	120
121	122	123
124	125	126
127	128	129
130	131	132
133	134	135
136	137	138
139	140	141
142	143	144
145	146	147
148	149	150
151	152	153
154	155	156
157	158	159
160	161	162
163	164	165
166	167	168
169	170	171
172	173	174
175	176	177
178	179	180
181	182	183
184	185	186
187	188	189
190	191	192
193	194	195
196	197	198
199	200	201

Eliminato rinvi e ridotto dimensioni	30-01-04	5684
MODIFICATION DESCRIPTION	Coordinates	Date
1	2	3
4	5	6
7	8	9
10	11	12
13	14	15
16	17	18
19	20	21
22	23	24
25	26	27
28	29	30
31	32	33
34	35	36
37	38	39
40	41	42
43	44	45
46	47	48
49	50	51
52	53	54
55	56	57
58	59	60
61	62	63
64	65	66
67	68	69
70	71	72
73	74	75
76	77	78
79	80	81
82	83	84
85	86	87
88	89	90
91	92	93
94	95	96
97	98	99
100	101	102
103	104	105
106	107	108
109	110	111
112	113	114
115	116	117
118	119	120





FLAT-CABLE

TERMINAL BOARD

S19/C

Modification Description	Coordinates	Date	Signature	Replaced N°	Replaced by N°


References	Date	Checked	Date	Approved	Date

**S19/C Board**  
**E100-30**  
**Inverter Control**

Code: 01.17.008.040  
 Doc: 9686

PG1/2





S19/C

PG2/2

References

Drawn: V. Retto 14-05-09  
 Checked: R. D. Dak 14-05-09  
 Approved: R. T. T. Dak 14-05-09

Code 01.17.008.040  
 Draw 9686


Modification Description Coordinates Date Signature Replaced N' Replaced by N'

Obtain design from your source. Reproduction or presentation without permission is prohibited. Manufacturer's name and address.









Modification Description	Coordinates	Date	Signature	Replaced N°	Replaced by N°
TH1 changed	4B	27-01-10	V. Ferrero	9675	

Drawn	Date	Checked	Date	Approved	Date
V. Ferrero	13-05-09	libe20081	27-01-10	R. Ferrero	27-01-10


**X-Ray Tube Stator Supply S21/B**  
 Cod 01.17.003.006  
 Draw 10192

References  
 Questo disegno non può essere prodotto  
 senza l'autorizzazione formale del  
 progettista o l'autorizzazione  
 esplicita del costruttore.  
 Riproduzione o presentazione in tutto o  
 in parte senza permesso scritto del  
 progettista è vietata.










Part Name	Value	Reference
IC9	68HC88T1	MEM93C88
IC1	CHURRES	MEM93C88
IC2C	74HC14	MEM93C88
OC1A	ILLD206	MEM93C88
OC1B	ILLD206	MEM93C88
IC9A	4093	MEM93C88
IC9B	4093	MEM93C88
IC9C	4093	MEM93C88
IC9D	4093	MEM93C88
IC9E	4093	MEM93C88
IC9F	4093	MEM93C88
IC9G	4093	MEM93C88
IC9H	4093	MEM93C88
IC9I	4093	MEM93C88
IC9J	4093	MEM93C88
IC9K	4093	MEM93C88
IC9L	4093	MEM93C88
IC9M	4093	MEM93C88
IC9N	4093	MEM93C88
IC9O	4093	MEM93C88
IC9P	4093	MEM93C88
IC9Q	4093	MEM93C88
IC9R	4093	MEM93C88
IC9S	4093	MEM93C88
IC9T	4093	MEM93C88
IC9U	4093	MEM93C88
IC9V	4093	MEM93C88
IC9W	4093	MEM93C88
IC9X	4093	MEM93C88
IC9Y	4093	MEM93C88
IC9Z	4093	MEM93C88





S22

PG3/5

Drawn:	Date:	Checked:	Date:	Approved:	Date:
V. Perico	29/05/09	R. Dolci	29/05/09	R. Trapletti	29/05/09

**References**

Questo disegno non può essere prodotto o reso noto a terzi senza autorizzazione scritta del costruttore.  
 Reproduction or Presentation to third persons Prohibited without manufacturer written Authorization

Code 01.17.010.004


Draw 9706

Modification Description	Coordinates	Date	Signature	Replaced N°	Replaced by N°









Drawn:	Date:	Checked:	Date:	Approved:	Date:
V. Perico	29/05/09	R. Dolet	29/05/09	R. Tripletti	29/05/09

**PGS/5**

Questo disegno non può essere prodotto  
 o reso noto a terzi senza l'autorizzazione  
 scritta del costruttore.  
 Reproduction or Presentation to third  
 persons prohibited without Manufacturer  
 written Authorization.



**Analogical  
 Input/Output**  
 Code 01.17.010.004  
 Draw **9706**

Modification Description	Coordinates	Date	Signature	Replaced N°	Replaced by N°









Modification Description	Coordinates	Date	Signature	Replaced N°	Replaced by N°

Drawn:	Date:	Checked:	Date:	Approved:	Date:
V. Perico	03-10-05	Ing. R. Sala	03-10-05	Ing. Fogliata	03-10-05

References	References


  

Code	Doc.
01.17.037.002	8202

PG1/2

S62/A Board





Modification Description	Coordinates	Date	Signature	Replaced N°	Replaced by N°

Drawn:	Date:	Checked:	Date:	Approved:	Date:
V. Perico	03-10-05	Ing. R. Sala	03-10-05	Ing. Fogliata	03-10-05

**References**

Questo disegno non può essere prodotto o reso noto a terzi senza l'autorizzazione scritta del costruttore. Reproduction or Presentation to third persons Prohibited without Manufacturer written Authorization.


**S62/A Board**

Code **01.17.037.002**

Doc. **8202**

**PG2/2**





\* Montaggio opzionale  
 Montare dissipatore su regolatori  
 Sotto gli opto prevedere cava di almeno  
 2mm per tutta la lunghezza dell'OC

S81

- MH1
- MH2
- MH3
- MH4

PGI/1

Modifications	Description	Coordinates	Date	Signature	Replaced N°	Replaced by N°


  

Dynamic	Date	Checked	Date	References
V. Perico	21/04/08	D. Minniti	10/04/09	10/04/09

Code	Code
01.17.081.001	9331





Drawn	M. Farina	Checked	C. Bonazzi	Approved	A. Mariani	Weight (kg)	Scale	
DATE	01-01-1999	DATE	01-01-1999	DATE	01-01-1999	Code Art.	01.17.005.005	


Questo disegno non può essere riprodotto o reso noto a terzi senza l'autorizzazione scritta del Costruttore. Reproduction or Presentation to Third Persons Prohibited Without Manufactured Written Authorization

TOPOGRAFIC  
S14 CARD  
LINE FILTER

MODIFICATION DESCRIPTION	Coordinates	Date	Signature	Replaced Draw N°	Replaced by Draw N°
1	2	3	4	5	6

Doc. N° DOC.	5040	Foglio N° Sheet	1/1	UNI	CEI	IEC
				PREVIOUS CODE	NEXT CODE	
Needed quantities for drawn groups				CARTOUCHÉ	A3-E	





Drawn		Checked		Approved		Weight (kg)	Scale			Needed quantities for drawn groups UNI CEI IEC PREVIOUS CODE NEXT CODE	
V. Perico		Ing. R. Dolci R. Trapletti				/	/			CARTOUCHE A3-E	
DATE	DATE	DATE	DATE	Code Art.							
13-05-09	13-05-09	13-05-09	13-05-09	01.17.001.029							
Questo disegno non può essere riprodotto o reso noto a terzi senza l'autorizzazione scritta del Costruttore.				Reproduction or Presentation to Third Persons Prohibited Without Manufacturer Written Authorization				Filament 400HZ Basic S15/B BOARD			
MODIFICATION DESCRIPTION		Coordinates	Date	Signature	Replaced Draw N°	Replaced by Draw N°					
1		2	3	4	5	6	7	8			



1	2	3	4
DESCRIZIONE MODIFICA	Coordinate	Data	Firma
Sost. il dis. N°	Sost. dal dis. N°		

Aggiornato cartiglio 30-01-04 5685

A

A

B


B

C

C

D

D




Pos.	Provenienza	Modello	Codice	Lung.(cm)	Colore	Ø Est. mm	Codice	Destinazione	Modello	Codice	
	Capicorda			Toll.+1cm					Capicorda		
ESTREMITA' DI COLLEGAMENTO				MEZZO TRASMISSIVO / CAVO				ESTREMITA' DI COLLEGAMENTO			
Disegnato	Controllato	Approvato	Peso (kg)	Scala				VEDI IMPLOSIONI			
V. Perico	M. Gotti	/	/	/							
DATA	DATA	DATA	Codice Art.	01.17.008.025			UNI				CEI
30-01-2004	30-01-2004	/	Code				COD. PRECEDENTE				
<b>international medical devices S.p.a.</b> via A. Moro 5/7 Scanzorosciate (Bg) ITALIA Tel. 035-668163/4 fax 035-668166 				<b>Monoblocco</b> <b>Scheda S17/A</b> <b>Monoblock</b> <b>S17/A Board</b>				Quant. necessaria per gruppi a dis.		CARTIGLIO A4-E	
								Dis. N° Draw		7547	

Questo disegno non può essere riprodotto o reso noto a terzi senza l'autorizzazione scritta della i.m.d Spa  
 Reproduction or Presentation to Third Persons Prohibited Without Written Authorization From i.m.d Spa









Drawn		Checked		Approved		Weight (kg)	Scale	Code Art.	01.17.006.003	Needed quantities for drawn groups		UNI	CEI	IEC	
V. Perico 12-05-09		Ing. Dolci 12-05-09		R. Trapletti 12-05-09						SCHEMA S20/A					CARTOUCHE A3-E
DATE		DATE		DATE		S20/A BOARD		SCHEDA S20/A		Doc. N° 9669		Foglio N° 1/1			
MODIFICATION DESCRIPTION		Coordinates		Date		Signature		Replaced Draw N°		Replaced by Draw N°		UNITS		PREVIOUS CODE	


Questo disegno non può essere riprodotto o reso noto a terzi senza l'autorizzazione scritta del Costruttore. Reproduction or Presentation to Third Persons Prohibited Without Manufacturer Written Authorization.

RA









Drawn	Checked	Approved	Weight (kg)	Scale
V. Perico	R. Dolci	R. Trapletti	/	1:1
DATE	DATE	DATE	Code Art.	
29-05-2009	29-05-2009	29-05-2009	01.17.010.004	


**S22UP Base Board**

MODIFICATION DESCRIPTION	Coordinates	Date	Signature	Replaced Draw N°	Replaced by Draw N°
1					
2					
3					
4					
5					
6					
7					

Doc. N° DOC. 9705	Foglio N° 1/1	UNI	CEI	IEC
		PREVIOUS CODE		
Needed quantities for drawn groups		CARTOUCHE A3-E		
NEXT CODE				

Questo disegno non può essere riprodotto o reso noto a terzi senza l'autorizzazione scritta del Costruttore  
 Reproduction or Presentation to Third Persons  
 Prohibited Without Manufacturer Written Authorization





Drawn V. Perico	Checked R. Sola	Approved /	Weight (kg) /	Scale /	
DATE 15-04-03	DATE 15-04-03	DATE /	Code Art. 01.17.005.009		

S25 BOARD  
SCHEDA S25


Needed quantities for drawn groups	01.02.005.071	CARTUCHE A3-E
	01.02.025.031	
	01.02.047.018	
		NEXT CODE
		PREVIOUS CODE
		UNI CEI IEC

Questo disegno non può essere riprodotto o reso noto a terzi senza l'autorizzazione scritta del Costruttore. Reproduction or Presentation to Third Persons Prohibited Without Manufacturer Written Authorization.

Doc. N° 7613 Foglio N° 1/1  
D.C.

MODIFICATION DESCRIPTION	Coordinates	Date	Signature	Replaced Draw N°	Replaced by Draw N°
1					
2					
3					
4					
5					
6					
7					
8					






CHANGED DISPLAY 20x1 BT12009 WITH									
THE NEW DISPLAY YC-2001A-SBYA-1-B									
MODIFICATION DESCRIPTION	Coordinates	Date	Signature	Replaced Draw N°	Replaced by Draw N°				

Drawn V. Perico	Checked R. Sala	Approved E. Fogliata	Weight (kg) /	Scale /		Needed quantities for drawn groups	UNI CEI IEC
DATE 03-10-05	DATE 03-10-05	DATE 03-10-05	Code Art. 01.17.037.002				
S62/A BOARD SCHEDA S62/A							
Questo disegno non può essere riprodotto o reso noto a terzi senza l'autorizzazione scritta del costruttore. Reproduction or Presentation to Third Persons Prohibited Without Written Manufacturer Authorization.							
						Doc. N° DDC: 8203	Foglio N° Sheet 1/1





**PG1/1**

**References**

<b>Drawn:</b>	<b>Date:</b>	<b>Checked:</b>	<b>Date:</b>	<b>Approved:</b>	<b>Date:</b>
V. Perico	27-01-09	D. Mainetti	10-04-09	R. Trapletti	10-04-09

Questo disegno non può essere prodotto o reso noto a terzi senza l'autorizzazione scritta del costruttore.  
 Reproduction or Presentation to third persons Prohibited without manufacturer written Authorization.

**S81 Driver Board layout**

**Code 01.17.081.001**

**Draw 9574**

<b>Modification Description</b>	<b>Coordinates</b>	<b>Date</b>	<b>Signature</b>	<b>Replaced N°</b>	<b>Replaced by N°</b>




115VAC  
POWER supply

## 1 THERMOMAGNETIC CIRCUIT-BREAKER

The 115V power supply circuit (SEE FIGURE 1) is protected by a thermomagnetic circuit-breaker (1); should this circuit-breaker trip, disassembly all the unit casings, as explained in the SERVICE MANUAL – MAINTENANCE SECTION, and reset the circuit-breaker.



**Any trip of the thermomagnetic circuit-breaker immediately after resetting can be an indication of incorrect operation of the system: contact the technical service assistance immediately.**



**FIGURE 1 - THERMOMAGNETIC CIRCUIT-BREAKER**

## 2 COMPONENTS REPLACEMENT

For the 115V power supply assembly dismantling and relative list of spare parts make reference to the EXPLODED DRAWINGS.

*NOTE: When you require spare part, it is necessary to tell the code of the piece and the **serial number of the unit** too.*

## 3 INDEX DRAWINGS


### ELECTRICAL

- 9250      Supply AUTOTRANSFORMER 115VAC


### Exploded

- 7856      115VAC POWER Supply ASSEMBLY

White Page



**MAINS**  
Gr. cavo  
Cod. 01.02.003.005  
(El. 100-15/30)



**GND** (01.02.047.018)  
**Morsetto di Terra su Gr. El. 100**  
**0V [SW1-3]** (01.02.047.018)

**References**  
Drawn: V. Perico 28-01-08  
Checked: R. Doldi 28-01-08  
Approved: E. Fogliata 28-01-08


**SUPPLY**  
**AUTOTRANSFORMER**  
**0V-115V/230V**  
**800 VA**

Code 01.02.040.003  
Draw **9250**

Questo disegno non può essere prodotto o reso noto a terzi senza l'autorizzazione scritta del costruttore.  
Reproduction or Presentation to third persons is prohibited without manufacturer written Authorization.

YG cable added				
Modification Description	Coordinates	Date	Signature	Replaced N°
		28-01-08	V. Perico	7662
				Replaced by N°





N°	DESCRIPTION	CODE	NOTE
<b>115V POWER SUPPLY ASSEMBLY - 15/30KW MODEL</b>			
F10	5x20 2A 250V delayed fuse	01.02.040.003 49.06.022.013	
K1	B6-30-01 220V contactor	49.06.015.016	
R1	22E 50W potentiometric resistor	49.06.011.094	
TC1	SM2070 N.C. bimetallic thermostat	49.06.006.003	
TR	115-230 800VA tor. autotransformer	01.05.001.023	

Drawing N° 01.02.040.003 - 007841