


# **HELIODENT**PLUS

# **Installation Requirements**



# Table of contents

1	Struc	cture of the document	4
	1.1	Identification of the danger levels	4
	1.2	Formats and symbols used	4
2 Sa	Safe	ty instructions	5
	2.1	Shielding of room	5
	2.2	Electromagnetic compatibility	5
	2.3	Modifications and extensions of the system	5
3 F	Prior	to installation	6
	3.1	Installation options	6
	3.2	On-site installation	9
4	Dime	ensions, technical data	12
	4.1	Dimensions with round support arm	12
		4.1.1 Dimensions of front view with all options	12
		4.1.2 Dimensions for 950 mm (37 3/8") support arm	13
		4.1.3 Dimensions for 700 mm (27 1/2") support arm	16
		4.1.4 Dimensions for 410 mm (16 1/8") support arm	19
		4.1.5 Dimensions mobile stand	22
	4.2	Dimensions with angular support arm	24
		4.2.1 Dimensions of front view with all options	24
		4.2.2 Dimensions for 910 mm (35 7/8") support arm	25
		4.2.3 Dimensions for 660 mm (26") support arm	28
		4.2.4 Dimensions for 370 mm (14 1/2") support arm	31
		4.2.5 Dimensions mobile stand	34
	4.3	Technical data	36
5	Elect	romagnetic compatibility	37
	5.1	Accessories	37
	5.2	Electromagnetic emission	37
	5.3	Immunity to interference	38
5.4	5.4	Working clearances	41

# Structure of the document

# 1.1 Identification of the danger levels

To prevent personal injury and material damage, please observe the warning and safety information provided in these operating instructions. Such information is highlighted as follows:

# DANGER

An imminent danger that could result in serious bodily injury or death.

### **WARNING**

A possibly dangerous situation that could result in serious bodily injury or death.

### **⚠** CAUTION

A possibly dangerous situation that could result in slight bodily injury.

#### NOTICE

A possibly harmful situation which could lead to damage of the product or an object in its environment.

#### **IMPORTANT**

Application instructions and other important information.

Tip: Information on making work easier.

# 1.2 Formats and symbols used

The formats and symbols used in this document have the following meaning:

✓ Prerequisite	Prompts you to do something.
1. First action step	
2. Second action step	
or	
Alternative action	
♥ Result	
➤ Individual action step	
See "Formats and symbols used [ → 4]"	Identifies a reference to another text passage and specifies its page number.
• List	Designates a list.
"Command / menu item"	Indicates commands, menu items or quotations.

# Safety instructions

# 2.1 Shielding of room

When using the HELIODENT Plus X-ray tube assembly, proper shielding of the room and operator position is essential.

It is the installer's responsibility to ensure that all local radiation regulations and safety measures are met.

# 2.2 Electromagnetic compatibility

The unit should not be operated in the immediate vicinity of other devices. If this proves to be unavoidable, the unit should be monitored to ensure that it is operating properly.

## 2.3 Modifications and extensions of the system

Modifications to this unit which might affect the safety of the system owner, patients or other persons are prohibited by law.

For reasons of product safety, this product may be operated only with original Sirona accessories or third-party accessories expressly approved by Sirona. The user assumes the risk of using non-approved accessories.

# Prior to installation

# 3.1 Installation options


#### Designations for release buttons and door contact

- Manual release S3
  - Coiled cable
- Release key on the control membrane S4
  - Directly connected to control board DX4
- Remote control release key S9
  - Integrated in remote control housing
- Door contact (safety circuit) S7

#### Installation option 1

Release in the treatment room without remote control


- Release
  - Manual release \$3



#### Installation option 2

Release in the treatment room with remote control


- Release
  - Manual release S3
    - or
  - Remote control release key S9



#### Installation option 3

Release in the treatment room with Remote Timer


- Release
  - Manual release S3
    - or
  - Release key on the control membrane S4




#### **NOTICE**

Length of cable supplied for Remote Timer approx. 10 meters (393") (must not be extended).

Conduit int. dia. at least 12 mm (1/2").



# 5 S<sub>a</sub> s3



#### Installation option 4

Release outside of the X-ray room with remote control

- Release
  - Manual release S3 or
  - Remote control release key S9

#### NOTICE

#### Installation prerequisites

Use of the remote control is permissible only if the yellow X-Ray LED is visible to the operating personnel during radiation release.

#### Installation option 5

Release outside of the X-ray room with Remote Timer

- Release
  - Manual release S3
    - or
  - Release key on the control membrane S4

#### **NOTICE**

Length of cable supplied for Remote Timer approx. 10 meters (393") (must not be extended).

Conduit int. dia. at least 12 mm (1/2").

#### Installation option 6


Release outside of the X-ray room **with** remote control, door contact safety circuit

- Door contact
  - Door contact \$7 wired to the wall adapter
- Release
  - Manual release S3
    - or
  - Remote control release key S9

#### **NOTICE**

#### Installation requirement

Use of the remote control is permissible only if the yellow X-Ray LED is visible to the operating personnel during radiation release.



#### Installation option 6.1

Release outside of the X-ray room **with** remote control, door contact safety circuit

- Door contact
  - Door contact **S7** wired to the remote control housing
- Release
  - Manual release S3

or

- Remote control release key S9

#### **NOTICE**

#### Installation requirement

Use of the remote control is permissible only if the yellow X-Ray LED is visible to the operating personnel during radiation release.

#### Installation option 7


Release outside of the X-ray room **with** Remote Timer, door contact safety circuit


- Door contact
  - Door contact S7
- Release
  - Release key on the control membrane S4

#### Installation option 8

Release outside of the X-ray room **with** Remote Timer, door contact safety circuit

- Door contact
  - Door contact \$7 wired to Remote Timer
- Release
  - Manual release \$3






### 3.2 On-site installation


# **CAUTION**

#### Observe wall properties

In installation situations, the technician is responsible for the assessment of wall properties and selecting the method of attaching the unit to the wall.



- K Tensile force per screw 3600 N (800 lbf) if L ≤ 700 mm (27 1/2") 4200N(950 lbf) if L ≤ 950 mm (37 3/8")
   L Length of support arm
   M Mounting plate (supplied)
   N Anchor plate
   O Threaded bolt M8
- The permissible tensile force of the selected attachment must at least equal the tensile force listed above.
- Matching wood screws for wooden beams are included in delivery.
- For all other wall structures, special wall anchors must be purchased from a selected dealer. The wall anchors and screws should be identical for every attachment point.
- Alternatively, an anchor plate can be used as a counter bearing. In this case, M8 threaded rods of the appropriate length for the wall (thickness of the wall + 2 x mounting plate thickness + attachment material) are required.



#### Re-use installation sites of old units

It is possible to conceal the installation site of an old unit when installing a HELIODENTPLUS.

- For the replacement of vertically mounted old units (e.g. HELIODENT DS, HELIODENT MD, Planmeca Intra<sup>1</sup>) an adapter plate is available for this purpose, REF 62 42 254.
- The drill holes of some vertically mounted units (e.g. Progeny Previa <sup>1</sup>, Gendex 765DC <sup>1</sup>) coincide with the dimensions of the drill holes of the Heliodent PLUS. No adapter plate is required.

#### NOTICE

Regardless of their prior use, the existing drill holes and wall plugs must comply with the installation regulations and must be checked by the person performing installation.

#### **NOTICE**

The different connection areas of the old units make it necessary to relocate the existing electrical connections (e.g. concealed installation) on-site.

#### **CAUTION**

The on-site electrical installation must be performed according to the valid regulations for medical electrical equipment (DIN VDE 0100-710).

- Cable for remote control or Remote Timer: Conduit ∅ int. min. 12mm (1/2"), requires an excess length of at least 0.25 m (10") at both ends.
- Power cable 3x1.5 mm<sup>2</sup> (AWG 16); required excess length for concealed installation: 0.25m (10").

#### CAUTION

Do not install the cables for Remote Timer and power cables in the same conduit.

#### NOTICE

The Heliodent<sup>Plus</sup> wall model is suitable for fixed installation only.

# DANGER

#### Perilous shock hazard!

#### Fixed connection!

Installing a mains plug instead of the specified fixed connection infringes international medical regulatory actions and is prohibited. In case of error, this puts patients, users, and other parties seriously at risk.

<sup>1.</sup> The product names mentioned may be copyrighted by their respective own-

# **CAUTION**

#### Observe the permissible nominal voltage range!

The unit can be connected to 120 V (1-phase connection) or to 200 - 240 V (1- or 2-phase connection),

for all other voltages a pre-transformer is required. As a ceiling model or device model,

the HELIODENTPLUS must

only be connected to 200 - 240 V (1 or 2-phase connection).

#### For the USA only


#### Power supply:

A separate three wire grounded circuit connected directly to the central distribution panel with an over-current protection rated for 20 amperes should be used.

# 4 Dimensions, technical data

# 4.1 Dimensions with round support arm

### 4.1.1 Dimensions of front view with all options




Α	Recommended installation height for the wall module
В	Cable bushing for network cable
С	Unit height
D	Ceiling height
Е	Wooden beam
F	Wall module cover
G	Cable bushing for remote control or Remote Timer

Н	Recommended installation height of remote control or Remote Timer
I	Remote control or Remote Timer
J	Cover for remote control or Remote Timer


# 4.1.2 Dimensions for 950 mm (37 3/8") support arm

#### Side view


Α	Recommended installation height: 1110 mm (43 3/4")
Т	X-ray tube assembly with standard tube, 200 mm (8") SSD




#### Top view




# Minimum dimensions for X-ray rooms with 950 mm (37 3/8") support arm







# 4.1.3 Dimensions for 700 mm (27 1/2") support arm Side view


Α	Recommended installation height: 1110 mm (43 3/4")
Т	X-ray tube assembly with standard tube, 200 mm (8") SSD





# Minimum dimensions for X-ray rooms with 700 mm (27 1/2") support arm







## 4.1.4 Dimensions for 410 mm (16 1/8") support arm


#### Side view


Α	Recommended installation height: 1110 mm (43 3/4")
Т	X-ray tube assembly with standard tube, 200 mm (8") SSD






# Minimum dimensions for X-ray rooms with 410 mm (16 1/8") support arm






### 4.1.5 Dimensions mobile stand

#### Front view




#### Top view



# 4.2 Dimensions with angular support arm

## 4.2.1 Dimensions of front view with all options




Α	Recommended installation height for the wall module
В	Cable bushing for network cable
С	Unit height
D	Ceiling height
Е	Wooden beam
F	Wall module cover
G	Cable bushing for remote control or Remote Timer
Н	Recommended installation height of remote control or Remote Timer

I	Remote control or Remote Timer
J	Cover for remote control or Remote Timer


# 4.2.2 Dimensions for 910 mm (35 7/8") support arm

#### Side view


Α	Recommended installation height: 1110 mm (43 3/4")
Т	X-ray tube assembly with standard tube, 200 mm (8") SSD



#### Top view



# Minimum dimensions for X-ray rooms with 910 mm (35 7/8") support arm





min. 1500 59"

# 4.2.3 Dimensions for 660 mm (26") support arm


#### Side view


Α	Recommended installation height: 1110 mm (43 3/4")
Т	X-ray tube assembly with standard tube, 200 mm (8") SSD





# Minimum dimensions for X-ray rooms with 660 mm (26") support arm







## 4.2.4 Dimensions for 370 mm (14 1/2") support arm


#### Side view


Α	Recommended installation height: 1110 mm (43 3/4")
Т	X-ray tube assembly with standard tube, 200 mm (8") SSD






# Minimum dimensions for X-ray rooms with 370 mm (14 1/2") support arm






# 4.2.5 Dimensions mobile stand Front view



### Top view



# 4.3 Technical data

Dimensions of the packaging	
HELIODENT <sup>PLUS</sup> Round support arm system	87 cm x 91 cm x 29 cm 34 1/4" x 35 7/8" x 11 1/2"
Angular support arm system	99 cm x 94 cm x 29 cm 39" x 37" x 11 1/2"
Mobile stand	120 cm x 80 cm x 40 cm 47 1/4" x 31 1/2" x 15 3/4"
Weight	incl. / without packaging
HELIODENT <sup>PLUS</sup> Round support arm system	31 kg / 24 kg 68.34 lb / 53 lb
Angular support arm system	37 kg / 30 kg 81.57 lb / 66.14
Mobile stand	128 kg / 105 kg 282.2 lb / 231.5 lb
Power supply connection	
Line voltage wall model ceiling model, unit model	120 V, 200 V - 240 V, 50 / 60 Hz 200 V - 240 V, 50 / 60 Hz
Line voltage variation	± 10%
Internal line resistance	At 120 V: 0.3 Ω for 200 V - 240 V: 0.8 Ω
Nominal current	At 120 V: 10 A for 200 V - 240 V: 6 A - 5 A
Fuse	16 A slow blow
Power consumption	≤ 1.2 kW

Transport and operating conditions:			
Transport and storage temperature:	-40°C - +70°C (-40°F - 158°F)		
Air humidity:	10% – 95%		
Operating conditions as specified in IEC 60601-1:	Ambient temperature +10 °C - +40 °C (50 °F - 104 °F)		
	Relative humidity: 30% – 75%		
Recommended operating temperature:	18°C - 35°C (64°F - 95°F) With room temperatures > 35°C (> 95°F) Sirona recommends the use of an air conditioning system.		
Operating altitude:	≤ 3000 m		

# 5 Electromagnetic compatibility

#### NOTICE

HELIODENT<sup>Plus</sup> complies with the requirements for electromagnetic compatibility (EMC) according to IEC 60601-1-2.

HELIODENT<sup>Plus</sup> is referred to in the following as "UNIT". Observance of the following information is necessary to ensure safe operation regarding EMC aspects.

#### 5.1 Accessories

Designation of the interface cables	Order no.
LIYCY 2x0.25mm <sup>2</sup> (AWG 24) remote cable L9 for remote release, 10m	62 42 064
LIYYC 8x0.22mm <sup>2</sup> (AWG 24) remote cable L6 for Remote Timer, 10 m	62 42 056
3x1.5mm <sup>2</sup> NYM	Commercially available

- The UNIT may be operated only with accessories and spare parts approved by Sirona. Unapproved accessories and spare parts may lead to an increased emission or to a reduced immunity to interference.
- The UNIT should not be operated immediately adjacent to other devices. If this proves to be unavoidable, the UNIT should be monitored to ensure that it is operating properly.

# 5.2 Electromagnetic emission

The **UNIT** is intended for operation in the electromagnetic environment specified below.

The customer or user of the **UNIT** should make sure that it is used in such an environment.

Emission measurement	Conformity	Electromagnetic environment - guidelines
RF emissions according to CISPR 11	Group 1	The <b>UNIT</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions according to CISPR 11	Class B	The <b>UNIT</b> is intended for use in all facilities,
Harmonics according to IEC 61000-3-2	Class A	including residential areas and in any facilities connected directly to a public power supply providing electricity to buildings used for residential
Voltage fluctuations / flicker according to IEC 61000-3-3	coincides	purposes.

# 5.3 Immunity to interference

The **UNIT** is intended for operation in the electromagnetic environment specified below.

The customer or user of the **UNIT** should make sure that it is used in such an environment.

Interference immunity tests	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be made of wood or concrete or finished with ceramic tiling. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst according to IEC 61000-4-4	± 1 kV for input and output lines ± 2 kV for power supply lines	± 1 kV for input and output lines ± 2 kV for power supply lines	The quality of the line power supply should be that of a typical commercial or hospital environment.
Surge voltages according to IEC 61000-4-5	± 1 kV differential mode voltage ± 2 kV common mode voltage	± 1 kV differential mode voltage ± 2 kV common mode voltage	The quality of the line power supply should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and variations of the power	Voltage dips: 0% U <sub>T</sub> with 1/2 period at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°		The quality of the line power supply should be that of a typical commercial or hospital environment.
supply according to IEC 61000-4-11	0% U <sub>T</sub> with 1 period and 70% U <sub>T</sub> with 25 periods at 50 Hz and 30 periods at 60 Hz each at 0°  Short interruptions: 0% U <sub>T</sub> with 250 periods at 50 Hz and 300 periods at 60 Hz		If the user of the <b>UNIT</b> requires it to continue functioning following interruptions of the power supply, it is recommended to have the <b>UNIT</b> powered by an uninterruptible power supply or a battery.
Magnetic field of power frequencies (50/60 Hz) according to IEC 61000-4-8	30 A/m 30 A/m		Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U <sub>T</sub> is the AC supply	voltage prior to application	on of the test level.	
			Portable and mobile radio equipment must not be used within the recommended working clearance from the UNIT and its cables, which is calculated based on the equation suitable for the relevant transmission frequency.  Recommended working clearance:

Interference immunity tests	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidelines
Conducted RF interference IEC 61000-4-6	3 V <sub>eff</sub> 150 kHz to 80 MHz <sup>1</sup>	3 V <sub>eff</sub>	d= [1.2] √P
Radiated RF interference IEC 61000-4-3	3 V/m 80 MHz to 800 MHz <sup>1</sup>	3 V <sub>eff</sub>	d= [1.2] √P at 80 MHz to 800 MHz
	3 V/m 800 MHz to 2.7 GHz <sup>1</sup>	3 V <sub>eff</sub>	$d= [2.3] \sqrt{p}$ at 800 MHz to 2.7 GHz
			with $P$ as the power rating of the transmitter in watts (W) according to the transmitter manufacturer's specifications and $d$ as recommended safety distance in meters (m).
			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>2</sup> should be less than the compliance level <sup>3</sup> in each frequency range.
			Interference is possible in the vicinity of equipment bearing the following
			graphic symbol.

- 1. The higher frequency range applies at 80 MHz and 800 MHz.
- 2. The field strengths of fixed transmitters, such as base stations of radiotelephones and mobile agricultural radio broadcast services, amateur radio stations, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. A site survey is recommended to assess the electromagnetic environment due to fixed RF transmitters. If the measured field strength in the location in which the UNIT is used exceeds the applicable RF compliance level above, the UNIT should be observed to verify normal operation. If unusual performance characteristics are observed, it may be necessary to take additional measures such as reorientation or repositioning of the UNIT.
- 3. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Test frequency (MHz)	Frequency band <sup>(a)</sup> (MHz)	Radio service (a)	Modulation (b)	Power max. (W)	Distance (m)	Immunity test level (V/m)
385	380 - 390	TETRA 400	Pulse modulation (b) 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM(c) ± 5 kHz stroke 1 kHz sinus	2	0.3	28
710 745 780	704 - 787	LTE band 13, 17	Pulse modulation (b) 217 Hz	0.2	0.3	9
810 870 930	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE band 5	Pulse modulation (b) 18 Hz	2	0.3	28
1720 1845 1970	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation (b) 217 Hz	2	0.3	28
2450	2400 - 2570	Bluetooth, WLAN, 802,11 b/g/n, RFID 2450, LTE band 7	Pulse modulation (b) 217 Hz	2	0.3	28
5240 5500 5785	5100 - 5800	WLAN 802.11 a/n	Pulse modulation (b) 217 Hz	0.2	0.3	9

#### Note:

If necessary, the distance between the transmitting antenna and the ME unit or ME system can be reduced to 1 m in order to achieve the immunity test level. The 1 m test distance conforms to IEC 61000-4-3.

- (a) For some radio services, only the frequencies for the radio link between the mobile communication device and the base station (uplink) have been recorded in the table.
- (b) The carrier must be modulated with a square-wave signal with 50% duty cycle.
- (c) Alternatively to the frequency modulation (FM), a pulse modulation with 50% duty cycle at 18 Hz can be used, as this, if not the actual modulation, would represent the worst case.

# 5.4 Working clearances

Recommended working clearances between portable and mobile HF communication devices and the UNIT The **UNIT** is intended for operation in an electromagnetic environment, where radiated RF interference is checked. The customer or the user of the **UNIT** can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **UNIT** - depending on the maximum output power of the communication device, as shown below.

Power rating of the transmitter	Working clearance according to transmission frequency [m]			
[W]	150 kHz - 80 MHz	80 MHz - 800 MHz	800 MHz - 2.5 GHz	
	d= [1,2] √P	d= [1,2] √P	d= [2,3] √P	
0.01	0.12	0.12	0.23	
0.1	0.38	0.38	0.73	
1	1.2	1.2	2.3	
10	3.8	3.8	7.3	
100	12	12	23	

The recommended safety distance  $\operatorname{\vec{d}}$  in meters (m) can be determined for transmitters, whose maximum power rating is not specified in the above table, using the equation that belongs to the corresponding column, wherein  $\operatorname{\vec{P}}$  is the maximum power rating of the transmitter in watts (W) according to the transmitter manufacturer.

#### Note 1

The higher frequency range applies at 80 MHz and 800 MHz.

#### Note 2

These guidelines may not apply in all cases. The propagation of electromagnetic waves is influenced by their absorption and reflection by buildings, objects and persons.

We reserve the right to make any alterations which may be required due to technical improvements.

© Sirona Dental Systems GmbH 2016 D3507.021.01.06.02 07.2016 Sprache: englisch Ä.-Nr.: 122 391

Printed in Germany

#### **Sirona Dental Systems GmbH**



Contact in the USA:

Sirona Dental, Inc. 4835 Sirona Drive Charlotte, NC 28273 USA

Order No

62 15 037 D3507