Neu ab:
New since:
Nouveau depuis:
Nuevo desde:

07.2012

sirona.
The Dental Company



ORTHOPHOS XG 5 / Ceph ORTHOPHOS XG 3

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Installation Requirements

Conditions d'installation

Requisitos de instalación

**ACHTUNG**

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**ATTENTION**

Proper shielding of room and operator position is essential.

Since these requirements vary from state to state it is the assembler's / installer's responsibility that all local radiation safety requirements are met.

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New since:

07.2012

Modification compared with last edition:
12.2010

Chapter or paragraph	page
Principle of On-site Installation.....	7

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List of Contents

Preparations

Possibilities of Installation	5
Principle of On-site Installation.....	7
Emergency Stop (if required by law)	9
On-site Installation for PC/Networks	11
For USA and Canada	12

Dimensions, technical data

ORTHOPHOS XG 5, ORTHOPHOS XG 3 1:20..	13
ORTHOPHOS XG 5 Ceph 1:20	15
Technical Data	16

Electromagnetic compatibility

Accessories	22
Electromagnetic emission	22
Immunity to interference.....	23
Working clearances	25

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
With remote control the emergency stop switch must be installed at a distance of < 1m from the control holder.

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
The media converter is necessary if no fiber optic network with SC connectors is present.

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For USA and Canada



Wire Size for Power Line

- The unit is designed to operate on a **nominal** 200 or 240VAC line.
Permitted line voltage variation $\pm 10\%$.
On request, the local Electrical Utility Company will perform a voltage regulation test to verify the line quality.
- The distributor box should be installed in the position as shown on page 6.
- To assure proper line quality, a separate three-core **grounded** power cable connected directly to the central distribution panel with an overcurrent circuit breaker rated for 25A must be used.
For an on-site installation with 14 AWG (3x2,5mm²) and an overcurrent circuit breaker rated for 20A, it is permissible to connect only the ORTHOPHOS XG or other such units that cause no danger to the patients or to the computer systems in case the automatic circuit breaker is activated.
- The line voltage drop in the power supply circuit from the central distribution panel to the distributor box depends on length and size of wire.
Measure the distance from the central distribution panel to the distributor box and select the correct wire size, see chart.

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


Technical Data

Dimensions packaging

ORTHOPHOS XG	199cm x 69cm x 122cm	(78 3/8" x 27 1/8" x 48")
Cephalometer	175cm x 78cm x 73cm	(68 7/8" x 30 3/4" x 28 3/4")
Floor stand	114cm x 105cm x 22cm	(56 3/4" x 41 3/8" x 8 5/8")

Weight

ORTHOPHOS XG	177kg / 99kg	including /without packaging (1kg=2.2lbs)
Cephalometer	40kg / 33kg	(390lb / 218lb)
Floor stand	50kg / 31kg	(88lb / 73lb)

Power supply

Line voltage	200V- 240V, 50 / 60Hz
Tolerance of line voltage	±10%
Power line resistance	max. 0,8Ω
Nominal current / Fuse	max. 12A / B 25A inert; with single connection: B 16A/20A inert
Power consumption	max. 2,0kW

Required transformer with 100V / 110V / 125V

Output	230V
Power	2kVA (permanent)
Maximal voltage breakdown with 10A ohmical load:	≤ 10%

Ambient conditions

Operating temperature	+10°C – +40°C (50°F – 104°F)
Transport and storage temperature	-10°C – +70°C (14°F – 158°F)
Air humidity	10% – 95%

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Electromagnetic compatibility

NOTE

The ORTHOPHOS XG 5 / Ceph, ORTHOPHOS XG 3 fulfills all requirements for electromagnetic compatibility (EMC) compliant with IEC 60601-1-2:2001.

The ORTHOPHOS XG 5 / Ceph, ORTHOPHOS XG 3 is referred to as "**UNIT**" in the following.

Observance of the following information is necessary to ensure safe operation regarding EMC aspects.

Accessories

Designation of interface cables	Supplier
PC as peripheral device.	Siemens Fujitsu
Remote cable L17/L117 XG, 15m	LEONI
Cable L25 OP-XG, 5m	EFB-Elektronik
Media converter	TTL-Network
LAN-cable Kat5, 3m	51 68 963 D3348
2nd protective ground wire, 1.5mm ²	58 72 648 D3285

- The **UNIT** may be operated only with accessories and spare parts approved by Sirona. Unapproved accessories and spare parts may lead to an increased emission of or a reduced immunity to interference.
- The **UNIT** should not be operated immediately adjacent to other devices. If this proves to be unavoidable, the **UNIT** should be monitored to check and make sure that it is used properly.

Electromagnetic emission

The **UNIT** is intended for operation in the electromagnetic environment specified below.

The customer or user of the **UNIT** should make sure that it is used in such an environment.

Emission measurement	Conformity	Electromagnetic environment guidelines
HF emission according to CISPR 11	Group 1	The UNIT uses HF energy only for its internal function. The HF emission is therefore very low, and it is improbable that nearby electronic devices might be disturbed.
HF emission according to CISPR 11	Class B	The UNIT is intended for use in all facilities, including residential areas and in any facilities connected directly to a public power supply providing electricity to buildings used for residential purposes.
Harmonics according to IEC 61000-3-2	Class A	
Voltage fluctuations / Flicker according to IEC 61000-3-3	compliant	

Immunity to interference

The **UNIT** is intended for operation in the electromagnetic environment specified below.

The customer or user of the **UNIT** should make sure that it is used in such an environment.

Immunity interference tests	IEC 60601-1-2 test level	Conformance level	Electromagnetic environment guidelines
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 6kV contact discharge ± 8 kV air discharge	± 6kV contact discharge ± 8kV air discharge	Floors should be made of wood or concrete or covered with ceramic tiling. If the floor surface consists of synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst according to IEC 61000-4-4	± 1kV for input and output lines ± 2kV power cables	± 1kV for input and output lines ± 2kV power cables	The quality of the supply voltage should conform to the typical business or hospital environment.
Surge voltages according to IEC 61000-4-5	± 1kV push-pull voltage ± 2kV push-pull voltage	± 1kV push-pull voltage ± 2kV push-pull voltage	The quality of the supply voltage should conform to the typical business or hospital environment.
Voltage dips, short interruptions and variations of the power supply according to IEC 61000-4-11	<5% U_T for ½ period (>95% dip of U_T) 40% U_T for 5 periods (60% dip of U_T) 70% U_T for 25 periods (30% dip of U_T) <5% U_T for 5sec. (>95% dip of U_T)	<5% U_T for ½ period (>95% dip of U_T) 40% U_T for 5 periods (60% dip of U_T) 70% U_T for 25 periods (30% dip of U_T) <5% U_T for 5sec. (>95% dip of U_T)	The quality of the supply voltage should correspond to the typical business or hospital environment. If the user of the UNIT requires it to continue functioning following interruptions of the power supply, it is recommended to have the UNIT powered by an uninterruptible power supply or a battery.
Magnetic field of power frequencies (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	The power frequency magnetic fields should correspond to the typical values found in the relevant business and hospital environment.
Remarks: U_T is the AC supply voltage prior to application of the test level.			
Conducted HF interference IEC 61000-4-6	$3V_{eff}$ 150 kHz to 80 MHz ^a	$3V_{eff}$	Portable and mobile radio equipment must not be used within the recommended working clearance from the UNIT and its cables, which is calculated based on the equation suitable for the relevant transmission frequency. Recommended working clearance: $d = [1, 2] \sqrt{P}$

Immunity interference tests	IEC 60601-1-2 test level	Conformance level	Electromagnetic environment guidelines
Radiated HF interference IEC 61000-4-3	<p>3V/m 80MHz to 800MHz^a</p> <p>3V/m 800MHz to 2.5GHz^a</p>	<p>$3V_{\text{eff}}$</p> <p>$3V_{\text{eff}}$</p>	<p>$d = [1, 2]\sqrt{P}$ at 80MHz to 800MHz</p> <p>$d = [2, 3]\sqrt{P}$ at 800MHz to 2.5GHz</p> <p>where P is the nominal transmitter output in watts (W) specified by the transmitter manufacturer and d is the recommended working clearance in meters (m).</p> <p>The field strength of stationary radio transmitters is based on a local investigation for all frequencies^b less than the conformance level for all frequencies^c.</p> <p>Interference is possible in the vicinity of equipment bearing the following graphic symbol.</p> 

- a. The higher frequency range applies at 80MHz and 800MHz.
- b. The field strength of stationary transmitters such as the base stations of radio telephones and land mobile services, amateur radio stations as well as AM and FM radio and television broadcasting stations cannot be accurately predetermined. An investigation of the location is recommended to determine the electromagnetic environment resulting from stationary HF transmitters. If the field strength measured at the **UNIT** location exceeds the conformance level specified above, the **UNIT** must be observed with respect to its normal operation at each application site. If unusual performance characteristics are observed, it may be necessary to take additional measures such as reorientation or repositioning of the **UNIT**.
- c. A frequency range of 150kHz to 80MHz results in a field strength of less than 3V/m.

Working clearances

Recommended working clearances between portable and mobile HF communication devices and the UNIT

The **UNIT** is intended for operation in an electromagnetic environment, where radiated HF interference is checked. The customer or the user of the **UNIT** can help prevent electromagnetic interference by duly observing the minimum distances between portable and/or mobile HF communication devices (transmitters) and the **UNIT**. These values may vary according to the output power of the relevant communication device as specified above.

Nominal transmitter output [W]	Working clearance according to transmission frequency [m]		
	150kHz to 80MHz	80MHz to 800MHz	800MHz to 2.5GHz
	$d = [1, 2]\sqrt{P}$	$d = [1, 2]\sqrt{P}$	$d = [2, 3]\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters whose maximum nominal output is not specified in the above table, the recommended working clearance d in meters (m) can be determined using the equation in the corresponding column, where P is the maximum nominal output of the transmitter in watts (W) specified by the transmitter manufacturer.

Annotation 1

The higher frequency range applies at 80 MHz and 800 MHz.

Annotation 2

These guidelines may not be applicable in all cases. The propagation of electromagnetic waves is influenced by their absorption and reflection by buildings, objects and persons.

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Sirona Dental Systems GmbH
Sirona Dental Systems LLC