R605DASM

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LIVELLO REVISIONE/REVISION LEVEL:

| This manual is a translation of the Italian text, which prevails in case of doubts. | | | |
|--------------------------------------------------------------------------------------------|------------------|-----------|--------------------------------|
| Rev. N | Data/Date | N | Descrizione/Description |
| A | 06.09.2007 | MRE/14/07 | Board updating |
| | | | Aggiornamento scheda |
| B | 06.11.2007 | MRE/29/07 | Board updating |
| | | | Aggiornamento scheda |

CONTENTS

| | |
|------------------------------------------------------------------------------------------|-----------------|
| livello revisione/revision level: | ----- iv |
| <u>A - AVVISI / FRONT MATTER</u> | |
| <u>B - DESCRIPTION</u> | |
| Characteristics: | ----- 2 |
| Versions: | ----- 3 |
| R605/170/DASM- | ----- 3 |
| R605/025/DASM- | ----- 4 |
| R605/027//DASM | ----- 5 |
| <u>C - SPECIFICHE/SPECIFICATIONS</u> | |
| Validation of specification data: | ----- 2 |
| Classification | ----- 2 |
| Operation environment: | ----- 2 |
| <u>D - SIMBOLI/SYMBOLS</u> | |
| <u>E - COMPATIBILITY WITH X-RAY TUBES</u> | |
| <u>F - MOUNTING THE COLLIMATOR TO THE X-RAY TUBE</u> | |
| Mechanical installation of the external interface unit: | ----- -3 |
| <u>G - COLLEGAMENTO ELETTRICO/ELECTRIC CONNECTION</u> | |
| Collegamento/Connection R605 DASM: | ----- 1 |
| <u>H - ELECTRICAL POWER CONNECTION OF THE EXTERNAL UNIT TO THE GENERAL SYSTEM</u> | |
| ----- | ----- 3 |
| Collegamento elettrico - ASR003 | ----- 4 |
| <u>I - ELECTRONIC SYSTEM, DESCRIPTION:</u> | |
| Electric connections to the board: | ----- -2 |
| J1 - Programmazione/Programming- | ----- 2 |
| J2 - Rete CAN-BUS / CAN-BUS network | ----- 2 |
| J3 - Fotocellule/Photocells | ----- 2 |
| J4 - Display | ----- 3 |
| J5 - Motore / Motor | ----- 3 |
| J6 - Alimentazione / Supply | ----- 3 |
| ----- | ----- 4 |
| CAN BUS INTERFACE | ----- 5 |
| ID - Identifications | ----- 6 |
| ID 700h (model) | ----- 6 |
| ID 701h (state) | ----- 6 |
| ID 702h (controls) | ----- 6 |
| ID 703h (serial number) | ----- 7 |
| ID 704h (drive mode)- | ----- 7 |

| | |
|--------------------------------|-----|
| ID 705h (drive speed)- - - - - | 8 |
| ID 706h (Options) - - - - - | -10 |
| Procedure - - - - - | -12 |

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L - SYSTEM CALIBRATION

M - ADJUSTMENTS

| | |
|---------------------------------------------------|-----------|
| Iris Motor Adjustment: - - - - - | -1 |
| Motor gear - intermediate gear: - - - - - | 1 |
| Intermediate gear- iris movement gear: - - - - - | 1 |
| Adjustment of Iris Photocell - - - - - | -1 |
| Adjustment of Shutter Motors - - - - - | -1 |
| Shutter gear - intermediate gear:- - - - - | 1 |
| Motor gear: - - - - - | 1 |
| Adjustment of Shutter photocells - - - - - | -2 |

N - OPERATION INSTRUCTIONS

| | |
|-------------------------|---|
| System setting- - - - - | 1 |
|-------------------------|---|

O - SUBSTITUTIONS, DISASSEMBLY, TRANSPORT

| | |
|----------------------------------------------------|-----------|
| Substitution of the Iris Motor - - - - - | -1 |
| Substitution of the Iris Photocell- - - - - | -1 |

| | |
|---------------------------------------------------|----------|
| Substitution of shutter Motors | 1 |
| Substitution of Photocell Shutters. | 2 |
| Disassembly | 2 |
| Transport and storage: | 2 |
| P - ROUTINE MAINTENANCE | |
| Cleaning recommendations | 1 |
| Recommended maintenance programme: | 1 |
| Q - TROUBLESHOOTING | |
| R - PARTI DI RICAMBIO / SPARE PARTS | |
| Etichetta/Label | 1 |
| R605/DASM | 1 |
| Esploso/Parts Breakdown: | 2 |
| R605DASM : | 2 |
| R605/025/DASM | 3 |
| R605/027/DASM | 4 |
| S - GENERAL | |
| Repairs | 1 |
| End of life disposal | 1 |
| Warranty | 1 |
| Components not covered by this warranty: | 2 |
| Safety/Responsibility | 3 |
| T - FIGURE E SCHEMI / FIGURES AND DIAGRAMS | |

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
B - DESCRIPTION

Compact radiological automatic collimating system for round and elliptic fields designed to operate with a mobile "C" arm Image Intensifier equipment.

R605DASM

The round and elliptical fields are defined as follows: The round field by 8 lead shutters located near the exit window and a brass cone near the x-ray focus; the elliptical field by the round field and two pairs of lead rectangular shutters located near the collimator entrance window.

Round and elliptical field shutters are controlled by 5 stepping motors.

The collimator features a microprocessor unit built into the collimator to control the 5 stepping-motors via external signal source with CanBus protocol.

The circuits return a CanBus protocol signal to indicate correct motor positioning.

The two pairs of lead rectangular shutters move jointly and both rotate $\pm 360^\circ$.

During rotation the elliptic field remains steady , with a tolerance of ± 10 mm, at a Focus Film Distance (FFD) of 100 cm.

The field dimension at FFD 100cm (39.4") with the collimator mounted at 80 mm (3.14") from the focus is 30 cm (12").


Characteristics:

- Radiation shielding 125 - 4 mA.
 - Inherent filtration absent
 - Round field variable from a minimum of 5cm to a maximum of 30 cm at 1m FFD.
- The radiological collimator, with round and elliptic x-ray fields, is designed and constructed for installation with **rotating or fixed anode x-ray tubes (EN 60601-1-3 par. 29.202.3)**. The motorised movement of its shutters permits adjustment of the x-ray field to the area under investigation.

NOTE:

STANDARD VERSION FORESEES 30 CM (12") AT A FOCUS DISTANCE OF 100CM, FOR LOWER OPENING IT IS NECESSARY TO SUBSTITUTE THE SECONDARY COLLIMATOR. SPECIFY THE MAXIMUM REQUIRED OPENING ON THE ORDER.

R605DASM


descrizione_605DASM.fm

Versions:


System can be provided in simplified and/or customised version, below described

R605/170/DASM

without external interface board located between the collimator and the external analogical components.

collimator\605\605_170_DASM


collimator\605\605_170_DASM

descrizione_605DASM.fm

R605/025/DASM

without external interface board located between the collimator and the external analogical components and without second pair of rectangular shutters.


descrizione_605DASM.fm

R605/027//DASM

without external interface board located between the collimator and the external analogical components and without second pair of rectangular shutters and with electronic circuit on the short side of the collimator .

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| | |
|-------------------------------------------|-----------------------------------|
| Collimator operation: | Mobile fluoroscopy systems |
| Dotazione Standard: | - Manuale Istruzioni |
| Items included with the collimator | - Instructions Manual |
| Peso/ Weight | 1,6 Kg |

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NB: the secondary collimator needs to be substitute if smaller fields, 9" and 7", are required. Please specify this requirement in your order.

Validation of specification data:

- Validation of minimum filtration of the radiation unit (x-ray tube, collimator and possible filters) must be performed on a completely installed system by whoever is responsible for the installation.
- Validation of light field luminosity and contrast is performed in-factory. See Chapter **COMPLIANCE VERIFICATION**.
- Validation of system x-ray leakage is to be performed following the installation of the system components under the responsibility of the person performing the installation.
- The correct installation of the system, electronic calibration and wiring layout of the same are the responsibility of the engineer responsible for the installation of the system.
- Protection of the cable is to be provided by the installation engineer.

Ralco is prepared to provide any information required regarding the validation methods described above.

Classification**WARNING**

GUIDES FOR ACCESSORIES PRIOR TO INSERTING ACCESSORIES IN THE GUIDES CHECK ON THE PERFECT FIT OF THE SUPPORT WITH THE MOUNTING SLOT ON THE COLLIMATOR (TOLERANCE MAX. $\pm 0,5$ MM). A FAULTY FIT COULD BE DANGEROUS AND IT COULD CAUSE THE FALL OF ACCESSORIES

EN 60601-1 par. 5

- Protection against electric hazards: "Class I" equipment
- Protection against direct and indirect contacts: Type B equipment with applied parts.
- Protection against water see page: "Common equipment"
- Safety of operation in the presence of inflammable anaesthetics with air or oxygen or nitrous oxide: Equipment not suited to application in the presence of inflammable anaesthetic mixtures containing air or oxygen or nitrous oxide.
- Operation conditions: Equipment for continuous operation at intermittent loads - See Chapter OPERATION INSTRUCTIONS.
- Should label data on the collimator not correspond to the specifications herein, inform Ralco of the non conformity.
- Verifications of the specifications are to be performed according to the indicated equipment standards.

Operation environment:

- Ambient temperature = from 10°C to 40°C
- Relative Humidity = from 30% to 75%
- Atm. Press. = from 700 to 1060 hPa.

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EN 60601-1-3 para. 29.205 - Focus-Skin Distance

Paragraphs 29.205.1 - 29.205.3 of the above standard deal with preventing the use of inappropriate focus/skin distances in order that the dose equivalent to the patient be kept as low as reasonably achievable.

29.205.1 - Radiation Equipment for fluoroscopy




X-ray devices for fluoroscopy must feature characteristics that, during fluoroscopy, will prevent focus-skin distances less than:

- 20 cm if the x-ray device is specifically for fluoroscopy
- 30 cm for other applications

Conformity is measured and checked visually. The x-ray system must comply with the norm; this means that the operator is responsible for the verification. A focus/skin spacer device is an option made available by Ralco.

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

E - COMPATIBILITY WITH X-RAY TUBES:




- Compatibility between the collimator and x-ray tube is determined by the mechanical possibility of mounting them - See the following **Table**.
- Collimator R605DASM has no inherent filtration. Check that the tube housing literature indicates congruent minimum inherent filtration (1mm or 0mm) and, that maximum radiation leakage is 30 mr/hour measured at one meter from the source when operating at its leakage technique factors (125 at 4 mA).
- Source values (tube housing-collimator) must not be less than 3mm Al for filtration and must never exceed 100 mR/hr for radiation leakage. (EN 60601-1-3 par. 29.201.7/29.201.6 e 21 CFR sub-chapter J, part 1020.30 (m) (1).
- The distance between the tube x-ray focus and the flange mounting plane (collimator upper plate) must be 80 mm (3.14"), tolerance +/- 1mm (0.04").

Table 1:

| Collimators | |
|----------------------|---------------------------------------------------------------------------------------------------------------------|
| R72 |  <p>general\vr72_cuffie.gif</p> |
| R103 R104 R105 |  <p>general\cuffie_104.gif</p> |

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



F - MOUNTING THE COLLIMATOR TO THE X-RAY TUBE

WARNING:

CAREFULLY FOLLOW THE MOUNTING INSTRUCTIONS AND MAKE SURE THAT THE COLLIMATOR IS CORRECTLY ASSEMBLED. INCORRECT MOUNTING COULD BE DANGEROUS: IT COULD CAUSE THE COLLIMATOR TO FALL OR TO OPERATE INACCURATELY

- a) Determine the distance from the focal spot to the tube port face from the X-ray tube housing literature
- b) Subtract the resulting distance from 80 mm. (3.15") and determine how many 1.5 mm (0,06") spacers combined with the thickness of the mounting flange will make up the difference (15mm for the flange). Allowable tolerance is 1 mm. (0.04")
- c) Carefully match the flange with the tube and make sure that the brass cone on the mounting flange is free from impediments on the side of the housing.
- d) Select four bolts of suitable thread (M4) and of such a length that they protrude through the flange and spacers far enough to engage at least 5 threads into the tube port face. Securely bolt the flange to the tube port face..


IMPORTANT

TO SAFEGUARD THE OPERATOR AND PATIENT AGAINST THE HAZARD OF A FALLING COLLIMATOR, THE FOLLOWING INDICATIONS ARE TO BE RESPECTED.

- e) Unscrew, but not remove, the three mounting Allen screws on the collimator metal ring. See **Figure 1** in this Chapter
- f) Fit the metal ring on the collimator over the mounting flange and tighten the screws.

Figura 1

collimator\ovs\montaggio_tubo_fig 1.gif

montag_coll-tubo.fm

MECHANICAL INSTALLATION OF THE EXTERNAL INTERFACE UNIT:


The external box is to be installed in such a way as to ensure connection of all the electric devices present in the radiological system.

We suggest that the box be installed at a safe distance from electromagnetic sources such as high voltage transformers.

Use the holes on the sides of the box to perform the installation.

Should it be necessary to install the back within a rack of the general system, remove the electronic board from the box prior to proceeding with the installation of the box itself.

Make sure that each portion of the rack that is to accommodate the box is adequately earthed.

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
WARNING

COLLIMATOR SUPPLY IS NOT PROTECTED BY A FUSE. PRIOR TO CONNECTING THE UNIT CHECK THAT SUPPLY IS PROTECTED BY AN EXTERNAL FUSE.
FUSE: 3,15 AMP

Collegamento/Connection R605 DASM:


collimatori\605\connett

Table 1: 9 pin (male) D-Sub CAN Bus PinOut

| Pin n. | Signal Names | Signal Description |
|--------|----------------|--------------------|
| 1 | | |
| 2 | CAN_L | Dominant Low |
| 3 | CAN_GND | Ground |
| 4 | | |
| 5 | | |
| 6 | GND (24 AC) | Ground |
| 7 | CAN_H | Dominant High |
| 8 | | |
| 9 | CAN_V+ (24 AC) | Power |

Table 2: R605/027/DASM

| Pin n. | Signal Names | Signal Description |
|--------|--------------|--------------------|
| 1 | 24 VAC | Supply |
| 2 | 24 VAC | Supply |
| 1 | CAN_L | Dominant Low |
| 2 | CAN_GND | Ground |
| 3 | CAN_H | Dominant High |

.images/colliamtori/605/col_ell_fig4

col_ell_fm

H - ELECTRICAL POWER CONNECTION OF THE EXTERNAL UNIT TO THE GENERAL SYSTEM

WARNING: THE SYSTEM MUST BE SUPPLIED BY 24VAC +/- 10%. ANY VOLTAGE IN EXCESS WILL DAMAGE STEPPER MOTOR CONTROL ELECTRONICS

:

WARNING

PRIOR TO CONNECTING THE COLLIMATOR CHECK THAT SUPPLY IS PROTECTED BY AN EXTERNAL FUSE 6,3A (SEE THE SPECIFICATIONS).


CABLES (24VAC) USED FOR COLLIMATOR CONNECTION MUST BE SUITABLE FOR PURPOSE AND COLLIMATOR CURRENT ABSORPTION.(MIN. SECTION 1.5MM²)

THE SYSTEM REQUIRES SCREENED CABLES WITH THE ADDITION OF FERRITE CORES TO AVOID INTERFERENCE OF ELECTROMAGNETIC DISTURBANCES WITH THE SERVER ELECTRONICS.

CABLES FROM THE SYSTEM TO THE EXTERNAL BOARD MUST BE SCREENED.

THE FERRITE CORES ARE TO BE MOUNTED ON THE CABLES INSIDE THE BOX AS SHOWN IN THE DRAWING BELOW. FERRITE CORES ARE USED TO ATTENUATE INDUCTION AND IRRADIATION OF ELECTROMAGNETIC INTERFERENCE.

WE SUGGEST A 279 FERRITE, ATTENUATION RANGE 100 MHZ TO 200 MHZ, SUITED TO YOUR CABLES; EG: MSFC-10 FERRITE BY RICHCO

elett_302_ACSdhhs.fm

Supply and signals to the collimator must be to 93/42CEE standards. Devices that supply to the collimator must therefore feature **double or reinforced insulation** as provided by the General Standard on Electro medical Equipment CEI62-5 +A2 (EN 60601-1)

The device is electrostatic sensitive consequently all the relating safety standards must be complied with.

All the system devices listed in the following tables (variable FFD table, variable FFD stand, fixed lateral FFD, vertical Bucky, lateral Bucky, etc.) are to be connected to the connectors of the external board.

WARNING: PREVENT METAL RESIDUE (SMALL CLIPPINGS OF CABLE WIRE) FROM DEPOSITING ON THE CIRCUIT TRACKS - THESE COULD DAMAGE THE BOARD OR CAUSE MALFUNCTIONS. PRIOR TO SUPPLYING THE SYSTEM, USE A VACUUM CLEANER TO CLEAN THE BOARD AND BOX. **Do NOT USE COMPRESSED AIR.**

Correct operation of the system is guaranteed if the return signals from the Bucky (if is foreseen) and FFD Potentiometres are not disturbed. Use tripolar screened cables with the addition of ferrite cores near the connections to the PC board (external unit)

The signal from the FFD variation control device must be analogic DC, 5VDC max

N.B. All transducers, Potter Bucky (if is foreseen) and FFD Potentiometres, fixed distance micro switches (if are foreseen), push buttons, etc. are analogic signals and are supplied by the board. The 5VDC power that supplies the Potentiometres is also used to supply the microprocessors that operate the collimator. The CAN signals are at 5V and would be damaged if they were to be connected to other circuits. It is therefore highly recommended to avoid reading electric signals shared with other electronic equipment; i.e. **do not use analogic signals generated by a unit external to the system** and make sure that the analogic signals are connected via screened cables.

Analogic signals from the main system need to be used (Volt direct current), these signals must never exceed 5V max. Any value in excess of this could damage the board.

Use appropriate cables to connect supply to the external unit (24V AC 5A - minimum supply transformer power 150 VA - external fuse 6A3) on contacts 83 and 84 (see the layout of board ASR003) and earth on contact 85. The metal box, besides being a mechanical protection, serves as electromagnetic shielding.

IMPORTANT:


THE 24VAC TO THE BOARD MUST BE SUPPLIED BY AN EXTERNAL TRANSFORMER WITH DOUBLE INSULATION ON THE SECONDARY WINDING. SECONDARY OUTPUTS MUST NOT BE EARTHED.

The earth connection of the collimator circuit is used to disperse noise or extra voltage coming from the collimator supply system.

The device that disperses noise is reversible, consequently noise instead of being eliminated is fed into the system power supply circuit. It is therefore essential that the earth outlet to which the system is connect be effectively the unipotential protection point of the system or the room.

NB: LEFT AND RIGHT SIDES HAVE BEEN DETERMINED BY
LOOKING AT THE COLLIMATOR FRONT PANEL.

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
I - ELECTRONIC SYSTEM, DESCRIPTION:

Electronic Board CAN-DC-STEP (MP11719) controls a single stepping-motor through CAN Bus interface. Shutter movements of collimator R605DASM are controlled by three motors; this means that the collimator includes three boards MP11719.

Motor positioning is controlled via an open loop; i.e. impulse count sent to the motor, a direct function of the motor shaft movement.

Stepping-motor operation provides a quasi-trapezoidal speed/time profile, characterized by a starting speed (V_{min}), acceleration time (T_R), running speed (V_{max}) and by deceleration time (T_R).

In the case of Dist travel lasting less than twice the ramp time T_R , running speed V_{max} would not be attained. In this case instant speed would reach its maximum value at mid-run and then reduce again to V_{min} speed at which point travel will stop.

Maximum current is set both by a trimmer situated on board and by CAN BUS setting over four levels between 50% and 100% of the value provided by the trimmer. Voltage in motor windings is controlled also in its damping stage and is adjustable over four levels.

ELECTRIC CONNECTIONS TO THE BOARD:**J1 - Programmazione/Programming**

| PIN | Descrizione | Description |
|------------|--------------------|--------------------|
| 1 | Alimentazione +5V | Supply + 5V |
| 2 | Busy | Busy |
| 3 | SCLK | SCLK |
| 4 | RX livello TTL | RX level TTL |
| 5 | CE | CE |
| 6 | n.c. | n.a. |
| 7 | Alimentazione GND | Supply GND |
| 8 | Reset | Reset |
| 9 | CNVSS | CNVSS |

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
images\elettrico\R605_27_CanDCStep.gif

CAN BUS INTERFACE

Board interfacing makes use of six CAN BUS identifications (stepping-motor)..

| ID | Tipo | byte | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | note |
|------|----------|------|-----|-----|-----|---------|------|------|------|------|--------------------|
| 700h | Remote | 8 | 'C | "S | " | "V | "1 | ". | "0 | "3 | 'Modello / release |
| 701h | Data out | 5 | P1 | P2 | P3 | P4 | ST 1 | --- | --- | --- | state |
| 702h | Data in | 1 | 00h | --- | --- | --- | --- | --- | --- | --- | Home |
| 702h | Data in | 5 | 01h | D1 | D2 | D3 | D4 | --- | --- | --- | Movim. assoluto |
| 702h | Data in | 5 | 02h | D1 | D2 | D3 | D4 | --- | --- | --- | Movim. relativo |
| 702h | Data in | 1 | FFh | --- | --- | --- | --- | --- | --- | --- | STOP |
| 703h | Remote | 8 | 'S | "/ | "n | 'X X | XX | XX | XX | XX | N. serie |
| 704h | Data in | 5 | M | K | T | D | L | --- | --- | --- | Drive mode |
| 705h | Data in | 7 | M | L1 | L2 | H1 | H2 | R1 | R2 | --- | Drive speed |
| 706h | Data in | 1 | O | --- | --- | --- | --- | --- | --- | --- | Opzioni |
| 707h | Data in | 8 | C1 | C2 | C3 | C4 | C5 | C6 | C7 | C8 | Display 1a riga |
| 708h | Data in | 8 | C9 | C10 | C11 | C1 2 | C1 3 | C1 4 | C1 5 | C1 6 | Display 1a riga |
| 709h | Data in | 8 | C17 | C18 | C19 | C2 0 | C2 1 | C2 2 | C2 3 | C2 4 | Display 2a riga |
| 70Ah | Data in | 8 | C25 | C26 | C27 | C2 8 | C2 9 | C3 0 | C3 1 | C3 2 | Display 2a riga |

P1,P2,P3,P4 = position at 32 bit P1 = bits 31..24

ST1 = status byte 1

D1,D2,D3,D4 = movement at 32 bit D1 = bits 31..24 absolute and relative

M = Mask, every bit determines a valid field

K = fraction of step

T = torque

D = decay

L = logical level per home

L1,L2= start frequency (Vmin) L1 = bits 15..8

H1,H2 = running frequency (Vmax) H1 = bits 15..8

R1,R2 = ramp time Tr R1 = bits 15..8

O = options

C1 ..C16 = characters on the first line of the display

C17.C32 = characters on the second line of the display

ID - Identifications

IDs and some other parameters may be varied by an appropriate programme and serial interface and/or using other systems.

ID 700h (model)

This ID, sent as remote request makes the board return its p/n, eg.: " CS V1.03".

ID 701h (state)

This frame is sent by the board independently every 250 mS. In a binary form it contains: current position determined on the base of counts and a status Byte coded as 5° byte

Transmission may also be forced by using the Force Update option.

5° byte

| | | | | | | | |
|------|---------|---------|---------|---------|--------------|-----------|--------|
| 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| Home | Riserv. | Riserv. | Riserv. | Riserv. | Pro- tez. | 0 degrees | operat |

- bit 0 is active during motor movement.
- bit 1 is active in correspondence to the controller 0 angle.
- bit 2 is active on start-up prior to the first stop control; it is activated in case of over current or excessive heat of the device.
- In these cases the board must be switch **off** and **on** again.
- bit 7 is active is input zero signal is high (e' attivo se il segnale di zero in ingresso e' a livello alto (as if the fork photocell were obscured).

ID 702h (controls)

This frame is used to send movement control to the board. In relation to, the contents of the first byte we have four different controls.

00h home

Performs a re-setting cycle of a given movement according to minimum speed. During this operation the system ceases to communicate.

01h absolute positioning

The subsequent 4 bytes represent the absolute position that needs to be reached in field 0 /1000000 i.e. 00 00 00 00 /00 0f 42 40. Any number exceeding this value is considered equal to the maximum value (1000000) of the domain.

The movement occurs at a quasi-trapezoidal speed profile; during the entire movement the status Byte of frame 701h signals that the system is moving and the position-related field is progressively updated.

02h relative positioning

The subsequent 4 bytes represent the relative movement that it has to be performed in field -100000/100000 i.e. FF F0 BD C0/00 0F 42 40. Any number that is entered is added/subtracted to the current position and verified that it falls in domain 0-1000000.

The movement occurs at a quasi-trapezoidal speed profile; during the entire movement the status Byte of frame 701h signals that the system is moving and the position-related field is progressively updated.

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FFh movement stop

This control stops the movement whatever the current position.

ID 703h (serial number)

This ID, sent as remote request makes the board respond with its S/n eg.: S/n_____.

ID 704h (drive mode)

This control sets some parameters related to the electromechanical system in question.

1°Byte (mask)

the scope of this byte is to select, wherever necessary, the parameters that require modification within the string. If the respective bit is set at 1, the corresponding value is transferred. It will, otherwise be ignored by the system.

The mask is necessary since a single string contains a number of parameters and only a third or whatever combination might want to be transferred as valid. In this case, parameters corresponding to non-selected positions (zero mask bits) will be ignored and may then assume any value.

1° Byte.

| | | | | | | | |
|---------|---------|---------|---------|---------|-------|--------|-------|
| 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| Riserv. | Riserv. | Riserv. | Riserv. | Home P. | Damp. | Torque | Frac. |

Example: if set at 0fh, all the ensuing fields are considered valid; if set at 02h the torque-related value only will be extracted from the string.

2°Byte (Step fraction)

The scope of this byte is to set the fraction required which may be chosen from:

00 =used only to cut-off supply to the motor. It makes sense only in the case of an complete step otherwise, actual position within the fraction itself would be lost.

The reason is that a motor stopped over a step fraction must be energised.

01 =invalid NEVER USE

02 =invalid NEVER USE

03 =entire step

04 =half step at low torque

05 =half step at high torque

06 =quarter of a step

07 =eighth of a step

08 =sixteenth of a step

NOTE:

THE EXCESSIVE MOTOR AND BOARD HEATING CAN CAUSE THE COLLIMATOR SHUTTERS TO BLOCK. THE OVERHEATING IS CAUSED BY A CONTINUOUS 100% SUPPLY OF HIGH VOLTAGE. IT IS ADVISED, WHEN THE MOTOR IS STOPPED, TO LOWER THE HIGH VOLTAGE RISE AS SHOWN BELOW.

3°Byte (Torque)

The scope of this byte is to set the maximum torque fraction required:

00 = 50% of the maximum value (160mA - 225mA)

01 = 70% of the maximum value (225mA - 315mA)

02 = 85% of the maximum value (270mA - 380mA)

03 = 100% of the maximum value (320mA - 450mA)

These values can be set by string sent to board even when it is working.

00 = 07x4 5 02 670067 67

01 = 07x4 5 02 670167 67

02 = 07x4 5 02 670267 67

03 = 07x4 5 02 670367 67

4°Byte (voltage decay)

The scope of this byte is to set phase voltage decay required:

00 = 12.5% low decay coefficient

01 = 37.5% medium decay coefficient

02 = 75% high decay coefficient

03 = 100% rapid decay coefficient

These values can be set by string sent to board even when it is working

00 = 07x4 5 04 676700 67

00 = 07x4 5 04 676701 67

02 = 07x4 5 04 676702 67

03 = 07x4 5 04 676703 67

Where x is a base value.

Base ID 720 724 x=2

5° Byte (Logical level of the home sensor)

The scope of this byte is to set the logical level for the input of the home sensor when the system is in the home position.

In the case of connection to a fork photocell, this value must be set at 00 for the photocell that is obscured at zero position (occupied slot) or set at 01 for the photocell that receives light at a zero position (free slot).

ID 705h (drive speed)

This controls sets some parameters related to the speed of the electromechanical system in question.

1° Byte (mask)

The scope of this byte is to select, wherever necessary, the parameters that need to be modified within the string. If the respective bit is set at 1, the corresponding value is transferred, otherwise it will be ignored by the system.

$$\frac{60 * \text{Freq}}{K * P}$$

The mask is necessary since a single string contains a number of parameters and only a third or whatever combination might want to be transferred as valid. In this case, parameters corresponding to non-selected positions (zero mask bits) will be ignored and may then assume any value

1° byte.

| | | | | | | | |
|---------|---------|---------|---------|---------|-------|-------|-------|
| 2 | 2 | 2 | 2 | 2 | 2 | 2 | 2 |
| Riserv. | Riserv. | Riserv. | Riserv. | Riserv. | Rampa | F.max | F.min |

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The scope of this byte is to set the minimum frequency corresponding to the motor start-up speed. Set value is in Hertz, minimum value is 125 Hz.

The corresponding angular speed is a function of the following equation:

R.p.m. =

where K is the step fraction and P the number of steps per revolution.

Example:

Motor 200 step/turn, set at a quarter step.

Start-up frequency set at 500 Hertz

Corresponding turns = 37.5 r.p.m.

since the stepping-stepping motor features a torque that reduces with the speed the a correct selection of this value makes the best of the motor characteristics.

4° e 5° Byte (Frequenza massima)

The scope of this byte is to set the maximum frequency corresponding to motor running speed. The set value is in Hertz, maximum value is 20000 Hertz.

The same considerations made for the minimum frequency apply. For the same reason, maximum frequency must not exceed the speed value at which the motor maintains a useful torque; a stall fault would otherwise occur resulting in loss of position.

6° e 7° Byte (Ramp Time)

The scope of this byte is to set motor acceleration and deceleration times in milliseconds which is the time value necessary to pass from minimum to maximum speed.

A too short a value could cause a stall fault; on the other hand an excessive value could enhance resonance phenomena, both electrical and mechanical, which could also cause stall faults.

A correct value allows a load/motor to accelerate and the system inertia allows speeds at which the system resonates without entering a stall condition.

ID 706h (Options)

This control serves to set some board operation options. These are PERMANENT options which means that they will still be in place even without power supply; they are normally used in the board installation phase.

It will be possible to change these options at least 100 times.

1°Byte (Option 1)

There are five fields:

1st field: LCD, refers to the possible use of a liquid-crystal display connected to the appropriate port.

Set at 1 the appropriate display is connected - set at 0 the display is not available

2ndfield: refers to the board-mounted LED operation modes:

- 00: Status display of Home sensor.
- 01: flashing at 1 Hz indicating board operation.
- 10: motor status; on = in movement , off = stopped
- 11: receipt of a valid CAN frame.

3rd field: possibility of moving upwards the addresses of the board output frames by a value equal to 0x80 when the field is set at 1.

4th field: consents the use of the second display line by suppressing the automatic value display. When set at 0 the value is displayed, when set at 1 the area may be used to send text via CAN BUS.

5th field:

Set at 1, this field forces transmission of ID 701h even in the presences of a movement control and movement actuation.

The above details refer to standard firmware installed on the board for the STEPPING-MOTOR VERSION.

Controls:

After having connected the external devices to the board and, set up stepping-motor voltage by means of trimmer R8, the board needs to be tested for CAN BUS net interrogation response.

The standard versions frame 701h, composed of 5 bytes, is sent automatically every 250 msec.

Receipt of this frame indicates correct CAN BUS net connection.

We start by sending **ID=700h DLC 8**; the board should reply with a string corresponding to the model.

We can then transmit control ID=703h DLC 8 to obtain the board serial number whenever necessary.

Frame 704h DLC = 5 is transmitted with the appropriate values bearing in mind that, if we want the board to consider the 4 parameters, the mask (1st byte) to be sent will be 0Fh.

This frame is not normally sent again

Frame 705h DLC =7 is then sent with values related to the two speeds (frequency) and acceleration time.

In this case also, the mask byte (07h) must be supplied.

Supposing we want a start frequency of 200 Hz and a running frequency of 1000Hz, acceleration time between the two frequencies of 300 mSec, we would send the following string: (besides ID e DLC) 07h 00h c8h 03h E8h 01h 2Ch.

We might want to perform a movement at a different running speed for some particular reason, eg.: running= 1500Hz. The string would be 02h 00h 00h 05h DCh 00h 00h.

The underlined digits can contain any value since, because of the mask 1st byte, they are not considered.

The same consideration applies to the variation of one or more parameters.

The correct procedure is now to send a re-set control to the motorised axis. To do this we send the following frame:

702h DLC = 1 with the 1° byte (type of movement) = 00h

Wait for the re-setting cycle to terminate

We can now, for example, send an absolute position control to get the actuator to move to position 5000. Following is the corresponding frame:

702h DLC = 5
01h 00h 00h 13h 88h

Frame 701h will now provide immediate information on actuator position (first 4 bytes) and the field relating to the 5th byte (movement status related) will remain at 1 until position is reached.

Should we now wish to move forwards by 400 points (relative movement):

702h DLC = 5
02h 00h 00h 01h 90h

In this case also we will monitor the status frame (701h) which will indicate movement status and ABSOLUTE axis position.

Let us now move by -1000 points; i.e. backwards with respect to the current position:



702h DLC = 5
02h FFh FFh FCh 18h

We might need to stop the movement so we'll send:

702h DLC = 1
FFh

The resulting position must be determined by the reading of frame 701h.

It is not normally necessary to perform a HOME procedure for the movements. It is sufficient to send an absolute movement control to position 0: this is faster than the re-setting process.

Procedure

With the board powered down, connect the serial adaptor to the connector

- Power up
- Use the increase and decrease buttons to read the desired value by 10h steps on the upward/downward scale.
- Remove the serial adaptor on connector J1.
- Switch the board off and remove the adjustment wiring.

The board is now ready to react to a new set of addresses.

A serial adapter may optionally be installed on connector J1 for a set of customised controls:

The board is must powered down: connect the serial adaptor on connector J1 and then to an ASCII terminal or PC featuring and appropriate terminal emulation programme) setting the following:

- Baud rate 9600
- 1 stop bit
- No parity
- Power up

Some parameters related to the board may now be read/changed as follows:

| Control | Example | Description | Min | Max | Notes |
|-------------|---------|-------------------------------|-------|-------|---------------------------------|
| Base_id | 0x200 | Basic address, board | 0x100 | 0x7F0 | |
| Update_rate | 200 | interval status transmission | 0 | 255 | if < 100 it is suppressed |
| Lcd | 1 | LCD in | 0 | 1 | See options, above |
| Led | 2 | LED operation | 0 | 3 | See options, above |
| Serial | Nessuno | Board S/N | 0 | 99999 | reserved |
| Can_prop | 1 | CAN propagation segment | 1 | 14 | (*) |
| Can_ph1 | 5 | CAN phase 1 segment | 1 | 14 | (*) |
| Can_ph2 | 5 | CAN phase 2 segment | 1 | 14 | (*) |
| Can_swj | 1 | CAN sync jump | 1 | 14 | (*) |
| Can_speed | 0 | CAN bit rate prescaler | 0 | 3 | (*) (**) |
| Offset | 1 | CAN data out offset | 0 | 1 | See options, above |
| Row_2 | 1 | Availability 2nd Display line | 0 | 1 | See options, above |
| Help / ? | - | Available control indications | - | - | Also indicates firmware version |

(*) random variation of these parameters could make the CAN-BUS inaccessible.

(**) legal values are:

0 = 500 Kbit/sec

1 = 250 Kbit/sec

2= NOT VALID

3 = 125 Kbit/sec

The control may be given in two ways: read only or setting.

An example of the read only is:

Base_id<cr>

The value of of this parameter will be viewed.

Or:

Base_id = 0x540

The value of this parameter will be varied. An inexistent or wrong control will return an Invalid command response. Wrong parameters could return an Invalid value response.

Warning!

These parameters are permanently memorized in the microprocessor FLASH memory and are retained even on powering down.

The FLASH memory has a finite number of operations consequently these settings are advisable only when strictly necessary. Values set beyond the limits indicated will not be memorised.

When setting is terminated, power down the board, remove the connections. On start up, the board will operate with the new settings.

Warning!

The optional bit related to the LCD display may be set at zero only also in versions designed for LCD display operation. Should a fork photocell type OPB840W51 need to be connected (J3), proceed as follows:

Pin 1 not connected

Pin 2 not connected

Pin 3 black

Pin 4 red

Pin 5 white

Pin 6 green

Resetting to 1 may be performed only in the factory.

J - ASR003 - EXTERNAL BOARD DESCRIPTION

Characteristics:

- Manual interface
- Potter Bucky Potentiometres interface.
- Vertical Focus-Film distance interface (FFD/SID) with potentiometers related to table and stand.
- Interface with 5 fixed lateral distances.
- Interface for general controls; e.g., open/close for shutters, tomography. See the specifications regarding this point described in this chapter.


Description:

Board ASR003 is an extension of the radiological system. It provides interfacing with the collimator which features board CAN DC STEP and analogic signals coming from other system devices.

It operates via CAN BUS with a Ralco proprietary protocol.

The user end features 32 digital inputs, 8 analogic inputs and 10 relay outputs.

So as to adapt to different devices, digital inputs are fixed for optoisolation so as to allow connection also with devices that supply output control voltages (typical range 12-32 Vdc).

Ingresso digitale/digital Input


collimator\acs\asr003_1.gif

Since the input signal sometimes comes from a voltage-free contact, the board features a number of pins to which an external supply may be connected. Alternatively these pin may be connect to a voltage incoming from board ASR003 itself by simply connecting pin M79 to pin M80 and pin M81 to pin M82.

The available voltage is not taken upstream of the voltage regulator but downstream of the fuse and rectifiers bridge.

The type of software applied determines the function of each input.

The analogic inputs are designed mainly for connection to 1Kohm potentiometers.

collimator\acs\asr003_2.gif

Ingresso analogico/analogic Input

In this case the input range will be reduced with respect to the digitalization capacity of the digital/ analogic convertor on the microprocessor. However, at the same time the board will be better protected against wrong Potentiometres connections.

Other protection elements have also been added since, here too, it is possible that the user does not employ Potentiometres but supplies an analogical signal directly.

The main purpose of these inputs is to furnish the collimator with the cassette sizes that may be present or, to supply information related to focus-film/focus- floor distances (depending on the configuration).

In general, non linearity introduced by the protection elements is contained and adequate to their functions. However, should a greater degree of precision be required (such as in the case of focus-film/focus-floor signals) we suggest that a voltage proportional to the case be supplied which is included between 0 and 6 Vdc so as to make the most of the convertor range.

Higher voltages that could occur from incorrect connections, up to 12V, will not damage the input.


the 10 output relays each have a voltage-free exchange contact and may consequently be adapted to a number of different configurations. Their use depends on the software applied.

The CanBus communication interface is opto-insulated and needs external supply with a direct voltage between 9 and 32Vdc.

In this case also, voltage for interface supply may be drawn from any one of the radiological system supply sources or from an external source.

This type of connection is recommended when there are problems related to the presence of differences in earth potential within the system. See the typical connections diagrams that follow in this section.

This board also supplies the go-ahead to programme some of the radiological system parameters. In this case a jumper is required on the board as shown in the following photograph.

collimator\acs\asr003_3.jpg

The white line and black arrow indicate where to position the jumper on board ASR003; this requires the connection of pins 2 & 4 of connector J3.

No other configuration must be employed otherwise permanent damage might be caused to the board itself.


The jumper must be removed according to indications provided by the collimator display.

J3 is a utility connector used for in-factory programming of board ASR003.

Some software versions could also cover the use of a serial RS232 port located on connector J2 for a possible OFF-LINE parameterization of the board itself.

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


collimatori\acs\elet_6_asr003.gif

K - ASR003 BOARD

ASR003 board is the central control unit for collimator R 605DASM.

Configuration

This board features ample configuration possibilities. To access configuration, connect to serial gate J2 via CAN communication protocol.

FFD CONFIGURATION

Stepless single: FFD value is obtained locally by readings of board analogic inputs. It may also be received in millimetres directly via CAN-Bus. The collimator continuously tracks FFD value for reproduces a format of constant dimensions on the receiver/plate.

FFD Measuring Source:

The value of a stepless FFD may be determined locally, readings of ASR003 inputs, or received via CAN-Bus in digital format from external equipment.

Focus Distance:

All the calculations performed by the collimator require an 80mm F-D (Focus Distance). If the mounting distance of the collimator is not 80mm, the correct value may be specified. The parameter may also be used to correct small defects in the measurement of projected formats.

Hystereris on FFD change

This is a small millimetric variation to the FFD which triggers format correction when it is detected by the analogic inputs of ASR003.

FFD Observation time

Through the analogic inputs, measured FFD is tracked until each variation becomes smaller than the hysteresis and the time limit has expired. The collimator starts to re-track FFD variation only when the variation exceeds the hysteresis window.

FORMAT CONFIGURATION

Format programmes:

8 format programmes are available. Each programme includes limit values for the open of Iris, Long, Cross filtering shutters and a rotation angle for the filtering shutters. Once a format is selected, the collimator positions the iris, long and cross and filtering shutters as indicated by the programme: It will now be possible to close the field but not open it further.

Format Selection Mode:

The format may be selected from 8 pre-set values but may also be modified. In the latter case the values are received via CAN-Bus (in millimetres or sexagesimal degrees) for Iris, Long, Cross.

Format Source

Allocation of one of the 8 pre-set formats may be local by activating one of the board digital inputs or remote via CAN-Bus.

Allocation of individual format measurements is remote only.

CONFIGURATION OF MANUAL MOVMENTS CONTROLS**Movments controls:**

A good number of the board digital inputs are dedicated to the control of manual movements: it is possible to open/close the projected field, rotate shutters, perform standard position of filters ect..ect.

Alternatively, these controls may be received by a remote unit via CAN-Bus.

FFD Mode

FFD value is obtained locally from readings of analogic inputs. It may also be received in millimetres via CAN-Bus. The collimator continuously tracks FFD value to reproduce a format of constant dimensions on the receiver/plate.

FFD Variation from 78 cm to 120 cm

OPERATION MODES**Axes Reset modes:**

When the boards are powered-up, the system automatically selects, through a memory, the mode of the previous exam and starts the axes resetting procedure, if there is no error (timeout zero set or comunication error). All axes linked to fixed format are positioned on the first format (Iris movement, long, cross, F1,RF2, F2, RF2)

Calibration Mode:

This mode serves to calibrate potentiometers for FFD measurement (whenever this configuration is applicable).

Insert jumper JP at any moment to access this mode.

Calibration proper starts when the configuration is in the following status:

- FFD acquisition in stepless mode (single or differential)
- FFD measurement acquisition by potentiometer readings through analogic inputs of ASR003.

Manual Mode:

This is the primary operation mode of the collimator.

Further to the configuration given to the system, the collimator determines set format projection (with correct millimetric and rotation values) at a given FFD.

The format is selected from 8 preset dimensions or is received via ModBus or CAN-Bus.

Manual movements are controlled and limited to the selected format. It is possible to perform adjustments on open/close and rotation by 1mm and 1° at a time

Note: The movements are at the operator's discretion and are not influenced by the FFD.

Error Mode:

An Error occurs upon system operation anomalies. The following tables serves to troubleshoot the causes:

| Error Type | Problem description | Solution |
|----------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Memory Access Error | The permanent memory has been damaged or its contents have been corrupted. Very rare occurrence. | Ask for technical assistance. |
| CAN DC STEP fails to respond to the reset control. | On power-up CAN DC STEP transmits the end stop search control to the ASR003 boards. If a confirmation is not produced, an error will set in. | A probable cause is faulty CAN connection or faulty supply to the CAN DC STEP boards. Check on these connections. |
| End stop search time-out. | CAN DC STEP has received control to reset but one or more movements have failed to receive the relative end stop within the allotted time - 10s for simple movements, 30 s for two-motor movements, | Observe the LEDs on CAN DC STEP. Flickering Leds indicate the end stop that has timed out: Check on electric and mechanical connections of the end stop; check that mechanisms are not blocked. |

| Error Type | Problem description | Solution |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|
| Communication Error on Service CAN | Once initialisation has terminated and movements have been re-set, ASR003 expects to receive data regularly from CAN DC STEP. An error will set in if this transmission does not take place and COMM relay will be activated. | Probable causes: interruption of CAN-Bus line of supply to boards CAN DC STEP. |

Default formats


| <i>Format number</i> | <i>Inch</i> | <i>Filter rotation 1</i> | <i>Filter rotation 2</i> |
|-----------------------------|--------------------|---------------------------------|---------------------------------|
| 1 | 12 | 90° | 0° |
| 2 | 9 | 90° | 0° |
| 3 | 7 | 90° | 0° |
| 4 | 6 | 90° | 0° |
| 5 | 4 | 90° | 0° |
| 6 | 3 | 90° | 0° |
| 7 | 8 | 90° | 0° |
| 8 | 2 | 90° | 0° |

Nota:Format programmes need do not contain position values for the strip of 4-position filters. These are always controlled manually and left to the operator's discretion.

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collimatori\acs\asr003_3.jpg

Calibration procedure

The calibration procedure needs to define the relationship between the stand and the variation of the potentiometer signal to determine position values for FFD calculation.

The relation between the movement and the signal must be linear. Power up the board and wait for collimator reset to be concluded.

At this point insert jumper between pin 2 and 4 of S3 to start the actual calibration.

istru_calib_605.fm

Attention: Calibration will not initiate if configuration is not correct.

During calibration ASR003 makes use of manual control pushbuttons and some relays as described in the following table:

| Element | Function description |
|-------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>Input</i> | |
| Iris close pushbutton (terminal M2) | The board performs the measurement, saves the results for table and/or stand position. |
| Iris open pushbutton (terminal M1) | Allows to return to previous calibration step and repeat the process if necessary. |
| <i>Output</i> | |
| MANUAL relay | Intermittent operation over 1 s periods to indicate that the board is on standby for the operator (depending on the calibration step). |
| ERROR relay | Activated if an error is detected. In case of an error the board sets to the point where the error occurred. Example: if at step 2 the table is lowered instead of being raised, an error is signalled and the board requires the correct table position instead of going to step 3). |

The procedure must be as follows:

1. Switch on collimator and wait for the correct setting of initial procedure and axes reset of all boards. Collimator will automatically achieve pre-defined format and will be ready to accept all external signals.
2. Insert the jumper between 2 and 4 pin of J3 programming connector of ASR003 control board. At this point all operations on the board will be inhibited and it will be possible to start the calibration procedure.
3. Correct insertion of jumper will be signalled by blinking led of the relay every second that controls MANUAL output.
4. At this point the system indicates the first calibration phase by activating output 1 and deactivating inputs 2 and 3.
5. Bring FFD to 78 cm and activate input 2 of the board (M39 IRIS CLOSE). When input pushbutton is released the average measurement of the signal coming from FFD potentiometer will be achieved on the system. Subsequently the signal will be compared to absolute minimum value in order to evaluate the measurement taken and saved.
6. At this point the system will deactivate output 1 and active input 2 in order to show that it is ready to obtain the second measurement. Bring FFD to 100 cm and repeat action by activating input 2 (M39 IRIS CLOSE).

7. The system will obtain the value and make the controls compared to previous values; it will then activate output 3 and deactivate output 2 in order to indicate it is ready for the third operation.
8. Operator must set FFD at 120 cm and re-activate input 2 (M39 IRIS CLOSE), the system usually follows up the control.
9. At the end of third operation, the operator can permanently save the obtained data by re-activating inputs 1, 2, 3.
10. At every step of measurement, operator can decide to return to the previous input 1 (M38 IRIS OPEN). In this way the data can be replaced without repeating the whole procedure. In other words, if 78 cm is achieved during the first step, because it is the first step, it will not be possible to return to the previous phase, but in the phase 100 cm it will be possible to return to 78 cm and from 120 cm phase it will be possible to go down to 100 cm

L - SYSTEM CALIBRATION

WARNING:

THE FOLLOWING PROCEDURES REQUIRE THAT X-RADIATION BE PRODUCED. TAKE ADEQUATE PRECAUTIONS TO SEE THAT NO PART OF THE HUMAN BODY IS EXPOSED TO X-RADIATION, DIRECT OR INDIRECT

- Produce X-rays and set the two rectangular shutter to maximum position.
- Widen the round field so that its projection is slightly smaller than the I.I. receptor surface.
- Observe the x-ray field edges: a hazy edge means that the primary cone, near the beam, is cutting the x-ray field instead of the iris shutters located at the collimator exit.
- This condition may be corrected by loosening the three hexagonal centring screws on the collimator attachment plane - see **Figure 1** in this chapter.
- Move the collimator so that the entire x-ray field edge is well defined.
- Tighten the three hexagonal screws.

NOTE: THE HEEL EFFECT WILL CAUSE THE FIELD TOWARD THE CATHODE TO BE SLIGHTLY LESS SHARP THAN ON THE OTHER THREE SIDES. THIS IS NORMAL AND CANNOT BE CORRECTED BY ADJUSTMENTS.

- Close the two pairs of lead rectangular shutters one at a time until the projection is aligned with the edges projected by the round field. If alignment is noted on one shutter only (while the other shutter juts out or is hidden by the round field), close the rectangular field further so that the other shutter is correctly positioned. This determines the amplitude of the correction.
- Should a correction become necessary, loosen the three hexagonal centring screws on the iris support plate - see **Figure 2** in this chapter.
- Place the iris collimator in the desired position and tighten the three screws.
- Rotate the rectangular shutters 360° to check on uniformity of alignment over the round field. If alignment is not uniform, repeat points 2,3, and 4.

Figura 1 R605DASM



Figura 2

M - ADJUSTMENTS

Adjustments are to be performed by technically qualified and authorized personnel only.

IRIS MOTOR ADJUSTMENT:

Motor gear - intermediate gear:

- Remove the two screws “A” and the hexagonal spacer “B” see **Figure 1** in this chapter.
- Turn the motor support plate to access the two motor fixing screws and loosen the two screws.
- Mount the hexagonal spacer “B” and the screw “A” - not tighten
- The plastic gears have a slightly oval shape. Find the critical point and adjust play on this point.
- Tighten the screws “A” and the hexagonal spacer “B”.


Intermediate gear- iris movement gear:

- Loosen the screws “A” and the hexagonal spacer “B” - see **Figure 1** in this chapter.
- Plastic gears are slight oval find the critical point and adjustment play accordingly.
- Tighten the screws “A” and the hexagonal “B”.

ADJUSTMENT OF IRIS PHOTOCELL

The serrated section of the gear ends with an indicator that serves to cut-off the photocell beam on completed iris closure.

- Loosen screws “F” to adjust the photocell - see **Figure 2** in this chapter.
- Move the photocell until a return signal of photocell cut-off is received.
- Tight the two screws “F”.

ADJUSTMENT OF SHUTTER MOTORS

Shutter gear - intermediate gear:

- Loosen the two screws “D” - see **Figure 3** of this chapter.
- Plastic gears are slight oval find the critical point and adjustment play accordingly.
- Tighten the screws “D”.


Motor gear:

- Loosen the two screws “C” - see **Figure 3** of this chapter.
- Plastic gears are slight oval find the critical point and adjustment play accordingly.
- Tighten the screws “C”.

ADJUSTMENT OF SHUTTER PHOTOCELLS


- Close the shutters completely, the indicator will appear on the corner of parallelogram, rotate the parallelogram until the indicator cuts-off the photocell beam. See **Figure 2 and 4** of this chapter.
- Adjust the parallelogram until the bar, with the indicator on its end, is parallel to the bent border of the collimator base.
- If necessary loosen the two screws “B” and adjust the photocell until a return signal indicating cell shut-off is obtained, Led on the edge of the collimator.

Figura 1

collimatori\605\reg_dasm_fig2.gif


Figura 2

collimatori\605\reg_dasm_fig4.gif

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
Figura 3

collimatori\605\reg_dasm_fig5.gif

Figura 4

collimatori\605\reg_dasm_fig3.gif

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N - OPERATION INSTRUCTIONS

System setting

- By supplying the system the collimator carries out Reset of all the axes.
Once Reset is ended the system sets at format 1
- The collimator is ready to operate.
- It is possible to limit the projected format by appropriate inputs or by selecting a new format. By varying FFD the system keeps the chosen format steady.

O - SUBSTITUTIONS, DISASSEMBLY, TRANSPORT

THE FOLLOWING OPERATIONS MUST BE PERFORMED BY AUTHORISED PERSONNEL WITH TECHNICAL TRAINING.

SUBSTITUTION OF THE IRIS MOTOR

- Disconnect supply
- Remove the cover
- Unscrew the screw "A" and the hexagonal spacer "B" - see **Figure 1** - chapter **ADJUSTMENTS**
- Rotate the support plate and unscrew the two motor fixing screws.
- Disconnect the board cable from connector J5V1.
- Substitute the motor with an identical item.
- Fasten the screw "A" and the hexagonal spacer "B" - not tighten.
- Re-mount the motor support plate.
- Adjust play between motor- and intermediate-gear - see **ADJUSTMENTS**
- Tighten the screw "A" and the hexagonal spacer "B"
- Connect the motor cable to the board.
- Re-mount the cover.

SUBSTITUTION OF THE IRIS PHOTOCELL

- Disconnect supply.
- Remove the cover
- Unscrew the two screws "F" - see **Figure 1** - chapter **ADJUSTMENTS**.
- Disconnect cables from the board.
- Remove the photocell support
- Substitute the photocell with an identical item.
- Connect the photocell cables to the board on connector J3V1 as follows:
 - Pin 1** not connected
 - Pin 2** not connected
 - Pin 3** black
 - Pin 4** red
 - Pin 5** white
 - Pin 6** green
- Remove the photocell support through the two screws "F"
- Adjust the photocell position - see **ADJUSTMENTS**
- Re-mount the cover.

SUBSTITUTION OF SHUTTER MOTORS

- Disconnect supply.
- Remove the cover.
- Unscrew the two screws "C" - see **Figure 3** - chapter **ADJUSTMENTS**.
- Disconnect cables from the board: J5V2 rotation - J5V3 travel
- Fasten the two screws "C" not tighten.
- Substitute the motor with an identical item.

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- Adjust play between motor- and intermediate-gear - see **ADJUSTMENTS**
- Tighten the two screws “C” - see **Figure 3** - chapter **ADJUSTMENTS**
- Connect the motor cable to the board.
- Re-mount the cover.

SUBSTITUTION OF PHOTOCELL SHUTTERS.


- Remove the cover
- Unscrew the screws “B” and “E” - see **Figure 2 and 4** - chapter **ADJUSTMENTS**
- Disconnect cables from the board
- Substitute the photocell with an identical item.
- Connect the photocell cables to the board on connector J3V1 as follows:
 - Pin 1** not connected
 - Pin 2** not connected
 - Pin 3** black
 - Pin 4** red
 - Pin 5** white
 - Pin 6** green
- Mount the photocell supports
- Adjust photocell position - see **Figure 3** - chapter **ADJUSTMENTS**
- Re-mount the cover

Disassembly

- Disconnect supply to the collimator.
- Remove the rear cover and disconnect the supply cables.
- Loosen the 4 fixing Allen screws on the upper part of the collimator - **care must be taken not to let the collimator fall.**

Transport and storage

- Suitable packing must be provided for the collimator.
- Place the collimator in a plastic bag to avoid packing material from entering the collimator.
- Use an appropriate box for transport, shipment or storage taking care to protect the collimator from rough handling. This will avoid damage to the collimator during transport shipment or storage.
- Limit Storage conditions:
 - Ambient Temperature = from -40°C to +70°C
 - Relative Humidity = from 10% to 95%
 - Atm. Pressure = from 500 a 1060 hPa

P - ROUTINE MAINTENANCE

To ensure constantly safe performance of the collimator and its compliance with applicable regulations, a maintenance program is indispensable.

It is the Owner's responsibility to supply or arrange for this service

Cleaning recommendations

If the collimator is not protected by an overall cover, the covering must be cleaned as prescribed by the sanitary regulations followed by the operator.

- Disconnect supply
- Use non abrasive cleaning products.

Care must be taken to prevent liquid from entering the collimator. **N.B. the collimator cover is not watertight.**

Do not re-apply power if inflammable liquids have leaked into the collimator.
See the following Maintenance Instructions.

Recommended maintenance programme:

Ralco suggests a yearly servicing programme. However shorter intervals are advisable when the collimator is subject to heavy workloads.

- Re-calibration of the collimator will be necessary whenever the x-ray tube is changed.
Calibration procedures must be performed as described in this manual.
- Check that the screws which serve to secure the collimator to the flange/tube adapter are correctly tightened.
- Remove the covers and panels from collimator. Inspect the moving parts for signs of wear or damage.
- Check the electric system and substitute parts that show wear.
- Check the Lexan panel and substitute if necessary.
- Clean the collimator with a soft cloth paying particular attention to the Lexan window. **Do not** use abrasive or inflammable cleaning products.
- Sparingly lubricate the moving parts using graphite oil.
- Wipe away all excess oil and remount the cover.

IMPORTANT: MAKE SURE TO TIGHTEN THE M6 ALLEN SCREWS SECURING THE CONTROL TABS.

APPROPRIATE TIGHTENING OF THE ALLEN SCREWS ENSURES SECURE MOUNTING OF THE COLLIMATOR. TIGHTENING FORCE USED MUST NOT EXCEED 0.50 NM.

Q - TROUBLESHOOTING

Should the Collimator become faulty do use it until it is completely repaired.

The use of a faulty collimator might impair the safety of the operator and patient.

Before returning the collimator to Ralco for repair, please make sure that it isn't one of the following problems to cause the fault.

| problem | cause | solution |
|-----------------------------------------|----------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
| The Shutters fail to move | The collimator is not controlled correctly | Check the board controlling the motor - See the instructions provided by the board manufacturer. |
| | The motor is faulty. | Substitute the motor - See Chapter SUBSTITUTIONS in this manual. |
| The shutters fail to position correctly | Faulty electronic control of the collimator | Check the board controlling the motor - See the instructions provided by the board manufacturer. |
| | Photocells are not positioned correctly for reset operations.. | Check photocell positioning - See Chapter CALIBRATIONS in this manual. |

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NOTE: when ordering spare parts, the customer is requested to specify the collimator model and serial number.

R 605 DASM

| Cod. | Descrizione | Description |
|--------|---------------------------------------------------|---------------------------------------------|
| RS 476 | Motore passo-passo iride e lamelle | stepper motor for iris and shutters |
| RS 786 | Cono ottone | Brass cone |
| RS 665 | Fotocellula | photocell |
| RS 477 | Flangia del montaggio, collimatore tubo | Mounting flange, collimator tube connection |
| RS 835 | Attacco tubo | Connection tube |
| RS 836 | Copertura | Housing |
| RS 837 | Scheda elettronica CAN DC STEP Iride | Can DC Step board, Iris |
| RS 838 | Scheda elettronica CAN DC STEP, movimento lamella | Can DC Step board, shutter movement |
| RS 839 | Scheda elettronica CAN DC STEP, rotazione lamella | Can DC Step board, shutter rotation |
| RS 840 | Scheda elettronica CAN DC STEP, movimento filtro | Can DC Step board, filter movement |
| RS 841 | Scheda elettronica CAN DC STEP, rotazione filtro | Can DC Step board, filter rotation |
| RS 847 | ASR003 PCB scheda esterna | ASR003 PCB outside collimator |

R 605/025/DASM

| Cod. | Descrizione | Description |
|-------------|---------------------------------------------------|---------------------------------------------|
| RS 476 | Motore passo-passo iride e lamelle | stepper motor for iris and shutters |
| RS 786 | Cono ottone | Brass cone |
| RS 665 | Fotocellula | photocell |
| RS 477 | Flangia del montaggio, collimatore tubo | Mounting flange, collimator tube connection |
| RS 835 | Attacco tubo | Connection tube |
| RS 836 | Copertura | Housing |
| RS 848 | Scheda elettronica CAN DC STEP Iride | Can DC Step board, Iris |
| RS 849 | Scheda elettronica CAN DC STEP, movimento lamella | Can DC Step board, shutter movement |
| RS 850 | Scheda elettronica CAN DC STEP, rotazione lamella | Can DC Step board, shutter rotation |

R 605/027/DASM

| Cod. | Descrizione | Description |
|-------------|---------------------------------------------------|---------------------------------------------|
| RS 476 | Motore passo-passo iride e lamelle | stepper motor for iris and shutters |
| RS 786 | Cono ottone | Brass cone |
| RS 665 | Fotocellula | photocell |
| RS 845 | Flangia del montaggio, collimatore tubo | Mounting flange, collimator tube connection |
| RS 835 | Attacco tubo | Connection tube |
| RS 846 | Copertura | Housing |
| RS 842 | Scheda elettronica CAN DC STEP Iride | Can DC Step board, Iris |
| RS 843 | Scheda elettronica CAN DC STEP, movimento lamella | Can DC Step board, shutter movement |
| RS 844 | Scheda elettronica CAN DC STEP, rotazione lamella | Can DC Step board, shutter rotation |

Etichetta/Label

R605/DASM


etichette\605_DASM.gif

ricambi_dasm.fm


Esploso/Parts Breakdown:

R605DASM :

esploso\605_170_dasm.gif


R605/025/DASM

esploso\605_DASM.gif

ricambi_dasm.fm

R605/027/DASM

esploso\R605_027_dasm.gif

S - GENERAL

REPAIRS

- Return the collimator to Ralco at the customer's expense if the unit is out of warranty.
- Provide the collimator with a detailed description (in Italian or English) of the functional problems and/or faults. It is important to indicate whether a repair or a complete overhaul is required.
- Our Quality Control will test the collimator.
- If the repair involved is extensive, Ralco will contact the customer to advise on the repair or possible substitution.

END OF LIFE DISPOSAL

Your collimator contains materials which can be recycled and reused. Specialised companies can recycle your product to increase the amount of reusable materials and to minimise the amount of materials to be disposed of.

The collimator does not contain polluting materials or products with the exception of the lead that composes the shutters - avoid direct contact with lead especially for prolonged periods.

It is recommended that you observe Local Laws regulating the disposal of your old set. Should this prove impossible, return the collimator to Ralco at the purchaser's expense and Ralco will take care of its correct disposal.

The product requires correct disposal at the end of its life-cycle and according to the current standards; it cannot be considered as normal waste. The unit must be disposed of through certified environmental management concerns or, should you need to replace the unit with new equipment, returned to Ralco. Please contact us if you require further information.

The product contains lead which can be highly contaminating if dispersed incorrectly. The following symbol signified that the products conforms to the environmental requirements of directives 202/95/EC, 2002/96/EC, 2003/108/EC; it must be disposed of correctly at the end of its life-cycle.

WARRANTY

Ralco undertakes to replace and repair any collimator part during a period of 24 months from the date of invoice and cover the labour costs involved.

The warranty applies provided the product has been handled properly in accordance with its operating instructions; presentation is required of the original invoice indicating the date of purchase, the model and serial number as well as other documents originally supplied with the set.

The warranty does not apply if:

- The documents have been altered in any way or made illegible;
- The model or production number on the product has been altered, deleted, removed or made illegible;

- Repairs or product modifications and alterations have been performed by unauthorised persons;
- Damage is caused by misuse or neglect, incorrect installation or accidental damage including but not limited to lighting, water or fire.
- Use of unoriginal spare parts and accessories.

In-warranty spares will be available only upon return to Ralco, at the customer's expense, of the parts considered to be faulty to allow assess the cause of the fault.

Components not covered by this warranty:

Consumable items such as lamps.

Items not produced by Ralco; these items will be accorded to warranty granted by the constructor:

- Motors - 1 year
- Potentiometres: 1 year
- Electronic boards: 1 year
- Laser 12 months

Ralco reserves the right to decide if the collimator is to be repaired or substituted.

Defective material is to be sent to:

RALCO SRL
VIA DEI TIGLI 13/G
20046 BIASSONO (MI) - ITALIA
FAX: ++39-039-2497.799
EMAIL: RALCO@RALCO.IT

SAFETY/RESPONSIBILITY

Ralco adheres to the directives governing manufacturers of electro-medical equipment:

Directive 93/42 CEE para.10 -Legislative Decree n° 46 para.10

Ralco shall not be held responsible when instructions provided in the present manual are not complied with.

Ralco shall not be held responsible if the collimator relates to one or several of the following instances:

- The unit is of Ralco construction to client specifications (no CE marking)
- The unit has been modified by the OEM or end user.
- The unit has been installed without respecting the instructions provided in this manual.
- The unit is used without respecting the instructions provided in this manual.
- The unit has not been subject to routine functional inspection.
- The unit has not been subject to routine maintenance.
- The unit has been repaired with unoriginal spare parts.
- Ralco shall decline all responsibility for any damage, direct or indirect, caused to persons or things by inappropriate accessories.

INFORMATION REGARDING ACCIDENTS THAT HAVE OCCURED WHILE USING THE RADIOLOGICAL COLLIMATOR MUST BE REPORTED IMMEDIATELY TO RALCO SRL.

T - FIGURE E SCHEMI / FIGURES AND DIAGRAMS

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
**NOTE: SHOULD THE ELECTRIC DIAGRAM
REQUESTED BY YOU DIFFER FROM THE STANDARD
UNIT SUPPLIED - PLEASE SEE THE ATTACHMENTS**

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WARNING
COLLIMATOR SUPPLY IS NOT PROTECTED BY A
FUSE. PRIOR TO CONNECTING THE UNIT
CHECK THAT SUPPLY IS PROTECTED BY AN
EXTERNAL FUSE.
FUSE: 3,15 AMP


Figura 1

605\coll_electfig2


figure_schemi_DASM.fm

Figura 2


figure_schemi_DASM.fm

Figura 3

figure_schemi_DASM.fm


Figura 4 - Cablaggio/Wiring DASM

figure_schemi_DASM.fm


elettrico\605_025DASM_fig1.gif

Figura 5

[illegible]

figure_schemi_DASM.fm


Figura 6

figure_schemi_DASM.fm

collimatori\elettrico\R605DASM.nif

Figura 7

figure_schemi_DASM.fm

Technical drawing of a mechanical assembly showing three components:


- ID.710 IRIS**: Located on the left side of the assembly.
- ID.720 FILTER**: Located in the top right section of the assembly.
- ID.730 ROTATION FILTER**: Located in the bottom right section of the assembly.

A small detail 'R' is shown next to the ID.720 FILTER component.

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