

Technical Service Bulletin # 5

Product:	Babytherm 8000 OC
Date:	November 30, 1999
Re:	Babytherm 8000 OC, product improvement Radiant Heater RH 600
Reference Doc:	Documents: <ul style="list-style-type: none">• "Best. LP Spannung", drawing no. 83 05 941 - 3• Conversion Instructions for RH600 (120V-devices)
Reason:	Reliability improvement of the On/Off switch.
Solution:	Exchange of the Voltage PCB. Fitting of an additional metal bracket below the solid state relay on the Relay PCB. New stand-by label for the power switch on the keypad.
Devices affected:	All RH 600 with following serial numbers: <ul style="list-style-type: none">• up to ARMC - 0101 (March 98),• ARMC-0103 to ARMC 0108 (March 98).• ARMC- 1010 (March 98).
When:	At next Service Call or Preventive Maintenance
Additional Info:	Modified units have a new stand-by label for the On/Off switch. This TSB is identical to TSB RH 600W no. 3 Check RH 600 after upgrade according to Test Certificate, see preventive maintenance RH 600.
Cost:	There is no charge for the upgrade kit.
Ordering Info:	Repair and upgrade kit RH 600 (120V) P/N 2M22196
Distribution:	Dräger Service Personnel and Authorized Service Organizations for NCS products.

If you have any questions, please contact Technical Support by phone at 1-800-543-5047 or by fax at 1-215-723-5935

Dräger Medical, Inc.
Technical Product Manager

Andreas Lenke