


PARTS LIST for the
Kodak X-Omat M43 and M43A PROCESSORS
and the
Kodak X-Omat Clinic 1 PROCESSOR



H130_0009DA



HEALTH SCIENCES DIVISION

PLEASE NOTE The information contained herein is based on the experience and knowledge relating to the subject matter gained by Eastman Kodak Company prior to publication.

No patent license is granted by this information.

Eastman Kodak Company reserves the right to change this information without notice, and makes no warranty, express or implied, with respect to this information. Kodak shall not be liable for any loss or damage, including consequential or special damages, resulting from any use of this information, even if loss or damage is caused by Kodak's negligence or other fault.



Possible damage from electrostatic discharge.



This manual is intended for use by qualified service personnel only.

Table of Contents


Description	Page
Covers and Panels	1-1
Roller Transport	2-1
Main Drive Assembly	3-1
Dryer Assembly	4-1
Plumbing	5-1
Electrical	6-1
Miscellaneous	7-1
Numerical Index	8-1
Alphabetical Index	9-1
Publication Change Notice Table	10-1

Section 1: Covers and Panels

List of Figures

Description	Page
Covers and Panels	1-2
Display Panels	1-4

Figure 1-1 Covers and Panels




H130_0168ECB
H130_0168EA


Figure 1-1

Item	Part No.	Description	Qty.	Notes
1	1C7748	Guide Shoe	1	
2	2B4910	Cover - Top, M43	1	
	1C7283	Cover - Top, M43A	1	
	2B4915	Cover - Top, Clinic 1	1	
3	696323	Label - Bin, M43 and M43A	2	
4	2B4907	Panel - Non-Drive Side Access, M43 and M43A	1	
	2B4912	Panel - Non-Drive Side Access, Clinic 1	1	
5	852805	Screw	2	
6	535004	Washer - Lock	2	
7	852799	Washer	2	
8	2B4916	Tray - Feed	1	
9	2B4908	Panel - Front Access, M43 and M43A	1	
	2B4913	Panel - Front Access, Clinic 1	1	
	955976	Label - Drain	1	Not shown.
10	918560	Alarm - Audio	1	LS1
11	459548	Plate - Floor	4	
12	1C4848	Screw - Leveling, Long	3	
	459229	Screw - Leveling, Short	1	Near the Replenishment Pump.
13	166331	Nut	8	
14	2B4906	Panel - Drive Side Access, M43 and M43A	1	
	2B4911	Panel - Drive Side Access, Clinic 1	1	
15	2B4909	Panel - Back Access, M43 and 43A	1	
	2B4914	Panel - Back Access, Clinic 1	1	
16	559793	Washer	1	
17	242038	O-Ring	1	Used to hold the Cover of the Electrical Box open.
18	166956	Wire Tie	1	
19	852777	Screw	9	
20	244917	Retainer	8	
21	852799	Washer	8	
22	535004	Washer - Lock	8	
23	852777	Screw	8	
24	963312	Magnet - Top Cover	1	For All Processors S/N 350 and Above, or with Mod 1 installed.

Figure 1-2 Display Panels
Top: M43 and M43A PROCESSOR
Bottom: Clinic 1 PROCESSOR



H130_0003ACF
H130_0003AA



H130_0018ACD
H130_0018AA

Figure 1–2


Item	Part No.	Description	Qty.	Notes
1	918554	Potentiometer - Dryer Heater	1	R3. Not shown.
	918560	Alarm - Audio	1	LS1. Not shown.
	914441	Circuit Breaker - CB1, M43A and Clinic 1	1	
	918552	Circuit Breaker - CB1, M43	1	
2	968762	Knob - Dryer Temperature Control	1	
3	963916	Display - Developer Temperature, M43 and M43A	1	
4	918566	Indicator Assembly	1	
5	955975	Label - Display, M43 and M43A	1	
6	963256	Label - Display, Clinic 1	1	

Section 2: Roller Transport

List of Figures

Description	Page
Complete Rack Assemblies	2-2
Developer and Fixer Rack Assemblies	2-4
Developer and Fixer Racks - Entrance Rollers	2-6
Developer and Fixer Racks - Exit Rollers	2-8
Developer and Fixer Rack - Other Rollers	2-10
Developer and Fixer Racks - Guide Shoe Assemblies	2-12
Wash Rack Assembly	2-14
Wash Rack - Entrance Rollers	2-16
Wash Rack - Exit Rollers	2-18
Wash Rack - Other Rollers	2-20

Figure 2-1 Complete Rack Assemblies



H130_0232DA

Figure 2-1

Item	Part No.	Description	Qty.	Notes
1	470086	Roller - Resilient	8	
2	1C0780	Roller - Rewet	2	
3	489693	Roller	14	
4	1C4022	Roller	4	
5	8B6568	Roller - Phenolic, 11/16 in.	1	
6	594237	Roller - Resilient	4	Includes Gear.
7	968796	Roller - Phenolic, Wash Rack	1	Includes Gear.
8	467016	Roller - Phenolic, 1 in.	1	
9	963235	Roller - Silicone, Wash Rack	2	Includes Shaft.
10	246979	Roller - Phenolic	3	
11	1C7690	Roller	1	

Figure 2-2 Developer and Fixer Rack Assemblies





Figure 2–2

Item	Part No.	Description	Qty.	Notes
1	1C0351	Rack - Fixer	1	Complete Rack Assembly.
	1C0352	Rack - Developer	1	Complete Rack Assembly including Cover and hardware.
	696333	Cover - Developer Rack	1	
2	544194	Screw	3	
	853096	Washer - Lock	3	
	886890	Washer	3	

Figure 2-3 Developer and Fixer Racks - Entrance Rollers




H130_0101ECA
H130_0101EA

Figure 2–3

Item	Part No.	Description	Qty.	Notes
1	696327	Sideplate - Non-Drive	1	Does not include the Manifold Box.
	1C0792	Manifold Box	1	Non-Drive Sideplate of the Developer Rack only. Not shown.
	852807	Screw	1	Not shown.
2	853100	Washer - Lock	1	Not shown.
	852611	Washer	1	Not shown.
	918503	Bearing	1	Developer Rack only.
3	535646	Pin	1	
4	696334	Gear	1	
5	1C0780	Roller - Rewet	1	
6	1C4022	Roller	1	
7	562228	Stud	4	
8	594273	Bearing - Spring Guide	2	
9	166326	E-Ring	6	
10	502464	Pin	1	
11	550774	Gear - Drive	1	
12	963274	Bushing - Thrust	1	
13	240059	Wire Tie - blue	1	Fixer Rack only.
14	590309	Spring - Short	2	
15	779119	Bearing	2	
16	489693	Roller	1	
17	1C0536	Gear - Snap, 1.0 in.	1	
18	1C0373	Gear - 1.0 in.	1	
19	184455	Pin	1	
20	470086	Roller - Resilient	1	
21	918626	Shaft - Drive, Long	1	Developer Rack only.
	1C4018	Shaft - Drive, Short	1	Fixer Rack only.
22	696326	Sideplate - Drive	1	

Figure 2-4 Developer and Fixer Racks - Exit Rollers




H130_0100ECA
H130_0100EA

Figure 2–4

Item	Part No.	Description	Qty.	Notes
	918503	Bearing	1	Fixer Rack only. Not shown.
1	166326	E-Ring	9	
2	470086	Roller - Resilient	2	Includes Gear.
3	918626	Shaft - Drive, Long	1	Fixer Rack only.
	1C4018	Shaft - Drive, Short	1	Developer Rack only.
4	696328	Shaft	1	
5	1C0537	Gear - Snap, 1.25 in.	1	
6	489693	Roller	1	
7	1C0536	Gear - Snap, 1.0 in.	1	
8	1C4022	Roller	1	
9	779119	Bearing	4	
10	246580	Spring	2	
11	511118	Spacer	2	
12	240058	Wire Tie - red	1	Developer Rack only.
13	963274	Bushing - Thrust	1	
14	550774	Gear - Drive	1	
15	502464	Pin	1	
16	562228	Stud	4	

Figure 2-5 Developer and Fixer Rack - Other Rollers




H130_0099ECA
H130_0099EA

Figure 2–5

Item	Part No.	Description	Qty.	Notes
1	696327	Sideplate - Non-Drive	1	Developer Rack does not include the Manifold Box 1C0792.
2	562228	Stud	10	
3	1C0534	Gear - 1.25 in.	2	
4	1C0373	Gear - 1.0 in.	1	
5	459343	Washer	3	
6	852682	Screw	3	
7	696326	Sideplate - Drive	1	
8	1C0537	Gear - Snap, 1.25 in.	5	
9	489693	Roller	5	
10	1C0536	Gear - Snap, 1.0 in.	1	

Figure 2-6 Developer and Fixer Racks - Guide Shoe Assemblies




H130_0098ECA
H130_0098EA

Figure 2–6

Item	Part No.	Description	Qty.	Notes
1	696332	Guide Shoe Assembly - Entrance, Upper	1	
2	1C0370	Guide Shoe Assembly - Entrance, Lower	1	
3	458219	Tie Rod	4	
4	852807	Screw	8	
5	853100	Washer - Lock	8	
6	474596	Bracket - Guide Shoe	8	
7	852941	Screw	16	
8	914703	Guide Shoe Assembly - Exit, Upper	1	
9	1C7448	Guide Shoe Assembly - Exit, Lower	1	

Figure 2-7 Wash Rack Assembly




H130_0044DA

Figure 2–7

Item	Part No.	Description	Qty.	Notes
	1C7722	Rack - Wash	1	Complete Rack Assembly

Figure 2-8 Wash Rack - Entrance Rollers



H130_0041ECA
H130_0041EA

Figure 2–8

Item	Part No.	Description	Qty.	Notes
1	779119	Bearing	4	
2	459344	Spring	2	
3	918504	Block - Bearing	2	
4	918506	Spacer	1	
5	166326	E-Ring	5	
6	963239	Spacer	1	
7	470086	Roller - Resilient	2	
8	918625	Shaft - Driven	1	
9	1C0369	Shaft - Drive	1	
10	963240	Spacer	1	
11	963274	Bushing - Thrust	1	
12	550774	Gear - Drive	1	
13	502464	Pin	1	
14	458384	Gear	1	
15	544081	Pin	1	
16	963238	Spacer	1	

Figure 2-9 Wash Rack - Exit Rollers





Figure 2–9

Item	Part No.	Description	Qty.	Notes
1	963237	Spacer	4	Small
2	779119	Bearing	4	
3	246580	Spring	2	
4	918504	Block - Bearing	2	Large with Tab
5	963241	Spacer	2	
6	166326	E-Ring	4	
7	963238	Spacer	2	Large
8	963235	Roller - Silicone, Wash Rack	1	
9	458385	Gear - Idler	2	
10	458386	Screw - Shoulder, Small	2	
11	458384	Gear	3	
12	1C0777	Sprocket	1	
13	1C0784	Chain	1	
14	1C0785	Screw - Shoulder	1	
15	544081	Pin	2	
16	1C7690	Roller	1	
17	1C7691	Shaft	1	
18	1C7664	Plate - Locking	2	
19	852871	Screw	2	
20	852799	Washer	2	

Figure 2-10 Wash Rack - Other Rollers



H130_0040ECA
H130_0040EA

Figure 2–10


Item	Part No.	Description	Qty.	Notes
1	459344	Spring	2	
2	852807	Screw	4	
3	853087	Washer - Lock	4	
4	963224	Sideplate, Non-Drive	1	
5	963230	Wash Tube - Upper	1	
6	963231	Wash Tube - Lower	1	
7	963229	Gasket	2	
8	458394	Tie Rod	2	
9	963233	Sideplate, Drive	1	
10	779119	Bearing	4	
11	786449	O-Ring	1	
12	166326	E-Ring	4	
13	1C0777	Sprocket	1	
14	852978	Pin	1	
15	8B6568	Roller - Phenolic, 11/16 in.	1	
16	1C0549	Shaft - Upper	1	Upper, drive roller with Gear and groove. Lower, driven roller with Gear and flange.
17	594237	Roller - Resilient	1	
18	246929	Roller - Phenolic	1	
19	1C0776	Shaft - Lower	1	
20	1C0548	Shaft - Tie Rod	1	
21	467016	Roller - Phenolic, 1 in.	1	
22	852805	Screw	2	
23	535004	Washer - Lock	2	
24	963239	Spacer	1	

Section 3: Main Drive Assembly

List of Figures

Description	Page
Main Drive Assembly	3-2

Figure 3-1 Main Drive Assembly



H130_0059DCB
H130_0059DA

Figure 3–1


Item	Part No.	Description	Qty.	Notes
1	918630	Shaft	1	Includes 2 Washers.
2	491893	Collar	1	
3	918511	Bearing - Thrust	1	
4	918512	Bearing	2	
5	853000	Pin	5	
6	748956	Worm	5	
7	1C4017	Fan	1	For M43A & Clinic 1 Processors. Includes Fan. For M43 Processor. Includes Fan.
8	1C7002	Motor - 115 V	1	
	1C7003	Motor - 230 V	1	
9	852799	Washer	2	
10	853004	Washer - Lock	2	
11	852805	Screw	2	
12	1C4012	Bracket - Drive Motor Mount	1	
13	852611	Washer	3	
14	853100	Washer - Lock	3	
15	852682	Screw	3	For M43 Processor. Includes Setscrew.
16	3B7037	Chain	1	
17	3B7039	Sprocket, 50 Hz	1	
	3B7038	Sprocket, 60 Hz	1	For M43A & Clinic 1 Processors. Includes Setscrew.
18	852558	Setscrew	2	For both 50 & 60 Hz. Includes Setscrew.
19	3B7040	Sprocket	1	

Section 4: Dryer Assembly

List of Figures

Description	Page
Dryer Assembly	4-2
Dryer Rack - Drive Rollers	4-4
Dryer Rack - Driven Rollers	4-6
Air Tubes	4-8
Dryer System	4-10

Figure 4-1 Dryer Assembly




H130_0048DA

Figure 4-1

Item	Part No.	Description	Qty.	Notes
	8B6741	Dryer Assembly	1	Rack Assembly Complete

Figure 4-2 Dryer Rack - Drive Rollers




H130_0037ECA
H130_0037EA

Figure 4-2

Item	Part No.	Description	Qty.	Notes
1	459344	Spring	6	Also shown on Page 4-6.
2	166326	E-Ring	6	
3	779119	Bearing	6	
4	246979	Roller - Phenolic	1	Exit Roller Center and Entrance Rollers
5	594237	Roller - Resilient	2	
6	143403	Nut - Lock	1	
7	963261	Bracket - Tensioner	1	
8	482271	Washer	1	
9	963260	Pulley - Idle	1	
10	990010	Screw - Shoulder	1	
11	963221	Belt - Drive	1	
12	5B6316	Screw	1	
13	852553	Setscrew	3	Includes 1 Setscrew each.
14	963220	Pulley	3	
15	960420	Washer	1	
16	458384	Gear	1	
17	544081	Pin	1	
18	963206	Sideplate - Drive	1	
19	963218	Shaft - Lower Dryer	3	
20	963207	Sideplate - Non-Drive	1	
21	458391	Shaft - Roller	1	
22	458949	Idler Roller Assembly	1	
23	1C7700	Spacer	1	

Figure 4-3 Dryer Rack - Driven Rollers




H130_0256ECA
H130_0256EA

Figure 4-3

Item	Part No.	Description	Qty.	Notes
1	459344	Spring	6	Also shown on Page 4-4.
2	166326	E-Ring	6	
3	779119	Bearing	6	
4	246979	Roller - Phenolic	2	Entrance and Center Roller Exit Roller
5	594237	Roller - Resilient	1	
6	963217	Shaft	3	
7	594237	Roller - Idler	1	
8	458391	Shaft - Idler	1	
9	1C7700	Spacer	1	

Figure 4-4 Air Tubes




H130_0035ECA
H130_0035EA

Figure 4-4

Item	Part No.	Description	Qty.	Notes
1	918663	Cover - Top, Dryer	1	
2	852682	Screw	8	
3	853100	Washer - Lock	8	
4	852611	Washer	8	
5	852807	Screw	8	
6	853100	Washer - Lock	8	
7	458394	Tie Rod	3	
8	963208	Tube - Air	8	
9	562609	Spring	8	
10	918664	Cover - Bottom, Dryer	1	

Figure 4-5 Dryer System



H130_0063ECB
H130_0063EA

Figure 4-5


Item	Part No.	Description	Qty.	Notes
1	963209	Plenum	1	
2	696339	Baffle	1	
3	1C7785	Gasket - Tank/Blower	1	
4	918516	Bracket - Motor	1	
	852830	Screw	2	Not shown.
	535004	Washer - Lock	2	Not shown.
	852799	Washer	2	Not shown.
5	853080	Nut	2	
6	433708	Washer - Star	4	
7	918544	Motor - Blower, 115 V	1	For M43A & Clinic 1 Processors.
	918570	Motor - Blower, 230 V	1	For M43 Processor.
8	165626	Screw	3	
9	853096	Washer - Lock	3	
10	1C7000	Heater - Dryer	1	
11	372056	Screw	2	
12	918550	Thermistor - Dryer	1	
13	5B6310	Bracket	1	
14	943375	Switch - Air Vane	1	Includes Sleevings.
15	852682	Screw	3	
16	853100	Washer - Lock	3	
17	696146	Harness - Air Vane Switch	1	
18	914599	Switch - DC Interlock S6	1	
19	963312	Magnet - Top Cover	1	For All Processors S/N 350 and Above, or with Mod 1 installed.
20	556032	Screw	2	
21	191098	Washer - Lock	4	
22	852619	Nut	2	
23	575822	Washer - Flat	2	

Section 5: Plumbing

List of Figures

Description	Page
Developer and Fixer Recirculation Pump	5-2
Replenishment Pump	5-4
Wash Pump	5-6
Drain Valve Assembly	5-8
Water Input Solenoid and Diverter Valve	5-10
Developer Thermistor	5-12
Developer Heater	5-14
Developer Plumbing Diagram	5-16
Fixer Plumbing Diagram	5-18
Wash Plumbing Diagram, M43 and M43A Processors	5-20
Wash Plumbing Diagram, Clinic 1 Processor	5-22

Figure 5-1 Developer and Fixer Recirculation Pump




H130_0062DCA
H130_0062DA

Figure 5–1

Item	Part No.	Description	Qty.	Notes
1	918543	Pump Assembly - Recirculation, 115 V	1	For M43A & Clinic 1 Processors.
1	918569	Pump Assembly - Recirculation, 230 V	1	For M43 Processor.
	619284	Screw	1	Not shown.
	535004	Washer - Lock	1	Not shown.
	852799	Washer	1	Not shown.
2	1C7009	Housing - Pump	1	
3	1C7010	Impeller	1	
4	473955	Shaft	1	
5	473954	Washer	1	
6	463363	O-Ring	1	
7	1C7008	Cover - Pump	1	
8	168843	Screw	7	

Figure 5-2 Replenishment Pump




H130_0072ECA
H130_0072EA

Figure 5–2

Item	Part No.	Description	Qty.	Notes
1	968766	Pump Assembly - Replenishment, 115 V	1	For M43A and Clinic 1 Processors. Also order 2 blue Wire Ties 240059 and 2 red Wire Ties 240058.
1	968767	Pump Assembly - Replenishment, 230 V	1	For M43 Processor. Also order 2 blue Wire Ties 240059 and 2 red Wire Ties 240058.
	619284	Screw	2	Not shown.
	535004	Washer - Lock	2	Not shown.
	852799	Washer	2	Not shown.
2	520390	Connector - Elbow	2	
3	520391	Nut - Connector	2	
4	520388	O-Ring	2	
5	607817	Valve Assembly	2	
6	1C7005	Body - Valve	1	
7	535639	O-Ring	2	
8	1C7006	Bellows	1	
9	594838	Cup - Displacement	1	
10	594839	Nut - Coupling	1	
11	1C7007	Rod - Connecting	1	
12	5B6289	Motor Assembly - 115 V	1	Includes Gear Box.
12	5B6290	Motor Assembly - 230 V	1	Includes Gear Box.

Figure 5-3 Wash Pump




H130_0075ECA
H130_0075EA

Figure 5–3

Item	Part No.	Description	Qty.	Notes
1	918545	Pump Assembly - Wash Recirculation, 115 V	1	For M43A and Clinic 1 Processors. Includes items 2 - 7.
1	918538	Pump Assembly - Wash Recirculation, 230 V	1	For M43 Processor. Includes items 2 - 7.
	619284	Screw	1	Not shown.
	535004	Washer - Lock	1	Not shown.
	852799	Washer	1	Not shown.
2	1C7011	Body - Pump	1	
3	1C7012	Bearing - Thrust	2	
4	1C7013	Shaft - Pump	1	
5	1C7014	Impeller and Magnet Assembly	1	
6	1C7015	O-Ring	1	
7	1C7016	Housing - Impeller, and Adapter	1	
8		Magnet - Cup	1	Part not available separately. Included with item 9 Wash Pump Motor.
9	5B6287	Motor - Wash Pump, 115 V	1	
9	5B6288	Motor - Wash Pump, 230 V	1	

Figure 5-4 Drain Valve Assembly




H130_0252DCA
H130_0252DA

Figure 5–4

Item	Part No.	Description	Qty.	Notes
1	5B6950	Drain Valve Complete	3	Includes items 2 and 4.
1	8B6945	Valve - Wash Drain, Orifice	1	Includes items 2 and 4.
2	5B6451	Body - Valve	3	
3	918515	Bracket	1	Includes Bracket only, not complete assembly.
4	5B6460	Drain Valve Stem Assembly	3	Used for Developer and Fixer Drain Valves and for Water Drain when in Low Water Usage Mode.
4	8B6948	Wash Valve Stem Assembly - with Orifice	1	Used for water when in Continuous Water Usage Mode.
5	968761	Knob - Drain Valve	3	Includes 2 Setscrews 852553.
6	852553	Setscrew	6	
	914908	Clamp - Thin Band	3	Not shown.

Figure 5-5 Water Input Solenoid and Diverter Valve



H130_0194ECA
H130_0194EA

Figure 5–5

Item	Part No.	Description	Qty.	Notes
1	918641	Solenoid - Water Input	1	L1
2	918937	Harness Assembly	3	Includes the Diode. One needed for the Water Input Solenoid, 2 for the Diverter Valve.
	186171	Screw	2	To hold the Diverter Valve to the Bracket. Not shown.
3	618502	Diode	3	
4	853080	Nut	1	
5	853100	Washer - Lock	1	
6	852611	Washer	1	
7	1C4010	Valve - Diverter	1	Includes Solenoids L2 and L3.
8	918631	Bracket	1	
9	619284	Screw	1	To hold Bracket to Base.
10	535004	Washer - Lock	1	
11	852799	Washer	1	
12	5B6293	Screen	1	
13	5B6294	C-Clip	1	
14	246800	Clamp	3	

Figure 5-6 Developer Thermistor




Figure 5–6

Item	Part No.	Description	Qty.	Notes
	918237	Plug	1	Not shown. For the Fixer Tank only. Also order O-Ring 453342.
	852943	Screw	1	Not shown.
	535004	Washer - Lock	1	Not shown.
1	853136	Nut	1	One needed for the Thermistor and one for the Plug.
2	853011	Washer - Lock	1	
3	453342	O-Ring	2	
4	960420	Washer	1	Also order O-Ring 453342.
5	918549	Thermistor - Developer	1	

Figure 5-7 Developer Heater





Figure 5–7

Item	Part No.	Description	Qty.	Notes
1	3B7276	Screw	3	
2	918548	Heater - Developer, 115 V	1	Also order O-Ring.
	918540	Heater - Developer, 230 V	1	Also order O-Ring.
3	918233	O-Ring - Heater	2	
	918238	Plug	1	Not shown. For Fixer Tank only. Also order O-Ring 918233.

Figure 5-8 Developer Plumbing Diagram




H130_0123ECB
H130_0123EA

Figure 5–8

Item	Part No.	Description	Qty.	Notes
1c	452990	Tubing - Bulk, 3/8-in. ID	6	Order by the foot. For use external to processor. Order by the foot. For use external to processor.
	760477	Tubing - Bulk, 3/4-in. ID, Drain		
	246808	Clamp - #10		
2	918674	Tubing - Formed	2	
3	1C0362	Tubing - Formed	1	
4	1C0361	Tubing - Formed	1	
5	3B7240	Tubing - 1/4-in. ID, 5 ft. length	3	
6	240058	Wire Tie - red		
7	242072	Adapter	1	
8c	246804	Clamp - #12	3	
9	1C0545	Orifice - Developer	1	
10	918248	Tubing - Formed	1	
11c	914908	Clamp - Narrow Band	1	
12	918543	Pump Assembly - Recirculation, 115 V	1	For M43A and Clinic 1 Processors. Also order 2 blue Wire Ties 240059 and 2 red Wire Ties 240058. For M43 Processor. Also order 2 blue Wire Ties 240059 and 2 red Wire Ties 240058.
13	918569	Pump Assembly - Recirculation, 230 V	1	
	968766	Pump Assembly - Replenishment, 115 V	1	
	968767	Pump Assembly - Replenishment, 230 V	1	

Figure 5-9 Fixer Plumbing Diagram




H130_0124ECB
H130_0124EA

Figure 5–9

Item	Part No.	Description	Qty.	Notes
1c	452990	Tubing - Bulk, 3/8-in. ID	6	Order by the foot. For use external to processor. Order by the foot. For use external to processor.
	760477	Tubing - Bulk, 3/4-in. ID, Drain		
	246808	Clamp - #10		
2c	246804	Clamp - #12	3	
3	963212	Tubing - Formed	1	
4	918242	Tubing - Formed	1	
5	3B7240	Tubing - 1/4-in. ID, 5 ft. length	3	
6	240059	Wire Tie - blue		
7	242072	Adapter	1	
8	918243	Tubing - Formed	1	
9	1C0546	Orifice - Fixer	1	
10	963213	Tubing - Formed	1	
11	918249	Tubing - Formed	1	
12c	914908	Clamp - Narrow Band	1	
13	918543	Pump Assembly - Recirculation, 115 V	1	
14	918569	Pump Assembly - Recirculation, 230 V	1	For M43A and Clinic 1 Processors. Also order 2 blue Wire Ties 240059 and 2 red Wire Ties 240058. For M43 Processor. Also order 2 blue Wire Ties 240059 and 2 red Wire Ties 240058.
	968766	Pump Assembly - Replenishment, 115 V	1	
	968767	Pump Assembly - Replenishment, 230 V	1	

Figure 5-10 Wash Plumbing Diagram, M43 and M43A Processors




H130_0121ECB
H130_0121EA

Figure 5–10

Item	Part No.	Description	Qty.	Notes
1	760477	Tubing - Bulk, 3/4-in. ID, Drain		Order by the foot.
	918231	Exchanger - Heat	2	
	658823	O-Ring	4	Not shown.
2c	918229	Plate	4	Not shown.
	3B7276	Screw	12	Not shown.
	246808	Clamp - #10	10	
3	1C0354	Tubing - Formed	1	
4	1C0358	Tubing - Formed	1	
5	1C0356	Tubing - Formed	1	
6	1C4010	Valve - Diverter	1	Includes Solenoids L2 and L3.
7c	246801	Clamp	10	
8	1C0355	Tubing - Formed	1	
9	918245	Tubing - Formed	1	
10	1C0359	Tubing - Formed	1	
11	918641	Solenoid - Water Input	1	
12c	246800	Clamp	2	
13	1C0367	Tubing - Formed	1	
14	963243	Connector - Y	1	
15	1C0353	Tubing - Formed	1	
16	963215	Tubing - Formed	1	
17	918545	Pump Assembly - Wash Recirculation, 115 V	1	M43A Processor
18	918538	Pump Assembly - Wash Recirculation, 230 V	1	M43 Processor
	1C0360	Tubing - Formed	1	
	1C0357	Tubing - Formed	1	
20	963242	Connector - T	1	
21	535004	Washer - Lock	1	
22	852805	Screw	1	
23	1C0783	Tube Assembly	1	Inside U.K.
23	8B6860	Tube Assembly	1	Outside U.K.
24c	914908	Clamp - Narrow Band	1	

Figure 5-11 Wash Plumbing Diagram, Clinic 1 Processor



H130_0122ECB
H130_0122EA

Figure 5-11


Item	Part No.	Description	Qty.	Notes
	658823	O-Ring	2	Used for Heat Exchanger.
	918229	Plate	2	Used for Heat Exchanger.
	3B7276	Screw	6	Used for Heat Exchanger.
	658823	O-Ring	2	Used in Fixer Tank with Plug. Clinic 1 only.
	696338	Plug	2	Used in Fixer Tank.
	796386	Screw	6	Used in Fixer Tank with Plug. Clinic 1 only.
1	918231	Exchanger - Heat	1	
1	918231	Exchanger - Heat	1	
2c	246808	Clamp - #10	7	
3	1C0363	Tubing - Formed	1	
4	1C0365	Tubing - Formed	1	
5	963242	Connector - T	1	
6c	246801	Clamp	5	
7	918245	Tubing - Formed	1	
8	1C0366	Tubing - Formed	1	
9	918641	Solenoid - Water Input	1	
10c	246800	Clamp	2	
11	1C0367	Tubing - Formed	1	
12	963215	Tubing - Formed	1	
13	1C0357	Tubing - Formed	1	
14	918545	Pump Assembly - Wash Recirculation, 115 V	1	
15	535004	Washer - Lock	1	
16	852805	Screw	1	
17	1C0783	Tube Assembly	1	Inside U.K.
17	8B6860	Tube Assembly	1	Outside U.K.
	760477	Tubing - Bulk, 3/4-in. ID, Drain		Order by the foot.
18	1C0547	Orifice - Wash	1	Clinic 1 only. Install at Pump output.

Section 6: Electrical

List of Figures

Description	Page
Electrical Box	6-2
AC Interlock Switch S5	6-4
DC Interlock Switch S6 for Processors Serial Number 350 and Above, or With Mod 1 Installed	6-6
Display Panel and Line Filter	6-8
Film Detector Assembly	6-10

Figure 6-1 Electrical Box





H130_0092ECA
H130_0092EA

Figure 6-1

Item	Part No.	Description	Qty.	Notes
	1C0791	Varistor - MOV Assembly	2	M43A and Clinic 1 Only. Not shown at TB2.
1	991743	Relay - Solid State, 20 A	3	Order new MOV.
2	968757	Standoff - Nylon, 1 in.	2	
3	943378	Varistor - MOV	5	
4	991740	Relay - Solid State, 10 A	2	Order new MOV.
5	8B6926	EEPROM - 100 Circuit Board, Main Program Version 1.80	1	U5, for M43 and M43A Processors S/N 349 and Below, and without Mod 1 installed.
	8B6925	EEPROM - 100 Circuit Board, Main Program, Version 1.80	1	U5, for Clinic 1 Processor S/N 349 and Below, and without Mod 1 installed.
	8B6926	EEPROM - 100 Circuit Board, Main Program, Version 1.90	1	U5, for M43 and M43A Processors S/N 350 and Above, or with Mod 1 installed.
	8B6925	EEPROM - 100 Circuit Board, Main Program, Version 1.90	1	U5, for Clinic 1 Processor S/N 350 and Above, or with Mod 1 installed.
6	335503	Microprocessor	1	U20.
7	292888	Jack - ESD	1	
8	696346	Circuit Board - 100	1	Includes the Microprocessor but not the EEPROM.
9	1C7282	Power Supply Assembly	1	Includes the Cable.
10	963923	Terminal Block - 6 Position	1	
	696148	Strip - Marker, 6 Position	1	
11	1C4009	Transformer - 20 V AC	1	
12	136963	Fuse - 3/4 A, 115 V, slow blow	1	For M43A and Clinic 1 Processors.
	337468	Fuse - 3/8 A, 230 V, slow blow	1	For M43 Processor.
13	943366	Holder - Fuse	2	
	477506	Carrier - Fuse	2	
14	126665	Fuse - 2 A, 115 V, slow blow	1	For M43A and Clinic 1 Processors.
	569165	Fuse - 1 A, 230 V, slow blow	1	For M43 Processor.
15	943381	Strap - EMI	1	
16	963914	Cover - Solid State Relay	1	
17	781202	Screw - Nylon	2	
18	154242	Screw	2	
19	433708	Washer - Star	2	
20	853080	Nut	2	

Figure 6-2 AC Interlock Switch S5




H130_0095DCA
H130_0095DA

Figure 6–2

Item	Part No.	Description	Qty.	Notes
1	746731	Switch - AC Interlock S5	1	
2	918234	Bracket	1	
3	853080	Nut	1	
4	433708	Washer - Star	2	
5	696144	Boot	1	
6	182915	Wire Tie	1	
7	182916	Wire Tie	1	
8	852799	Washer	2	
9	535004	Washer - Lock	2	
10	852683	Screw	2	Order TL-6240. Not shown.
		Connector Kit	1	

**Figure 6-3 DC Interlock Switch S6 for Processors
Serial Number 350 and Above, or With Mod 1 Installed**




H130_0228GCA
H130_0228GA

Figure 6–3

Item	Part No.	Description	Qty.	Notes
1	914599	Switch - DC Interlock S6	1	For All Processors S/N 350 and Above, or with Mod 1 installed.
2	963312	Magnet - Top Cover	1	
3	852619	Nut	2	
4	556032	Screw	2	
5	191098	Washer - Lock	4	
6	575822	Washer - Flat	2	
7	5B6310	Bracket	1	

Figure 6-4 Display Panel and Line Filter




H130_0083DCA
H130_0083DA

Figure 6-4

Item	Part No.	Description	Qty.	Notes
1	918537	Filter - Line	1	FL1
2	853080	Nut	1	Used in recent serial numbers. Screw and Flat Washer, not shown, used in early serial numbers.
3	583100	Washer - Lock	1	
4	968762	Knob - Dryer Temperature Control	1	
5	918566	Indicator Assembly	1	Includes DS1, DS2, DS3, 3 Jam Nuts, and Connector.
6	963916	Display - Developer Temperature, M43 and M43A	1	Includes Bracket and Nuts.
7	918559	Strip - Marker	1	
8	598808	Grounding Lug	1	
9	852799	Washer	1	
10	535004	Washer - Lock	1	
11	852683	Screw	1	
12	852976	Screw	1	
13	853096	Washer - Lock	1	
14	968763	Terminal Block	1	TB1
15	918560	Alarm - Audio	1	LS1
16	575822	Washer - Flat	1	
17	191098	Washer - Lock	1	
18	575821	Screw	1	
19	443700	Washer - Star	2	
20	853080	Nut	1	
21	914441	Circuit Breaker - CB1, M43A and Clinic 1	1	
	918552	Circuit Breaker - CB1, M43	1	
22	918554	Potentiometer - Dryer Heater	1	R3. Includes Nut and Washer.

Figure 6-5 Film Detector Assembly



H130_0094DCB
H130_0094DA

Figure 6-5


Item	Part No.	Description	Qty.	Notes
1	918557	Sensor - Film	3	Includes 3 Film Sensors.
2	918643	Film Detector Assembly	1	
3	852805	Screw	2	
4	535004	Washer - Lock	2	
5	852799	Washer	2	
6	575821	Screw	6	
7	852619	Nut	6	

Section 7: Miscellaneous

List of Figures

Description	Page
Mounting Stand and Seismic Brackets	7-2
Lighttight Feed Tray Kit	7-4
Kits	7-6
Replenishment Tanks	7-8
Accessories	7-10

Figure 7-1 Mounting Stand and Seismic Brackets




H130_0091DCB
H130_0091DA

Figure 7-1

Item	Part No.	Description	Qty.	Notes
1	261413	Seismic Bracket Kit Mounting Stand - <i>Kodak X-Omat</i> M35, M43, and Clinic 1	1 1	Includes items 2 - 6. CAT 808-1176
2	261898	Floor Plate	4	
3	240104	Bracket	4	
4	261929	Anchor	8	
5	188623	Screw - Stainless Steel	8	
6	173824	Washer - Stainless Steel	8	
7	166331	Nut	4	
8	1C4849	Spacer	4	Provided with Stand CAT 808-1176.
9	1C4848	Screw - Leveling, Long	3	
	459229	Screw - Leveling, Short	1	Near the Replenishment Pump.

Figure 7-2 Lighttight Feed Tray Kit




H130_0120HA

Figure 7-2

Item	Part No.	Description	Qty.	Notes
		Lighttight Feed Tray Kit - <i>Kodak X-Omat</i> M43	1	CAT 188-0335
		Gasket - Foam	1	Not shown.
		Block - Foam	2	Not shown.
		Plate - Face	1	Not shown.
		Screws	4	Not shown.
		Wire Tie	2	Not shown.
		Seal - Light, Lid	1	Not shown.
		Seal - Light, Feed Tray	2	Not shown.
		Center Punch	1	Not shown.
		Gasket - Front Access Panel	1	Not shown.
		Plate - Nut	1	Not shown.

Figure 7-3 Kits




H130_0225DCA
H130_0225DA

Figure 7-3

Item	Part No.	Description	Qty.	Notes
1	261413	Vent Duct Adapter Kit	1	Includes Plate, Screws, and Washers. CAT 143-4943.
		Lighttight Feed Tray Kit - <i>Kodak X-Omat M43</i>	1	CAT 188-0335. See the figure on Page 7-4.
		Seismic Bracket Kit	1	See the figure on Page 7-2.
	264503	Ventilation Fan Kit - 110 V	1	<i>Kodak</i> Auxiliary. Not shown. To be used with item 1, if necessary for correct ventilation.
		Transformer/Circuit Breaker Kit	1	CAT 167-4340 for M43A and Clinic 1 Processor sites not within the 104 - 127 V AC range.
		Transformer/Circuit Breaker Kit	1	CAT 171-0292 for M43 Processor sites not within the 207 - 253 V AC range.
		M43 Through-the-Wall Kit	1	Not shown. CAT 871-3109. Includes: Foam Block, Foam Gasket, Hose Clamps #2 and #12, $\frac{3}{8}$ - and $\frac{3}{4}$ -in. opaque Hose, $\frac{3}{8}$ - and $\frac{3}{4}$ -in. Barbed Fittings, Eyebolts, Nuts, & Washers.
		Clinic 1 Installation Kit	1	Not shown. CAT 863-2754. Includes: Drain Tubing, Drip Tray, Graduated Cylinder, Leveling Foot Washers, Power Cord, Clamps, Splash Guard, & Strainer Assemblies.

Figure 7-4 Replenishment Tanks




H112_0010ACB
H112_0010AC

Figure 7-4

Item	Part No.	Description	Qty.	Notes
1	472261	Strainer Assembly	2	Contains Strainer 457995 and Gasket 613114.
2		Replenishment Tanks	2	Includes Bulkhead Fitting 492578.
		7 Gallon - Rectangular		CAT 150-0537.
		14 Gallon - Round		Developer, CAT 151-1740.
		14 Gallon - Round		Fixer, CAT 151-1765.
		30 Gallon - Round		Developer, CAT 102-2987.
		30 Gallon - Round		Fixer, CAT 102-2961.
		55 Gallon - Round		Developer, CAT 151-1823.
		55 Gallon - Round		Fixer, CAT 151-1849.
3	492578	Fitting - Bulkhead	2	
3	246800	Clamp	2	

Figure 7-5 Accessories



H130_0169DCA
H130_0169DA

Figure 7-5

Item	Part No.	Description	Qty.	Notes
	760477	Tubing - Bulk, 3/4-in. ID, Drain		Order by the foot.
	760476	Tubing - Bulk, 1/2-in. ID		Order by the foot.
	452990	Tubing - Bulk, 3/8-in. ID		Order by the foot.
1	918648	Cylinder - Graduated, 500 mL	1	Prepack 2-Roll Pack.
2	650656	Tubing - Drain, 18 ft., 3/4-in. ID		
3	650938	Gasket	1	
4	613122	Cover - Evaporation	2	
5	1C4011	Tray - Drip	1	
6	1C4019	Guard - Splash	1	

Section 8: Numerical Index

Numerical Index

Part No.	Description	Figure No.
126665	Fuse - 2 A, 115 V, slow blow	6-1
136963	Fuse - 3/4 A, 115 V, slow blow	6-1
143403	Nut - Lock	4-2
154242	Screw	6-1
165626	Screw	4-5
166326	E-Ring	2-3, 2-4, 2-8, 2-9, 2-10, 4-2, 4-3
166331	Nut	1-1, 7-1
166956	Wire Tie	1-1
168843	Screw	5-1
173824	Washer - Stainless Steel	7-1
182915	Wire Tie	6-2
182916	Wire Tie	6-2
184455	Pin	2-3
186171	Screw	5-5
188623	Screw - Stainless Steel	7-1
191098	Washer - Lock	4-5, 6-3, 6-4
1C0351	Rack - Fixer	2-2
1C0352	Rack - Developer	2-2
1C0353	Tubing - Formed	5-10
1C0354	Tubing - Formed	5-10
1C0355	Tubing - Formed	5-10
1C0356	Tubing - Formed	5-10
1C0357	Tubing - Formed	5-10, 5-11
1C0358	Tubing - Formed	5-10
1C0359	Tubing - Formed	5-10
1C0360	Tubing - Formed	5-10
1C0361	Tubing - Formed	5-8
1C0362	Tubing - Formed	5-8
1C0363	Tubing - Formed	5-11
1C0365	Tubing - Formed	5-11
1C0366	Tubing - Formed	5-11
1C0367	Tubing - Formed	5-10, 5-11
1C0369	Shaft - Drive	2-8
1C0370	Guide Shoe Assembly - Entrance, Lower	2-6
1C0373	Gear - 1.0 in.	2-3, 2-5
1C0534	Gear - 1.25 in.	2-5
1C0536	Gear - Snap, 1.0 in.	2-3, 2-4, 2-5
1C0537	Gear - Snap, 1.25 in.	2-4, 2-5
1C0545	Orifice - Developer	5-8
1C0546	Orifice - Fixer	5-9
1C0547	Orifice - Wash	5-11

Numerical Index

Part No.	Description	Figure No.
1C0548	Shaft - Tie Rod	2-10
1C0549	Shaft - Upper	2-10
1C0776	Shaft - Lower	2-10
1C0777	Sprocket	2-9, 2-10
1C0780	Roller - Rewet	2-1, 2-3
1C0783	Tube Assembly	5-10, 5-11
1C0784	Chain	2-9
1C0785	Screw - Shoulder	2-9
1C0791	Varistor - MOV Assembly	6-1
1C0792	Manifold Box	2-3
1C4009	Transformer - 20 V AC	6-1
1C4010	Valve - Diverter	5-5, 5-10
1C4011	Tray - Drip	7-5
1C4012	Bracket - Drive Motor Mount	3-1
1C4017	Fan	3-1
1C4018	Shaft - Drive, Short	2-3, 2-4
1C4019	Guard - Splash	7-5
1C4022	Roller	2-1, 2-3, 2-4
1C4848	Screw - Leveling, Long	1-1, 7-1
1C4849	Spacer	7-1
1C7000	Heater - Dryer	4-5
1C7002	Motor - 115 V	3-1
1C7003	Motor - 230 V	3-1
1C7005	Body - Valve	5-2
1C7006	Bellows	5-2
1C7007	Rod - Connecting	5-2
1C7008	Cover - Pump	5-1
1C7009	Housing - Pump	5-1
1C7010	Impeller	5-1
1C7011	Body - Pump	5-3
1C7012	Bearing - Thrust	5-3
1C7013	Shaft - Pump	5-3
1C7014	Impeller and Magnet Assembly	5-3
1C7015	O-Ring	5-3
1C7016	Housing - Impeller, and Adapter	5-3
1C7282	Power Supply Assembly	6-1
1C7283	Cover - Top, M43A	1-1
1C7448	Guide Shoe Assembly - Exit, Lower	2-6
1C7664	Plate - Locking	2-9
1C7690	Roller	2-1, 2-9
1C7691	Shaft	2-9
1C7700	Spacer	4-2, 4-3
1C7722	Rack - Wash	2-7
1C7748	Guide Shoe	1-1
1C7785	Gasket - Tank/Blower	4-5

Numerical Index

Part No.	Description	Figure No.
240058	Wire Tie - red	2-4, 5-8
240059	Wire Tie - blue	2-3, 5-9
240104	Bracket	7-1
242038	O-Ring	1-1
242072	Adapter	5-8, 5-9
244917	Retainer	1-1
246580	Spring	2-4, 2-9
246800	Clamp	5-5, 5-10, 5-11, 7-4
246801	Clamp	5-10, 5-11
246804	Clamp - #12	5-8, 5-9
246808	Clamp - #10	5-8, 5-9, 5-10, 5-11
246929	Roller - Phenolic	2-10
246979	Roller - Phenolic	2-1, 4-2, 4-3
261413	Seismic Bracket Kit	7-1, 7-3
261898	Floor Plate	7-1
261929	Anchor	7-1
264503	Ventilation Fan Kit - 110 V	7-3
292888	Jack - ESD	6-1
2B4906	Panel - Drive Side Access, M43 and M43A	1-1
2B4907	Panel - Non-Drive Side Access, M43 and M43A	1-1
2B4908	Panel - Front Access, M43 and M43A	1-1
2B4909	Panel - Back Access, M43 and 43A	1-1
2B4910	Cover - Top, M43	1-1
2B4911	Panel - Drive Side Access, Clinic 1	1-1
2B4912	Panel - Non-Drive Side Access, Clinic 1	1-1
2B4913	Panel - Front Access, Clinic 1	1-1
2B4914	Panel - Back Access, Clinic 1	1-1
2B4915	Cover - Top, Clinic 1	1-1
2B4916	Tray - Feed	1-1
335503	Microprocessor	6-1
337468	Fuse - 3/8 A, 230 V, slow blow	6-1
372056	Screw	4-5
3B7037	Chain	3-1
3B7038	Sprocket, 60 Hz	3-1
3B7039	Sprocket, 50 Hz	3-1
3B7040	Sprocket	3-1
3B7240	Tubing - 1/4-in. ID, 5 ft. length	5-8, 5-9
3B7276	Screw	5-7, 5-10, 5-11
433708	Washer - Star	4-5, 6-1, 6-2
443700	Washer - Star	6-4
452990	Tubing - Bulk, 3/8-in. ID	5-8, 5-9, 7-5
453342	O-Ring	5-6
458219	Tie Rod	2-6

Numerical Index

Part No.	Description	Figure No.
458384	Gear	2-8, 2-9, 4-2
458385	Gear - Idler	2-9
458386	Screw - Shoulder, Small	2-9
458391	Shaft - Idler	4-3
458391	Shaft - Roller	4-2
458394	Tie Rod	2-10, 4-4
458949	Idler Roller Assembly	4-2
459229	Screw - Leveling, Short	1-1, 7-1
459343	Washer	2-5
459344	Spring	2-8, 2-10, 4-2, 4-3
459548	Plate - Floor	1-1
463363	O-Ring	5-1
467016	Roller - Phenolic, 1 in.	2-1, 2-10
470086	Roller - Resilient	2-1, 2-3, 2-4, 2-8
472261	Strainer Assembly	7-4
473954	Washer	5-1
473955	Shaft	5-1
474596	Bracket - Guide Shoe	2-6
477506	Carrier - Fuse	6-1
482271	Washer	4-2
489693	Roller	2-1, 2-3, 2-4, 2-5
491893	Collar	3-1
492578	Fitting - Bulkhead	7-4
502464	Pin	2-3, 2-4, 2-8
511118	Spacer	2-4
520388	O-Ring	5-2
520390	Connector - Elbow	5-2
520391	Nut - Connector	5-2
535004	Washer - Lock	1-1, 2-10, 4-5, 5-1, 5-2, 5-3, 5-5, 5-6, 5-10, 5-11, 6-2, 6-4, 6-5
535639	O-Ring	5-2
535646	Pin	2-3
544081	Pin	2-8, 2-9, 4-2
544194	Screw	2-2
550774	Gear - Drive	2-3, 2-4, 2-8
556032	Screw	4-5, 6-3
559793	Washer	1-1
562228	Stud	2-3, 2-4, 2-5
562609	Spring	4-4
569165	Fuse - 1 A, 230 V, slow blow	6-1

Numerical Index

Part No.	Description	Figure No.
575821	Screw	6-4, 6-5
575822	Washer - Flat	4-5, 6-3, 6-4
583100	Washer - Lock	6-4
590309	Spring - Short	2-3
594237	Roller - Idler	4-3
594237	Roller - Resilient	2-1, 2-10, 4-2, 4-3
594273	Bearing - Spring Guide	2-3
594838	Cup - Displacement	5-2
594839	Nut - Coupling	5-2
598808	Grounding Lug	6-4
5B6287	Motor - Wash Pump, 115 V	5-3
5B6288	Motor - Wash Pump, 230 V	5-3
5B6289	Motor Assembly - 115 V	5-2
5B6290	Motor Assembly - 230 V	5-2
5B6293	Screen	5-5
5B6294	C-Clip	5-5
5B6310	Bracket	4-5, 6-3
5B6316	Screw	4-2
5B6451	Body - Valve	5-4
5B6460	Drain Valve Stem Assembly	5-4
5B6950	Drain Valve Complete	5-4
607817	Valve Assembly	5-2
613122	Cover - Evaporation	7-5
618502	Diode	5-5
619284	Screw	5-1, 5-2, 5-3, 5-5
650656	Tubing - Drain, 18 ft., 3/4-in. ID	7-5
650938	Gasket	7-5
658823	O-Ring	5-10, 5-11
696144	Boot	6-2
696146	Harness - Air Vane Switch	4-5
696148	Strip - Marker, 6 Position	6-1
696323	Label - Bin, M43 and M43A	1-1
696326	Sideplate - Drive	2-3, 2-5
696327	Sideplate - Non-Drive	2-3, 2-5
696328	Shaft	2-4
696332	Guide Shoe Assembly - Entrance, Upper	2-6
696333	Cover - Developer Rack	2-2
696334	Gear	2-3
696338	Plug	5-11
696339	Baffle	4-5
696346	Circuit Board - 100	6-1
746731	Switch - AC Interlock S5	6-2
748956	Worm	3-1

Numerical Index

Part No.	Description	Figure No.
760476	Tubing - Bulk, 1/2-in. ID	7-5
760477	Tubing - Bulk, 3/4-in. ID, Drain	5-8, 5-9, 5-10, 5-11, 7-5
779119	Bearing	2-3, 2-4, 2-8, 2-9, 2-10, 4-2, 4-3
781202	Screw - Nylon	6-1
786449	O-Ring	2-10
796386	Screw	5-11
852553	Setscrew	4-2, 5-4
852558	Setscrew	3-1
852611	Washer	2-3, 3-1, 4-4, 5-5
852619	Nut	4-5, 6-3, 6-5
852682	Screw	2-5, 3-1, 4-4, 4-5
852683	Screw	6-2, 6-4
852777	Screw	1-1
852799	Washer	1-1, 2-9, 3-1, 4-5, 5-1, 5-2, 5-3, 5-5, 6-2, 6-4, 6-5
852805	Screw	1-1, 2-10, 3-1, 5-10, 5-11, 6-5
852807	Screw	2-3, 2-6, 2-10, 4-4
852830	Screw	4-5
852871	Screw	2-9
852941	Screw	2-6
852943	Screw	5-6
852976	Screw	6-4
852978	Pin	2-10
853000	Pin	3-1
853004	Washer - Lock	3-1
853011	Washer - Lock	5-6
853080	Nut	4-5, 5-5, 6-1, 6-2, 6-4
853087	Washer - Lock	2-10
853096	Washer - Lock	2-2, 4-5, 6-4
853100	Washer - Lock	2-3, 2-6, 3-1, 4-4, 4-5, 5-5
853136	Nut	5-6
886890	Washer	2-2
8B6568	Roller - Phenolic, 11/16 in.	2-1, 2-10
8B6741	Dryer Assembly	4-1
8B6860	Tube Assembly	5-10, 5-11
8B6925	EEPROM - 100 Circuit Board, Main Program, Version 1.80	6-1

Numerical Index

Part No.	Description	Figure No.
8B6925	EEPROM - 100 Circuit Board, Main Program, Version 1.90	6-1
8B6926	EEPROM - 100 Circuit Board, Main Program Version 1.80	6-1
8B6926	EEPROM - 100 Circuit Board, Main Program, Version 1.90	6-1
8B6945	Valve - Wash Drain, Orifice	5-4
8B6948	Wash Valve Stem Assembly - with Orifice	5-4
914441	Circuit Breaker - CB1, M43A and Clinic 1	1-2, 6-4
914599	Switch - DC Interlock S6	4-5, 6-3
914703	Guide Shoe Assembly - Exit, Upper	2-6
914908	Clamp - Narrow Band	5-8, 5-9, 5-10
914908	Clamp - Thin Band	5-4
918229	Plate	5-10, 5-11
918231	Exchanger - Heat	5-10, 5-11
918233	O-Ring - Heater	5-7
918234	Bracket	6-2
918237	Plug	5-6
918238	Plug	5-7
918242	Tubing - Formed	5-9
918243	Tubing - Formed	5-9
918245	Tubing - Formed	5-10, 5-11
918248	Tubing - Formed	5-8
918249	Tubing - Formed	5-9
918503	Bearing	2-3, 2-4
918504	Block - Bearing	2-8, 2-9
918506	Spacer	2-8
918511	Bearing - Thrust	3-1
918512	Bearing	3-1
918515	Bracket	5-4
918516	Bracket - Motor	4-5
918537	Filter - Line	6-4
918538	Pump Assembly - Wash Recirculation, 230 V	5-3, 5-10
918540	Heater - Developer, 230 V	5-7
918543	Pump Assembly - Recirculation, 115 V	5-1, 5-8, 5-9
918544	Motor - Blower, 115 V	4-5
918545	Pump Assembly - Wash Recirculation, 115 V	5-3, 5-10, 5-11
918548	Heater - Developer, 115 V	5-7
918549	Thermistor - Developer	5-6
918550	Thermistor - Dryer	4-5
918552	Circuit Breaker - CB1, M43	1-2, 6-4
918554	Potentiometer - Dryer Heater	1-2, 6-4
918557	Sensor - Film	6-5
918559	Strip - Marker	6-4
918560	Alarm - Audio	1-1, 1-2, 6-4
918566	Indicator Assembly	1-2, 6-4
918569	Pump Assembly - Recirculation, 230 V	5-1, 5-8, 5-9
918570	Motor - Blower, 230 V	4-5

Numerical Index

Part No.	Description	Figure No.
918625	Shaft - Driven	2-8
918626	Shaft - Drive, Long	2-3, 2-4
918630	Shaft	3-1
918631	Bracket	5-5
918641	Solenoid - Water Input	5-5, 5-10, 5-11
918643	Film Detector Assembly	6-5
918648	Cylinder - Graduated, 500 mL	7-5
918663	Cover - Top, Dryer	4-4
918664	Cover - Bottom, Dryer	4-4
918674	Tubing - Formed	5-8
918937	Harness Assembly	5-5
943366	Holder - Fuse	6-1
943375	Switch - Air Vane	4-5
943378	Varistor - MOV	6-1
943381	Strap - EMI	6-1
955975	Label - Display, M43 and M43A	1-2
955976	Label - Drain	1-1
960420	Washer	4-2, 5-6
963206	Sideplate - Drive	4-2
963207	Sideplate - Non-Drive	4-2
963208	Tube - Air	4-4
963209	Plenum	4-5
963212	Tubing - Formed	5-9
963213	Tubing - Formed	5-9
963215	Tubing - Formed	5-10, 5-11
963217	Shaft	4-3
963218	Shaft - Lower Dryer	4-2
963220	Pulley	4-2
963221	Belt - Drive	4-2
963224	Sideplate, Non-Drive	2-10
963229	Gasket	2-10
963230	Wash Tube - Upper	2-10
963231	Wash Tube - Lower	2-10
963233	Sideplate, Drive	2-10
963235	Roller - Silicone, Wash Rack	2-1, 2-9
963237	Spacer	2-9
963238	Spacer	2-8, 2-9
963239	Spacer	2-8, 2-10
963240	Spacer	2-8
963241	Spacer	2-9
963242	Connector - T	5-10, 5-11
963243	Connector - Y	5-10
963256	Label - Display, Clinic 1	1-2
963260	Pulley - Idle	4-2
963261	Bracket - Tensioner	4-2

Numerical Index

Part No.	Description	Figure No.
963274	Bushing - Thrust	2-3, 2-4, 2-8
963312	Magnet - Top Cover	1-1, 4-5, 6-3
963914	Cover - Solid State Relay	6-1
963916	Display - Developer Temperature, M43 and M43A	1-2, 6-4
963923	Terminal Block - 6 Position	6-1
968757	Standoff - Nylon, 1 in.	6-1
968761	Knob - Drain Valve	5-4
968762	Knob - Dryer Temperature Control	1-2, 6-4
968763	Terminal Block	6-4
968766	Pump Assembly - Replenishment, 115 V	5-2
968766	Pump Assembly - Replenishment, 115 V	5-8, 5-9
968767	Pump Assembly - Replenishment, 230 V	5-2
968767	Pump Assembly - Replenishment, 230 V	5-8, 5-9
968796	Roller - Phenolic, Wash Rack	2-1
990010	Screw - Shoulder	4-2
991740	Relay - Solid State, 10 A	6-1
991743	Relay - Solid State, 20 A	6-1

Section 9: Alphabetical Index

Alphabetical Index

Part No.	Description	Figure No.
242072	Adapter	5-8, 5-9
918560	Alarm - Audio	1-1, 1-2, 6-4
261929	Anchor	7-1
696339	Baffle	4-5
779119	Bearing	2-3, 2-4, 2-8, 2-9, 2-10, 4-2, 4-3
918503	Bearing	2-3, 2-4
918512	Bearing	3-1
594273	Bearing - Spring Guide	2-3
1C7012	Bearing - Thrust	5-3
918511	Bearing - Thrust	3-1
1C7006	Bellows	5-2
963221	Belt - Drive	4-2
918504	Block - Bearing	2-8, 2-9
1C7011	Body - Pump	5-3
1C7005	Body - Valve	5-2
5B6451	Body - Valve	5-4
696144	Boot	6-2
240104	Bracket	7-1
5B6310	Bracket	4-5, 6-3
918234	Bracket	6-2
918515	Bracket	5-4
918631	Bracket	5-5
1C4012	Bracket - Drive Motor Mount	3-1
474596	Bracket - Guide Shoe	2-6
918516	Bracket - Motor	4-5
963261	Bracket - Tensioner	4-2
963274	Bushing - Thrust	2-3, 2-4, 2-8
5B6294	C-Clip	5-5
477506	Carrier - Fuse	6-1
1C0784	Chain	2-9
3B7037	Chain	3-1
696346	Circuit Board - 100	6-1
918552	Circuit Breaker - CB1, M43	1-2, 6-4
914441	Circuit Breaker - CB1, M43A and Clinic 1	1-2, 6-4
246800	Clamp	5-5, 5-10, 5-11, 7-4
246801	Clamp	5-10, 5-11
246808	Clamp - #10	5-8, 5-9, 5-10, 5-11
246804	Clamp - #12	5-8, 5-9
914908	Clamp - Narrow Band	5-8, 5-9, 5-10

Alphabetical Index

Part No.	Description	Figure No.
914908	Clamp - Thin Band	5-4
491893	Collar	3-1
520390	Connector - Elbow	5-2
963242	Connector - T	5-10, 5-11
963243	Connector - Y	5-10
696333	Cover - Developer Rack	2-2
613122	Cover - Evaporation	7-5
1C7008	Cover - Pump	5-1
963914	Cover - Solid State Relay	6-1
2B4915	Cover - Top, Clinic 1	1-1
918663	Cover - Top, Dryer	4-4
2B4910	Cover - Top, M43	1-1
1C7283	Cover - Top, M43A	1-1
594838	Cup - Displacement	5-2
918648	Cylinder - Graduated, 500 mL	7-5
618502	Diode	5-5
963916	Display - Developer Temperature, M43 and M43A	1-2, 6-4
5B6950	Drain Valve Complete	5-4
5B6460	Drain Valve Stem Assembly	5-4
8B6741	Dryer Assembly	4-1
166326	E-Ring	2-3, 2-4, 2-8, 2-9, 2-10, 4-2, 4-3
8B6926	EEPROM - 100 Circuit Board, Main Program Version 1.80	6-1
8B6925	EEPROM - 100 Circuit Board, Main Program, Version 1.80	6-1
8B6925	EEPROM - 100 Circuit Board, Main Program, Version 1.90	6-1
8B6926	EEPROM - 100 Circuit Board, Main Program, Version 1.90	6-1
918231	Exchanger - Heat	5-10, 5-11
1C4017	Fan	3-1
918643	Film Detector Assembly	6-5
918537	Filter - Line	6-4
492578	Fitting - Bulkhead	7-4
261898	Floor Plate	7-1
569165	Fuse - 1 A, 230 V, slow blow	6-1
126665	Fuse - 2 A, 115 V, slow blow	6-1
136963	Fuse - 3/4 A, 115 V, slow blow	6-1
337468	Fuse - 3/8 A, 230 V, slow blow	6-1
650938	Gasket	7-5
963229	Gasket	2-10
1C7785	Gasket - Tank/Blower	4-5
458384	Gear	2-8, 2-9, 4-2
696334	Gear	2-3
1C0373	Gear - 1.0 in.	2-3, 2-5
1C0534	Gear - 1.25 in.	2-5
550774	Gear - Drive	2-3, 2-4, 2-8

Alphabetical Index

Part No.	Description	Figure No.
458385	Gear - Idler	2-9
1C0536	Gear - Snap, 1.0 in.	2-3, 2-4, 2-5
1C0537	Gear - Snap, 1.25 in.	2-4, 2-5
598808	Grounding Lug	6-4
1C4019	Guard - Splash	7-5
1C7748	Guide Shoe	1-1
1C0370	Guide Shoe Assembly - Entrance, Lower	2-6
696332	Guide Shoe Assembly - Entrance, Upper	2-6
1C7448	Guide Shoe Assembly - Exit, Lower	2-6
914703	Guide Shoe Assembly - Exit, Upper	2-6
696146	Harness - Air Vane Switch	4-5
918937	Harness Assembly	5-5
918548	Heater - Developer, 115 V	5-7
918540	Heater - Developer, 230 V	5-7
1C7000	Heater - Dryer	4-5
943366	Holder - Fuse	6-1
1C7016	Housing - Impeller, and Adapter	5-3
1C7009	Housing - Pump	5-1
458949	Idler Roller Assembly	4-2
1C7010	Impeller	5-1
1C7014	Impeller and Magnet Assembly	5-3
918566	Indicator Assembly	1-2, 6-4
292888	Jack - ESD	6-1
968761	Knob - Drain Valve	5-4
968762	Knob - Dryer Temperature Control	1-2, 6-4
696323	Label - Bin, M43 and M43A	1-1
963256	Label - Display, Clinic 1	1-2
955975	Label - Display, M43 and M43A	1-2
955976	Label - Drain	1-1
963312	Magnet - Top Cover	1-1, 4-5, 6-3
1C0792	Manifold Box	2-3
335503	Microprocessor	6-1
1C7002	Motor - 115 V	3-1
1C7003	Motor - 230 V	3-1
918544	Motor - Blower, 115 V	4-5
918570	Motor - Blower, 230 V	4-5
5B6287	Motor - Wash Pump, 115 V	5-3
5B6288	Motor - Wash Pump, 230 V	5-3
5B6289	Motor Assembly - 115 V	5-2
5B6290	Motor Assembly - 230 V	5-2
166331	Nut	1-1, 7-1
852619	Nut	4-5, 6-3, 6-5
853080	Nut	4-5, 5-5, 6-1, 6-2, 6-4
853136	Nut	5-6

Alphabetical Index

Part No.	Description	Figure No.
520391	Nut - Connector	5-2
594839	Nut - Coupling	5-2
143403	Nut - Lock	4-2
1C7015	O-Ring	5-3
242038	O-Ring	1-1
453342	O-Ring	5-6
463363	O-Ring	5-1
520388	O-Ring	5-2
535639	O-Ring	5-2
658823	O-Ring	5-10, 5-11
786449	O-Ring	2-10
918233	O-Ring - Heater	5-7
1C0545	Orifice - Developer	5-8
1C0546	Orifice - Fixer	5-9
1C0547	Orifice - Wash	5-11
2B4914	Panel - Back Access, Clinic 1	1-1
2B4909	Panel - Back Access, M43 and 43A	1-1
2B4911	Panel - Drive Side Access, Clinic 1	1-1
2B4906	Panel - Drive Side Access, M43 and M43A	1-1
2B4913	Panel - Front Access, Clinic 1	1-1
2B4908	Panel - Front Access, M43 and M43A	1-1
2B4912	Panel - Non-Drive Side Access, Clinic 1	1-1
2B4907	Panel - Non-Drive Side Access, M43 and M43A	1-1
184455	Pin	2-3
502464	Pin	2-3, 2-4, 2-8
535646	Pin	2-3
544081	Pin	2-8, 2-9, 4-2
852978	Pin	2-10
853000	Pin	3-1
918229	Plate	5-10, 5-11
459548	Plate - Floor	1-1
1C7664	Plate - Locking	2-9
963209	Plenum	4-5
696338	Plug	5-11
918237	Plug	5-6
918238	Plug	5-7
918554	Potentiometer - Dryer Heater	1-2, 6-4
1C7282	Power Supply Assembly	6-1
963220	Pulley	4-2
963260	Pulley - Idle	4-2
918543	Pump Assembly - Recirculation, 115 V	5-1, 5-8, 5-9
918569	Pump Assembly - Recirculation, 230 V	5-1, 5-8, 5-9
968766	Pump Assembly - Replenishment, 115 V	5-2
968766	Pump Assembly - Replenishment, 115 V	5-8, 5-9
968767	Pump Assembly - Replenishment, 230 V	5-2

Alphabetical Index

Part No.	Description	Figure No.
968767	Pump Assembly - Replenishment, 230 V	5-8, 5-9
918545	Pump Assembly - Wash Recirculation, 115 V	5-3, 5-10, 5-11
918538	Pump Assembly - Wash Recirculation, 230 V	5-3, 5-10
1C0352	Rack - Developer	2-2
1C0351	Rack - Fixer	2-2
1C7722	Rack - Wash	2-7
991740	Relay - Solid State, 10 A	6-1
991743	Relay - Solid State, 20 A	6-1
244917	Retainer	1-1
1C7007	Rod - Connecting	5-2
1C4022	Roller	2-1, 2-3, 2-4
1C7690	Roller	2-1, 2-9
489693	Roller	2-1, 2-3, 2-4, 2-5
594237	Roller - Idler	4-3
246929	Roller - Phenolic	2-10
246979	Roller - Phenolic	2-1, 4-2, 4-3
467016	Roller - Phenolic, 1 in.	2-1, 2-10
8B6568	Roller - Phenolic, 11/16 in.	2-1, 2-10
968796	Roller - Phenolic, Wash Rack	2-1
470086	Roller - Resilient	2-1, 2-3, 2-4, 2-8
594237	Roller - Resilient	2-1, 2-10, 4-2, 4-3
1C0780	Roller - Rewet	2-1, 2-3
963235	Roller - Silicone, Wash Rack	2-1, 2-9
5B6293	Screen	5-5
154242	Screw	6-1
165626	Screw	4-5
168843	Screw	5-1
186171	Screw	5-5
372056	Screw	4-5
3B7276	Screw	5-7, 5-10, 5-11
544194	Screw	2-2
556032	Screw	4-5, 6-3
575821	Screw	6-4, 6-5
5B6316	Screw	4-2
619284	Screw	5-1, 5-2, 5-3, 5-5
796386	Screw	5-11
852682	Screw	2-5, 3-1, 4-4, 4-5
852683	Screw	6-2, 6-4
852777	Screw	1-1
852805	Screw	1-1, 2-10, 3-1, 5-10, 5-11, 6-5

Alphabetical Index

Part No.	Description	Figure No.
852807	Screw	2-3, 2-6, 2-10, 4-4
852830	Screw	4-5
852871	Screw	2-9
852941	Screw	2-6
852943	Screw	5-6
852976	Screw	6-4
1C4848	Screw - Leveling, Long	1-1, 7-1
459229	Screw - Leveling, Short	1-1, 7-1
781202	Screw - Nylon	6-1
1C0785	Screw - Shoulder	2-9
990010	Screw - Shoulder	4-2
458386	Screw - Shoulder, Small	2-9
188623	Screw - Stainless Steel	7-1
261413	Seismic Bracket Kit	7-1, 7-3
918557	Sensor - Film	6-5
852553	Setscrew	4-2, 5-4
852558	Setscrew	3-1
1C7691	Shaft	2-9
473955	Shaft	5-1
696328	Shaft	2-4
918630	Shaft	3-1
963217	Shaft	4-3
1C0369	Shaft - Drive	2-8
918626	Shaft - Drive, Long	2-3, 2-4
1C4018	Shaft - Drive, Short	2-3, 2-4
918625	Shaft - Driven	2-8
458391	Shaft - Idler	4-3
1C0776	Shaft - Lower	2-10
963218	Shaft - Lower Dryer	4-2
1C7013	Shaft - Pump	5-3
458391	Shaft - Roller	4-2
1C0548	Shaft - Tie Rod	2-10
1C0549	Shaft - Upper	2-10
696326	Sideplate - Drive	2-3, 2-5
963206	Sideplate - Drive	4-2
696327	Sideplate - Non-Drive	2-3, 2-5
963207	Sideplate - Non-Drive	4-2
963233	Sideplate, Drive	2-10
963224	Sideplate, Non-Drive	2-10
918641	Solenoid - Water Input	5-5, 5-10, 5-11
1C4849	Spacer	7-1
1C7700	Spacer	4-2, 4-3
511118	Spacer	2-4
918506	Spacer	2-8

Alphabetical Index

Part No.	Description	Figure No.
963237	Spacer	2-9
963238	Spacer	2-8, 2-9
963239	Spacer	2-8, 2-10
963240	Spacer	2-8
963241	Spacer	2-9
246580	Spring	2-4, 2-9
459344	Spring	2-8, 2-10, 4-2, 4-3
590309	Spring - Short	2-3
1C0777	Sprocket	2-9, 2-10
3B7040	Sprocket	3-1
3B7039	Sprocket, 50 Hz	3-1
3B7038	Sprocket, 60 Hz	3-1
968757	Standoff - Nylon, 1 in.	6-1
472261	Strainer Assembly	7-4
943381	Strap - EMI	6-1
918559	Strip - Marker	6-4
696148	Strip - Marker, 6 Position	6-1
562228	Stud	2-3, 2-4, 2-5
746731	Switch - AC Interlock S5	6-2
943375	Switch - Air Vane	4-5
914599	Switch - DC Interlock S6	4-5, 6-3
968763	Terminal Block	6-4
963923	Terminal Block - 6 Position	6-1
918549	Thermistor - Developer	5-6
918550	Thermistor - Dryer	4-5
458219	Tie Rod	2-6
458394	Tie Rod	2-10, 4-4
1C4009	Transformer - 20 V AC	6-1
1C4011	Tray - Drip	7-5
2B4916	Tray - Feed	1-1
963208	Tube - Air	4-4
1C0783	Tube Assembly	5-10, 5-11
8B6860	Tube Assembly	5-10, 5-11
3B7240	Tubing - 1/4-in. ID, 5 ft. length	5-8, 5-9
760476	Tubing - Bulk, 1/2-in. ID	7-5
760477	Tubing - Bulk, 3/4-in. ID, Drain	5-8, 5-9, 5-10, 5-11, 7-5
452990	Tubing - Bulk, 3/8-in. ID	5-8, 5-9, 7-5
650656	Tubing - Drain, 18 ft., 3/4-in. ID	7-5
1C0353	Tubing - Formed	5-10
1C0354	Tubing - Formed	5-10
1C0355	Tubing - Formed	5-10
1C0356	Tubing - Formed	5-10
1C0357	Tubing - Formed	5-10, 5-11

Alphabetical Index

Part No.	Description	Figure No.
1C0358	Tubing - Formed	5-10
1C0359	Tubing - Formed	5-10
1C0360	Tubing - Formed	5-10
1C0361	Tubing - Formed	5-8
1C0362	Tubing - Formed	5-8
1C0363	Tubing - Formed	5-11
1C0365	Tubing - Formed	5-11
1C0366	Tubing - Formed	5-11
1C0367	Tubing - Formed	5-10, 5-11
918242	Tubing - Formed	5-9
918243	Tubing - Formed	5-9
918245	Tubing - Formed	5-10, 5-11
918248	Tubing - Formed	5-8
918249	Tubing - Formed	5-9
918674	Tubing - Formed	5-8
963212	Tubing - Formed	5-9
963213	Tubing - Formed	5-9
963215	Tubing - Formed	5-10, 5-11
1C4010	Valve - Diverter	5-5, 5-10
8B6945	Valve - Wash Drain, Orifice	5-4
607817	Valve Assembly	5-2
943378	Varistor - MOV	6-1
1C0791	Varistor - MOV Assembly	6-1
264503	Ventilation Fan Kit - 110 V	7-3
963231	Wash Tube - Lower	2-10
963230	Wash Tube - Upper	2-10
8B6948	Wash Valve Stem Assembly - with Orifice	5-4
459343	Washer	2-5
473954	Washer	5-1
482271	Washer	4-2
559793	Washer	1-1
852611	Washer	2-3, 3-1, 4-4, 5-5
852799	Washer	1-1, 2-9, 3-1, 4-5, 5-1, 5-2, 5-3, 5-5, 6-2, 6-4, 6-5
886890	Washer	2-2
960420	Washer	4-2, 5-6
575822	Washer - Flat	4-5, 6-3, 6-4
191098	Washer - Lock	4-5, 6-3, 6-4
535004	Washer - Lock	1-1, 2-10, 4-5, 5-1, 5-2, 5-3, 5-5, 5-6, 5-10, 5-11, 6-2, 6-4, 6-5

Alphabetical Index

Part No.	Description	Figure No.
583100	Washer - Lock	6-4
853004	Washer - Lock	3-1
853011	Washer - Lock	5-6
853087	Washer - Lock	2-10
853096	Washer - Lock	2-2, 4-5, 6-4
853100	Washer - Lock	2-3, 2-6, 3-1, 4-4, 4-5, 5-5
173824	Washer - Stainless Steel	7-1
433708	Washer - Star	4-5, 6-1, 6-2
443700	Washer - Star	6-4
166956	Wire Tie	1-1
182915	Wire Tie	6-2
182916	Wire Tie	6-2
240059	Wire Tie - blue	2-3, 5-9
240058	Wire Tie - red	2-4, 5-8
748956	Worm	3-1

Section 10: Publication Change Notice Table

Print Date	Pub. No.	ECO No.	Affected Pages	File Name	Notes
Jan. 1994	981091	2622-064	All Pages	3229pl_b.txt	1st Printing
March 1994	981091	2622-082	Front Cover, 3, 4, 2-3, 3-5, 3-18, 3-19, 5-3, 6-9, 6-13, 6-21, 7-3, 9-1 through 10-7, 11-1, Back Cover	3229pl_b_082.txt	
April 1994	981091	2622-085	Front Cover, 2, 1-2, 2-2, 2-3, 5-8, 5-9, 7-2 through 7-4, 9-1 through 10-7, 11-1, Back Cover	3229pl_b_085.txt	Added DC Interlock Switch S6
May 1994	981091	2622-098	Front Cover, 3-10, 3-11, 11-1, Back Cover	3229pl_b_098.txt	Update the Print Date
Oct. 1994	981091	2622-117	All Pages	3229pl_b_117.txt	Update and Revise the Manual to Bring it to Final Release Status
Sept. 1995	981091	2622-157	All Pages	ip3229_1_157.book	Revision due to product changes; reprinted entire manual due to new publication software used for documentation

Kodak and X-Omat are trademarks.

